*Breast Biopsy System*

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**BAIRD**

**BIOPSY SYSTEMS**

A Business Unit of  
Bard Peripheral Vascular, Inc.

## Instructions for Use

**Caution:** Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

### A. Device Description

The FINESSE™ ULTRA Breast Biopsy System is composed of a driver and a probe and allows for obtaining and storing multiple samples with a single probe insertion. The components of the system are designed to operate safely when used together for diagnostic sampling of tissue during a biopsy procedure. The driver is a self-contained, handheld, reusable electro-mechanical vacuum-assisted biopsy device (Figure 1) that utilizes a rechargeable lithium-ion battery and comes complete with AC power accessories (Figure 2) and carrying case. The biopsy probe is composed of a hollow outer cutting cannula, an inner needle with a notched sample chamber, and an integrated coaxial cannula. The outer cannula and integrated coaxial cannula have sharpened cutting edges and the inner needle has a sharpened trocar tip. The driver utilizes DC operated motors to create a vacuum while simultaneously moving the outer cannula backward to expose the sample chamber to draw tissue in. The outer cannula is released and moves forward over the sample chamber to cut the tissue. The tissue sample is transported via the sample chamber to the sample container. The sample chamber is returned to the original position and is available to acquire another sample. The integrated coaxial cannula may be detached after the biopsy and remain in the breast to retain a track to the biopsy site when placing a tissue marker.

Figure 1



Figure 2

### B. Indications for Use

The FINESSE™ ULTRA Breast Biopsy System is indicated to obtain tissue samples from the breast or axillary lymph nodes for diagnostic analysis of breast abnormalities. The instrument is intended to provide breast tissue for histologic examination with partial or complete removal of the imaged abnormality.

The extent of histologic abnormality cannot be reliably determined from its mammographic appearance. Therefore, the extent of removal of the imaged evidence of an abnormality does not predict the extent of removal of a histologic abnormality e.g., malignancy. When the sampled abnormality is not histologically benign, it is essential that the tissue margins be examined for completeness of removal using standard surgical procedures.

### C. Contraindications

1. The FINESSE™ ULTRA Breast Biopsy System is for diagnostic use only, NOT for therapeutic use.
2. The FINESSE™ ULTRA Breast Biopsy System is contraindicated for those patients where, in the physician's judgment, there is an increased risk of complications associated with percutaneous removal of tissue samples.

### D. Warnings

1. Patients who may have a bleeding disorder or who are receiving anticoagulant therapy may be at increased risk of complications.
2. As with any biopsy instrument, there is a potential for infection.
3. The driver should not be used in a MRI suite.
4. The driver charger should not be used in an operating room.

5. The driver is not classified as an AP or APG classified device. The driver is not suitable for use in the presence of flammable anesthetic.
6. The driver must only be used with FINESSE™ ULTRA Biopsy Probes.
7. Do not reuse FINESSE™ ULTRA Biopsy Probes. Reusing the probe bears the risk of cross-patient contamination as biopsy probes – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the probe for an indeterminable period of time. The residue of biological material can promote the contamination of the probe with pyrogens or microorganisms which may lead to infectious complications.
8. Do not resterilize FINESSE™ ULTRA Biopsy Probes. After resterilization, the sterility of the probe is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilizing the probe increases the probability that it will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.
9. All breast biopsies should be performed under imaging guidance to confirm the probe position relative to the target region to be sampled and to help mitigate the occurrence of a false negative biopsy.

## E. Precautions

1. Do not use the FINESSE™ ULTRA Biopsy Probe without the integrated coaxial cannula. The integrated coaxial cannula may be removed after the biopsy to retain a track to the biopsy site when placing a tissue marker.
2. The FINESSE™ ULTRA Breast Biopsy System should only be used by a physician trained in its indicated use, limitations, and possible complications of percutaneous needle techniques.
3. Do not demonstrate sample acquisition sequence in air prior to performing a biopsy procedure. Damage may occur to the needle or cannula tip.

## F. Potential Complications

1. Potential complications are those associated with any percutaneous removal/biopsy technique for tissue collection. Potential complications are limited to the region surrounding the biopsy site and include hematoma, hemorrhage, infection, a non-healing wound, pain and tissue adherence to the biopsy probe while removing it from the breast.
2. As per routine biopsy procedures, it may be necessary to cut tissue adhering to the biopsy probe while removing it from the breast.

## G. Equipment Required

The following equipment is required for a biopsy procedure:

- Appropriate imaging modality and accessories
- FINESSE™ ULTRA Biopsy Driver
- FINESSE™ ULTRA Biopsy Probe
- Surgical gloves and drapes
- Local anesthetic
- Scalpel
- Other equipment as necessary


## H. Directions for Use

### Driver Button Operations (Figure 3):

Two buttons operate the driver (**PRIME/PIERCE** and **SAMPLE**). Lights indicate the operational status of the driver. The driver has an optional piercing mode that retracts the probe 20-22mm and fires the probe forward 20-22mm to assist in penetrating dense breast tissue.



- A continuously illuminated green light indicates an operational step that can be started using the corresponding button on the key pad.
- A continuously illuminated orange light indicates the device is currently engaged. When engaged, the device should not be moved or changed in any way.
- When no light is illuminated, no operation/action can be performed using the corresponding button on the key pad.
- Flashing orange lights indicate the driver has experienced an error that may be reset by the user. To clear the flashing orange lights error see Troubleshooting.

- Flashing red lights indicate the driver has experienced an error that cannot be reset with the probe in place. To clear the flashing red lights error, see Troubleshooting.
- One of two operational buttons can be selected when both the “PRIME/PIERCE” and the “SAMPLE” indicator lights are continuously illuminated green.
- The “PRIME/PIERCE” button initiates the optional piercing sequence. When pressed, the “PRIME/PIERCE” indicator light continuously illuminates orange indicating the probe is being retracted approximately 20-22mm inside the driver, the “SAMPLE” indicator light is off and the “SAMPLE” button is deactivated. Once the probe is retracted, the “PRIME/PIERCE” indicator light continuously illuminates green indicating the probe is primed, the “SAMPLE” indicator light is off and the “SAMPLE” button is deactivated. Pressing the “PRIME/PIERCE” button again advances the probe approximately 20-22mm into the breast. Using the appropriate imaging guidance, the probe's post-pierce position can be verified. After the piercing sequence is complete, the “PRIME/PIERCE” and “SAMPLE” indicator lights continuously illuminate green indicating another piercing sequence may be performed or a sample may be acquired.
- The “SAMPLE” button initiates the sampling sequence. When pressed, the “SAMPLE” indicator light continuously illuminates orange indicating the sampling sequence is in operation, the “PRIME/PIERCE” indicator light is turned off and the “PRIME/PIERCE” button is deactivated. During the sampling process, the vacuum is turned on, the probe's cutting cannula is pulled back exposing the sample chamber, tissue is drawn into the sample chamber, the cutting cannula is released to cut the tissue, the tissue sample is transported to the sample container and the system resets for the next sampling or piercing sequence.

#### Loading the Probe into the Driver

1. Ensure that the driver has been fully charged. The driver should be charged overnight after daily use.
2. Remove the probe from the package using aseptic technique by grasping the probe from below. Ensure that the probe remains sterile.
3. Insert the probe into the driver (Figure 4).
4. Remove the retaining pin (Figure 5).
5. Slide the probe carriage back ensuring the tab securely locks into place (Figure 6).

**NOTE:** The driver does not contain an on-off switch. The driver automatically recognizes the probe type, senses when the probe has been properly inserted and performs an initialization to synchronize the driver and probe. The driver sounds one beep and the “PRIME/PIERCE” and “SAMPLE” indicator lights display green indicating the probe is properly inserted. During the initialization the motor is audible.

6. After initialization is complete, remove the protective sleeve from the probe.

**NOTE:** If the driver is left undisturbed for 120 seconds after loading a probe, the driver will enter a sleep mode to conserve battery charge. During sleep mode, the “PRIME/PIERCE” and “SAMPLE”

indicator lights will turn off and the buttons will be deactivated. Once the driver is moved, a motion sensor is activated and the driver automatically turns on. The driver will not enter sleep mode if the probe is primed.

### Performing a Biopsy

1. Prepare the biopsy site using appropriate aseptic techniques and adequate local anesthesia. Appropriate precautions should be followed for any percutaneous procedure (personal protective gear such as gloves and glasses/shields, etc. should be worn).
2. Make a small incision or puncture the skin with a scalpel.
3. Using imaging guidance, insert the probe through the incision and position the tip to the appropriate depth:
  - If the lesion is going to be pierced: Position the tip of the probe proximal to the edge of the lesion.
  - If the lesion is not going to be pierced: Position the sampling chamber to the center of the target region. (The center of the sampling chamber is approximately 17mm from the tip of the probe.)
4. If piercing is desired, prior to inserting the probe into the breast press the **"PRIME/PIERCE"** button to retract the probe 20-22mm. After inserting the probe press the **"PRIME/PIERCE"** button to advance the probe 20-22mm into the region of interest prior to taking the first tissue sample.

**NOTE:** Image guidance should be used to confirm probe position relative to the target region to be sampled.

**NOTE:** Do not leave the system unattended with the probe in the primed state.

**NOTE:** During the initial positioning of the probe, it is recommended to advance the probe into the tissue using the optional piercing mode. During subsequent samplings, the physician may wish to forego this step.

**NOTE:** The probe may be primed for piercing while inside the patient if additional piercing shots are desired to access dense tissue or lesions.

5. Begin tissue sampling by pressing the **"SAMPLE"** button. The sampling process proceeds automatically: A vacuum is created, the outer cannula of the probe is automatically retracted, tissue is drawn by vacuum into the sampling chamber and the outer cannula is released cutting the tissue. The tissue sample is then transported to the sample container. During this process the **"SAMPLE"** indicator light continuously illuminates orange and the motor will be audible. The sample container illuminates to aid in visualizing the tissue sample. After the tissue sample is deposited into the sample container the system resets to acquire another sample if desired.

**NOTE:** Do not press on the sample container or remove the sample container while acquiring tissue samples. If the system detects an error related to the sample container it will illuminate the **SAMPLE** indicator flashing orange and sound 2 beeps. Remove the sample container from the probe, clear any tissue underneath the container or from the probe where the container is seated, insert the container back into the probe. The **SAMPLE** light will flash green to indicate that the sample container is properly inserted. Press the **SAMPLE** button to resume. The system will show the **SAMPLE** light continuous green. If error is not resolved see Troubleshooting.


6. To acquire multiple samples, repeat Step 5 above.

**NOTE:** Collecting multiple tissue samples may help minimize the risk of a false negative biopsy.

7. After obtaining the final tissue sample, remove the probe from the breast and care for the incision site as appropriate.


**NOTE:** Prior to removing the probe from the breast, the integrated coaxial cannula may be detached and remain in the breast to retain a track to the biopsy site while placing a tissue marker (Figure 7).

Figure 7

8. Remove the tissue samples from the sample container by twisting the lid counter clockwise. If additional samples are required, ensure the sample container lid is in place while acquiring samples.
- NOTE:** The sample container may be removed completely from the probe by pulling up on the lid and may be removed before or after the probe is removed from the driver. If the probe is loaded into the driver and the sample container is removed from the probe, the "PRIME/PIERCE" and "SAMPLE" indicator lights will flash orange. After the sample container is re-inserted into the probe, the "PRIME/PIERCE" and "SAMPLE" indicator lights will continuously illuminate green.
9. Remove the probe from the driver by pressing down on the locking tab, pushing the carriage completely forward and pulling the probe straight from the driver (Figure 8).

Figure 8

**NOTE:** Once a probe has been removed from the driver following a procedure, the probe cannot be re-inserted. If additional tissue samples are required, insert a new probe into the driver.

**NOTE:** Driver will initiate mechanical reset following probe removal. Do not insert a new probe prior to completion of this reset. Five beeps and "PRIME/PIERCE" and "SAMPLE" indicators flashing orange, will warn the user if a probe is inserted before the reset is completed.

10. Refer to the specific FINESSE™ ULTRA Biopsy Probe Instructions for Use for additional important information.

**NOTE:** The probe, coaxial cannula and sample container are for single use only. After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

## I. Cleaning and Maintenance

- After each use, clean the driver thoroughly with a germicidal disposable wipe such as Sani-Cloth® or a damp cloth moistened with a 10% bleach solution. Other than routine cleaning, no other user level inspection or preventive maintenance is required. However, it is recommended to return the device to an authorized Bard Service & Repair Facility once per year for inspection and preventive maintenance.
- Do not spray the driver with any fluids. Spraying the driver may cause it to malfunction and will void the warranty.
- Do not submerge the driver in fluids. Submerging the driver in fluids may cause it to malfunction and will void the warranty.
- Sterilization and exposure to fluids may damage the electrical components of the device. If the driver is inappropriately cleaned it may cause the driver to malfunction and will void the warranty.
- Do not autoclave. Do not heat over 54°C (129°F).

## J. Driver Charging


The driver utilizes a rechargeable lithium-ion battery. Completely charge the driver prior to first use. The battery must be adequately charged prior to using the driver for the first time as well as prior to beginning each procedure. The battery charge level may be checked by pressing the "SAMPLE" button while the driver is not loaded with a biopsy probe and is not in the charging stand (Figure 2). If the battery is sufficiently charged to perform a biopsy procedure, the battery indicator light will continuously illuminate green for 5 seconds. If the battery is not sufficiently charged to perform a biopsy procedure, the battery indicator light will continuously illuminate red for 5 seconds and the driver must be charged. If the battery is completely discharged, the battery light will not illuminate and the driver must be charged.

To charge the driver, follow these instructions:

1. Ensure the appropriate adapter plug is connected to the AC power adapter.
2. Connect the AC adapter cord to the charging stand.

3. Plug the AC power adapter into a power source. If the AC power adapter has current, the green power light on the rear of the charging stand will immediately illuminate green and remain illuminated.
4. Place the driver in the charging stand by first placing the front of the driver on the stand and then rotating the driver onto the charging stand as shown in Figure 9. If the battery is not sufficiently charged to perform a biopsy, the battery indicator light will continuously illuminate red. If the battery is sufficiently charged to perform a biopsy, but is not fully charged, the battery indicator light will flash green.

Figure 9

5. Once the battery is fully charged, the battery indicator light will continuously illuminate green. A completely exhausted battery will require approximately 120 minutes to fully charge. After every use or when the driver is not in use, the driver should remain in the charging stand to keep the battery fully charged. A regulator restricts overloading the battery.


**NOTE:** Use only with supplied AC power accessories. The power cord is to be used for mains disconnection.

**NOTE:** If during a biopsy procedure the battery indicator light illuminates orange, the battery has sufficient charge to complete the biopsy procedure, but the driver must be charged prior to starting another biopsy procedure.

**NOTE:** The battery may only be replaced by an authorized Bard Service and Repair facility.

**NOTE:** Do not pull up or out on the driver when removing it from the charging stand. Rotate the driver on and off of stand pivoting on bottom of charging stand as shown in Figure 10. Improper loading of driver on and off of charging stand may result in damage to the driver or charging stand.

Figure 10

## K. Troubleshooting

### Driver Keypad Light and Audible Alarm Conditions

PP = PRIME / PIERCE button, S=SAMPLE button

LED	Color / Status	Beeps	Condition	Action Required
No lights		None	No probe loaded or driver in sleep mode if not moved for 120 seconds.	<p>Move driver or place driver on charging stand. PP and S lights will turn solid green.</p> <ul style="list-style-type: none"> <li>If the indicator lights do not turn on after loading a probe into the driver, ensure the probe cover and locking tab are securely in place and ensure the probe is new and has never been used before.</li> <li>If the above conditions are satisfied, remove the probe and check to see if the driver has sufficient charge by pressing and holding the S button. If the Battery indicator light indicates sufficient charge, replace the probe with a new probe.</li> <li>If the indicator lights do not turn on after properly inserting a new probe, contact Bard Service and Repair.</li> </ul> <p>See section J, Driver Charging, for Battery light and charging stand instructions.</p>
PP and S	Green / Solid	1 beep (only when loading probe)	Probe loaded in driver.	Driver ready to perform biopsy. Option to either "Prime" or "Sample."
PP and S	Orange / Flashing alternating	None	Sample container removed.	Ensure sample container is installed in probe.
PP and S	Orange / Flashing simultaneously	2 beeps	Tissue sample stuck in sample chamber or probe chamber is open and probe is stuck in driver.	<p>If procedure is not yet complete, detach coaxial from probe to maintain location in breast. Remove probe from the breast.</p> <ul style="list-style-type: none"> <li>If sample chamber is closed, press the PP and S buttons simultaneously for at least 3 seconds until the chamber opens. Remove any tissue from sample chamber and then press and hold the PP and S buttons simultaneously for at least 3 seconds. If both the PP and S lights are solid green, the biopsy can be continued. Insert the probe and reattach the coaxial to continue the biopsy.</li> <li>If the sample chamber is open, remove any tissue from the chamber and then press and hold the PP and S buttons simultaneously for at least 3 seconds. If both the PP and S lights are solid green, the biopsy can be continued.</li> </ul> <p>If the PP and S lights are flashing red after holding the PP and S buttons simultaneously for at least 3 seconds, discard the biopsy probe, reset the driver, and replace with a new probe if additional samples are required.</p>
PP and S	Orange / Flashing simultaneously	Multiple beeps	Keypad button held down after completing prime, pierce, or sample acquisition sequence.	Release the keypad button. Either PP or PP and S button(s) will turn solid green, and the biopsy procedure can be continued.
PP and S	Orange / Flashing Alternating	5 beeps	Probe inserted before driver has reset.	Remove probe from driver, allow driver to reset and reload probe after driver has reset.

LED	Color / Status	Beeps	Condition	Action Required
PP and S	Red / Flashing alternating	3 beeps	Probe or driver failure.	Remove and discard probe. Reset driver by pressing the PP and S button simultaneously for at least 3 seconds until the motors engage. If reset does not occur and lights remain flashing red, place driver on charging stand and wait for driver to reset (PP and S lights will cycle from red to orange to green to off). If the indicator lights do not turn off or the probe cannot be removed, contact Bard Service and Repair.
PP	Green / Solid	None	Driver primed (probe retracted).	Press "PP" when in desired location to pierce target.
PP	Orange / Solid	None	Device is performing either the prime or pierce function.	Wait for driver to complete function and indicator light to turn solid green before moving or repositioning.
S	Orange / Solid	None	Device is engaged in sample acquisition sequence.	Wait for driver to complete sample acquisition sequence and indicator light to turn solid green before moving or repositioning.
S	Orange / Flashing	2 beeps	Sample cup obstructed.	Remove sample cup. Inspect cup base and below cup and remove any tissue outside cup. Replace sample cup. S light will flash green when obstruction is resolved. Press and hold S button for 3 seconds to reset probe. If S light remains flashing orange, discard the biopsy probe, reset the driver, and replace with a new probe if additional samples are required.
S	Green / Flashing	None	Sample cup has been reinserted (observed after S orange / flashing error).	Press and hold sample button for at least 3 seconds to reset probe. S will turn solid orange while reset occurs. When reset is complete both PP and S will be solid green and the biopsy can be continued.

## L. Specifications

Electrical Conformity	This medical equipment has passed all required testing for electric shock, fire and mechanical hazards in accordance with UL 60601-1, IEC/EN 60601-1, CAN/CSA C22.2 No 601-1.
Electromechanical Compatibility	Portable and mobile radio frequency communications equipment (cell phones) should not be used at close distances. Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Classification	Class II, internally powered Type BF equipment.
Environmental Conditions	Transport and store in a cool dry place between -20°C to +50°C, 10% to 90% relative humidity and 500 to 1060 hPa atmospheric pressure.
Ingress of Water	Not protected against the ingress of water. IPX0.
Flammability	Equipment not suitable for use in the presence of flammable anesthetic.
Mode of Operation	Continuous.
Dimensions	24.0cm L X 4.5cm W X 5.0cm H
Weight	355g
AC Power Accessories	Use only with supplied AC power accessories: Power Input: 250 mA Input Voltage Range: 100 – 240 VAC Input Frequency Range: 50 – 60 Hz Output Voltage: 5 V Output Current Max: 1600 mA Max Watts: 10.0 W

### M. How Supplied

- The FINESSE™ ULTRA Biopsy Driver is supplied nonsterile and should be cleaned prior to each use.
- The FINESSE™ ULTRA Biopsy Probes are sold separately. The probes are supplied sterile for single patient use only.

### N. Warranty

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

**TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.**

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

Extended warranty service agreements are available. Consult a Bard representative for details regarding terms and conditions.

An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

Assembled in Poland. Charger and Power Adapters made in Germany.

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**Copyright ©2010 C. R. Bard, Inc. All rights reserved.**
Copyright ©2010 C. R. Bard, Inc. Tous droits réservés.
Copyright ©2010 C. R. Bard, Inc. Alle Rechte vorbehalten.
Copyright ©2010 C. R. Bard, Inc. Tutti i diritti riservati.
Copyright ©2010 C. R. Bard, Inc. Todos los derechos reservados.
Copyright ©2010 C.R. Bard, Inc. Alle rechten voorbehouden.