

Journal of Visualized Experiments

Laser capture microdissection of pure populations of neurons or single neurons for gene expression analysis after traumatic brain injury

--Manuscript Draft--

Manuscript Number:	JoVE50308R2
Article Type:	Invited Methods Article - JoVE Produced Video
Corresponding Author:	Helen Lee Hellmich, Ph.D. University of Texas Medical Branch Galveston, TX UNITED STATES
First Author:	Deborah R. Boone, BS
Order of Authors:	Deborah R. Boone, BS
	Stacy L. Sell, Ph.D
	Helen Lee Hellmich, Ph.D.



Helen L. Hellmich, PhD
Associate Professor

August 13, 2012

Dr. Angela Messmer-Blust
Science Editor
Journal of Visualized Experiments

Dear Dr. Messmer-Blust,

On behalf of my coauthors, I am submitting our manuscript "*Laser capture microdissection of pure populations of neurons or single neurons for gene expression analysis after traumatic brain injury*" for consideration in JoVE's online video journal.

We have used this laser capture microdissection (LCM) technique for over nine years to investigate the molecular mechanisms of brain injury in identified populations of brain cells. LCM has allowed us to gain valuable insight into cell survival and cell death mechanisms after traumatic brain injury hitherto not available through conventional microdissection and gene expression analysis techniques. JoVE's unique video format is perfect for demonstrating the fine art of laser microdissection which is difficult to visualize and replicate from merely published protocols.

Thank you for consideration of our manuscript.

Regards,



Helen Lee Hellmich, Ph.D
Associate Professor
Department of Anesthesiology

Laser capture microdissection of enriched populations of neurons or single neurons for gene expression analysis after traumatic brain injury

AUTHORS:

Deborah R. Boone

Department of Anesthesiology
University of Texas Medical Branch
601 Harborside Dr, Bldg #21
Galveston TX, USA
drkenned@utmb.edu

Stacy L. Sell

Department of Anesthesiology
University of Texas Medical Branch
601 Harborside Dr, Bldg #21
Galveston TX, USA
slsell@utmb.edu

Helen L. Hellmich (corresponding author)

Department of Anesthesiology
University of Texas Medical Branch
601 Harborside Dr, Bldg #21
Galveston TX, USA
hhellmic@utmb.edu
409-772-4216
Fax #: 409-772-6409

AUTHOR'S AFFILIATION: University of Texas Medical Branch

KEYWORDS: Laser capture microdissection; traumatic brain injury, hippocampus, Fluoro-Jade, gene expression analysis

SHORT ABSTRACT: We describe how to use laser capture microdissection (LCM) to obtain enriched populations of hippocampal neurons or single neurons from frozen sections of the injured rat brain for subsequent gene expression analysis using quantitative real time PCR and/or whole-genome microarrays.

LONG ABSTRACT: Long-term cognitive disability after TBI is associated with injury-induced neurodegeneration in the hippocampus—a region in the medial temporal lobe that is critical for learning, memory and executive function.^{1,2} Hence our studies focus on gene expression analysis of specific neuronal populations in distinct subregions of the hippocampus. The technique of laser capture microdissection (LCM), introduced in 1996 by Emmert-Buck, et al.,³ has allowed for significant advances in gene expression analysis of single cells and enriched populations of cells from heterogeneous tissues such as the mammalian brain that contains thousands of functional cell types.⁴ We use LCM and a well established rat model of traumatic brain injury

(TBI) to investigate the molecular mechanisms that underlie the pathogenesis of TBI. Following fluid-percussion TBI, brains are removed at pre-determined times post-injury, immediately frozen on dry ice, and prepared for sectioning in a cryostat. The rat brains can be embedded in OCT and sectioned immediately, or stored several months at -80°C before sectioning for laser capture. Additionally, we use LCM to study the effects of TBI on circadian rhythms. For this, we capture neurons from the suprachiasmatic nuclei that contain the master clock of the mammalian brain. Here, we demonstrate the use of LCM to obtain single identified neurons (injured and degenerating, Fluoro-Jade-positive, or uninjured, Fluoro-Jade-negative) and enriched populations of hippocampal neurons for subsequent gene expression analysis by real time PCR and/or whole-genome microarrays. These LCM-enabled studies have revealed that the selective vulnerability of anatomically distinct regions of the rat hippocampus are reflected in the different gene expression profiles of different populations of neurons obtained by LCM from these distinct regions. The results from our single-cell studies, where we compare the transcriptional profiles of dying and adjacent surviving hippocampal neurons, suggest the existence of a cell survival rheostat that regulates cell death and survival after TBI.

INTRODUCTION: Gene expression analysis of heterogeneous tissues has always been problematic; this is particularly true in the mammalian brain, which has approximately 5000 different cell types.⁴ Prior to the development of the laser capture microdissection (LCM) technique, genomic studies of the effects of TBI *in vivo* were based on analysis of gene expression in a mixed population of brain cells comprised, not only of different neuronal cell types, but also of supporting glial and immunomodulatory cells. The resulting complex gene expression profiles obtained from these heterogeneous tissues, and the often conflicting patterns of injury-induced cellular signals, may be one explanation for the failure in human clinical trials of therapeutic strategies proven effective in pre-clinical studies of TBI.⁵

To obtain a clear understanding of injury-induced gene expression in vulnerable populations of neurons from the rat hippocampus, we adopted the technique of LCM, first reported by Emmert-Buck et al.³ Subsequently, we modified and optimized this microdissection technique for efficient capture of enriched populations of neurons and single neurons for mRNA profiling using quantitative real time PCR and microarray analysis. To maintain the integrity of mRNA in frozen sections of brain tissue for genomic analysis, we modified existing protocols for fixing, staining and rapid capture of neurons from frozen rat brain sections. For identification and isolation of injured and dying hippocampal neurons, we also optimized the Fluoro-Jade staining technique for LCM. Fluoro-Jade does not distinguish between apoptotic and necrotic cell death. Thus, all types of degenerating neurons can be detected by this stain.^{6,7}

Here, we describe the protocol used in our laboratory to obtain pools of single dying or surviving neurons as well as swaths of enriched populations of distinct neuronal cell types (i.e., CA1-CA3 neurons for gene expression analysis after TBI). The procedure for fluid percussion brain injury performed in our laboratory is described in detail in Shimamura et al.,⁸ and is very similar to the lateral fluid-percussion injury protocol for mice published in JoVE by Alder et al.⁹ Since the LCM technique has been shown to have minimal or no effect on the integrity of DNA, RNA and protein in tissues, this is an excellent tool for molecular and protein analysis of defined cell types.

PROCEDURE:

1.) Surgical procedures and Fluid percussion TBI

- 1.1) All animal experiments are first approved by the Institutional Animal Care and Use Committee of the University of Texas Medical Branch, Galveston, Texas and the National Institutes of Health Guide for the Care and Use of Laboratory Animals (8th Edition, National Research Council).
- 1.2) Rats (adult male Sprague-Dawley rats, 400-500g obtained from vendor Charles Rivers, Portland, Maine) are housed two per cage and provided food and water ad libitum in a vivarium with these constant conditions: light cycle (600 h to 1800 h) temperature (21°C to 23°C), and humidity (40% to 50%) one week prior to use.
- 1.3) Anesthetize rats with 4% isoflurane, intubate, and mechanically ventilate (NEMI Scientific; New England Medical Instruments, Medway, MA) the rats with 1.5–2.0% isoflurane in oxygen: air (70:30) and prepare them for parasagittal fluid-percussion injury as previously described.^{8,10}
- 1.4) Sacrifice rats at the appropriate time point after injury depending on the experimental design. Quickly remove brains, freeze immediately on dry ice, and store at -80°C in a 50 ml tube or proceed immediately to embed in OCT for frozen sectioning.

2.) Sectioning and Staining of Rat Brain. (Note: the sectioning, staining and LCM procedures can all be filmed, although only the LCM part needs to be shown in detail. We suggest that 2.1-2-10 need not be filmed in their entirety but a snapshot of the sectioning and staining procedures would be sufficient. The subsequent RNA isolation procedures will just be described.)

- 2.1) Brain tissue that is not used immediately can be kept at -80°C for up to one month if kept at a constant temperature. Brains that are frozen at -80°C, then thawed to -20 °C for sectioning, and then re-frozen do not yield quality RNA after the second thawing. Once the brain is thawed, mounted in OCT and sectioned, slides should be stained within 24 hours and used for LCM within 30 minutes to one hour of staining. This will ensure good quality RNA. Following LCM, captured cells on LCM caps can be stored in lysis buffer at -80 °C for up to one week, but RNA should be isolated in a timely manner to ensure the highest quality.
- 2.2) Prior to sectioning the brain, wipe the cryostat with RNase-Zap and clean the brushes with ETOH (replace the disposable blade between each brain).
- 2.3) Retrieve the brain in the 50 ml tube from the -80°C freezer; place it into the cryostat at a temperature of -22°C and thaw in tube for approximately 10 minutes. Remove the brain from the tube and place onto the stage on gauze, ventral side up

- 2.4) Using a razor blade, slice the brain to remove the posterior portion of the brain just rostral to the cerebellum and the anterior portion at the optic chiasm. Fill a cryomold with OCT mounting medium (Tissue Tek), and place the brain into the mold with the anterior side down. Allow the brain tissue to freeze in the mounting medium until it turns white (approximately 10 minutes).
- 2.5) Freeze the specimen disc (Tissue Tek) onto the brain with OCT. Remove the brain from the mold. Insert the brain into the specimen head and tighten the screws.
- 2.6) Insert a disposable, low profile blade (Fisher Scientific) into the knife holder and tighten the lever down.
- 2.7) Begin slicing the brain at 20 microns to remove the outer layer of OCT. Once the hippocampal region is reached, set the micron dial to 10. Collect coronal serial sections by placing a glass slide or plus-glass slide onto the tissue section (Fisher Scientific). Slides are preserved at -20°C in an RNase-free staining rack until sectioning is complete.
- 2.8) To remove all RNases from glassware where tissue sections are processed, wipe down all staining containers and graduated cylinders with Eliminase (Fisher Scientific) and rinse in Milli Q water. Prepare all solutions with RNase-free water and filter the cresyl violet (Sigma-Aldrich) and Fluoro-Jade (Histo-chem) stain with a 0.2 micron filter prior to use.
- 2.9) Thaw brain sections at room temperature for 30 seconds and fix in 75% ETOH (1min).
- 2.10) For LCM of single injured neurons after fixation, rinse slides in RNase-free water (1 min), counterstain with 1% cresyl violet (15-20 sec), rinse in RNase-free water (2×30 sec), stain with Fluoro-Jade (4 min), rinse in RNase-free water (3×1 min), dehydrate with 95% ETOH **made from RNase-free water** (30 sec), 100% ETOH (30 sec), and xylene (2×3 min).
- 2.11) For LCM of swaths of neurons after fixation, rinse slides in RNase-free water (1 min), stain with 1% cresyl violet (1 min), rinse in RNase-free water for (3×1 min), dehydrate with 95% ETOH (2×30 sec), 100% ETOH (2×30 sec), and xylene (2×3 min).
- 2.12) Air dry all sections in a fume hood for 15 minutes prior to LCM. Slides can be stored, section sides up, in a slide box lined with dessicant, if they need to be transferred from room to room. However, optimal results are obtained if LCM is performed immediately after sections are dry. LCM should be limited from 30 minutes to 1 hour. In the next section, we describe first how to capture continuous swaths of cells, which is easiest to master for those who are learning this technique. Then we describe how to precisely capture single neurons. In our procedure, we identify dying neurons by their affinity for Fluoro-Jade—a marker of degenerating neurons. Fluoro-Jade-negative cells are presumed to be surviving neurons.

3.) Laser Capture Microdissection (LCM)

LCM of swaths of enriched populations of hippocampal neurons

- 3.1.) LCM is performed using a PixCell Iie microscope with an infrared diode laser (Life Technologies).
- 3.2.) Adjust the center joystick to the vertical position. Rotate and lock the 4× objective into position.
- 3.3.) Place a slide at the center of the microscope stage and manually adjust until the region of the hippocampus is in view. Use coarse and fine focus adjustments to bring the image into focus.
- 3.4.) Activate the vacuum chuck by pressing the vacuum button on the controller. Rotate and lock the 10× objective into place. Focus again with coarse and fine adjustments.
- 3.5.) Load the CapSure Macro LCM caps (Life Technologies) into the CapSure cassette module and position a cap at the load line position. Rotate the Cap Placement arm and position around the cap. Raise the Cap Placement arm to remove the CapSure cap from the cassette module.
- 3.6.) Rotate the Cap Placement arm over the slide and lower the arm down over the slide; this will place the cap onto the slide. Adjust the fine focus and move the joystick to check if cells are inside the black circle on the cap. All cell captures should be inside this black circle.
- 3.7.) Set the laser parameters. First, press the laser-enable button on the controller. The laser target spot will be visible as a pink dot in the field of view on the computer monitor.
- 3.8.) Set the laser spot size to small (7.5 μm) to focus the laser and adjust the power and duration to 65 – 75 mW, for 1.0-2.0 ms as needed for optimal cell capture.
- 3.9.) To test the laser, use the joystick to move the laser spot to an area where there are no cells to pick up. Fire the laser using the thumb switch. Use the joystick to move the laser spot away from the melted polymer spot and check to see if a visible dark ring surrounds a clear area where the laser was fired.
- 3.10) To capture cells, use the joystick to move the laser spot over the area of cells to capture. Fire the laser using the thumb switch. Move the joystick while firing the laser to capture a large area of cells. When the laser is fired it melts the thermoplastic polymer film on the cap to the surface of the target cells. The polymer absorbs the laser radiation, thus not altering the RNA, DNA or protein ensuring the integrity of the sample for future molecular applications.
- 3.11) After firing all cells of interest, lift the Cap Placement arm. The captured cells will separate from the tissue section and adhere to the CapSure cap. The remaining tissue and

missing cells will be visible in the field of view. The cells on the cap can also be viewed by placing the cap on an empty portion of the slide.

- 3.12) Lift the Cap Placement arm and rotate it to the Cap Unload Station. Lower the cap into the station and rotate the Cap Placement arm back into its previous position.
- 3.13) Slide the CapSure insertion tool along the platform and onto the cap. Lift the tool from the platform. The cap will remain attached. Lightly touch the cap to the CapSure clean pad to clean away unwanted tissue.
- 3.14) Place the cap into a 0.5 ml RNase-free tube filled with 100 ul of lysis buffer obtained from the RNAqueous Micro-Isolation Kit. Immediately vortex the cap for 30 seconds to lyse the cells. Incubate the samples at 42°C for 30 minutes and freeze at -80°C until RNA isolation.

LCM of single FJ+ neurons in the hippocampus.

- 3.15) To capture single cells, load CapSure HS LCM caps into the CapSure cassette module and position a cap at the load line position. Rotate the Cap Placement arm and position it around the cap. Raise the Cap Placement arm to remove the CapSure cap from the cassette module.
- 3.16) Set the laser spot size to small (7.5 μm) to focus the laser and adjust the laser power and duration to 65 -75 mW for 0.45 to 0.50 ms for optimal cell capture of single cells.
- 3.17) Rotate the Cap Placement arm over the slide and lower the arm down toward the slide. This will place the cap onto the slide. The Capsure HS cap surface sits elevated from the sample during LCM.
- 3.18) Use the joystick to move the laser over single FJ+ cells. All cells being captured must be inside the small black ring on the cap. Fire the laser using the thumb switch.
- 3.19) When all cells within the black circle area have been fired with the laser, lift the Cap Placement arm. The captured cells will separate from the tissue section and adhere to the CapSure HS cap.
- 3.20) Place another slide on the microscope stage and collect FJ+ cells until the HS cap is full. Place the cap into a 0.5 ml RNase-free tube filled with 40-50 ul of lysis buffer. Immediately vortex the cap for 30 seconds to lyse the cells. Incubate the samples at 42°C for 30 minutes and freeze at -80°C until RNA isolation.

4.) Total RNA isolation from LCM samples (this section will just be described in the video presentation)

- 4.1.) Total RNA from cells is isolated using the RNAqueous-Micro Kit (Ambion) following the manufacturer's protocols. All reagents and wash solutions are provided in this kit. Pre-wet

the micro filter cartridge assembly by adding 30 ul of lysis solution buffer to the filter. After 5 minutes, centrifuge the filter for 30 seconds at $10,000 \times g$.

- 4.2.) Add 3 ul of LCM additive to the sample lysate, pipette to mix, and centrifuge.
- 4.3.) Add 52 ul of 100% ETOH (1/2 volume) to the sample lysate; pipette to mix, and transfer the lysate to the filter column.
- 4.4) Centrifuge the filter column at $10,000 \times g$ for 1 min to bind the RNA to the column. If there are multiple caps for one sample, spin each lysate through the same column separately.
- 4.5) Complete the RNA isolation procedure as described in the RNAqueous Micro kit.
- 4.6) Perform DNase treatment and DNase inactivation of LCM RNA samples as described in the kit, transfer the RNA to a RNase-free tube and store at -80°C .

5.) RNA Quantification and Downstream applications

Assessment of RNA quality and quantification is performed with an Agilent Bioanalyzer using reagents from the Pico Kit (Agilent Technologies, Santa Clara, CA) following the manufacturer's protocol

- 5.1.) Total RNA is analyzed on an Agilent Bioanalyzer 2100 (Agilent Technologies). The software evaluates the concentration and the integrity of the RNA sample by assigning an RNA Integrity Number (RIN) to each sample ranging from 1 to 10. 1-5 indicates degraded RNA and 7-10 indicates good quality RNA.
- 5.2.) Non-degraded intact RNA will be used for Real-Time PCR; individual Taqman assays or RT2 profiler arrays using Sybr-green chemistry. The RNA can also be used for mRNA microarray analysis and microRNA microarray analysis.

REPRESENTATIVE RESULTS:

In our first LCM study, we were able to show that functionally distinct subregions (CA1 and CA3) of the rat hippocampus have gene expression profiles that reflect their vulnerability to TBI.⁸ Figure 1A-C shows laser capture of hippocampal pyramidal neurons (CA1-CA3), dentate gyrus neurons and SCN neurons before, and after LCM. Clean captures of pyramidal, granule or SCN neurons, respectively, are shown on the macro caps. Figure 2A-B illustrates dying, Fluoro-Jade positive pyramidal neurons and surviving, Fluoro-Jade negative pyramidal neurons before and after LCM. The clean capture of the dying or surviving neurons is shown by viewing the captured cells on the LCM caps.

After LCM, the total RNA can be used for a diverse array of molecular studies including quantitative real-time PCR analysis of gene expression using TaqMan or SYBR green

chemistries (Figure 3A), small focused PCR arrays with SYBR green probes (Figure 3B) or whole genome microarray studies (Figure 3C).

DISCUSSION:

This LCM technique has enabled us to make significant advances in understanding the molecular mechanisms of TBI.^{8,10,11,12,13} We were the first investigators to utilize this technique to demonstrate that distinct subregions of the rat hippocampus have gene expression profiles that correlate with their selective vulnerability to injury. **Our studies showed that we could quantify gene expression, by qPCR, from as few as 10 laser captured cells and that we could perform genome-wide microarray analysis from as few as 600 captured cells.** LCM would be a valuable tool in similar studies of other brain regions that are known to be implicated in diverse neurological and neuropsychiatric disorders. For example, studies using LCM have shed light into the pathogenesis of Parkinson's disease in dopamine neurons of the substantia nigra,¹⁴ and aided genome-wide profiling studies of the nucleus accumbens, which is involved in reward circuits implicated in substance abuse disorders.¹⁵

Neuronal heterogeneity is reflected at the genome level¹⁶ and may contribute to the lack of success of experimental treatments for TBI in clinical trials. Thus, our goal is to utilize this technique to investigate the critical elements influencing neuronal survival after TBI. Our recent genome-wide profiling study of dying and adjacent surviving hippocampal neurons after TBI suggested that a cellular rheostat that reflects the ratio of the expression levels of cell survival to cell death genes regulates cell fate after TBI. These ongoing studies will contribute to the design and development of pharmacotherapeutic strategies that can positively influence the cell survival rheostat. Furthermore, we currently use LCM to track and monitor the effects of potential therapeutic drug treatments in hippocampal neurons after TBI. Thus, our studies demonstrate that given careful RNase-free handling techniques and some modifications of existing protocols, it is possible to obtain high quality RNA samples from LCM for accurate quantitative gene expression analysis.

However there are a few pitfalls associated with LCM techniques. For example, collecting only neuronal cells without any microglia contamination can be almost impossible. In the pyramidal cell layer of the hippocampus it has been estimated that 95% of the cells are neurons with 90% of this population to be pyramidal cells and 10% interneurons, leaving a small percentage of glial and other cell types.^{8,17,18} In our studies we use Fluoro-jade a fluorescent dye that specifically labels only injured neurons in brain tissue. A different stain that is a marker for GFAP would be necessary to stain microglia¹⁹. Collecting only the Fluoro-jade stained single neuronal bodies ensures a more homogenous population of neurons. Another common problem that may require troubleshooting is that a circle is not defined when the laser is fired. If this occurs, it is necessary to check the laser settings and adjust the power and duration as necessary. Also it is important to ensure that the cap is placed flat on the tissue and seated correctly in the arm. Referring to the LCM technical manual or calling technical support will sometimes be necessary. Following LCM and RNA isolation, the quality of the RNA should always be assessed using a bioanalyzer prior to any gene expression analysis.

There are several types of laser capture instruments currently on the market including laser cutting (Life Technologies, Leica Microsystems) and laser catapulting (Zeiss) instruments. Our LCM system works well for capturing small numbers of cells. Other high-throughput and more automated systems may be more suitable for obtaining larger numbers of cells for genomic and particularly proteomic analysis. Indeed, LCM using these automated systems has great potential for analysis of protein expression in identified cells or enriched populations of cells.²⁰ Furthermore, LCM can facilitate gene expression analysis of immunolabeled cells,²¹ allowing us to investigate gene expression in defined cell types regardless of the complexity in the most heterogeneous tissues. Thus, LCM is an excellent tool for cutting-edge molecular studies of single or **enriched** populations of cells.

ACKNOWLEDGMENTS:

This work is funded by R01 NS052532 (to HLH), the Moody Foundation, and the Department of Anesthesiology. We thank Laurie Bolding, Christine Courteau Butler and Christy Perry for their editorial assistance, and Christy Perry for layout and excellent production of all figures, tables and artwork.

DISCLOSURES:

- The authors declare that they have no conflicts of interest.

REAGENTS AND EQUIPMENT:

Name of the reagent	Company	Catalogue number
PixCell Ite Laser Capture Microscope	Life Technologies (Arcturus)	
Agilent 2100 Bioanalyzer	Agilent Technologies	
RNAqueous Micro- Kit	Life Technologies (Ambion)	AM1931
Agilent 6000 Pico Kit	Agilent Technologies	5067-1513
Capsure LCM caps	Applied Biosystems	LCM0211/LCM0214
Fluoro-Jade	Histo-chem Inc.	#1FJ
Cresyl Violet Acetate	Sigma-aldrich	C5042-10G
Xylene (Histological grade)	Fisher Scientific	X3P-1GAL
Ethanol 200 proof (Histological grade)	UTMB pharmacy	
RNase free barrier pipette tips	Fisher Scientific	2123635, 212364, 212361, 21-236-2A

Arcturus Laser Capture Microdissection system (Life Technologies)
http://tools.invitrogen.com/content/sfs/brochures/cms_086283.pdf

Zeiss Laser Microdissection The PALM family
http://microscopy.zeiss.com/microscopy/en_de/downloads/brochure-downloads.html?catalog=laser_microdissection

Leica Microsystems
http://www.leica-microsystems.com/fileadmin/downloads/Leica%20DM6000%20B/Brochures/Leica_DM4000-6000_B-Brochure_en.pdf

FIGURE LEGENDS

Figure 1. Laser capture microdissection of swaths of defined cell populations from the rat brain. Frozen 10µm coronal sections of the rat brain are fixed in ethanol, stained with a nissl stain (cresyl violet) and prepared for LCM. Shown are pictures of the rat hippocampus and suprachiasmatic nucleus before and after LCM and the captured cells as visualized on the macro caps. A. Pyramidal neurons from the CA1-CA3 subfields of the rat hippocampus. B. Granule cells in the hippocampal dentate gyrus. C. Bilateral suprachiasmatic nuclei located on either side of the third ventricle and located above the optic chiasma.


Figure 2. Laser capture microdissection of single neurons from rat hippocampus. Shown are A. Degenerating, Fluoro-Jade-positive (stained) hippocampal CA3 neurons and B. Surviving, Fluoro-Jade-negative (unstained) CA3 neurons before and after LCM. Captured cells are visualized.


Figure 3. Gene expression analysis of total RNA isolated from laser captured neurons. A. Quantitative real time PCR data of circadian clock gene expression in SCN. B Heatmap of gene expression in dying and surviving hippocampal neurons showing expression of apoptosis-related genes using focused PCR arrays. C. Agilent whole-genome microarray analysis of gene expression in hippocampal CA1-CA3 neurons from rats that received sham injury, TBI or TBI plus a recombinant adeno-associated siRNA virus designed to knockdown expression of genes induced by brain injury (neuronal nitric oxide synthase and glutathione peroxidase-1).

Reference List


1. Squire, L. R., Stark, C. E., Clark, R. E. The medial temporal lobe. *Annual Review of Neuroscience* **27** 279-306 (2004).
2. Bast, T. Toward an integrative perspective on hippocampal function: from the rapid encoding of experience to adaptive behavior. *Rev Neurosci* **18** (3-4), 253-281 (2007).
3. Emmert-Buck, M. R., Bonner, R. F. *et al.* Laser capture microdissection. *Science* **274** (5289), 998-1001 (1996).
4. Bota, M., Dong, H. W., Swanson, L. W. From gene networks to brain networks. *Nat Neurosci* **6** (8), 795-799 (2003).
5. Schouten, J. W. Neuroprotection in traumatic brain injury: a complex struggle against the biology of nature. *Curr Opin Crit Care* **13** (2), 134-142 (2007).
6. Schmued, L. C., Albertson, C., Slikker, W., Jr. Fluoro-Jade: a novel fluorochrome for the sensitive and reliable histochemical localization of neuronal degeneration. *Brain Research* **751** (1), 37-46 (1997).
7. Ye, X., Carp, R. I., Schmued, L. C., Scallet, A. C. Fluoro-Jade and silver methods: application to the neuropathology of scrapie, a transmissible spongiform encephalopathy. *Brain Res Brain Res Protoc* **8** (2), 104-112 (2001).
8. Shimamura, M., Garcia, J. M., Prough, D. S., Hellmich, H. L. Laser capture microdissection and analysis of amplified antisense RNA from distinct cell populations of the young and aged rat brain: effect of traumatic brain injury on hippocampal gene expression. *Mol Brain Res* **17** (1), 47-61 (2004).
9. Alder, J., Fujioka, W., Lifshitz, J., Crockett, D. P., Thakker-Varia, S. Lateral fluid percussion: model of traumatic brain injury in mice. *J Vis Exp*(54), (2011).
10. Shimamura, M., Garcia, J. M. *et al.* Analysis of long-term gene expression in neurons of the hippocampal subfields following traumatic brain injury in rats. *Neuroscience* **131** (1), 87-97 (2005).
11. Shah, S. A., Prough, D. S., Garcia, J. M., DeWitt, D. S., Hellmich, H. L. Molecular correlates of age-specific responses to traumatic brain injury in mice. *Exp Gerontol* **2006 Nov**; **41**(11): **1201-5**. *Epub* **2006 Sep 15** (2006).
12. Hellmich, H. L., Garcia, J. M. *et al.* Traumatic brain injury and hemorrhagic hypotension suppress neuroprotective gene expression in injured hippocampal neurons. *Anesthesiology* **102** (4), 806-814 (2005).


13. Rojo, D. R., Prough, D. S. *et al.* Influence of stochastic gene expression on the cell survival rheostat after traumatic brain injury. *PLoS One* **6** (8), e2311-(2011).
14. Elstner, M., Morris, C. M. *et al.* Expression analysis of dopaminergic neurons in Parkinson's disease and aging links transcriptional dysregulation of energy metabolism to cell death. *Acta Neuropathol* **122** (1), 75-86 (2011).
15. Chen, H., Liu, Z. *et al.* Genome-Wide Gene Expression Profiling of Nucleus Accumbens Neurons Projecting to Ventral Pallidum Using both Microarray and Transcriptome Sequencing. *Front Neurosci* **5** 98-(2011).
16. Lein, E. S., Hawrylycz, M. J. *et al.* Genome-wide atlas of gene expression in the adult mouse brain. *Nature* **445** (7124), 168-176 (2007).
17. Amaral, D. G., Ishizuka, N., Claiborne, B. Neurons, numbers and the hippocampal network. *Prog Brain Res* **83** 1-11 (1990).
18. Zeng, Y. C., Bongrani, S. *et al.* Effect of long-term treatment with L-deprenyl on the age-dependent microanatomical changes in the rat hippocampus. *Mech Ageing Dev* **79** (2-3), 169-185 (1995).
19. Schmued, L. C., Hopkins, K. J. Fluoro-Jade: novel fluorochromes for detecting toxicant-induced neuronal degeneration. *Toxicol Pathol* **28** (1), 91-99 (2000).
20. Banks, R. E., Dunn, M. J. *et al.* The potential use of laser capture microdissection to selectively obtain distinct populations of cells for proteomic analysis-- preliminary findings. *Electrophoresis* **20** (4-5), 689-700 (1999).
21. Fend, F., Emmert-Buck, M. R. *et al.* Immuno-LCM: laser capture microdissection of immunostained frozen sections for mRNA analysis. *Am J Pathol* **154** (1), 61-66 (1999).





*Figure3
[Click here to download high resolution image](#)












Animated/Video Figure 04
[Click here to download high resolution image](#)




Animated/Video Figure 05
[Click here to download high resolution image](#)






















CAPTURED CELLS





Animated/Video Figure 13


[Click here to download high resolution image](#)








Animated/Video Figure 16
[Click here to download high resolution image](#)



Animated/Video Figure 17
[Click here to download high resolution image](#)



*Table of Reagents and Equipment
[Click here to download Table of Reagents/ Materials Used: Reagents and Equipment.xlsx](#)

Name of the reagent	Company	Catalogue number
PixCell Iie Laser Capture	Life Technologies	
Agilent 2100 Bioanalyzer	Agilent Technologies	
RNAqueous Micro- Kit	Life Technologies	AM1931
Agilent 6000 Pico Kit	Agilent Technologies	5067-1513
Capsure LCM caps	Applied Biosystems	LCM0211/LCM0214
Fluoro-Jade	Histo-chem Inc.	#1FJ
Cresyl Violet Acetate	Sigma-aldrich	C5042-10G
Xylene (Histological grade)	Fisher Scientific	X3P-1GAL
Ethanol 200 proof (Histological	UTMB pharmacy	
RNase free barrier pipette tips	Fisher Scientific	2123635, 212364, 212361, 21-236-

ARTICLE AND VIDEO LICENSE AGREEMENT

Title of Article:

Author(s):

Item 1 (check one box): The Author elects to have the Materials be made available (as described at <http://www.jove.com/publish>) via: ☒ Standard Access ☐ Open Access

Item 2 (check one box):

- ☒ The Author is NOT a United States government employee.
- ☐ The Author is a United States government employee and the Materials were prepared in the course of his or her duties as a United States government employee.
- ☐ The Author is a United States government employee but the Materials were NOT prepared in the course of his or her duties as a United States government employee.

ARTICLE AND VIDEO LICENSE AGREEMENT

1. **Defined Terms.** As used in this Article and Video License Agreement, the following terms shall have the following meanings: **"Agreement"** means this Article and Video License Agreement; **"Article"** means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein; **"Author"** means the author who is a signatory to this Agreement; **"Collective Work"** means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; **"CRC License"** means the Creative Commons Attribution-Non Commercial-No Derivs 3.0 Unported Agreement, the terms and conditions of which can be found at: <http://creativecommons.org/licenses/by-nc-nd/3.0/legalcode>; **"Derivative Work"** means a work based upon the Materials or upon the Materials and other pre-existing works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; **"Institution"** means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; **"JoVE"** means MyJoVE Corporation, a Massachusetts corporation and the publisher of *The Journal of Visualized Experiments*; **"Materials"** means the Article and / or the Video; **"Parties"** means the Author and JoVE; **"Video"** means any video(s) made by the Author, alone or in conjunction with any other parties, or by JoVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion of the Article, and in which the Author may or may not appear.

2. **Background.** The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JoVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.

3. **Grant of Rights in Article.** In consideration of JoVE agreeing to publish the Article, the Author hereby grants to JoVE, subject to **Sections 4 and 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the "Open Access" box has been checked in **Item 1** above, JoVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.

4. **Retention of Rights in Article.** Notwithstanding the exclusive license granted to JoVE in **Section 3** above, the

ARTICLE AND VIDEO LICENSE AGREEMENT

Author shall, with respect to the Article, retain the non-exclusive right to use all or part of the Article for the non-commercial purpose of giving lectures, presentations or teaching classes, and to post a copy of the Article on the Institution's website or the Author's personal website, in each case provided that a link to the Article on the JoVE website is provided and notice of JoVE's copyright in the Article is included. All non-copyright intellectual property rights in and to the Article, such as patent rights, shall remain with the Author.

5. Grant of Rights in Video – Standard Access. This **Section 5** applies if the "Standard Access" box has been checked in **Item 1** above or if no box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby acknowledges and agrees that, Subject to **Section 7** below, JoVE is and shall be the sole and exclusive owner of all rights of any nature, including, without limitation, all copyrights, in and to the Video. To the extent that, by law, the Author is deemed, now or at any time in the future, to have any rights of any nature in or to the Video, the Author hereby disclaims all such rights and transfers all such rights to JoVE.

6. Grant of Rights in Video – Open Access. This **Section 6** applies only if the "Open Access" box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby grants to JoVE, subject to **Section 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Video in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Video into other languages, create adaptations, summaries or extracts of the Video or other Derivative Works or Collective Works based on all or any portion of the Video and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. For any Video to which this Section 6 is applicable, JoVE and the Author hereby grant to the public all such rights in the Video as provided in, but subject to all limitations and requirements set forth in, the CRC License.

7. Government Employees. If the Author is a United States government employee and the Article was prepared in the course of his or her duties as a United States government employee, as indicated in **Item 2** above, and any of the licenses or grants granted by the Author hereunder exceed the scope of the 17 U.S.C. 403, then the rights granted hereunder shall be limited to the maximum rights permitted under such statute. In such case, all provisions contained herein that are not in conflict with such statute shall remain in full force and effect, and all provisions contained herein that do so conflict

shall be deemed to be amended so as to provide to JoVE the maximum rights permissible within such statute.

8. Likeness, Privacy, Personality. The Author hereby grants JoVE the right to use the Author's name, voice, likeness, picture, photograph, image, biography and performance in any way, commercial or otherwise, in connection with the Materials and the sale, promotion and distribution thereof. The Author hereby waives any and all rights he or she may have, relating to his or her appearance in the Video or otherwise relating to the Materials, under all applicable privacy, likeness, personality or similar laws.

9. Author Warranties. The Author represents and warrants that the Article is original, that it has not been published, that the copyright interest is owned by the Author (or, if more than one author is listed at the beginning of this Agreement, by such authors collectively) and has not been assigned, licensed, or otherwise transferred to any other party. The Author represents and warrants that the author(s) listed at the top of this Agreement are the only authors of the Materials. If more than one author is listed at the top of this Agreement and if any such author has not entered into a separate Article and Video License Agreement with JoVE relating to the Materials, the Author represents and warrants that the Author has been authorized by each of the other such authors to execute this Agreement on his or her behalf and to bind him or her with respect to the terms of this Agreement as if each of them had been a party hereto as an Author. The Author warrants that the use, reproduction, distribution, public or private performance or display, and/or modification of all or any portion of the Materials does not and will not violate, infringe and/or misappropriate the patent, trademark, intellectual property or other rights of any third party. The Author represents and warrants that it has and will continue to comply with all government, institutional and other regulations, including, without limitation all institutional, laboratory, hospital, ethical, human and animal treatment, privacy, and all other rules, regulations, laws, procedures or guidelines, applicable to the Materials, and that all research involving human and animal subjects has been approved by the Author's relevant institutional review board.

10. JoVE Discretion. If the Author requests the assistance of JoVE in producing the Video in the Author's facility, the Author shall ensure that the presence of JoVE employees, agents or independent contractors is in accordance with the relevant regulations of the Author's institution. If more than one author is listed at the beginning of this Agreement, JoVE may, in its sole discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JoVE reserves the right, in its absolute and sole discretion and without giving any reason therefore, to accept or decline any work submitted to JoVE. JoVE and its employees, agents and independent contractors shall have full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JoVE has sole discretion as to the method of making and publishing the Materials, including,

ARTICLE AND VIDEO LICENSE AGREEMENT

without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

11. **Indemnification.** The Author agrees to indemnify JoVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JoVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JoVE, making of videos by JoVE, or publication in JoVE or elsewhere by JoVE. The Author shall be responsible for, and shall hold JoVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to the making of a video by JoVE its employees, agents or independent contractors. All sterilization, cleanliness or decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's expense. All indemnifications provided herein shall include JoVE's attorney's fees and costs related to said losses or

damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JoVE, its employees, agents or independent contractors.

12. **Fees.** To cover the cost incurred for publication, JoVE must receive payment before production and publication the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JoVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.

13. **Transfer, Governing Law.** This Agreement may be assigned by JoVE and shall inure to the benefits of any of JoVE's successors and assignees. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to be one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement required per submission.

AUTHOR:

Name:

Helen L. Hellmich

Department:

Department of Anesthesiology

Institution:

University of Texas Medical Branch

Article Title:

Laser capture microdissection of pure populations of neurons or single neurons for gene expression analysis after traumatic brain injury

Signature:



Date:

August 13, 2012

Please submit a signed and dated copy of this license by one of the following three methods:

- 1) Upload a scanned copy as a PDF to the JoVE submission site upon manuscript submission (preferred);
- 2) Fax the document to +1.866.381.2236; or
- 3) Mail the document to JoVE / Attn: JoVE Editorial / 17 Sellers St / Cambridge, MA 02139

For questions, please email editorial@jove.com or call +1.617.945.9051.

MS # (internal use):

Dr. Angela Messmer-Blust, Ph.D
Science Editor
Journal of Visualized Experiments

15 October, 2012

MS: JoVE50308R1

“Laser capture microdissection of pure populations of neurons or single neurons for gene expression analysis after traumatic brain injury”

Dear Dr. Messmer-Blust,

Thank you for the opportunity to submit our revised manuscript to *Journal of Visualized Experiments*. We appreciate the insightful comments and suggestions of the reviewer and have revised our manuscript accordingly. Our response to the reviewer's comments is below.

Reviewer:

Major Concerns:

1. "...there should be some discussion of how one would confirm the purity of the sample for neurons with minimal microglia contamination..."

Response: We agree that when collecting a swath of neurons in the CA1-CA3 region of the hippocampus, it is impossible not to obtain some microglial contamination. However, greater than 95% of the cells in the pyramidal cell layer of the CA1 region are neurons [1]. Thus, microglia contamination should be minimal. Nevertheless, to be sure, we can obtain an assessment of glial contamination in the samples using qPCR to determine the expression levels of a glial marker such as GFAP. In our experiments we use Fluoro-Jade staining, which stains only injured neurons, allowing the identification and subsequent laser capture of single neuronal bodies with minimal to no microglia contamination. Uninjured neurons are identified by staining with Cresyl Violet, a Nissl stain that is compatible with preservation of RNA integrity and is also neuron specific.[2] We have revised the discussion to include an explanation of how to confirm the purity of a sample.

2. "...there is limited discussion of the pitfalls of trying to do this procedure, hints for how to troubleshoot, how to confirm quality, the minimal number of captured cells, any difficulties that may arise with using "injury" tissue.

Response: We agree that it is important to discuss the problems and pitfalls of LCM procedures. In particular, collecting single neurons is tedious and requires many subjects; cells from multiple animals need to be pooled to accumulate enough tissue for one sample. Furthermore to ensure good quality RNA, LCM sessions are limited to 30 min to 1 hour per, cap and the quality of the RNA is measured using an Agilent Bioanalyzer (Agilent Technologies, Santa Clara CA), discussed in Sections 5.1 and 5.2.

Minor Concerns:

1. Section 2. - The length of time the tissue can be held at -80°C needs to be mentioned.

Response: Tissue storage and timing is, indeed, an important issue. In our experience, brain tissue can be stored at -80°C for up to one month prior to sectioning. Once the brain is thawed, mounted in OCT and sliced, we promptly (within 24 h) stain the sections and use them for LCM within 30 minutes to one hour of staining. This ensures good quality RNA. We found that brain tissue that has been sectioned, then re-frozen and used for LCM at a later time yields lower quality RNA than whole brain tissue that has been frozen at -80°C for months prior to sectioning-- a technique which can still yield high quality RNA. We have revised and addressed in section 2.1 the comment pertaining to the length of time tissue can be stored at -80°C prior to sectioning and LCM.

2.” ... having one technician cutting the tissue staining and then giving to a second person doing the capture significantly decreased the number of samples required for microarray analysis.”

Response: We agree that it would be useful, when possible, to have one technician sectioning and staining while another technician performs LCM, as it would reduce the time required to process all of the sections. However, this protocol makes it possible for the procedure to be carried out by a single technician. Not all investigators have the number of research personnel needed for complex procedures such as these and in our laboratory one person is routinely able to obtain sufficient numbers of neurons for gene expression analysis using the methods described.

3.” ... it would be of benefit to mention the timing aspect of tissue collection...mention of how long the collected tissue can be stored prior to RNA isolation would be helpful.

Response: Collected tissue in lysis buffer can be stored at -80C for a week, but RNA should be isolated in a timely manner to insure highest quality. See Section 2.1.

4. Was the cryostat blade and the area on which the brain slice will rest after cutting washed with RNase free water? Or some other treatment? Is the 95% alcohol made from 100% with RNase free water?

Response: The cryostat is cleaned with ETOH and RNase-zap before sectioning. Disposable low profile blades are changed in between each tissue sample. All staining solutions are prepared with RNase free water.

5. If the authors have experience with other RNA isolation kits, such as a lipid-rich tissue kit, or others it would be helpful to mention this somewhere for what worked and what did not.

Response: The kit used here, RNAqueous-Micro Kit (Ambion), is specifically designed to isolate RNA from small populations of cells (1 – 10). At the time of our development of LCM techniques, it was the only kit available for this purpose and was particularly suited to LCM. There are other kits on the market, but we do not have any experience with them; this kit has worked quite well.

6. What is the minimum number of cells required for analysis? Thus, what is the approximate yield needed for each of the types of analysis.

Response: The minimum number of cells required for analysis depends on the downstream application. For microarray analysis, we have obtained excellent data from as few as 600 cells; for qPCR we have performed gene expression analysis from as few as 10 cells. In our first gene expression studies of neurons obtained by LCM, we were able to linearly amplify total RNA from as few as 5 cells, (3 captures from 10µm sections is estimated to be equivalent to one cell since the neurons are approximately 30 µm in size) and obtain excellent results using ribonuclease protection assay.[3] The discussion has been revised accordingly.

Regards,

Helen Lee Hellmich, Ph.D
Associate Professor
Department of Anesthesiology
University of Texas Medical Branch

Reference List

1. Zeng YC, Bongrani S, Bronzetti E, Cadel S, Ricci A, Valsecchi B, Amenta F (1995) Effect of long-term treatment with L-deprenyl on the age-dependent microanatomical changes in the rat hippocampus. *Mech Ageing Dev* 79: 169-185.
2. Shimamura M, Garcia JM, Prough DS, DeWitt DS, Uchida T, Shah SA, Avila MA, Hellmich HL (2005) Analysis of long-term gene expression in neurons of the hippocampal subfields following traumatic brain injury in rats. *Neuroscience* 131: 87-97.
3. Hellmich HL, Garcia JM, Shimamura M, Shah SA, Avila MA, Uchida T, Parsley MA, Capra BA, Eidson KA, Kennedy DR, Winston JH, DeWitt DS, Prough DS (2005) Traumatic brain injury and hemorrhagic hypotension suppress neuroprotective gene expression in injured hippocampal neurons. *Anesthesiology* 102: 806-814.