# **Journal of Visualized Experiments**

# Isolation of viable multicellular glands from tissue of the carnivorous plant, Nepenthes --Manuscript Draft--

Manuscript Number:	JoVE50993R2
Full Title:	Isolation of viable multicellular glands from tissue of the carnivorous plant, Nepenthes
Article Type:	Methods Article - JoVE Produced Video
Corresponding Author:	Axel Mithöfer Max Planck Institute Jena, Thuringia GERMANY
Corresponding Author Secondary Information:	
Corresponding Author E-Mail:	amithoefer@ice.mpg.de
Corresponding Author's Institution:	Max Planck Institute
Corresponding Author's Secondary Institution:	
First Author:	Sandy Rottloff
First Author Secondary Information:	
Other Authors:	Sandy Rottloff
	Ute Müller
	Roland Kilper
Order of Authors Secondary Information:	
Abstract:	Many plants possess specialized structures that are involved in the production and secretion of specific low molecular weight compounds and proteins. These structures are almost always localized on plant surfaces. Among them are nectaries or glandular trichomes. The secreted compounds are often employed in interactions with the biotic environment, for example as attractants for pollinators or deterrents against herbivores. Glands that are unique in several aspects can be found in carnivorous plants. In so-called pitcher plants of the genus Nepenthes, bi-functional glands inside the pitfall-trap on the one hand secrete the digestive fluid, including all enzymes necessary for prey digestion, and on the other hand take-up the released nutrients. Thus, these glands represent an ideal, specialized tissue predestinated to study the underlying molecular, biochemical, and physiological mechanisms of protein secretion and nutrient uptake in plants. Moreover, generally the biosynthesis of secondary compounds produced by many plants equipped with glandular structures could be investigated directly in glands.  In order to work on such specialized structures, they need to be isolated efficiently, fast, metabolically active, and without contamination with other tissues. Therefore, a mechanical micropreparation technique was developed and applied for studies on Nepenthes digestion fluid. Here, a protocol is presented that was used to successfully prepare single bi-functional glands from Nepenthes traps, based on a mechanized microsampling platform. The glands could be isolated and directly used further for gene expression analysis by PCR techniques after preparation of RNA.
Author Comments:	
Additional Information:	

Cover Letter

Dear Editor,

Enclosed we submit our manuscript "Isolation of viable multicellular glands from tissue of the carnivorous plant, *Nepenthes*" by Sandy Rottloff et al. to the Journal of Visualized Experiments.

In this paper we describe an efficient and rapid method to prepare single glands from plant tissues using microdissection. As model plant we chose the carnivorous species Nepenthes because here these glands are involved in both, enzyme secretion as well as nutrient uptake, and represent an ideal tissue to study plant cell biology. Our approach allows a fast mechanical preparation of almost intact glands which can be directly used for subsequent experiments, and serves as an alternative to laser capture dissection. This might launch new strategies to isolate and study glandular structures in plants in general.

All authors have seen and approved the paper, the manuscript has not been, and will not be published elsewhere. SR and UM did almost all of the experiments; AM and RK designed the experiments. SR and AM wrote the manuscript.

This preparation of this paper has been assisted by E. Sheeley.

We thank you for your kind attention.

Yours sincerely,

Axel Mithöfer

#### Suggested referees:

Günther Heubl; University of Munich; heubl@lrz.uni-muenchen.de

Veronika. Mayer; University of Vienna; veronika. mayer @univie.ac. at

Dirk Hoelscher; Max Planck Institute CE; hoelscher@ice.mpg.de

Irene Lichtscheidl; University of Vienna; irene.lichtscheidl@univie.ac.at

Harald Meimberg; Technical University Munich; meimberg@wzw.tum.de

Angelos Kanellis; University Thessaloniki; kanellis@pharm.auth.gr

## Isolation of viable multicellular glands from tissue of the carnivorous plant, Nepenthes

#### **Authors:**

Sandy Rottloff, Axel Mithöfer, Ute Müller, Roland Kilper

#### Authors: institution(s)/affiliation(s) for each author:

Rottloff, Sandy Laboratoire Agronomie et Environnement Université de Lorraine

Mithöfer, Axel
Department of Bioorganic Chemistry
Max Planck Institute for Chemical Ecology

Müller, Ute aura optik

Kilper, Roland aura optik

Corresponding author: Dr. Axel Mithöfer.

Telephone: +49 (0) 3641 57 1263

Fax: +49 (0)3641 57 1202

# **Keywords:**

Plant science, micropreparation, mechanical dissection, glands, carnivory, Nepenthes, PCR.

# **Short Abstract:**

Plants glands are specialized structures responsible for the biosynthesis and secretion of many compounds involved in interactions with the biotic environment. To enable studies on their molecular and biochemical features, a mechanical micropreparation technique was established in order to isolate single metabolically active glands, here from the carnivorous plant *Nepenthes*.

#### Long Abstract:

Many plants possess specialized structures that are involved in the production and secretion of specific low molecular weight compounds and proteins. These structures are almost always localized on plant surfaces. Among them are nectaries or glandular trichomes. The secreted compounds are often employed in interactions with the biotic environment, for example as attractants for pollinators or deterrents against herbivores.

Glands that are unique in several aspects can be found in carnivorous plants. In so-called pitcher plants of the genus *Nepenthes*, bi-functional glands inside the pitfall-trap on the one hand secrete the digestive fluid, including all enzymes necessary for prey digestion, and on the other hand take-up the released nutrients. Thus, these glands represent an ideal, specialized tissue predestinated to study the underlying molecular, biochemical, and

physiological mechanisms of protein secretion and nutrient uptake in plants. Moreover, generally the biosynthesis of secondary compounds produced by many plants equipped with glandular structures could be investigated directly in glands.

In order to work on such specialized structures, they need to be isolated efficiently, fast, metabolically active, and without contamination with other tissues. Therefore, a mechanical micropreparation technique was developed and applied for studies on *Nepenthes* digestion fluid. Here, a protocol is presented that was used to successfully prepare single bi-functional glands from *Nepenthes* traps, based on a mechanized microsampling platform. The glands could be isolated and directly used further for gene expression analysis by PCR techniques after preparation of RNA.

#### Introduction:

Plants produce a wide huge variety of compounds with different industrial and pharmaceutical applications, ranging from small molecular weight molecules to polymers. Such compounds are frequently produced and secreted by highly specialized structures, many of which are localized at the surface of the plants such as glandular trichomes or nectaries, or internally such as idioblats or latex and resin ducts, respectively. However, relatively little is known about the biology of these specialized secretory structures, although all types of those glands display many features that are indicative of active metabolism. In particular, glands of carnivorous plants such as *Nepenthes* can be seen as model systems in plant cell biology<sup>1</sup>.

Carnivorous plants are the object of several studies since 1875 when Charles Darwin's 'Insectivorous Plants' was published<sup>2</sup>. However, only in the last few years research on the molecular level has been done and still our knowledge is limited. For example, the protein composition of the digestion fluid in the pitcher plants of the genus *Nepenthes* is still not completely known. Only recently several hydrolytic proteins have been identified and a few of the corresponding genes as well<sup>3</sup>. In the pitfall-traps of *Nepenthes* (Figure 1), multicellular glands are localized at the inner bottom of the pitchers. These glands do both production and secretion of the digestive enzymes into the pitcher fluid and absorption of dissolved nutrients from the pitcher fluid<sup>4</sup>. Hence, these bi-functional glands represent a unique micro-tissue with a central position and function in the carnivory of *Nepenthes* that is worth to be studied in detail.

# [Place Figure 1 here]

In plant sciences, many experimental approaches rely on 'bulk material' because no defined tissue is investigated but entire organs. As a consequence, their homogenization and analysis can only provide results that are diluted and proportionately with respect to individual tissues and cells<sup>5</sup>. To solve this problem, micropreparation techniques, for example the so-called laser capture microdissection (LCM), were developed and successfully used to harvest specific plant tissues and analyze their individual contents<sup>6</sup>. However, LCM is often limited if unstable and/or fast degrading cell components are targeted such as RNA. In such a case, the need of previous tissue fixation that is time consuming and often accompanied by degradation of the biomolecules, is of disadvantage. Moreover, due to the high water content of many plant tissues and the strength of cell walls, LCM often failed to prepare straight from fresh tissue<sup>7</sup>. Any approach to investigate cellular components that

are present in secretory glands from genes, mRNA, and proteins to secondary compounds and their biosynthesis needs alternative techniques for the particular preparation. Thus, specific techniques to isolate single cells or even multicellular glands that remain intact after preparation are necessary.

Here a method is presented, employing a mechanized microsampling platform that features a direct microscopic visualization. This technique is both rapid and efficient in the preparation of individual glands from the pitchers of *Nepenthes* species. These multicellular glands are bi-functional, *i.e.* they are involved in secretion as well as in uptake processes, they are sessile, and not exposed cell layers are impregnated to form an endodermis. Entire glands can be isolated without neighboring tissue and deposited in a well-defined manner into a PCR tube. This technique is applied to study specific genes as well as gene expression in the glands' tissues.

#### **Protocol Text:**

# 1. Plant material preparation

- 1.1 All *Nepenthes* species were grown either in the greenhouses of the MPI for Chemical Ecology or the Botanical Garden in Jena, Germany.
- 1.2 For the preparation of secreting tissue harvest a pitcher and remove the upper one third including lid, peristome, and slippery zone. The pitcher can be used either freshly or frozen in liquid nitrogen and stored at -80 °C.

# 2. Microdissection of pitcher tissue

- 2.1. The microsampling platform consists of a stereomicroscope and a micromanipulator, which is directly attached to the microscope body and is operated by a 3D joystick. Here, an early prototype of aureka was used as documented (Figure 2A).
- 2.2 For the dissection of secretory glands from *Nepenthes* pitchers use a microforceps and attach it to the manipulator (Figure 2A, B).
- 2.3 For RNA extraction, treat all instruments (razor blades, tweezers, micromanipulator, microforceps, object slide) with RNaseZAP.
- 2.4 Prepare the collection tubes by adding 9  $\mu$ l nuclease-free water and 1  $\mu$ l RNase inhibitor in a PCR tube for RNA extraction. For direct PCR from tissue samples use only nuclease-free water (10  $\mu$ l).
- 2.4.1 Take the fresh or frozen pitcher and remove a 1 to 2 cm<sup>2</sup> piece of tissue containing the secretory glands by using a razor blade and tweezers. Wash the tissue with water; add 0.1% diethylpyrocarbonate (DEPC) if RNA will be extracted.
- 2.4.2 Fix the plant material with double-faced adhesive tape on an object slide (the gland covered surface should face upwards, Figure 2C, D) and secure the slide with adhesive tape on a universal stage, which should be applied with a moveable table.

- 2.4.3 For extracting nucleic acids the preparation should be done on ice. Therefore take a flat Styrofoam box, place the universal stage containing the sample inside and fill the box with ice. Afterwards place the box under the microscope.
- 2.4.4 For preparation the positioning of the glands is very important. The secretory glands are partial covered by a protective hood (Figure 3) and only accessible by the forceps from one side. Furthermore the angle of the sample to the forceps has to be adjusted by tilting the table of the universal stage approximately 135°, depending on the *Nepenthes* species.
- 2.4.5 The microdissection is done under visual control. For gland isolation, use a non-UV filter, to limit the degradation of DNA and RNA.
- 2.4.6 Move micromanipulator and forceps via the control panel to the selected gland. When reached the gland, place each arm of the forceps on either side of the gland. Remove the gland from the remaining tissue by closing the forceps (Figure 4A).
- 2.4.7 Direct the gland attached to the forceps away from the object and to the collection tube. Bring the very tips of forceps in contact with the solution in the tube (Figure 4B). The release of the gland by opening the forceps has to be well observed, because the glands tend to stick on the tool tip. The tissue can be stored at -80 °C until PCR or RNA extraction.

# 3. Direct PCR amplification from Nepenthes tissue

- 3.1 To prove the efficiency of the preparation method, a PCR approach directly from the gland tissue without a DNA extraction step can be performed.
- 3.2 For the PCR amplification use gene specific primers; here encoding for the sequence of a thaumatin-like protein (TLP) (forward: 5′-CAATGAGCCAATTCATAAAATTCATTG-3′, reverse: 5′-CAGTTATACTTTAAGGGCAAAACACAAC-3′).
- 3.2.1 Collect 1 to 5 glands in a PCR tube containing 10  $\mu$ l of nuclease-free water. For a 50  $\mu$ l reaction add 5  $\mu$ l of 10x Taq buffer, 1  $\mu$ l 10 mM deoxynucleoside triphosphates (dNTPs), 3  $\mu$ l 25 mM MgCl<sub>2</sub>, 0.4  $\mu$ l (1.25 U) Taq DNA polymerase, 1  $\mu$ l (100 pmol) of each primer, and 28.6  $\mu$ l of nuclease-free water.
- 3.2.2 Carry out PCR reactions in a thermocycler with the following program: initial denaturation 94 °C for 3 min; 35 cycles of denaturation (94 °C for 30 s), annealing (59 °C for 30 s), and extension (72 °C for 1 min); and a final extension step of 10 min at 72 °C.
- 3.2.3 Prove the amplification of the PCR product by using agarose gel electrophoresis. Therefor use 1.5% (w/v) agarose and dissolve it in 1x TAE (40 mM Tris base, 20 mM glacial acetic acid, 1 mM EDTA, pH 8) buffer. Heat the solution in a microwave until the agarose is dissolved. To 100 ml agarose solution, add 1  $\mu$ l ethidium bromide (10  $\mu$ g/ $\mu$ l) for downstream visualization. Run the gel in 1x TAE buffer at 300 V, 115 mA for approximately 20 min in an electrophoresis chamber. The correct size of the TLP amplicon must be at 690 bp.

# 4. RNA extraction from small tissue samples and reverse transcription

- 4.1 Perform the following steps on ice. For the extraction of RNA two extraction methods were used, a Micro Kit for the extraction of total RNA and a mRNA isolation kit applied with magnetic beads.
- 4.2 Collect at least 30 secretory glands by micropreparation in a PCR tube containing 1  $\mu$ l RNase inhibitor and 9  $\mu$ l of nuclease-free water. Afterwards lyse the tissue in a thermocycler for 10 min at 75 °C.
- 4.3 Extraction of total RNA.
- 4.3.1 After preparation add 90  $\mu$ l of lysis solution to the sample and incubate for 10 min at 75 °C. Afterwards load the reaction mix on the column (provided with the kit) and proceed with the washing steps as described in the manual of the manufacturer.
- 4.3.2 Elute the RNA two times with 10  $\mu$ l Elution Solution respectively and continue with the DNase I treatment.
- 4.3.3 To 20  $\mu$ l extracted RNA add 2  $\mu$ l of 10x DNase I buffer and 1  $\mu$ l DNase I (provided with RNA extraction kit) and incubate 1 h at 37 °C. Stop the reaction by adding 2  $\mu$ l of DNase Inactivation Reagent and incubate for 2 min at room temperature. After a centrifugation step for 1.5 min at maximum speed, the RNA can be transferred to a fresh RNase-free tube. Determine the RNA concentration of each sample using a spectrophotometer. The samples should be stored at -80 °C.
- 4.4 Extraction of mRNA using magnetic beads.
- 4.4.1 Adjust the sample volume to 100  $\mu$ l with lysis/binding buffer (supplied in the kit) and incubate for 10 min at 75 °C. Add 20  $\mu$ l of magnetic beads supplied with an Oligo (dT)<sub>25</sub> overhang and mix the reaction for 5 min at room temperature. The mRNA binds by the poly(A) tail to the Oligo (dT)<sub>25</sub> overhang of the beads.
- 4.4.2 Perform the washing steps of the mRNA bound to the magnetic beads by following the manufacturer's protocol. Solutions can be removed by using a magnet. Resuspend the beads in 20  $\mu$ l 10 mM Tris-HCl buffer (pH 7.5).
- 4.4.3 Remove contaminations of genomic DNA by treating the sample with 2 U (1  $\mu$ l) of DNase (2  $\mu$ l 10x DNase buffer and 17  $\mu$ l nuclease-free water) for 1 h at 37 °C. Afterwards the enzyme solution is removed by the magnet and the beads are resuspended in 20  $\mu$ l of 10 mM Tris-HCl buffer. It is not possible to measure the RNA concentration since it is attached to the beads. The samples should be stored at -80 °C.
- 4.5 For the reverse transcription use a cDNA Synthesis kit. Reactions were performed by following the manufacturer's protocol.
- 4.5.1 For the reverse transcription of total RNA add 1  $\mu$ l (100 pmol) Oligo (dT)<sub>20</sub> primer or 1  $\mu$ l (0.4  $\mu$ g) Random Hexamers, 1  $\mu$ l RNase inhibitor, and 1  $\mu$ l 10 mM dNTP Mix to 1  $\mu$ g of

RNA and adjust the volume to 14  $\mu$ l with nuclease-free water. Incubate at 65 °C for 5 min and afterwards place the samples on ice. Add 1  $\mu$ l 0.1 M DTT, 1  $\mu$ l reverse transcriptase, and 4  $\mu$ l of 5x RT buffer. The reactions are incubated for 55 min at 50 °C, followed by an inactivation step at 70 °C for 10 min. After reverse transcription samples can be stored at -20 °C.

- 4.5.2 For mRNA bound to magnetic beads, replace the Tris-HCl buffer with the first reaction mix (1  $\mu$ l RNase inhibitor, 1  $\mu$ l 10 mM dNTP Mix, and 12  $\mu$ l nuclease-free water) and incubate for 5 min at 65 °C. Because of the Oligo (dT)<sub>25</sub> overhang of the beads it is not necessary to add an Oligo (dT)<sub>20</sub> primer or Random Hexamers to the samples. Afterwards follow the protocol above (4.5.1).
- 4.5.3 Measure the cDNA concentration with a spectrophotometer. Store the samples at 20 °C.

#### 5. Real-time PCR

- 5.1 For the analyses with Real-time PCR use total RNA extracted with the Micro Kit and a SYBR Green QPCR Master Mix. To ensure an equal amount of transcripts, adjust the cDNA concentration of each sample to the same value after the reverse transcription.
- 5.2 In addition to the samples insert control reactions. Use a non-template control to screen for contamination in the reagents and to check the quality of your primers. To verify the absence of genomic DNA run a non-RT control. Both reactions should be negative (no fluorescence signal) in the Real-time PCR.
- 5.3 To guaranty an equal transcript level between the different samples use at least one reference gene. Here, designed primer pairs encoding for an actin (forward: 5′-CTCTTAACCCCAAAGCAAACAGG-3′, reverse: 5′-GTGAGAGAACAGCCTGGATG-3′) and an 18SrRNA (forward: 5′-CTTGATTCTATGGGTGGTGGTG-3′, reverse: GTTAGCAGGCTGAGGTCTC-3′) were used. As a target gene the TLP (forward: 5′- GTGGGCAACAATTGGACCCAG-3′, reverse: 5′-CATTCGATCCGGAAAGGCTAC-3′) was chosen. The size of each amplified fragment was about 100 bp.
- 5.4 For a 25  $\mu$ l reaction add 12.5  $\mu$ l 2x SYBR Green QPCR Master Mix, cDNA ( $\geq$  20 ng), 400 nM of each primer, and 30 nM of the passive reference dye ROX (provided with kit). Adjust the sample volume to 25  $\mu$ l with nuclease-free water.
- 5.5 Reactions were performed in a QPCR thermocycler with the following program: initial denaturation at 95 °C for 10 min to active the polymerase; 40 cycles of denaturation (95 °C for 30 s), annealing (61 °C for 60 s), and extension (72 °C for 1 min); followed by the dissociation curve (95 °C for 60 s, 55 °C for 30 s, 95 °C for 30 s). The dissociation curve is necessary to determine the absence of primer dimers, which can be detected by an additional peak.
- 5.6 Determine the efficiency (E) of each PCR reaction by using the method of Liu and Saint<sup>8</sup>. The efficiency can be deflected from the linear part of the amplification curve. Here the

relative amplicon number (Rn, A) of one cycle point (CPA) is linked to the amount of amplicons (Rn, B) of another cycle (CPB).

$$E = \left(\frac{Rn, A}{Rn, B}\right)^{\frac{1}{CPA - CPB}} - 1.$$

5.7 Analyze the Real-time PCR results with the Relative Expression Software Tool (REST 2009, <a href="http://www.gene-quantification.de/rest-2009.html">http://www.gene-quantification.de/rest-2009.html</a>) or a similar software. Use the program to calculate the expression ratio of a gene in a treated sample compared to a control.

# **Representative Results:**

For the preparation of single glands from *Nepenthes* pitcher tissue, a mechanized microdissection technique was developed and applied. The micropreparation platform needs to be equipped with microforceps, attached to the micromanipulator, which in turn are operated by a control panel. Individual glands, visible under the microscope, could be targeted with the forceps that grasped the gland tissue and removed it mechanically (Figure 4A).

The glandular tissue was transferred directly into a PCR tube for subsequent analysis (Figure 4B). The success of gland removal could be observed because it uncovered the plant's vascular system below the glands (Figure 4A).

In order to further prove the quality and efficiency of the preparation method and to make sure that this is a useful method, a PCR approach was performed. As a candidate gene a thaumatin-like protein (TLP) from Nepenthes gracilis (NCBI Acc. No DQ352144) was selected and amplified successfully (Figure 4C). However, not only single genes could be detected but, in addition, a direct isolation of both total- and mRNA was effectively performed and used for cDNA synthesis. By using TLP again, the amplification of that cDNA by PCR was demonstrated (Figure 4D). Furthermore, the cDNA was introduced in various Real-time PCR experiments. Therefore glands and fragments of the epidermal tissue, which is surrounding the glands, were isolated by the micropreparation technique from closed and open pitchers of N. alata and N. mirabilis. In a Real-time PCR reaction the TLP expression was analyzed. Using this method it was shown that the TLP transcript level is higher in the secretory glands compared to the epidermal tissue. After prey trapping (simulated by feeding the plant with Drosophila melanogaster, fruit flies) the epidermal tissue exhibited a higher TLP expression (Figure 4D-1). In addition, the inducibility of gene expression by the presence of prey was analyzed. Here, the TLP transcript level of glands and epidermal tissue in a control (non-fed) and treated (addition of fruit flies) pitcher of N. mirabilis was compared (Figure 4D-2). The upregulation of TLP expression in both tissue types indicates that the production of pitcher fluid proteins can indeed be induced by prey.

#### **Tables and Figures:**

**Figure 1:** *Nepenthes alata* pitcher. The lower part of the pitcher contains the digestive fluid and glands on the inside.

**Figure 2:** Microsampling platform. (A) aureka platform consisting of a stereomicroscope and a micromanipulator, equipped with microforceps. (B) Detailed view of the microforceps; in

the yellow oval, the universal stage with a moveable table to present the sample can been seen. (C) The same as in (B) with an object slide presenting the plant tissue. (D) Inner tissue of a *Nepenthes* pitcher spotted with glands. For photography, an UV-filter with an excitation wavelength of 365 nm was used.

**Figure 3:** Individual glands of different *Nepenthes* species, covered by a protective hood. Note the strong autofluorescence of the cell walls in the glands of *N. thorelii* (1), *N. superba* (2), *N. mirabilis* (7), and *N. superpum* (9), which is not present in the glands of the other species (3, 4, 5, 6, and 8).

Figure 4: Micropreparation of individual Nepenthes glands. (A) Forceps taking gland tissue, thereby uncovering the vascular system below the gland. (B) Direct transfer of an isolated gland into a PCR tube. For photo-documentation the usage of an UV-filter was necessary (excitation at 365 nm). (C) Individual gene-expression analysis. Here, five glands were combined and prepared by lysis for 20 min at 70 °C for mRNA extraction and reverse transcription of the mRNA. Genomic DNA was digested with DNase. The cDNA was used for the subsequent PCR-reaction to amplify the gene sequence of interest, indicated by the arrow (lane 1; M: DNA-ladder). (D) Analyses of TLP gene expression by Real-time PCR. Secretory glands and fragments of the glands surrounding epidermal tissue of closed and open pitchers (treated with Drosophila melanogaster) were microdissected. Total RNA was extracted from all samples. After digestion of genomic DNA and reverse transcription, the cDNA was used in various Real-time PCR experiments. (1) The TLP expression of glands (control) and epidermal tissue (treatment) of closed pitchers from N. alata (a) and N. mirabilis (b) was compared and showed a higher TLP transcript level in the secretory glands. After feeding D. melanogaster to an open N. mirabilis (c) pitcher the TLP transcript level was 3-fold higher in the epidermal tissue compared to the glands. (2) Analysis of TLP expression of two open N. mirabilis pitchers after feeding with D. melanogaster. The level of TLP transcripts in glands (d) and epidermal tissue (e) of the control (non-fed) and the treated (fed with D. melanogaster) pitcher was compared. The TLP-mRNA level in both tissue types was higher after feeding with fruit flies, which indicates the inducibility of the gene by prey; (expression ratio R(E) (R(E) < 1: downregulation, R(E) > 1: upregulation); \*\*\*, indicates a significant difference between groups compared, P(H1) < 0.001).

# **Discussion:**

For a long time, studies of plant secretory structures, in particular trichomes have been performed essentially on anatomy and were purely descriptive. There is a large body of publications of this type as summarized in reviews by<sup>9-11</sup>. Our biochemical and molecular knowledge of plant secretory structures is still limited in spite of some early remarkable work e.g. on mint glandular trichomes showing that glandular trichomes are the site of biosynthesis of monoterpenoids which are secreted into the subcuticular space<sup>12</sup>. Thus, it is necessary to create tools and to develop techniques for the isolation of such specialized organs in order to study them independently. Direct access to these glands would offer the possibility to directly investigate and understand both physiological and molecular processes in secretory systems of plants. The function of genes playing a key role in the biosynthesis and the transport of the secreted compounds could be addressed as well as the question how development of such highly specialized organs is regulated.

It should be emphasized that the production and secretion capacities of many secreting glands identified them as prime targets for metabolic engineering in plants with potential applications in a wide range of areas including human health or durable crop protection<sup>13,14</sup>. In addition, secretions by these different structures play important roles in the interactions of the plants with their abiotic and biotic environment in particular with insects. Secreted compounds can serve as attractants to pollinating insects, as in the case of volatile monoterpenes such as linalool produced by glandular trichomes, and sugars secreted by flower nectaries<sup>15</sup>. Resistance against insect pests is another important aspect as illustrated by various reports on this topic<sup>16,17</sup>.

In order to isolate single multicellular glands from the pitcher tissue of carnivorous *Nepenthes* species a fast and efficient micropreparation technique was developed and employed. This method is an alternative to the widespread but, in terms of applications, limited technique of laser capture dissection. The advantage of the new methods is due to the fact that no previous fixation of plant material is necessary and, instead, the freshly harvested plant tissue can be used immediately. Experiments where PCR studies of individual genes were performed<sup>18</sup>, as well as successful RNA extraction combined with subsequent gene expression analysis showed that the new preparation method is highly efficient. Thus, its application is suitable for many other molecular and biochemical studies.

However, also some drawbacks exit using the mechanical micropreparation. In particular, it often happens that the prepared samples stick to the forceps and it is hard to stripe off the preparations. That will cause a rapid drying-out and, as a consequence, a loss of the samples. Thus, experience with the technique is important to be fast and precise in any step, starting from positioning of the forceps to transferring the sample into the target tubes. Automation of such steps might be helpful. In addition, for any work with very small tissue samples it is necessary to provide a (micro)-environment that guarantees non-destructive conditions. That means, for example, high humidity, low temperatures, and low UV-irradiation in order to keep subcellular structures (DNA, RNA, proteins) intact.

The technique presented appears to be adaptable to a wide range of plants that carry multicellular glands, offering the ability to study the different secretory systems. However, preparation of glands or other microstructures from living tissues is only the first step providing the samples of interest for further analyses. The subsequent analytical techniques that will be applied in the particular studies have different requirements in terms of the nature and quality of the sample. This might require individual methodical developments for different samples such as newly designed forceps or the use of anti-adhesive materials. In the long term, this technique will create a series of subsequent approaches which are relevant to understand plant secretory structures, including deep sequencing of cDNAs, proteomics or transcription promoter analyses specific to glandular cells could be performed in various plant species, not only carnivorous plants such as *Nepenthes* or *Drosera* but also, for example, from *Artemisia*, *Mentha* or *Cistus* or many others. In general, applications beyond plant gland physiology are easily conceivable such as single cell analysis in biological and medicinal sciences or even forensic.

**Acknowledgments:** We thank the greenhouse team at the MPI-CE at the Botanical Garden of the Friedrich Schiller University, Jena, for cultivating plants and Wilhelm Boland and the Max Planck Society for support.

**Disclosures:** We have nothing to disclose.

#### **References:**

- [1] Adlassnig, W., Peroutka, M., Land, I. & Lichtscheidl, I.K. Glands of carnivorous plants as a model system in cell biological research. *Acta Botanica Gallica*. **152**, 111-124 (2005).
- [2] Darwin, C. Insectivorous plants. London: John Murray (1875).
- [3] Mithöfer A. Carnivorous pitcher plants: Insights in an old topic. *Phytochemistry* **72,** 1678-1682 (2011).
- [4] Owen Jr. T.P., Lennon, K.A., Santo, M.J. & Anderson, A.N. Pathways for nutrient transport in the pitchers of the carnivorous plant *Nepenthes alata. Annals of Botany* **84**, 459-466 (1999).
- [5] Brandt, S.P. Microgenomics: gene expression analysis at the tissue-specific and single-cell levels. *Journal of Experimental Botany* **56**, 495-505 (2005).
- [6] Hölscher, D. & Schneider, B. Application of laser-assisted microdissection for tissue and cell-specific analysis of RNA, proteins, and metabolites. *Progress in Botany* **69**, 141-167 (2008).
- [7] Rottloff, S., Müller, U., Kilper, R., Mithöfer, A. Micropreparation of single secretory glands from the carnivorous plant *Nepenthes*. *Analytical Biochemistry* **394,** 135-137 (2009).
- [8] Liu, W.H. & Saint, D.H. A new quantitative method of real time reverse transcription polymerase chain reaction assay based on simulation of polymerase chain reaction kinetics. *Analytical Biochemistry* **302**, 52-59 (2002).
- [9] Werker, E. Trichome diversity and development. *Advances in Botanical Research Incorporating Advances in Plant Pathology* **31,** 1-35 (2000).
- [10] Fahn, A. Structure and function of secretory cells. *Advances in Botanical Research Incorporating Advances in Plant Pathology* **31,** 37-75 (2000).
- [11] Juniper, B.E., Robins, R.J. & Joel, D.M. The carnivorous plants. London: Academic Press (1989).
- [12] Gershenzon, J., McCaskill, D., Rajaonarivony, J., Mihaliak, C., Karp., F. & Croteau, R. Isolation of secretory-cells from plant glandular trichomes and their use in biosynthetic-studies of monoterpenes and other gland products. *Analytical Biochemistry* **200**, 130-138 (1992).
- [13] Rontein, D., Onillon, S., et al. CYP725A4 from yew catalyzes complex structural rearrangement of taxa-4(5),11(12)-diene into the cyclic ether 5(12)-oxa-3(11)-cyclotaxane. *Journal of Biological Chemistry* **283**, 6067-75 (2008).
- [14] Aharoni, A., Jongsma, M.A. & Bouwmeester, H.J. Volatile science? Engineering of terpenoid biosynthesis in plants. *Trends in Plant Science* **10**, 594-602 (2005).
- [15] Mithöfer, A., Boland, W., Maffei, M.E. Chemical ecology of plant-insect interactions. *Annual Plant Reviews* **34**, 261-291 (2009).
- [16] Kennedy, G.G. Tomato, pests, parasitoids, and predators: Tritrophic interactions involving the genus Lycopersicon. Annual Review of Entomology 48, 51-72 (2003).


- [17] Mithöfer, A. & Boland, W. Plant defense against herbivores: Chemical aspects. *Annual Review of Plant Biology* **63**, 431-450 (2012).
- [18] Rottloff, S., Stieber, R., Maischak, H., Turini, F.G., Heubl, G. & Mithöfer, A. Functional characterization of a class III acid endochitinase from the traps of the carnivorous pitcher plant genus, Nepenthes. *Journal of Experimental Botany* **62**, 4649-4647 (2011).


Fig 1 Click here to download high resolution image




Figure 1

Fig 2 Click here to download high resolution image










Figure 2

Fig 3 Click here to download high resolution image




Fig 4 Click here to download high resolution image




Figure 4

Name of Material/ Equipment	Company	Catalog Number		
RNAqueous-Micro Kit	Life Technologies GmbH, Germany	AM1931		
Dynabeads mRNA DIRECT Micro Kit	Life Technologies GmbH, Germany	61011		
DynaMag-2 magnet	Life Technologies GmbH, Germany	12321D		
SuperScript CellsDirect cDNA Synthesis System	Life Technologies GmbH, Germany	18080-200		
RNaseOUT	Life Technologies GmbH, Germany	10777-019		
Random Hexamers (0.4 μg/μl)	QIAGEN GmbH, Germany	79236		
10 mM dNTP Mix, molecular biology grade	Fermentas GmbH, Germany	R0192		
Taq DNA polymerase, (recombinant)	Fermentas GmbH, Germany	EP0402		
TURBO DNase	Life Technologies GmbH, Germany	AM2238		
Agarose	Carl Roth GmbH Co. Germany	T846.3		
GeneRuler Low Range DNA ladder, ready-to-				
use	Fermentas GmbH, Germany	SM1193		
Brilliant II SYBR Green QPCR Master Mix	Agilent Technologies, Inc., USA	600828		
RNaseZAP	Sigma-Aldrich, Taufkirchen, Germany	R2020		
Razor blade	Carl Roth GmbH & Co., Germany	CK08.1		
tweezers	Carl Roth GmbH & Co., Germany	2801.1 / 2855.1		
GE NanoVue spectrophotometer	GE Healthcare Europe GmbH, Germany	28-9569-62		
Mastercycler gradient	Eppendorf AG, Hamburg, Germany	5331 000.010		
Mx3000P Real-Time PCR System	Stratagene, USA	401403		
Agagel Standard Horizontal Gel Electrophoresis				
chamber	Biometra	Discontinued		
Micromanipulator	aura optik gmbh	aureka		
Stereomicroscope	Carl Zeiss Microscopy, Jena, Germany	SteREO Lumar.V12		
Lumar Filter set 01	Carl Zeiss, Microscopy, Jena, Germany	485001-0000-000		
Lumar Filter set 09	Carl Zeiss, Microscopy, Jena, Germany	485009-0000-000		

# **Comments/Description**

Comments/ Description
Isolation of total RNA, 50 samples
5 ml Dynabeads Oligo (dT) <sub>25</sub>
25 reactions
Rnase inhibitor, 5000 units
100 μΙ
1 ml
Supplied with 10x <i>Taq</i> buffer (with KCl or
(NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub> ) and 25 mM MgCl <sub>2</sub>
1000 units
50 μg
Single kit
Equipped with a microforceps
Equipped with an AxioCam
UV filter
non-UV filter



# ARTICLE AND VIDEO LICENSE AGREEMENT

Title of Article:	Isolation of viable multicellular glands from tirsue of the connicovous plant,
Author(s):	ROTTLOFF, MITHÖFER, MÜLLER, WILPER
Item 1 (check one http://www	box): The Author elects to have the Materials be made available (as described at .jove.com/publish) via: Standard Access  Open Access
Item 2 (check one b	ox):
The Aut	hor is NOT a United States government employee.
	thor is a United States government employee and the Materials were prepared in the sor her duties as a United States government employee.
	thor is a United States government employee but the Materials were NOT prepared in the sor her duties as a United States government employee.

# ARTICLE AND VIDEO LICENSE AGREEMENT

1. Defined Terms. As used in this Article and Video License Agreement, the following terms shall have the following meanings: "Agreement" means this Article and Video License Agreement; "Article" means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein; "Author" means the author who is a signatory to this Agreement; "Collective Work" means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; "CRC License" means the Creative Commons Attribution-Non Commercial-No Derivs 3.0 Unported Agreement, the terms and conditions of which can be found http://creativecommons.org/licenses/by-ncnd/3.0/legalcode; "Derivative Work" means a work based upon the Materials or upon the Materials and other preexisting works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; "Institution" means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; "JoVE" means MyJove Corporation, a Massachusetts corporation and the publisher of The Journal of Visualized Experiments; "Materials" means the Article and / or the Video; "Parties" means the Author and JoVE; "Video" means any video(s) made by the Author, alone or in conjunction with any other parties, or by JoVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion of the Article, and in which the Author may or may not appear.

- 2. <u>Background</u>. The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JoVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.
- 3. Grant of Rights in Article. In consideration of JoVE agreeing to publish the Article, the Author hereby grants to JoVE, subject to Sections 4 and 7 below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the "Open Access" box has been checked in Item 1 above, JoVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.
- 4. Retention of Rights in Article. Notwithstanding the exclusive license granted to JoVE in Section 3 above, the



# ARTICLE AND VIDEO LICENSE AGREEMENT

Author shall, with respect to the Article, retain the non-exclusive right to use all or part of the Article for the non-commercial purpose of giving lectures, presentations or teaching classes, and to post a copy of the Article on the Institution's website or the Author's personal website, in each case provided that a link to the Article on the JoVE website is provided and notice of JoVE's copyright in the Article is included. All non-copyright intellectual property rights in and to the Article, such as patent rights, shall remain with the Author.

- 5. Grant of Rights in Video Standard Access. This Section 5 applies if the "Standard Access" box has been checked in Item 1 above or if no box has been checked in Item 1 above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby acknowledges and agrees that, Subject to Section 7 below, JoVE is and shall be the sole and exclusive owner of all rights of any nature, including, without limitation, all copyrights, in and to the Video. To the extent that, by law, the Author is deemed, now or at any time in the future, to have any rights of any nature in or to the Video, the Author hereby disclaims all such rights and transfers all such rights to JoVE.
- 6. Grant of Rights in Video Open Access. This Section 6 applies only if the "Open Access" box has been checked in Item 1 above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby grants to JoVE, subject to Section 7 below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Video in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Video into other languages, create adaptations, summaries or extracts of the Video or other Derivative Works or Collective Works based on all or any portion of the Video and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. For any Video to which this Section 6 is applicable, JoVE and the Author hereby grant to the public all such rights in the Video as provided in, but subject to all limitations and requirements set forth in, the CRC License.
- 7. Government Employees. If the Author is a United States government employee and the Article was prepared in the course of his or her duties as a United States government employee, as indicated in Item 2 above, and any of the licenses or grants granted by the Author hereunder exceed the scope of the 17 U.S.C. 403, then the rights granted hereunder shall be limited to the maximum rights permitted under such statute. In such case, all provisions contained herein that are not in conflict with such statute shall remain in full force and effect, and all provisions contained herein that do so conflict

shall be deemed to be amended so as to provide to JoVE the maximum rights permissible within such statute.

- 8. <u>Likeness</u>, <u>Privacy</u>, <u>Personality</u>. The Author hereby grants JoVE the right to use the Author's name, voice, likeness, picture, photograph, image, biography and performance in any way, commercial or otherwise, in connection with the Materials and the sale, promotion and distribution thereof. The Author hereby waives any and all rights he or she may have, relating to his or her appearance in the Video or otherwise relating to the Materials, under all applicable privacy, likeness, personality or similar laws.
- 9. Author Warranties. The Author represents and warrants that the Article is original, that it has not been published, that the copyright interest is owned by the Author (or, if more than one author is listed at the beginning of this Agreement, by such authors collectively) and has not been assigned, licensed, or otherwise transferred to any other party. The Author represents and warrants that the author(s) listed at the top of this Agreement are the only authors of the Materials. If more than one author is listed at the top of this Agreement and if any such author has not entered into a separate Article and Video License Agreement with JoVE relating to the Materials, the Author represents and warrants that the Author has been authorized by each of the other such authors to execute this Agreement on his or her behalf and to bind him or her with respect to the terms of this Agreement as if each of them had been a party hereto as an Author. The Author warrants that the use, reproduction, distribution, public or private performance or display, and/or modification of all or any portion of the Materials does not and will not violate, infringe and/or misappropriate the patent, trademark, intellectual property or other rights of any third party. The Author represents and warrants that it has and will continue to comply with all government, institutional and other regulations, including, without limitation all institutional, laboratory, hospital, ethical, human and animal treatment, privacy, and all other rules, regulations, laws, procedures or guidelines, applicable to the Materials, and that all research involving human and animal subjects has been approved by the Author's relevant institutional review board.
- 10. JoVE Discretion. If the Author requests the assistance of JoVE in producing the Video in the Author's facility, the Author shall ensure that the presence of JoVE employees, agents or independent contractors is in accordance with the relevant regulations of the Author's institution. If more than one author is listed at the beginning of this Agreement, JoVE may, in its sole discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JoVE reserves the right, in its absolute and sole discretion and without giving any reason therefore, to accept or decline any work submitted to JoVE. JoVE and its employees, agents and independent contractors shall have full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JoVE has sole discretion as to the method of making and publishing the Materials, including,



# ARTICLE AND VIDEO LICENSE AGREEMENT

without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

11. Indemnification. The Author agrees to indemnify JoVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JoVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JoVE, making of videos by JoVE, or publication in JoVE or elsewhere by JoVE. The Author shall be responsible for, and shall hold JoVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to the making of a video by JoVE its employees, agents or independent contractors. All sterilization, cleanliness or decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's expense. All indemnifications provided herein shall include JoVE's attorney's fees and costs related to said losses or damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JoVE, its employees, agents or independent contractors.

- 12. Fees. To cover the cost incurred for publication, JoVE must receive payment before production and publication the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JoVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.
- 13. Transfer, Governing Law. This Agreement may be assigned by JoVE and shall inure to the benefits of any of JoVE's successors and assignees. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to me one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement required per submission.

<b>AUTHOR:</b>	
Name:	HITHOFER, Axel
Department:	Bioorganic Chemishy
Institution:	MAX PLANCK INSTITUTE for Chemical Ecology
Article Title:	Inolation of viable would rellular glands from time of the
	X.0 0-18.1
Signature:	Sul lu Cot Date: Haver, 26, 2013

Please submit a signed and dated copy of this license by one of the following three methods:

- 1) Upload a scanned copy as a PDF to the JoVE submission site upon manuscript submission (preferred);
- 2) Fax the document to +1.866.381.2236; or
- 3) Mail the document to JoVE / Attn: JoVE Editorial / 17 Sellers St / Cambridge, MA 02139

For questions, please email editorial@jove.com or call +1.617.945.9051.

	 	***************************************	22000	 Annual Ann	 	WA A-1000	***************************************	
MS # (internal use):								

Response letter (JoVE50993R1):

Dear Editor,

Thank you for the comments of the referees on our manuscript. In the new version we addressed almost all points they raised.

All changes have been made using the track change function. Please find the responses to the reviewers' comments below:

#### Reviewer#1

**RESPONSE:** We corrected all typos or phrasing that could be misunderstood. In addition, we explained in more detail the nature of the glands in *Nepenthes* (last paragraph of the Introduction) as requested. We also included an additional reference for the attractants (#15; Mithöfer et al., 2009).

#### Reviewer#2

**RESPONSE:** We changed that particular sentence and, in addition, went through the whole manuscript to improve the language.

#### Reviewer#3

**RESPONSE:** We corrected "°C" and "%" throughout the manuscript. We also weakened the statement that "laser capture microdissection" technique is too limited for many approaches. Concerning the referee's comments about more recent publications on secretory structures and a deep discussion on "laser capture microdissection" versus the mechanical microdissection technique that is described in our paper, we think that such a discussion is far beyond the scope of this article.

## Reviewer#4

**RESPONSE:** See response to reviewers 1 & 2. For us and the other three referees a complete restructuring of the manuscript seems not to be necessary.

Sincerely,

Axel Mithöfer