

Procedure - Synergy 2 (Com3)

Select steps

Actions

- Read
- Set Temperature
- Shake
- Dispense

Kinetic

- Start Kinetic
- Monitor Well

Pause

- Delay
- Plate Out/In
- Stop/Resume

Process Mode

- Well Mode
- Plate Mode


Other

- Comment
- Options

Plate Type: 96 WELL PLATE

Cuvette

Select wells: Per step At runtime

Description	Comments
 Dispense 10 µL using dispenser 1	

Shake Step

Intensity: Medium

Duration: 0:30 MM:SS

Continuous Shake

OK Cancel Help

Validate OK Cancel Help