

Journal of Visualized Experiments

Measuring protein stability in living zebrafish embryos using Fluorescence Decay After Photoconversion (FDAP) --Manuscript Draft--

Manuscript Number:	JoVE52266R2
Full Title:	Measuring protein stability in living zebrafish embryos using Fluorescence Decay After Photoconversion (FDAP)
Article Type:	Invited Methods Article - JoVE Produced Video
Keywords:	Fluorescence Decay After Photoconversion (FDAP); protein stability; clearance rate; confocal microscopy; photoconvertible protein; Dendra2; embryo; zebrafish
Manuscript Classifications:	2.1.50.150.900.493.200.244.828: Zebrafish; 5.5.595.395: Microscopy, Confocal; 7.2.111.700: Protein Stability
Corresponding Author:	Patrick Müller Friedrich Miescher Laboratory of the Max Planck Society Tübingen, Baden-Württemberg GERMANY
Corresponding Author Secondary Information:	
Corresponding Author E-Mail:	patrick.mueller@tuebingen.mpg.de
Corresponding Author's Institution:	Friedrich Miescher Laboratory of the Max Planck Society
Corresponding Author's Secondary Institution:	
First Author:	Katherine Rogers
First Author Secondary Information:	
Other Authors:	Katherine Rogers Alexander Bläßle Alexander Schier
Order of Authors Secondary Information:	
Abstract:	Protein stability influences many aspects of biology, and measuring the clearance kinetics of proteins can provide important insights into biological systems. In FDAP experiments, the clearance of proteins within living organisms can be measured. A protein of interest is tagged with a photoconvertible fluorescent protein, expressed in vivo and photoconverted, and the decrease in the photoconverted signal over time is monitored. The data is then fitted with an appropriate clearance model to determine the protein half-life. Importantly, the clearance kinetics of protein populations in different compartments of the organism can be examined separately by applying compartmental masks. This approach has been used to determine the intra- and extracellular half-lives of secreted signaling proteins during zebrafish development. Here, we describe a protocol for FDAP experiments in zebrafish embryos. It should be possible to use FDAP to determine the clearance kinetics of any taggable protein in any optically accessible organism.
Author Comments:	
Additional Information:	
Question	Response
If this article needs to be "in-press" by a certain date to satisfy grant requirements, please indicate the date below and explain in your cover letter.	
If this article needs to be filmed by a	

certain date to due to
author/equipment/lab availability, please
indicate the date below and explain in
your cover letter.



Dr. Nam Nguyen
JoVE Editorial Board
1 Alewife Center, Suite 200
Cambridge, MA 02140

Dr. Patrick Müller
Research Group Leader
Systems Biology of
Development
Tel.: +49-(0)7071-601-815
Fax: +49-(0)7071-601-801
pmueller@tuebingen.mpg.de

FDAP invited manuscript

June 24, 2014

Dear Dr. Nguyen,

Thank you for your editorial comments on our manuscript “Measuring protein stability in living zebrafish embryos using Fluorescence Decay After Photoconversion (FDAP)”. Below we address your comments.

Editorial comments

1. Please take this opportunity to thoroughly proofread the manuscript to ensure that there are no spelling or grammar issues. The JoVE editor will not copy-edit your manuscript and any errors in the submitted revision may be present in the published version.

2. Please ensure that **all text** in the protocol section is written in the imperative tense as if telling someone how to do the technique (e.g., “Do this,” “Ensure that,” etc.). The actions should be described in the imperative tense in complete sentences wherever possible. Avoid usage of phrases such as “could be,” “should be,” and “would be” throughout the Protocol. Any text that cannot be written in the imperative tense may be added as a “Note.” However, notes should be concise and used sparingly.

We have changed several sentences in Section 5 to the imperative mood. We also eliminated “should” and “would” where possible and removed the note at the beginning of Section 5.

For example: Step 5 could be moved to the Representative Results.

Multiple reviewers requested further explanation of control experiments, and reviewer #3 specifically asked for the control experiments to appear in the numbered protocol. We therefore created Section 5 to address these important concerns.

3. The Protocol should contain only action items that direct the reader to do something. Please move the discussion about the protocol to the Discussion.

We have changed several sentences to the imperative mood (see response to point #2 above) and include brief notes to clarify subsequent instructions.

4. The Figure Legends should include a title and a short description of the data presented in the Figure and relevant symbols. The Discussion of the Figures should be placed in the Representative Results. Details of the methodology should not be in the Figure Legends, but rather the Protocol.

We have shortened the figure legends where possible.

Yours sincerely,

Patrick Müller

Title:
Measuring protein stability in living zebrafish embryos using Fluorescence Decay After Photoconversion (FDAP)

Authors:
Katherine W. Rogers
Department of Molecular and Cellular Biology
Harvard University
Cambridge, USA
kwrogers@fas.harvard.edu

Alexander Bläßle
Systems Biology of Development Group
Friedrich Miescher Laboratory of the Max Planck Society
Tübingen, Germany
alexander.blaessle@tuebingen.mpg.de

Alexander F. Schier
Department of Molecular and Cellular Biology
Harvard University
Cambridge, USA
schier@fas.harvard.edu

Patrick Müller
Systems Biology of Development Group
Friedrich Miescher Laboratory of the Max Planck Society
Tübingen, Germany
patrick.mueller@tuebingen.mpg.de

Corresponding author: Patrick Müller

Keywords:
Fluorescence Decay After Photoconversion (FDAP), protein stability, clearance rate, confocal microscopy, photoconvertible protein, Dendra2, embryo, zebrafish

Short abstract:
Protein levels in cells and tissues are often tightly regulated by the balance of protein production and clearance. Using Fluorescence Decay After Photoconversion (FDAP), the clearance kinetics of proteins can be experimentally measured *in vivo*.

Long abstract:
Protein stability influences many aspects of biology, and measuring the clearance kinetics of proteins can provide important insights into biological systems. In FDAP experiments, the clearance of proteins within living organisms can be measured. A protein of interest is tagged

with a photoconvertible fluorescent protein, expressed *in vivo* and photoconverted, and the decrease in the photoconverted signal over time is monitored. The data is then fitted with an appropriate clearance model to determine the protein half-life. Importantly, the clearance kinetics of protein populations in different compartments of the organism can be examined separately by applying compartmental masks. This approach has been used to determine the intra- and extracellular half-lives of secreted signaling proteins during zebrafish development. Here, we describe a protocol for FDAP experiments in zebrafish embryos. It should be possible to use FDAP to determine the clearance kinetics of any taggable protein in any optically accessible organism.

Introduction:

The levels of proteins in cells and organisms are determined by their rates of production and clearance. Protein half-lives can range from minutes to days¹⁻⁴. In many biological systems, the stabilization or clearance of key proteins has important effects on cellular activity. Modulation of intracellular protein stability is required for cell cycle progression^{5,6}, developmental signaling⁷⁻⁹, apoptosis¹⁰, and normal function and maintenance of neurons^{11,12}. Extracellular protein stability affects the distribution and availability of secreted proteins, such as morphogens^{13,14}, within a tissue.

Over the last few decades, protein stability has mainly been assessed in cell culture using radioactive pulse-labeling or cycloheximide chase experiments¹⁵. In such pulse-chase experiments, cells are either transiently exposed to a “pulse” of radioactive amino acid precursors that are incorporated into newly synthesized proteins, or they are exposed to cycloheximide, which inhibits protein synthesis. Cultured cells are then collected at different time points, and either immunoprecipitation followed by autoradiography (in radioactive pulse-chase experiments) or western blotting (in cycloheximide experiments) is used to quantify the clearance of protein over time.

Conventional protein stability assays have several shortcomings. First, proteins in these assays are often not expressed in their endogenous environments, but rather in tissue culture and sometimes in cells from different species. For proteins whose stability is context-dependent, this approach is problematic. Second, it is not possible to follow protein clearance in individual cells or organisms over time, and the data from these assays reflects an average of different populations of cells at different time points. Since individual cells may have started with different amounts of protein, may have taken up the radioactive label or cycloheximide at different times, or may have different clearance kinetics, such aggregate data could be misleading. Finally, in the case of cycloheximide chase experiments, addition of the protein synthesis inhibitor may have unintended physiological effects that could artificially alter protein stability¹⁶⁻¹⁸. These shortcomings can be avoided by using Fluorescence Decay After Photoconversion (FDAP), a technique that utilizes photoconvertible proteins to measure protein clearance dynamically in living organisms¹⁹⁻²⁵ (see **Discussion** for limitations of the FDAP technique).

Photoconvertible proteins are fluorescent proteins whose excitation and emission properties

change after exposure to specific wavelengths of light²⁶. One commonly used variant is Dendra2, a “green-to-red” photoconvertible protein that initially has excitation and emission properties similar to green fluorescent proteins, but after exposure to UV light—“photoconversion”—its excitation/emission properties become similar to those of red fluorescent proteins^{23,27}. Importantly, new protein produced after photoconversion will not have the same excitation/emission properties as the photoconverted protein, allowing decoupling of production and clearance upon photoconversion and observation of only a pool of photoconverted protein. Tagging proteins of interest with photoconvertible proteins thus provides a convenient way to pulse-label proteins in intact, optically accessible living organisms.

In FDAP assays (**Figure 1A**), proteins of interest are tagged with a photoconvertible protein and expressed in living organisms (**Figure 1B**). The fusion proteins are photoconverted, and the decrease in photoconverted signal over time is monitored by fluorescence microscopy (**Figure 1C**). The data is then fitted with an appropriate model to determine the half-life of the fusion protein (**Figure 1D**).

The FDAP assay described here was designed to determine the extracellular half-lives of secreted signaling proteins in zebrafish embryos during early embryogenesis¹⁹. However, this approach can be adapted to any transparent model organism that tolerates live imaging, and could be used to monitor the clearance of any taggable intracellular or extracellular protein. Variations of the technique described here have been performed in cultured cells^{20,23} and *Drosophila*²² and mouse²¹ embryos.

Protocol:

1. Generating a photoconvertible fusion construct and injecting dechorionated zebrafish embryos

1.1) Generate a functional construct containing the protein of interest fused to a green-to-red photoconvertible protein (see **Discussion**), then use *in vitro* transcription to generate capped mRNA encoding the fusion protein as in Müller *et al.*, 2012¹⁹.

1.2) Use pronase to remove the chorions from about 30 zebrafish embryos at the one-cell stage. Alternatively, manually dechorionate embryos using forceps²⁸.

Note: Embryos must be dechorionated for subsequent imaging. If desired, embryos can also be dechorionated just prior to imaging.

1.2.1) Make a 5 mg/ml stock solution of pronase from *Streptomyces griseus* in standard zebrafish embryo medium¹⁹. Rock the solution gently at room temperature for 10 min to allow the protease to dissolve. Aliquot 2 ml into microcentrifuge tubes and freeze at -20 °C.

1.2.2) Transfer one-cell stage embryos to a 5 cm diameter glass or agarose-coated plastic petri dish containing ~8 ml embryo medium. Add 2 ml of thawed pronase stock solution to the dish and incubate at room temperature for 5 min.

1.2.3) Avoid exposing embryos to air or plastic, as contact with either will cause dechorionated embryos to rupture. Fill a 200 ml glass beaker with embryo medium. Transfer the embryos to the beaker by tilting the petri dish while submerging it in the medium.

1.2.4) After the embryos have settled to the bottom of the beaker, pour out most of the embryo medium, then pour fresh embryo medium into the beaker. The mild swirling of the medium pouring into the beaker causes embryos to lose their weakened chorions.

1.2.5) Repeat step 1.2.4.

1.3) Transfer the dechorionated embryos to an agarose-coated injection dish²⁹ using a glass Pasteur pipette with a flamed tip. Flaming the pipette tip prevents jagged edges from injuring embryos.

1.4) Co-inject the mRNA and a 3 kDa Alexa488-dextran conjugate^{29,30} (**Figure 1B**; see **Discussion** for suggested mRNA and Alexa488-dextran amounts). Inject directly into the center of the cell (not the yolk) to ensure even distribution of mRNA and fluorescent dye once cleavage commences.

Note: The Alexa488 signal will be used during data analysis to generate compartmental masks in order to distinguish between intracellular and extracellular fluorescence.

1.5) Transfer injected embryos to a 1 - 2% agarose-coated well of a six-well plastic dish filled with embryo medium. Incubate in the dark at 28 °C until embryos have reached late sphere stage³¹ (approximately five hours post fertilization). Check embryos every one to two hours under a stereomicroscope and remove any debris generated by embryos that have died.

2. Mounting zebrafish embryos for photoconversion and imaging on an inverted confocal microscope

2.1) Use a stereomicroscope to identify one to five healthy embryos, and use a glass Pasteur pipette with a flamed tip to remove them from the dish.

2.2) Gently eject the embryos into a microcentrifuge tube containing ~1 ml of melted 1% low melting point agarose in 1x Danieau's embryo medium (see **Table of Materials/Equipment**) (**Figure 2A**).

Note: Agarose should have a temperature between 40 and 42 °C; higher temperatures could damage the embryos.

2.3) Draw the embryos back into the pipette along with some agarose. Gently eject the agarose and embryos onto the cover glass of a glass-bottom dish (**Figure 2B**). Ensure that the thickness of the cover glass is compatible with the objective on the confocal microscope.

2.3.1) Re-use the glass pipette if desired. To clean the residual agarose out of the pipette and prevent clogging, quickly pipette embryo medium up and down. Place a 15 ml tube filled with ~5 ml of embryo medium next to the stereomicroscope for this purpose.

2.4) Use a metal probe to position the embryos so that the animal pole (blastoderm) faces the cover glass. Work quickly since the agarose will solidify in 20 - 30 s. Use the stereomicroscope to monitor the embryos' positions and readjust as necessary until the agarose hardens.

2.5) Repeat steps 2.1 – 2.4 until the desired number of embryos has been mounted.

Note: In a typical experiment, four agarose drops containing four or five embryos each will fit easily on the cover glass. About 16 embryos can be imaged during a single ideal experiment (**Figure 2C**).

2.6) When the agarose has solidified, fill the glass-bottom dish with 1x Danieau's embryo medium.

3. Photoconverting and measuring the decrease of the photoconverted signal

Note: A 25x or 40x water objective is appropriate for the size and refractive index of zebrafish embryos. It is best to use immersion oil with the same refractive index as water rather than actual water, since water will evaporate during the course of the five-hour experiment. Ensure that the immersion oil is designed to be used with a water (not oil) objective.

3.1) Place a large drop of immersion oil on the objective to ensure that the oil film between the objective and cover glass will not break as the stage moves to different embryo positions during imaging. Securely place the glass-bottom dish onto the stage so that the dish will not shift when the stage moves. If possible, use a heated stage at 28 °C, the optimal temperature for zebrafish development.

3.2) Define each embryo's position in the confocal microscope's software package. Adjust the z-depth for each embryo, and attempt to target roughly the same plane in each embryo.

Note: About 30 μm from the animal pole is a good depth since at this depth the enveloping layer of the embryo can be avoided, imaging area is maximized, and light scattering is minimal. A single optical slice with a thickness of ~3.3 μm provides sufficient data; there is no need to acquire a z-stack (see **Section 5**).

3.3) Collect two signals during the experiment: the “green” signal from the Alexa488-dextran conjugate—which will be used during data analysis to isolate extracellular and intracellular fluorescence—and the “red” signal from the fusion protein after it is photoconverted.

3.3.1) Excite Alexa488 using a 488 nm laser, and collect emitted fluorescence between ~500 and 540 nm.

Note: After photoconversion, many green-to-red photoconvertible proteins (*e.g.*, Dendra2) can be excited with a 543 nm laser and emit fluorescence between ~550 and 650 nm. Adjust as necessary based on the photoconvertible protein used.

3.4) Acquire “pre-photoconversion” images, and configure the confocal microscope’s software to image each of the previously defined positions (from step 3.2) with the appropriate imaging conditions every 10 or 20 min for a five-hour time course (see **Section 5** and **Discussion**).

3.5) To photoconvert the fusion protein, switch to a 10x objective and expose groups of embryos to UV light from a mercury arc lamp with a ~300-400 nm excitation filter at 100% output for 2 min. Shift the focus along the z-axis to promote uniform photoconversion (see **Section 5**). Ensure that the immersion oil does not drip onto the 10x objective during photoconversion.

Note: The shifting of focus during photoconversion could be automated in order to avoid variability among experimenters.

3.6) Switch back to the 25x or 40x objective immediately after photoconversion. Ensure that the previously defined positions from step 3.2 are still accurate. If the dish shifted during photoconversion, re-define the positions of the embryos.

3.6.1) Start the program created in step 3.4 and allow imaging to continue for five hours. Note the time elapsed between photoconversion and the start of imaging for each embryo.

3.7) Occasionally check on the experiment. Monitor the level of Danieau’s medium and add more if necessary. Restart the software if it has stalled.

3.8) In order to determine the background fluorescence values that will be used during data analysis to estimate the asymptote of an exponentially decreasing model, include some uninjected embryos in the experiment. To determine the instrument noise, which will also be used during subsequent data analysis, acquire an image in the absence of a sample.

4. Analyzing the data using PyFDAP

4.1) Visually inspect the time course data sets from each embryo. Discard data sets from embryos that died during imaging, that shifted significantly, that have very low levels of

photoconverted signal, or that contain regions of cells that look unusual and have stopped moving and dividing (typical of injured or sick embryos).

Note: Occasionally, bubbles in the immersion oil or other artifacts will appear in one or two images in an otherwise usable data set. Note any images that contain artifacts; they will be discarded later, and the remaining time points from such a data set can still be analyzed.

4.2) Use the Python-based software package PyFDAP to analyze the FDAP data. PyFDAP calculates half-lives by determining the average intracellular and extracellular red fluorescence intensity in each image and fitting the data with an exponentially decreasing function (**Figure 3**).

4.2.1) Download PyFDAP (see **Table of Materials/Equipment**).

4.2.2) Use PyFDAP to separate intracellular and extracellular photoconverted signal (**Figure 3A,B**). Use the Alexa488 signal, which is strictly intracellular, to create an intracellular mask. Apply this mask to the corresponding red channel image to prevent intracellular pixels from being considered when calculating average extracellular intensity. To measure average intracellular intensity, invert the mask.

4.2.3) In PyFDAP, display the masked images generated in step 4.2.2. Visually inspect these images and discard data sets in which masks do not accurately distinguish intracellular from extracellular space (this should be rare; note that cell membranes are included in images in which extracellular space has been masked, but they could be removed by altering the thresholding algorithm or by introducing a membrane mask (*e.g.*, using membrane-CFP)). Also discard any single images containing artifacts (*e.g.*, bubbles in the immersion oil) identified in step 4.1.

4.2.4) Use PyFDAP to calculate average extracellular and intracellular fluorescence intensities for each image. PyFDAP calculates these averages by summing the intensities of pixels that fall outside of the mask and dividing by the total number of pixels summed.

4.2.5) Fit the fluorescence data (**Figure 3C**) with the following exponential function:

$$c(t) = c_0 e^{-kt} + y_0$$

where t is time post-photoconversion, $c(t)$ is intensity at a given value of t , c_0 is the intensity at $t = 0$, k is the clearance rate constant, and y_0 is the asymptote that the function approaches as fluorescence decreases (**Figure 1D**). y_0 can be constrained based on the measurements in step 3.8¹⁹.

4.2.6) Use PyFDAP to calculate the extracellular and intracellular protein half-lives (τ) from the clearance rate constants (k) using the following relationship:

$$\text{Half-life } \tau = \frac{\ln(2)}{k}$$

5. Control experiments to assess photobleaching, inadvertent photoconversion, and photoconversion uniformity

5.1) Assessing photobleaching

Note: Photobleaching could cause an artifactual decrease in fluorescence intensity that reflects the bleaching properties of the fluorescent protein in addition to the clearance of the protein of interest.

5.1.1) To assess possible photobleaching, perform one set of FDAP experiments with 10 min intervals between imaging and a second set with 20 min intervals between imaging (**Figure 4**). Analyze the data from both sets of experiments using PyFDAP as described in Section 4.

5.1.2) Compare the half-lives from the 10 and 20 min interval experiments. Longer half-lives from 20 min interval experiments indicate significant photobleaching. If the half-lives from both experiments are identical, photobleaching is not a significant concern.

5.1.3) Alternatively, assess photobleaching by acquiring a series of ~30 images immediately after photoconversion. A significant decrease in fluorescence intensity indicates significant photobleaching.

5.1.4) If photobleaching is detected, use lower laser power, decrease imaging time, or consider using a more photostable photoconvertible protein³².

5.2) Assessing inadvertent photoconversion

Note: Dendra2 can be photoconverted using 488 nm illumination²⁷. When exciting Alexa488 with the 488 nm laser as described in step 3.3.1, inadvertent photoconversion and therefore an artifactual increase in the apparent half-life of the protein of interest is possible. However, we and others³³ have found that 488 nm illumination is an inefficient method of photoconversion in zebrafish embryos.

5.2.1) Use the control experiment described in step 5.1.1 to detect inadvertent photoconversion. Compare the half-lives from the 10 and 20 min interval experiments. Shorter half-lives from 20 min interval experiments indicate significant inadvertent photoconversion. If the half-lives from both experiments are identical, inadvertent photoconversion is not a significant concern.

5.2.2) If inadvertent photoconversion is detected, use a lower 488 nm laser power and shorter imaging times to avoid inadvertently photoconverting Dendra2.

5.3) Assessing photoconversion uniformity

Note: If photoconversion is biased toward the animal pole of the embryo, the decrease in fluorescence will be influenced by protein diffusion or cell movement into deeper planes (**Figure 5A**).

5.3.1) To determine whether photoconversion is uniform, express a secreted photoconvertible protein (for experiments with extracellular fusion proteins) or a cytoplasmic photoconvertible protein (for experiments with intracellular fusion proteins). Photoconvert as usual, then acquire a z-stack encompassing most of the blastoderm every 20 min for 80 min.

5.3.2) If photoconversion is biased toward the animal pole, the fluorescence intensity in deeper planes will increase over time due to diffusion or cell movement (**Figure 5B**). If non-uniform photoconversion is detected, focus deeper into the embryos during photoconversion.

Representative Results:

FDAP has been used to determine the half-lives of extracellular signaling proteins in zebrafish embryos¹⁹. One of these proteins, Squint, induces expression of mesendodermal genes during embryogenesis³⁴. Squint-Dendra2 activates expression of mesendodermal genes at levels similar to untagged Squint, as demonstrated by qRT-PCR and *in situ* hybridization assays¹⁹. Embryos were co-injected with Alexa488-dextran and mRNA encoding Squint-Dendra2 and subjected to the FDAP assay. A decrease in the extracellular photoconverted signal intensity over time is evident (**Figure 4A**). Extracellular intensity profiles from 23 embryos were generated using PyFDAP. The resulting data was fitted in PyFDAP with a first-order clearance kinetics model, and an average clearance rate constant k of $1.00 \times 10^{-4}/s$, corresponding to an average half-life τ of 116 min, was determined. Similar intensity profiles and clearance rate constants were obtained when the intervals between imaging were 10 or 20 min, suggesting that photobleaching or inadvertent photoconversion did not contribute significantly to intensity changes (**Figure 4B**).

Figure 1: Fluorescence Decay After Photoconversion (FDAP) overview. A) Workflow of an FDAP experiment. B) Injection of mRNA and a fluorescent dye into a zebrafish embryo at the one-cell stage. Protein is produced from the mRNA as the embryo develops over about five hours prior to imaging. The dye labels cells (green circles). C) The fusion protein is photoconverted using a UV pulse, and the decrease in the intensity of the photoconverted signal over time is monitored. D) The data is fitted with an exponentially decreasing function to obtain clearance rate constants (k) and half-lives (τ).

Figure 2: Mounting zebrafish embryos for FDAP experiments. A) Zebrafish embryos (blastoderm = white, yolk = black) are transferred from embryo medium (blue) into melted agarose (yellow). B) Embryos and agarose are placed onto the cover glass of a glass-bottom dish. Embryos are then manually positioned so that the animal pole faces the cover glass. A cross-section of a glass-bottom dish is shown. C) Schematic overview of a glass-bottom dish with several agarose drops containing four embryos each (view looking down into the dish).

Figure 3: Data analysis using PyFDAP. A) PyFDAP uses the Otsu thresholding algorithm³⁵ to generate intra- and extracellular masks from the intracellular Alexa488 signal (green). B) Photoconverted signal (red) from an embryo expressing a secreted Dendra2 fusion protein (Squint-Dendra2¹⁹). Average extra- and intracellular fluorescence intensities were calculated using the masks shown in (A). The space outside of the embryo was excluded from these calculations by discarding pixels outside of the yellow circle. C) PyFDAP screenshot displaying extracellular intensity data from an FDAP experiment (black circles) fitted with an exponentially decreasing function (red dashed line). The extracellular half-life is indicated by the red arrow.

Figure 4: Representative FDAP results. A) A representative embryo expressing secreted Squint-Dendra2 just prior to photoconversion (far left) and 27, 87, and 287 min post-photoconversion. B) To control for photobleaching and inadvertent photoconversion (see **Section 5**), experiments with 10 or 20 min intervals were performed (data from Müller *et al.*, 2012¹⁹). In PyFDAP, extracellular intensity profiles were generated, fitted with exponentially decreasing functions, and normalized by subtracting the fitted y_0 value from each data point and dividing by the fitted c_0 value. Data from the 10 min (black, $n = 11$) and 20 min (blue, $n = 12$) interval experiments were then averaged, respectively. Error bars indicate standard deviation.

Figure 5: Assessing photoconversion uniformity. A) Non-uniform photoconversion of an extracellular protein can lead to an erroneously short apparent half-life if photoconverted protein diffuses into deeper planes over time. B) To determine whether photoconversion is uniform, a z-stack covering most of the blastoderm is acquired at several time points post-photoconversion. Fluorescence intensity will increase in deeper planes over time if photoconversion was non-uniform (note that light scattering causes deeper planes to appear dimmer than higher planes).

Discussion:

The success of an FDAP experiment relies on the generation of a functional photoconvertible fusion protein. Tagging a protein can affect its biological activity and/or biophysical properties, including its localization, solubility, and stability³⁶⁻⁴¹. Be prepared to test the activity of several different fusion constructs in order to find one that is active. We have found that changing the position of the photoconvertible protein relative to the protein of interest or using longer linkers (*e.g.*, using the amino acid sequence LGDPPVAT¹⁹) can enhance the activity of the fusion protein. In the case of signaling proteins, the activity of the fusion protein can be determined by testing its ability to induce expression of target genes¹⁹. qRT-PCR or *in situ* hybridization provide good readouts of target gene expression¹⁹. Note that the protocol described here is designed to determine the stability of proteins in the early zebrafish embryo and would require modification to assess protein stability in other contexts.

The green-to-red photoconvertible protein Dendra2²⁷ has been used successfully in zebrafish FDAP experiments¹⁹ (**Figure 4**), but other options are available^{26,42}. To avoid potential artifacts due to aggregation of the fusion protein, choose a monomeric photoconvertible protein. Photoactivatable proteins can also be used in FDAP assays^{20-22,26}.

Before performing the FDAP experiment, several aspects of the protocol need to be optimized. Inject different amounts of mRNA to determine the lowest amount that provides useable signal after photoconversion; ~50 pg mRNA is a good starting amount. In order to generate meaningful compartmental masks (**Figure 3**), the Alexa488 signal must be bright enough to compete against the signal from the non-photoconverted fusion protein that is constantly produced from the injected mRNA; inject between 0.2 and 4 ng of Alexa488-dextran to find the optimal amount of fluorescent dye. Find the optimal post-conversion imaging conditions and use the same conditions for all experiments with a given construct. Use good quantitative imaging practices⁴³, and choose an appropriate dynamic range to avoid saturated pixels in the red channel. Determine the appropriate imaging interval for each fusion protein. Proteins with very short half-lives may require more frequent imaging over a shorter total time period. Establish the optimal photoconversion technique based on the organism and photoconvertible protein used. We describe one robust photoconversion method using a mercury arc lamp in step 3.5, but Dendra2 can also be photoconverted with a 405³³ or 488 nm laser²⁷.

One limitation of this FDAP protocol is that overexpression of the protein of interest is required. Overexpression could affect protein stability, for instance, through abnormal expression of other genes that modify the protein's clearance kinetics^{14,44}. If the protein of interest is a signaling molecule, consider performing experiments in the presence of a signaling inhibitor to determine whether blocking expression of target genes affects protein stability. In the future, it may be possible to generate transgenic embryos expressing photoconvertible fusions under the control of endogenous expression elements⁴⁵⁻⁵⁰. If transgenic embryos produce sufficient signal, FDAP experiments in a non-overexpression context are conceivable, perhaps using light-sheet microscopy to observe fluorescence decrease in the entire embryo⁵¹.

Possible further applications of FDAP include investigating the mechanisms that regulate protein stability by examining the effects of different perturbations (*e.g.*, expression of putative proteases) or protein modifications (*e.g.*, phosphorylation) on stability. For example, the factors controlling the extracellular stability of Squint are currently unknown. Many secreted developmental signals are internalized by cells⁵²⁻⁵⁵, which could contribute to the clearance of Squint and other ligands from the extracellular space. FDAP experiments in which internalization is blocked might provide information about mechanisms controlling extracellular protein clearance.

In contrast to conventional assays for measuring protein stability¹⁵, FDAP offers a microscopy-based alternative in which the clearance of proteins can be monitored over time within living organisms. Similar methods have been used to monitor protein clearance in model systems other than zebrafish embryos²⁰⁻²³. This FDAP protocol has the potential to be adapted to determine the half-life of any taggable protein in biological systems where live imaging is feasible.

Acknowledgements:

The authors would like to thank Jeffrey Farrell, James Gagnon, and Jennifer Bergmann for comments on the manuscript. KWR was supported by the National Science Foundation Graduate Research Fellowship Program during the development of the FDAP assay. This work was supported by grants from the NIH to AFS and by grants from the German Research Foundation (Emmy Noether Program), the Max Planck Society, and the Human Frontier Science Program (Career Development Award) to PM.

Disclosures:

The authors have no conflicts to disclose.


References:

1. Schwanhäusser, B., *et al.* Global quantification of mammalian gene expression control. *Nature*. **473** (7347), 337–342, doi: 10.1038/nature10098 (2011).
2. Boisvert, F. M., *et al.* A Quantitative Spatial Proteomics Analysis of Proteome Turnover in Human Cells. *Molecular & Cellular Proteomics*. **11** (3), M111.011429–M111.011429, doi: 10.1074/mcp.M111.011429 (2012).
3. Belle, A., Tanay, A., Bitincka, L., Shamir, R. & O'Shea, E. K. Quantification of protein half-lives in the budding yeast proteome. *Proceedings of the National Academy of Sciences of the United States of America*. **103** (35), 13004–13009, doi: 10.1073/pnas.0605420103 (2006).
4. Eden, E., *et al.* Proteome Half-Life Dynamics in Living Human Cells. *Science*. **331** (6018), 764–768, doi: 10.1126/science.1199784 (2011).
5. Parry, D. H. & O'Farrell, P. H. The schedule of destruction of three mitotic cyclins can dictate the timing of events during exit from mitosis. *Current Biology*. **11** (9), 671–683 (2001).
6. Holloway, S. L., Glotzer, M., King, R. W. & Murray, A. W. Anaphase is initiated by proteolysis rather than by the inactivation of maturation-promoting factor. *Cell*. **73** (7), 1393–1402 (1993).
7. Dharmasiri, N. & Estelle, M. Auxin signaling and regulated protein degradation. *Trends Plant Sci*. **9** (6), 302–308, doi: 10.1016/j.tplants.2004.04.003 (2004).
8. MacDonald, B. T., Tamai, K. & He, X. Wnt/beta-catenin signaling: components, mechanisms, and diseases. *Developmental Cell*. **17** (1), 9–26, doi: 10.1016/j.devcel.2009.06.016 (2009).
9. Chen, X., *et al.* Processing and turnover of the Hedgehog protein in the endoplasmic reticulum. *The Journal of Cell Biology*. **192** (5), 825–838, doi: 10.1083/jcb.201008090 (2011).
10. Elmore, S. Apoptosis: A Review of Programmed Cell Death. *Toxicologic Pathology*. **35** (4), 495–516, doi: 10.1080/01926230701320337 (2007).
11. Yi, J. J. & Ehlers, M. D. Emerging Roles for Ubiquitin and Protein Degradation in Neuronal Function. *Pharmacological Reviews*. **59** (1), 14–39, doi: 10.1124/pr.59.1.4 (2007).
12. Rubinsztein, D. C. The roles of intracellular protein-degradation pathways in neurodegeneration. *Nature*. **443** (7113), 780–786, doi: 10.1038/nature05291 (2006).
13. Müller, P. & Schier, A. F. Extracellular Movement of Signaling Molecules. *Developmental Cell*. **21** (1), 145–158, doi: 10.1016/j.devcel.2011.06.001 (2011).

14. Eldar, A., Rosin, D., Shilo, B.-Z. & Barkai, N. Self-enhanced ligand degradation underlies robustness of morphogen gradients. *Developmental Cell*. **5** (4), 635–646 (2003).
15. Zhou, P. Determining protein half-lives. *Methods In Molecular Biology (Clifton, N.J.)*. **284**, 67–77, doi: 10.1385/1-59259-816-1:067 (2004).
16. Woodside, K. H. Effects of cycloheximide on protein degradation and gluconeogenesis in the perfused rat liver. *Biochim Biophys Acta*. **421** (1), 70–79 (1976).
17. Schimke, R. T. & Doyle, D. Control of enzyme levels in animal tissues. *Annual Review of Biochemistry*. **39**, 929–976, doi: 10.1146/annurev.bi.39.070170.004433 (1970).
18. Kenney, F. T. Turnover of rat liver tyrosine transaminase: stabilization after inhibition of protein synthesis. *Science*. **156** (3774), 525–528 (1967).
19. Müller, P., *et al.* Differential Diffusivity of Nodal and Lefty Underlies a Reaction-Diffusion Patterning System. *Science*. **336** (6082), 721–724, doi: 10.1126/science.1221920 (2012).
20. Kiuchi, T., Nagai, T., Ohashi, K. & Mizuno, K. Measurements of spatiotemporal changes in G-actin concentration reveal its effect on stimulus-induced actin assembly and lamellipodium extension. *The Journal of Cell Biology*. **193** (2), 365–380, doi: 10.1083/jcb.201101035.dv (2011).
21. Plachta, N., Bollenbach, T., Pease, S., Fraser, S. E. & Pantazis, P. Oct4 kinetics predict cell lineage patterning in the early mammalian embryo. *Nature Cell Biology*. **13** (2), 117–123, doi: 10.1038/ncb2154 (2011).
22. Drocco, J. A., Grimm, O., Tank, D. W. & Wieschaus, E. Measurement and Perturbation of Morphogen Lifetime: Effects on Gradient Shape. *Biophys J*. **101** (8), 1807–1815, doi: 10.1016/j.bpj.2011.07.025 (2011).
23. Zhang, L., *et al.* Method for real-time monitoring of protein degradation at the single cell level. *BioTechniques*. **42** (4), 446–450, doi: 10.2144/000112453 (2007).
24. Miyawaki, A. Proteins on the move: insights gained from fluorescent protein technologies. *Nat Rev Mol Cell Biol*. **12** (10), 656–668, doi: 10.1038/nrm3199 (2011).
25. Pantazis, P. & Supatto, W. Advances in whole-embryo imaging: a quantitative transition is underway. *Nat Rev Mol Cell Biol*. **15** (5), 327–339, doi: 10.1038/nrm3786 (2014).
26. Lukyanov, K. A., Chudakov, D. M., Lukyanov, S. & Verkhusha, V. V. Innovation: Photoactivatable fluorescent proteins. *Nat Rev Mol Cell Biol*. **6** (11), 885–891, doi: 10.1038/nrm1741 (2005).
27. Gurskaya, N. G., *et al.* Engineering of a monomeric green-to-red photoactivatable fluorescent protein induced by blue light. *Nature Biotechnology*. **24** (4), 461–465, doi: 10.1038/nbt1191 (2006).
28. Zou, J. & Wei, X. Transplantation of GFP-expressing Blastomeres for Live Imaging of Retinal and Brain Development in Chimeric Zebrafish Embryos. *Journal of Visualized Experiments*. **41**, doi: 10.3791/1924 (2010).
29. Yuan, S. & Sun, Z. Microinjection of mRNA and Morpholino Antisense Oligonucleotides in Zebrafish Embryos. *Journal of Visualized Experiments*. **27**, doi: 10.3791/1113 (2009).
30. Rosen, J. N., Sweeney, M. F. & Mably, J. D. Microinjection of Zebrafish Embryos to Analyze Gene Function. *Journal of Visualized Experiments*. **25**, doi: 10.3791/1115 (2009).
31. Kimmel, C. B., Ballard, W. W., Kimmel, S. R., Ullmann, B. & Schilling, T. F. Stages of embryonic development of the zebrafish. *Developmental Dynamics*. **203** (3), 253–310, doi: 10.1002/aja.1002030302 (2004).

32. McKinney, S. A., Murphy, C. S., Hazelwood, K. L., Davidson, M. W. & Looger, L. L. A bright and photostable photoconvertible fluorescent protein. *Nature Methods*. **6** (2), 131–133, doi: 10.1038/nmeth.1296 (2009).
33. Dempsey, W. P., Qin, H. & Pantazis, P. In Vivo Cell Tracking Using PHOTO Zebrafish. *Methods in Molecular Biology*. **1148** (Chapter 14), 217–228, doi: 10.1007/978-1-4939-0470-9_14 (2014).
34. Schier, A. F. Nodal Morphogens. *Cold Spring Harbor Perspectives in Biology*. **1** (5), a003459–a003459, doi: 10.1101/cshperspect.a003459 (2009).
35. Otsu, N. A threshold selection method from gray-level histograms. *IEEE Trans Syst Man Cybern*. **9** (1), 62–66. doi: 10.1109/TSMC.1979.4310076 (1979).
36. Pédelacq, J.-D., Cabantous, S., Tran, T., Terwilliger, T. C. & Waldo, G. S. Engineering and characterization of a superfolder green fluorescent protein. *Nature Biotechnology*. **24** (1), 79–88, doi: 10.1038/nbt1172 (2005).
37. Swulius, M. T. & Jensen, G. J. The Helical MreB Cytoskeleton in Escherichia coli MC1000/pLE7 Is an Artifact of the N-Terminal Yellow Fluorescent Protein Tag. *Journal of Bacteriology*. **194** (23), 6382–6386, doi: 10.1128/JB.00505-12 (2012).
38. Landgraf, D., Okumus, B., Chien, P., Baker, T. A. & Paulsson, J. Segregation of molecules at cell division reveals native protein localization. *Nature Methods*. **9** (5), 480–482, doi: 10.1038/nmeth.1955 (2012).
39. Stadler, C., *et al.* Immunofluorescence and fluorescent-protein tagging show high correlation for protein localization in mammalian cells. *Nature Methods*. **10** (4), 315–323, doi: 10.1038/nmeth.2377 (2013).
40. Quattrocchio, F. M., Spelt, C. & Koes, R. Transgenes and protein localization: myths and legends. *Trends Plant Sci*. **18** (9), 473–476, doi: 10.1016/j.tplants.2013.07.003 (2013).
41. Morimoto, Y. V., Kojima, S., Namba, K. & Minamino, T. M153R Mutation in a pH-Sensitive Green Fluorescent Protein Stabilizes Its Fusion Proteins. *PLoS ONE*. **6** (5), e19598, doi: 10.1371/journal.pone.0019598.t002 (2011).
42. Shaner, N. C., Steinbach, P. A. & Tsien, R. Y. A guide to choosing fluorescent proteins. *Nature Methods*. **2** (12), 905–909, doi: 10.1038/nmeth819 (2005).
43. Waters, J. C. Accuracy and precision in quantitative fluorescence microscopy. *The Journal of Cell Biology*. **185** (7), 1135–1148, doi: 10.1021/ar040136s (2009).
44. Moll, U. M. & Petrenko, O. The MDM2-p53 interaction. *Mol Cancer Res*. **1** (14), 1001–1008 (2003).
45. Auer, T. O., Duroure, K., De Cian, A., Concordet, J. P. & Del Bene, F. Highly efficient CRISPR/Cas9-mediated knock-in in zebrafish by homology-independent DNA repair. *Genome Research*. **24** (1), 142–153, doi: 10.1101/gr.161638.113 (2014).
46. Bedell, V. M., *et al.* In vivo genome editing using a high-efficiency TALEN system. *Nature*. **490** (7422), 114–118, doi: 10.1038/nature11537 (2013).
47. Hruscha, A., *et al.* Efficient CRISPR/Cas9 genome editing with low off-target effects in zebrafish. *Development*. **140** (24), 4982–4987, doi: 10.1242/dev.099085 (2013).
48. Hwang, W. Y., *et al.* Heritable and Precise Zebrafish Genome Editing Using a CRISPR-Cas System. *PLoS ONE*. **8** (7), e68708, doi: 10.1371/journal.pone.0068708.s005 (2013).
49. Hwang, W. Y., *et al.* Efficient genome editing in zebrafish using a CRISPR-Cas system. *Nature Biotechnology*. **31** (3), 227–229, doi: 10.1038/nbt.2501 (2013).

50. Zu, Y., *et al.* TALEN-mediated precise genome modification by homologous recombination in zebrafish. *Nature Methods*. **10** (4), 329–331, doi: 10.1038/nmeth.2374 (2013).
51. Keller, P. J., Schmidt, A. D., Wittbrodt, J. & Stelzer, E. H. K. Reconstruction of Zebrafish Early Embryonic Development by Scanned Light Sheet Microscopy. *Science*. **322** (5904), 1065–1069, doi: 10.1126/science.1162493 (2008).
52. Blanchet, M. H., *et al.* Cripto Localizes Nodal at the Limiting Membrane of Early Endosomes. *Science Signaling*. **1** (45), ra13–ra13, doi: 10.1126/scisignal.1165027 (2008).
53. Jullien, J. & Gurdon, J. Morphogen gradient interpretation by a regulated trafficking step during ligand-receptor transduction. *Genes & Development*. **19** (22), 2682–2694, doi: 10.1101/gad.341605 (2005).
54. Incardona, J. P., *et al.* Receptor-mediated endocytosis of soluble and membrane-tethered Sonic hedgehog by Patched-1. *Proceedings of the National Academy of Sciences of the United States of America*. **97** (22), 12044–12049, doi: 10.1073/pnas.220251997 (2000).
55. Scholpp, S. & Brand, M. Endocytosis Controls Spreading and Effective Signaling Range of Fgf8 Protein. *Current Biology*. **14** (20), 1834–1841, doi: 10.1016/j.cub.2004.09.084 (2004).









Figure 3

[Click here to download Figure: FDAP_JoVE_figure3.pdf](#)







Name of Material / Equipment
PyFDAP (download from the following website: http://people.tuebingen.mpg.de/mueller-lab)
mMessage mMachine Sp6 Transcription Kit (Life Technologies, AM1340)
Alexa488-dextran conjugate, 3 kDa (Life Technologies, D34682)
6-well plastic dish (BD Falcon)
Embryo medium
Protease from <i>Streptomyces griseus</i> (Sigma, P5147)
5 cm diameter glass petri dish
200 ml glass beaker
Microinjection apparatus
Stereomicroscope
1x Danieau's medium
UltraPure low melting point agarose (Invitrogen, 16520-100)
Glass Pasteur pipette (Kimble Chase (via Fisher), 63A53WT)
Metal probe
Glass bottom dishes (MatTek, P35G-1.5-14-C)
15 ml tube (BD Falcon) filled with ~5 ml embryo medium
Inverted laser scanning confocal microscope
Heated stage
Confocal software capable of time-lapse imaging
25x or 40x water objective
10x air objective
Immersion oil

Comments / Description

Install and operate using the instructions provided on the PyFDAP website; PyFDAP is compatible with Linux, Mac, and Windows operating

To generate capped mRNA for injection into embryos

Co-inject with mRNA to create intracellular and extracellular masks

Incubate embryos in agarose-coated wells until ready for mounting

250 mg/l Instant Ocean salt, 1 mg/l methylene blue in reverse osmosis water adjusted to pH 7 with NaHCO_3

Make a 5 mg/ml stock and use at 1 mg/ml to dechorionate embryos at the one-cell stage

For embryo dechoronation

For embryo dechoronation

For injection of mRNA and dye into embryos at the one-cell stage

For injecting and mounting embryos

Dilute low melting point agarose and perform imaging in this medium; recipe: 0.2 mM filtered solution of 58 mM NaCl, 0.7 mM KCl, 0.4 mM MgCl_2

For mounting embryos; use at a concentration of 1% in Danieau's medium: add 200 mg to 20 ml Danieau's medium, microwave until dissolve

For mounting embryos; flame the tip to prevent jagged edges from injuring embryos

For positioning embryos during mounting

Use the appropriate cover glass thickness for your objective; part number listed here is for cover glass No. 1.5

For rinsing residual agarose from the Pasteur pipette

A mercury arc lamp, 488 nm laser, 543 nm laser, and the appropriate filter sets are required

To maintain embryos at the optimal temperature of 28 °C during the experiment

Must be able to define multiple positions and automatically image them at defined intervals

Objective for imaging

Objective for photoconversion

Immersion oil with the same refractive index as water

1gSO₄, 0.3 mM CaCl₂, 5 mM HEPES pH 7.2

d, then aliquot 1 ml into microcentrifuge tubes; aliquots can be stored at 4 °C, re-melted at 70 °C, and cooled to 40 - 42 °C when ready

/ to use



1 Alewife Center #200
Cambridge, MA 02140
tel. 617.945.9051
www.jove.com

ARTICLE AND VIDEO LICENSE AGREEMENT

Title of Article:

Measuring protein stability in living zebrafish embryos using Fluorescence Decrease After Photoconversion (FDAP)

Author(s):

Katherine W. Rogers, Alexander Bläßle, Alexander F. Schier, and Patrick Müller

Item 1 (check one box): The Author elects to have the Materials be made available (as described at

<http://www.jove.com/publish>) via: ☒ Standard Access ☐ Open Access

Item 2 (check one box):

- ☒ The Author is NOT a United States government employee.
- ☐ The Author is a United States government employee and the Materials were prepared in the course of his or her duties as a United States government employee.
- ☐ The Author is a United States government employee but the Materials were NOT prepared in the course of his or her duties as a United States government employee.

ARTICLE AND VIDEO LICENSE AGREEMENT

1. Defined Terms. As used in this Article and Video License Agreement, the following terms shall have the following meanings: “**Agreement**” means this Article and Video License Agreement; “**Article**” means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein; “**Author**” means the author who is a signatory to this Agreement; “**Collective Work**” means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; “**CRC License**” means the Creative Commons Attribution-Non Commercial-No Derivs 3.0 Unported Agreement, the terms and conditions of which can be found at: <http://creativecommons.org/licenses/by-nc-nd/3.0/legalcode>; “**Derivative Work**” means a work based upon the Materials or upon the Materials and other pre-existing works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; “**Institution**” means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; “**JoVE**” means MyJoVE Corporation, a Massachusetts corporation and the publisher of *The Journal of Visualized Experiments*; “**Materials**” means the Article and / or the Video; “**Parties**” means the Author and JoVE; “**Video**” means any video(s) made by the Author, alone or in conjunction with any other parties, or by JoVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion of the Article, and in which the Author may or may not appear.

2. Background. The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JoVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.

3. Grant of Rights in Article. In consideration of JoVE agreeing to publish the Article, the Author hereby grants to JoVE, subject to **Sections 4** and **7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the “Open Access” box has been checked in **Item 1** above, JoVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.

ARTICLE AND VIDEO LICENSE AGREEMENT

4. Retention of Rights in Article. Notwithstanding the exclusive license granted to JoVE in **Section 3** above, the Author shall, with respect to the Article, retain the non-exclusive right to use all or part of the Article for the non-commercial purpose of giving lectures, presentations or teaching classes, and to post a copy of the Article on the Institution's website or the Author's personal website, in each case provided that a link to the Article on the JoVE website is provided and notice of JoVE's copyright in the Article is included. All non-copyright intellectual property rights in and to the Article, such as patent rights, shall remain with the Author.

5. Grant of Rights in Video – Standard Access. This **Section 5** applies if the "Standard Access" box has been checked in **Item 1** above or if no box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby acknowledges and agrees that, Subject to **Section 7** below, JoVE is and shall be the sole and exclusive owner of all rights of any nature, including, without limitation, all copyrights, in and to the Video. To the extent that, by law, the Author is deemed, now or at any time in the future, to have any rights of any nature in or to the Video, the Author hereby disclaims all such rights and transfers all such rights to JoVE.

6. Grant of Rights in Video – Open Access. This **Section 6** applies only if the "Open Access" box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby grants to JoVE, subject to **Section 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Video in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Video into other languages, create adaptations, summaries or extracts of the Video or other Derivative Works or Collective Works based on all or any portion of the Video and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. For any Video to which this Section 6 is applicable, JoVE and the Author hereby grant to the public all such rights in the Video as provided in, but subject to all limitations and requirements set forth in, the CRC License.

7. Government Employees. If the Author is a United States government employee and the Article was prepared in the course of his or her duties as a United States government employee, as indicated in **Item 2** above, and any of the licenses or grants granted by the Author hereunder exceed the scope of the 17 U.S.C. 403, then the rights granted hereunder shall be limited to the maximum rights permitted under such

statute. In such case, all provisions contained herein that are not in conflict with such statute shall remain in full force and effect, and all provisions contained herein that do so conflict shall be deemed to be amended so as to provide to JoVE the maximum rights permissible within such statute.

8. Likeness, Privacy, Personality. The Author hereby grants JoVE the right to use the Author's name, voice, likeness, picture, photograph, image, biography and performance in any way, commercial or otherwise, in connection with the Materials and the sale, promotion and distribution thereof. The Author hereby waives any and all rights he or she may have, relating to his or her appearance in the Video or otherwise relating to the Materials, under all applicable privacy, likeness, personality or similar laws.

9. Author Warranties. The Author represents and warrants that the Article is original, that it has not been published, that the copyright interest is owned by the Author (or, if more than one author is listed at the beginning of this Agreement, by such authors collectively) and has not been assigned, licensed, or otherwise transferred to any other party. The Author represents and warrants that the author(s) listed at the top of this Agreement are the only authors of the Materials. If more than one author is listed at the top of this Agreement and if any such author has not entered into a separate Article and Video License Agreement with JoVE relating to the Materials, the Author represents and warrants that the Author has been authorized by each of the other such authors to execute this Agreement on his or her behalf and to bind him or her with respect to the terms of this Agreement as if each of them had been a party hereto as an Author. The Author warrants that the use, reproduction, distribution, public or private performance or display, and/or modification of all or any portion of the Materials does not and will not violate, infringe and/or misappropriate the patent, trademark, intellectual property or other rights of any third party. The Author represents and warrants that it has and will continue to comply with all government, institutional and other regulations, including, without limitation all institutional, laboratory, hospital, ethical, human and animal treatment, privacy, and all other rules, regulations, laws, procedures or guidelines, applicable to the Materials, and that all research involving human and animal subjects has been approved by the Author's relevant institutional review board.

10. JoVE Discretion. If the Author requests the assistance of JoVE in producing the Video in the Author's facility, the Author shall ensure that the presence of JoVE employees, agents or independent contractors is in accordance with the relevant regulations of the Author's institution. If more than one author is listed at the beginning of this Agreement, JoVE may, in its sole discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JoVE reserves the right, in its absolute and sole discretion and without giving any reason therefore, to accept or decline any work submitted to JoVE. JoVE and its employees, agents and independent contractors shall have

ARTICLE AND VIDEO LICENSE AGREEMENT

full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JoVE has sole discretion as to the method of making and publishing the Materials, including, without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

11. **Indemnification.** The Author agrees to indemnify JoVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JoVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JoVE, making of videos by JoVE, or publication in JoVE or elsewhere by JoVE. The Author shall be responsible for, and shall hold JoVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to the making of a video by JoVE its employees, agents or independent contractors. All sterilization, cleanliness or decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's


expense. All indemnifications provided herein shall include JoVE's attorney's fees and costs related to said losses or damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JoVE, its employees, agents or independent contractors.

12. **Fees.** To cover the cost incurred for publication, JoVE must receive payment before production and publication the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JoVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.

13. **Transfer, Governing Law.** This Agreement may be assigned by JoVE and shall inure to the benefits of any of JoVE's successors and assignees. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to be one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement required per submission.

CORRESPONDING AUTHOR:

Name:	Patrick Müller	
Department:	Systems Biology of Development Group	
Institution:	Friedrich Miescher Laboratory of the Max Planck Society	
Article Title:	Measuring protein stability in living zebrafish embryos using Fluorescence Decrease After Photoconversion (FDAP)	
Signature:		Date: April 30, 2014

Please submit a signed and dated copy of this license by one of the following three methods:

- 1) Upload a scanned copy of the document as a pdf on the JoVE submission site;
- 2) Fax the document to +1.866.381.2236;
- 3) Mail the document to JoVE / Attn: JoVE Editorial / 1 Alewife Center #200 / Cambridge, MA 02139

For questions, please email submissions@jove.com or call +1.617.945.9051



Dr. Nam Nguyen
JoVE Editorial Board
1 Alewife Center, Suite 200
Cambridge, MA 02140

Dr. Patrick Müller
Research Group Leader
Systems Biology of
Development
Tel.: +49-(0)7071-601-815
Fax: +49-(0)7071-601-801
pmueller@tuebingen.mpg.de

FDAP invited manuscript

June 24, 2014

Dear Dr. Nguyen,

Thank you for your editorial comments on our manuscript “Measuring protein stability in living zebrafish embryos using Fluorescence Decay After Photoconversion (FDAP)”. Below we address your comments.

Editorial comments

1. Please take this opportunity to thoroughly proofread the manuscript to ensure that there are no spelling or grammar issues. The JoVE editor will not copy-edit your manuscript and any errors in the submitted revision may be present in the published version.

2. Please ensure that **all text** in the protocol section is written in the imperative tense as if telling someone how to do the technique (e.g., “Do this,” “Ensure that,” etc.). The actions should be described in the imperative tense in complete sentences wherever possible. Avoid usage of phrases such as “could be,” “should be,” and “would be” throughout the Protocol. Any text that cannot be written in the imperative tense may be added as a “Note.” However, notes should be concise and used sparingly.

We have changed several sentences in Section 5 to the imperative mood. We also eliminated “should” and “would” where possible and removed the note at the beginning of Section 5.

For example: Step 5 could be moved to the Representative Results.

Multiple reviewers requested further explanation of control experiments, and reviewer #3 specifically asked for the control experiments to appear in the numbered protocol. We therefore created Section 5 to address these important concerns.

3. The Protocol should contain only action items that direct the reader to do something. Please move the discussion about the protocol to the Discussion.

We have changed several sentences to the imperative mood (see response to point #2 above) and include brief notes to clarify subsequent instructions.

4. The Figure Legends should include a title and a short description of the data presented in the Figure and relevant symbols. The Discussion of the Figures should be placed in the Representative Results. Details of the methodology should not be in the Figure Legends, but rather the Protocol.

We have shortened the figure legends where possible.

Yours sincerely,

Patrick Müller