

Journal of Visualized Experiments

Sedimentation Equilibrium of a Small Oligomer-forming Membrane Protein: effect of Histidine Protonation on Pentameric Stability.

--Manuscript Draft--

Manuscript Number:	JoVE52404R4
Full Title:	Sedimentation Equilibrium of a Small Oligomer-forming Membrane Protein: effect of Histidine Protonation on Pentameric Stability.
Article Type:	Methods Article - JoVE Produced Video
Keywords:	Analytical ultracentrifugation; sedimentation equilibrium; molecular weight; membrane proteins; small hydrophobic; respiratory syncytial virus; detergent; density matching; oligomeric size; histidine protonation; oligomer stability
Manuscript Classifications:	3.2.782.580.600: Paramyxoviridae Infections; 3.8: Respiratory Tract Diseases; 4.12.125.142.308: Histidine; 4.12.776.543: Membrane Proteins; 4.12.776.543.550.425: Ion Channels; 4.27.720.877.265: Detergents; 5.5.181.724: Ultracentrifugation; 5.5.196.941.336: Centrifugation, Density Gradient; 8.1.158.344: Biophysics
Corresponding Author:	Jaume Torres Nanyang Technological University Singapore, Singapore SINGAPORE
Corresponding Author Secondary Information:	
Corresponding Author E-Mail:	jtorres@ntu.edu.sg
Corresponding Author's Institution:	Nanyang Technological University
Corresponding Author's Secondary Institution:	
First Author:	Wahyu Surya
First Author Secondary Information:	
Other Authors:	Wahyu Surya
Order of Authors Secondary Information:	
Abstract:	<p>Analytical ultracentrifugation (AUC) can be used to study reversible interactions between macromolecules over a wide range of interaction strengths and under physiological conditions. This makes AUC a method of choice to quantitatively assess stoichiometry and thermodynamics of homo- and hetero-association that are transient and reversible in biochemical processes. In the modality of sedimentation equilibrium (SE), a balance between diffusion and sedimentation provides a profile as a function of radial distance that depends on a specific association model. Herein, a detailed SE protocol is described to determine the size and monomer-monomer association energy of a small membrane protein oligomer using an analytical ultracentrifuge. AUC-ES is label-free, only based on physical principles, and can be used on both water soluble and membrane proteins. An example is shown of the latter, the small hydrophobic (SH) protein in the human respiratory syncytial virus (hRSV), a 65-amino acid polypeptide with a single α-helical transmembrane (TM) domain that forms pentameric ion channels. NMR-based structural data shows that SH protein has two protonatable His residues in its transmembrane domain that are oriented facing the lumen of the channel. SE experiments have been designed to determine how pH affects the association constant and the oligomeric size of SH protein. While the pentameric form was preserved in all cases, its association constant was reduced at low pH. These data are in agreement with a similar pH dependency observed for SH channel activity, consistent with a luminal orientation of the two His residues in SH protein. The latter may experience electrostatic repulsion and reduced oligomer stability at low pH. In summary, this method is applicable whenever quantitative information on subtle protein-protein association changes in physiological conditions have to be measured.</p>
Author Comments:	

Additional Information:	
Question	Response
If this article needs to be "in-press" by a certain date to satisfy grant requirements, please indicate the date below and explain in your cover letter.	Aug 01, 2014
If this article needs to be filmed by a certain date to due to author/equipment/lab availability, please indicate the date below and explain in your cover letter.	

Sedimentation Equilibrium of a Small Oligomer-forming Membrane Protein: effect of Histidine Protonation on Pentameric Stability.

AUTHORS:

SURYA, Wahyu
School of Biological Sciences
Nanyang Technological University
Singapore, Singapore
wsurya1@e.ntu.edu.sg

TORRES, Jaume
School of Biological Sciences
Nanyang Technological University
Singapore, Singapore
jtorres@ntu.edu.sg

CORRESPONDING AUTHOR:

TORRES, Jaume
School of Biological Sciences
Nanyang Technological University
Singapore, Singapore
jtorres@ntu.edu.sg

KEYWORDS:

Analytical ultracentrifugation, sedimentation equilibrium, molecular weight, membrane proteins, small hydrophobic, respiratory syncytial virus, detergent, density matching, oligomeric size, histidine protonation, oligomer stability

SHORT ABSTRACT:

Sedimentation equilibrium (SE) can be used to study protein-protein interactions in a physiological environment. This manuscript describes the use of this technique to determine the effect of pH on the stability of a homo-pentamer formed by the small hydrophobic (SH) protein encoded by the human syncytial respiratory virus (hRSV).

LONG ABSTRACT:

Analytical ultracentrifugation (AUC) can be used to study reversible interactions between macromolecules over a wide range of interaction strengths and under physiological conditions. This makes AUC a method of choice to quantitatively assess stoichiometry and thermodynamics of homo- and hetero-association that are transient and reversible in biochemical processes. In the modality of sedimentation equilibrium (SE), a balance between diffusion and sedimentation provides a profile as a function of radial distance that depends on a specific association model.

Herein, a detailed SE protocol is described to determine the size and monomer-monomer association energy of a small membrane protein oligomer using an analytical ultracentrifuge. AUC-ES is label-free, only based on physical principles, and can be used on both water soluble and membrane proteins. An example is shown of the latter, the small hydrophobic (SH) protein in the human respiratory syncytial virus (hRSV), a 65-amino acid polypeptide with a single α -helical transmembrane (TM) domain that forms pentameric ion channels. NMR-based structural data shows that SH protein has two protonatable His residues in its transmembrane domain that are oriented facing the lumen of the channel. SE experiments have been designed to determine how pH affects association constant and the oligomeric size of SH protein. While the pentameric form was preserved in all cases, its association constant was reduced at low pH. These data are in agreement with a similar pH dependency observed for SH channel activity, consistent with a luminal orientation of the two His residues in SH protein. The latter may experience electrostatic repulsion and reduced oligomer stability at low pH. In summary, this method is applicable whenever quantitative information on subtle protein-protein association changes in physiological conditions have to be measured.

INTRODUCTION:

Analytical ultracentrifugation¹⁻⁵ is one of the most important methods to study interactions of macromolecules under physiological conditions, being accessible to both weak and strong interactions. The method is label-free and uses light absorption or interference, and even fluorescence optical systems can be used to access concentration ranges over several orders of magnitude⁶.

This method is especially useful since most biochemical processes depend on reversible interactions. The stoichiometry and strength of these interactions have to be quantitatively characterized to understand biological processes, and a number of methods exist for this purpose^{7,8}. However, transient interactions are difficult to study⁹.

The choice of a method to characterize macromolecular interactions depends on its static or dynamic nature. In the first case, sedimentation velocity (SV) is used, where the rate of radial transport is measured and complexes are fractionated on the basis of differences in buoyant mass and shape.

In contrast, dynamic associations that are reversible on the time scale of the experiment cannot be physically separated. In this case, self- or hetero-interactions leading to non-covalent interactions are in an equilibrium that depends on the total protein concentration. These dynamic interactions can be studied by both sedimentation equilibrium (SE) and sedimentation velocity (SV)¹⁰. However, the first method is simpler to perform and is described here. In SE, centrifugation is performed at a sufficiently low speed so that an equilibrium is reached between diffusion and sedimentation. At this point, the equilibrium profile of an optical signal (UV-VIS) as a function of radial distance, can be analyzed using pre-set thermodynamic models for associations¹¹.

In the present paper, a sedimentation equilibrium study is presented of the self-association of a viral membrane protein that forms ion channels. Because of its hydrophobicity, the experiment is run in presence of detergent, and in this case the density of solvent has to be matched to that of the detergent. However, the protocol described would be identical in the case of a water soluble protein, except that no solvent density matching would be required.

The protein used is encoded in the human respiratory syncytial virus (hRSV), an enveloped pneumovirus in the paramyxoviridae family that causes lower respiratory tract disease in infants, elderly and immunocompromised populations worldwide¹². Up to 64 million reported cases of hRSV infection and 160,000 deaths occur each year.

The hRSV genome transcribes 11 proteins, including the three membrane proteins F, G, and small hydrophobic (SH). SH protein is involved in the pathogenesis of RSV infection. RSV lacking the SH gene (RSV Δ SH) was viable, caused formation of syncytia and grew as well as the wild-type (WT) virus¹³⁻¹⁶. However, RSV Δ SH virus replicated 10-fold less efficiently than the WT in the upper respiratory tract^{15,16}. Also, RSV Δ SH virus was attenuated in *in vivo* mouse and chimpanzee models^{13,17}.

The SH protein is a 64 (RSV subgroup A) or 65 (RSV subgroup B) amino acids long type II integral membrane protein that accumulates mostly at the membranes of the Golgi compartment¹⁸. SH protein has a single predicted α -helical transmembrane (TM) domain¹⁹ which is highly conserved^{20,21}. The C- and N-terminal extramembrane domains are oriented lumenally/extracellularly and cytoplasmically, respectively.

Both synthetic TM domain (residues 18-43) and full length SH protein have been shown to form homopentamers in a variety of detergents. The homopentameric form is responsible for channel activity in planar lipid bilayers^{22,23}. The correct orientation of the TM monomers in the lipid bilayer was first determined using site specific infrared dichroism²³, which showed His-22 to be in a luminal, close to inter-helical, orientation. The same TM domain orientation was confirmed by NMR studies that reconstructed the pentameric α -helical bundle of the full-length protein in dodecylphosphocholine (DPC) micelles²². In this 'micelle' model, a single α -helical TM domain was flanked N-terminally by an α -helix, and C-terminally by an extended β -hairpin. The two protonatable residues of SH protein, His-22 and His-51, are located in the TM domain (lumenally oriented), and at the tip of the extramembrane C-terminal β hairpin (far from the channel pore), respectively. In a bicellar environment, however, the TM α -helix extends up to His-51, and both His residues are accessible to the lumen of the channel²⁴. The channel structure adopts a funnel-like architecture²², where the narrower region (Ser-29 to Cys-45)²² is lined with hydrophobic side chains (Ile-32, Ile-36, Ile-40 and Leu-44), and Ile-36 defines the narrowest point in the channel lumen. His-22 is located at the largest opening of this funnel, whereas His-51 is at the tip of the smallest opening.

In the present paper, analytical centrifugation in a sedimentation equilibrium mode has been used to determine if His protonation affects the stability of the SH protein pentamer. In this case, SH protein was solubilized in C14-betaine detergent, which has been used previously to

show that SH protein forms pentameric oligomers²².

PROTOCOL:

This protocol is based on the following resources, which are to be referred for more details and special considerations^{3,25-28}.

1. Density matching of detergent micelles with $^2\text{H}_2\text{O}$

Note: The density of the buffer solution needs to be matched to the density of the detergent micelles. Common density-adjusting agents include $^2\text{H}_2\text{O}$, H_2^{18}O , $^2\text{H}_2^{18}\text{O}$, glycerol and sucrose²⁹. H_2^{18}O has the same density as $^2\text{H}_2\text{O}$ and may be a better choice if deuteration of exchangeable protons in the protein is not desired. In this procedure, the density of 3-(N,N-dimethylmyristylammonio)-propanesulfonate (C14SB) detergent in 50 mM Tris pH 7.3, 100 mM NaCl will be matched with $^2\text{H}_2\text{O}$. As an initial guess the following concentrations of $^2\text{H}_2\text{O}$ will be used: 10, 30, and 50% v/v.

1.1. Sample preparation

1.1.1. Prepare the following stock solutions and filter sterilize through a 0.2 μm syringe filter: 50 mL 500 mM Tris pH 7.3 and 1 M NaCl (10X buffer solution); 1 mL 250 mM C14SB (50X detergent solution).

1.1.2. Prepare 200 μL of **sample** solution by mixing 20 μL 10X buffer solution, 4 μL 50X detergent solution, 20 μL $^2\text{H}_2\text{O}$ (99.9%), and 156 μL deionized H_2O . Prepare also 200 μL of **reference** solution by mixing 20 μL 10X buffer solution, 20 μL $^2\text{H}_2\text{O}$ (99.9%), and 160 μL deionized H_2O .

1.1.3. Repeat step 1.1.2 for the other $^2\text{H}_2\text{O}$ concentrations, i.e. 30% and 50%, adjusting the $^2\text{H}_2\text{O}$ and H_2O amounts appropriately.

1.2. Assembly of 6-channel AUC cells and sample loading into the cells

Note: There are two types of AUC cell depending on the sample loading method. Cells without external fill has to be loaded prior to sealing the cell, whereas external-fill cells can be loaded after the cells are sealed. Assembly of an external-fill AUC cell has been described previously³. In this protocol, the assembly of a **6-channel AUC cell without external fill** is described. The main difference is that it has screw rings on both sides which need to be tightened separately, and it doesn't need housing plugs (Fig. 1). The difference in assembly steps are highlighted below.

[Place Figure 1 here]

1.2.1. Prepare two window assemblies for each AUC cell with **sapphire window** instead of

quartz window (Fig. 2). Place the window gasket into the window holder. Slightly bend the window liner and place it into the window holder such that the gap is formed opposite to the window holder keyway. Place the sapphire window inside the window liner, aligning the mark with the window holder keyway.

Note: The quartz window is compressible and thus will produce more light refraction at high speed^{28,30}. Therefore for interference measurements above 30000 rpm, such as in this density matching experiment, sapphire windows are used. A sapphire window is heavier than quartz window and has an "X" etched onto its side.

[Place Figure 2 here]

1.2.2. Place the cell housing with the **part number upside-down**. With the keyways aligned with the housing key, slide into the cell housing **firstly a 6-sector centerpiece with beveled side down, followed by one window assembly** with the window facing down (Fig. 1, left).

1.2.3. Lightly coat the screw ring threads and screw ring washer with spinkote. Place a screw ring washer on top of the window assembly. Install the screw ring into the window housing with the word "OUT" facing outside. Hand-tighten the screw ring by using the cell aligning tool.

1.2.4. Using the torque wrench, tighten the screw ring to **only 60 inch-pounds**.

1.2.5. Place the cell with the **part number upright and positioned at 12-noon**. Load 120 μ l reference into the left rows and 110 μ l sample into the right rows. Ensure that each sample and reference is correctly paired. Note: Exact sample volume is not critical, but the reference needs to have slightly more volume than the sample (5-10 μ l) so that the sample meniscus will be distinct.

1.2.6. Carefully slide into the cell housing one window assembly with the window facing down (Fig. 1, right). Take care not to disturb the cell excessively and spill the contents.

1.2.7. **Repeat step 1.2.33 and 1.2.44**, tightening the second screw ring to **120 inch-pounds**. Invert the cell and re-tighten the first screw ring to **120 inch-pounds**.

1.2.8. Load the cells into the rotor, install the rotor into the centrifuge and install the monochromator according to the manufacturer's instructions²⁸. Note: Details on this step can also be found in this reference³.

1.3. Setting up interference measurement

1.3.1. Start up the user interface software for the AUC instrument and perform laser setup and radial calibration for each cell at 3000 rpm, according to the manufacturer's instructions, which will be briefly summarized in the following step.

1.3.1.1 After sufficient vacuum has been reached (< 100 microns), run the centrifuge at 3000 rpm. Preview the interference pattern in the user interface software and adjust the laser parameters to obtain the highest contrast.

1.3.2. Create a new set up file (File | New File) specifying “Equilibrium” and “Interference” measurement. Set up a sedimentation equilibrium method (“Method” button) to run at 45000 rpm or the highest speed anticipated for protein samples, whichever is higher, with run temperature at 20°C and collect 1 scan every 15 minutes. Monitor equilibrium progress by opening the data files in HeteroAnalysis and selecting “Match” function after at least 12 hours (approximately overnight, Fig. 3).

Note: a similar function is also available in SEDFIT (Options | Loading Options | Test Approach to Equilibrium).

[Place Figure 3 here]

1.4. Data analysis

1.4.1. For each set of samples, plot the slope of the radial distribution profile against the $^2\text{H}_2\text{O}$ concentration. Note: The distribution will be a very shallow exponential that approaches linearity. X-axis intercept corresponds to the matching $^2\text{H}_2\text{O}$ concentration.

1.4.2. For more accurate results, perform the experiment in several replicates. Alternatively, repeat the experiment with a narrower range of $^2\text{H}_2\text{O}$ concentration.

2. Sedimentation equilibrium of SH in C14SB micelles

2.1. Run parameters

2.1.1. Calculate buffer density and viscosity, protein partial specific volume and centrifugation speed by using SEDNTERP. To calculate buffer density and viscosity, select Compute in “Buffer Data Select” section and enter the buffer components accordingly, including the D_2O concentration.

2.1.1.1 To calculate protein partial specific volume, select Compute in “V-bar” section and enter the protein amino acid sequence. Specify the highest expected oligomeric size in “Make an oligomer from this monomer: N=” field, in this case $N = 5$. Calculate the speed by entering values in the RPM field on the main window until $\sigma \approx 1$; this is a rule of thumb to ensure a good exponential shape of the radial distribution profile³¹.

Note: The values calculated for this experiment was as follows: $\rho = 1.03839 \text{ g/ml}$, $\eta = 1.0267 \text{ cP}$, $\bar{v} = 0.7569 \text{ ml/g}$, $\omega_1 = 16000 \text{ rpm}$.

2.1.2. Calculate subsequent speeds to follow $\omega_n^2/\omega_{n-1}^2 \approx 1.5$ to ensure sufficient difference

between the distribution profile at one speed and the next³¹. Note: this can also be done from the function "Estimate equilibrium rotor speeds function" in SEDFIT, which takes into account the solution column (filling volume).

2.2. Sample preparations

2.2.1. Prepare 1 ml **reference** solution with 5 mM C14SB and 32.3% ²H₂O as determined from density matching experiment (section 1), by mixing 100 µl 10X buffer solution (step 1.1.1), 20 µl 50X detergent solution (step 1.1.1), 323 µl ²H₂O (99.9%) and 527 µl deionized H₂O.

2.2.2. Dissolve lyophilized, HPLC-purified SH peptides (expression and purification described previously³²) in appropriate solvent such as methanol or 50% v/v aqueous acetonitrile. Measure A₂₈₀ of the dissolved peptides in a microlitre-scale UV/Vis spectrophotometer and aliquot for three samples to give A_{280, 12mm} = 0.3, 0.5, and 0.8 (A_{280, 10mm} = 0.25, 0.417, and 0.67) each when diluted to 130 µl. Lyophilize the samples overnight and resuspend in 130 µl reference solution (step 2.2.1) to give the **sample** solutions.

Note: SH protein can be detected from UV/Vis absorbance at 280 nm because it contains Trp and Tyr residues. Proteins without aromatic residues can be detected by tagging them with a suitable chromophore, using a Trp-containing mutant, or by using interference measurements instead of absorbance.

2.2.3. Follow steps in section 1.2 to assemble a 6-channel AUC cell **with quartz windows**. Load the highest concentration sample (A_{280, 12mm} = 0.8) in the channel nearest to the rotor center and lowest concentration sample (A_{280, 12mm} = 0.3) furthest from the rotor center.

2.3. Setting up absorbance measurements

2.3.1. Create a new set up file (File | New File) specifying "Equilibrium" and "Absorbance" measurement. Specify 280 nm as the detector wavelength.

2.3.2. Perform radial calibration at 3000 rpm by checking "Radial calibration before first scan" in Scan Options, specifying data collection at low resolution, e.g. with Radial step size = 0.01 cm, Replicates = 3 (low resolution, fast), and executing a Single Scan. After the scan is completed, uncheck the option.

2.3.3. Set up a sedimentation equilibrium method ("Method" button) to run at the first speed calculated in step 2.1.2, at 20°C, and collect 1 scan every 30 minutes. In each cells' "Detail", specify data collection at low resolution as in step 2.3.2. Monitor equilibrium progress by opening the data files in HeteroAnalysis and selecting "Match" function after at least 18 hours (approximately overnight, Fig. 4).

Note: attainment of equilibrium may take substantially longer time for the first speed, whereas subsequent speeds will take less time.

[Place Figure 4 here]

2.3.4. Once equilibrium has been reached, collect a single scan at high resolution, e.g. with Radial step size = 0.001 cm, Replicates = 10 (high resolution, slow).

2.3.5. After the scan is completed, repeat step 2.3.3 and 2.3.4 for the next speed.

2.3.6. Optionally, when the time required to reach equilibrium for each speed is known (calculated or from experience), set up the sedimentation equilibrium method scan to include all the speeds calculated in step 2.1.2 and to collect 1 scan after the equilibration time for each speed. In this case, specify high resolution data collection in each cells' "Detail".

2.4. Data analysis in SEDFIT and SEDPHAT

Note: For further details and considerations in data analysis reader is referred to the following website: www.analyticalultracentrifugation.com.

2.4.1. Open high resolution scans in SEDFIT (Data | load sedimentation equilibrium data) and split the data into 3 channels (these correspond to different concentrations for each sample; Options | Loading Options | Save 6-channel Raw Data in 3 Subsets).

2.4.2. Re-open data files that belong to the same sample and same concentration but different speeds in SEDFIT. Adjust the meniscus (vertical red line), cell bottom (vertical blue line) and fitting limits (vertical green lines), and export the data for use in SEDPHAT (Data | Export Data to SEDPHAT). Input the parameters calculated in step 2.1.1, as well as rotor type and centerpiece type as requested. Repeat this step for every sample and concentration.

2.4.3. Open all data from the same sample (all concentrations and speeds) in SEDPHAT and fill in Experiment Parameters; an example is shown in Fig. 5.

Note: when D₂O is added in the buffer, deuteration of exchangeable protons could significantly alter the molecular weight of the protein, especially for water-soluble proteins. Membrane proteins, especially small ones like SH protein, are less affected because membrane-embedded regions are protected from exchange. To correct this, input the "buffer D mol fraction".

Note: At this step it is recommended to save the edited dataset separately by selecting menu Data | Copy All Data And Save As New Config.

2.4.4. Select a Model and fill in Global Parameters for that model. Note: As an example, the "Monomer-n-Mer Self-Association" model and its parameters are shown in Fig. 6.

[Place Figure 5 here]

[Place Figure 6 here]

2.4.5. Run a Global Fit by selecting menu Fit | Global Fit and wait until the fit converges. Note down (or take a screenshot of) the fitting results, especially the global reduced chi-square and log K_a values. Extract other data such as fit data and fitting residuals from menu Copy and Display | Display Thermodynamic Information.

2.4.6. Return to Global Parameters and check M(1) to fit the monomer molecular weight and repeat step 2.4.5. Note down the fitted molecular weight and global reduced chi-square.

2.4.7. Repeat step 2.4.4 to 2.4.6 for each model to be tested and compare the fit quality of each model by comparing the reduced global chi-square value as well as fitting residuals.

Note: Small and random fitting residuals generally indicates a good fit, and the model that fits best would have the smallest global reduced chi-square. The fitted monomer molecular weight and its chi-square value should not differ substantially from that of the fixed (theoretical) molecular weight.

2.4.8. Calculate the confidence interval for the obtained log K_a by value by first selecting Statistics | Critical chi-square for error surface projections and inputting the desired confidence interval. Next, go to Statistics | Generate 1-dimensional error surface projection and unselect log K_a in Global Parameters dialog to obtain the chi-square values for log K_a .

Note: readers are advised to consult the following sources (<http://www.analyticalultracentrifugation.com/sedphat/statistics.htm>) for more details on the method ³³ as well as illustration of this method ³⁴.

REPRESENTATIVE RESULTS:

The radial distribution profile of C14SB detergent micelles in 50 mM Tris, 100 mM NaCl pH 7.3 forms a very shallow exponential that could be fitted to a linear model (Fig. 7A). The slope of this distribution is inversely correlated to D₂O concentration (Fig. 7B). The point where the slope is zero, i.e. the matching D₂O concentration, was found to be 32.3%.

[Place Figure 7 here]

The same experiment was repeated for different buffers: 50 mM sodium phosphate, 100 mM NaCl pH 5.5 and 50 mM phosphate-citrate, 100 mM NaCl pH 3 to obtain matching D₂O concentrations of 30.3% and 41.0%, respectively.

Samples of SH wild type (wt) in detergent were exposed to pH 3, 5.5 and 7.3 (total 6 samples), followed by centrifugation at 15,000, 19,000, 23,000, 28,000, 34,000 and 42,000 rpm. Data obtained at lower speed could not be reliably fitted, possibly because equilibrium had not been attained, therefore data obtained from the four higher speeds were used. At pH 7.3, SH WT was

found to form pentamers (Fig. 8 and Table 1) with apparent $\log K_a = 21.35$ (Table 2). The association constant did not change at pH 5.5, but it was significantly reduced at pH 3. This is consistent with previous reports²⁴, which reported a decrease of conductance at lower pH, with a pK of 4.5.

[Place Figure 8 here]

[Place Table 1 here]

Figure Legends:

Figure 1. Exploded view of a 6-channel AUC cell without external fill. This figure has been modified from Beckman Coulter An-50 Ti and An-60 Ti Analytical Rotor, Cells, and Counterbalance user manual³⁵.

Figure 2. Exploded view of the window assembly. This figure has been modified from Beckman Coulter An-50 Ti and An-60 Ti Analytical Rotor, Cells, and Counterbalance user manual³⁵.

Figure 3. Result from HeteroAnalysis Match function. The Match function can be used to monitor equilibrium progress by comparing RMSD between successive scans and the last scan. This example shows attainment of equilibrium after 8 hours as indicated by RMSD values asymptotic to X-axis.

Figure 4. Results from HeteroAnalysis Match function. The first and second speeds (top left and right) appear to have reached equilibrium, but it is better to wait a few more hours to be sure. In comparison, the third and fourth speeds (bottom left and right) have clearly reached equilibrium in a shorter time.

Figure 5. An example on how to fill in Experimental Parameters.

Figure 6. An example on how to fill in Global Parameters for Monomer-n-Mer Self-Association model.

Figure 7. Density matching of C14SB with D₂O. (A) Radial distribution profile formed by C14SB detergent micelles in 50 mM Tris, 100 mM NaCl pH 7.3 with 25, 30, and 35% D₂O. The data was separately fitted to a linear model (black). (B) The slopes were plotted against D₂O concentration (black squares) and fitted to a linear model (red line). The matching D₂O concentration is indicated (red arrow).

Figure 8. Fitting of SH in C14SB to monomer-pentamer model. Radial distribution profile of SH in C14SB at pH 7 (open circles) was found to fit best to the monomer-pentamer self-association model (black solid line). Fitting residual is shown below.

Table 1. Comparison of global reduced chi-square values of different monomer-n-mer models

Table 2. Comparison of apparent log K_a values at different pH.

DISCUSSION:

This paper provides an experimental protocol for sample preparation and analysis of oligomerization of a small membrane protein in detergent using equilibrium sedimentation. The protocol described is equally valid –and simpler– for soluble proteins, as the density matching step is not required. Indeed, the system is constituted by a mixture of detergent and protein. To conduct sedimentation studies, the detergent must be invisible to the gravitational field so that it does not contribute to the particle flotation. Thus, the density of the detergent has to be carefully matched by addition of D₂O to the buffer, with the limitation that if the detergent used is too dense, e.g., SDS, not even 100 % D₂O can match it. The density matching step is not required when working with a water soluble protein, since the sample does not have detergent. In our application, a small viral protein that forms ion channels has been used. Like all viroporins, SH protein has an α -helical transmembrane domain, and therefore has to be studied in a detergent that will not disrupt its native oligomeric size.

Prior to these studies, therefore, suitable detergents have to be screened. For example, in previous SE studies SH protein formed pentameric oligomers in DPC, C8E5 and C14-betaine, and pentamers were also observed during electrophoresis in the mild detergent PFO²². All these detergents are suitable for SE, since their density can be matched by addition of D₂O. Samples should also be tested by sedimentation velocity (SV), which provides the minimum number of species present³⁶.

Also, unsuitable detergents will produce multiple species in SV, and multiple bands in PFO electrophoresis, indicative of multiple non-specific association. Useful information from SE is obtained when the predominant number of species present in the system is not higher than two; in the case of SH protein, the data was fitted to an equilibrium between monomers and pentamers.

It is important to note that the model chosen to fit the SE data, should also be as simple as possible, i.e., a single species of unknown molecular weight should be tried first, followed by reversible equilibria between monomer and oligomers of increasing size. More complex models have associated a higher ambiguity. Also, small populations of other smaller or larger oligomers cannot be detected, and the model will just indicate what are the predominant species present. In the case of viroporins, this is especially evident because oligomeric size changes depending on subtle experimental conditions, e.g., detergent used, pH, protein concentration, centrifugation speed or mutations, eg, in Hepatitis C Virus p7³⁷, influenza A M2³⁸ or SH protein²⁴.

The results obtained for WT SH protein clearly show that the association constant of the pentamer is reduced at pH 3. These results match channel activity measurements obtained with SH protein in synthetic lipid bilayers, where conductance was markedly reduced at pH 3 (pK ~ 4.5), whereas it remained constant between pH 7 and 5²⁴. Based on these results, it may be

possible that SH channel activity is somewhat regulated by lower pH when present in native biological membranes. Although in the Golgi lumen the pH is only one unit below that of the cytoplasm, intravesicular pH drops along the endocytic pathway from pH 6.0–6.5 in early endosomes to pH 4.5–5.5 in late endosomes and lysosomes³⁹. In the infected cell, the pK for conductance changes or pentameric stability may be higher than those obtained *in vitro* herein and previously²⁴, therefore pH could indeed play a role in modulating channel activity during the life cycle of the virus. Mutation of protonatable His residues in the context of the infected cell is an interesting avenue for future experiments⁴⁰.

Finally, the study of SH protein is facilitated by the presence of Trp in its sequence, facilitating UV absorption measurements. However, knowledge of the structure of the oligomer can enable introduction of Trp at insensitive locations, e.g., lipid –or solvent- facing parts of the protein in membrane – or soluble- proteins. Alternatively the protein may be tagged with a Visible-absorbing or fluorescent label.

In summary, this protocol describes the application of SE to determine oligomeric size and association constants of a defined viral channel, when experimental parameters are changed. In this case, pH has been varied to study the effect on stability of His protonation, but many other hypothesis can be tested, such as effects of mutations on the structural integrity of these oligomers.

ACKNOWLEDGEMENTS:

This work has been funded by the National Research Foundation grant NRF-CRP4-2008-02 (J.T.) and Tier 1 grant RG 51/13.

DISCLOSURES:

The authors have nothing to disclose.

REFERENCES:

- 1 Laue, T. M. & Stafford, W. F., 3rd. Modern applications of analytical ultracentrifugation. *Annu. Rev. Biophys. Biomol. Struct.* **28**, 75-100 (1999).
- 2 Lebowitz, J., Lewis, M. S. & Schuck, P. Modern analytical ultracentrifugation in protein science: A tutorial review. *Protein Sci.* **11**, 2067-2079, doi:Doi 10.1110/Ps.0207702 (2002).
- 3 Balbo, A., Zhao, H., Brown, P. H. & Schuck, P. Assembly, Loading, and Alignment of an Analytical Ultracentrifuge Sample Cell. e1530, doi:doi:10.3791/1530 (2009).
- 4 Rivas, G., Stafford, W. & Minton, A. P. Characterization of heterologous protein-protein interactions using analytical ultracentrifugation. *Methods-a Companion to Methods in Enzymology* **19**, 194-212, doi:DOI 10.1006/meth.1999.0851 (1999).
- 5 Howlett, G. J., Minton, A. P. & Rivas, G. Analytical ultracentrifugation for association and assembly the study of protein. *Curr. Opin. Chem. Biol.* **10**, 430-436, doi:DOI

- 10.1016/j.cbpa.2006.08.017 (2006).
- 6 MacGregor, I. K., Anderson, A. L. & Laue, T. M. Fluorescence detection for the XLI analytical ultracentrifuge. *Biophys. Chem.* **108**, 165-185, doi:DOI 10.1016/j.bpc.2003.10.018 (2004).
- 7 Phizicky, E. M. & Fields, S. Protein-Protein Interactions - Methods for Detection and Analysis. *Microbiol. Rev.* **59**, 94-123 (1995).
- 8 Alexandrov, A. *et al.* A facile method for high-throughput co-expression of protein pairs. *Mol. Cell. Proteomics* **3**, 934-938, doi:DOI 10.1074/mcp.T400008-MCP200 (2004).
- 9 Nooren, I. M. A. & Thornton, J. M. Structural characterisation and functional significance of transient protein-protein interactions. *J. Mol. Biol.* **325**, 991-1018, doi:DOI 10.1016/S0022-2836(02)01281-0 (2003).
- 10 Ebel, C. Sedimentation velocity to characterize surfactants and solubilized membrane proteins. *Methods* **54**, 56-66, doi:DOI 10.1016/j.ymeth.2010.11.003 (2011).
- 11 Minton, A. P. Quantitative characterization of reversible macromolecular associations via sedimentation equilibrium: an introduction. *Exp. Mol. Med.* **32**, 1-5 (2000).
- 12 Dowell, S. F. *et al.* Respiratory syncytial virus is an important cause of community-acquired lower respiratory infection among hospitalized adults. *J. Infect. Dis.* **174**, 456-462 (1996).
- 13 Bukreyev, A., Whitehead, S. S., Murphy, B. R. & Collins, P. L. Recombinant respiratory syncytial virus from which the entire SH gene has been deleted grows efficiently in cell culture and exhibits site-specific attenuation in the respiratory tract of the mouse. *J. Virol.* **71**, 8973-8982 (1997).
- 14 Fuentes, S., Tran, K. C., Luthra, P., Teng, M. N. & He, B. Function of the respiratory syncytial virus small hydrophobic protein. *J. Virol.* **81**, 8361-8366 (2007).
- 15 Jin, H. *et al.* Recombinant respiratory syncytial viruses with deletions in the NS1, NS2, SH, and M2-2 genes are attenuated in vitro and in vivo. *Virology* **273**, 210-218, doi:DOI 10.1006/viro.2000.0393 (2000).
- 16 Karron, R. A. *et al.* Respiratory syncytial virus (RSV) SH and G proteins are not essential for viral replication in vitro: clinical evaluation and molecular characterization of a cold-passaged, attenuated RSV subgroup B mutant. *Proc. Nat. Acad. Sci. USA* **94**, 13961-13966 (1997).
- 17 Whitehead, S. S. *et al.* Recombinant respiratory syncytial virus bearing a deletion of either the NS2 or SH gene is attenuated in chimpanzees. *J. Virol.* **73**, 3438-3442 (1999).
- 18 Rixon, H. W. *et al.* The small hydrophobic (SH) protein accumulates within lipid-raft structures of the Golgi complex during respiratory syncytial virus infection. *J. Gen. Virol.* **85**, 1153-1165 (2004).
- 19 Collins, P. L. & Mottet, G. Membrane orientation and oligomerization of the small hydrophobic protein of human respiratory syncytial virus. *J. Gen. Virol.* **74**, 1445-1450 (1993).
- 20 Collins, P. L., Olmsted, R. A. & Johnson, P. R. The small hydrophobic protein of human respiratory syncytial virus: comparison between antigenic subgroups A and B. *J. Gen. Virol.* **71**, 1571-1576 (1990).
- 21 Chen, M. D., Vazquez, M., Buonocore, L. & Kahn, J. S. Conservation of the respiratory syncytial virus SH gene. *J. Infect. Dis.* **182**, 1228-1233 (2000).

- 22 Gan, S. W. *et al.* The small hydrophobic protein of the human respiratory syncytial virus forms pentameric ion channels. *J. Biol. Chem.* **287**, 24671-24689 (2012).
- 23 Gan, S. W., Ng, L., Xin, L., Gong, X. & Torres, J. Structure and ion channel activity of the human respiratory syncytial virus (hRSV) small hydrophobic protein transmembrane domain. *Protein Sci.* **17**, 813-820 (2008).
- 24 Li, Y. *et al.* Inhibition of the Human Respiratory Syncytial Virus Small Hydrophobic Protein and Structural variations in a bicelle environment. *J. Virol.* **In press** (2014).
- 25 Burgess, N. K., Stanley, A. M. & Fleming, K. G. Determination of membrane protein molecular weights and association equilibrium constants using sedimentation equilibrium and sedimentation velocity. *Meth. Cell. Biol.* **84**, 181-211 (2008).
- 26 Cole, J. L., Lary, J. W., P. Moody, T. & Laue, T. M. Analytical Ultracentrifugation: Sedimentation Velocity and Sedimentation Equilibrium. *Meth. Cell. Biol.* **84**, 143-179 (2008).
- 27 Fleming, K. G. Determination of membrane protein molecular weight using sedimentation equilibrium analytical ultracentrifugation. *Curr. Protoc. Prot. Sci.* **53**, 7.12:17.12.11–17.12.13., doi:10.1002/0471140864.ps0712s53 (2008).
- 28 Beckman Coulter, I. (2005).
- 29 Mayer, G. *et al.* Studying membrane proteins in detergent solution by analytical ultracentrifugation: Different methods for density matching. *Prog. Colloid Polym. Sci.* **113**, 176-181 (1999).
- 30 Laue, T. in *Current Protocols in Protein Science* Ch. 20.3, 20.23.21–20.23.13 (John Wiley & Sons, Inc., 2001).
- 31 Burgess, N. K., Stanley, A. M. & Fleming, K. G. Vol. 84 181-211 (2008).
- 32 Gan, S. W. *et al.* The Small Hydrophobic Protein Of The Human Respiratory Syncytial Virus Forms Pentameric Ion Channels. *J. Biol. Chem.* **287**, 24671-24689 (2012).
- 33 Bevington, P. R. & Robinson, D. K. *Data reduction and error analysis for the physical sciences*. Vol. 336 (McGraw-Hill New York, 1969).
- 34 Schuck, P., Radu, C. G. & Ward, E. S. Sedimentation equilibrium analysis of recombinant mouse FcRn with murine IgG1. *Molecular Immunology* **36**, 1117-1125, doi:http://dx.doi.org/10.1016/S0161-5890(99)00093-0 (1999).
- 35 www.beckmancoulter.com.
<(https://www.beckmancoulter.com/wsrportal/techdocs?docname=LXLA-TB-003).>
- 36 Gan, S. W., Vararattanavech, A., Nordin, N., Eshaghi, S. & Torres, J. A cost-effective method for simultaneous homo-oligomeric size determination and monodispersity conditions for membrane proteins. *Anal. Biochem.* **416**, 100-106 (2011).
- 37 Montserret, R. *et al.* NMR structure and ion channel activity of the p7 protein from hepatitis C virus. *J. Biol. Chem.* **285**, 31446-31461 (2010).
- 38 Stouffer, A. L., DeGrado, W. F. & Lear, J. D. Analytical Ultracentrifugation Studies of the Influenza M2 Homotetramerization Equilibrium in Detergent Solutions. *Progr Colloid Polym Sci* **131**, 108-115 (2006).
- 39 Sorkin, A. & von Zastrow, M. Signal transduction and endocytosis: Close encounters of many kinds. *Nat. Rev. Mol. Cell Biol.* **3**, 600-614, doi:Doi 10.1038/Nrm883 (2002).
- 40 Gan, S. W., Ng, L., Lin, X., Gong, X. & Torres, J. Structure and ion channel activity of the human respiratory syncytial virus (hRSV) small hydrophobic protein transmembrane

domain. *Protein science : a publication of the Protein Society* **17**, 813-820, doi:10.1110/ps.073366208 (2008).

Figure 1
[Click here to download high resolution image](#)




Figure 2
[Click here to download high resolution image](#)




Figure 3
[Click here to download high resolution image](#)




Figure 4
[Click here to download high resolution image](#)







Figure 5

[Click here to download high resolution image](#)

Experimental Parameters

(1) ABSORBANCE data for MULTISPEED EQ

Commentinner section of 23000 RA3

☒ active

noise0.0050

☐ σ_{rot}
(N1/Nx)

Centripiece6

Pathlength1.200000

Rotor type8

☒ TI noise rotor stretch
mass cons. w. rotor stretch

v-bar (ml/g)0.7569

buffer density (g/ml)1.038390

Temperature20.0

☒ fit baseline

☒ fit Rd Noise

☐ fit TI Noise

☐ Meniscus5.8399

☒ Bottom6.0998

☐ redirect men./bot
1

For Associating Systems:

☐ extinction coefficient A8370.0000

☐ extinction coefficient B0.0000

☐ extinction coefficient C0.0000

☐ redirect at A
1

☐ redirect at B
1

☐ redirect at C
1

☐ use for signs of MC sims:0.0000

☒ use local mod

☒ correct for HD Exchange: buffer D mol fraction0.0000

Cancel

OK

Figure 6
[Click here to download high resolution image](#)

A + A + ... <--> [nA]

Monomer

☒ vary local concs
c(xp1) (uM)
ext.coeff (xp1)
☐ M(1)
☐ s(1)

☐ Mass Conservation
☐ add non-part. species

Complex [nA]

n
☒ log(Ka)
☐ log(k-)
☐ s (n)

☐ log10[ks/(mg/ml)] nonid sec
☐
☐ with 2nd virial coeff. MB [ml.
☐
☐ with hydrodyn kc

incompetent fraction

n
☐ fraction

Cancel

OK

Figure 7
[Click here to download high resolution image](#)






Figure 8
[Click here to download high resolution image](#)



Model (n-mer)	Chi-square (fixed MW)	Chi-square (fitted MW)	Fitted MW
3	15.444	1.0492	13477 Da
4	3.8094	1.0469	9889 Da
5	1.0499	1.0497	7822 Da
6	2.5994	1.547	6504 Da
7	6.1667	1.2112	5743 Da

pH	Lower limit (1σ)	Log Ka	Upper limit (1σ)
3	17.432	17.576	17.737
5.5	20.064	20.419	20.839
7.3	20.687	21.052	21.492

Name of Material/ Equipment	Company
3-(N,N-dimethylmyristylammonio)propanesulfonate	Sigma
Deuterium oxide 99.8%	Cambridge Isotope
An-50 Ti Rotor, Analytical, 8-Place	Beckman Coulter
An-60 Ti Rotor, Analytical, 4-Place	Beckman Coulter
Cell housing	Beckman Coulter
12 mm six-channel centerpiece, epon charcoal-filled	Beckman Coulter
Window holder	Beckman Coulter
Window gasket	Beckman Coulter
Window liner	Beckman Coulter
Sapphire window	Beckman Coulter
Quartz window	Beckman Coulter
Screw-ring washer	Beckman Coulter
Screw ring	Beckman Coulter
Spinkote	Beckman Coulter
Torque stand assembly	Beckman Coulter
Counterbalance	Beckman Coulter
Cell alignment tool	Beckman Coulter
SEDNTERP	
HeteroAnalysis	
SEDFIT	
SEDPHAT	

Catalog Number	Comments/Description
T0807	
DLM-4-99.8	
363782	
361964	
334784	
331376	
305037	
327021	
362329	
307177	
301730	
362328	
301922	
306812	
361318	
360219	
362340	
	http://bitcwiki.sr.unh.edu/index.php/Main_Page
	http://www.biotech.uconn.edu/auf/?i=aufftp
	http://www.analyticalultracentrifugation.com/sedfit.htm
	http://www.analyticalultracentrifugation.com/sedphat/default.htm

ARTICLE AND VIDEO LICENSE AGREEMENT

Title of Article:

Author(s):

HISTIDINE PROTONATION AFFECTS PENTAMERIC STABILITY OF A VIROPROLIN...

MAHYU SURYA AND JAUME TORRES

Item 1 (check one box): The Author elects to have the Materials be made available (as described at

http://www.jove.com/publish) via: ☒ Standard Access ☐ Open Access

Item 2 (check one box):

☒ The Author is NOT a United States government employee.

☐ The Author is a United States government employee and the Materials were prepared in the course of his or her duties as a United States government employee.

☐ The Author is a United States government employee but the Materials were NOT prepared in the course of his or her duties as a United States government employee.

ARTICLE AND VIDEO LICENSE AGREEMENT

1. Defined Terms. As used in this Article and Video license Agreement, the following terms shall have the following meanings: "Agreement" means this Article and Video license Agreement; "Article" means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein; "Author" means the author who is a signatory to this Agreement; "Collective Work" means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; "CRC License" means the Creative Commons Attribution-Non Commercial-No Derivs 3.0 Unported Agreement, the terms and conditions of which can be found at: <http://creativecommons.org/licenses/by-nc-nd/3.0/legalcode>; "Derivative Work" means a work based upon the Materials or upon the Materials and other pre-existing works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; "Institution" means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; "JOVE" means MyJove Corporation, a Massachusetts corporation and the publisher of *The Journal of Visualized Experiments*; "Materials" means the Article and / or the Video; "Parties" means the Author and JOVE; "Video" means any video(s) made by the Author, alone or in conjunction with any other parties, or by JOVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion of the Article, and in which the Author may or may not appear.

2. Background. The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JOVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.

3. Grant of Rights in Article. In consideration of JOVE agreeing to publish the Article, the Author hereby grants to JOVE, subject to **Sections 4 and 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the "Open Access" box has been checked in **Item 1** above, JOVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.

ARTICLE AND VIDEO LICENSE AGREEMENT

4. Retention of Rights in Article. Notwithstanding the exclusive license granted to JOVE in **Section 3** above, the Author shall, with respect to the Article, retain the non-exclusive right to use all or part of the Article for the non-commercial purpose of giving lectures, presentations or teaching classes, and to post a copy of the Article on the Institution's website or the Author's personal website, in each case provided that a link to the Article on the JOVE website is provided and notice of JOVE's copyright in the Article is included. All non-copyright intellectual property rights in and to the Article, such as patent rights, shall remain with the Author.

5. Grant of Rights in Video – Standard Access. This **Section 5** applies if the “Standard Access” box has been checked in **Item 1** above or if no box has been checked in **Item 1** above. In consideration of JOVE agreeing to produce, display or otherwise assist with the Video, the Author hereby acknowledges and agrees that, Subject to **Section 7** below, JOVE is and shall be the sole and exclusive owner of all rights of any nature, including, without limitation, all copyrights, in and to the Video. To the extent that, by law, the Author is deemed, now or at any time in the future, to have any rights of any nature in or to the Video, the Author hereby disclaims all such rights and transfers all such rights to JOVE.

6. Grant of Rights in Video – Open Access. This **Section 6** applies only if the “Open Access” box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby grants to JoVE, subject to **Section 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Video in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Video into other languages, create adaptations, summaries or extracts of the Video or other Derivative Works or Collective Works based on all or any portion of the Video and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. For any Video to which this Section 6 is applicable, JoVE and the Author hereby grant to the public all such rights in the Video as provided in, but subject to all limitations and requirements set forth in, the CRC License.

7. Government Employees. If the Author is a United States government employee and the Article was prepared in the course of his or her duties as a United States government employee, as indicated in Item 2 above, and any of the licenses or grants granted by the Author hereunder exceed the scope of the 17 U.S.C. 403, then the rights granted hereunder shall be limited to the maximum rights permitted under such

statute. In such case, all provisions contained herein that are not in conflict with such statute shall remain in full force and effect, and all provisions contained herein that do so conflict shall be deemed to be amended so as to provide to JOVE the maximum rights permissible within such statute.

8. Likeness, Privacy, Personality. The Author hereby grants JOVE the right to use the Author's name, voice, likeness, picture, photograph, image, biography and performance in any way, commercial or otherwise, in connection with the Materials and the sale, promotion and distribution thereof. The Author hereby waives any and all rights he or she may have, relating to his or her appearance in the Video or otherwise relating to the Materials, under all applicable privacy, likeness, personality or similar laws.

9. Author Warranties. The Author represents and warrants that the Article is original, that it has not been published, that the copyright interest is owned by the Author (or, if more than one author is listed at the beginning of this Agreement, by such authors collectively) and has not been assigned, licensed, or otherwise transferred to any other party. The Author represents and warrants that the author(s) listed at the top of this Agreement are the only authors of the Materials. If more than one author is listed at the top of this Agreement and if any such author has not entered into a separate Article and Video License Agreement with JOVE relating to the Materials, the Author represents and warrants that the Author has been authorized by each of the other such authors to execute this Agreement on his or her behalf and to bind him or her with respect to the terms of this Agreement as if each of them had been a party hereto as an Author. The Author warrants that the use, reproduction, distribution, public or private performance or display, and/or modification of all or any portion of the Materials does not and will not violate, infringe and/or misappropriate the patent, trademark, intellectual property or other rights of any third party. The Author represents and warrants that it has and will continue to comply with all government, institutional and other regulations, including, without limitation all institutional, laboratory, hospital, ethical, human and animal treatment, privacy, and all other rules, regulations, laws, procedures or guidelines, applicable to the Materials, and that all research involving human and animal subjects has been approved by the Author's relevant institutional review board.

10. JOVE Discretion. If the Author requests the assistance of JOVE in producing the Video in the Author's facility, the Author shall ensure that the presence of JOVE employees, agents or independent contractors is in accordance with the relevant regulations of the Author's institution. If more than one author is listed at the beginning of this Agreement, JOVE may, in its sole discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JOVE reserves the right, in its absolute and sole discretion and without giving any reason therefore, to accept or decline any work submitted to JOVE. JOVE and its employees, agents and independent contractors shall have

ARTICLE AND VIDEO LICENSE AGREEMENT

full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JOVE has sole discretion as to the method of making and publishing the Materials, including, without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

expense. All indemnifications provided herein shall include JOVE's attorney's fees and costs related to said losses or damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JOVE, its employees, agents or independent contractors.

11. Indemnification. The Author agrees to indemnify JOVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JOVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JOVE, making of videos by JOVE, or publication in JOVE or elsewhere by JOVE. The Author shall be responsible for, and shall hold JOVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to the making of a video by JOVE its employees, agents or independent contractors. All sterilization, cleanliness or decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's

12. Fees. To cover the cost incurred for publication, JOVE must receive payment before production and publication the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JOVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.

13. Transfer, Governing Law. This Agreement may be assigned by JOVE and shall inure to the benefits of any of JOVE's successors and assigns. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to me one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement required per submission.


CORRESPONDING AUTHOR:

Name: JAUME TORRES

Department: SCHOOL OF BIOLOGICAL SCIENCES

Institution: NANYANG TECHNOLOGICAL UNIVERSITY (SINGAPORE)

Article Title: HISTIDINE PRODUCTION AFFECTS PENTAMERIC STABILITY OF A VIROPROTEIN...

Signature: 

Date: 11. JUNE. 14

Please submit a signed and dated copy of this license by one of the following three methods:

- 1) Upload a scanned copy of the document as a pdf on the JOVE submission site;
- 2) Fax the document to +1.866.381.2236;
- 3) Mail the document to JOVE / Attn: JOVE Editorial / 1 Alewife Center #200 / Cambridge, MA 02139

For questions, please email submissions@jove.com or call +1.617.945.9051

Editorial comments:

Dear Editor,

Please see reply to the requested changes in blue.

Best regards

Jaume Torres

#####

Dear Dr Torres,

Your manuscript JoVE52404R3 'Sedimentation Equilibrium of a Small Oligomer-forming Membrane Protein: effect of Histidine Protonation on Pentameric Stability.' has been editorially reviewed and the following comments need to be addressed.

Please use the "track-changes" function in Microsoft Word or change the text color to identify all of your manuscript edits. When you have revised your submission, please also upload a list of changes, where you respond to each of the comments individually, in a separate document at the same time as you submit your revised manuscript.

Editorial comments:

Thank you for the Peer Review revisions that were made to the manuscript. The following editorial comments need to be addressed in order to finalize your manuscript.

- 1) All of your previous revisions have been incorporated into the most recent version of the manuscript, attached to this email. Please use the attached updated manuscript for all future revisions.

- 2) The length of your Short Abstract exceeds our 50 word limit. Please revise the Short Abstract so that it clearly states the goal of the protocol within 50 words. For example, “This protocol/manuscript describes...” Being a methods journal, the main focus of the JoVE article should be on the technique/method being presented. As written, the Short Abstract reads too much like an original research article. Please revise it so that the focus is on the methods.

Corrected

3) Please remove the embedded Tables from the manuscript. All tables should be uploaded to the Editorial Manager site in the form of Excel files. A description of the table should be included with the Figure legends.

Corrected

4) Editor removed the embedded figures from the manuscript as figures should be uploaded only to the JoVE submission site.

The figures were already deleted. One had to just accept the changes.

5) In the Table of Materials, please include the references for the software packages used as part of the main table. For e.g. "SEDNTERP" should be added to the "Name of Material/Equipment" column of the Table and the websites should be added to the "Comments/Description" column.

Corrected

6) Please provide higher resolution figures for Figures 3 and 4.

Provided

7) Figures 1 and 2 have been reproduced from a previous publication. If you are re-using figures from a previous publication, you must obtain explicit permission to re-use the figure from the previous publisher (this can be in the form of a letter from an editor or a link to the editorial policies that allows you to re-publish the figure). Please upload the text of the re-print permission (may be copied and pasted from an email/website) as a Word document to the Editorial Manager site in the "Supplemental files (as requested by JoVE)" section. Please also cite the figure appropriately in the figure legend, i.e. "This figure has been modified from [citation]."

Message from Beckman is provided, uploaded as Supplemental file. Reference to the manual is included (as a website, because it is an online document).

Your revision is due by **Nov 12, 2014**.

To submit a revision, go to the [JoVE submission site](#) and log in as an author. You will see a menu item called 'Submission Needing Revision'. You will find your submission record there.

Sincerely,

Sephorah Zaman, Ph.D.
Science Editor

[JoVE](#)

1 Alewife Center, Suite 200, Cambridge, MA 02140

Supplemental File (as requested by JoVE)

[Click here to download Supplemental File \(as requested by JoVE\): Randy_Pawlovich_reply\(Beckman\).tif](#)

Supplemental File (as requested by JoVE)_email from Beckman

[Click here to download Supplemental File \(as requested by JoVE\): Email_exchange_Beckman.pdf](#)