

Journal of Visualized Experiments

Labeling of Neuronal Morphology using Custom Diolistic Techniques

--Manuscript Draft--

Manuscript Number:	JoVE55574R4
Full Title:	Labeling of Neuronal Morphology using Custom Diolistic Techniques
Article Type:	Invited Methods Article - JoVE Produced Video
Keywords:	neuroscience; Neuronal Labeling; Diolistic; Neuronal Morphology; Synaptic Density; Gene Gun.
Manuscript Classifications:	1.8: Nervous System; 1.8.663: Neurons; 1.8.663.256: Dendrites; 1.8.663.256.200: Dendritic Spines; 1.8.850: Synapses; 5.1.370.225.750: Histological Techniques; 5.5.200.750.600.670: Staining and Labeling; 8.1.158.610: Neurosciences
Corresponding Author:	Lyon Hough, Ph.D. Missouri State University Springfield, Missouri UNITED STATES
Corresponding Author Secondary Information:	
Corresponding Author E-Mail:	LyonHough@missouristate.edu
Corresponding Author's Institution:	Missouri State University
Corresponding Author's Secondary Institution:	
First Author:	Lyon Hough, Ph.D.
First Author Secondary Information:	
Other Authors:	Michael E. Brown, Ph.D.
Order of Authors Secondary Information:	
Abstract:	<p>Diolistic labeling is increasingly utilized in neuroscience as a highly efficient, reproducible method for visualization and analysis of neuronal morphology. The use of lipophilic carbocyanine dyes, combined with particle-mediated biolistic delivery allows for non-toxic fluorescent labeling of multiple neurons with their dendritic arbors and spines in both living and fixed tissue. Since first described, this novel labeling method has been modified and adapted to fit a variety of research goals and laboratory settings. Diolistic labeling has traditionally relied on the use of a commercially available hand-held gene gun for the propulsion of coated micro-particles into tissue sections. Recently, laboratory built biolistic devices have been developed which allow for increased availability and customization for this method. Here, we discuss one such custom biolistic device and provide a detailed protocol for its use in Diolistic labeling. While decreasing the associated costs, the laboratory built device also overcomes many of the obstacles normally experienced with traditional diolistics to allow for reliable and reproducible neuronal labeling. The versatility of this method allows for its adaptation to a variety of laboratory settings and neuroscience related research goals.</p>
Author Comments:	<p>All editor and reviewer comments have either been addressed or rebutted in the "Hough 2016 Rebuttal to editorial and peer reviewer comments" file. Thank you for your consideration.</p> <p>The manuscript has been altered to adhere to editor suggestions and copy edited as requested.</p>
Additional Information:	
Question	Response
If this article needs to be "in-press" by a certain date, please indicate the date below and explain in your cover letter.	

TITLE:

Labeling Neuronal Morphology Using Custom Diolistic Techniques

AUTHORS:

Hough, Lyon H

Department of Biomedical Sciences

Missouri State University

Springfield, MO, USA

LyonHough@Missouristate.edu

Brown, Michael E.

Department of Pathology

Saint Louis University School of Medicine

Saint Louis, MO, USA

mbrown65@slu.edu

CORRESPONDING AUTHOR:

Lyon H. Hough

KEYWORDS:

Neuroscience, Labeling, Diolistic, Neuronal Morphology, Synaptic Density, Gene Gun.

SHORT ABSTRACT:

Here, we present a protocol for custom diolistic labeling. The customization of this fluorescent neuronal labeling method provides a modifiable technique that can be adapted to a wide variety of research goals and applications in the analysis of neuronal morphology.

LONG ABSTRACT:

Diolistic labeling is increasingly utilized in neuroscience as a highly efficient and reproducible method for the visualization and analysis of neuronal morphology. The use of lipophilic carbocyanine dyes, combined with particle-mediated biolistic delivery, allows for the non-toxic fluorescent labeling of multiple neurons, including their dendritic arbors and spines, in both living and fixed tissue. Since first described, this novel labeling method has been modified and adapted to fit a variety of research goals and laboratory settings. Diolistic labeling has traditionally relied on the use of a commercially available, hand-held gene gun for the propulsion of coated micro-particles into tissue sections. Recently, laboratory-built biolistic devices have been developed and allow for the increased availability and customization of this method. Here, we discuss one such custom biolistic device and provide a detailed protocol for its use in diolistic labeling. In addition to decreasing the associated costs, the laboratory-built device also overcomes many of the obstacles normally experienced with traditional diolistics, allowing for reliable and reproducible neuronal labeling. The versatility of this method allows for its adaptation to a variety of laboratory settings and neuroscience-related research goals.

INTRODUCTION:

45 Three-dimensional morphological reconstructions of individual neurons and their dendritic
46 arbors have served as the bases for analyzing the structure-function relationships within the
47 nervous system¹⁻³. For over a century, the main method for these studies consisted of various
48 modifications of the famed Golgi staining procedure⁴, which has proved invaluable in
49 developing our modern understanding of the nervous system^{1,2}. However, this method and the
50 various modifications of the silver impregnation of a percentage of the neuronal population are
51 not without drawbacks. In particular, the small population of stained cells, while an advantage
52 in morphological studies, elicits the lingering concern of selection bias^{2,5}. The Golgi method also
53 suffers from a limited compatibility with modern immunolabeling techniques and confocal
54 microscopy. Additionally, data analysis following the Golgi method of staining may be
55 negatively influenced by inconsistent neuronal impregnation, sectioning artifacts, and the
56 overlapping of fine dendritic processes that leads to an indistinguishable morphology.

57
58 More recently, the use of neuronal transfection and electroporation methods have allowed for
59 neuronal morphological labeling that circumvent the limitations experienced with Golgi
60 staining^{6,7}. While the administration of dye into individual cells through microinjection using
61 intracellular or patch pipettes generates excellent single-cell labeling, the technique is
62 technically demanding and may be vulnerable to sampling bias^{6,8-10}. Neuronal transfection
63 relies on the introduction of DNA constructs into target cells and tissues through a variety of
64 methods. One such method of DNA transfer is known as “Biolistic” delivery, in which a gene gun
65 utilizes a pressurized release of gas to propel DNA-coated micro-particles into tissue, crossing
66 the plasma membrane to target the cells³. This method does not require the same level of
67 technical expertise needed for conventional intracellular injections¹¹. Furthermore, the random
68 sampling used in biolistic delivery is ideal for full-scale quantitative analysis. However, the
69 reliance on DNA transfection for fluorescence expression limits the method to use on living
70 tissues in prepared cultures.

71
72 A more recent advancement from the biolistic approach to the morphological labeling of live or
73 fixed neurons was first reported by Gan *et al.* in 2000¹² and is known as diolistic labeling.
74 Diolistic labeling utilizes the lipophilic fluorescent dye dialkylcarbocyanine (Dil). While Dil has
75 traditionally been used in anterograde and retrograde neuronal tracing, it has also proven to be
76 effective in the fluorescent labeling of the neuronal cell membrane. Utilizing a diolistic
77 approach¹³, the ballistic delivery of Dil-coated micro-carriers to fixed or cultured tissue slices
78 allows the Dil to be incorporated into the cellular membrane. This occurs through lateral
79 diffusion in living or fixed tissue, illuminating the cellular morphology. When ballistically
80 delivered to tissue sections through a single pulse of high-purity gas (*i.e.*, helium or nitrogen),
81 individual Dil-coated tungsten particles are embedded in various neurons of the tissue. The
82 micro-carriers pass into the soma of the neuron, while the Dil is captured in the neuronal
83 membrane. Thus, the Dil is allowed to diffuse along the cellular membrane of a single neuron,
84 fluorescently labeling the fine neuronal architecture of dendritic branches and spines¹³. Dil-
85 labeled neurons can be observed through high-resolution imaging, such as confocal or two-
86 photon microscopy, and can be digitally reconstructed in precise detail. The quantification and
87 classification of dendritic branching and dendritic spines can be accomplished with appropriate
88 software packages^{14,15}.

Since its development, diolistic labeling has relied on the use of a commercially available gene gun, designed for biolistic transfection, for the propulsion of dye-coated micro-particles into tissue. However, these handheld devices have several drawbacks in the proposed staining method. First, the handheld nature of the device can alter the exact angle of delivery of the particles, causing inconsistent dye patterning in the tissue. Next, the standard tubular barrel of the gun causes increased particle density in the middle of the field and produces a burst of gas strong enough to damage the superficial layer of fixed tissue. Finally, the cost of commercially available devices and the associated materials may preclude some laboratories from using this type of methodology. In an effort to circumvent the aforementioned obstacles, Bridgman *et al.*¹¹ designed and constructed a custom-built device for use in biolistic applications (Figure 1). With protocol modifications, the device has been optimized for the diolistic labeling of fixed tissues. The device consists of a solenoid valve triggered by a relay switch to fire for 50 ms and is a modification of an original model constructed by Dr. David Kirk at Washington University in Saint Louis. The main features of the device that contributes to its reliability and reproducibility include a precisely timed trigger for the solenoid valve, a precise height adjustment system, a narrow baffled barrel, and a small pore size filter¹¹. Since the gun is mounted on a fixed base, the angle of delivery remains the same throughout all procedures, while the baffled barrel (constructed based on the design described by O'Brien)¹⁴ limits the amount of pressurized gas that contacts the tissue. Particle carriers (cartridges) are supported by a plastic ring that fits inside a modified filter holder (cartridge holder)¹¹. The particle carriers are cut-off yellow (200 µL) pipette tips. The device components listed here allow for a similar device to be constructed and operated with commercially available materials for a lower cost than traditional gene gun systems.

PROTOCOL:

All protocols using live animals must first be reviewed and approved by an Institutional Animal Care and Use Committee (IACUC) and must follow officially approved procedures for the care and use of laboratory animals. This protocol has been approved by the IACUC at Missouri State University.

Note: The following protocol will detail the methods used in the diolistic labeling of fixed brain tissue obtained from rats prepared using transcardial perfusion. Previous studies¹⁶ have utilized similar methods of labeling but without the transcardial perfusion and from various animal species.

1. Dil/Tungsten-Coated Bead Preparation

Note: The Dil (1-1'-Diocetadecyl- 3,3,3',3'-tetra methylindole carbocyanine perchlorate)-coated beads should be prepared at least 72 h in advance of tissue preparation to allow for adequate drying and optimal labeling. Prior reports¹ have shown increased success using CM-Dil when performing additional immunolabeling.

1.1. Place a clean glass slide on waxed weigh paper under the fume hood.

1.2. Add 90 mg of 1.3- μ m tungsten powder to a glass slide and dice thoroughly using a clean razor blade. Spread the fine powder evenly over the slide in a thin layer free of clumps.

1.3. Dissolve 2 mg of DiI with 75 μ L of methylene chloride in a 200- μ L microfuge tube under the fume hood. Seal the tube quickly.

CAUTION: Methylene chloride is caustic and evaporates extremely quickly; care should be taken to work under the fume hood and to minimize the time uncovered.

1.4. Vortex the microfuge tube to ensure that the DiI dissolves completely.

1.5. Slowly release the DiI/methylene chloride mixture over the surface of the 90 mg of evenly-spread, 1.3- μ m tungsten particles on the glass slide, ensuring that the DiI is applied evenly across the tungsten particles. The pipette tip can be used to ensure the even mixing of the tungsten particles and the DiI solution.

1.6. Allow the methylene chloride to dry under the hood until the bead mixture turns a light gray. Using a single-edged razor blade, scrape the DiI-coated tungsten particles loose from the slide. Continue to chop the DiI/tungsten particles evenly across the glass slide until the particles appear as a fine powder.

1.7. Scrape the coated tungsten beads carefully onto weigh paper and transfer them to a 15-mL conical tube containing 10 mg/mL PVP in distilled water. Vortex vigorously for 30 s.

1.8. Sonicate the DiI/PVP solution for 10 min in a water bath at room temperature.

2. Cartridge Preparation

Note: Cartridges are made using standard 200- μ L yellow pipette tips. Low-retention pipette tips should not be used in this protocol, as they may interfere with particle adherence to the inner wall of the tip. The protocol outlined here utilizes standard polypropylene research-grade pipette tips throughout the testing.

2.1. Immediately following sonication, vortex the DiI/PVP solution vigorously for 10-20 s.

2.2. Pull 100 μ L of the DiI/PVP solution from the middle of the conical tube into a mounted pipette tip cartridge.

2.3. Hold the pipette in a horizontal position, parallel to the ground, and rotate it slowly over 360°.

Note: Some particles suspended in the solution should attach to the pipette tip inner wall in an even fashion. Rotating too quickly or too slowly will affect the adherence of particles to the wall.

2.4. Slowly release the PVP solution and unbound particles back into the conical vial and gently disconnect the pipette tip from the pipette while continuing a slight axial rotation. Lay the pipette tip horizontally on the benchtop and avoid any unnecessary disturbance.

2.5. Re-vortex the original Dil/PVP solution for 10-20 s and repeat the procedure (steps 2.1-2.4) with the next cartridge. Repeat the procedure until the Dil/PVP solution appears relatively clear (typically ~25 cartridges are prepared).

2.6. Gently remove the cartridge/pipette loading end and 10 mm of the tip using a single-edge razor blade, revealing ~20 mm of the middle of the cartridge/pipette. Determine the successful loading of the cartridge by confirming the presence of a light haze of tungsten beads adhering to the inside wall of the cartridge, without the presence of large aggregations or clumps of beads. Discard any unsuccessful cartridges.

2.7. Protect the cartridges from light and store them with a desiccant for a minimum of 72 h and a maximum of 1 month if stored at 4 °C.

3. Tissue Preparation

Note: Through repeated trials, it has been found that tissue fixation by the transcardial perfusion of a low-concentration (1.5%) aldehyde prepared with 25 mM PBS greatly enhances the results of neuronal labeling. Higher concentrations (4%) were determined to over-fixate the tissues and result in incomplete labeling. Lower fixative levels resulted in excessive dye leakage from cellular membranes.

3.1. Anesthetize the animal and confirm deep anesthesia through an IACUC-approved method. Open the abdominal and thoracic walls, exposing the peritoneal and thoracic cavities.

3.2. Transcardially perfuse through the left ventricle with ice-cold 25 mM PBS (pH 7.2) until a clear perfusate consistently returns through the opened right atrium; the clearing of the liver (*i.e.*, a color change from dark red to pale brown) also indicates a successful exsanguination.

3.3. Transcardially perfuse with ice-cold 1.5% paraformaldehyde until the neck and hindquarters are rigid.

Note: The times and amounts of fixative for the perfusion may vary between labs and between the species, age, and size of the animals.

3.4. Carefully remove the brain by exposing the skull and removing the calvarium with rongeurs. Take care to avoid any unnecessary or direct handling of the brain tissue.

3.5. Postfix the intact brain for 2 h in 1.5% paraformaldehyde.

3.6. Prepare regions of interest or whole-brain specimens under a dissecting scope, taking

care to prevent tissue drying.

3.7. Prepare a 4% solution of low-melting point agarose in PBS. Heat the agarose until molten with constant stirring to avoid excessive bubble formation.

3.8. Pour the molten agarose into an appropriately sized histology embedding mold for the brain tissue and keep it on a 40 °C warming plate. Orient the brain tissue in the desired plane for sectioning and place it in the molten agarose. Remove the agarose from the heat and allow it to cool to room temperature before handling.

3.8.1. Avoid placing it in a refrigerated location to speed cooling, as this will negatively affect the embedding and may cause the agarose to break away from the tissue during sectioning.

3.9. Remove the cooled block from the mold and trim the excess agarose from the tissue, making sure that all areas of the brain tissue are surrounded by at least 5 mm of agarose.

3.10. Secure the tissue block to a mounting block using a cyanoacrylate glue and allow it to dry for a few minutes.

3.11. While drying, prepare culture plates by filling each well with 2 mL of PBS.

3.12. Place the mounting block onto the vibratome stage and secure it. Lower the stage completely. Attach a clean blade, free of any preserving oils.

3.13. Fill the vibratome chamber with cold 25 mM PBS to the specified level so that the blade is covered.

Note: Only PBS should be used to fill the vibratome well. Lower concentrations or pure water will cause the tissue to swell.

3.14. Trim the block in 300- μ m increments until the desired section depth is reached.

3.15. Vibratome-section the brain into 150- to 300- μ m serial sections throughout the region of interest. A high amplitude (8.5 Hz) and low speeds (0.5) are recommended for sectioning brain tissue but should be tested on the vibratome in each laboratory.

3.16. Using a paintbrush, remove the sliced tissues from the PBS vibratome bath and transfer them into the tissue culture dishes containing PBS. Store them until the DiI labeling is to be performed.

Note: Prior experiments have found that tissues may be stored at 4 °C for up to one week prior to labeling without significant impact on tissue labeling^{1,17}.

4. Diolistic Cellular Labeling

265 4.1. Carefully remove previously prepared Dil/tungsten-coated cartridges from the storage
266 container and load them into individual cartridge holders.

267
268 4.2. Load each individual cartridge holder into the device and set the gas pressure to 120 psi
269 using the regulator.

270
271 Note: The device presented here uses high-purity compressed nitrogen as a propellant, but
272 helium may be used instead.

273
274 4.3. Place the culture plate under the device barrel and carefully pull off the PBS using a
275 pipette. Take care to avoid damage to tissue, and ensure that the tissue lays flat on the bottom
276 of the culture plate. Adjust the position using a paint brush, if needed.

277
278 Note: Fluid left covering the tissue can cause beads to ricochet from the surface, negatively
279 impacting tissue labeling.

280
281 4.4. Place a 3.0- μ m microplate filter insert over the tissue and lower the barrel to a set point
282 15 mm above the tissue surface. Avoid handling the filter inserts with ungloved hands, as dirt
283 and oils may clog the filter.

284
285 4.5. Activate the device using a relay switch that allows a 50-ms burst of gas to pass through
286 the device and cartridge.

287
288 4.6. After the shot, immediately raise the barrel and remove the filter insert. Return the PBS
289 to the well, being careful not to damage the tissue. Take care to ensure that the tissue surface,
290 vertical at the time of labeling, remains vertical for the remainder of the procedure; this is
291 especially important for thicker tissue sections.

292
293 4.7. Load a new cartridge into the device and repeat the process (steps 4.1-4.6) for each
294 tissue section in the culture plate.

295
296 4.8. Once completed, protect the entire culture plate from light and allow it to incubate for
297 24 h at room temperature.

298
299 Note: Through testing, it was determined that room temperature (20-23.5 °C) was ideal to
300 allow the diffusion of the dye throughout all lipid neuronal membranes of all tested neuronal
301 populations.

302
303 4.9. Post-fix the sections for 1 h in 2 mL of 4% paraformaldehyde; this step has proven to be
304 useful in minimizing "overstaining," or the leaching of dye out of neuronal membranes. Wash
305 briefly with PBS and store in fresh PBS at 4 °C after post-fixation.

306
307 Note: Tissues can remain stored at 4 °C for up to a week before mounting but should be
308 mounted and imaged as soon as possible for best results.

4.10. Mount tissue slices on coated or charged slides and coverslip with an appropriate antifade agent (5% n-propyl-gallate).

Note: If time allows for imaging in the 24-72 h post-labeling, wet-mounting the slide in PBS and sealing the edges with nail polish or another sealant provides a cost-effect and time-sensitive mounting method. Wet-mounted tissues should be stored in humid chamber at 4 °C when not undergoing microscopy.

5. Confocal Microscopy and Analysis

5.1. Image the mounted slides with any confocal microscope equipped with a motorized z-stage and lasers between 543 and 561 nm.

Note: The use of epifluorescence microscopy equipped with TRITC filters will permit the visualization of the Dil-impregnated neurons and can be an important tool in locating areas of interest. However, the use of epifluorescence microscopy in the quantitative three-dimensional morphological analysis of neurons is discouraged. Confocal microscopy allows for increased resolution in images of Dil-impregnated tissues and for optical sectioning through the cells of interest. Digital reconstructions of the optical sections can be assembled into a three-dimensional model for the analysis of morphology¹. The representative results presented in this report were obtained using sampling parameters detailed by Staffened and Meisel, 2011¹. It is recommended that sampling parameters and experimentally/statistically appropriate cell numbers and dendritic segments be established before imaging. There are a variety of commercially available software extensions that can be used for image analysis and morphological reconstructions. These product choices should be based on individual laboratory needs.

REPRESENTATIVE RESULTS:

Using the procedure presented here, custom diolistic labeling has been used to characterize alterations in neuronal morphology of the lateral/dentate nucleus of the cerebellum. Here, we show representative labeling results from developing rat pups exposed to abnormally high levels of the serotonergic agonist 5-methoxytryptamine (5-MT), both pre-and postnatally¹⁸. With the aid of quantitative software, dendritic branching morphology (Figure 2), architecture (Figure 3), and synaptic connections (Figure 4) were analyzed.

Successful diolistic labeling can be rapidly confirmed through cursory visualization with bright-field and epifluorescent microscopes equipped with red (TRITC) filter (Figure 5). Areas of interest (Figure 5A and B) and potential cells for imaging can be identified (Figures 5C, D, and E), though morphological details are often obscured prior to confocal imaging. Unsuccessful labeling (Figure 6) can be caused by improper bead and cartridge preparation or by storage, and care should be taken during steps 1 and 2 of this protocol. Over-labeling (Figure 6) is commonly caused by the aggregations of coated beads in one area of the cartridge (step 2.6) and can be minimized by using the proper technique or by decreasing the bead concentration in the slurry solution (step 1.7). Following confocal image acquisition, max projection z-stacks

can be quickly obtained through the use of the open-source image processing resource, ImageJ (NIH), and used to evaluate the criteria for acceptable labeling (Figure 7).

Figure 1: Custom biolistic device.

Key features include a baffled barrel to improve particle delivery, a knob for height adjustment, and a fixed base. The firing mechanism consists of a solenoid valve activated by a relay switch that allows a burst of nitrogen gas to pass through the cartridge and filter assembly and into the tissue specimens.

Figure 2: Dil-labeled neuron of the lateral/dentate nucleus of the rat cerebellum.

A Z-projection of Dil-labeled neurons are used to illustrate the dendritic morphological analysis. The inset image illustrates the three-dimensional computer tracing models used for the quantification of the morphology. A morphometric software package was used to contour the dendritic processes of the digitized 2-D images and to develop three-dimensional models.

Figure 3: Sholl analysis and a comparison of dendritic branching patterns.

Significant differences in branching morphology were observed in experimental (5-MT) groups compared to control (saline) populations. (**) denotes a significance value of $p \leq 0.001$, (*) denotes a significance value of $p \leq 0.05$. This figure has been modified from Hough and Segal 2016¹⁸.

Figure 4: Analysis of the total dendritic spine numbers and dendritic spine density.

The mean total number of spines per neuron per group (A)¹⁸ and the mean dendritic spine density per group (B)¹⁸. Representative Dil-labeled dendritic segments used for analysis of control (C) and experimental (D) groups. (**) denotes a significance value of $p \leq 0.001$, (*) denotes a significance value of $p \leq 0.05$. This figure has been modified from Hough and Segal 2016¹⁸.

Figure 5: Bright-field and epifluorescence overview of tissue sections and an assessment of cytoarchitectural boundaries¹⁸.

Composite image (A) of low-magnification bright-field photomicrographs. The lateral/dentate nucleus of the cerebellum (LN) (B) (solid arrow) can be viewed at the lateral aspect of the deep cerebellar nuclei, surrounded laterally by white matter. The medial aspect of the LN abuts the dorsolateral hump (DL) (dashed arrow) of the anterior interposed (AI) nuclei (dotted arrow), with the superior cerebellar peduncle (arrowhead) extending from the hilus of the LN. Epifluorescence photomicrographs at low (10x) (C), medium (20x) (D), and high (60x) (E) magnifications. A single LN neuron (arrow) (C and D) is illuminated with the TRITC filter and can be further imaged using the confocal microscope. High magnification (E) resolves some dendritic segments that extend beyond the focal plane (arrowhead), illustrating the need for confocal imaging at a specific z-depth. The LN is outlined in white (C). This figure has been modified from Hough and Segal 2016¹⁸.

Figure 6: Representative over-labeling.

Composite image (A) of low-magnification bright-field photomicrographs representative of

failed (over-)labeling with tungsten microparticles (arrow). Dil staining and over-labeling can be visually appreciated by the presence of the pink (arrowhead) coloration outside of a single cell. The epifluorescence photomicrograph (B) shows the over-labeled areas (arrows), with indistinguishable cellular structure.

Figure 7: Retrograde versus anterograde labeling.

A laser scanning confocal photomicrograph z-project of LN neurons, representative of retrograde (A) and anterograde (B) Dil labeling. Note the differential appearance of the cell soma (white arrow) with (A) appearing shadowed and incomplete compared to the soma pictured in (B). Care should be taken to avoid including retrograde-labeled neurons in the quantitative analysis if possible.

DISCUSSION:

Here, we demonstrate a method of custom diolistic labeling to quantitatively analyze neuronal morphology and synaptic connectivity. The versatility of this method allows for its adaptation to a variety of laboratory settings and research goals. While the results presented here exhibit its use in rat neural tissue, other studies have used diolistic labeling to investigate diverse species through various neuroscience-related applications¹⁶. The method is relatively fast, as it takes 1-2 days from tissue fixation to image acquisition and has shown promise in its ability to be combined with other, more classical labeling approaches, such as immunostaining^{16,19}. With a laboratory-designed and -built device, researchers can take advantage of low-cost materials, which may allow for the method to be incorporated into budget-conscious investigations. Material availability may also assist with increasing the number of tissue sections labeled, in turn increasing the cell numbers used in the analysis while decreasing the effects of failed labeling on a particular study.

While diolistic labeling has been shown to be a reliable and reproducible method of dendritic labeling, failures in labeling have been reported. To mitigate these failures, it is important to limit variables in cartridge preparation and storage while outlining clear selection criteria for the inclusion of stained neurons for analysis. Neurons suitable for morphological analysis should: (1) be contained within a region of interest, as defined by cytoarchitectural boundary criteria; (2) demonstrate complete fluorescent labeling of well-defined, tapered dendritic endings; and (3) show primary and secondary dendritic trees with clear branching points. To select individual dendritic segments for imaging and spine analysis, the segment should satisfy the following criteria: (1) the entire dendritic segment, as assessed by tapered, well-defined endings, should fall within the optical depth of the tissue section; (2) dendritic segments should be either parallel or at acute angles to the plane of the section to allow for visualization; and (3) segments should not be obscured by the overlap of other dendritic segments at the same z-depth, which may interfere with the quantification of dendritic spines^{18,20}. By setting criteria such as these, the impact of sampling bias and variations in Dil labeling can be minimized.

The most critical step in this protocol is the preparation of the cartridges. Properly prepared cartridges should only exhibit a light haze of beads adhering to the inner wall of the cartridge. Large aggregations of beads lead to the overstaining of the tissue. Additionally, the cartridges

should be adequately stored: protected from light, under desiccation, and at a cool temperature. High temperatures can cause the dye to separate from the tungsten bead, while light will decrease the fluorescence of the dye. Harsh handling of the dried cartridges can cause the “dust” of coated beads to fall from the cartridge prior to shooting, resulting in sparse labeling. Current limitations of this method exist in the truncated timetable from labeling to imaging. While an antifade mounting medium can preserve the labeling for extended times, these times pale in comparison to those of other neuronal labeling methods, such as Golgi staining. Experimentation on prolonging tissue labeling by using submersion in antifreeze solutions and storage in sub-zero temperatures is currently underway, with positive preliminary results. Secondary limitations in the method concern the ability to precisely target neuronal populations in a large tissue sample. The device is currently designed to provide a sparse labeling pattern, uniformly spread over a tissue. This approach limits the ability of an investigator to selectively target a population of cells with any accuracy. Future advancements in the firing mechanism and micro-barrel design may facilitate adjustments and refinements in aiming and selection in the future.

The diolistic approach to fluorescently label neurons in their entirety has many advantages when compared to more traditional techniques. The rapid nature of the method, the passive diffusion of the dyes, and the compatibility with confocal microscopy provide an unparalleled opportunity to examine features of neuronal cells at a high spatial resolution in a complex, three-dimensional tissue environment. Diolistic labeling should be viewed as a modifiable technique that can be customized to individual research goals and applications.

ACKNOWLEDGMENTS:

The authors are grateful for the technical assistance provided by Dr. Paul Bridgman of Washington University School of Medicine. Appreciation is also given to Dr. David Kirk and Dr. Michael Nonet of Washington University School of Medicine, who have made custom biolistic devices for their respective laboratories.

DISCLOSURES:

The authors have nothing to disclose.

REFERENCES:

- 1 Staffend, N. A. & Meisel, R. L. DiOlistic labeling in fixed brain slices: phenotype, morphology, and dendritic spines. *Curr Protoc Neurosci*. **Chapter 2** (2011).
- 2 Staffend, N. A. & Meisel, R. L. DiOlistic Labeling of Neurons in Tissue Slices: A Qualitative and Quantitative Analysis of Methodological Variations. *Front Neuroanat*. **5** 14-14 (2011).
- 3 Wu, C.-C., Reilly, J. F., Young, W. G., Morrison, J. H. & Bloom, F. E. High-throughput morphometric analysis of individual neurons. *Cereb Cortex*. **14** (5), 543-554 (2004).
- 4 Fairen, A. Pioneering a golden age of cerebral microcircuits: the births of the combined Golgi-electron microscope methods. *Neuroscience*. **136** (3), 607-614 (2005).
- 5 Spergel, D. J., Kruth, U., Shimshek, D. R., Sprengel, R. & Seeburg, P. H. Using reporter genes to label selected neuronal populations in transgenic mice for gene promoter, anatomical, and physiological studies. *Prog Neurobiol*. **63** (6), 673-686 (2001).

- 6 Buhl, E. H. & Lubke, J. Intracellular lucifer yellow injection in fixed brain slices combined with retrograde tracing, light and electron microscopy. *Neuroscience*. **28** (1), 3-16 (1989).
- 7 Spacek, J. Dynamics of the Golgi method: a time-lapse study of the early stages of impregnation in single sections. *J Neurocytol*. **18** (1), 27-38 (1989).
- 8 de Lima, A. D., Voigt, T. & Morrison, J. H. Morphology of the cells within the inferior temporal gyrus that project to the prefrontal cortex in the macaque monkey. *J Comp Neurol*. **296** (1), 159-172 (1990).
- 9 Schmidt, M., Lehnert, G., Baker, R. G. & Hoffmann, K. P. Dendritic morphology of projection neurons in the cat pretectum. *J Comp Neurol*. **369** (4), 520-532 (1996).
- 10 Surkis, A., Taylor, B., Peskin, C. S. & Leonard, C. S. Quantitative morphology of physiologically identified and intracellularly labeled neurons from the guinea-pig laterodorsal tegmental nucleus in vitro. *Neuroscience*. **74** (2), 375-392 (1996).
- 11 Bridgman, P. C., Brown, M. E. & Balan, I. Biolistic transfection. *Methods Cell Biol*. **71** 353-368 (2003).
- 12 Gan, W. B., Grutzendler, J., Wong, W. T., Wong, R. O. & Lichtman, J. W. Multicolor "DiOlistic" labeling of the nervous system using lipophilic dye combinations. *Neuron*. **27** (2), 219-225 (2000).
- 13 Gan, W. B., Grutzendler, J., Wong, R. O. & Lichtman, J. W. Ballistic delivery of dyes for structural and functional studies of the nervous system. *Cold Spring Harb Protoc*. **2009** (4), pdb.prot5202, doi:10.1101/pdb.prot5202, (2009).
- 14 O'Brien, J. A., Holt, M., Whiteside, G., Lummis, S. C. & Hastings, M. H. Modifications to the hand-held Gene Gun: improvements for in vitro biolistic transfection of organotypic neuronal tissue. *J Neurosci Methods*. **112** (1), 57-64 (2001).
- 15 Noterdaeme, M., Mildenerger, K., Minow, F. & Amorosa, H. Evaluation of neuromotor deficits in children with autism and children with a specific speech and language disorder. *Eur Child Adolesc Psychiatry*. **11** (5), 219-225 (2002).
- 16 Seabold, G. K., Daunais, J. B., Rau, A., Grant, K. A. & Alvarez, V. A. DiOLISTIC labeling of neurons from rodent and non-human primate brain slices. *J Vis Exp*. (41), (2010).
- 17 Forlano, P. M. & Woolley, C. S. Quantitative analysis of pre- and postsynaptic sex differences in the nucleus accumbens. *J Comp Neurol*. **518** (8), 1330-1348, doi:10.1002/cne.22279, (2010).
- 18 Hough, L. H. & Segal, S. Effects of Developmental Hyperserotonemia on the Morphology of Rat Dentate Nuclear Neurons. *Neuroscience*. doi:10.1016/j.neuroscience.2016.02.021, (2016).
- 19 Lee, P. R., Brady, D. L., Shapiro, R. A., Dorsa, D. M. & Koenig, J. I. Social interaction deficits caused by chronic phencyclidine administration are reversed by oxytocin. *Neuropsychopharmacology*. **30** (10), 1883-1894 (2005).
- 20 Rasmuson, S. *et al*. Environmental enrichment selectively increases 5-HT1A receptor mRNA expression and binding in the rat hippocampus. *Brain Res Mol Brain Res*. **53** (1-2), 285-290 (1998).

Figure 1





Figure 2



Dendritic Intersections




Figure 4

[Click here to download Figure Figure 4.TIF](#)




Figure 5




Figure 6

[Click here to download Figure Figure 6.TIF](#)





Figure 7

[Click here to download Figure Figure 7.TIF](#)



Name of Reagent/ Equipment

Granular Paraformaldehyde
Sodium Phosphate Monobasic Anhydrous
Sodium Phosphate Dibasic Anhydrous
Type 1-A, Low EEO Gel Agarose,
Polyvinylpyrrolidone
Dichloromethane (Methylene Chloride)
Yellow Pipet tips 1000pk
Corning Transwell Multiple Well Plate with Permeable Polyester Membrane Inserts
Carbocyanine fluorescent DiI 100mg
ProLong Gold Antifade
1.3- μ m Tungsten particles 6G
Millipore Swinnex Filter Holder
Neurolucida Neuron Tracing Software

Company	Catalog Number	Comments/Description
Sigma-Aldrich	P6148-5KG	
Sigma-Aldrich	S2554-500G	
Sigma-Aldrich	71640-1KG	
Sigma Aldrich	A0169-250G	
Sigma Aldrich	77627-100G	
Sigma Aldrich	D65100-1L	
Fisher Scientific	2681151	
Fisher Scientific	07-200-155	
Invitrogen	D-282	
Invitrogen	P36930	
BioRad	165-2269	
EDMmillipore	SX0001300	
MicroBrightField	n/a	



1 Alewife Center #200
Cambridge, MA 02140
Tel. 617.945.9051
www.jove.com

ARTICLE AND VIDEO LICENSE AGREEMENT

Title of Article:

Labeling of Neuronal Morphology using Custom Diolistic Techniques

Author(s):

Hough, Lyon H.; Brown, Michael E.

Item 1 (check one box): The Author elects to have the Materials be made available (as described at <http://www.jove.com/publish>) via: ☒ Standard Access ☐ Open Access

Item 2 (check one box):

- ☐ The Author is NOT a United States government employee.
- ☒ The Author is a United States government employee and the Materials were prepared in the course of his or her duties as a United States government employee.
- ☐ The Author is a United States government employee but the Materials were NOT prepared in the course of his or her duties as a United States government employee.

ARTICLE AND VIDEO LICENSE AGREEMENT

1. **Defined Terms.** As used in this Article and Video License Agreement, the following terms shall have the following meanings: "**Agreement**" means this Article and Video License Agreement; "**Article**" means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein; "**Author**" means the author who is a signatory to this Agreement; "**Collective Work**" means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; "**CRC License**" means the Creative Commons Attribution-Non Commercial-No Derivs 3.0 Unported Agreement, the terms and conditions of which can be found at: <http://creativecommons.org/licenses/by-nc-nd/3.0/legalcode>; "**Derivative Work**" means a work based upon the Materials or upon the Materials and other pre-existing works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; "**Institution**" means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; "**JoVE**" means MyJoVE Corporation, a Massachusetts corporation and the publisher of *The Journal of Visualized Experiments*; "**Materials**" means the Article and / or the Video; "**Parties**" means the Author and JoVE; "**Video**" means any video(s) made by the Author, alone or in conjunction with any other parties, or by JoVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion of the Article, and in which the Author may or may not appear.

2. **Background.** The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JoVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.

3. **Grant of Rights in Article.** In consideration of JoVE agreeing to publish the Article, the Author hereby grants to JoVE, subject to **Sections 4 and 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the "Open Access" box has been checked in **Item 1** above, JoVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.

ARTICLE AND VIDEO LICENSE AGREEMENT

4. Retention of Rights in Article. Notwithstanding the exclusive license granted to JoVE in **Section 3** above, the Author shall, with respect to the Article, retain the non-exclusive right to use all or part of the Article for the non-commercial purpose of giving lectures, presentations or teaching classes, and to post a copy of the Article on the Institution's website or the Author's personal website, in each case provided that a link to the Article on the JoVE website is provided and notice of JoVE's copyright in the Article is included. All non-copyright intellectual property rights in and to the Article, such as patent rights, shall remain with the Author.

5. Grant of Rights in Video – Standard Access. This **Section 5** applies if the "Standard Access" box has been checked in **Item 1** above or if no box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby acknowledges and agrees that, Subject to **Section 7** below, JoVE is and shall be the sole and exclusive owner of all rights of any nature, including, without limitation, all copyrights, in and to the Video. To the extent that, by law, the Author is deemed, now or at any time in the future, to have any rights of any nature in or to the Video, the Author hereby disclaims all such rights and transfers all such rights to JoVE.

6. Grant of Rights in Video – Open Access. This **Section 6** applies only if the "Open Access" box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby grants to JoVE, subject to **Section 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Video in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Video into other languages, create adaptations, summaries or extracts of the Video or other Derivative Works or Collective Works based on all or any portion of the Video and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. For any Video to which this **Section 6** is applicable, JoVE and the Author hereby grant to the public all such rights in the Video as provided in, but subject to all limitations and requirements set forth in, the CRC License.

7. Government Employees. If the Author is a United States government employee and the Article was prepared in the course of his or her duties as a United States government employee, as indicated in **Item 2** above, and any of the licenses or grants granted by the Author hereunder exceed the scope of the 17 U.S.C. 403, then the rights granted hereunder shall be limited to the maximum rights permitted under such

statute. In such case, all provisions contained herein that are not in conflict with such statute shall remain in full force and effect, and all provisions contained herein that do so conflict shall be deemed to be amended so as to provide to JoVE the maximum rights permissible within such statute.

8. Likeness, Privacy, Personality. The Author hereby grants JoVE the right to use the Author's name, voice, likeness, picture, photograph, image, biography and performance in any way, commercial or otherwise, in connection with the Materials and the sale, promotion and distribution thereof. The Author hereby waives any and all rights he or she may have, relating to his or her appearance in the Video or otherwise relating to the Materials, under all applicable privacy, likeness, personality or similar laws.

9. Author Warranties. The Author represents and warrants that the Article is original, that it has not been published, that the copyright interest is owned by the Author (or, if more than one author is listed at the beginning of this Agreement, by such authors collectively) and has not been assigned, licensed, or otherwise transferred to any other party. The Author represents and warrants that the author(s) listed at the top of this Agreement are the only authors of the Materials. If more than one author is listed at the top of this Agreement and if any such author has not entered into a separate Article and Video License Agreement with JoVE relating to the Materials, the Author represents and warrants that the Author has been authorized by each of the other such authors to execute this Agreement on his or her behalf and to bind him or her with respect to the terms of this Agreement as if each of them had been a party hereto as an Author. The Author warrants that the use, reproduction, distribution, public or private performance or display, and/or modification of all or any portion of the Materials does not and will not violate, infringe and/or misappropriate the patent, trademark, intellectual property or other rights of any third party. The Author represents and warrants that it has and will continue to comply with all government, institutional and other regulations, including, without limitation all institutional, laboratory, hospital, ethical, human and animal treatment, privacy, and all other rules, regulations, laws, procedures or guidelines, applicable to the Materials, and that all research involving human and animal subjects has been approved by the Author's relevant institutional review board.

10. JoVE Discretion. If the Author requests the assistance of JoVE in producing the Video in the Author's facility, the Author shall ensure that the presence of JoVE employees, agents or independent contractors is in accordance with the relevant regulations of the Author's institution. If more than one author is listed at the beginning of this Agreement, JoVE may, in its sole discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JoVE reserves the right, in its absolute and sole discretion and without giving any reason therefore, to accept or decline any work submitted to JoVE. JoVE and its employees, agents and independent contractors shall have

ARTICLE AND VIDEO LICENSE AGREEMENT

full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JoVE has sole discretion as to the method of making and publishing the Materials, including, without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

11. Indemnification. The Author agrees to indemnify JoVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JoVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JoVE, making of videos by JoVE, or publication in JoVE or elsewhere by JoVE. The Author shall be responsible for, and shall hold JoVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to the making of a video by JoVE its employees, agents or independent contractors. All sterilization, cleanliness or decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's

expense. All indemnifications provided herein shall include JoVE's attorney's fees and costs related to said losses or damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JoVE, its employees, agents or independent contractors.

12. Fees. To cover the cost incurred for publication, JoVE must receive payment before production and publication the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JoVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.

13. Transfer, Governing Law. This Agreement may be assigned by JoVE and shall inure to the benefits of any of JoVE's successors and assignees. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to be one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement required per submission.

CORRESPONDING AUTHOR:

Name:

Department:

Institution:

Article Title:

Signature:

Date:

Please submit a signed and dated copy of this license by one of the following three methods:

- 1) Upload a scanned copy of the document as a pdf on the JoVE submission site;
- 2) Fax the document to +1.866.381.2236;
- 3) Mail the document to JoVE / Attn: JoVE Editorial / 1 Alewife Center #200 / Cambridge, MA 02139

For questions, please email submissions@jove.com or call +1.617.945.9051

Rebuttal to editorial and peer review comments

TITLE:

Labeling of Neuronal Morphology using Custom Diolistic Techniques

AUTHORS:

Hough, Lyon H

Department of Biomedical Sciences

Missouri State University

Springfield, MO, USA

LyonHough@Missouristate.edu

Brown, Michael E.

Department of Pathology

Saint Louis University School of Medicine

Saint Louis, MO, USA

mbrown65@slu.edu

Reviewers' comments:

Reviewer #1:

Manuscript Summary:

The JoVE article by Hough and Brown describes the protocol for diolistic labeling of neurons from fixed tissue. It is This article will be a valuable resource for other researchers, particularly because it provides information to improve upon the existing equipment so that it is adapted for this specific application. The BioRad gene gun was not intended for diolistic labeling, and the authors explain how their equipment will improve random labeling of neurons. The authors clearly explain the improvements of their system over existing techniques and over Golgi staining. I look forward to watching the video. I think that it will be important to provide visual details for the assembly of the equipment and for the production of the Dil coated bullet cartridges.

Major Concerns:

There are no major concerns.

Minor Concerns:

It will be important to provide sufficient detail in the video to demonstrate how the equipment is assembled and how the bullets are placed in the chamber.

Author response to Reviewer comments:

Thank you for your comments and thorough review of our manuscript. We agree that a properly narrated video will facilitate the understanding of equipment set-up and cartridge preparation. To insure this, the authors have...

Reviewer #2:

Manuscript Summary:

In this paper "Labeling of Neuronal Morphology using Custom Diolistic Techniques" Lyon Hough and Michael E. Brown describe a detailed protocol for lipophilic dye delivery transfection (diolistics) using a custom built device. The main purpose of this paper is to present a "low-cost" diolistic protocol with high reliability and reproducibility. The efficacy of this method is supported by referring to recently published work "Effects of developmental hyperserotonemia on the morphology of rat dentate nuclear neurons".

Major Concerns:

- 1.) Although the authors tackled in principle a very interesting technical issue for neuroscience, there are a number of excellent reviews already available the same very topic.
- 2.) The major challenge of this MS is the description of how to utilize an alternative gene-gun apparatus that however has not been developed by these authors and that has already been described by the original authors (Bridgman et al).
- 3.) However, in this MS there is no sufficient technical description for building up such device. Moreover, another interesting aspect is the use of pipettes tips as cartridges, however also in this case no sufficient description was given : "Particle carriers (cartridges) are supported by a 108 plastic ring that fits inside a modified Millipore filter holder (cartridge holder)¹²". How should the modified millipore filter be modified in order to obtain a cartridge holder?
- 4.) The description of neuronal imaging after labelling is very well reported but also not novel. Maybe the authors could have focused their attention to something new such as for instances (as they mention in the MS) the long-time preservation of Dil fluorescence, something that the community is looking for. Or, combining immunofluorescence for synaptic proteins with DiOlistic.
- 5.) Moreover, although the protocol steps are clearly explained in "protocol section", the reproducibility of this methodology might be hampered by inaccurate description of reagents and materials. In fact, some information are missing (i.e. 3.0 um filters, methylene chloride, pipette tips) and should be added in the table of Material/Equipment. Moreover, it is reasonable that pipette tips retention index might affect loading of Dil-/PVP in cartridge and then compromise Diolistic labeling. In the "protocol section" it is unclear if aldehyde are dissolved in PBS and, differently from the above mentioned paper PBS 25mM is used for all procedures.
- 6.) Finally, in the "discussion section" data obtained with this protocol should be compared to that obtained using standard gene-gun or other methodology (e.g. Golgi staining) in the same regions.

Minor Concerns:

- 7.) There are several typos throughout the MS (e.g.: Cartilages should be Cartridges)

Author response to Reviewer comments:

- 1.) The authors recognize the availability of reviews and even an available JoVE video covering the topic of Diolistic labeling. However, the goal of this manuscript is to provide an alternative approach to the methods previously reported. We feel the strength of the discussed method comes from its adaptability to a variety of laboratory settings and research aims, utilizing more readily available equipment than that previously covered. Additionally, previously reported conflicting information concerning reagents, concentrations, and techniques have been clearly addressed in this updated protocol.
- 2.) The device previously reported by Bridgeman and Brown was developed, built and tested by Michael E. Brown, coauthor of this manuscript. This same device was modified and protocol optimized by Lyon H. Hough (author) for its use in diolistic labeling and published results in

“Effects of developmental hyperserotonemia on the morphology of rat dentate nuclear neurons”. Additional devices have been constructed by Hough utilizing the same basic parameters and available equipment.

- 3.) While the construction of the device was not covered at length in this manuscript, the authors feel sufficient information and description is given in the manuscript, and will be covered in the video, for readers to undertake construction a similar device with available resources. The authors agree that all modifications to products such as the Millipore filter holders and pipet tips need to have sufficient explanation to allow for reproducibility. We feel these modifications can be best described in a narrated video. Therefore we have highlighted sections of the manuscript protocol to be included in the video demonstration which specifically cover these areas.
- 4.) The authors have elaborated the statement on line 449 to read “Experimentation on prolonging tissue labeling utilizing submersion in antifreeze solutions and storage in sub-zero temperatures is currently underway with positive preliminary results.” We feel the novelty of the manuscript comes from its description of the preparation and labeling protocol, and supports pervious descriptions of neuronal labeling. The goal in this section is to demonstrate that a reader can adapt the method presented here and produce neuronal labeling results previously only possible with commercialized equipment and proprietary reagents.
- 5.) Thank you for the suggestion, microfilter inserts, pipette tips, and methylene chloride used in this procedure have been added to the included equipment table. The authors have added the statement “Low-retention pipette tips should not be used in this protocol as it may interfere with particle adherence to the inner wall of the tip. The protocol outlined here has utilized standard polypropylene research grade pipette tips (Fisher Scientific) throughout testing.” To line 159. The statement “prepared with 25mM PBS” has been added to line 192.
- 6.) While the authors agree that the comparison between Golgi, standard Diolistic, and our customized method would be an interesting prospect, we find that to be outside the goals of this particular work. The goal of this manuscript is to inform the reader of a reproducible protocol which has been adapted from previous works and to provide representative results which can be expected if said protocol is followed. The authors feel the representative e results and previously published works provide sufficient evidence as to the efficacy of the method.
- 7.) Thank you for your critical review of our manuscript. Typographical errors have been addressed.

27 December 2016

JoVE Editorial Board,
1 Alewife Center, Suite 200,
Cambridge, MA 02140

Dear JoVE Editorial Board,
All editorial suggestions have been addressed.

Thank you

Lyon H. Hough, Ph.D.
Assistant Professor, Biomedical Sciences,
Missouri State University,
Professional Building, 409
Springfield, MO 65897
LyonHough@missouristate.edu
Tel: 417 836-6485

Editor comments:

1. Before we can accept the manuscript, some of the language has to be revised to be original. There are a few lines in the Introduction (second and third paragraph) and Section 5 of the protocol need to be revised. Please see the attached Ithenticate report. We cannot accept previously published text as it poses copyright issues in publication. As section 5 is not highlighted for filming, one option would be to drastically reduce this section and cite a reference instead.

Author response: The sections have been reworded so as to not infringe on copywrite restrictions. Citations continue in the document as needed.

2. Please revise the numbering of the references in order of appearance. Currently, the first cited reference in number 2 in the Introduction. Please do not use references in the short and long abstract.

Author response: References have been reordered as necessary

3. Is it Diolistic, diolistic, or DiOlistic? Also, there continues to be many typos throughout the manuscript. "while and advantage"; "to the well carful not the damage the tissue" Please copy-edit the manuscript.


Author response: The term has been published several different ways in the past. The authors believe the use of diolistic labeling to be the most effective as it closely corresponds to the term biolistic. The term does not warrant capitalization on any grammatical terms and should not be used in quotation as it is no longer deemed a novel method named only by the current author. We have used the un-capitalized diolistic throughout the manuscript except where it begins a sentence.

4. Please remove the commercial branding in the manuscript. Please use the Materials Table to identify the materials and vendor instead: Fisher Scientific, Millipore, Neurolucida, etc.

Author response: commercial branding has been removed. We have elected to leave ImageJ (NIH) in the paper due to its open access nature.

5. For Figure 5, please describe Panels D and E. Is there a typo here? C is 10X, D is 20X, and E is 60X? This is not explained
- 6.

Author response: There was a typo found here which has been removed. The authors find the addition to clarify any confusion.



[Click here to access/download](#)

Supplemental File (as requested by JoVE)
RightsLink Printable License.pdf

