

Journal of Visualized Experiments

On the preparation and testing of fuel cell catalysts using the thin film rotating disc electrode method --Manuscript Draft--

Article Type:	Invited Methods Article - JoVE Produced Video
Manuscript Number:	JoVE57105R2
Full Title:	On the preparation and testing of fuel cell catalysts using the thin film rotating disc electrode method
Keywords:	Oxygen reduction reaction; rotating disk electrode; polymer electrolyte fuel cells; electrocatalysts; high surface area catalysts; catalyst inks; catalyst thin films
Corresponding Author:	Matthias Arenz Universitat Bern Bern, Bern SWITZERLAND
Corresponding Author's Institution:	Universitat Bern
Corresponding Author E-Mail:	matthias.arenz@dcb.unibe.ch
First Author:	Masanori Inaba
Other Authors:	Masanori Inaba
	Jonathan Quinson
	Jan Rudolf Bucher
Author Comments:	Dear Editor, we changed the manuscript according the the comments and suggestions. We submit one manuscript with the answers and changes tracked - as recommended as well as a clean manuscript with all changes accepted and with all comments removed. Best regards, Matthias Arenz
Additional Information:	
Question	Response
If this article needs to be "in-press" by a certain date, please indicate the date below and explain in your cover letter.	

Prof. Dr. Matthias Arenz

Department of Chemistry and Biochemistry
University of Bern
Freiestrasse 3
CH - 3012 Bern
Tel +41 31 631 53 84
Email: matthias.arenz@dcb.unibe.ch



^b
UNIVERSITÄT
BERN

Bern, 2nd August 2017

Dear Editor,

Please find attached the manuscript entitled

“On the preparation and testing of fuel cell catalysts using the thin film rotating disc electrode method”

by Inaba *et al.* for publication in Jove. All authors have seen and approved the final version of the manuscript.

Fuel cell technology is a promising piece for future energy schemes relying on renewable energy sources. They can efficiently re-convert the energy carrier hydrogen, which is produced by water electrolysis using renewable energy source like wind or solar power.

In our work, we present a tutorial for the preparation and testing of a Pt/C fuel cell catalyst using the thin-film rotating disk electrode technique. Although this method is widely applied in the scientific community, there is still a need for common procedures and best practice as several research groups outlined in recent years. We believe that Jove provides an ideal platform for this purpose.

Sincerely yours,

Matthias Arenz

TITLE:

On the Preparation and Testing of Fuel Cell Catalysts Using the Thin Film Rotating Disk Electrode Method

AUTHORS AND AFFILIATIONS:

Masanori Inaba¹, Jonathan Quinson¹, Jan Bucher², Matthias Arenz^{1,2}

¹Nano-Science Center, Department of Chemistry, University of Copenhagen, Copenhagen, Denmark

²Department of Chemistry and Biochemistry, University of Bern, Bern, Switzerland

EMAIL ADDRESSES:

Masanori Inaba (masanori@chem.ku.dk)

Jonathan Quinson (jonathan.quinson@chem.ku.dk)

Jan Bucher (jan.bucher@dcb.unibe.ch)

CORRESPONDING AUTHOR:

Matthias Arenz (Matthias.arenz@dcb.unibe.ch)

Tel: +41-31-631-53-84

KEYWORDS:

Oxygen Reduction Reaction, Rotating Disk Electrode, Polymer Electrolyte Fuel Cells, Electrocatalysts, High Surface Area Catalysts, Catalyst Inks, Catalyst Thin Films

SUMMARY:

Preparing and testing Pt/C fuel cell catalysts is subject to continuous discussion in the scientific community with respect to reproducibility and best practice. With the presented work, we intend to present a step-by-step tutorial to make and test Pt/C catalysts, which can serve as benchmark for novel catalyst systems.

ABSTRACT:

We present a step-by-step tutorial to prepare proton exchange membrane fuel cell (PEMFC) catalysts, consisting of Pt nanoparticles (NPs) supported on a high surface area carbon, and to test their performance in thin film rotating disk electrode (TF-RDE) measurements. The TF-RDE methodology is widely used for catalyst screening; nevertheless, the measured performance sometimes considerably differs among research groups. These uncertainties impede the advancement of new catalyst materials and, consequently, several authors discussed possible best practice methods and the importance of benchmarking.

The visual tutorial highlights possible pitfalls in the TF-RDE testing of Pt/C catalysts. A synthesis and testing protocol to assess standard Pt/C catalysts is introduced that can be used together with polycrystalline Pt disks as benchmark catalysts. In particular, this study highlights how the properties of the catalyst film on the glassy carbon (GC) electrode influence the measured performance in TF-RDE testing. To obtain thin, homogeneous catalyst films, not only the catalyst

preparation, but also the ink deposition and drying procedures are essential. It is demonstrated that an adjustment of the ink's pH might be necessary, and how simple control measurements can be used to check film quality. Once reproducible TF-RDE measurements are obtained, determining the Pt loading on the catalyst support (expressed as Pt wt%) and the electrochemical surface area is necessary to normalize the determined reaction rates to either surface area or Pt mass. For the surface area determination, so-called CO stripping, or the determination of the hydrogen underpotential deposition (H_{upd}) charge, are standard. For the determination of the Pt loading, a straightforward and cheap procedure using digestion in *aqua regia* with subsequent conversion of Pt(IV) to Pt(II) and UV-vis measurements is introduced.

INTRODUCTION:

The increasing use of renewable energy sources like wind or solar power require the ability to store and re-convert large amounts of electrical energy. In this respect, the application of chemical energy carriers, such as hydrogen, is a promising pathway¹. Water electrolysis can produce hydrogen from excess electrical energy, whereas PEMFCs efficiently re-convert the hydrogen into electricity. However, for large scale applications, such as automobiles, the precious group metal (PGM) content in PEMFCs needs to be significantly reduced to lower the costs. In the scientific literature, several catalyst systems² are presented, which are reported to exhibit the potential to meet these requirements based on TF-RDE half-cell testing³.

TF-RDE measurements became extremely popular in PEMFC research and are now used as a standard method to compare different PEMFC catalysts. However, for similar catalysts, significant differences in oxygen reduction reaction (ORR) activity or electrochemical surface areas (ECSA) are often reported by different research groups⁴. As a consequence, several research groups worked on improving experimental procedures and defining best practice procedures of how to test PEMFC catalysts⁵. At first, factors like catalyst loading, contaminations from the reference electrode, iR compensation, etc. were highlighted⁶⁻⁸; whereas, in recent years, the influence of the catalyst film's properties on performance came into focus^{9,10}. For example, it could be shown that both ORR activity and ECSA depend on the catalyst ink composition, which in turn influences the homogeneity of the catalyst film on the glassy carbon (GC) electrode used as substrate. The observed increase in the measured Pt surface area upon the addition of isopropanol to the catalyst ink is especially surprising at first sight, but these results point towards the importance of carbon support wetting. Another important (and correlated) factor is mass transport. The TF-RDE is limited to low current densities due to the low gas solubility in the electrolyte¹¹⁻¹³. Therefore, it is often assumed that in such extremely thin catalyst layers (*ca.* 1 μm), reactant mass transport plays a limited role. Nevertheless, in very recent work, it was shown that adjustment of the ink pH, leads to a considerable increase in ORR activity in home-made catalysts¹⁴. These examples highlight that in TF-RDE measurements, careful controls are necessary, and that depending on the catalyst, it might be difficult to define a single standard testing protocol/recipe.

Therefore, in the presented work, we discuss a step-by-step procedure for synthesizing and testing a standard Pt/C fuel cell catalyst. The procedure includes catalyst synthesis, its characterization, and the ink preparation and adjustment as well as the electrochemical TF-RDE

measurements. The purpose of the procedure is to enhance consciousness concerning certain sources of errors and experimental pitfalls, as well as to provide a recipe for a possible Pt/C benchmark catalyst.

PROTOCOL:

1. Synthesis of 50 wt% Pt/C Catalysts

1.1. Colloidal synthesis of Pt nanoparticles (NPs)

1.1.1. To obtain a colloidal suspension of 2 nm Pt NPs, mix 4 mL of a solution of NaOH at 0.4 M in ethylene glycol (EG) with 4 mL of $\text{H}_2\text{PtCl}_6 \cdot 6\text{H}_2\text{O}$ at 40 mM in EG in a microwave reaction vessel.

1.1.2. Heat the mixture for 3 min at 160 °C with the microwave reactor (dynamic mode, heating power: 100 W, stirring: low).

1.2. Immobilization of Pt NPs onto carbon supports

1.2.1. Add 30 mL of 1 M HCl solution to 7.3 mL of the colloidal Pt NPs suspension for precipitation. Centrifuge the mixture at a relative centrifugal force of 2,900 x g (5,000 rotations per minute (rpm)) for 5 min and discard the supernatant solvent. Repeat these washing/centrifugation steps twice.

Note: All chemicals should be collected and disposed of according to the rules at the local institution.

1.2.2. Re-disperse the Pt NPs in 7 mL of acetone. Disperse 27.5 mg of carbon black in 1 mL of acetone. Mix both acetone dispersions and evaporate the acetone in a rotary evaporator (40 °C, 200 mbar) or Schlenk line (40 °C, 200 mbar) until completely dried.

1.2.3. Dry the obtained Pt/C catalyst powder in an oven (120 °C) over night.

1.2.4. Add deionized (DI) water to the Pt/C catalyst powder and sonicate (35 kHz, 160 Watt) it in an ultrasonic bath filled with cold water (<5 °C) for 3 min. Filter out and wash the catalyst with DI water on filter paper (4-7 µm). Dry the catalyst powder in a vacuum oven (100-120 °C, 10 kPa).

Note: This procedure leads to 55 mg of a Pt/C catalyst with a Pt wt% of 50 and a Pt particle size of ca. 2 nm. To prepare larger amounts of Pt/C catalyst, scale up the respective amounts in a ratio of 1:1. By varying the amount of Pt NPs to carbon black, the Pt loading can be adjusted between 10 and 70-80 wt% without significant NP agglomeration.

1.2.5. To make a catalyst ink, mix 6.3 mg of 50 wt% Pt/C catalyst powder with 6 mL of DI water and 2 mL of isopropanol (IPA). Add a small amount of 1 M KOH solution (~10 µL) to adjust the pH of the ink to be around 10. Place the glass vial containing the mixture in an ultrasonic bath filled

with cold water (<5 °C) and sonicate (35 kHz, 160 Watt) it for 15 min.

Note: Do not add IPA to the catalyst first. If alcohol is added to a dried Pt/C catalyst powder, the catalyst catches fire.

2. Characterization of Pt NPs and Pt/C Catalysts

2.1. Transmission electron microscope (TEM) characterization of Pt NPs

2.1.1. Add 30 mL of 1 M HCl solution to 7.3 mL of the colloidal Pt NPs suspension for precipitation. Centrifuge the mixture at a relative centrifugal force of 2,900 x g (5,000 rotations per minute (rpm)) for 5 min and discard the supernatant solvent. Repeat these washing/centrifugation steps twice.

2.1.2. Take a drop (*ca.* 100 µL) of the Pt NP suspension and re-disperse it in 1 mL of acetone.

2.1.3. Put a drop of the diluted, re-dispersed NPs on a TEM grid (with a normal carbon film) placed on a filter paper to absorb excess solvent and let the grid dry in air.

2.1.4. Analyze the particle size distribution by imaging at least 5 different areas of the grid, using at least 3 different magnifications (*e.g.*, 100,000X, 300,000X, 400,000X). Measure the size (and related size distribution) of at least 300 NPs from the micrographs.

2.2. TEM characterization of Pt/C catalyst

2.2.1. Dilute the catalyst ink (described above) twice with 100% IPA. Put a drop of the diluted catalyst ink on a TEM grid (with a holey or lacey carbon film) placed on a filter paper to absorb the excess of solvent and let the grid dry in air.

2.2.2. Analyze the particle size distribution by imaging at least 5 different areas of the grid, using at least 3 different magnifications (*e.g.*, 100,000X, 300,000X, 400,000X). Measure the size (and related size distribution) of at least 300 NPs from the micrographs.

Note: TEM images with better contrast are obtained by observing the Pt/C catalyst particles, which do not overlap with a carbon film of the grid. This is the advantage of using TEM grids with a holey or lacey carbon film here.

2.3. Digestion of Pt NPs by *aqua regia*

2.3.1. Dry the Pt/C catalyst in a vacuum oven (80 °C, 10 kPa) overnight. Weigh the catalyst powder in a ceramic crucible.

2.3.2. Heat the catalyst powder in the crucible of a muffle furnace (air, 900 °C, 30 min) to burn off the carbon support. Put a lid on the crucible to prevent the catalyst powder from flying during

the heating.

2.3.3. Add 4 mL of *aqua regia* (mixture of 30% HCl and 65% HNO₃ in a volume ratio of 3:1) into the crucible and heat the mixture on a hot plate at *ca.* 80 °C for 2 h.

CAUTION: *Aqua regia* is an aggressive acid and should be handled with care.

2.3.4. Adjust the final volume to 10 mL with DI water.

2.3.5. Determine the concentration of Pt in the *aqua regia* solution by an inductively coupled plasma-mass spectrometry (ICP-MS)¹⁵ or an ultraviolet-visible spectroscopy (UV-Vis, described below).

2.4. Determination of Pt loading via UV-Vis

2.4.1. Mix 1 mL of the *aqua regia* sample with 0.25 mL of 1 M SnCl₂ prepared in 4 M HCl, and 0.75 mL of 2 M HCl in a quartz cuvette (the concentration of Pt in the *aqua regia* solution should be around 5 ppm). Stir the mixture well in the cuvette using a small magnetic stir bar.

2.4.2. Measure the UV-vis spectrum of a mixture of 0.25 mL of 1 M SnCl₂ in 4 M HCl and 1.75 mL of 2 M HCl as a baseline.

2.4.3. Measure the UV-Vis spectrum of the sample mixture (700 nm – 350 nm).

2.4.4. Add 5 µL of Pt standard solution (1,000 ppm) to the sample solution, stir it well, and measure the UV-Vis spectrum. Repeat the addition of the Pt standard solution and the measurement 4 times.

2.4.5. Subtract the absorbance at 680 nm (baseline) from the absorbance at 402 nm (peak top), and plot it against the concentration of the added Pt. Determine the concentration of Pt in the *aqua regia* sample from the x-intercept of the calibration curve.

3. Rotating Disc Electrode (RDE) Measurement of Pt/C Catalysts

3.1. Characterization of the catalyst ink via light scattering measurements

3.1.1. Measure the pH of the catalyst ink with a pH meter.

3.1.2. Dilute the catalyst ink 50 times by adding IPA/H₂O (1:3, v:v) mixed solvent. Add acid or base to keep the pH of the catalyst ink constant.

3.1.3. Measure the zeta potential of the 50X diluted sample using the electrophoretic light scattering (ELS)¹⁶ mode of the light scattering device.

3.1.4. Measure the size of aggregates in the diluted sample using dynamic light scattering (DLS)¹⁷ of the light scattering device.

3.2. Catalyst thin film fabrication

3.2.1. Polish a glassy carbon (GC) disk electrode (5 mm in diameter) to mirror finish using alumina paste (first 0.3 μm , then 0.05 μm), and clean it ultrasonically in DI water.

3.2.2. Pipette a 5 μL aliquot of the catalyst ink onto the GC electrode leading to a Pt loading of 10 $\mu\text{g}_{\text{Pt}} \text{cm}^{-2}$. Dry the ink under Ar gas flow humidified with a mixture of IPA and DI water (17:3, v:v) in a bubbler.

Note: The drying condition can be optimized by changing the ratio between IPA and DI water in the bubbler.

3.2.3. Confirm that the surface of the GC electrode is uniformly covered with the catalyst thin film by a CCD camera.

3.3. Electrochemical measurements

3.3.1. Employ an electrochemical glass or polytetrafluoroethylene (PTFE) cell based on a three-compartment configuration in all electrochemical measurements. Use the modified GC disk electrode, platinum mesh, and a trapped hydrogen/saturated calomel electrode (SCE) as working electrode, a counter electrode, and a reference electrode, respectively.

Note: A sub-compartment of the cell for the reference electrode needs to be separated from a main compartment by a perfluorinated membrane (*e.g.*, nafion membrane) when a SCE is used to avoid the diffusion of chloride ions.

3.3.2. Soak the cell and components in mixed acid ($\text{H}_2\text{SO}_4:\text{HNO}_3 = 1:1$, v:v) overnight. Subsequently, rinse the cell and other components thoroughly with DI water, and boil them in DI water twice.

Note: Clean the cell with the mixed acid at least once a month (depending on the frequency of use).

CAUTION: Handle the mixed acid with care. Between electrochemical experiments, boil and keep the cell in DI water covered by a lid.

3.3.3. Pour 0.1 M HClO_4 electrolyte solution into the electrochemical cell. Place the 3 electrodes (the working electrode, the counter electrode, and the reference electrode) in each compartment of the cell, and connect the electrodes to the potentiostat.

3.3.4. Measure the resistance in the electrolyte between the working electrode and the Luggin

capillary (iR-drop; $\sim 30\ \Omega$), and compensate it during the whole measurement process by using the analog positive feedback scheme of the potentiostat. The resulting effective solution resistance is less than $3\ \Omega$.

Note: If the potentiostat does not have an analog feedback scheme, apply a software based iR-drop compensation.

3.3.5. Calibrate the reference electrode. Purge the electrolyte with H_2 gas and scan the potential at a scan rate of $50\ \text{mV s}^{-1}$. Determine the potential where the current density is equal to $0\ \text{mA cm}^{-2}$ as $0\ \text{V}$ versus reversible hydrogen electrode (RHE).

3.3.6. Purge the electrolyte with Ar gas (*ca.* $0.3\ \text{L min}^{-1}$). Clean the catalysts by potential cycles between 0.05 and $1.20\ \text{V}$ at a scan rate of $0.05\ \text{V s}^{-1}$, until a stable cyclic voltammogram is observed (~ 50 cycles).

3.3.7. Perform background measurements.

3.3.7.1. Record a cyclic voltammogram (CV) as a background measurement. Purge the electrolyte with Ar gas (*ca.* $0.3\ \text{L min}^{-1}$). Sweep the potential between 0.05 and $1.10\ \text{V}$ at a scan rate of $0.05\ \text{V s}^{-1}$.

3.3.8. Perform ORR activity measurement.

3.3.8.1. Record a cyclic voltammogram in the O_2 -purged electrolyte (*ca.* $0.3\ \text{L min}^{-1}$). Sweep the potential between 0.05 and $1.10\ \text{V}$ at a scan rate of $0.05\ \text{V s}^{-1}$ and a rotation speed of $1,600\ \text{rpm}$.

3.3.9. Perform CO stripping voltammetry.

3.3.9.1. Hold the working electrode at $0.05\ \text{V}$ during a CO purge (*ca.* $0.3\ \text{L min}^{-1}$) through the electrolyte for $5\ \text{min}$ followed by an Ar purge (*ca.* $0.3\ \text{L min}^{-1}$) for $10\ \text{min}$. Subsequently sweep the potential from $0.05\ \text{V}$ to $1.10\ \text{V}$ at a scan rate of $0.05\ \text{V s}^{-1}$.

CAUTION: Check if the use of CO gas is allowed; if so, this work must be done in a fume hood and safety measures need to be followed.

3.3.10. After the electrochemical measurements, transfer the catalyst films on the GC electrode onto a tissue paper by pressing the electrodes on the paper wetted with water and photograph the transferred catalyst films.

3.4. Data analysis

3.4.1. Calculate the electrochemical surface area (ECSA) of the catalyst from the charge Q_{CO} of the CO oxidation peak, taking into account the scan speed ν , the Pt loading L , and using a conversion coefficient of $390\ \mu\text{C cm}^{-2}_{\text{Pt}}$ (two electron process)⁶.

Note: Alternatively, the ECSA can be calculated from the H_{upd} (Q_H) area of a cyclic voltammogram recorded in Ar-purged electrolyte using a conversion coefficient of $195 \mu\text{C cm}^{-2}_{\text{Pt}}$ (one electron process).

$$Q_{\text{CO,H}} = \int \frac{j dE}{v}$$

$$\text{ECSA} = \frac{Q_{\text{CO}}}{390 \mu\text{C cm}_{\text{Pt}}^{-2} 10 \mu\text{g}_{\text{Pt}} \text{cm}_{\text{geo}}^{-2}}$$

3.4.2. Correct the nonfaradaic background by subtracting the cyclic voltammogram recorded in Ar-purged electrolyte from the voltammogram recorded in O_2 -purged electrolyte. Determine the ORR kinetic current density, j_k , using the anodic scan of the corrected voltammogram and Koutecky-Levich equation:

$$\frac{1}{j} = \frac{1}{j_k} + \frac{1}{j_d}$$

Note: Here j and j_d are representative of the measured current density and diffusion limiting current density, respectively.

REPRESENTATIVE RESULTS:

A dark brown, colloidal suspension of Pt NPs is obtained from the protocol section 1.1 (**Figure 1**, left). The colloidal suspension is so stable that it can be stored for more than 1 month without any precipitation. Larger Pt NPs are synthesized by decreasing the NaOH concentration¹⁸. However, the colloidal suspension becomes less stable by decreasing the NaOH concentration. As an extreme example, completely agglomerated Pt NPs are obtained when the Pt precursor is heated in EG without NaOH (**Figure 1**, right).

[Place Figure 1 here]

Well-dispersed Pt/C catalyst is obtained from the protocol section 1.2, as shown in **Figure 2a**. When a less stable colloidal suspension, *e.g.*, of larger Pt NPs is used, significant agglomeration of Pt NPs on the carbon support might be observed (**Figure 2b**).

[Place Figure 2 here]

Representative results of the UV-Vis measurement for the determination of the Pt concentration in *aqua regia* (protocol 2.4) are shown in **Figure 3a**. When SnCl_2 is added to the *aqua regia* sample, the Pt in the *aqua regia* is reduced from Pt(IV) to Pt(II), leading to a yellow colored solution. **Figure 3b** is the calibration curve obtained from the spectra in **Figure 3a**. From this calibration curve, the Pt concentration in the sample mixture is determined to be 3.68 ppm.

[Place Figure 3 here]

Figures 4a-b are representative examples of a homogeneous catalyst thin film fabricated on a GC electrode by following the protocol 3.2. The whole surface of the GC electrode is covered by the catalyst film uniformly without any significant agglomeration. **Figure 4c-d** are typical inhomogeneous catalyst thin films fabricated by drying the catalyst ink in air. The catalyst agglomerates on the edge of the GC electrode and forms a so-called coffee ring¹⁹.

[Place Figure 4 here]

Figure 5 is an example of a cyclic voltammogram measured in H₂ saturated electrolyte for calibration of the reference electrode potential against RHE (protocol 3.3.5). The average of the potentials where the current density is 0 mA cm⁻² both in the positive going scan and the negative going scan is defined as '0 V vs. RHE'.

[Place Figure 5 here]

Shown in **Figure 6** are cyclic voltammograms in Ar saturated electrolyte and linear sweep voltammograms (LSV) in O₂ saturated electrolyte for the 50 wt% Pt/Vulcan catalyst obtained from the RDE measurement (protocol 3.3.7 and 3.3.8). After sufficient cleaning cycles in the protocol step 3.4.6, the red CV in **Figure 6a** is obtained. When the catalyst is not sufficiently cleaned, the shoulder of the Pt oxidation peak around 0.8 V is less sharp than that of the well-cleaned catalyst (**Figure 6a**, gray CV). The shape of the LSV in the O₂ saturated electrolyte is highly sensitive to the quality of the catalyst thin film (**Figure 6b**). When the catalyst thin film is homogeneous like in **Figure 4a**, the O₂ diffusion limiting current density (below 0.8 V) is observed around -6 mA cm⁻², and the shoulder of the LSV curve around 0.8 V is sharp (**Figure 6b**, red). On the other hand, the O₂ diffusion limiting current density is smaller and the shoulder of the LSV curve is less sharp when the catalyst thin film is nonhomogeneous, as in **Figure 4c**, or the surface of the GC electrode is not fully covered with the catalyst thin film.

[Place Figure 6 here]

Shown in **Figure 7** are a cyclic voltammogram recorded in Ar-saturated electrolyte with the H_{upd} area indicated (**a**) and a CO stripping voltammogram with the stripping charge indicated (**b**). Both measurements were obtained for a 50 wt% Pt/Vulcan catalyst. From the area under the peaks, the Pt surface area (ECSA) is calculated, see step 3.4.1.

[Place Figure 7 here]

Results of the RDE measurement of the 50 wt% Pt/Vulcan catalyst (ECSA, ORR specific activity and mass activity at 0.9 V_{RHE}) are summarized in **Table 1**. The results for the commercial Pt/C catalyst with similar Pt particle size and similar Pt loading are shown in the table as well.

[Place Table 1 here]

FIGURE AND TABLE LEGENDS:

Figure 1: Photograph of the colloidal suspensions of Pt NPs in EG. (Left) Well-dispersed Pt NPs were obtained by following the protocol 1.1. (Right) Agglomerated Pt NPs obtained by heating the Pt precursor in EG without NaOH.

Figure 2: TEM images of 50 wt.% Pt/Vulcan catalysts. (a) 2 nm Pt NPs are well-dispersed on the carbon supports. (b) 4 nm Pt NPs are agglomerated on the carbon supports.

Figure 3: Representative results of the UV-Vis measurement for determination of Pt concentration in aqua regia. (a) UV-vis spectra of *aqua regia* sample with SnCl_2 before and after adding the Pt standard solution (5–20 μL). (b) Calibration curve obtained by plotting differences between the absorbance at 402 nm and the absorbance at 680 nm against the concentration of the added Pt. The Pt concentration in the sample solution is determined from the x-intercept of the calibration curve.

Figure 4: Photographs and SEM images of homogeneous and inhomogeneous catalyst thin film fabricated on a GC electrode. (a,b) Homogeneous catalyst thin film fabricated by following the protocol. (c,d) Inhomogeneous catalyst thin film with a “coffee ring” fabricated by drying the catalyst ink in air. (a) and (c) are whole pictures of the catalyst thin films. (b) and (d) are SEM images of the catalyst thin films taken at the edge of the GC disks.

Figure 5: Cyclic voltammogram in H_2 saturated electrolyte for the calibration of the potential the reference electrode against RHE (50 mV s^{-1} , 1600 rpm). The average of the potentials where the current density is 0 mA cm^{-2} both in the positive going scan and the negative going scan is defined as 0 V vs. reversible hydrogen electrode (RHE) (see dotted red lines).

Figure 6: Representative examples for “good” and “bad” voltammograms for the 50 wt% Pt/Vulcan catalyst. (a) Cyclic voltammograms in Ar saturated electrolyte (50 mV s^{-1}). (b) LSVs in O_2 saturated electrolyte (1,600 rpm, 50 mV/s , positive going scan, background subtracted, iR compensated).

Figure 7: Representative cyclic voltammogram (a) and CO stripping voltammogram (b) for the 50 wt% Pt/Vulcan catalyst (50 mV s^{-1}). The H_{upd} charge and the CO stripping charge are highlighted in sky blue and pink, respectively.

Figure 8: Photographs of the catalyst thin film fabricated on GC electrodes with various drying conditions. The ratio between isopropanol (IPA) and H_2O in the bubbler was changed to optimize drying conditions. (a) 100% IPA, (b) 90% IPA / 10% H_2O , (c) 80% IPA / 20% H_2O , (d) 70% IPA / 30% H_2O . (c) features the best drying conditions in this case, as the most homogeneous catalyst thin film is fabrication on the GC electrode.

Figure 9: Photographs of the catalyst inks with different pH: (left) pH \approx 4, (middle) pH \approx 7, (right) pH \approx 10. The photographs are taken 1 week after the sonication. The dispersion of the acidic ink is much less stable than the neutral and the alkaline ink.

Figure 10: Photographs of the catalyst films transferred on a paper (a-c) and SEM images of the catalyst films on GC electrodes (d-f). (a,d) Catalyst film fabricated from acidic (pH \approx 4) ink, (b,e) Catalyst film fabricated from neutral (pH \approx 7) ink, and (c,f) Catalyst film fabricated from alkaline (pH \approx 10) ink.

Figure 11: Influence of iR-compensation on the measurement. (a,b) Screenshots of the interface of “EC4 DAQ” program with and without iR-compensation. (c) LSVs measured in O₂ saturated electrolyte with iR-compensation (red) and without iR-compensation (gray) (50 mV s⁻¹, 1,600 rpm, back ground subtracted).

Figure 12: Overview of the experimental steps. Corresponding characterization methods and controllable parameters for each experimental step in this protocol are shown. TEM = transmission electron microscopy, SAXS = small angle X-ray scattering, EXAFS = extended X-ray absorption fine structure, XANES = X-ray absorption fine structure, ICP-MS = inductively coupled plasma mass spectrometry, DLS = dynamic light scattering, CCD = charge-coupled device, SEM = scanning electron microscopy.

Table 1: Summary of results of the TF-RDE measurements for the 50 wt% Pt/Vulcan catalyst and the commercial 46 wt% Pt/C catalyst. The measured values are shown with standard deviation. ECSA = electrochemical surface area, ORR = oxygen reduction reaction, TTK = Tanaka Kikinzoku Kogyo.

DISCUSSION:

It is well known that the ORR activity of high surface area catalysts measured by the TF-RDE technique highly depends on the uniformity of the catalyst thin film^{6,9,10,14,19}. Several groups have reported fabrication methods for homogeneous catalyst thin films on GC electrodes and researchers should carefully optimize their drying method when entering this research area²⁰. Using rotational methods¹⁹ allows more flexibility in coating of thicker catalyst films, whereas stationary methods have the advantage that multiple electrodes can be prepared at the same time. As an example of a stationary drying method, Shinozaki *et al.* recently reported that uniform catalyst thin films are fabricated by drying the catalyst ink in an IPA atmosphere¹⁰. The fabrication of the catalyst thin films in this protocol is based on their method. However, the catalyst inks are dried in a humidified gas flow instead of a stationary atmosphere (protocol 3.3.2). The advantage of this modified method is that the drying condition can easily be adjusted by changing the ratio between IPA and water in the bubbler. **Figure 8** indicates that the uniformity of the catalyst thin film is optimized by changing the condition of the humidification.

[Place Figure 8 here]

We also found that the pH of the catalyst inks is an important parameter to be optimized to

obtain homogeneous catalyst thin films, as reported in the literature¹⁴. **Figure 9** demonstrates the varying stability of catalyst inks with different pH. Since HCl is used for the washing of Pt NPs in the protocol 1.3.3, the catalyst ink usually becomes acidic without the addition of KOH. The acidic catalyst ink is not particularly stable and most of the Pt/C catalyst particles settle to the bottom 1 week after the sonication. The neutral ink is more stable than the acidic ink, even though some precipitates are seen in the bottom. The alkaline ink is the most stable and no precipitate is seen 1 week after the sonication. This pH dependence of the catalyst ink stability is explained by the magnitude of the zeta potential, which becomes larger with increasing pH¹⁴.

[Place Figure 9 here]

Not only the stability of the catalyst inks, but also the uniformity of the obtained catalyst thin film depends on the pH of the catalyst ink. Highly agglomerated catalyst thin films are obtained when acidic inks are used (**Figure 10a,d**). Although there is no significant difference visible to the eye between the catalyst thin films fabricated from the neutral ink and from the alkaline ink (**Figure 10 b,c**), SEM images reveal that there are some agglomerates in the catalyst thin film obtained from the neutral ink, whereas no significant agglomerates are seen in the catalyst thin film obtained from the alkaline ink (**Figure 10e,f**).

[Place Figure 10 here]

In addition to the catalyst film quality, the scan rate⁴ and the compensation of the cell resistance influences the ORR activity determination⁸. In our measurement setup, the cell resistance without iR-compensation is usually around 30 Ω . It is decreased to be less than 3 Ω by using the iR-compensation of the potentiostat (**Figure 11a,b**). The LSVs measured in O₂ saturated electrolyte with and without iR-compensation are compared in **Figure 11c**. A noticeable potential shift due to the cell resistance is seen at around 0.9 V_{RHE}, where ORR activity is evaluated. As a scan rate for the LSV, we chose 50 mV s⁻¹, while other groups⁵ prefer 20 mV s⁻¹ and membrane electrode assemblies (MEAs) are often tested under steady state conditions⁴. In general, it can be stated for TF-RDE studies that the lower the scan rate, the more the measurement becomes susceptible to possible contaminations. In MEA testing, considerably higher currents are applied. In particular for automotive applications, high power performance is particularly interesting²¹. Scanning the potential would lead to significant errors if no online iR compensation is applied.

Due to these differences, direct predictions of MEA performance based on TF-RDE measurements should be taken with caution. The TF-RDE should be seen as a fast method to screen test the intrinsic ORR activity of PEMFC catalysts, rather than an alternative to MEA testing.

[Place Figure 11 here]

The overall procedure is summarized in **Figure 12**. In addition to the discussed standard characterization methods, the obtained colloidal Pt NPs suspension and Pt/C catalyst can also be investigated by more advanced methods such as small angle X-ray scattering (SAXS)²² or X-ray absorption spectroscopy (XAS)²³.

[Place Figure 12 here]

ACKNOWLEDGMENTS:

J.Q. and M.A. acknowledge support from the Villum Foundation in form of a block stipend. M.I. and M.A. acknowledge support from Toyota Central R&D Labs., Inc. J.Q. has received funding from the European Union's Horizon 2020 research and innovation program under the Marie Skłodowska-Curie grant agreement No 703366.

DISCLOSURES:


The authors have nothing to disclose.



REFERENCES:



- 1 Stephens, I. E. L., Rossmeisl, J. & Chorkendorff, I. Toward sustainable fuel cells. *Science*. **354** (6318), 1378-1379 (2016).
- 2 Chen, C. *et al.* Highly Crystalline Multimetallic Nanoframes with Three-Dimensional Electrocatalytic Surfaces. *Science*. **343** (6177), 1339-1343 (2014).
- 3 Schmidt, T. J. *et al.* Characterization of high-surface area electrocatalysts using a rotating disk electrode configuration. *J. Electrochem. Soc.* **145** (7), 2354-2358 (1998).
- 4 Gasteiger, H. A., Kocha, S. S., Sompalli, B. & Wagner, F. T. Activity benchmarks and requirements for Pt, Pt-alloy, and non-Pt oxygen reduction catalysts for PEMFCs. *Appl. Catal. B-Environm.* **56** (1-2), 9-35 (2005).
- 5 Kocha, S. S. *et al.* Best Practices and Testing Protocols for Benchmarking ORR Activities of Fuel Cell Electrocatalysts Using Rotating Disk Electrode. *Electrocatalysis*. **8** (4), 366-374 (2017).
- 6 Mayrhofer, K. J. J. *et al.* Measurement of oxygen reduction activities via the rotating disc electrode method: From Pt model surfaces to carbon-supported high surface area catalysts. *Electrochim. Acta* **53** (7), 3181-3188 (2008).
- 7 Mayrhofer, K. J. J., Ashton, S. J., Kreuzer, J. & Arenz, M. An Electrochemical Cell Configuration Incorporating an Ion Conducting Membrane Separator between Reference and Working Electrode. *Intern. J. Electrochem. Sci.* **4** (1), 1-8 (2009).
- 8 Nesselberger, M. *et al.* The Particle Size Effect on the Oxygen Reduction Reaction Activity of Pt Catalysts: Influence of Electrolyte and Relation to Single Crystal Models. *J. Am. Chem. Soc.* **133** (43), 17428-17433 (2011).
- 9 Ke, K., Hiroshima, K., Kamitaka, Y., Hatanaka, T. & Morimoto, Y. An accurate evaluation for the activity of nano-sized electrocatalysts by a thin-film rotating disk electrode: Oxygen reduction on Pt/C. *Electrochim. Acta* **72** 120-128 (2012).
- 10 Shinozaki, K., Zack, J. W., Richards, R. M., Pivovar, B. S. & Kocha, S. S. Oxygen Reduction Reaction Measurements on Platinum Electrocatalysts Utilizing Rotating Disk Electrode Technique. *J. Electrochem. Soc.* **162** (10), F1144-F1158 (2015).
- 11 Zalis, C. M., Sharman, J., Wright, E. & Kucernak, A. R. Properties of the hydrogen oxidation reaction on Pt/C catalysts at optimised high mass transport conditions and its relevance to the anode reaction in PEMFCs and cathode reactions in electrolyzers. *Electrochim. Acta* **176**, 763-776 (2015).

- 12 Wiberg, G. K. H., Fleige, M. J. & Arenz, M. Design and test of a flexible electrochemical setup for measurements in aqueous electrolyte solutions at elevated temperature and pressure. *Rev. Sci. Instr.* **85** (8), 085105 (2014).
- 13 Wiberg, G. K. H., Fleige, M. & Arenz, M. Gas diffusion electrode setup for catalyst testing in concentrated phosphoric acid at elevated temperatures. *R Rev. Sci. Instr.* **86** (2), 024102 (2015).
- 14 Inaba, M., Quinson, J. & Arenz, M. pH matters: The influence of the catalyst ink on the oxygen reduction activity determined in thin film rotating disk electrode measurements. *J. Power Sources* **353** 19-27 (2017).
- 15 Cherevko, S. *et al.* Dissolution of Platinum in the Operational Range of Fuel Cells. *ChemElectroChem.* **2** (10), 1471-1478 (2015).
- 16 Ware, B. R. Electrophoretic light scattering. *Adv. Colloid Interface Sci.* **4** (1), 1-44, (1974).
- 17 Berne, B. J. & Pecora, R. *Dynamic Light Scattering, with application to Chemistry, Biology and Physics*. (Dover Publication Inc, 2000).
- 18 Baranova, E. A., Bock, C., Ilin, D., Wang, D. & MacDougall, B. Infrared spectroscopy on size-controlled synthesized Pt-based nano-catalysts. *Surf. Sci.* **600** (17), 3502-3511 (2006).
- 19 Garsany, Y., Singer, I. L. & Swider-Lyons, K. E. Impact of film drying procedures on RDE characterization of Pt/VC electrocatalysts. *J. of Electroanal. Chem.* **662** (2), 396-406 (2011).
- 20 Garsany, Y., Baturina, O. A., Swider-Lyons, K. E. & Kocha, S. S. Experimental Methods for Quantifying the Activity of Platinum Electrocatalysts for the Oxygen Reduction Reaction. *Anal. Chem.* **82** (15) (2010).
- 21 Kongkanand, A. & Mathias, M. F. The Priority and Challenge of High-Power Performance of Low-Platinum Proton-Exchange Membrane Fuel Cells. *J. Phys. Chem. Lett.* **7** (7), 1127-1137 (2016).
- 22 Speder, J. *et al.* Pt based PEMFC catalysts prepared from colloidal particle suspensions - a toolbox for model studies. *Phys. Chem. Chem. Phys.* **15** (10), 3602-3608 (2013).
- 23 Neumann, S. *et al.* Nanoparticles in a box: a concept to isolate, store and re-use colloidal surfactant-free precious metal nanoparticles. *J. Mat. Chem. A.* **5** (13), 6140-6145 (2017).

Figure 1








(a)



(c)



(b)



(d)







Figure 5




(a)




(b)



(a)



(b)









Figure 9







(a)




(b)



(c)





	Pt particle size	ECSA	ORR specific activity	ORR mass activity
	[nm]	[m ² g ⁻¹ _{Pt}]	[μA cm ⁻² _{Pt}]	[A g ⁻¹ _{Pt}]
50 wt.% Pt/Vulcan	2	102 ± 3	852 ± 66	879 ± 82
46 wt.% Pt/C (TKK)	2-3	93 ± 3	738 ± 30	683 ± 31

Name of Material/ Equipment	Company	Catalog Number	Comments/Description
Ethylene glycol	Sigma-Aldrich	293237-1L	99.8%
Sodium hydroxide	Fisher Chemical	AC134070010	98%
Hexachloroplatinic(IV) acid hexahydrate	Alfa Aesar	11051	99.9%
Hydrochloric acid	Merck	1003171000	37%
Nitric acid	Any		60%
Sulfuric acid	Any		96%
Perchloric acid	Merck	1005170250	70%, Suprapur
Potassium hydroxide hydrate	Merck	1050020500	99.995%, Suprapur
Tin(II) chloride dihydrate	Sigma-Aldrich	31669-100G	>98%
Platinum standard for AAS	Sigma-Aldrich	47037-100ML	1000 mg/L \pm 4 mg/L
Acetone	Alfa Aesar	30698	99.5+%
Isopropanol	Alfa Aesar	36644	99.5%
Carbon black	Cabot	Vulcan XC 72R	
46 wt.% Pt/C	Tanaka Kikinzoku Kogyo	TEC10E50E	
Ar gas	Air Liquide		99.999%
O ₂ gas	Air Liquide		99.999%
CO gas	Air Liquide		99.97%
H ₂ gas	Air Liquide		99.999%
Microwave reactor	CEM	Discover SP	
Centrifuge	Corning	6759	
Rotary evaporator	KNF	RC600	
Ultrasonic bath	Qualilab	USR 54 H	35 kHz, 160/320 W
Filter paper	Albet	DP595 055	Retention rate 4-7 μ m, ϕ 55 mm
Crucible	VWR	459-0202	12 mL
Lid for crucible	VWR	459-0216	ϕ 35 mm
pH meter	VWR	Symphony SP70P	
Glass electrode (for pH meter)	Mettler Toledo	InLab Routine	
Light scattering	Anton Paar	Litesizer	
TEM microscope	FEI	Tecani Spirit	
TEM grid	Quantifoil	N1-C73nCu20-01	Classic carbon film, Cu 200 mesh
TEM grid	Quantifoil		Classic holey carbon film, Cu 200 mesh
TEM grid	Quantifoil	N1-C74nCu20-01	Classic lacey carbon film, Cu 200 mesh
UV-vis spectrophotometer	Varian	CARY 5E	
Saturated calomel electrode	Schott	B3510+	
Perfluorinated membrane	Fuel Cell Store	591239	Nafion 117
RDE rotator	Raiometer Analytical	EDI101	
Potentiostat	Nordic Electrochemistry	ECi 200	
Control and data acquisition program	Nordic Electrochemistry	EC4 DAQ	
Data analysis program	Nordic Electrochemistry	EC4 VIEW	



1 Alewife Center #200
 Cambridge, MA 02140
 tel. 617.945.9051
www.jove.com

ARTICLE AND VIDEO LICENSE AGREEMENT

Title of Article:

On the preparation and testing of fuel cell catalysts using the thin film rotating disc electrode method

Author(s):

Masanori Inaba, Jonathan Quinson, Jan Bucher, Matthias Arenz

Item 1 (check one box): The Author elects to have the Materials be made available (as described at <http://www.jove.com/author>) via: ☐ Standard Access ☒ Open Access

Item 2 (check one box):

- ☒ The Author is NOT a United States government employee.
- ☐ The Author is a United States government employee and the Materials were prepared in the course of his or her duties as a United States government employee.
- ☐ The Author is a United States government employee but the Materials were NOT prepared in the course of his or her duties as a United States government employee.

ARTICLE AND VIDEO LICENSE AGREEMENT

1. **Defined Terms.** As used in this Article and Video License Agreement, the following terms shall have the following meanings: **"Agreement"** means this Article and Video License Agreement; **"Article"** means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein; **"Author"** means the author who is a signatory to this Agreement; **"Collective Work"** means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; **"CRC License"** means the Creative Commons Attribution-Non Commercial-No Derivs 3.0 Unported Agreement, the terms and conditions of which can be found at: <http://creativecommons.org/licenses/by-nc-nd/3.0/legalcode>; **"Derivative Work"** means a work based upon the Materials or upon the Materials and other pre-existing works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; **"Institution"** means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; **"JoVE"** means MyJoVE Corporation, a Massachusetts corporation and the publisher of *The Journal of Visualized Experiments*; **"Materials"** means the Article and / or the Video; **"Parties"** means the Author and JoVE; **"Video"** means any video(s) made by the Author, alone or in conjunction with any other parties, or by JoVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion of the Article, and in which the Author may or may not appear.

2. **Background.** The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JoVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.

3. **Grant of Rights in Article.** In consideration of JoVE agreeing to publish the Article, the Author hereby grants to JoVE, subject to **Sections 4** and **7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the "Open Access" box has been checked in **Item 1** above, JoVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.

ARTICLE AND VIDEO LICENSE AGREEMENT

4. Retention of Rights in Article. Notwithstanding the exclusive license granted to JoVE in **Section 3** above, the Author shall, with respect to the Article, retain the non-exclusive right to use all or part of the Article for the non-commercial purpose of giving lectures, presentations or teaching classes, and to post a copy of the Article on the Institution's website or the Author's personal website, in each case provided that a link to the Article on the JoVE website is provided and notice of JoVE's copyright in the Article is included. All non-copyright intellectual property rights in and to the Article, such as patent rights, shall remain with the Author.

5. Grant of Rights in Video – Standard Access. This **Section 5** applies if the "Standard Access" box has been checked in **Item 1** above or if no box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby acknowledges and agrees that, Subject to **Section 7** below, JoVE is and shall be the sole and exclusive owner of all rights of any nature, including, without limitation, all copyrights, in and to the Video. To the extent that, by law, the Author is deemed, now or at any time in the future, to have any rights of any nature in or to the Video, the Author hereby disclaims all such rights and transfers all such rights to JoVE.

6. Grant of Rights in Video – Open Access. This **Section 6** applies only if the "Open Access" box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby grants to JoVE, subject to **Section 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Video in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Video into other languages, create adaptations, summaries or extracts of the Video or other Derivative Works or Collective Works based on all or any portion of the Video and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. For any Video to which this Section 6 is applicable, JoVE and the Author hereby grant to the public all such rights in the Video as provided in, but subject to all limitations and requirements set forth in, the CRC License.

7. Government Employees. If the Author is a United States government employee and the Article was prepared in the course of his or her duties as a United States government employee, as indicated in **Item 2** above, and any of the licenses or grants granted by the Author hereunder exceed the scope of the 17 U.S.C. 403, then the rights granted hereunder shall be limited to the maximum rights permitted under such

statute. In such case, all provisions contained herein that are not in conflict with such statute shall remain in full force and effect, and all provisions contained herein that do so conflict shall be deemed to be amended so as to provide to JoVE the maximum rights permissible within such statute.

8. Likeness, Privacy, Personality. The Author hereby grants JoVE the right to use the Author's name, voice, likeness, picture, photograph, image, biography and performance in any way, commercial or otherwise, in connection with the Materials and the sale, promotion and distribution thereof. The Author hereby waives any and all rights he or she may have, relating to his or her appearance in the Video or otherwise relating to the Materials, under all applicable privacy, likeness, personality or similar laws.

9. Author Warranties. The Author represents and warrants that the Article is original, that it has not been published, that the copyright interest is owned by the Author (or, if more than one author is listed at the beginning of this Agreement, by such authors collectively) and has not been assigned, licensed, or otherwise transferred to any other party. The Author represents and warrants that the author(s) listed at the top of this Agreement are the only authors of the Materials. If more than one author is listed at the top of this Agreement and if any such author has not entered into a separate Article and Video License Agreement with JoVE relating to the Materials, the Author represents and warrants that the Author has been authorized by each of the other such authors to execute this Agreement on his or her behalf and to bind him or her with respect to the terms of this Agreement as if each of them had been a party hereto as an Author. The Author warrants that the use, reproduction, distribution, public or private performance or display, and/or modification of all or any portion of the Materials does not and will not violate, infringe and/or misappropriate the patent, trademark, intellectual property or other rights of any third party. The Author represents and warrants that it has and will continue to comply with all government, institutional and other regulations, including, without limitation all institutional, laboratory, hospital, ethical, human and animal treatment, privacy, and all other rules, regulations, laws, procedures or guidelines, applicable to the Materials, and that all research involving human and animal subjects has been approved by the Author's relevant institutional review board.

10. JoVE Discretion. If the Author requests the assistance of JoVE in producing the Video in the Author's facility, the Author shall ensure that the presence of JoVE employees, agents or independent contractors is in accordance with the relevant regulations of the Author's institution. If more than one author is listed at the beginning of this Agreement, JoVE may, in its sole discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JoVE reserves the right, in its absolute and sole discretion and without giving any reason therefore, to accept or decline any work submitted to JoVE. JoVE and its employees, agents and independent contractors shall have

ARTICLE AND VIDEO LICENSE AGREEMENT

full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JoVE has sole discretion as to the method of making and publishing the Materials, including, without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

11. **Indemnification.** The Author agrees to indemnify JoVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JoVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JoVE, making of videos by JoVE, or publication in JoVE or elsewhere by JoVE. The Author shall be responsible for, and shall hold JoVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to the making of a video by JoVE its employees, agents or independent contractors. All sterilization, cleanliness or decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's


expense. All indemnifications provided herein shall include JoVE's attorney's fees and costs related to said losses or damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JoVE, its employees, agents or independent contractors.

12. **Fees.** To cover the cost incurred for publication, JoVE must receive payment before production and publication the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JoVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.

13. **Transfer, Governing Law.** This Agreement may be assigned by JoVE and shall inure to the benefits of any of JoVE's successors and assignees. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to be one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement required per submission.

CORRESPONDING AUTHOR:

Name:	<u>Matthias Arenz</u>	
Department:	<u>Chemistry and Biochemistry</u>	
Institution:	<u>University of Bern</u>	
Article Title:	<u>On the preparation and testing of fuel cell catalysts using the thin film rotating disc electrode method</u>	
Signature:		Date: <u>2nd August 2017</u>

Please submit a signed and dated copy of this license by one of the following three methods:

- 1) Upload a scanned copy of the document as a pdf on the JoVE submission site;
- 2) Fax the document to +1.866.381.2236;
- 3) Mail the document to JoVE / Attn: JoVE Editorial / 1 Alewife Center #200 / Cambridge, MA 02139

For questions, please email submissions@jove.com or call +1.617.945.9051

Prof. Dr. Matthias Arenz

Department of Chemistry and Biochemistry
University of Bern
Freiestrasse 3
CH - 3012 Bern
Tel +41 31 631 53 84
Email: matthias.arenz@dcb.unibe.ch

u^b

u^b
**UNIVERSITÄT
BERN**

Bern, 21st November 2017

Dear Editor,

Please find attached the revised version of the manuscript entitled

“On the preparation and testing of fuel cell catalysts using the thin film rotating disc electrode method”

by Inaba *et al.* for publication in Jove.

We addressed all points from the Editor and revised the manuscript accordingly. We submit one manuscript with all changes visible and one with all changes accepted.

Sincerely yours,

Matthias Arenz