

# Journal of Visualized Experiments

## Quantitative Analysis of Vacuum Induction Melting by Laser-induced Breakdown Spectroscopy

--Manuscript Draft--

<b>Article Type:</b>	Invited Methods Article - JoVE Produced Video
<b>Manuscript Number:</b>	JoVE57903R3
<b>Full Title:</b>	Quantitative Analysis of Vacuum Induction Melting by Laser-induced Breakdown Spectroscopy
<b>Keywords:</b>	Laser-induced breakdown spectroscopy (LIBS), Plasma spectroscopy, Elemental composition, On-line process monitoring, Vacuum induction melting (VIM), Liquid steel, Molten alloy, Relative standard deviation (RSD), Limit of detection (LOD)
<b>Corresponding Author:</b>	Tianzhuo Zhao Academy of Opto-Electronics Beijing, Beijing CHINA
<b>Corresponding Author's Institution:</b>	Academy of Opto-Electronics
<b>Corresponding Author E-Mail:</b>	zhaotianzhuo@aoe.ac.cn
<b>Order of Authors:</b>	Tianzhuo Zhao Zhongwei Fan Fuqiang Lian Hong Xiao Shuzhen Nie Yang Liu Xin Li Qixiu Zhong Sining Sun
<b>Additional Information:</b>	
<b>Question</b>	<b>Response</b>
Please indicate whether this article will be Standard Access or Open Access.	Standard Access (US\$2,400)
Please indicate the <b>city, state/province, and country</b> where this article will be <b>filmed</b> . Please do not use abbreviations.	Zhongkeheguang Applied Laser Technology Institute Company, Ltd., Low-Carbon Industry Base G2, Huaming High-tech Zone, Hongcheng Road, Tianjin, China

**TITLE:**

Quantitative Analysis of Vacuum Induction Melting by Laser-induced Breakdown Spectroscopy

**AUTHORS & AFFILIATIONS:**

Tianzhuo Zhao<sup>1,2,3</sup>, Zhongwei Fan<sup>1,2,3</sup>, Fuqiang Lian<sup>1,3</sup>, Hong Xiao<sup>1,3</sup>, Shuzhen Nie<sup>1,3</sup>, Yang Liu<sup>1,3</sup>,  
Xin Li<sup>1,2</sup>, Qixiu Zhong<sup>1,2</sup>, Sining Sun<sup>4</sup>

<sup>1</sup>*The Academy of Opto-electronics, Chinese Academy of Sciences, Beijing, China*

<sup>2</sup>*University of Chinese Academy of Sciences, Beijing, China*

<sup>3</sup>*National Engineering Research Center for DPSSL, Beijing, China*

<sup>4</sup>*Beijing GK Laser Technology Co., Ltd., Beijing, China*

*Corresponding Author:*

Tianzhuo Zhao

[zhaotianzhuo@aoe.ac.cn](mailto:zhaotianzhuo@aoe.ac.cn)

Zhongwei Fan

[fanzhongwei@aoe.ac.cn](mailto:fanzhongwei@aoe.ac.cn)

*Email Addresses of Co-authors:*

Fuqiang Lian ([lianfuqiang@aoe.ac.cn](mailto:lianfuqiang@aoe.ac.cn))

Hong Xiao ([xiaohongava@163.com](mailto:xiaohongava@163.com))

Shuzhen Nie ([nieshuzhen@aoe.ac.cn](mailto:nieshuzhen@aoe.ac.cn))

Yang Liu ([lyang@aoe.ac.cn](mailto:lyang@aoe.ac.cn))

Xin Li ([lixin815@mails.ucas.edu.cn](mailto:lixin815@mails.ucas.edu.cn))

Qixiu Zhong ([zhongqiuxiu16@mails.ucas.edu.cn](mailto:zhongqiuxiu16@mails.ucas.edu.cn))

Sining Sun ([sunsining@gklaser.com](mailto:sunsining@gklaser.com))

**KEYWORDS:**

Laser-induced breakdown spectroscopy (LIBS), Plasma spectroscopy, Elemental composition, On-line process monitoring, Vacuum induction melting (VIM), Liquid steel, Molten alloy, Relative standard deviation (RSD), Limit of detection (LOD)

**SUMMARY:**

During vacuum induction melting, laser-induced breakdown spectroscopy is used to perform real-time quantitative analysis of the main-ingredient elements of a molten alloy.

**ABSTRACT:**

Vacuum induction melting is a popular method for refining high purity metal and alloys. Traditionally, standard process control in metallurgy involves several steps, include drawing samples, cooling, cutting, transport to the laboratory, and analysis. The whole analysis process requires more than 30 minutes, which hinders on-line process control. Laser-induced breakdown spectroscopy is an excellent on-line analysis method that can satisfy the requirements of vacuum induction melting because it is fast and noncontact and does not

require sample preparation. The experimental facility uses a lamp-pumped Q-switched laser to ablate melted liquid steel with an output energy of 80 mJ, a frequency of 5 Hz, a FWHM pulse width of 20 ns, and a working wavelength of 1064 nm. A multi-channel linear charge coupled device (CCD) spectrometer is used to measure the emission spectrum in real time, with a spectral range from 190 to 600 nm and a resolution of 0.06 nm at a wavelength of 200 nm. The protocol includes several steps: standard alloy sample preparation and an ingredient test, smelting of standard samples and determination of the laser breakdown spectrum, and construction of the elements concentration quantitative analysis curve of each element. To realize the concentration analysis of unknown samples, the spectrum of a sample also needs to be measured and disposed with the same process. The composition of all main elements in the melted alloy can be quantitatively analyzed with an internal standard method. The calibration curve shows that the limit of detection of most metal elements ranges from 20-250 ppm. The concentration of elements, such as Ti, Mo, Nb, V, and Cu, can be lower than 100 ppm, and the concentrations of Cr, Al, Co, Fe, Mn, C, and Si range from 100-200 ppm. The  $R^2$  of some calibration curves can exceed 0.94.

## INTRODUCTION:

Due to its unique features, such as remote sensing, fast analysis, and no need for sample preparation, laser-induced breakdown spectroscopy (LIBS) offers unique capabilities for on-line concentration determination<sup>1-3</sup>. Although the use of the LIBS technique in different fields has been investigated<sup>4-6</sup>, a considerable attempt to develop its capabilities in industrial applications is ongoing.

Analysis of molten material contents during the course of industrial processes can effectively improve the product quality, which is a promising development direction of LIBS. Experimental findings have been reported about the application of LIBS in the industrial field, such as findings about argon oxygen liquid steel<sup>7-11</sup>, molten aluminium alloy<sup>12</sup>, molten salt<sup>13</sup>, and molten silicon<sup>14</sup>. The majority of these materials exist in the environment of air or an assistant gas. However, vacuum induction melting (VIM) is another good application field of LIBS to realize processing control. A VIM furnace can realize smelting at temperatures higher than 1700 °C for alloy refining; it is the most popular method for refining high-purity metal and alloys such as iron-base or nickel-base alloys, high purity alloys, and clean magnetic alloys. During the course of melting, the pressure in a furnace is always in the region of 1-10 Pa, and the composition of air in the furnace mainly includes the air absorbed on the sample or the inner wall of the furnace and some vaporous oxide or nitride metal. These working situations induce quite different LIBS measurement situations for smelting in air. Here, we report an experimental investigation of the analysis of molten alloy during the course of VIM by LIBS.

An optical window is added to a furnace for laser ablation and radiant light detection. A silica glass with a diameter of 80 mm serves as the window. An emitting laser and gathering of radiant light employ the same window; it is a co-axial optical structure that focuses on the same point. The working focal length is approximately 1.8 m, and the focusing length of the experimental setup can be adjusted from 1.5 to 2.5 m.

Based on the practicality of industrial online analysis, precision, repeatability and stability is more important than the low limit of detection (LOD) during molten alloy ingredient analysis. The technical route of a four-channel linear CCD spectrometer is chosen, the spectral range of the spectrometer ranges from 190 to 600 nm, the resolution is 0.06 nm, and the wavelength is 200 nm. A laser diode pumped Q-switched laser (constructed in house) is used to ablate molten alloy, with an output energy of 100 mJ, a frequency of 5 Hz, an FWHM pulse width of 20 ns, and a working wavelength of 1064 nm. The remaining part will present the VIM LIBS-analysing process and live measurement, followed by an introduction of the data processing results.

## **PROTOCOL:**

### **1. Preparation of Standard Samples**

NOTE: This step is not essential.

1.1. Prepare raw material (**Table 1**). To make a 100 kg of sample #1, add 12.82 kg of Cr, 3.39 kg of Mo, 4.79 kg of Al, 1.00 kg of Ti, 0.60 kg of Cu, and approximately 77.4 kg of Ni to the crucible. During the melting process, some elements will be burned. The final ingredient is determined by the melting temperature, melting duration, and other working parameters. The ingredient test shows the quantity of each element inside the alloys.

1.2. Perform vacuum induction melting at approximately 1700 °C for approximately 45 min for each standard sample. The furnace used to make standard samples can melt approximately 100 kg alloys each time for 2 sets of standard samples.

1.3. Pour all molten liquid steel into a sticky mold to make standard samples, and naturally cool for at least 4 h. The size of the standard samples is determined by the furnace in the experiment. Use rod-shaped standard samples in the experiments with a rod diameter of 100 mm. The shape of the crucible in the furnace is a frustum cone with a cup-like container. The diameter of the rim is 150 mm, the bottom is 100 mm, and the depth is 200 mm. Chunked standard samples can also be employed.

1.4. Use an electric saw to cut standard samples for the experiment. The length of the standard rod sample is determined by the furnace. Use a length of 150 mm for the experimental smelting system. Melt a piece of the samples for each experiment.

1.5. Repeat these steps and make all standard samples. Ten samples are employed in this experiment.

### **2. Test Ingredient of Standard Alloy Samples**

2.1. Use a chemical analysis method to test the composition of all standard alloy samples. Test all elements in each sample.

NOTE: We strongly recommend sending these samples to an authority organization to perform the analysis. These samples are sent to the Central Iron & Steel Research Institute of China for the ingredients test. The test results for these samples are listed in **Table 2**.

### **3. Smelt Samples**

3.1. Check the security of the smelting system, which includes the power supply, availability of each pump, vacuum hold ability of the experimental furnace, cooling water, and current.

3.2. Place the standard samples into the smelting system. To ensure that a small amount of each sample is ablated by the laser, use small samples for the experiment. Due to the size of the crucible in the furnace, smelt approximately 10 kg samples each time.

3.3. Open the vacuum pump until the pressure is lower than 0.1 Pa. Use 2 levels of pump to make the vacuum: mechanical pump and diffusion pump. The mechanical pump can reach approximately 1 Pa in 15 min, and the diffusion pump can reach 0.01 Pa after 40 min.

3.4. Melt samples. Increase the furnace working current to approximately 130 A; this parameter is determined by the ingredients of the samples and the size of the furnace. A standard sample requires approximately 15 min to become molten. Due to oxidation or nitridation, the ingredients of liquid steel slowly change during the course of smelting.

3.4.1. To ensure precision of the experiment, determine the spectrum within 15 min after the standard samples are molten.

### **4. Determine Laser Breakdown Spectra of Standard Samples**

4.1. Check availability of laser focusing and spectrum gathering system, laser generator, and spectrometer.

4.2. Set the spectrometer and laser generator to work synchronously. Use a spectrometer output synchronization signal and laser passive working method in the system. The method of laser generator synchronization signal or method of synchronizer output synchronization signal can also be employed to control the laser generator and spectrometer.

4.3. Open the laser generator and spectrometer; prepare to generate the pulse laser. The pulse width is 20 ns, the frequency is 5 Hz, and the energy of each pulse is 90 mJ.

4.4. Use spectrum deposit software to trigger the laser output and gather the spectrum. Set the integration time of the spectrometer to 10 ms, and each laser pulse generates a frame of the spectrum. If the integration time is too short, the spectrum signal intensity will be too weak. If the integration time is too long, more background signals will be gathered.

4.5. Adjust the laser focusing position, and effectively ablate the sample. Optimize the focusing position until the strongest spectrum signal is obtained. This process is used to adjust the focusing point. The spectrometer numerical divided signal intensity ranges from 0 to 65535. In most cases, the intensity of a signal should exceed 15% of the saturation signal, which indicates that the highest peak intensity should exceed 10000. If the signal intensity is too small, the quantitative analysis will have low precision.

4.6. Optimize the delay time. Choose the delay after bremsstrahlung, and the strength of the signal need with the optimized delay time should be sufficient.

4.7. Use the spectrometer to gather a spectrum for the analysis. Gather 20 frames of the spectrum, and obtain the average for LIBS analysis.

4.8. Shut off the working current of the furnace and cool the samples. Solidification of the samples requires approximately 15 min.

4.9. Inject nitrogen into the experimental furnace to break the vacuum.

4.10. Open the lid of the experimental furnace and remove the solidification samples.

4.11. Repeat step 3.3 to 4.10 until all samples are measured.

## **5. Construct Calibration Curve of Quantitative Analysis**

5.1. Spectrum pretreatment

5.1.1. Background correction. Delete the background effect caused by braking radiation. The method of baseline correction is employed in the experiment.

5.1.2. Spectrum-peak searching. Use a two-order derivative method to identify peaks of each element; local minimum points are weighted.

5.1.3. Spectrum fitting. Apply a Lorentz spectra overlay to the selected peaks to prevent self-corrosion or overlap. The spectral peak intensity, stretch status, and center wavelength is obtained by a fitting algorithm.

5.2. Import the chemical ingredient analysis results of all standard samples.

5.3. Construct the calibration curve.

5.3.1. Choose an inner relative standard wavelength. The main element spectral lines are always selected.

5.3.2. Choose a calibration wavelength. Select from the NIST spectrum database<sup>15</sup>.

5.3.3. Fit the curve. Use linear fitting or quadratic fitting.

5.4. Achieve analysis precision. Calculate the fitting factor and relative standard error after fitting. A program is used to automatically select the best relative standard wavelength and calibration wavelength from the wavelength base of NIST<sup>15</sup>.

## **6. Elemental Composition Analysis of Molten Alloy**

NOTE: The experimental setup has been divided to two parts, namely, the detector head and the control cabinet, as shown in **Figure 1**. The same laser and spectrometer parameters, moulting, and spectrum gathering process employed in the previous process are utilized to ensure accurate quantitative analysis results.

6.1. Put the unknown sample into the smelting system.

6.2. Vacuum the experimental system.

6.3. Increase the smelting current until the sample is molten. The melting temperature is approximately 1700 °C, and the melting time is approximately 45 min.

6.4. Open the laser generator and realize the pulse laser output. Use the following laser parameters: pulse width is 20 ns, frequency is 5 Hz, and energy of each pulse is 90 mJ.

6.5. Open the spectrometer and the spectrum deposit software to determine the spectrum. Employ the same spectrometer with a spectral range from 190 to 600 nm and resolution of 0.06 nm at a wavelength of 200 nm. The integration time of the spectrometer is 10 ms. The spectrometer is used to trigger the laser and determine the spectrum.

6.6. Adjust laser focusing position. Optimize focusing position until the strongest spectrum signal is attained; the value of the highest peak should exceed 10000.

6.7. Determine laser breakdown spectrum. Each laser pulse generates a frame of the spectrum; 20 frames of the spectrum are obtained and averaged for the analysis.

6.8. Spectrum pretreatment. Perform background correction, such as deleting the background effect caused by braking radiation, as mentioned in 5.1.3, to perform spectrum fitting.

6.9. Elemental concentration calculation. Perform analysis elemental concentration by the internal standard method from the calibration curve.

## **REPRESENTATIVE RESULTS:**

Ten nickel-based alloy samples (#1-#10) are used to construct internal-standard calibration curves. The compositions of all samples are listed in **Table 1**. The elemental concentrations of these samples are orthogonally designed to avoid signal interference. The concentration of each element in all samples is measured with chemical analysis methods.

Nickel is the internal standard element. The calibration curves of Cu, Ti, Mo, Al, and Cr are constructed. **Figure 2** to **Figure 6** show the calibration results. In these figures, the X-axis represents the concentration of the calibrated elements, and the Y-axis represents the relative signal intensity ratio of the calibrated element after the disposal process of background correction and peak fitting. The error bar of each point in these figures shows the fluctuation range of the signal strength with twenty frame measurements. The calibration parameters of these elements are listed in **Table 3** to **Table 7**. The linear curve fitting results, including the residual sum of squares, Pearson's  $r$ , and the linear fitting coefficient  $R^2$ , are shown from **Figure 2** to **Figure 6**. The intercept and slope of the coefficient of determination are also shown in these figures. The calibration curves show a near-linear relationship between the concentration of the element and the peak intensity. The spectral lines used for each element were introduced in the legend of these figures. These lines are searched by a method of filtration. All signal peaks are filtered by the signal intensity, the central of wavelength, and the Lorenz fitting effect. These selected peaks are chosen by a permutation-combination analysis of the fitting factor  $R^2$ .

According to the standard of International Union of Pure and Applied Chemistry (IUPAC),  $3\sigma$  limit of detection (LOD) of Cu, Ti, Mo, Al and Cr are calculated and listed in **Table 8**. Other elements, such as Si, C, and Nb, are analyzed. The RSD ranges from 4-6%, and the  $R^2$  exceeds 0.93. The precision can be improved if a better relative standard is employed.

#### FIGURE AND TABLE LEGENDS:

**Figure 1:** Experimental setup of quantitative analysis in the process of vacuum induction melting by laser-induced breakdown spectroscopy.

**Figure 2:** Calibration curves of Cu. Internal standard lines include Cu: 224.70 nm, Ni: 241.61 nm and 233.75 nm.

**Figure 3:** Calibration curves of Ti. Internal standard lines include Ti: 444.38 nm and 337.22 nm, Ni: 445.90 nm and 313.41 nm.

**Figure 4:** Calibration curves of Mo. Internal standard lines include Mo: 342.23 nm, 346.02 nm, and 277.44 nm, Ni: 440.16 nm and 336.68 nm.

**Figure 5:** Calibration curves of Al. Internal standard lines include Al: 272.31 nm, 231.22 nm, and 334.85 nm, Ni: 221.65 nm, 332.23 nm, and 440.16 nm.

**Figure 6:** Calibration curves of Cr. Internal standard lines include Cr: 286.51 nm, 302.67 nm and 342.12 nm, Ni: 224.27 nm, 233.75 nm, and 350.08 nm.



**Table 1:** Raw material ingredients in the experiment.

**Table 2:** Standard nickel-based alloy samples ingredient measured results.

**Table 3:** Calibration data of Cu.

**Table 4:** Calibration data of Ti.

**Table 5:** Calibration data of Mo.

**Table 6:** Calibration data of Al.

**Table 7:** Calibration data of Cr.

**Table 8:** Limit of detection of Cu, Ti, Mo, Al, and Cr.

## **DISCUSSION:**

For elemental analysis, popular methods are X-ray fluorescence (XRF), spark discharge optical emission spectrometry (SD-OES), atomic absorption spectroscopy (AAS), and inductive couple plasma (ICP). These methods are mainly suited for a laboratory and industrial online application for molten alloys, which is determined by the characters of these technologies, is difficult. XRF uses X-rays to shock samples, and SD-OES makes sparks on the samples. The working distance of these two methods are always in the range of several centimetres. AAS and ICP yields liquid or powder samples, which requires several tens of minutes for preparation. These methods are not suitable for high-temperature samples or measurements from a distance of several metres. Compared with these analysis methods, LIBS has the advantages of long-distance analysis, fast analysis, and the need to prepare samples. LIBS is the only good method for realizing melting alloys ingredient online analysis.

The protocol includes three critical steps: using a laser to burn the molten alloy, using a spectrometer to determine the spectrum of the plasma, and quantitatively analyzing the elemental composition with the calibration curve. Preparation of the samples with gradient components and construction of the calibration curve to demonstrate the relation between the laser breakdown spectrum intensity and the elemental content are preparative steps.

Use of the LIBS to analyze the elemental composition of molten alloy has some limitations. The precision of the quantitative analysis is the most important problem. The precision of LIBS is expected to improve by an order of magnitude. The gas pressure, surface status of the samples, and focusing precision have a distinct influence on the precision; however, compensation of these errors is difficult<sup>1,2,6</sup>.

Use of the LIBS system for on-line analysis of elemental composition during vacuum melting is proven by experiments. The experimental results have shown that the plasma spectrum can be

determined in a typical industrial vacuum melting furnace situation. The calibration results show that the major components of molten alloys can be quantitatively analysed.

#### ACKNOWLEDGMENTS:

This study was financially supported by the National Key Scientific Instrument and Equipment Development Projects (Grant No. 2014YQ120351), the Youth Innovation Promotion Association of CAS (Grant No. 2014136), and the China Innovative Talent Promotion Plans for Innovation Team in Priority Fields (Grant No. 2014RA4051).


#### DISCLOSURES:


The authors have nothing to disclose.


#### REFERENCES:


1. Radziemski L., Cremers D., A brief history of laser-induced breakdown spectroscopy: From the concept of atoms to LIBS 2012, *Spectrochimica Acta Part B: Atomic Spectroscopy*, **87**, 3-10, doi:10.1016/j.sab.2013.05.013 (2013).
2. El Haddad J., Canioni L., Bousquet B., Good practices in LIBS analysis: Review and advices, *Spectrochimica Acta Part B: Atomic Spectroscopy*, **101**, 171-182, doi:10.1016/j.sab.2014.08.039 (2014).
3. Mueller M., Gornushkin I.B., Florek S., Mory D., Panne U., Approach to Detection in Laser-Induced Breakdown Spectroscopy, *Analytical Chemistry*, **79** (12), 4419-4426, doi:10.1021/ac0621470 (2007).
4. Noll R., Fricke-Begemann C., Brunk M., Connemann S., Meinhardt C., Schsrn M., Sturm V., Makowe J., Gehlen C., Laser-induced breakdown spectroscopy expands into industrial applications, *Spectrochimica Acta Part B: Atomic Spectroscopy*, **93**, 41-51, doi:10.1016/j.sab.2014.02.001 (2014).
5. Leon R., David C., A brief history of laser-induced breakdown spectroscopy: From the concept of atoms to LIBS 2012, *Spectrochimica Acta Part B: Atomic Spectroscopy*, **87**, 3-10 (2013).
6. El Haddad J., Canioni L., Bousquet B., Good practices in LIBS analysis: Review and advices, *Spectrochimica Acta Part B: Atomic Spectroscopy*, **101**, 171-182, doi:10.1016/j.sab.2014.08.039 (2014).
7. Gonzaga B.F., Pasquini C., A compact and low cost laser induced breakdown spectroscopic system: Application for simultaneous determination of chromium and nickel in steel using multivariate calibration, *Spectrochimica Acta Part B: Atomic Spectroscopy*, **69**, 20-24, doi:10.1016/j.sab.2012.02.007 (2012).
8. Peter L., Sturm V., Noll R., Liquid steel analysis with laser-induced breakdown spectrometry in the vacuum ultraviolet *Applied Optics*, **42** (30), 6199-6204, doi:10.1364/AO.42.006199 (2003).
9. Hubmer G., Kitzberger R., Mörwald K., Application of LIBS to the in-line process control of liquid high-alloy steel under pressure *Analytical and Bioanalytical Chemistry*, **385** (2), 219-224, doi:10.1007/s00216-006-0321-9 (2006).


10. Sun L.X., Yu H.B., Automatic estimation of varying continuum background emission in laser-induced breakdown spectroscopy *Spectrochimica Acta Part B: Atomic Spectroscopy*, **64**, 278-287, doi:10.1016/j.sab.2009.02.010 (2009).
11. Lin X.M., Chang P.H., Chen G.H., Lin J.J., Liu R.X., Yang H., Effect of melting iron-based alloy temperature on carbon content observed in laser-induced breakdown spectroscopy *Plasma Science & Technology*, **17** (11), 933-937, doi:10.1088/1009-0630/17/11/08 (2015).
12. Rai A.K., Yueh F.Y., Singh J.P., Laser-induced breakdown spectroscopy of molten aluminum alloy *Applied Optics*, **42** (12), 2078-2084, doi:10.1364/AO.42.002078 (2003).
13. Hanson C., Phongikaroon S., Scott J.R., Temperature effect on laser-induced breakdown spectroscopy spectra of molten and solid salts *Spectrochimica Acta Part B: Atomic Spectroscopy*, **97**, 79-85, doi:10.1016/j.sab.2014.04.012 (2014).
14. Darwiche S., Benrabbah R., Benmansour M., Morvan D., Impurity detection in solid and molten silicon by laser induced breakdown spectroscopy *Spectrochimica Acta Part B: Atomic Spectroscopy*, **74**, 115-118, doi:10.1016/j.sab.2012.06.020 (2012).
15. P.J. Linstrom and W.G. Mallard, Eds., NIST Chemistry WebBook, NIST Standard Reference Database Number 69, National Institute of Standards and Technology, Gaithersburg MD, 20899, doi:10.18434/T4D303 (2018).
















						unit: %
Sample \ Element	Cu	Ti	Mo	Al	Cr	Ni
1#	0.60	1.00	3.39	4.79	12.82	rest
2#	0.80	0.39	5.79	7.92	10.08	rest
3#	0.89	0.26	5.16	5.44	13.49	rest
4#	0.71	1.41	5.43	5.96	11.75	rest
5#	1.00	1.23	3.02	5.05	14.69	rest
6#	0.35	1.56	4.17	7.05	12.33	rest
7#	0.53	1.80	4.87	6.61	14.08	rest
8#	0.45	2.00	4.51	4.27	11.20	rest
9#	0.28	0.64	3.62	6.35	15.40	rest
10#	0.17	0.82	3.98	7.56	10.59	rest

						unit:%
Element Sample	Cu	Ti	Mo	Al	Cr	Ni
1#	0.60	1.11	3.42	5.33	12.62	rest
2#	0.78	0.41	5.85	8.64	9.77	rest
3#	0.87	0.28	5.15	6.12	13.11	rest
4#	0.69	1.56	5.56	6.34	11.50	rest
5#	0.93	1.35	3.01	5.63	14.57	rest
6#	0.36	1.75	4.23	7.64	11.97	rest
7#	0.54	2.04	4.94	7.21	13.50	rest
8#	0.46	2.18	4.59	4.61	11.04	rest
9#	0.30	0.65	3.70	6.74	15.11	rest
10#	0.19	0.92	4.04	8.12	10.18	rest

Calibration data of Cu			
Concentration (%)	Signal Intensity		
	Average	Maxium	Minium
0.6	0.1742	0.0136	-0.0105
0.77	0.1901	0.0168	-0.0227
0.87	0.1994	0.0222	-0.0195
0.68	0.1898	0.011	-0.01
0.92	0.2149	0.011	-0.0099
0.36	0.1551	0.0163	-0.0146
0.53	0.1838	0.0133	-0.0154
0.46	0.1623	0.0082	-0.0136
0.29	0.1552	0.0146	-0.0174
0.18	0.1343	0.0145	-0.0156

Calibration data of Ti			
Concentration (%)	Signal Intensity		
	Average	Maxium	Minium
1.11	1.1755	0.0539	-0.0891
0.412	0.8338	0.0205	-0.0232
0.284	0.7254	0.0309	-0.0401
1.56	1.4204	0.0761	-0.0466
1.35	1.3082	0.0296	-0.0587
1.75	1.5039	0.0363	-0.066
2.04	1.6569	0.093	-0.0653
2.18	1.7069	0.1017	-0.084
0.652	0.952	0.0667	-0.0605
0.921	1.0868	0.032	-0.032

Calibration data of Mo			
Concentration (%)	Signal Intensity		
	Average	Maxium	Minium
3.42	0.4184	0.0196	-0.0188
5.85	0.5508	0.0114	-0.0149
5.15	0.5156	0.0244	-0.0356
5.56	0.54	0.0537	-0.0582
3.01	0.4174	0.0164	-0.0202
4.23	0.4714	0.0404	-0.0544
4.94	0.5074	0.0206	-0.0272
4.59	0.4636	0.0117	-0.0093
3.7	0.4518	0.04	-0.0432
4.04	0.4847	0.0157	-0.0141

Calibration data of Al			
Concentration (%)	Signal Intensity		
	Average	Maxium	Minium
5.33	2.8048	0.1672	-0.1334
8.64	3.9482	0.2182	-0.2832
6.12	3.1027	0.1528	-0.1383
6.34	3.1552	0.1401	-0.1593
5.63	2.8836	0.1025	-0.0698
7.64	3.6064	0.1523	-0.1258
7.21	3.5292	0.2307	-0.1389
4.61	2.4939	0.0421	-0.05
6.74	3.2185	0.0868	-0.087
8.12	3.777	0.1284	-0.1071

Calibration data of Cr			
Concentration (%)	Signal Intensity		
	Average	Maxium	Minium
12.62	2.5908	0.0923	-0.0764
9.77	2.2203	0.095	-0.0859
13.11	2.7817	0.1338	-0.0874
11.5	2.4624	0.1285	-0.1273
14.57	3.0377	0.0465	-0.0527
11.97	2.5475	0.1476	-0.1173
13.5	2.8318	0.1108	-0.1147
11.04	2.3723	0.0352	-0.0563
15.11	3.094	0.144	-0.1524
10.18	2.2301	0.0747	-0.0604



	Cu	Ti	Mo	Al	Cr
LOD(ppm)	37	19	98	258	166

Name of Material/ Equipment	Company	Catalog Number	Comments/Description
Laser source	Gklaser Co.,Ltd.		
Molten alloy to be measured			
Smelting furnace	Tianyu Co.,Ltd.		
Spectrometer	Avantes		
standard samples			Well known of its composition



1 Alewife Center #200  
Cambridge, MA 02140  
tel. 617.945.9051  
www.jove.com

## ARTICLE AND VIDEO LICENSE AGREEMENT

Title of Article:

Realize Elements Concentration Analysis for Vacuum Induction Melting by Laser-Induced Breakdown Spectroscopy

Author(s):

Tianzhuo Zhao, Zhongwei Fan, Fuqiang Lian, Yang Liu, Hong Xiao, Xin Li, Qixiu Zhong, Pu Wang, Shuzhen Nie

Item 1 (check one box): The Author elects to have the Materials be made available (as described at <http://www.jove.com/author>) via: ☒ Standard Access ☐ Open Access

Item 2 (check one box):

- ☒ The Author is NOT a United States government employee.
- ☐ The Author is a United States government employee and the Materials were prepared in the course of his or her duties as a United States government employee.
- ☐ The Author is a United States government employee but the Materials were NOT prepared in the course of his or her duties as a United States government employee.

### ARTICLE AND VIDEO LICENSE AGREEMENT

1. **Defined Terms.** As used in this Article and Video License Agreement, the following terms shall have the following meanings: "**Agreement**" means this Article and Video License Agreement; "**Article**" means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein; "**Author**" means the author who is a signatory to this Agreement; "**Collective Work**" means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; "**CRC License**" means the Creative Commons Attribution-Non Commercial-No Derivs 3.0 Unported Agreement, the terms and conditions of which can be found at: <http://creativecommons.org/licenses/by-nc-nd/3.0/legalcode>; "**Derivative Work**" means a work based upon the Materials or upon the Materials and other pre-existing works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; "**Institution**" means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; "**JoVE**" means MyJoVE Corporation, a Massachusetts corporation and the publisher of *The Journal of Visualized Experiments*; "**Materials**" means the Article and / or the Video; "**Parties**" means the Author and JoVE; "**Video**" means any video(s) made by the Author, alone or in conjunction with any other parties, or by JoVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion of the Article, and in which the Author may or may not appear.

2. **Background.** The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JoVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.

3. **Grant of Rights in Article.** In consideration of JoVE agreeing to publish the Article, the Author hereby grants to JoVE, subject to Sections 4 and 7 below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the "Open Access" box has been checked in Item 1 above, JoVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.



## ARTICLE AND VIDEO LICENSE AGREEMENT

4. **Retention of Rights in Article.** Notwithstanding the exclusive license granted to JoVE in **Section 3** above, the Author shall, with respect to the Article, retain the non-exclusive right to use all or part of the Article for the non-commercial purpose of giving lectures, presentations or teaching classes, and to post a copy of the Article on the Institution's website or the Author's personal website, in each case provided that a link to the Article on the JoVE website is provided and notice of JoVE's copyright in the Article is included. All non-copyright intellectual property rights in and to the Article, such as patent rights, shall remain with the Author.
5. **Grant of Rights in Video – Standard Access.** This **Section 5** applies if the "Standard Access" box has been checked in **Item 1** above or if no box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby acknowledges and agrees that, Subject to **Section 7** below, JoVE is and shall be the sole and exclusive owner of all rights of any nature, including, without limitation, all copyrights, in and to the Video. To the extent that, by law, the Author is deemed, now or at any time in the future, to have any rights of any nature in or to the Video, the Author hereby disclaims all such rights and transfers all such rights to JoVE.
6. **Grant of Rights in Video – Open Access.** This **Section 6** applies only if the "Open Access" box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby grants to JoVE, subject to **Section 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Video in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Video into other languages, create adaptations, summaries or extracts of the Video or other Derivative Works or Collective Works based on all or any portion of the Video and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. For any Video to which this Section 6 is applicable, JoVE and the Author hereby grant to the public all such rights in the Video as provided in, but subject to all limitations and requirements set forth in, the CRC License.
7. **Government Employees.** If the Author is a United States government employee and the Article was prepared in the course of his or her duties as a United States government employee, as indicated in **Item 2** above, and any of the licenses or grants granted by the Author hereunder exceed the scope of the 17 U.S.C. 403, then the rights granted hereunder shall be limited to the maximum rights permitted under such statute. In such case, all provisions contained herein that are not in conflict with such statute shall remain in full force and effect, and all provisions contained herein that do so conflict shall be deemed to be amended so as to provide to JoVE the maximum rights permissible within such statute.
8. **Likeness, Privacy, Personality.** The Author hereby grants JoVE the right to use the Author's name, voice, likeness, picture, photograph, image, biography and performance in any way, commercial or otherwise, in connection with the Materials and the sale, promotion and distribution thereof. The Author hereby waives any and all rights he or she may have, relating to his or her appearance in the Video or otherwise relating to the Materials, under all applicable privacy, likeness, personality or similar laws.
9. **Author Warranties.** The Author represents and warrants that the Article is original, that it has not been published, that the copyright interest is owned by the Author (or, if more than one author is listed at the beginning of this Agreement, by such authors collectively) and has not been assigned, licensed, or otherwise transferred to any other party. The Author represents and warrants that the author(s) listed at the top of this Agreement are the only authors of the Materials. If more than one author is listed at the top of this Agreement and if any such author has not entered into a separate Article and Video License Agreement with JoVE relating to the Materials, the Author represents and warrants that the Author has been authorized by each of the other such authors to execute this Agreement on his or her behalf and to bind him or her with respect to the terms of this Agreement as if each of them had been a party hereto as an Author. The Author warrants that the use, reproduction, distribution, public or private performance or display, and/or modification of all or any portion of the Materials does not and will not violate, infringe and/or misappropriate the patent, trademark, intellectual property or other rights of any third party. The Author represents and warrants that it has and will continue to comply with all government, institutional and other regulations, including, without limitation all institutional, laboratory, hospital, ethical, human and animal treatment, privacy, and all other rules, regulations, laws, procedures or guidelines, applicable to the Materials, and that all research involving human and animal subjects has been approved by the Author's relevant institutional review board.
10. **JoVE Discretion.** If the Author requests the assistance of JoVE in producing the Video in the Author's facility, the Author shall ensure that the presence of JoVE employees, agents or independent contractors is in accordance with the relevant regulations of the Author's institution. If more than one author is listed at the beginning of this Agreement, JoVE may, in its sole discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JoVE reserves the right, in its absolute and sole discretion and without giving any reason therefore, to accept or decline any work submitted to JoVE. JoVE and its employees, agents and independent contractors shall have



## ARTICLE AND VIDEO LICENSE AGREEMENT

full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JoVE has sole discretion as to the method of making and publishing the Materials, including, without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

11. **Indemnification.** The Author agrees to indemnify JoVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JoVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JoVE, making of videos by JoVE, or publication in JoVE or elsewhere by JoVE. The Author shall be responsible for, and shall hold JoVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to the making of a video by JoVE its employees, agents or independent contractors. All sterilization, cleanliness or decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's


expense. All indemnifications provided herein shall include JoVE's attorney's fees and costs related to said losses or damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JoVE, its employees, agents or independent contractors.

12. **Fees.** To cover the cost incurred for publication, JoVE must receive payment before production and publication the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JoVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.

13. **Transfer, Governing Law.** This Agreement may be assigned by JoVE and shall inure to the benefits of any of JoVE's successors and assignees. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to be one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement required per submission.

### CORRESPONDING AUTHOR:

Name:	Tianzhuo Zhao	
Department:	Laser Engineering Center	
Institution:	Academy of Opto-electronics, Chinese Academy of Sciences	
Article Title:	Realize Elements Concentration Analysis for Vacuum Induction Melting by Laser-Induced Breakdown Spectroscopy	
Signature:		Date: 2018.01.23

Please submit a signed and dated copy of this license by one of the following three methods:

- 1) Upload a scanned copy of the document as a pdf on the JoVE submission site;
- 2) Fax the document to +1.866.381.2236;
- 3) Mail the document to JoVE / Attn: JoVE Editorial / 1 Alewife Center #200 / Cambridge, MA 02139

For questions, please email [submissions@jove.com](mailto:submissions@jove.com) or call +1.617.945.9051



Dear editor:

We wish to thank you for the time and effort you have spent reviewing our paper.

Motivated by your comments, we have reconsidered the content of our work and tried to fix all the problems mentioned. In particular, this revised manuscript of our resubmitted letter has been improved mainly as follows:

1. American Journal Experts (AJE) edited language of the manuscript.
2. Each point of comments is revised, and we introduce how the manuscript is revised after this letter.

I will be very appreciated getting your advices soon. My telephone is 86-10-82178672, cell phone is +86 15010891828, and Email is zhaotianzhuo@aoe.ac.cn.

Yours,  
Sincerely  
Tianzhuo Zhao

How much is poured into the mold?

All the molten liquid steel. Revised in line 110 of the manuscript.

How is this testing done? What is actually done?

We send these samples to Central Iron & Steel Research Institute of China for test. It has official Authorization about ingredient test of alloys.

We do not know the test process exactly, and it is not important to the experiment of this paper.

You have chosen this for filming. Please provide more details here as need specific values in order to film.

Details parameter of spectrometer and laser has been added in the part of filming. The length of this part is increased from 122 words to 260 words.

Are these Figure numbers correct?

It has been correct now.

Where is this mentioned in the text?

It has been mentioned in Line 224 now.

Please spell out all journal names.

It is OK now.