Journal of Visualized Experiments Mouse distraction osteogenesis model --Manuscript Draft--

Article Type:	Invited Methods Article - JoVE Produced Video		
Manuscript Number:	JoVE57925R2		
Full Title:	Mouse distraction osteogenesis model		
Keywords:	distraction osteogenesis; Bone; Angiogenesis; endothelial progenitor cells; bone marrow stromal cells; callus		
Corresponding Author:	Fujio Masahito		
	JAPAN		
Corresponding Author's Institution:			
Corresponding Author E-Mail:	m-fujio@med.nagoya-u.ac.jp		
Order of Authors:	Masahito Fujio		
	Yusuke Osawa		
	Masaki Matsushita		
	Kota Ogisu		
	Shuhei Tsuchiya		
	Hiroshi Kitoh		
	Hideharu Hibi		
Additional Information:			
Question	Response		
Please indicate whether this article will be Standard Access or Open Access.	Standard Access (US\$2,400)		
Please indicate the city, state/province, and country where this article will be filmed . Please do not use abbreviations.	Nagoya University Graduate School of Medicine 65 Tsurumai-cho, Showa-ku, Nagoya 466-8550 Japan		

KEYWORDS:

distraction osteogenesis, bone, angiogenesis, endothelial progenitor cells, bone marrow stromal cells, callus

SHORT ABSTRACT:

We present a mouse tibial distraction osteogenesis model developed using a custom-made distractor. The use of a mouse as an analysis target is advantageous for advancing research.

LONG ABSTRACT:

Distraction osteogenesis (DO) is a surgical procedure that involves skeletal tissue regeneration without cell transplantation. A DO model consists of the following three phases: the latency phase after osteotomy and placement of the external distractor; the distraction phase, wherein the separated bone ends are gradually and continuously distracted; and the consolidation phase. This custom-made distractor used for DO is comprised of two incomplete acrylic resin rings and an expansion screw. The process was initiated by making a mold with silicone impression material and then creating the custom-made distractor. Dental resin was poured into the formwork made of silicone impression material, and it was allowed to polymerize to create the incomplete resin rings required for the custom-made distractor. These rings were fixed with an expansion screw using transparent resin. The custom-made distractor created via this approach was attached to the tibia of mice. The tibia was fixed to the device using one pair of 25-gauge needles proximally, one pair of 27-gauge needles distally, and acrylic resin. After a latency period of 5 days, distraction was initiated at a rate of 0.2 mm/12 h. The lengthening was continued for 8 days, resulting in a total gap of 3.2 mm. The mice were sacrificed 4 weeks after distraction. Bone formation in the distraction gap was confirmed using both radiography and histology.

INTRODUCTION:

Distraction osteogenesis (DO) is an established treatment method for a variety of skeletal disorders, such as limb length discrepancies, bone defects, and limb deformities¹. This unique treatment strategy is based on the "tension–stress principle" proposed by Ilizarov. The method requires several days for latency, several weeks for active distraction, and several months for consolidation until the mature bone is formed².

The local hypoxic conditions due to blockage of blood flow^{3,4} and mechanical stimulation^{5,6} are particularly important in the healing process of DO. Hypoxia-induced angiogenesis carries oxygen, nutrients, soluble factors and cells necessary for tissue repair locally through the blood flow. Mechanical stimulation by extension operation causes biological reactions such as differentiation of mesenchymal stem cells, bone formation, calcification and remodeling. Serial DO treatment allows the formation of not only hard tissues but also soft tissues, including nerves, muscles, blood vessels, and skin tissues, without the need for stem cell transplantation. Therefore, a DO model is considered to be an excellent model for analyzing the regeneration of various tissues.

Rabbits and dogs are the most widely used animals in basic research for DO; however, there are few analysis tools available for these animals. The use of a mouse DO model facilitates a more detailed analysis. It is particularly suitable for experiments using knockout mice. However, when using a mouse as an experimental animal, an extension device should be created. Here, we present a mouse tibial DO model developed using a custom-made distractor created using a dental laboratory tool and technique, which has been used in a previous study.

94 95

PROTOCOL:

96 97

All experiments were carried out in accordance with protocols approved by the Animal Care and Use Committee of our institution.

98 99 100

1. Preparation of a Mold for Creating the Custom-Made Distractor

101102

1.1. Make two incomplete rings (outer diameter, 20 mm; inner diameter, 10 mm), which are a part of the distractor, with one sheet of paraffin wax (145 mm x 74 mm) using an Evans wax carver.

104 105 106

107

103

1.1.1. Make 4 of the same pieces. Use a wax spatula heated with a gas burner. Stack the four rings for a thickness of 5 mm. Provide space for the expansion screw and needle, 8 mm x 2 mm and 5 mm x 2 mm respectively with the Evans wax carver (Figures 1A, 1B).

108 109

1.2. Embed wax patterns in the silicone impression material and create a mold for the resin rings (**Figure 1C**).

112

1.3. After curing the silicone impression material (approximately 5 min at room temperature), remove the wax patterns.

115116

2. Production of the Custom-Made Distractor

117

2.1. Apply petroleum jelly thinly and evenly to the silicone mold. Mix 3 g of polymerization dental resin and immediately pour it into the silicone mold. Set for 5 min at room temperature (RT).

121

2.2. Remove the polymerized resin ring from the mold and polish it with a carbide bar and a
 dental micromotor to remove the burr (about 1 min). Two resin rings are required for each
 custom-made distractor.

125

2.2.1. Repeat the procedure to create the required number of resin rings (**Figure 1D**).

127

Leave a gap of about 2 mm to fix the resin rings with utility wax and secure the outside
 of the expansion screw completely with transparent resin. Set the resin for 5 min at RT. When
 polymerization is complete, remove the expansion screw knob (Figure 1E).

3. **Surgical Protocol**

132 133

134 3.1. Use 8-week-old male mice.

135

136 3.2. Induce anesthesia with an intraperitoneal injection of medetomidine hydrochloride at 137 0.3 mg/g, midazolam at 4 mg/kg, and butorphanol tartrate at 5 mg/kg of body weight.

138

139 3.2.1. Carefully shave and disinfect the surgical area with 10% iodine solution, and then, 140 administer 0.5% lidocaine hydrochloride at the right lower limb.

141 142

143

144

Make a longitudinal skin incision (approximately 15 mm in length) at the right lower leg with a No. 15 scalpel. Bluntly separate the underlying muscles, taking care not to remove all of the periosteum. Approach from the outside to easily reach the fibula. Cut the fibula with scissors.

145 146 147

148

149

Grasp the ankle with narrow thin forceps, and use a 27-gauge needle to make a hole in the bone about 5 mm from the heel. The needle should penetrate the skin, bone, and skin in that order. When the needle penetrates, cut the tip and root with a nipper so that the needle is 150 about 15 mm.

151

152 3.4.1. Make another hole in the same manner, about 2-3 mm proximal. Hold around the ankle 153 and pass two 25-gauge needles under the knee in the same manner. After the needles 154 penetrate, cut the tips and roots with a nipper (Figure 1F).

155 156

3.5. Place the custom-made distractor such that it is parallel to the extension direction. Fix the needles and device with enough polymerization dental resin to fill the grooves of the device. Wait for the polymerization to complete (approximately 5 min) (Figure 1G).

158 159 160

157

3.6. Be careful not to damage the surrounding tissues. Cut the middle of the diaphysis of the tibia using a very thin cutting disc while applying a saline solution (1-2 mL).

161 162 163

Close the wound with a 4–0 nylon suture. 3.7.

164 165

166

Give subcutaneous injections of buprenorphine (0.1 mg/kg) for analgesia immediately after operation. Continue buprenorphine every 12 hours through postoperative day 7 and as needed thereafter.

167 168 169

Distraction Protocol

170

171 Note: There are various reports on the latency period and distraction rate, but here, 172 representative protocols are shown.

- 4.1. After a latency period of 5 days, start distraction at a rate of 0.2 mm/12 h. For extension, use the pin attached to the rapid expansion screw. Move the pin in the direction of the yellow arrow attached to the extension screw (**Figure 1E**).
- 4.2. Anesthesia is not necessary during extension. Hold the tail with the little finger and
 palm, and fix the extension device with the forefinger and thumb.
- 181 4.3. Perform extension (0.2 mm for 1/4 turn). Continue lengthening for 8 days, which will result in a total gap of 3.2 mm.

Analysis

186 5.1. Radiography analysis: After completion of extension, evaluate bone regeneration with computed tomography under general anesthesia by using 2.0% isoflurane.

Note: If care has been taken with regard to the position of the extension device, it is possible to evaluate it with simple radiography even with the device attached.

5.2. Histological analysis: Remove the apparatus carefully so as not to cause a fracture distraction site during sampling. Fix the sample in 10% neutral buffered formalin for 24 hours at room temperature. Decalcify with Morse's solution (10% sodium citrate, 20% formic acid) overnight at 4 °C.

REPRESENTATIVE RESULTS:

Figures 1A and 1B present incomplete rings (outer diameter, 20 mm; inner diameter, 10 mm; thickness, 5 mm) with paraffin wax. Two wax patterns were embedded in silicone impression material, and a mold for the resin rings (Figure 1C) was formed. Polymerized resin was immediately poured into this mold, and resin rings were obtained (Figure 1D). A custom-made distractor was created by combining two resin rings and an expansion screw (Figure 1E). Figure 2A presents typical radiography findings at 4 weeks after distraction in the model. Newly formed bone was observed at the distraction site. Figure 2B presents the hematoxylin and eosin staining results. A newly formed bone bridge was observed, and the newly formed bone could be easily distinguished from the native bone. There was no intraoperative fracture or postoperative infection in the 5 mice used in this experiment.

FIGURE AND TABLE LEGENDS:

Figure 1. Custom-made distractor. (A, B) The approximate dimensions of the ring are as follows: outer diameter, 20 mm; inner diameter, 10 mm; thickness, 5 mm. The spaces for the expansion screw (white arrowhead) and needle (black arrowhead) are secured. (C) Wax pattern embedded in the silicone impression material. (D) Completed resin rings. (E) The two resin rings and the expansion screw are secured with resin. (F) Four needles are inserted into the proximal and distal metaphysis of the tibia. (G) Distraction device placed in the mouse limb.

Page 4 of 6 revised October 2016

Figure 2. Radiological and histological findings. (A) Radiologically, newly formed bone is observed 4 weeks after distraction. (B) Histologically, the gap is filled with newly formed bone (bar = $500 \mu m$).

DISCUSSION:

When a large animal is used as an experimental model, a ready-made extension device can be used, and it is easy to obtain good fixation and assess the extension operation itself and the extension amount. However, when a mouse is used as an experimental model, it is necessary to develop some or all of the equipment. Isefuku *et al.* and Tay *et al.* made the device and created a mouse model^{7,8}. Carvahjo *et al.* adopted a method involving fixation of a ready-made expansion device (track distractor: KLS Martin) with a ligature wire⁹. We could also create a mouse DO model using existing products for extension devices and devising a fixation method. This approach would make it easy to manufacture the extension device, perform extension, and assess the extension amount. Moreover, it is expected that a tibial segmental defect model might be created using the same device. Additionally, the length of the defect can be adjusted in steps of 0.2 mm, and thus, the device is considered to be versatile.

In such an experiment, it is technically difficult to thread the needle through the tibia. In particular, the diameter of the distal part of the tibia is small, and some training is required to pass two needles. If sufficient bone is not noted, it is necessary to check whether good fixation is obtained. In order to obtain good fixation, it might be better to pass the needles at an angle and with a gap as wide as possible between the needles^{7,8}. Additionally, as the rate of healing and the histology during healing would change according to the method of osteotomy, it is necessary to perform the procedure in a certain manner.

The size of the experimental animal can be considered as a limiting factor of this experiment. There are no problems when the mouse is large; however, if it is too small, surgery itself may not be possible. Previously, we used 4-week-old mice and were able to perform experiments without any problems¹⁰. With regard to the evaluation period, extension devices could be installed without problems until 42 days after cutting the bone. However, as further observations were not made, additional research is necessary in the future. The occurrence of complications might be an issue when advancing experiments. Fracture during surgery might occur in approximately 1 out of 30 mice. Similarly, postoperative wound infection might occur in approximately 1 out of 30 mice.

In conclusion, we presented a mouse tibial DO model developed using a custom-made distractor.

ACKNOWLEDGMENTS:

The authors thank Ms. Makiko Kato for providing encouragement to complete this study. We also thank the Division of Experimental Animals and Medical Research Engineering, Nagoya University Graduate School of Medicine, for the housing of mice.

DISCLOSURES:


Page 5 of 6 revised October 2016


The authors have nothing to disclose.


261262263

REFERENCES:

- 1. Watson, J.T. Distraction osteogenesis. *The Journal of the American Academy of Orthopaedic Surgeons.* **14**, S168 174, (2006).
- 2. Ilizarov, G.A. The tension-stress effect on the genesis and growth of tissues: Part II. The
- influence of the rate and frequency of distraction. *Clinical Orthopaedics and Related Research*.
- 268 **239**, 263 285, (1989).
- 3. Wan, C., et al. Activation of the hypoxia-inducible factor-1 alpha pathway accelerates
- 270 bone regeneration. Proceedings of the National Academy of Sciences of the United States of
- 271 *America*. **105** (2), 686 691, (2008).
- 272 4. Fujio, M., et al. Stromal cell-derived factor-1 enhances distraction osteogenesis-
- 273 mediated skeletal tissue regeneration through the recruitment of endothelial precursors. *Bone.*
- **49** (4), 693-700, (2011).
- 5. Tong, L., et al. Focal adhesion kinase expression during mandibular distraction
- osteogenesis: evidence for mechanotransduction. Plastic and reconstructive surgery. 111 (1),
- 277 211-222, (2003).
- 278 6. Rhee, S.T., El-Bassiony, L, Buchman, S.R. Extracellular signal- related kinase and bone
- 279 morphogenetic protein expression during distraction osteogenesis of the mandible:
- in vivo evidence of mechanotransduction mechanism for differentiation and osteogenesis by
- mesenchymal precursor cells. *Plastic and reconstructive surgery.* **117** (7), 2243-2249, (2006).
- 7. Isefuku, S., Joyner, C.J., Simpson, A.H. A murine model of distraction osteogenesis. *Bone.*
- 283 **27** (5), 661 665, (2000).
- 284 8. Tay, B.K., Le, A.X., Gould, S.E., Helms, J.A. Histochemical and molecular analyses of
- distraction osteogenesis in a mouse model. *Journal of Orthopaedic Research*. **16** (5), 636 642,
- 286 (1998).
- 287 9. Carvalho, R.S., et al. The role of angiogenesis in a murine tibial model of distraction
- 288 osteogenesis. *Bone*. **34** (5), 849 861, (2004).
- 289 10. Osawa, Y., et al. Activated FGFR3 promotes bone formation via accelerating
- endochondral ossification in mouse model of distraction osteogenesis. *Bone.* **105**, 42 49,
- 291 10.1016/j.bone.2017.05.016 (2017).







Name of Material/ Equipment	Company	Catalog Number	Comments/Description	
Paraffin wax	YAMAHACHI DENTAL MFG. CO.	_	For preparation a mold for resin rings	
Labocone putty	GC Corporation	_	For preparation a mold for resin rings	
Utility wax	GC Corporation	_	For preparation a mold for resin rings	
Expansion screw	Ortho Dentaurum	600-301-30	Component of custom-made distractor	
			Immediate polymerization resin	
			Component of custom-made	
Unifast III	GC Corporation	_	distractor	
			Transparent resin	
			Component of custom-made	
Ortho Crystal	NISSIN	_	distractor	
25-gauge needle	TERUMO	NN-2516R	For custom-made distractor	
27-gauge needle	TERUMO	NN-2719S	For custom-made distractor	
	Chubu Kagaku Shizai			
ICR mouse	Corporation	_	Experimental animal	
Somnopentyl	Kyoritsu Seiyaku		Pentobarbital sodium salt	
	FUJIFILM Wako Pure Chemical			
Isoflurane	Corporation	099-06571	Isoflurane inhalation solution	



1 Alewife Center #200 Cambridge, MA 02140 tel. 617.945.9051 www.jove.com

ARTICLE AND VIDEO LICENSE AGREEMENT

Title of Article:	Mouse	distraction	osteogenesis	model	
Author(s):	Masahito Fujio	, Yusuke Osawa	, Masaki Matsushi	Ta, Kota Ogisu, Shuhei Tsuchi	ira, Hiroshi Kitoh.
				,	Hideharu Hibi
Item 1 (check one	box): The A	uthor elects t	o have the Mate	erials be made availabl	e (as described at
http://www.j	ove.com/aut	hor) via: 🔽	Standard Access	Open Access	
Item 2 (check one box	:):				
The Auth	or is NOT a U	nited States go	overnment employ	ree.	
		-	ernment employe ates government e	e and the Materials we employee.	re prepared in the
			nment employee l ates government e	out the Materials were N employee.	OT prepared in the

ARTICLE AND VIDEO LICENSE AGREEMENT

- 1. Defined Terms. As used in this Article and Video License Agreement, the following terms shall have the following meanings: "Agreement" means this Article and Video License Agreement; "Article" means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein; "Author" means the author who is a signatory to this Agreement; "Collective Work" means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; "CRC License" means the Creative Commons Attribution-Non Commercial-No Derivs 3.0 Unported Agreement, the terms and conditions of which can be found http://creativecommons.org/licenses/by-ncnd/3.0/legalcode; "Derivative Work" means a work based upon the Materials or upon the Materials and other preexisting works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; "Institution" means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; "JoVE" means MyJove Corporation, a Massachusetts corporation and the publisher of The Journal of Visualized Experiments; "Materials" means the Article and / or the Video; "Parties" means the Author and JoVE; "Video" means any video(s) made by the Author, alone or in conjunction with any other parties, or by JoVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion of the Article, and in which the Author may or may not appear.
- 2. <u>Background</u>. The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JoVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.
- 3. Grant of Rights in Article. In consideration of JoVE agreeing to publish the Article, the Author hereby grants to JoVE, subject to Sections 4 and 7 below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish. reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats. whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the "Open Access" box has been checked in Item 1 above, JoVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.



1 Alewife Center #200 Cambridge, MA 02140 tel. 617.945.9051 www.jove.com

ARTICLE AND VIDEO LICENSE AGREEMENT

- 4. Retention of Rights in Article. Notwithstanding the exclusive license granted to JoVE in Section 3 above, the Author shall, with respect to the Article, retain the non-exclusive right to use all or part of the Article for the non-commercial purpose of giving lectures, presentations or teaching classes, and to post a copy of the Article on the Institution's website or the Author's personal website, in each case provided that a link to the Article on the JoVE website is provided and notice of JoVE's copyright in the Article is included. All non-copyright intellectual property rights in and to the Article, such as patent rights, shall remain with the Author.
- 5. Grant of Rights in Video Standard Access. This Section 5 applies if the "Standard Access" box has been checked in Item 1 above or if no box has been checked in Item 1 above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby acknowledges and agrees that, Subject to Section 7 below, JoVE is and shall be the sole and exclusive owner of all rights of any nature, including, without limitation, all copyrights, in and to the Video. To the extent that, by law, the Author is deemed, now or at any time in the future, to have any rights of any nature in or to the Video, the Author hereby disclaims all such rights and transfers all such rights to JoVE.
- 6. Grant of Rights in Video Open Access. This Section 6 applies only if the "Open Access" box has been checked in Item 1 above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby grants to JoVE, subject to Section 7 below, the exclusive. royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Video in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Video into other languages, create adaptations, summaries or extracts of the Video or other Derivative Works or Collective Works based on all or any portion of the Video and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. For any Video to which this Section 6 is applicable, JoVE and the Author hereby grant to the public all such rights in the Video as provided in, but subject to all limitations and requirements set forth in, the CRC License.
- 7. Government Employees. If the Author is a United States government employee and the Article was prepared in the course of his or her duties as a United States government employee, as indicated in Item 2 above, and any of the licenses or grants granted by the Author hereunder exceed the scope of the 17 U.S.C. 403, then the rights granted hereunder shall be limited to the maximum rights permitted under such

- statute. In such case, all provisions contained herein that are not in conflict with such statute shall remain in full force and effect, and all provisions contained herein that do so conflict shall be deemed to be amended so as to provide to JoVE the maximum rights permissible within such statute.
- 8. <u>Likeness</u>, <u>Privacy</u>, <u>Personality</u>. The Author hereby grants JoVE the right to use the Author's name, voice, likeness, picture, photograph, image, biography and performance in any way, commercial or otherwise, in connection with the Materials and the sale, promotion and distribution thereof. The Author hereby waives any and all rights he or she may have, relating to his or her appearance in the Video or otherwise relating to the Materials, under all applicable privacy, likeness, personality or similar laws.
- 9. Author Warranties. The Author represents and warrants that the Article is original, that it has not been published, that the copyright interest is owned by the Author (or, if more than one author is listed at the beginning of this Agreement, by such authors collectively) and has not been assigned, licensed, or otherwise transferred to any other party. The Author represents and warrants that the author(s) listed at the top of this Agreement are the only authors of the Materials. If more than one author is listed at the top of this Agreement and if any such author has not entered into a separate Article and Video License Agreement with JoVE relating to the Materials, the Author represents and warrants that the Author has been authorized by each of the other such authors to execute this Agreement on his or her behalf and to bind him or her with respect to the terms of this Agreement as if each of them had been a party hereto as an Author. The Author warrants that the use, reproduction, distribution, public or private performance or display, and/or modification of all or any portion of the Materials does not and will not violate, infringe and/or misappropriate the patent, trademark, intellectual property or other rights of any third party. The Author represents and warrants that it has and will continue to comply with all government, institutional and other regulations, including, without limitation all institutional, laboratory, hospital, ethical, human and animal treatment. privacy, and all other rules, regulations, laws, procedures or guidelines, applicable to the Materials, and that all research involving human and animal subjects has been approved by the Author's relevant institutional review board.
- 10. JoVE Discretion. If the Author requests the assistance of JoVE in producing the Video in the Author's facility, the Author shall ensure that the presence of JoVE employees, agents or independent contractors is in accordance with the relevant regulations of the Author's institution. If more than one author is listed at the beginning of this Agreement, JoVE may, in its sole discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JoVE reserves the right, in its absolute and sole discretion and without giving any reason therefore, to accept or decline any work submitted to JoVE. JoVE and its employees, agents and independent contractors shall have



1 Alewife Center #200 Cambridge, MA 02140 tel. 617.945.9051 www.jove.com

ARTICLE AND VIDEO LICENSE AGREEMENT

full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JoVE has sole discretion as to the method of making and publishing the Materials, including, without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

11. Indemnification. The Author agrees to indemnify JoVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JoVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JoVE, making of videos by JoVE, or publication in JOVE or elsewhere by JoVE. The Author shall be responsible for, and shall hold JoVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to the making of a video by JoVE its employees, agents or independent contractors. All sterilization, cleanliness or decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's

expense. All indemnifications provided herein shall include JoVE's attorney's fees and costs related to said losses or damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JoVE, its employees, agents or independent contractors.

- 12. <u>Fees</u>. To cover the cost incurred for publication, JoVE must receive payment before production and publication the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JoVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.
- 13. <u>Transfer, Governing Law.</u> This Agreement may be assigned by JoVE and shall inure to the benefits of any of JoVE's successors and assignees. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to me one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement required per submission.

CORRESPONDING AUTHOR:

Name:	Masahito Fujio
Department:	Oral and Makillofacial Surgery
Institution:	Nagoya University Graduate School of Medicine.
Article Title:	Mouse distraction Ostegenesis model
Signature:	Masahiro Fyjio Date: 26,01,2018

Please submit a signed and dated copy of this license by one of the following three methods:

- 1) Upload a scanned copy of the document as a pfd on the JoVE submission site;
- 2) Fax the document to +1.866.381.2236;
- 3) Mail the document to JoVE / Attn: JoVE Editorial / 1 Alewife Center #200 / Cambridge, MA 02139

For questions, please email submissions@jove.com or call +1.617.945.9051

Editorial comments:

- 1. Please expand the Introduction to be 3-6 paragraphs with sufficient references.
- The authors expand the Introduction to be 3 paragraphs with references.
- 2. Additional details are needed in the manuscript. Please see the comments in the attached manuscript.
- The authors rewrote them basically according to the instructions of the editor. Also, I think that lack of information is sufficiently supplemented with video.

3. Veterinary Review:

The use of the needles is an unique method for providing external fixation for the fracture that was created (see attached picture). Thus, the use of the needles don't specifically concern me. However, the lack of analgesia is a major welfare concern and I am quite surprised that this was not brought up by the reviewers. Fracture models are considered to cause significant pain and only using lidocaine on the skin is unacceptable. The mice should have received pre-emptive analgesics prior to the surgery and been maintained on appropriate analgesics for the length of the study.

- We apologize for the incomplete description. I described the current anesthesia that improved the method that was done in the past, consideration for distress.