Journal of Visualized Experiments

A novel micro agar salt bridge electrode for analyzing proton turnover rate of recombinant membrane proteins --Manuscript Draft--

Article Type:	Invited Methods Article - JoVE Produced Video	
Manuscript Number:	JoVE58552R1	
Full Title:	A novel micro agar salt bridge electrode for analyzing proton turnover rate of recombinant membrane proteins	
Keywords:	Electrophysiology, Electrogenic transport, diffusion potential, Salt-bridge, protein turnover number, Nernst potential, uncoupling proteins, membrane transporter	
Corresponding Author:	E. Dr. Pohl University of Veterinary Medicine Vienna Vienna, AUSTRIA	
Corresponding Author's Institution:	University of Veterinary Medicine Vienna	
Corresponding Author E-Mail:	Elena.Pohl@vetmeduni.ac.at	
Order of Authors:	Jürgen Kreiter	
	E. Dr. Pohl	
Additional Information:		
Question	Response	
Please indicate whether this article will be Standard Access or Open Access.	Standard Access (US\$2,400)	
Please indicate the city, state/province, and country where this article will be filmed . Please do not use abbreviations.	A-1210 Veterinaerplatz, 1 Vienna, Austria	



Dr Ronald Myers Science editor JoVE

Institute for Physiology, Pathophysiology und Biophysics

Elena E. Pohl, MD University Professor

Tel. +43 1 25077 4570

Tel. +43 1 25077 4571 (Secretary)

Fax +43 1 25077 4595 Email: elena.pohl@vetmeduni.ac.at

Homepage: http://www.vetmeduni.ac.at/en/biophysics/

Vienna, 18/07/18

Dear Editor,

We resubmit our manuscript with a modified title: "A novel micro agar salt bridge electrode for analysing proton turnover rate of recombinant membrane proteins" by Elena E. Pohl and Jürgen Kreiter, which now addresses all of the points that were raised by the reviewers.

In the course of the manuscript revision we additionally changed the term "Liquid Junction Potential" to "Diffusion Potential", which is more precise for the described setup.

We very much appreciate all of the constructive feedback that our original submission received. After having addressed the issues that had been raised, we feel that the quality of the paper has improved and we hope you will also agree.

Sincerely yours,

TITLE:

2 A Micro-agar Salt Bridge Electrode for Analyzing the Proton Turnover Rate of Recombinant

Membrane Proteins

4 5

1

3

AUTHORS & AFFILIATIONS:

6 Jürgen Kreiter¹, Elena E. Pohl¹

7 8

¹Department of Physiology and Biophysics, University of Veterinary Medicine, Vienna, Austria

9 10

Corresponding Authors:

Elena E. Pohl (Elena.Pohl@vetmeduni.ac.at 11

12 Tel: +43-1-25077-4570

13

14 Jürgen Kreiter (Juergen.Kreiter@vetmeduni.ac.at)

15 Tel: +43-1-25077-4503

16 17

KEYWORDS:

Electrophysiology, electrogenic transport, diffusion potential, protein turnover number, Nernst potential, uncoupling proteins

19 20 21

22

23

24

18

SUMMARY:

In electrophysiological measurements, the presence of a diffusion potential disturbs the precise measurement of the reverse potential by altering the electrode potential. Using a micro-agar salt bridge, the impact of the diffusion potential is minimized, which allows a more precise measurement of substrate turnover numbers of reconstituted recombinant membrane proteins.

25 26 27

28

29

30

31

32

33

34

35

36

37

38

39

40

41

42

43

44

ABSTRACT:

To date, more than 50% of all pharmacological drugs target the transport kinetics of membrane proteins. The electrophysiological characterization of membrane carrier proteins reconstituted in lipid bilayer membranes is a powerful but delicate method for the assessment of their physicochemical and pharmacological properties. The substrate turnover number is a unique parameter that allows the comparison of the activity of different membrane proteins. In an electrogenic transport, the gradient of the translocated substrate creates a membrane potential that directly correlates to the substrate turnover rate of the protein. By using silver chloride electrodes, a diffusion potential, also called liquid junction potential, is induced, which alters electrode potential and significantly disturbs precise membrane potential measurements. Diffusion potential can be minimized by a salt bridge, which balances electrode potential. In this article, a micro-agar salt bridge is designed to improve the electrophysiological set-up, which uses micropipettes for the membrane formation. The salt solution is filled into a microcapillary pipette tip, stabilized by the addition of agarose, and can be easily mounted to a standard electrode. The electrode potential of a micro-salt bridge electrode is more stable compared to a standard electrode. The implementation of this system stabilizes electrode potential and allows more precise measurements of membrane potential generated by a pH gradient. Using this system, the proton turnover rates of the mitochondrial carriers UCP1 and UCP3 are reinvestigated and compared to earlier measurements.

INTRODUCTION:

Membrane proteins are targeted by up to 60% of all known pharmaceutical drugs¹. Electrophysiological measurements of membrane proteins are a powerful but delicate tool to analyze the electrogenic transport of substrates mediated by membrane carrier proteins. The modulation of the transmembrane current by the application of constant voltages or voltage ramps allows assessing the pharmacological and physical properties of the carriers, for instance, the activation and inhibition by substrates or the transport kinetics. Of special interest is the substrate turnover number, which displays the amount of substrate that is translocated by a membrane protein per time unit. It is a major parameter when comparing the kinetics of various membrane proteins. Establishing a concentration gradient of the charged substrate across the membrane generates an electromotive force from which the turnover number of the substrate is deduced.

Using an AgCl electrode, the presence of a chloride-free buffer creates a diffusion potential that alters electrode potential and leads to a shift in current-voltage measurements². Although always present, it is negligible for standard conductance and capacity measurements since these parameters are either dependent on the slope of the current-voltage recording (conductance) or are the difference of a single recording (capacity), which cancels the potential. However, the recording of the reverse potential, which is created by the transport of substrate, can be significantly disturbed by the diffusion potential. Thus, for exact measurements of the reverse potential, the electrode potentials have to be kept constant.

The diffusion potential can be minimized by two methods: (i) in the presence of a bilayer membrane, a substrate concentration has to be increased on one side of the membrane^{3,4}, or (ii) a salt bridge balances the electrode potential⁵. The first method is highly dependent on the stability of the measurements. The membrane has to survive for several minutes, from the addition of substrate under stirring until the substrate is almost equally distributed in the solution. If the membrane ruptures in between, the substrate gradient is altered by the free exchange of charged molecules, and measurements turn inaccurate. The latter method balances the diffusion potential but is limited by the size of the set-up. Implementing a small but functioning salt bridge in a micro range to an electrophysiological set-up is challenging⁶. For the latter method, the salt solution is filled into a microcapillary tip and stabilized by the addition of agarose to prevent diffusion of the salt solution to the buffer solution.

In this protocol, a straightforward production of a micro-agar salt bridge and implementation into an electrophysiological set-up based on the pipette set-up⁷ is described. A microcapillary tip is adjusted both to contain a 3 M KCl solution with 1 mol% (w/v) agarose and to bridge an AgCl electrode and buffer solution. The advantage of the micro-salt bridge is displayed by time recordings of the electrode potential shift and the more precise measurements of membrane potential at various pH gradients. In the model system of recombinant proteins reconstituted into liposomes, the turnover rates of mitochondrial carriers UCP1 and UCP3 produced under similar conditions are reinvestigated and compared to previous results^{3,8}.

PROTOCOL: 1. Production of Recombinant Uncoupling Proteins (UCPs) and the Formation of Planar Bilayer **Membranes** 1.1. Produce recombinant UCP1 and UCP3 as described by Rupprecht et al.⁹ and Hilse et al.¹⁰ 1.2. Form planar bilayer membranes on the tips of conventional dispensable plastic pipettes as described by Beck et al.⁷ 2. Preparation of the Micro-agar Salt Bridge Electrode 2.1. Adjust a microcapillary pipette tip (see **Table of Materials**) to the appropriate length. 2.1.1. Mark the position on an empty tip at which the buffer-containing tip will be attached and use a sliding caliper to measure the length of the micro-agar salt bridge electrode. CAUTION: The microcapillary tip must be long enough to enter the buffer solution for the measurements. 2.1.2. Cut the microcapillary tip with a sharp knife or blade and clean the cut surface with ethanol and water. 2.2. Prepare the AgCl-coated electrode. 2.2.1. Cut off an Ag wire of approximately 8 cm in length and clean it with ethanol and water. 2.2.2. Take a piece of sandpaper and smooth down the surface of 1-cm length at the end of the wire, which should be in contact with the salt solution. 2.2.3. Dip the smoothed end into a 3 M KCl solution and coat it electrochemically with chloride by using a DC supply at 1.5 V for 10 s. 2.2.4. Clean the electrode with water, dry it, and adjust the length of the AgCl electrode from the uncoated side so that it penetrates the microcapillary tip as deep as possible. Note: The protocol can be paused here. 2.3. Prepare a 3 M KCl salt solution with 1% (v/v) agarose. 2.3.1. Weigh out 4.47 g of KCl and dissolve it in 20 mL of water by stirring it in a flask. 2.3.2. Remove the stirrer, weigh out 0.2 g of agarose, and add it to the flask.

2.3.3. Heat up the solution to 100 °C to melt the agarose and prevent it from clotting. CAUTION: The flask will be very hot. Do not touch it with bare hands. Use gloves for handling. 2.4. Fill the microcapillary tip with the agarose salt solution. CAUTION: The solution is hot. Protect hands and work carefully to avoid splashes. 2.4.1. If the agarose starts clotting, heat up the salt solution to fully melt the agarose again. 2.4.2. Soak 10 µL of the agar salt solution into the microcapillary tip. Soak it slowly and carefully to avoid air bubbles in the tip. 2.4.3. Remove the tip from the pipette and push in the AgCl electrode from the broader side of the tip. Ensure that the electrode penetrates the salt solution. 2.4.4. Cool down the electrode to room temperature and plug it into the amplifier. 2.5. Prepare the buffer for the measurements. 2.5.1. Weigh out 0.710 g of Na₂SO₄, 0.195 g of MES, 0.121 g of TRIS, and 0.023 g of EGTA, then add 100 mL of distilled water to a beaker and stir the solution. 2.5.2. Check for the pH value of the buffer by using a pH electrode and adjust the pH value to 7.32.

2.5.3. Check that the reference electrode and the agar salt bridge electrode are in electrical contact.

2.5.3.1. Add 1 mL of buffer to a plastic container.

2.5.3.2. Dip in the reference electrode and the agar salt bridge electrode and check for a signal response. If the signal response is correct, proceed to step 2.7.

2.6. If there is a false signal response or no electrical contact at all, perform the following troubleshooting.

2.6.1. Check if the electrode is in contact with the salt solution and push the electrode into the solution.

Note: If the solution is already too sticky, remove the micro-agar salt bridge and prepare a new microcapillary tip.

- 177 2.6.2. Check if there are air bubbles within the salt solution. If yes, prepare a new microcapillary 178 179 180 2.6.3. Check if the salt solution is in contact with the buffer solution. If not, then cut off another 181 1 mm of the tube end of the tip. 182 183 CAUTION: Be sure that the tip is still long enough to penetrate the buffer-containing tip. If there 184 is still no contact, prepare a new microcapillary tip. 185 186 2.6.4. If none of these steps help, prepare a new microcapillary tip. 187 188 2.7. For storage, dip the agar salt bridge electrode into a 3 M KCl solution. 189 190 Note: The protocol can be stopped here. For a pause overnight, store the electrode in 3 M KCl 191 salt solution at 4 °C. 192 193 2.8. Prepare the buffer-containing plastic tip. 194 195 Note: If the electrode was stored overnight, take it out and let it heat up to room temperature 196 for 30 min. 197 2.8.1. Take a microcapillary tip and bend the tube, 2 cm from the edge of the narrow part, 198 199 approximately 90 degrees using a heating wire. 200 201 2.8.2. Use a very sharp knife or blade and cut off the tube around 5 mm from the bend. 202 203 2.8.3. Clean the surface at the end with ethanol and water and measure the diameter of the hole 204 of the tip. From this, calculate the area of the surface. 205 206 2.8.4. Coat the surface of the tip with 85:15 (v:v) hexane:hexadecane. 207
- 208 2.8.4.1. Pipette 3 μL of the solvent and remove it from the tip.

209

211

213

216

218

- 2.8.4.2. Wait 1 min so that all of the residual solvent in the tip has evaporated.
- 2.8.5. Fill the measurement tip with 3 μL of buffer and plug it into the salt bridge electrode.
- Note: Check again if the salt bridge electrode and the reference electrode are in electrical contact by performing step 2.5.3.
- 2.8.6. If there is no electrical contact, check for the following:
- 2.8.6.1. Check whether the salt bridge electrode is in contact with the buffer solution. If not, then
 either increase the volume of the buffer in the tip or prepare a longer microcapillary tip.

2.8.6.2. If there is an air bubble, remove the tip from the micro-salt bridge electrode, refill the buffer in the tip, and plug it into the electrode again.

3. Measurement of the Electrical Parameters of the Membrane Reconstituted with Recombinant Protein

3.1. Apply a triangular alternating voltage signal with maximal voltage $U_{max} = 50$ mV and $\Delta T_{ramp} = 50$ ms, which creates a rectangular alternating current response. From the mean values of the positive and negative currents (I₊ and I₋), calculate the capacity of the membrane according to the following formula:

$$C = \frac{(I_{+} - I_{-}) \cdot \Delta T_{ramp}}{2 \cdot U_{max}}$$

3.2. Apply a voltage ramp ranging from -50 mV to + 50 mV and record the current. Fit a linear function to the data—the slope is the conductance—and calculate the x-axis intersection point of the fit.

3.3. Discharge the solution in the plastic tip and fill a new one with a buffer containing an increased substrate concentration.

3.3.1. If no membrane is formed within the first 20 - 30 s, discharge the volume and refill it. This guarantees that the concentration gradient across the membrane does not significantly change during the membrane formation.

246 3.3.2. After the formation of a membrane, perform steps 3.1 and 3.2 again to verify proper membrane formation and to get the x-axis intersection point.

4. Calculation of the Substrate Turnover Rate

Note: See previous work for details^{3,7}.

4.1. Estimate the amount of protein in the membrane from the molecular mass of lipid and protein (M_{Lipid} and $M_{Protein}$), the area of the membrane and of one lipid head group ($A_{Membrane}$ and A_{Lipid}), and the mass ratio of the protein per lipid (r).

$$N_{Protein} = \frac{M_{Lipid} \cdot A_{Membrane}}{M_{Protein} \cdot A_{Lipid}} \cdot r$$

4.2. Calculate the Nernst potential for the transported substrate. R is the gas constant, T the temperature, z the charge of the transported substrate, F Faraday's constant, and c_1 and c_2 the concentrations of the substrate of both sides of the membrane.

$$\Phi_{Nernst} = \frac{RT}{zF} \cdot \ln \left(\frac{c_1}{c_2} \right)$$

4.3. From the current-voltage recordings, take the reverse potential calculated with the difference of the x-axis intersection points from the linear fits in the presence and absence of the substrate gradient.

4.4. Calculate the proportion of substrate conductance, $G_{Substrate}$, to the total membrane conductance, G_{total} , by the ratio of reverse potential to Nernst potential¹¹.

$$G_{Substrate} = \frac{\Phi_{reverse}}{\Phi_{Nernst}} \cdot G_{total}$$

4.5. Calculate the substrate turnover number $\Delta N_{Substrate}$ per time unit ΔT from the substrate conductance (G_{Substrate}), the applied voltage U, and the charge of the substrate z:

$$\frac{\Delta N_{Substrate}}{\Delta T} = \frac{G_{Substrate} \cdot U}{z \cdot e}$$

4.6. From the ratio of transported substrate per time to the number of proteins, calculate the
 turnover rate κ.

$$\kappa = \frac{\Delta N_{Substrate}/\Delta T}{N_{Protein}}$$

REPRESENTATIVE RESULTS:

To verify the minimization of the diffusion potential, the stability of the current-voltage measurements of an intact membrane was measured. In **Figure 2A**, representative-current voltage recordings are depicted in the presence (white dots) and absence (black dots) of a pH gradient. According to the Nernst equation, the pH gradient induces a shift in voltage. From the x-axis intersection point of a linear fit to the data, the potential shift is calculated. In order to test both electrodes, the shift in x-axis intersection point was analyzed for a standard AgCl electrode (**Figure 2B**; white dots) and a micro-agar salt bridge (**Figure 2B**; black dots). A voltage ramp was recorded ten times in a row and the mean shift in the x-axis is depicted against time. Whereas the agar salt bridge electrode had a maximum shift of less than 5 mV even after 300 seconds of measurement, the standard electrode varied up to 30 mV in unpredictable, random, behavior.

Next, both electrodes were tested at different pH gradients (**Figure 3A**). For the standard electrode, a pH gradient of 0.35 and 1.0 (white dots) was generated; for the agar salt bridge electrode, pH gradients of 0.35, 0.7, and 1.0 (black dots). The shift in potential was analyzed in three independent measurements. In contrast to a gradient of 0.35, where the measured shift only varies slightly, the voltage shift significantly alters at the pH gradient of 1.0 in the absence of the micro-agar salt bridge. From a linear fit to the data, the slope of the function is 26.4 ± 2.3 mV/ Δ pH for the standard electrode and 50.1 ± 4.6 mV/ Δ pH for the micro-agar salt bridge electrode. According to the Nernst equation, the calculated potential shift is 60.7 mV/ Δ pH at T =

303 32 °C.

Using the micro-agar salt bridge, the proton turnover number, κ , of mitochondrial UCP1 and UCP3 was measured and compared to previous measurements (**Figure 3B**). Similar to **Figure 3A**, $\Delta pH = 1.0$ was generated and the reverse potential was measured. The amount of protein in the membrane was estimated according to the formula provided in section 4 of the protocol, with a protein to lipid ratio of $4 \mu g/(mg \text{ of lipid})$, a molecular mass of 33,000, and 750 Da for the protein and lipid, a membrane surface area of $3.53 \times 10^{-4} \text{ cm}^2$, and an area per lipid of $7.8 \times 10^{-15} \text{ cm}^2$. The obtained κ was $5.56 \pm 0.38 \text{ s}^{-1}$ and $4.10 \pm 0.71 \text{ s}^{-1}$ for UCP1 and UCP3, respectively (**Figure 3B**).

FIGURE AND TABLE LEGENDS:

Figure 1: The electrophysiological set-up with the micro-agar salt bridge. (A) This panel shows a sketch of the set-up. The microcapillary tip containing the agar salt solution (depicted in orange) is placed between the electrode (black) and the buffer-containing tip (blue). The membrane is formed at the surface at the end of the buffer-containing tip (indicated by the arrow). (B) This panel shows an image of the electrophysiological set-up with the implementation of the microagar salt bridge. The arrows point to the electrode, the micro-agar salt bridge, the reference electrode, and the container with buffer solution.

Figure 2: The comparison of a micro-agar salt bridge and a standard AgCl electrode. (A) This panel shows a representative current-voltage measurement in the presence (grey dots) and absence (white dots) of a pH gradient of 1. The lines represent a linear fit to the data, from which conductance and x-axis intersection values are obtained. The voltage shift is evaluated by the difference of the intersection values of both recordings. (B) This panel shows the shift of membrane potential of a standard AgCl electrode (white dots) to a micro-agar salt bridge electrode (black dots) in time. Ten current-voltage measurements were recorded in a row and the mean voltage shift by diffusion potential is plotted against time. In all experiments, the membrane was made of 45:45:10 mol% DOPC:DOPE:CL reconstituted with 15 mol% arachidonic acid at a concentration of 1.5 mg/mL. The buffer contained 50 mM Na₂SO₄, 10 mM MES, 10 mM TRIS, and 0.6 mM EGTA at pH = 7.34 and T = 32 °C.

Figure 3: Proton turnover number of UCP1 and UCP3 calculated from the reverse potential in the presence of a pH gradient. (A) This panel shows the shift in potential of a UCP1-containing membrane of various pH gradients of a standard AgCl electrode (white dots) and a micro-agar salt bridge electrode (black dots). The lines represent a linear fit to the data. (B) This panel shows the proton turnover number of UCP1 (first data set) and UCP3 (second data set) as calculated from the ratio of voltage shift to Nernst potential according to the formulas in section 4 of the protocol. The first bar of each set represents turnover rates measured with the micro-agar salt bridge. The second bar of each data set represents previous measurements using a standard AgCl electrode. Values for UCP1 and UCP3 were taken from Urbankova *et al.*³ and Macher *et al.*⁸. In all measurements, the membrane was made of 45:45:10 mol% DOPC:DOPE:CL reconstituted with 15 mol% AA and UCP1/UCP3. The concentration of lipid and protein was 1.5 mg/mL and 4 μ g/mg of lipid, respectively. The buffer contained 50 mM Na₂SO₄, 10 mM MES, 10 mM TRIS, and 0.6 mM

EGTA at pH = 7.34 and T = 32 °C. The pH of the buffer for the gradient measurements was increased to 7.66, 8.00, or 8.33 by adding TRIS and was changed in the solution-containing pipette. The values are the mean \pm standard deviation of three independent measurements.

DISCUSSION:

 The implementation of the micro-agar salt bridge with the electrode minimizes its diffusion potential and allows more precise measurements of membrane potential generated by a pH gradient. In the presence of various transmembrane pH gradients, the potential shift of both electrodes was acceptable at $\Delta pH = 0.35$ when comparing to the theoretical value of the Nernst equilibrium potential ($\Phi_{Nernst} = 23.8$ mV for $\Delta pH = 0.4$). However, at more physiological pH gradients, as for instance in mitochondria between the matrix and intermembrane space, the standard AgCl electrode failed to precisely measure the potential shift at $\Delta pH = 1$ (**Figure 3A**). The electrode bridged with micro-agar salt delivered the values which were much more comparable to the theory.

Diffusion potential may also occur at the AgCl reference electrode if the buffer solution is changed during the experiment. Chloride-free buffer solution was used in the experiments since uncoupling proteins were suggested to transport chloride ions, and the pH was adjusted using Tris or MES. The electrode potential, in the absence of a substantial concentration of chloride, primarily depends on chloride impurities in the buffer solution. As its composition is unchanged during the experiments, it will simply result in a constant offset potential. However, for the measurement of an absolute potential difference between the two electrodes, a simple agar salt-bridge system (Ag/AgCl 3 M KCl) could also be used for the reference electrode.

A micro-agar salt bridge balances the diffusion potential by an equilibration of the electrode potential. In order to stabilize the salt solution, 1% (w/v) agarose was added to prevent the mixing of the salt solution with the buffer solution. The salt ions K^+ and Cl^- have similar mobilities in liquid and balance the electrode potential. To properly install the salt bridge, the agar salt solution has to be sufficiently heated up to fill the microcapillary tip without any air bubbles and to cover the AgCl electrode. Before further usage, electrical contact between the salt bridge electrode and the reference electrode has to be checked. Depending on the time the salt bridge is used, the salt solution has to be sufficiently gelled to prevent any mixing of the salt solution with the buffer. This is especially critical if K^+ or Cl^- transporters are investigated. The salt bridge was used for a very short time and the elution of agarose is negligible in this time range. A higher concentration of agarose of up to 5%, or of agar (3% - 5%), allows using the salt bridge for a longer period of time^{6,12}.

This method allows determining the transport kinetics of a membrane transporter (i) with low turnover rates and (ii) of mitochondrial proteins of the inner membrane, which can hardly be investigated in standard patch clamp set-ups¹³. Its precision is mainly dependent on the reverse potential measurement, which accuracy is decreased at a low total membrane conductance and small concentration gradients which induce a membrane potential below the noise of recording.

Using this set-up, the turnover rates of UCP1 and UCP3 as produced under the same conditions

were measured. Due to the higher pH gradient, the obtained rates seem to be more precise and unperturbed by artefacts resulting from the minor electrode potential shift. It can be used to further analyze and compare mitochondrial membrane transporters produced under similar conditions.

395 396

ACKNOWLEDGMENTS:

This work was supported by the Austrian Research Fund (P31559-B20 to E.E.P.). The authors thank Sarah Bardakji for the excellent technical assistance in the production and reconstitution of mouse UCP1 and UCP3 into proteoliposomes.

400 401

DISCLOSURES:

The authors have nothing to disclose.

402 403 404

REFERENCES:

405 1. Terstappen, G. C., Reggiani, A. In silico research in drug discovery. *Trends in Pharmacological* 406 *Sciences.* **22** (1), 23-26 (2001).

407

2. Raynauld, J. P., Laviolette, J. R. The silver-silver chloride electrode: A possible generator of offset voltages and currents. *Journal of Neuroscience Methods.* **19** (3), 249-255 (1987).

410

3. Urbankova, E., Voltchenko, A., Pohl, P., Jezek, P., Pohl, E. E. Transport kinetics of uncoupling proteins. Analysis of UCP1 reconstituted in planar lipid bilayers. *Journal of Biological Chemistry*. **278** (35), 32497-32500 (2003).

414

4. Beck, V. *et al.* Polyunsaturated fatty acids activate human uncoupling proteins 1 and 2 in planar lipid bilayers. *The FASEB Journal.* **21** (4), 1137-1144 (2007).

417

5. Shao, X. M., Feldman, J. L. Micro-agar salt bridge in patch-clamp electrode holder stabilizes electrode potentials. *Journal of Neuroscience Methods*. **159** (1), 108-115 (2007).

420

421 6. Kleene, S. J. A simple intrapipette salt bridge. *Journal of Neuroscience Methods.* **46** (1), 11-16 (1993).

423

7. Beck, V. et al. A new automated technique for the reconstitution of hydrophobic proteins
 into planar bilayer membranes. Studies of human recombinant uncoupling protein 1.
 Biochimica et Biophysica Acta-Bioenergetics. 1757 (5-6), 474-479 (2006).

427

428 8. Macher, G. *et al.* Inhibition of mitochondrial UCP1 and UCP3 by purine nucleotides and phosphate. *Biochimica et Biophysica Acta*. **1860** (3), 664-672 (2018).

430

9. Rupprecht, A. *et al.* Role of the transmembrane potential in the membrane proton leak.
Biophysical Journal. 98 (8), 1503-1511 (2010).

433

10. Hilse, K. E. et al. The expression of UCP3 directly correlates to UCP1 abundance in brown


436
 437 11. Fuks, B., Homble, F. Permeability and electrical properties of planar lipid membranes from


adipose tissue. Biochimica et Biophysica Acta. 1857 (1), 72-78 (2016).


435

443

- 438 thylakoid lipids. *Biophysical Journal*. **66** (5), 1404-1414 (1994).
- 12. Barry, P. H., Lewis, T. M., Moorhouse, A. J. An optimised 3 M KCl salt-bridge technique used
 to measure and validate theoretical liquid junction potential values in patch-clamping and
 electrophysiology. *European Biophysics Journal.* 42 (8), 631-646 (2013).
- Huang, D., Li, J. The feasibility and limitation of patch-clamp recordings from neonatal rat
 cardiac ventricular slices. *In Vitro Cellular & Developmental Biology.* 47 (4), 269-272 (2011).







Name of Material/ Equipment

Company

Microloader tips Eppendorf

Ethanol 99% AustrAlco Österr. Agrar-Alkohol Handelsges.m.b.H

Kaliumchlorid Carl Roth GmbH + Co. Kg

DC supply Voltcraft

Agarose Standard Carl Roth GmbH + Co. Kg

Patch Clamp Amplifier Heka

Sample tube Carl Roth GmbH + Co. Kg

Na2SO4 Carl Roth GmbH + Co. Kg

MES Carl Roth GmbH + Co. Kg

TRIS Carl Roth GmbH + Co. Kg

EGTA Carl Roth GmbH + Co. Kg

Hexane Sigma-Aldrich

Hexadecane Sigma-Aldrich

Heating wire Voltcraft

Catalog Number

Comments/Description

5242956.003

Microcapillary pipette tip

AAAH-5020-07025-230317

6781.3

V10/CPG 1940 -01

3810.2

5863.1

8560.3

4256.2

AE15.2

3054.1

296090-100ML

296317-100ML

USPS-2250



1 Alewife Center #200 Cambridge, MA 02140 tel. 617.945.9051 www.jove.com

ARTICLE AND VIDEO LICENSE AGREEMENT

Title of Article:		, A novel micro agar salt bridge electrode for analyzing proton turnover rate of recombina	
		nt membrane proteins.	
Author(nor(s): Elena E. Pohl and Jürgen Kreiter		
ltem 1	(check one	box): The Author elects to have the Materials be made available (as described at	
	http://www.j	ove.com/author) via: X tandard Access	
	X The Auth	or is NOT a United States government employee.	
	The Author is a United States government employee and the Materials were prepared in the course of his or her duties as a United States government employee.		
		or is a United States government employee but the Materials were NOT prepared in the or her duties as a United States government employee.	

ARTICLE AND VIDEO LICENSE AGREEMENT

- 1. Defined Terms. As used in this Article and Video License Agreement, the following terms shall have the following meanings: "Agreement" means this Article and Video License Agreement; "Article" means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein: "Author" means the author who is a signatory to this Agreement; "Collective Work" means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; "CRC License" means the Creative Commons Attribution-Non Commercial-No Derivs 3.0 Unported Agreement, the terms and conditions of which can be found http://creativecommons.org/licenses/by-ncnd/3.0/legalcode; "Derivative Work" means a work based upon the Materials or upon the Materials and other preexisting works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; "Institution" means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; "JoVE" means MyJove Corporation, a Massachusetts corporation and the publisher of The Journal of Visualized Experiments; "Materials" means the Article and / or the Video; "Parties" means the Author and JoVE; "Video" means any video(s) made by the Author, alone or in conjunction with any other parties, or by JoVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion of the Article, and in which the Author may or may not appear.
- 2. <u>Background</u>. The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JoVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.
- 3. Grant of Rights in Article. In consideration of JoVE agreeing to publish the Article, the Author hereby grants to JoVE, subject to Sections 4 and 7 below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the "Open Access" box has been checked in Item 1 above, JoVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.



1 Alewife Center #200 Cambridge, MA 02140 tel. 617.945.9051 www.jove.com

ARTICLE AND VIDEO LICENSE AGREEMENT

- 4. Retention of Rights in Article. Notwithstanding the exclusive license granted to JoVE in Section 3 above, the Author shall, with respect to the Article, retain the non-exclusive right to use all or part of the Article for the non-commercial purpose of giving lectures, presentations or teaching classes, and to post a copy of the Article on the Institution's website or the Author's personal website, in each case provided that a link to the Article on the JoVE website is provided and notice of JoVE's copyright in the Article is included. All non-copyright intellectual property rights in and to the Article, such as patent rights, shall remain with the Author.
- 5. <u>Grant of Rights in Video Standard Access</u>. This **Section 5** applies if the "Standard Access" box has been checked in **Item 1** above or if no box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby acknowledges and agrees that, Subject to **Section 7** below, JoVE is and shall be the sole and exclusive owner of all rights of any nature, including, without limitation, all copyrights, in and to the Video. To the extent that, by law, the Author is deemed, now or at any time in the future, to have any rights of any nature in or to the Video, the Author hereby disclaims all such rights and transfers all such rights to JoVE.
- 6. Grant of Rights in Video Open Access. This Section 6 applies only if the "Open Access" box has been checked in Item 1 above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby grants to JoVE, subject to Section 7 below, the exclusive, royalty-free perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Video in all forms. formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Video into other languages, create adaptations, summaries or extracts of the Video or other Derivative Works or Collective Works based on all or any portion of the Video and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. For any Video to which this Section 6 is applicable, JoVE and the Author hereby grant to the public all such rights in the Video as provided in, but subject to all limitations and requirements set forth in, the CRC License.
- 7. <u>Government Employees</u>. If the Author is a United States government employee and the Article was prepared in the course of his or her duties as a United States government employee, as indicated in **Item 2** above, and any of the licenses or grants granted by the Author hereunder exceed the scope of the 17 U.S.C. 403, then the rights granted hereunder shall be limited to the maximum rights permitted under such

- statute. In such case, all provisions contained herein that are not in conflict with such statute shall remain in full force and effect, and all provisions contained herein that do so conflict shall be deemed to be amended so as to provide to JoVE the maximum rights permissible within such statute.
- 8. Likeness, Privacy, Personality. The Author hereby grants JoVE the right to use the Author's name, voice, likeness, picture, photograph, image, biography and performance in any way, commercial or otherwise, in connection with the Materials and the sale, promotion and distribution thereof. The Author hereby waives any and all rights he or she may have, relating to his or her appearance in the Video or otherwise relating to the Materials, under all applicable privacy, likeness, personality or similar laws.
- 9. Author Warranties. The Author represents and warrants that the Article is original, that it has not been published, that the copyright interest is owned by the Author (or, if more than one author is listed at the beginning of this Agreement, by such authors collectively) and has not been assigned, licensed, or otherwise transferred to any other party. The Author represents and warrants that the author(s) listed at the top of this Agreement are the only authors of the Materials. If more than one author is listed at the top of this Agreement and if any such author has not entered into a separate Article and Video License Agreement with JoVE relating to the Materials, the Author represents and warrants that the Author has been authorized by each of the other such authors to execute this Agreement on his or her behalf and to bind him or her with respect to the terms of this Agreement as if each of them had been a party hereto as an Author. The Author warrants that the use, reproduction, distribution, public or private performance or display, and/or modification of all or any portion of the Materials does not and will not violate, infringe and/or misappropriate the patent, trademark, intellectual property or other rights of any third party. The Author represents and warrants that it has and will continue to comply with all government, institutional and other regulations, including, without limitation all institutional, laboratory, hospital, ethical, human and animal treatment, privacy, and all other rules, regulations, laws, procedures or guidelines, applicable to the Materials, and that all research involving human and animal subjects has been approved by the Author's relevant institutional review board.
- 10. JoVE Discretion. If the Author requests the assistance of JoVE in producing the Video in the Author's facility, the Author shall ensure that the presence of JoVE employees, agents or independent contractors is in accordance with the relevant regulations of the Author's institution. If more than one author is listed at the beginning of this Agreement, JoVE may, in its sole discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JoVE reserves the right, in its absolute and sole discretion and without giving any reason therefore, to accept or decline any work submitted to JoVE. JoVE and its employees, agents and independent contractors shall have



1 Alewife Center #200 Cambridge, MA 02140 tel. 617.945.9051 www.jove.com

ARTICLE AND VIDEO LICENSE AGREEMENT

full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JoVE has sole discretion as to the method of making and publishing the Materials, including, without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

11. Indemnification. The Author agrees to indemnify JoVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JoVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JoVE, making of videos by JoVE, or publication in JoVE or elsewhere by JoVE. The Author shall be responsible for, and shall hold JoVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to the making of a video by JoVE its employees, agents or independent contractors. All sterilization, cleanliness or decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's

expense. All indemnifications provided herein shall include JoVE's attorney's fees and costs related to said losses or damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JoVE, its employees, agents or independent contractors.

- 12. Fees. To cover the cost incurred for publication, JoVE must receive payment before production and publication the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JoVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.
- 13. <u>Transfer, Governing Law.</u> This Agreement may be assigned by JoVE and shall inure to the benefits of any of JoVE's successors and assignees. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to me one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement required per submission.

CORRESPONDING AUTHOR:

Name:

Department:

Institution:

Article Title:

A novel micro agar salt bridge electrode for analysing proton turnover rate of recombinant membrane proteins.

Signature:

Date:

Please submit a signed and dated copy of this license by one of the following three methods:

- 1) Upload a scanned copy of the document as a pfd on the JoVE submission site;
- 2) Fax the document to +1.866.381.2236;
- 3) Mail the document to JoVE / Attn: JoVE Editorial / 1 Alewife Center #200 / Cambridge, MA 02139

For questions, please email submissions@jove.com or call +1.617.945.9051

Response to Reviewers

We thank the reviewers for their helpful and constructive comments. Below we have addressed point by point all issues raised by the reviewers. The newly introduced/modified text is underlined in the resubmitted manuscript.

Reviewer #1:

The manuscript by Kreiter et al. deals with the elaboration of the electrophysiological technique for the reconstitution and functional characterization of the reconstituted carrier proteins. The results are sound, the manuscript is well written and quite interesting.

We thank the reviewer for the positive feedback.

1) I have however, several comments which can potentially improve the impact of the work. 1). I did not find the description of the reference electrode preparation in the manuscript. It seems that the authors did not use it. I think that it is necessary to test the effect of the agar bridge presence on the reference electrode on the measurements of reverse potentials. Since the micro agar salt bridge improved substantially the results of the measurements, the presence of the agar bridge on the reference electrode may also be helpful.

We thank the reviewer for addressing this point. Indeed, if substrates would be added directly to the buffer reservoir containing the reference electrode, the usage of the salt bridge would balance substrate-related diffusion potential. However, we used buffer solution which is chloride-free to exclude the putative possibility that uncoupling proteins transport chloride ions. For the adjusting of pH we used Tris or MES. The electrode potential in the absence of substantial concentration of chloride is primarily dependent on chloride impurities in the buffer solution. As its composition is unchanged during the experiments it will simply result in a constant offset potential. However, for a different situation where it is required to obtain an absolute potential difference between the two electrodes, a simple agar salt-bridge system (Ag/AgCl 3 M KCl) could also be used for the reference electrode.

We made a corresponding comment in the revised manuscript (underlined).

2). The manuscript would benefit if the authors describe in more details the background and the meaning of the equations which transform the measurable parameter, a reverse potential, into another parameter, i.e. proton turnover rate. For example, it is hard to understand the

relationship between the proton turnover rate and delta_T which is a time parameter of the triangular alternating voltage applied.

We thank the reviewer to mention this point. The more detailed deviation of the conversion of the measurable parameter, reverse potential, to the turnover number is described in Urbankova et al (2003) JBC and Beck et al (2007) FASEB. We adjusted the symbols to ΔT_{ramp} to indicate the time constant of the voltage ramp in step 3 in order to distinguish it from the time unit of the substrate turnover rate.

3). I suggest to run control measurements with a protonophore (i.e. CCCP), K+-carrier (i.e. nonactin), and H+/K+-channel (gramicidin) which can validate the approach.

We appreciate the suggestion of the reviewer. Our main goal in this study was to develop a salt bridge for the rapid and stable determination of the substrate turnover number and we therefore did not focus on reverse potential measurements and conversions to substrate turnover rates per se, as it was already established in Urbankova et al. 2003. Proposed experiments will substantially confirm the method, but are out of the scope for this work.

Reviewer #2:

This manuscript clearly demonstrates that using such a 'novel micro agar salt bridge electrode configuration' gives vastly superior results compared to the large inaccuracies that resulting from the use of a simple Ag/AgCl electrode as the test electrode for analysing proton turnover rate of membrane proteins reconstituted in the lipid bilayer membrane. The problem with the use of a simple Ag/AgCl electrode results from the fact that its electrical potential in a salt solution is highly dependent on the Cl- activity in the solution (e.g., Eqs. 1 & 2 in Raynauld & Laviolette, 1987 paper) and that in the absence of a significant Cl- concentration, the electrode potentials are poorly defined and unstable. This creates a major problem at the internal membrane interfacial solution in the absence of a significant Cl-, and where the solution composition is being altered. The problem is essentially overcome with the use of this novel 3M KCl micro-agar electrode, as described in this JoVE manuscript.

Thank you for the positive comment.

There are some relatively minor suggestions, comments and questions, which the authors should take note of:

1. With respect to the Ag/AgCl electrode in the large external Cl--free buffer solution, there should also be some small comment indicating that while its electrode potential is poorly defined and primarily dependent on chloride impurities in that solution, as its composition is unchanged during the experiments, it will simply result in a constant offset potential, during an experiment. However, for a different situation where it was required to obtain an absolute potential difference between the two electrodes, a simpler Ag/AgCl 3 M KCl agar salt-bridge system could also be used for the reference electrode.

We appreciate the advice of the reviewer and introduced the following comment to the text of the article.

Since uncoupling proteins were suggested to transport chloride ions, chlorid- free buffer solution was used and the pH was adjusted using Tris or MES. The electrode potential in the absence of substantial concentration of chloride is primarily dependent on chloride impurities in the buffer solution. As its composition is unchanged during the experiments it will simply result in a constant offset potential. However, for a different situation where it is required to obtain an absolute potential difference between the two electrodes, a simple agar salt-bridge system (Ag/AgCl 3 M KCl) could also be used for the reference electrode.

2. With respect to the micro agar salt bridge on Line 77, I was surprised that the agarose concentration used was only 1%, as the more standard concentration in most salt-bridges 3-4% (w/v) agarose is used (Barry et al., 2013, Eur Biophys J 42-631-646; cp Kleene, 1993, who even used 5% agarose for his salt bridges) and I wondered if the salt-bridge in this JoVE manuscript is sufficiently gelled to retain the agar in the salt bridge for a reasonably long time.

We thank the reviewer for addressing this point. Indeed, the concentration of agarose in the salt bridge is lower than compared to the before mentioned references. Since our measurements are relatively short (less than 5 minutes) and we prepare new salt bridges for each new measurement, the usage of 1 % (w/v) agarose is enough to gel the salt bridge in this short period. Generally we suggest using agar in a concentration range of 1 - 5% (w/v), depending on the time scale of the

measurements and how often the salt bridge is reused. We also introduced the recommended citations.

3. I found the statement that the 'salt ions K+ and Cl- have equal mobility in liquid' (Lines 316f) was strictly inaccurate and misleading with respect to the role of 3 M KCl at such junctions. The relative mobilities of the Cl- and K+ ions uCl/uK=1.04, and 3 M KCl microelectrodes or saltbridges are useful for minimizing liquid junction potentials (LJPs), because of the high concentration of the 3 M KCl and because the mobilities of the Cl- and Cl- ions are similar Cl- Cl-

We apologize for the inaccuracy regarding the mobilities of the ion species. The values are of course not equal but similar. We have corrected the phrases in the manuscript (underlined) to avoid incorrect statements.

4. Fig. 1a. With respect to the precise construction of the Microload pipette, it is difficult to see from Fig. 1A, the precise relationship between the chlorided part of the Ag wire and the part of it in contact with the 3 M KCl agar-Ag/AgCl region (compare Fig. 1 in Shao & Feldman). It could also help to replace 'Agar - Salt' by 'Agar - 3 M KCl' in Fig. 1A).

We have made changes to Fig. 1 A so that the contact side of the chloride part of the Ag wire and the 3 M KCl agarose solution is outlined more clearly.

5. Section 3.1 (Lines 191-195). It would be more accurate to define ΔT as "...and the duration of each ramp as ΔT , calculate..." especially, as 'time constant' tends to have a well-defined usage in electrophysiology and electricity. I also take it that Umax = 100 mV?

We have defined the parameters of the voltage ramp - the maximum voltage U_{max} (50 mV) and the time of the ramp ΔT_{ramp} - more clearly.