

Journal of Visualized Experiments

Large-scale Top-down Proteomics Using Capillary Zone Electrophoresis Tandem Mass Spectrometry --Manuscript Draft--

Article Type:	Invited Methods Article - JoVE Produced Video
Manuscript Number:	JoVE58644R1
Full Title:	Large-scale Top-down Proteomics Using Capillary Zone Electrophoresis Tandem Mass Spectrometry
Keywords:	top-down proteomics, capillary zone electrophoresis, electrospray ionization, tandem mass spectrometry, proteoform, Escherichia coli
Corresponding Author:	Liangliang Sun Michigan State University East Lansing, UNITED STATES
Corresponding Author's Institution:	Michigan State University
Corresponding Author E-Mail:	lsun@chemistry.msu.edu
Order of Authors:	Liangliang Sun Elijah N. McCool Rachele Lubeckyj Xiaojing Shen Qiang Kou Xiaowen Liu
Additional Information:	
Question	Response
Please indicate whether this article will be Standard Access or Open Access.	Standard Access (US\$2,400)
Please indicate the city, state/province, and country where this article will be filmed . Please do not use abbreviations.	East Lansing, Michigan, United States

MICHIGAN STATE UNIVERSITY

Cover Letter

Dear Benjamin,

On behalf of myself and all the co-authors, Elijah N. McCool, Rachele A. Lubeckyj, Xiaojing Shen, Qiang Kou, Xiaowen Liu, I am submitting the revised manuscript titled "Large-scale Top-down Proteomics Using Capillary Zone Electrophoresis Tandem Mass Spectrometry" to *JoVE*.

In the revised manuscript, we made changes based on the editor's and reviewers' comments. Since the changes are so thorough, we did not highlight the changes. We also provided a "responses to the comments" file containing point-by-point response to the comments.



We greatly appreciate your kind consideration.

Sincerely,

Liangliang Sun

**DEPARTMENT OF
CHEMISTRY**

College of Natural Science

Michigan State University
578 S. Shaw Lane • Room 485
East Lansing, MI
48824-1322
517-353-1086
FAX: 517-353-1793

Assistant Professor at Department of Chemistry, Michigan State University
578 S Shaw Ln, East Lansing, MI 48824
Email: lsun@chemistry.msu.edu Phone: 517-353-0498

TITLE:**Large-scale Top-down Proteomics Using Capillary Zone Electrophoresis Tandem Mass Spectrometry****AUTHORS & AFFILIATIONS:**

Elijah N. McCool¹ *, Rachele Lubecky¹ *, Xiaojing Shen¹, Qiang Kou², Xiaowen Liu^{2,3}, Liangliang Sun¹

¹Department of Chemistry, Michigan State University, East Lansing, MI, United States of America

²Department of BioHealth Informatics, Indiana University-Purdue University Indianapolis, Indianapolis, IN, United States of America

³Center for Computational Biology and Bioinformatics, Indiana University School of Medicine, Indianapolis, IN, United States of America

* These authors equally contributed.

Corresponding Author:

Liangliang Sun (lsun@chemistry.msu.edu)

E-mail Addresses of the Co-authors:

Elijah N. McCool (mccool@chemistry.msu.edu)

Rachele Lubecky (lubecky1@chemistry.msu.edu)

Xiaojing Shen (shenxia7@chemistry.msu.edu)

Qiang Kou (qkou@uimail.iu.edu)

Xiaowen Liu (xwliu@iupui.edu)

KEYWORDS:

Top-down proteomics, capillary zone electrophoresis, electrospray ionization, tandem mass spectrometry, proteoform, *Escherichia coli*

SUMMARY:

A detailed protocol is described for the separation, identification, and characterization of proteoforms in protein samples using capillary zone electrophoresis-electrospray ionization-tandem mass spectrometry (CZE-ESI-MS/MS). The protocol can be used for the high-resolution characterization of proteoforms in simple protein samples and the large-scale identification of proteoforms in complex proteome samples.

ABSTRACT:

Capillary zone electrophoresis-electrospray ionization-tandem mass spectrometry (CZE-ESI-MS/MS) has been recognized as a useful tool for top-down proteomics that aims to characterize proteoforms in complex proteomes. However, the application of CZE-MS/MS for large-scale top-down proteomics has been impeded by the low sample-loading capacity and narrow separation window of CZE. Here, a protocol is described using CZE-MS/MS with a microliter-scale sample-loading volume and a 90-min separation window for large-scale top-down proteomics. The CZE-

MS/MS platform is based on a linear polyacrylamide (LPA)-coated separation capillary with extremely low electroosmotic flow, a dynamic pH-junction-based online sample concentration method with a high efficiency for protein stacking, an electro-kinetically pumped sheath flow CE-MS interface with extremely high sensitivity, and an ion trap mass spectrometer with high mass resolution and scan speed. The platform can be used for the high-resolution characterization of simple intact protein samples and the large-scale characterization of proteoforms in various complex proteomes. As an example, a highly efficient separation of a standard protein mixture and a highly sensitive detection of many impurities using the platform is demonstrated. As another example, this platform can produce over 500 proteoform and 190 protein identifications from an *Escherichia coli* proteome in a single CZE-MS/MS run.

INTRODUCTION:

Top-down proteomics (TDP) aims for the large-scale characterization of proteoforms within a proteome. TDP relies on the effective liquid-phase separation of intact proteins before electrospray ionization-tandem mass spectrometry (ESI-MS/MS) analysis due to the high complexity and large concentration dynamic range of the proteome¹⁻⁵. Capillary zone electrophoresis (CZE) is a powerful technique for the separation of biomolecules based on their size-to-charge ratios⁶. CZE is relatively simple, requiring only an open tubular-fused silica capillary, a background electrolyte (BGE), and a power supply. A sample of intact proteins can be loaded into the capillary using pressure or voltage, and separation is initiated by immersing both ends of the capillary in the BGE and applying a high voltage. CZE can approach ultra-high separation efficiency (> one million theoretical plates) for the separation of biomolecules⁷. CZE-MS has a drastically higher sensitivity than widely used reversed-phase liquid chromatography (RPLC)-MS for the analysis of intact proteins⁸. Although CZE-MS has a great potential for large-scale top-down proteomics, its wide application in proteomics has been impeded by several issues, including a low sample-loading capacity and narrow separation window. The typical sample loading volume in CZE is about 1% of the total capillary volume, which usually corresponds to less than 100 nL⁹⁻¹¹. The separation window of CZE is usually less than 30 min due to the strong electroosmotic flow (EOF)^{9,10}. These issues limit the CZE-MS/MS for the identification of a large number of proteoforms and low abundant proteoforms from a complex proteome.

Much effort has been made to improve the sample loading volume of CZE *via* online sample concentration methods (*e.g.*, solid-phase microextraction [SPME]^{12,13}, field-enhanced sample stacking [FESS]^{9,11,14}, and dynamic pH junction¹⁵⁻¹⁸). FESS and dynamic pH junction are simpler than SPME, only requiring a significant difference between the sample buffer and the BGE in conductivity and pH. FESS employs a sample buffer with much lower conductivity than the BGE, leading to a stacking of analytes on the boundary between the sample zone and the BGE zone in the capillary. Dynamic pH junction utilizes a basic sample plug (*e.g.*, 50 mM ammonium bicarbonate, pH 8) and an acidic BGE (*e.g.*, 5% [v/v] acetic acid, pH 2.4) on both sides of the sample plug. Upon application of a high positive voltage at the injection end of the capillary, titration of the basic sample plug occurs, focusing the analytes into a tight plug before undergoing a CZE separation. Recently, the Sun group systematically compared FESS and dynamic pH junction for the online stacking of intact proteins, demonstrating that dynamic pH junction could produce

much better performance than FESS for the online concentration of intact proteins when the sample injection volume was 25% of the total capillary volume¹⁹.

Neutrally coated separation capillaries (*e.g.*, linear polyacrylamide [LPA]) have been employed to reduce the EOF in the capillary, slowing down the CZE separation and widening the separation window^{20,21}. Recently, the Dovichi group developed a simple procedure for the preparation of stable LPA coating on the inner wall of capillaries, utilizing ammonium persulfate (APS) as the initiator and temperature (50 °C) for free radical production and polymerization²². Very recently, the Sun group employed the LPA-coated separation capillary and the dynamic pH junction method for the CZE separation of intact proteins, reaching a microliter-scale sample loading volume and a 90-min separation window¹⁹. This CZE system opens the door to using CZE-MS/MS for large-scale top-down proteomics.

CZE-MS requires a highly robust and sensitive interface to couple CZE to MS. Three CE-MS interfaces have been well developed and commercialized in the history of CE-MS, and they are the co-axial sheath-flow interface²³, the sheathless interface using a porous tip as the ESI emitter²⁴, and the electro-kinetically pumped sheath flow interface^{25,26}. The electro-kinetically pumped sheath-flow-interface-based CZE-MS/MS has reached a low zeptomole peptide detection limit⁹, over 10,000 peptide identifications (IDs) from the HeLa cell proteome in a single run¹⁴, a fast characterization of intact proteins¹¹, and highly stable and reproducible analyses of biomolecules²⁶. Recently, the LPA-coated separation capillary, the dynamic pH junction method, and the electro-kinetically pumped sheath flow interface were used for large-scale top-down proteomics of an *Escherichia coli* (*E. coli*) proteome^{19,27}. The CZE-MS/MS platform approached over 500 proteoform IDs in a single run¹⁹ and nearly 6,000 proteoform IDs *via* coupling with size-exclusion chromatography (S)-RPLC fractionation²⁷. The results clearly show the capability of CZE-MS/MS for large-scale top-down proteomics.

Herein, a detailed procedure of using CZE-MS/MS for large-scale top-down proteomics is described. The CZE-MS/MS system employs the LPA-coated capillary to reduce the EOF in the capillary, the dynamic pH junction method for the online concentration of proteins, the electro-kinetically pumped sheath flow interface for coupling CZE to MS, an orbitrap mass spectrometer for the collection of MS and MS/MS spectra of proteins, and a TopPIC (TOP-Down Mass Spectrometry-Based Proteoform Identification and Characterization) software for proteoform ID *via* database search.

PROTOCOL:

1. Preparation of LPA Coating on the Inner Wall of the Separation Capillary

1.1. Pretreatment of the capillary

1.1.1. Flush a fused silica capillary (120 cm in length, 50 µm in inner diameter [i.d.], 360 µm in outer diameter [o.d.]) successively with 500 µL of 1 M sodium hydroxide, deionized water, 1 M hydrochloric acid, deionized water, and LC-MS grade methanol using a syringe pump.

1.1.2. Dry the capillary with nitrogen gas (10 psi, ≥ 12 h) and fill the capillary with 50% (v/v) 3-(trimethoxysilyl)propyl methacrylate in methanol using a syringe pump. Seal both ends of the capillary with silica rubber and incubate it at room temperature for at least 24 h.

Note: During this step, the inner wall of the capillary is functionalized with C=C. Longer incubation results in better capillary coating due to a more complete reaction.

1.1.3. Cut a small portion (~ 5 mm) of the capillary from both ends with a cleaving stone. Rinse the capillary with methanol (500 μ L) using a syringe pump to clean up the unreacted reagents. Dry the capillary with nitrogen gas (10 psi, ≥ 12 h).

1.2. Preparation of the LPA coating

Note: This procedure is based on Xhu *et al.*²² with minor modifications.

1.2.1. Prepare an acrylamide solution (40 mg of acrylamide in 1 mL of water) and ammonium persulfate (APS) solution (5% [w/v] in water).

1.2.2. Add 2 - 3 μ L of the APS solution to 500 μ L of the acrylamide solution, vortex the mixture, and degas it with nitrogen gas for 5 min to remove the oxygen in the solution.

1.2.3. Load the mixture into the pretreated capillary using a vacuum, seal both ends of the capillary with silica rubber, and incubate it in a water bath at 50 $^{\circ}$ C for 40 min.

1.2.4. Remove a small portion (~ 5 mm) of the capillary from both ends with a cleaving stone. Push the unreacted solution out of the capillary with water (200 μ L), using the syringe pump.

Note: Ensure the polymer pushed out of the capillary is an agarose gel-like consistency. A longer incubation step (up to ~ 45 - 50 min) can result in a better capillary coating. With longer reaction periods, the capillary may become blocked and high pressure is required to push out the polymer with water. An HPLC pump can be used for this purpose.

2. Etching of the Capillary with Hydrofluoric Acid

CAUTION: Use appropriate safety procedures while handling hydrofluoric acid (HF) solutions. All the HF-related operations need to be done in a chemical hood. Before any HF-related operation, make sure that 2.5% calcium gluconate gel is available for use in the case of exposure. Double gloves are required, a typical nitrile glove inside and a heavy neoprene glove outside. Wear a lab coat and chemical safety goggles. After the HF operations, keep liquid and solid hazardous waste separate. The liquid HF waste must be neutralized immediately with a high-concentration sodium hydroxide solution for temporary storage before waste pick-up. The solid HF waste needs to be temporarily stored in a plastic container that is lined with two thick plastic one-gallon Ziploc bags and a lid. Both the solid and liquid waste must be labeled properly.

2.1. Burn a portion of the capillary using a gentle flame (*e.g.*, pocket lighter) to remove the polyimide outer-coating (1 cm in length) around 4 cm away from one end of the capillary. Gently clean the burnt portion of the capillary by wiping to remove the polyimide coating completely.

Note: Proceed with caution, as the portion of the capillary without the polyimide coating will be a little fragile.

2.2. Drill a small hole at the end of a 200- μ L tube around the same size as the capillary outer diameter to sufficiently hold the capillary in place once it is threaded through this hole. Thread the end of the capillary that is close to the burnt portion through the hole until the burnt portion is in the tube.

Note: It is recommended to test the size of the hole with the end of the capillary away from the burnt portion to ensure the correct size, as the burnt portion of the capillary is fragile.

2.3. Add \sim 150 μ L of HF (48% - 51% solution in water) to the 200- μ L tube so that the HF solution is about halfway up the burnt portion on the capillary. Incubate the capillary in the HF solution at room temperature for 90 - 100 min.

Note: It is recommended that another hole is created in the lid of the tube and some small object (*e.g.*, a small pipette tip) is used to hold up the tube while the incubation is taking place. The pipette tip can be used to puncture foam or some other solid platform from which the reaction chamber can hang from, creating a safe environment. Place paper towels below the tube containing HF and follow proper procedures in case of a spill.

2.4. Remove the capillary from the tube and wash the exterior with deionized water to remove any residual HF.

Note: Ensure that the outer diameter of the capillary at the narrowest portion of the burnt part is now smaller than 100 μ m, as it should be.

2.5. Cut the etched part of the capillary at the middle of the narrowest portion using a cleaving stone to produce a separation capillary with less than 100 μ m in o.d. at one end. Cut the non-etched end of the capillary to get a 1-m separation capillary.

Note: Dispose of HF waste according to the institutionally prescribed protocol.

3. Preparation of the Samples

3.1. Preparation of a standard protein mixture

3.1.1. Prepare a standard protein mixture containing cytochrome c (Cyto.c, 12 kDa, 0.1 mg/mL), lysozyme (14.3 kDa, 0.1 mg/mL), β -casein (24 kDa, 0.4 mg/mL), myoglobin (16.9 kDa, 0.1 mg/mL),

carbonic anhydrase (CA, 29 kDa, 0.5 mg/mL), and bovine serum albumin (BSA, 66.5 kDa, 1.0 mg/mL) in LC-MS grade water as a stock solution.

Note: The stock solution can be aliquoted and stored at -80 °C for use.

3.1.2. Dilute the stock solution by a factor of 10 with a 50 mM ammonium bicarbonate solution (pH 8.0) in LC-MS grade water for CZE-MS analysis.

Note: Filter the ammonium bicarbonate solution before use with a membrane filter (*e.g.*, 0.2 µm of cellulose nitrate membrane and 50 mm in diameter).

3.2. Preparation of an *E. coli* sample

3.2.1. Culture *E. coli* (K-12 MG1655) cells in LB medium at 37 °C while shaking at 225 rpm until the OD600 value approaches 0.7. Collect *E. coli* cells *via* centrifugation (3283 x *g*, 10 min) and wash them multiple times with phosphate-buffered saline (PBS) to remove the medium.

3.2.2. Suspend the *E. coli* cells in the lysis buffer containing 8 M urea, protease inhibitors, and 100 mM ammonium bicarbonate (pH 8.0) and ultrasonicate on ice for 15 min using an ultrasonic cell disruptor with a cup horn for complete cell lysis.

3.2.2.1. Set the **Duty Cycle (%)** to 50 and set the **Output Control** to 7 for the ultrasonic cell disruptor.

3.2.3. Centrifuge the lysate at 18,000 x *g* for 10 min. Collect the supernatant and discard the pellet. Use a small aliquot of the supernatant for protein concentration measurement with the bicinchoninic acid (BCA) assay.

Note: The lysate needs to be diluted at least by a factor of three to reduce the urea concentration lower than 3 M in order to become compatible with the BCA assay. The BCA assay is performed based on the manufacturer's procedure.

3.2.4. Mix 1 mg of *E. coli* proteins with cold acetone with a volume ratio of 1:4 and keep the mixture at -20 °C overnight to precipitate the proteins.

Note: Perform all experiments using acetone in a chemical hood.

3.2.5. Spin down precipitated proteins (12,000 x *g*, 5 min) and remove the supernatant. Wash the pellet with cold acetone again (same volume as before) and spin down again.

3.2.6. Remove the supernatant and allow the pellet to dry in the chemical hood for a couple of minutes.

Note: Do not overdry the protein pellet.

3.2.7. Dissolve the precipitated *E. coli* proteins (1 mg) in 200 μ L of 8 M urea in 100 mM ammonium bicarbonate (pH 8.0).

3.2.8. Denature, reduce, and alkylate the sample.

3.2.8.1. Incubate the sample at 37 °C for 30 min to denature the proteins.

3.2.8.2. Reduce with dithiothreitol (DTT) by adding 2 μ L of 1 M DTT solution (in 100 mM ammonium bicarbonate) into the sample and incubating the sample at 37 °C for 30 min.

3.2.8.3. Alkylate the proteins with iodoacetamide (IAA) by adding 6 μ L of 1 M IAA solution (in 100 mM ammonium bicarbonate) to the sample and incubate the sample for 20 min at room temperature in the dark.

3.2.8.4. Quench the excess IAA with DTT by adding 2 μ L of 1 M DTT solution and incubate for 5 min at room temperature. Acidify the sample with formic acid (FA) to get a final FA concentration of 1% (v/v).

Note: Reduction and alkylation are performed to aid the unfolding of proteins for downstream fragmentation. Those steps will lose the information of disulfide bonds. If the goal is to study the disulfide bonds on the proteins, avoid those steps. Remember to handle all concentrated acid solutions in a chemical hood.

3.2.9. Desalt the sample using a C4 trap column (*e.g.*, 4 mm of i.d., 10 mm in length, packed with 3- μ m particles having 300-Å pores) with an HPLC system. Use a UV detector at a wavelength of 254 nm for detection.

3.2.9.1. Set the mobile phase flow rate to 1 mL/min. Activate the column by flushing it with 80% (v/v) acetonitrile (ACN), 0.1% FA in water for 10 min, and equilibrate by flushing it with 2% (v/v) ACN, 0.1% FA in water for 10 min.

3.2.9.2. Load 500- μ g proteins onto the column and desalt by flushing them with 2% (v/v) ACN, 0.1% FA in water for 10 min at a 1 mL/min flow rate.

3.2.9.3. Elute the proteins with 80% (v/v) ACN, 0.1% FA in water for 3 min at a 1 mL/min flow rate. Collect the eluate and lyophilize it with a vacuum concentrator.

Note: The C4 trap column can be used to desalt large amounts of proteins but not low micrograms of proteins due to possible sample loss. For limited protein materials, consider using pipette tips with C4 media for protein desalting. Be careful not to overdry the proteins when lyophilizing the sample.

3.2.10. Dissolve the 500 µg of proteins (assuming no sample loss during the sample preparation) in 250 µL of 50 mM ammonium bicarbonate (pH 8.0) for CZE-MS experiments.

4. Set-up of the CZE-MS/MS System and Analysis of the Samples

4.1. Set-up of the CZE system

4.1.1. Prepare a BGE containing 5% or 10% (v/v) acetic acid (pH ~2.4 or ~2.2) in LC-MS grade water. Put ~1.5 mL of the BGE into a BGE vial that matches with the buffer tray of the CZE autosampler.

Note: Make fresh BGE every week and change the BGE solution in the vial every day to avoid any contamination.

4.1.2. Add a sample solution (5 - 200 µL) to an insert tube and put the insert tube into a sample vial that matches with the sample tray of the CZE autosampler.

4.1.3. Load the sample, the BGE, and the separation capillary into the CZE autosampler.

Note: The language used in this protocol most accurately reflects the use of one particular autosampler (see **Table of Materials**), but the principles can be applied to other autosamplers.

4.1.3.1. Select **Sample load** in the manual control page of the CZE autosampler. Load the BGE vial into the buffer tray of the autosampler, noting its position in the buffer tray. Load the sample vial into the sample tray of the autosampler, noting its position in the sample tray.

Note: The instrument can be operated in manual control by clicking on the picture of the instrument for **Device Monitor** on the main instrument page and then selecting **Manual** under **Direct Control**.

4.1.3.2. Load the separation capillary into the CZE autosampler, using the non-etched end of the capillary. Put the autosampler to the **Alignment** position using manual control. Thread the non-etched end of the capillary into a hole that is used to hold the separation capillary until it cannot be physically threaded further.

Note: In this case, the end of the capillary is almost at the same height as the end of the electrode, and at least 50 µL of the sample is required in the sample vial to make sure the end of the capillary is immersed in the sample during injection. If the sample volume is low (*i.e.*, 5 µL), the height of the injection end of the capillary needs to be adjusted as described in step 4.1.3.2.1.

4.1.3.2.1. Adjust the height of the injection end of the capillary if the sample volume is lower than 50 µL (*i.e.*, 5 µL). Switch the autosampler to **Standby** using manual control. Type the sample position in the system and move the capillary to the sample. Push the injection end of the capillary further down to reach the bottom of the sample vial.

Note: In this case, the sample injection can be performed with a sample of only 5 μ L in the vial.

4.1.3.3. Continuing to use the manual control, flush the capillary with the BGE for 20 min, using a pressure of 20 psi.

4.2. Set-up of the CZE-MS interface

4.2.1. Pull a borosilicate glass capillary (1 mm in o.d., 0.75 mm in i.d., 10 cm long) into two electrospray emitters with an orifice of 20 - 40 μ m in o.d.

Note: Keep the length of the electrospray emitter as 4 - 5 cm and cut it with a cleaving stone if necessary. Dispose of all glass waste in a sharps container. The settings to get 20- to 40- μ m tips are as follows. Set the heat to 498, the pull to 5, the velocity to 10, the delay to 150, and the pressure to 200. Check the size of the emitters with a microscope before use. Adjust the parameters slightly if necessary.

4.2.2. Mount a commercial electro-kinetically pumped sheath flow interface (see **Table of Materials**) to the front of the mass spectrometer. Fill the sheath buffer reservoir with a buffer containing 10% (v/v) methanol and 0.2% (v/v) FA in LC-MS grade water.

4.2.3. Flush the **T** in the interface with the sheath buffer *via* applying pressure manually with a syringe. Thread the 1-mm-o.d. end of an electrospray emitter through a sleeve tubing and connect the emitter with one port of the **T** *via* a fitting. Simply flush the **T** manually again with a syringe to fill the emitter with the sheath buffer.

4.2.4. Adjust the distance between the orifice of the emitter and the entrance of the mass spectrometer to \sim 2 mm with the help of the camera that came with the interface. Apply a 2- to 2.2-kV voltage at the sheath buffer vial for electrospray. Adjust the spray voltage to reach a stable electrospray.

Note: If no electrospray or an unstable electrospray is observed, check the interface to make sure there are no large bubbles in the emitter, in the **T**, and in the tubing that used to connect the sheath buffer vial and the **T**. The distance between the orifice of the emitter and the mass spectrometer entrance can be roughly estimated by comparing it with the 1-mm o.d. of the emitter. The spray voltage depends on the size of the emitter. For 20- to 40- μ m emitters, 2 - 2.2 kV is good enough. Remember to always be careful when using high voltages.

4.2.5. Turn off the spray voltage and gently thread the etched end of the separation capillary through the **T** into the emitter until it cannot be pushed further. Apply a low pressure (*i.e.*, 5 psi) at the injection end of the capillary during this process to make sure there are no bubbles in the capillary.

Note: The capillary is connected to the **T** *via* a sleeve tubing and a fitting. Adjust the position of the etched end of the capillary in the emitter with the help of the camera to within 500 µm. The distance can be roughly estimated by comparing it with the 1-mm o.d. of the emitter.

4.2.6. Stop the low pressure at the injection end of the capillary. Flush the emitter a little bit with the sheath buffer. Again, apply 2 - 2.2 kV of spray voltage to test the spray.

4.3. Set-up of a CZE method for sample injection, separation, capillary flushing, and triggering of the mass spectrometer for data acquisition

4.3.1. Select **New method** under the **File** drop-down menu on the main instrument screen to initiate a new method.

4.3.2. Select **Inlet** as SV/EV, **Tray** as Sample, **Pressure** as 5 psi, **KV** as 0, and **Duration** as 95 s for a sample injection.

Note: The sample is injected into the capillary *via* applying pressure. The sample injection volume can be calculated based on the pressure and injection time using Poiseuille's Law. For example, 5 psi for a 95-s sample injection corresponds to about 500 nL of sample-loading volume for a 1-m-long separation capillary (50-µm i.d.).

4.3.3. Select parameters for the CZE separation and set up parameters for triggering the data acquisition.

4.3.3.1. Set **Inlet** as InV, **Tray** as Buffer, **Pressure** as 0 psi, **KV** as 30, and **Duration** as 4,200 s for the separation of the standard protein mixture sample. Set **Inlet** as InV, **Tray** as Buffer, **Pressure** as 0 psi, **KV** as 20, and **Duration** as 6,600 s for the separation of the *E. coli* proteome sample.

Note: The voltage applied for the separation can be adjusted based on the sample complexity. For simple protein samples, 30 kV can be applied to speed up the analysis. For complex samples, 20 kV is applied to slow down the separation and acquire more MS/MS spectra for proteoform IDs.

4.3.3.2. Set up the parameters for triggering the data acquisition under **Timed Events**. Set **Time (s)** as 0.0 and **Type** as Relay 1 for **Activate** in step 1; Set **Time (s)** as 1.0 and **Type** as Relay 1 for **Deactivate** in step 2.

4.3.4. Select **Inlet** as InV, **Tray** as Buffer, **Pressure** as 10 psi, **KV** as 30, and **Duration** as 600 s for capillary flushing.

4.3.5. Save the method files for the CZE-MS and MS/MS experiments.

4.4. Set-up of the MS and MS/MS

4.4.1. Set up the MS and MS/MS parameters for intact protein analysis using a quadrupole-ion trap mass spectrometer (see **Table of Materials**).

Note: The positive ion mode and higher energy collisional dissociation (HCD) for fragmentation are employed.

4.4.1.1. Adjust the tune file settings: turn on **Intact protein mode** and use a trapping pressure of 0.2. Set the ion transfer capillary temperature to 320 °C and the s-lens RF level to 55.

4.4.1.2. Build the MS and MS/MS method file.

4.4.1.2.1. For full MS, set the number of microscans to 3, the resolution to 240,000 (at m/z 200), the automatic gain control (AGC) target value to 1E6, the maximum injection time to 50 ms, and the scan range to 600 - 2000 m/z .

4.4.1.3. For MS/MS, use data-dependent acquisition (DDA) for the eight most intense ions in a full MS spectrum. Set the isolation window to 4 m/z and the normalized collisional energy (NCE) to 20%. Set the number of microscans to 1, the resolution to 120,000 (at m/z 200), the AGC target value to 1E5, and the maximum injection time to 200 ms.

4.4.1.4. Set the intensity threshold for triggering fragmentation to 1E5 and set the ions with a charge state higher than 5 to be isolated for fragmentation. Turn on the **exclude isotopes** and set the dynamic exclusion to 30 s.

Note: The number of microscans for MS/MS can be increased to 3. In this case, a Top3 DDA method can be used in order to reduce the cycle time.

4.4.2. Set up a wired connection between the CE autosampler and the mass spectrometer for the automatic triggering of the data acquisition. Connect one end of a wire to the **Relay 1 Contact A** and **Relay 1 Contact B** on the back of the CE autosampler. Connect the other end of the wire to **Start In -** and **Start In +** on the side of the mass spectrometer.

4.5. CZE-MS/MS experiment and analysis of the samples

4.5.1. Apply 30 kV using the CE manual control at the sample injection end and 2 - 2.2 kV at the sheath buffer vial using the external power supply to test the separation capillary and electrospray.

Note: The separation current, in this case, is around 8 - 9 μ A if 5% (v/v) acetic acid is used as the BGE.

4.5.2. Set up a data acquisition sequence on the mass spectrometer computer by selecting a new sequence.

4.5.2.1. Select **Unknown** as the sample type and specify a file name and path. Indicate the MS method to be used for each sample run under **Instrument Method**.

Note: Since there is a separate computer for controlling the CE autosampler, the sample position does not need to be specified here.

4.5.3. Set up a CZE separation sequence on the CE computer to indicate the buffer position, sample position, and the CZE method. To start a new sequence, select **Sequence** under the **Analysis** drop-down menu on the main instrument page.

4.5.4. Start the data acquisition method on the mass spectrometer computer first by selecting **Run sequence** (or **Run sample** for single runs). Then, start the CE sequence on the CE autosampler computer.

Note: When the separation voltage is on after the sample injection, a signal is sent from the CE autosampler to the mass spectrometer for triggering the data acquisition. The separation current will decrease at the beginning during the separation due to the dynamic pH junction sample stacking. Then, the current recovers gradually.

5. Database Search of the Collected Raw Files with the TopPIC Software

5.1. Transfer the .raw files, including the MS and MS/MS data, from the mass spectrometer computer to a computer dedicated to database searching.

Note: These .raw files may be large and require a powerful computer in order to reach a fast database search.

5.2. Download ProteoWizard and TopPIC suite onto the computer dedicated to database searching.

Note: ProteoWizard and TopPIC suite are open-source software and can be found online (ProteoWizard: <http://proteowizard.sourceforge.net>; TopPIC suite: <http://proteomics.informatics.iupui.edu/software/toppic/>). UniProt databases are used for database searches and can be found on the UniProt website.

5.3. Convert the .raw files to .mzML files with msconvert, an MS file format conversion tool in ProteoWizard²⁸. In the GUI of msconvert (MSConvertGUI.exe), click the button **Browse** and select the .raw file to be converted. Select mzML as the output format and keep all other options at their default values. Click the button **Start** at the bottom right of the GUI.

5.4. Convert the .mzML files to .msalign files with the TopFD tool in TopPIC suite²⁹.

5.4.1. In the GUI of TopFD (topfd_gui.exe), click the button **File** and select the .mzML file created in the previous step. Keep the default values of the parameters and, then, click the button **Start** at the bottom right of the GUI.

Note: TopFD converts precursor and fragment isotopic clusters to monoisotopic masses and identifies possible proteoform features in MS1 data by combining precursor isotopic clusters with similar monoisotopic masses and close migration times. Three files will be generated from TopFD, including an ms1.msalign file, an ms2.msalign file, and a .feature file. The ms1.msalign and ms2.msalign files contain deconvoluted monoisotopic masses of MS1 and MS/MS spectra, respectively. The .feature file contains possible proteoform features.

5.5. Perform a database search with TopPIC (version 1.1.3) in TopPIC suite³⁰.

5.5.1. In the TopPIC GUI (toppic_gui.exe), click on **Database file** and select an appropriate database for searching.

5.5.2. Click on **Spectrum file** and select the ms2.msalign file generated in step 5.4.1 as the input. Select the **MS1 feature file** and choose the .feature file generated in step 5.4.1.

5.5.3. Set cysteine carbamidomethylation (C57) as a fixed modification due to the IAA treatment.

Note: A text file can also be created with various fixed modifications by a text editor. Then, the text file can be selected by using the **File** button next to **Fixed modifications** in the TopPIC GUI.

5.5.4. Select the **Decoy database** feature and, under cutoff settings, select **FDR** (false discovery rate) in the drop-down menu next to **Spectrum level**. Set the FDR to 0.01 at the spectrum level.

Note: TopPIC provides multiple options for filtering the results: spectrum-level FDR, proteoform-level FDR, or both. The FDR is evaluated based on the target-decoy approach³¹.

5.5.5. Leave **Generating function** unselected and set the **Error tolerance (ppm)** to 15.

5.5.6. Under **Advanced Parameters**, select 2 for the maximum number of mass shifts in the drop-down menu. Set the **Maximum mass shift (Da)** of unknown modifications to 500 Da. Leave all other parameters at their default values and click the button **Start** at the bottom right of the GUI.

Note: The TopPIC software generates two resulting text files. The file with the suffix .OUTPUT_TABLE contains the list of identified proteoform-spectrum matches (PrSMs), and the one with a suffix .FORM_OUTPUT_TABLE contains the list of identified proteoforms. For each PrSM, the software provides general information on the match, the observed fragmentation pattern, matched fragment ions, and detected mass shifts.

REPRESENTATIVE RESULTS:

Figure 1 shows a diagram of the dynamic pH-junction-based CZE-ESI-MS system used in the experiment. A long plug of the sample in a basic buffer is injected into an LPA-coated separation capillary filled with an acidic BGE. After applying high voltages I and II, the analytes in the sample zone will be concentrated *via* the dynamic pH junction method. To evaluate the performance of the CZE-MS system, a standard protein mixture (cytochrome c, lysozyme, β -casein, myoglobin, CA, and BSA) is typically analyzed. The representative electropherogram for the standard protein mixture is shown in **Figure 2A**. The standard protein mixture is typically run at least in duplicate to evaluate the separation efficiency and the reproducibility of the system. The separation efficiency can be evaluated with the number of theoretical plates of some proteins, as shown in **Figure 2B**. The reproducibility can be evaluated by the relative standard deviations of protein intensity and migration time. **Figure 3A** shows a zoomed-in view of an electropherogram of the *E. coli* protein sample analyzed by the dynamic pH-junction-based CZE-MS/MS. The normalized-level (NL) protein intensity should be on the scale of 10^8 if 1 μ g of *E. coli* proteins are loaded for the analysis with a quadrupole ion trap mass spectrometer. A zoomed-in view of the electropherogram can be used to assess the separation window of the system. In this case, the separation window is 80 - 90 min. **Figure 3B** shows an example PrSM, including the general corresponding proteoform information, protein sequence, observed fragmentation pattern, and modifications. The very low E-Value ($2.11E-48$) and Spectral FDR (0) suggest the high confidence of the proteoform ID. The high number of matched fragment ions (60) further indicates the high confidence of the ID. The observed fragmentation pattern shows that the fragmentation of the proteoform is highly efficient covering the termini and middle part of the proteoform. The N-terminal cleavage of three amino acids (MTM) is determined through the database search.

FIGURE LEGENDS:

Figure 1: Diagram of the dynamic pH-junction-based CZE-ESI-MS system. The sample is dissolved in 50 mM ammonium bicarbonate, pH 8. The BGE is 5% or 10% (v/v) acetic acid, pH ~ 2.4 or ~ 2.2 . A 1-m LPA-coated capillary is used for separation. An electro-kinetically pumped sheath flow interface is used to couple CZE to MS. A quadrupole ion trap mass spectrometer is used. The high voltage I (HV I) is provided by the power supply integrated into the CE autosampler for separation. High voltage II (HV II) is provided by a separate power supply for electrospray. The mass spectrometer is grounded. The distance between the orifice of the emitter and the entrance of the mass spectrometer is ~ 2 mm. The distance between the end of the etched capillary in the emitter and the orifice of the emitter is less than 500 μ m. The size of the orifice of the electrospray emitter is 20 - 40 μ m.

Figure 2: Data of a standard protein mixture analyzed by the dynamic pH-junction-based CZE-MS. (A) This panel shows a base peak electropherogram. (B) This panel shows the number of theoretical plates (N) of three proteins. The sample injection volume was 500 nL. HV I was 30 kV for separation, and HV II was 2.2 kV for electrospray. The BGE was 5% (v/v) acetic acid, pH 2.4. The sample was in 50 mM ammonium bicarbonate (pH 8) and contains cytochrome c (0.01 mg/mL), lysozyme (0.01 mg/mL), β -casein (0.04 mg/mL), myoglobin (0.01 mg/mL), carbonic anhydrase (CA, 0.05 mg/mL), and bovine serum albumin (BSA, 0.1 mg/mL). No MS/MS spectra were acquired for the standard protein mixture sample.

Figure 3: Data of the *E.coli* proteome analyzed by the dynamic pH-junction-based CZE-MS/MS.

(A) This panel shows a zoomed-in view of the based peak electropherogram of the *E. coli* protein sample. **(B)** This panel shows an example PrSM identified through the TopPIC database search of the acquired MS/MS spectra. The general corresponding proteoform information, protein sequence, observed fragmentation pattern, and modifications are presented. The matched fragment ions from individual PrSM can also be viewed if needed. HV I was 20 kV for separation, and HV II was 2.2 kV for electrospray. The BGE was 10% (v/v) acetic acid, pH ~2.2. The sample was in 50 mM ammonium bicarbonate (pH 8) and the protein concentration was 2 mg/mL. The sample injection volume was 500 nL. A top8 DDA method was used to acquire the data.

DISCUSSION:

Here we provide a detailed protocol to use CZE-MS/MS for the high-resolution characterization of proteoforms in simple protein samples and for the large-scale identification of proteoforms in complex proteome samples. A diagram of the CZE-ESI-MS/MS system is shown in **Figure 1**. There are four critical steps in the protocol. First, the preparation of high-quality LPA coating on the inner wall of the separation capillary is extremely important. An LPA-coated separation capillary can reduce the EOF in the capillary, widen the separation window of CZE, and reduce protein adsorption on its inner wall¹⁹. The polymerization reaction for making the LPA coating is performed *via* filling the capillary with a mixture of acrylamide (monomer) and ammonium persulfate (polymerization initiator), followed by heating the capillary in a water bath to trigger the reaction²². The timing in the water bath is critical. Too short a reaction time will lead to an incomplete reaction and poor coating. We typically allow the reaction to go up to 50 minutes, allowing the reaction to proceed for a longer period for a better coating. With longer reaction periods, the capillary may become blocked and an HPLC pump might need to be used to push out the polymer inside of the capillary with water. One LPA-coated capillary can be continuously used for about one week, based on the literature data and our experience²².

The second critical step is coupling CZE to MS with the electro-kinetically pumped sheath flow CE-MS interface^{25,26}. The distance between the end of the separation capillary in the ESI emitter and the emitter orifice affects the sensitivity of CZE-MS significantly, and a shorter distance produces a higher sensitivity²⁵. The end of the separation capillary needs to be etched with HF, to reduce its outer diameter from 360 μm to less than 100 μm , allowing the end of the capillary to be pushed close to the emitter orifice. The distance between the emitter orifice and mass spectrometer entrance has been optimized to reach the best sensitivity and good stability²⁶. We typically keep the distance around 2 mm.

The third critical step is the CZE-MS and MS/MS analysis with the dynamic pH-junction-based online sample stacking¹⁹. The pH of the sample buffer and the BGE need to be significantly different in order to ensure efficient sample stacking. The protein concentration in the sample needs to be appropriate. If the protein concentration is too high, the proteins can be precipitated in the capillary during sample stacking. The protein concentrations of the samples used for **Figures 2 and 3** can be considered as typical examples. The last critical step is the database search with TopPIC for proteoform IDs³⁰. Multiple steps are required for file conversions before the

database search with TopPIC can be performed. A proper FDR filter is necessary to ensure the quality of the identified proteoforms. We typically use a 1% spectrum-level FDR to filter the database search results. We note that the top-down proteomics initiative is a very good resource for top-down proteomics methods and data analysis software^{32,33}.

We must note that some parts of the protocol may need to be modified slightly and some troubleshooting may be required. The time of polymerization in the water bath for the preparation of the LPA coating may vary at different labs. The time for HF etching of the separation capillary may need to be slightly modified due to a different temperature in different labs. When setting up the CE-MS interface, make sure the electrospray is stable before threading the separation capillary into the electrospray emitter. If no electrospray or an unstable electrospray is observed, check the interface to make sure there are no big bubbles in the emitter, in the T, or in the tubing that used to connect the sheath buffer vial and the T. In addition, the distance between the orifice of the emitter and the mass spectrometer entrance can affect the electrospray stability and is usually about 2 mm. The electrospray voltage can also influence the spray stability. For 20- to 40- μ m emitters, 2 - 2.2 kV is usually good enough. As described in the previous paragraph, the protein concentration in the sample needs to be appropriate. If the protein concentration is too high, the proteins can be precipitated in the capillary during sample stacking. Pay attention to the current profile on the CE autosampler computer. If the current is zero at the very beginning of the CZE run, it means that a plug of air is injected into the capillary during the sample injection step. Make sure the volume of the sample in the vial is large enough. If the current is normal at the beginning and suddenly becomes zero during the run, it usually means that the protein concentration is too high.

CZE-MS/MS based on this protocol enables the high-resolution characterization of simple protein samples and the large-scale top-down proteomics of complex proteomes. As an example, the CZE-MS approached the high-resolution characterization of a standard protein mixture containing six proteins¹⁹. As shown in **Figure 2**, the CZE-MS clearly separated the six proteins with a high separation efficiency, revealed three forms of β -casein, and detected many impurities. As another example, single-shot CZE-MS/MS reproducibly yielded over 500 proteoform IDs and 190 protein IDs from the *E. coli* proteome, using only 1 μ g of *E. coli* proteins with a 1% spectrum-level FDR¹⁹. The CZE-MS/MS improved the number of proteoform IDs from a complex proteome by at least three times, compared with previous single-shot CZE-MS/MS studies. As shown in **Figure 3A**, the CZE-MS/MS simultaneously reached a 500-nL sample-loading volume and a 90-min separation window for the analysis of the *E. coli* proteome¹⁹. The TopPIC software can provide comprehensive information about the identified proteoforms (**Figure 3B**). The CZE-MS/MS system provides the proteomics community with a useful tool for large-scale and high-resolution top-down proteomics.

CZE-MS/MS still has some limitations for deep top-down proteomics. Although we demonstrated that the CZE-MS/MS can reach over 500 proteoform IDs from the *E. coli* proteome in a single run, the number of proteoform IDs from single-shot CZE-MS/MS is only roughly 50% of that from state-of-the-art RPLC-MS/MS^{34,35}. In this protocol, only HCD is used for protein fragmentation, leading to limited fragmentation coverages of identified proteoforms. Several improvements can

be made to the CZE-MS/MS protocol to improve both the number and quality of the proteoform IDs from single-shot CZE-MS/MS. First, increasing the length of the separation capillary should provide a wider separation window, leading to more proteoform IDs. Second, using a mass spectrometer with a much higher resolving power, faster scan rate, and multiple fragmentation methods (*e.g.*, HCD,³⁶ electron transfer dissociation [ETD],³⁷ and ultraviolet photodissociation [UVPD]³⁸) will certainly improve both the number of proteoform IDs and the fragmentation coverages of identified proteoforms.

ACKNOWLEDGMENTS:

The authors thank Heedeok Hong's group at the Department of Chemistry, Michigan State University, for kindly providing the *Escherichia coli* cells for the experiments. The authors thank the support from the National Institute of General Medical Sciences, the National Institutes of Health (NIH) through Grant R01GM118470 (to X. Liu) and Grant R01GM125991 (to L. Sun and X. Liu).

DISCLOSURES:

The authors have nothing to disclose.

REFERENCES:

1. Aebersold, R. *et al.* How many proteoforms are there? *Nature Chemical Biology*. **14** (3), 206-214 (2018).
2. Tran, J. C. *et al.* Mapping intact protein isoforms in discovery mode using top-down proteomics. *Nature*, **480** (7376), 254-258 (2011).
3. Catherman, A. D. *et al.* Large-scale Top-down Proteomics of the Human Proteome: Membrane Proteins, Mitochondria, and Senescence. *Molecular and Cellular Proteomics*. **12** (12), 3465-3473 (2013).
4. Cai, W. *et al.* Top-Down Proteomics of Large Proteins up to 223 kDa Enabled by Serial Size Exclusion Chromatography Strategy. *Analytical Chemistry*. **89** (10), 5467-5475 (2017).
5. Toby, T. K., Fornelli, L., Kelleher, N. L. Progress in Top-Down Proteomics and the Analysis of Proteoforms. *Annual Review of Analytical Chemistry*. **9** (1), 499-519 (2016).
6. Jorgenson, J. W., Lukacs, K. D. Capillary Zone Electrophoresis. *Science*. **222** (4621), 266-272 (1983).
7. Keithley, R. B. Capillary electrophoresis with three-color fluorescence detection for the analysis of glycosphingolipid metabolism. *Analyst*. **138** (1), 164-70 (2013).
8. Han, X. *et al.* In-Line Separation by Capillary Electrophoresis Prior to Analysis by Top-Down Mass Spectrometry Enables Sensitive Characterization of Protein Complexes. *Journal of Proteome Research*. **13** (12), 6078-6086 (2014).

9. Sun, L., Zhu, G., Zhao, Y., Yan, X., Mou, S., Dovichi, N. J. Ultrasensitive and Fast Bottom-up Analysis of Femtogram Amounts of Complex Proteome Digests. *Angewandte Chemie International Edition*. **52** (51), 13661-13664 (2013).
10. Choi, S. B., Lombard-Banek, C., Muñoz-LLancao, P., Manzini, M. C., Nemes, P. Enhanced Peptide Detection Toward Single-Neuron Proteomics by Reversed-Phase Fractionation Capillary Electrophoresis Mass Spectrometry. *Journal of the American Society for Mass Spectrometry*. **29** (5), 913-922 (2018).
11. Sun, L., Knierman, M. D., Zhu, G., Dovichi, N. J. Fast Top-Down Intact Protein Characterization with Capillary Zone Electrophoresis-Electrospray Ionization Tandem Mass Spectrometry. *Analytical Chemistry*. **85** (12) 5989-5995 (2013).
12. Wang, Y., Fonslow, B. R., Wong, C. C., Nakorchevsky, A., Yates 3rd, J.R. Improving the comprehensiveness and sensitivity of sheathless capillary electrophoresis-tandem mass spectrometry for proteomic analysis. *Analytical Chemistry*. **84** (20), 8505-8513 (2012).
13. Zhang, Z., Yan, X., Sun, L., Zhu, G., Dovichi, N.J. Detachable strong cation exchange monolith, integrated with capillary zone electrophoresis and coupled with pH gradient elution, produces improved sensitivity and numbers of peptide identifications during bottom-up analysis of complex proteomes. *Analytical Chemistry*. **87** (8), 4572-4577 (2015).
14. Sun, L. *et al.* Over 10,000 peptide identifications from the HeLa proteome by using single-shot capillary zone electrophoresis combined with tandem mass spectrometry. *Angewandte Chemie International Edition in English*. **53** (50), 13931-13933 (2014).
15. Aebersold, R., Morrison, H. D. Analysis of dilute peptide samples by capillary zone electrophoresis. *Journal of Chromatography*. **516** (1), 79-88 (1990).
16. Britz-McKibbin, P., Chen, D. D. Y. Selective Focusing of Catecholamines and Weakly Acidic Compounds by Capillary Electrophoresis Using a Dynamic pH Junction. *Analytical Chemistry*. **72** (6), 1242-1252 (2000).
17. Zhao, Y., Sun, L., Zhu, G., Dovichi, N. J. Coupling Capillary Zone Electrophoresis to a Q Exactive HF Mass Spectrometer for Top-down Proteomics: 580 Proteoform Identifications from Yeast. *Journal of Proteome Research*. **15** (10), 3679-3685 (2016).
18. Zhu, G., Sun, L., Yan, X., Dovichi, N. J. Bottom-Up Proteomics of *Escherichia coli* Using Dynamic pH Junction Preconcentration and Capillary Zone Electrophoresis-Electrospray Ionization-Tandem Mass Spectrometry. *Analytical Chemistry*. **86** (13), 6331-6336 (2014).
19. Lubeckyj, R. A., McCool, E. N., Shen, X., Kou, Q., Liu, X., Sun, L. Single-Shot Top-Down Proteomics with Capillary Zone Electrophoresis-Electrospray Ionization-Tandem Mass

Spectrometry for Identification of Nearly 600 *Escherichia coli* Proteoforms. *Analytical Chemistry*. **89** (22), 12059-12067 (2017).

20. Busnel, J. M. High capacity capillary electrophoresis-electrospray ionization mass spectrometry: coupling a porous sheathless interface with transient-isotachopheresis. *Analytical Chemistry*. **82** (22), 9476-9483 (2010).

21. Zhang, Z., Peuchen, E. H., Dovichi, N. J. Surface-Confined Aqueous Reversible Addition-Fragmentation Chain Transfer (SCRAFT) Polymerization Method for Preparation of Coated Capillary Leads to over 10 000 Peptides Identified from 25 ng HeLa Digest by Using Capillary Zone Electrophoresis-Tandem Mass Spectrometry. *Analytical Chemistry*. **89** (12), 6774-6780 (2017).

22. Zhu, G., Sun, L., Dovichi, N.J. Thermally-initiated free radical polymerization for reproducible production of stable linear polyacrylamide coated capillaries, and their application to proteomic analysis using capillary zone electrophoresis-mass spectrometry. *Talanta*. **146**, 839-843 (2016).

23. Smith, R. D., Barinaga, C. J., Udseth, H. R. Improved Electrospray Ionization Interface for Capillary Zone Electrophoresis-Mass Spectrometry. *Analytical Chemistry*. **60** (16), 1948-1952 (1988).

24. Moini, M. Simplifying CE-MS Operation. 2. Interfacing Low-Flow Separation Techniques to Mass Spectrometry Using a Porous Tip. *Analytical Chemistry*. **79** (11), 4241-4246 (2007).

25. Wojcik, R., Dada, O. O., Sadilek, M., Dovichi, N. J. Simplified capillary electrophoresis nanospray sheath-flow interface for high efficiency and sensitive peptide analysis. *Rapid Communications in Mass Spectrometry*. **24** (17), 2554-2560 (2010).

26. Sun, L., Zhu, G., Zhang, Z., Mou, S., Dovichi, N. J. Third-Generation Electrokinetically Pumped Sheath-Flow Nanospray Interface with Improved Stability and Sensitivity for Automated Capillary Zone Electrophoresis-Mass Spectrometry Analysis of Complex Proteome Digests. *Journal of Proteome Research*. **14** (5), 2312-2321 (2015).

27. McCool, E. N. *et al.* Deep Top-Down Proteomics Using Capillary Zone Electrophoresis-Tandem Mass Spectrometry: Identification of 5700 Proteoforms from the *Escherichia coli* Proteome. *Analytical Chemistry*. **90** (9), 5529-5533 (2018).


28. Kessner, D., Chambers, M., Burke, R., Agus, D., Mallick, P. ProteoWizard: open source software for rapid proteomics tools development. *Bioinformatics*. **24** (21), 2534-2536 (2008).


29. TopPIC Suite. <http://proteomics.informatics.iupui.edu/software/topfd/>. June 8, 2017 (2017).

30. Kou, Q., Xun, L., Liu, X. TopPIC : a software tool for top-down mass spectrometry-based proteoform identification and characterization. *Bioinformatics*. **32** (22), 3495-3497 (2016).

31. Elias, J. E., Gygi, S. P. Target-decoy search strategy for increased confidence in large-scale protein identifications by mass spectrometry. *Nature Methods*. **4** (3), 207-214 (2007).
32. Methods – Consortium for Top Down Proteomics. <http://www.topdownproteomics.org/resources/methods/> (2017-2018).
33. Software – Consortium for Top Down Proteomics. <http://www.topdownproteomics.org/resources/software/> (2017-2018).
34. Ansong, C. *et al.* Top-down proteomics reveals a unique protein S-thiolation switch in *Salmonella Typhimurium* in response to infection-like conditions. *Proceedings of the National Academy of Sciences of the United States of America*. **110**(25), 10153-10158 (2013).
35. Shen, Y. *et al.* High-resolution ultrahigh-pressure long column reversed-phase liquid chromatography for top-down proteomics *Journal of Chromatography A*. **1498**, 99-110 (2017).
36. Olsen, J. V., Macek, B., Lange, O., Makarov, A., Horning, S., Mann, M. Higher-energy C-trap dissociation for peptide modification analysis. *Nature Methods*. **4** (9), 709-712 (2007).
37. Syka, J. E., Coon, J. J., Schroeder, M. J., Shabanowitz, J., Hunt, D. F. Peptide and protein sequence analysis by electron transfer dissociation mass spectrometry. *Proceedings of the National Academy of Sciences of the United States of America*. **101** (26), 9528-9533 (2004).
38. Shaw, J. B. *et al.* Complete Protein Characterization Using Top-Down Mass Spectrometry and Ultraviolet Photodissociation. *Journal of the American Chemical Society*. **135** (34), 12646-12651 (2013).


Figure 1





(B)

Peak	Protein	MW (Da)	N
d	myoglobin	17000	90000
e	CA	29000	45000
f3	β -casein	24000	32000



Name of Material/ Equipment	Company	Catalog Number
Fused silica capillary	Polymicro Technologies	1068150017
Sodium hydroxide pellets	Macron Fine Chemicals	7708-10
LC-MS grade water	Fisher Scientific	W6-1
Hydrochloric acid	Fisher Scientific	SA48-1
Methanol	Fisher Scientific	A456-4
3-(Trimethoxysilyl)propyl methacrylate	Sigma-Aldrich	M6514
Hydrofluoric acid	Acros Organics	423805000
Acrylamide	Acros Organics	164855000
Ammonium persulfate	Sigma-Aldrich	A3678
lysozyme	Sigma-Aldrich	L6876
Cytochrome C	Sigma-Aldrich	C7752
Myoglobin	Sigma-Aldrich	M1882
β -casein	Sigma-Aldrich	C6905
Carbonic anhydrase	Sigma-Aldrich	C3934
Bovine serum albumin	Sigma-Aldrich	A2153
Urea	Alfa Aesar	36428-36
DL-Dithiothreitol	Sigma-Aldrich	D0632
Iodoacetamide	Fisher Scientific	AC122270250
Formic Acid	Fisher Scientific	A117-50
C4 trap column	Sepax Technologies	110043-4001C
Acetonitrile	Fisher Scientific	A998SK-4
Ammonium bicarbonate	Sigma-Aldrich	1066-33-7
Nalgene rapid-flow filters	Thermo Scientific	126-0020
E. coli cells		
Dulbecco's phosphate-buffered saline	Sigma-Aldrich	D8537
BCA assay	Thermo Scientific	23250
Acetone	Fisher Scientific	A11-1
HPLC system for protein desalting	Agilent	1260 Infinity II
Acetic Acid	Fisher Scientific	A38-212
CE autosampler	CMP Scientific	ECE-001
Electro-kinetically pumped sheath flow interface	CMP Scientific	
Q Exactive HF Hybrid Quadrupole-Orbitrap Mass Spectrometer	Thermo Fisher Scientific	
Sutter flaming/brown micropipette puller	Sutter Instruments	P-1000
Ultrasonic cell disruptor for cell lysis	Branson	101063196

Vaccum concentrator

Thermo Fisher
Scientific

SPD131DDA-115

Comments/Description

50 μm i.d. 360 μm o.d.

Corrosive

Corrosive

Toxic, Health Hazard

Moisture and heat sensitive

Highly toxic

Toxic, health hazard

Health hazard, Oxidizer

Health Hazard

Health Hazard

Corrosive, Health Hazard

3 μm particles, 300 Å pores, 4.0 mm i.d. 10 mm long

Toxic, Oxidizer

0.2 μm CN membrane, and 50 mm diameter

K-12 MG1655

Model S-250A



1 Alewife Center #200
 Cambridge, MA 02140
 tel. 617.945.9051
www.jove.com

ARTICLE AND VIDEO LICENSE AGREEMENT

Title of Article:	Large-scale Top-down Proteomics Using Capillary Zone Electrophoresis Tandem Mass Spectrometry
Author(s):	Elijah N. McCool , Rachele Lubeckyj , Xiaojing Shen, Qiang Kou, Xiaowen Liu, Liangliang Sun

Item 1: The Author elects to have the Materials be made available (as described at <http://www.jove.com/publish>) via:

☒ Standard Access ☐ Open Access

Item 2: Please select one of the following items:

- ☒ The Author is **NOT** a United States government employee.
- ☐ The Author is a United States government employee and the Materials were prepared in the course of his or her duties as a United States government employee.
- ☐ The Author is a United States government employee but the Materials were NOT prepared in the course of his or her duties as a United States government employee.

ARTICLE AND VIDEO LICENSE AGREEMENT

1. **Defined Terms.** As used in this Article and Video License Agreement, the following terms shall have the following meanings: “**Agreement**” means this Article and Video License Agreement; “**Article**” means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein; “**Author**” means the author who is a signatory to this Agreement; “**Collective Work**” means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; “**CRC License**” means the Creative Commons Attribution-Non Commercial-No Derivs 3.0 Unported Agreement, the terms and conditions of which can be found at: <http://creativecommons.org/licenses/by-nc-nd/3.0/legalcode>; “**Derivative Work**” means a work based upon the Materials or upon the Materials and other pre-existing works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; “**Institution**” means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; “**JoVE**” means MyJoVE Corporation, a Massachusetts corporation and the publisher of The Journal of Visualized Experiments; “**Materials**” means the Article and / or the Video; “**Parties**” means the Author and JoVE; “**Video**” means any video(s) made by the Author, alone or in conjunction with any other parties, or by JoVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion

of the Article, and in which the Author may or may not appear.

2. **Background.** The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JoVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.

3. **Grant of Rights in Article.** In consideration of JoVE agreeing to publish the Article, the Author hereby grants to JoVE, subject to **Sections 4** and **7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the “Open Access” box has been checked in **Item 1** above, JoVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.

ARTICLE AND VIDEO LICENSE AGREEMENT

4. **Retention of Rights in Article.** Notwithstanding the exclusive license granted to JoVE in **Section 3** above, the Author shall, with respect to the Article, retain the non-exclusive right to use all or part of the Article for the non-commercial purpose of giving lectures, presentations or teaching classes, and to post a copy of the Article on the Institution's website or the Author's personal website, in each case provided that a link to the Article on the JoVE website is provided and notice of JoVE's copyright in the Article is included. All non-copyright intellectual property rights in and to the Article, such as patent rights, shall remain with the Author.

5. **Grant of Rights in Video – Standard Access.** This **Section 5** applies if the "Standard Access" box has been checked in **Item 1** above or if no box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby acknowledges and agrees that, Subject to **Section 7** below, JoVE is and shall be the sole and exclusive owner of all rights of any nature, including, without limitation, all copyrights, in and to the Video. To the extent that, by law, the Author is deemed, now or at any time in the future, to have any rights of any nature in or to the Video, the Author hereby disclaims all such rights and transfers all such rights to JoVE.

6. **Grant of Rights in Video – Open Access.** This **Section 6** applies only if the "Open Access" box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby grants to JoVE, subject to **Section 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Video in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Video into other languages, create adaptations, summaries or extracts of the Video or other Derivative Works or Collective Works based on all or any portion of the Video and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. For any Video to which this **Section 6** is applicable, JoVE and the Author hereby grant to the public all such rights in the Video as provided in, but subject to all limitations and requirements set forth in, the CRC License.

7. **Government Employees.** If the Author is a United States government employee and the Article was prepared in the course of his or her duties as a United States government employee, as indicated in **Item 2** above, and any of the licenses or grants granted by the Author hereunder exceed the scope of the 17 U.S.C. 403, then the rights granted hereunder shall be limited to the maximum

rights permitted under such statute. In such case, all provisions contained herein that are not in conflict with such statute shall remain in full force and effect, and all provisions contained herein that do so conflict shall be deemed to be amended so as to provide to JoVE the maximum rights permissible within such statute.

8. **Protection of the Work.** The Author(s) authorize JoVE to take steps in the Author(s) name and on their behalf if JoVE believes some third party could be infringing or might infringe the copyright of either the Author's Article and/or Video.

9. **Likeness, Privacy, Personality.** The Author hereby grants JoVE the right to use the Author's name, voice, likeness, picture, photograph, image, biography and performance in any way, commercial or otherwise, in connection with the Materials and the sale, promotion and distribution thereof. The Author hereby waives any and all rights he or she may have, relating to his or her appearance in the Video or otherwise relating to the Materials, under all applicable privacy, likeness, personality or similar laws.

10. **Author Warranties.** The Author represents and warrants that the Article is original, that it has not been published, that the copyright interest is owned by the Author (or, if more than one author is listed at the beginning of this Agreement, by such authors collectively) and has not been assigned, licensed, or otherwise transferred to any other party. The Author represents and warrants that the author(s) listed at the top of this Agreement are the only authors of the Materials. If more than one author is listed at the top of this Agreement and if any such author has not entered into a separate Article and Video License Agreement with JoVE relating to the Materials, the Author represents and warrants that the Author has been authorized by each of the other such authors to execute this Agreement on his or her behalf and to bind him or her with respect to the terms of this Agreement as if each of them had been a party hereto as an Author. The Author warrants that the use, reproduction, distribution, public or private performance or display, and/or modification of all or any portion of the Materials does not and will not violate, infringe and/or misappropriate the patent, trademark, intellectual property or other rights of any third party. The Author represents and warrants that it has and will continue to comply with all government, institutional and other regulations, including, without limitation all institutional, laboratory, hospital, ethical, human and animal treatment, privacy, and all other rules, regulations, laws, procedures or guidelines, applicable to the Materials, and that all research involving human and animal subjects has been approved by the Author's relevant institutional review board.

11. **JoVE Discretion.** If the Author requests the assistance of JoVE in producing the Video in the Author's facility, the Author shall ensure that the presence of JoVE employees, agents or independent contractors is in accordance with the relevant regulations of the Author's institution. If more than one author is listed at the beginning of this Agreement, JoVE may, in its sole

ARTICLE AND VIDEO LICENSE AGREEMENT

discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JoVE reserves the right, in its absolute and sole discretion and without giving any reason therefore, to accept or decline any work submitted to JoVE. JoVE and its employees, agents and independent contractors shall have full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JoVE has sole discretion as to the method of making and publishing the Materials, including, without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

12. **Indemnification.** The Author agrees to indemnify JoVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JoVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JoVE, making of videos by JoVE, or publication in JoVE or elsewhere by JoVE. The Author shall be responsible for, and shall hold JoVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to

the making of a video by JoVE its employees, agents or independent contractors. All sterilization, cleanliness or decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's expense. All indemnifications provided herein shall include JoVE's attorney's fees and costs related to said losses or damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JoVE, its employees, agents or independent contractors.

13. **Fees.** To cover the cost incurred for publication, JoVE must receive payment before production and publication the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JoVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.

14. **Transfer, Governing Law.** This Agreement may be assigned by JoVE and shall inure to the benefits of any of JoVE's successors and assignees. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to be one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement is required per submission.

CORRESPONDING AUTHOR

Name:

Liangliang Sun

Department:

Chemistry

Institution:

Michigan State University

Title:

Assistant Professor

Signature:



Date:

06-09-2018

Please submit a **signed** and **dated** copy of this license by one of the following three methods:

1. Upload an electronic version on the JoVE submission site
2. Fax the document to +1.866.381.2236
3. Mail the document to JoVE / Attn: JoVE Editorial / 1 Alewife Center #200 / Cambridge, MA 02140

Responses to the comments

We thank the comments from the editor and the reviewers. Those comments really helped us improve the manuscript. We made changes in the revised manuscript based on those comments. Since many changes have been made based on the comments, we did not highlight the changes in the revised manuscript. We provided responses to the comments one by one.

Editorial comments:

Changes to be made by the Author(s) regarding the written manuscript:

1. Please take this opportunity to thoroughly proofread the manuscript to ensure that there are no spelling or grammar issues.

Reply: We have thoroughly reviewed the manuscript for these issues.

2. Please label/number the institutional affiliation of each author sequentially.

Reply: Institutional affiliation labels have been added to the manuscript.

3. Please expand the Short Abstract to briefly describe the applications of this protocol.

Reply: we made the change based on the comment.

4. Please convert centrifuge speeds to centrifugal force (x g) instead of revolutions per minute (rpm).

Reply: We made the change in the revised manuscript.

5. Please revise the protocol text to avoid the use of any personal pronouns (e.g., "we", "you", "our" etc.).

Reply: we made the changes.

6. JoVE cannot publish manuscripts containing commercial language. This includes trademark symbols (™), registered symbols (®), and company names before an instrument or reagent. Please remove all commercial language from your manuscript and use generic terms instead. All commercial products should be sufficiently referenced in the Table of Materials and Reagents. For example: Kimwipe, Eppendorf, Nalgene Rapid-Flow, Branson Sonifier, VWR Scientific, Agilent, etc.

Reply: All language in the protocol text has been changed to avoid using commercial language.

7. Please revise the protocol (e.g., 1.1.2, 1.2.3, 3.1.1-3.2.4, etc.) to contain only action items that direct the reader to do something (e.g., "Do this," "Ensure that," etc.). The actions should be described in the imperative tense in complete sentences wherever possible. Avoid usage of phrases such as "could be," "should be," and "would be" throughout the Protocol. Any text that cannot be written in the imperative tense may be added as a "Note." Please include all safety

procedures and use of hoods, etc. Please move the discussion about the protocol to the Discussion.

Reply: Tense has been corrected in the protocol and certain phrases have been avoided when possible. Safety has been added to appropriate steps.

8. The Protocol should be made up almost entirely of discrete steps without large paragraphs of text between sections. Please simplify the Protocol so that individual steps contain only 2-3 actions per step and a maximum of 4 sentences per step. Use sub-steps as necessary. Please move the discussion about the protocol to the Discussion.

Reply: Individual steps have been shortened to contain less than or equal to 4 sentences.

9. Please note that your protocol will be used to generate the script for the video and must contain everything that you would like shown in the video. Please add more specific details (e.g. button clicks for software actions, numerical values for settings, etc.) to your protocol steps. There should be enough detail in each step to supplement the actions seen in the video so that viewers can easily replicate the protocol. Some examples:

1.1.3: How is methanol (500 μ L) injected, using a syringe pump?

2.2.4: What is used to cut?

4.3.2.1-4.3.2.3, 5.2.3: Please describe how these are done. Software must have a GUI (graphical user interface) and software steps must be more explicitly explained ('click', 'select', etc.).

Reply: Changes have been made in the revised manuscript based on the comments. Details have been added where appropriate while trying to keep language as general as possible for other instrument users.

10-12. After you have made all the recommended changes to your protocol (listed above), please highlight 2.75 pages or less of the Protocol (including headings and spacing) that identifies the essential steps of the protocol for the video, i.e., the steps that should be visualized to tell the most cohesive story of the Protocol. Please highlight complete sentences (not parts of sentences). Please ensure that the highlighted part of the step includes at least one action that is written in imperative tense. Please include all relevant details that are required to perform the step in the highlighting. For example: If step 2.5 is highlighted for filming and the details of how to perform the step are given in steps 2.5.1 and 2.5.2, then the sub-steps where the details are provided must be highlighted.

Reply: we highlighted the essential steps of the protocol in the revised manuscript based on the comment.

13. Please revise to explain the Representative Results in the context of the technique you have described, e.g., how do these results show the technique, suggestions about how to analyze the outcome, etc. The paragraph text should refer to all of the figures. However for figures showing the experimental set-up, please reference them in the Protocol.

Reply: The representative results have been revised and should have a better explanation of the technique for the readers.

14. As we are a methods journal, please revise the Discussion to explicitly cover the following in detail in 3-6 paragraphs with citations:

- a) Critical steps within the protocol
- b) Any modifications and troubleshooting of the technique
- c) Any limitations of the technique
- d) The significance with respect to existing methods
- e) Any future applications of the technique

Reply: we made the changes based on the comment.

Reviewers' comments:

Reviewer #1:

Manuscript Summary:

This is a well-written, well-referenced and timely manuscript describing the use of CZE-MS for proteoform analysis. Proteoform analysis is an emerging area of proteomics and the use of CZE for proteoform identification is under-developed. This article should have a strong impact on the field. The descriptions contain sufficient detail such that others should easily be able to easily reproduce the results.

Major Concerns:

None.

Minor Concerns:

-The introduction is especially well-written. The protocol contains several minor grammatical errors that should be corrected by the copy editors.

Reply: Thank you for the nice comment. We went through the manuscript carefully and fixed some grammatical errors.

-The quality of the figures in the review copy were low. These figures should be substituted with high-quality/high-resolution graphs prior to final publication.

Reply: We provided high-resolution figures in the revision.

-It would have been nice for the authors to show the improvement that injecting a larger amount of sample allows (e.g. by reporting the number of proteoforms observed in the same set-up when injecting 1/10th the amount). However, I don't view this as a necessary experiment prior to acceptance. Simply would have better validated the claims made. This protocol will be highly valuable and is the major contribution.

Reply: In our paper "Single-Shot Top-Down Proteomics with Capillary Zone Electrophoresis-Electrospray Ionization-Tandem Mass Spectrometry for Identification of Nearly 600 Escherichia coli Proteoforms. Analytical Chemistry. 89 (22), 12059-12067 (2017)", we showed the difference in proteoform identifications from 100 nL injection volume to 1 µL injection volume.

Reviewer #2:

The manuscript: Large-scale Top-down Proteomics Using Capillary Zone Electrophoresis Tandem Mass Spectrometry from Liangliang et al present an extended description of the approach used to perform top down analysis of intact proteins using CE-MS with neutral coated capillary, taken from previously published on Analytical Chemistry 89 (22), 12059-12067 (2017) and Analytical Chemistry. 90 (9), 5529-5533 (2018).

I found the text part of the material well written and organized and I look forward to see the video part of the contribution, being particularly interested on the fabrication of the sheath flow interface of the described setup.

One minor note:

As the article may interest also people approaching this field I suggest to include references on the top down proteomics initiative as from this resource are accessible settings and method that would extend the compatibility of the method to other MS platforms:

<http://www.topdownproteomics.org/resources/methods/>

And also include a repository concerning tools for TD data analysis

<http://www.topdownproteomics.org/resources/software/>

Reply: Thank you for the nice comments. We included those two references in the revised manuscript.

Reviewer #3:

Manuscript Summary:

The authors have submitted a JOVE article to describe their CE-MS setup for intact protein analysis. The topic is of interest to a general audience. The items below are provided to improve to the readability of the article. At times, the language gets very lax and is not technical.

Major Concerns:

-line 93: How is CE more sensitive than HPLC? Sensitivity refers to the ability to detect changes of a signal. Therefore, CE is not generally considered a detector.

Reply: As we described in the manuscript, “CZE-MS has drastically higher sensitivity than widely used reversed-phase liquid chromatography (RPLC)-MS for analysis of intact proteins.” We did not say “CE is more sensitive than HPLC”. When a same amount of protein material is loaded into the CZE capillary and the RPLC column, CZE can deliver much narrower protein peaks to ESI-MS for detection with much lower sample loss than RPLC. The narrower peaks and lower sample loss from CZE are due to the lack of stationary phase with a large surface area. Therefore, CZE-MS is more sensitive than HPLC-MS for protein detection.

-line 155: function should be functionalized?

Reply: we made the change in the revised manuscript.

-line 176: A more detailed safety description should be provided for using HF including the antedote

Reply: Thank you for the great comment. We added the detailed safety description in the revised manuscript.

-line 228: Settings are missing for this sonication. Is this a bath or a probe sonicator? Needs details.

Reply: We made the change based on the comment.

-line 252: what units? wt/v or v/v?

Reply: We made the change.

-line 255: ".... the goal is to approach the disulfide bonds..." I have no idea what this means
Reply:

Reply: we made the change.

-line 305: Fix is a generic term for repair. Do you mean to mount?

Reply: Thank you for the comment. We made the change.

-line 383: I would refrain from "obviously" it is negative voice. Such an article as this is a teaching tool.

Reply: we made the change.

-Figure 1: please refer as a diagram as opposed to a sketch

Reply: We made the change.

-Figure 3: Why is the top of the electropherogram cut off? It should be scaled.

Reply: There are very high abundant peaks in the electropherogram. If we scaled to 100%, the low abundant peaks cannot be seen clearly. Therefore, we used a zoomed-in electropherogram to show the low abundant peaks better.