

Journal of Visualized Experiments

Identification of coding and non-coding RNA classes expressed in swine whole blood --Manuscript Draft--

Article Type:	Invited Methods Article - JoVE Produced Video
Manuscript Number:	JoVE58689R2
Full Title:	Identification of coding and non-coding RNA classes expressed in swine whole blood
Keywords:	mRNA; Non-coding RNA; RNA-Seq; Whole blood; Epigenetic; Globin reduction
Corresponding Author:	Laura C Miller UNITED STATES
Corresponding Author's Institution:	
Corresponding Author E-Mail:	laura.miller@ars.usda.gov
Order of Authors:	Damarius S. Fleming Sarah J. Anderson Laura C Miller
Additional Information:	
Question	Response
Please indicate whether this article will be Standard Access or Open Access.	Standard Access (US\$2,400)
Please indicate the city, state/province, and country where this article will be filmed . Please do not use abbreviations.	Building 20, 1920 Dayton Ave, Ames IA 50010

TITLE:

Identification of Coding and Non-coding RNA Classes Expressed in Swine Whole Blood

AUTHORS AND AFFILIATIONS:

Damarius S. Fleming^{1,2}, Sarah J. Anderson¹, Laura C. Miller¹

¹Virus and Prion Research Unit, National Animal Disease Center, USDA, Agricultural Research Service, Ames, IA, USA

²ORAU/ORISE Oak Ridge, TN, USA

Corresponding Author:

Laura C. Miller (laura.miller@ars.usda.gov)

Email Addresses of Co-authors:

Damarius Fleming (damarius.fleming@ars.usda.gov)

Sarah Anderson (sarah.anderson@ars.usda.gov)

KEYWORDS:

mRNA, Non-coding RNA, RNA-Seq, Whole blood, Epigenetic, Globin reduction

SUMMARY:

Here, we present a protocol optimized for the processing of coding (mRNA) and non-coding RNAs (ncRNA) globin reduced RNA-seq libraries from a single whole blood sample.

ABSTRACT:

The advent of innovative and increasingly powerful next generation sequencing techniques has opened new avenues into the ability to examine the underlying gene expression related to biological processes of interest. These innovations not only allow researchers to observe expression from the mRNA sequences that code for genes that effect cellular function, but also the non-coding RNA (ncRNA) molecules that remain untranslated, but still have regulatory functions. Although researchers have the ability to observe both mRNA and ncRNA expression, it has been customary for a study to focus on one or the other. However, when studies are interested in both mRNA and ncRNA expression, many times they use separate samples to examine either coding or non-coding RNAs due to the difference in library preparations. This can lead to the need for more samples which can increase time, consumables, and animal stress. Additionally, it may cause researchers to decide to prepare samples for only one analysis, usually the mRNA, limiting the number of biological questions that can be investigated. However, ncRNAs span multiple classes with regulatory roles that effect mRNA expression. Because ncRNA are important to fundamental biologic processes and disorder of these processes in during infection, they may, therefore, make attractive targets for therapeutics. This manuscript demonstrates a modified protocol for the generation mRNA and non-coding RNA expression libraries, including viral RNA, from a single sample of whole blood. Optimization of this protocol, improved RNA purity, increased ligation for recovery of methylated RNAs, and omitting size selection, to allow capture of more RNA species.

INTRODUCTION:

Next generation sequencing (NGS) has emerged as a powerful tool for the investigation of the changes that occur at the genomic level of biological organisms. Sample preparation for NGS methods can be varied depending on the organism, tissue type, and more importantly the questions the researchers are keen to address. Many studies have turned to NGS as a means of studying the differences in gene expression between states such as healthy and sick individuals¹⁻⁴. The sequencing take place on a whole genome basis and allows a researcher to capture the most, if not all, of the genomic information for a particular genetic marker at a time point.

The most common markers of expression observed are the messenger RNAs (mRNAs). The most used procedures for prepping libraries for RNA-seq are optimized for the recovery of mRNA molecules through the use of a series of purifications, fragmentations, and ligations^{5,6}. However, the decision on how a protocol is to be performed relies heavily on the sample type and the questions being posed about said sample. In most cases total RNA is extracted; yet, not all RNA molecules are of interest and in cases such as mRNA expression studies overly abundant RNA species, like ribosomal RNAs (rRNA) need to be removed to increase the number of detectable transcripts associated with the mRNAs. The most popular and widely used method for removing the abundant rRNA molecules is the reduction of polyadenylated RNA transcripts referred to as polyA depletion⁷. This approach works well for the analysis of mRNA expression as it does not affect the mRNA transcripts. However, in studies that are interested in non-coding or viral RNAs, polyA depletion also removes these molecules.

Many studies choose to focus on the RNA sequence library preparation to examine either mRNA expression (coding) or a particular class of small or large non-coding RNA. Although there are other procedures⁸ like ours that allow for the dual sample preparation, many studies prepare libraries from separate samples for separate studies when available. For a study like ours, this would normally require multiple blood samples increasing time, consumables, and animal stress. The goal of our study was to be able to use whole blood from animals to identify and quantify the different classes of both mRNA and non-coding RNA expressed between healthy and highly pathogenic porcine reproductive and respiratory syndrome virus (HP-PRRSV) challenged pigs^{9,10} despite having only a single whole blood sample (2.5 mL) from each pig. In order to do this, we needed to optimize the typical extraction and library creation protocols to generate the proper data to allow for analysis of both mRNA and non-coding RNA (ncRNA) expression¹¹ from a single sample.

This prompted a need for a protocol that allowed for mRNA and non-coding RNA analysis because the available standard kits and methods for RNA-extraction and library creation were intended chiefly for mRNA and use a poly-A depletion step¹². This step would have made it impossible to recover non-coding RNA or viral transcripts from the sample. Therefore, an optimized method was needed that allowed for total RNA extraction without sample polyA depletion. The method presented in this manuscript has been optimized to allow for the use of whole blood as a sample type and to build sequencing libraries for both mRNA and ncRNAs of

small and large sizes. The method has been optimized to allow for the analysis of all detectable non-coding RNAs as well as retain viral RNAs for later investigation¹³. In all, our optimized library preparation protocol allows for the investigation of multiple RNA molecules from a single whole blood sample.

The overall goal behind the use of this method was to develop a process that allowed for the collection of both non-coding RNA and mRNA from one sample of whole blood. This allows us to have mRNA, ncRNA, and viral RNA for each animal in our study sourced from a single sample⁹. This, ultimately, allows for more scientific discovery without additional animal costs and gives a more complete picture of the expression of each individual sample. The described method allows for the examination of the regulators of gene expression as well as allowing for completion of correlative studies comparing both mRNA and non-coding RNA expression using a single whole blood sample. Our study used this protocol to examine the changes in gene expression and possible epigenetic regulators in virally infected 9-week old male commercial pigs.

PROTOCOL:

Animal protocols were approved by the National Animal Disease Center (USDA-NADC-ARS) Animal Care and Use Committee.

1. Collection of Swine Blood Samples

1.1 Collect blood samples into RNA tubes. Collect ~2.5 mL or more if larger collection tubes are available.

2. Processing of Swine Blood Samples

2.1 Centrifuge the blood tubes at 5020 x g for 10 min at room temperature (15-25 °C). If processing frozen samples incubate tube at room temperature for a **minimum** of 2 h prior to centrifugation.

2.2 Remove supernatant and add 8 mL of RNase-free water to the pellet. Close and vortex the pellet until it is visibly dissolved. Centrifuge sample tubes at 5020 x g for 10 min at room temperature to recover the pellet. Discard all supernatant and preserve the pellet.

3. Organic Extraction for Total RNA and Small RNA (miRNA Isolation Kit)

3.1 Begin total RNA extraction by pipetting 300 µL of lysis binding buffer to the pellet from the step 2.2.

3.2 Vortex and transfer the mixture to a new labelled 1.5 mL centrifuge tube. Add 30 µL of homogenate additive from the kit. Vortex the tube and place on the ice for 10 min.

3.3 Remove the tube and add 300 μ L of the acid-phenol: chloroform reagent from the kit.
Vortex the tube to mix. Centrifuge at 10,000 x g for 5 min at room temperature.

3.4 Carefully remove aqueous phase (300-350 μ L) to a fresh tube. Note the volume for the next step.

4. Total RNA Isolation Procedure

4.1 Based on the amount of aqueous recovery (300-350 μ L) add 1.25x volume of 100% (~375 μ L) ethanol to aqueous phase. Mix the sample using a pipette.

4.2 Set-up new collection tubes containing a filter cartridge for each sample. Pipette ~ 675 μ L of the lysate/ethanol mixture onto the filter cartridge.

Note: Do not add > 700 μ L to filter cartridge at one time. For larger volumes apply in succession.

4.3 Centrifuge briefly (~15-20 s) at 10,000 x g to pass liquid through the filter. Do not spin harder than this.

4.4 Discard the flow through and repeat the centrifugation with the remaining lysate/ethanol mixture until it has all been applied. Retain the same filter cartridge and collection tube for the next step.

4.5 Add 700 μ L of Wash Solution 1 from the kit to the filter cartridge and centrifuge briefly (~10 s) to pull through the filter. Discard the flow through and retain the same filter cartridge and collection tube.

4.6 Add 500 μ L of Wash Solution 2/3. Centrifuge to draw the liquid through the filter cartridge. Discard the flow through. Repeat wash step.

4.7 In the same tube, spin the filter cartridge an additional 60 s to remove any residual liquid from filter. Transfer the filter cartridge to a fresh collection tube.

4.8 Add 100 μ L of pre-heated (95 $^{\circ}$ C) nuclease-free water to the center of the filter cartridge. Spin for approximately 20-30 s at the tabletop centrifuge max speed.

Note: RNA is contained in the eluate and can now be further processed or stored at -20 $^{\circ}$ C or below. Enrichment for small RNAs was not performed.

5. Globin Reduction (based on a protocol optimized for porcine whole blood samples)^{14,15}

Note: Globin reduction is performed so that libraries are not overpopulated with reads mapping to globin genes, which would lower the number of reads available to map to other genes of greater interest^{14,15}.

5.1 Hybridization with globin reduction oligos

5.1.1 Denature the RNA (6 µg RNA sample in maximum 7 µL volume) by pipetting the extracted sample into a 0.2 mL thin-walled nuclease-free reaction tube and placing in a thermal cycler at 70 °C for 2 min. It is key to ice tubes immediately after the first denature step for optimal RNA quality. No DNase treatment is needed.

5.1.2 While tubes cool prepare a 400 µL of 10x globin reduction oligo mix: 100 µL each of two *HBA* oligos (5'-GATCTCCGAGGCTCCAGCTTAACGGT-3', and 5'-TCAACGATCAGGAGGTCAGGGTGCAA-3') at 30 µM, two *HBB* oligos (5'-AGGGGAAGTCTAGTGGTACTTGTGGGC-3', and 5'-GGTTCAGAGGAAAAAGGGCTCCTCCT-3') at 120 µM per reaction, yielding a final concentration of 7.5 µM *HBA* oligos and 30 µM *HBB* oligos. Prepare 10x oligo hybridization buffer: 100 mM Tris-HCL, pH 7.6; 200 mM KCl.

5.1.3 Prepare the hybridization mix: 6 µg of RNA sample (maximum 7 µL volume), 2 µL of the 400 µL of 10x globin reduction oligo mix (final concentration 2X), 1 µL of 10x oligo hybridization buffer (final concentration 1x). Add nuclease-free water to a final volume of 10 µL.

5.1.4 Set thermocycler at 70 °C for 5 min. Cool immediately to 4 °C and proceed to RNase H digestion.

5.2 RNase H digestion

5.2.1 Dilute 10x RNase H (10 U/ µL) to 1x RNase H with 1x RNase H buffer.

Note: The RNase buffer comes at 10x and will need to be diluted with 1x RNase H buffer to 1x prior to use.

5.2.2 Prepare RNase H reaction mix by combining: 2 µL of 10x RNase buffer, 1 µL of RNase inhibitor in 2 µL of 1x RNase H, and 5 µL of nuclease-free water to a total volume 10 µL.

5.2.3 Mix thoroughly the Globin Reduction hybridization samples with 10 µL of the RNase H reaction mix and digest at 37 °C for 10 min. Cool to 4 °C.

5.2.4 Stop digestion by addition of 1.0 µL of 0.5 M EDTA to each sample and proceed immediately to the cleanup step.

5.3 RNase H-Treated Total RNA Cleanup.

5.3.1 Purify the RNase H treated RNA using a silica-membrane-based elution purification cleanup kit according to manufacturer's instructions. Premix buffers: For the mild washing buffer, add 44 mL of 100% ethanol.

5.3.2 Do not transfer the sample to a new tube. Add 80 μ L of RNase-free water and 350 μ L of lysis buffer. Add 250 μ L of 100% ethanol to the diluted RNA and mix well by pipetting.

Note: Do not centrifuge. Proceed immediately to step 5.3.3.

5.3.3 Now transfer the 700 μ L sample to an elution filter cartridge placed in a 2 mL collection tube to collect the flow through. Centrifuge for 15 s at $\geq 8,000 \times g$. Discard the flow-through.

5.3.4 Repeat this by placing the elution filter cartridge into a new 2 mL collection tube. Add 500 μ L of the mild washing buffer to the filter cartridge and centrifuge for 15 s at $\geq 8,000 \times g$ in order to wash the filter cartridge membrane. Discard the flow-through. Reuse the collection tube in step 5.3.5.

5.3.5 Now using the same sample tube, add another 500 μ L of 80% ethanol to the filter cartridge. This time centrifuge the tube for 2 min at $\geq 8,000 \times g$. Collect the elution spin column for the next step and discard both the flow-through and collection tube.

5.3.7 Put the elution filter cartridge from the last step into a new 2 mL collection tube. Leave the lid open on filter cartridge, and centrifuge at full speed for 5 min to dry the spin column membrane and prevent ethanol carry over. Discard the flow-through and collection tube.

Note: To avoid damage angle the lids to point in a direction opposite to that of the rotor.

5.3.8. Take the dried filter cartridge and place into a new 1.5 mL collection tube. Add 14 μ L of RNase-free water to filter cartridge membrane making sure to add the RNase-free water directly to the center. Centrifuge for 60 s at full speed to elute the RNA.

5.4 Assess quality of globin-reduced RNA samples (step 6). Proceed to mRNA sample preparation (step 7) and small RNA library preparation (step 8)^{16,17}.

Note: Globin-reduced RNA samples can now be stored at -20 °C, however storage at -80 °C is recommended for long-term preservation.

6. Assessment of RNA

6.1 Quantify RNA concentration using a spectrophotometer. Examine the wavelength ratio of 260 and 280 nm. Ratios of ~ 2 or greater are considered pure for RNA and lower values indicate some contamination. The instrument uses this ratio to determine RNA concentration as ng/ μ L.

6.2 Assess RNA quality by using 1 μ L (100 ng) of sample on an appropriate chip. Final product should be a RIN of \sim 2 or higher and a peak at \sim 280 nt for mRNA single read libraries; small RNA libraries peaks at 143 correspond to miRNAs.

7. Stranded Total RNA Sample Preparation for the mRNA and long ncRNA libraries.¹⁶

Note: Bring elution buffer and rRNA removal beads to room temperature. Pre-label 0.2ml thin-walled PCR tubes (plates may also be used). Protocol steps based on manufacturer's instructions¹⁶.

7.1.1 Start with 4 μ L of globin-reduced RNA. Add 6 μ L of nuclease-free water, 5 μ L of rRNA binding buffer, and 5 μ L of rRNA removal mix per tube (1st tube). Pipette to mix and re-cap. Place in thermocycler with 100 °C heated-lid for 5 min at 68° C. Remove and leave at room temperature for 60 s.

7.1.2 Resuspend rRNA removal beads by vortex; add 35 μ L of beads to new (2nd) PCR tube and pipette sample from the 1st tubes onto beads in the 2nd tubes. Incubate 2nd PCR tube at room temperature for 3 min. Place on the magnetic stand for 7 min.

Note: Mix each sample thoroughly by pipetting for optimal rRNA depletion. Raising/Lowering tube in the stand may help expedite separation process.

7.1.4 Transfer the supernatant from the 2nd PCR tube to matching 3rd PCR tube and place on to the magnetic stand for a minimum of 60 s. Repeat the transfer into a new PCR tube only if the beads are not moved to the tube sides.

7.1.5 Mix Sample Purification beads by vortex; add 99 μ L to each 3rd PCR tube. Allow to sit at room temperature for 15 min. Place on magnetic stand for 5 min. Ensure beads moves to sides. Pipette to discard the supernatant.

7.1.6 Keep 3rd tube set on magnetic stand. Add 200 μ L of 70% ethanol whilst being careful not to jostle the beads. Allow to sit for at least 30 s and pipette to discard the supernatant. Repeat step.

7.1.7 Allow sample to dry at room temperature for 15 min on the magnetic stand. Centrifuge the kit elution buffer for 5 s at 600 x g. Remove tube and add 11 μ L of elution buffer mixing thoroughly. Incubate for 2 min on the bench then at least 5 min on the magnetic stand at room temperature.

7.1.8 Transfer 8.5 μ L of supernatant from the 3rd to a new (4th) PCR tube. Add 8.5 μ L of Elute/Prime/Fragment High mix from the kit. Mix well. Cap and place in a thermocycler with pre-heated lid for 8 min at 94 °C, hold at 4 °C. Remove and centrifuge briefly.

7.2 Synthesize First Strand cDNA

305

306 7.2.1 Allow first strand synthesis mix from kit to come to room temperature and centrifuge at
307 600 × g for 5 s. Transfer 50 µL of reverse transcriptase into the first strand synthesis mix.

308 Reverse transcriptase can be added to first strand synthesis at a ratio of 1:9.

309

310 7.2.2 Pipette 8 µL of the combined reverse transcriptase/first strand synthesis mix into the 4th
311 tube with sample. Cap and centrifuge 600 x g for 5 s. Place into thermocycler with pre-heated
312 lid set at 100 °C. Run for: 10 min at 25 °C, 15 min at 42 °C, 15 min at 70 °C, then rest at 4 °C.

313 Final volume is ~25 µL per well. Move immediately to next step.

314

315 7.3 Synthesize Second Strand cDNA

316

317 7.3.1 Bring End Repair Reagent (ERR) and Second Strand Mix (SSM) to room temperature, and
318 centrifuge at 600 x g for 5 s. Mix ERR with resuspension buffer at 1:50 dilution. Uncap 4th tube
319 and add 5 µL of diluted ERR and 20 µL of SSM to each; Pipette to mix well. Final volume ~50 µL
320 ds cDNA.

321

322 7.3.2 Cap and incubate in thermocycler for 1 hour at 16 °C. When the cycle is complete, remove
323 caps and allow to come to room temperature on bench top.

324

325 7.4 ds cDNA clean-up step

326

327 7.4.1 Begin by mixing the Solid Phase Reversible Immobilization (SPRI) paramagnetic beads by
328 vortex. Transfer 90 µL of the beads to the sample in the 4th tube (ds cDNA) and mix. Final
329 volume is 140 µL. Allow tubes to sit at room temperature for 15 min. incubate on magnetic
330 stand for ~5 min. Remove ~135 µL of supernatant from each well. There should 5 µL left in each
331 well.

332

333 7.4.2 Leave the tube on the magnetic stand and wash by adding 200 µL of 80% ethanol. Allow
334 to sit for 30 s then remove and discard supernatant. Repeat wash step. Leave tubes on the
335 magnetic stand for 15 min to allow to dry.

336

337 7.4.3 Centrifuge the resuspension buffer at 600 x g for 5 s after coming to room temperature.
338 Remove tubes from stand. Transfer 17.5 µL of resuspension buffer to the tube and mix. Let
339 tubes incubate on the bench top for 2 min then move to magnetic stand for 5 min.

340

341 7.4.4 Pipette 15 µL of the supernatant containing the ds cDNA samples to new set of (5th) 0.2
342 mL thin-walled tubes. Move to next step promptly or seal and store at -15 °C to -25 °C for no
343 more than 7 days.

344

345 7.5 Adenylate 3' Ends.

346

347 7.5.1 Pipette 2.5 µL of room temperature resuspension buffer into the sample tube then add
348 12.5 µL of thawed A-Tailing mix. Mix completely by pipetting.

Note: No A-tailing control used.

7.5.2 Cap and incubate in a thermocycler with 100 °C pre-heated lid. Run at 37 °C for 30 min, then at 70 °C for 5 min with heated lid. Allow the thermocycler to rest at 4 °C.

7.6 Ligate index adaptors

7.6.1 Bring the RNA Adaptor tubes and the Stop Ligation Buffer mix to room temperature. For each, centrifuge at 600 x g for 5 s. Leave ligation mix in freezer until ready for use. Sample pooling arrangement for indexing should be known.

7.6.2 Add 2.5 µL of resuspension buffer and 2.5 µL of ligation mix to each sample tube. Now pipette in 2.5 µL of the proper RNA Adaptor into each sample tube. Adaptor selection should be performed by manufacturer's instructions for the chosen kit.

7.6.3 Recap and mix by centrifugation for 1 min at 280 x g. incubate in a thermal cycler for 10 min at 30 °C. Add 5 µL of Stop Ligation Buffer to sample to halt the reaction and mix well.

7.6.4 Begin clean-up by mixing SPRI paramagnetic beads by vortex for 60 s. Transfer 42 µL of beads to each tube. Mix thoroughly then incubate for 15 min on bench top.

7.6.5 Move tubes to magnetic stand and leave until liquid turns clear (~5 min). Once clear remove and discard 79.5 µL of supernatant. Leave tubes on the magnetic stand and wash by adding 200 µL of 80% ethanol. Allow to sit for at least 30 s then remove and discard supernatant. Repeat wash step. Leave on magnetic stand for additional ~15 min to allow for drying. Do not disrupt beads during wash steps.

7.6.6 Add 52.5 µL of resuspension buffer to the sample tube and mix until the beads are completely re-suspended. Incubate for 2 min on bench top. Move sample tube back to magnetic stand until liquid is clear (~5 min).

7.6.7 Leave on stand and gently pipette 50 µL of supernatant to new (6th) tubes. Add 50 µL of vortexed SPRI beads. Allow to incubate on bench top for 15 min.

7.6.8 Move to magnetic stand and let plate stay until liquid turns clear (~5 min). Discard 95 µL of supernatant leaving 5 µL in each tube.

7.6.9 Leave tubes on magnetic stand and wash by adding 200 µL of 80% ethanol. Allow tube to sit for 30 s then remove and discard supernatant. Repeat wash. Let tubes dry on magnetic stand for 15 min.

7.6.10 Add in 22.5 μL of resuspension buffer and mix until beads are suspended. Incubate on bench for 2 min then move to magnetic stand until liquid turns clear (~5 min). Pipette 20 μL of sample into new PCR tubes (7th).

7.7 Bring needed kit components to room temperature. Transfer 5 μL of PCR primer mix and 25 μL of PCR master mix from kit to sample (7th PCR tube) containing indexed samples. Mix sample and cap tube.

7.7.1 Place tubes into thermal cycler with pre-heated lid. Run program: 98 °C for 30 s, then 15 cycles of 98 °C for 10 s, 60 °C for 30 s, 72 °C for 30 s, 72 °C for 5 min. Hold at 4 °C.

7.8 Clean up samples by adding 50 μL of sufficiently mixed SPRI beads to sample tube and pipette to combine. Allow the reaction to incubate on the bench top for 15 min.

7.8.1 Move the tube to a magnetic stand and incubate until liquid turns clear (~5 min). Discard 95 μL of supernatant leaving 5 μL in each tube.

7.8.2 Leave tube on magnetic stand and wash by adding 200 μL of 80% ethanol. Allow to sit for at least 30 s and then remove and discard the supernatant. Repeat this wash step. Leave on magnetic stand for 15 min to allow to dry.

7.8.3 Add in 32.5 μL of resuspension buffer and mix until beads are completely re-suspended. Incubate on bench for 2 min then move to magnetic stand until liquid turns clear (~5 min). Pipette 30 μL of sample from each tube into a new PCR tube.

7.8.4 Assess DNA quantity and quality by repeating steps 6.1 and 6.2 with the appropriate chip. Move on to pooling step (step 9).

8. Small RNA Library Preparation for the sncRNAs.¹⁷

Note: Protocol steps based on manufacturer's instructions¹⁷.

8.1 Start small RNA library prep with ~220 ng - ~1.1 $\mu\text{g}/\mu\text{L}$ of globin-reduced RNA samples (4 μL).

8.2 Adaptor ligation

8.2.1 Ligate the 3'-adaptor by mixing thoroughly 1 μL of adaptor, 2 μL of nuclease-free water, and 4 μL of globin-reduced RNA sample. Incubate at 70 °C for 2 min and then immediately place on ice.

8.2.2 Add 10 μL of 2x ligation buffer, and 3 μL of ligation enzyme mix, mix and incubate 16 °C for 18 h.

Note: Increase sample incubation time to 18 hours at a lower temperature of 16 °C is recommended for studies interested in methylated RNA species, such as piwi RNAs. The longer incubation time allows for greater ligation efficiency due to the class of modifications during the 3' adaptor ligation phase.

8.2.3 Add 1 µL of reverse transcription primer in 4.5 µL of nuclease-free water to the ligation mixture to prevent adaptor-dimer formation of excess 3' adaptor.

8.2.4 Incubate in a preheated thermal cycler programmed for 5 min at 75 °C, 15 min at 37 °C, 15 min at 25 °C, and hold at 4°C.

8.2.5 Denature 1 µL of pre-diluted 5'-adaptor per sample in a thermal cycler at 70 °C for 2 min and then immediately place on ice.

8.2.6 Add 1 µL of denatured 5'-adaptor, 1 µL of 10x 5' ligation buffer and 2.5 µL of 5' ligation enzyme. Mix and incubate 25 °C for 1 h.

8.3 cDNA synthesis and amplification

8.3.1 Mix thoroughly by pipetting 30 µL of 5'/3'-adapter-ligated RNA, 8 µL of first-strand buffer, 1 µL of RNase inhibitor and 1 µL of reverse transcriptase. Incubate the mixture at 50 °C for 60 min. Proceed immediately to PCR amplification or heat inactivate the reaction at 70 °C for 15 min and store at -20 °C.

8.3.2 PCR amplify the 40 µL of Reverse Transcriptase reaction by adding 50 µL of PCR master mix, 2.5 µL of RNA PCR primer, add 2.5 µL of designated RNA PCR Primer Index to sample, and nuclease-free water to a total volume of 100 µL. Mix by pipetting. Run PCR in the thermal cycler using the following cycling conditions: initial denaturation at 94 °C for 30 s; 11 cycles of 94 °C for 15 s, annealing at 62 °C for 30 s, and extension at 70 °C for 15 s; followed by a final extension at 70 °C for 5 min; samples can be held at 4 °C in the cycler.

Note: Indexes are added for purpose of sample pooling.

8.4 Sample clean-up

8.4.1 From the DNA cleanup kit add 500 µL of binding buffer (5 M Gu-HCl, 30% isopropanol) to the 100 µL of PCR amplified sample to enable efficient binding to the spin-column membrane.

Note: If the binding buffer has pH indicator included check that the color of the mixture is yellow.

8.4.2 Pipette 600 µL mixture of sample and binding buffer into a filter cartridge inside a 2 mL collection tube, spin for 30-60 s at 17,900 x g in a tabletop centrifuge. Discard flow-through.

8.4.3 Wash the filter cartridge into the same collection tube with 0.75 mL of kit supplied wash buffer (10 mM Tris-HCl pH 7.5, 80% ethanol) by spinning for 30-60 s at 17,900 x g in a tabletop centrifuge and discarding the flow-through. Spin the filter cartridge dry for an additional 60 s.

8.4.4 Place the filter cartridge in a clean 1.5 mL collection tube. Add 30 µL of elution buffer (10 mM Tris-HCl pH 8.5), let the column stand for 60 s, and then spin for 60 s at 17,900 x g in a tabletop centrifuge.

8.5 Assess sample quality by repeating steps 6.1 and 6.2 with the appropriate DNA chip. Move to pooling step (step 9).

9. Sample pooling for sequencing

9.1 Pool barcoded and QC'ed samples by setting up and labeling new tubes (or a 96 well PCR plate) to contain either the mRNA or ncRNA samples. Transfer 13 µL of each barcoded 10nM library to corresponding tubes (or wells in the new plate) per manufacturer's instructions^{16,17}. Submit pooled samples for sequencing, be certain to choose similar read lengths if samples are to be run on the same chip.

REPRESENTATIVE RESULTS:

The representative samples in our study are the globin and ribo-depleted whole blood samples. The representative outcome of the protocol consists of a globin depleted library sample with an RNA integrity number (RIN) above 7 (**Figure 1a**) and 260/280 nm concentration ratios at or above 2 (**Figure 1b and 1c**). Validation of the sample outcome was performed using spectrophotometer to give the final concentration of each library and chip-based electrophoresis to give RIN number along with a graph of peaks that show which molecules (mRNA or ncRNA) were captured based on insert size in the library sample prior to pooling and sequencing. For the current study focus was placed on the mRNA and small ncRNAs only. For the mRNA libraries the representative result is an electropherogram peak at ~280 bp (**Figure 1d**). For the small ncRNAs representative results consist of a range of peaks from ~100-400 bp (**Figure 1e**), with peaks at ~143 and ~153 bp correspond to miRNAs and piRNAs, respectively. Our sample results showed that our optimized technique resulted in libraries with RIN numbers that ranged from 6.5 (sub-optimal) to 9.2 (above optimal). This proved to be an improvement compared to other studies that used globin depletion methods and were only able to achieve RIN numbers at or near 6. The chip-based electrophoresis results also showed that from one single blood sample it is possible to achieve RNA molecule capture of peaks representing both mRNA and ncRNA (**Figures 1d and 1e**), and sample insert sizes that covered both small and large non-coding RNA molecules. These results are representative of the optimal RIN scores and insert sizes needed to ensure quality transcript reads can be sequenced from the prepared libraries for nearly all NGS RNA-seq platforms.

FIGURE LEGENDS:

Figure 1: Representative results for mRNA and non-coding RNA expression libraries from single samples of porcine whole blood. (a) Electrophoresis file run summary representation of

RIN numbers pre-globin reduction and post-globin reduction. Representative electropherograms (b) pre-globin reduction and (c) post-globin reduction of a single porcine whole blood sample. (d) Representative electropherograms of globin and ribo-depleted whole blood mRNA library samples prior to pooling and sequencing. (e) Representative electropherograms of globin depleted sncRNA whole blood libraries of the same samples featured in panel d prior to pooling and sequencing

DISCUSSION:

The first critical step in the protocol that made it optimized included the added globin depletion steps, which made it possible to get quality reads from whole blood samples. One of the largest limitations on using whole blood in sequencing studies are the high numbers of reads in the sample that will map to globin molecules and reduce the reads that could map to other molecules of interest¹⁸. Therefore, in optimizing the protocol for our sample type, we needed to incorporate a globin depletion step to ensure the highest possible mRNA and non-coding RNA capture through sequencing. All of the samples were globin depleted to account for high levels of globin transcripts using porcine specific hemoglobin A and B (HBA and HBB) oligonucleotides based on the procedure from Choi *et al.*, 2014¹⁴. Another critical step was the use of a ribo-depletion extraction kit which allowed for the ability remove unwanted RNA molecules from our samples while simultaneously retaining polyadenylated non-coding and viral RNA molecules of experimental interest within our libraries for sequencing¹⁹. Additionally, we were able to work with single samples to create both coding (mRNA) and non-coding (small and long) RNA libraries by removing the small RNA enrichment and size selection steps. By doing this, we were able to maximize the RNA aliquots for pooling and to guarantee we sent enough volume to allow for sequencing. Optimization of all of the manufacturer's protocols were done to increase mRNA and non-coding RNA recovery for downstream library creation. We also optimized the non-coding RNA library preparation for the small non-coding RNA portion of our downstream analysis by adding additional time to the PCR incubation per the manufacturers' instructions to increase the ligation efficiency of methylated RNAs and the likelihood of recovering piwi-RNAs.

The limitations of the optimized protocol are rooted in the lack of specificity that also makes it ideal for multiple analyses from a single sample. By removing the enrichment and size selection portions of the non-coding RNA library preparation, we are limiting novel small RNA discovery²⁰ to balance it against acquiring both small and long non-coding RNA molecules. Hence, the optimized protocol gives a good overview of the types of non-coding expression but, requires extra processing during the analytical stage of mapping to gain more specific information on the classes and sizes of non-coding RNA expression captured in the sample. Another limitation to the method is due to the incorporation of the globin depletion steps, which will lower the overall RIN numbers. We were able to troubleshoot this by first examining the RIN numbers before and after globin depletion to understand how large of a reduction in RNA integrity we would experience. This yielded results that showed we could experience a drop of ~1-2 points. To account for this drop we further optimized the globin depletion protocol to 6 µg of sample instead of the recommended 10 µg. This helped to improve our RIN numbers, as well as, conserve sample. Additionally, we also used thin-walled microamp reaction tubes and iced the

tubes immediately after the first denaturing step. This allowed us to quench the reaction quicker, allowing for improved RIN numbers on the globin depleted samples.

The modified protocol used in this study has several advantages over the base library preparation methods used in other studies where only mRNA or non-coding RNA are studied separately. The changes we made for optimization allowed for efficient use of whole blood as our sample type. This is advantageous to future studies because whole blood as a sample type has less steps for collection and processing than the leukocyte portion of blood. Additionally, whole blood also has the added advantages of allowing researchers to examine systemic responses and can be repeatedly collected from the animal. However, there is a caveat to the use of whole blood samples in that the amount of blood collected is dependent on the age and size of the organism in question. Smaller amounts of blood will lead to lower RNA yields, however the method presented here will allow dual library creation from at least 2.5 mL of whole blood. This allowed for us to use one sample to investigate both mRNA and non-coding expression and correlate both to an individual at a snapshot in time, effectively allowing us to collect more information than traditional library preparation methods. Also, by performing the ribo-depletion we made it possible to reduce transcripts that could deplete sequencing reads, while still retaining the non-coding and viral transcripts that can be lost during traditional library creation to study mRNA. In this way our optimized protocol triples the amount of information that can be gained from a single sample. Other significant changes we made to the method to facilitate the capture of more RNA species information were: changing the type of PCR plate to quickly stop the reaction, which improved RNA purity yield; a longer thermocycler incubation to increase ligation for possible recovery of methylated RNAs; and not employing a size selection gel, to allow for all detectable ncRNAs to be identified regardless of length. For the downstream analysis this allowed for the capture of multiple non-coding RNAs between 18nt-200nt in length.

By employing this optimized method, our group was able to put forward a protocol that can be applied to whole blood transcriptomic analyses and allow for mRNA and non-coding RNA from a single sample.

ACKNOWLEDGMENTS:

This work was mainly supported by the by USDA NIFA AFRI 2013-67015-21236, and in part by USDA NIFA AFRI 2015-67015-23216. This study was supported in part by an appointment to the Agricultural Research Service Research Participation Program administered by the Oak Ridge Institute for Science and Education (ORISE) through an interagency agreement between the US Department of Energy (DOE) and the US Department of Agriculture. ORISE is managed by Oak Ridge Associated Universities under DOE contract no. DE-AC05-06OR2310.

We would like to thank Dr. Kay Faaberg for the HP-PRRSV infectious clones, Dr. Susan Brockmeier for her help with animals involved in the experiment, and Sue Ohlendorf for secretarial assistance in preparation of the manuscript.

DISCLOSURES:

The authors have nothing to disclose.

Mention of trade names or commercial products in this article is solely for the purpose of providing specific information and does not imply recommendation or endorsement by the U.S. Department of Agriculture. USDA is an equal opportunity provider and employer.

REFERENCES:


- 1 Finotello, F. & Di Camillo, B. Measuring differential gene expression with RNA-seq: challenges and strategies for data analysis. *Briefings in Functional Genomics*. **14** (2), 130-142, 10.1093/bfpg/elu035 (2015).
- 2 Coble, D. J. *et al.* RNA-seq analysis of broiler liver transcriptome reveals novel responses to high ambient temperature. *BMC Genomics*. **15** 1084, 10.1186/1471-2164-15-1084 (2014).
- 3 Koltes, J. E. *et al.* Identification of a putative quantitative trait nucleotide in guanylate binding protein 5 for host response to PRRS virus infection. *BMC Genomics*. **16** 412, 10.1186/s12864-015-1635-9 (2015).
- 4 Miller, L. C. *et al.* Comparative analysis of signature genes in PRRSV-infected porcine monocyte-derived cells to different stimuli. *PLoS One*. **12** (7), e0181256, 10.1371/journal.pone.0181256 (2017).
- 5 Head, S. R. *et al.* Library construction for next-generation sequencing: overviews and challenges. *Biotechniques*. **56** (2), 61-64, 66, 68, passim, 10.2144/000114133 (2014).
- 6 Wang, Z., Gerstein, M. & Snyder, M. RNA-Seq: a revolutionary tool for transcriptomics. *Nature reviews. Genetics*. **10** (1), 57-63, 10.1038/nrg2484 (2009).
- 7 Dominic, O. N., Heike, G. & Martin, S. Ribosomal RNA Depletion for Efficient Use of RNA-Seq Capacity. *Current Protocols in Molecular Biology*. **103** (1), 4.19.11-14.19.18, 10.1002/0471142727.mb0419s103 (2013).
- 8 Pollet, S. *FR-AgEncode: a French pilot project to enrich the annotation of livestock genomes- RNA Extraction Protocol*, <ftp://ftp.faang.ebi.ac.uk/ftp/protocols/assays/INRA_SOP_RNA_extraction_20160504.pdf> (2015).
- 9 Fleming, D. S. & Miller, L. C. Identification of small non-coding RNA classes expressed in swine whole blood during HP-PRRSV infection. *Virology*. **517** 56-61, 10.1016/j.virol.2018.01.027 (2018).
- 10 Fleming, D. S. & Miller, L. C. Small non-coding RNA expression status in animals faced with highly pathogenic porcine reproductive and respiratory syndrome virus (HP-PRRSV). *Proceedings of the World Congress on Genetics Applied to Livestock Production. Species - Porcine 2* (11) (2018).
- 11 Bivens Nathan, J. & Zhou, M. RNA-Seq Library Construction Methods for Transcriptome Analysis. *Current Protocols in Plant Biology*. **1** (1), 197-215, 10.1002/cppb.20019 (2016).
- 12 Cui, P. *et al.* A comparison between ribo-minus RNA-sequencing and polyA-selected RNA-sequencing. *Genomics*. **96** (5), 259-265, <https://doi.org/10.1016/j.ygeno.2010.07.010> (2010).
- 13 Guo, Y. *et al.* RNAseq by Total RNA Library Identifies Additional RNAs Compared to Poly(A) RNA Library. *BioMed Research International*. **2015** 9, 10.1155/2015/862130 (2015).
- 14 Choi, I. *et al.* Increasing gene discovery and coverage using RNA-seq of globin RNA reduced porcine blood samples. *BMC Genomics*. **15** 954, 10.1186/1471-2164-15-954 (2014).

- 15 Wu, K., Miyada, G., Martin, J. & Finkelstein, D. *Globin reduction protocol: A method for*
16 *processing whole blood RNA samples for improved array results.*,
17 <http://media.affymetrix.com:80/support/technical/technotes/blood2_technote.pdf> (2007).
18 Illumina. TruSeq® Stranded Total RNA Sample Preparation Guide.
19 New England Biolabs Inc. *Protocol for use with NEBNext Multiplex Small RNA Library*
20 *Prep Kit for Illumina (Index Primers 1-48) (E7560)*, <[https://www.neb.com/~media/Catalog/All-](https://www.neb.com/~media/Catalog/All-Protocols/758F75A29CDE4D03954E1EF75E78EA3D/Content/manualE7560_Figure_1.jpg)
21 [Protocols/758F75A29CDE4D03954E1EF75E78EA3D/Content/manualE7560_Figure_1.jpg](https://www.neb.com/~media/Catalog/All-Protocols/758F75A29CDE4D03954E1EF75E78EA3D/Content/manualE7560_Figure_1.jpg)> (
22 Shin, H. *et al.* Variation in RNA-Seq transcriptome profiles of peripheral whole blood
23 from healthy individuals with and without globin depletion. *PLoS One*. **9** (3), e91041,
24 10.1371/journal.pone.0091041 (2014).
25 Costa, V., Angelini, C., De Feis, I. & Ciccodicola, A. Uncovering the complexity of
26 transcriptomes with RNA-Seq. *Journal of Biomedicine and Biotechnology*. **2010** 853916,
27 10.1155/2010/853916 (2010).
28 Martens-Uzunova, E. S., Olvedy, M. & Jenster, G. Beyond microRNA--novel RNAs derived
29 from small non-coding RNA and their implication in cancer. *Cancer Letters*. **340** (2), 201-211,
30 10.1016/j.canlet.2012.11.058 (2013).


a.

Sample Name	Sample Comment	Status	Result Label	Result Color
Pre-globin reduction		✓	RIN: 8.90	
Pre-globin reduction		✓	RIN: 9.10	
Pre-globin reduction		✓	RIN: 9.50	
Pre-globin reduction		✓	RIN: 8.70	
Post-globin reduction		✓	RIN: 8.40	
Post-globin reduction		✓	RIN: 7.20	
Post-globin reduction		✓	RIN: 7.20	
Post-globin reduction		✓	RIN: 6.80	
Post-globin reduction		✓	RIN: 8.10	
Post-globin reduction		✓	RIN: 7.60	
Post-globin reduction		✓	RIN: 6.30	
Post-globin reduction		✓	RIN: 7.10	
Ladder		✓	All Other Samples	

b.



c.



RNA Area:

219.8

RNA Integrity Number (RIN):

9.2 (8.02.08)

RNA Concentration:

286 ng/ul

Result Flagging Color:


rRNA Ratio (28s / 18s):


1.4

Result Flagging Label:

RIN: 9.20

Name	Start Size [nt]	End Size [nt]	Area	% of total Area
18S	677	1,135	43.0	19.6
28S	1,536	3,937	61.7	28.1

d.

e.

Materials	Company	Catalog Number	Comments
PAXgene Tubes	PreAnalytix	762165	
Molecular Biology Grade Water	ThermoFisher	10977-015	
mirVana miRNA Isolation Kit	ThermoFisher	AM1560	
Rneasy MinElute Clean Up Kit	QIAGEN	74204	
100% Ethanol	Decon Labs, Inc.	2716	
0.2 mL thin-walled tubes	ThermoFisher	98010540	
1.5 mL RNase/DNase - free tubes	Any supplier		
Veriti 96-well Thermocycler	ThermoFisher	4375786R	
Globin Reduction Oligo (α 1)	Any supplier		Sequence GAT CTC CGA GGC TCC AGC TTA ACG GT
Globin Reduction Oligo (α 2)	Any supplier		Sequence TCA ACG ATC AGG AGG TCA GGG TGC AA
Globin Reduction Oligo (β 1)	Any supplier		Sequence AGG GGA ACT TAG TGG TAC TTG TGG GT
Globin Reduction Oligo (β 2)	Any supplier		Sequence GGT TCA GAG GAA AAA GGG CTC CTC CT
10X Oligo Hybridization Buffer			
-Tris-HCl, pH 7.6	Fisher Scientific	BP1757-100	
-KCl	Millipore Sigma	60142-100ML-F	
10X RNase H Buffer			
-Tris-HCl, pH 7.6	Fisher Scientific	BP1757-100	
-DTT	ThermoFisher	Y00147	
-MgCl ₂	Promega	A351B	
-Molecular Biology Grade Water	ThermoFisher	10977-015	
RNase H	ThermoFisher	AM2292	
SUPERase-IN	ThermoFisher	AM2694	Rnase inhibitor
EDTA	Millipore Sigma	E7889	
Microcentrifuge	Any supplier		
2100 Electrophoresis BioAnalyzer Instrument	Agilent Technologies	G2938C	
Agilent RNA 6000 Nano Kit	Agilent Technologies	5067-1511	
Agilent High Sensitivity DNA Kit	Agilent Technologies	5067-4626	
TruSeq Stranded Total RNA Library Prep Kit with Ribo-Zero	Illumina	RS-122-2201	mRNA kit; Human/Mouse/Rat Set A (48 samples, 12 index) Available on-line
TruSeq Stranded Total RNA Sample Preparation Guide	Illumina		
RNAClean XP Beads	BeckmanCoulter	A63987	
AMPure XP Beads	BeckmanCoulter	A63880	
MicroAmp Optical 8-tube Strip	ThermoFisher	N8010580	0.2 ml thin-walled tubes
MicroAmp Optical 8-tube Strip Cap	ThermoFisher	N801-0535	
RNase/DNase - free reagent reservoirs	Any supplier		
SuperScript II Reverse Transcriptase	ThermoFisher	18064-014	
MicroAmp Optical 96 well plates	ThermoFisher	N8010560	These were used in place of .3mL plates as needed
MicroAmp Optical adhesive film	ThermoFisher	4311971	
NEBNext Multiplex Small RNA Library Prep Set for Illumina® (Set 1)	New England Biolabs	E73005	small RNA kit
NEBNext Multiplex Small RNA Library Prep Set for Illumina® (Set 2)	New England Biolabs	E75805	small RNA kit
QIAQuick PCR Purification Kit	QIAGEN	28104	
96S Super Magnet Plate	ALPAQUA	A001322	



1 Alewife Center #200
Cambridge, MA 02140
tel. 617.945.9051
www.jove.com

ARTICLE AND VIDEO LICENSE AGREEMENT

Title of Article:

Identification of coding and non-coding RNA classes expressed in swine whole blood

Author(s):

Damarius S. Fleming, Sarah J. Anderson, Laura C. Miller

Item 1: The Author elects to have the Materials be made available (as described at <http://www.jove.com/publish>) via:

☒

Standard Access

☐

Open Access

Item 2: Please select one of the following items:

☐

The Author is **NOT** a United States government employee.

☒

The Author is a United States government employee and the Materials were prepared in the course of his or her duties as a United States government employee.

☐

The Author is a United States government employee but the Materials were NOT prepared in the course of his or her duties as a United States government employee.

ARTICLE AND VIDEO LICENSE AGREEMENT

1. **Defined Terms.** As used in this Article and Video License Agreement, the following terms shall have the following meanings: **"Agreement"** means this Article and Video License Agreement; **"Article"** means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein; **"Author"** means the author who is a signatory to this Agreement; **"Collective Work"** means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; **"CRC License"** means the Creative Commons Attribution-Non Commercial-No Derivs 3.0 Unported Agreement, the terms and conditions of which can be found at: <http://creativecommons.org/licenses/by-nc-nd/3.0/legalcode>; **"Derivative Work"** means a work based upon the Materials or upon the Materials and other pre-existing works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; **"Institution"** means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; **"JoVE"** means MyJoVE Corporation, a Massachusetts corporation and the publisher of The Journal of Visualized Experiments; **"Materials"** means the Article and / or the Video; **"Parties"** means the Author and JoVE; **"Video"** means any video(s) made by the Author, alone or in conjunction with any other parties, or by JoVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion

of the Article, and in which the Author may or may not appear.

2. **Background.** The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JoVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.

3. **Grant of Rights in Article.** In consideration of JoVE agreeing to publish the Article, the Author hereby grants to JoVE, subject to **Sections 4** and **7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the "Open Access" box has been checked in **Item 1** above, JoVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.

ARTICLE AND VIDEO LICENSE AGREEMENT

4. **Retention of Rights in Article.** Notwithstanding the exclusive license granted to JoVE in **Section 3** above, the Author shall, with respect to the Article, retain the non-exclusive right to use all or part of the Article for the non-commercial purpose of giving lectures, presentations or teaching classes, and to post a copy of the Article on the Institution's website or the Author's personal website, in each case provided that a link to the Article on the JoVE website is provided and notice of JoVE's copyright in the Article is included. All non-copyright intellectual property rights in and to the Article, such as patent rights, shall remain with the Author.

5. **Grant of Rights in Video – Standard Access.** This **Section 5** applies if the "Standard Access" box has been checked in **Item 1** above or if no box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby acknowledges and agrees that, Subject to **Section 7** below, JoVE is and shall be the sole and exclusive owner of all rights of any nature, including, without limitation, all copyrights, in and to the Video. To the extent that, by law, the Author is deemed, now or at any time in the future, to have any rights of any nature in or to the Video, the Author hereby disclaims all such rights and transfers all such rights to JoVE.

6. **Grant of Rights in Video – Open Access.** This **Section 6** applies only if the "Open Access" box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby grants to JoVE, subject to **Section 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Video in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Video into other languages, create adaptations, summaries or extracts of the Video or other Derivative Works or Collective Works based on all or any portion of the Video and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. For any Video to which this **Section 6** is applicable, JoVE and the Author hereby grant to the public all such rights in the Video as provided in, but subject to all limitations and requirements set forth in, the CRC License.

7. **Government Employees.** If the Author is a United States government employee and the Article was prepared in the course of his or her duties as a United States government employee, as indicated in **Item 2** above, and any of the licenses or grants granted by the Author hereunder exceed the scope of the 17 U.S.C. 403, then the rights granted hereunder shall be limited to the maximum

rights permitted under such statute. In such case, all provisions contained herein that are not in conflict with such statute shall remain in full force and effect, and all provisions contained herein that do so conflict shall be deemed to be amended so as to provide to JoVE the maximum rights permissible within such statute.

8. **Protection of the Work.** The Author(s) authorize JoVE to take steps in the Author(s) name and on their behalf if JoVE believes some third party could be infringing or might infringe the copyright of either the Author's Article and/or Video.

9. **Likeness, Privacy, Personality.** The Author hereby grants JoVE the right to use the Author's name, voice, likeness, picture, photograph, image, biography and performance in any way, commercial or otherwise, in connection with the Materials and the sale, promotion and distribution thereof. The Author hereby waives any and all rights he or she may have, relating to his or her appearance in the Video or otherwise relating to the Materials, under all applicable privacy, likeness, personality or similar laws.

10. **Author Warranties.** The Author represents and warrants that the Article is original, that it has not been published, that the copyright interest is owned by the Author (or, if more than one author is listed at the beginning of this Agreement, by such authors collectively) and has not been assigned, licensed, or otherwise transferred to any other party. The Author represents and warrants that the author(s) listed at the top of this Agreement are the only authors of the Materials. If more than one author is listed at the top of this Agreement and if any such author has not entered into a separate Article and Video License Agreement with JoVE relating to the Materials, the Author represents and warrants that the Author has been authorized by each of the other such authors to execute this Agreement on his or her behalf and to bind him or her with respect to the terms of this Agreement as if each of them had been a party hereto as an Author. The Author warrants that the use, reproduction, distribution, public or private performance or display, and/or modification of all or any portion of the Materials does not and will not violate, infringe and/or misappropriate the patent, trademark, intellectual property or other rights of any third party. The Author represents and warrants that it has and will continue to comply with all government, institutional and other regulations, including, without limitation all institutional, laboratory, hospital, ethical, human and animal treatment, privacy, and all other rules, regulations, laws, procedures or guidelines, applicable to the Materials, and that all research involving human and animal subjects has been approved by the Author's relevant institutional review board.

11. **JoVE Discretion.** If the Author requests the assistance of JoVE in producing the Video in the Author's facility, the Author shall ensure that the presence of JoVE employees, agents or independent contractors is in accordance with the relevant regulations of the Author's institution. If more than one author is listed at the beginning of this Agreement, JoVE may, in its sole

ARTICLE AND VIDEO LICENSE AGREEMENT

discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JoVE reserves the right, in its absolute and sole discretion and without giving any reason therefore, to accept or decline any work submitted to JoVE. JoVE and its employees, agents and independent contractors shall have full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JoVE has sole discretion as to the method of making and publishing the Materials, including, without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

12. **Indemnification.** The Author agrees to indemnify JoVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JoVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JoVE, making of videos by JoVE, or publication in JoVE or elsewhere by JoVE. The Author shall be responsible for, and shall hold JoVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to

the making of a video by JoVE its employees, agents or independent contractors. All sterilization, cleanliness or decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's expense. All indemnifications provided herein shall include JoVE's attorney's fees and costs related to said losses or damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JoVE, its employees, agents or independent contractors.

13. **Fees.** To cover the cost incurred for publication, JoVE must receive payment before production and publication the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JoVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.

14. **Transfer, Governing Law.** This Agreement may be assigned by JoVE and shall inure to the benefits of any of JoVE's successors and assignees. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to be one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement is required per submission.

CORRESPONDING AUTHOR

Name:

Laura C. Miller

Department:

Virus and Prion Research Unit

Institution:

National Animal Disease Center, USDA, Agricultural Research Service

Title:

Research Microbiologist

Signature:

Laura Miller

Date:

06/26/2018

Please submit a **signed** and **dated** copy of this license by one of the following three methods:

1. Upload an electronic version on the JoVE submission site
2. Fax the document to +1.866.381.2236
3. Mail the document to JoVE / Attn: JoVE Editorial / 1 Alewife Center #200 / Cambridge, MA 02140

JoVE58689R1: Identification of coding and non-coding RNA classes expressed in swine whole blood

Editorial comments:

1. The editor has formatted the manuscript to match the journal's style. Please retain the same.
The formatting of the manuscript by the editor has been retained.

2. Please address all the specific comments marked in the manuscript.
All the specific comments marked in the manuscript have been addressed and are included in the Rebuttal Letter_Replies document following the Rebuttal Letter.

3. Please ensure that the highlight steps protocol reflects the title of the manuscript.
Highlighted text now reflects title. Non-highlighted protocol sections are performed per manufacturer's instructions.

4. After formatting please ensure that the protocol is no more than 10 pages long and the highlight is no more than 2.75 pages both including the headings and spacing.
A clean copy of the revised manuscript after formatting is included. The protocol itself is 9 pages total and the highlighted portion is ~2.75 pages.

5. Figures: Please either combine the panels or number each panel as individual figure. Also please remove the legend from the figure and provide a detail description of each figure in the figure legend. Each Figure Legend should include a title and a short description of the data presented in the Figure and relevant symbols.
The figures have been combined into a single Figure 1 with 5 panels a, b, c, d, and e. Legend has been removed from the figure itself and the Figure Legend includes a title and a short description of the data presented in the Figure panels.

TITLE:

Identification of Coding and Non-coding RNA Classes Expressed in Swine Whole Blood

AUTHORS AND AFFILIATIONS:

Damarius S. Fleming^{1,2}, Sarah J. Anderson¹, Laura C. Miller¹

¹Virus and Prion Research Unit, National Animal Disease Center, USDA, Agricultural Research Service, Ames, IA, USA

²ORAU/ORISE Oak Ridge, TN, USA

Corresponding Author:

Laura C. Miller (laura.miller@ars.usda.gov)

Tel: (515)-337-6794

Email Addresses of Co-authors:

Damarius Fleming (damarius.fleming@ars.usda.gov)

Sarah Anderson (sarah.anderson@ars.usda.gov)

KEYWORDS:

mRNA, Non-coding RNA, RNA-Seq, Whole blood, Epigenetic, Globin reduction

SUMMARY:

Here, we present a protocol optimized for the processing of coding (mRNA) and non-coding RNAs (ncRNA) globin reduced RNA-seq libraries from a single whole blood sample.

ABSTRACT:

The advent of innovative and increasingly powerful next generation sequencing techniques has opened new avenues into the ability to examine the underlying gene expression related to biological processes of interest. These innovations not only allow researchers to observe expression from the mRNA sequences that code for genes that effect cellular function, but also the non-coding RNA (ncRNA) molecules that remain untranslated, but still have regulatory functions. Although researchers have the ability to observe both mRNA and ncRNA expression, it has been customary for a study to focus on one or the other. However, when studies are interested in both mRNA and ncRNA expression, many times they use separate samples to examine either coding or non-coding RNAs due to the difference in library preparations. This can lead to the need for more samples which can increase time, consumables, and animal stress. Additionally, it may cause researchers to decide to prepare samples for only one analysis, usually the mRNA, limiting the number of biological questions that can be investigated. However, ncRNAs span multiple classes with regulatory roles that effect mRNA expression. Because ncRNA are important to fundamental biologic processes and disorder of these processes in during infection, they may, therefore, make attractive targets for therapeutics. This manuscript demonstrates a modified protocol for the generation mRNA and non-coding RNA expression libraries, including viral RNA, from a single sample of whole blood. Optimization of this protocol, improved RNA purity, increased ligation for recovery of methylated RNAs, and omitting size selection, to allow capture of more RNA species.

Commented [A1]: The highlighted protocol only shows how to prepare, purify and sequence coding RNA. Non-coding RNA prep is not highlighted. Please ensure that the highlight is in alignment with the title of the protocol.

Also please remember that no more than 2.75 pages can be highlighted. This includes heading and spacings as well.

Commented [DSF2]: Response to Editor ques #1: Highlighted text now reflects title. Non-highlighted protocol sections are performed per manufacturer's instructions.

INTRODUCTION:

Next generation sequencing (NGS) has emerged as a powerful tool for the investigation of the changes that occur at the genomic level of biological organisms. Sample preparation for NGS methods can be varied depending on the organism, tissue type, and more importantly the questions the researchers are keen to address. Many studies have turned to NGS as a means of studying the differences in gene expression between states such as healthy and sick individuals¹⁻⁴. The sequencing take place on a whole genome basis and allows a researcher to capture the most, if not all, of the genomic information for a particular genetic marker at a time point.

The most common markers of expression observed are the messenger RNAs (mRNAs). The most used procedures for prepping libraries for RNA-seq are optimized for the recovery of mRNA molecules through the use of a series of purifications, fragmentations, and ligations^{5,6}. However, the decision on how a protocol is to be performed relies heavily on the sample type and the questions being posed about said sample. In most cases total RNA is extracted; yet, not all RNA molecules are of interest and in cases such as mRNA expression studies overly abundant RNA species, like ribosomal RNAs (rRNA) need to be removed to increase the number of detectable transcripts associated with the mRNAs. The most popular and widely used method for removing the abundant rRNA molecules is the reduction of polyadenylated RNA transcripts referred to as polyA depletion⁷. This approach works well for the analysis of mRNA expression as it does not affect the mRNA transcripts. However, in studies that are interested in non-coding or viral RNAs, polyA depletion also removes these molecules.

Many studies choose to focus on the RNA sequence library preparation to examine either mRNA expression (coding) or a particular class of small or large non-coding RNA. Although there are other procedures⁸ like ours that allow for the dual sample preparation, many studies prepare libraries from separate samples for separate studies when available. For a study like ours, this would normally require multiple blood samples increasing time, consumables, and animal stress. The goal of our study was to be able to use whole blood from animals to identify and quantify the different classes of both mRNA and non-coding RNA expressed between healthy and highly pathogenic porcine reproductive and respiratory syndrome virus (HP-PRRSV) challenged pigs^{9,10} despite having only a single whole blood sample (2.5 ml) from each pig. In order to do this, we needed to optimize the typical extraction and library creation protocols to generate the proper data to allow for analysis of both mRNA and non-coding RNA (ncRNA) expression¹¹ from a single sample.

This prompted a need for a protocol that allowed for mRNA and non-coding RNA analysis because the available standard kits and methods for RNA-extraction and library creation were intended chiefly for mRNA and use a poly-A depletion step¹². This step would have made it impossible to recover non-coding RNA or viral transcripts from the sample. Therefore, an optimized method was needed that allowed for total RNA extraction without sample polyA depletion. The method presented in this manuscript has been optimized to allow for the use of whole blood as a sample type and to build sequencing libraries for both mRNA and ncRNAs of small and large sizes. The method has been optimized to allow for the analysis of all detectable non-coding RNAs as well as retain viral RNAs for later investigation¹³. In all, our optimized library preparation protocol allows

Commented [A3]: Volume?

Commented [DSF4]: Response to Editor ques #2

for the investigation of multiple RNA molecules from a single whole blood sample.

The overall goal behind the use of this method was to develop a process that allowed for the collection of both non-coding RNA and mRNA from one sample of whole blood. This allows us to have mRNA, ncRNA, and viral RNA for each animal in our study sourced from a single sample⁹. This, ultimately, allows for more scientific discovery without additional animal costs and gives a more complete picture of the expression of each individual sample. The described method allows for the examination of the regulators of gene expression as well as allowing for completion of correlative studies comparing both mRNA and non-coding RNA expression using a single whole blood sample. Our study used this protocol to examine the changes in gene expression and possible epigenetic regulators in virally infected 9-week old male commercial pigs.

PROTOCOL:

Animal protocols were approved by the National Animal Disease Center (USDA-NADC-ARS) Animal Care and Use Committee.

1. Collection of Swine Blood Samples

1.1 Collect blood samples into RNA tubes. Collect ~2.5 mL or more if larger collection tubes are available.

2. Processing of Swine Blood Samples

2.1 Centrifuge the blood tubes at 5020 x g for 10 min at room temperature (15-25° C). If processing frozen samples incubate tube at room temperature for a **minimum** of 2 h prior to centrifugation.

2.2 Remove supernatant and add 8 mL RNase-free water to the pellet. Close and vortex the pellet until it is visibly dissolved. Centrifuge sample tubes at 5020 x g for 10 min at room temperature to recover the pellet. Discard all supernatant and preserve the pellet.

3. Organic Extraction for Total RNA and Small RNA (miRNA Isolation Kit)

3.1 Begin total RNA extraction by pipetting 300 µL of lysis binding buffer to the pellet from the step 2.2.

3.2 Vortex and transfer the mixture to a new labelled 1.5 mL centrifuge tube. Add 30 µL homogenate additive from kit. Vortex the tube and place on the ice for 10 min.

3.3 Remove the tube and add 300 µL of the Acid-Phenol: Chloroform reagent from the kit. Vortex the tube to mix. Centrifuge at 10,000 x g for 5 min at room temperature.

Commented [DSF5]: Response to Editor ques #3:
Institute added

Commented [A6]: Of which Institute?

132 3.4 Carefully remove aqueous phase (300-350 μ L) to a fresh tube. Note the volume for the next
133 step.

134
135 **4. Total RNA Isolation Procedure**

136
137 4.1 Based on the amount of aqueous recovery (300-350 μ L) add 1.25x volume of 100% (~375 μ L)
138 ethanol to aqueous phase. Mix the sample using a pipette.

139
140 4.2 Set-up new collection tubes containing a filter cartridge for each sample. Pipette ~ 675 μ L of
141 the lysate/ethanol mixture onto the filter cartridge.

142 **Note:** Do not add > 700 μ L to filter cartridge at one time. For larger volumes apply in succession.

143
144 4.3 Centrifuge briefly (~15-20 s) at 10,000 x g to pass liquid through the filter. Do not spin harder
145 than this.

146
147 4.4 Discard the flow through and repeat the centrifugation with the remaining lysate/ethanol
148 mixture until it has all been applied. Retain the same filter cartridge and collection tube for the
149 next step.

150
151 4.5 Add 700 μ L Wash Solution 1 from kit to filter cartridge and centrifuge briefly (~10 s) to pull
152 through the filter. Discard the flow through and retain the same filter cartridge and collection
153 tube.

154
155 4.6 Add 500 μ L Wash Solution 2/3. Centrifuge to draw the liquid through the filter cartridge.
156 Discard the flow through. Repeat wash step.

157
158 4.7 In the same tube, spin the filter cartridge an additional 60 s to remove any residual liquid
159 from filter. Transfer the filter cartridge to a fresh collection tube.

160
161 4.8 Add 100 μ L of pre-heated (95 $^{\circ}$ C) nuclease-free water to the center of the filter cartridge. Spin
162 for approximately 20-30 s at the tabletop centrifuge max speed.

163 **Note:** RNA is contained in the eluate and can now be further processed or stored at -20 $^{\circ}$ C or
164 below. Enrichment for small RNAs was not performed.

165
166 **5. Globin Reduction (based on a protocol optimized for porcine whole blood samples)^{14,15}**

167
168 **Note:** Globin reduction is performed so that libraries are not overpopulated with reads mapping
169 to globin genes, which would lower the number of reads available to map to other genes of
170 greater interest^{14,15}.

171
172 **5.1 Hybridization with globin reduction oligos**

173
174 5.1.1 Denature the RNA (6 μ g RNA sample in maximum 7 μ L volume) by pipetting the extracted
175 sample into a 0.2 ml thin-walled nuclease-free reaction tube and placing in a thermal cycler at 70

Commented [A7]: Include a note stating why globin reduction needs to be performed in one-two lines with citation.

Commented [DSF8]: Response to Editor ques #6
Inline response to editor comments on globin reduction.

Commented [DSF9]: Response to Editor ques #7:
Commercial term removed.

Commented [A10]: Please do not use commercial term.
Please use generic term instead.

176 °C for 2 min. It is key to ice tubes immediately after the first denature step for optimal RNA
177 quality. No DNase treatment is needed.

178
179 5.1.2 While tubes cool prepare a 400 µL 10X globin reduction oligo mix; 100 µL each of two HBA
180 oligos (5'-GATCTCCGAGGCTCCAGCTTAACGGT-3', and 5'-TCAACGATCAGGAGGTCAGGGTGCAA-3')
181 at 30 µM, two HBB oligos (5'-AGGGGAAGCTTAGTGGTACTTGTGGGC-3', and 5'-
182 GGTTCAAGAGGAAAAAGGGCTCCTCCT-3') at 120 µM per reaction, yielding a final concentration of
183 7.5 µM HBA oligos and 30 µM HBB oligos. Prepare 10X oligo hybridization buffer: 100 mM Tris-
184 HCL, pH 7.6; 200 mM KCl.

185
186 5.1.3 Prepare the hybridization mix: 6 µg RNA sample (maximum 7 µL volume), 2 µL of the 400
187 µL 10X globin reduction oligo mix (final concentration 2X), 1 µL 10X oligo hybridization buffer
188 (final concentration 1X). Add nuclease-free water to a final volume of 10 µL.

189
190 5.1.4 Set thermocycler at 70°C for 5 min. Cool immediately to 4 °C and proceed to RNase H
191 digestion.

192
193 5.2 RNase H digestion

194
195 5.2.1 Dilute 10X RNase H (10 U/ µL) to 1X RNase H with 1X RNase H buffer.

196 Note: The RNase buffer comes at 10X and will need to be diluted with 1X RNase H buffer to 1X
197 prior to use.

198
199 5.2.2 Prepare RNase H reaction mix by combining: 2 µL 10X RNase buffer, 1 µL RNase inhibitor In
200 2 µL 1X RNase H, 5 µL nuclease-free water to a total volume 10 µL.

201
202 5.2.3 Mix thoroughly the Globin Reduction hybridization samples with 10 µL of the RNase H
203 reaction mix and digest at 37°C for 10 min. Cool to 4°C.

204
205 5.2.4 Stop digestion by addition of 1.0 µL 0.5 M EDTA to each sample and proceed immediately
206 to the cleanup step.

207
208 5.3 RNase H-Treated Total RNA Cleanup.

209
210 5.3.1 Purify the RNase H treated RNA using a silica-membrane-based elution purification cleanup
211 kit according to manufacturer's instructions. Premix buffers: For the mild washing buffer, add 44
212 ml of 100% ethanol.

213
214 5.3.2 Do not transfer the sample to a new tube. Add 80 µL of RNase-free water and 350 µL lysis
215 buffer. Add 250 µL of 100% ethanol to the diluted RNA and mix well by pipetting.

216 Note: Do not centrifuge. Proceed immediately to step 5.3.3.

217
218 5.3.3 Now transfer the 700 µL sample to an elution filter cartridge placed in a 2 ml collection tube
219 to collect the flow through. Centrifuge for 15 s at ≥ 8,000 x g. Discard the flow-through.

Commented [A11]: It is critical to place the tube on ice.. ?

Commented [DSF12]: Response to Editor ques #8: Yes, this step is vital.

Commented [A13]: Does this oligo binds to the RNA as well since it contains T instead of U?

Commented [DSF14]: Response to Editor ques #9: Oligos bind globin poly-A sites for removal during RT-PCR.

Commented [DSF15]: Response to Editor ques #10: Hybridize removed.

Commented [DSF16]: Response to Editor ques #11

Commented [A17]: Diluted with?

Commented [A18]: We cannot have commercial language in the manuscript. Please use generic terms instead.

Commented [DSF19]: Response to Editor ques #12: Replaced in text.

Commented [DSF20]: Response to Editor ques #13: Line altered to address question. Not necessary to check because EDTA addition quenches the reaction (digestion).

220
221 5.3.4 Repeat this by placing the elution filter cartridge into a new 2 ml collection tube. Add 500
222 μL of the mild washing buffer to the filter cartridge and centrifuge for 15 s at $\geq 8,000 \times g$ in order
223 to wash the filter cartridge membrane. Discard the flow-through. Reuse the collection tube in
224 step 5.3.5.

225
226 5.3.5 Now using the same sample tube, add another 500 μL of 80% ethanol to the filter cartridge.
227 This time centrifuge the tube for 2 min at $\geq 8,000 \times g$. Collect the elution spin column for the next
228 step and discard both the flow-through and collection tube.

229
230 5.3.7 Put the elution filter cartridge from the last step into a new 2 ml collection tube. Leave the
231 lid open on filter cartridge, and centrifuge at full speed for 5 min to dry the spin column
232 membrane and prevent ethanol carry over. Discard the flow-through and collection tube.
233 **Note:** To avoid damage angle the lids to point in a direction opposite to that of the rotor.

234
235 5.3.8. Take the dried filter cartridge and place into a new 1.5 ml collection tube. Add 14 μL RNase-
236 free water to filter cartridge membrane making sure to add the RNase-free water directly to the
237 center. Centrifuge for 60 s at full speed to elute the RNA.

238
239 5.4 Assess quality of globin-reduced RNA samples (step 6). Proceed to mRNA sample preparation
240 (step 7) and small RNA library preparation (step 8)^{16,17}.

241 Note: Globin-reduced RNA samples can now be stored at -20°C , however storage at -80°C is
242 recommended for long-term preservation.

243
244 **6. Assessment of RNA**

245
246 6.1 Quantify RNA concentration using a spectrophotometer. Examine the wavelength ratio of
247 260 and 280 nm. Ratios of ~ 2 or greater are considered pure for RNA and lower values indicate
248 some contamination. The instrument uses this ratio to determine RNA concentration as $\text{ng}/\mu\text{L}$.

249
250 6.2 Assess RNA quality by using 1 μL (100 ng) of sample on an appropriate chip. Final product
251 should be a RIN of ~ 2 or higher and a peak at ~ 280 nt for mRNA single read libraries; small RNA
252 libraries peaks at 143 correspond to miRNAs.

253
254 **7. Stranded Total RNA Sample Preparation for the mRNA and long ncRNA libraries.**¹⁶

255
256 **Note:** Bring elution buffer and rRNA removal beads to room temperature. Pre-label 0.2ml thin-
257 walled PCR tubes (plates may also be used). Protocol steps based on manufacturer's
258 instructions¹⁶.

259
260 **7.1.1 Start with 4 μL of globin-reduced RNA.** Add: 6 μL of nuclease-free water, 5 μL rRNA binding
261 buffer, 5 μL rRNA removal mix per tube (1st tube). Pipette to mix and re-cap. Place in thermocycler
262 with 100°C heated-lid for 5 min at 68°C . Remove and leave at room temperature for 60 s.

263

Commented [A21]: Please make two action steps.

Commented [DSF22]: Response to Editor ques #14: 7.1.1 now note and action step.

264 7.1.2 Resuspend rRNA removal beads by vortex; add 35 μ L of beads to new (2nd) PCR tube and
265 pipette sample from the 1st tubes onto beads in the 2nd tubes. Incubate 2nd PCR tube at room
266 temperature for 3 min. Place on the magnetic stand for 7 min.

267 **Note:** Mix each sample thoroughly by pipetting for optimal rRNA depletion. Raising/lowering
268 tube in the stand may help expedite separation process.

269
270 7.1.4 Transfer the supernatant from the 2nd PCR tube to matching 3rd PCR tube and place on to
271 the magnetic stand for a minimum of 60 s. Repeat the transfer into a new PCR tube only if the
272 beads are not moved to the tube sides.

273
274 7.1.5 Mix Sample Purification beads by vortex; add 99 μ L to each 3rd PCR tube. Allow to sit at
275 room temperature for 15 min. Place on magnetic stand for 5 mins. Ensure beads moves to sides.
276 Pipette to discard the supernatant.

277
278 7.1.6 Keep 3rd tube set on magnetic stand. Add 200 μ L of 70% ethanol whilst being careful not to
279 jostle the beads. Allow to sit for at least 30 s and pipette to discard the supernatant. Repeat step.

280
281 7.1.7 Allow sample to dry at room temperature for 15 min on the magnetic stand. Centrifuge the
282 kit elution buffer for 5 s at 600 x g. Remove tube and add 11 μ L of elution buffer mixing
283 thoroughly. Incubate for 2 min on the bench then at least 5 min on the magnetic stand at room
284 temperature.

285
286 7.1.8 Transfer 8.5 μ L of supernatant from the 3rd to a new (4th) PCR tube. Add 8.5 μ L of
287 Elute/Prime/Fragment High mix from the kit. Mix well. Cap and place in a thermocycler with pre-
288 heated lid for 8 min at 94 °C, hold at 4 °C. Remove and centrifuge briefly.

289
290 7.2 Synthesize First Strand cDNA

291
292 7.2.1 Allow first strand synthesis mix from kit to come to room temperature and centrifuge at
293 600 x g for 5 s. Transfer 50 μ L of reverse transcriptase into the first strand synthesis mix. reverse
294 transcriptase can be added to first strand synthesis at a ratio of 1:9.

295
296 7.2.2 Pipette 8 μ L of the combined reverse transcriptase / first strand synthesis mix into the 4th
297 tube with sample. Cap and centrifuge 600 x g for 5 s. Place into thermocycler with pre-heated lid
298 set at 100 °C. Run for: 10 min at 25 °C, 15 min at 42 °C, 15 min at 70 °C, then rest at 4 °C. Final
299 volume is ~25 μ L per well. Move immediately to next step.

300
301 7.3 Synthesize Second Strand cDNA

302
303 7.3.1 Bring End Repair Reagent (ERR) and Second Strand Mix (SSM) to room temperature, and
304 centrifuge at 600 x g for 5 s. Mix ERR with resuspension buffer at 1:50 dilution. Uncap 4th tube
305 and add 5 μ L of diluted ERR and 20 μ L of SSM to each; Pipette to mix well. Final volume ~50 μ L
306 ds cDNA.

307

Commented [DSF23]: Response to Editor ques #15: Line removed.

Commented [A24]: We cannot publish commercial terms. Please use generic term instead.

Commented [DSF25]: Response to Editor ques #16: commercial language removed.

Commented [A26]: We cannot publish commercial terms. Please use generic term instead.

Commented [DSF27]: Response to Editor ques #17: commercial language removed.

308 7.3.2 Cap and incubate in thermocycler for 1 hour at 16 °C. When the cycle is complete, remove
309 caps and allow to come to room temperature on bench top.

310
311 7.4 ds cDNA clean-up step

312
313 7.4.1 Begin by mixing the Solid Phase Reversible Immobilization (SPRI) paramagnetic beads by
314 vortex. Transfer 90 µL of the beads to the sample in the 4th tube (ds cDNA) and mix. Final volume
315 is 140 µL. Allow tubes to sit at room temperature for 15 min. incubate on magnetic stand for ~5
316 min. Remove ~135 µL of supernatant from each well. There should 5 µL left in each well.

317
318 7.4.2 Leave the tube on the magnetic stand and wash by adding 200 µL of 80% ethanol. Allow to
319 sit for 30 s then remove and discard supernatant. Repeat wash step. Leave tubes on the magnetic
320 stand for 15 min to allow to dry.

321
322 7.4.3 Centrifuge the resuspension buffer at 600 x g for 5 s after coming to room temperature.
323 Remove tubes from stand. Transfer 17.5 µL of resuspension buffer to the tube and mix. Let tubes
324 incubate on the bench top for 2 min then move to magnetic stand for 5 min.

325
326 7.4.4 Pipette 15 µL supernatant containing the ds cDNA samples to new set of (5th) 0.2 ml thin-
327 walled tubes. Move to next step promptly or seal and store at -15°C to -25°C for no more than 7
328 days.

329
330 7.5 Adenylate 3' Ends.

331
332 7.5.1 Pipette 2.5 µL of room temperature resuspension buffer into the sample tube then add 12.5
333 µL thawed A-Tailing mix. Mix completely by pipetting. Note: No A-tailing control used.

334
335 7.5.2 Cap and incubate in a thermocycler with 100 °C pre-heated lid. Run at 37 °C for 30 min, then
336 at 70 °C for 5 min with heated lid. Allow the thermocycler to rest at 4 °C.

337
338 7.6 Ligate index adaptors

339
340 7.6.1 Bring the RNA Adaptor tubes and the Stop Ligation Buffer mix to room temperature. For
341 each, centrifuge at 600 x g for 5 s. Leave ligation mix in freezer until ready for use. Sample pooling
342 arrangement for indexing should be known.

343
344 7.6.2 Add 2.5 µL of resuspension buffer and 2.5 µL of ligation mix to each sample tube. Now
345 pipette in 2.5 µL of the proper RNA Adaptor into each sample tube. Adaptor selection should be
346 performed by manufacturer's instructions for the chosen kit.

347
348 7.6.3 Recap and mix by centrifugation for 1 min at 280 x g. incubate in a thermal cycler for 10 min
349 at 30 °C. Add 5 µL of Stop Ligation Buffer to sample to halt the reaction and mix well.

350
351 7.6.4 Begin clean-up by mixing SPRI paramagnetic beads by vortex for 60 s. Transfer 42 µL of

352 beads to each tube. Mix thoroughly then incubate for 15 min on bench top.
353
354 7.6.5 Move tubes to magnetic stand and leave until liquid turns clear (~5 min). Once clear remove
355 and discard 79.5 µL of supernatant. Leave tubes on the magnetic stand and wash by adding 200
356 µL of 80% ethanol. Allow to sit for at least 30 s then remove and discard supernatant. Repeat
357 wash step. Leave on magnetic stand for additional ~15 min to allow for drying. Do not disrupt
358 beads during wash steps.
359
360 7.6.6 Add 52.5 µL of resuspension buffer to the sample tube and mix until the beads are
361 completely re-suspended. Incubate for 2 min on bench top. Move sample tube back to magnetic
362 stand until liquid is clear (~5 min).
363
364 7.6.7 Leave on stand and gently pipette 50 µL of supernatant to new (6th) tubes. Add 50 µL of
365 vortexed SPRI beads. Allow to incubate on bench top for 15 min.
366
367 7.6.8 Move to magnetic stand and let plate stay until liquid turns clear (~5 min). Discard 95 µL of
368 supernatant leaving 5 µL in each tube.
369
370 7.6.9 Leave tubes on magnetic stand and wash by adding 200 µL of 80% ethanol. Allow tube to
371 sit for 30 s then remove and discard supernatant. Repeat wash. Let tubes dry on magnetic stand
372 for 15 min.
373
374 7.6.10 Add in 22.5 µL of resuspension buffer and mix until beads are suspended. Incubate on
375 bench for 2 min then move to magnetic stand until liquid turns clear (~5 min). Pipette 20 µL of
376 sample into new PCR tubes (7th).
377
378 7.7 Bring needed kit components to room temperature. Transfer 5 µL of PCR primer mix and 25
379 µL of PCR master mix from kit to sample (7th PCR tube) containing indexed samples. Mix sample
380 and cap tube.
381
382 7.7.1 Place tubes into thermal cycler with pre-heated lid. Run program: 98 °C for 30 s, then 15
383 cycles of 98 °C for 10 s, 60 °C for 30 s, 72 °C for 30 s, 72 °C for 5 min. Hold at 4 °C.
384
385 7.8 Clean up samples by adding 50 µL of sufficiently mixed SPRI beads to sample tube and pipette
386 to combine. Allow the reaction to incubate on the bench top for 15 min.
387
388 7.8.1 Move the tube to a magnetic stand and incubate until liquid turns clear (~5 min). Discard
389 95 µL of supernatant leaving 5 µL in each tube.
390
391 7.8.2 Leave tube on magnetic stand and wash by adding 200 µL of 80% ethanol. Allow to sit for
392 at least 30 s and then remove and discard the supernatant. Repeat this wash step. Leave on
393 magnetic stand for 15 min to allow to dry.
394
395 7.8.3 Add in 32.5 µL of resuspension buffer and mix until beads are completely re-suspended.

Commented [DSF28]: Response to Editor ques #18: PCR amplification step is used to enrich DNA fragments that have ligated adapters.

Incubate on bench for 2 min then move to magnetic stand until liquid turns clear (~5 min). Pipette 30 µL of sample from each tube into a new PCR tube.

7.8.4 Assess DNA quantity and quality by repeating steps 6.1 and 6.2 with the appropriate chip. Move on to pooling step (step 9).

8. Small RNA Library Preparation for the sncRNAs.¹⁷

Note: Protocol steps based on manufacturer's instructions¹⁷.

8.1 Start small RNA library prep with ~220 ng - ~1.1 µg/ µL of globin-reduced RNA samples (4 µL).

8.2 Adaptor ligation

8.2.1 Ligate the 3'-adaptor by mixing thoroughly 1 µL adaptor, 2 µL nuclease-free water, and 4 µL of globin-reduced RNA sample. Incubate at 70 °C for 2 min and then immediately place on ice.

8.2.2 Add 10 µL 2X ligation buffer, and 3 µL ligation enzyme mix, mix and incubate 16 °C for 18 h. Note: Increase sample incubation time to 18 hours at a lower temperature of 16°C is recommended for studies interested in methylated RNA species, such as piwi RNAs. The longer incubation time allows for greater ligation efficiency due to the class of modifications during the 3' adaptor ligation phase.

8.2.3 Add 1 µL reverse transcription primer in 4.5 µL nuclease-free water to the ligation mixture to prevent adaptor-dimer formation of excess 3' adaptor.

8.2.4 Incubate in a preheated thermal cycler programmed for 5 min at 75 °C, 15 min at 37 °C, 15 min at 25 °C, and hold at 4°C.

8.2.5 Denature 1 µL of pre-diluted 5'-adaptor per sample in a thermal cycler at 70 °C for 2 min and then immediately place on ice.

8.2.6 Add 1 µL denatured 5'-adaptor, 1 µL 10X 5' ligation buffer and 2.5 µL 5' ligation enzyme. Mix and incubate 25 °C for 1 h.

8.3 cDNA synthesis and amplification

8.3.1 Mix thoroughly by pipetting 30 µL of 5'/3'-adapter-ligated RNA, 8 µL first-strand buffer, 1 µL RNase inhibitor and 1 µL of reverse transcriptase. Incubate the mixture at 50 °C for 60 min. Proceed immediately to PCR amplification or heat inactivate the reaction at 70 °C for 15 min and store at -20 °C.

8.3.2 PCR amplify the 40 µL Reverse Transcriptase reaction by adding 50 µL PCR master mix, 2.5 µL of RNA PCR primer, add 2.5 µL of designated RNA PCR Primer Index to sample, and nuclease-

Commented [DSF29]: Response to Editor ques #19: Line shortened to express only the action in the step.

440 free water to a total volume of 100 μ L. Mix by pipetting. Run PCR in the thermal cycler using the
441 following cycling conditions: initial denaturation at 94°C for 30 s; 11 cycles of 94°C for 15 s,
442 annealing at 62 °C for 30 s, and extension at 70 °C for 15 s; followed by a final extension at 70 °C
443 for 5 min; samples can be held at 4 °C in the cycler.

444 Note: Indexes are added for purpose of sample pooling.

445

446 8.4 Sample clean-up

447

448 8.4.1 From the DNA cleanup kit add 500 μ L binding buffer (5 M Gu-HCl, 30% isopropanol) to the
449 100 μ L of PCR amplified sample to enable efficient binding to the spin-column membrane.

450 Note: If the binding buffer has pH indicator included check that the color of the mixture is yellow.

451

452 8.4.2 Pipette 600 μ L mixture of sample and binding buffer into a filter cartridge inside a 2 mL
453 collection tube, spin for 30-60 s at 17,900 x g in a tabletop centrifuge. Discard flow-through.

454

455 8.4.3 Wash the filter cartridge into the same collection tube with 0.75 mL kit supplied wash buffer
456 (10mM Tris-HCl pH 7.5, 80% ethanol) by spinning for 30-60 s at 17,900 x g in a tabletop centrifuge
457 and discarding the flow –through. Spin the filter cartridge dry for an additional 60 s.

458

459 8.4.4 Place the filter cartridge in a clean 1.5 mL collection tube. Add 30 μ L elution buffer (10 mM
460 Tris-HCl pH 8.5), let the column stand for 60 s, and then spin for 60 s at 17,900 x g in a tabletop
461 centrifuge.

462

463 8.5 Assess sample quality by repeating steps 6.1 and 6.2 with the appropriate DNA chip. Move to
464 pooling step (step 9).

465

466 9. Sample pooling for sequencing

467

468 9.1 Pool barcoded and QC'ed samples by setting up and labeling new tubes (or a 96 well PCR
469 plate) to contain either the mRNA or ncRNA samples. Transfer 13 μ L of each barcoded 10nM
470 library to corresponding tubes (or wells in the new plate) per manufacturer's instructions^{16,17}.
471 Submit pooled samples for sequencing, be certain to choose similar read lengths if samples are
472 to be run on the same chip.

473

474 REPRESENTATIVE RESULTS:

475 The representative samples in our study are the globin and ribo-depleted whole blood samples.
476 The representative outcome of the protocol consists of a globin depleted library sample with an
477 RNA integrity number (RIN) above 7 (Figure 1a) and 260/280 nm concentration ratios at or above
478 2 (Figure 1b and 1c). Validation of the sample outcome was performed using spectrophotometer
479 to give the final concentration of each library and chip-based electrophoresis to give RIN number
480 along with a graph of peaks that show which molecules (mRNA or ncRNA) were captured based
481 on insert size in the library sample prior to pooling and sequencing. For the current study focus
482 was placed on the mRNA and small ncRNAs only. For the mRNA libraries the representative result
483 is an electropherogram peak at ~280 bp (Figure 1d). For the small ncRNAs representative results

Commented [A30]: Please use imperative tense.

Commented [DSF31]: Response to Editor ques #20:
Imperative tense inserted.

Commented [DSF32]: Response to Editor ques #21:
Original line #s 520-522 removed.

Commented [A33]: Please combine all the panels in one
figure. Alternatively

Commented [DSF34]: Response to Editor ques #22: All
panels combined into 1 figure.

consist of a range of peaks from ~100-400 bp (**Figure 1e**), with peaks at ~143 and ~153 bp correspond to miRNAs and piRNAs, respectively. Our sample results showed that our optimized technique resulted in libraries with RIN numbers that ranged from 6.5 (sub-optimal) to 9.2 (above optimal). This proved to be an improvement compared to other studies that used globin depletion methods and were only able to achieve RIN numbers at or near 6. The chip-based electrophoresis results also showed that from one single blood sample it is possible to achieve RNA molecule capture of peaks representing both mRNA and ncRNA (**Figures 1d and 1e**), and sample insert sizes that covered both small and large non-coding RNA molecules. These results are representative of the optimal RIN scores and insert sizes needed to ensure quality transcript reads can be sequenced from the prepared libraries for nearly all NGS RNA-seq platforms.

FIGURE LEGENDS:

Figure 1: Representative results for mRNA and non-coding RNA expression libraries from single samples of porcine whole blood. (a) Electrophoresis file run summary representation of RIN numbers pre-globin reduction and post-globin reduction. Representative electropherograms (b) pre-globin reduction and (c) post-globin reduction of a single porcine whole blood sample. (d) Representative electropherograms of globin and ribo-depleted whole blood mRNA library samples prior to pooling and sequencing. (e) Representative electropherograms of globin depleted sncRNA whole blood libraries of the same samples featured in panel d prior to pooling and sequencing.

DISCUSSION:

The first critical step in the protocol that made it optimized included the added globin depletion steps, which made it possible to get quality reads from whole blood samples. One of the largest limitations on using whole blood in sequencing studies are the high numbers of reads in the sample that will map to globin molecules and reduce the reads that could map to other molecules of interest¹⁸. Therefore, in optimizing the protocol for our sample type, we needed to incorporate a globin depletion step to ensure the highest possible mRNA and non-coding RNA capture through sequencing. All of the samples were globin depleted to account for high levels of globin transcripts using porcine specific hemoglobin A and B (HBA and HBB) oligonucleotides based on the procedure from Choi *et al.*, 2014¹⁴. Another critical step was the use of a ribo-depletion extraction kit which allowed for the ability remove unwanted RNA molecules from our samples while simultaneously retaining polyadenylated non-coding and viral RNA molecules of experimental interest within our libraries for sequencing¹⁹. Additionally, we were able to work with single samples to create both coding (mRNA) and non-coding (small and long) RNA libraries by removing the small RNA enrichment and size selection steps. By doing this, we were able to maximize the RNA aliquots for pooling and to guarantee we sent enough volume to allow for sequencing. Optimization of all of the manufacturer's protocols were done to increase mRNA and non-coding RNA recovery for downstream library creation. We also optimized the non-coding RNA library preparation for the small non-coding RNA portion of our downstream analysis by adding additional time to the PCR incubation per the manufacturers' instructions to increase the ligation efficiency of methylated RNAs and the likelihood of recovering piwi-RNAs.

Commented [DSF35]: Response to Editor ques #23: All panels combined into 1 figure.

Commented [DSF36]: Response to Editor ques #24: Legends made to be more descriptive.

Commented [A37]: Figure legends should have one liner title and description of all the panels in detail and description of the all the abbreviations.

Commented [DSF38]: Response to Editor ques #25: personal pronoun removed.

The limitations of the optimized protocol are rooted in the lack of specificity that also makes it ideal for multiple analyses from a single sample. By removing the enrichment and size selection portions of the non-coding RNA library preparation, we are limiting novel small RNA discovery²⁰ to balance it against acquiring both small and long non-coding RNA molecules. Hence, the optimized protocol gives a good overview of the types of non-coding expression but, requires extra processing during the analytical stage of mapping to gain more specific information on the classes and sizes of non-coding RNA expression captured in the sample. Another limitation to the method is due to the incorporation of the globin depletion steps, which will lower the overall RIN numbers. We were able to troubleshoot this by first examining the RIN numbers before and after globin depletion to understand how large of a reduction in RNA integrity we would experience. This yielded results that showed we could experience a drop of ~1-2 points. To account for this drop we further optimized the globin depletion protocol to 6 µg of sample instead of the recommended 10 µg. This helped to improve our RIN numbers, as well as, conserve sample. Additionally, we also used thin-walled microamp reaction tubes and iced the tubes immediately after the first denaturing step. This allowed us to quench the reaction quicker, allowing for improved RIN numbers on the globin depleted samples.

The modified protocol used in this study had several advantages over the base library preparation methods used in other studies where only mRNA or non-coding RNA are studied separately. The changes we made for optimization allowed for efficient use of whole blood as our sample type. This is advantageous to future studies because whole blood as a sample type has less steps for collection and processing than the leukocyte portion of blood. Additionally, whole blood also has the added advantages of allowing researchers to examine systemic responses and can be repeatedly collected from the animal. However, there is a caveat to the use of whole blood samples in that the amount of blood you are able to collect is dependent on the age and size of the organism in question. Smaller amounts of blood will lead to lower RNA yields, however the method presented here will allow dual library creation from at least 2.5ml of whole blood. This allowed for us to use one sample to investigate both mRNA and non-coding expression and correlate both to an individual at a snapshot in time, effectively allowing us to collect more information than traditional library preparation methods. Also, by performing the ribo-depletion we made it possible to reduce transcripts that could deplete sequencing reads, while still retaining the non-coding and viral transcripts that can be lost during traditional library creation to study mRNA. In this way our optimized protocol triples the amount of information that can be gained from a single sample. Other significant changes we made to the method to facilitate the capture of more RNA species information were: changing the type of PCR plate to quickly stop the reaction, which improved RNA purity yield; a longer thermocycler incubation to increase ligation for possible recovery of methylated RNAs; and not employing a size selection gel, to allow for all detectable ncRNAs to identified regardless of length. For the downstream analysis this allowed for the capture of multiple non-coding RNAs between 18nt-200nt in length.

By employing this optimized method, our group was able to put forward a protocol that can be applied to whole blood transcriptomic analyses and allow for mRNA and non-coding RNA from a single sample.

ACKNOWLEDGMENTS:

This work was mainly supported by the by USDA NIFA AFRI 2013-67015-21236, and in part by USDA NIFA AFRI 2015-67015-23216. This study was supported in part by an appointment to the Agricultural Research Service Research Participation Program administered by the Oak Ridge Institute for Science and Education (ORISE) through an interagency agreement between the US Department of Energy (DOE) and the US Department of Agriculture. ORISE is managed by Oak Ridge Associated Universities under DOE contract no. DE-AC05-06OR2310.

We would like to thank Dr. Kay Faaberg for the HP-PRRSV infectious clones, Dr. Susan Brockmeier for her help with animals involved in the experiment, and Sue Ohlendorf for secretarial assistance in preparation of the manuscript.

DISCLOSURES:

The authors have nothing to disclose.

Mention of trade names or commercial products in this article is solely for the purpose of providing specific information and does not imply recommendation or endorsement by the U.S. Department of Agriculture. USDA is an equal opportunity provider and employer.

REFERENCES:

- 1 Finotello, F. & Di Camillo, B. Measuring differential gene expression with RNA-seq: challenges and strategies for data analysis. *Briefings in Functional Genomics*. **14** (2), 130-142, doi:10.1093/bfpg/elu035, (2015).
- 2 Coble, D. J. *et al.* RNA-seq analysis of broiler liver transcriptome reveals novel responses to high ambient temperature. *BMC Genomics*. **15** 1084, doi:10.1186/1471-2164-15-1084, (2014).
- 3 Koltes, J. E. *et al.* Identification of a putative quantitative trait nucleotide in guanylate binding protein 5 for host response to PRRS virus infection. *BMC Genomics*. **16** 412, doi:10.1186/s12864-015-1635-9, (2015).
- 4 Miller, L. C. *et al.* Comparative analysis of signature genes in PRRSV-infected porcine monocyte-derived cells to different stimuli. *PLoS One*. **12** (7), e0181256, doi:10.1371/journal.pone.0181256, (2017).
- 5 Head, S. R. *et al.* Library construction for next-generation sequencing: overviews and challenges. *Biotechniques*. **56** (2), 61-64, 66, 68, passim, doi:10.2144/000114133, (2014).
- 6 Wang, Z., Gerstein, M. & Snyder, M. RNA-Seq: a revolutionary tool for transcriptomics. *Nature reviews. Genetics*. **10** (1), 57-63, doi:10.1038/nrg2484, (2009).
- 7 Dominic, O. N., Heike, G. & Martin, S. Ribosomal RNA Depletion for Efficient Use of RNA-Seq Capacity. *Current Protocols in Molecular Biology*. **103** (1), 4.19.11-14.19.18, doi:doi:10.1002/0471142727.mb0419s103, (2013).
- 8 Pollet, S. *FR-AgEncode: a French pilot project to enrich the annotation of livestock genomes- RNA Extraction Protocol*, <[ftp://ftp.faang.ebi.ac.uk/ftp/protocols/assays/INRA SOP RNA_extraction_20160504.pdf](ftp://ftp.faang.ebi.ac.uk/ftp/protocols/assays/INRA_SOP_RNA_extraction_20160504.pdf)> (2015).
- 9 Fleming, D. S. & Miller, L. C. Identification of small non-coding RNA classes expressed in

Commented [DSF39]: Response to Editor ques #26: This language is required on all ARS manuscripts and must be stated under disclosures regardless of if they appear in the article or table of materials.

Commented [A40]: We cannot have commercial names in the manuscript. Please use generic terms and refer to the commercial term in the table of materials.

swine whole blood during HP-PRRSV infection. *Virology*. **517** 56-61, doi:10.1016/j.virol.2018.01.027, (2018).

10 Fleming, D. S. & Miller, L. C. Small non-coding RNA expression status in animals faced with highly pathogenic porcine reproductive and respiratory syndrome virus (HP-PRRSV). *Proceedings of the World Congress on Genetics Applied to Livestock Production. Species - Porcine* **2**, (11), (2018).

11 Bivens Nathan, J. & Zhou, M. RNA-Seq Library Construction Methods for Transcriptome Analysis. *Current Protocols in Plant Biology*. **1** (1), 197-215, doi:10.1002/cppb.20019, (2016).

12 Cui, P. *et al.* A comparison between ribo-minus RNA-sequencing and polyA-selected RNA-sequencing. *Genomics*. **96** (5), 259-265, doi:<https://doi.org/10.1016/j.ygeno.2010.07.010>, (2010).

13 Guo, Y. *et al.* RNAseq by Total RNA Library Identifies Additional RNAs Compared to Poly(A) RNA Library. *BioMed Research International*. **2015** **9**, doi:10.1155/2015/862130, (2015).

14 Choi, I. *et al.* Increasing gene discovery and coverage using RNA-seq of globin RNA reduced porcine blood samples. *BMC Genomics*. **15** 954, doi:10.1186/1471-2164-15-954, (2014).

15 Wu, K., Miyada, G., Martin, J. & Finkelstein, D. *Globin reduction protocol: A method for processing whole blood RNA samples for improved array results.*, <http://media.affymetrix.com:80/support/technical/technotes/blood2_technote.pdf> (2007).

16 Illumina. TruSeq® Stranded Total RNA Sample Preparation Guide.

17 New England Biolabs Inc. *Protocol for use with NEBNext Multiplex Small RNA Library Prep Kit for Illumina (Index Primers 1-48) (E7560)*, <https://www.neb.com/~media/Catalog/All-Protocols/758F75A29CDE4D03954E1EF75E78EA3D/Content/manualE7560_Figure_1.jpg> (

18 Shin, H. *et al.* Variation in RNA-Seq transcriptome profiles of peripheral whole blood from healthy individuals with and without globin depletion. *PLoS One*. **9** (3), e91041, doi:10.1371/journal.pone.0091041, (2014).

19 Costa, V., Angelini, C., De Feis, I. & Ciccociola, A. Uncovering the complexity of transcriptomes with RNA-Seq. *Journal of Biomedicine and Biotechnology*. **2010** 853916, doi:10.1155/2010/853916, (2010).

20 Martens-Uzunova, E. S., Olvedy, M. & Jenster, G. Beyond microRNA--novel RNAs derived from small non-coding RNA and their implication in cancer. *Cancer Letters*. **340** (2), 201-211, doi:10.1016/j.canlet.2012.11.058, (2013).