Journal of Visualized Experiments

Analysis and Visualization of Mouse and Human Antibody Repertoires by Next-generation Sequencing --Manuscript Draft--

Article Type:	Invited Methods Article - JoVE Produced Video		
Manuscript Number:	JoVE58804R2		
Full Title:	Analysis and Visualization of Mouse and Human Antibody Repertoires by Next- generation Sequencing		
Keywords:	Antibody repertoire; next-generation sequencing; NGS; RNAseq; immune response; antigen recognition; V-D-J recombination; Immunoglobulins		
Corresponding Author:	Kazuo Ohnishi National Instutute of Infectious Diseases Shinjuku, Tokyo JAPAN		
Corresponding Author's Institution:	National Instutute of Infectious Diseases		
Corresponding Author E-Mail:	ohnishik@nih.go.jp		
Order of Authors:	Lin Sun		
	Naoko Kono		
	Hiroyuki Toh		
	Hanbing Xue		
	Kaori Sano		
	Tadaki Suzuki		
	Akira Ainai		
	Yasuko Orba		
	Junya Yamagishi		
	Hideki Hasegawa		
	Yoshimasa Takahashi		
	Shigeyuki Itamura		
	Kazuo Ohnishi		
Additional Information:			
Question	Response		
Please indicate whether this article will be Standard Access or Open Access.	Open Access (US\$4,200)		
Please indicate the city, state/province, and country where this article will be filmed . Please do not use abbreviations.	Department of Immunology, National Institute of Infectious Diseases, Tokyo 162-8640, Japan		

Department of Immunology National Institute of Infectious Diseases 1-23-1 Toyama, Shinjuku-ku, Tokyo 162-8640, JAPAN

Phone: +81-3-4582-2711 (direct), +81-3-5285-1111 (op) Fax: +81-3-5285-1156

E-mail: ohnishik@nih.go.jp HP: http://www.nih.go.jp/niid/



国立感染症研究所

September 27, 2018

Dr. Vineeta Bajaj Review Editor JoVE

Dear Editor:

I wish to submit a revised manuscript for publication in JoVE titled "Analysis and Visualization of Mouse and Human Antibody Repertoires by Next-generation Sequencing." The paper was coauthored by Lin Sun, Naoko Kono, Hiroyuki Toh, Hanbing Xue, Kaori Sano, Tadaki Suzuki, Akira Ainai, Yasuko Orba, Junya Yamagishi, Hideki Hasegawa, Yoshimasa Takahashi, and Shigeyuki Itamura.

We would like to thank the editorial comments to improve the quality of the manuscript. Our point-by-point responses are uploaded as rebuttal documents and all changes to the manuscript are tracked.

We also uploaded revised files of Table_1 and Table_of_Materials.

We also uploaded the file for copyright permission.

We hope that the revised version of this manuscript is now suitable for publication in JoVE.

We look forward to hearing from you at your earliest convenience.

Sincerely,

Kazuo Ohnishi, Ph.D.

Department of Immunology, National Institute of Infectious Diseases,

Tokyo 162-8640, Japan Email: ohnishik@nih.go.jp Phone: +81-3-5285-1111

Fax: +81-3-5285-1156

1 TITLE: 2 Identification of Mouse and Human Antibody Repertoires by Next-Generation Sequencing 3 4 **AUTHORS & AFFILIATIONS:** Lin Sun^{1,*}, Naoko Kono^{2,*}, Hiroyuki Toh³, Hanbing Xue¹, Kaori Sano^{4,5}, Tadaki Suzuki⁴, Akira 5 Ainai⁴, Yasuko Orba⁶, Junya Yamagishi^{7,8}, Hideki Hasegawa^{4,5}, Yoshimasa Takahashi⁹, 6 7 Shigeyuki Itamura², Kazuo Ohnishi^{9,10} 8 9 ¹Graduate School of Life and Environmental Sciences, University of Tsukuba, Ibaraki, Japan 10 ²Center for Influenza Virus Research, National Institute of Infectious Diseases, Tokyo, Japan ³School of Science and Technology, Kwansei Gakuin University, Hyogo, Japan 11 12 ⁴Department of Pathology, National Institute of Infectious Diseases, Tokyo, Japan 13 ⁵Division of Infectious Diseases Pathology, Department of Global Infectious Diseases, Tohoku 14 University Graduate School of Medicine, Miyagi, Japan 15 ⁶Division of Molecular Pathobiology, Research Center for Zoonosis Control, Hokkaido 16 University, Sapporo, Japan. 17 ⁷Division of Collaboration and Education, Research Center for Zoonosis Control, Hokkaido 18 University, Sapporo, Japan. 19 ⁸Global Station for Zoonosis Control, GI-CoRE, Hokkaido University, Sapporo, Japan. 20 ⁹Department of Immunology, National Institute of Infectious Diseases, Tokyo, Japan 21 ¹⁰Faculty of Life and Environmental Sciences, University of Tsukuba, Ibaraki, Japan 22 23 *These authors contributed equally to this work. 24 25 **Corresponding Author:** 26 Kazuo Ohnishi (ohnishik@nih.go.jp) 27 28 **Email Addresses of Co-authors:** 29 Lin Sun (sunlin@niid.go.jp) 30 Naoko Kono (nao2.26torak@gmail.com) 31 Hiroyuki Toh (tohhir@kwansei.ac.jp) 32 Hanbing Xue (aloisxue@nih.go.jp) (ka--sano@nih.go.jp) 33 Kaori Sano 34 (tksuzuki@nih.go.jp) Tadaki Suzuki 35 Akira Ainai (ainai@nih.go.jp) 36 Yasuko Orba (orbay@czc.hokudai.ac.jp) 37 Junya Yamagishi (junya@czc.hokudai.ac.jp) 38 Hideki Hasegawa (hasegawa@nih.go.jp) 39 Yoshimasa Takahashi (ytakahas@niid.go.jp) 40 Shigeyuki Itamura (sitamura@nih.go.jp) 41

KEYWORDS:

43 Antibody Repertoire, Next-Generation Sequencing, NGS, RNA-seq, Immune Response,

44 Antigen Recognition, V-D-J Recombination, Immunoglobulin

46 **SUMMARY:**

42

45

47 Here, we describe protocols for the analysis and visualization of the structure and constitution

of whole antibody repertoires. This involves the acquisition of vast sequences of antibody RNA using next-generation sequencing.

LONG ABSTRACT:

The immense adaptability of antigen recognition by antibodies is the basis of the acquired immune system. Despite our understanding of the molecular mechanisms underlying the production of the vast repertoire of antibodies by the acquired immune systems, it has not yet been possible to arrive at a global view of a complete antibody repertoire. In particular, B cell repertoires have been regarded as a black box because of their astronomical number of antibody clones. However, next-generation sequencing technologies are enabling breakthroughs to increase our understanding of the B cell repertoire. In this report, we describe a simple and efficient method to visualize and analyze whole individual mouse and human antibody repertoires. From the immune organs, representatively from spleen in mice and peripheral blood mononuclear cells in humans, total RNA was prepared, reverse transcribed, and amplified using the 5'-RACE method. Using a universal forward primer and antisense primers for the antibody class-specific constant domains, antibody mRNAs were uniformly amplified in proportions reflecting their frequencies in the antibody populations. The amplicons were sequenced by next-generation sequencing (NGS), yielding more than 10⁵ antibody sequences per immunological sample. We describe the protocols for antibody sequence analyses including V(D)J-gene-segment annotation, a bird's-eye view of the antibody repertoire, and our computational methods.

INTRODUCTION:

The antibody system is one of the fundamentals of the acquired immune system. It is highly potent against invading pathogens due to its vast diversity, fine antigen recognition specificity, and the clonal expansion of antigen-specific B cells. The repertoire of antibody-producing B cells is estimated to be more than 10^{15} in a single individual¹. This immense diversity is generated with the help of VDJ gene recombination in the immunoglobulin genetic loci². Description of the entire B cell repertoires and their dynamic changes in response to antigen-immunization is therefore challenging, but essential for a complete understanding of the antibody response against invading pathogens.

Because of their astronomical diversity, B cell repertoires have been regarded as a black box; however, the advent of NGS technology has enabled breakthroughs to an enhanced understanding of their complexity^{3,4}. Whole antibody repertoires have been successfully analyzed, firstly in zebrafish⁵, then mice⁶, and humans^{6,7}. Although NGS has now become a powerful tool in the study of the adaptive immune response, basic analyses of the commonalities and differences in antibody repertoires among individual animals are lacking.

In mice, it was reported that the IgM repertoires are almost identical between individuals, whereas those of IgG1 and IgG2c are substantially different between individuals⁸. In addition to V-gene usage profile, the observed frequency of VDJ-profile in naive peripheral B cells is highly similar between individuals⁸. The analysis of the amino acid sequences of the VDJ-region also showed the occurrence of the same junctional sequences in different mice much more frequently than previously thought⁸. These results indicate that the mechanisms for the antibody repertoire formation can be deterministic rather than stochastic^{5,8,9}. The process of antibody repertoire development in mice has also been successfully analyzed using NGS to

Page 1 of 6 revised October 2016

95 further highlight the potential of NGS to uncover the antibody immune system in detail¹⁰.

96 97

In this report, we describe a simple and efficient method to visualize and analyze an antibody repertoire at a global level.

98 99 100

PROTOCOL:

101102

103

104

105

106

All animal experiments were performed according to institutional guidelines and with the approval of the National Institute of Infectious Diseases Animal Care and Use Committee. Sampling of PBMCs from healthy adult volunteers, used as the representative result in this report, was performed with the approval of the Ethics Committee of the National Institute of Infectious Diseases, Tokyo, Japan, and written informed consent was obtained from each participant using an ethics committee-approved form.

107108109

1. Primer Design

110

1.1. Design a universal forward primer to cDNA to amplify the immunoglobulin mRNA without
 bias from PCR primers, as used in the 5'-RACE^{11,12} and SMART-PCR¹³ techniques.

113114

1.2. For the immunoglobulin VH gene amplification, design the immunoglobulin class-specific sequences in the constant region as reverse primers^{8,14} (**Figure 1A**).

115116117

[Please insert Table 1 here]

118119

NOTE: Multiplex tag sequences can be added to any of these primers to label the library molecules from different sample sources. Sequences for nested PCR can also be added, according to the manual of the kit used¹⁵.

121 122

120

2. Nucleic Acid Isolation from Immune Cells and Tissues

123124125

NOTE: The procedure given below is for extracting nucleic acids from the mouse spleen. However, it is applicable to other immune tissues and human cells such as lymph nodes or peripheral blood mononuclear cells (PBMCs) (Figure 1B).

127128

126

2.1. Dissect the tissue, e.g., spleen from an 8-week-old C57BL/6 mouse and pass it through a stainless-steel mesh (200 to 400 μm) with 2 mL of PBS buffer to obtain dispersed cells.
 Transfer the cell suspension to a 2.0 mL microcentrifuge tube, and centrifuge for 5 min at 600

132 $\times g$ and 4 °C. Discard the supernatant.

133134

2.2. Add 800 μL of ACK lysing buffer (150 mM NH₄Cl, 1 mM KHCO₃, 0.1 mM Na₂EDTA, pH 7.2) to the pellet, and incubate on ice for 2 min to lyse red blood cells in the tissue.

135136

2.3. Wash the tissue cells with 2 mL of PBS 3x, followed by centrifugation for 5 min at $600 \times g$ and 4 °C.

139

2.4. Add 800 μL of phenol/guanidine isothiocyanate reagent to the pellet, vortex thoroughly,
 and incubate at about 25 °C for 5 min.

Page 2 of 6

2.5. Add chloroform (200 μ L), shake manually for 15 s, and then incubate for 2 min at about 25 °C.

2.6. Separate the phases by centrifugation for 15 min at 12,000 × g and 25 °C and transfer the upper aqueous phase to a fresh tube.

2.7. Add one volume of 70% ethanol, vortex briefly and apply it to the silica spin column.

2.8. Elute the RNA with 30–100 μL of water.

2.9. Quantitate initial RNA concentration using a fluorometer (**Table of Materials**).

2.10. Store the purified RNA at -80 °C.

3. cDNA Synthesis and PCR Amplification

NOTE: The method described below is based on the 5'-RACE^{11,12} and SMART-PCR techniques¹³. The details and optimization of the reaction are described in the manual of the kit ¹⁵. The starting materials for mouse immunoglobulin are the sample from step 2.10. The starting materials for human immunoglobulin are the sample from human tissues, ex. PBMC, treated as described in steps 2.3 to 2.10.

3.1 Synthesize the first-strand cDNA from 2 to 10 μ g of total RNA template using 5′-RACE CDS primer (oligo-dT-containing) and SMART-PCR oligonucleotide (**Table of Materials**) according to the manufacturer's instructions¹⁵.

3.2.1 For the mouse immunoglobulin, PCR-amplify cDNA with high-fidelity DNA Polymerase using the universal forward primer and immunoglobulin class-specific reverse primers (**Table 1**). Set the thermal cycling conditions as: 94 °C for 2 min, then 40 cycles of 94 °C for 30 s, 59 °C for 30 s, and 72 °C for 30 s, followed by a final extension step at 72 °C for 5 min.

NOTE: Typical experiments amplify IgM, IgG1, IgG2c, Ig κ and Ig λ immunoglobulin classes to look at the naive, Th1-dependent and Th2-dependent B cells (**Figure 3**).

3.2.2 For the human immunoglobulin, perform the 1st PCR using the universal forward primer and immunoglobulin class-specific reverse primers (**Table 1**) with tag sequences. Include the index sequences for each sample by 2nd PCR using index sequence primers. Use the following PCR conditions and the Taq polymerase: 94 °C for 2 min, 21 cycles (1st PCR) or 32 cycles (2nd PCR) at 94 °C for 30 s, 59 °C for 30 s, 72 °C for 30 s.

NOTE: Typical experiments amplify IgM, IgD, IgG (IgG1, IgG2, IgG3 and IgG4), IgA (IgA1 and IgA2), IgE, Igκ and Igλ immunoglobulin classes to look at all B cell populations (**Figure 4**).

3.3 Electrophorese the PCR products on an agarose gel and purify 600 to 800 bp fragments using a silica membrane spin-column.

189	3.3.1 Electrophorese the sample from 3.2.1 or 3.2.2 on 2% agarose gel.
190	
191	3.3.2 Visualize the DNA bands on UV-transilluminator and excise the gel-slice containing the
192	broad band between 600 to 800 bp.
193	
194	3.3.3 Add 10 µL of membrane binding solution per 10 mg of gel slice. Mix and incubate at 50–
195	65 °C until the gel slice is completely dissolved.
196	, ,
197	3.3.4 Transfer the gel solution on silica membrane spin-column. Wash once with washing
198	buffer and elute DNA with 50 μL of nuclease-free water (Table of Materials).
199	in the state of th
200	3.4 Quantify the purified amplicons with a fluorometer and pool amplicons from each
201	immunoglobulin class in equal amounts for NGS sequencing.
202	minuments and some equal announce for the sequencing.
203	NOTE: Typically, 2-10 μg amplicon DNA was recovered for each immunoglobulin class. Mix
204	each sample solution equally in DNA amount to give rise 50 μ L solution containing 10-20 ng
205	DNA/μL.
206	3.5 Determine the size and concentration of libraries using a micro-capillary based
207	
208	electrophoresis with DNA sizing chip (Table of Materials). Store the libraries at - 20 °C.
209	4 NGC Companying of Librarian
210	4. NGS Sequencing of Libraries
211	4.4. Consider a Conside Charles of the later and a series of the second consideration
212	4.1 Generate a SampleSheet.cvs for the sequencing run specifying sample name, index
	information and inclused to platein forth files and .
213	information and instruct to obtain .fastq files only.
214	
214 215	information and instruct to obtain .fastq files only. 4.2 Thaw the reagent cartridge (Table of Materials) and the libraries.
214215216	4.2 Thaw the reagent cartridge (Table of Materials) and the libraries.
214 215 216 217	
214 215 216 217 218	4.2 Thaw the reagent cartridge (Table of Materials) and the libraries.4.3 Make 0.2 N NaOH and dilute the libraries to obtain the desired molar concentration.
214 215 216 217 218 219	 4.2 Thaw the reagent cartridge (Table of Materials) and the libraries. 4.3 Make 0.2 N NaOH and dilute the libraries to obtain the desired molar concentration. 4.4 Rinse and dry the flow cell. Add 600 μL of diluted and denatured library solution into the
214 215 216 217 218 219 220	4.2 Thaw the reagent cartridge (Table of Materials) and the libraries.4.3 Make 0.2 N NaOH and dilute the libraries to obtain the desired molar concentration.
214 215 216 217 218 219 220 221	 4.2 Thaw the reagent cartridge (Table of Materials) and the libraries. 4.3 Make 0.2 N NaOH and dilute the libraries to obtain the desired molar concentration. 4.4 Rinse and dry the flow cell. Add 600 μL of diluted and denatured library solution into the well of the reagent cartridge.
214 215 216 217 218 219 220 221 222	 4.2 Thaw the reagent cartridge (Table of Materials) and the libraries. 4.3 Make 0.2 N NaOH and dilute the libraries to obtain the desired molar concentration. 4.4 Rinse and dry the flow cell. Add 600 μL of diluted and denatured library solution into the
214 215 216 217 218 219 220 221 222 223	 4.2 Thaw the reagent cartridge (Table of Materials) and the libraries. 4.3 Make 0.2 N NaOH and dilute the libraries to obtain the desired molar concentration. 4.4 Rinse and dry the flow cell. Add 600 μL of diluted and denatured library solution into the well of the reagent cartridge. 4.5 Start the sequencing run.
214 215 216 217 218 219 220 221 222 223 224	 4.2 Thaw the reagent cartridge (Table of Materials) and the libraries. 4.3 Make 0.2 N NaOH and dilute the libraries to obtain the desired molar concentration. 4.4 Rinse and dry the flow cell. Add 600 μL of diluted and denatured library solution into the well of the reagent cartridge.
214 215 216 217 218 219 220 221 222 223 224 225	 4.2 Thaw the reagent cartridge (Table of Materials) and the libraries. 4.3 Make 0.2 N NaOH and dilute the libraries to obtain the desired molar concentration. 4.4 Rinse and dry the flow cell. Add 600 μL of diluted and denatured library solution into the well of the reagent cartridge. 4.5 Start the sequencing run. 5. Quality Control of NGS Data
214 215 216 217 218 219 220 221 222 223 224 225 226	 4.2 Thaw the reagent cartridge (Table of Materials) and the libraries. 4.3 Make 0.2 N NaOH and dilute the libraries to obtain the desired molar concentration. 4.4 Rinse and dry the flow cell. Add 600 μL of diluted and denatured library solution into the well of the reagent cartridge. 4.5 Start the sequencing run.
214 215 216 217 218 219 220 221 222 223 224 225 226 227	 4.2 Thaw the reagent cartridge (Table of Materials) and the libraries. 4.3 Make 0.2 N NaOH and dilute the libraries to obtain the desired molar concentration. 4.4 Rinse and dry the flow cell. Add 600 μL of diluted and denatured library solution into the well of the reagent cartridge. 4.5 Start the sequencing run. 5. Quality Control of NGS Data 5.1. Perform the quality control of FASTQ data using the "FASTX-Toolkit"¹⁶.
214 215 216 217 218 219 220 221 222 223 224 225 226	 4.2 Thaw the reagent cartridge (Table of Materials) and the libraries. 4.3 Make 0.2 N NaOH and dilute the libraries to obtain the desired molar concentration. 4.4 Rinse and dry the flow cell. Add 600 μL of diluted and denatured library solution into the well of the reagent cartridge. 4.5 Start the sequencing run. 5. Quality Control of NGS Data 5.1. Perform the quality control of FASTQ data using the "FASTX-Toolkit" 16. NOTE: A basic example of the parameter settings used is as follows:
214 215 216 217 218 219 220 221 222 223 224 225 226 227	 4.2 Thaw the reagent cartridge (Table of Materials) and the libraries. 4.3 Make 0.2 N NaOH and dilute the libraries to obtain the desired molar concentration. 4.4 Rinse and dry the flow cell. Add 600 μL of diluted and denatured library solution into the well of the reagent cartridge. 4.5 Start the sequencing run. 5. Quality Control of NGS Data 5.1. Perform the quality control of FASTQ data using the "FASTX-Toolkit"¹⁶. NOTE: A basic example of the parameter settings used is as follows: fastq_quality_trimmer -v -t 20 -l 200 -i [InFilename.fastq] -o [InFilename.fastq]
214 215 216 217 218 219 220 221 222 223 224 225 226 227 228	 4.2 Thaw the reagent cartridge (Table of Materials) and the libraries. 4.3 Make 0.2 N NaOH and dilute the libraries to obtain the desired molar concentration. 4.4 Rinse and dry the flow cell. Add 600 μL of diluted and denatured library solution into the well of the reagent cartridge. 4.5 Start the sequencing run. 5. Quality Control of NGS Data 5.1. Perform the quality control of FASTQ data using the "FASTX-Toolkit" 16. NOTE: A basic example of the parameter settings used is as follows:
214 215 216 217 218 219 220 221 222 223 224 225 226 227 228 229	 4.2 Thaw the reagent cartridge (Table of Materials) and the libraries. 4.3 Make 0.2 N NaOH and dilute the libraries to obtain the desired molar concentration. 4.4 Rinse and dry the flow cell. Add 600 μL of diluted and denatured library solution into the well of the reagent cartridge. 4.5 Start the sequencing run. 5. Quality Control of NGS Data 5.1. Perform the quality control of FASTQ data using the "FASTX-Toolkit"¹⁶. NOTE: A basic example of the parameter settings used is as follows: fastq_quality_trimmer -v -t 20 -l 200 -i [InFilename.fastq] -o [InFilename.fastq]
214 215 216 217 218 219 220 221 222 223 224 225 226 227 228 229 230	 4.2 Thaw the reagent cartridge (Table of Materials) and the libraries. 4.3 Make 0.2 N NaOH and dilute the libraries to obtain the desired molar concentration. 4.4 Rinse and dry the flow cell. Add 600 μL of diluted and denatured library solution into the well of the reagent cartridge. 4.5 Start the sequencing run. 5. Quality Control of NGS Data 5.1. Perform the quality control of FASTQ data using the "FASTX-Toolkit"¹⁶. NOTE: A basic example of the parameter settings used is as follows: fastq_quality_trimmer -v -t 20 -l 200 -i [InFilename.fastq] -o [InFilename.fastq] fastq_quality_filter -v -q 20 -p 80 -i [InFilename.fastq] -o [InFilename.fastq]
214 215 216 217 218 219 220 221 222 223 224 225 226 227 228 229 230 231	 4.2 Thaw the reagent cartridge (Table of Materials) and the libraries. 4.3 Make 0.2 N NaOH and dilute the libraries to obtain the desired molar concentration. 4.4 Rinse and dry the flow cell. Add 600 μL of diluted and denatured library solution into the well of the reagent cartridge. 4.5 Start the sequencing run. 5. Quality Control of NGS Data 5.1. Perform the quality control of FASTQ data using the "FASTX-Toolkit"¹⁶. NOTE: A basic example of the parameter settings used is as follows: fastq_quality_trimmer -v -t 20 -l 200 -i [InFilename.fastq] -o [InFilename.fastq] fastq_quality_filter -v -q 20 -p 80 -i [InFilename.fastq] -o [InFilename.fastq]
214 215 216 217 218 219 220 221 222 223 224 225 226 227 228 229 230 231 232	 4.2 Thaw the reagent cartridge (Table of Materials) and the libraries. 4.3 Make 0.2 N NaOH and dilute the libraries to obtain the desired molar concentration. 4.4 Rinse and dry the flow cell. Add 600 μL of diluted and denatured library solution into the well of the reagent cartridge. 4.5 Start the sequencing run. 5. Quality Control of NGS Data 5.1. Perform the quality control of FASTQ data using the "FASTX-Toolkit"¹⁶. NOTE: A basic example of the parameter settings used is as follows: fastq_quality_trimmer -v -t 20 -l 200 -i [InFilename.fastq] -o [InFilename.fastq] fastq_quality_filter -v -q 20 -p 80 -i [InFilename.fastq] -o [InFilename.fastq] fastx_reverse_complement -v -i [InFilename.fastq] -o [InFilename.fastq]

6. Extraction and Analysis of Immunoglobulin Sequences from .fna Data

- NOTE: The example programs were implemented in a UNIX environment. Please use them as
- an example references because performance may depend on the operating system and
- hardware environment. The authors do not accept any liability for errors or omissions. The
- programming languages, Perl¹⁷, R¹⁸, and required modules need to be installed according to
- the instructions on the cited websites. the IgBLAST program need to be installed according to
- 243 the instructions on the appropriate website^{19,20}.
- 244

236

- 245 6.1 Download the following examples of in-house programs for repertoire analyses from
- 246 https://github.com/KzPipeLine/KzPipeLine:
- 247 03_PipeLine_Mouse.zip; A set of example programs for the analyses of mouse antibody
- 248 sequences.
- 249 05 PipeLine Human.zip; A set of example programs for the analyses of human antibody
- 250 sequences.
- 251
- 252 6.2 Extract the antibody reads in the sequence data: Extract the immunoglobulin (Ig)
- sequences of each Ig-class from the data (.fna) by a Perl program that searches the signature
- sequences in each immunoglobulin constant region (**Table 2**).
- 255
- 256 6.2.1 For the mouse immunoglobulin heavy chain (IgH) genes, extract the reads by the
- 257 following command:
- \$ perl 01 KzMFTlgCmgggaNtdVer3 Kz160607.pl [Input filename] [Output filename (suffix)]
- 259
- 260 6.2.2 For the mouse immunoglobulin light chain (IgL) genes extract the reads by the following
- 261 command:
- \$ perl 01_KzMFTCkltNtdVer1_170810.pl [Input filename] [Output filename (suffix)]
- 263
- 264 6.2.3 For the human immunoglobulin heavy chain (IgH) genes, extract the reads by the
- 265 following command:
- \$\footnote{\text{y}} \text{ perl } 01_KzMfHulgHCmgadeNtdVer1_Kz180312.pl [Input filename] [Output filename]
- 267 (suffix)]
- 268
- 269 6.2.4 For the human immunoglobulin light chain (IgL) genes, extract the reads by the following
- 270 command:
- 271 \$ perl 01 KzMfHulgCkltNtd 180316.pl [Input filename] [Output filename (suffix)]
- 272
- 273 6.3 Annotate and check the productivity of V(D)J gene recombination:
- 274
- NOTE: The method described below utilizes standalone IgBLAST¹⁹ for the annotation of V(D)J
- 276 gene segments in the sequence. Set the database for the V(D)J genes and the parameter
- 277 settings for IgBLAST as described²⁰.
- 278
- 279 6.3.1 Annotate the mouse immunoglobulin heavy chain (IgH) genes by the following
- 280 command:
- 281 \$ igblastn -germline db V \$IGDATA/ImtgMouselghV NtdDb.txt -germline db J
- 282 \$IGDATA/ImtgMouseIghJ NtdDb.txt -germline db D \$IGDATA/ImtgMouseIghD NtdDb.txt -

- 283 organism mouse -domain_system imgt -query ./\$InFile -auxiliary_data
- 284 \$IGDATA/optional file/mouse gl.aux -show translation -outfmt 7 >> ./\$OutName

285

- 286 6.3.2 Annotate the mouse immunoglobulin light chain (IgL) genes by the following command:
- 287 \$ igblastn -germline_db_V \$IGDATA/ImtgMouseIgkV_NtdDb.txt -germline_db_J
- 288 \$IGDATA/ImtgMouselgkJ NtdDb.txt -germline db D \$IGDATA/ImtgMouselghD NtdDb.txt -
- 289 organism mouse -domain_system imgt -query ./\$InFile -auxiliary_data
- 290 \$IGDATA/optional file/mouse gl.aux -show translation -outfmt 7 >> ./\$OutName

291

- 292 6.3.3 Annotate the human immunoglobulin heavy chain (IgH) genes by the following 293 command:
- 294 \$ igblastn -germline_db_V \$IGDATA/ImtgHumanIghV_NtdDb.txt -germline_db_J
- 295 \$IGDATA/ImtgHumanlghJ_NtdDb.txt -germline_db_D \$IGDATA/ImtgHumanlghD_NtdDb.txt -
- 296 organism human -domain_system imgt -query ./\$InFile -auxiliary_data
- 297 \$IGDATA/optional file/Human gl.aux -show translation -outfmt 7 >> ./\$OutName

298

- 299 6.3.4 Annotate the human immunoglobulin light chain (IgL) genes by the following command:
- 300 \$ igblastn -germline_db_V \$IGDATA/ImtgHumanIgkV_NtdDb.txt -germline_db_J
- 301 \$IGDATA/ImtgHumanlgkJ_NtdDb.txt -germline_db_D \$IGDATA/ImtgHumanlghD_NtdDb.txt -
- 302 organism human -domain system imgt -query ./\$InFile -auxiliary data
- 303 \$IGDATA/optional_file/human_gl.aux -show_translation -outfmt 7 >> ./\$OutName

304

305 6.4 Visualize the global feature of an antibody repertoire.

306

- 307 6.4.1 Visualize the mouse IgH repertoire by the following command:
- 308 \$. 00a1_3DView_MolgH_Kz180406.sh

309

- NOTE: The input file is filename.fna (sequence data), preferably the output file from 6.2.1.
- 311 This file needs to be placed in a lower directory (folder) named "filename". In line 50 of the
- shell script, assign a "filename" for Para_4.

313

- 314 6.4.2 Visualize the human IgH repertoire by the following command:
- 315 \$. 00a1_3DView_HulgH_Kz180411.sh

316

- NOTE: The input file is filename.fna sequence data, preferably the output file of 6.2.3. This
- 318 file needs to be placed in the lower directory (folder) that name is "filename". In line 46 of
- 319 the shell script, assign a "filename" for Para 4.

320

- 321 6.4.3 Visualize the mouse IgL repertoire by the following command:
- 322 \$.00 2DViewS MolgL Kz180406.sh

323

- NOTE: With this pipeline, $Ig\kappa$ and $Ig\lambda$ are processed concomitantly. The input file is
- 325 filename.fna (sequence data), preferably the output file from 6.2.2. This file needs to be
- 326 placed in the lower directory (folder) named "filename". In line 53 of the shell script, assign a
- 327 "filename" for Para_4. The output file's name, ending with "_lgKlCount.txtDim2Rpm.txt"
- gives the coordinates for a two-dimensional bar graph (Figure 3, IgL).

330 6.4.4 Visualize the human IgL repertoire by the following command:

\$.00 2DView HulgL Kz180319.sh

NOTE: With this pipeline, $Ig\kappa$ and $Ig\lambda$ are processed concomitantly. The input file is filename.fna sequence data, preferably the output file of 6.2.4. This file needs to be placed in a lower directory (folder) named "filename". In line 53 of the shell script, assign "filename" for Para_4. The output file name ending with "_IgKlCount.txtDim2Rpm.txt" gives the coordinates for a two-dimensional bar graph (**Figure 4**, IgL).

REPRESENTATIVE RESULTS:

Antibody repertoires of mouse

A perspective of a murine antibody repertoire as a whole can be obtained from cells or tissues such as the spleen, bone marrow, lymph node, or blood. **Figure 3** shows representative results of IgM, IgG1, IgG2c, and immunoglobulin light chain (IgL) repertoires from a naïve mouse spleen. The summary of the read numbers is shown in **Table 3**. For example, 166,175/475,144 reads contained IgM-specific signature sequence (**Table 2**) and 133,371/166,175 reads were VDJ-productive inferred by IgBLAST¹⁹.

Figure 3 shows a repertoire profile of VDJ-rearrangement by 3D-VDJ-plot, in which the size of each ball represents the relative number of reads; in other words, the number of antibody mRNAs in whole B cells. The 3-D mesh consists of 110 IGHV, 12 IGHD, and 4 IGHJ, which are aligned to reflect their order on the chromosome. In addition, the genes ambiguously assigned by IgBLAST were collected separately in the last position for each IGHV, IGHD and IGHJ line, giving rise to 7,215 nodes in the cuboid.

Also, shown in **Figure 3** is a 2D-VJ-plot showing the profile of VJ-rearrangement in the IgL repertoire. The length of each bar on this plot represents the relative number of reads. The x-axis represents 101 IGLV κ and 3 IGLV λ genes, and the y-axis represents 4 IGLJ κ and 3 IGLJ λ genes. The unannotated V- and J-genes are represented on the right borderline.

The complementarity-determining region 3 (CDR3) sequences of these productive reads, which give rise to the majority of antigen-binding specificity, are given in IgBLAST outputs. The CDR3 sequences can be analyzed statistically, including biological or technical replicates, as described previously^{8,10}.

Human antibody repertoires

A perspective of a human antibody repertoire as a whole can be analyzed from various tissues including peripheral blood mononuclear cells (PBMCs) or pathological tissues. **Figure 4** shows representative results of IgM, total IgG (IgG1, IgG2, IgG3, and IgG4), total IgA (IgA1 and IgA2), IgD, IgE and IgL repertoires from normal PBMCs. A summary of the read numbers is shown in **Table 3**. For example, 90,238/1,582,754 reads contained IgM-specific signature sequence and 67,896/90,238 reads were VDJ-productive.

The repertoire profile of VDJ rearrangement is shown on a 3D-VDJ-plot in which the size of each ball represents the relative number of reads; in other words, the number of antibody mRNAs from whole PBMCs (**Figure 4**). The 3-D mesh consists of 56 IGHV, 27 IGHD, and 6 IGHJ,

aligned in the order they appear on the chromosome. In addition, genes ambiguously assigned by IgBLAST are represented separately in the last position for each IGHV, IGHD and IGHJ line, giving rise to 11,172 nodes in the cuboid.

The profile of VJ-rearrangement in the IgL repertoire is depicted in a 2D-VJ-plot in which the length of each bar represents the relative number of reads (**Figure 4**). The x-axis represents 41 IGLV κ and 32 IGLV λ genes, and the y-axis represents 5 IGLJ κ and 5 IGLJ λ genes. The unannotated V- and J-genes are represented on the right borderline.

The human CDR3 sequences are given in IgBLAST outputs and can be analyzed statistically as described previously^{8,10}.

FIGURE AND TABLE LEGENDS:

Table 1: Primer sequences for PCR-amplification of immunoglobulins

Table 2: Summary of the immunoglobulin signature sequences

Table 3: Summary of the read numbers in the experiments

Figure 1: Schematic representation of sequencing strategy for analyzing antibody repertoires in individual mice. (A) Total RNA from the immune cells or tissues was reverse-transcribed and PCR-amplified using the universal forward primer and immunoglobulin class-specific reverse primers. The amplicons from each immunoglobulin class were pooled and rendered for next-generation sequencing. (B) The biological replicates such as spleens from C57BL/6 mice were treated as follows: total RNAs were purified from spleen samples, and cDNAs were amplified by 5'-RACE using the universal primer and antibody class-specific primer. They were then rendered for next-generation sequencing with labeling primers for individual mice. Parts of the figure are adapted from with permission.

Figure 2: Schematic of data-processing flowchart for analyzing antibody repertoires in individual mice. Amplicon reads obtained after next-generation sequencing were processed as follows: (1) read sequences were checked for the presence of antibody class-specific signature sequences; (2) sequences were examined for the V, D, and J gene fragments using IMGT/HighV-Quest and/or IgBLAST; (3) the sequences containing a productive VDJ junction were collected; and (4) these sequences were used for the analysis of overall repertoire features, CDR3, etc.

Figure 3: Global data visualization for mouse antibody repertoires. The overall repertoire profiles of each antibody class were visualized by 3D-VDJ-plot. The x-axis represents 110 IGHV genes ordered as on the chromosome. The y- and z-axis represents 12 IGHD and 4 IGHJ genes, respectively. The volume of spheres on each node represents the number of reads. Red spheres: un-annotated V, D, and J genes. The IgL read distributions are shown on a 2D-VJ-plot in which the length of each bar represents the relative number of reads. The x-axis represents $101 \times IGLV\kappa$ and $3 \times IGLV\lambda$ genes, and the y-axis represents $4 \times IGLJ\kappa$ and $3 \times IGLJ\lambda$ genes. The un-annotated V and J genes are represented on the right borderline.

Figure 4: Global data visualization for human antibody repertoires. The overall repertoire profiles of each antibody class were visualized by 3D-VDJ-plot. The x-axis represents 56 IGHV genes ordered as on the chromosome. The y- and z-axis represents 27 IGHD and 6 IGHJ genes, respectively. The volume of spheres on each node represents the number of reads. Red spheres: un-annotated V, D, and J genes. The IgL reads are arrayed on the 2D-VJ-plot in which the length of each bar represents the relative number of the reads. The x-axis represents 41 x IGLV κ and 32 x IGLV λ genes, and the y-axis represents 5 x IGLJ κ and 5 x IGLJ λ genes. The unannotated V- and J-genes are represented on the right borderline.

DISCUSSION:

The method described here utilizes NGS for antibody RNA amplified using the 5'-RACE method. In contrast to methods that use degenerate 5'-V_H gene primers, mRNAs of each antibody class are amplified evenly using universal forward primers. In addition, the use of antisense primers specific for the constant-region 1 (CH1) of the antibody gene enables repertoire profiling of specific immunoglobulin classes. This is very beneficial for dissecting the class-specific antibody response, as well as for comparing naive and immunized repertoires^{8,9}.

A most likely pitfall of the method is a paucity of amplified immunoglobulin messages. The depth of antibody repertoire obtained by this protocol substantially depends on the PCR amplification described in steps 3.1 and 3.2. If the repertoire depth is not properly obtained, changing the ratios of template cDNA and primers in steps 3.2.1 or 3.2.2 is strongly recommended.

Generally, approximately 20% of the antibody reads produced by NGS are ambiguous sequences²¹. Even with established "correction methods", 5-10% remain ambiguous³. We, therefore, analyzed the sequence and filtered raw reads containing signature sequences corresponding to immunoglobulin constant regions ($C\mu H1$, $C\gamma 1H1$, $C\gamma 2cH1$, etc.). Hence the analysis of somatic hyper-mutations needs the careful examinations.

One of the limitations of this method is that immunoglobulin heavy and light chain pair is unable to be inferred. Hence the repertoire view obtained by this method is not holistic. However, it is possible to approximate the top-ranking pairs by statistical analysis of the data¹⁰. Also, a novel method to sequence the immunoglobulin pairs was reported recently^{3,4}.

The immunoglobulin sequences in the output .fna data were extracted based on the presence of immunoglobulin gene signature sequences. The V, D, and J gene segments were then annotated and the productivity of V(D)J rearrangements were assessed. The complementarity-determining region 3 (CDR3) sequences were also annotated. These systematic examinations of immunoglobulin sequences in .fna data were usefully provided by the IMGT/HighV-QUEST server²²⁻²⁴. However, building an automated processing pipeline has the merit to analyze the big experimental data. The pipeline customized for each purpose is possible to set up by using the standalone IgBLAST protocol¹⁹. This approach needs basic programming literacy but is very useful for detailed analyses of the immunoglobulin system. The pipelines described are the examples of the customized protocol (**Figure 2**).

The number of antibody reads is proportional to the amount of antibody RNAs in the sample,

reflecting the antibody constituents of the antibody system at given time points^{5,8,25}. The method described here gives a bird's eye view of the V(D)J constitution of an antibody repertoire using R programs^{8,18,26}.

 The global view of IgM antibody repertoires of individual naive mice revealed a highly conserved VDJ-profile as compared to those of IgG1 or IgG2c⁸. It was reported that VDJ combinations of immature zebrafish are highly stereotyped⁹. In contrast, human VDJ combinations are reported to be highly skewed⁶. The highly conserved deterministic VDJ-profiles in naive B cells are probably generated either by skewed VDJ-rearrangements or negative selection with auto-antigens presented in the body. For example, IGHV11-2 is expressed preferentially in the fetal IgM repertoire²⁷ and this predominance is attributed to the autoreactivity of IGHV11-2 against senescent erythrocytes²⁷. Interestingly, IGHV11-2 was also the most common major repertoire in our previously published analysis of naive IgM⁸.

 The method described here is useful for deciphering antigen-responsive antibody repertoires by inclusively analyzing the antibody-repertoire space generated in individual bodies, avoiding inadvertent omission of key antibody repertoires^{8,10}. This method also allows the examination of detailed antibody network dynamism, which would facilitate accelerated discovery of protective antibodies against newly emerging pathogens.

ACKNOWLEDGMENTS:

This work was supported by a grant from AMED under Grant Number JP18fk0108011 (KO and SI) and JP18fm0208002 (TS, KO, and YO), and a Grant-in-Aid from the Ministry of Education, Culture, Sports, Science and Technology (15K15159) to KO. We thank Sayuri Yamaguchi and Satoko Sasaki for the valuable technical assistance. We would like to thank Editage (www.editage.jp) for English language editing.

DISCLOSURES:

The authors have no conflicts of interest to disclose.


REFERENCES:


- Schroeder, H. W., Jr. Similarity and divergence in the development and expression of the mouse and human antibody repertoires. *Developmental & Comparative Immunology.* **30** (1-2), 119-135, doi:10.1016/j.dci.2005.06.006 (2006).
- 505 2 Tonegawa, S. Somatic generation of antibody diversity. *Nature.* **302** (5909), 575-581 506 (1983).
- Georgiou, G. *et al.* The promise and challenge of high-throughput sequencing of the antibody repertoire. *Nature Biotechnology.* **32** (2), 158-168, doi:10.1038/nbt.2782 (2014).
- Lees, W. D. & Shepherd, A. J. Studying Antibody Repertoires with Next-Generation Sequencing. *Methods in Molecular Biology.* **1526,** 257-270, doi:10.1007/978-1-4939-6613-4 15 (2017).
- 5 Weinstein, J. A., Jiang, N., White, R. A., 3rd, Fisher, D. S. & Quake, S. R. High-throughput 513 sequencing of the zebrafish antibody repertoire. *Science*. **324** (5928), 807-810, 514 doi:10.1126/science.1170020 (2009).
- Arnaout, R. *et al.* High-resolution description of antibody heavy-chain repertoires in humans. *PLoS One.* **6** (8), e22365, doi:10.1371/journal.pone.0022365 (2011).
- 517 7 Boyd, S. D. *et al.* Individual variation in the germline Ig gene repertoire inferred from


- variable region gene rearrangements. Journal of Immunology. 184 (12), 6986-6992,
- 519 doi:10.4049/jimmunol.1000445 (2010).
- 520 8 Kono, N. et al. Deciphering antigen-responding antibody repertoires by using next-
- 521 generation sequencing and confirming them through antibody-gene synthesis. *Biochemical*
- *and Biophysical Research Communications.* **487** (2), 300-306, doi:10.1016/j.bbrc.2017.04.054
- 523 (2017).
- 524 9 Jiang, N. et al. Determinism and stochasticity during maturation of the zebrafish
- antibody repertoire. Proceedings of the National Academy of Sciences of the United States of
- 526 America. **108** (13), 5348-5353, doi:10.1073/pnas.1014277108 (2011).
- 527 10 Sun, L. et al. Distorted antibody repertoire developed in the absence of pre-B cell
- receptor formation. Biochemical and Biophysical Research Communications. 495 (1), 1411-
- 529 1417, doi:10.1016/j.bbrc.2017.11.171 (2018).
- 530 11 Olivarius, S., Plessy, C. & Carninci, P. High-throughput verification of transcriptional
- starting sites by Deep-RACE. *Biotechniques.* **46** (2), 130-132, doi:10.2144/000113066 (2009).
- 532 12 Yeku, O. & Frohman, M. A. Rapid amplification of cDNA ends (RACE). Methods in
- 533 *Molecular Biology.* **703,** 107-122, doi:10.1007/978-1-59745-248-9_8 (2011).
- 534 13 Zhu, Y. Y., Machleder, E. M., Chenchik, A., Li, R. & Siebert, P. D. Reverse transcriptase
- 535 template switching: a SMART approach for full-length cDNA library construction.
- 536 *Biotechniques.* **30** (4), 892-897 (2001).
- 537 14 Vollmers, C., Sit, R. V., Weinstein, J. A., Dekker, C. L. & Quake, S. R. Genetic
- 538 measurement of memory B-cell recall using antibody repertoire sequencing. Proceedings of
- the National Academy of Sciences of the United States of America. 110 (33), 13463-13468,
- 540 doi:10.1073/pnas.1312146110 (2013).
- 541 15 *SMARTer RACE 5'/3' Kit User Manual (634858, 634859)* (2018).
- 542 16 FASTX-Toolkit, http://hannonlab.cshl.edu/fastx_toolkit/ (2018).
- 543 17 *Perl*, < https://perldoc.perl.org/> (2018).
- 544 18 R: A language and environment for statistical computing. (R Foundation for Statistical
- 545 Computing, Vienna, Austria., 2016).
- 546 19 Ye, J., Ma, N., Madden, T. L. & Ostell, J. M. IgBLAST: an immunoglobulin variable
- domain sequence analysis tool. Nucleic Acids Research. 41 (Web Server issue), W34-40,
- 548 doi:10.1093/nar/gkt382 (2013).
- 549 20 IgBLAST, https://www.ncbi.nlm.nih.gov/igblast/faq.html-standalone (2018).
- 550 21 Prabakaran, P., Streaker, E., Chen, W. & Dimitrov, D. S. 454 antibody sequencing error
- characterization and correction. *BMC Research Notes.* **4,** 404, doi:10.1186/1756-0500-4-404
- 552 (2011).
- 553 22 Lefranc, M. P. et al. IMGT, the international ImMunoGeneTics information system.
- 554 *Nucleic Acids Research.* **37** (Database issue), D1006-1012, doi:10.1093/nar/gkn838 (2009).
- Alamyar, E., Duroux, P., Lefranc, M. P. & Giudicelli, V. IMGT((R)) tools for the
- 556 nucleotide analysis of immunoglobulin (IG) and T cell receptor (TR) V-(D)-J repertoires,
- 557 polymorphisms, and IG mutations: IMGT/V-QUEST and IMGT/HighV-QUEST for NGS. *Methods*
- *in Molecular Biology.* **882** 569-604, doi:10.1007/978-1-61779-842-9 32 (2012).
- 559 24 *IMGT/HighV-QUEST*, < http://www.imgt.org/HighV-QUEST/login.action> (2018).
- 560 25 Glanville, J. et al. Precise determination of the diversity of a combinatorial antibody
- library gives insight into the human immunoglobulin repertoire. *Proceedings of the National*
- 562 Academy of Sciences of the United States of America. 106 (48), 20216-20221,
- 563 doi:10.1073/pnas.0909775106 (2009).
- rgl: 3D Visualization Using OpenGL. R package version 0.95.1247 (2015).


Hardy, R. R., Wei, C. J. & Hayakawa, K. Selection during development of VH11+ B cells: a model for natural autoantibody-producing CD5+ B cells. *Immunological Reviews.* **197** 60-74 (2004).

Page 12 of 6 revised October 2016









Universal forward primer 5'- AAGCAGTGGTATCAACGCAGAGT-3'

Reverse primers for the mouse immunoglobulins (Ref.8)

IgM_CH1:5'- CACCAGATTCTTATCAGACAGGGGGCTCTC -3'IgG1_CH1:5'- CATCCCAGGGTCACCATGGAGTTAGTTTGG -3'IgG2c_CH1:5'- GTACCTCCACACACAGGGGCCAGTGGATAG -3'IgG3_CH1:5'-ATGTGTCACTGCAGCCAGGGACCAAGGGA-3'IgA_CH1:5'-GAATCAGGCAGCCGATTATCACGGGATCAC-3'Igk_CH1:5'- GCTCACTGGATGGTGGGAAGATGGATACAG -3'Igλ_CH1:5'- CTBGAGCTCYTCAGRGGAAGGTGGAAACA -3'

Reverse primers for the human immunoglobulins (Ref.14)

5'- GGGAATTCTCACAGGAGACG -3' IgM CH1: IgG_CH1: 5'- AAGACCGATGGGCCCTTG -3' IgD_CH1: 5'- GGGTGTCTGCACCCTGATA -3' IgA_CH1: 5'- GAAGACCTTGGGGCTGGT -3' IgE1 CH1: 5'- GAAGACGGATGGGCTCTGT -3' IgE2 CH1: 5'- TTGCAGCAGCGGGTCAAGGG -3' 5'- TGCTCATCAGATGGCGGGAAGAT -3' Igk_CH1: Igλ_CH1: 5'- AGAGGAGGGCGGGAACAGAGTGA -3'

Immunoglobulin class	Sense	Antisense
Mouse immunoglobulin he	eavy chains (C57BI /6)	
IgM	AGTCAGTCCTTCCCAAATGTC	GACATTTGGGAAGGACTGACT
lgG1	AAAACGACACCCCATCTGTC	GACAGATGGGGGTGTCGTTTT
(IgG1 variant)	AAAACAACACCCCCATCAGTC	GACTGATGGGGGTGTTTTT
lgG2c	AAAACAACAGCCCCATCGGTC	GACCGATGGGGCTGTTGTTTT
lgG3	GTGATCCCGTGATAATCGGCT	AGCCGATTATCACGGGATCAC
IgA	TCCCTTGGTCCCTGGCTGCAG	TCCCTTGGTCCCTGGCTGCAG
Mouse immunoglobulin lig		
lgк	CTGTATCCATCTTCCCACCATCCAGTGAGC	GCTCACTGGATGGTGGGAAGATGGATACAG
lgλ1	TGTTTCCACCTTCCTCTGAAGAGCTCGAG	CTCGAGCTCTTCAGAGGAAGGTGGAAACA
lgλ2	TGTTTCCACCTTCCTCTGAGGAGCTCAAG	CTTGAGCTCCTCAGAGGAAGGTGGAAACA
lgλ3	TGTTTCCACCTTCCCCTGAGGAGCTCCAG	CTGGAGCTCCTCAGGGGAAGGTGGAAACA
lgλ4	TGTTCCCACCTTCCTCTGAAGAGCTCAAG	CTTGAGCTCTTCAGAGGAAGGTGGGAACA
Human immunoglobulin he	eavy chains	
IgM	GGGAGTGCATCCGCCCAAC	GTTGGGGCGGATGCACTCCC
IgG	GCTTCCACCAAGGGCCCATC	GATGGGCCCTTGGTGGAAGC
IgA	GCATCCCGACCAGCCCCAA	GACCGATGGGGCTGTTGTTTT
IgD	GCACCCACCAAGGCTCCGGA	TCCGGAGCCTTGGTGGGTGC
lgE	GCCTCCACACAGAGCCCATC	GATGGGCTCTGTGTGGAGGC
Human immunoglobulin lig	ght chains	
lgк	ACTGTGGCTGCACCATCTGC	GCAGATGGTGCAGCCACAGT
lgλ1,2,6	GTCACTCTGTTCCCGCCCTC	GAGGGCGGAACAGAGTGAC
Igλ3,7	GTCACTCTGTTCCCACCCTC	GAGGGTGGGAACAGAGTGAC

Mouse IgH	Total reads	IgM	lgG1	IgG2c
Input	475,144			_
IgC-containing		166,175	229,671	36,628
VDJ-productive		133,371	196,583	31,446
				<u> </u>

Mouse IgL	Total reads	IgKappa	IgLambda
Input	527,668		_
IgC-containing		178,948	21,446
VJ-productive		160,924	16,988

Human IgH	Total reads	lgM	IgG	lgA	lgD	lgE
Input	1,582,754					
IgC-containing		90,238	5,298	94,061	75,549	2,932
VDJ-productive		67,896	2,775	78,203	56,495	3

Human IgL	Total reads	IgKappa	IgLambda
Input	1,582,754		
IgC-containing		120,316	64,148
VJ-productive		97,169	52,324

Name of Material/ Equipment	Company	Catalog Number	Comments/Description
0.2 mL Strip Tubes	Thermo Fisher Scientific	AB0452	120 strips
100 bp DNA Ladder	ТОУОВО	DNA-035	0.5 mL
2100 Bioanalyzer Systems	Agilent Technologies	G2939BA /2100	
Acetic Acid	Wako	017-00256	500 mL
Agarose, NuSieve GTG	Lonza	50084	
Ammonium Chloride	Wako	017-02995	500 g
Chloroform	Wako	038-02606	500 mL
Dulbecco's PBS (-)"Nissui"	NISSUI	08192	
Ethylenediamine-N,N,N',N'-tetraacetic Acid Di	sc Wako	345-01865	500 g
Falcon 40 µm Cell Strainer	Falcon	352340	50/Case
ling lock tube 1.7 mL	BM EQUIPMENT	BM-15	
ling lock tube 2.0 mL	BM EQUIPMENT	BM-20	
MiSeq Reagent Kit v2	illumina	MS-102-2003	500 cycles
MiSeq System	illumina	SY-410-1003	
NanoDrop 2000c Spectrophotometer	Thermo Fisher Scientific		
Potassium Hydrogen Carbonate	Wako	166-03275	500 g
PureLink RNA Mini Kit	life technologies	12183018A	
Qubit 3.0 Fluorometer	Thermo Fisher Scientific	Q33216	
Qubit dsDNA HS Assay Kit	Thermo Fisher Scientific	Q32854	500 assays
SMARTer RACE 5'/3' Kit	Clontech	634858	
TaKaRa Ex Taq Hot Start Version	Takara Bio Inc.	RR006A	
Trizma base	Sigma	T6066	1 kg
TRIzol Reagent	AmbionThermo Fisher Scientific	15596026	100 mL
Ultra Clear qPCR Caps	Thermo Fisher Scientific	AB0866	120 strips
UltraPure Ethidium Bromide	Thermo Fisher Scientific	15585011	
Wizard SV Gel and PCR Clean-Up System	Promega	A9282	



ARTICLE AND VIDEO LICENSE AGREEMENT

Title of Article:

Holistic Visualization and Analysis of Mouse and Human Antibody Repertoires by Next-generation Sequencing

Author(s):

Lin Sun, Naoko Kono, Hiroyuki Toh, Hanbing Xue, Kaori Sano, Tadaki Suzuki, Akira Ainai, Yasuko Orba, Junya Yamagishi, Hideki Hasegawa, Yoshimasa Takahashi, Shiqeyuki Itamura, Kazuo Ohnishi

Item 1: The Author elects to have the Materials be made available (as described at http://www.jove.com/publish) via:

Standard Access

V Open Access

Item 2: Please select one of the following items:

V The Author is NOT a United States government employee.

The Author is a United States government employee and the Materials were prepared in the course of his or her duties as a United States government employee.

The Author is a United States government employee but the Materials were NOT prepared in the course of his or her duties as a United States government employee.

ARTICLE AND VIDEO LICENSE AGREEMENT

- 1. Defined Terms. As used in this Article and Video License Agreement, the following terms shall have the following meanings: "Agreement" means this Article and Video License Agreement; "Article" means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein; "Author" means the author who is a signatory to this Agreement; "Collective Work" means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; "CRC License" means the Creative Commons Attribution-Non Commercial-No Derivs 3.0 Unported Agreement, the terms and conditions of which can be found at: http://creativecommons.org/licenses/by-nc-
- nd/3.0/legalcode; "Derivative Work" means a work based upon the Materials or upon the Materials and other preexisting works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; "Institution" means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; "JoVE" means MyJove Corporation, a Massachusetts corporation and the publisher of The Journal of Visualized Experiments; "Materials" means the Article and / or the Video; "Parties" means the Author and JoVE; "Video" means any video(s) made by the Author, alone or in conjunction with any other parties, or by JoVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion

- of the Article, and in which the Author may or may not appear.
- 2. **Background.** The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JoVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.
- Grant of Rights in Article. In consideration of JoVE agreeing to publish the Article, the Author hereby grants to JoVE, subject to Sections 4 and 7 below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and(c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the "Open Access" box has been checked in Item 1 above, JoVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.



ARTICLE AND VIDEO LICENSE AGREEMENT

- 4. **Retention of Rights in Article.** Notwithstanding the exclusive license granted to JoVE in **Section 3** above, the Author shall, with respect to the Article, retain the non-exclusive right to use all or part of the Article for the non-commercial purpose of giving lectures, presentations or teaching classes, and to post a copy of the Article on the Institution's website or the Author's personal website, in each case provided that a link to the Article on the JoVE website is provided and notice of JoVE's copyright in the Article is included. All non-copyright intellectual property rights in and to the Article, such as patent rights, shall remain with the Author.
- 5. **Grant of Rights in Video Standard Access.** This **Section 5** applies if the "Standard Access" box has been checked in **Item 1** above or if no box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby acknowledges and agrees that, Subject to **Section 7** below, JoVE is and shall be the sole and exclusive owner of all rights of any nature, including, without limitation, all copyrights, in and to the Video. To the extent that, by law, the Author is deemed, now or at any time in the future, to have any rights of any nature in or to the Video, the Author hereby disclaims all such rights and transfers all such rights to JoVE.
- 6. Grant of Rights in Video - Open Access. This Section 6 applies only if the "Open Access" box has been checked in Item 1 above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby grants to JoVE, subject to Section 7 below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Video in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Video into other languages, create adaptations, summaries or extracts of the Video or other Derivative Works or Collective Works based on all or any portion of the Video and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. For any Video to which this **Section 6** is applicable, JoVE and the Author hereby grant to the public all such rights in the Video as provided in, but subject to all limitations and requirements set forth in, the CRC License.
- 7. **Government Employees.** If the Author is a United States government employee and the Article was prepared in the course of his or her duties as a United States government employee, as indicated in **Item 2** above, and any of the licenses or grants granted by the Author hereunder exceed the scope of the 17 U.S.C. 403, then the rights granted hereunder shall be limited to the maximum

- rights permitted under such statute. In such case, all provisions contained herein that are not in conflict with such statute shall remain in full force and effect, and all provisions contained herein that do so conflict shall be deemed to be amended so as to provide to JoVE the maximum rights permissible within such statute.
- 8. **Protection of the Work.** The Author(s) authorize JoVE to take steps in the Author(s) name and on their behalf if JoVE believes some third party could be infringing or might infringe the copyright of either the Author's Article and/or Video.
- 9. **Likeness, Privacy, Personality.** The Author hereby grants JoVE the right to use the Author's name, voice, likeness, picture, photograph, image, biography and performance in any way, commercial or otherwise, in connection with the Materials and the sale, promotion and distribution thereof. The Author hereby waives any and all rights he or she may have, relating to his or her appearance in the Video or otherwise relating to the Materials, under all applicable privacy, likeness, personality or similar laws.
- Author Warranties. The Author represents and warrants that the Article is original, that it has not been published, that the copyright interest is owned by the Author (or, if more than one author is listed at the beginning of this Agreement, by such authors collectively) and has not been assigned, licensed, or otherwise transferred to any other party. The Author represents and warrants that the author(s) listed at the top of this Agreement are the only authors of the Materials. If more than one author is listed at the top of this Agreement and if any such author has not entered into a separate Article and Video License Agreement with JoVE relating to the Materials, the Author represents and warrants that the Author has been authorized by each of the other such authors to execute this Agreement on his or her behalf and to bind him or her with respect to the terms of this Agreement as if each of them had been a party hereto as an Author. The Author warrants that the use, reproduction, distribution, public or private performance or display, and/or modification of all or any portion of the Materials does not and will not violate, infringe and/or misappropriate the patent, trademark, intellectual property or other rights of any third party. The Author represents and warrants that it has and will continue to comply with all government, institutional and other regulations, including, without limitation all institutional, laboratory, hospital, ethical, human and animal treatment, privacy, and all other rules, regulations, laws, procedures or guidelines, applicable to the Materials, and that all research involving human and animal subjects has been approved by the Author's relevant institutional review board.
- 11. **JoVE Discretion.** If the Author requests the assistance of JoVE in producing the Video in the Author's facility, the Author shall ensure that the presence of JoVE employees, agents or independent contractors is in accordance with the relevant regulations of the Author's institution. If more than one author is listed at the beginning of this Agreement, JoVE may, in its sole



ARTICLE AND VIDEO LICENSE AGREEMENT

discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JoVE reserves the right, in its absolute and sole discretion and without giving any reason therefore, to accept or decline any work submitted to JoVE. JoVE and its employees, agents and independent contractors shall have full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JoVE has sole discretion as to the method of making and publishing the Materials, including, without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

Indemnification. The Author agrees to indemnify JoVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JoVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JoVE, making of videos by JoVE, or publication in JoVE or elsewhere by JoVE. The Author shall be responsible for, and shall hold JoVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to the making of a video by JoVE its employees, agents or independent contractors. All sterilization, cleanliness or decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's expense. All indemnifications provided herein shall include JoVE's attorney's fees and costs related to said losses or damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JoVE, its employees, agents or independent contractors.

- 13. **Fees.** To cover the cost incurred for publication, JoVE must receive payment before production and publication the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JoVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.
- 14. **Transfer, Governing Law.** This Agreement may be assigned by JoVE and shall inure to the benefits of any of JoVE's successors and assignees. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to me one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement is required per submission.

CORRESPONDING AUTHOR

Name:

Kazuo OHNISHI

Department:

Department of Immunology

Institution:

National Institute of Infectious Diseases

Title:

Senior Scientist, Ph.D.

Signature:

Date:

July 19, 2018

Please submit a **signed** and **dated** copy of this license by one of the following three methods:

- 1. Upload an electronic version on the JoVE submission site
- 2. Fax the document to +1.866.381.2236
- 3. Mail the document to JoVE / Attn: JoVE Editorial / 1 Alewife Center #200 / Cambridge, MA 02140

Responses to the Editorial comments:

Editorial comments:

1. The editor has formatted the manuscript to match the journal's style. Please retain the same.

Response: We retained the format made in the manuscript.

2. Please address specific comments marked in the manuscript.

Response: We addressed the specific comments as follows:

Comment[1]: We agree with this comment.

Comment[2]: We agree with this comment.

Comment[3]: We agree with this comment.

Comment[4]: The word 'female' was removed. The word 'pass' was replaced by 'smash'.

Comment[5]: The word 'in the tissue' was added.

Comment[6]: The word 'tissue' was added.

Comment[7]: A step describing RNA quantitation was added as step 2.9.

Comment[8]: We explained about the human sample in Note, line 169-171.

Comment[9]: The words 'Illumina index primers' were changed to 'index sequence primers'.

Comment[10]: The step was expanded to explain this procedure (3.3.1 - 3.3.4).

Comment[11]: The number of immunoglobulin classes amplified and the yields were explained in the Notes added.

Comment[12]: The step was explained further.

Comment[13]: The explanation was expanded.

Comment[14]: The step was explained further.

Comment[15]: We understand this comment.

Comment[16]: Table 2 (signature sequences) was referred.

Comment[17]: IgBLAST was referred.

Comment[18]: Some troubleshooting and limitations were added.

3. We cannot film scripting steps. We can film steps involving graphical user interface,

button clicks, knob turns, discrete experimental steps only.

- # Response: We understand this and changed the highlighted part in the manuscript.
- 4. if any of the figures are reused from previous publications, please obtain explicit copyright permission and upload it to your editorial manager account.
- # Response: We obtained the explicit copyright permission for the Figure 1 and 2.
- 5. Please reword lines 444 and line 446-447 as it matches with the previously published literature.
- # Response: The lines 444 and line 446-447 were reworded.

ELSEVIER LICENSE TERMS AND CONDITIONS

Sep 25, 2018

This Agreement between Prof. Kazuo OHNISHI ("You") and Elsevier ("Elsevier") consists of your license details and the terms and conditions provided by Elsevier and Copyright Clearance Center.

License Number 4435791402085 License date Sep 25, 2018

Licensed Content Publisher Elsevier

Licensed Content Publication Biochemical and Biophysical Research Communications

Licensed Content Title Deciphering antigen-responding antibody repertoires by using next-

generation sequencing and confirming them through antibody-gene

synthesis

Licensed Content Author Naoko Kono, Lin Sun, Hiroyuki Toh, Takeyuki Shimizu, Hanbing

Xue, Osamu Numata, Manabu Ato, Kazuo Ohnishi, Shigeyuki Itamura

Licensed Content Date May 27, 2017

Licensed Content Volume 487
Licensed Content Issue 2
Licensed Content Pages 7
Start Page 300

Type of Use reuse in a journal/magazine

306

Requestor type author of new work

Intended publisher of new

work

End Page

Other

Portion figures/tables/illustrations

Number of figures/tables

/illustrations

2

Format both print and electronic

Are you the author of this

Elsevier article?

Yes

Will you be translating? No

Original figure numbers Figures 1A and 1C

Title of the article Extensive Analysis and Visualization of Mouse and Human Antibody

Repertoires by Next-generation Sequencing

Publication new article is in Journal of Visualized Experiments

Publisher of the new article Journal of Visualized Experiments

Author of new article Lin Sun, Naoko Kono, Hiroyuki Toh, Hanbing Xue, Kaori Sano, Tadaki

Suzuki, Akira Ainai, Yasuko Orba, Junya Yamagishi, Hideki

Hasegawa, Yoshimasa Takahashi, Shigeyuki Itamura, Kazuo Ohnishi

Expected publication date Nov 2018

Estimated size of new article 5

(number of pages)

Requestor Location Prof. Kazuo OHNISHI

NIID

1-23-1 Toyama

Shinjuku, Tokyo 162-8640

Japan

Attn: Prof. Kazuo OHNISHI

Publisher Tax ID JP00022
Total 0 JPY

Terms and Conditions

INTRODUCTION

1. The publisher for this copyrighted material is Elsevier. By clicking "accept" in connection with completing this licensing transaction, you agree that the following terms and conditions apply to this transaction (along with the Billing and Payment terms and conditions established by Copyright Clearance Center, Inc. ("CCC"), at the time that you opened your Rightslink account and that are available at any time at http://myaccount.copyright.com).

GENERAL TERMS

- 2. Elsevier hereby grants you permission to reproduce the aforementioned material subject to the terms and conditions indicated.
- 3. Acknowledgement: If any part of the material to be used (for example, figures) has appeared in our publication with credit or acknowledgement to another source, permission must also be sought from that source. If such permission is not obtained then that material may not be included in your publication/copies. Suitable acknowledgement to the source must be made, either as a footnote or in a reference list at the end of your publication, as follows:
- "Reprinted from Publication title, Vol /edition number, Author(s), Title of article / title of chapter, Pages No., Copyright (Year), with permission from Elsevier [OR APPLICABLE SOCIETY COPYRIGHT OWNER]." Also Lancet special credit "Reprinted from The Lancet, Vol. number, Author(s), Title of article, Pages No., Copyright (Year), with permission from Elsevier."
- 4. Reproduction of this material is confined to the purpose and/or media for which permission is hereby given.
- 5. Altering/Modifying Material: Not Permitted. However figures and illustrations may be altered/adapted minimally to serve your work. Any other abbreviations, additions, deletions and/or any other alterations shall

be made only with prior written authorization of Elsevier Ltd. (Please contact Elsevier at permissions@elsevier.com). No modifications can be made to any Lancet figures/tables and they must be reproduced in full.
6. If the permission fee for the requested use of our material is waived in this instance, please be advised that your future requests for Elsevier materials may attract a fee.

- 7. Reservation of Rights: Publisher reserves all rights not specifically granted in the combination of (i) the license details provided by you and accepted in the course of this licensing transaction, (ii) these terms and conditions and (iii) CCC's Billing and Payment terms and conditions.
- 8. License Contingent Upon Payment: While you may exercise the rights licensed immediately upon issuance of the license at the end of the licensing process for the transaction, provided that you have disclosed complete and accurate details of your proposed use, no license is finally effective unless and until full payment is received from you (either by publisher or by CCC) as provided in CCC's Billing and Payment terms and conditions. If full payment is not received on a timely basis, then any license preliminarily granted shall be deemed automatically revoked and shall be void as if never granted. Further, in the event that you breach any of these terms and conditions or any of CCC's Billing and Payment terms and conditions, the license is automatically revoked and shall be void as if never granted. Use of materials as described in a revoked license, as well as any use of the materials beyond the scope of an unrevoked license, may constitute copyright infringement and publisher reserves the right to take any and all action to protect its copyright in the materials.
- 9. Warranties: Publisher makes no representations or warranties with respect to the licensed material.
- 10. Indemnity: You hereby indemnify and agree to hold harmless publisher and CCC, and their respective officers, directors, employees and agents, from and against any and all claims arising out of your use of the licensed material other than as specifically authorized pursuant to this license.
- 11. No Transfer of License: This license is personal to you and may not be sublicensed, assigned, or transferred by you to any other person without publisher's written permission.
- 12. No Amendment Except in Writing: This license may not be amended except in a writing signed by both parties (or, in the case of publisher, by CCC on publisher's behalf).
- 13. Objection to Contrary Terms: Publisher hereby objects to any terms

contained in any purchase order, acknowledgment, check endorsement or other writing prepared by you, which terms are inconsistent with these terms and conditions or CCC's Billing and Payment terms and conditions. These terms and conditions, together with CCC's Billing and Payment terms and conditions (which are incorporated herein), comprise the entire agreement between you and publisher (and CCC) concerning this licensing transaction. In the event of any conflict between your obligations established by these terms and conditions and those established by CCC's Billing and Payment terms and conditions, these terms and conditions shall control.

14. Revocation: Elsevier or Copyright Clearance Center may deny the permissions described in this License at their sole discretion, for any reason or no reason, with a full refund payable to you. Notice of such denial will be made using the contact information provided by you. Failure to receive such notice will not alter or invalidate the denial. In no event will Elsevier or Copyright Clearance Center be responsible or liable for any costs, expenses or damage incurred by you as a result of a denial of your permission request, other than a refund of the amount(s) paid by you to Elsevier and/or Copyright Clearance Center for denied permissions.

LIMITED LICENSE

The following terms and conditions apply only to specific license types: 15. **Translation**: This permission is granted for non-exclusive world **English** rights only unless your license was granted for translation rights. If you licensed translation rights you may only translate this content into the languages you requested. A professional translator must perform all translations and reproduce the content word for word preserving the integrity of the article.

16. Posting licensed content on any Website: The following terms and conditions apply as follows: Licensing material from an Elsevier journal: All content posted to the web site must maintain the copyright information line on the bottom of each image; A hyper-text must be included to the Homepage of the journal from which you are licensing at http://www.sciencedirect.com/science/journal/xxxxx or the Elsevier homepage for books at http://www.elsevier.com; Central Storage: This license does not include permission for a scanned version of the material to be stored in a central repository such as that provided by Heron/XanEdu.

Licensing material from an Elsevier book: A hyper-text link must be

included to the Elsevier homepage at http://www.elsevier.com. All content posted to the web site must maintain the copyright information line on the bottom of each image.

Posting licensed content on Electronic reserve: In addition to the above the following clauses are applicable: The web site must be password-protected and made available only to bona fide students registered on a relevant course. This permission is granted for 1 year only. You may obtain a new license for future website posting.

17. **For journal authors:** the following clauses are applicable in addition to the above:

Preprints:

A preprint is an author's own write-up of research results and analysis, it has not been peer-reviewed, nor has it had any other value added to it by a publisher (such as formatting, copyright, technical enhancement etc.). Authors can share their preprints anywhere at any time. Preprints should not be added to or enhanced in any way in order to appear more like, or to substitute for, the final versions of articles however authors can update their preprints on arXiv or RePEc with their Accepted Author Manuscript (see below).

If accepted for publication, we encourage authors to link from the preprint to their formal publication via its DOI. Millions of researchers have access to the formal publications on ScienceDirect, and so links will help users to find, access, cite and use the best available version. Please note that Cell Press, The Lancet and some society-owned have different preprint policies. Information on these policies is available on the journal homepage. **Accepted Author Manuscripts:** An accepted author manuscript is the manuscript of an article that has been accepted for publication and which typically includes author-incorporated changes suggested during submission, peer review and editor-author communications. Authors can share their accepted author manuscript:

immediately

- via their non-commercial person homepage or blog
- by updating a preprint in arXiv or RePEc with the accepted manuscript
- via their research institute or institutional repository for internal institutional uses or as part of an invitation-only research collaboration work-group

- directly by providing copies to their students or to research collaborators for their personal use
- for private scholarly sharing as part of an invitation-only work group on commercial sites with which Elsevier has an agreement
- After the embargo period
 - via non-commercial hosting platforms such as their institutional repository
 - o via commercial sites with which Elsevier has an agreement

In all cases accepted manuscripts should:

- link to the formal publication via its DOI
- bear a CC-BY-NC-ND license this is easy to do
- if aggregated with other manuscripts, for example in a repository or other site, be shared in alignment with our hosting policy not be added to or enhanced in any way to appear more like, or to substitute for, the published journal article.

Published journal article (JPA): A published journal article (PJA) is the definitive final record of published research that appears or will appear in the journal and embodies all value-adding publishing activities including peer review co-ordination, copy-editing, formatting, (if relevant) pagination and online enrichment.

Policies for sharing publishing journal articles differ for subscription and gold open access articles:

<u>Subscription Articles:</u> If you are an author, please share a link to your article rather than the full-text. Millions of researchers have access to the formal publications on ScienceDirect, and so links will help your users to find, access, cite, and use the best available version.

Theses and dissertations which contain embedded PJAs as part of the formal submission can be posted publicly by the awarding institution with DOI links back to the formal publications on ScienceDirect.

If you are affiliated with a library that subscribes to ScienceDirect you have additional private sharing rights for others' research accessed under that agreement. This includes use for classroom teaching and internal training at the institution (including use in course packs and courseware programs), and inclusion of the article for grant funding purposes.

<u>Gold Open Access Articles:</u> May be shared according to the authorselected end-user license and should contain a <u>CrossMark logo</u>, the end

user license, and a DOI link to the formal publication on ScienceDirect. Please refer to Elsevier's posting policy for further information.

- 18. For book authors the following clauses are applicable in addition to the above: Authors are permitted to place a brief summary of their work online only. You are not allowed to download and post the published electronic version of your chapter, nor may you scan the printed edition to create an electronic version. Posting to a repository: Authors are permitted to post a summary of their chapter only in their institution's repository.
- 19. **Thesis/Dissertation**: If your license is for use in a thesis/dissertation your thesis may be submitted to your institution in either print or electronic form. Should your thesis be published commercially, please reapply for permission. These requirements include permission for the Library and Archives of Canada to supply single copies, on demand, of the complete thesis and include permission for Proquest/UMI to supply single copies, on demand, of the complete thesis. Should your thesis be published commercially, please reapply for permission. Theses and dissertations which contain embedded PJAs as part of the formal submission can be posted publicly by the awarding institution with DOI links back to the formal publications on ScienceDirect.

Elsevier Open Access Terms and Conditions

You can publish open access with Elsevier in hundreds of open access journals or in nearly 2000 established subscription journals that support open access publishing. Permitted third party re-use of these open access articles is defined by the author's choice of Creative Commons user license. See our open access license policy for more information.

Terms & Conditions applicable to all Open Access articles published with Elsevier:

Any reuse of the article must not represent the author as endorsing the adaptation of the article nor should the article be modified in such a way as to damage the author's honour or reputation. If any changes have been made, such changes must be clearly indicated.

The author(s) must be appropriately credited and we ask that you include the end user license and a DOI link to the formal publication on ScienceDirect.

If any part of the material to be used (for example, figures) has appeared in our publication with credit or acknowledgement to another source it is

the responsibility of the user to ensure their reuse complies with the terms and conditions determined by the rights holder.

Additional Terms & Conditions applicable to each Creative Commons user license:

CC BY: The CC-BY license allows users to copy, to create extracts, abstracts and new works from the Article, to alter and revise the Article and to make commercial use of the Article (including reuse and/or resale of the Article by commercial entities), provided the user gives appropriate credit (with a link to the formal publication through the relevant DOI), provides a link to the license, indicates if changes were made and the licensor is not represented as endorsing the use made of the work. The full details of the license are available at http://creativecommons.org/licenses/by/4.0.

CC BY NC SA: The CC BY-NC-SA license allows users to copy, to create extracts, abstracts and new works from the Article, to alter and revise the Article, provided this is not done for commercial purposes, and that the user gives appropriate credit (with a link to the formal publication through the relevant DOI), provides a link to the license, indicates if changes were made and the licensor is not represented as endorsing the use made of the work. Further, any new works must be made available on the same conditions. The full details of the license are available at http://creativecommons.org/licenses/by-nc-sa/4.0.

CC BY NC ND: The CC BY-NC-ND license allows users to copy and distribute the Article, provided this is not done for commercial purposes and further does not permit distribution of the Article if it is changed or edited in any way, and provided the user gives appropriate credit (with a link to the formal publication through the relevant DOI), provides a link to the license, and that the licensor is not represented as endorsing the use made of the work. The full details of the license are available at http://creativecommons.org/licenses/by-nc-nd/4.0. Any commercial reuse of Open Access articles published with a CC BY NC SA or CC BY NC ND license requires permission from Elsevier and will be subject to a fee. Commercial reuse includes:

- Associating advertising with the full text of the Article
- Charging fees for document delivery or access
- Article aggregation
- Systematic distribution via e-mail lists or share buttons

Posting or linking by commercial companies for use by customers of those companies.

20. Other Conditions:

v1.9

Questions? $\underline{\text{customercare@copyright.com}}$ or +1-855-239-3415 (toll free in the US) or +1-978-646-2777.