

Journal of Visualized Experiments

Single liposome measurements for the study of proton-pumping membrane enzymes using a combination of electrochemistry and fluorescent microscopy --Manuscript Draft--

Article Type:	Invited Methods Article - JoVE Produced Video
Manuscript Number:	JoVE58896R2
Full Title:	Single liposome measurements for the study of proton-pumping membrane enzymes using a combination of electrochemistry and fluorescent microscopy
Keywords:	cytochrome bo3; single enzyme; proteoliposomes; pH-sensitive dye; proton translocation; bioelectrochemistry
Corresponding Author:	Lars Jeuken, Ph.D. UNITED KINGDOM
Corresponding Author's Institution:	
Corresponding Author E-Mail:	L.J.C.Jeuken@leeds.ac.uk
Order of Authors:	levgen Mazurenko Nikos Hatzakis Lars J C Jeuken, Ph.D.
Additional Information:	
Question	Response
Please indicate whether this article will be Standard Access or Open Access.	Open Access (US\$4,200)
Please indicate the city, state/province, and country where this article will be filmed . Please do not use abbreviations.	Leeds, United Kingdom

School of Biomedical Sciences



UNIVERSITY OF LEEDS

Dear Jialan Zhang, Editor,

Please see included our second revision of the JoVE article "Single liposome measurements for the study of proton-pumping membrane enzymes using a combination of electrochemistry and fluorescent microscopy". As mentioned in the cover letter of our first submission, our method can be applied to a number of proton and ion-pumping membrane proteins and hence will be of interest to a large biochemical community working on these enzymes. As the technique is novel and multidisciplinary, it requires a specialized journal like JoVE to present the method in sufficient detail to allow colleagues around the world to take up this method in their future endeavours.

In this revision cycle, the editors requested a number of minor changes. Comments by the editors were made directly into Word document. In a separately uploaded document, we have answered the editors' questions using the 'Reply to comments' feature in Word. In the document, we have used tracked-changes to indicate any changes made in the text. We include an 'Answer to comments' document to summarize the changes we have made.

Finally, we have explained before that the commercial software package 'Matlab' is integral to the Methods paper and that the submitted software code is specifically written for Matlab and hence this code cannot be used with any other software without completely recoding the script. Matlab is a frequently used package and known throughout the community. Nonetheless, we were again requested to remove references to the software package 'Matlab' (in spite of being used in JoVE before). We have now removed 'Matlab' from the manuscript, but note that the latter will make the manuscript LESS user friendly as a reader will have to confirm in the Table of Materials that the software code is written in Matlab.

We hope that we hereby have addressed all editors' comments. We look forward to hearing from you.

Warm regards,

Lars Jeuken
Yours sincerely,



Lars Jeuken
School of Biomedical Sciences
University of Leeds
LS2 9JT Leeds
United Kingdom

TITLE:

Single Liposome Measurements for the Study of Proton-Pumping Membrane Enzymes Using Electrochemistry and Fluorescent Microscopy

AUTHORS AND AFFILIATIONS:

Ievgen Mazurenko,¹ Nikos S. Hatzakis,² Lars J. C. Jeuken¹

¹ School of Biomedical Sciences & the Astbury Centre for Structural Molecular Biology, University of Leeds, Leeds, U.K.

²Department of Chemistry and Nano-Science Center, University of Copenhagen, Copenhagen, Denmark

Corresponding Author:

Lars J. C. Jeuken (l.j.c.jeuken@leeds.ac.uk)

Email Addresses of Co-Authors:

Nikos S. Hatzakis (hatzakis@nano.ku.dk)

Ievgen Mazurenko (i.mazurenko@leeds.ac.uk)

KEYWORDS:

cytochrome *bo*₃; single enzyme; proteoliposomes; pH-sensitive dye; proton translocation; bioelectrochemistry

SUMMARY:

Here, we present a protocol to study the molecular mechanism of proton translocation across lipid membranes of single liposomes, using cytochrome *bo*₃ as an example. Combining electrochemistry and fluorescence microscopy, pH changes in the lumen of single vesicles, containing single or multiple enzyme, can be detected and analyzed individually.

ABSTRACT:

Proton-pumping enzymes of electron transfer chains couple redox reactions to proton translocation across the membrane, creating a proton-motive force used for ATP production. The amphiphilic nature of membrane proteins requires particular attention to their handling, and reconstitution into the natural lipid environment is indispensable when studying membrane transport processes like proton translocation. Here, we detail a method that has been used for the investigation of the proton-pumping mechanism of membrane redox enzymes, taking cytochrome *bo*₃ from *Escherichia coli* as an example. A combination of electrochemistry and fluorescence microscopy is used to control the redox state of the quinone pool and monitor pH changes in the lumen. Due to the spatial resolution of fluorescent microscopy, hundreds of liposomes can be measured simultaneously while the enzyme content can be scaled down to a single enzyme or transporter per liposome. The respective single enzyme analysis can reveal patterns in the enzyme functional dynamics that might be otherwise hidden by the behavior of the whole population. We include a description of a script for automated image analysis.

INTRODUCTION:

Information about enzyme mechanisms and kinetics is usually obtained on the ensemble or macroscale level with enzyme population in the thousands to millions of molecules, where measurements represent a statistical average. It is known, however, that complex macromolecules such as enzymes may demonstrate heterogeneity in their behavior and molecular mechanisms observed at the ensemble level are not necessarily valid for every molecule. Such deviations on the individual molecule scale have been extensively confirmed by studies of single enzymes with a variety of methods emerging during the last two decades¹. Notably, fluorescence detection of individual enzyme activity has been used to investigate heterogeneity of enzymes activity^{2,3} or discover the so-called memory effect (periods of high enzymes activity succeeded by periods of low activity and *vice versa*)^{4,5}.

Many single enzyme studies require that the enzymes are immobilized on the surface or spatially fixed in another way to remain sufficiently long in the field of view for continuous observation. Enzyme encapsulation into liposomes has been shown to enable enzyme immobilization while preventing any negative impact due the surface-enzyme or protein-protein interactions^{6,7}. In addition, liposomes offer a unique possibility to study single membrane proteins in their natural lipid bilayer environment⁸⁻¹⁰.

A class of membrane proteins, transporters, exercises a directional translocation of substances across the cell membrane, a behavior that can only be studied when proteins are reconstituted into the lipid bilayers (*e.g.*, liposomes)¹¹⁻¹³. For example, proton translocation, exhibited by several enzymes of prokaryotic and eukaryotic electron transport chains, plays an important role in cellular respiration by creating a proton-motive force used for ATP synthesis. In this case, the proton pumping activity is coupled to the electron transfer, although the detailed mechanism of this process often remains elusive.

Recently, we demonstrated the possibility to couple fluorescent detection with electrochemistry to study proton pumping activity of single enzymes of the terminal ubiquinol oxidase of *Escherichia coli* (cytochrome *bo*₃) reconstituted in the liposomes¹⁴. This was achieved by encapsulation of a pH-sensitive membrane-impermeable fluorescent dye into the lumen of liposomes prepared from *E. coli* polar lipids (**Figure 1A**). The protein amount was optimized so that most liposomes either contained no or only one reconstituted enzyme molecule (according to Poisson distribution). The two substrates of cytochrome *bo*₃ were provided by adding ubiquinone to the lipid mix that formed the liposomes and (ambient) oxygen in solution. The liposomes are then sparsely adsorbed on a semi-transparent ultra-smooth gold electrode, covered with a self-assembled monolayer of 6-mercaptohexanol. Finally, the electrode is mounted on the bottom of a simple spectroelectrochemical cell (**Figure 1B**). Electrochemical control of the quinone pool redox state allows one to flexibly trigger or to stop the enzymatic reaction at any moment, while the pH-sensitive dye is used to monitor pH changes inside the lumen of the liposomes as a result of proton translocation by the enzymes. By using the fluorescence intensity of a second, lipid-bound fluorescent dye, the size and volume of individual liposomes can be determined and thus the quantification of enzyme proton pumping activity. Using this technique, we notably found that cytochrome *bo*₃ molecules are able to enter into a

spontaneous leak state that rapidly dissipates the proton motive force. The goal of this article is to introduce the technique of single liposome measurements in detail.

PROTOCOL:

1. Preparation of an *E. coli* Polar Lipids Mix

NOTE: *E. coli* lipids used for the liposomes' preparation should be aliquoted and thoroughly mixed with ubiquinone-10 (enzyme substrate) and long-wavelength fluorescent dye-labeled lipids (for liposomes size determination) prior to the reconstitution.

1.1. Using a glass syringe, transfer 200 μ L of chloroform stock of lipid polar extract from *Escherichia coli* (25 mg/mL) into glass vials to make 5 mg aliquots.

1.2. Add 50 μ L of 1 mg/mL ubiquinone-10 (UQ-10; in chloroform) to the lipids to make the final ratio of UQ-10: lipids 1:100 (1% w/w).

1.3. Add 20 μ L of 1 mg/mL (0.4% w/w) of a long-wavelength fluorescent dye-labeled lipid (FDLL) to the lipids/UQ-10 mix.

1.4. Homogenize the chloroform solution by short vortexing and evaporate most of the chloroform under a gentle nitrogen or argon flow. Remove the chloroform traces entirely by further evaporation under vacuum for at least 1 h.

NOTE: Lipid aliquots can be stored under an inert atmosphere at -20 °C for several months.

2. Reconstitution of Cytochrome *bo*₃

NOTE: For the purification of cytochrome *bo*₃ from *E. coli*, follow the protocol from Rumbley *et al.*¹⁵ To ensure high purity of natively-folded enzyme samples, add size-exclusion chromatography after the affinity purification step described by Rumbley *et al.*¹⁵

2.1. Add 312.5 μ L of 40 mM MOPS-KOH/60 mM K₂SO₄, pH 7.4, to one aliquot of lipids/UQ-10/FDLL dry mix (5 mg, step 1.4) and re-suspend using vortex mixing followed by 2 min of treatment in an ultrasonic bath.

2.2. Add 125 μ L of 25 mM 8-hydroxypyrene-1,3,6-trisulphonic acid (HPTS), a pH-sensitive fluorescent dye that needs to be encapsulated inside the liposomes.

2.3. Add 137.5 μ L of 250 mM n-octyl β -D-glucopyranoside (OGP) surfactant, mix using vortex and sonicate in an ultrasonic water bath for 10 min to ensure all lipids are solubilized into surfactant micelles. Transfer the dispersion into a 1.5 mL plastic tube.

NOTE: Cloudy suspension of lipids should become transparent after the solubilization with

surfactant.

2.4. Add the required amount of cytochrome *bo*₃ (see the note below) and add ultrapure water to make a total volume of 50 µL (cytochrome *bo*₃ solution plus water). Incubate at 4 °C for 10 min on a roller mixer.

NOTE: Typical amount for single enzyme conditions is 0.1 – 0.2% (w/w protein-to-lipid, *i.e.* 5 – 10 µg of protein), although it can be increased till 1 – 2% (50 – 100 µg of protein) if the goal is to observe only the electrochemical activity (see below). As a negative control, liposomes without cytochrome *bo*₃ can be prepared.

2.5. Weigh 2x 50 mg and 2x 100 mg of polystyrene microbeads into four 1.5 mL-tube caps and close with paraffin film to prevent drying.

NOTE: Before use, polystyrene microbeads should be washed with methanol, water and stored in water according to the manufacturer's manual.

2.6. Add the 1st 50 mg of polystyrene microbeads into the reconstitution mixture (step 2.4) by putting the cap with polystyrene microbeads on the 1.5 mL-plastic tube with dispersion and performing a short spin for a couple of seconds. Incubate at 4 °C on a roller mixer for polystyrene microbeads to adsorb the surfactant for 30 min.

2.6.1. Repeat the additions of polystyrene microbeads and incubations as follows: add 50 mg of microbeads for 60 min of incubation; add 100 mg of microbeads for 60 min of incubation; and add 100 mg of microbeads for 120 min of incubation.

NOTE: The solution above the settled microbeads will turn translucent during step 2.6 as proteoliposomes are formed.

2.7. Separate the proteoliposome solution from the polystyrene microbeads using a micropipette with a thin tip. Dilute the dispersion in 90 mL of 20 mM MOPS-KOH/30 mM K₂SO₄, pH 7.4 (MOPS buffer) and transfer in a Ti45 ultracentrifuge tube.

2.8. Ultracentrifuge the dispersion using Type 45 Ti rotor at 125,000 x g (at *r*_{max}) for 1 h to pellet the proteoliposomes.

NOTE: Smaller centrifuge tubes can be used although the dilution in large buffer volume helps to reduce the concentration of the non-encapsulated HPTS in the final suspension.

2.9. Discard the supernatant, rinse the pellets with 20 mM MOPS-KOH/30 mM K₂SO₄, pH 7.4 buffer (without resuspending the pellet). Then re-suspend the proteoliposomes in 500 µL of MOPS buffer by pipetting it back and forth with thin tip micropipette. Then transfer to a 1.5 mL plastic tube.

2.10. Centrifuge the suspension for 5 min at 12,000 x g to remove the debris. Transfer the supernatant (reconstituted proteoliposomes) into a new vial.

2.11. Store the reconstituted proteoliposomes dispersion at 4 °C overnight and use within 2 days.

3. Fabrication of Semi-Transparent Gold Electrodes

NOTE: The smooth gold surface is obtained by a template stripping method of a 30 nm-thick layer of 99.99% gold from an atomically smooth silicon wafer. The small thickness of the gold layer is important since it must be semi-transparent to permit fluorescence observation. Details of gold evaporation (physical vapor deposition, PVD) can be find elsewhere¹⁶ and only template-stripping is covered here. Alternatively, ultra-smooth gold chips can be purchased elsewhere (see **Table of Materials**).

3.1. Glue up to 9 glass cover slips (0.17 mm thick) onto the evaporated gold surface using bi-component low-fluorescence epoxy. Cure the glue at 80 °C for 4 h.

3.2. Just before modification with a self-assembled monolayer (step 4), detach the glass cover slips from the silicon wafers with a blade. Due to the thinness of the cover slips, take care when detaching the cover slips to not crack or break the glass slides.

4. Modification of the Gold Surface with Self-Assembled Monolayer (SAM)

4.1. Prepare 5 mL of water solution of 1 mM 6-mercaptohexanol (6MH).

4.2. Dip freshly detached gold-coated cover slips (step 3.1) into 6MH solution and leave at 20 – 25 °C overnight (>16 h) to form the SAM.

NOTE: Thiol solutions have an unpleasant smell, so a closed vessel should be used when incubating the gold-coated cover slip.

4.3. The next day, remove the gold-coated cover slip from the 6MH solution, wash briefly with water or methanol and then with isopropanol. Dry under a gentle gas flow.

5. Electrochemical Testing of Proteoliposomes Activity

NOTE: The electrochemical activity of enzyme is first verified on a closely packed liposomes layer (step 5) and lower vesicle coverages are used in single vesicles experiment to measure pH changes in the lumen of the liposomes (step 6).

5.1. Assemble the gold-coated cover slip in a spectro-electrochemical cell (see **Figure 1**). Make contact to the gold with a flat wire outside an area defined by a rubber O-ring.

5.2. Add 2 mL of the electrolyte buffer solution and place the reference and auxiliary electrodes in the cell.

NOTE: As further discussed in the Discussion section, references electrodes without chloride are preferred to prevent formation of Au(I)Cl during electrochemistry. Here, a Hg/Hg₂SO₄ (sat. K₂SO₄) reference electrode is used and potentials are given versus Standard Hydrogen Electrode (SHE) using 0.658 V vs SHE for the Hg/Hg₂SO₄ (sat. K₂SO₄).

5.3. Run electrochemical impedance spectroscopy (0.1 Hz – 100 kHz) at open cell potential (OCP) potential to assess the quality of SAM. Convert impedance values to admittance and divide by $2\pi\omega$ to plot a Cole-Cole plot, where ω is the frequency (see the following references^{16–18} for details on converting and interpreting impedance spectra).

NOTE: Compact and dense SAMs of 6MH should give a close to semi-circle Cole-Cole plot and give capacitance values in the range 2.5–3.0 $\mu\text{F}/\text{cm}^2$ (**Figure 2A**). If significant deviation from semi-circle shape or out-of-range capacitance values are obtained, change the electrode.

5.4. Run blank cyclic voltammograms (CVs) with scan rates 100 and 10 mV/s in the potential region -0.3 – 0.8 V. A typical CV is shown on (**Figure 2B**, dashed line).

NOTE: This should demonstrate almost pure capacitive behavior and an absence of significant faradaic current, even under ambient oxygen conditions as used here.

5.5. Add proteoliposomes (0.5 mg/mL final lipids concentration, 1–2% (w/w) ratio of cytochrome *bo*₃ to lipid) to the electrochemical cell and mix slightly with a pipette. Wait until the adsorption of proteoliposomes on the electrode surface is finished (30–60 min at room temperature).

NOTE: Cyclic voltammograms (CVs) at 10 mV/s can be run during the proteoliposomes adsorption to follow the process. The adsorption is finished when consecutive CVs stop changing. Information about CV techniques can be found in the following textbooks.^{19, 20}

5.6. Wash the cell by changing the buffer solution at least 10 times but avoid leaving the electrode surface completely dry.

5.7. Run the electrochemical impedance spectroscopy at OCP (**Figure 2A**, blue line) to confirm the SAM on the gold electrode remains unchanged and CVs with scan rates 10 and 100 mV/s to observe catalytic ubiquinol oxidation (and oxygen reduction) by cytochrome *bo*₃ at onset potentials of electrochemical quinone reduction about 0 V vs SHE) (**Figure 2B**).

6. Detection of Enzymatic Proton Pumping by Fluorescence Microscopy

6.1. Modify the gold electrode as in step 5 but using 100x less proteoliposomes compared to step 5.5 (*i.e.*, 5 $\mu\text{g}/\text{mL}$). For single enzyme studies, reduce the cytochrome *bo*₃ to lipid ratio to

0.1-0.2% (w/w).

NOTE: Liposomes will sparsely adsorb on the electrode surface enabling single vesicle monitoring by fluorescence microscopy. Under these single-enzyme conditions, the amount of enzyme immobilized on the electrode surface is insufficient for the observation of a catalytic current.

6.2. Place the electrochemical cell on the oil objective (60X) of an inverted fluorescence microscope with a drop of immersion oil. Using appropriate filters for FDLL fluorescence, focus on the electrode surface. Single liposomes should appear as bright spots at the diffraction limit of the microscope/objective. Take an image of FDLL fluorescence.

6.3. Switch to the one of HPTS fluorescence filter sets on the microscope to verify that HPTS fluorescence is clearly visible and distinguishable from the background (at the chosen exposure time, see step 6.4). Increase the light intensity of the lamp source if it is not the case.

6.4. Program the microscope software to perform a timed image acquisition by alternating two HPTS filter sets (menu **Applications | Define/Run ND Acquisition**). Set the delay between image acquisitions at minimum. In this experiment, use a 1 s exposure and 0.3 s delay (due to the turret movement).

NOTE: The ratio between fluorescence intensities at these two channels will be later converted to pH inside the liposomes at each time point. The duration of the acquisition can vary according to the expected pH change rate; 5 min is used in this article.

6.5. Adjust the settings of the potentiostat to change the potential during the image acquisition. For example, in this experiment, use the following sequence: 0 – 60 s: no potential applied (*i.e.* OCP); 60 – 180 s: -0.2 V (vs SHE); 180 – 300 s: 0.4 V (vs SHE). In this case, only the applied potential in the second phase (-0.2 V vs SHE) is sufficient to efficiently reduce the quinone pool.

6.6. Run simultaneously the timed images acquisition (microscope) and the potential sequence (potentiostat) by manually starting both measurements at the same time.

NOTE: The experiment can be repeated on the same electrode several time by moving the microscope stage to a different area on the surface. The delay between the acquisitions should be at least 5-10 min to insure complete pH equilibration of liposomes on the surface. Different durations and potential patterns can be applied depending on the need, although imaging time is limited due to photobleaching of HPTS during acquisition.

7. Analysis of Fluorescence Images

NOTE: A typical experiment produces a set of images with a time step (*e.g.*, 2.6 s) for each of the two channels, *i.e.*, (duration * 2 / 2.6) images. An example of such image set recorded during a single enzyme experiment can be accessed via the Research Data Leeds repository²¹. An image

treatment consists of several steps by using Fiji (ImageJ) and high-level mathematical analysis programming language software (henceforth referred to as scripting software, see **Table of Materials** for details).

7.1. Use Fiji to separate time lapse file, align images and save as separate channels and timeframes.

7.1.1. Open time lapse file using the Bioformats importer provided within Fiji (macro command: run("Bio-Formats Importer")).

7.1.2. Use the plugin StackReg to align each frame to the first one to account for possible stage movements or thermal drift occurred during the acquisition (macro command: run("StackReg ", "transformation=Translation")).

7.1.3. Save separate uncompressed image files for each channel and each time step in TIFF format into a single folder.

7.1.4. Extract the exact time stamps for each frame from the time lapse file metadata and save them as a CSV-file in the same folder as the images. Alternatively, extract time stamps using the acquisition software and save manually into a CSV-file using a spreadsheet software.

NOTE: The folder with images and time stamps is now ready to be analyzed by the scripting software. Steps 7.1.1. – 7.1.4. can be automatized using a Fiji script (a script written in Python for batch processing of time lapse files is provided, use “**Ctrl-Shift-N**” (in Windows), then **File | Open** to load the script).

7.2. Use scripting software for automatized processing of the images. Load the provided code to the software and click **Run to end**. When prompted, select the folder containing the images from step 7.1.

NOTE: A script for analysis is provided as mlx-format live script, with extensive comments, to identify single liposomes, fit them to 2D-Gaussian function, filter them and quantify the pH values at each point of time. The following sub-steps are executed automatically by the script.

7.2.1. Load all time frames images for a single experiment into the memory.

7.2.2. Average all images for a single channel (select the channel with the highest fluorescence intensity) and use the averaged image to identify all maximums that might correspond to the single liposomes.

7.2.3. Fit the identified maximums on the averaged image to the 2D-Gaussian function and save the resulted fit parameters for each liposome (see **Figure 4A** for fitted liposome example).

7.2.4. Filter the maximums according to the expected single liposomes criteria, such as size,

circularity and intensity. Reject liposomes that are poorly fitted due to low signal to noise, closely neighboring liposome or being too close to the image edge.

7.2.5. Load time stamps from the external file and fit each filtered liposome to the 2D-gaussian function on each time frame image separately.

NOTE: The use of parallel programming and multicore CPU can significantly enhance the calculation speed at this stage.

7.2.6. Filter the liposomes again according to the similar criteria as in step 7.2.4 but applied to each time frame.

7.2.7. On each time step, define the fluorescence intensity ratio of a liposome as the ratio of volumes enclosed by fitted 2D-Gaussian function at two channels. Calculate the pH values from the ratio of intensities determined from the two HPTS channels and using a calibration curve (step 8). Plot resulting pH-time curves for the liposomes and observe their pH change when the potential is applied.

8. Performing a Calibration Curve of HPTS Fluorescence

NOTE: To convert HPTS fluorescence ratio to an intravesicular pH, a calibration curve must be first established that would take into account particular conditions of experiments such as gold transmittance, filters quality, *etc.* This step has to be performed only once or twice and the calibration data can be used as long as the setup and measurement parameters remain the same in step 6.

8.1. Prepare an electrode with sparsely adsorbed liposomes (without cytochrome *bo*₃ as described in Steps 5 and 6.1).

8.2. Add 2 μL of 0.1 mg/mL gramicidin solution in ethanol to 2 mL of buffer to create the concentration 100 ng/mL.

8.3. Capture two fluorescence images for the two HPTS channels.

8.4. Change the pH of the cell by addition of small aliquots of 1 M HCl or 1 M H₂SO₄. Measure the pH of the buffer with a standard pH meter and capture two fluorescence images for the two HPTS channels.

8.5. Repeat steps 8.4 for the pH range 6 to 9.

8.6. Using the algorithm from 7.2, fit, filter and calculate the average HPTS fluorescence ratio of individual liposomes at every pH. Take the HPTS ratio average over all liposomes.

8.7. Fit the resulting pH-ratio dependence to the following equation:

$$pH = pK_a + \log \left(\frac{Ratio - R_a}{R_b - Ratio} \right), \text{ where } pK_a, R_a \text{ and } R_b \text{ are fitting parameters.}$$

8.8. Use pK_a , R_a and R_b to convert the HPTS ratio of individual liposomes in step 7.2.7. to pH values.

REPRESENTATIVE RESULTS:

The quality of the gold-modified cover slip (the electrode) with a SAM of 6MH is checked before each experiment with electrochemical impedance spectroscopy. **Figure 2A** shows representative Cole-Cole plots measured using electrochemical impedance spectroscopy before and after liposomes are adsorbed. If the quality of SAM is sufficient, impedance spectroscopy should demonstrate an almost pure capacitive behavior resulting in a semi-circle Cole-Cole plot. The diameter of the semi-circle in a Cole-Cole plot equals the double layer capacitance of the electrode surface, which should be in the 2.5-3.0 $\mu\text{F}/\text{cm}^2$ range. Note that the capacitance should not change significantly upon liposomes adsorption, although a slight increase in capacitance might be observed (high liposomes coverage, **Figure 2A**, bottom) or a slight shift in impedance can be observed at low frequencies (low liposome coverages, **Figure 2A**, top).

Electrochemistry can be used to test the catalytic activity of cytochrome bo_3 in the proteoliposomes, where catalytic currents measured with cyclic voltammetry reflect the ubiquinol-oxygen oxidoreductase activity. However, significant catalytic currents can only be detected at high quantities of adsorbed proteoliposomes and a closely packed layer of proteoliposomes on the gold-modified cover slip is needed. **Figure 2B** demonstrates representative cyclic voltammograms (CVs) of the electrode before and after liposomes adsorption with either low or high liposome coverage. No faradaic current is observed on the blank CV because oxygen reduction by the bare gold electrode is blocked by the SAM. When the surface is saturated with liposomes with high cytochrome bo_3 content (1.3% (w/w) in this case), a clear catalytic wave due to oxygen reduction is observed with an onset potential of 0 V, *i.e.*, the potential of ubiquinone reduction (**Figure 2B**, blue line). The ubiquinone pool acts both as the natural substrate for the enzyme and as an electron mediator, transferring electrons from the enzyme to the electrode surface. We note that under high liposome coverage, no distinction of individual vesicles is possible by microscopy and lower coverage is needed for single liposome studies. At low liposome coverage (but high protein to lipid ratio), the catalytic current is significantly reduced, barely distinguishable from background (**Figure 2B**, red line). Note that under single enzyme conditions (low protein to lipid ratio), the catalytic current is even lower and cannot be measured reliably.

Figure 3 shows fluorescence images of liposomes adsorbed on the electrodes at three different coverages. All images were taken in identical light and exposure conditions and their brightness was adjusted equally to enable direct comparison. The dye-containing-liposomes are visible on the images as bright spots. The central part of the image was photobleached for a couple of minutes to reveal background fluorescence level (we note that FDLL is relatively photostable and is not photobleached completely in **Figure 3**). The images on two HPTS channels are superimposable, where the ratio between the two channels (410/535 and 470/535) corresponds

to pH 7.4 used in this experiment. A larger number of liposomes are visible with the FDLL channel, which indicates the presence of liposomes that have no HPTS encapsulated. The difference between the HPTS and FDLL channels is more pronounced at higher coverages, possibly because at high liposome coverage, liposomes are more likely to burst or fuse on the surface.

The image analysis requires the alignment of frames (against the first frame) using the ImageJ, plugin StackReg. The alignment is indispensable in most situations since a slight thermal drift of the sample usually occurs within the timescale of the experiment, changing the liposome coordinates. All further analysis is performed using a scripting software code. This code performs automatic liposome identification. As the size of liposomes are below the Abbe diffraction limit, their fluorescence is seen as a point spread function much larger than the actual liposome diameter (**Figure 4A**). The code fits the fluorescence intensity of each vesicle on two channels to a 2D-Gaussian function (**Figure 4A**) and calculates their volumetric intensity ratio that can be converted to pH using a calibration curve. By performing these actions on all time frames, the pH is obtained inside every vesicle within timescale of the experiment (**Movie 1**). **Figure 4B** shows the medians of all vesicle pH changes within a single image when cytochrome *bo₃* content is 1.3%. An increase in pH (proton pumping) is clearly visible when a potential between -0.1 and -0.3 V is applied, but not for 0 V since the latter is not sufficient to reduce the quinone pool. When cytochrome *bo₃* content is much lower (protein-to-lipid ratio of 0.1%, **Figure 4C**), the median curves become almost indistinguishable from those of empty liposomes (**Figure 4D**). The difference becomes evident when considering individual pH traces of random liposomes (**Figure 5, 6 and 7**). While liposomes without cytochrome *bo₃* display no significant pH changes with respect to noise (**Figure 7**), a selection of liposomes show an increase (dominantly) or a decrease of pH when a potential is applied that activates cytochrome *bo₃* (**Figures 6**, grey zone). The obvious difference in pH-traces between liposomes is consistent with the low protein-to-lipid ratio, where cytochrome *bo₃* is only present in a small subset of liposomes activity. The prevalence of a pH increase over decrease is explained by the reconstitution method that favors an “outward” orientation of enzyme molecules (*ca.* 75%). The fraction of liposomes that display a pH change increases when the cytochrome *bo₃* to lipid ratio is higher (**Figure 5**). Note also that some liposomes stop proton pumping and enter into proton dissipation mode before the end of the potential application. We attribute this behavior to the cytochrome *bo₃* molecules entering a “leak state”, which allows protons to flow back into the liposome lumen¹⁴.

The further analysis of the pH traces including their fitting and determination of proton pumping/leaking rates can be done using a script we have published previously^{14,22} and can be obtained from the Research Data Leeds Repository^{23,24}.

FIGURE AND TABLE LEGENDS:

Figure 1: Principle of the method. (A) Principle of single enzyme activity monitoring and (B) the scheme of the experimental setup used in this work with a cutaway to demonstrate an internal part of the cell.

Figure 2: Electrochemical response of liposomes. (A) Cole-Cole plots measured at OCP (0.22 V vs SHE; black squared line) before and after liposomes (1.3% w/w cytochrome *bo₃*) adsorption from

solutions at concentrations of 5 $\mu\text{g/mL}$ (red circle line) and 500 $\mu\text{g/mL}$ (blue circle line). (B) Cyclic voltammograms at 100 mV/s (left panel) and 10 mV/s (right panel) of Au-SAM (black dashed line) before and after liposomes (1.3% w/w cytochrome bo_3) adsorption from solutions at concentrations of 5 $\mu\text{g/mL}$ (red line) and 500 $\mu\text{g/mL}$ (blue line).

Figure 3: Fluorescence microscopy of liposomes. Fluorescence images of liposomes at different surface coverages: surfaces incubated for 30 min with a liposomes solution of 5 $\mu\text{g/mL}$ (top row), 20 $\mu\text{g/mL}$ (middle row) and 500 $\mu\text{g/mL}$ (bottom row). Left column corresponds to 470/535 nm excitation/emission filter setup (1st HPTS channel); middle column – 410/535 nm (2nd HPTS channel); right column – 560/645 nm (FDLL channel). The images were taken at magnification 90X (60X objective and 1.5X magnification by the microscope prior to the camera), 1 s exposure and similar light intensity. Scale bars correspond to 50 μm .

Figure 4: Liposome identification and pH change. (A) 3D-view of a part of fluorescence image containing a typical single liposome. Its fluorescence intensity is seen as a point spread function (color surface) that is fitted to a 2D-Gaussian function (black mesh). (B-D) pH, displayed as the median of all vesicles within single image area (typically several hundred) at different applied potentials. From 0-60 s, no potential is applied (OCP). Between 60 and 180 s (grey zone) different potentials were applied: 0 V (black), -0.1 V (red), -0.2 V (green), -0.3 V (blue). Between 180 and 300 s, 0.4 V vs. SHE was applied. The cytochrome bo_3 to lipid ratios were (B) 1.3%, (C) 0.1%, (D) 0%. The traces are offset for clarity.

Figure 5: pH traces of liposomes containing 1.3% of enzyme. Traces of pH change of 72 single liposomes, randomly selected, containing cytochrome bo_3 (1.3% w/w) measured and analyzed during one experiment. From 0-60 s, no potential is applied (OCP); between 60 and 180 s (grey zone) -0.2 V vs. SHE; between 180 and 300 s, 0.4 V vs. SHE was applied.

Figure 6: pH traces of liposomes containing 0.1% of enzyme. Traces of pH change of 72 single liposomes, randomly selected, containing cytochrome bo_3 (0.1% w/w) measured and analyzed during one experiment. From 0-60 s, no potential is applied (OCP); between 60 and 180 s (grey zone) -0.2 V vs. SHE; between 180 and 300 s, 0.4 V vs. SHE was applied.

Figure 7: pH traces of liposomes without enzyme. Traces of pH change of 72 single liposomes, randomly selected, without cytochrome bo_3 measured and analyzed during one experiment. From 0-60 s: no potential is applied (OCP); between 60 and 180 s (grey zone) -0.2 V vs. SHE; between 180 and 300 s, 0.4 V vs. SHE was applied.

Movie 1: Animation of single liposome pH change. (top panel) Change of a liposome fluorescence on two HPTS channels (shown as 3D-surface plot of the corresponding area) during 300 s of the experiment. The reductive potential was applied between 60 s and 180 s. (bottom panel) Corresponding plot of volumetric intensity ratio of two HPTS channels versus time. The zone of reductive potential application is shaded in grey.

DISCUSSION:

The method described is suitable to study proton pumping by respiratory membrane proteins that can be reconstituted into liposomes and are able to exchange electrons with the quinone pool. Proton pumping activity can be monitored at the single-enzyme level using pH-sensitive (ratiometric) dyes encapsulated in the liposome lumen (**Figure 1A**).

The method relies on the ability of ubiquinone (or other quinones), incorporated into the lipid bilayer, to exchange electrons with electrodes modified with a SAM¹⁸. The properties of the gold electrode, modified with a SAM, are very specific. Liposomes need to adsorb onto the SAM and the SAM needs to be thin enough to enable rapid electrochemical oxidation and reduction of the quinone pool in the liposomes. Importantly, the interaction between the electrode and liposome needs to be weak enough not to impair the integrity of the lipid membrane. Furthermore, depending on the details of the experimental system, it is desirable if the SAM prevents gold-catalyzed side reactions, such as oxygen reduction. Finally, the SAM needs to be stable within the potential window used. The use of ultra-flat gold electrodes modified with high-quality SAM of 6MH fulfills these requirements in the case of cytochrome *bo*₃. Stripping gold electrodes from atomically-flat silicon wafers creates an ultra-smooth gold surface with a near ideal SAM as indicated by the impedance spectra. Please note that we form the SAM from 6MH in water. We previously reported on SAMs with 6MH in isopropanol¹⁶, but these SAMs exhibit higher double-layer capacitance values and impedance spectra that indicate a more heterogeneous surface (*i.e.*, more defects) with a lower resistance. Furthermore, when SAMs from 6MH are prepared from a water solution, less background oxygen reduction (due to gold catalyzed oxygen reduction) is observed compared to SAMs formed from an ethanol/isopropanol solution. All these observations indicate that SAMs from 6MH in water solutions have a higher quality with less defects. Higher-quality SAMs can also be formed from thiol-compounds with longer alkane chains (*e.g.*, 1-hydroxy-decane-thiol), but the electron transfer kinetics become slower as the SAM thickness increases.

During (spectro)electrochemistry, chloride ions should be avoided in the buffer solution since they may chemisorb on the gold surface at high potentials, possibly deteriorating the SAM quality. For the same reason, a chloride-free reference electrode should be preferred.

Another important point is the background permeability of the liposomes to protons, which should be low enough to allow proton accumulation as a result of enzymatic proton translocation on the timescale of experiment. The exact lipid composition may influence this permeability. In our work, we use polar *E. coli* lipid extracts supplemented with ubiquinone-10. Proton permeability has been shown to depend strongly on the fatty-acid chain length²⁵, although this has been contradicted in studies with more complex lipid compositions²⁶. Interestingly, ubiquinone has recently been shown to enhance membrane stability²⁷. Finally, residual detergents from the protein reconstitution procedure can influence proton leakage and hence it is important to remove detergents as completely as possible using hydrophobic polystyrene microbeads, dilution followed by centrifugation and thorough rinsing of the electrochemical cell after the proteoliposomes are adsorbed on surface. If doubtful, liposome permeability can be verified by following the lumen pH change (via HPTS) in a response to an external pH jump²⁸.

Single enzyme measurement require unilamellar liposomes and smaller liposomes are expected to show faster or higher pH changes as the lumen's volume decreases. To achieve this, polystyrene microbeads are added gradually in small amounts leading to slow rate of liposomes formation. In such conditions, small liposomes of homogeneous size are formed with an average diameter of 70 nm and a polydispersity index of 0.24 in our case¹⁴. HPTS also adsorbs on the polystyrene microbeads, reducing the capacity of polystyrene microbeads to adsorb detergent. Hence, excess polystyrene microbeads should be added to compensate for its loss. Treatment with polystyrene microbeads was observed to reduce the HPTS concentration in the lipid-protein suspension by about a quarter (from 5 mM to 3 - 4 mM). However, the actual HPTS concentration in the liposome's lumen might be higher since the liposomes are formed early on the treatment with polystyrene microbeads. Instead of HPTS, other pH-sensitive dyes can be used. However, it is important that the dye is membrane-impermeable and ratiometric. The latter avoids experimental errors due to photobleaching and varying amounts of encapsulated dye.

Simple calculations can be made to estimate whether, stochastically, most of the non-empty liposomes contain only one enzyme molecule. Assuming an average liposome diameter of 70 nm, a lipid molecular weight of 750 Da and a lipid area of 0.65 nm² gives us a liposome mass of 5.2×10^{-17} g, *i.e.*, 9.6×10^{13} liposomes per preparation (5 mg of lipid). At 0.1% of cytochrome *bo*₃ (MW = 144 kDa), 2×10^{13} enzyme molecules are present, corresponding to a 0.2 protein:liposome ratio. Using a Poisson distribution, one may further calculate that the probability to find one enzyme molecule (17%) per liposome is ten times higher than probability to find more than one molecule (1.7%). More precise calculations can be made if the losses of enzyme and lipids during the reconstitution determined from corresponding assays are taken into account. The latter can be estimated from a protein assay and by measuring FDLL fluorescence of prepared liposomes.

Once the protocol is established and single enzyme traces are recorded, further modifications of the method is feasible depending on the aim. One might think about an enzyme inhibitor addition or introduction of an ionophore into the system. Care should be taken to verify that addition of solvents used to prepare ionophore or inhibitor stock do not influence the state of SAM or permeability of the liposomes.

The technique as described here is limited to a particular group of enzymes that are both membrane proton transporters and quinone-converting. The equipment and experimental conditions as used here were optimized for the enzymatic activity of cytochrome *bo*₃. Here, the time resolution is 2-3 seconds, determined by the exposure time and the time it takes for the turret to change filters. The duration of the experiment is limited by photobleaching of fluorescent dye and, thus, by the light intensity. We have previously found average proton translocation rate of cytochrome *bo*₃ to be 73 ± 2.2 protons/s using this technique, although activities down to 20 protons/s were detected. To use these techniques for enzymes with either more or less activity, image acquisition parameters need to be adapted (*i.e.*, light intensity, exposure time and duration of experiment). In the future, this method can be extended towards other electron transport driving proton pumping enzymes, an example being mitochondrial complex I. Other enzymes might require different reconstitution protocols, and this would need to be optimized for each different transporter. Proton pumps that are not quinone-converting

enzymes can also be studied, *e.g.*, ATP-driven²⁹, although in this case proton translocation cannot be triggered electrochemically, but addition of an initiator (*e.g.*, ATP) is required. In the latter case, there is no need to use a gold-modified cover slip. Moreover, provided that a suitable membrane-impermeable and ion-sensitive fluorescent dye is used, this method can be extended to other ionic pumps, *e.g.*, sodium and potassium.

ACKNOWLEDGMENTS:

The authors acknowledge the BBSRC (BB/P005454/1) for financial support. NH was funded by the VILLUM Foundation Young Investigator Program.

DISCLOSURES:

The authors have nothing to disclose.

REFERENCES:


1. Claessen, V.I. *et al.* Single-Biomolecule Kinetics: The Art of Studying a Single Enzyme. *Annual Review of Analytical Chemistry*. **3** (1), 319–340, 10.1146/annurev.anchem.111808.073638 (2010).
2. Rojek, M.J., Walt, D.R. Observing Single Enzyme Molecules Interconvert between Activity States upon Heating. *PLoS ONE*. **9** (1), e86224, 10.1371/journal.pone.0086224 (2014).
3. Velonia, K. *et al.* Single-Enzyme Kinetics of CALB-Catalyzed Hydrolysis. *Angewandte Chemie International Edition*. **44** (4), 560–564, 10.1002/anie.200460625 (2005).
4. Lu, H.P., Xun, L., Xie, X.S. Single-molecule enzymatic dynamics. *Science (New York, N.Y.)*. **282** (5395), 1877–82, 10.1126/science.282.5395.1877 (1998).
5. Engelkamp, H., Hatzakis, N.S., Hofkens, J., De Schryver, F.C., Nolte, R.J.M., Rowan, A.E. Do enzymes sleep and work? *Chemical Communications*. (9), 935–940, 10.1039/b516013h (2006).
6. Boukobza, E., Sonnenfeld, A., Haran, G. Immobilization in Surface-Tethered Lipid Vesicles as a New Tool for Single Biomolecule Spectroscopy. *The Journal of Physical Chemistry B*. **105** (48), 12165–12170, 10.1021/jp012016x (2001).
7. Hsin, T.-M., Yeung, E.S. Single-Molecule Reactions in Liposomes. *Angewandte Chemie International Edition*. **46** (42), 8032–8035, 10.1002/anie.200702348 (2007).
8. García-Sáez, A.J., Schwille, P. Single molecule techniques for the study of membrane proteins. *Applied Microbiology and Biotechnology*. **76** (2), 257–266, 10.1007/s00253-007-1007-8 (2007).
9. Onoue, Y. *et al.* A giant liposome for single-molecule observation of conformational changes in membrane proteins. *Biochimica et Biophysica Acta (BBA) - Biomembranes*. **1788** (6), 1332–1340, 10.1016/J.BBAMEM.2009.01.015 (2009).
10. Jefferson, R.E., Min, D., Corin, K., Wang, J.Y., Bowie, J.U. Applications of Single-Molecule Methods to Membrane Protein Folding Studies. *Journal of Molecular Biology*. **430** (4), 424–437, 10.1016/J.JMB.2017.05.021 (2018).
11. Rigaud, J.L., Pitard, B., Levy, D. Reconstitution of membrane proteins into liposomes: application to energy-transducing membrane proteins. *BBA - Bioenergetics*. **1231** (3), 223–246, 10.1016/0005-2728(95)00091-V (1995).
12. Jesorka, A., Orwar, O. Liposomes: Technologies and Analytical Applications. *Annual Review of Analytical Chemistry*. **1** (1), 801–832, 10.1146/annurev.anchem.1.031207.112747

- 660 (2008).
- 661 13. Seddon, A.M., Curnow, P., Booth, P.J. Membrane proteins, lipids and detergents: Not just
662 a soap opera. *Biochimica et Biophysica Acta - Biomembranes*. **1666** (1–2), 105–117,
663 10.1016/j.bbamem.2004.04.011 (2004).
- 664 14. Li, M. *et al.* Single Enzyme Experiments Reveal a Long-Lifetime Proton Leak State in a
665 Heme-Copper Oxidase. *Journal of the American Chemical Society*. **137** (51), 16055–16063,
666 10.1021/jacs.5b08798 (2015).
- 667 15. Rumbley, J.N., Furlong Nickels, E., Gennis, R.B. One-step purification of histidine-tagged
668 cytochrome bo₃ from *Escherichia coli* and demonstration that associated quinone is not required
669 for the structural integrity of the oxidase. *Biochimica et Biophysica Acta (BBA) - Protein Structure
670 and Molecular Enzymology*. **1340** (1), 131–142, 10.1016/S0167-4838(97)00036-8 (1997).
- 671 16. Jeuken, L.J.C. *et al.* Phase separation in mixed self-assembled monolayers and its effect
672 on biomimetic membranes. *Sensors and Actuators, B: Chemical*. **124** (2), 501–509,
673 10.1016/j.snb.2007.01.014 (2007).
- 674 17. Barsoukov, E., Macdonald, J.R. *Impedance Spectroscopy Theory, Experiment, and
675 Applications*. Chapters 1-2. John Wiley & Sons. (2005).
- 676 18. Jeuken, L.J.C., Connell, S.D., Henderson, P.J.F., Gennis, R.B., Evans, S.D., Bushby, R.J. Redox
677 Enzymes in Tethered Membranes. *Journal of the American Chemical Society*. **128** (5), 1711–1716,
678 10.1021/ja056972u (2006).
- 679 19. Compton, R.G., Banks, C.E. *Understanding voltammetry*. Chapter 4. World Scientific
680 Publishing Europe Ltd. Singapore. (2018).
- 681 20. Girault, H.H. *Analytical and Physical Electrochemistry*. Chapter 10. EPFL Press. Lausanne.
682 (2004).
- 683 21. Mazurenko, I., Jeuken, L.J.C. Timelapse (movie) of single vesicles fluorescence change
684 upon potential application. 10.5518/442 (2018).
- 685 22. Li, M., Khan, S., Rong, H., Tuma, R., Hatzakis, N.S., Jeuken, L.J.C. Effects of membrane
686 curvature and pH on proton pumping activity of single cytochrome bo₃enzymes. *Biochimica et
687 Biophysica Acta - Bioenergetics*. **1858** (9), 763–770, 10.1016/j.bbabbio.2017.06.003 (2017).
- 688 23. Li, M., Tuma, R., Jeuken, L.J.C. MATLAB code for the analysis of proton transport activity
689 in single liposomes. 10.5518/150 (2017).
- 690 24. Li, M., Tuma, R., Jeuken, L.J.C. Time-resolved fluorescence microscopic data (traces) of
691 individual lipid vesicles with proton transport activity. 10.5518/151 (2017).
- 692 25. Paula, S., Volkov, A.G., Van Hoek, A.N., Haines, T.H., Deamer, D.W. Permeation of protons,
693 potassium ions, and small polar molecules through phospholipid bilayers as a function of
694 membrane thickness. *Biophysical Journal*. **70** (1), 339–348, 10.1016/S0006-3495(96)79575-9
695 (1996).
- 696 26. Brookes, P.S., Hulbert, A.J., Brand, M.D. The proton permeability of liposomes made from
697 mitochondria) inner membrane phospholipids: No effect of fatty acid composition. *Biochimica et
698 Biophysica Acta (BBA) - Biomembranes*. **1330** (2), 157–164, 10.1016/S0005-2736(97)00160-0
699 (1997).
- 700 27. Eriksson, E.K., Agmo Hernández, V., Edwards, K. Effect of ubiquinone-10 on the stability
701 of biomimetic membranes of relevance for the inner mitochondrial membrane. *Biochimica et
702 Biophysica Acta - Biomembranes*. **1860** (5), 1205–1215, 10.1016/j.bbamem.2018.02.015 (2018).
- 703 28. Seneviratne, R. *et al.* A reconstitution method for integral membrane proteins in hybrid

704 lipid-polymer vesicles for enhanced functional durability. *Methods.* 1–8,
705 10.1016/j.ymeth.2018.01.021 (2018).
706 29. Veshaguri, S. *et al.* Direct observation of proton pumping by a eukaryote P-type ATPase.
707 *Science.* **351** (6280), 1–6 (2016).
708

Figure 1

[Click here to access/download;Figure;Figure1.TIF](#)







Figure 3

[Click here to access/download;Figure;Figure3.TIF](#)




Figure 4

[Click here to access/download;Figure;Figure4.TIF](#)







Figure 5




Figure 6

[Click here to access/download;Figure;Figure6.TIF](#)








Figure 7

[Click here to access/download;Figure;Figure7.TIF](#)





Click here to access/download
Video or Animated Figure
Movie_200ms.mp4



Name of Material/ Equipment	Company	Catalog Number
6-Mercapto-1-hexanol (6MH)	Sigma	451088
8-hydroxypyrene-1,3,6-trisulphonic acid (HPTS)	BioChemika	56360
Aluminium holder (Electrochemical cell)	Custom-made	
Auxiliary electrode		
Chloroform	VWR Chemicals	83627
E.coli polar lipids	Avanti	100600C
Epoxy	EPO-TEK	301-2FL
Fiji (ImageJ 1.52d)		
Filter cube ("ATTO633")	Chroma Technology Corporation	
Filter cube ("HPTS1")	Chroma Technology Corporation	
Filter cube ("HPTS2")	Chroma Technology Corporation	
Filter cube ("Texas Red")	Chroma Technology Corporation	
Fluorescent dye-labelled lipids (FDLL)	ThermoFisher Scientific	T1395MP
Fluorescent dye-labelled lipids (FDLL) (alternative)	ATTO-TEC	AD 633-161
Gel filtration column	GE Healthcare	28-9893-33
Glass coverslips	VWR International	631-0172
Glass syringe	Hamilton	1725 RNR
Glass vials	Scientific Glass Laboratories Ltd	T101/V1
Gold	GoodFellow	
Gramacidin	Sigma	G5002
Mercury sulfate reference electrode	Radiometer (Hash)	E21M012
Microcentrifuge	Eppendorf	Minispin NL040
Microscope	Nikon	Eclipse Ti
Microscope Camera	Andor	Zyla 5.5 sCMOS
Microscope Lamp	Nikon	Intensilight C-HGFI
NIS-Elements AR 5.0.2	Nikon	
n-Octyl β -D-glucopyranoside	Melford Laboratories	B2007
Nova 1.10	Metrohm	
Objective	Nikon	Plan Apo λ 60x/1.4 oil
OriginPro 2017	OriginLab	
O-ring (Electrochemical cell)	Orinoko	
Plastic tubes	Eppendorf	3810X

Polystyrene microbeads	Bio-RAD	152-3920
Potentiostat	Metrohm Autolab	PGSTAT 128N
Potentiostat	CH Instruments	CHI604C
Scripting software	Matlab	R2017a
Silicon wafers	IDB Technologies LTD	Si-C2 (N<100>P)
Teflon cell (Electrochemical cell)	Custom-made	
Temple-Stripped Ultra-Flat Gold Surfaces	Platypus Technologies	AU.1000.SWTSG
Thin micropipette tips	Sarstedt	70.1190.100
Ubiquinone-10	Sigma	C-9538
Ultracentrifuge	Beckman-Coulter	L-80XP
Ultrasonic bath	Fisher Scientific	FB15063

[illegible]

[illegible]

[illegible]

[illegible]

ARTICLE AND VIDEO LICENSE AGREEMENT - UK

Title of Article:

Single liposome measurements for the study of proton-pumping membrane enzymes using a combination of electrochemistry.....

Author(s):

Ievgen Maruzenko and Lars J.C. Jeuken.

Item 1: The Author elects to have the Materials be made available (as described at <http://www.jove.com/publish>) via:

☐ Standard Access

☒ Open Access

Item 2: Please select one of the following items:

☒ The Author is **NOT** a United States government employee.

☐ The Author is a United States government employee and the Materials were prepared in the course of his or her duties as a United States government employee.

☐ The Author is a United States government employee but the Materials were NOT prepared in the course of his or her duties as a United States government employee.

ARTICLE AND VIDEO LICENSE AGREEMENT

1. **Defined Terms.** As used in this Article and Video License Agreement, the following terms shall have the following meanings: "**Agreement**" means this Article and Video License Agreement; "**Article**" means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein; "**Author**" means the author who is a signatory to this Agreement; "**Collective Work**" means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; "**CRC License**" means the Creative Commons Attribution 3.0 Agreement (also known as CC-BY), the terms and conditions of which can be found at: <http://creativecommons.org/licenses/by/3.0/us/legalcode>; "**Derivative Work**" means a work based upon the Materials or upon the Materials and other pre-existing works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; "**Institution**" means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; "**JoVE**" means MyJoVE Corporation, a Massachusetts corporation and the publisher of The Journal of Visualized Experiments; "**Materials**" means the Article and / or the Video; "**Parties**" means the Author and JoVE; "**Video**" means any video(s) made by the Author, alone or in conjunction with any other parties, or by JoVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion

of the Article, and in which the Author may or may not appear.

2. **Background.** The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JoVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.

3. **Grant of Rights in Article.** In consideration of JoVE agreeing to publish the Article, the Author hereby grants to JoVE, subject to **Sections 4 and 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the "Open Access" box has been checked in **Item 1** above, JoVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.

4. **Retention of Rights in Article.** Notwithstanding the exclusive license granted to JoVE in **Section 3** above, the Author shall, with respect to the Article, retain the non-exclusive right to use all or part of the Article for the non-commercial purpose of giving lectures, presentations or teaching classes, and to post a copy of the Article on the Institution's website or the Author's personal website, in each case provided that a link to the Article on the JoVE website is provided and notice of JoVE's copyright in the Article is included. All non-copyright intellectual property rights in and to the Article, such as patent rights, shall remain with the Author.

5. **Grant of Rights in Video - Standard Access.** This **Section 5** applies if the "Standard Access" box has been checked in **Item 1** above or if no box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby acknowledges and agrees that, subject to **Section 7** below, JoVE is and shall be the sole and exclusive owner of all rights of any nature, including, without limitation, all copyrights, in and to the Video. To the extent that, by law, the Author is deemed, now or at any time in the future, to have any rights of any nature in or to the Video, the Author hereby disclaims all such rights and transfers all such rights to JoVE.

6. **Grant of Rights in Video - Open Access.** This **Section 6** applies only if the "Open Access" box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby grants to JoVE, subject to **Section 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Video in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Video into other languages, create adaptations, summaries or extracts of the Video or other Derivative Works or Collective Works based on all or any portion of the Video and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats.

7. **Government Employees.** If the Author is a United States government employee and the Article was prepared in the course of his or her duties as a United States government employee, as indicated in **Item 2** above, and any of the licenses or grants granted by the Author hereunder exceed the scope of the 17 U.S.C. 403, then the rights granted hereunder shall be limited to the maximum rights permitted under such statute. In such case, all provisions contained herein that are not in conflict with such statute shall remain in full force and effect, and all provisions contained herein that do so conflict shall be

deemed to be amended so as to provide to JoVE the maximum rights permissible within such statute.

8. **Protection of the work.** The Author(s) authorize JoVE to take steps in the Author(s) name and on their behalf if JoVE believes some third party could be infringing or might infringe the copyright of either the Author's Article and/or Video.

9. **Likeness, Privacy, Personality.** The Author hereby grants JoVE the right to use the Author's name, voice, likeness, picture, photograph, image, biography and performance in any way, commercial or otherwise, in connection with the Materials and the sale, promotion and distribution thereof. The Author hereby waives any and all rights he or she may have, relating to his or her appearance in the Video or otherwise relating to the Materials, under all applicable privacy, likeness, personality or similar laws.

10. **Author Warranties.** The Author represents and warrants that the Article is original, that it has not been published, that the copyright interest is owned by the Author (or, if more than one author is listed at the beginning of this Agreement, by such authors collectively) and has not been assigned, licensed, or otherwise transferred to any other party. The Author represents and warrants that the author(s) listed at the top of this Agreement are the only authors of the Materials. If more than one author is listed at the top of this Agreement and if any such author has not entered into a separate Article and Video License Agreement with JoVE relating to the Materials, the Author represents and warrants that the Author has been authorized by each of the other such authors to execute this Agreement on his or her behalf and to bind him or her with respect to the terms of this Agreement as if each of them had been a party hereto as an Author. The Author warrants that the use, reproduction, distribution, public or private performance or display, and/or modification of all or any portion of the Materials does not and will not violate, infringe and/or misappropriate the patent, trademark, intellectual property or other rights of any third party. The Author represents and warrants that it has and will continue to comply with all government, institutional and other regulations, including, without limitation all institutional, laboratory, hospital, ethical, human and animal treatment, privacy, and all other rules, regulations, laws, procedures or guidelines, applicable to the Materials, and that all research involving human and animal subjects has been approved by the Author's relevant institutional review board.

11. **JoVE Discretion.** If the Author requests the assistance of JoVE in producing the Video in the Author's facility, the Author shall ensure that the presence of JoVE employees, agents or independent contractors is in accordance with the relevant regulations of the Author's institution. If more than one author is listed at the beginning of this Agreement, JoVE may, in its sole discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JoVE reserves the right, in its absolute and sole

discretion and without giving any reason therefore, to accept or decline any work submitted to JoVE. JoVE and its employees, agents and independent contractors shall have full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JoVE has sole discretion as to the method of making and publishing the Materials, including, without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

12. **Indemnification.** The Author agrees to indemnify JoVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JoVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JoVE, making of videos by JoVE, or publication in JoVE or elsewhere by JoVE. The Author shall be responsible for, and shall hold JoVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to the making of a video by JoVE its employees, agents or independent contractors. All sterilization, cleanliness or

decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's expense. All indemnifications provided herein shall include JoVE's attorney's fees and costs related to said losses or damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JoVE, its employees, agents or independent contractors.

13. **Fees.** To cover the cost incurred for publication, JoVE must receive payment before production and publication the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JoVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.

14. **Transfer, Governing Law.** This Agreement may be assigned by JoVE and shall inure to the benefits of any of JoVE's successors and assignees. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to be one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement is required per submission.

CORRESPONDING AUTHOR

Name:

Lars J.C. Jeulien

Department:

School of Biomedical Science

Institution:

University of Leeds

Title:

Professor of Molecular Biophysics

Signature:



Date:

2 Aug. 2018

Please submit a **signed** and **dated** copy of this license by one of the following three methods:

1. Upload an electronic version on the JoVE submission site
2. Fax the document to +1.866.381.2236
3. Mail the document to JoVE / Attn: JoVE Editorial / 1 Alewife Center #200 / Cambridge, MA 02140

Editorial comments:**1. The editor has formatted the manuscript to match the journal's style. Please retain the same.**

We have kept the formatting.

2. Please address all the specific comments marked in the manuscript.

We have addressed all the comments mentioned in the manuscript. Please, see the replies to these comments as well as tracked changes in the manuscript.

3. For the units please use the format mg/mL and not mg/ml or mg ml⁻¹.

We have changed all the units according to suggestions.

4. Please combine some action steps so there are 2-3 actions per step in order to have a 2.75-page limit highlight.

We have merged steps 2.5.1-2.5.4 into a single step 2.6 and now adhere to the 2.75-page highlight limit.

5. Figure 1: Please remove the word teflon.

We have replaced the word Teflon for Polytetrafluoroethylene.

6. Other figures: Please format the units as suggested in comment 3.

We have changed all the units in the figures



Click here to access/download
Supplemental Coding Files
Main_Matlab_LiveScript.mlx






Click here to access/download
Supplemental Coding Files
Main_Matlab_Script.m






Click here to access/download
Supplemental Coding Files
Python_Script_for_Fiji.py








Click here to access/download
Supplemental Coding Files
Batch.m





Click here to access/download
Supplemental Coding Files
FastPeakFind.m





[Click here to access/download](#)
Supplemental Coding Files
FastPeakFind_license.txt






Click here to access/download
Supplemental Coding Files
fCalculatepH.m








Click here to access/download
Supplemental Coding Files
fCalculateThreshold.m








Click here to access/download
Supplemental Coding Files
fConvertVesicles.m








[Click here to access/download](#)
Supplemental Coding Files
fCreateVesicles.m








Click here to access/download
Supplemental Coding Files
fFilterTraces.m







[Click here to access/download](#)
Supplemental Coding Files
fGauss2D.m






Click here to access/download
Supplemental Coding Files
fGauss2DFit.m





Click here to access/download
Supplemental Coding Files
fGauss2DFitFrames.m






Click here to access/download
Supplemental Coding Files
fGetATTO.m







[Click here to access/download](#)
Supplemental Coding Files
fGetMedian.m







[Click here to access/download](#)
Supplemental Coding Files
fSaveTraces.m





Click here to access/download
Supplemental Coding Files
fSaveVesicles.m






[Click here to access/download](#)
Supplemental Coding Files
fSeeRandomTraces.m








Click here to access/download
Supplemental Coding Files
fSeeRandomVesicles.m







[Click here to access/download](#)
Supplemental Coding Files
fSeeVesicle.m





Click here to access/download
Supplemental Coding Files
fTifLoad.m






Click here to access/download
Supplemental Coding Files
fVesicleFilter.m






[Click here to access/download](#)



Supplemental Coding Files
fVesicleFilterFrames.m






[Click here to access/download](#)

Supplemental Coding Files
fWriteGif.m





[Click here to access/download](#)
Supplemental Coding Files
uigetfile_n_dir.m