

# Journal of Visualized Experiments

## Live Cell Analysis of Shear Stress on *Pseudomonas aeruginosa* using the Automated Higher-Throughput Microfluidic System

--Manuscript Draft--

<b>Article Type:</b>	Methods Article - JoVE Produced Video
<b>Manuscript Number:</b>	JoVE58926R2
<b>Full Title:</b>	Live Cell Analysis of Shear Stress on <i>Pseudomonas aeruginosa</i> using the Automated Higher-Throughput Microfluidic System
<b>Keywords:</b>	<i>Pseudomonas aeruginosa</i> ; enhanced green fluorescent protein; Microfluidics; shear; biofilm; minimal media; Fluorescence microscopy; bright-field microscopy; biofilm coverage; Image analysis; Montage software
<b>Corresponding Author:</b>	Arin L Sutlief Doane University Crete, NE UNITED STATES
<b>Corresponding Author's Institution:</b>	Doane University
<b>Corresponding Author E-Mail:</b>	arin.sutlief@doane.edu
<b>Order of Authors:</b>	Arin L Sutlief Helena Valquier-Flynn Christina Wilson Marco Perez Hunter Kleinschmidt Brett J. Schofield Elizabeth Delmain Andrea E. Holmes Christopher D. Wentworth
<b>Additional Information:</b>	
<b>Question</b>	<b>Response</b>
Please indicate whether this article will be Standard Access or Open Access.	Open Access (US\$4,200)
Please indicate the <b>city, state/province, and country</b> where this article will be <b>filmed</b> . Please do not use abbreviations.	68333

**TITLE:**

**Live Cell Analysis of Shear Stress on *Pseudomonas aeruginosa* using an Automated Medium-Throughput Microfluidic System**

**AUTHORS AND AFFILIATIONS:**

Arin L. Sutlief<sup>1\*</sup>, Helena Valquier-Flynn<sup>1</sup>, Christina Wilson<sup>1</sup>, Marco Perez<sup>1</sup>, Hunter Kleinschmidt<sup>2</sup>, Brett J. Schofield<sup>2</sup>, Elizabeth Delmain<sup>3</sup>, Andrea E. Holmes<sup>1</sup>, Christopher D. Wentworth<sup>4\*</sup>

<sup>1</sup>Department of Chemistry, Doane University, Crete, NE, USA

<sup>2</sup>Department of Biology, Doane University, Crete, NE, USA

<sup>3</sup>Department of Pathology and Microbiology, University of Nebraska Medical Center  
Omaha, NE, USA

<sup>4</sup>Department of Physics and Engineering, Crete, NE, USA

**Corresponding Authors:**

Christopher Wentworth (chris.wentworth@doane.edu)

Arin L. Sutlief (arin.sutlief@doane.edu)

Tel: 1-(402)-826-8325

**Co-authors:**

Helena Valquier-Flynn (helena.valquierflynn@doane.edu)

Christina Wilson (christina.wilson@doane.edu)

Marco Perez (marco.perez@doane.edu)

Hunter Kleinschmidt (hunter.kleinschmidt@doane.edu)

Brett Schofield (brett.schofield@doane.edu)

Elizabeth Delmain (elizabeth.delmain@unmc.edu)

Andrea E. Holmes (andrea.holmes@doane.edu)

**KEYWORDS:**

*Pseudomonas aeruginosa*, enhanced green fluorescent protein, microfluidics, shear, biofilm, minimal media, fluorescence microscopy, bright-field microscopy, biofilm coverage, image analysis

**SUMMARY:**

Here, we describe the use of a medium-throughput microfluidic bioreactor coupled with a fluorescent microscope for the analysis of shear stress effects on *Pseudomonas aeruginosa* biofilms expressing green fluorescent proteins, including instrument set up, the determination of biofilm coverage, growth rate, and morphological properties.

**ABSTRACT:**

A medium-throughput microfluidic *in vitro* bioreactor coupled with fluorescence microscopy has been used to study bacterial biofilm growth and morphology, including *Pseudomonas aeruginosa* (*P. aeruginosa*). Here, we will describe how the system can be used to study the growth kinetics

and the morphological properties such as the surface roughness and textural entropy of *P. aeruginosa* strain PA01 that expresses an enhanced green fluorescent protein (PA01-EGFP). A detailed protocol will describe how to grow and seed PA01-EGFP cultures, how to set up the microscope and autorun, and conduct the image analysis to determine growth rate and morphological properties using a variety of shear forces that are controlled by the microfluidic device. This article will provide a detailed description of a technique to improve the study of PA01-EGFP biofilms which eventually can be applied towards other strains of bacteria, fungi, or algae biofilms using the microfluidic platform.

## INTRODUCTION:

Here, we will demonstrate a method to measure the effect of shear stress on the formation of fluorescent *Pseudomonas aeruginosa* (*P. aeruginosa*) PA01 biofilms using an automated medium-throughput microfluidic system.

Biofilms are communities of microorganisms, such as bacteria, organized by an extracellular polymeric substance that are attached to a support, and are typically found at the interface between a liquid and a solid surface<sup>1</sup>. These biofilm communities can be beneficial to the environment, such as improving water quality in water supply lines and in bioremediation of recalcitrant compounds<sup>2,3</sup>. However, biofilms can also be highly harmful to human health with undesirable consequences. For example, medical devices, such as hip and knee implants, are one type of surface where biofilm accumulation has been a challenge and causes severe medical complications<sup>4,5</sup>. Biofilms can also enter natural water systems, such as rivers and lakes, and infiltrate water supply pipes leading to bacteria contamination in drinking water resulting in infections<sup>6-8</sup>. Biofilms formed in marine environments adhere to ships and other man-made substrates and present a major economic and environmental problem as increased friction leads to increase fuel consumption<sup>9,10</sup>. Antimicrobial coatings, such as Tributyltin, have been developed to prevent these problems but are toxic to marine life<sup>11</sup>.

*P. aeruginosa* is a Gram-negative bacterium with high thriving capabilities in a variety of environmental and nutritional conditions<sup>12</sup>. *P. aeruginosa* is a common cause of community- and hospital-acquired infections and found to be closely associated to injuries, such as severe burns, and immunocompromised hosts, such as in cystic fibrosis (CF)<sup>5,12,13</sup>, AIDS, and cancer patients<sup>5,13</sup>. The formation of *P. aeruginosa* biofilms has been most seriously connected to CF, where chronic lung infections are the leading cause of death for this disease<sup>5</sup>.

A reference strain of *P. aeruginosa*, PA01, is used in this report and is genetically modified to express an enhanced green fluorescent protein (PA-EGFP). EGFP represents a mutant form of GFP with greater fluorescence properties which allows for *in situ* biofilm analysis using fluorescence microscopy<sup>14-16</sup>. This type of fluorescence analysis is advantageous for the study of biofilms because GFP does not interfere significantly with cell growth and function<sup>17</sup>. For example, *Escherichia coli* cells that were tagged with GFP grew well and continuously without having suffered any toxic effects in comparison to the control bacteria<sup>17</sup>. Other reports substantiate this claim<sup>18-20</sup>. Furthermore, the use of a fluorescent reporter such as EGFP is quick and simple, yet only live cells will be measured because dead cells quickly cease to fluoresce<sup>21</sup>.

Biofilms can grow under various environmental conditions including those with different flow rates. For example, films can grow in high shear stress, such as in rivers, where the high water flow conditions lead to a greater microbial diversity<sup>22</sup>. Contrarily, standing water in ponds or oral biofilms experience a much lower shear force<sup>23</sup>. In addition to flow rate, there are other factors that influence biofilm adhesion, including surface roughness and hydrophobicity, media composition, and even the bacterial cell surface<sup>1,4,7,24</sup>. Conditions can also cause variation in the spatial structure or morphology of a biofilm. This includes environmental conditions such as shear stress exerted by a moving fluid or gradients in nutrient availability and biological factors such as the species present in the system, motility of cells, and the specific proteins present in the extracellular polymeric substance<sup>25–27</sup>. Under some conditions, the biofilm will be lawn-like (smooth and flat), while under other conditions the biofilm will be rough, fluffy, or even mushroom-like<sup>28</sup>. While the qualitative difference between biofilm lawns and mushroom structures can be clearly seen in microscopic images, understanding the relationship between film structure and biological processes within the film requires systematic and quantitative methods of describing the morphology. Morphological properties suggested for study by researchers include porosity, fractal dimension, diffusional length, microcolony area at the substratum, microcolony volume, roughness coefficient, and textural entropy<sup>29,30</sup>.

Bioreactors are used in the study of biofilms to mimic real life conditions<sup>31</sup>. Drip flow reactors (DFR) represent a low-shear environment where nutrients in media flow slowly across cells that are attaching to a surface over time to form a biofilm with a high cell density<sup>32</sup>. CDC reactors are bioreactors that create a high shear stress fluid environment by control of a stirring rod that is continuously rotating within the media filled tank<sup>33</sup>. These types of bioreactors are simple to set up, but they are limited in scope because of the relatively low sample size, the high consumption of media, large amounts of biohazardous waste produced from media drip flows ranging from 125  $\mu$ L/minute for drip flow reactors to more than 1 mL/min for CDC reactors, and the need to autoclave large amounts of glassware and waste media<sup>34</sup>. Biofilms do not grow evenly across the surface in a drip flow reactor because the low shear of media causes trailing along the larger conglomerates of *P. aeruginosa* bacteria therefore, the biofilm growth is not very smooth and the uneven samples cannot be easily analyzed with fluorescence microscopy<sup>35,36</sup>.

Some common bioreactor limitations are overcome by using a medium-throughput microfluidic bioreactor, where only milliliters of media are required, and the reaction plates are small and readily disposable after autoclaving<sup>37</sup>. Furthermore, depending on the number of wells, many replications can be performed in just one reactor run, which provides sufficient amount of data to conduct meaningful statistical analysis. In **Figure 1**, the different components of the microfluidic-microscopy system that allow for controlled conditions, including temperature and flow rate<sup>38–40</sup>, are shown. The bioreactor is coupled with fluorescence microscopy to visualize the fluorescence of the EGFP tag in PA01 under applied low through high shear conditions that will mimic more realistic scenarios that are encountered in the environment or in the biomedical field.

[Place Figure 1 here]



An excerpt of the microfluidic plate is shown in **Figure 2**. The most commonly used plates consist of 48 wells. One experiment requires one inlet and one outlet wells, a total of 2 wells. This allows for 24 simultaneous experiments that can be performed with various experimental conditions, such as bacterial strains, antimicrobial treatments, and media varied from channel to channel and controlled shear flow for each column of six channels. The experimental temperature is also controlled with one temperature setting throughout the plate. The microfluidic channels show that each channel has a serpentine region to provide sufficient back pressure and controlled shear.

[Place Figure 2 here]

A step by step guide is provided to allow users of microfluidic bioreactors that are coupled with fluorescence microscopy to conduct novel biofilm experiments using different shear environments. This method will allow for the expansion of experiments involving other microorganisms besides bacteria, such as fungi and algae, that have medical and environmental applications<sup>41–43</sup>. The detailed approach describes how to culture PA01-EGFP, inoculate a 48-well plate and set up the microfluidic device and software, set up the fluorescent microscope, and demonstrate the software analysis to obtain the biofilm coverage, growth rate, and morphological properties such as surface roughness.

## **PROTOCOL:**

### **1. Media Preparation**

1.1. Prepare minimal media (MM) with 0.25% glucose. To make 1 L of MM with 0.25% glucose, add 200 mL of sterile M9 salt solution, 2 mL of sterile 1 M  $\text{MgSO}_4$ , 100  $\mu\text{L}$  of sterile 1 M  $\text{CaCl}_2$ , and 12.5 mL of sterile 20% (w/v) glucose to water ( $\text{dH}_2\text{O}$ ) at a final volume of 1 L.

1.2. Transfer required media to a sterile bottle using sterile techniques. Prepare the amount required depending on the number of channels being used. Typically, each channel requires 200  $\mu\text{L}$  for priming, 300  $\mu\text{L}$  for seeding and 1300  $\mu\text{L}$  for the 24 h experiment.

1.3. Place the bottle in an incubator or water bath set at the experimental temperature. The media should be at the experiment temperature before use to avoid bubbles forming in the microchannels.

### **2. Preparation of an Overnight and Experimental Culture of PA01-EGFP.**

2.1. Place 15 mL of experimental media into a sterile sidearm flask and inoculate with one or two colonies from an agar plate streaked with PA01-EGFP. Expand the culture for 12-16 h on an incubator/shaker table at 37 °C and 180-220 rpm.

2.2. Measure the  $\text{OD}_{600}$  of the overnight culture. When the  $\text{OD}_{600}$  is greater than 0.80, dilute the

overnight culture to the final OD of 0.8 using fresh MM. Place the experimental culture in an incubator or water bath at 37 °C until needed for seeding the microfluidic channels. Check the OD<sub>600</sub> again immediately before seeding to ensure that it has not changed significantly from the target OD<sub>600</sub>, 0.8 in this case.

### 3. Equipment Startup

3.1. Set up the system station according to the user guide, similar to the setup in **Figure 1**. To avoid an error in the connection of the instrument to the software, turn on the instrument in the following order:

PC Workstation

Fluorescence Module. Make sure the fluorescence shutter is ON (blue light by shutter button on)

Hardware Controllers

Imaging system Controller (see **Table of Materials**)

CCD Camera

Imaging Station (Microscope)

NOTE: The temperature of the heated plate should be adjusted to the desired experimental temperature.

3.2. Start the control application and enter the plate number located on the label on the side of the plate.

NOTE: After control application starts there will be two separate application windows, one for the software that controls the microscope/imaging software and one for the control module that controls the pump and the well-plate interface. This is shown in **Figure S1**, where the “Multi Dimensional Acquisition” menu is open on the top and the Control Module is set for an autorun on the bottom.

### 4. Priming and Seeding the Microfluidic Plate

NOTE: The priming, seeding, bacterial attachment and growth are illustrated in **Figure 3**.

4.1. Remove the 48-well microfluidic plate from the packaging making sure to not touch the glass surface at the bottom of the plate. Clean the glass slide at the bottom of the plate well with a lens tissue, a lint free cloth, or a low-lint wipe.

4.2. To prime the microfluidic channels, pipette 200 µL of 37 °C MM into the output well, being careful to avoid bubbles. Place the plate into the plate stage and wipe the interface with ethanol, allowing to dry, before sealing it onto the plate stage.

4.3. In **Manual Mode** on the Control Module, set **Fluid** as LB at 37 °C and **Max Shear** at 5.00 dyne/cm<sup>2</sup>. Click the output wells to activate flow from the output to the input well, priming the channels. After the 5 min priming, pause the flow to prepare for seeding. Carefully remove the

plate from the stage and pipette residual media from the output but not removing any of the media from the inner circle that leads to the microfluidic channels.

4.4. To seed the experimental channels, first pipette 300  $\mu\text{L}$  of MM into the input well followed by pipetting 300  $\mu\text{L}$  of the bacterial culture into the output well into the output well. Place the plate back into the plate stage making sure to wipe the interface before placing it on the plate.

4.5. On the Control Module, focus on a single channel using the live camera feed after placing the plate stage onto the microscope stage. While visually monitoring by the live feed, resume flow at 1.00-2.00 dyne/cm<sup>2</sup> for approximately 2-4 s to allow cells to enter the experimental channel but not in the serpentine channels. Leave the plate on the temperature controlled stage for 1 h to allow for cell attachment. During the 1 h incubation, the Control Module and Montage software can be set up for the autorun (Step 5).

NOTE: The amount of time needed for seeding will vary with media and organism, so it should be monitored closely by the live feed until optimized and used as a general time frame. Seeding throughout the plate may vary which would require more time of applied flow to certain columns of channels for complete seeding.

4.6. After the attachment period, gently remove the plate from the stage and pipette the bacteria from the output first avoiding disturbing the channel. With a new pipette tip, remove the media from the input wells.

[Place Figure 3 here]

## 5. Setting the Software

5.1. In the software, open “**Multiple Dimensional Acquisition**” to control the microscope image acquisition. Under the “**Main**” menu, select “**Timelapse**”, “**Multiple Stage Positions**”, and “**Multiple Wavelengths**” (Figure S1A).

5.1.1. Setting up the **Saving** settings, create a simple **Base Name**, making sure that **Increment base name if file exists** is checked. Click **Select Directory** to select the folder in which all files will be saved. Include essential details of the experiment in the **Description**.

5.1.2. Under the **Timelapse** tab, adjust the **Duration** of the experimental time to 24 h. To control how often images are acquired, set the **Time Interval** so images are acquired every 5 min throughout the experiment. The **Number of Time Points** automatically adjusts with the set parameters.

5.1.3. Set stage positions using the live camera feed in the bright field microscope, allowing for correct placement and focus. Starting with the 10X objective, focus on the center of the channel, located above or below the channel numbers engraved on the plate. Switch to the 20X objective, finding the optimal viewing area and focal plane within the channel. Add the position to the list

or replace the preset channel position with the new settings.

5.1.4. Under the **Wavelengths** menu, set the number of wavelengths to 3. For wavelength 1 (W1), select the FITC 100% CAM with a 10 ms exposure time. For wavelength 2 (W2), select the Brightfield 50% CAM 50% VIS with a minimal exposure time of 3 ms. For wavelength 3 (W3), set as the **All Closed** filter so no light remains on the last channel between acquisition times.

5.2. When setting up the Control Module, check to make sure that the manual run has been stopped.

5.2.1. In the **Edit Autorun** menu under the **Protocol Setup** tab, set a new protocol for 24 h of duration time, flowing in the forward direction at the desired shear rate. Shear rate is varied to study the effect of biofilm growth with shear rate. Click **Add** and **Save As** the protocol.

5.2.2. Under the **Sequence Setup** tab, create a new sequence by selecting **LB@37degrees** as the default fluid for all channels. Under **Step Iteration** 1, for channels 1-12, select the protocol with the desired shear rate and **enable** all channels. For channels 13-24, select the protocol with the second desired shear rate and **enable** all channels. Select **Apply** and **Save As** the sequence.

5.2.3. In the **Autorun** menu, select the saved sequence to be used for the autorun.

## 6. Timed Biofilm Growth Experiment Setup

6.1. Pipette a maximum of 1300  $\mu$ L of sterile MM into the input well of the microfluidic plate. Place the plate back onto the plate stage and wipe the interface with ethanol, allowing ethanol to dry completely, before sealing the plate. Recheck the saved stage positions after placing the plate back onto the plate stage to ensure the saved positions are still centered within the viewing channels.

6.2. Placing the plate stage onto the microscope stage, and be sure that the protocol(s) and sequence(s) are set up correctly. Select **Start** to start the autorun quickly followed by clicking **Acquire** to start the microscope image collection.

6.3. Adjust the focus and placement of all stage positions between image acquisitions. Select **Pause** and, using the live image mode in the bright field wavelength, select **Go To** to view each stage positions. Set new settings by selecting **Set to Current**. Click **Resume** before the next scheduled acquisition time.

## 7. Review and Analyze Image Sequences after the Timed Biofilm Growth Experiment

7.1. To review each channel following the 24 h autorun, in the software, open **Review Multidimensional Data**, located under **Analysis Tools**. Click on **Select Base File | Select Directory** to navigate to the folder with image sequences for the experiment of interest. In the list of **Data Sets** in the folder that appears, select the base name of the data of interest.

7.2. Once the data of interest is selected, click **View** to review this data. Select a wavelength to review that data under **Wavelength** window (**Figure S2**).

7.3. Review the image sequence for each channel by selecting the channel number under **Stage Position** and using the video controls to analyze the data for the 289 time points (assuming that images are obtained every 5 min over a 24 h time frame). Take note of useable growth data.

NOTE: In each channel, watch for the development of air bubbles and clogs that will disturb the biofilm growth within the channel, affecting the data. However, if these occur later in the run, data prior to these circumstances may be found useful.

7.4. After reviewing and selecting the desired data for analysis, create stacks of the images for each channel. Under the **File** menu, select **Open Special | Build Stack | Numbered Names**. In the **Build Stack** window, select the **Select First Image** button then select the first image for the stack in the file folder; select the **Select Last Image** button and then select the file corresponding to the last image in the stack. Click on **OK** to open the stack with all the images of the channel in chronological order. Stacks can be saved under the **File** menu, **Save As** as tiff files.

NOTE: Do not create stacks with the thumb images. These files are not very useful and can be deleted to save room on the hard drive. Also, image files are saved as the following:

BASE NAME\_WAVELENGTH\_sCHANNEL#\_tIMAGE#.

For example, 07062018\_W1FITC 100\_s12\_t112 is for the channel 12, time point 112, image in the first wavelength-FITC 100% camera data with the base name "07062018".

7.5. Finally, select **Save As** from the **File** drop-down menu to save the sequence as a tiff stack. Stacks can also be overlaid and saved as a movie. This is addressed in step 9.

7.6. Before quantifying the % surface coverage, calibrate the image distances. Under **Analysis Tools**, click on **Calibrate Distances**. Select the appropriate calibration measurement and click **Apply to All Open Images**.

7.7. While the stack is on the first image of the sequence, click on the **Threshold** button. Use **Auto Threshold for Light Objects** to threshold for fluorescent signal (FITC wavelength) and **Auto Threshold for Dark Objects** to threshold in bright field. Adjust threshold to the coverage that represents the coverage of the image.

7.7.1. For fluorescent thresholds, utilize a channel with non-fluorescent cells to establish any background signal and set the minimum threshold value to exclude signal detected in channels containing these non-fluorescent cells to ensure the signal measurement does not overestimate the area of fluorescence. This may be different from run to run so it is a good idea to use at least one, if not more, channel(s) for non-fluorescent controls.

7.7.2. For bright field thresholds, if all cells are not covered by the orange threshold signature, adjust the maximum threshold value either via the sliding toolbar or using the **Threshold Image** (found under the **Analysis Tools**) so that all cells are covered but any background is excluded.

NOTE: While these threshold values ideally would be used for all images within a stack and for all channels, threshold values selected for one image/channel may not be appropriate for other images or channels, therefore the user may need to adjust the threshold range periodically. For this reason, it is a good idea to always record the maximum and minimum threshold values used should the data need to be reviewed later.

7.8. To quantify the coverage, under **Analysis Tools**, click **Show Region Statistics**. Make sure that **Use Threshold** is checked, the **Entire Image** is selected, and under **Obtained Data**, make sure the desired measurements/settings are selected (threshold area, average, standard deviation, min, max, and % threshold area). Click on **Open Log** and make sure DDE file is selected before clicking **OK**. Select **Microsoft Excel** and click **OK**. In the spreadsheet that will open, click **Log Data** to automatically record the measurements of the image being analyzed.

7.9. Leaving this spreadsheet open, use the same sheet to collect all image analysis data for a single stack. In the stack, go to the fourth image and threshold to the optimal settings. Click on **Log Data** in the **Show Region Statistics** screen and the values are logged into the spreadsheet. Repeat this for analysis of images every 15 min.

## 8. Other Analysis Including Morphology and Surface Coverage Measure Using Python scripts from the Biofilm Morphology Suite

8.1. Install a distribution of Python 3.6 that includes standard scientific modules by using a standard scientific Python distribution such as Anaconda, available at <https://www.anaconda.com/download>.

8.2. Obtain the Biofilm Morphology Suite from GitHub in a browser by navigating to <https://github.com/cdwentworth/Biofilm-Morphology-Suite.git>, then select **Clone or download** (green button) and select **Download Zip**. Unzip the file and move the code folders to a working directory.

8.3. Measure percent thresholded coverage by copying the image stack into the “**Coverage**” folder. Use a terminal window to navigate to that folder and then execute the script with the command

```
python bfCoverage.py tiffStackName.tif
```

from the command line interface, where *tiffStackName.tif* is the name of the file containing the tiff stack of images. A text file containing the coverage measurements at each time point will be created with the name *tiffStackName.txt* and a graphics file with a plot of coverage as a function

of time will be created with the name *tiffStackName.png*.

8.4. Use the same procedure specified in step 8.3 to measure the other morphological characteristics: biofilm accumulation, roughness coefficient, and textural entropy. Use the bfAcc.py script for the accumulation measurement, the roughCoef.py script for the roughness coefficient measurement, and the TEvsTime.py script for the textural entropy measurement.

## 9. Other Software Applications – Making Overlays and Movies

9.1. If multiple wavelengths are used, create overlay stacks including these wavelengths. Open all stacks of interest and calibrate the distance as done in step 7.6.

9.1.1. Under **Analysis Tools**, select **Overlay Images**. Set the sources as the stacks of interest. In the FITC stack, click on the rainbow circle and select green color filter to be added in order to highlight this wavelength.

9.1.2. Use **bal** to adjust the overlay so that the FITC stack is apparent in the images. After all adjustments have been made, select **All Planes** for both stacks and click **Apply**. Save the overlays in the same way as stacks described in step 7.5 or as a movie described in step 9.2.

9.2. To save stacks or overlays as a timelapse movie, under **MM Standard** and **Stack**, click on **Make Movie**. Select the source as the stack or overlay that is desired. Click **Save** and save as a **Microsoft Video 1** with a quality set between 70-80.

## REPRESENTATIVE RESULTS:

**Figure 4a** demonstrates the percent thresholded area over time from a 24 h run at shear flows of 0.2, 0.5, and 1.0 dyne/cm<sup>2</sup>. The biofilm coverage or percent threshold surface area (*C* [%]) was different for all three shear settings. The biofilm coverage was fastest at 1.0 dyne/cm<sup>2</sup> shear, where the threshold area increased from 2-5% to 100% after 200 min of growth and reached a stationary phase after 400 min. At 0.5 dyne/cm<sup>2</sup>, the biofilm coverage was delayed and began to increase at 400 min reaching 100% coverage after 800 min. The lowest shear at 0.2 dyne/cm<sup>2</sup> clearly demonstrated the slowest increase in biofilm coverage, where percent coverage started to increase at 500 min but never reached beyond the 65% threshold area. These results indicated shear had a direct influence in biofilm surface coverage. The higher shear seemed to be a more optimal condition for biofilm growth, possibly due to the fact that the media provided the bacteria with more nutrients so that the biofilm can proliferate quicker.

[Place Figure 4 here]

**Figure 4b** is consistent with the findings shown in **Figure 4a**. The biofilm accumulation (*N*[rel]) was measured based on the assumption that the GFP signal at a point in an image is proportional to the live cell density at that position. It demonstrates that total relative measurement biofilm accumulation increased as a function of time using 0.2, 0.5, and 1.0 dyne/cm<sup>2</sup> over a 24 h time period, and the growth rate declined from highest shear stress to lowest shear stress. There is a

clear time period giving exponential growth from which a quantitative growth rate can be calculated.

**Figure 4c** demonstrates the roughness coefficient ( $R_a$ ) at shear flows of 0.2, 0.5, and 1.0 dyne/cm<sup>2</sup>. Surface roughness, quantified by the roughness coefficient, measures the variance in the thickness profile of the film. The formal definition is

$$R = \frac{1}{N} \sum_{i=1}^N \frac{|T_a - T_i|}{T_a}$$

where  $T_i$  is the  $i$ -th thickness measurement,  $T_a$  is the average thickness, and  $N$  is the number of thickness measurement<sup>30</sup>. The procedure described in this investigation measures the thickness associated with living cells. A flat arrangement of cells will yield a roughness coefficient of zero while significant variations of thickness from the average will yield a roughness coefficient greater than one. Similar to the shear influence on the growth rate and the percent threshold coverage of the film, the biofilms exhibited different topographies over time. Overall,  $R_a$  decreased over time for all shear conditions indicating that all surfaces became smoother. However, in comparison to the lowest shear of 0.2 dyne/cm, the higher shear settings of 0.5 and 1.0 dyne/cm<sup>2</sup> resulted in a smoother surface over time indicating that a faster shear flow contributed to smoother and more even surface, and the highest shear of 1.0 dyne/cm<sup>2</sup> reaching below 0.2  $R_a$ .

The smoothness, regularity, or coarseness of a surface can also be expressed in textural entropy (TE). TE is a property used in image analysis to measure the degree of randomness in a two-dimensional image. Its calculation is based on the gray level co-occurrence matrix, defined by Haralick *et al.*, which looks at whether pixel values at one location are correlated with pixel values at another location<sup>44</sup>. A high degree of correlation will lead to low entropy. **Figure 4d** depicts the TE at shear flows of 0.2, 0.5, and 1.0 dyne/cm<sup>2</sup>. The TE increased over time for all shear conditions but the highest shear stress at 1.0 dyne/cm reached the maximum TE (1.0) earlier than the lower shear stresses at 900 min. The lowest shear of 0.2 dyne/cm<sup>2</sup> had the lowest TE (0.8) reaching its maximum after 1000 min. However, the intermediate shear stress of 0.5 dyne/cm<sup>2</sup> reached its maximum TE (1.2) much later than the high or low shear stress conditions.

The roughness coefficient and TE measure different features. While a flat film would have a low roughness coefficient and low entropy, a film with significant variation in thickness would have a high roughness coefficient but could still be low entropy if the variation is sinusoidal rather than random. In this case,  $R_a$  decreased with increased shear stress and time while the TE trends cannot be directly related to shear stress applied to biofilm formation over time.

#### FIGURE AND TABLE LEGENDS:

**Figure 1: Individual components of the Microfluidic System.** The individual components are listed from left to right: 1. CCD Camera, 2. High Resolution Inverted Microscope with Automated Stage, Automated Fluorescence Module, and Autofocus Module, 3. Plate Stage, 4: Imaging



system Interface, 5: Manual Microscope Stage Control, 6: Vapor Trap, 7: Imaging system Controller (including Temperature Controller), 8: Hardware Controllers, 9: Fluorescence Controller, 10: Uninterruptible Power Supply, 11: External Hard Drive for Image Storage, 12: PC Workstation.

**Figure 2: Visualization of the microfluidic channels and viewing window.** Two inlet and outlet wells with the microfluidic channels connecting them are highlighted with red and green dyes. The dye makes visible a serpentine region in each channel which creates sufficient back pressure and controlled shear during fluid flow. Each viewing channel (within the red circle) can be imaged with desired wavelengths. Shown are bright field (top) and fluorescent (bottom) microscopy images of a single channel with a PA01-EGFP biofilm using a 20X objective. Scale bar = 80  $\mu\text{m}$ .

**Figure 3: Experimental overview of priming, seeding, and attachment of PA01-EGFP in the microfluidic channels.** The priming, seeding, and attachment is outlined. The first step of priming requires fresh media introduced to the output. The seeding entails equal volumes of media and bacterial culture in the input and output wells, respectively. The culture should not pass the viewing segment of the experimental channel (red line) to avoid clogging the serpentine channels. After the incubation period is over, fresh media continuously flows from the input, to the viewing chamber and into the outlet. This initiates the attachment and growth of the bacterial biofilm.

**Figure 4: Percent threshold area, total biofilm accumulation, roughness coefficient, and textural entropy.** a) Percent threshold area ( $C$  [%]) over time using 0.2, 0.5, and 1.0  $\text{dyne/cm}^2$  over a 24 h time period. Data was obtained from one channel for each shear condition. b) Total biofilm accumulation (relative measure) as a function of time using 0.2, 0.5, and 1.0  $\text{dyne/cm}^2$  over a 24 h time period. The black lines are least-squares fits to an exponential model. Data was obtained from one channel for each shear condition. c) Roughness coefficient of PA01-EGFP using shear stress values of 0.2, 0.5, and 1.0  $\text{dyne/cm}^2$  monitored over 24 h. Data was obtained from one channel for each shear condition. d) Textural entropy of PA01-EGFP using shear stress values of 0.2, 0.5, and 1.0  $\text{dyne/cm}^2$ . Data was obtained from one channel for each shear condition (Valquier-Flynn, H., Sutlief, A.L., Wentworth, C.D., 2018).

**Figure S1: Image capture of the control software windows.** Software window with **Multi Dimensional Acquisition** menu open (Top). Control Module set for an AutoRun (Bottom).

**Figure S2: Image capture of the software for reviewing data.** Application window after selecting the Data Set of interest using the **Review Multi Dimensional Data** tool from the **Analysis Tools** menu.

## DISCUSSION:

The microfluidic system and image analysis procedure discussed here focuses on the execution of microfluidic biofilm experiments to determine morphological properties that do not require the full three-dimensional information typically found from confocal microscope studies. These included microcolony substratum coverage (percent coverage), surface roughness measured by

the roughness coefficient, and textural entropy. A method for estimating total relative biofilm cell accumulation is also presented from which growth rates in log phase can be calculated.

There are several small but significant steps that should be highlighted in this method. Wiping the interface with alcohol helps avoid contamination of other bacteria in going from experiment to experiment but also from well to well in one experiment. Priming and seeding are also very important because priming allows the user to determine which channels are allowing the media to flow through the channels without any disturbances or clogging. Channels should also not be disturbed (*i.e.* they should be constantly full of media) after priming to increase the chances of a successful experiment with no air bubbles or clogging. The seeding step can be varied according the type of bacteria and so should be optimized for cell attachment. For example, if cells do not seem to attach, surface modifications may have to occur on the microfluidic plate prior to seeding or longer incubations times may be needed. It is also critical to make sure that the microscope is correctly set up to obtain an in-focus image and should be monitored periodically throughout the experiment to assure that quality images are obtained. If the focus is off, the microscope can and should be adjusted as the experiment continues. During the image acquisition time, the last wavelength should be set to **All Closed** in order to avoid the exposure of filters and lighting to only one channel during the waiting time that occurs in between the image acquisitions. Also, the image analysis that determines the % surface coverage was designed in house because the Montage software manual did not explicitly describe the procedure. Furthermore, in order to expand on the image analysis and determine other characteristics such as surface roughness, *etc.*, open source Python-based code<sup>45</sup> was developed in house and shared on the gitHub repository. There are also limitations on how much data can be stored and managed on the local hard drive, so an external hard drive or online data sharing is needed such as CyVerse<sup>46</sup>.

Conventional bioreactors, such as the CDC reactor and the drip flow reactor<sup>34</sup>, require a lot media, provide fewer sample sizes, and demand a high amount of sterilization of equipment. In contrast, the advantages of this higher throughput platform include the ability to control shear, flow rates, and the assumption that the *in vitro* experiments closely resemble *in vivo* conditions. Disadvantages of the system include the multiple accessories and software that require a meticulous setup that must be performed in the correct order of events. Furthermore, the manual that is provided for the equipment does not fully explain each step of the experiments, and the software commands, and consequently, many mistakes occur during the experiments, including clogging of channels, lack of growth or attachment, or the lack high quality microscopy images or movies. The instrument itself and consumables, such as the microfluidic plates, are also relatively expensive with a price tag of over \$200 per plate and are not reusable. Thus, while the technique lends powerful results, the technical expertise required for its use is relatively high and requires repeated training by experts in the field. This report attempts to resolve this issue by providing a guide to new users of these bioreactors to study biofilms characteristics.

The microfluidic system, which is able to perform cellular analysis, has gained considerable attention for various scientific modalities, such as in microbiology, immunology, hematology, oncology, and stem cell research. More specifically, the technology has resulted in many publications describing topics that are highly relevant to medical applications<sup>37,47</sup>, including

microbial oral adhesion<sup>48</sup>, determining the effects of biosurfactants on *Pseudomonas aeruginosa* and *Staphylococcus aureus*<sup>49,50</sup>, host pathogen interactions in *E. coli*<sup>51</sup>, *Streptococcus* adherence<sup>52</sup>, and treatment of cystic fibrosis<sup>53</sup>. Given the fact that this microfluidic system is very versatile, it is anticipated that more and more systems will be distributed throughout the world.

Some specific protocol steps should be carefully considered. Media can be diluted to 50% with dH<sub>2</sub>O to help prevent bubbles and clogging but was not required in this case. The specific value of OD<sub>600</sub> used for seeding must be determined using trial runs of a growth experiment to see what works best for the particular set of conditions used. Bubbles in the wells prior to sealing can lead to bubbles in the microfluidic channels and should be removed by either popped or sucked out with a pipette tip. It is important to keep bacteria out of the small serpentine channels. By having equal volumes of media in the input and output during the seeding process, flow due to pressure from the liquid volume will be controlled so flow is only due to the applied pressure from the system. The calibration distances should be set up during installation by the company representative. These settings are specific for each camera.

There are several challenges that occur when finding the most representative threshold for an image. Setting maximum threshold values can be difficult if the average pixel intensity across regions of the background are not consistent caused by either selecting a stage position that is not at the center of the channel or from debris on the plate. Under the **MM Standard**, click on **Process**, and then select **Background and Shading Correction** tool to correction for these inconsistencies. However, this tool is generally only helpful if the user has taken images of the channels before seeding that they can use as reference images. Or, if background/shading reference images are unavailable, the user will need to use their judgment to set a threshold value that covers the most cell area without including background for the entire image. Alternatively, select representative areas to measure that exclude regions of inconsistency by click **Rectangular Region**, **Ellipse Region**, or **Trace Region** to select a region and select **Active Region** rather than **Entire Image** on the **Show Region Statistics** window (under **Analysis Tools**). If a representative region is utilized for thresholding the bright field image, the same region should be used for measurement of the corresponding FITC image. It is helpful to record the specific spatial statistics (**Left**, **Top**, **Width**, **Height**, **Area**, **Perimeter**) associated with that representative region so the same region will be found and measured on the corresponding FITC image.

To prevent a buildup of data on the hard drive that will cause the computer to slow down, an external hard drive can be purchased for data storage. Another option for data storage and facilitating data sharing is the CyVerse bioinformatics platform. Create an account on the CyVerse system by going to <http://www.cyverse.org/> Once logged in, launch the Discovery Environment then select "Log into CyVerse". Select "**Data**" and navigate to your the folder. If the image stack is on the local computer then select "**Upload**" then "**Simple Upload from Desktop**". Find the image stack file and select for upload. The file or a folder can be shared with collaborators if they have a CyVerse account and are granted permission. Sharing of the data folder to the general public requires that metadata be added for each file using CyVerse approved standards. This procedure will not be discussed here because this is not within the scope of this work.

**ACKNOWLEDGMENTS:**

This work was made possible by grants from the National Institute for General Medical Science (NIGMS) (5P20GM103427), a component of the National Institutes of Health (NIH)

**DISCLOSURES:**

The authors have nothing to disclose.

**REFERENCES:**

1. Donlan, R.M. *et al.* Model System for Growing and Quantifying *Streptococcus pneumoniae* Biofilms In Situ and in Real Time. *Applied and Environmental Microbiology*. **70** (8), 4980–4988, doi: 10.1128/AEM.70.8.4980-4988.2004 (2004).
2. Lührig, K. *et al.* Bacterial Community Analysis of Drinking Water Biofilms in Southern Sweden. *Microbes and Environments*. **30** (1), 99–107, doi: 10.1264/jsme2.ME14123 (2015).
3. Singh, R., Paul, D., Jain, R. Biofilms: implications in bioremediation. *Trends in Microbiology*. **14** (9), 389–397, doi: 10.1016/j.tim.2006.07.001 (2006).
4. Veerachamy. *et al.* Bacterial adherence and biofilm formation on medical implants: A review. *Proceedings of the Institution of Mechanical Engineers, Part H: Journal of Engineering in Medicine*. **228** (10), 1083–1099, doi: 10.1177/0954411914556137 (2014).
5. Lyczak, J.B., Cannon, C.L., Pier, G.B. Establishment of *Pseudomonas aeruginosa* infection: lessons from a versatile opportunist. *Microbes and Infection*. **2** (9), 1051–1060, doi: 10.1016/S1286-4579(00)01259-4 (2000).
6. Percival, S.L., Knapp, J.S., Edyvean, R., Wales, D.S. Biofilm Development On Stainless Steel In Mains Water. *Water Research*. **32** (1), 243–253, doi: 10.1016/S0043-1354(97)00132-2 (1998).
7. Percival, S.L., Knapp, J.S., Wales, D.S., Edyvean, R.G.J. The effect of turbulent flow and surface roughness on biofilm formation in drinking water. *Journal of Industrial Microbiology and Biotechnology*. **22** (3), 152–159, doi: 10.1038/sj.jim.2900622 (1999).
8. Ling, F., Whitaker, R., LeChevallier, M.W., Liu, W.-T. Drinking water microbiome assembly induced by water stagnation. *The ISME Journal*. **1**, doi: 10.1038/s41396-018-0101-5 (2018).
9. Moss, B. Water pollution by agriculture. *Philosophical Transactions of the Royal Society B: Biological Sciences*. **363** (1491), 659–666, doi: 10.1098/rstb.2007.2176 (2008).
10. Schultz, M.P., Bendick, J.A., Holm, E.R., Hertel, W.M. Economic impact of biofouling on a naval surface ship. *Biofouling*. **27** (1), 87–98, doi: 10.1080/08927014.2010.542809 (2011).
11. Tornero, V., Hanke, G. Chemical contaminants entering the marine environment from sea-based sources: A review with a focus on European seas. *Marine Pollution Bulletin*. **112** (1–2), 17–38, doi: 10.1016/j.marpolbul.2016.06.091 (2016).
12. LaBauve Annette E., Wargo Matthew J. Growth and Laboratory Maintenance of *Pseudomonas aeruginosa*. *Current Protocols in Microbiology*. **25** (1), 6E.1.1-6E.1.8, doi: 10.1002/9780471729259.mc06e01s25 (2012).
13. He, J. *et al.* The broad host range pathogen *Pseudomonas aeruginosa* strain PA14 carries two pathogenicity islands harboring plant and animal virulence genes. *Proceedings of the National Academy of Sciences of the United States of America*. **101** (8), 2530–2535, doi: 10.1073/pnas.0304622101 (2004).
14. Shaner, N.C., Steinbach, P.A., Tsien, R.Y. A guide to choosing fluorescent proteins. *Nature*

- Methods*. **2** (12), 905–909, doi: 10.1038/nmeth819 (2005).
15. Zhang, G., Gurtu, V., Kain, S.R. An Enhanced Green Fluorescent Protein Allows Sensitive Detection of Gene Transfer in Mammalian Cells. *Biochemical and Biophysical Research Communications*. **227** (3), 707–711, doi: 10.1006/bbrc.1996.1573 (1996).
16. Patterson, G.H., Knobel, S.M., Sharif, W.D., Kain, S.R., Piston, D.W. Use of the green fluorescent protein and its mutants in quantitative fluorescence microscopy. *Biophysical Journal*. **73** (5), 2782–2790, doi: 10.1016/S0006-3495(97)78307-3 (1997).
17. Chalfie, M., Tu, Y., Euskirchen, G., Ward, W.W., Prasher, D.C. Green fluorescent protein as a marker for gene expression. *Science*. **263** (5148), 802–805, doi: 10.1126/science.8303295 (1994).
18. Heim, R., Prasher, D., Tsien, R. Wavelength mutations and posttranslational autoxidation of green fluorescent protein. *Proceedings of the National Academy of Sciences of the United States of America*. **91**, 12501–12504 (1994).
19. Bakke, R., Kommedal, R., Kalvenes, S. Quantification of biofilm accumulation by an optical approach. *Journal of Microbiological Methods*. **44** (1), 13–26, doi: 10.1016/S0167-7012(00)00236-0 (2001).
20. Heydorn, A. *et al.* Statistical analysis of *Pseudomonas aeruginosa* biofilm development: impact of mutations in genes involved in twitching motility, cell-to-cell signaling, and stationary-phase sigma factor expression. *Applied and Environmental Microbiology*. **68** (4), 2008–2017, doi: 10.1128/AEM.68.4.2008-2017.2002 (2002).
21. Wilson, E. *et al.* Using Fluorescence Intensity of Enhanced Green Fluorescent Protein to Quantify *Pseudomonas aeruginosa*. *Chemosensors*. **6** (21), doi: doi:10.3390/chemosensors6020021.
22. Fang, H., Chen, Y., Huang, L., He, G. Analysis of biofilm bacterial communities under different shear stresses using size-fractionated sediment. *Scientific Reports*. **7** (1), 1299, doi: 10.1038/s41598-017-01446-4 (2017).
23. Saunders, K.A., Greenman, J. The formation of mixed culture biofilms of oral species along a gradient of shear stress. *Journal of Applied Microbiology*. **89** (4), 564–572, doi: 10.1046/j.1365-2672.2000.01148.x (2001).
24. Valquier-Flynn, H., Wilson, C.L., Holmes, A.E., Wentworth, C.D. Growth Rate of *Pseudomonas aeruginosa* Biofilms on Slippery Butyl Methacrylate-Co-Ethylene Dimethacrylate (BMA-EDMA), Glass and Polycarbonate Surfaces. *Journal of Biotechnology & Biomaterials*. **7** (4), doi: 10.4172/2155-952X.1000274 (2017).
25. Trulear, M.G., Characklis, W.G. Dynamics of biofilm processes. *Journal (Water Pollution Control Federation)*. **54** (9), 1288–1301, doi: 10.2307/25041684 (1982).
26. Machineni, L., Rajapantul, A., Nandamuri, V., Pawar, P.D. Influence of Nutrient Availability and Quorum Sensing on the Formation of Metabolically Inactive Microcolonies Within Structurally Heterogeneous Bacterial Biofilms: An Individual-Based 3D Cellular Automata Model. *Bulletin of Mathematical Biology*. **79** (3), 594–618, doi: 10.1007/s11538-017-0246-9 (2017).
27. Jiao, Y. *et al.* Identification of Biofilm Matrix-Associated Proteins from an Acid Mine Drainage Microbial Community. *Applied and Environmental Microbiology*. **77** (15), 5230–5237, doi: 10.1128/AEM.03005-10 (2011).
28. Flemming, H.C., Wingender, J. The Biofilm Matrix. *Nature Reviews Microbiology*. **8** (9), 632–


- 633, doi: 10.1038/nrmicro2415 (2010).
29. Yang, X., Beyenal, H., Harkin, G., Lewandowski, Z. Quantifying biofilm structure using image analysis. *Journal of Microbiological Methods*. **39** (2), 109–119, doi: 10.1016/S0167-7012(99)00097-4 (2000).
30. Heydorn, A. *et al.* Quantification of biofilm structures by the novel computer program COMSTAT. *Microbiology*. **146** (10), 2395–2407, doi: 10.1099/00221287-146-10-2395 (2000).
31. Azeredo, J. *et al.* Critical review on biofilm methods. *Critical Reviews in Microbiology*. **43** (3), 313–351, doi: 10.1080/1040841X.2016.1208146 (2017).
32. Goeres, D.M. *et al.* A method for growing a biofilm under low shear at the air liquid interface using the drip flow biofilm reactor. *Nature Protocols*. **4** (3), 783–788 (2009).
33. Goeres, D.M. *et al.* Statistical assessment of a laboratory method for growing biofilms. *Microbiology*. **151**, 757–762, doi: 10.1099/mic.0.27709-0 (2005).
34. Swartz, K., Stephenson, R., Hernandez, M., Jambang, N., Boles, B.R. The use of Drip Flow and Rotating Disk Reactors for Staphylococcus aureus biofilm analysis. *Journal of Visual Experiments*. **46**, e2470, doi: 10.3791/2470 (2010) (2010).
35. Werner, E. *et al.* Stratified Growth in Pseudomonas aeruginosa Biofilms. *Applied and Environmental Microbiology*. **70** (10), 6188–6196, doi: 10.1128/AEM.70.10.6188-6196.2004 (2004).
36. Wilson, C. *et al.* The Quantitative Assessment of Pseudomonas aeruginosa (PA)14 Biofilm Surface Coverage on Slippery Liquid Infused Polymer Surfaces (SLIPS). *International Journal of Nanotechnology in Medicine & Engineering*. **3** (3), 35–42 (2018).
37. Samarian, D.S., Jakubovics, N.S., Luo, T.L., Rickard, A.H. Use of a High-throughput In Vitro Microfluidic System to Develop Oral Multi-species Biofilms. *Journal of Visualized Experiments*. (94), 52467, doi: 10.3791/52467 (2014).
38. Rusconi, R., Garren, M., Stocker, R. Microfluidics Expanding the Frontiers of Microbial Ecology. *Annual Review of Biophysics*. **43**, 65–91, doi: 10.1146/annurev-biophys-051013-022916 (2014).
39. Kim J, Park H-D, Chung S. Microfluidic Approaches to Bacterial Biofilm Formation. *Molecules*. **17** (8), 9818–9834 (2012).
40. Mosier AP, Cady NC. Analysis of bacterial surface interactions using microfluidic systems. *Science Progress*. **94** (4), 431–50 (2011).
41. Nobile, C.J., Johnson, A.D. Candida albicans Biofilms and Human Disease. *Annual Review of Microbiology*. **69**, 71–92, doi: 10.1146/annurev-micro-091014-104330 (2015).
42. Chandra, J., Kuhn, D.M., Mukherjee, P.K., Hoyer, L.L., McCormick, T., Ghannoum, M.A. Biofilm Formation by the Fungal Pathogen Candida albicans: Development, Architecture, and Drug Resistance. *Journal of Bacteriology*. **183** (18), 5385–5394, doi: 10.1128/JB.183.18.5385-5394.2001 (2001).
43. Vasconcelos, M.A. *et al.* Effect of Algae and Plant Lectins on Planktonic Growth and Biofilm Formation in Clinically Relevant Bacteria and Yeasts. *BioMed Research International*. **2014**, 9 (2014).
44. Haralick, R.M., Shanmugam, K., Dinstein, I. Textural Features for Image Classification. *IEEE Transactions on Systems, Man, and Cybernetics*. **SMC-3** (6), 610–621, doi: 10.1109/TSMC.1973.4309314 (1973).
45. Wentworth, C.D. *Biofilm Morphology Suite*. at <<https://github.com/cdwentworth/Biofilm->

- Morphology-Suite.git>. (2018).
46. Cyverse. at <<https://user.cyverse.org/services/mine>>. (2018).
47. Benoit, M.R., Conant, C.G., Ionescu-Zanetti, C., Schwartz, M., Matin, A. New Device for High-Throughput Viability Screening of Flow Biofilms. *Applied and Environmental Microbiology*. **76** (13), 4136–4142, doi: 10.1128/AEM.03065-09 (2010).
48. Ding, A.M., Palmer, R.J., Cisar, J.O., Kolenbrander, P.E. Shear-Enhanced Oral Microbial Adhesion. *Applied and Environmental Microbiology*. **76** (4), 1294–1297, doi: 10.1128/AEM.02083-09 (2010).
49. Diaz De Rienzo, M.A., Stevenson, P.S., Marchant, R., Banat, I.M. Effect of biosurfactants on *Pseudomonas aeruginosa* and *Staphylococcus aureus* biofilms in a BioFlux channel. *Applied Microbiology and Biotechnology*. **100** (13), 5773–5779, doi: 10.1007/s00253-016-7310-5 (2016).
50. Moormeier, D.E., Bayles, K.W. *Staphylococcus aureus* Biofilm: A Complex Developmental Organism. *Molecular Microbiology*. **104** (3), 365–376, doi: 10.1111/mmi.13634 (2017).
51. Tremblay, Y.D.N., Vogeleer, P., Jacques, M., Harel, J. High-Throughput Microfluidic Method To Study Biofilm Formation and Host-Pathogen Interactions in Pathogenic *Escherichia coli*. *Applied and Environmental Microbiology*. **81** (8), 2827–2840, doi: 10.1128/AEM.04208-14 (2015).
52. Nobbs, A.H., Lamont, R.J., Jenkinson, H.F. *Streptococcus* Adherence and Colonization. *Microbiology and Molecular Biology Reviews*. **73** (3), 407–450, doi: 10.1128/MMBR.00014-09 (2009).
53. Díez-Aguilar, M., Morosini, M.I., Köksal, E., Oliver, A., Ekkelenkamp, M., Cantón, R. Use of Calgary and Microfluidic BioFlux Systems To Test the Activity of Fosfomycin and Tobramycin Alone and in Combination against Cystic Fibrosis *Pseudomonas aeruginosa* Biofilms. *Antimicrobial Agents and Chemotherapy*. **62** (1), e01650-17, doi: 10.1128/AAC.01650-17 (2018).



Figure 1

[Click here to access/download;Figure;Figure1.pdf](#)



1

2

3

4

5

6

7

8

9

10

11

12






Figure 3









Figure 4





Click here to access/download  
**Video or Animated Figure**  
Figure4\_edited.svg

Name of Material/ Equipment	Company	Catalog Number	Comments/Description
Ammonium Chloride, ACS	VWR	BDH9208-500G	Part of the minimal media composition
	Fluxion		
BioFlux 1000 48 Well Low Shear Plate	Biosciences	910-0047	
BioFlux 1000Z Microfluidic Imaging System	Fluxion		
	Biosciences	BF 1000Z	
Calcium Chloride Dihydrate, ACS Grade	VWR	97061-904	Part of the minimal media composition
Dextrose, Anhydrous, ACS	VWR	BDH9230-500G	Part of the minimal media composition
Magnesium Sulfate ACS Grade	VWR	EM-MX0070-1	Part of the minimal media composition
Potassium Phosphate Monobasic, ACS Grade	VWR	BDH9268-500G	Part of the minimal media composition
Pseudomonas Aeruginosa GFP	ATCC	15692GFP	Pseudomonas aeruginosa bacterial strain PA01 with GFP modification used for this study.
Sodium Chloride, ACS	VWR	BDH9286-500G	Part of the minimal media composition
Sodium Phosphate, Monobasic, Anhydrous, Reagent Grade	VWR	97061-942	Part of the minimal media composition



1 Alewife Center #200  
Cambridge, MA 02140  
tel. 617.945.9051  
[www.jove.com](http://www.jove.com)

## ARTICLE AND VIDEO LICENSE AGREEMENT

Title of Article:	Live Cell Analysis of Shear Stress on <i>Pseudomonas aeruginosa</i> using the Automated High-Throughput BioFlux 1000Z System
Author(s):	Arin L. Sutlief, Helena Valquier-Flynn, Christina Wilson, Marco Perez, Hunter, Kleinschmidt, Brett Schofield, Elizabeth Delmain, Andrea E. Holmes, Christopher D. Wentworth

Item 1: The Author elects to have the Materials be made available (as described at <http://www.jove.com/publish>) via:

☐ Standard Access ☒ Open Access

Item 2: Please select one of the following items:

- ☒ The Author is **NOT** a United States government employee.
- ☐ The Author is a United States government employee and the Materials were prepared in the course of his or her duties as a United States government employee.
- ☐ The Author is a United States government employee but the Materials were NOT prepared in the course of his or her duties as a United States government employee.

### ARTICLE AND VIDEO LICENSE AGREEMENT

1. **Defined Terms.** As used in this Article and Video License Agreement, the following terms shall have the following meanings: **"Agreement"** means this Article and Video License Agreement; **"Article"** means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein; **"Author"** means the author who is a signatory to this Agreement; **"Collective Work"** means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; **"CRC License"** means the Creative Commons Attribution-Non Commercial-No Derivs 3.0 Unported Agreement, the terms and conditions of which can be found at: <http://creativecommons.org/licenses/by-nc-nd/3.0/legalcode>; **"Derivative Work"** means a work based upon the Materials or upon the Materials and other pre-existing works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; **"Institution"** means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; **"JoVE"** means MyJoVE Corporation, a Massachusetts corporation and the publisher of The Journal of Visualized Experiments; **"Materials"** means the Article and / or the Video; **"Parties"** means the Author and JoVE; **"Video"** means any video(s) made by the Author, alone or in conjunction with any other parties, or by JoVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion

of the Article, and in which the Author may or may not appear.

2. **Background.** The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JoVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.

3. **Grant of Rights in Article.** In consideration of JoVE agreeing to publish the Article, the Author hereby grants to JoVE, subject to **Sections 4** and **7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the "Open Access" box has been checked in **Item 1** above, JoVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.

## ARTICLE AND VIDEO LICENSE AGREEMENT

4. **Retention of Rights in Article.** Notwithstanding the exclusive license granted to JoVE in **Section 3** above, the Author shall, with respect to the Article, retain the non-exclusive right to use all or part of the Article for the non-commercial purpose of giving lectures, presentations or teaching classes, and to post a copy of the Article on the Institution's website or the Author's personal website, in each case provided that a link to the Article on the JoVE website is provided and notice of JoVE's copyright in the Article is included. All non-copyright intellectual property rights in and to the Article, such as patent rights, shall remain with the Author.

5. **Grant of Rights in Video – Standard Access.** This **Section 5** applies if the "Standard Access" box has been checked in **Item 1** above or if no box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby acknowledges and agrees that, Subject to **Section 7** below, JoVE is and shall be the sole and exclusive owner of all rights of any nature, including, without limitation, all copyrights, in and to the Video. To the extent that, by law, the Author is deemed, now or at any time in the future, to have any rights of any nature in or to the Video, the Author hereby disclaims all such rights and transfers all such rights to JoVE.

6. **Grant of Rights in Video – Open Access.** This **Section 6** applies only if the "Open Access" box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby grants to JoVE, subject to **Section 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Video in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Video into other languages, create adaptations, summaries or extracts of the Video or other Derivative Works or Collective Works based on all or any portion of the Video and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. For any Video to which this **Section 6** is applicable, JoVE and the Author hereby grant to the public all such rights in the Video as provided in, but subject to all limitations and requirements set forth in, the CRC License.

7. **Government Employees.** If the Author is a United States government employee and the Article was prepared in the course of his or her duties as a United States government employee, as indicated in **Item 2** above, and any of the licenses or grants granted by the Author hereunder exceed the scope of the 17 U.S.C. 403, then the rights granted hereunder shall be limited to the maximum

rights permitted under such statute. In such case, all provisions contained herein that are not in conflict with such statute shall remain in full force and effect, and all provisions contained herein that do so conflict shall be deemed to be amended so as to provide to JoVE the maximum rights permissible within such statute.

8. **Protection of the Work.** The Author(s) authorize JoVE to take steps in the Author(s) name and on their behalf if JoVE believes some third party could be infringing or might infringe the copyright of either the Author's Article and/or Video.

9. **Likeness, Privacy, Personality.** The Author hereby grants JoVE the right to use the Author's name, voice, likeness, picture, photograph, image, biography and performance in any way, commercial or otherwise, in connection with the Materials and the sale, promotion and distribution thereof. The Author hereby waives any and all rights he or she may have, relating to his or her appearance in the Video or otherwise relating to the Materials, under all applicable privacy, likeness, personality or similar laws.

10. **Author Warranties.** The Author represents and warrants that the Article is original, that it has not been published, that the copyright interest is owned by the Author (or, if more than one author is listed at the beginning of this Agreement, by such authors collectively) and has not been assigned, licensed, or otherwise transferred to any other party. The Author represents and warrants that the author(s) listed at the top of this Agreement are the only authors of the Materials. If more than one author is listed at the top of this Agreement and if any such author has not entered into a separate Article and Video License Agreement with JoVE relating to the Materials, the Author represents and warrants that the Author has been authorized by each of the other such authors to execute this Agreement on his or her behalf and to bind him or her with respect to the terms of this Agreement as if each of them had been a party hereto as an Author. The Author warrants that the use, reproduction, distribution, public or private performance or display, and/or modification of all or any portion of the Materials does not and will not violate, infringe and/or misappropriate the patent, trademark, intellectual property or other rights of any third party. The Author represents and warrants that it has and will continue to comply with all government, institutional and other regulations, including, without limitation all institutional, laboratory, hospital, ethical, human and animal treatment, privacy, and all other rules, regulations, laws, procedures or guidelines, applicable to the Materials, and that all research involving human and animal subjects has been approved by the Author's relevant institutional review board.

11. **JoVE Discretion.** If the Author requests the assistance of JoVE in producing the Video in the Author's facility, the Author shall ensure that the presence of JoVE employees, agents or independent contractors is in accordance with the relevant regulations of the Author's institution. If more than one author is listed at the beginning of this Agreement, JoVE may, in its sole



## ARTICLE AND VIDEO LICENSE AGREEMENT

discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JoVE reserves the right, in its absolute and sole discretion and without giving any reason therefore, to accept or decline any work submitted to JoVE. JoVE and its employees, agents and independent contractors shall have full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JoVE has sole discretion as to the method of making and publishing the Materials, including, without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

12. **Indemnification.** The Author agrees to indemnify JoVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JoVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JoVE, making of videos by JoVE, or publication in JoVE or elsewhere by JoVE. The Author shall be responsible for, and shall hold JoVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to

the making of a video by JoVE its employees, agents or independent contractors. All sterilization, cleanliness or decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's expense. All indemnifications provided herein shall include JoVE's attorney's fees and costs related to said losses or damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JoVE, its employees, agents or independent contractors.

13. **Fees.** To cover the cost incurred for publication, JoVE must receive payment before production and publication the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JoVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.

14. **Transfer, Governing Law.** This Agreement may be assigned by JoVE and shall inure to the benefits of any of JoVE's successors and assignees. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to be one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement is required per submission.

### CORRESPONDING AUTHOR

Name:	Arin L. Sutlief	
Department:	Chemistry	
Institution:	Doane University	
Title:	Post Doctoral Research Fellow	
Signature:		Date: 08/08/2018

Please submit a **signed** and **dated** copy of this license by one of the following three methods:

1. Upload an electronic version on the JoVE submission site
2. Fax the document to +1.866.381.2236
3. Mail the document to JoVE / Attn: JoVE Editorial / 1 Alewife Center #200 / Cambridge, MA 02140



Dear Dr. Steindel,

We would like to thank you and the reviewers for the insightful comments and suggestions. We have implemented revisions, including the remaking and reordering of the existing images, review of the protocol, and we tried to address every single comment the reviewers pointed out. Please see our revisions based on the reviewer's comments as detailed below, with the reviewers comments in black text and our corrections shown in red. The manuscript has also been edited throughout to correct grammatical or spelling errors. All changes are highlighted in red font color within the revised manuscript.

Sincerely,

Arin Sutlief

Editorial comments:

1. Please take this opportunity to thoroughly proofread the manuscript to ensure that there are no spelling or grammar issues. The JoVE editor will not copy-edit your manuscript and any errors in the submitted revision may be present in the published version.

**The article was proofread to avoid any spelling or grammar issues.**

2. There are missing articles throughout.

**Missing articles were added throughout the manuscript where appropriate**

3. Please submit the figures as a vector image file to ensure high resolution throughout production: (.svg, .eps, .ai). If submitting as a .tif or .psd, please ensure that the image is 1920 pixels x 1080 pixels or 300dpi.

**Files have been converted to .svg files.**

4. Figures 3 and 5 can be made supplemental. Please note that there is a video here.

**These were revised as described.**

5. Figures 6-9 can be combined into one figure with 4 panels.

**Images were reformatted as described.**

6. Please revise the table of the essential supplies, reagents, and equipment. The table should include the name, company, and catalog number of all relevant materials in separate columns in an xls/xlsx file.

**The table of essential supplies, reagents, and equipment was revised as described.**

7. JoVE cannot publish manuscripts containing commercial language. This includes trademark symbols (™), registered symbols (®), and company names before an instrument or reagent. Please remove all commercial language from your manuscript and use generic terms instead. All commercial products should be sufficiently referenced in the Table of Materials and Reagents.

For example: BioFlux Montage, etc.

**The single mention of the company Fluxion Biosciences was removed. BioFlux is not the company and BioFlux Montage was used in previously published articles about the BioFlux including in JoVE.**

8. Please ensure that all text in the protocol section is written in the imperative tense as if telling someone how to do the technique (e.g., "Do this," "Ensure that," etc.). The actions should be described in the imperative tense in complete sentences wherever possible. Avoid usage of phrases such as "could be," "should be," and "would be" throughout the Protocol. Any text that cannot be written in the imperative tense may be added as a "Note." However, notes should be concise and used sparingly. Please include all safety procedures and use of hoods, etc.

**The protocol section was revised to contain only the imperative mood, except for a small number of "notes".**

9. The Protocol should contain only action items that direct the reader to do something. Please move the discussion about the protocol to the Discussion.

The protocol has been reviewed to reword or remove all non-action items. The discussion has been modified to include items of priority removed from the protocol.

10. Please obtain explicit copyright permission to reuse any figures from a previous publication. Explicit permission can be expressed in the form of a letter from the editor or a link to the editorial policy that allows re-prints. Please upload this information as a .doc or .docx file to your Editorial Manager account. The Figure must be cited appropriately in the Figure Legend, i.e. "This figure has been modified from [citation]."

All figures were originally made by the authors.

Reviewers' comments:

Reviewer #1:

Manuscript Summary:

The authors describe the use of the commercially available microfluidics platform from Fluxion, the Bioflux Z1000 system. They describe the use of this system for analysis of PAO1 biofilm formation. The described protocol is reasonable complete and should be useful for readers in the present state. The cited literature is not very complete in my eyes as there have been more studies describing the use of this system that are not mentioned. e.g., this system has been used with *Candida albicans*.

Line 59: I suggest changing "...the method..." to "...a method..."

Revised as described.

Line 61 and general: I have some doubts whether the Bioflux system is truly a high throughput system. I suggest the term medium throughput.

Taking into consideration this comment, all mentions of the throughput of the BioFlux system was described as a "higher throughput".

Line 69. There is no doubt that microbial biofilms are present in all natural streams, rivers and lakes. Estimates suggest that near 80% of all microbial life resides in biofilms. Biofilms are therefore not an exceptional community, just exceptionally understudied...

We are unsure of the meaning of this comment since there is no mention of an exceptional community in the article but this section was reviewed and reworded.

Line 76: I object to the abbreviation of *Pseudomonas aeruginosa* as PA. There is an agreement on the use of species names for a long time already. Please use *P. aeruginosa* as the standard abbreviation throughout (with exception of first mention).

Throughout the article, this was revised as described.

Line 83: the term prototype strain is strange. Please use reference strain.

Revised as described.

Line 86: ...This type of fluorescence analysis is advantageous for the study of biofilms because GFP does not interfere with cell growth and function... This is debatable as all heterologous expressed proteins, even presence of "extra" DNA poses a burden on the cell. The protein has to be synthesized at the expense of energy. I agree that the effects are generally marginal and probably not relevant for the described experiments, yet the above statement is scientifically incorrect.

This was revised to "...GFP does not interfere significantly with cell growth and function..."

Line 102: Please replace...film... with ...biofilm...

Revised as described.

Line 121: Biofilms in nature and medicine are notoriously heterogeneous. It is therefore relevant

to include this heterogeneity in model systems, even though it does make the job harder...

This phrase was reworded for clarity "Biofilms do not grow evenly across the surface in a drip flow reactor because the low shear of media causes trailing along the larger conglomerates of *P. aeruginosa* bacteria therefore the biofilm growth is not very smooth and the uneven samples cannot be easily analyzed."

Line 125: A typical biofilm experiment in the Bioflux requires at least 0.5-1 ml medium

This was revised to read as "...where only milliliters of media are required..."

Line 130: Please replace fluorescence spectroscopy with microscopy. Also, do you measure GFP or visualize. I think this is an important drawback of microscopic methods, you can't really get objective quantification of GFP fluorescence. It should be made clear that the use of the bioflux method allows for semi-quantitative measurements at best.

This section was reworded to "It is coupled with fluorescence microscopy to visualize the fluorescence of the EGFP tag in PA01..."

Line 138: Important to mention is that temperature will be constant over the entire plate and not a variable within an experiment. Also, shear and thus flow is managed per column limiting the variation to 4 different values within a single experiment.

The section was reworded to the following to enhance its description: "This allows for 24 simultaneous experiments that can be performed with various experimental conditions, such as bacterial strains, antimicrobial treatments, and media varied from channel to channel and controlled shear flow for each column of six channels. The experimental temperature is also controlled with one temperature setting throughout the plate."

Line 147: Please cite literature on the use of the bioflux for studying *Candida albicans* and other species. I would not state...It is hoped that... However it is the aim of this protocol/report/manuscript...

This section was reworded and three references were added: This report will allow for the expansion of experiments involving other microorganism besides bacteria, such as fungi and algae that have medical and environmental applications.<sup>41-43</sup>

Line 155: I disagree with this statement. We always use partial plates without any problems. It is better to use partial plates in light of observed slight differences in flow within a column depending on the position of the well (near air connection on top or far away from it).

This statement has been removed in order to stay neutral to the company's claims.

Line 176: Does this strain require Antibiotic pressure to maintain eGFP? If so, please mention this.

No chemicals or antibiotic is required for the expression of eGFP and a reduction in eGFP signal has not been observed within the strain.

Line 186: The order of PC, Microscope and Bioflux is not really important, if the CCD is not on when the software is started, it will mention this.

In order to make this a more neutral statement the wording was adjusted to "In order to avoid an error in the connection of the instrument to the software, it is advised that the order of turning on the instrument is performed in the following order:...."

Section 4.2 and further: airbubbles will always form, Our experience is to always make sure to remove any bubbles prior to sealing the plate. Maybe this is important to mention, also how to do that...

To address this we added the following wording "Bubbles in the wells prior to sealing can lead to bubbles in the microfluidic channels and should be removed by either popped or sucked out with a pipette tip."

Also, we always prime from the input well. This is not really an issue with the 48-well plate, however when using the 2 inlet system (24-well plate) this is critical to prevent air to be caught in the system.

It was recommended by Fluxion by representatives and the manual to prime from the output well. In this case, the procedure described is a system for a 48-well plate system.

Specialty plates such as the 24-well plate may require modification to the set-up process. The following sentence at the end of introduction section was modified to make sure it was clear we were describing a method for a 48-well plate system: "The detailed approach describes how to culture PA01-EGFP, inoculate a 48-well plate and set up the microfluidic device and software, set up the fluorescent microscope, and demonstrate the software analysis to obtain the biofilm coverage, growth rate, and morphological properties such as surface roughness."

Section 4.5: it is critical to monitor this by eye due to slight differences throughout the length of the column. It is better to not give a time for the flow as this might vary depending on temp. medium and organism used. Also, maybe it is good to mention that all channels need to be checked individually to ensure all of them are inoculated. The microscope only visualizes 1 or 2 channels.

Given that this is a paper describing a protocol for a specific experimental set-up with a single experimental temperature and organism, a time frame would be appropriate for the protocol to act as a frame of reference. The phrase (below) was reworded to reflect that the live camera feed should be used but a note (below) was included to caution that the time frame would be optimal for our conditions only and across the plate there may be variation which would require additional flow time for certain columns of channels.

"On the Montage Control Module, focus on a single channel using the live camera feed after placing the plate stage onto the microscope stage. Visually monitoring by the live feed, resume flow at 1.00-2.00 dyne/cm<sup>2</sup> for approximately two to four s to allow cells to enter the experimental channel but not in the serpentine channels."

"Note: The amount of time needed for seeding will vary with media and organism so it should be monitored closely by the live feed until optimized and used as a general time frame. Seeding throughout the plate may vary which would require more time of applied flow to certain columns of channels for complete seeding."

Section 6.1: We were advised never to use more than 1 mL of liquid per well...

Another group that regularly use the BioFlux advised us that 1.3 mL maximum could be put into the well as long as certain cautions were made such as the closer the media is to the cover the more slight movements could cause bacteria to come into contact with the cover. This would contaminant the cover which could be significant in the case of organisms not susceptible to ethanol killing. The amount of volume can be adjusted based on desired experimental time/flow rate but would be advised to keep it a consistent variable.

Reviewer #2:

Manuscript Summary:

This manuscript provides a clear procedure for the effective use of the BioFlux device to study *P. aeruginosa* PAO1 biofilms in more detail and in new ways that hasn't to my knowledge been described previously. This procedure will be extremely helpful to other researchers interested in using the BioFlux microfluidic device to study biofilms produced by *P. aeruginosa* as well as other bacterial species. The use of the enhanced GFP is very important as this eliminates the need for another fluorescent stain, such as propidium iodide (which is not that effective in determining dead cells in *P. aeruginosa*) to distinguish between live and dead cells in the biofilm. Also, the methodology for determining textural entropy and roughness coefficient are fascinating as is the concept that the amount of shear flow applied affects the formation of the biofilms. This will help other researchers maximize or optimize the conditions they may wish to use depending on what they are testing. Finally, the procedure described is detailed and clear and the figures help enhance the procedures and the results that are obtained.

Major Concerns:

Lines 180-182 - I found the OD used to be quite high, but maybe this is an optimal OD600

based on the number of hours the biofilms are allowed to develop. We had problems in the past and had to use PAO1 at a lower OD600 (around 0.3) or we would get over-development by 8 hours. Also I am wondering about holding the cells at 37 degrees Celcius, I am assuming this would keep them warm as I know we have to warm our media and such before adding to the device, but do you think the OD600 will increase if the cells are left at 37 too long before adding to the plate?

The following was added in the discussion under "protocol considerations" to clarify the choice of OD600 used for seeding:

"The specific value of OD<sub>600</sub> used for seeding must be determined using trial runs of a growth experiment to see what works best for the particular set of conditions used."

The following sentence was added to the protocol to mitigate the concern about holding the cells at 37 [C] while waiting for the seeding procedure:

"The OD<sub>600</sub> should be checked again immediately before seeding to insure that it has not changed significantly from the target OD<sub>600</sub>, 0.8, in this case."

Figure 2 - the images for the biofilms could be a bit more resolved. I understand it is hard to see individual cells normally, but a bit more definition of the cells would be best. Also it would be nice to see some information on how many cells are living versus dead on the fluorescent image, if this can be distinguished. Finally, what objective was used on the microscope, are these images taken from a 10X lens?

The image has been remade with biofilm images that are a little clearer.

The GFP marker is known to dissipate upon the death of GFP modified *P. aeruginosa* and therefore, acts as a method to disguise live versus dead cells.

Step 5.3.1 states that the objective used for the experiment is a 20x lens. Wording was added to Figure 2 caption further clarification. "Shown are bright field (top) and fluorescent (bottom) microscopy images of a single channel with a PA01-EGFP biofilm using a 20x objective."

Minor Concerns:

Lines 184-306 - why are these sections (3-6) highlighted in yellow?

Sections of the manuscript are highlighted as instructed by JoVE requirements to indicate the sections the authors find important to address in the complimentary video.

Lines 647-648 - should this read biofilm characteristics?

The indicated lines does not include any biofilm characteristics so we are not sure what is being addressed with this comment. We reviewed the section to catch any possible errors.

Reviewer #3:

Manuscript Summary:

The manuscript entitled "Live Cell Analysis of Shear Stress on *Pseudomonas aeruginosa* using an Automated High-Throughput Microfluidic System" by Arin Sutlief and colleagues provides a comprehensive guidance to study microbial biofilms by time-lapse microscopy using the BioFlux 1000Z microfluidic system. The manuscript is clearly written and provides all necessary technical details required for successful and reproducible application of the described methodology. The subject described in the manuscript fits well within the scope of the journal and the methodology will undoubtedly be interesting to a broad audience. I have only a few minor comments:

-Lines 65-66 - This sentence requires rephrasing as it is a bit unclear in the present form. Specifically:

"...be beneficial to the bacteria and their environment...." - what kind of environments the authors are referring to?

"....they can also be highly harmful with undesirable consequences" - harmful to what?

This section has been reworded and references were added for clarity. "These biofilm communities can be beneficial to the environment, such as improving water quality in water supply lines and in bioremediation of recalcitrant compounds.<sup>2, 3</sup> However, biofilms can also be highly harmful to human health with undesirable consequences. For example, medical devices, such as hip and knee implants, are one type of surface where biofilm accumulation has been a challenge and causes severe medical complications.<sup>4, 5</sup>

-Lines 77-78 - "...community-acquired and hospital-acquired..." - can be changed to ... community- and hospital-acquired...


Revised as described.

-Line 84 - "EGFP is a mutant of GFP..." - can be changed to -... EGFP represents a mutant form of GFP...

Revised as described.

-Line 93 - Can be rephrased to -.... Biofilms can grow under various environmental conditions including those with different flow rates.






File Edit Devices BioFlux Workflows Data Log Journal Analysis Tools MM Standard Help

Move Stage to Absolute Position Sample Reload Adjustment Control Panel Multi Dimensional Acquisition

Threshold Image Integrated Morphometry Analysis

Locator Rectangular Region Ellipse Region Trace Region Single Line Multi-Line Traced Line Auto-Trace Region

Icons: [File Explorer] [Save] [Print] [Zoom In] [Zoom Out] [Fit] [Reset] [Undo] [Redo] [Help]



### Review Multi Dimensional Data

Select Base File... H:\Wentworth2018\03302018\033018.ND Open Sequential...

Wavelengths:

- ☒ FITC 100% ...
- ☐ Brightfield 5...
- ☐ All Closed

Stage Position: 15, 1 Time [Navigation Icons] Timepoints: 1 to 289

☐ Enable Montage # of Thumbnails: 1 x 1 [Navigation Icons] Apply

Selections [X's] Display Event marks

Full Image Active Region

☐ Color Composite Red: <None> Green: <None> Blue: <None> Gray: <None>

Reset Image Displays Run Journal Loop... Close