

Journal of Visualized Experiments

High-throughput and Comprehensive Drug Surveillance Using Multisegment Injection-capillary Electrophoresis Mass Spectrometry

--Manuscript Draft--

Article Type:	Invited Methods Article - JoVE Produced Video
Manuscript Number:	JoVE58986R2
Full Title:	High-throughput and Comprehensive Drug Surveillance Using Multisegment Injection-capillary Electrophoresis Mass Spectrometry
Keywords:	Drug testing, drug surveillance, drugs of abuse, drug metabolites, urine, capillary electrophoresis, mass spectrometry, clinical depression, chronic pain, high throughput screening
Corresponding Author:	Philip Britz-McKibbin McMaster University Hamilton, ON CANADA
Corresponding Author's Institution:	McMaster University
Corresponding Author E-Mail:	britz@mcmaster.ca
Order of Authors:	Philip Britz-McKibbin Meera Shanmuganathan Sabrina Macklai Clara Barrenas Cárdenas Zachary Kroezen William Zizek Howard Lee Marcus Kim
Additional Information:	
Question	Response
Please indicate whether this article will be Standard Access or Open Access.	Open Access (US\$4,200)
Please indicate the city, state/province, and country where this article will be filmed . Please do not use abbreviations.	Hamilton, ON, Canada

CHEMISTRY and CHEMICAL BIOLOGY



27 September, 2018

Dr. Phillip Steindel
Review Editor, JoVE

Re: Manuscript Submission: "High Throughput and Comprehensive Drug Surveillance using Multisegment Injection-Capillary Electrophoresis-Mass Spectrometry (JoVE58986)"

Dear Dr. Steindel,

We wish to re-submit our revised manuscript to JoVE as an article (JoVE58986) that describes an innovative method for *high throughput and comprehensive surveillance of drugs of abuse* (DoA) without complicated sample handling when using multisegment injection-capillary electrophoresis (MSI-CE-MS). This is urgently needed given a *global opioid drug crisis and its impact on public health*, including an alarming number of overdose deaths attributed to synthetic opioids and other addictive psychoactive medications. We have included a highlighted version of our manuscript indicating all changes made, as well as a detailed response to reviewers/editors critiques in our submission.

Currently, a two-tiered strategy is used for urine drug testing based on a primary screen using immunoassays followed by confirmatory testing by more specific yet lower throughput GC/MS and increasingly LC-MS/MS methods. However, this screening algorithm is prone to both false positives and false negatives due to the cross-reactivity and poor sensitivity of antibody reagents for many classes of DoA in human urine. Additionally, conventional GC-MS and LC-MS methods have limited sample throughput for resolving a chemically diverse range of DoA and their metabolites due to times required for elution and column conditioning. Herein we report an improved method for drug surveillance using MSI-CE-MS that is an extension to our recently published work (*Anal. Chem.* **2017**, 89,11853). In this contribution, we focus on describing a **detailed protocol** for routine drug screening by MSI-CE-MS using a **modified methodology** to improve resolution of complex drug mixtures (isobars/isomers) in authentic urine samples from a cohort of **clinically depressed patients**. We also demonstrate unambiguous identification and quantification of diverse classes of DoA as required for confirming patient drug adherence and potential drug misuse/substitution based on their known prescription.

In summary, we anticipate that our manuscript will appeal to a broad audience while offering a powerful new solution for systematic drug surveillance for an alarming array of illicit drugs, prescription medications and over-the-counter drugs prone to abuse/misuse and inadvertent overdose. We hope that publishing in JoVE will allow for better translation of our methodology to other research laboratories in academia, as well as commercial forensic toxicology/therapeutic drug monitoring laboratories.

Sincerely,

Dr. Philip Britz-McKibbin
Department of Chemistry and Chemical Biology
McMaster University
Hamilton, ON, Canada

TITLE:

High-throughput and Comprehensive Drug Surveillance Using Multisegment Injection-capillary Electrophoresis Mass Spectrometry

AUTHORS & AFFILIATIONS:

Meera Shanmuganathan¹, Sabrina Macklai¹, Clara Barrenas Cárdenas¹, Zachary Kroezen¹, Marcus Kim², William Zizek³, Howard Lee³, Philip Britz-McKibbin¹

¹Department of Chemistry and Chemical Biology, McMaster University, Hamilton, ON, Canada

²Agilent Technologies Inc., Mississauga, ON, Canada

³Seroclinix Corporation, Mississauga, ON, Canada

Corresponding Author:

Philip Britz-McKibbin (britz@mcmaster.ca)

Email Addresses of Co-authors:

Meera Shanmuganathan (shanmm2@mcmaster.ca)

Sabrina Macklai (macklais@mcmaster.ca)

Clara Barrenas Cárdenas (c.barrena1@usp.ceu.es)

Marcus Kim (marcus.kim@agilent.com)

Zachary Kroezen (kroezep@mcmaster.ca)

Howard Lee (hlee@seroclinix.com)

William Zizek (bzizek@seroclinix.com)

KEYWORDS:

Drug testing, drug surveillance, drugs of abuse, drug metabolites, urine, capillary electrophoresis, mass spectrometry, clinical depression, chronic pain, high-throughput screening

SUMMARY:

Here we describe a high-throughput method for comprehensive drug surveillance that allows for the improved resolution and detection of large panels of drugs of abuse and their metabolites, with quality control based on multisegment injection-capillary electrophoresis-mass spectrometry.

ABSTRACT:

New analytical methods are urgently needed to enable high-throughput, yet comprehensive drug screening, given an alarming opioid and prescription drug crisis in public health. Conventional urine drug testing based on a two-tier immunoassay screen followed by a gas chromatography-tandem mass spectrometry (GC-MS/MS) or liquid chromatography-tandem mass spectrometry (LC-MS/MS) method are expensive and prone to bias while being limited to targeted panels of known drugs of abuse (DoA). Herein, we outline an improved method for drug surveillance that allows for the resolution and detection of an expanded panel of DoA and their metabolites when using multisegment injection-capillary electrophoresis-mass spectrometry (MSI-CE-MS). Multiplexed separations of ten urine samples with a quality control by CE (< 3 min/sample) in

conjunction with full-scan data acquisition using a time-of-flight mass spectrometer (TOF-MS) under positive ion mode detection allows for the identification and quantification of DoA above recommended cut-off levels. An excellent resolution of drug isomers and isobars, including background interferences, are achieved when using MSI-CE-MS with an electrokinetic spacer between sample segments, where accurate mass/molecular formula together with the comigration of a matching deuterated internal standard and the detection of one or more biotransformed metabolites facilitate DoA identification over a wider detection window. Additionally, urine samples can be analyzed directly without enzyme deconjugation for the rapid screening without complicated sample workup. MSI-CE-MS enables the surveillance of a broad spectrum of DoA that is required for the treatment monitoring of high-risk patients, including confirming prescribed drug adherence, revealing illicit drug use/substitution, and evaluating optimal dosage regimes as required for new advances in precision medicine.

INTRODUCTION:

An alarming increase in the misuse of and the addiction to opioids for the management of chronic pain represents a growing menace to public health, with over 70,000 drug overdose deaths in the USA estimated in 2017¹. Similarly, various other psychotropic medications are also widely prescribed to children and young adults for the treatment of anxiety, depression, and mental health issues². As a result, urine drug testing methods, widely developed for the workplace and forensic toxicology, is of vital importance to the therapeutic monitoring of prescribed medications that are prone to tolerance and dependence^{3,4}. This is needed to ensure adherence, optimal therapeutic efficacy, and patient safety, while revealing potential substitution, including illicit or nonprescribed drug misuse. Currently, urine drug tests for DoA rely on a two-tiered approach, comprising of an initial competitive immunoassay screen *via* point-of-care devices or laboratory analyzers, followed by a confirmatory test with greater specificity based on GC-MS/MS and, increasingly, LC-MS/MS⁵. However, immunoassays are prone to the bias, as antibody reagents bind nonspecifically to various drug classes to generate a presumptive screen-positive result, which prevents a reliable identification and quantification of specific drugs or complex drug mixtures⁶. In this context, more accurate urine drug tests are urgently needed given the exorbitant costs of comprehensive polydrug screens,⁷ including designer drugs and synthetic urine products that elude conventional targeted assays.

High-efficiency separations in conjunction with high-resolution MS (HRMS) using TOF or orbitrap mass analyzers have been proposed as a nontargeted strategy for drug surveillance in an era of polypharmacy and expanding panels of DoA^{8,9}. However, conventional LC separations are slow (> 15 min) due to the long elution times for the resolution of chemically diverse classes of DoA and their metabolites using gradient elution programs, which limits the sample throughput for routine drug screening. Alternatively, direct analysis methods based on the desorption ionization (DESI)¹⁰ and laser diode thermal desorption (LDTD) allow for faster drug screening without separation¹¹. However, these ambient ionization methods are prone to isobaric/isomeric interferences when analyzing complex urine specimens, thus requiring independent confirmatory testing.

Our laboratory has recently developed a multiplexed separation platform to increase sample

throughput while retaining the resolution and data fidelity of a high-efficiency separation based on MSI-CE-MS¹². Novel data workflows can be designed in MSI-CE-MS for the nontargeted metabolite profiling (*i.e.*, metabolomics) of volume-restricted biospecimens, whereas quality controls (QC) allow for the batch correction as required for large-scale population studies¹³. In this case, separations are performed using an isocratic buffer with solute ionization occurring under steady-state solvent conditions when using a conventional sheath liquid interface, which allows for serial injections of 10 or more samples within a single run for the rapid yet selective screening of diverse classes of DoA and their metabolites¹⁴. The focus of this study is to further validate MSI-CE-MS for the direct analysis of authentic urine samples from a representative cohort of clinically depressed patients using a “dilute-and-shoot” method that avoids the need for enzyme hydrolysis¹⁵. Additionally, the implementation of an electrokinetic spacer between serial sample plugs in MSI-CE-MS¹⁶ was performed to further improve the resolution of various drug isomers/isobars, background urinary interferences, and/or cross-interferences when screening large panels of DoA. The absolute quantification of DoA in urine specimens when using matching deuterated internal standards (d-IS) is demonstrated when using MSI-CE-MS. This approach also facilitates drug identification, as well as deducing the correct sample position of screen-positive cases when compared to a drug panel mixture at the recommended screening cut-off level that also functions as an internal reference/QC within the same run.

PROTOCOL:

Blinded urine samples were kindly provided by Dr. Zainab Samaan from the Moods Disorder Clinic at St. Joseph’s Hospital (Hamilton, ON, Canada), whose study was approved by the Hamilton Integrated Research Ethics Board.

1. Reagents and Sample Solutions Preparation

1.1. Background electrolyte preparation

1.1.1. Prepare 50 mL of background electrolyte (BGE) biweekly, 1 M formic acid, pH 1.8, with 15% v/v acetonitrile as an organic modifier.

1.1.2. Aliquot 1.9 mL of concentrated formic acid stock (26.5 M) into a 50 mL volumetric flask and add 7.5 mL of acetonitrile and 40.6 mL of deionized water to bring the total volume to 50 mL. Then, sonicate the solution for 15 min and transfer it to a sterilized bottle with tight sealing.

1.2. Sheath liquid preparation

1.2.1. Prepare 200 mL of sheath liquid solution weekly, comprising 0.1% formic acid in 60:40 (MeOH:H₂O).

1.2.2. Aliquot 200 µL of formic acid (26.5 M stock) into a bottle containing 120 mL of MeOH and 80 mL of H₂O to make sheath liquid and degas it for 15 min.

1.2.3. Next, add 10 μL of each of the following reference ions, purine and hexakis(2,2,3,3-tetrafluoropropoxy)phosphazine (HP-921), into the sheath liquid to provide constant mass signals at m/z 121.05087 and m/z 922.00979, respectively.

NOTE: Those reference ions allow for real-time mass correction and also monitor for potential matrix-induced ion suppression or enhancement effects during separation.

1.3. Standard preparation

1.3.1. Prepare a standard solution containing an 84-drug panel mixture at the 3x clinical screening cut-off level (L6)¹⁴ in 1 mL of methanol, using the drug standards purchased from a chemical supplier. Prepare separately a mixture of 48 matching deuterated internal drug standard (d-ISs) at a 10-fold higher concentration than the 84-drug mixture in 1 mL of methanol. In addition, prepare a stock solution containing 200 μM of 4-fluoro-*L*-phenylalanine (F-Phe) and 3-chloro-*L*-tyrosine (Cl-Tyr) in 1 mL of deionized H_2O .

1.3.2. Perform a fivefold dilution of the above-mentioned 84-drug mixture (20 μL) and the 48 d-ISs (20 μL), as well as a tenfold dilution of 4-fluorophenylalanine (F-Phe, 10 μL) and 3-chlorotyrosine (Cl-Tyr, 10 μL) in a certified synthetic urine matrix (40 μL) to make a total volume of 100 μL in a 250 μL centrifuge tube.

1.4. External calibration curve preparation

1.4.1. Prepare five-point external calibration curves (as described below) in quadruplicate ($n = 4$) over a 20-fold dynamic range using the 84-drug standard mixture and 48 of their corresponding d-ISs made in step 1.3.1. For example, for making a five-point external calibration curve, ensure diluting the 84-drug mixture (L6, step 1.3.1) to 2-, 4-, 10-, 20-, and 40-fold, whereas it is necessary to prepare a fivefold dilution of matching d-ISs (20 μL) and a tenfold dilution of 4-F-Phe (10 μL) and Cl-Tyr (10 μL) as an additional IS in a synthetic urine matrix to make a total volume of 100 μL . In cases when a matching d-IS is not available for a specific drug, use 4-F-Phe as a surrogate IS for data normalization.

NOTE: Avoid the inclusion of additional d-ISs that cross-interfere (isobaric) with other DoA within the panel if they are not fully resolved by the CE separation.

1.5. Urine sample preparation

1.5.1. Thaw deidentified morning urine samples from a representative cohort of ten clinically depressed patients with a known prescription drug history. Store the urine samples at -80°C after collection until they have to be thawed for analysis.

NOTE: Avoid multiple freeze-thaw cycles of urine or a delayed storage at room temperature following the collection due to their effect on the chemical stability of certain DoA metabolites and their conjugates.

175
176 1.5.2. Vortex the deidentified morning urine samples for 30 s and, then, centrifuge them for 1
177 min at 14,000 x *g* for sedimentation. After that, aliquot 10 µL of the above-mentioned processed
178 urine, 10 µL of d-IS, and 5 µL of F-Phe/Cl-Tyr and 25 µL of deionized H₂O, and vortex for 1 min
179 (repeat these steps in triplicate to assess technical precision). Transfer a 20 µL aliquot of the
180 above-described mixture to a polypropylene vial for analysis.

181 182 **2. Setup of the CE-TOF-MS System**

183 184 **2.1. Uncoated fused-silica capillary conditioning parameters**

185
186 2.1.1. Make sure that all standards, external calibration curves, and urine samples prepared in
187 sections 1.3 - 1.5 are separated using an uncoated open-tubular fused silica polyimide-coated
188 capillary with an internal diameter of 50 µm, an outer diameter of 360 µm, and a total capillary
189 length of 120 cm.

190
191 NOTE: A diamond cutter tool is used to ensure consistent capillaries that are further inspected to
192 ensure a flush and smooth cut.

193
194 2.1.2. Remove about 7 mm of the polyimide coating from both the distal capillary ends, using a
195 capillary window maker to reduce sample carry-over and to prevent potential polyimide swelling
196 when in contact with organic solvent.

197
198 2.1.3. Make sure to protrude about 2 mm of the capillary outlet out of the CE sprayer needle and
199 install the capillary inlet in the CE cartridge by looping two 360° turns. Then, carefully install the
200 CE cartridge in the CE and place the CE sprayer out of the ion source when conditioning the
201 capillary.

202
203 NOTE: Cut the capillary outlet consistently as it is crucial for ensuring a stable electrospray current
204 when coupled to MS. Two different sprayers were used in this work: a CE sprayer was used for
205 the sample analysis, and an LC sprayer was used for the mass calibration.

206
207 2.1.4. Make sure to perform the capillary conditioning with the CE sprayer not placed in the ion
208 source. Put the LC sprayer in the ion source.

209
210 NOTE: This allows for mass calibration while avoiding any contamination of the coaxial sheath
211 liquid interface that is subsequently coupled to the TOF-MS.

212
213 2.1.5. Condition new capillaries by selecting the **flush function** (940 mbar) on the vendor's
214 software used to control the instrument and set the duration of the flush at 30 min each for four
215 different solvents in the following order: methanol, 1M NaOH, deionized water, and BGE. Make
216 sure to put the isocratic pump at standby during the capillary conditioning period.

2.1.6. Wipe the CE-MS sprayer with a tissue wipe soaked in MeOH/H₂O to remove any residual salt deposits.

2.1.7. Perform the daily preventative maintenance of the CE-MS system before analyzing the samples. Use methanol to wipe down of CE electrode (at inlet) but use an isopropanol and water (50:50) mixture to clean the CE-MS interface, which is important to avoid salt build-up and limit potential sample carry-over.

2.1.8. Remove the LC sprayer out of the ion source and position the cleaned sprayer with a newly conditioned capillary in the coaxial sheath liquid interface for CE-MS.

2.1.9. Turn the isocratic pump on and apply a voltage of 30 kV for 15 min to ensure that a stable CE current profile is achieved prior to the urine analysis.

2.2. Injection and separation conditions using CE

2.2.1. Open the vendor software used for controlling the instrument to set the CE parameters and select the **preconditioning step**. Set **flush function** to 600 s and specify a BGE vial position.

NOTE: An external water chiller may need to be connected to CE systems that lack sample tray cooling features, as this is required to prevent sample evaporation when analyzing large batches of volume-restricted samples.

2.2.2. Set the **injection function** to hydrodynamically inject the samples at 100 mbar for 5 s and electrokinetically inject BGE spacers at 30 kV for 75 s.

NOTE: This process is repeated with each sample injection (from a different vial position) followed by a BGE spacer (from the same vial position) that is performed automatically within a method program on the software for a total of 11 discrete samples analyzed in the same run when using an MSI format. Analyze 10 urine samples in random order after an 84-panel drug mixture at the screening cut-off level is injected as the first sample position for each MSI-CE-MS run.

2.2.3. Set the **applied voltage** to be 30 kV, the **cartridge temperature** to be 25 °C, and the **total runtime** to be 45 min. Using the timetable, apply a pressure gradient of 2 mbar/min from 0 min to 45 min during the separation.

NOTE: A gradient pressure is applied during the separation to allow for the elution of slow-migrating drugs, which includes some weakly, basic benzodiazepines and neutral/acidic drugs (e.g., acetaminophen, barbiturates). However, most basic DoA and their metabolites migrate as cations within 25 min, including amphetamines, opioids, and other classes of alkaloids.

2.2.4. After the sample acquisition is completed, make sure to rinse the capillary at a low pressure (50 mbar) with BGE until the next day. Otherwise, rinse the capillary for 600 s at a high pressure

(900 mbar) with water and then for 600 s with air, and store it in a sprayer holder until the next use.

2.3. Isocratic pump and sheath liquid

2.3.1. Use an Infinity Isocratic Pump and an Infinity Degasser to deliver the sheath liquid (60:40 MeOH:H₂O with 0.1% v/v formic acid) at a rate of 10 µL/min to the CE-MS sprayer.

2.4. TOF-MS settings

2.4.1. Ensure that a TOF-MS, with a coaxial sheath-liquid electrospray ion source and heated nitrogen gas, is equipped with a CE unit.

2.4.2. Operate the TOF-MS under **positive ion detection** that spanned a mass of m/z 50 - 1,700, with a data acquisition rate of 500 ms/spectrum. Ensure that both the profile and the centroid data are stored in a ".d" file format.

2.4.3. Set the **electrospray ionization conditions** to be: Vcap and nozzle voltage at 2,000 V, the nebulizer gas at 10 psi, and the drying gas delivered at 8 L/min at 300 °C, with a sheath gas flow of 3.5 L/min at 195 °C. In addition, set the MS voltage settings of the fragmentor, skimmer, and Oct1 RF to 120, 65, and 750 V, respectively.

2.4.4. Perform the instrument control and data acquisition using the vendor's software (**Table of Materials**).

2.5. Data analysis

2.5.1. Analyze the MSI-CE-MS data using the vendor's software, including processing, data smoothing, and the integration of extracted ion electropherograms (EIEs) for representative DoA and their metabolites.

2.5.2. Open the software and set the following parameters.

2.5.2.1. Under **chromatograms**, select the **extraction data format**, and **chromatogram and mass spectral data format** to **Profile** mode.

2.5.2.2. Select **Quadratic/Cubic Savitzky Golay** under the **smoothing function** and **function width** to 15 points.

2.5.2.3. Then, click **Integrate (MS)**, set the integrator to be **Agile**, and select the maximum number of peaks limit to the largest 11 under **peak filters**.

2.5.2.4. Click **View**, click **Integration peak list**, and select the **peak number**, **retention time (RT)**, **peak area**, **peak height**, and **signal-to-noise ratio (SNR)**.

2.5.2.5. Save these parameters under a unique method name. Apply this method to the process and interpret each data set.

Note: A summary table with peak number, the RT, the peak area and height, and the SNR displays on the right-hand side. Perform data normalization for measured RT and integrated peak area for each DoA to either a matching d-IS or (if not available) to F-Phe to improve analytical precision.

3. Analysis of the Urine Samples and External Calibration Curves

3.1. Different serial injection configurations used in MSI-CE-MS

3.1.1. For the first serial sample injection configuration, used to demonstrate drug isomer/isobar resolution (**Figure 1A**), ensure that five repeated injections are performed on a mixture of 84 drug standards with d-IS and 4-F-Phe/Cl-Tyr, followed by a sixth injection, which is a blank sample in synthetic urine, and five additional injections of the standard drug mixture.

3.1.2. For the second serial sample injection configuration, used to demonstrate the detection of DoA in individual urine samples from patients with a known prescription (**Figure 2A**), ensure that the first sample plug is a mixture of 84-drug standard mixture at a 1x cut-off screening level (positive control) followed by a randomized injection of 10 diluted urine samples from patients: #293, #88, #309, #43, #281, #64, #221, #208, #183, and #50.

3.1.3. For the third serial sample injection configuration, used to illustrate the acquisition of external calibration curves in a single run, ensure that calibrant solutions for DoA are prepared over a 20-fold linear dynamic range corresponding to 0.25x, 0.5x, 1x, 2.5x, and 5x cut-off concentration levels, together with matching d-IS and F-Phe/Cl-Tyr as additional IS in a synthetic urine matrix ($n = 4$).

REPRESENTATIVE RESULTS:

MSI-CE-MS enables the serial injection of ten or more discrete samples within a single run, which greatly enhances throughput (<3 min/sample) without complicated instrumental modifications, column-switching programs, or costly infrastructure investments (**Figure 1A**). An alternating series of hydrodynamic injections of the sample and electrokinetic spacer of BGE is performed within an unmodified fused-silica capillary, where zonal electrophoretic separations of ions occur under strongly acidic electrolyte conditions (pH 1.8). Solute ionization also occurs under steady-state conditions. In this case a coaxial sheath liquid, with a mass calibrant, is used as an interface for CE-MS under the positive ion mode detection with minimal ion suppression or enhancement effects as monitored by the mass calibrant ions signal. TOF represents a robust yet cost-effective HRMS with a fast data acquisition ideally suited for the nontargeted screening/drug surveillance applications when using MSI-CE-MS. For example, impressive resolution is achieved for various isobaric/isomeric DoA and their metabolites, including structural isomers of three opioid structural isomers, namely norhydrocodone, hydromorphone, and morphine (**Figure 1B**). In this

case, 30 resolved peaks from 10 independent samples of a drug mixture are detected without sample carry-over. A synthetic urine blank/negative urine control is also included within the serial injection series (sample position #6). Also, two other isobaric opioids, namely 6-acetylmorphine (an inactive heroin metabolite) and naloxone (an opioid receptor antagonist used for the emergency treatment of opioid overdose) are fully resolved as 20 distinct peaks with full-scan data acquisition (**Figure 1C**). Similarly, two amphetamine positional isomers are fully resolved in MSI-CE-MS to distinguish illicit methamphetamine abuse from the potential misuse of phentermine, a prescribed stimulant that is used as an appetite suppressant for weight loss.

[Insert **Figure 1**]

To screen DoA, a serial injection configuration in MSI-CE-MS used an 84-drug panel mixture at the recommended screening cut-off concentration levels (as the first injection position for all runs) followed by a randomized analysis of 10 representative urine samples from clinically depressed patients with a known prescription history. For instance, a positive screening test result for methadone (m/z 310.2165) from patient #208 (**Figure 2A**) is deduced by a detection of a large signal peak (originating from injection #9) which comigrates with methadone-d3 with a low mass error (<5 ppm). No other signals are detected in any other urine samples within the same run. Methadone concentration exceeds 13x the recommended cut-off limit when comparing its measured ion response ratio in the sample (injection #9) with the reference drug mixture/QC (injection #1) and corrected by a fourfold urine dilution factor. Thus, this result confirms the patient's adherence to methadone maintenance therapy. Evidence of nonprescribed amphetamine intake (**Figure 2B**) is shown by elevated amphetamine levels (m/z 136.1121) which was only detected in one patient (patient #50, in injection #11) within the MSI-CE-MS run. The measured concentration slightly exceeded the recommended cut-off levels (1.3x). This comigrates with amphetamine-d5 and had a low mass error with a top-ranked molecular formula match. A positive test result for the antidepressant venlafaxine (m/z 278.2115), a selective serotonin-norepinephrine reuptake inhibitor, is also demonstrated in patient #281 (injection #6). In this case, the concentration exceeds the recommended cut-off levels (15x). This is identified by its accurate mass- or molecular formula-match, together with comigrating venlafaxine-d6 (**Figure 2C**). The latter criterion is not met for an unknown isobar that is also detected within the EIE trace, which highlights the need for caution when relying solely on its accurate mass. Also, a definitive detection of prescribed pregabalin, which is prescribed for the treatment of neuropathic pain, as well as of generalized anxiety, is also demonstrated for patient #309, based on its grossly elevated concentration (64x) above the recommended cut-off levels (**Figure 2D**). It is, however, not detected in the nine other patient urine samples analyzed within the same run. Similar to the other screen-positive cases, injection #4, comigrates with pregabalin-d6, which is included in all urine samples analyzed by MSI-CE-MS.

[Insert **Figure 2**]

The absolute quantification of DoA and their metabolites is also achieved by MSI-CE-MS based on external calibration curves, using reference calibrant standards for DoA that are acquired rapidly within a single run. For instance, the serial dilution of an 84-drug calibrant mixture

together with a matching d-IS at a fixed concentration provides a reliable quantitative analysis in complex urine samples. This compensates for the potential matrix-induced ion suppression or enhancement but also for variations in the injection volume in the capillary between the samples. In cases when a matching d-IS is commercially unavailable or cost-prohibitive for certain DoA, a surrogate IS is used for data normalization, such as a d-IS from within the same drug class or a synthetic compound not found in urine (F-Phe), as demonstrated previously¹⁴. Representative EIEs and external calibration curves for oxycodone are shown in **Figure 3A, B.** and for citalopram in **Figure 3C, D.** These are widely prescribed analgesics and antidepressants, respectively, with abuse potential. Relative ion response ratios to their matching d-IS are measured. In this case, a serial injection configuration in MSI-CE-MS, comprising five different drug calibrants, are analyzed in duplicate within a single run together with a synthetic urine blank. Overall, good linearity ($R^2 > 0.990$) over a 20-fold concentration range was achieved with adequate sensitivity for the detection of the majority of DoA and their metabolites (*i.e.*, cationic alkaloids) within the 84-drug panel¹⁴. In all cases, drug metabolites are detected as their protonated molecular ion $[MH^+]$ above their screening cut-off limits (>50 ng/mL), with the exception of certain acidic/neutral drugs that have a poor ionization efficiency under the positive ion mode, such as cannabinoids (*e.g.*, THC-COOH), barbiturates (*e.g.*, secobarbital), and carbamates (*e.g.*, carisoprodol).

[Insert **Figure 3**]

FIGURE AND TABLE LEGENDS:

Figure 1: A series of extracted ion electropherograms (EIE) for representative DoA isomers/isobars that are resolved by MSI-CE-MS by analyzing ten samples and a blank within a single run. (A) Schematic of MSI-CE-MS that depicts the serial injection configuration used for an 84-DoA panel in synthetic urine. This multiplexed separation method uses an alternating series of hydrodynamic injections of 11 discrete samples and blank with an electrokinetic injection of a buffer to initiate the zonal electrophoretic separation of ions, followed by a full-scan data acquisition by TOF-MS with positive ion mode detection. **(B)** Three isobaric functional group opioid isomers (m/z 286.1438), separated by CE, comprising 30 resolved peaks from 10 discrete sample injections, including norhydrocodone, hydromorphone, and morphine. **(C)** Two isobaric opioid drugs and their metabolites (m/z 328.1543), separated by CE, comprising 20 resolved peaks from 10 discrete sample injections, including 6-acetylmorphine (heroin metabolite) and naloxone. **(D)** Two isobaric functional group amphetamine isomers, separated by CE, comprising 20 resolved peaks from 10 discrete sample injections, including methamphetamine and phentermine. In all cases, negative urine controls/blanks introduced at the sixth sample position in MSI-CE-MS had negligible evidence of sample carry-over.

Figure 2: A series of extracted ion electropherograms (EIE) for representative screen-positive urine drug test results from a cohort of 10 clinically depressed patients as confirmed when using MSI-CE-MS, which includes an 84-drug panel at recommended cut-off level, injected as the first sample injection position serving as internal reference/QC. (A) EIE overlay corresponding to methadone-d3, which comigrates with methadone, highlighting that only one

urine sample (injection position #9) has elevated methadone concentrations far exceeding the cut-off level. (B) EIE overlay corresponding to amphetamine-d3, which comigrates with amphetamine, highlighting that only one urine sample (injection position #11) has elevated concentrations above cut-off levels. (C) EIE overlay corresponding to venlafaxine-d6, which comigrates with venlafaxine, highlighting that only one urine sample (injection position #6) has elevated concentrations above cut-off levels. (D) EIE overlay corresponding to pregabalin-d6, which comigrates with pregabalin, highlighting that only one urine sample (injection position #4) has grossly elevated concentrations exceeding cut-off levels. All urine samples were analyzed directly after a fivefold dilution in deionized water together with the addition of a matching d-IS (when available). A screen-positive result in MSI-CE-MS corresponds to a drug that comigrates with d-ISs with a low mass error (<5 ppm) and the correct molecular formula, whose concentration exceeds the cut-off that is analyzed within the same run as internal reference/QC.

Figure 3: A series of extracted ion electropherograms (EIE) for the quantification of the representative DoA. This is done by generating external calibration curves based on their relative ion response ratio with a matching d-IS over a 20-fold linear dynamic range. (A) Duplicate injection of a five-point calibration curve for oxycodone calibrants together with oxycodone-d3 and a blank. (B) External calibration curve for oxycodone, following linear regression to derive sensitivity (slope) and linearity (R^2), with error bars representing $\pm 1\sigma$ ($n = 4$). (C) Duplicate injection of a five-point calibration curve for citalopram calibrants, together with citalopram-d6 and a blank. (D) External calibration curve for citalopram, following linear regression to derive sensitivity (slope) and linearity (R^2), with error bars representing ± 1 SD ($n = 4$).

DISCUSSION:

Conventional chromatographic separations typically rely on a single sample injection per run, which is then followed by a gradient elution for resolving complex drug mixtures and column reconditioning. These requirements fundamentally limit sample throughput and duty cycle even when using optimal column-switching programs. In this context, high-volume urine drug analysis for the workplace, toxicology, or therapeutic monitoring applications by GC-MS/MS and, increasingly, LC-MS/MS are thus performed in parallel but used preferentially as a second-tier or confirmatory test when necessary (*i.e.*, legal or medical context). This is due to the much higher capital investments and operational costs, as well as complications with data analysis when comparing samples analyzed across multiple instrumental platforms within an accredited laboratory. As a result, immunoassays still remain the primary method for routine drug screening despite being prone to false-positives and false-negatives among many drug classes, with limited access to antibody reagents for emerging designer drugs. Previous studies have demonstrated that multiplexed separations based on MSI-CE-MS offer a simple solution to enhance the sample throughput up to one order of magnitude. Additionally, this approach allows for the design of novel data workflows for biomarker discovery with high data fidelity based on the implementation of effective batch correction and quality control¹²⁻¹⁴. However, the introduction of 10 or more serial hydrodynamic injections in MSI-CE-MS shortens the effective capillary length required to maintain high-efficiency separations, which may compromise the selectivity when analyzing DoA and their metabolites in human urine.

Herein, we introduced an electrokinetic injection of the BGE following each hydrodynamic sample injection, such that the separation takes advantage of the full capillary length (120 cm), thus improving the resolution of important drug isobars/isomers when using full-scan data acquisition by TOF-MS (**Figure 1A**). In comparison to a recent report¹⁴, better resolution is achieved for several important isobaric/isomeric DoA and their metabolites, which is critical when screening for large drug panels. For instance, the resolution of several important structural/functional group isomers was realized, including three analogous opioid drugs/metabolites (**Figure 1B**), two unrelated opioid isobars (**Figure 1C**), and two methamphetamine positional isomers (**Figure 1D**). In all cases, sample carry-over from serial injection of different calibrant solutions within the same run is not significant as confirmed by a negative urine control/blank. Additionally, improved resolution is also realized for several cross-interferences involving certain d-ISs with isobaric drugs within the panel, such as cotinine-d3 and 3,4-methylenedioxyamphetamine (MDA), EDDP-d3 and imipramine, norfentanyl-d5 and ketamine, codeine-d6 and sertraline, and cocaine-d3 and zolpidem. This outcome is extended to other isobaric interferences within the drug panel that were better resolved in this study (*e.g.*, noroxycodone/oxymorphone, normeperidine/methylphenidate), as well as major background urinary interferences (*e.g.*, in-source fragment ions of creatinine with amphetamine)¹⁴. Indeed, access to two orthogonal parameters to satisfy the putative identification of a specific DoA is critical to reducing false-positives due to isobaric interferences, namely accurate mass combined with comigration with d-ISs. Also, the detection of one or more biotransformed metabolites of the parent drug within the same sample, such as hydroxylated, demethylated, or intact glucuronide drug conjugate(s), adds further confidence toward drug identification while expanding the window for detection.

The recommended cut-off levels for urine drug screening vary widely for different classes of DoA (ranging from 50 to 1,000 ng/mL) depending on their pharmacokinetics, toxicity, and background interferences, in order to reduce method bias based on guidelines from the Substance Abuse and Mental Health Services Administration (SAMHSA)¹⁴. The potential for MSI-CE-MS to detect and identify diverse classes of DoA directly in urine with minimal sample pretreatment was applied to a group of clinically depressed patients with a known prescription record. The definitive identification of methadone (prescribed), amphetamine (nonprescribed/illicit), venlafaxine (prescribed), and pregabalin (prescribed) in diluted, yet nonhydrolyzed, urine samples was demonstrated when using MSI-CE-MS. This was based on a direct comparison of a drug detected in a specific injection position relative to the 84-drug mixture introduced in the first sample position at the recommended screening cut-off level which serves as an internal reference/QC and positive control (**Figure 2**). Also, absolute drug quantification is feasible when using external calibration curves (**Figure 3**) based on the ion response ratio measured for a drug relative to its d-IS. Unlike most chromatographic separations, there is no deuterium effect impacting migration time differences between a d-IS and its nondeuterated drug since they possess analogous electrophoretic mobilities in free solution in CE. Indeed, the migration behavior of DoA is accurately modeled in CE, based on their fundamental physicochemical properties/chemical structure, namely molecular volume and effective charge (pK_a)¹⁴. Since many drugs also undergo significant secondary metabolism prior to excretion in urine (*e.g.*, morphine glucuronide), screening cut-off levels require adjustment as their measured concentrations are lower than

anticipated when compared to methods that utilize enzyme hydrolysis for total drug detection. A major benefit of “dilute-and-shoot” urine drug testing, besides reducing costs/time, sample handling, and a potential bias or batch variations due to incomplete enzyme hydrolysis, is that presumptive screen-positive cases are further confirmed by the detection of one or more related drug metabolites within the same sample. This also provides deeper insights into drug pharmacokinetics and optimum dosage requirements for individual patients while improving the detection for “fast” metabolizers or prescribed/illicit drugs with short half-lives. Additionally, comigration with a matching d-IS plays two important functions for reliable drug screening when using MSI-CE-MS—namely, it identifies the exact sample injection position (*i.e.*, patient #) while also correcting for differences in ion responses/injection volumes for improved precision and accuracy.

Future work aims to develop customized software tools to facilitate automated data processing of multiplexed separations coupled to HRMS as required for high-volume urine drug testing with QC/QA. Rigorous validation of MSI-CE-MS for the broad-spectrum screening of DoA will also be investigated among a larger cohort of high-risk patients in order to objectively evaluate prescribed drug adherence and potential misuse/substitution that may compromise treatment efficacy, patient safety, and/or psychiatric evaluation/diagnosis. A complementary analysis of acidic/anionic classes of DoA and their metabolites will also be performed by MSI-CE-MS under alkaline conditions with negative ion mode detection as required for the comprehensive screening of natural/synthetic cannabinoids. This is important given the looming public health implications of the legalization of recreational marijuana across Canada and several US states. A major advantage of full-scan data acquisition by TOF-MS is that retrospective analysis of samples can be performed even when urine specimens are no longer available for follow-up testing, whereas other lifestyle or dietary exposures can be assessed to better understand differential responses to drug therapy. In summary, a rapid yet accurate drug surveillance method by MSI-CE-MS offers significant advantages to conventional targeted immunoassays, as well as direct infusion/ambient ionization-MS/MS methods that are prone to interferences/bias when resolving expanded panels of DoA and their metabolites in complex biological samples at incremental costs.

ACKNOWLEDGMENTS:

P.B.M. wishes to acknowledge funding support from the Natural Sciences and Engineering Research Council of Canada, the Canada Foundation for Innovation, Genome Canada and McMaster University. The authors thank Howard Lee at Seroclinx Corporation and Dr. Marcus Kim from Agilent Technologies for their insightful discussions. Additionally, the authors acknowledge Dr. Zainab Samaan from the Department of Psychiatry and Behavioral Sciences at McMaster University and the Moods Disorder Clinic at St. Joseph’s Hospital for access to the deidentified patient urine samples used in this study.

DISCLOSURES:


The authors disclose a US patent (PCT/CA2014/050454) on MSI-CE-MS as a multiplexed screening platform and data workflow for chemical analysis.

REFERENCES:


1. Products - Vital Statistics Rapid Release - Provisional Drug Overdose Data. <https://www.cdc.gov/nchs/nvss/vsrr/drug-overdose-data.htm>. September 27, 2018 (2018).
2. Olfson, M., King, M., Schoenbaum, M. Treatment of Young People With Antipsychotic Medications in the United States. *JAMA Psychiatry*. **72**, 867-874 (2015).
3. Moeller, K. E., Kissack, J. C., Atayee, R. S., Lee, K. C. Clinical Interpretation of Urine Drug Tests: What Clinicians Need to Know About Urine Drug Screens. *Mayo Clinic Proceedings*. **92**, 774-796 (2017).
4. Levy, S.; Siqueira, L. M., Committee on Substance Abuse. Testing for Drugs of Abuse in Children and Adolescents. *American Academy of Pediatrics*. **133**, e1798 (2015).
5. Pesce, M., Mikel, C., West, C. A tale of two drug testing technologies: GC-MS and LC-MS/MS. *American Society of Interventional Pain Physicians*. **13**, 91-92 (2010).
6. Saitman, A., Park, H-D., Fitzgerald, R. L. False-Positive Interferences of Common Urine Drug Screen Immunoassays: A Review. *Journal of Analytical Toxicology*. **38**, 387-396 (2014).
7. Segal, D. In Pursuit of Liquid Gold. *The New York Times*. <https://www.nytimes.com/interactive/2017/12/27/business/urine-test-cost.html>. September 27, 2018 (2017).
8. Wu, A. H. *et al.* Role of liquid chromatography-high resolution mass spectrometry (LC-HR/MS) in clinical toxicology. *Clinical Toxicology*. **50**, 733-742 (2012).
9. Gualé, F. *et al.* Validation of LC-TOF-MS Screening for Drugs, Metabolites, and Collateral Compounds in Forensic Toxicology Specimens. *Journal of Analytical Toxicology*. **37**, 17-24 (2013).
10. Kaupilla, T. J. *et al.* Rapid analysis of metabolites and drugs of abuse from urine samples by desorption electrospray ionization-mass spectrometry. *The Analyst*. **132**, 868-875 (2007).
11. Bynum, N. D., Moore, K. N., Grabenauer, M. Evaluation of Laser Diode Thermal Desorption-Tandem Mass Spectrometry (LDTD-MS-MS) in Forensic Toxicology. *Journal of Analytical Toxicology*. **38**, 528-535 (2014).
12. Kuehnbaum, N.L., Kormendi, A., Britz-McKibbin, P. Multisegment Injection-Capillary Electrophoresis-Mass Spectrometry: A High-Throughput Platform for Metabolomics with High Data Fidelity. *Analytical Chemistry*. **85**, 10664-10669 (2013).
13. Nori de Macedo, A. *et al.* The Sweat Metabolome of Screen-Positive Cystic Fibrosis Infants: Revealing Mechanisms Beyond Impaired Chloride Transport. *ACS Central Science*. **3**, 904-913 (2017).
14. DiBattista, A., Rampersaud, D., Lee, H., Kim, M., Britz-McKibbin P. High Throughput Screening Method for Systematic Surveillance of Drugs of Abuse by Multisegment Injection-Capillary Electrophoresis-Mass Spectrometry. *Analytical Chemistry*. **89**, 11853-11861 (2017).
15. Cao, Z., Kaleta, E., Wang, P. Simultaneous Quantitation of 78 Drugs and Metabolites in Urine with a Dilute-And-Shoot LC-MS-MS Assay. *Journal of Analytical Toxicology*. **29**, 335-346 (2015).
16. Drouin, N., Rudaz, S., Schappler, J. New Supported Liquid Membrane for Electromembrane Extraction of Polar Basic Endogenous Metabolites. *Journal of Pharmaceutical and Biomedical Analysis*. **159**, 53-59 (2018).

(A)


Serial 11 Sample Injection Series in MSI-CE-MS




(B)




(C)




(D)




(A)




(B)







(C)



(D)





Name of Material/ Equipment	Company	Catalog Number
7100 Capillary Electrophoresis System	Agilent Technologies Inc.	G7100A
6230 Series Time-of-Flight Mass Spectrometer	Agilent Technologies Inc.	G6230B
CE-ESI-MS Sprayer Kit	Agilent Technologies Inc.	G1603A
1260 Infinity Isocratic Pump and Degasser	Agilent Technologies Inc.	G1310B
MassHunter Workstation Data Acquisition Software (B.06.01)	Agilent Technologies Inc.	--
MassHunter Qualitative Analysis Software (B.06.01)	Agilent Technologies Inc.	--
Shortix Capillary Cutter	Agilent Technologies Inc.	5813-4620
Capillary Window Maker	Microsolv Inc.	07200-S
Flexible Fused-silica Capillary Tubing	Polymicro Technologies Inc.	TSP05375
Drug standards, deuterated internal standards, synthetic urine matrix (SURINE)	Cerilliant Inc.	Miscellaneous

Comments/Description

CE instrument used for separation of drug mixtures, desalting and anotation

HRMS mass analyzer used for drug detection and anotation

CE/MS coaxial sheath liquid interface and capillary cassette

Isocratic pump to deliver sheath liquid/mass calibrant

Software used for control of CE-MS system

Software used for processing of CE-MS data

Cutting tool with diamond blade used to cut capillaries

Burner with 7 mm window size to remove polyimide coating from CE capillary

Standard polyimide coated fused-silica capillary for CE separation (50 micron ID; 360 micron OD)

Certified drugs of abuse reference standards (86 drug panel) with 48 deuterated internal standards and negative urine control (Surine)

ARTICLE AND VIDEO LICENSE AGREEMENT

Title of Article:	High Throughput and Comprehensive Drug Surveillance using Multisegment Injection-Capillary Electrophoresis-Mass Spectrometry
Author(s):	Meera Shanmuganathan, Sabrina Macklai, Clara Barrenas Cárdenas, Zachary Kroezen, William Zizek, Marcus Kim, Howard Lee and Philip Britz-McKibbin*

Item 1: The Author elects to have the Materials be made available (as described at <http://www.jove.com/publish>) via:

☐

Standard Access

☒

Open Access

Item 2: Please select one of the following items:

☒

The Author is **NOT** a United States government employee.

☐

The Author is a United States government employee and the Materials were prepared in the course of his or her duties as a United States government employee.

☐

The Author is a United States government employee but the Materials were NOT prepared in the course of his or her duties as a United States government employee.

ARTICLE AND VIDEO LICENSE AGREEMENT

1. **Defined Terms.** As used in this Article and Video License Agreement, the following terms shall have the following meanings: **"Agreement"** means this Article and Video License Agreement; **"Article"** means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein; **"Author"** means the author who is a signatory to this Agreement; **"Collective Work"** means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; **"CRC License"** means the Creative Commons Attribution-Non Commercial-No Derivs 3.0 Unported Agreement, the terms and conditions of which can be found at: <http://creativecommons.org/licenses/by-nc-nd/3.0/legalcode>; **"Derivative Work"** means a work based upon the Materials or upon the Materials and other pre-existing works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; **"Institution"** means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; **"JoVE"** means MyJoVE Corporation, a Massachusetts corporation and the publisher of The Journal of Visualized Experiments; **"Materials"** means the Article and / or the Video; **"Parties"** means the Author and JoVE; **"Video"** means any video(s) made by the Author, alone or in conjunction with any other parties, or by JoVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion

of the Article, and in which the Author may or may not appear.

2. **Background.** The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JoVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.

3. **Grant of Rights in Article.** In consideration of JoVE agreeing to publish the Article, the Author hereby grants to JoVE, subject to **Sections 4** and **7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the "Open Access" box has been checked in **Item 1** above, JoVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.

ARTICLE AND VIDEO LICENSE AGREEMENT

4. **Retention of Rights in Article.** Notwithstanding the exclusive license granted to JoVE in **Section 3** above, the Author shall, with respect to the Article, retain the non-exclusive right to use all or part of the Article for the non-commercial purpose of giving lectures, presentations or teaching classes, and to post a copy of the Article on the Institution's website or the Author's personal website, in each case provided that a link to the Article on the JoVE website is provided and notice of JoVE's copyright in the Article is included. All non-copyright intellectual property rights in and to the Article, such as patent rights, shall remain with the Author.

5. **Grant of Rights in Video – Standard Access.** This **Section 5** applies if the "Standard Access" box has been checked in **Item 1** above or if no box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby acknowledges and agrees that, Subject to **Section 7** below, JoVE is and shall be the sole and exclusive owner of all rights of any nature, including, without limitation, all copyrights, in and to the Video. To the extent that, by law, the Author is deemed, now or at any time in the future, to have any rights of any nature in or to the Video, the Author hereby disclaims all such rights and transfers all such rights to JoVE.

6. **Grant of Rights in Video – Open Access.** This **Section 6** applies only if the "Open Access" box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby grants to JoVE, subject to **Section 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Video in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Video into other languages, create adaptations, summaries or extracts of the Video or other Derivative Works or Collective Works based on all or any portion of the Video and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. For any Video to which this **Section 6** is applicable, JoVE and the Author hereby grant to the public all such rights in the Video as provided in, but subject to all limitations and requirements set forth in, the CRC License.

7. **Government Employees.** If the Author is a United States government employee and the Article was prepared in the course of his or her duties as a United States government employee, as indicated in **Item 2** above, and any of the licenses or grants granted by the Author hereunder exceed the scope of the 17 U.S.C. 403, then the rights granted hereunder shall be limited to the maximum

rights permitted under such statute. In such case, all provisions contained herein that are not in conflict with such statute shall remain in full force and effect, and all provisions contained herein that do so conflict shall be deemed to be amended so as to provide to JoVE the maximum rights permissible within such statute.

8. **Protection of the Work.** The Author(s) authorize JoVE to take steps in the Author(s) name and on their behalf if JoVE believes some third party could be infringing or might infringe the copyright of either the Author's Article and/or Video.

9. **Likeness, Privacy, Personality.** The Author hereby grants JoVE the right to use the Author's name, voice, likeness, picture, photograph, image, biography and performance in any way, commercial or otherwise, in connection with the Materials and the sale, promotion and distribution thereof. The Author hereby waives any and all rights he or she may have, relating to his or her appearance in the Video or otherwise relating to the Materials, under all applicable privacy, likeness, personality or similar laws.

10. **Author Warranties.** The Author represents and warrants that the Article is original, that it has not been published, that the copyright interest is owned by the Author (or, if more than one author is listed at the beginning of this Agreement, by such authors collectively) and has not been assigned, licensed, or otherwise transferred to any other party. The Author represents and warrants that the author(s) listed at the top of this Agreement are the only authors of the Materials. If more than one author is listed at the top of this Agreement and if any such author has not entered into a separate Article and Video License Agreement with JoVE relating to the Materials, the Author represents and warrants that the Author has been authorized by each of the other such authors to execute this Agreement on his or her behalf and to bind him or her with respect to the terms of this Agreement as if each of them had been a party hereto as an Author. The Author warrants that the use, reproduction, distribution, public or private performance or display, and/or modification of all or any portion of the Materials does not and will not violate, infringe and/or misappropriate the patent, trademark, intellectual property or other rights of any third party. The Author represents and warrants that it has and will continue to comply with all government, institutional and other regulations, including, without limitation all institutional, laboratory, hospital, ethical, human and animal treatment, privacy, and all other rules, regulations, laws, procedures or guidelines, applicable to the Materials, and that all research involving human and animal subjects has been approved by the Author's relevant institutional review board.

11. **JoVE Discretion.** If the Author requests the assistance of JoVE in producing the Video in the Author's facility, the Author shall ensure that the presence of JoVE employees, agents or independent contractors is in accordance with the relevant regulations of the Author's institution. If more than one author is listed at the beginning of this Agreement, JoVE may, in its sole

ARTICLE AND VIDEO LICENSE AGREEMENT

discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JoVE reserves the right, in its absolute and sole discretion and without giving any reason therefore, to accept or decline any work submitted to JoVE. JoVE and its employees, agents and independent contractors shall have full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JoVE has sole discretion as to the method of making and publishing the Materials, including, without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

12. **Indemnification.** The Author agrees to indemnify JoVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JoVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JoVE, making of videos by JoVE, or publication in JoVE or elsewhere by JoVE. The Author shall be responsible for, and shall hold JoVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to


the making of a video by JoVE its employees, agents or independent contractors. All sterilization, cleanliness or decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's expense. All indemnifications provided herein shall include JoVE's attorney's fees and costs related to said losses or damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JoVE, its employees, agents or independent contractors.

13. **Fees.** To cover the cost incurred for publication, JoVE must receive payment before production and publication the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JoVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.

14. **Transfer, Governing Law.** This Agreement may be assigned by JoVE and shall inure to the benefits of any of JoVE's successors and assignees. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to be one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement is required per submission.

CORRESPONDING AUTHOR

Name:	Philip Britz-McKibbin		
Department:	Department of Chemistry and Chemical Biology		
Institution:	McMaster University		
Title:	High Throughput and Comprehensive Drug Surveillance using Multisegment Injection-Capillary Electrophoresis-Mass Spectrometry		
Signature:		Date:	Aug. 23, 2018

Please submit a **signed** and **dated** copy of this license by one of the following three methods:

1. Upload an electronic version on the JoVE submission site
2. Fax the document to +1.866.381.2236
3. Mail the document to JoVE / Attn: JoVE Editorial / 1 Alewife Center #200 / Cambridge, MA 02140



www.chemistry.mcmaster.ca

Department of Chemistry, McMaster University
1280 Main Street West, Hamilton, ON Canada L8S 4M1
(905) 525-9140 ext. 24504 (voice); (905) 522-2509 (fax)

27 September, 2018

Dr. Phillip Steindel
Review Editor, JoVE

Re: Revised Manuscript: "High Throughput and Comprehensive Drug Surveillance using Multisegment Injection-Capillary Electrophoresis-Mass Spectrometry (JoVE58986)"

Dear Dr. Steindel,

Please note that the editorial and reviewer's comments are in *italic*, whereas our responses are in normal font. All changes made to the revised manuscript based on the reviewer comments are highlighted using track changes to identify all of the edits.

Editorial Request for Formatting Changes:

We thank the editor for highlighting the changes required for our manuscript and the changes are highlighted (track change format) in the revised version of the manuscript.

-Please take this opportunity to thoroughly proofread the manuscript to ensure that there are no spelling or grammar issues.

-Please revise lines 157-159, 189-190, and 214-215 to avoid previously published text.

-Please define all abbreviations before use.

We have updated these changes in the revised manuscript

-JoVE cannot publish manuscripts containing commercial language. This includes trademark symbols (™), registered symbols (®), and company names before an instrument or reagent. Please remove all commercial language from your manuscript and use generic terms instead. All commercial products should be sufficiently referenced in the Table of Materials and Reagents. You may use the generic term followed by "(see Table of Materials)" to draw the readers' attention to specific commercial names. Examples of commercial sounding language in your manuscript are: Surine™, Cerilliant, Round Rock, Polymicro Technologies Inc., MicroSolv, Agilent, Wavemetric Inc., Igor Pro 5.0, Kim wipe, Agilent MassHunter, Oxycontin™, Celexa™, Effexor™, etc.

We have removed these names from the revised manuscript, and we have included these commercial trade names or suppliers only in the Table of Materials as requested.

-Please adjust the numbering of the Protocol to follow the JoVE Instructions for Authors. For example, 1 should be followed by 1.1 and then 1.1.1 and 1.1.2 if necessary. Please refrain from using bullets, dashes, or indentations.

We have updated the numbering of the protocol based on the recommendations above.

-Please revise the protocol to contain only action items that direct the reader to do something (e.g., “Do this,” “Ensure that,” etc.). The actions should be described in the imperative tense in complete sentences wherever possible. Avoid usage of phrases such as “could be,” “should be,” and “would be” throughout the Protocol. Any text that cannot be written in the imperative tense may be added as a “Note.” Please include all safety procedures and use of hoods, etc. However, notes should be used sparingly and actions should be described in the imperative tense wherever possible.

-Line 105: The meaning of this statement is unclear. If details are presented in lines 107-110, please combine these two steps.

We have combined lines 107-110 with line 105

-Line 109: Do you mean that sonicate and degas the solution at the same time? If not, please specify the time needed for each procedure.

We have updated the statement to “sonicate the solution for 15 min”

-Line 112: The meaning of this statement is unclear. If details are presented in lines 113-119, please combine these two steps.

We have combined lines 113-119 with line 112

-Lines 121-128: Please describe how to prepare standards in the imperative tense.

-Lines 130-142: Please describe how to prepare calibration curves in the imperative tense.

-Line 165: Please describe how capillary conditioning is actually done in the imperative tense.

We have described the protocols in the imperative tense.

-Line 183: What is the flow rate of BGE?

The BGE flow rate is approximately 10 nL/min but it is determined by the magnitude of the electroosmotic flow and the applied pressure gradient. As a result, it is not a fixed flow rate setting like in conventional chromatography.

- Line 201: Where is sheath liquid delivered?

We have updated the requested information. It is being delivered to the CE-MS interface to provide a make-up flow for stable spray formation.

-Lines 205-206: Please move the equipment information to the Table of Materials.

Yes, we have moved the equipment information to the Table of Materials.

-Lines 218-225: Please describe how to actually perform data analysis. To be filmed, software must have a GUI (graphical user interface) and software steps must be more explicitly explained ('click', 'select', etc.). Please add more specific details (e.g. button clicks for software actions, numerical values for settings, etc.) to your protocol steps.

We have described data analysis section in details as requested.

-Lines 230-244: Please describe in the imperative tense the specific actions being performed here.

We have described the specific action in the imperative tense.

-Please include single-line spaces between all paragraphs, headings, steps, etc.

-After you have made all the recommended changes to your protocol (listed above), please highlight 2.75 pages or less of the Protocol (including headings and spacing) that identifies the essential steps of the protocol for the video, i.e., the steps that should be visualized to tell the most cohesive story of the Protocol.

-Please highlight complete sentences (not parts of sentences). Please ensure that the highlighted part of the step includes at least one action that is written in imperative tense.

-Please include all relevant details that are required to perform the step in the highlighting. For example: If step 2.5 is highlighted for filming and the details of how to perform the step are given in steps 2.5.1 and 2.5.2, then the sub-steps where the details are provided must be highlighted.

We have updated the recommended changes and highlighted the section for filming.

-Discussion: As we are a methods journal, please also discuss critical steps within the protocol, any modifications and troubleshooting of the technique, and any limitations of the technique.

We have addressed the above recommendations.

-References: Please do not abbreviate journal titles.

We have made the requested changes in the revised manuscript.

-Reviewer #1:

Comments to the Author

This publication described the protocol for a high throughput and sensitive screening method for identification and quantification drugs of abuse in urine. Preparation of BGE and solution are clearly described.

-However, I think more information should be provided regarding compounds identification and quantification when no d-IS is available.

We thank **reviewer #1** for the critical comments concerning drug identification and quantification when matching deuterated internal standards are unavailable. Supporting information was provided in the long abstract and discussion sections. For instance, “detection of one or more bio-transformed metabolites of parent drug adds further confidence towards unambiguous identification” and “when a matching D-IS is unavailable for certain DoA, a surrogate IS is used for data normalization using a synthetic chemical standard not found in urine (F-Phe) as described previously [ref. 14]”. The latter point is repeated both in results and discussion section of the revised manuscript, which we have also previously demonstrated in our previous work.

-Do you think this approach could also be used for instance for highly polar substances, such as pesticides screening?

Indeed, our methodology is applicable for high throughout screening of other classes of polar/ionic chemicals in complex sample matrices, such as anionic pesticides (e.g., glyphosate, glufosinate, ethephon etc.).

-I think this manuscript could gain in clarity with a table describing all the steps of the injection process.

We have updated the description of the serial injection process then in the revised manuscript, but have kept the format requested by the journal for the protocol. This process will be much clearer as well during the recoding session of our protocol by video.

-I think description of the procedure for Matrix effect identification and correction should be discussed. (line 309-311), especially when no d-IS is available.

Please refer to our responses above. Additionally, we have also mentioned in the discussion and protocol sections that the mass calibrant ions added to the sheath liquid also serve as an effective way to monitor for potential ion suppression or enhancement effects in “real-time” during separation, which was also described extensively in our previous manuscript [ref. 14].

-As shown highlighted in figure 1B, you have a diminution of the isobaric resolution, probably due to the gradient pressure used (2mbar/min). How do you manage pic integration to get reliable quantification in such circumstance?

Indeed, Figure 1B is an extreme (yet impressive) example of simultaneous resolution of three isomeric opioids (equimolar concentrations of standards in synthetic urine) from 10 independent samples injected within a single run when using MSI-CE-MS. It is true that baseline resolution is not fully achieved for the later migrating peaks due to increased longitudinal diffusion; however, this is still sufficient for reliable drug screening purposes which often requires semi-quantitative performance to detect drugs above a minimum threshold/cutoff level. Moreover, detection of all three opioid isomers in all urine samples analyzed is also not likely to occur under most circumstances. The reviewer is reminded that drug screening by immunoassays has no ability to discriminate among such complex mixtures of opioid isomers with adequate selectivity.

-Reference 1 and 7: for internet references, please add the visiting date.

-Reference 10, the complete title is "Rapid analysis of metabolites and drugs of abuse from urine samples by desorption electrospray ionization-mass spectrometry" instead of "Rapid analysis of metabolites and drugs of abuse from urine samples"

-Reference 13, the right title should be "The Sweat Metabolome of Screen-Positive Cystic Fibrosis Infants: Revealing Mechanisms beyond Impaired Chloride Transport". Also correct the referenced page of the publication (p904-913)

-Reference 14, the right title should be "High Throughput Screening Method for Systematic Surveillance of Drugs of Abuse by Multisegment Injection-Capillary Electrophoresis-Mass Spectrometry" instead of "A high throughput screening method for systematic surveillance of drugs of abuse by multisegment injection capillary electrophoresis-mass spectrometry"

We have updated all the references cited above as recommended by the reviewer.

-Line 392: this work is not the first introduction of MSI-CE-MS based on electrokinetic plug. This approach was already described and used in another publication (DOI: 10.1016/j.jpba.2018.06.029). This should be cited and reference added in the list.

We have updated the recently published reference [ref 16] in the introduction. In fact, as this manuscript was primarily focused on optimization of electroextraction preconcentration, we were not aware of the use of MSI-CE-MS as it is only briefly described in the methodology section. Also, there were no explicit comments on the rationale of using the electrokinetic spacer plugs for sample, which we have elaborated in the current manuscript.

-Figure 2C is referenced twice (lines 290 and 298). According to your manuscript the referenced figure in line 298 should be Figure 2D.

We have updated the figure number correctly as requested


-Reviewer #2:

-This is an excellent manuscript showing the utility of multi-segment injection CE-MS for the quantitative analysis of various drugs of abuse in patient urine samples. The experimental procedures are described clearly and in detail. Given the opioid epidemic or opioids crisis in various countries, I expect that the proposed analytical tool will have a great impact in the screening of DoAs.

We thank reviewer # 2 for the positive comments about our manuscript.

We thank the editor and both external reviewers for their constructive feedback, and hope that our revised manuscript is deemed acceptable for publication.

With best regards,

A handwritten signature in blue ink, appearing to read "Philip Britz-McKibbin".

Philip Britz-McKibbin
Professor and Cystic Fibrosis Canada Researcher
Department of Chemistry and Chemical Biology
McMaster University
E-mail: britz@mcmaster.ca