

Journal of Visualized Experiments

Use of single chain MHC technology to investigate co-agonism in human CD8+ T cell activation --Manuscript Draft--

Article Type:	Methods Article - JoVE Produced Video
Manuscript Number:	JoVE59126R1
Full Title:	Use of single chain MHC technology to investigate co-agonism in human CD8+ T cell activation
Keywords:	Single chain MHC, human T cell activation, T cell antigen recognition, self peptide, stable cell line, engineered antigen presenting cells, PBMC isolation, human T cell degranulation assay, human T cell cytokine production
Corresponding Author:	Joanna Brzostek National University of Singapore Singapore , Singapore SINGAPORE
Corresponding Author's Institution:	National University of Singapore
Corresponding Author E-Mail:	micjmb@nus.edu.sg
Order of Authors:	Joanna Brzostek Xiang Zhao Maryam Hamidinia Joanna Ai Ling Choo Chien Tei Too Zi Zong Ho Ee Chee Ren Antonio Bertoletti Paul A. MacAry Keith G. Gould Nicholas R.J. Gascoigne
Additional Information:	
Question	Response
Please indicate whether this article will be Standard Access or Open Access.	Open Access (US\$4,200)
Please indicate the city, state/province, and country where this article will be filmed . Please do not use abbreviations.	Singapore

September 17th, 2018

The Editors
Journal of Visualized Experiments


Dear Editors,

We wish to submit an original article entitled "Use of single chain MHC technology to investigate co-agonism in human CD8⁺ T cell activation" for consideration by Journal of Visualized Experiments.

We have recently shown that non-stimulatory peptide MHC class I complexes can enhance human CD8⁺ T cell activation through process termed co-agonism (Zhao et al., 2018 Nonstimulatory peptide-MHC enhances human T-cell antigen-specific responses by amplifying proximal TCR signaling. Nature Communications 13:2716. doi: 10.1038/s41467-018-05288-0). Our methods manuscript describes an experimental system to investigate co-agonism during human CD8⁺ T cell activation by expressing human MHC class I molecules presenting pre-determined peptides as single polypeptides (single chain MHC) in a xenogeneic cell line. The protocol describes: cell line transfection with single chain MHC constructs, generation of stable cell lines, culture of hepatitis B virus-specific human CD8⁺ T cells and T cell activation experiments simultaneously quantifying cytokine production and degranulation. The presented methods and the single chain MHC class I antigen presentation system can be used for research on different aspects of CD8⁺ T cell activation in human T cell systems with known peptide MHC specificity.

We confirm that this work is original and has not been published elsewhere, nor is it currently under consideration for publication elsewhere.

Sincerely yours,



Joanna Brzostek, Ph.D.
Senior Research Fellow
Prof. N. Gascoigne laboratory
Department of Microbiology and Immunology
Yong Loo Lin School of Medicine,
National University of Singapore
micjmb@nus.edu.sg

TITLE:

Use of Single Chain MHC Technology to Investigate Co-Agonism in Human CD8+ T Cell Activation

AUTHORS AND AFFILIATIONS:

Xiang Zhao^{1*}, Maryam Hamidinia¹, Joanna Ai Ling Choo¹, Chien Tei Too^{1,2}, Zi Zong Ho^{3*}, Ee Chee Ren⁴, Antonio Bertoletti³, Paul A. MacAry^{1,2,5}, Keith G. Gould⁶, Joanna Brzostek¹, and Nicholas R.J. Gascoigne^{1,2,5}

¹Department of Microbiology and Immunology, Yong Loo Lin School of Medicine, National University of Singapore, Singapore

²Immunology Programme, Life Sciences Institute, National University of Singapore, Singapore

³Emerging Infectious Diseases Program, Duke-NUS Graduate Medical School, Singapore

⁴Singapore Immunology Network, A*STAR, Singapore

⁵NUS Graduate School for Integrative Sciences and Engineering (NGS), National University of Singapore, National University of Singapore

⁶Department of Immunology, Wright-Fleming Institute, Imperial College London, United Kingdom

*Current Addresses: XZ: Department of Molecular and Cellular Physiology, Stanford University School of Medicine, Stanford, CA, USA; ZZH: Lion TCR, Singapore.

Xiang Zhao
xiangzh@stanford.edu

Maryam Hamidinia
e0081746@u.nus.edu

Joanna Ai Ling Choo
joannachooal@gmail.com

Chien Tei Too
mictct@nus.edu.sg

Zi Zong Ho
zack.ho@liontcr.com

Ee Chee Ren
ren_ee_chee@immunol.a-star.edu.sg

Antonio Bertoletti
antonio@duke-nus.edu.sg

Paul A. MacAry

micpam@nus.edu.sg
Keith G. Gould
k.gould@imperial.ac.uk

Corresponding authors:
Joanna Brzostek
micjmb@nus.edu.sg

Nicholas R.J. Gascoigne
micnrjg@nus.edu.sg

KEYWORDS:

Single chain MHC, human T cell activation, T cell antigen recognition, self peptide, stable cell line, engineered antigen presenting cells, PBMC isolation, human T cell degranulation assay, human T cell cytokine production

SUMMARY:

This protocol describes the use of single chain MHC class I complexes to investigate molecular interactions in human CD8+ T cell activation: generation of engineered antigen presenting cells expressing single chain constructs, culture of human CD8+ T cell clone and T cell activation experiments.

ABSTRACT:

Non-stimulatory self peptide MHC (pMHC) complexes do not induce T cell activation and effector functions, but can enhance T cell responses to agonist pMHC, through a process termed co-agonism. This protocol describes an experimental system to investigate co-agonism during human CD8+ T cell activation by expressing human MHC class I molecules presenting pre-determined peptides as single polypeptides (single chain MHC) in a xenogeneic cell line. We expressed single chain MHCs under conditions where low levels of agonist single chain p-MHC complexes and high levels of non-stimulatory single chain p-MHC complexes were expressed. Use of this experimental system allowed us to compare CD8+ T cell responses to agonist pMHC in the presence or absence of non-stimulatory pMHC. The protocol describes cell line transfection with single chain MHC constructs, generation of stable cell lines, culture of hepatitis B virus-specific human CD8+ T cells and T cell activation experiments simultaneously quantifying cytokine production and degranulation. The presented methods can be used for research on different aspects of CD8+ T cell activation in human T cell systems with known peptide MHC specificity.

INTRODUCTION:

MHC class I (MHC-I) molecules present short (8-10 amino acids) peptides derived from proteins synthesized by each cell. MHC-I molecules on healthy cells present self peptides derived from endogenous proteins. Viral infection leads to presentation of non-self virus-derived peptides on MHC-I, but even infected cells present non-self peptides in the presence of large amounts of self peptide-MHC (pMHC). MHC-I is recognized by the T cell receptor (TCR) and the CD8 co-

receptor on T cells. TCR binding affinity to peptide pMHC is highly dependent on the sequence of presented peptides, whereas CD8 binding is peptide-independent. TCR binding to non-self agonist pMHC complex leads to T cell activation and effector functions. Non-stimulatory self pMHC complexes do not induce T cell activation and effector functions, but can enhance T cell responses to agonist pMHC, through a process termed co-agonism¹⁻³. Co-agonism has been observed in mouse CD4⁺ T cells⁴⁻⁶ as well as in mouse and human CD8⁺ T cells⁷⁻¹³. Under physiological conditions, T cells need to recognize limited quantities of agonist pMHC on antigen presenting cells (APCs) co-presenting high levels of non-stimulatory pMHC complexes, suggesting that co-agonism contributes to optimal T cell responses in vivo. Total cell surface MHC class I levels are frequently downregulated during viral infections¹⁴ and on tumors¹⁵. This immune evasion strategy decreases agonist pMHC presentation, but also reduces the amount of co-agonist pMHC I available for T cell activation enhancement. Moreover, high cell surface expression levels of the MHC class I molecule HLA-C have been shown to correlate with improved HIV control¹⁶, suggesting a critical role for total MHC class I expression in T cell activation. Despite the physiological significance of co-agonism, its molecular mechanism is still incompletely understood. This is largely due to technical challenges of experimentally controlling the amount of presented co-agonist pMHC while keeping agonist pMHC constant.

We developed an experimental system to investigate co-agonism during human CD8⁺ T cell activation by expressing human MHC class I molecules presenting pre-determined peptide as single polypeptide (single chain MHC-I, **Figure 1A**) in a xenogeneic cell line^{9,10}. Single chain (sc) agonist pMHC was expressed under control of tetracycline-inducible promoter in the hamster-derived T-REx CHO cell line expressing tetracycline repressor; this allows very low “leaky” expression of sc agonist pMHC in the absence of tetracycline and high sc agonist pMHC expression after addition of tetracycline (**Figure 1B,C**). The sc agonist pMHC-expressing T-REx CHO cells were then supertransfected with non-stimulatory, co-agonist sc pMHC-I under control of a constitutively active promoter (**Figure 1B,C**). Use of this experimental system allowed us to compare CD8⁺ T cell responses to agonist pMHC in the presence or absence of high levels of non-stimulatory pMHC. Critically, since peptide and MHC-I heavy chain are covalently linked in the scMHC-I format, this allows introductions of mutations into MHC molecules presenting pre-determined peptides. Use of scMHC-I technology thus allowed us to specifically modulate TCR or CD8-binding properties of agonist or co-agonist pMHC.

Co-agonism in CD8⁺ T cell activation has been observed using purified MHC-I complexes presenting agonist and co-agonist peptides in the form of heterodimers^{11,17} or presented on quantum dots^{12,13}. Early work on co-agonism in CD8⁺ T cell activation in the context of agonist pMHC recognition on APC used TAP2-deficient cell lines as APCs. TAP is required for peptide transport from the cytoplasm to the endoplasmic reticulum, and TAP-deficiency severely reduces the pool of available MHC class I-binding peptides. In the absence of peptides, the nascent complex of MHC class I heavy chain and β_2 -microglobulin is unstable, resulting in very low MHC-I cell surface expression. Initially, comparison of T cells stimulated with antigen loaded on TAP-sufficient and -deficient mouse RMA and RMA-S cell lines, respectively, as APCs, did not reveal any evidence for co-agonism during mouse T cell activation¹⁸. However, use of this experimental system did not allow precise control over the amount of agonist pMHC

presented independently of non-stimulatory self pMHC. Moreover, these two cell lines may also differ in expression of other molecules that modulate T cell activation, such as ligands to T cell co-stimulatory or inhibitory receptors, or adhesion molecules.

Subsequently, we developed a protocol for loading TAP2-deficient RMA-S cells with exogenous peptides to allow presentation of a fixed amount of agonist pMHC in the presence or absence of non-stimulatory pMHC. This was achieved through the following: incubation of TAP2-deficient RMA-S cells at lower temperatures (28 °C, instead of the usual 37 °C) to stabilize empty MHC class I heavy chain/ β_2 microglobulin complexes¹⁹; incubation at 28 °C with exogenous agonist peptide; incubation at 28 °C in the presence or absence of exogenous non-stimulatory peptide; and incubation at 37 °C to reduce cell surface expression of empty MHC class I complexes^{8,9}. Use of this experimental set-up revealed that the presence of non-stimulatory pMHC enhanced activation of mouse thymocytes, naïve peripheral CD8⁺ T cells and CTLs⁸. Mechanistically, the presence of non-stimulatory pMHC can induce CD8 recruitment to the T cell:APC immune synapse even in the absence of agonist pMHC, and non-stimulatory pMHC can enhance TCR interaction with CD8 coreceptor^{7,8}. However, this experimental system does not allow testing the relative contributions of TCR and CD8 coreceptor binding to agonist and co-agonist pMHC in mediating the activation enhancement, as any modifications altering TCR or CD8 binding to MHC class I would affect all MHC molecules, regardless of the presented peptide. We therefore used MHC class I single chain technology to overcome this problem.

The single chain (sc) MHC class I format has been used before to improve antigenic pMHC presentation for therapeutic strategies, and to answer fundamental questions on mechanisms of TCR and NK cell receptor activation and function^{20,21}. scMHC-I consists of peptide, β_2 -microglobulin and MHC-I heavy chain joined by two flexible serine/glycine-rich linkers (**Figure 1A**), several different linker sequences have been developed and tested²². scMHC-I can fold independently of the TAP complex and chaperone proteins required for conventional pMHC complex assembly²⁰. scMHC-I have been shown to fold correctly, as determined using conformation-specific antibody staining²² and by solving sc structure with X-ray crystallography²³. The basic scMHC-I design has been modified to improve binding of low affinity peptides^{24,25}.

This protocol describes the use of human scMHC-I to investigate co-agonism during human CTL clone activation. The protocol has been optimized for human CTL clone specific for an epitope of Hepatitis B virus (HBV). HBV is a severe healthcare burden worldwide, as there are 350 million people infected with HBV and chronic HBV infection can lead to development of hepatocellular carcinoma (HCC). HCC cells resulting from HBV infection can usually present HBV-derived epitopes and can be recognized by specific human T cells. HBV and HCC clearance relies on the human T cell responses²⁶, but HCC patients often suffer from T cell exhaustion and reduced T cell responses²⁷. Introduction of HBV-specific TCR into T cells of HBV has been tried as a potential immunotherapy strategy to eliminate chronic HBV infection²⁸. However, it is unknown if this method will be effective in a clinical setting, due to the low levels of surface MHC-I in human hepatocytes²⁹. Therefore, understanding the mechanism of coagonism may provide alternative perspectives for immunotherapy of HBV, possibly by enhancing the

presentation of endogenous pMHC to induce co-agonism-mediated T cell responses. The protocol covers all the necessary steps for research on human co-agonism: transfection of T-REx CHO cells with agonist and co-agonist sc pMHC, generation of stable T-REx CHO cell lines, culture of HBV-specific human CTL CD8⁺ T cell line and CTL activation experiments quantifying cytokine production and degranulation in response to T-REx CHO cells. Although we focus on using the sc MHC class I technology to investigate co-agonism during HBV T cell activation, it must be noted that the presented methods can easily be adapted for research on other aspects of T cell activation in human T cell systems with known pMHC-I specificity.

This protocol uses a human CD8⁺ CTL line specific for HBV-derived peptide E183-91 (FLLTRILT: "E183")³⁰ presented on HLA-A2. sc HLA molecules consisting of a signal sequence from human β 2-microglobulin, a HLA-A2 binding peptide of interest, a glycine/serine linker, a human β 2-microglobulin, a glycine/serine linker and a A2 heavy chain can be commercially synthesized (**Figure 1A**). Plasmids encoding scA2 constructs with agonist E183 peptide, or coagonist HIV GAG (SLYNTVATL)³¹ peptides are available from the authors upon request. The peptide sequence can be changed using site directed-mutagenesis. A tetracycline-inducible vector pcDNA5/TO was used to express the agonist single chain HLA-A2 with E183 peptide (sc E183-HLA-A2). In the presence of tetracycline repressor, the pcDNA5/TO plasmid allows only very low, "leaky" expression of the protein of interest; high expression can be induced by the addition of tetracycline (**Figure 1B,C**). Use of a tetracycline-inducible expression system allows control over the amount of cell surface agonist pMHC expression. A constitutive expression plasmid pcDNA3.1 was used to express coagonist scA2 with GAG peptide (sc GAG-HLA-A2). The tetracycline-regulated expression (T-REx) CHO cell line (hamster cell line, no endogenous human MHC) was used to express agonist and co-agonist sc-HLA-A2 constructs. This experimental system allows precise control over pMHC presentation (**Figure 1C**): no pMHC presentation (untransfected T-REx), a high amount of co-agonist pMHC presentation (constitutive sc GAG-HLA-A2 expression), a low amount of agonist pMHC presentation (sc E183-HLA-A2 in the absence of tetracycline), a high amount of agonist pMHC presentation (sc E183-HLA-A2 in the presence of tetracycline), a low amount of agonist pMHC presentation and high expression of co-agonist pMHC (sc E183-HLA-A2 in the absence of tetracycline, and constitutive sc GAG-HLA-A2 expression). Use of this experimental system allows precise control of cell surface amounts of agonist and coagonist pMHC.

PROTOCOL:

NOTE: All steps involving use of live, non-fixed cells should be performed in a biosafety hood, and any biohazardous waste should be discarded according to the local health and safety regulations.

1 CHO cell transfection and generation of stable CHO cell lines.

1.1 CHO cell transfection

NOTE: This section describes CHO cell transfection using polyethylenimine (PEI). However, alternative methods can be used, as CHO cells are relatively easy to transfect using commercially available lipid transfection reagents.

1.1.1 Culture CHO cells in complete F12 (cF12, F12 supplemented with 5% FBS, 100 U/mL penicillin, 100 µg/mL streptomycin and 10 µg/mL blasticidin to maintain expression of the tetracycline repressor). Passage the cells every 3-4 days at a 1:10 ratio.

1.1.2 One day before transfection, prepare the CHO cell suspension. Wash CHO cells in the flask using PBS (10 mL for PBS for a T75 flask, 5 mL PBS for a T25 flask). Add 0.05% Trypsin/0.02% EDTA in PBS (2 mL for a T75 flask, 1 mL for a T25 flask) to the flask and incubate for 5-10 min at room temperature (RT). Using a light microscope, check that the cells have detached, and stop trypsinization by adding cF12 to the flask (8 mL for a T75 flask, 4 mL for a T25 flask).

1.1.3 Transfer the CHO cell suspension to a 15 mL tube, and centrifuge at 300-400 x g for 5 min at 4 °C or RT. Remove supernatant, re-suspend in cF12 (10 mL for a T75 flask, 5 mL for a T25 flask) and count using a hemocytometer.

1.1.4 Adjust cell concentration to 60,000-100,000 CHO cells/mL. Add 5 mL of CHO cell suspension (300,000-500,000 CHO cells) per well in a 6 well plate. Incubate overnight at 37 °C, 5% CO₂.

1.1.5 On the day of the transfection, add 400 µL of serum-free F12 media or low serum media to a 1.5 mL tube. Add 3 µg of plasmid to 400 µL of media. Mix by pipetting. Add 6 µg of PEI (6 µL of 1 mg/mL PEI solution) to the media/DNA mix, and vortex immediately for 10 s.

1.1.3 Incubate the media/DNA/PEI mix at RT for 10-15 min. During this incubation, replace CHO cell media in the 6 well plate with 5 mL of fresh cF12 media.

1.1.4 Add the media/DNA/PEI solution to CHO cells. Add the solution dropwise and evenly over the well. Incubate overnight at 37 °C, 5% CO₂.

1.2 Generation of stable CHO cell lines

1.2.1 One day after transfection, remove media from the wells and add 5 mL of cF12 containing the appropriate selection antibiotics to each well. Use 0.3 mg/mL hygromycin for pcDNA5TO-based plasmids, and 1 mg/mL geneticin for pcDNA3-based plasmids.

1.2.2 If the cells reach more than 90% confluency at 24 h post transfection, trypsinize the cells. Wash wells with 1 mL of PBS per well, add 1 mL of 0.05% trypsin/0.02% EDTA in PBS and incubate for 5-10 min at RT. Using a light microscope, check that the cells have detached, and stop trypsinization by adding 3-4 mL of cF12 per well.

1.2.3 Transfer the CHO cell suspension to a 15 mL tube, and centrifuge at 300-400 x *g* for 5 min at 4 °C or RT. Remove the supernatant and re-suspend in 10 mL of cF12. Add 5 mL of cell suspension per well in 6 well plate. Use two wells to maximize yield of the transfected cells.

1.2.4 Incubate the 6 well plates at 37 °C, 5% CO₂. Monitor cell death and the growth of resistant colonies using a light microscope. Remove media and replace it with fresh media containing the appropriate antibiotics after 4-5 days, or when significant cell death is observed.

1.2.5 Once the resistant cells reach more than 70% confluence, trypsinize the cells as described in step 1.2.2 and transfer to a T25 flask for expansion. The antibiotic selection takes 1-2 weeks.

1.2.6 Once the antibiotic-resistant CHO cells in the T25 flask reach 70-80% confluence, perform antibody staining to verify cell surface expression of the constructs of interest. One day before antibody staining, trypsinize the cells as described in step 1.2.2 and adjust the cell concentration to 200,000 cells/mL, and add 1 mL of cell suspension to 12 well plate. For testing of tetracycline-inducible constructs, add 50-60 ng/mL of tetracycline. Incubate overnight at 37 °C, 5% CO₂.

1.2.7 Use cell scraper to detach CHO cells from bottom of the well. This method is preferable to trypsinization, as it reduces the risk of epitope damage due to proteolytic cleavage.

1.2.8 Transfer the 1 mL of media containing the CHO cells to a FACS tube. Spin down at 300-400 x *g*, 4 °C for 5 min. Re-suspend in 100 µL of anti-HLA antibody diluted in FWB. Incubate 30 min at 4 °C on ice. Perform this staining with anti-A2 antibody to detect total A2 cell surface expression levels, and TCR-like antibody specific for A2 in complex with E183 peptide to detect agonist pMHC-I cell surface expression levels.

1.2.9 After verifying MHC expression, FACS-sort the positive cells, and maintain the sorted cells in media with antibiotics to reduce transgene loss. For CHO cells expressing agonist pMHC, single cell sorting is recommended, to ensure comparable agonist pMHC expression with and without supertransfection with co-agonist pMHC. 2-3 weeks are required for CHO cell growth following the single cell sorting.

2) HBV-specific CTL culture

2.1 Preparation of PBMCs as feeders for HBV-specific human CTL re-stimulation

NOTE: For the A2-E183-specific CTL clone, re-stimulation with freshly isolated PBMCs, rather than thawed PBMCs, gives the best results.

CAUTION: Obtain all the necessary ethical approvals before recruiting volunteers for blood donations. For the purpose of this work, PBMC were collected from healthy volunteers under the protocol approved by NUS IRB. Informed consent was obtained from all donors. Follow all

the necessary precautions when working with human blood to reduce risk of transmission of blood borne diseases.

2.1.1 Before starting, set the centrifuge temperature at 19-20 °C and pre-warm the density gradient (e.g., Ficoll) to RT.

2.1.2 Add 15 mL of room temperature density gradient to 50 mL tube. Centrifuge at 400 x *g* for 5 min at 19-20 °C. This is to ensure there is no density gradient on the side of the tube.

2.1.3 Collect 20 mL of blood from a healthy volunteer donor through venipuncture. Add 15 mL of PBS to 20 mL of blood in a 50 mL tube. Mix by pipetting up and down.

2.1.4 Slowly, add the 35 mL of blood diluted in PBS on top of the 15 mL of density gradient from step 2.1.2. Make sure that the blood and the density gradient do not mix. Centrifuge the tubes at 19 °C, 1,200 x *g* for 20 min with the brakes off.

2.1.5 Remove approximately 70% of the plasma layer. Carefully aspirate the layer containing the PBMC (the cloudy layer on top of the clear density gradient layer) using a Pasteur pipette and place it in a new 50 mL tube.

2.1.6 Add PBS to PBMCs to obtain a 50 mL total volume. Centrifuge tube at 4 °C, 300-400 x *g* for 10 min. Carefully remove the supernatant using the vacuum aspirator and sterile glass Pasteur pipette, or by decanting.

2.1.7 Re-suspend the cells in 50 mL of PBS. Centrifuge tube at 4 °C, 300 x *g* for 10 min. Carefully remove the supernatant using the vacuum aspirator and sterile glass Pasteur pipette, or by decanting.

2.1.8 Re-suspend cells in 2 mL of complete human T cell media (serum free media supplemented with 2% human serum). Count the cells using hemocytometer. This protocol gives on average $1-2 \times 10^6$ PBMCs per 1 mL of blood used.

2.2 CTL re-stimulation using PBMCs as feeder cells

2.2.1 Irradiate PBMCs at 30 Gy to inhibit proliferation in response to subsequent stimulation.

CAUTION: Ensure adequate training and compliance with the laboratory safety requirements when using the irradiator.

2.2.2 Re-suspend the irradiated PBMCs in complete human T cell media at 2×10^6 cells/mL. Add cytokines and lectin from *Phaseolus vulgaris* (PHA) for the CTL re-stimulation at 2x concentration into the PBMC suspension. The final cytokine and PHA concentrations are: 20 U/mL recombinant human IL-2, 10 ng/mL recombinant IL-7, 10 ng/mL recombinant human IL-15 and 1.5 µg/mL PHA.

2.2.3 Thaw the frozen CTL vial in a 37 °C water bath. Transfer the thawed cells into a 15 mL tube with 10 mL of complete human T cell media. Spin down at 300-400 x *g* at 4 °C for 5 min. Carefully remove the supernatant using the vacuum aspirator and sterile glass Pasteur pipette, or by decanting.

2.2.4 Re-suspend CTLs in 2 mL of complete human T cell media. Count the cells using hemocytometer. Add complete human T cell media to obtain cell concentration of 10⁶ cells/mL.

2.2.5 In a 24 well plate, add 0.5 mL of irradiated PBMCs (10⁶ cells) with 2x cytokines and PHA (step 2.2.2), and 0.5 mL of CTLs (0.5 x 10⁶ cells) per well. The number of wells used depends on the total number of CTLs or PBMCs available. Incubate at 37 °C, 5% CO₂.

2.2.7 After 4-5 days, when CTL blasts are visible and the media becomes yellow, transfer the cultures into 6 well plates, and supplement the initial 1 mL volume with 4 mL of complete human T cell media with 20 U/mL IL-2, 10 ng/mL IL-7 and 10 ng/mL IL-15 (no PHA).

2.2.8 After 2-3 days, transfer the cultures into a T75 flask, add 40 mL of complete human T cell media with cytokines (step 2.2.2, no PHA) and culture in an upright position.

2.2.9 Maintain CTLs at a concentration of 1-2 x 10⁶ cells/mL by adding fresh media supplemented with cytokines (step 2.2.2, no PHA) at a volume equal to that of the existing culture volume at every other day. CTLs can be used for experiments 7-14 days after re-stimulation.

3) T cell activation experiments

NOTE: A typical experiment requires several different T-RExCHO lines: untransfected T-REx CHO cells (negative control), T-REx CHO cells expressing non-stimulatory pMHC only (scA2-GAG; negative control), T-REx CHO cells expressing low amounts of agonist pMHC-I (scA2-ENV without tetracycline, experimental sample), T-REx CHO cells expressing low amounts of agonist pMHC and high amounts of non-stimulatory co-agonist pMHC (scA2-ENV and scA2-GAG, without tetracycline, experimental sample), and T-REx CHO cells expressing high amounts of agonist pMHC-I scA2-ENV with tetracycline, positive control), as shown on **Figure 1C,D**.

3.1 Plating CHO cells as antigen presenting cells

3.1.1 Culture CHO cells in T75 flask, as described in step 1.1.1. Use cells at 70-90% confluence for best results. One day before the T cell activation experiment, trypsinize the cells as described in section 1.1.2. After trypsinization, transfer the cell suspension into a 15 mL tube, centrifuge at 400 x *g* for 5 min at 4 °C or RT, remove supernatant by pipetting or decanting.

3.1.2 Re-suspend the cell pellet in 5 mL of cF12. Count CHO cells using a hemocytometer.

3.1.3 Adjust the cell concentration to 2×10^5 /mL using cF12. Add 50-60 ng/mL of tetracycline to the agonist only samples to induce high amounts of agonist pMHC-I expression as positive control (**Figure 1C**). CHO cells will settle down in 15 mL tubes, so make sure to mix them by pipetting or vortexing before plating.

3.1.4 For T cell activation assay, add 100 μ L of the CHO cell suspension per well to U-bottom 96 well plate, and use triplicates per CHO cell line. For CHO cell anti-MHC staining, add 1 mL of cell suspension to 12 well plate, and use triplicates per cell line. Incubate overnight at 37 °C, 5% CO₂.

3.2 T cell activation assay

NOTE: This T cell activation assay allows simultaneous quantification of T cell degranulation (CD107a staining, a measure of T cell cytotoxicity) and cytokine production.

3.2.1 The day before the experiment, centrifuge CTLs cells at 400 x *g* at 4 °C for 5 min, count using hemocytometer and re-suspend the cells at 10^6 cells/mL in complete human T cell media without cytokines or PHA. Incubate at 37 °C, 5% CO₂. The purpose of this rest step is to reduce CTL sensitivity, as we have observed maximum activation in response to low amounts of antigen without this step.

3.2.2 On the day of the experiment, centrifuge T cells at 400 x *g* at 4 °C for 5 min, remove the supernatant by decanting or pipetting, and re-suspend the cells in 2 mL complete serum-free media (without cytokines or PHA). Count live T cells using a hemocytometer, and adjust T cell concentration to 10^6 live cells/mL using complete human T cell media (without cytokines or PHA).

3.2.3 Add anti-human CD107a antibody at a 1:100 dilution (final concentration of 2.5 μ g/mL), add brefeldin A at a 1:1000 dilution (final concentration of 1 μ g/mL).

3.2.4 Remove media from the CHO cells in the 96 well plate by aspirating with glass Pasteur pipette or by flicking the plate. Per well, add 200 μ L of the T cell suspension (0.2×10^6 cells) containing anti-CD107a and brefeldin A. Co-culture T cells with CHO cells for 3-4 h.

3.2.5 At the end of the co-culture, perform T cell surface antigen staining. Spin down the 96 well plate at 300-400 x *g*, 5 min at 4 °C, and remove the supernatant by aspirating or flicking. Add 2.5 μ L of CD3 and 2.5 μ L of CD8 antibodies in 50 μ L of 0.5% BSA in PBS (flow wash buffer, FWB) per well for cell surface antigen staining. Incubate samples on ice for 30 min in the dark.

3.2.6 Wash samples by adding 150 μ L of FWB per well. Spin down the 96 well plate at 300-400 x *g*, 5 min at 4 °C, and remove the supernatant by aspirating or flicking.

3.2.7 Perform fixation/permeabilization steps using a fixation/permeabilization kit. Add 100 μ L of fixation/permeabilization solution (from the fixation/permeabilization kit) per well, mix by pipetting. Incubate on ice for 20 min.

3.2.8 Centrifuge at 300-400 x *g* and 4 °C for 5 min, and remove the supernatant by aspirating or flicking. Add 200 μ L of 1x perm/wash buffer (dilute the 10x perm/wash stock buffer from the kit prior to use) per well. Repeat this step.

3.2.9 Re-suspend the cell pellet in 50 μ L of 1x perm/wash buffer containing anti-IFN- γ antibody at 0.3 μ g/mL. Incubate on ice for 30 min in dark.

3.2.10 Add 150 μ L of 1x perm/wash buffer. Centrifuge at 300-400 x *g* for at 4 °C for 5 min, and remove the supernatant by aspirating or flicking. Add 200 μ L of 1x perm/wash buffer, centrifuge 300-400 x *g* for at 4 °C for 5 min, and remove the supernatant by aspirating or flicking.

3.2.11 Re-suspend each sample in 200 μ L of FWB. Proceed with flow cytometric analysis.

3.3 CHO cell MHC class I staining

NOTE: It is critical to verify scMHC-I cell surface levels on cells used on each experiment, as cell lines may lose transgene expression. When preparing 96 well plate for T cell stimulation, set up 12 well plate for CHO cell staining, as outlined in Section 3.1.4.

3.3.1 Use a cell scraper to detach CHO cells from the bottom of the well. This method is preferable to trypsinization, as it reduces the risk of epitope damage due to proteolytic cleavage.

3.3.2 Transfer the 1 mL of media containing the CHO cells to the FACS tube. Spin down at 300-400 x *g* at 4 °C for 5 min. Re-suspend in 100 μ L of anti-HLA antibody diluted in FWB. Incubate for 30 min at 4 °C on ice.

3.3.2.1. Perform this staining with anti-A2 antibody to detect total A2 cell surface expression levels, and TCR-like antibody specific for A2 in complex with E183 peptide³² to detect agonist pMHC-I cell surface expression levels.

3.3.3 Add 1 mL of FWB to each CHO sample, spin down at 300-400 x *g*, 4 °C, 5 min, and then re-suspend the sample in 0.3 mL of FWB. Proceed with flow cytometry analysis.

REPRESENTATIVE RESULTS:

The sc MHC class I technology used here allows precise control of cell surface expression of MHC class I with known, covalently linked peptides, to induce T cell activation. We have generated an engineered xenogeneic APC system (**Figure 1A,B,C**) that allows presentation of low levels of agonist pMHC in the presence or absence of co-agonist pMHC (**Figure 1D,E**).

Critically, we have used a TCR-like antibody specific for HLA-A2 presenting E183 peptide to show that presentation of the agonist sc E183-HLA-A2 is not altered by co-expression of co-agonist sc GAG-HLA-A2 (**Figure 1E**). However, it must be noted that the E183-HLA-A2 specific TCR-like antibody shows some binding to HLA-A2 presenting GAG peptide (**Figure 1E**). Similar engineered antigen presenting cell systems can be constructed using different peptides or MHC molecules to investigate molecular requirements for TCR and/or coreceptor interactions in human T cells with different specificities.

We used these engineered APCs to stimulate an E183-HLA-A2-specific human CD8⁺ T cell clone. Three negative controls were used: T cells only, untransfected T-REx CHO cells, and T-REx CHO cells expressing high amounts of co-agonist sc GAG-HLA-A2 to test for potential T cell activation by any molecules expressed on T-REx CHO cells (untransfected T-REx CHO cells as compared to T cell only), and for T cell activation by sc-GAG-HLA-A2 (T-REx CHO cells expressing high levels of sc-GAG-HLA-A2 as compared to untransfected T-REx CHO cells). Untransfected T-REx CHO cells or T-REx CHO cells expressing sc-GAG-HLA-A2 did not induce activation of E183-HLA-A2-specific CD8⁺ CTL clone, as determined by quantifying IFN- γ production and expression of degranulation marker CD107a as a measure of T cell cytotoxicity (**Figure 2A,B**). Expression of high levels of agonist sc E183-HLA-A2, used as positive control, induced very efficient IFN- γ production and degranulation (**Figure 2A,B**), demonstrating that the sc HLA-A2 construct can be recognized by a specific TCR to activate T cell effector functions. Expression of low levels of sc E183-HLA-A2 induced lower levels of IFN- γ production and degranulation, and this was enhanced by presence of co-agonist sc GAG-HLA-A2 (**Figure 2A,B**). It must be noted that very high percentages of CD107a⁺ T cells were observed even in response to low levels of antigenic pMHC (**Figure 2B**), consistent with previous reports of relatively low requirements for the amount of agonist pMHC for induction of T cell cytotoxicity³³. The CD107a MFI, indicating the amount of degranulation on a single cell basis, provides a better dynamic range (**Figure 2B**). The T cell activation assay used allows simultaneous quantification of two major effector functions: cytokine production and cytotoxicity, and provides very sensitive readout with good dynamic range.

We used the single chain technology to test whether TCR binding to co-agonist pMHC is required for co-agonist pMHC-dependent activation enhancement. We mutated valine at position 152 at the edge of the HLA-A2 peptide binding groove (**Figure 3B**) to glutamic acid to create V152E mutant, as this mutation has been reported to abolish TCR binding to HLA-A2³⁴. We then tested if the V152E mutation can abolish TCR binding for the E183-HLA-A2 CTL clone used, by introducing this mutation into agonist sc E183-HLA-A2 and expressing the mutant agonist sc construct in T-REx CHO cells. Wild type, but not V152E, sc E183-HLA-A2 induced activation of the E183-HLA-A2-specific CTL clone used (**Figure 3A**). It must be noted that any TCR-binding mutation on MHC class I needs to be separately tested for each individual TCR/T cell clone used, as the effect of mutation will depend on the precise molecular interactions between TCR and pMHC complex, and these will differ between different TCRs. For example, we tested eight mutations, either previously reported TCR-binding mutations or mutations predicted to disrupt TCR binding without abrogating peptide binding into HLA-A2. Of these, V152E was the only mutation that abolished activation of the E183-specific T cell clone used¹⁰.

We then introduced the V152E mutation into the co-agonist sc GAG-HLA-A2, and co-expressed WT or V152E sc GAG-HLA-A2 with low levels of agonist sc E183-HLA-A2. Both WT and V152E sc GAG-HLA-A2 enhanced IFN- γ production and degranulation (**Figure 3C,D**), suggesting that TCR binding to co-agonist pMHC is not required for co-agonist pMHC-dependent CD8⁺ T cell activation enhancement. Single chain MHC class I technology, as presented here, can be used to modify TCR and/or coreceptor binding to MHC class I with known peptide to investigate the molecular requirements for T cell antigen recognition and activation.

FIGURE AND TABLE LEGENDS:

Figure 1. Single chain HLA-A2 constructs used to investigate co-agonist in human CD8⁺ T cell activation.

(A) Schematic diagram of a scHLA-A2 construct, consisting of covalently signal sequence from fused human β_2 microglobulin, peptide of choice, flexible glycine-serine linker (GGGGSGGGGS GGGGS), human β_2 microglobulin, flexible glycine-serine linker and HLA-A2 heavy chain.

(B) Schematic diagram of the plasmids used to generate engineered antigen presenting cells to investigate co-agonism in human CD8⁺ T cell activation: tetracycline repressor, agonist sc E183-HLA-A2 under control of tetracycline-inducible promoter, non-stimulatory co-agonist sc GAG-HLA-A2 under control of constitutively active promoter.

(C) Schematic diagram of the engineered antigen presenting cells used to investigate co-agonism: T-REx CHO cells (no human MHC expression), T-REx CHO cells expressing constitutively high levels of sc GAG-HLA-A2 (non-stimulatory co-agonist pMHC), T-REx CHO cells expressing low “leaky” levels of sc E183-HLA-A2 in the absence of tetracycline (low level of agonist peptide), T-REx CHO cells expressing high levels of sc E183-HLA-A2 in the presence of tetracycline (high level of agonist peptide), T-REx CHO cells expressing low levels of sc E183-HLA-A2 and high levels of sc GAG-HLA-A2 (low level of agonist peptide and high level of co-agonist peptide).

(D) Cell surface staining of the engineered antigen presenting cell lines using anti-HLA-A2 antibody. The numbers shown denote MFI values for MHC stain in the live cell gate. Data representative of 5 independent experiments (E) Cell surface staining of the engineered antigen presenting cell lines using TCR-like antibody specific against HLA-A2 presenting E183 peptide. The numbers shown denote MFI values for MHC stain in the live cell gate. Data representative of 5 independent experiments.

Figure 2. Presence of co-agonist pMHC enhanced cytokine production and degranulation of E183-HLA-A2-specific human CD8⁺ CTL clone.

(A) Intracellular IFN- γ staining of HLA-A2-specific human CD8⁺ CTL clone after 3h stimulation with the indicated engineered antigen presenting cells. The numbers shown denote the percentage positive for IFN- γ staining. Data representative of 3 independent experiments, each performed with technical triplicates.

(B) Cell surface CD107a staining of HLA-A2-specific human CD8⁺ CTL clone after 3h stimulation with the indicated engineered antigen presenting cells. The numbers shown denote the percentage positive for CD107a staining or the CD107a MFI for the CD8⁺ T cell population. Data representative of 3 independent experiments, each performed with technical triplicates.

Figure 3. Co-agonist – mediated T cell activation enhancement does not require TCR binding to the co-agonist pMHC complex.

(A) V152E mutation in sc E183-HLA-A2 abolishes activation of the E183-HLA-A2-specific human CD8+ CTL clone. T-REx CHO cells expressing high levels (tetracycline induction) of the WT or V152E mutant sc E183-HLA-A2 were used to stimulate E183-HLA-A2-specific human CD8+ CTL clone for 3h. T cell activation was quantified using intracellular IFN- γ staining. The numbers shown denote the percentage positive for IFN- γ staining. Data representative of 3 independent experiments, each performed with technical triplicates.

(B) Location of V152E mutation within the peptide binding groove of HLA-A2.

(C) V152E mutation on sc GAG-HLA-A2 co-agonist pMHC complex does not abolish co-agonist-dependent activation enhancement. Intracellular IFN- γ staining of HLA-A2-specific human CD8+ CTL clone after 4h stimulation with the indicated engineered antigen presenting cells. The numbers shown denote the percentage positive for IFN- γ staining. Data representative of 3 independent experiments, each performed with technical triplicates.

(D) V152E mutation on sc GAG-HLA-A2 co-agonist pMHC complex does not abolish co-agonist-dependent activation enhancement. Cell surface CD107a staining of HLA-A2-specific human CD8+ CTL clone after 3h stimulation with the indicated engineered antigen presenting cells. The numbers shown denote the percentage positive for CD107a staining or the CD107a MFI for the CD8+ T cell population. Data representative of 3 independent experiments, each performed with technical triplicates.

DISCUSSION

The presented protocol provides a robust tool for investigation of molecular interactions during human CD8+ T cell activation. A critical step for investigation of co-agonism during human T cell activation is control over agonist pMHC cell surface expression to ensure equal agonist pMHC presentation on the APCs with or without co-agonist pMHC expression. In our experimental system, this is achieved by using a single cell clone expressing low amounts of agonist pMHC, followed by supertransfection with a co-agonist pMHC construct and agonist pMHC quantification using TCR-like antibody^{10,32}. As agonist sc pMHC expression can change over time, it is critical to quantify agonist pMHC surface expression on APCs used in each experiment using TCR-like antibodies. If no specific TCR-like antibody is available, agonist scMHC-I can be expressed as fluorescent protein fusion, and its cell surface expression can then be quantified using a sensitive microscopy method, such as total internal reflection fluorescence microscopy^{35,36}. Moreover, cell surface expression of co-agonist pMHC decreases over time, due to methylation-dependent and –independent silencing of the CMV promoter³⁷, and should be monitored using antibody staining and flow cytometry analysis, with repeated FACS-sorting when needed. We have cultured CHO cells for up to 5 weeks with relatively limited loss of sc pMHC expression. We highly recommend freezing the CHO cells soon after selection/sorting to ensure a reliable stock of cells with known sc MHC expression.

The several steps of the presented protocols can be modified, depending on the availability of equipment, or depending on the specific research aims. For example, PBMC proliferation potential can be inhibited (Step 3.2.1) using alternative methods that do not require a gamma

(or X-) irradiator, such as treatment with mitomycin C³⁸. A non-enzymatic cell dissociation buffers containing metal chelators³⁹ can be used instead of cell scraping in step 3.3.1. Moreover, the antigen presenting system can be modified to make it more suitable for human CTL clones expressing different TCRs. CHO cells express hamster MHC class I molecules, and although these are irrelevant for the CTL line studied here (**Figure 2A** and **Figure 3A**, we observed no T cell responses to untransfected CHO cells), there is a possibility that other human TCRs may show some xenoreactivity. In that case, hamster MHC class I expression can be reduced by knocking out hamster β 2-microglobulin or TAP using CRISPR/Cas9; or a human APC line can be used after CRISPR/Cas9-mediated β 2-microglobulin or TAP knock out.

The experimental system presented here allows precise control of TCR and CD8 interactions with MHC molecules presenting pre-defined peptides. For example, we have previously shown that mutations abolishing CD8 binding⁴⁰ to coagonist pMHC eliminated coagonist-mediated activation enhancement in murine and human CD8+ T cells^{9,10}, whereas abrogating TCR interaction with co-agonist pMHC had no effect on coagonist-mediated activation enhancement¹⁰. However, the use of single chain MHC constructs has several limitations. For example, it is time and labor-consuming to test multiple co-agonist peptide sequences using this experimental system, as this involves generation of multiple plasmids encoding sc MHC with different peptides, and subsequent transfections and cell sorting. Previously, we used the murine TAP-deficient cell line RMA-S to test multiple co-agonist peptides in mouse T cell activation^{8,9}, but this approach was not successful with the TAP-deficient human T2 cell line¹⁰, most likely due to presentation of TAP-independent peptides⁴¹. Another limitation of our experimental system is that it does not provide a fast method for altering the amount of co-agonist pMHC molecules for the purpose of determining the critical amount of co-agonist pMHC required to trigger T cell activation. Moreover, the tetracycline-inducible system used does not allow titration of agonist pMHC quantity at the single cell level by changing the tetracycline concentration, as varying tetracycline concentration alters the proportion of HLA-A2 positive cells, rather than HLA-A2 levels on per cell basis. An alternative approach that we are currently investigating is to use supported lipid bilayers⁴², previously used to investigate co-agonism during murine CD4+ T cell activation⁴ or beads presenting fixed amount of agonist pMHC in the presence of varying amounts of co-agonist pMHC. When used in conjunction with UV-cleavable peptide technology^{4,43}, these alternative experimental systems should allow fast testing of multiple co-agonist peptide sequences as well as different amounts of co-agonist pMHC.

The sc MHC technology can be also applicable to investigation of co-agonism in human or murine CD4 T cells, as single chain MHC class II molecules presenting pre-determined peptides have been successfully expressed on mammalian cell surfaces⁴⁴. Moreover, the use of sc MHC technology can be extended to investigation of non-classical MHC molecules such as HLA-E. Sc HLA-E has been described before, and its expression has been shown to inhibit human T cell activation⁴⁵. Another non-classical MHC molecule of interest is CD1, which presents lipids instead of peptides⁴⁶. Sc constructs consisting of CD1 heavy chain and β 2-microglobulin have been shown to be expressed on the cell surface and to bind relevant lipid ligands⁴⁷. Single chain

MHC technology has great potential as a research tool for investigating molecular interactions during T and NK cell activation.

ACKNOWLEDGMENTS:

This research was supported by the Singapore Ministry of Health's National Medical Research Council under its CBRG/0064/2014 to N.R.J.G., by CREATE under its R571-002-012-592 to P.A.M., National Research Foundation Investigatorship R571-000-272-281 to P.A.M., and by a Singapore Translational Research (STaR) Investigator Award (NMRC/STaR/013/2012) to A.B. We are grateful to Dr. Paul Hutchinson and Mr. Teo Guo Hui from NUS Immunology Programme Flow Core facility for the expert cell sorting. We are grateful to Elijah Chen for critical reading of the manuscript.

DISCLOSURES:

The authors declare no competing interests.

REFERENCES:

- 1 Gascoigne, N. R. Do T cells need endogenous peptides for activation? *Nature Reviews. Immunology*. **8** (11), 895-900, doi:10.1038/nri2431 (2008).
- 2 Gascoigne, N. R., Zal, T., Yachi, P. P. & Hoerter, J. A. Co-receptors and recognition of self at the immunological synapse. *Current Topics in Microbiology and Immunology*. **340** 171-189, doi:10.1007/978-3-642-03858-7_9 (2010).
- 3 Krogsgaard, M., Juang, J. & Davis, M. M. A role for "self" in T-cell activation. *Seminars in Immunology*. **19** (4), 236-244, doi:10.1016/j.smim.2007.04.003 (2007).
- 4 Krogsgaard, M. et al. Agonist/endogenous peptide-MHC heterodimers drive T cell activation and sensitivity. *Nature*. **434** (7030), 238-243, doi:10.1038/nature03391 (2005).
- 5 Ebert, P. J., Jiang, S., Xie, J., Li, Q. J. & Davis, M. M. An endogenous positively selecting peptide enhances mature T cell responses and becomes an autoantigen in the absence of microRNA miR-181a. *Nature Immunology*. **10** (11), 1162-1169, doi:10.1038/ni.1797 (2009).
- 6 Lo, W. L. et al. An endogenous peptide positively selects and augments the activation and survival of peripheral CD4+ T cells. *Nature Immunology*. **10** (11), 1155-1161, doi:10.1038/ni.1796 (2009).
- 7 Yachi, P. P., Ampudia, J., Gascoigne, N. R. & Zal, T. Nonstimulatory peptides contribute to antigen-induced CD8-T cell receptor interaction at the immunological synapse. *Nature Immunology*. **6** (8), 785-792, doi:10.1038/ni1220 (2005).
- 8 Yachi, P. P., Lotz, C., Ampudia, J. & Gascoigne, N. R. T cell activation enhancement by endogenous pMHC acts for both weak and strong agonists but varies with differentiation state. *The Journal of Experimental Medicine*. **204** (11), 2747-2757, doi:10.1084/jem.20062610 (2007).
- 9 Hoerter, J. A. et al. Coreceptor affinity for MHC defines peptide specificity requirements for TCR interaction with coagonist peptide-MHC. *The Journal of Experimental Medicine*. **210** (9), 1807-1821, doi:10.1084/jem.20122528 (2013).
- 10 Zhao, X. et al. Nonstimulatory peptide-MHC enhances human T-cell antigen-specific

responses by amplifying proximal TCR signaling. *Nature Communications*. **9** (1), 2716, doi:10.1038/s41467-018-05288-0 (2018).

11 Juang, J. et al. Peptide-MHC heterodimers show that thymic positive selection requires a more restricted set of self-peptides than negative selection. *The Journal of Experimental Medicine*. **207** (6), 1223-1234, doi:10.1084/jem.20092170 (2010).

12 Anikeeva, N. et al. Quantum dot/peptide-MHC biosensors reveal strong CD8-dependent cooperation between self and viral antigens that augment the T cell response. *Proceedings of the National Academy of Sciences of the United States of America*. **103** (45), 16846-16851, doi:10.1073/pnas.0607771103 (2006).

13 Anikeeva, N., Gakamsky, D., Scholler, J. & Sykulev, Y. Evidence that the density of self peptide-MHC ligands regulates T-cell receptor signaling. *PloS One*. **7** (8), e41466, doi:10.1371/journal.pone.0041466 (2012).

14 Schuren, A. B., Costa, A. I. & Wiertz, E. J. Recent advances in viral evasion of the MHC Class I processing pathway. *Current Opinion in Immunology*. **40** 43-50, doi:10.1016/j.coi.2016.02.007 (2016).

15 Bubenik, J. MHC class I down regulation, tumour escape from immune surveillance and design of therapeutic strategies. *Folia Biologica*. **51** (1), 1-2 (2005).

16 Apps, R. et al. Influence of HLA-C expression level on HIV control. *Science (New York, N.Y.)*. **340** (6128), 87-91, doi:10.1126/science.1232685 (2013).

17 Cebecauer, M. et al. CD8+ cytotoxic T lymphocyte activation by soluble major histocompatibility complex-peptide dimers. *The Journal of Biological Chemistry*. **280** (25), 23820-23828, doi:10.1074/jbc.M500654200 (2005).

18 Sporri, R. & Reis e Sousa, C. Self peptide/MHC class I complexes have a negligible effect on the response of some CD8+ T cells to foreign antigen. *European Journal of Immunology*. **32** (11), 3161-3170, doi:10.1002/1521-4141(200211)32:11<3161::aid-immu3161>3.0.co;2-6 (2002).

19 Ljunggren, H. G. et al. Empty MHC class I molecules come out in the cold. *Nature*. **346** (6283), 476-480, doi:10.1038/346476a0 (1990).

20 Kotsiou, E., Brzostek, J. & Gould, K. G. Properties and applications of single-chain major histocompatibility complex class I molecules. *Antioxidants & Redox Signaling*. **15** (3), 645-655, doi:10.1089/ars.2010.3694 (2011).

21 Hansen, T. H., Connolly, J. M., Gould, K. G. & Fremont, D. H. Basic and translational applications of engineered MHC class I proteins. *Trends in Immunology*. **31** (10), 363-369, doi:10.1016/j.it.2010.07.003 (2010).

22 Yu, Y. Y., Netuschil, N., Lybarger, L., Connolly, J. M. & Hansen, T. H. Cutting edge: single-chain trimers of MHC class I molecules form stable structures that potently stimulate antigen-specific T cells and B cells. *Journal of Immunology (Baltimore, Md. : 1950)*. **168** (7), 3145-3149 (2002).

23 Mitaksov, V. et al. Structural engineering of pMHC reagents for T cell vaccines and diagnostics. *Chemistry & Biology*. **14** (8), 909-922, doi:10.1016/j.chembiol.2007.07.010 (2007).

24 Truscott, S. M. et al. Disulfide bond engineering to trap peptides in the MHC class I binding groove. *Journal of Immunology (Baltimore, Md. : 1950)*. **178** (10), 6280-6289 (2007).

743 25 Truscott, S. M. et al. Human major histocompatibility complex (MHC) class I molecules
744 with disulfide traps secure disease-related antigenic peptides and exclude competitor
745 peptides. *The Journal of Biological Chemistry*. **283** (12), 7480-7490,
746 doi:10.1074/jbc.M709935200 (2008).

747 26 Bertoletti, A. & Ferrari, C. Innate and adaptive immune responses in chronic hepatitis B
748 virus infections: towards restoration of immune control of viral infection. *Postgraduate*
749 *Medical Journal*. **89** (1051), 294-304, doi:10.1136/postgradmedj-2011-301073rep
750 (2013).

751 27 Boni, C. et al. Characterization of hepatitis B virus (HBV)-specific T-cell dysfunction in
752 chronic HBV infection. *Journal of Virology*. **81** (8), 4215-4225, doi:10.1128/jvi.02844-06
753 (2007).

754 28 Gehring, A. J. et al. Engineering virus-specific T cells that target HBV infected
755 hepatocytes and hepatocellular carcinoma cell lines. *Journal of Hepatology*. **55** (1), 103-
756 110, doi:10.1016/j.jhep.2010.10.025 (2011).

757 29 Gehring, A. J. et al. The level of viral antigen presented by hepatocytes influences CD8 T-
758 cell function. *Journal of Virology*. **81** (6), 2940-2949, doi:10.1128/jvi.02415-06 (2007).

759 30 Bertoni, R. et al. Human histocompatibility leukocyte antigen-binding supermotifs
760 predict broadly cross-reactive cytotoxic T lymphocyte responses in patients with acute
761 hepatitis. *The Journal of Clinical Investigation*. **100** (3), 503-513, doi:10.1172/jci119559
762 (1997).

763 31 Goulder, P. J. et al. Patterns of immunodominance in HIV-1-specific cytotoxic T
764 lymphocyte responses in two human histocompatibility leukocyte antigens (HLA)-
765 identical siblings with HLA-A*0201 are influenced by epitope mutation. *The Journal of*
766 *Experimental Medicine*. **185** (8), 1423-1433 (1997).

767 32 Sastry, K. S. et al. Targeting hepatitis B virus-infected cells with a T-cell receptor-like
768 antibody. *Journal of Virology*. **85** (5), 1935-1942, doi:10.1128/jvi.01990-10 (2011).

769 33 van den Berg, H. A. et al. Cellular-level versus receptor-level response threshold
770 hierarchies in T-cell activation. *Frontiers in Immunology*. **4** 250,
771 doi:10.3389/fimmu.2013.00250 (2013).

772 34 Moots, R. J., Matsui, M., Pazmany, L., McMichael, A. J. & Frelinger, J. A. A cluster of
773 mutations in HLA-A2 alpha 2 helix abolishes peptide recognition by T cells.
774 *Immunogenetics*. **34** (3), 141-148 (1991).

775 35 Crites, T. J., Chen, L. & Varma, R. A TIRF microscopy technique for real-time,
776 simultaneous imaging of the TCR and its associated signaling proteins. *Journal of*
777 *Visualized Experiments*. (61), doi:10.3791/3892 (2012).

778 36 Axmann, M., Schutz, G. J. & Huppa, J. B. Single Molecule Fluorescence Microscopy on
779 Planar Supported Bilayers. *Journal of Visualized Experiments*. (105), e53158,
780 doi:10.3791/53158 (2015).

781 37 Mariati, Koh, E. Y., Yeo, J. H., Ho, S. C. & Yang, Y. Toward stable gene expression in CHO
782 cells. *Bioengineered*. **5** (5), 340-345, doi:10.4161/bioe.32111 (2014).

783 38 Matsushita, S., Tanaka, Y., Matsuoka, T. & Nakashima, T. Clonal expansion of freshly
784 isolated CD4T cells by randomized peptides and identification of peptide ligands using
785 combinatorial peptide libraries. *European Journal of Immunology*. **31** (8), 2395-2402,
786 doi:10.1002/1521-4141(200108)31:8<2395::aid-immu2395>3.0.co;2-0 (2001).

787 39 Heng, B. C., Cowan, C. M. & Basu, S. Comparison of enzymatic and non-enzymatic means
788 of dissociating adherent monolayers of mesenchymal stem cells. *Biological Procedures*
789 *Online*. **11** 161-169, doi:10.1007/s12575-009-9001-4 (2009).

790 40 Smith, K. et al. Sensory adaptation in naive peripheral CD4 T cells. *The Journal of*
791 *Experimental Medicine*. **194** (9), 1253-1261 (2001).

792 41 Wei, M. L. & Cresswell, P. HLA-A2 molecules in an antigen-processing mutant cell
793 contain signal sequence-derived peptides. *Nature*. **356** (6368), 443-446,
794 doi:10.1038/356443a0 (1992).

795 42 Vardhana, S. & Dustin, M. Supported Planar Bilayers for the Formation of Study of
796 Immunological Synapses and Kinapse. *Journal of Visualized Experiments*. (19),
797 doi:10.3791/947 (2008).

798 43 Bakker, A. H. et al. Conditional MHC class I ligands and peptide exchange technology for
799 the human MHC gene products HLA-A1, -A3, -A11, and -B7. *Proceedings of the National*
800 *Academy of Sciences of the United States of America*. **105** (10), 3825-3830,
801 doi:10.1073/pnas.0709717105 (2008).

802 44 Rhode, P. R. et al. Single-chain MHC class II molecules induce T cell activation and
803 apoptosis. *Journal of Immunology (Baltimore, Md. : 1950)*. **157** (11), 4885-4891 (1996).


804 45 Crew, M. D., Cannon, M. J., Phanavanh, B. & Garcia-Borges, C. N. An HLA-E single chain
805 trimer inhibits human NK cell reactivity towards porcine cells. *Molecular Immunology*.
806 **42** (10), 1205-1214, doi:10.1016/j.molimm.2004.11.013 (2005).

807 46 Van Kaer, L., Wu, L. & Joyce, S. Mechanisms and Consequences of Antigen Presentation
808 by CD1. *Trends in Immunology*. **37** (11), 738-754, doi:10.1016/j.it.2016.08.011 (2016).


809 47 Im, J. S. et al. Direct measurement of antigen binding properties of CD1 proteins using
810 fluorescent lipid probes. *The Journal of Biological Chemistry*. **279** (1), 299-310,
811 doi:10.1074/jbc.M308803200 (2004).

812


A




B




C









D



E



A**B**



Name of Material/ Equipment	Company	Catalog Number	Comments/Description
Aim-V (serum free media)	Gibco	12055091	
blasticidin	Gibco	R21001	
Cytofix/Cytoperm (fixation/permealibisation kit)	BD	554714	
F12	Gibco	10565-018	
Fetal bovine serum (FBS)	HyClone	SH3007103	
geneticin	Gibco	10131035	
GolgiPlug (Brefeldin A)	BD	555029	
Human serum	Sigma-Aldrich	H4522	
hygromycin	Gibco	10687010	
Lectin from Phaseolus vulgaris (PHA)	Sigma-Aldrich	L4144	
Opti-MEM (low serum media)	Bibco	31985070	
10 X PBS	Vivantis	PB0344 – 1L	
Penicillin/Streptomycin	HyClone	SV30010	
polyethylenimine (PEI)	Sigma-Aldrich	408727	
recombinant human IL2	R&D Systems	202-IL	
recombinant human IL7	R&D Systems	207-IL	
recombinant human IL15	R&D Systems	247-ILB	
tetracycline	Sigma-Aldrich	T7660	
T-REx CHO cell line	Thermo Fisher Scientific	R71807	
10 x Trypsin/EDTA	Gibco	15400054	
Antibodies			
CD3	BD	347347	Clone SK7, RRID:AB_400287
CD8 α	BD	563795	Clone RPA-T8, RRID:AB_2722501
CD107a	BD	566261	Clone H4A3, RRID:AB_2739639
HLA-A2	eBioscience	17-9876-41	Clone BB7.2, RRID:AB_11151522
IFN- γ	eBioscience	12-7319-41	Clone 4S.B3, RRID:AB_1311250



1 Alewife Center #200
Cambridge, MA 02140
tel. 617.945.9051
www.jove.com

ARTICLE AND VIDEO LICENSE AGREEMENT - UK

Title of Article:	Use of single chain MHC technology to investigate co-agonism in human cost T cell activation
Author(s):	Zhao, Hamidiya, Choo, Too, Ho, Ren, Bertolotti, MacAry, Gould, Buostek and Gascoigne

Item 1: The Author elects to have the Materials be made available (as described at <http://www.jove.com/publish>) via:

☐ Standard Access

☒ Open Access

Item 2: Please select one of the following items:

☒ The Author is **NOT** a United States government employee.

☐ The Author is a United States government employee and the Materials were prepared in the course of his or her duties as a United States government employee.

☐ The Author is a United States government employee but the Materials were NOT prepared in the course of his or her duties as a United States government employee.

ARTICLE AND VIDEO LICENSE AGREEMENT

1. **Defined Terms.** As used in this Article and Video License Agreement, the following terms shall have the following meanings: "**Agreement**" means this Article and Video License Agreement; "**Article**" means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein; "**Author**" means the author who is a signatory to this Agreement; "**Collective Work**" means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; "**CRC License**" means the Creative Commons Attribution 3.0 Agreement (also known as CC-BY), the terms and conditions of which can be found at: <http://creativecommons.org/licenses/by/3.0/us/legalcode>; "**Derivative Work**" means a work based upon the Materials or upon the Materials and other pre-existing works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; "**Institution**" means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; "**JoVE**" means MyJoVE Corporation, a Massachusetts corporation and the publisher of The Journal of Visualized Experiments; "**Materials**" means the Article and / or the Video; "**Parties**" means the Author and JoVE; "**Video**" means any video(s) made by the Author, alone or in conjunction with any other parties, or by JoVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion

of the Article, and in which the Author may or may not appear.

2. **Background.** The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JoVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.

3. **Grant of Rights in Article.** In consideration of JoVE agreeing to publish the Article, the Author hereby grants to JoVE, subject to **Sections 4 and 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the "Open Access" box has been checked in **Item 1** above, JoVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.

612542.6 For questions, please contact us at submissions@jove.com or +1.617.945.9051.

4. **Retention of Rights in Article.** Notwithstanding the exclusive license granted to JoVE in **Section 3** above, the Author shall, with respect to the Article, retain the non-exclusive right to use all or part of the Article for the non-commercial purpose of giving lectures, presentations or teaching classes, and to post a copy of the Article on the Institution's website or the Author's personal website, in each case provided that a link to the Article on the JoVE website is provided and notice of JoVE's copyright in the Article is included. All non-copyright intellectual property rights in and to the Article, such as patent rights, shall remain with the Author.

5. **Grant of Rights in Video - Standard Access.** This **Section 5** applies if the "Standard Access" box has been checked in **Item 1** above or if no box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby acknowledges and agrees that, subject to **Section 7** below, JoVE is and shall be the sole and exclusive owner of all rights of any nature, including, without limitation, all copyrights, in and to the Video. To the extent that, by law, the Author is deemed, now or at any time in the future, to have any rights of any nature in or to the Video, the Author hereby disclaims all such rights and transfers all such rights to JoVE.

6. **Grant of Rights in Video - Open Access.** This **Section 6** applies only if the "Open Access" box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby grants to JoVE, subject to **Section 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Video in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Video into other languages, create adaptations, summaries or extracts of the Video or other Derivative Works or Collective Works based on all or any portion of the Video and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats.

7. **Government Employees.** If the Author is a United States government employee and the Article was prepared in the course of his or her duties as a United States government employee, as indicated in **Item 2** above, and any of the licenses or grants granted by the Author hereunder exceed the scope of the 17 U.S.C. 403, then the rights granted hereunder shall be limited to the maximum rights permitted under such statute. In such case, all provisions contained herein that are not in conflict with such statute shall remain in full force and effect, and all provisions contained herein that do so conflict shall be

deemed to be amended so as to provide to JoVE the maximum rights permissible within such statute.

8. **Protection of the work.** The Author(s) authorize JoVE to take steps in the Author(s) name and on their behalf if JoVE believes some third party could be infringing or might infringe the copyright of either the Author's Article and/or Video.

9. **Likeness, Privacy, Personality.** The Author hereby grants JoVE the right to use the Author's name, voice, likeness, picture, photograph, image, biography and performance in any way, commercial or otherwise, in connection with the Materials and the sale, promotion and distribution thereof. The Author hereby waives any and all rights he or she may have, relating to his or her appearance in the Video or otherwise relating to the Materials, under all applicable privacy, likeness, personality or similar laws.

10. **Author Warranties.** The Author represents and warrants that the Article is original, that it has not been published, that the copyright interest is owned by the Author (or, if more than one author is listed at the beginning of this Agreement, by such authors collectively) and has not been assigned, licensed, or otherwise transferred to any other party. The Author represents and warrants that the author(s) listed at the top of this Agreement are the only authors of the Materials. If more than one author is listed at the top of this Agreement and if any such author has not entered into a separate Article and Video License Agreement with JoVE relating to the Materials, the Author represents and warrants that the Author has been authorized by each of the other such authors to execute this Agreement on his or her behalf and to bind him or her with respect to the terms of this Agreement as if each of them had been a party hereto as an Author. The Author warrants that the use, reproduction, distribution, public or private performance or display, and/or modification of all or any portion of the Materials does not and will not violate, infringe and/or misappropriate the patent, trademark, intellectual property or other rights of any third party. The Author represents and warrants that it has and will continue to comply with all government, institutional and other regulations, including, without limitation all institutional, laboratory, hospital, ethical, human and animal treatment, privacy, and all other rules, regulations, laws, procedures or guidelines, applicable to the Materials, and that all research involving human and animal subjects has been approved by the Author's relevant institutional review board.

11. **JoVE Discretion.** If the Author requests the assistance of JoVE in producing the Video in the Author's facility, the Author shall ensure that the presence of JoVE employees, agents or independent contractors is in accordance with the relevant regulations of the Author's institution. If more than one author is listed at the beginning of this Agreement, JoVE may, in its sole discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JoVE reserves the right, in its absolute and sole

discretion and without giving any reason therefore, to accept or decline any work submitted to JoVE. JoVE and its employees, agents and independent contractors shall have full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JoVE has sole discretion as to the method of making and publishing the Materials, including, without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

12. **Indemnification.** The Author agrees to indemnify JoVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JoVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JoVE, making of videos by JoVE, or publication in JoVE or elsewhere by JoVE. The Author shall be responsible for, and shall hold JoVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to the making of a video by JoVE its employees, agents or independent contractors. All sterilization, cleanliness or

decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's expense. All indemnifications provided herein shall include JoVE's attorney's fees and costs related to said losses or damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JoVE, its employees, agents or independent contractors.

13. **Fees.** To cover the cost incurred for publication, JoVE must receive payment before production and publication the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JoVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.

14. **Transfer, Governing Law.** This Agreement may be assigned by JoVE and shall inure to the benefits of any of JoVE's successors and assignees. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to be one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement is required per submission.

CORRESPONDING AUTHOR

Name:

JOANNA BRZOSTEK

Department:

MICROBIOLOGY AND IMMUNOLOGY

Institution:

NATIONAL UNIVERSITY OF SINGAPORE

Title:

SENIOR RESEARCH FELLOW

Signature:

J. Brzostek

Date:

18-9-2018

Please submit a **signed** and **dated** copy of this license by one of the following three methods:

1. Upload an electronic version on the JoVE submission site
2. Fax the document to +1.866.381.2236
3. Mail the document to JoVE / Attn: JoVE Editorial / 1 Alewife Center #200 / Cambridge, MA 02140

We are grateful to the reviewers and the editor for useful suggestions and corrections that have helped us to improve the manuscript. The manuscript has been substantially revised and edited. Our responses to specific points raised can be found below, highlighted in yellow.

Responses to editor's comments:

1. Please take this opportunity to thoroughly proofread the manuscript to ensure that there are no spelling or grammar issues.

The manuscript has been proofread to ensure that there are no spelling or grammar issues.

2. Please provide an email address for each author.

This has been provided.

3. Please use SI abbreviations for all units: L, mL, μ L, h, min, s, etc.

SI abbreviations are used throughout the manuscript text.

4. Lines 165-199, etc.: Please revise the protocol to contain only action items that direct the reader to do something (e.g., "Do this," "Ensure that," etc.). The actions should be described in the imperative tense in complete sentences wherever possible. Avoid usage of phrases such as "could be," "should be," and "would be" throughout the Protocol. Any text that cannot be written in the imperative tense may be added as a "Note." Please include all safety procedures and use of hoods, etc. However, notes should be used sparingly and actions should be described in the imperative tense wherever possible. Please move the discussion about the protocol to the Discussion.

The protocol has been substantially revised, with the descriptive sections 1 and 2 moved to the introduction. Information on safety procedures has been included throughout the manuscript.

5. The Protocol should be made up almost entirely of discrete steps without large paragraphs of text between sections. Please simplify the Protocol so that individual steps contain only 2-3 actions per step and a maximum of 4 sentences per step. Use sub-steps as necessary. Please move the discussion about the protocol to the Discussion.

The protocol has been substantially revised according to these guidelines.

6. Please remove commercial language from the manuscript: Opti-MEM, Cytofix/Cytoperm, etc.

This has been corrected.

7. Line 198: Please provide the composition of the culture media for T-Rex CHO cells.

This has been provided. This section of the protocol has been substantially edited for better clarity.

8. 2.2.3: Please specify what is considered to be high confluence (80%? 90%?). Please describe how to trypsinise the cells and specify the concentration of trypsin and reaction conditions.

This information is provided in the revised version of the manuscript.

9. 2.2.4: Please describe how to perform antibody staining.

This information is provided in the revised version of the manuscript.

10. 4.1.1: Please specify the conditions for culturing CHO cells.

This information is provided in the revised version of the manuscript.

11. 4.1.2: Please split into several substeps.

This section of the protocol has been substantially edited for better clarity.

12. After you have made all the recommended changes to your protocol (listed above), please highlight 2.75 pages or less of the Protocol (including headings and spacing) that identifies the essential steps of the protocol for the video, i.e., the steps that should be visualized to tell the most cohesive story of the Protocol.

The relevant sections are highlighted in yellow.

13. Please highlight complete sentences (not parts of sentences). Please ensure that the highlighted part of the step includes at least one action that is written in imperative tense.

14. Please include all relevant details that are required to perform the step in the highlighting. For example: If step 2.5 is highlighted for filming and the details of how to perform the step are given in steps 2.5.1 and 2.5.2, then the sub-steps where the details are provided must be highlighted.

15. Please note that JoVE requires Representative Results and Discussion to be two separate sections. Representative Results sections should explain the results in the context of the technique you have described, e.g., how do these results show the technique, suggestions about how to analyze the outcome, etc. The paragraph text should refer to all of the figures. However for figures showing the experimental set-up, please reference them in the Protocol. Data from both successful and sub-optimal experiments can be included.

The revised version of the manuscript contains Representative Results and Discussion.

16. JoVE articles are focused on the methods and the protocol, thus the discussion should be similarly focused. Please add a Discussion section to explicitly cover the following in detail in 3-6 paragraphs with citations:

- a) Critical steps within the protocol
- b) Any modifications and troubleshooting of the technique
- c) Any limitations of the technique
- d) The significance with respect to existing methods
- e) Any future applications of the technique

The revised version of the manuscript contains Discussion.

17. References: Please do not abbreviate journal titles.

The references have been re-formatted.

18. Table of Equipment and Materials: Please provide lot numbers and RRIDs of antibodies, if available. Please sort the items in alphabetical order according to the Name of Material/ Equipment.

The table has been re-formatted.

Reviewers' comments:

Reviewer #1:

Manuscript Summary:

* Are the title and abstract appropriate for this methods article?

Yes

* Are there any other potential applications for the method/protocol the authors could discuss?

Yes - the authors could mention that linked peptide technology could also be used to measure coagonism by MHC class II for CD4, and study the potential influence of non-classical MHC class I molecules such as HLA-E or even CD1.

This has been included in the discussion of the revised manuscript.

* Are all the materials and equipment needed listed in the table? (Please note that any basic lab materials or equipment do not need to be listed, e.g. pipettes.)

Yes

* Do you think the steps listed in the procedure would lead to the described outcome?

Yes

* Are the steps listed in the procedure clearly explained?

Yes- note there is an error on line 203, streptomycin is 100µg/ml not 10mg/ml! Also there is a discrepancy between the timing of co-culture in the procedure (4h) and Fig 3 (3h for A, B, and D). Line 355 should be amended to 3-4h to avoid confusion.

Thank you for spotting the mistakes. This has been corrected in the revised manuscript.

* Are any important steps missing from the procedure?

No

* Are appropriate controls suggested?

Yes, although one is absent- specifically mutation of scGAG to abolish CD8 binding.

This has been added to the discussion in the revised manuscript.

* Are all the critical steps highlighted?

Yes

* Is there any additional information that would be useful to include?

Yes, some alternatives for key reagents/instruments would be useful. eg if no gamma irradiator is available another method for inactivating the PBMCs (such as mitomycin C or UV irradiation) could be

given. Similarly, if no specific "TCR-like" antibody is available, an alternative such as using an sc-peptide-HLAB7 construct and mAb BB7.1 could be discussed. Likewise, in section 4.3.1 the use of a non-enzymatic tissue dissociation medium as an alternative to scraping should be discussed. Also, as the CMV promoter is particularly sensitive to methylation dependent inactivation, some indication of how long the clones based on pCDNA3 can be maintained in culture should be given, and/or alternatives such as a constitutive vector with an IRES driven drug selection cassette mentioned. Finally, it should be noted that the CHO cells will still express Hamster MHC class I molecules. Although these appear irrelevant for the system studied, the possibility that other human TCRs might show xenoreactivity should be mentioned, and the use of CRISPR/Cas9 technology to eliminate the hamster beta 2 microglobulin given as an alternative in the event of a background signal is detected.

We are grateful to the reviewer for these excellent suggestions, which have been included in the discussion in the revised manuscript.

* Are the anticipated results reasonable, and if so, are they useful to readers?

Yes

* Are any important references missing and are the included references useful?

Yes, at least one reference discussing mutations to the alpha 3 domain that abolish CD8 binding (eg Purbhoo et al J. Biol Chem Vol. 276, pp. 32786-32792, 2001) should be given. Given that the V152E scGAG was still active as a co-agonist (although there was no control using a GAG-specific T cell showing that this mutation was also effective for this pMHC complex- an important point given the author's comments about the ineffectiveness of 7 other mutations on E183) the use of a D227K/T228A mutant scGAG that would be predicted to eliminate/reduce its co-agonism would also be an important control.

This has been added to the discussion in the revised manuscript.

Minor point - I found the yellow circles in Fig 1C difficult to see - I would recommend filling the "cells" with gray to increase the contrast, or substituting with a different color such as cyan.

Thank you for the suggestion. The figure has been re-formatted.

Reviewer #2:

Manuscript Summary:

This article by Brzostek et al. presents the development of an important protocol to develop single chain (sc) MHC class I to investigate molecular interaction in human CD8+T cell activation. The article focuses on sc MHC class I technology to investigate co-agonism during HBV T cell activation. Though not-stimulatory self-peptide MHC (pMHC) complexes do not induce T cell activation and effector functions, it can enhance T cell responses to agonist pMHC. This 'co-agonism' process may be useful in cases where there is only low level of surface MHC-1 presented in the cell. Therefore understanding the mechanism of co-agonism may provide an alternative perspective for immunotherapy, possibly by enhancing the presentation of endogenous pMHC to induce co-agonism-mediated T cell responses. In addition the method developed in this paper can be easily adapted for research on different aspects of CD8+ T cell activation in human T cell systems with known pMHC-I specificity.

Previous research showed that when a fixed amount of agonist pMHC are allowed in the presence or absence of non-stimulatory pMHC, the presence of non-stimulatory pMHC enhanced activation of mouse thymocytes, naïve peripheral CD8⁺ T cells and CTLs. However, this experimental system does not allow testing the relative contributions of TCR and CD8 coreceptor binding to agonist and co-agonist pMHC in mediating the activation enhancement, as any modifications altering TCR or CD8 binding to MHC class I would affect all MHC molecules, regardless of the presented peptide. To overcome this shortcoming, the authors developed a clever approach by using MHC class I single chain technology (scMHC I). The scMHC-I consists of peptide, β 2-microglobulin and MHC-I heavy chain joined by two flexible linkers. This protocol describes use of human scMHC-I to investigate co-agonism during human CTL clone activation.

They used the tetracycline-regulated expression (T-Rex) CHO cell line to express agonist and co-agonist constructs. The advantage of using this is in the presence of tetracycline repressor, pcDNA5/TO plasmid allows only very low, "leaky" expression of the protein of interest, high expression levels can be induced by addition of tetracycline. This way they precisely control the pMHC presentation in the cell surface to induce T cell activation. This engineered xenogeneic APC system allows presentation of low levels of agonist pMHC in the presence or absence of co-agonist pMHC. They tested potential CD8⁺T cell activation in engineered APCs in samples and controls: three negative controls (T cell only, untransfected T-Rex CHO cells and T-Rex CHO cells expressing high levels of co-agonist), two samples (low level of agonist in the absence and in the presence of co-agonist), and one positive control (high level of agonist). Simultaneous quantification of the two major effector functions, cytokine production and cytotoxicity using T cell activation assay showed expression of low level of agonist induced lower levels of IFN- γ production and degranulation and these outcomes were enhanced by the presence of co-agonist. They also used mutations in TCR to show that TCR binding to co-agonist pMHC is not required for co-agonist pMHC dependent CD8⁺T cell activation.

Major Concerns:

I think this is a very exciting, innovative protocol to study general T cell activation. It is clear and well written. The sample and controls were designed with a clear concept and the data presented here validate their claim without ambiguity. The statistical analysis of their data can be improved, as it is not clear how many times each experiment was repeated and what the sample size is for the functional analysis.

This information has been added to the revised version of the manuscript.

It is also not clear from the manuscript the critical amount of co-agonist needed to trigger T cell activation when low level of agonist present in the cell. Is it possible to vary the amount of tetracycline to regulate agonist antigen expression? Or is it possible to vary the amount of leaky agonist for a fixed amount of co-agonist? Can they measure quantitatively the relative contribution of agonist and co-agonist, which trigger T cell activation as in the case of high amount of agonist.

These are great questions, and we have addressed them in the discussion in the revised manuscript.

Minor Concerns:

1. A brief description of why they used the tetracycline-regulated expression?

This information has been added to the introduction in the revised manuscript.

2. The protocol steps 4.1, 4.1.1. and 4.1.2 seem to miss the continuity (in line 312, 315 and 319).

This section of the protocol has been substantially edited to improve clarity.

3. The centrifugation time is missing in line 317, 340 and 345.

Thank you for spotting this. This has been corrected.

4. The axis labels in the flow histograms are small and x-axis values in Fig 1D and 1E are missing.

Thank you for spotting this. This has been corrected.