

Show/Hide MAN Settings

Select a Sample to Subtract

Adjust Protein Dendrogram

📄 Download Current Network Data

Hint 1: Use mouse to select parts of the tree and display the MAN of corresponding samples.

Hint 2: Use mouse to click & drag parts (nodes) of the MAN if it appears congested.

Note 1: For publication-quality networks click "Download Current Network Data". This will download a .csv file of the currently-displayed network that can be easily imported into Gephi or Cytoscape.

Suggestions for Reporting MAN Analysis:

This MAN was created by analyzing 3 samples, retaining peaks with a signal to noise ratio above 4 , and occurring in greater than 70 % of replicate spectra. Peaks occurring below 200 m/z or above 2000 m/z were removed from the analysis.

