

IDBac

127.0.0.1:5125

IDBac

Introduction

Starting With Raw Data

Work With Previous Experiments

Create an IDBac experiment

- Click here to convert Bruker files
- Click here to convert mzML/mzXML files
- Click here to convert txt files
- Click here to convert Microtyper files

Starting with mzML or mzXML Data:

1: Enter a name for this new experiment

This will become a filename, non-valid characters will be removed.

Hint: Intead of a space, use "_".

Enter Experiment Name Here

Once entered, the filename-friendly version of the entered name will appear here once.

This will be the version of your experiment name that is saved.

2: Click to select the location of your mzML files

Raw Data Folder

No Folder Selected

Samples will be named according to the file name of the provided files

4: Click "Process Data" to begin spectra conversion.

Process Data