

IDBac

127.0.0.1:5125

☆

C

Transfer samples from previous experiments to new/other experiments.

Select an existing experiment to copy samples from:

example

Location of experiment file:

[1] "C:/Users/CMC/Documents/GitHub/IDBac_App/example.sqlite"

New experiment name:

Enter new experiment name here

Move samples between boxes by clicking the strain's name and then an arrow. Samples in the right box will be used for analysis.

Available Samples

172-1
172-10
172-11
172-7
Matrix

Selected Samples

➔

➞

Add to new Experiment