

# Journal of Visualized Experiments

## High Throughput Traction Force Microscopy using PDMS Reveals Dose-Dependent Effects of TGF- $\beta$ on the Epithelial-to-Mesenchymal Transition

--Manuscript Draft--

|  |   |
|--|---|
| <b>Article Type:</b>   | Invited Methods Article - Author Produced Video   |
| <b>Manuscript Number:</b>  | JoVE59364R1   |
| <b>Full Title:</b>   | High Throughput Traction Force Microscopy using PDMS Reveals Dose-Dependent Effects of TGF- $\beta$ on the Epithelial-to-Mesenchymal Transition   |
| <b>Keywords:</b>   | PDMS (Polydimethylsiloxane); Traction Force Microscopy; High throughput; Contractile Force Screening; mechanobiology; cancer; Silicone Elastomers; Epithelial-Mesenchymal Transition; EMT; Traction; biophysics |
| <b>Corresponding Author:</b>   | Allen Ehrlicher<br>McGill University<br>Montreal, QC CANADA   |
| <b>Corresponding Author's Institution:</b>                                   | McGill University   |
| <b>Corresponding Author E-Mail:</b>  | aje.mcgill@gmail.com  |
| <b>Order of Authors:</b>   | Haruka Yoshie<br>Newsha Koushki<br>Clayton Molter<br>Peter M. Siegel<br>Ramaswamy Krishnan<br>Allen J. Ehrlicher  |
| <b>Additional Information:</b>   |   |
| <b>Question</b>  | <b>Response</b>   |
| Please indicate whether this article will be Standard Access or Open Access. | Standard Access (US\$1200)  |



**Allen Ehrlicher**  
**Canada Research Chair in Active Biological Mechanics**  
**Associate Professor**  
**Department of Bioengineering**  
**Faculty of Engineering**  
**McGill University**

**Associate Member, Departments of:**  
**Department of Mechanical Engineering,**  
**Department of Biomedical Engineering,**  
**Goodman Cancer Research Centre,**  
**Department of Anatomy and Cell Biology**

McGill University  
358 McConnell Engineering Building  
3480 University Street  
Montreal, Quebec, Canada H3A 2A7

Université McGill  
270 Pavillon McConnell  
3840, rue Université  
Montréal, Québec, Canada H3A 2A7

Tel/Tél : (514) 714-8239  
Fax/Télécopieur: (514) 398-7379  
[allen.ehrlicher@mcgill.ca](mailto:allen.ehrlicher@mcgill.ca)  
<http://bam.lab.mcgill.ca>

February 27, 2019

Dear Dr DSouza:

Please find attached our revised manuscript; the editorial and reviewer comments were constructive and improved the manuscript, and we have addressed all points raised in the text, table, and Figures. These are addressed line by line in the rebuttal letter.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Allen J Ehrlicher", with a stylized, flowing script.

Allen J Ehrlicher

**TITLE:**

High Throughput Traction Force Microscopy using PDMS Reveals Dose-Dependent Effects of Transforming Growth Factor- $\beta$  on the Epithelial-to-Mesenchymal Transition

**AUTHORS AND AFFILIATIONS:**

Haruka Yoshie<sup>1</sup>, Newsha Koushki<sup>1</sup>, Clayton Molter<sup>1</sup>, Peter M. Siegel<sup>2,3</sup>, Ramaswamy Krishnan<sup>4</sup>, Allen J. Ehrlicher<sup>1,2</sup>

<sup>1</sup>Department of Bioengineering, McGill University, Montreal, QC, Canada

<sup>2</sup>Goodman Cancer Research Centre, McGill University, Montreal, QC, Canada

<sup>3</sup>Department of Medicine, McGill University, Montreal, QC, Canada

<sup>4</sup>Department of Emergency Medicine, Beth Israel Deaconess Medical Center, Boston, MA, USA

**Corresponding Author:**

Allen J. Ehrlicher (allen.ehrlicher@mcgill.ca)

**Email Addresses of Co-authors:**

Haruka Yoshie ([haruka.yoshie@mail.mcgill.ca](mailto:haruka.yoshie@mail.mcgill.ca))

Newsha Koushki ([newsha.koushki@mail.mcgill.ca](mailto:newsha.koushki@mail.mcgill.ca))

Clayton Molter ([molter@ualberta.ca](mailto:molter@ualberta.ca))

Peter M. Siegel ([peter.siegel@mcgill.ca](mailto:peter.siegel@mcgill.ca))

Ramaswamy Krishnan ([rkrishn2@bidmc.harvard.edu](mailto:rkrishn2@bidmc.harvard.edu))

**KEYWORDS:**

PDMS, polydimethylsiloxane, traction force microscopy, high throughput, contractile force screening, epithelial-to-mesenchymal transition, mechanobiology, cancer

**SUMMARY:**

We present a high throughput traction force assay fabricated with silicone rubber (PDMS). This novel assay is suitable for studying physical changes in cell contractility during various biological and biomedical processes and diseases. We demonstrate this method's utility by measuring a TGF- $\beta$  dependent increase in contractility during the epithelial-to-mesenchymal transition.

**ABSTRACT:**

Cellular contractility is essential in diverse aspects of biology, driving processes that range from motility and division, to tissue contraction and mechanical stability, and represents a core element of multi-cellular animal life. In adherent cells, acto-myosin contraction is seen in traction forces that cells exert on their substrate. Dysregulation of cellular contractility appears in a myriad of pathologies, making contractility a promising target in diverse diagnostic approaches using biophysics as a metric. Moreover, novel therapeutic strategies can be based on correcting the apparent malfunction of cell contractility. These applications; however, require direct quantification of these forces.

We have developed silicone elastomer-based traction force microscopy (TFM) in a parallelized multiwell format. Our use of a silicone rubber, specifically polydimethylsiloxane (PDMS), rather than the commonly employed hydrogel polyacrylamide (PAA) enables us to make robust and inert substrates with indefinite shelf-lives requiring no specialized storage conditions. Unlike pillar-PDMS based approaches that have a modulus in the GPa range, the PDMS used here is very compliant, ranging from approximately 0.4 kPa to 100 kPa. We create a high-throughput platform for TFM by partitioning these large monolithic substrates spatially into biochemically independent wells, creating a multiwell platform for traction force screening that is compatible with existing multiwell systems.

In this manuscript, we use this multiwell traction force system to examine the Epithelial to Mesenchymal Transition (EMT); we induce EMT in NMuMG cells by exposing them to TGF- $\beta$ , and to quantify the biophysical changes during EMT. We measure the contractility as a function of concentration and duration of TGF- $\beta$  exposure. Our findings here demonstrate the utility of parallelized TFM in the context of disease biophysics.

## INTRODUCTION:

Acto-myosin contractility is an essential element of active cell mechanics, impacting cell behaviors from motility and proliferation to stem cell differentiation. In tissues, contractility drives activity from polar separation in embryogenesis, to airway constriction and cardiac activity. Critically, to generate tension, cells must first adhere to their extracellular environment. In doing so, this contractility generates traction forces on their surroundings. Traction Force Microscopy (TFM) has emerged in a multitude of forms as a way to quantify these forces from diverse cells under different conditions.

The field of TFM has seen an exceptional breadth of innovation and application, and the results have paved the way for new perspectives in biology, which incorporate mechanics and physical forces. Starting with wrinkling silicone substrates<sup>1</sup>, researchers have applied various techniques to measure cell traction forces. These approaches have been continuously improved and have now reached a level of resolution on the order of several microns<sup>2</sup>. However, one principal problem has emerged, which is the difficulty in creating substrates of suitably low moduli using the available silicones. To circumvent this problem, polyacrylamide was adopted as a replacement due to the ease of creating substrates on the order of 1–20 kPa<sup>3</sup>. We recently implemented very compliant silicones in TFM<sup>4</sup>, allowing us to fabricate the same range of moduli as polyacrylamide, but with the advantages of inert and robust silicone.

TFM approaches have enabled valuable mechano-biological discoveries, however, a persistent shortcoming is their complexity, often restricting their use to researchers in the engineering or physical sciences disciplines. This is due, in large part, to the detailed calibrations and challenging calculations that are required to quantify contractility. Another significant challenge is that TFM methods are largely low-throughput and therefore ill-suited to study many different conditions or populations simultaneously<sup>5</sup>. This has presented a bottleneck, which has hampered transfer of TFM from a specialist biophysics setting into broader biological sciences and pharmacology applications.



We have recently developed a multi-well format TFM plate, which allows researchers to parallelize their TFM measurements for faster quantification of contractility metrics, while exploring the impact of different compounds and also using less reagents<sup>4</sup>. This methodology has broad utility in diverse mechanobiology studies, from evaluating the effects of compounds on cellular activity, to quantifying the contractile changes in differentiation or disease.

One area of biomedical research that will benefit greatly from TFM is the study of how physical cues impact the malignant phenotypes of cancer cells. Metastasis, responsible for 90% of cancer-related deaths, is characterized by cancerous cells leaving their original tumor site and colonizing a secondary site. For cells to migrate through tissue and pass in and out of the vascular system, they must radically change their shapes to squeeze through these physical barriers while generating substantial forces to pull their way along extracellular matrix or move between other cells. These forces are transmitted to the substrate through focal adhesion interactions<sup>2-3</sup>, and can be quantified using TFM. While cancers are biochemically exceptionally diverse, with an expanding repertoire of known mutations and protein changes, some common physical changes have been observed; in a variety of cancers, including breast, prostate, and lung cancers, metastatic cells have been shown to exert 2–3 times the traction forces of non-metastatic cells<sup>6-8</sup>. These results suggest that there may be a strong correlation between metastatic progression and the traction forces exerted by cells; however, the detailed time-dependent changes in contractility are difficult to examine.

The epithelial-to-mesenchymal transition (EMT) is a process whereby cells reduce adherens- and tight-junction mediated cell-cell adhesion, becoming more migratory and invasive. In addition to physiological functions that include wound healing and developmental processes, EMT is also a process exploited during metastasis, making it a useful model system to study this process. Using TGF- $\beta$ , we can induce the EMT in ErbB2-transformed murine mammary epithelial cells (NMuMG-derived)<sup>9</sup> to directly quantify the physical changes during this transformation, and characterize the time and dose-dependent effects of TGF- $\beta$  on EMT and cell contractility. In this article, we demonstrate the utility of this approach by measuring the changes in contractility during an induced EMT.

## PROTOCOL:

NOTE: The following protocol will guide researchers in fabricating and using the multi-well TFM dish shown in **Figure 1**.

### 1. Preparation of PDMS silicone substrates

1.1. Preparation of PDMS silicone rubber mixture based on a composite mixture of two commercially available kits.

1.1.1. Add Part A and Part B of PDMS kit (e.g., GEL-8100, see the **Table of Materials**) in a 1:1 weight ratio into the 50 mL tube.

NOTE: The mixture is mixed on a rotator at a speed slow enough for the mixture to flow back and forth during revolution to ensure complete mixing.

1.1.2. Add the required amount of curing agent for the desired modulus of the substrate.

NOTE: The amount of curing agent to be added to the mixture depends on the desired modulus of the substrate and may typically range from 0.1% to 1.8%. Refer to **Table 1** and **Figure 3** for a guide to specific crosslinker concentrations and resulting moduli.

1.1.3. Mix the formulation on the rotator for 30–45 min; ensure the rotation is slow enough for thorough mixing.

1.2. Bottom layer: coating PDMS substrates on the glass slide

1.2.1. Place the custom-built chuck illustrated in **Figure 2** on the spin-coater. Clean the glass with ethanol or isopropanol, and dry with a lint-free wipe. Place the glass slide in the chuck and turn on the vacuum to hold the slide in place.

1.2.2. Apply uncured PDMS approximately 1 cm from the edges and work in towards the center. Apply enough (3–4 mL) PDMS to ensure the whole surface will be covered.

NOTE: To ensure that the PDMS is evenly spread on the surface of the glass, a pipette tip may be helpful to spread the PDMS mixture from the center to the edges.

1.2.3. Spin the glass with the PDMS mixture on a spin-coater with the following protocol.

1.2.3.1. To spread the uncured PDMS on the slide, accelerate at 50 rpm/s from 0 to 200 rpm; hold at 200 rpm for 1 min.

1.2.3.2. To achieve a 100  $\mu\text{m}$  thick PDMS layer, accelerate at 50 rpm/s to 300 rpm and hold for 1 min at 300 rpm. Different desired thicknesses other than 100  $\mu\text{m}$  will require specific speed values.

1.2.3.3. To remove, decelerate at 50 rpm/s to 0 rpm. Disable vacuum and remove the coated slide, taking care not to touch the coated surface.

NOTE: It is important to include the acceleration and deceleration steps to ensure a smooth continuous surface.

CAUTION: To ensure that the sample does not fly off the chuck, the custom-made holder should be used to hold the slide, and not rely simply on vacuum and the existing flat chuck. The details and specifications of this holder are given in **Figure 2**.

1.2.4. Place the spin-coated sample in the oven at the manufacturer recommended temperature (100 °C) for 2 h.

NOTE: The surface of the oven where the sample is placed should be solid (i.e., not a wire rack) and level surface to ensure the uniform heating and thickness of the sample. A ceramic or steel plate makes an ideal surface.

### 1.3. Top bead layer

1.3.1. Add bead solution in the appropriate ratio to the remaining PDMS mixture.

NOTE: This ratio depends on the concentration of the stock bead solution and the desired bead density on the sample. Typical final values are  $9.2 \times 10^{11}$  beads/mL and 0.05–0.2 beads/ $\mu\text{m}^2$ , and an excess of beads is generally preferable to an inadequate amount.

1.3.2. Mix the bead solution with the uncured PDMS. This may be accomplished by placing the tube on a rotator for approximately 30 min, vortexing for 1–2 min, or sonication for 30 min. These methods may be combined.

NOTE: In our application, we find that sonication is effective in breaking bead aggregates, and rotation is effective in mixing. Synthesized fiduciary beads may aggregate in storage. Prior to use, they may be resuspended in hexane and sonicated. If there are significant large aggregates, one can filter the bead suspension through a 5  $\mu\text{m}$  syringe filter. This filtration step is optional; it helps to coat the sample with monodispersed beads, but a significant fraction of beads may be lost in the filter.

1.3.3. Take out the slide from the oven, allow to cool to room temperature, and place it on the spin-coater.

1.3.4. Add 3–4 mL of the bead and uncured PDMS mixture onto the surface of the coated sample.

NOTE: The mixture with beads added is less viscous due to the hexane. Make sure not to touch the surface of the substrate as it may damage the already-coated PDMS substrate. Additionally, the bead mixture may initially not wet the surface; take care that the mixture does not immediately flow off the cured PDMS surface.

1.3.5. Spin the sample with the following protocol.

1.3.5.1. To spread the bead and uncured PDMS mixture, accelerate at 100 rpm/s from 0 to 500 rpm; hold at 500 rpm for 1 min.

1.3.5.2. To achieve a thin layer of bead-embedded PDMS ( $\sim 1 \mu\text{m}$ ), accelerate at 200 rpm/s from 500 to 5000 rpm; hold at 5000 rpm for 10 s.

1.3.5.3. To remove, decelerate at 100 rpm/s to 0 rpm. Disable vacuum and remove coated slide, taking care not to touch the coated surface.

1.3.6. Place the spin-coated sample in the oven at 100 °C for 1 h.

NOTE: Temperatures above 100 °C or durations longer than 1 h can reduce the bead fluorescence. Make sure that the oven temperature is set to 100 °C and not higher.

1.3.7. The protocol can be paused here. To store the sample, cover the surface to avoid dust and light exposure. Make sure that nothing touches the surface. The sample is shelf-stable at room temperature indefinitely.

## 1.4. Assembling the plate

1.4.1. Add Part A (base) and Part B (curing agent, e.g., Sylgard) of a PDMS elastomer kit in a 1:1 weight ratio into the 50 mL tube.

1.4.2. Mix the mixture on the rotator for 30–45 min.

NOTE: For one plate, mix 5 mL of base with 0.5 mL of curing agent. Up to 1 mL of hexane can be added to reduce the viscosity of the mixture.

1.4.3. Apply the mixture to the bottom of the divider and spread the mixture.

NOTE: The divider should be placed upside-down.

1.4.4. Lay the substrate onto the divider upside-down.

1.4.5. Place the sample upside-down in the oven at 65 °C for 2 h.

1.4.6. Take out samples from the oven and clean the bottom of the glass with 70% ethanol or isopropanol to remove any PDMS residue.

NOTE: The protocol can be paused here. To store the sample, place a lid on the substrate and wrap the device in aluminum foil to avoid exposure to the light. The sample is shelf-stable at room temperature indefinitely. The divider utilized in this method is in a 96-well format; however, researchers may employ other formats (384-well, 2-well, 4-well, 8-well, etc.) depending on desired experiment setups and availability of dividing structures. Some further optimization may be required.

## 2. Surface functionalization

2.1. Dissolve 80 µL of a Sulfo-SANPAH aliquot in 40 mL of 0.1 M HEPES buffer (pH 7–9).

NOTE: Prepare Sulfo-SANPAH aliquots by dissolving 100 mg of Sulfo-SANPAH powder in 2 mL of sterile dimethyl sulfoxide (DMSO). Prepare 0.1 M HEPES buffer by diluting 50 mL of HEPES in 450 mL of sterile deionized water and filter through a 0.22  $\mu$ m pore filter.

2.2. Add 200  $\mu$ L of diluted Sulfo-SANPAH solution to each well of the 96-well plate.

2.3. Expose the plate to UV (300–460nm) light at appropriate distance and duration.

NOTE: After UV exposure, the color of the solution should be darker. UV exposure distance and duration depend on the UV lamp power. In our application, we expose for 10-15 min at a distance of 5 cm.

2.4. Remove the Sulfo-SANPAH solution from the wells and add 200  $\mu$ L of 5  $\mu$ g/mL fibronectin solution to each well.

NOTE: Researcher-specified protein can be used for surface coating. Some commonly used proteins are collagen, fibronectin, and laminin. We have found Sulfo-SANPAH to be the most effective method, and plasma cleaning while sometimes employing in PDMS is discouraged as it creates a silicon dioxide layer and visibly damages the surface.

2.5. Incubate the plate at 4  $^{\circ}$ C overnight.

NOTE: Different incubation methods can be applied depending on the protein used for coating.

2.6. Remove the fibronectin solution and wash each well with phosphate-buffered saline (PBS) twice.

2.7. Place a lid on the sample.

2.8. Add 200  $\mu$ L of PBS to each well.

NOTE: The protocol can be paused here. The fibronectin-coated samples can be stored at 4  $^{\circ}$ C for up to 2 weeks.

### 3. UV sterilization

3.1. Sterilize the sample under UV light in a biological safety cabinet for 30 min.

NOTE: Longer UV exposure times may negatively impact the bead fluorescence. All subsequent steps must be performed under sterile conditions.

### 4. Cell culture

4.1. Remove PBS from each well and add 200  $\mu$ L of cells suspended in culture media to each well.

4.1.1. Plate the cells at the desired cell density. Cell density depends on the desired experiment. For single cell studies, cells should be minimum of 50  $\mu\text{m}$  apart and cells near the edges of the imaging window should not be included in the TFM measurement. For monolayer cells, the imaging window should have the viewing field covered with a confluent layer of cells.

4.1.2. Prepare the complete growth culture media for NMuMG cells by supplementing DMEM with 5% FBS, 10 mM HEPES, 10  $\mu\text{g}/\text{mL}$  insulin, 1% penicillin-streptomycin, 1 mM L-glutamine, and 0.5  $\mu\text{g}/\text{mL}$  amphotericin B.

4.1.3. Prepare insulin stock by reconstituting in acidified water (2.5 mL of glacial acetic acid in 130 mL of deionized water) to a concentration of 10 mg/mL. Store the stock solution at 4  $^{\circ}\text{C}$ . Wait until the solution is clear, and then filter through a 0.22  $\mu\text{m}$  pore filter.

## 4.2. TGF- $\beta$ addition

4.2.1. To prepare of TGF- $\beta$  stock, dissolve 2  $\mu\text{g}$  of TGF- $\beta$  in 100  $\mu\text{L}$  of 10 mM citric acid (pH 3.0) and filter sterilize with 0.22  $\mu\text{m}$  pores. Vortex the tube and aliquot into the desired volumes. Store the aliquots at -80  $^{\circ}\text{C}$ .

4.2.2. Add 1.5  $\mu\text{L}$  of TGF- $\beta$  stock solution to 10 mL of the complete cell culture media to constitute the cell culture media with the final TGF- $\beta$  concentration of 3 ng/mL.

NOTE: To make 10 mM pH 3.0 citric acid, dilute the acid in water and adjust pH to 3.0 by adding HCl.

## 5. Data acquisition

5.1. For each position, acquire at least one image of fiduciary particles and cells. Focus on the bead layer.

NOTE: Pixel size should be optimized based on size of the fiduciary particles and image processing method being used. In this application, the authors use a 10x 0.4 NA objective, and images are acquired with 1024 x 1024 resolution, with 455 nm/pixel. In general, it is helpful to retain a resolution of at least approximately 1–5 pixels per bead; here, beads are polydisperse and have an individual size of 300–500 nm. It is critical that the fluorescent fiduciary beads be in focus for images to be used for TFM calculations. The focus and imaging quality of the beads should be prioritized over imaging the cells themselves. There should be no cross-talk between different channels, particularly any fluorescence not from the fiduciary beads which appears in the imaging spectra of the beads.

5.2. Once all the positions of interest have been recorded, add detachment solution to each well to acquire force-free reference images of the fiduciary particles.

5.2.1. To prepare cell detachment solution, mix an aqueous solution containing 2% TritonX-100, 50 mM sodium azide, and 500 mM potassium hydroxide.

NOTE: The above is provided as an example of an effective detachment solution. Different detachment solutions at researchers' discretion may be used to detach the cells off the substrate surface.

## 6. Image analysis

6.1. Perform appropriate image analysis as desired.

NOTE: Analysis software was developed in-house. Image analysis may be done with custom-made software or software available online.

## 7. Bead synthesis

NOTE: The following protocol is based on the synthesis method described by Klein et al.<sup>10</sup>.

7.1. Under a fume hood, prepare the three-neck flask with a water-cooled reflux condenser.

CAUTION: Set up a synthesis in a well-ventilated chemical fume hood.

7.2. Add 0.5 mL of PDMS stabilizer and fluorophore to the flask.

7.3. Equip one neck with a rubber septum with a nitrogen inlet needle and an outlet needle and equip the other neck with a rubber septum for adding reagent with a syringe.

7.4. Add 100 mL of anhydrous hexane in 250 mL to the flask and add a small magnetic stir bar.

7.5. Place the flask in the mineral oil bath at 75 °C and purge it with nitrogen gas for 1 h.

7.6. Add 6 mL of methylmethacrylate to a 25 mL round bottom flask.

7.7. Add 0.100 g of 2,2'-azobisisobutyronitrile (AIBN) to the round bottom flask and purge the mixture with nitrogen gas for 1 h.

NOTE: Flush methylmethacrylate through preppacked column to remove inhibitors before use. Add methylmethacrylate and AIBN mixture to the three-neck flask.

7.8. The solution initially becomes cloudy and turns milky. Let the reaction run for 3 h after the solution becomes cloudy.

7.9. After 3 h, place the flask in an ice water bath.

7.10. Vacuum-filter the solution with coarse filter paper.

7.11. Centrifuge the filtrate and re-suspend the particles in hexane.

NOTE: The volume of the hexane to be added depends on the desired concentration of the bead solution. Sonication facilitates re-dispersion of the bead particles in hexane solution. Beads produced by the authors have polydisperse diameters of approximately 300–500nm. Due to the use of cross-correlation pattern tracking to determine displacements, monodisperse beads are not required.

## 8. Rheology measurement protocol

NOTE: Rheology is not required for every researcher or experiment, but is necessary to quantify the moduli for new formulations of PDMS. In this protocol, we employ a shear rheometer to measure the effects of crosslinker, frequency, and strain on moduli of PDMS samples. Depending on the available tools and expertise, moduli may also be measured using many other mechanical analysis approaches. Additionally, researchers using this protocol may elect to use our published moduli presented in **Table 1**, **Figure 3** and **Figure 4**.

8.1. Use a rheometer with a 25 mm diameter parallel plate geometry. Other geometries may be used.

8.2. Initialize the system and calibrate the device and measuring system (parallel plate,  $d = 25$  mm). After measuring the zero gap, begin loading the PDMS sample.

8.3. As soon as the PDMS elastomer and crosslinking agent have been mixed, pipet the mixture onto the bottom plate of rheometer.

8.3.1. Move the spindle down to completely contact the top of the PDMS sample.

8.3.2. Carefully trim the loaded sample excess from the bottom plate.

8.4. In a strain sweep test, apply increasing strains for each composition with different crosslinking density to ensure the polymer structure remains in the linear viscoelastic regime during all shear measurements.

NOTE: Strain values relevant for cell studies are typically in the range of 0.1–10%. We have found PDMS to be linear up to approximately 100% strain.

8.5. Measure the dynamic shear storage modulus ( $G'$ ), and loss modulus ( $G''$ ) of the PDMS network in a time sweep test with frequency of 1 Hz and oscillatory shear strain of 0.5% at 100 °C.



8.6. To determine the viscoelasticity and time dependency of the final PDMS network, apply a frequency sweep test with frequency ranging 0.1–100 Hz and oscillatory shear strain of 0.5%.

#### REPRESENTATIVE RESULTS:

Before addition of TGF- $\beta$ , a confluent monolayer of cells has a cobblestone like shape and is tightly packed. Upon TGF- $\beta$  treatment, cells become more elongated in morphology, enlarging the cell area and acquiring a more mesenchymal phenotype. Utilizing the multi-well device fabricated with soft PDMS elastomers, the physical properties of cells in a total 17 different conditions were studied. The cells were treated with four different TGF- $\beta$  concentrations (0.5, 1, 2, and 4 ng/mL) and four different incubation times (12, 24, 48, 96 h), and these results are summarized in **Figure 5**. The cells treated with TGF- $\beta$  applied larger traction stresses and strain energies than the cells cultured without TGF- $\beta$ . Cells incubated with TGF- $\beta$  for 96 h showed the largest traction stresses and strain energies. The cells applied larger stresses and strain energies when treated with higher concentration of TGF- $\beta$ . The difference in tractions and strain energies were more distinct at longer incubation times.

The surface of the substrate needs to be smooth and uniformly coated with ligands, such as fibronectin or collagen. Scratched surface and/or non-uniform coating of ligands may lead to improper cell attachment, resulting in inaccurate traction measurement. **Figure 6** shows the localized traction stresses due to the non-uniform substrate surface.

#### FIGURE AND TABLE LEGENDS:

**Figure 1: Overview of multi-well plate fabrication.** (A) A custom-cut glass slide is the starting point. (B) The glass slide is coated with a thick ( $\sim 100$   $\mu$ m) layer of PDMS. (C) A layer of fiduciary beads (shown in green) are then spin-coated in a  $\sim 1$   $\mu$ m thick layer on top of the previous layer. (D) The multi-well divider is carefully placed on top of the fiduciary bead layer. (E) The complete multi-well plate is assembled and ready for use or storage.

**Figure 2: Custom chuck for spin-coater.** To prevent the bottom glass (where PDMS is coated) from flying off of the standard spin-coater chuck, a custom-made holder is placed onto the chuck. (A) A custom-made holder for the spin-coater chuck. (B) Engineering drawing of the holder with all the dimensions.

**Figure 3: Substrate stiffness with changing crosslinker concentration (wt%).** PDMS mixtures are prepared as described, with the additional crosslinker percentage increasing the elastic modulus until a max of 100 kPa at 1.8 weight percent (wt%) crosslinker. As a guide, the modulus as a function of crosslinker is fit linearly, which can be used by researchers to formulate their own mixture at a particular desired modulus within this range.

**Figure 4: Rheology of PDMS elastomer containing 0, 0.15, 0.36 and 1 wt% of crosslink agents at a temperature of 100 °C.** Triangle symbols indicate loss modulus ( $G''$ ) and square symbols indicate storage modulus ( $G'$ ). (A) Oscillatory time sweep at frequency of 1 Hz and shear strain of

0.5% during gelation. (B) Oscillatory frequency sweep of PDMS network at shear strain of 0.5%. (C) Strain sweep of PDMS network at frequency of 1 Hz. All data points were acquired in triplicate.

**Figure 5: Representative results utilizing traction stress device with multiwell format.** Monolayer of cells in different conditions were cultured in a multi-well device to measure contractility and strain energy. (A) Graph of traction stresses with increasing TGF- $\beta$  incubation time for different TGF- $\beta$  concentrations. Traction stresses increased with increasing TGF- $\beta$  concentrations and incubation time. All data are statistically significant with respect to control (i.e., no TGF- $\beta$ ) except the following: 12 h - 0.5 ng, 12 h - 1 ng, 48 h - 1 ng. (B) Graph of strain energy with increasing TGF- $\beta$  incubation time for different TGF- $\beta$  concentrations. Strain energies increased with increasing TGF- $\beta$  concentrations and incubation time. Sample sizes range from  $n = 7$  to  $n = 15$ . All data are statistically significant with respect to control (i.e., no TGF- $\beta$ ) except the following: 12 h - 0.5 ng, 12 h - 1 ng, 24 h - 0.5 ng, 24 h - 1 ng, 48 h - 1 ng. Statistical significance was determined using the Kruskal Wallis test, which does not assume a normal distribution.

**Figure 6: Representative results from a sub-optimal (A, B) and a satisfactory (C, D) experiment.** (A) Reflection image of the substrate surface. Undesired contaminants, which can include dust or fibers, or material defects such as scratches are present on the substrate surface. (B) Stress map of cell contractility: due to the non-uniform surface, irregular and inaccurate traction stresses are observed around the objects. (C) Reflection image of the substrate surface. No apparent contaminants are visible. (D) Stress map of cell contractility: no artifact discontinuities are visible.

**Table 1: Substrate stiffness with changing crosslinker concentration (wt%).** By changing additional crosslinker concentrations, PDMS substrates of the desired Young's moduli are fabricated. PDMS shear moduli were measured with a rheometer and Young's moduli were determined. For each data point, 3 independent preparations were done and the standard deviation is given.

## DISCUSSION:

For the success of this method, it is critical to have a uniformly coated sample with a constant thickness of approximately 100  $\mu\text{m}$ . The modulus should be carefully chosen to examine the physical significance of the biological system of interest. When fabricating a top layer, the concentration of the fiducial fluorescent particles should be optimized for accurate analysis of displacement and traction stress. Analyzing isolated single cells requires a denser fiduciary layer than measuring confluent monolayers. Additionally, the surface of the substrate should have stable and uniform coating of adhesion molecules such as collagen, fibronectin, and laminin to ensure proper attachment of the cells to the substrate surface. Particular care must be taken when attaching the multiwell divider; it should be placed without sliding to prevent sealing PDMS glue from being smeared on the culture surface, and the outer edge must be carefully aligned with the glass dimensions to prevent any leaks on the border wells.

To prevent wells from being leaky, an adequate amount of PDMS glue needs to be applied to the well-divider to hold it in place. However, use of excessive amounts will cover cell culture surface. Sufficient volume of PDMS needs to be added during the spin-coating procedure. When curing, the oven and racks need to be level. Inadequate PDMS or an unlevelled oven may lead to uneven PDMS thickness.

Bead density needs to be optimized for proper analysis. When the bead density is not adequate, concentration of the beads in fiduciary PDMS mixture needs to be increased. Additionally, an adequate amount of the mixture needs to be added for spin coating. To ensure proper fluorescent signals of the fiduciary bead particles, curing condition should be carefully monitored. Longer times and higher temperatures than the ones specified in this protocol may degrade the fluorescence. Improper protein coating on the surface of the sample may lead to poor cell adhesion. To prevent this, the surface needs to be cleaned and the expiration of reagents such as Sulfo-SANPAH and ligands needs to be checked. UV exposure time and strength should be optimized. Immunofluorescence or other fluorescent protein construct may be tested to ensure uniform ligand coating.

The device fabricated with this method and set of PDMS formulations is only applicable to substrates with Young's moduli ranging from 0.4 kPa to 100 kPa. Other moduli may be possible using alternative formulations, however, the current range spans a large set of physiologically relevant values. While this method is specifically for multi-well fabrication, it may also be used to prepare individual slides or other slide geometries, and in these instances further user troubleshooting may be necessary. In this embodiment, a relatively thick (1 mm) glass slide is employed, precluding common high numerical aperture (NA) objectives on an inverted microscope. While it is technically feasible to construct these multi-well TFM dishes using thinner glass, the fragility of these in processing and assembly may result in a high rate of breakage, and can run contrary to the higher-throughput approach. Furthermore, surface coating may be investigated further for wider applications of the device.

Utilizing a multiwell format, high-throughput experiments are possible. Multiwell format enabled us to examine the physical changes of NMuMG cells during EMT with TGF- $\beta$  treatment with combinations of different concentrations and different incubation times in a single dish. Fabrication of traction stress devices with hydrogels such as polyacrylamide gel have several limitations. Dominant ingredient being water, hydrogels are sensitive to salt concentrations (osmolarity), temperature, humidity, and pH changes. Changes in any of these properties can lead to the change in mechanical properties of the material, requiring extra care in using samples and interpreting results. Furthermore, having water, hydrogels are more susceptible to infections, requiring extra cautions. In contrast to water-based hydrogels, silicone-based PDMS is stable and inert. Once fabricated, it can be stored at room temperature indefinitely without requiring any special handling and storing condition.

The impermeable nature of PDMS allows us to fabricate a monolithic substrate, ensuring that all wells are truly identical in substrate thickness and composition, while allowing unique biochemistry to be applied in the media, or different cells to be utilized in each well. A hydrogel-

based surface allows diffusive factors to permeate and travel through it, necessitating adding the hydrogel to wells that are otherwise partitioned to prevent cross-contamination.

This method allows researchers to study cell contractility, a basic and ubiquitous aspect of cell biology, in a high-throughput manner, enabling more efficient and reproducible data acquisition. This helps to translate Traction Force Microscopy into Contractile Force Screening, allowing researchers to truly use contractility as a metric for cell activity, or efficacy of a compound. As a methodology, this has broad applications across biophysical and biomedical sciences, from understanding the physics of cells, to characterizing and testing pharmaceutical compounds against standardized cell types, or perhaps highly specialized individual cells for personalized medicine.

#### ACKNOWLEDGMENTS:

The authors thank Tom Kodger, Michael Landry, and Christopher J. Barrett for assistance with bead synthesis. A.J.E. acknowledges Natural Sciences and Engineering Research Council grants RGPIN/05843-2014 and EQPEQ/472339-2015, Canadian Institutes of Health Research grant no. 143327, Canadian Cancer Society grant no. 703930, and Canadian Foundation for Innovation Project #32749. R. Krishnan acknowledges National Institutes of Health grant no. R21HL123522 and R01HL136209. H.Y. was supported by Fonds de recherche Santé Québec, and Fonds de recherche Nature et Technologies Québec. The authors thank Johanan Idicula and Zixin He for assistance in preparing the video component of this manuscript.

#### DISCLOSURES:

AJE and RK have interest in Live Cell Technologies, a company which fabricates materials described in this article.

#### REFERENCES:

1. Harris, A., Wild, P., Stopak, D. Silicone rubber substrata: a new wrinkle in the study of cell locomotion. *Science* **208** (4440), 177–179 (1980).
2. Oliver, T., Dembo, M., Jacobson, K. Traction forces in locomoting cells. *Cell Motil Cytoskeleton* **31** (3), 225–240 (1995).
3. Dembo, M., Wang, Y.-L. Stresses at the Cell-to-Substrate Interface during Locomotion of Fibroblasts. *Biophysical Journal* **76** (4), 2307–2316 (1999).
4. Yoshie, H., Koushki, N., et al. Traction Force Screening Enabled by Compliant PDMS Elastomers. *Biophysical Journal* **114** (9), 2194–2199 (2018).
5. Park, C. Y., Zhou, E. H., et al. High-throughput screening for modulators of cellular contractile force. *Integrative biology : quantitative biosciences from nano to macro* **7** (10), 1318–24 (2015).
6. Kraning-Rush, C. M., Califano, J. P., Reinhart-King, C. A. Cellular traction stresses increase with increasing metastatic potential. *PLoS ONE* **7** (2), e32572 (2012).
7. Physical Sciences - Oncology Centers Network, Agus, D. B., et al. A physical sciences network characterization of non-tumorigenic and metastatic cells. *Scientific Reports* **3** (1), 1449 (2013).
8. Guo, M., Ehrlicher, A. J., et al. Probing the Stochastic, Motor-Driven Properties of the

612 Cytoplasm Using Force Spectrum Microscopy. *Cell* **158** (4), 822–832 (2014).  
613 9. Ngan, E., Northey, J. J., Brown, C. M., Ursini-Siegel, J., Siegel, P. M. A complex containing  
614 LPP and alpha-actinin mediates TGFbeta-induced migration and invasion of ErbB2-expressing  
615 breast cancer cells. *Journal of Cell Science* **126** (Pt 9), 1981–1991 (2013).  
616 10. Klein, S. M., Manoharan, V. N., Pine, D. J., Lange, F. F. Preparation of monodisperse PMMA  
617 microspheres in nonpolar solvents by dispersion polymerization with a macromonomeric  
618 stabilizer. *Colloid & Polymer Science* **282** (1), 7–13 (2003).  
619








**A****B****C****D****E**



Figure 2

A

[Click here to access/download;Figure;Fig2\\_Chuck\\_complete.pdf](#)



2




1

B

B

B

A



All dimensions are in millimeters (mm)

Author: Clayton Molter

TITLE:

Custom Chuck to Locate and Support Slides during Spin Coating

SIZE

DWG. NO.

REV


A

1




A



SCALE: 1:2

SHEET 1 OF 1



















| Crosslinker concentration (wt%) | G' (Pa) | Standard Deviation | E (Pa) | Standard Deviation |
|---------------------------------|---------|--------------------|--------|--------------------|
| 0                               | 0.135   | 0.014              | 0.405  | 0.043              |
| 0.15                            | 1       | 0.1                | 3      | 0.35               |
| 0.36                            | 4.027   | 0.245              | 12.081 | 0.73               |
| 0.75                            | 16.01   | 0.49               | 48.03  | 1.2                |
| 1                               | 18.44   | 0.989              | 55.32  | 1.94               |
| 1.5                             | 27.638  | 0.93               | 82.91  | 1.64               |
| 1.8                             | 33.986  | 0.88               | 101.94 | 1.088              |
| 2                               | 33.36   | 0.67               | 100.08 | 1.1                |

| Name of Material/ Equipment                                   | Company                    | Plate                   | Catalog Number |
|---|----------------------------|-------------------------|----------------|
| GEL-8100  | Nusil Technology           | GEL-8100                |                |
| Dow Corning Sylgard 184 Silicone Encapsulant Clear 0.5 kg Kit | Ellsworth Adhesives        | 184 SIL ELAST KIT 0.5KG |                |
| Custom Cut Glass  | Hausser Scientific Company |                         |                |
| Target 2TM Nylon Syringe Filter                               | ThermoFisher Scientific    | F2513-4                 |                |
| 96-well Stripwell Egg Crate Strip Holder                      | Corning                    | 2572                    |                |
| Polystyrene Universal Microplate Lid With Corner Notch        | Corning                    | 3099                    |                |
| Ethyl alcohol   | Greenfield Global          | P016EA95                |                |
| 2-Propanol  | Sigma-Aldrich              | 190764                  |                |
|   |                            | Surface Coating         |                |
| Sulfo-SANPAH Crosslinker                                      | Proteochem                 | c1111-100mg             |                |
| Fibronectin bovine plasma                                     | Sigma-Aldrich              | F1141-1MG               |                |
| PBS, 1X   | Wisent                     | 319-005-CL              |                |
| DMSO  | Sigma-Aldrich              |                         | 472301         |
|   |                            | Cell Culture            |                |
| DMEM, 1X  | Wisent                     | 319-005-CL              |                |
| FBS (Fetal Bovine Serum)                                      | Wisent                     | 080-150                 |                |
| HEPES   | Wisent                     | 330-050-EL              |                |
| Human Insulin Recombinant                                     | Wisent                     | 511-016-CM              |                |
| Penicillin-Streptomycin Solution                              | Wisent                     | 450-201-EL              |                |
| L-Glutamine solution  | Wisent                     | 609-065-EL              |                |

|   |                      |                       |        |
|---|----------------------|-----------------------|--------|
| Amphotericine B   | Wisent               | 450-105-QL            |        |
| Recombinant Human TGF- $\beta$ 1  | Peprotech            | 100-21                |        |
| Acetic acid   | Sigma-Aldrich        |                       | 537020 |
| Citric acid   | Sigma-Aldrich        | 251275                |        |
| NMuMG   | ATCC                 | CRL-1636              |        |
| Sodium azide  | Fisher Scientific    | AC190385000           |        |
| Potassium hydroxide   | Sigma-Aldrich        |                       | 221473 |
| TritonX-100   | Sigma-Aldrich        | X100                  |        |
| <b>Bead Synthesis</b>   |                      |                       |        |
| 1,1'-Dioctadecyl-3,3,3',3'-<br>tetramethylindocarbocyanine perchlorate<br>(DiI) | Sigma-Aldrich        | 468495-100MG          |        |
| Methyl methacrylate   | Sigma-Aldrich        | M55909-500ML          |        |
| Inhibitor Remover   | Sigma-Aldrich        | 306312-1EA            |        |
| Methacryloxypropyl Terminated<br>Polydimethylsiloxane                           | Gelest               | DMS-R31 (25,000g/mol) |        |
| 2,2'-Azobis(2-methylpropionitrile) (AIBN)                                       | Sigma-Aldrich        | 441090-25G            |        |
| Hexane  | Sigma-Aldrich        | 296090-2L             |        |
| Hexane, mixture of isomers  | Sigma-Aldrich        | 227064-1L             |        |
| Whatman qualitative filter paper, Grade 1                                       | Sigma-Aldrich        | WHA1001055            |        |
| <b>Equipment</b>  |                      |                       |        |
| Laurell WS-650Mz-23NPPB   | Laurell Technologies |                       |        |
| UVP Handheld UV Lamp Model UVGL-58  | VWR                  | 21474-622             |        |
| Rheometer   | Anton Paar           | MCR 302 WESP          |        |

## Comments/Description

High Purity Dielectric, Soft Silicone Gel kit

curing agent

109.6mm± x 72.8mm± x 1mm thickness

0.95

ACS reagent, ≥99.5%

solution, sterile-filtered, BioReagent,  
suitable for cell culture

pH 7.4, without calcium and magnesium

4.5g/L glucose, with L-glutamine, sodium  
pyruvate and phenol red

Premium Quality, Endotoxin <1,

Hemoglobin <25

1M, free acid

USP grade

100 X, sterile filtered for cell culture

200mM solution, sterile filtered for cell  
culture

250µg/ml, sterile filtered for cell culture

HEK293 Derived

Glacial, ≥99.85%

ACS reagent, ≥99.5%

Mouse Mammary Gland Cell Line

99%, extra pure, ACROS Organics

ACS reagent, ≥85%, pellets

laboratory grade

97%

contains ≤30 ppm MEHQ as inhibitor, 99%

Prepacked column for removing

hydroquinone and monomethyl ether

hydroquinone

Polydimethylsiloxane stabilizer,

25,000g/mol, 1,000 cSt

98%

anhydrous, 95%

anhydrous, ≥99%

circles, diam. 55 mm,





1 Alewife Center #200  
Cambridge, MA 02140  
tel. 617.945.9051  
www.jove.com

## ARTICLE AND VIDEO LICENSE AGREEMENT

Title of Article:

High Throughput Traction Force Microscopy using PDMS Reveals Dose-Dependent Effects of Transforming Growth Factor- $\beta$  on the Epithelial-to-Mesenchymal Transition

Author(s):

Haruka Yoshie, Newsha Koushki, Clayton Molter, Peter M. Siegel, Ramaswamy Krishnan, Allen J. Ehrlicher

Item 1 (check one box): The Author elects to have the Materials be made available (as described at

<http://www.jove.com/author>) via: ☒ Standard Access ☐ Open Access

Item 2 (check one box):

- ☒ The Author is NOT a United States government employee.
- ☐ The Author is a United States government employee and the Materials were prepared in the course of his or her duties as a United States government employee.
- ☐ The Author is a United States government employee but the Materials were NOT prepared in the course of his or her duties as a United States government employee.

### ARTICLE AND VIDEO LICENSE AGREEMENT

1. **Defined Terms.** As used in this Article and Video License Agreement, the following terms shall have the following meanings: “**Agreement**” means this Article and Video License Agreement; “**Article**” means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein; “**Author**” means the author who is a signatory to this Agreement; “**Collective Work**” means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; “**CRC License**” means the Creative Commons Attribution-Non Commercial-No Derivs 3.0 Unported Agreement, the terms and conditions of which can be found at: <http://creativecommons.org/licenses/by-nc-nd/3.0/legalcode>; “**Derivative Work**” means a work based upon the Materials or upon the Materials and other pre-existing works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; “**Institution**” means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; “**JoVE**” means MyJoVE Corporation, a Massachusetts corporation and the publisher of *The Journal of Visualized Experiments*; “**Materials**” means the Article and / or the Video; “**Parties**” means the Author and JoVE; “**Video**” means any video(s) made by the Author, alone or in conjunction with any other parties, or by JoVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion of the Article, and in which the Author may or may not appear.

2. **Background.** The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JoVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.

3. **Grant of Rights in Article.** In consideration of JoVE agreeing to publish the Article, the Author hereby grants to JoVE, subject to **Sections 4 and 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the “Open Access” box has been checked in **Item 1** above, JoVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.

## ARTICLE AND VIDEO LICENSE AGREEMENT

4. **Retention of Rights in Article.** Notwithstanding the exclusive license granted to JoVE in **Section 3** above, the Author shall, with respect to the Article, retain the non-exclusive right to use all or part of the Article for the non-commercial purpose of giving lectures, presentations or teaching classes, and to post a copy of the Article on the Institution's website or the Author's personal website, in each case provided that a link to the Article on the JoVE website is provided and notice of JoVE's copyright in the Article is included. All non-copyright intellectual property rights in and to the Article, such as patent rights, shall remain with the Author.

5. **Grant of Rights in Video – Standard Access.** This **Section 5** applies if the "Standard Access" box has been checked in **Item 1** above or if no box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby acknowledges and agrees that, Subject to **Section 7** below, JoVE is and shall be the sole and exclusive owner of all rights of any nature, including, without limitation, all copyrights, in and to the Video. To the extent that, by law, the Author is deemed, now or at any time in the future, to have any rights of any nature in or to the Video, the Author hereby disclaims all such rights and transfers all such rights to JoVE.

6. **Grant of Rights in Video – Open Access.** This **Section 6** applies only if the "Open Access" box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby grants to JoVE, subject to **Section 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Video in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Video into other languages, create adaptations, summaries or extracts of the Video or other Derivative Works or Collective Works based on all or any portion of the Video and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. For any Video to which this Section 6 is applicable, JoVE and the Author hereby grant to the public all such rights in the Video as provided in, but subject to all limitations and requirements set forth in, the CRC License.

7. **Government Employees.** If the Author is a United States government employee and the Article was prepared in the course of his or her duties as a United States government employee, as indicated in **Item 2** above, and any of the licenses or grants granted by the Author hereunder exceed the scope of the 17 U.S.C. 403, then the rights granted hereunder shall be limited to the maximum rights permitted under such

statute. In such case, all provisions contained herein that are not in conflict with such statute shall remain in full force and effect, and all provisions contained herein that do so conflict shall be deemed to be amended so as to provide to JoVE the maximum rights permissible within such statute.

8. **Likeness, Privacy, Personality.** The Author hereby grants JoVE the right to use the Author's name, voice, likeness, picture, photograph, image, biography and performance in any way, commercial or otherwise, in connection with the Materials and the sale, promotion and distribution thereof. The Author hereby waives any and all rights he or she may have, relating to his or her appearance in the Video or otherwise relating to the Materials, under all applicable privacy, likeness, personality or similar laws.

9. **Author Warranties.** The Author represents and warrants that the Article is original, that it has not been published, that the copyright interest is owned by the Author (or, if more than one author is listed at the beginning of this Agreement, by such authors collectively) and has not been assigned, licensed, or otherwise transferred to any other party. The Author represents and warrants that the author(s) listed at the top of this Agreement are the only authors of the Materials. If more than one author is listed at the top of this Agreement and if any such author has not entered into a separate Article and Video License Agreement with JoVE relating to the Materials, the Author represents and warrants that the Author has been authorized by each of the other such authors to execute this Agreement on his or her behalf and to bind him or her with respect to the terms of this Agreement as if each of them had been a party hereto as an Author. The Author warrants that the use, reproduction, distribution, public or private performance or display, and/or modification of all or any portion of the Materials does not and will not violate, infringe and/or misappropriate the patent, trademark, intellectual property or other rights of any third party. The Author represents and warrants that it has and will continue to comply with all government, institutional and other regulations, including, without limitation all institutional, laboratory, hospital, ethical, human and animal treatment, privacy, and all other rules, regulations, laws, procedures or guidelines, applicable to the Materials, and that all research involving human and animal subjects has been approved by the Author's relevant institutional review board.

10. **JoVE Discretion.** If the Author requests the assistance of JoVE in producing the Video in the Author's facility, the Author shall ensure that the presence of JoVE employees, agents or independent contractors is in accordance with the relevant regulations of the Author's institution. If more than one author is listed at the beginning of this Agreement, JoVE may, in its sole discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JoVE reserves the right, in its absolute and sole discretion and without giving any reason therefore, to accept or decline any work submitted to JoVE. JoVE and its employees, agents and independent contractors shall have

## ARTICLE AND VIDEO LICENSE AGREEMENT

full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JoVE has sole discretion as to the method of making and publishing the Materials, including, without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

11. **Indemnification.** The Author agrees to indemnify JoVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JoVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JoVE, making of videos by JoVE, or publication in JoVE or elsewhere by JoVE. The Author shall be responsible for, and shall hold JoVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to the making of a video by JoVE its employees, agents or independent contractors. All sterilization, cleanliness or decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's

expense. All indemnifications provided herein shall include JoVE's attorney's fees and costs related to said losses or damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JoVE, its employees, agents or independent contractors.

12. **Fees.** To cover the cost incurred for publication, JoVE must receive payment before production and publication the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JoVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.

13. **Transfer, Governing Law.** This Agreement may be assigned by JoVE and shall inure to the benefits of any of JoVE's successors and assignees. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to me one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement required per submission.

### CORRESPONDING AUTHOR:

Name:

Allen Ehrlicher

Department:

Bioengineering


Institution:

McGill University

Article Title:

High Throughput Traction Force Microscopy using PDMS Reveals Dose-Dependent Effects of Transforming Growth Factor- $\beta$  on the Epithelial-to-Mesenchymal Transition

Signature:



Date:

Nov 5 2018

Please submit a signed and dated copy of this license by one of the following three methods:

- 1) Upload a scanned copy of the document as a pdf on the JoVE submission site;
- 2) Fax the document to +1.866.381.2236;
- 3) Mail the document to JoVE / Attn: JoVE Editorial / 1 Alewife Center #200 / Cambridge, MA 02139

For questions, please email [submissions@jove.com](mailto:submissions@jove.com) or call +1.617.945.9051



## Rebuttal letter.

**Page 1 addresses editorial comments, and pages 2 & 3 address reviewer comments.**

### Editorial Comments:

- Please take this opportunity to thoroughly proofread the manuscript to ensure that there are no spelling or grammatical errors.

**We have checked the manuscript again and not found any spelling or grammatical errors.**

- **Protocol Detail:** Please add more specific details (e.g. button clicks for software actions, numerical values for settings, etc) to your protocol steps. Some examples of missing items are listed below. Please double check all your protocol steps and ensure that the manuscript matches the video.

**We believe that button clicks and software actions are not constructive for the reader, as these vary substantially between versions and systems, whereas the physical aspects do not. We have added details to acquisition specifics that should allow methodology, platform, and system independence.**

- 1) 5.1: how is the image acquired? Mention magnification, type of microscopy and any relevant settings.

**The method of image acquisition is not key, simply that it is fluorescence microscopy, of which confocal microscopy is employed and subset thereof. We have added this in the protocol, and provided details about minimum resolution required.**

- 2) Please ensure that all specific details from the video are mentioned in the manuscript (the manuscript should be able to serve as a standalone article). Examples of missing items: 05:44-06:09, 06:18-07:20, 07:32-07:45, etc.

**We have corrected these in the manuscript.**

- 3) Please be consistent with naming between the manuscript text and video (For eg, in section 1.2 you use "silicone" in the manuscript and PDMS in the video.

**Silicone is more general term for PDMS, but we have now specified PDMS throughout.**

- 4) 1.3.2: details provided here don't match the video.

**We have corrected these in the manuscript.**

- If your figures and tables are original and not published previously or you have already obtained figure permissions, please ignore this comment. If you are re-using figures from a previous publication, you must obtain explicit permission to re-use the figure from the previous publisher (this can be in the form of a letter from an editor or a link to the editorial policies that allows you to re-publish the figure). Please upload the text of the re-print permission (may be copied and pasted from an email/website) as a Word document to the Editorial Manager site in the "Supplemental files (as requested by JoVE)" section. Please also cite the figure appropriately in the figure legend, i.e. "This figure has been modified from [citation]."

**Figure 1 is similar to that published in our previous manuscript (Yoshie et al Biophysical Journal 2018), however, it was newly created by the authors for this JoVE manuscript, and thus is not a reproduction, or modification of a reproduction.**

- **Video Comment: at 13:31** - Please remove the '[www.jove.com](http://www.jove.com)' from the bottom of the title card, since the video wasn't produced by JoVE.

**We have removed this.**

### Comments from Peer-Reviewers:

**Reviewer #1:**

*The authors present a high throughput traction force assay fabricated with silicone rubber. In the last two decades, the key role of cell mechanics in health and disease (e.g. cancer) has spurred the development of several tools to measure forces at single cell level and at optical resolution. However, such tools generally lack any true high throughput capability and have been mainly used by specialist biophysics lab rather than biology/pharmacology groups. The present contractile force screening assay highly contributes to overcoming these limitations and its power is exemplified by quantifying changes in the mechanical response of epithelial cells after induction of EMT through TGF- $\beta$ .*

We thank the reviewer for their positive appraisal.

The authors are encouraged to address the following points:

*\*Line 57: "Unlike pillar-PDMS based approaches that have a modulus in the GPa range". It is not clear to which approaches the authors refer to. The common PDMS micropillar/micropost arrays are prepared using 10:1 base:curing agent ratio (e.g. Sylgard 184) that confers a modulus in the MPa range for the bulk array. By properly tuning the geometry (pillar length and diameter), the effective stiffness that cells sense on top of the pillars is in the 10-150 kPa range (e.g. as described in Ghibaudo et al. Soft Matter 2008 - Ladoux group). Therefore the authors are encouraged to clarify their sentence.*

We regret the typo of GPa, which should indeed read as MPa; this has been corrected, and we thank the reviewer for pointing this out. We are indeed referring to approaches such as that referenced by the reviewer, where the pillar geometry results in a given stiffness, from a MPa modulus. As any technique has its merits and drawbacks, and the discontinuous nature of pillar-based TFM has both, we are simply differentiating our approach from those to prevent confusion regarding the use of PDMS in TFM, as virtually all PDMS applications in TFM employ Sylgard 184.

*\*While it is clear from the video, the protocol reads confusing on which PDMS kit is being used. Readers that are mainly familiar with the Sylgard 184 kit (base and curing agent) might get confused when reading Part A and Part B (of GEL-8100?) followed by the addition of curing agent (from Sylgard 184?).*

We agree with the reviewer and had omitted the specifications based on our understanding of JoVE's request to leave out product names, however, these are specified in the list of materials. Nevertheless, to prevent confusion, we have specified these names in the protocol, and that the general approach is a combination of two commercially available kits.

*\*The glass coverslip is 1 mm thick. This prevents the use of high NA oil immersion objectives for high resolution imaging. The achievable resolution is clearly enough to calculate the force fields, whereas insufficient to resolve subcellular structures such as focal adhesions. In the "Limitations of the method" section, the authors should address whether it is possible to use thinner coverslips (#0, 1, 1.5).*

We have added a comment in the text. While it is not the goal of this protocol to resolve subcellular structures, and so it is somewhat beyond the scope to address this, we have made multi-well plates with coverslips, as it is technically possible, but we do not pursue it as it is problematic in this approach for several reasons, and thus not recommended for a high-throughput general contractility assay for which this is designed. First, the general fragility of coverslips, particularly in a large format, results in a large percent of breakage in fabrication. Second, immersion objectives require immersion fluid, and using this over a large area again is problematic. TFM with PDMS on coverslips is performed by the authors and others, but not employing this current multiwell high-throughput methodology for the above reasons.

*\*Statistical tests of data in Figure 5 are necessary to support the significant change in the mechanical response as claimed in the manuscript.*

We have quantified statistical significance using the Kruskal Wallis test, which does not assume a normal distribution, and have added this information in Figure 5.

*\*Figure 6 should contain also a representative image of an optimal experiment.*  
We have added this as Figure 6C&D.

*\*To achieve surface coating, would it be possible to use UV/Ozone treatment instead of Sulfo-SANPAH? Please address this point in the manuscript.*

While we believe that there are likely many methods to functionalize the surface, in this protocol we report sulfo-Sanpah as it has been the most reliable method. Plasma cleaning is discouraged as it creates a silicon dioxide layer and visibly damages the surface. We have added a note to this effect in the text.

**Reviewer #2:**

*Manuscript Summary:*

*The protocol is very clear and easy to follow with the video and the manuscript. Overall, I think it is very useful protocol to perform high-throughput TFM. I like the strategy presented for making the plates; however, it's unclear how traction stress measurements will be affected by the multi-well plate format. It would be helpful to include information about data acquisition, including limitations in resolution and drift correction. It's unclear how accurate the measurements are, directly comparing single well traction to multi-well plate format might clear this up.*

We thank the reviewer for their comments.

*Major Concerns:*

*The paper doesn't include some essential information. Questions below address these points.*

- 1. What is the size range for the synthesized beads that are being used?*
- 2. What is the objective magnification? If air vs oil/water objective will provide any distinct advantage?*
- 3. What is the UV exposure time used for the reported experiment to activate sulfo-sanpah?*
- 4. Table 1 is missing from the document.*
- 5. Will the thickness of the glass plate interfere with image resolution at higher magnifications. Discussion of image resolution seems important.*
- 6. Is the final PDMS coated plate able to be autoclaved? Sterility would be a concern unless the whole process is performed in a sterile environment.*
- 7. Data analysis would be problematic if there is high amount of stage drift between initial acquisition and cell detachment, as each point would have to be corrected for. However, with adequate software this drift should be able to be corrected post acquisition. Some discussion of difficulties in multiwell plate analysis would be helpful.*

**Author replies:**

1. Beads are polydisperse and approximately 300nm- this has been added to the manuscript
2. A 10x air objective is used, this has been added to the manuscript. Most immersion objectives are not suitable due to a) the small working distance offered, requiring the substrate and plate to be thinner increasing fragility and bending, b) the necessity of immersion fluid being applied over a large surface causes problems with dripping. Employing immersion optics and higher magnification could provide up to approximately 3-fold improvements in resolution, but would not result in qualitative changes, and is better suited for lower throughput sub-cellular studies.
3. We have added a note in the text- in our application we have used 10-15 min at a distance of 5 cm. These values may vary depending on effective illumination power.
4. We regret the omission. Thank you- Table 1 has been added.
5. The thickness of the glass does not interfere with 10x objective usage. We have added notes in the manuscript on image resolution and size.
6. The PDMS may be autoclaved, however, this damages the fluorescence of the beads, and catastrophically damages the well-divider employed here. The plate is recommended to be sterilized under biosafety cabinet UV light as described in Section 3.1
7. Dedrift software accommodates for this, provided that the drift distance is not overly large; in our usage, 90 microns seems to be a practical limit. This limit may be increased by imaging a larger field of view, and increasing the search window size.