

Journal of Visualized Experiments

Measurement of Interaction Force Between Droplet and Super-hydrophobic Substrate by Optical Lever Method --Manuscript Draft--

Article Type:	Methods Article - JoVE Produced Video
Manuscript Number:	JoVE59539R2
Full Title:	Measurement of Interaction Force Between Droplet and Super-hydrophobic Substrate by Optical Lever Method
Keywords:	Super-hydrophobic, interaction force, cantilever, optical lever, calibration, electrostatic force
Corresponding Author:	Yelong Zheng Tianjin University Tianjin, Tianjin CHINA
Corresponding Author's Institution:	Tianjin University
Corresponding Author E-Mail:	zhengyelongby@163.com
Order of Authors:	Shuya Zhuang Meirong Zhao Zhiyi Wang Lele Zhang Yinguo Huang Yelong Zheng
Additional Information:	
Question	Response
Please indicate whether this article will be Standard Access or Open Access.	Standard Access (US\$2,400)
Please indicate the city, state/province, and country where this article will be filmed . Please do not use abbreviations.	Tianjin, China

Dear editor of *JoVE*,

We *submit* a manuscript titled “Measurement of Interaction Force Between Droplet and Super-hydrophobic Substrate by Optical Lever Method” to *JoVE*.

The paper is dedicated to investigating the interaction between droplets and super-hydrophobic substrates in the air. The protocol includes calibrating the measurement system and measuring the interaction force at super-hydrophobic substrates with different grid fractions.

Thanks a lot for considering our manuscript to be published on *JoVE*.

Yours sincerely!

Dr. Yelong Zheng

State Key Laboratory of Precision Measuring Technology and Instruments

Tianjin University

Tianjin 300072, China.

Tel: +86 15122346153

Email: zhengyelongby@163.com

Dec. 12. 2018

TITLE:

Measuring the Interaction Force Between a Droplet and a Super-hydrophobic Substrate by the Optical Lever Method

AUTHORS AND AFFILIATIONS:

Shuya Zhuang, Meirong Zhao, Zhiyi Wang, Lele Zhang, Yinguo Huang, Yelong Zheng

State Key Laboratory of Precision Measuring Technology and Instruments, Tianjin University, Tianjin, China.

Corresponding Author:

Yelong Zheng (zhengyelongby@163.com)

Email Addresses of Co-authors:

Shuya Zhuang (syzhuang@tju.edu.cn)

Meirong Zhao (meirongzhao@tju.edu.cn)

Zhiyi Wang (zhiyi_wang@qq.com)

Lele Zhang (1147485882@qq.com)

Yinguo Huang (hyg@tju.edu.cn)

KEYWORDS:

Super-hydrophobic, interaction force, cantilever, optical lever, calibration, electrostatic force

SUMMARY:

The protocol aims to investigate the interaction between droplets and super-hydrophobic substrates in the air. This includes calibrating the measurement system and measuring the interaction force at super-hydrophobic substrates with different grid fractions.

ABSTRACT:

This goal of the paper is to investigate the interaction force between droplets and super-hydrophobic substrates in the air. A measurement system based on an optical lever method is designed. A millimetric cantilever is used as a force sensitive component in the measurement system. Firstly, the force sensitivity of the optical lever is calibrated using electrostatic force, which is the critical step in measuring interaction force. Secondly, three super-hydrophobic substrates with different grid fractions are prepared with nanoparticles and copper grids. Finally, the interaction forces between droplets and super-hydrophobic substrates with different grid fractions are measured by the system. This method can be used to measure the force on the scale of sub-micronewton with a resolution on the scale of nanonewton. The in-depth study of the contact process of droplets and super-hydrophobic structures can help to improve the production efficiency in coating, film and printing. The force measurement system designed in this paper can also be used in other fields of microforce measurement.

INTRODUCTION:

The contact between a droplet and a super-hydrophobic surface is very common in daily life and

industrial production: water droplets sliding from the surface of lotus leaf^{1,2}, and a water strider traveling rapidly over the water³⁻⁶. A super-hydrophobic coating on the exterior surface of a ship can help reduce the corrosion degree of the ship and reduce the resistance of the navigation⁷⁻¹⁰. There is great value for industrial production and bionics research in studying the contact process between a droplet and a super-hydrophobic surface.

To observe the spreading process of droplets on a solid surface, Biance used a high-speed camera to photograph the contact process and found that the duration of the inertial regime is mainly fixed by the drop size¹¹. Eddi photographed the contact process between the droplet and the transparent plate from the bottom and side using a high-speed camera, which comprehensively revealed the variation of the contact radius of the viscous droplet with time¹². Paulsen combined an electrical method with high-speed camera observation, thus reducing the response time to 10 ns^{13,14}.

Atomic force microscopy (AFM) has also been used to measure the interaction force between the droplet/bubble and solid surfaces. Vakarelski used an AFM cantilever to measure the interaction forces between two small bubbles (approximately 80-140 μm) in aqueous solution during controlled collisions on the scale of micrometers to nanometers¹⁵. Shi used a combination of AFM and reflection interference contrast microscopy (RICM) to simultaneously measure the interaction force and the spatiotemporal evolution of the thin water film between an air bubble and mica surfaces of different hydrophobicity^{16,17}.

However, since commercial cantilevers used in AFM are too small, the laser spot irradiated on the cantilever would be submerged by droplets or bubbles. The AFM has difficulties in measuring the interaction force between droplets and droplets/substrates in the air.

In this paper, a measurement system based on an optical lever method is designed to measure the interaction force between droplets and super-hydrophobic substrates. The force sensitivity of the optical lever (S_{OL}) is calibrated by electrostatic force¹⁸, and then the interaction forces between droplets and different super-hydrophobic substrates are measured by the measurement system.

The schematic diagram of the measurement system is shown in **Figure 1**. The laser and position sensitive detector (PSD) constitute the optical lever system. A millimetric silicon cantilever is used as a sensitive component in the system. The substrate is fixed on the nanopositioning z-stage, which can move in vertical direction. When the substrate approaches the droplet, the interaction force causes the cantilever to bend. Thus, the position of the laser spot on PSD will change, and the output voltage of PSD will change. The output voltage of PSD V_p is proportional to the interaction force F_i , as shown in Eq. (1).

$$F_i = S_{OL} V_p \quad (1)$$

In order to acquire the interaction force, S_{OL} must be calibrated first. The electrostatic force is used as the standard force in the calibration of S_{OL} . As shown in **Figure 2**, the cantilever and the electrode make up a parallel plate capacitor, which could generate electrostatic force in a vertical

direction. The electrostatic force F_{es} is determined by the voltage of the DC power supply V_s , as shown in Eq. (2)¹⁹⁻²¹.

$$F_{es} = \frac{1}{2} V_s^2 \frac{dC}{dz} \quad (2)$$

where C is the capacitance of the parallel plate capacitor, z is the displacement of the cantilever free end, and dC/dz is called capacitance gradient. The capacitance could be measured by the capacitance bridge. The mathematical relationship between C and z can be fitted by a quadratic polynomial, as shown in Eq. (3).

$$C = Qz^2 + Pz + CT \quad (3)$$

where Q , P and CT are the coefficients of the quadratic term, the primary term and the constant term respectively. Therefore, the electrostatic force F_{es} can be expressed as Eq. (4).

$$F_{es} = \frac{1}{2} V_s^2 (2Qz + P) \quad (4)$$

Since the overlap area of two plates of the capacitor is very small, the elastic force acted on the cantilever can be expressed as Eq. (5), according to Hooke's law:

$$F_i = kz \quad (5)$$

where k is the stiffness of the cantilever.

When the elastic force and electrostatic force applied on the cantilever are equal (i.e., $F_i = F_{es}$), the cantilever is in equilibrium. Eq. (6) can be derived from Eqs. (1), (2) and (5):

$$2Q \frac{S_{OL}}{k} V_P V_s^2 + P V_s^2 = 2S_{OL} V_P \quad (6)$$

In order to reduce the uncertainty of calibration results, a difference method is used to calculate S_{OL} . The results of two experiments are taken as V_{s1} , V_{p1} and V_{s2} , V_{p2} , and are substituted into Eq. (6):

$$\begin{cases} 2Q \frac{S_{OL}}{k} V_{P1} V_{s1}^2 + P V_{s1}^2 = 2S_{OL} V_{P1} \\ 2Q \frac{S_{OL}}{k} V_{P2} V_{s2}^2 + P V_{s2}^2 = 2S_{OL} V_{P2} \end{cases} \quad (7)$$

Transforming the equations and subtracting the lower equation from the upper equation in Eq. (7), the parameters Q and k are eliminated. Then the calibration formula of S_{OL} is obtained, as shown in Eq. (8):

$$2S_{OL} \left(\frac{1}{V_{s1}^2} - \frac{1}{V_{s2}^2} \right) = P \left(\frac{1}{V_{P1}} - \frac{1}{V_{P2}} \right) \quad (8)$$

Performing a series of experiments, the curve is drawn with $P(1/V_{p1}-1/V_{p2})$ as the ordinate and $2(1/V_{s1}^2-1/V_{s2}^2)$ as the abscissa. The slope of the curve is S_{OL} .

After obtaining S_{OL} , the electrode will be replaced by different super-hydrophobic substrates. The interaction forces between droplets and super-hydrophobic substrates will be measured by the system shown in **Figure 1**.

129 130 **PROTOCOL:**

131 132 **1. Assembly of the S_{OL} calibration system**

133
134 1.1 Assemble the S_{OL} calibration system according to the schematic diagram shown in **Figure**
135 **2.**

136
137 1.2 Fix the laser to a support, making the angle between the laser and the horizontal direction
138 be 45° .

139
140 1.3 Fix the PSD to another support, making the PSD perpendicular to the laser. Connect the
141 PSD to the data acquisition device and the data acquisition device to the computer.

142
143 NOTE: These angles are determined by the experimenter's visual measurement and need not be
144 exactly 45° or 90° .

145
146 1.4 Fix the wider end of the cantilever to a holding device while the other end is suspended.
147 Fix the holding device to a two-dimensional high-precision displacement stage.

148
149 NOTE: The dimensions of the cantilever are shown in **Figure 3.**

150
151 1.5 Fix the plate electrode to the nanopositioning z-stage by a clamping device.

152
153 NOTE: The nanopositioning z-stage can bring the electrode to move along the z-axis with a
154 displacement resolution of 1 nm.

155
156 1.6 Connect the positive pole of the capacitive bridge with the cantilever and the negative
157 pole with the plate electrode.

158
159 1.7 Install a high-speed camera, whose line of sight is perpendicular to the cantilever.

160
161 1.8 Adjust the position of the plate electrode, making the vertical distance between the plate
162 electrode and the cantilever be about $100\ \mu\text{m}$, and the overlap length about 0.5 mm.

163
164 NOTE: These distances are checked by image processing.

165 166 **2. Measurement of the capacitance gradient**

167
168 2.1 Use the computer to control the capacitance bridge to collect the capacitance changes
169 between the plate electrode and the cantilever in real time. Set the sampling rate to 0.5 Hz.

170
171 2.2 Control the nanopositioning z-stage by computer to drive the plate electrode to step
172 upward with a step of $10\ \mu\text{m}$ and a step number of 6 and stay for 10 s after each movement.

- 173
- 174 2.3 Change the movement direction of the plate electrode to downward, and repeat step 2.2.
- 175
- 176 2.4 Determine the relationship between the capacitance and displacement of the plate
- 177 electrode in the measurement result, and obtain the value of P according to the Eq. (3).
- 178
- 179 2.5 Repeat steps 2.1–2.4 5x and calculate the average value of P .
- 180

181 3. Calibration of the optical lever

182

183 3.1 Disconnect the connection between the capacitive bridge and the cantilever/plate

184 electrode.

185

186 3.2 Connect the positive pole of the DC power supply with the cantilever and the negative

187 pole with the plate electrode.

188

189 3.3 Adjust the relative position between the laser, PSD and cantilever to make the laser

190 reflected on PSD by cantilever.

191

192 NOTE: The laser spot is a circle about 2 mm in diameter.

193

194 3.4 Control the DC power supply by computer to apply voltage varying with time on the

195 parallel plate capacitor. At the same time, collect the output voltage of PSD in real time by the

196 data acquisition device.

197

198 3.4.1 Set the sampling rate of the data acquisition device to 1,000 Hz.

199

200 3.4.2 Set the initial voltage of the DC power supply to 0 V and hold for 5 s.

201

202 3.4.3 Increase the voltage by 25 V and hold for 5 s.

203

204 3.4.4 Repeat step 3.4.2 4x until the voltage increases to 125 V.

205

206 3.4.5 Decrease the voltage by 25 V and hold for 5 s.

207

208 3.4.6 Repeat step 3.4.4 4x until the voltage decreases to 0 V.

209

210 3.5 Determine the relationship about the output voltage of PSD and the DC supply voltage in

211 the measurement result, and obtain the value of S_{OL} according to Eq. (8).

212

213 3.6 Repeat steps 3.4–3.5 5x and calculate the average value of S_{OL} .

214

215 4. Preparation of super-hydrophobic substrates

216

4.1 Prepare three circular copper grids with the same diameter of 3 mm and different grid fractions. Their grid fractions are 46.18%, 51.39% and 58.79% respectively.

NOTE: These copper grids are commercial products that were purchased.

4.2 Spray the nanoparticles onto three copper grids to obtain super-hydrophobic substrates with micro structure and hydrophobicity.

4.2.1 Spray the base coat onto the copper grid.

4.2.2 Spray the top coat onto the copper grid when the first coat is dry.

NOTE: The nanoparticles are packed in a can with a spray head. The nanoparticles will be sprayed by pressing the spray head when used.

4.3 Glue the copper grids on the side of cylinders with a diameter of 3 mm to obtain a surface super-hydrophobic structure with a curvature of $1/3 \text{ mm}^{-1}$.

5. Measurement of the interaction force between droplets and super-hydrophobic substrates

5.1 Disconnect the connection between the DC power supply and the cantilever/plate electrode. Remove the plate electrode from the nanopositioning z-stage.

5.2 Fix a plate support to the nanopositioning z-stage.

5.3 Install a high-speed camera, whose line of sight is perpendicular to the cantilever.

5.4 Suspend a droplet on the lower surface of the free end of the cantilever.

5.4.1 Place a super-hydrophobic structure with a contact angle of almost 180° on the plate support.

5.4.2 Place a $2 \mu\text{L}$ droplet on the super-hydrophobic structure using a micropipettor.

5.4.3 Control the nanopositioning z-stage using software (e.g., PIMikroMove) to drive the droplet to move upward.

5.4.3.1 In the dialog box, set the velocity to $10 \mu\text{m/s}$.

5.4.3.2 Click the **Forward** button and the droplet begins to move upwards.

5.4.3.3 Click the **Stop** button when the droplet contacts with the free end of the cantilever.

5.4.4 Stay for 1 or 2 s, and then control the nanopositioning z-stage to drive the slide glass away

from the cantilever.

NOTE: Since the silicon cantilever is hydrophilic, the droplet is suspended on the lower surface of the free end of the cantilever, forming a hemispherical droplet with a diameter of about 0.5 mm.

5.4.5 Remove the super-hydrophobic structure with a contact angle of almost 180° from the plate support.

5.5 Place the super-hydrophobic substrate with a grid fraction of 46.18% on the plate support.

5.6 Adjust the position of the plate support, making the vertical distance between the super-hydrophobic substrate and the hemispherical droplet be about 100 μm .

NOTE: The distance is checked by image processing.

5.7 Turn on the PSD, laser, and high-speed camera.

5.8 Control the data acquisition device by computer to collect the output voltage of PSD in real time. Set the sampling rate to 100 kHz.

5.9 Set the velocity to 10 $\mu\text{m/s}$ in the software, and then click the **Forward** button, so that the super-hydrophobic substrate moves gradually closer to the droplet.

5.10 Click the **Stop** button when the super-hydrophobic substrate and the droplet contact.

5.11 Set the velocity to 10 $\mu\text{m/s}$ in the software, and then click the **Back** button to drive the super-hydrophobic substrate to move downward.

5.12 Click the **Stop** button when the super-hydrophobic substrate is separated from the droplet.

5.13 Draw the curve of the output voltage of PSD varying with time.

5.14 Repeat steps 5.3–5.11 using super-hydrophobic substrates with grid fractions of 51.39% and 58.79%.

5.15 Analyze the relationship between the interaction force and the grid fraction of super-hydrophobic substrate.

REPRESENTATIVE RESULTS:

The displacement of the plate electrode and the corresponding capacitance between the cantilever and the electrode measured in one experiment are shown in **Table 1**. The relationship between capacitance C and displacement z is fitted by quadratic polynomial using the polyfit

function in MATLAB, as shown in **Figure 4**. The first order coefficient P can be obtained by the fitting function. The final value of P is 0.2799 pF/mm, which is the average calculated from ten experimental results.

The supply voltage and the corresponding output voltage of PSD in one experiment are shown in **Table 2**. The relationship between the output voltage of PSD V_p and the supply voltage V_s is fitted by linear function, as shown in **Figure 5**, where V_{p1} and V_{p2} are the differences between the measured values and V_{p0} (the initial value of the output voltage of PSD at $V_s=0$). S_{OL} can be obtained by the slope of the curve in **Figure 5**. The final value of S_{OL} is 8.847 $\mu\text{N/V}$, which is the average calculated from twelve experimental results.

The curves of measured interaction forces between droplets and super-hydrophobic substrates varying with time are shown in **Figure 6**. The interaction forces are calculated by Eq. (1), where V_p are the differences between the measured output voltages of PSD and initial output voltages of PSD.

Before point A, the substrate has not contacted with the droplet, so the interaction force is 0. In stage AB, the distance between the substrate and the droplet is very small. Due to the influence of aerodynamic, there will be a repulsive force between the substrate and the droplet, which shows a rising curve in the figure. Point B is the critical point at which the substrate and droplet begin to contact. After point B, the interaction force between them becomes attractive force. In stage BC, the droplet gradually wets the super-hydrophobic substrate under the action of capillary force. The cantilever will bend downward in this stage, showing a decreasing curve in the figure. At point C, the system reaches equilibrium again, and the cantilever begins to oscillate in the equilibrium position.

As shown in **Figure 6**, the interaction force between the droplet and the substrate decreases with the increase of the grid fraction. The reason is that the contact between droplets and super-hydrophobic substrates is a process of releasing energy. The hydrophobicity of the substrate is positively correlated with the grid fraction. The stronger the hydrophobicity is, the less energy released during contact, so the contact force is smaller.

During the experiment we found that the repulsive force only exists in the contact process between the droplet and the substrate with a grid fraction of 46.18%. With the increase of hydrophobicity, the surface energy of the substrate decreases. When the repulsive force cannot reach the resolution of the system, it is difficult to measure the repulsive force.

The magnitude of the force is in a direct relationship with the droplet volume. A supplementary experiment was conducted to illustrate the relationship between the interaction force and droplet volume. Three droplets of different sizes were used in the contact experiment, as shown in **Figure 7**. The volume of droplet (a), (b), and (c) are 0.0135 μL , 0.0087 μL , and 0.0073 μL , respectively. In the experiment, the volume of droplets is measured by the change of PSD output voltage. The PSD output voltages before and after suspension of the droplet by cantilever are measured, and their difference V_d is multiplied by S_{OL} to obtain the gravity of droplet. The volume

value of the droplet is converted by gravity. The substrate with a grid fraction of 51.39% is used in the experiments. The measured interaction force between the three droplets and the substrate are shown in **Figure 8**. It is obvious that the interaction force increases with the increase of droplet volume.

FIGURE AND TABLE LEGENDS:

Figure 1: The interaction force measurement system. The schematic diagram of the measurement system based on optical lever method for measuring the interaction between droplets and super-hydrophobic substrates.

Figure 2: The calibration system of S_{OL} . The schematic diagram of the system designed for calibrating the force sensitivity of optical lever using electrostatic force.

Figure 3: The dimensions of the millimetric cantilever. Top view and elevation view of the millimetric cantilever.

Figure 4: The relationship between capacitance and displacement. The quadratic polynomial fitting curve of capacitance C and displacement z in the calibration experiment.

Figure 5: The relationship between the output voltage of PSD and the supply voltage. The linear fitting curve of $P(1/V_{p1}-V_{p2})$ and $2(1/V_{s1}^2-V_{s2}^2)$ according to Eq. (8).

Figure 6: The measurement results of interaction force. The interaction forces between droplets and super-hydrophobic substrates with different grid fractions in the air. Different colors represent different grid fractions.

Figure 7: Images of three drops of different volumes used in the experiment. The volume of droplets (a), (b), and (c) are $0.0135\mu\text{L}$, $0.0087\mu\text{L}$, and $0.0073\mu\text{L}$ respectively.

Figure 8: Interaction force between the three droplets of different volumes and the substrate. Different colors represent different droplets.

Table 1: The calibration results of capacitance gradient. The displacement of the plate electrode and the corresponding capacitance between the cantilever and the electrode in one experiment.

Table 2: The calibration results of S_{OL} . The supply voltage and the corresponding output voltage of PSD in one experiment.

DISCUSSION:

In this protocol, a measurement system based on optical lever method is assembled and calibrated, which is designed for measuring the interaction force between the droplets and super-hydrophobic substrates. Among the all steps, it is critical to calibrate S_{OL} using electrostatic force. The results of the calibration experiment verify Eq. (8): $P(1/V_{p1}-1/V_{p2})$ is proportional to $2(1/V_{s1}^2-1/V_{s2}^2)$ and make it possible to obtain the value of the force to be measured through the

output voltage of the PSD. Through the experiment of measuring the interaction force between the droplets and super-hydrophobic substrates of different hydrophobicity, the interaction force decreases with the increase of hydrophobic capacity, which verifies the relationship between hydrophobicity and surface energy of the substrates.

The force measurement method based on a millimetric silicon cantilever is an important complement to the traditional methods. Compared with the high-speed camera method, the optical lever method can accurately measure the force on the nanonewton scale. AFM is usually used to measure the interaction force between micron-scale objects, while the system designed in this paper can be applied on objects of millimeter-scale in the air. This method can be used to measure the force in the scale of sub-micronewtons, and its resolution can reach the scale of nanonewtons.

The method proposed for measuring interaction force in this paper is limited to a small measurement range. A large force will lead to plastic deformation or even break of the silicon cantilever, which will cause incorrect results. In addition, because the principle of this experiment is to measure the interaction force between the droplet and the substrate under the balance of the elastic force of the cantilever and the interaction force, the cantilever can only measure the quasi-static force, but not dynamic force.

The in-depth study of the contact process of droplets and super-hydrophobic structures can help people to improve production efficiency in coating, film, printing and other industrial production. As a generalized adhesion force measurement technique, the substrates in the system can be replaced with substrates made of other materials. For example, a super-hydrophobic substrate with multi-stage microstructures which is made of PDMS (polydimethylsiloxane) can be used. The force measurement system based on optical lever method can also be used in other fields of micro force measurement, such as the interaction force during the coalescence of two droplets and the interaction force between super-hydrophobic substrates and droplets of different surface tension.

ACKNOWLEDGMENTS:

The authors thank the Tianjin Natural Science Foundation (No. 18JCQNJC04800), Tribology Science Fund of State Key Laboratory of Tribology (No. SKLTKF17B18) and National Natural Science Foundation of China (Grant No. 51805367) for their support.


DISCLOSURES:


The author has nothing to disclose.


REFERENCES:


1. Guo, Z., Liu, W. Biomimic from the superhydrophobic plant leaves in nature: Binary structure and unitary structure. *Plant Science*. **172** (6), 1103-1112 (2007).
2. Koch, K., Bhushan, B., Barthlott, W. Diversity of structure, morphology and wetting of plant surfaces. *Soft Matter*. **4** (10), 1943-1963 (2008).


3. Gao, X., Jiang, L. Biophysics: Water-repellent legs of water striders. *Nature*. **432** (7013), 36 (2004).
4. Yin, W., Zheng, Y.L., Lu, H.Y. Three-dimensional topographies of water surface dimples formed by superhydrophobic water strider legs. *Applied Physics Letters*. **109** (16), 163701 (2016).
5. Zheng, Y. et al. Elegant Shadow Making Tiny Force Visible for Water-Walking Arthropods and Updated Archimedes' Principle. *Langmuir*. **32** (41), 10522-10528 (2016).
6. Hu, D.L., Chan, B, Bush, J.W.M. The hydrodynamics of water strider locomotion. *Nature*. **424** (6949), 663-666 (2003).
7. Jiang, C.G., Xin, S.C., Wu, C.W. Drag reduction of a miniature boat with superhydrophobic grille bottom. *AIP Advances*. **1** (3), 032148 (2011).
8. Guo, Z., Liang, J., Fang, J., Guo, B., Liu, W. A Novel Approach to the Robust Ti6Al4V-Based Superhydrophobic Surface with Crater-like Structure. *Advanced Engineering Materials*. **9** (4), 316-321 (2007).
9. Guo, Z., Liu, W., Su, B. A stable lotus-leaf-like water-repellent copper. *Applied Physics Letters*, **92** (6), 063104 (2008).
10. Liu, K., Zhang, M., Zhai, J., Wang, J., Jiang, L. Bioinspired construction of Mg-Li alloys surfaces with stable superhydrophobicity and improved corrosion resistance. *Applied Physics Letters*. **92** (18), 61 (2008).
11. Biance, A.L., Clanet, C., Quéré, D. First steps in the spreading of a liquid droplet. *Physical Review E Statistical Nonlinear & Soft Matter Physics*. **69** (1), 016301 (2004).
12. Eddi, A., Winkels, K.G., Snoeijer, J.H. Short time dynamics of viscous drop spreading. *Physics of Fluids*. **25**(1), 77-177 (2017).
13. Paulsen, J.D., Burton, J.C., Nagel, S.R. Viscous to Inertial Crossover in Liquid Drop Coalescence. *Physical Review Letters*. **106** (11), 114501 (2011).
14. Paulsen, J.D. Approach and coalescence of liquid drops in air. *Physical Review E*. **88** (6), 063010 (2013).
15. Vakarelski, I.U. et al. Bubble colloidal AFM probes formed from ultrasonically generated bubbles. *Langmuir*. **24** (3), 603-605 (2008).
16. Shi, C. et al. Measuring forces and spatiotemporal evolution of thin water films between an air bubble and solid surfaces of different hydrophobicity. *ACS Nano*. **9** (1), 95-104 (2015).
17. Shi, C., Chan, D.Y.C., Liu, Q., Zeng, H. Probing the Hydrophobic Interaction between Air Bubbles and Partially Hydrophobic Surfaces Using Atomic Force Microscopy. *The Journal of Physical Chemistry C*. **118** (43), 25000-25008 (2014).
18. Chung, K., Shaw, G.A., Pratt, J.R. Accurate noncontact calibration of colloidal probe sensitivities in atomic force microscopy. *Review of Scientific Instruments*. **80** (6), 930 (2009).
19. Zheng, Y. et al. Improving environmental noise suppression for micronewton force sensing based on electrostatic by injecting air damping. *Review of Scientific Instruments*. **85** (5), 055002 (2014).
20. Song, L. et al. Highly sensitive, precise, and traceable measurement of force. *Instrumentation Science & Technology*. **44** (4), 15 (2016).
21. Zheng, Y. et al. The multi-position calibration of the stiffness for atomic-force microscope cantilevers based on vibration. *Measurement Science & Technology*. **26** (5), 055001 (2012).



















Time(s)	0	10	20	30	40	50	60
Displacement(μm)	0	10	20	30	40	50	60
Capacitance(pF)	2.399	2.402	2.406	2.411	2.416	2.422	2.429
Time(s)	70	80	90	100	110	120	
Displacement(μm)	50	40	30	20	10	0	
Capacitance(pF)	2.422	2.416	2.411	2.407	2.403	2.400	

Time(s)	0	5	10	15	20	25
Supply Voltage(V)	0	25	50	75	100	125
Output voltage of PSD(V)	-3.5757	-3.5656	-3.5327	-3.4797	-3.3775	-3.1733
Time(s)	30	35	40	45	50	
Supply Voltage(V)	100	75	50	25	0	
Output voltage of PSD(V)	-3.3765	-3.4786	-3.5321	-3.5644	-3.5755	

Name of Material/ Equipment	Company	Catalog Number	Comments/Description
Camera	Shenzhen Andonstar Tech Co., Ltd	digital microscope A1	Frame rate: 30 frames/sec; Focal distance: 5 mm - 30 mm
Capacitive bridge	Andeen-Hagerling	AH2550A	The capacitive bridge is used to measure the capacitance between the cantilever and the plate electrode.
Data acquisition device	National Instruments	USB-4431	The data acquisition device is used to read the output voltage data.
DC power supply	Keithley	2410	Voltage range: $\pm 5\text{ }\mu\text{V}$; Accuracy: 0.012%
Grid	Electron Microscopy China	AGH100, AGH150, AGH300	The grid fractions of AGH100, AGH150 and AGH300 are 46.18%, 51.39% and 58.79% respectively
Laser	Shenzhen Infrared Laser Technology Co., Ltd.	HW650AD100-10BD	Laser wavelength: 650 nm
Nanoparticle	Rust-Oleum	274232	NeverWet Multi-Surface Liquid Repelling Treatment is a revolutionary super hydrophobic coating.
Nanopositioning z- stage	Physik Instrumente	P622.ZCD	Travel ranges 50 μm to 250 μm (350 μm open loop); Resolution to 0.1 nm; Linearity error only 0.02%
Position sensitive detector	Hamamatsu Photonics K.K.	S1880	The two-dimensional PSD is used to translate optical signals into electrical signals.



1 Alewife Center #200
Cambridge, MA 02140
tel. 617.945.9051
www.jove.com

ARTICLE AND VIDEO LICENSE AGREEMENT

Title of Article:

Measurement of Interaction Force Between Droplet and Super-hydrophobic Substrate by Optical Lever Method

Author(s):

Shuya Zhuang, Meirong Zhao, Zhiyi Wang, Lele Zhang, Yinguo Huang, Yelang Zheng

Item 1: The Author elects to have the Materials be made available (as described at <http://www.jove.com/publish>) via:

☒ Standard Access

☐ Open Access

Item 2: Please select one of the following items:

☒ The Author is **NOT** a United States government employee.

☐ The Author is a United States government employee and the Materials were prepared in the course of his or her duties as a United States government employee.

☐ The Author is a United States government employee but the Materials were NOT prepared in the course of his or her duties as a United States government employee.

ARTICLE AND VIDEO LICENSE AGREEMENT

1. **Defined Terms.** As used in this Article and Video License Agreement, the following terms shall have the following meanings: **"Agreement"** means this Article and Video License Agreement; **"Article"** means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein; **"Author"** means the author who is a signatory to this Agreement; **"Collective Work"** means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; **"CRC License"** means the Creative Commons Attribution-Non Commercial-No Derivs 3.0 Unported Agreement, the terms and conditions of which can be found at: <http://creativecommons.org/licenses/by-nc-nd/3.0/legalcode>; **"Derivative Work"** means a work based upon the Materials or upon the Materials and other pre-existing works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; **"Institution"** means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; **"JoVE"** means MyJoVE Corporation, a Massachusetts corporation and the publisher of The Journal of Visualized Experiments; **"Materials"** means the Article and / or the Video; **"Parties"** means the Author and JoVE; **"Video"** means any video(s) made by the Author, alone or in conjunction with any other parties, or by JoVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion

of the Article, and in which the Author may or may not appear.

2. **Background.** The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JoVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.

3. **Grant of Rights in Article.** In consideration of JoVE agreeing to publish the Article, the Author hereby grants to JoVE, subject to **Sections 4 and 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the "Open Access" box has been checked in **Item 1** above, JoVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.

ARTICLE AND VIDEO LICENSE AGREEMENT

4. **Retention of Rights in Article.** Notwithstanding the exclusive license granted to JoVE in **Section 3** above, the Author shall, with respect to the Article, retain the non-exclusive right to use all or part of the Article for the non-commercial purpose of giving lectures, presentations or teaching classes, and to post a copy of the Article on the Institution's website or the Author's personal website, in each case provided that a link to the Article on the JoVE website is provided and notice of JoVE's copyright in the Article is included. All non-copyright intellectual property rights in and to the Article, such as patent rights, shall remain with the Author.

5. **Grant of Rights in Video – Standard Access.** This **Section 5** applies if the "Standard Access" box has been checked in **Item 1** above or if no box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby acknowledges and agrees that, Subject to **Section 7** below, JoVE is and shall be the sole and exclusive owner of all rights of any nature, including, without limitation, all copyrights, in and to the Video. To the extent that, by law, the Author is deemed, now or at any time in the future, to have any rights of any nature in or to the Video, the Author hereby disclaims all such rights and transfers all such rights to JoVE.

6. **Grant of Rights in Video – Open Access.** This **Section 6** applies only if the "Open Access" box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby grants to JoVE, subject to **Section 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Video in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Video into other languages, create adaptations, summaries or extracts of the Video or other Derivative Works or Collective Works based on all or any portion of the Video and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. For any Video to which this **Section 6** is applicable, JoVE and the Author hereby grant to the public all such rights in the Video as provided in, but subject to all limitations and requirements set forth in, the CRC License.

7. **Government Employees.** If the Author is a United States government employee and the Article was prepared in the course of his or her duties as a United States government employee, as indicated in **Item 2** above, and any of the licenses or grants granted by the Author hereunder exceed the scope of the 17 U.S.C. 403, then the rights granted hereunder shall be limited to the maximum

rights permitted under such statute. In such case, all provisions contained herein that are not in conflict with such statute shall remain in full force and effect, and all provisions contained herein that do so conflict shall be deemed to be amended so as to provide to JoVE the maximum rights permissible within such statute.

8. **Protection of the Work.** The Author(s) authorize JoVE to take steps in the Author(s) name and on their behalf if JoVE believes some third party could be infringing or might infringe the copyright of either the Author's Article and/or Video.

9. **Likeness, Privacy, Personality.** The Author hereby grants JoVE the right to use the Author's name, voice, likeness, picture, photograph, image, biography and performance in any way, commercial or otherwise, in connection with the Materials and the sale, promotion and distribution thereof. The Author hereby waives any and all rights he or she may have, relating to his or her appearance in the Video or otherwise relating to the Materials, under all applicable privacy, likeness, personality or similar laws.

10. **Author Warranties.** The Author represents and warrants that the Article is original, that it has not been published, that the copyright interest is owned by the Author (or, if more than one author is listed at the beginning of this Agreement, by such authors collectively) and has not been assigned, licensed, or otherwise transferred to any other party. The Author represents and warrants that the author(s) listed at the top of this Agreement are the only authors of the Materials. If more than one author is listed at the top of this Agreement and if any such author has not entered into a separate Article and Video License Agreement with JoVE relating to the Materials, the Author represents and warrants that the Author has been authorized by each of the other such authors to execute this Agreement on his or her behalf and to bind him or her with respect to the terms of this Agreement as if each of them had been a party hereto as an Author. The Author warrants that the use, reproduction, distribution, public or private performance or display, and/or modification of all or any portion of the Materials does not and will not violate, infringe and/or misappropriate the patent, trademark, intellectual property or other rights of any third party. The Author represents and warrants that it has and will continue to comply with all government, institutional and other regulations, including, without limitation all institutional, laboratory, hospital, ethical, human and animal treatment, privacy, and all other rules, regulations, laws, procedures or guidelines, applicable to the Materials, and that all research involving human and animal subjects has been approved by the Author's relevant institutional review board.

11. **JoVE Discretion.** If the Author requests the assistance of JoVE in producing the Video in the Author's facility, the Author shall ensure that the presence of JoVE employees, agents or independent contractors is in accordance with the relevant regulations of the Author's institution. If more than one author is listed at the beginning of this Agreement, JoVE may, in its sole

ARTICLE AND VIDEO LICENSE AGREEMENT

discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JoVE reserves the right, in its absolute and sole discretion and without giving any reason therefore, to accept or decline any work submitted to JoVE. JoVE and its employees, agents and independent contractors shall have full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JoVE has sole discretion as to the method of making and publishing the Materials, including, without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

12. **Indemnification.** The Author agrees to indemnify JoVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JoVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JoVE, making of videos by JoVE, or publication in JoVE or elsewhere by JoVE. The Author shall be responsible for, and shall hold JoVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to

the making of a video by JoVE its employees, agents or independent contractors. All sterilization, cleanliness or decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's expense. All indemnifications provided herein shall include JoVE's attorney's fees and costs related to said losses or damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JoVE, its employees, agents or independent contractors.

13. **Fees.** To cover the cost incurred for publication, JoVE must receive payment before production and publication of the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JoVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.

14. **Transfer, Governing Law.** This Agreement may be assigned by JoVE and shall inure to the benefits of any of JoVE's successors and assignees. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to me one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement is required per submission.

CORRESPONDING AUTHOR

Name:

Yelong Zheng

Department:

State Key Laboratory of Precision Measuring Technology and Instruments

Institution:

Tianjin University

Title:

Measurement of Interaction Force Between Droplet and Super-hydrophobic Substrate by Optical
Lever Method

Signature:

Yelong Zheng

Date:

2018.12.11

Please submit a **signed** and **dated** copy of this license by one of the following three methods:

1. Upload an electronic version on the JoVE submission site
2. Fax the document to +1.866.381.2236
3. Mail the document to JoVE / Attn: JoVE Editorial / 1 Alewife Center #200 / Cambridge, MA 02140

Editorial comments:

1. Could you include some of your response to Reviewer 1 (from the rebuttal) in the manuscript itself (including Figures)? Other readers may have similar concerns about the broad applicability of this procedure. If you do so, please clarify how you estimate the volume of droplets (area measurement?).

Re: Thanks for your comment. The response to Reviewer 1 has been included in the manuscript itself now, and the method to estimate the volume of droplets are clarified too. The new content is shown in blue.

2. It's a little difficult to tell what you want to film now; please re-highlight the protocol with only yellow.

Re: Thanks for your comment. The protocol to be filmed is re-highlighted with yellow now.