

Journal of Visualized Experiments

Isolation, culture, characterization, and differentiation of human muscle progenitor cells from the skeletal muscle biopsy procedure --Manuscript Draft--

Article Type:	Invited Methods Article - JoVE Produced Video
Manuscript Number:	JoVE59580R2
Full Title:	Isolation, culture, characterization, and differentiation of human muscle progenitor cells from the skeletal muscle biopsy procedure
Keywords:	Human adult stem cells; muscle progenitor cells; skeletal muscle biopsy; cell culture; satellite cell; skeletal muscle
Corresponding Author:	Anna Thalacker-Mercer UNITED STATES
Corresponding Author's Institution:	
Corresponding Author E-Mail:	aet74@cornell.edu
Order of Authors:	Brandon J. Gheller Jamie Blum Sharon Soueid-Baumgarten Erica Bender Benjamin D. Cosgrove Anna Thalacker-Mercer
Additional Information:	
Question	Response
Please indicate whether this article will be Standard Access or Open Access.	Standard Access (US\$2,400)
Please indicate the city, state/province, and country where this article will be filmed . Please do not use abbreviations.	Ithaca, New York, USA

TITLE:

Isolation, Culture, Characterization, and Differentiation of Human Muscle Progenitor Cells from the Skeletal Muscle Biopsy Procedure

AUTHORS AND AFFILIATIONS:

Brandon J. Gheller^{1,*}, Jamie Blum^{1,*}, Sharon Soueid-Baumgarten², Erica Bender¹, Benjamin D. Cosgrove², Anna Thalacker-Mercer¹

¹Division of Nutritional Sciences, Cornell University, Ithaca, NY, United States

²Meinig School of Biomedical Engineering, Cornell University, Ithaca, NY, United States

*These authors contributed equally.

Corresponding author:

Anna Thalacker-Mercer (aet74@cornell.edu)

Email addresses of co-authors:

Brandon J. Gheller (bjg223@cornell.edu)

Jamie Blum (jeb462@cornell.edu)

Sharon Soueid-Baumgarten (syb8@cornell.edu)

Erica Bender (eb572@cornell.edu)

Benjamin D. Cosgrove (bdc68@cornell.edu)

KEYWORDS:

human adult stem cells, muscle progenitor cells, skeletal muscle biopsy, cell culture, satellite cell, skeletal muscle

SUMMARY:

We present techniques for isolating, culturing, characterizing, and differentiating human primary muscle progenitor cells (*hMPCs*) obtained from skeletal muscle biopsy tissue. *hMPCs* obtained and characterized through these methods can be used to subsequently address research questions related to human myogenesis and skeletal muscle regeneration.

ABSTRACT:

The use of primary human tissue and cells is ideal for the investigation of biological and physiological processes such as the skeletal muscle regenerative process. There are recognized challenges to working with human primary adult stem cells, particularly human muscle progenitor cells (*hMPCs*) derived from skeletal muscle biopsies, including low cell yield from collected tissue and a large degree of donor heterogeneity of growth and death parameters among cultures. While incorporating heterogeneity into experimental design requires a larger sample size to detect significant effects, it also allows us to identify mechanisms that underlie variability in *hMPC* expansion capacity, and thus allows us to better understand heterogeneity in skeletal muscle regeneration. Novel mechanisms that distinguish the expansion capacity of

cultures have the potential to lead to the development of therapies to improve skeletal muscle regeneration.

INTRODUCTION:

Skeletal muscle is the largest organ system in the human body, accounting for 30–40% of whole body mass¹. In addition to its well-recognized role in locomotion, skeletal muscle maintains body temperature and posture, and plays a central role in whole body nutrient homeostasis. Research involving human participants, animals, and cell culture models are all valuable to address questions pertaining to skeletal muscle biology and regeneration. Isolation and culture of human primary muscle progenitor cells (*hMPCs*) provides a robust model that allows for cell culture techniques and manipulations to be applied to human samples. An advantage of using *hMPCs* is that they retain the genetic and metabolic phenotype from each donor^{2,3}. Maintenance of the donor phenotype allows researchers to examine inter-individual variation in the myogenic process. For example, we have employed our *hMPC* characterization method to identify age- and sex-related differences in *hMPC* population expansion capacity⁴.

The purpose of this protocol is to detail techniques to isolate, culture, characterize, and differentiate human primary muscle progenitor cells (*hMPCs*) from skeletal muscle biopsy tissue. Building on previous work that described *hMPCs* and identified potential cell surface makers for *hMPC* isolation^{5,6}, this protocol fills a critical gap in knowledge by linking the isolation to the characterization of *hMPCs*. Further, the detailed step-by-step instructions included in this protocol make *hMPC* isolation and characterization accessible to a broad scientific audience, including those with limited prior experience with *hMPCs*. Our protocol is among the first to describe use of an imaging cytometer to track cell populations. Newly designed imaging cytometers are state-of-the-art, high-throughput, and microplate-based, enabling live cell imaging, cell counting, and multichannel fluorescence analysis of all cells in each well of a culture vessel within minutes. This system allows for rapid quantification of dynamic changes in proliferation and viability of an entire cell population with only minimal disruption to the culture. For example, we are able to perform objective measures of confluence on successive days *in vitro* to determine growth kinetics of each culture derived from different donors. Many protocols in the literature, particularly those involving the differentiation of *MPCs*, require cells to reach a defined level of confluence before initiating differentiation or treatment⁷. Our method allows for objective determination of the confluence of each well in a culture vessel allowing researchers to initiate treatment in an unbiased, non-subjective manner.

In the past, a major limitation of using primary *hMPCs* was low yields that limit the number of cells available for experiments. We and others have shown the yield of *MPCs* from skeletal muscle biopsy tissue is 1–15 *MPCs* per milligram of tissue (**Figure 1**)⁸. Because our protocol allows for four passages of the cells prior to purification with fluorescence activated cell sorting (FACS), our cryopreserved *hMPC* yields, derived from small amounts of biopsy tissue (50–100 mg), are sufficient to address research aims where multiple experiments are required. Our FACS protocol produces a ~80% pure (Pax7 positive) *MPC* population, thus our protocol is optimized for both yield and purity.

PROTOCOL:

This protocol was approved by the Institutional Review Board at Cornell University. All participants were screened for underlying health conditions and gave informed consent.

1. Obtaining human muscle tissue via the skeletal muscle biopsy

1.1. Identify the vastus lateralis muscle via palpation, anatomical landmarks, and active contraction of the muscle.

1.2. Palpate the vastus lateralis by having the participant tighten their quadriceps muscle and find the belly of the muscle, about 1/3 of the way from the patella to the hip, and just ventrolateral to the femur.

1.3. Use a paper measuring tape to mark an area with a surgical marker about 4–6 inches from the top of the patella representing the belly of the muscle.

1.4. Shave any hair near the marked area and then cleanse with a surgical grade antiseptic that is used before medical procedures (For example: 2% chlorhexidine gluconate (CHG)/ 70% isopropyl alcohol (IPA)) .

1.5. Anesthetize the biopsy site with 1% lidocaine.

NOTE: Avoid injecting the lidocaine into the skeletal muscle.

1.6. Make a 4–5 mm incision pathway at the marked site.

1.7. After the area is anesthetized, push a Bergstrom needle with trocar through the fascia into the belly of the muscle (approximately 2–3 cm).

1.8. Draw muscle tissue into the Bergstrom needle by withdrawing the cutting trocar 2–3 cm while an assistant applies suction by withdrawing an attached 60 mL syringe plunger.

1.9. Close the cutting chamber.

1.10. Rotate the base of the biopsy needle a quarter turn and repeat the sequence of opening the needle, applying suction, and closing the needle.

NOTE: Each time the needle is closed, a cut of muscle tissue should be obtained inside the Bergstrom needle. During sample collection, the needle will remain in the thigh in the same vertical position while it is rotated longitudinally through 360° to obtain each sample. This sequence can be repeated 3–4 times.

1.11. Expel the collected biopsy tissue from the needle into a sterile nonadherent dressing pad.

1.12. Examine the tissue and remove any visible adipose and fibrotic tissues using sterile scalpel blades or tweezers.

1.13. Place a portion (50–100 mg) of tissue into a sterile tube containing CO₂-independent nutrient medium (e.g., Hibernate A).

NOTE: The weight of the tissue is determined as follows. First place the capped tube containing the nutrient medium on a scale and tare. Next place the tissue into the tube, recap the tube, and then weigh the tube containing the tissue in the nutrient medium.

1.14. Store the tube containing the tissue and nutrient medium at 4 °C until processing.

NOTE: Biopsy tissue should be processed within 48 hours of collection.

2. Isolating human muscle progenitor cells from biopsy tissue

2.1. Mince biopsy tissue.

2.1.1. Transfer the muscle tissue from the nutrient medium into a sterile Petri dish in a biological safety cabinet or similar sterile workplace.

2.1.2. Use sterile scalpel blades to mince tissue into ~1 mm³ pieces.

2.2. Wash minced biopsy tissue to remove residual debris.

2.2.1. Transfer the muscle tissue into a 15 mL conical tube containing 10 mL of phosphate-buffered saline (PBS) and allow the tissue to settle via gravity.

2.2.2. Aspirate the supernatant being careful not to disturb the muscle tissue.

2.2.3. Add 10 mL of PBS to the washed tissue and allow the muscle tissue to resettle via gravity.

2.2.4. Aspirate the supernatant being careful not to disturb the muscle tissue.

2.2.5. Add 10 mL of low glucose Dulbecco's modified Eagle medium (DMEM) and allow the tissue to resettle via gravity.

2.2.6. Aspirate the supernatant being careful not to disturb the muscle tissue.

2.3. Incubate the tissue in digest media.

2.3.1. Resuspend the muscle tissue in 3 mL of digest media (2 mg/mL collagenase D in low glucose DMEM).

NOTE: Collagenase D stock solutions can be prepared at 20 mg/mL in low glucose DMEM and stored at -80 °C, then diluted 1:10 in low glucose DMEM at time of use.

2.3.2. Incubate the tube containing the muscle tissue in a water bath at 37 °C for 30 min.

2.3.3. Triturate the muscle tissue suspension every 10 min using a 10 mL serological pipette.

2.3.4. Add 83 µL of additional digest media and 24 µL of dispase solution (0.76 mM calcium acetate, 1.39 mM sodium acetate, 13.91 mg/mL dispase II in low glucose DMEM). Return the suspension to the 37 °C water bath.

2.3.5. Triturate the muscle tissue suspension every 5–10 min with a wide bore pipet tip.

2.3.6. Continue trituration until a uniform slurry is achieved, and for no longer than 1.5 h to prevent over digestion.

NOTE: Sufficient mincing of the tissue and the enzymatic activity of collagenase and dispase can impact sample digestion. A uniform slurry about to pass through a normal pipette tip within minimal occlusion after 1.5 h. If a uniform slurry is not achieved after 1.5 h, double check that the tissue mincing is sufficient, and tissue is minced to the appropriate size (step 2.1.2). Additionally, the activity of the enzymes should be tested as lot to lot variation in enzymatic activity does occur.

2.4. Pellet cells and freeze in recovery media.

2.4.1. Add 6 mL of growth media (GM; 79% Ham's F12, 20% fetal bovine serum, 1% penicillin/streptomycin, 5 ng/mL basic fibroblast growth factor [bFGF], pH 7.4, passed through a 0.22 µm filter) to the muscle slurry.

2.4.2. Pipette the suspension through a 70 µm cell strainer into a 50 mL conical tube. Rinse the strainer with an additional 4 mL of GM.

NOTE: Sample can be transferred back into a 15 mL conical tube for ease of visualizing the pellet. A larger wash volume can be used if subsequent cell yields are low.

2.4.3. Centrifuge at 300 x g for 5 min at room temperature (RT).

2.4.4. Aspirate the supernatant being sure not to disrupt the pellet. Flick the tube to resuspend the pellet in the remaining media. If cells are not going to be cultured immediately proceed to step 2.4.5 for instructions on long-term storage, otherwise proceed to step 3.1.6.

NOTE: If immediate culture is desired, GM in cell culture plates should be pre-incubated (steps 3.1.1 and 3.1.2) before biopsy processing begins.

2.4.5. Add 1.5 mL of recovery cell culture freezing media (**Table of Materials**) to the recovered pellet.

2.4.6. Transfer the slurry into a cryovial. Place the cryovial in an isopropanol-controlled rate chiller overnight at -80 °C. Store cryovial at -80 °C until ready for culture.

NOTE: Cryovials containing slurry can be kept in liquid nitrogen for long term storage (greater than 1 month). Using an isopropanol-controlled rate chiller allows cells to be reproducibly frozen at a rate of 1 °C/ min resulting in optimal cell recovery rates as evidenced by viability after freeze and increased likelihood of attachment to cell culture vessels⁹. Alternatively, a beaker containing RT isopropanol, placed in a -80 °C freezer may be used.

3. Culturing human muscle progenitor cells

3.1. Thaw the muscle slurry.

3.1.1. Prior to thawing and seeding *hMPCs*, pre-equilibrate cell culture plates in GM. Add 250 µL of GM to 4 wells of a collagen-coated (collagen I), 24-well cell culture plate. Fill the remaining wells with 500 µL of sterile base media (i.e., F12) to prevent evaporation during long-term culture.

NOTE: For biopsies between 50 g and 100 g, four wells of a 24-well plate are appropriate. If the biopsy tissue is less than 50 g, two wells may be more appropriate. Residual debris and low cell number prevents the use of a cell counting technique at this step; therefore, the exact seeding density is not determined and varies slightly from biopsy to biopsy.

3.1.2. Incubate cell culture vessels with GM at 37 °C in 5% CO₂ for 1 h prior to seeding.

3.1.3. Remove cryovials from the -80 °C freezer and place in a 37 °C water bath; the liquid in the tube should be completely submerged in water.

3.1.4. When thawing is almost complete (~5 min), transfer the cryovial to a sterile laminar flow hood. Transfer the muscle slurry out of the cryovial tube and add to a 15 mL conical tube containing 6 mL of pre-warmed GM at a rate of 1 mL/min.

3.1.5. Centrifuge the muscle slurry at 300 x *g* for 5 min at RT.

NOTE: A swinging bucket centrifuge can make it easier to visualize small pellets at this step but is not essential.

3.1.6. Aspirate the supernatant being sure not to disrupt the pellet. Flick the tube to resuspend the pellet in the remaining media and resuspend in 1 mL of GM.

3.1.7. Transfer 250 μ L of the resuspended pellet (now considered an *hMPC* suspension) to each of the 4 wells of the pre-incubated 24-well cell culture plate.

3.2. Culture and passage *hMPCs*.

3.2.1. Maintain *hMPC* cultures in GM at 37 °C in 5% CO₂.

NOTE: Stocks of GM without the bFGF (i.e., F12, antibiotics, FBS) are prepared in batches and kept at 4 °C for up to two weeks. bFGF is added to GM on the day of use. Media containing bFGF can be stored for 48 h at 4 °C.

3.2.2. Twenty-four hours after isolation, aspirate the GM, gently wash the culture vessel 2x with pre-warmed GM, and add fresh GM.

NOTE: Twenty-four hours should provide *hMPCs* adequate time to attach; no *hMPCs* should be removed after washing. This step will also remove remaining debris that was generated during the cell harvest, making the vessel more amenable to accurate confluence scanning using an imaging cytometer (section 5 below). After this initial media change, GM is changed every 48 h.

3.2.3. Passage *hMPCs* when 70% confluence is achieved or when they have remained on the same culture dish for 10 days, whichever occurs first.

NOTE: 70% confluence is approximately 55,000 cells/cm².

3.2.4. To passage, add 250 μ L of pre-warmed trypsin to each well of a 24-well cell culture plate and incubate for ~5 min.

3.2.5. Tap the cell culture vessel on a firm surface to detach *hMPCs*. A light microscope can be used to verify that the *hMPCs* have detached.

3.2.6. Transfer the trypsin/*hMPC* suspensions to 5 mL of GM in a 15 mL conical tube.

3.2.7. Centrifuge the *hMPC* suspension at 300 x *g* for 5 min at RT.

3.2.8. Remove *hMPC* suspensions from the centrifuge and place on ice in the sterile laminar flow hood.

NOTE: Keeping *hMPCs* on ice during passing results in less cell aggregation.

3.2.9. Aspirate supernatant from *hMPCs* and gently resuspend pellet in 1 mL of GM.

3.2.10. Count cells using a hemocytometer or an automated cell counter.

NOTE: A 1:5 dilution of cell suspension to cell counting buffer is generally appropriate.

3.2.11. Seed *hMPCs* onto collagen-coated culture dishes containing pre-warmed GM at a density of 3,500 cells per cm². A combination of 10 cm plates and 24-well plates is ideal. The 24-well plates can be scanned daily on an imaging cytometer to monitor growth.

3.2.12. Follow the same procedure (starting at step 3.2.4.) for subsequent passages.

NOTE: Cell health and purity can be monitored across passages using a marker of cellular senescence (e.g., β -galactosidase) and immunostaining for Pax7.

3.3. Cryopreserve the *hMPCs*.

3.3.1. Cryopreserve excess cells for later use to make the culture volume more manageable.

3.3.2. For cryopreservation, follow the trypsinization procedure described in steps 3.2.4.-3.2.9.

3.3.3. While the cell suspension is being centrifuged, prepare a mixture of 20% (e.g., 400 μ L) dimethyl sulfoxide (DMSO) and 80% (e.g., 800 μ L) GM. Pipette up and down 20x to ensure adequate mixing. Leave the mixture to rest on ice.

3.3.4. Based on cell count, dilute the *hMPC* suspension to 2×10^6 /mL in GM.

3.3.5. Combine the *hMPC* suspension and the 20% sterile DMSO/80% GM mixture 50/50. The final cryopreservation media is a combination of DMSO (10%) and GM (90%) containing 1×10^6 *hMPCs*/mL.

3.3.6. Place 1 mL of the *hMPC* suspension prepared in step 3.3.5 into as many cryovials as the initial cell count allows. Then place aliquoted *hMPCs* in an isopropanol-controlled rate chiller in a -80 °C freezer overnight.

NOTE: At this step, typically between 5 and 15 million cells are obtained, 30–60% of which are identified as *hMPCs* by FACS.

4. Isolating Pax7⁺ human muscle progenitor cells via flow cytometry

4.1. To prepare *hMPCs* from live cultures (step 3.2.11), trypsinize enough culture vessels for 1×10^6 – 2×10^6 total cells. Prepare cell suspensions as described in steps 3.2.4–3.2.9.

4.2. Prepare *hMPCs* from cryopreserved cultures (step 2.4.6 or 3.3.6).

4.2.1. Select 1–2 tubes containing $\sim 1 \times 10^6$ passage four *hMPCs* for each donor.

NOTE: Using passage four *hMPCs* allows for adequate cell yield while still maintaining proliferative potential until passage six (**Figure 2**).

4.2.2. Rapidly thaw cells by removing them from the -80 °C freezer and placing them in a 37 °C water bath with the liquid in the tube completely submerged.

4.2.3. When thawing is almost complete (~5 min), transfer *hMPCs* to a sterile laminar flow hood and transfer the contents of the cell suspension into 6 mL of pre-warmed GM in a 15 mL conical tube at a rate of 1 mL/min.

4.2.4. Prepare cell suspensions as described in steps 3.2.7–3.2.9.

4.2.5. Resuspend each pellet in 3 mL of FACS buffer (500 mL PBS, 12.5 mL normal goat serum, 2 mL 0.5 M EDTA, pH 7.4, passed through a 0.22 µm filter).

4.2.6. Centrifuge the cell suspension at 300 x *g* for 5 min at RT.

4.2.7. Aspirate supernatant from *hMPCs* and gently resuspend pellet in 250 µL of FACS buffer; place on ice.

4.2.8. Place 100 µL of resuspended cells in a 1.7 mL microcentrifuge tube and leave on ice to be used as live/dead controls (step 4.5.1).

4.3. Stain *hMPCs* for cell surface markers.

4.3.1. Prepare the following antibody cocktail (amounts shown are for each culture, scale-up appropriately): 5 µL CD56 (neural cell adhesion molecule [NCAM]; PE-Cy7-conjugated) and 5 µL CD29 (β1-integrin; AlexFluor488-conjugated). The final antibody dilution is 1:50.

4.3.2. Add 10 µL of the antibody cocktail to each *hMPC* suspension and incubate on ice for 30 min in the dark.

4.3.3. Add 10 mL of FACS buffer to each suspension/antibody cocktail mixture.

4.3.4. Centrifuge each tube at 300 x *g* for 5 min at RT.

4.3.5. Resuspend the pellet in 300 µL of FACS buffer and transfer to a 5 mL round bottom polystyrene tube.

4.3.6. Keep the suspension on ice and protected from light.

4.4. Prepare the compensation beads.

NOTE: Positive controls for each antibody and a negative control (unstained beads) should be prepared using compensation beads. Compensation bead controls should be run on the flow cytometer prior to the samples to be sorted. These controls allow for determination of how much

the fluorochrome from each antibody can be detected in the other's channel so that compensation can be applied.

4.4.1. Vortex compensation beads.

4.4.2. In three 1.7 mL microcentrifuge tubes, add one drop (50 μ L) of compensation beads.

4.4.3. Add 1 μ L of each antibody to the appropriate tube, or none to the negative control.

4.4.4. Incubate in the dark for 15 min.

4.4.5. Add 1.5 mL of FACS buffer to each tube and vortex.

4.4.6. Centrifuge each tube at 300 x *g* in a benchtop centrifuge for 5 min at RT.

4.4.7. Aspirate supernatant and resuspend the pellet in 200 μ L of FACS buffer.

NOTE: The pellet may be very difficult to visualize; care should be taken when aspirating the supernatant.

4.4.8. Place the suspension in a 5 mL round bottom polystyrene tube and then keep the tube on ice and protected from light.

4.5. Prepare live and dead controls.

4.5.1. Divide the 100 μ L of *h*MPCs from step 4.2.8 into two 1.7 mL microcentrifuge tubes (one for the live control and one for the dead control).

4.5.2. Place the dead control tube in a heating block pre-warmed to 70 °C for >10 min. During this time, keep the live control on ice.

4.5.3. Centrifuge each tube at 300 x *g* in a benchtop centrifuge for 5 min at RT.

4.5.4. Aspirate supernatant and resuspend the pellet in 150 μ L of FACS buffer.

4.5.5. Combine the live/dead controls into a single 1.7 mL microcentrifuge tube.

4.6. Perform viability staining.

4.6.1. Add 1 μ L of 7-aminoactinomycin D (7-AAD) to all tubes containing cells to be sorted and the live/dead control.

4.6.2. Add GM to collagen-coated 24-well (250 μ L/well) and 10 cm (8 mL/well) plates and allow to equilibrate in a cell culture incubator at 37 °C in 5% CO₂.

4.7. Sort *hMPCs* via flow cytometry.

4.7.1. Sort *hMPCs* on a flow cytometer with 100 μm nozzle at a pressure of 20 psi. Determine live *hMPCs* based on side and forward scatter as well as 7-AAD negativity (sorting parameters and a representative sort can be seen in **Figure 3**).

4.7.2. Cells positive for PE-Cy7 (780/60 nm) and AlexFluor488 (530/30 nm) represent $\text{CD56}^+/\text{CD29}^+$ cells and thus are Pax7 positive cells (*hMPCs*). Sort double positive cells into collection tubes containing a 300 μL cushion of FACS buffer.

NOTE: The purity of sorted cultures can be determined by immunostaining for Pax7 followed by downstream flow cytometry analysis (**Figure 4**).

4.7.3. Discard all other cells.

4.7.4. Transfer the sorted cells to a 15 mL conical tube.

4.7.5. Spin the *hMPC* suspension at $300 \times g$ for 5 min at RT.

4.7.6. Carefully aspirate the supernatant and resuspend in 1 mL of GM.

4.7.7. Count cells using a hemocytometer or an automated cell counter.

NOTE: While the number of cells may be inferred as the number of events that are positively sorted, this count has been found to overestimate the number of cells by 10–40% compared to traditional cell counting methods.

4.7.8. Seed the *hMPCs* as outlined in step 3.2.11.

NOTE: Verification of the purity of *hMPC* population may be determined by immunostaining for Pax7 and analysis via flow cytometry.

5. Characterizing *hMPC* cultures using an imaging cytometer

5.1. Use the imaging cytometer (**Table of Materials**) to perform confluence scans.

5.1.1. Select **Create New Scan** and then enter/select the **Plate Category**, **Plate Profile** and **Plate ID**. The imaging cytometer can accommodate 96, 24, and 6 well plates.

5.1.2. Under **Application**, select **Confluence > Confluence 1**.

5.1.3. Select a well to view using the navigator on the right of the screen.

5.1.4. Find a plane of focus where the cells appear white compared to the background and select **register manual** or allow the imaging cytometer to select the plane of focus by selecting **register auto**.

NOTE: For the 96 well plates suggested in the **Table of Materials**, an approximate position is 4.3 units.

5.1.5. Use the mouse to highlight all wells that need to be scanned. The entire well area will automatically be scanned. Select **Start Scan**.

5.1.6. Select analysis settings optimal for plate and cell type.

NOTE: For *hMPCs*, the following settings are optimal: intensity = 6, saturated intensity = 0, diameter = 8, minimum thickness = 3, min cluster size = 500, and precision = high.

5.1.7. Select **Start Analysis**.

NOTE: Visually inspect images to ensure that the user and imaging cytometer agree regarding the perimeter of cells. If a large discrepancy exists, rescan the plate using a different plane of focus or different analysis settings. Confluence scanning can be used to compare growth between cultures, time the administration of a treatment, or to monitor cells and decide when to initiate differentiation. If differentiation is desired, details are provided in section 6.

5.2. Use the imaging cytometer to count cells.

5.2.1. Prepare a staining solution containing 12 mL of Ham's F12, 6 µL of Hoechst 33342 and 24 µL of propidium iodide.

5.2.2. Aspirate media and replace with staining solution, incubate at 37 °C for 15 min.

5.2.3. Aspirate staining solution and replace with Ham's F12.

NOTE: It is important to use a media with low levels of phenol red to obtain clear images. Alternatively, PBS may be used.

5.2.4. Load the plate onto the imaging cytometer and under **Application**, select **Cell Viability > Live + Dead + Total**. The live channel will be a bright field image, the dead channel will be the propidium iodide, and the total channel will be the Hoechst 33342.

5.2.5. Set the **Live** channel to bright field and click **Register Auto** under **Focus Setup**.

5.2.6. Under the **Total** channel, set the channel to blue. Use **Find focus** to bring the nuclei into focus or use the up and down arrows to manually set the focus. Click **Set offset** to select the focus.

5.2.7. Under the **Dead** channel, set the channel to red. Use **Find focus** to bring the dead cells into focus or use the up and down arrows to manually set the focus. Click **Set offset** to select the focus.

5.2.8. Select **Start Scan** to start the scan.

5.2.9. Select analysis parameters that are appropriate for the plate and cell type. Select **Start Analysis** to start the analysis. After the analysis is complete, visually verify that the analysis is appropriately counting cells (e.g., for the total channel check that every nucleus is being recognized as a nucleus and that the imaging cytometer is not counting any spots in the background as nuclei). If the scan and analysis are satisfactory, return the plate to the incubator, otherwise rescan or modify analysis parameters.

6. Differentiating *hMPCs* to myotubes

6.1. Determine the *hMPC* confluence using the protocol in section 5.1. To differentiate *hMPCs* into myotubes, it is ideal for cells to be 80–85% confluent as determined by the confluence function on the imaging cytometer (see section 5). If *hMPCs* from more than one unique donor are used, confluence scanning allows cells from different donors to begin differentiation at the time they reach confluence.

NOTE: At 80% confluence, there are approximately 66,000 cells/cm².

6.2. Wash cells 2x with differentiation media (97% Ham's F12, 2% equine serum, 1% penicillin/streptomycin, pH 7.4, passed through a 0.22 µm filter).

6.3. Maintain cells in differentiation media for as long as the experiment dictates while switching media daily.

6.4. Directly assess myotube formation by immunostaining for embryonic myosin heavy chain, which is typically detectable after 7–10 days of differentiation.

NOTE: The time required to form myotubes may differ slightly depending on donor and number of cells when differentiation was initiated; therefore, it is recommended that cultures from each donor be monitored individually.

REPRESENTATIVE RESULTS:

Representative flow cytometry results of *hMPC* isolation from human muscle tissue can be viewed in **Figure 1**. *hMPCs* can be identified by first gating events based on side scatter and forward scatter to eliminate dead cells or debris, followed by selecting only cells which are negative for 7-AAD and therefore are viable. Selection of cells positive for both the cell surface markers CD56 and CD29 represents the *hMPC* population. A biopsy of 60 mg only provides approximately 75–250 *hMPCs*.

A representative confluence scan is shown in **Figure 2A**. The green outlines in **Figure 2B** were generated by the imaging cytometer and show how the imaging cytometer determined confluence based on the analysis settings selected. In both images shown, the confluence determined by the imaging cytometer agrees with the confluence visually determined by the user. However, these images also highlight that the imaging cytometer is not perfect. For example, the red arrow highlights an imperfection in the plate which is being counted as cells. If a large number of these imperfection are present on the plate, the resulting confluence will not accurately represent the confluence of the cells. **Figure 2C** shows a representation of all three channels used to count cells (brightfield, blue and red) and the merged image. **Figure 2D** shows heterogeneity between donors. Discovering and accurately characterizing this heterogeneity is a key application of the imaging cytometer. **Figure 2E** shows that confluence measurements and nuclei counts from unique donors are highly correlated.

Figure 3 provides a guideline for how to identify *hMPCs* based on the FACS procedure described in this protocol. Gating based on forward and side scatter allows for separation of viable cells from debris (**Figure 3A**). A comparison of forward scatter by height to forward scatter by area was used to distinguish events representing single cells (the boxed area in **Figure 3B**). Viable cells are marked by a lack of incorporation of the viability stain 7-AAD (**Figure 3C**). Finally, cells positive for both CD56 and CD29 are identified as *hMPCs* (**Figure 3D**). To validate the FACS procedure described in this protocol the selection of Pax7 positive cells can be determined by immunostaining cells with a Pax7-specific antibody and measuring expression on a flow cytometer. **Figure 4** compares the number of *hMPCs* as determine by Pax7 immunostaining (**Figure 4A**) vs. the FACS procedure (**Figure 4B**) detailed in this protocol from the same population of cells. **Figure 5** uses Pax7 immunostaining and analysis via flow cytometry to show the enrichment of Pax7 expressing cells in the total population after FACS (**Figure 5A** compared to **Figure 5B**) as well as the maintenance of the number of Pax7 expressing cells after passaging (**Figure 5B** compared to **Figure 5C**).

hMPCs can be differentiated to form myotubes by following section 6. To determine whether isolated *hMPCs* maintain myogenic capacity cells can be immunostained with an antibody specific for embryonic myosin heavy chain and visualized using fluorescent microscopy. Representative images of embryonic myosin heavy chain positive myotubes can be viewed in **Figure 6**.

FIGURE LEGENDS:

Figure 1: Representative flow cytometry images for *hMPCs* sorted directly out of muscle biopsy tissue. (A) Side scatter area (y-axis) vs. forward scatter area (x-axis) gating strategy used to identify all cells in a donor sample. (B) 7-AAD viability stain (y-axis) vs. forward scatter (x-axis) to identify live cells. (C) Negative selection markers (y-axis) vs. the negative selection marker CD34. (D) Negative selection marker CD34 (y-axis) vs. the positive selection marker CD29. (E) Positive selection marker CD56 (y-axis) vs. the positive selection marker CD29 (x-axis). (F) Yields from three different skeletal muscle biopsies. FSC, forward scatter; SSC, side scatter; *hMPCs*, human muscle progenitor cells.

Figure 2: Proliferative potential of *hMPCs* derived from human donors is maintained after 6 passages. (A) Representative confluence scan. (B) Representative confluence scan with green outlines showing how the imaging cytometer determined confluence. The red arrow shows an imperfection on the plate. (C) Representative brightfield, Hoechst 33342 (blue), and propidium iodide (red) staining and a merged image of these three channels. Scale bars indicate 1mm. (D) Nuclei counts from 6 unique donors highlights *hMPC* heterogeneity. (E) Confluence and nuclei count are highly correlated.

Figure 3: Representative flow cytometry sorting parameters for *hMPC* isolation. (A) Side scatter area (y-axis) vs. forward scatter area (x-axis) gating strategy used to identify all cells in a donor sample. (B) Forward scatter height (y-axis) vs. forward scatter area (x-axis) to disqualify doublets from the sorting population. (C) Forward scatter height (y-axis) vs. 7-AAD incorporation to identify viable cells. (D) PE-Cy7 (CD56) vs. AF488 (CD29) staining with Q2 representing double positive cells (*hMPCs*). Color represents the density of events.

Figure 4: Comparison of CD29/CD56 positivity vs. Pax7 positivity in passage 4 *hMPCs* from the same donor. (A) Count (y-axis) vs. Pax7 expression (x-axis). (B) Positive selection marker CD56 (y-axis) vs. positive selection marker CD29 (x-axis). NCAM, neural cell adhesion molecule; *hMPC*, human muscle progenitor cell.

Figure 5: Pax7 positivity is maintained in *hMPCs* derived from the same donor over multiple passages. (A) Pax7 expression (x-axis) of cells which had been passaged 4 times prior to FACS sorting. (B) Pax7 expression (x-axis) of *hMPCs* 1 passage after FACS sorting for CD29 and CD56 positivity (5 total passages). (C) Pax7 expression (x-axis) of *hMPCs* 2 passages after FACS sorting for CD29 and CD56 positivity (6 total passages).

Figure 6: Staining of *hMPC* derived myotubes for embryonic myosin heavy chain. Representative microscopic images of differentiated *hMPCs* co-stained with DNA stain (Hoechst 33342, blue) and embryonic myosin heavy chain (green) (n = 3).

DISCUSSION:

Primary *hMPCs* are an important research model used to understand skeletal muscle biology and the regenerative process. Additionally, *hMPCs* have the potential to be used for therapy. However, there are recognized challenges in using primary *hMPCs* for both research and therapy, including the limited understanding of cells derived from humans¹⁰. There is also a large degree of variation in the expansion capacity among donor cultures, which limits the potential for use of *hMPCs* and can affect research results¹¹. For example, it has recently been demonstrated that *hMPCs* isolated using the method described in this protocol produced cells that varied in expansion capacity and this was driven by transcriptional profile that did not align with sex or age¹¹.

Here, we present an efficient and high yield method for isolating *hMPCs* in sufficient quantities for a number of downstream applications including differentiation into myotubes, thereby

modeling the myogenic process in vitro. Our method relies on FACS to isolate CD56⁺/CD29⁺ cells, which have previously been identified as representing an enriched Pax7 expressing population¹². Other equally valid cell sorting strategies have also been described^{5,6}.

The most important aspect of our *hMPC* isolation and culture method is the careful monitoring of individual cultures. The donor-based heterogeneity includes number of *hMPCs* obtained from each biopsy, their time to initiate proliferation in vitro, and their population expansion rate. Therefore, if differences in myotube formation after a treatment is the question of interest individual, *hMPC* cultures should be switched to treatment based on a confluence level determined *a priori* and assessed objectively using the imaging cytometer.

Our method to isolate and culture *hMPCs* has been optimized to obtain yields of cells in the millions for in vitro experiments from relatively little starting biopsy tissue (50–100 mg). The major benefit of this approach is that multiple assays and experiments can be performed on *hMPCs* from the same donor, allowing for maintenance of donor phenotype across experiments while introducing little inter-experiment variability. Our method differs from recently published methods that focus on extracting the purest population of Pax7 expressing cells directly out of biopsy tissue where the focus is xenotransplantation⁸. The challenge with these previous methods however, is that they require a starting tissue weight of greater than one gram. Therefore, they are not practical for research involving human skeletal muscle biopsies. One limitation of our method is that the potential for xenotransplantation has not been assessed. However, this type of procedure was not one of the original intended downstream applications. Future directions of this method may include assessing the engraftment capacity of *hMPCs*, after undergoing our isolation procedure, to determine whether there is potential for use of this procedure in MPC transplantation studies.

ACKNOWLEDGMENTS:

The authors thank the Cornell University, Biotechnology Resource Center Imaging Facility for their help with fluorescence activated cell sorting. We also thank Molly Gheller for her help with participant recruitment and Erica Bender for conducting the skeletal muscle biopsies. Finally, we thank the participants for their time and participation in the study. This work was supported by the National Institute on Aging of the National Institutes of Health under Award Number R01AG058630 (to B.D.C. and A.E.T.), by a Glenn Foundation for Medical Research and American Federation for Aging Research Grant for Junior Faculty (to B.D.C.), and by the President's Council for Cornell Women (to A.E.T.).


DISCLOSURES:

The authors have nothing to disclose.

REFERENCES:

1. Janssen, I., Heymsfield, S. B., Wang, Z., Ross, R. Skeletal muscle mass and distribution in 468 men and women aged 18–88 yr. *Journal of Applied Physiology*. **89** (1), 81–88 (2000).
2. D'Souza, D. M. et al. Decreased satellite cell number and function in humans and mice with type 1 diabetes is the result of altered notch signaling. *Diabetes*. **65** (10), 3053-3061 (2016).

3. Scheele, C. et al. Satellite cells derived from Obese humans with type 2 diabetes and differentiated into Myocytes in vitro exhibit abnormal response to IL-6. *PLoS One*. **7** (6), e39657 (2012).
4. Riddle, E. S., Bender, E. L., Thalacker-Mercer, A. E. Expansion capacity of human muscle progenitor cells differs by age, sex, and metabolic fuel preference. *American Journal of Physiology: Cell Physiology*. **315** (5), C643-C652 (2018).
5. Charville, G. W. et al. Ex vivo expansion and in vivo self-renewal of human muscle stem cells. *Stem Cell Reports*. **5** (4), 621-632 (2015).
6. Alexander, M. S. et al. CD82 Is a Marker for Prospective Isolation of Human Muscle Satellite Cells and Is Linked to Muscular Dystrophies. *Cell Stem Cell*. **19** (6), 800–807 (2016).
7. Thalacker-Mercer, A. et al. Cluster analysis reveals differential transcript profiles associated with resistance training-induced human skeletal muscle hypertrophy. *Physiological Genomics*. **45** (12), 499–507 (2013).
8. Garcia, S. M. et al. High-Yield Purification, Preservation, and Serial Transplantation of Human Satellite Cells. *Stem Cell Reports*. **10** (3), 1160–1174 (2018).
9. Yokoyama, W. M., Thompson, M. L., Ehrhardt, R. O. Cryopreservation and thawing of cells. *Current Protocols in Immunology*. Appendix 3, 3G (2012).
10. Bareja, A., Billin, A. N. Satellite cell therapy - from mice to men. *Skeletal Muscle*. **3** (1), 2 (2013).
11. Riddle, E. S., Bender, E. L., Thalacker-Mercer, A. Transcript profile distinguishes variability in human myogenic progenitor cell expansion capacity. *Physiological Genomics*. **50** (10), 817–827 (2018).
12. Xu, X. et al. Human Satellite Cell Transplantation and Regeneration from Diverse Skeletal Muscles. *Stem Cell Reports*. **5** (3), 419–434 (2015).






Figure 3

[Click here to access/download;Figure;Figure 3.tiff](#)




Figure 4

[Click here to access/download;Figure;Figure 4.tiff](#)










Figure 5

[Click here to access/download;Figure;Figure 5.tiff](#)






DONOR 1



DONOR 2



DONOR 3

Name of Material/ Equipment	Company	Catalog Number
0.25% Trypsin, 2.21 mM EDTA	Corning	25-053-CI
10 cm cell culture plate	VWR	664160
15 mL Falcon tube	Falcon	352196
24 well cell culture plate	Grenier Bio-One	662 160
7-AAD Viability Staining Solution	eBioscience	00-6993-50
Alexa Fluor 488 anti-human CD29, Clone: TS2/16	BioLegend	303016
Black 96-well cell culture plate	Grenier Bio-One	655079
Celigo S	Nexcelcom Bioscience	
Cell Strainer	VWR	352350
Collagen Type I (Rat Tail)	Corning	354236
Collagenase D	Roche	11 088 882 001
Dimethyl Sulfoxide	VWR	WN182
Dispase II	Sigma Life Sciences	D4693
Dulbecco's Modified Eagle Medium Low Glucose powder	Gibco	31600-034
Dulbecco's Phosphate Buffered Saline	Gibco	21600-010
EDTA Disodium Salt Dihydrate	J.T. Baker	4040-01
Fetal Bovine Serum	VWR	89510-186
Ham's F12	Gibco	21700-026
Heat Inactivated Equine Serum	Gibco	26-050-070
Hemocytometer	iNCyto	DHC-N0105
Hibernate A	Gibco	A1247501
Hoechst 33342, trihydrochloride, trihydrate	Life Technologies	H21492
Isopropanol	Fisher Scientific	A416P-4
Moxi buffer	Orflo	MXA006
Moxi Cassettes	Orflo	MXC002
Moxi z Mini Automated Cell Counter	Orflo	
Mr. Frosty Freezing Container	Thermo Fisher Scientific	5100-0001
Normal Goat Serum (10%)	Thermo Fisher Scientific	50062Z
PE-Cy7 Mouse Anti-human CD56 , Clone: B159	BD Pharmingen	557747
Penicillin/Streptomycin 100X Solution	Corning	30-002-CI
Propidium iodide	Thermo Fisher Scientific	P3566

Recombinant Human basic fibroblast growth factor	Promega	G5071
Recovery Cell Culture Freezing Medium	Gibco	12648-010
Sodium Bicarbonate	Fisher Scientific	S233-3
Sterile Round Bottom 5 mL tubes	VWR	60818-565
UltraComp eBeads	eBioscience	01-2222-42

Comments/Description
Trypsin used for removing adherent <i>h</i> MPCs from cell culture vessels
Plates used for culturing <i>h</i> MPCs
15 mL conical tubes used throughout the <i>h</i> MPC isolation and culturing protocols
Plates used for culturing <i>h</i> MPCs
Viability stain for identifying living cells during FACS sorting
Conjugated antibody for FACS
96-well cell culture plate ideal for fluorescent imaging using the Celigo S
Imaging cytometer used to track <i>h</i> MPC cultures
Cell strainer to eliminate large pieces of debris during muscle biopsy processing
Collagen for coating cell culture plates
Used for degradation of collagen and other connective tissue in the skeletal muscle biopsy tissue
Used for cryopreservation of <i>h</i> MPCs
A protease used for enzymatic digestion of skeletal muscle biopsy tissue
Low glucose DMEM for muscle biopsy processing
PBS for muscle biopsy processing
Required for FACS buffer
Fetal bovine serum used for <i>h</i> MPC growth media
Base media for <i>h</i> MPCs
Horse serum used to make <i>h</i> MPC differentiation media
Used to count cells
Media for preserving skeletal muscle biopsy tissue
DNA stain for identifying all cells using the Celigo S
Used for controlled rate freezing of <i>h</i> MPCs
Buffer for automated cell counter
Cassettes for automated cell counter
Automated cell counter
Commerically available controlled rate cell freezing container
Goat serum used in FACS buffer
Conjugated antibody for FACS
Antibiotics added to culture media
DNA stain for identifying dead cells using the Celigo S

Supplement in <i>h</i> MPC growth media to prevent spontaneous differentiation
Media used to cryopreserve muscle biopsy slurries
Added to Ham's F12
Tubes used for FACS
Compensation beads for calibrating flow FACS settings



1 Alewife Center #208
Cambridge, MA 02140
Tel: 617.945.9051
www.jove.com

ARTICLE AND VIDEO LICENSE AGREEMENT

Title of Article: Isolate, culture, characterize, and differentiate human muscle progenitor cells from the skeletal muscle biopsy procedure
Author(s): Brandon T. O'Brien, Jamie Blum, Sharon Seneit-Baumgarten, Erica Bender, Benjamin D. Cosgrove, Anna Thalacker-Merk

Item 1: The Author elects to have the Materials be made available (as described at <http://www.jove.com/publish>) via:

☒ Standard Access

☐ Open Access

Item 2: Please select one of the following items:

☒ The Author is **NOT** a United States government employee.

☐ The Author is a United States government employee and the Materials were prepared in the course of his or her duties as a United States government employee.

☐ The Author is a United States government employee but the Materials were NOT prepared in the course of his or her duties as a United States government employee.

ARTICLE AND VIDEO LICENSE AGREEMENT

1. **Defined Terms.** As used in this Article and Video License Agreement, the following terms shall have the following meanings: **"Agreement"** means this Article and Video License Agreement; **"Article"** means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein; **"Author"** means the author who is a signatory to this Agreement; **"Collective Work"** means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; **"CRC License"** means the Creative Commons Attribution-Non Commercial-No Derivs 3.0 Unported Agreement, the terms and conditions of which can be found at: <http://creativecommons.org/licenses/by-nc-nd/3.0/legalcode>; **"Derivative Work"** means a work based upon the Materials or upon the Materials and other pre-existing works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; **"Institution"** means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; **"JoVE"** means MyJoVE Corporation, a Massachusetts corporation and the publisher of The Journal of Visualized Experiments; **"Materials"** means the Article and / or the Video; **"Parties"** means the Author and JoVE; **"Video"** means any video(s) made by the Author, alone or in conjunction with any other parties, or by JoVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion

of the Article, and in which the Author may or may not appear.

2. **Background.** The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JoVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.

3. **Grant of Rights in Article.** In consideration of JoVE agreeing to publish the Article, the Author hereby grants to JoVE, subject to **Sections 4 and 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the "Open Access" box has been checked in **Item 1** above, JoVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.

ARTICLE AND VIDEO LICENSE AGREEMENT

discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JoVE reserves the right, in its absolute and sole discretion and without giving any reason therefore, to accept or decline any work submitted to JoVE. JoVE and its employees, agents and independent contractors shall have full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JoVE has sole discretion as to the method of making and publishing the Materials, including, without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

12. **Indemnification.** The Author agrees to indemnify JoVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JoVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JoVE, making of videos by JoVE, or publication in JoVE or elsewhere by JoVE. The Author shall be responsible for, and shall hold JoVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to

the making of a video by JoVE its employees, agents or independent contractors. All sterilization, cleanliness or decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's expense. All indemnifications provided herein shall include JoVE's attorney's fees and costs related to said losses or damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JoVE, its employees, agents or independent contractors.

13. **Fees.** To cover the cost incurred for publication, JoVE must receive payment before production and publication the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JoVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.

14. **Transfer, Governing Law.** This Agreement may be assigned by JoVE and shall inure to the benefits of any of JoVE's successors and assignees. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to me one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement is required per submission.

CORRESPONDING AUTHOR

Name:

Anna Thalacker-Merker

Department:

Division of Nutritional Sciences


Institution:

Cornell University

Title:

Assistant Professor

Signature:



Date:

12.18.2018

Please submit a **signed** and **dated** copy of this license by one of the following three methods:

1. Upload an electronic version on the JoVE submission site
2. Fax the document to +1.866.381.2236
3. Mail the document to JoVE / Attn: JoVE Editorial / 1 Alewife Center #200 / Cambridge, MA 02140

Dear Dr. Cao,

We appreciate the opportunity to revise our manuscript (JoVE submission JoVE59580R1 Isolation, culture, characterization, and differentiation of human muscle progenitor cells from the skeletal muscle biopsy procedure). We have addressed all of your comments and believe the quality of the manuscript is improved.

1. Please note that the editor has formatted the manuscript to match the journal's style. Please retain the same. The updated manuscript is attached and please use this version to incorporate the changes that are requested. Please email me your revised manuscript and figures at your earliest convenience.

We have maintained the formatting of the manuscript and uploaded it via the online portal.

2. Please address specific comments marked in the attached manuscript. Please turn on Track Changes to keep track of the changes you make to the manuscript.

We have addressed all comments using Track Changes.

3. Figure 2A and 2B: Please change the scale bar unit “mM” to “mm” and include a space between the number and its unit (i.e., 1 mm).

We have made the suggested change.

4. Figure 6: Please change the scale bar unit “ μm ”.

We have made the suggested changed.

5. Please ensure that all uploaded figures are at least 300 dpi.

We have remade figures that were not 300 dpi.