




↓ Dilution

10^{-1}

10^{-2}

10^{-3}

10^{-4}

10^{-5}

Focus Forming Units/mL =

(mean foci/well) x (dilution factor)

mL inoculum

$$\frac{(17+16) \times 10^5}{0.1}$$

0.1

Plate: D:\Plates\James\1

Double-click well to: ☒ Recount / QC individual well

☐ Set to zero ☐ Set to 'Rejected' ☐ Accept count

☐ Set to TNTC ☐ Set to 'Not Done'

Can be applied to individual wells, rows or columns.

