

Journal of Visualized Experiments

Quantifying subcellular ubiquitin-proteasome activity in the rodent brain

--Manuscript Draft--

Article Type:	Invited Methods Article - JoVE Produced Video
Manuscript Number:	JoVE59695R2
Full Title:	Quantifying subcellular ubiquitin-proteasome activity in the rodent brain
Keywords:	Ubiquitin; Proteasome; Nucleus; synapse; Cytoplasm; memory; Hippocampus; Amygdala
Corresponding Author:	Timothy Jarome, Ph.D. Virginia Polytechnic Institute and State University Blacksburg, VA UNITED STATES
Corresponding Author's Institution:	Virginia Polytechnic Institute and State University
Corresponding Author E-Mail:	tjarome@vt.edu
Order of Authors:	Taylor McFadden Rishi K. Devulapalli Timothy Jarome, Ph.D.
Additional Information:	
Question	Response
Please indicate whether this article will be Standard Access or Open Access.	Standard Access (US\$2,400)
Please indicate the city, state/province, and country where this article will be filmed . Please do not use abbreviations.	Blacksburg, VA, USA

TITLE:

Quantifying Subcellular Ubiquitin-Proteasome Activity in the Rodent Brain

AUTHORS AND AFFILIATION:

Taylor McFadden¹, Rishi K. Devulapalli² & Timothy J. Jarome^{1,2}

¹Department of Animal and Poultry Sciences, Virginia Polytechnic Institute and State University, Blacksburg, VA, USA

²School of Neuroscience, Virginia Polytechnic Institute and State University, Blacksburg, VA, USA

Corresponding author:

Timothy J. Jarome

Email Address: tjarome@vt.edu

KEYWORDS:

ubiquitin, proteasome, nucleus, synapse, cytoplasm, memory, hippocampus, amygdala

SHORT ABSTRACT:

This protocol is designed to efficiently quantify ubiquitin-proteasome system (UPS) activity in different cellular compartments of the rodent brain. Users are able to examine UPS functioning in nuclear, cytoplasmic and synaptic fractions in the same animal, reducing the amount of time and number of animals needed to perform these complex analyses.

LONG ABSTRACT

The ubiquitin-proteasome system is a key regulator of protein degradation and a variety of other cellular processes in eukaryotes. In the brain, increases in ubiquitin-proteasome activity are critical for synaptic plasticity and memory formation and aberrant changes in this system are associated with a variety of neurological, neurodegenerative and psychiatric disorders. One of the issues in studying ubiquitin-proteasome functioning in the brain is that it is present in all cellular compartments, in which the protein targets, functional role and mechanisms of regulation can vary widely. As a result, the ability to directly compare brain ubiquitin protein targeting and proteasome catalytic activity in different subcellular compartments within the same animal is critical for fully understanding how the UPS contributes to synaptic plasticity, memory and disease. The method described here allows collection of nuclear, cytoplasmic and crude synaptic fractions from the same rodent (rat) brain, followed by simultaneous quantification of proteasome catalytic activity (indirectly, providing activity of the proteasome core only) and linkage-specific ubiquitin protein tagging. Thus, the method can be used to directly compare subcellular changes in ubiquitin-proteasome activity in different brain regions in the same animal during synaptic plasticity, memory formation and different disease states. This method can also be used to assess the subcellular distribution and function of other proteins within the same animal.

INTRODUCTION

The ubiquitin-proteasome system (UPS) is a complex network of interconnected protein structures and ligases that controls the degradation of most short-lived proteins in cells¹. In this

system, proteins are marked for degradation or other cellular processes/fates by the small modifier ubiquitin. A target protein can acquire 1-7 ubiquitin modifications, which can link together at one of seven lysine (K) sites (K6, K11, K27, K29, K33, K48 and K63) or the N-terminal methionine (M1; as known as linear) in the previous ubiquitin². Some of these polyubiquitin tags are degradation-specific (K48)³, while others are largely independent of the protein degradation process (M1)⁴⁻⁶. Thus, the protein ubiquitination process is incredibly complex and the ability to quantify changes in a specific polyubiquitin tag is critical for ultimately understanding the role of that given modification in cellular functioning. Further complicating the study of this system, the proteasome, which is the catalytic structure of the UPS⁷, both degrades proteins but can also be involved in other non-proteolytic processes^{8,9}. Not surprisingly then, since its initial discovery, normal and aberrant ubiquitin-proteasome activity has been implicated in long-term memory formation and a variety of disease states, including many neurological, neurodegenerative and psychiatric disorders^{10,11}. As a result, methods which can effectively and efficiently quantify UPS activity in the brain are critical for ultimately understanding how this system is dysregulated in disease states and the eventual development of treatment options targeting ubiquitin and/or proteasome functioning.

There are a number of issues in quantifying ubiquitin-proteasome activity in brain tissue from rats and mice, which are the most common model systems used to study UPS function, including 1) the diversity of ubiquitin modifications, and 2) distribution and differential regulation of UPS functioning across subcellular compartments¹²⁻¹⁴. For example, many of the early demonstrations of ubiquitin-proteasome function in the brain during memory formation used whole cell lysates and indicated time-dependent increases in both protein ubiquitination and proteasome activity¹⁵⁻²⁰. However, we recently found that ubiquitin-proteasome activity varied widely across subcellular compartments in response to learning, with simultaneous increases in some regions and decreases in others, a pattern that differs significantly from what was previously reported in whole cell lysates²¹. This is consistent with the limitation of a whole cell approach, as it cannot dissociate the contribution of changes in UPS activity across different subcellular compartments. Though more recent studies have employed synaptic fraction protocols to study the UPS specifically at synapses in response to learning²²⁻²⁴, the methods used occlude the ability to measure nuclear and cytoplasmic ubiquitin-proteasome changes in the same animal. This results in an unnecessary need to repeat experiments multiple times, collecting a different subcellular fraction in each. Not only does this result in a greater loss of animal lives, but it eliminates the ability to directly compare UPS activity across different subcellular compartments in response to a given event or during a specific disease state. Considering that protein targets of ubiquitin and the proteasome vary widely throughout the cell, understanding how ubiquitin-proteasome signaling differs in distinct subcellular compartments is critical for identifying the functional role of the UPS in the brain during memory formation and neurological, neurodegenerative and psychiatric disorders.

To address this need, we recently developed a procedure in which nuclear, cytoplasmic and synaptic fractions could be collected for a given brain region from the same animal²¹. Additionally, to account for the limited amount of protein that can be obtained from collecting multiple subcellular fractions from the same sample, we optimized previously established

protocols to assay in vitro proteasome activity and linkage-specific protein ubiquitination in lysed cells collected from rodent brain tissue. Using this protocol, we were able to collect and directly compare learning-dependent changes in proteasome activity, K48, K63, M1 and overall protein polyubiquitination levels in the nucleus and cytoplasm and at synapses in the lateral amygdala of rats. Here, we describe in detail our procedure (**Figure 1**), which could significantly improve our understanding of how the UPS is involved in long-term memory formation and various disease states. However, it should be noted that the in vitro proteasome activity discussed in our protocol, while widely used, does not directly measure the activity of complete 26S proteasome complexes. Rather, this assay measures the activity of the 20S core, meaning it can only serve as a proxy to understand the activity of the core itself as opposed to the entire 26S proteasome complex.

PROTOCOL

All procedures including animal subjects have been approved by the Virginia Polytechnic Institute and State University Institutional Animal Care and Use Committee (IACUC).

1. Collection and dissection of rodent brain tissue

NOTE: This protocol can be applied to a variety of brain regions and used with various tissue collection procedures. Below is the procedure used in our lab for subcellular of rat brain tissue, using 8-9 week old male Sprague Dawley rats. In order to process all cellular compartments in the same animal, section 1.3. must be followed regardless of the tissue collection procedure used.

1.1. Extract the rat brain and place into a cryogenic cup pre-chilled on dry ice. Flash-freeze the brain on dry ice, or use liquid nitrogen if available, and transfer to a -80 °C freezer. Brains can be dissected the same day or at a later time.

1.2. Remove the frozen brain from the cryogenic cup and place into a rat brain matrix chilled with dry ice. Each slot on the matrix corresponds to 0.5 mm, which can be used to determine the approximate location of the region of interest (ROI).

1.2.1. Using a Rat Brain Atlas, insert a razor blade into the matrix immediately in front (anterior) of the predicted start of ROI. Next, insert a razor blade immediately at the predicted end (posterior) of ROI.

1.2.2. Remove the slice between the razor blades using a scalpel and place onto a microscope slide chilled on dry ice. Typically, sections are 2-3 mm thick but this varies based on the ROI.

1.3. Using a scalpel, dissect out the ROI one hemisphere at a time. Place each hemisphere into separate 1.5 mL microcentrifuge tubes pre-chilled on dry ice. Frozen tissue can immediately be used for subcellular fractionation or processed at a later time if stored at -80 °C.

2. Nuclear and cytoplasmic extraction

NOTE: This protocol uses premade stock solutions of common lab chemicals, including 0.1 M HEPES, 1 M MgCl₂, 1 M Dithiothreitol (DTT), 0.5 M ethylenediaminetetraacetic acid (EDTA), 5 M NaCl, 10% NP-40 (IGEPAL) and 50% glycerol. If endpoint is used in proteasome activity assays, glycerol and ATP can be added to all buffers to help prevent disassembly of proteasome complexes during lysis.

2.1. Prepare the Lysis Buffer. Add 8.63 mL of ultrapure (deionized and distilled) water to a sterile 15 mL conical tube.

2.1.1. To the conical tube, add 1,000 µL of 0.1 M HEPES, 15 µL of 1 M MgCl₂, 100 µL of 1 M DTT and 50 µL of 10% NP-40.

2.1.2. Add 100 µL of protease inhibitor cocktail and 100 µL of phosphatase inhibitor cocktail. Briefly vortex the solution to mix and place on wet ice to chill. Note that the solution may have a yellow tint from the phosphatase inhibitor.

NOTE: The use of protease inhibitors is essential to preserve protein and ensure the specificity of the in vitro proteasome activity assay. However, these inhibitors may also result in a slight reduction in proteasome activity, meaning that the activity measured on the in vitro assay may be an underestimate of actual activity levels.

2.2. Prepare the Extraction Buffer. Add 5.925 mL of ultrapure (deionized and distilled) water to a sterile 15 mL conical tube.

2.2.1. To the conical tube, add 2,000 µL of 0.1 M HEPES, 1250 µL of 50% glycerol, 15 µL of 1 M MgCl₂, 5 µL of 1 M DTT, 5 µL of 0.5 M EDTA and 600 µL of 5 M NaCl.

2.2.2. Add 100 µL of protease inhibitor cocktail and 100 µL of phosphatase inhibitor cocktail. Briefly vortex the solution to mix and place on wet ice to chill. Note that the solution may have a yellow tint from the phosphatase inhibitor.

2.3. Remove the 1.5 mL centrifuge microtube containing one hemisphere of the ROI from the -80 °C freezer. Ensure that the hemisphere used is counterbalanced across extraction conditions for each experimental group. For example, in a two-group experiment, use the left hemisphere for the first two animals in each group and the right hemisphere for the next two samples, etc.

2.4. Using a scalpel, transfer the frozen brain tissue to a 2 mL glass Teflon homogenizer. Add 500 µL of Lysis buffer to the Teflon tube.

2.4.1. Using pestle B, homogenize the same tissue using 15 strokes until no visible amount of solid material are present. Use a turning motion (clockwise or counter-clockwise) during each stroke.

2.5. Using a 1,000 μ L pipette, transfer the homogenized sample to a new 1.5 mL microcentrifuge tube. Place the tube on wet ice and incubate for 30 min.

2.6. Place tube in microcentrifuge and spin for 10 min at 845 x *g* and 4 °C. After completion, carefully remove the supernatant by pipetting and place into a new 1.5 mL microcentrifuge tube; this is the Cytoplasmic Fraction, which can be stored on ice or at 4 °C until Section 2.9.

2.7. Add 50 μ L of Extraction buffer to the resulting pellet and resuspend by pipetting. Do not vortex the pellet. Place the tube containing the resuspended pellet on ice and incubate for 30 min.

2.8. Place the tube in microcentrifuge and spin for 20 min at 21,456 x *g*, 4 °C. After completion, carefully remove the supernatant by pipetting and place into a new 1.5 mL microcentrifuge tube; this is the Nuclear Fraction. The pellet can now be discarded.

2.9. Measure protein concentration of the Cytoplasmic and Nuclear extractions using the Dc protein assay (according to the manufacturer's instructions), or any equivalent assay for samples collected using non-ionic detergents. Proceed immediately to proteasome activity assay (Section 4) or Western blotting (Section 5). Alternatively, samples can be stored at -80 °C until needed.

3. Synaptic fraction collection

NOTE: This protocol uses premade stock solutions of common lab chemicals, including 1 M Tris (pH 7.5), 0.5 M EDTA, 5 M NaCl and 10% SDS. If endpoint is used in proteasome activity assays, glycerol and ATP can be added to all buffers to help prevent disassembly of proteasome complexes during lysis

3.1. Prepare the TEVP buffer with 320 mM sucrose. Add 60 mL of ultrapure (deionized and distilled) water to a clean 100 mL beaker.

3.1.1. To the beaker, add 1,300 μ L of 1 M Tris (pH 7.5), 260 μ L of 0.5 M EDTA, 100 μ L of protease inhibitor cocktail, 100 μ L of phosphatase inhibitor cocktail and 10.944 g of sucrose. Mix with a stir bar until sucrose is completely dissolved.

3.1.2. Bring pH to ~7.4 by adding hydrochloric acid (HCl) dropwise to the solution using a pipette.

3.1.3. Transfer the solution to a 100 mL volumetric flask. Bring volume to 100 mL by adding ultrapure water. Chill the final solution on ice until needed.

3.2. Prepare the Homogenization Buffer. Add 60 mL of ultrapure (deionized and distilled) water to a clean 100 mL beaker.

3.2.1. To the beaker, add 5,000 μL of 1 M Tris (pH 7.5), 3,000 μL of 5 M NaCl, 100 μL of protease inhibitor cocktail, 100 μL of phosphatase inhibitor cocktail and 2,000 μL of 10% SDS. Mix with a stir bar.

NOTE: Ionic detergents can interfere with proteasome activity assays by resulting in denaturing of the proteasome complex. The volume of SDS could be lowered to help preserve proteasome function if desired.

3.2.2. Bring pH to ~ 7.4 by adding HCl dropwise to the solution using a pipette.

3.2.3. Transfer the solution to a 100 mL volumetric flask. Bring the volume to 100 mL by adding ultrapure water. Store the final solution at room temperature; do not chill as SDS will precipitate out of solution.

3.3. Remove the 1.5 mL centrifuge containing one hemisphere of the ROI from the $-80\text{ }^{\circ}\text{C}$ freezer. Ensure that the hemisphere used is counterbalanced across extraction conditions for each experimental group. For example, in a two-group experiment, use the left hemisphere for the first two animals in each group and the right hemisphere for the next two samples, etc.

3.4. Using a scalpel, transfer the frozen brain tissue to a 2 mL glass Teflon homogenizer. Add 500 μL of TEVP buffer to the Teflon tube.

3.4.1. Using pestle B, homogenize the same tissue using 15 strokes until no visible transfers of solid material are present. Use a turning motion (clockwise or counter-clockwise) during each stroke.

3.5. Using a 1,000 μL pipette, transfer the homogenized sample to a new 1.5 mL microcentrifuge tube. Centrifuge the sample at $1,000 \times g$ for 10 min, $4\text{ }^{\circ}\text{C}$.

3.6. Collect the supernatant and transfer to a new 1.5 mL microcentrifuge tube using a 1,000 μL pipette. Centrifuge the sample at $10,000 \times g$ for 10 min, $4\text{ }^{\circ}\text{C}$. The original pellet (P1) contains nuclei and the large debris and can be discarded.

3.7. Transfer the supernatant to a new 1.5 mL microcentrifuge tube. This is a Cytosolic Fraction. Add 50 μL of Homogenization buffer to the pellet (P2) and resuspend by pipetting until no solid material is visible.

3.8. Centrifuge the sample at $20,000 \times g$ for 10 min, $4\text{ }^{\circ}\text{C}$. Transfer the supernatant to a new 1.5 mL microcentrifuge tube; this is the Crude Synaptosomal Membrane (Synaptic) Fraction. The pellet can be discarded.

3.9. Measure the protein concentration of the Synaptic fraction using the Dc protein assay (according to the manufacturer's instructions), or any equivalent assay for samples collected

using ionic detergents. Proceed immediately to proteasome activity assay (Section 4) or Western blotting (Section 5). Alternatively, samples can be stored at -80 °C until needed.

4. Proteasome activity assay

NOTE: Proteasome activity can be measured in homogenized brain tissue using a slightly modified version of the 20S Proteasome Activity Kit. This assay does not directly measure the activity of complete 26S proteasome complexes. Rather, it measures the activity of the 20S core, meaning it can only serve as a proxy to understand the activity of the core itself as opposed to the entire 26S proteasome complex. The success of this assay declines with repeated freeze-thaw cycles and/or increasing levels of detergents, particularly ionic, and requires the use of a plate reader with a 360/460 (excitation/emission) filter set and heating capabilities up to 37 °C.

4.1. Plate reader settings: Pre-warm to 37 °C and hold through the run.

4.1.1. Set excitation to 360 and emission to 460. If the 96 well plate used is clear, set the optics position to **Bottom**. If a dark/black 96 well plate is used, set optics position to **Top**.

4.1.2. Set older plate reader models to **Auto-Gain** under the fluorescent read options; newer models are preset for this. Program a kinetic run with a time of 2 h, scanning (reading) every 30 min.

4.2. Reconstitute the 10x assay buffer provided in the kit with 13.5 mL of ultrapure water. Add 14 µL of 100 mM ATP to the now 1x buffer; this significantly enhances proteasome activity in the samples and improves assay reliability. The final 20S Assay buffer can be stored on ice or at 4 °C until needed and is stable for several months.

4.3. Reconstitute the AMC standard provided in kit with 100 µL of DMSO. Perform this step in the dark or under low light conditions, as the standard is light sensitive.

4.4. Create a stand curve of AMC using the reconstituted standard, in the dark or under low light conditions.

4.4.1. In separate 0.5 mL microcentrifuge tubes, add 16, 8, 6.4, 3.2, 1.6, 0.8, 0.4 and 0 µL of the AMC standard, which corresponds to 20, 10, 8, 4, 2, 1, 0.5 and 0 µM AMC concentration.

4.4.2. To these tubes, in the same order, add 84, 92, 93.6, 96.8, 98.4, 99.2, 99.6 and 100 µL of 20S Assay buffer. This creates a series of high to low AMC concentrations, which will be used for plate reader calibration and analysis of proteasome activity in the homogenized samples.

4.4.3. Store all diluted standards on ice in the dark until needed.

4.5. Reconstitute the proteasome substrate (Suc-LLVY-AMC) provided in kit with 65 µL of DMSO. Perform this step in the dark or under low light conditions, as the substrate is light

sensitive. Create a 1:20 dilution of the proteasome substrate in a new 1.5 mL microcentrifuge tube using 20S Assay buffer. Store the diluted substrate on ice in the dark until needed.

NOTE: For example, if the plate will have 10 samples and 1 blank, you will need enough diluted substrate for 22 wells (with duplicates) at 10 μ L per well. This equates to 220 μ L need + 30 μ L for pipetting errors, requiring 12.5 μ L of substrate and 237.5 μ L of 20S Assay buffer.

4.6. Thaw desired samples (if frozen) and add a normalized amount to a 96 well plate. Run each sample in duplicates. The amount of sample needed varies based on tissue preparation. Generally, 10-20 μ g is sufficient for any subcellular fraction.

4.7. Bring the sample well volume to 80 μ L with ultrapure water. The amount added depends on the volume of sample added. For example, if sample 1 was 4.5 μ L of protein and sample 2 was 8.7 μ L, then the amount of water needed will be 75.5 μ L and 71.3 μ L, respectively. In two separate wells, add 80 μ L of water alone; these will be the Assay Blanks.

NOTE: In order to limit changes in protein volume due to pipetting errors, protein concentrations can be normalized to the least concentrated sample. This will allow use of the same sample volume in all conditions.

4.8. Add 10 μ L of 20S Assay buffer to each well, including Assay Blanks. A repeater/automated pipette is recommended here to ensure consistent assay volume across wells.

4.9. **Optional:** At this stage, introduce in vitro manipulations if desired; this will require that each sample has an additional 2 wells per treatment, including the vehicle. If so, add 5-10 μ L of drug/compound of interest to 2 of the sample wells and an equivalent volume of control/vehicle to another 2 wells. Place the plate onto the pre-warmed plate reader or into a 37 $^{\circ}$ C incubator for 30 min.

4.10. Turn off the lights or enter a dark room. Add all 100 μ L of diluted AMC standards to a new well; each standard will have a single well.

4.11. In the dark, add 10 μ L of diluted proteasome substrate to wells containing sample and Assay Blanks, but not to the AMC standard. A repeater/automated pipette is recommended here to ensure consistent assay volume across wells.

4.12. Place the plate into the plate reader and start the kinetic run.

NOTE: The plate does not need to be under constant agitation during the kinetic run, however, the user can choose to if desired

4.13. At end of the kinetic run, export raw 360/460 fluorescent values to Microsoft Excel.

4.13.1. Average together duplicate wells for each standard, each sample and the Assay Blank for all 5 scans. Raw fluorescent values should increase across scans for the samples but remain stable (or decrease slightly) for standards and Assay Blanks.

4.13.2. Take the highest AMC standard well average and divide by the known concentration (20 μ M). Divide this value by the sample concentration used in the assay to get standardized AMC value. For each sample and Assay Blank average, divide by the standardized AMC value.

4.13.3. Take this final value and divide by the sample concentration used in the assay to get the normalized value for each sample and the Assay Blank. Do this for all 5 scans.

4.13.4. Subtract the normalized Assay Blank value from each normalized sample value for all 5 scans.

5. Quantification of linkage-specific protein ubiquitination

5.1. Quantification of diverse polyubiquitin tags in different subcellular fractions collected from rodent brain tissue should be done using a variety of standard Western blotting protocols in combination with unique, linkage-specific polyubiquitin antibodies.

5.1.1. For denaturing, mix the normalized samples with an equal volume of Laemmli sample buffer supplemented with β -mercaptoethanol to a percentage of 5% by volume, per the manufacturer's instruction.

5.1.2. For quantification of all monoubiquitination and polyubiquitination modifications regardless of linkage, use a pan-ubiquitin.

5.1.3. To recognize all polyubiquitinated proteins, use an ubiquitin antibody that does not cross-react with monoubiquitination.

5.1.4. For linkage-specific polyubiquitination, use antibodies that can detect Lysine-27, Lysine-48, Lysine-63 and linear (M1) polyubiquitination.

5.2. To prevent cross-contamination between developments, which could result in a false positive or interfere with imaging of a different ubiquitin modification, strip membranes between developments using 0.1 M NaOH for 10 min.

5.2.1. Wash the membranes in TBS with 0.1% tween twice for 10 min and reblock (using whatever blocking agent is preferred in the Western blot protocol used).

5.2.2. Incubate the membrane with the secondary antibody and redevelop in order to confirm that the membrane was stripped of primary and secondary antibodies properly.

5.2.3. **Recommended:** For successful development of different ubiquitin antibodies without contamination, use fluorescent or near-infrared imaging systems. This often times can prevent carryover between antibodies.

5.3. Some ubiquitin Western blot images will provide columns of distinct bands (such as M1) while others produce a smear-like pattern with few or no clear lines (common with K48). For quantification of imaged ubiquitin Western blots, draw a box around the column that extends the entire molecular standards ladder.

5.3.1. Adjust the box up (or down) if ubiquitin staining does extend through the entire ladder; this is common for Lysine-48 modifications and varies widely across subcellular compartments.

5.3.2. Subtract out the background, which is calculated as the mean optical density of the background immediately surrounding the column on all sides.

REPRESENTATIVE RESULTS

Using the procedure described here, nuclear, cytoplasmic and synaptic fractions were collected from the lateral amygdala of the rat brain (**Figure 1**). Purity of the individual fractions were confirmed via Western blotting, probing with antibodies against proteins that should be enriched or depleted in the lysate. In the first hemisphere where a crude synaptic fraction was collected, postsynaptic density protein 95 (PSD95) was present in the synaptic, but not nuclear, fraction, with lower levels in the cytoplasm (**Figure 2A**). This is consistent with previous work demonstrating that the synaptic fraction preparation isolates both presynaptic and postsynaptic components²⁵. Conversely, the nuclear protein histone H3 was present in the nuclear, but not the synaptic, fraction, with lower levels in the cytoplasm (**Figure 2B**). The presence of PSD95 and H3 in the cytoplasm is consistent with their cytoplasmic translation. The cytoplasmic protein β -tubulin was present in our cytoplasmic fraction, but was largely absent from the nuclear lysate (**Figure 2C**), with lower levels in the synaptic region. This suggests that we were able to produce a nuclear fraction that was largely absent of cytoplasmic proteins. The presence of tubulin in the synaptic region is consistent with previous studies²⁶. All three fractions showed similar levels of the housing keeping protein β -actin (**Figure 2D**), which was used as a loading control. Collectively, these results confirm that the purity of the nuclear, cytoplasmic and synaptic fractions collected from a single rat lateral amygdala.

Next, all lysates were confirmed for functional proteasome activity using the described modified version of the in vitro 20S proteasome activity assay. In all lysates, success of the assay was defined as an increase in raw fluorescent units (RFU) detected from the first scan (0 min) to the fifth/final scan (120 min). For all of these analyses, 10 μ M AMC was used as the highest standard for normalization of the raw fluorescent units (RFU). In the crude synaptic fraction, RFU peaked at scan 5 (**Figure 3A**), resulting in a normalized RFU of just under 0.1 (**Figure 3B**). In the cytoplasmic fraction, RFU increased across scans (**Figure 3C**) with a final normalized RFU of \sim 1.6 (**Figure 3D**). The nuclear fraction also displayed time-dependent changes in RFU (**Figure 3E**), with a final normalized RFU of 0.3 (**Figure 3F**). The differences in proteasome activity across compartments likely reflects the availability of proteasomes in the fraction, which are generally

most abundant in the cytoplasm and nucleus²⁷, compartments which have the highest level of activity in our preparation. The lowest level of activity in the synaptic region is consistent with it being the only lysate that was collected using ionic detergents, which can reduce proteasome activity due to the harsher denaturing condition. Importantly, RFU did not increase across time in the Assay Blanks or in lysates (synaptic) incubated with the highly specific and potent proteasome inhibitor clasto-lactacystin- β -lactone (β lac; **Figure 3G**), displaying final normalized RFU levels of 0.01 and 0.001, respectively (**Figure 3H**). This suggests that the observed change in RFU was due specifically to activity of the proteasome and not other proteases. Furthermore, when analyzed across experimental conditions, there was an increase in nuclear, but not cytoplasmic, proteasome activity in the lateral amygdala following learning, which occurred simultaneously with a decrease in synaptic proteasome activity in comparison to control animals (**Figure 4**). Collectively, these results confirm that proteasome activity could be accurately measured in all three subcellular fractions collected from a single rat lateral amygdala.

One of the advantages of the proteasome activity assay is that in vitro manipulations can be introduced to the samples immediately prior to addition of the proteasome substrate, which allows identification of specific molecules which can regulate the proteasome in that particular subcellular fraction. As an example of this, the role of CaMKII (calcium/calmodulin dependent protein kinase II) was assessed in the synaptic fraction, the harshest denatured lysate collected, since CaMKII is thought to regulate proteasome function at synapses. A 30 min incubation with the CaMKII inhibitor myr-AIP (myristolayted autocamide-2 related inhibitory peptide) resulted in a significantly diminished increase in proteasome activity on the assay, which only reached levels that were ~61% of the untreated controls (**Figure 5A**). Conversely, the same manipulation applied to the cytoplasmic did not result in a change in proteasome activity from vehicle treated wells (**Figure 5B**). These results confirm that proteasome activity can be further manipulated in vitro and that this manipulation can have different effects depending on the subcellular fraction collected.

In addition to quantifying proteasome activity, the described protocol can be used to measure subcellular differences in diverse ubiquitin modifications using Western blot procedures. It is important to note that the ubiquitin tags that can be detected are limited by the availability of linkage-specific antibodies, which currently include K48, K63 and M1 for rats (note: a K27 antibody is available, though did not produce a detectable image in any lateral amygdala fraction or whole cell lysates under a variety of conditions). Overall polyubiquitination, degradation-independent linear/M1 and K63 ubiquitination and degradation-specific K48 ubiquitination were detected in all subcellular fractions. Importantly, when analyzed across different experimental conditions, there was an increase in overall (**Figure 6A**), linear (**Figure 6B**), K63 (**Figure 6C**) and K48 (**Figure 6D**) polyubiquitination in the lateral amygdala nuclear fraction following learning in comparison to control animals. At the same time, in the cytoplasmic region, overall polyubiquitination decreased and K48 ubiquitination increased following learning, while synaptic K63 ubiquitination was reduced. Collectively, these results indicate that within the same animal subcellular differences in linkage-specific protein ubiquitination can be accurately detected.

FIGURE LEGENDS

Figure 1. Schematic for subcellular fractionation of rat brain tissue. Rodent brain is extracted, brain region dissected and hemispheres split. Using a series of buffers and centrifuging steps, nuclear and cytoplasmic fractions are collected from one hemisphere while a crude synaptic fraction is collected from the other. Both fractions are then used for proteasome activity assays and Western blotting, to assay protein polyubiquitination levels.

Figure 2. Confirmation of crude synaptic, nuclear and cytoplasmic fraction purity. (A) The synaptic protein postsynaptic density protein 95 (PSD95; 1:1000) was present in the synaptic, but not nuclear fraction with lower levels in the cytoplasm. (B) Histone H3 (1:500) was present in the nuclear, but not synaptic, fraction with lower expression in the cytoplasm. (C) β -Tubulin (1:1000) was present in the cytoplasmic, but largely absent from the nuclear, fraction with lower expression in the synaptic lysate. (D) Housekeeping protein β -actin (1:1000) was present in all subcellular compartments. Areas indicate the expected size of the target protein. This figure has been modified from Orsi, S. A. et al.²¹.

Figure 3. Quantification of proteasome activity in nuclear, cytoplasmic and synaptic fractions collected from the lateral amygdala of the same animal. During the in vitro proteasome activity assay, relative fluorescent units (RFU) detected increased from the beginning (Scan 1) to the end (Scan 5) of the assay in the synaptic (A), Cytoplasmic (C) and nuclear (E) fractions. Quantification of this change from baseline indicated a normalized RFU value (relative to 10 μ M AMC) of 0.1 in the synaptic (B), 1.6 in the cytoplasmic (D) and 0.3 in the nuclear (F) fractions. The proteasome inhibitor β lac prevented RFUs from changing across the session (G-H). This figure has been modified from Orsi, S. A. et al.²¹.

Figure 4. Subcellular differences in proteasome activity in the lateral amygdala of the same animal. An increase in nuclear proteasome activity was detected in trained (fear conditioned) animals relative to controls, which corresponded to a decrease in activity within the synaptic fraction. Cytoplasmic proteasome activity remained at baseline. * $P < 0.05$ from Control. This figure has been modified from Orsi, S. A. et al.²¹.

Figure 5. In vitro manipulation of proteasome activity in collected synaptic and cytoplasmic fractions. In vitro manipulation of CaMKII signaling via the inhibitor AIP reduced proteasome activity in the synaptic (A), but not the cytoplasmic (B), fraction from the rat lateral amygdala. This figure has been modified from Jarome, T. J. et al.²³.

Figure 6. Subcellular differences in linkage-specific protein ubiquitination in the lateral amygdala of the same animal following learning. (A) There was an increase in overall ubiquitination in the nuclear fraction following learning, which correlated with a decrease in the cytoplasmic fraction. (B) There was an increase in linear ubiquitination in the nuclear, but not cytoplasmic or synaptic, fraction following learning. (C) There was an increase in K63 ubiquitination in the nuclear fraction following learning, which correlated with a decrease in the synaptic fraction. (D) K48 ubiquitination increased in the nuclear and cytoplasmic, but not synaptic, fraction following learning. * $P < 0.05$ from Control. All obtained ubiquitin optical

densities were normalized to β -actin levels (lower representative Western blot images in A), which was used as a loading control. This figure has been modified from Orsi, S. A. et al.²¹.

DISCUSSION

Here, we demonstrate an efficient method for quantifying changes in ubiquitin-proteasome activity across different subcellular compartments in the same animal. Currently, most attempts at measuring subcellular changes in activity of the ubiquitin-proteasome system have been limited to a single compartment per sample, resulting in the need to repeat experiments. This leads to significant costs and loss of animal life. Our protocol alleviates this problem by splitting hemispheres, allowing different cellular fractions to be collected from each hemisphere of the same animal. Using this protocol, we were able to show for the first time that in the same animal ubiquitin-proteasome activity differentially changes in nuclear, cytoplasmic and synaptic fractions in response to learning²¹.

The main limitation of the protocol described here is it is reliant upon the amount of brain tissue (sample) obtained. For example, as outlined above, this protocol requires splitting hemispheres of a given brain region. However, this may not always be possible, such as in the case of certain areas that are only present in one hemisphere. In these cases, the protocol could be modified by first homogenizing the entire brain region in the TEVP buffer used in the synaptic fraction step (Section 3.1), since this buffer is free of all denaturing agents. The sample can then be split into two equal parts by volume. The first part can be used for the synaptic fraction, following the protocol as described. For the second half of the sample, the non-ionic detergent NP-40 can be added to a final concentration of 0.05%, followed by centrifuging as described in Section 2.6. This will allow the separation of cytoplasmic proteins into the supernatant and nuclear proteins into the pellet, which can be further isolated following the remaining steps in Section 2. Another concern in tissue quantity is that some brain regions are very small in size, such as the prelimbic cortex. However, in these cases, the above protocol can still be used by reducing the volumes of the buffers used, which would have to be determined empirically based on the size of the brain region collected. Thus, this protocol can be amended to even the more difficult brain regions in which less tissue is available.

One of the major advantages of the protocol we outline here is that it uses common laboratory equipment and reagents that can be found in most facilities, allowing this methodology to be amendable to those even on a restricted budget or with limited resources. Furthermore, while we outline this protocol as a way of measuring subcellular changes in ubiquitin-proteasome signaling, this methodology can also be applied to any other protein or cellular process in which understanding the cellular localization and function is important. Thus, this protocol could have broad applications to understanding the subcellular functions of specific proteins or complexes during learning and memory or different disease states.

ACKNOWLEDGMENTS

This work was supported by startup funds from the College of Agricultural and Life Sciences and the College of Science at Virginia Tech. T.M. is supported by the George Washington Carver Program at Virginia Tech.


DISCLOSURES


The authors have nothing to disclose.


REFERENCES


- 1 Hershko, A. & Ciechanover, A. The ubiquitin system. *Annu Rev Biochem.* **67** 425-479, (1998).
- 2 Akutsu, M., Dikic, I. & Bremm, A. Ubiquitin chain diversity at a glance. *Journal of Cell Science.* **129** (5), 875-880, (2016).
- 3 Ravid, T. & Hochstrasser, M. Diversity of degradation signals in the ubiquitin-proteasome system. *Nature Reviews Molecular Cell Biology.* **9** (9), 679-690, (2008).
- 4 Erpapazoglou, Z., Walker, O. & Hagenauer-Tsapis, R. Versatile roles of k63-linked ubiquitin chains in trafficking. *Cells.* **3** (4), 1027-1088, (2014).
- 5 Iwai, K., Fujita, H. & Sasaki, Y. Linear ubiquitin chains: NF-kappaB signalling, cell death and beyond. *Nature Review Molecular Cell Biology.* **15** (8), 503-508, (2014).
- 6 Rieser, E., Cordier, S. M. & Walczak, H. Linear ubiquitination: a newly discovered regulator of cell signalling. *Trends in Biochemical Sciences.* **38** (2), 94-102, (2013).
- 7 Collins, G. A. & Goldberg, A. L. The Logic of the 26S Proteasome. *Cell.* **169** (5), 792-806, (2017).
- 8 Bhat, K. P. et al. The 19S proteasome ATPase Sug1 plays a critical role in regulating MHC class II transcription. *Molecular Immunology.* **45** (8), 2214-2224, (2008).
- 9 Ezhkova, E. & Tansey, W. P. Proteasomal ATPases link ubiquitylation of histone H2B to methylation of histone H3. *Molecular Cell.* **13** (3), 435-442, (2004).
- 10 Jarome, T. J. & Helmstetter, F. J. The ubiquitin-proteasome system as a critical regulator of synaptic plasticity and long-term memory formation. *Neurobiology of Learning and Memory.* **105** 107-116, (2013).
- 11 Hegde, A. N. & Upadhy, S. C. The ubiquitin-proteasome pathway in health and disease of the nervous system. *Trends in Neuroscience.* **30** (11), 587-595, (2007).
- 12 Adori, C. et al. Subcellular distribution of components of the ubiquitin-proteasome system in non-diseased human and rat brain. *Journal of Histochemistry and Cytochemistry.* **54** (2), 263-267, (2006).
- 13 Upadhy, S. C., Ding, L., Smith, T. K. & Hegde, A. N. Differential regulation of proteasome activity in the nucleus and the synaptic terminals. *Neurochemistry International.* **48** (4), 296-305, (2006).
- 14 Enenkel, C., Lehmann, A. & Kloetzel, P. M. Subcellular distribution of proteasomes implicates a major location of protein degradation in the nuclear envelope-ER network in yeast. *EMBO Journal.* **17** (21), 6144-6154, (1998).
- 15 Jarome, T. J., Kwapis, J. L., Ruenzel, W. L. & Helmstetter, F. J. CaMKII, but not protein kinase A, regulates Rpt6 phosphorylation and proteasome activity during the formation of long-term memories. *Frontiers in Behavioral Neuroscience.* **7** 115, (2013).
- 16 Jarome, T. J., Werner, C. T., Kwapis, J. L. & Helmstetter, F. J. Activity dependent protein degradation is critical for the formation and stability of fear memory in the amygdala. *PLoS One.* **6** (9), e24349, (2011).



- 17 Lopez-Salon, M. et al. The ubiquitin-proteasome cascade is required for mammalian long-term memory formation. *European Journal of Neuroscience*. **14** (11), 1820-1826, (2001).
- 18 Reis, D. S., Jarome, T. J. & Helmstetter, F. J. Memory formation for trace fear conditioning requires ubiquitin-proteasome mediated protein degradation in the prefrontal cortex. *Frontiers in Behavior Neuroscience*. **7** 150, (2013).
- 19 Rosenberg, T., Elkobi, A. & Rosenblum, K. mAChR-dependent decrease in proteasome activity in the gustatory cortex is necessary for novel taste learning. *Neurobiology of Learning and Memory*. **135** 115-124, (2016).
- 20 Rosenberg, T., Elkobi, A., Dieterich, D. C. & Rosenblum, K. NMDAR-dependent proteasome activity in the gustatory cortex is necessary for conditioned taste aversion. *Neurobiology of Learning and Memory*. **130** 7-16, (2016).
- 21 Orsi, S. A. et al. Distinct subcellular changes in proteasome activity and linkage-specific protein polyubiquitination in the amygdala during the consolidation and reconsolidation of a fear memory. *Neurobiology of Learning and Memory*. **157** 1-11, (2019).
- 22 Cullen, P. K., Ferrara, N. C., Pullins, S. E. & Helmstetter, F. J. Context memory formation requires activity-dependent protein degradation in the hippocampus. *Learning and Memory*. **24** (11), 589-596, (2017).
- 23 Jarome, T. J., Ferrara, N. C., Kwapis, J. L. & Helmstetter, F. J. CaMKII regulates proteasome phosphorylation and activity and promotes memory destabilization following retrieval. *Neurobiology of Learning and Memory*. **128** 103-109, (2016).
- 24 Lee, S. H. et al. Synaptic protein degradation underlies destabilization of retrieved fear memory. *Science*. **319** (5867), 1253-1256, (2008).
- 25 Dunah, A. W. & Standaert, D. G. Dopamine D1 receptor-dependent trafficking of striatal NMDA glutamate receptors to the postsynaptic membrane. *Journal of Neuroscience*. **21** (15), 5546-5558, (2001).
- 26 Kelly, P. T. & Cotman, C. W. Synaptic proteins. Characterization of tubulin and actin and identification of a distinct postsynaptic density polypeptide. *Journal of Cell Biology*. **79** (1), 173-183, (1978).
- 27 Mengual, E., Arizti, P., Rodrigo, J., Gimenez-Amaya, J. M. & Castano, J. G. Immunohistochemical distribution and electron microscopic subcellular localization of the proteasome in the rat CNS. *Journal of Neuroscience*. **16** (20), 6331-6341, (1996).














Name of Material/ Equipment	Company	Catalog Number
0.5M EDTA	Fisher	15575020
20S Proteasome Activity Kit	Millipore Sigma	APT280
ATP	Fisher	FERR1441
Beta-actin antibody	Cell signaling	4967S
Beta-tubulin antibody	Cell signaling	2128T
BioTek Synergy H1 plate reader	BioTek	VATECHH1MT3
B-mercaptoethanol	Fisher	ICN19024280
clasto lactacystin β -lactone	Millipore Sigma	L7035
Cryogenic cup	Fisher	033377B
DMSO	DMSO	D8418
DTT	Millipore Sigma	D0632
Glycerol	Millipore Sigma	G5516
H3 antibody	Abcam	ab1791
HEPES	Millipore Sigma	H3375
Hydrochloric acid	Fisher	SA48
IGEPAL (NP-40)	Millipore Sigma	I3021
K48 Ubiquitin Antibody	Abcam	ab140601
K63 Ubiquitin Antibody	Abcam	ab179434
KCl	Millipore Sigma	P9541
KONTES tissue grinder	VWR	KT885300-0002
Laemmli sample buffer	Bio-rad	161-0737

Linear Ubiquitin Antibody	Life Sensors	AB-0130-0100
MgCl	Millipore Sigma	442611
Microcentrifuge	Eppendorf	2231000213
myr-AIP	Enzo Life Sciences	BML-P212-0500
NaCl	Millipore Sigma	S3014
Odyssey Fc Imaging System	LiCor	2800-02
Phosphatase Inhibitor	Millipore Sigma	524625
Precision Plus Protein Standard	Bio-rad	161-0373
Protease Inhibitor	Millipore Sigma	P8340
PSD95 antibody	Cell signaling	3450T
SDS	Millipore Sigma	L3771
Sodium hydroxide	Fisher	SS255
Sucrose	Millipore Sigma	S0389
TBS	Alfa Aesar	J62938
Tris	Millipore Sigma	T1503
Tween-20	Fisher	BP337-100
Ubiquitin Antibody	Enzo Life Sciences	BML-PW8810

Comments/Description

Various other vendors

Other vendors carry different versions

Various other vendors

Various other vendors

Various other vendors

Other vendors carry different versions

Various other vendors

Various other vendors

Various other vendors

Various other vendors

Various other vendors

Various other vendors

Various other vendors

Various other vendors

Various other vendors

Various other vendors

Various other vendors

Various other vendors

Various other vendors

Various other vendors

Various other vendors

Only M1 antibody

Various other vendors

Various other manufacturers/models

Carried by Millipore-Sigma

Various other vendors

Other vendors carry different versions

Various other vendors

Various other vendors

Various other vendors

Various other vendors

Various other vendors

Various other vendors

Various other vendors

Various other vendors

Various other vendors

Various other vendors

Various other vendors



1 Alewife Center #200
Cambridge, MA 02140
tel. 617.945.9051
www.jove.com

ARTICLE AND VIDEO LICENSE AGREEMENT

Title of Article:

Quantifying Subcellular Ubiquitin-proteasome activity in the rodent brain

Author(s):

Taylor McFadden, Rishi K. Devulapalli, Timothy J. Jarome

Item 1: The Author elects to have the Materials be made available (as described at <http://www.jove.com/publish>) via:

☒ Standard Access

☐ Open Access

Item 2: Please select one of the following items:

☒ The Author is **NOT** a United States government employee.

☐ The Author is a United States government employee and the Materials were prepared in the course of his or her duties as a United States government employee.

☐ The Author is a United States government employee but the Materials were NOT prepared in the course of his or her duties as a United States government employee.

ARTICLE AND VIDEO LICENSE AGREEMENT

1. **Defined Terms.** As used in this Article and Video License Agreement, the following terms shall have the following meanings: "**Agreement**" means this Article and Video License Agreement; "**Article**" means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein; "**Author**" means the author who is a signatory to this Agreement; "**Collective Work**" means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; "**CRC License**" means the Creative Commons Attribution-Non Commercial-No Derivs 3.0 Unported Agreement, the terms and conditions of which can be found at: <http://creativecommons.org/licenses/by-nc-nd/3.0/legalcode>; "**Derivative Work**" means a work based upon the Materials or upon the Materials and other pre-existing works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; "**Institution**" means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; "**JoVE**" means MyJoVE Corporation, a Massachusetts corporation and the publisher of The Journal of Visualized Experiments; "**Materials**" means the Article and / or the Video; "**Parties**" means the Author and JoVE; "**Video**" means any video(s) made by the Author, alone or in conjunction with any other parties, or by JoVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion

of the Article, and in which the Author may or may not appear.

2. **Background.** The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JoVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.

3. **Grant of Rights in Article.** In consideration of JoVE agreeing to publish the Article, the Author hereby grants to JoVE, subject to Sections 4 and 7 below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the "Open Access" box has been checked in Item 1 above, JoVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.

612542.6 For questions, please contact us at submissions@jove.com or +1.617.945.9051.

ARTICLE AND VIDEO LICENSE AGREEMENT

4. **Retention of Rights in Article.** Notwithstanding the exclusive license granted to JoVE in Section 3 above, the Author shall, with respect to the Article, retain the non-exclusive right to use all or part of the Article for the non-commercial purpose of giving lectures, presentations or teaching classes, and to post a copy of the Article on the Institution's website or the Author's personal website, in each case provided that a link to the Article on the JoVE website is provided and notice of JoVE's copyright in the Article is included. All non-copyright intellectual property rights in and to the Article, such as patent rights, shall remain with the Author.

5. **Grant of Rights in Video – Standard Access.** This Section 5 applies if the "Standard Access" box has been checked in Item 1 above or if no box has been checked in Item 1 above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby acknowledges and agrees that, Subject to Section 7 below, JoVE is and shall be the sole and exclusive owner of all rights of any nature, including, without limitation, all copyrights, in and to the Video. To the extent that, by law, the Author is deemed, now or at any time in the future, to have any rights of any nature in or to the Video, the Author hereby disclaims all such rights and transfers all such rights to JoVE.

6. **Grant of Rights in Video – Open Access.** This Section 6 applies only if the "Open Access" box has been checked in Item 1 above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby grants to JoVE, subject to Section 7 below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Video in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Video into other languages, create adaptations, summaries or extracts of the Video or other Derivative Works or Collective Works based on all or any portion of the Video and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. For any Video to which this Section 6 is applicable, JoVE and the Author hereby grant to the public all such rights in the Video as provided in, but subject to all limitations and requirements set forth in, the CRC License.

7. **Government Employees.** If the Author is a United States government employee and the Article was prepared in the course of his or her duties as a United States government employee, as indicated in Item 2 above, and any of the licenses or grants granted by the Author hereunder exceed the scope of the 17 U.S.C. 403, then the rights granted hereunder shall be limited to the maximum

rights permitted under such statute. In such case, all provisions contained herein that are not in conflict with such statute shall remain in full force and effect, and all provisions contained herein that do so conflict shall be deemed to be amended so as to provide to JoVE the maximum rights permissible within such statute.

8. **Protection of the Work.** The Author(s) authorize JoVE to take steps in the Author(s) name and on their behalf if JoVE believes some third party could be infringing or might infringe the copyright of either the Author's Article and/or Video.

9. **Likeness, Privacy, Personality.** The Author hereby grants JoVE the right to use the Author's name, voice, likeness, picture, photograph, image, biography and performance in any way, commercial or otherwise, in connection with the Materials and the sale, promotion and distribution thereof. The Author hereby waives any and all rights he or she may have, relating to his or her appearance in the Video or otherwise relating to the Materials, under all applicable privacy, likeness, personality or similar laws.

10. **Author Warranties.** The Author represents and warrants that the Article is original, that it has not been published, that the copyright interest is owned by the Author (or, if more than one author is listed at the beginning of this Agreement, by such authors collectively) and has not been assigned, licensed, or otherwise transferred to any other party. The Author represents and warrants that the author(s) listed at the top of this Agreement are the only authors of the Materials. If more than one author is listed at the top of this Agreement and if any such author has not entered into a separate Article and Video License Agreement with JoVE relating to the Materials, the Author represents and warrants that the Author has been authorized by each of the other such authors to execute this Agreement on his or her behalf and to bind him or her with respect to the terms of this Agreement as if each of them had been a party hereto as an Author. The Author warrants that the use, reproduction, distribution, public or private performance or display, and/or modification of all or any portion of the Materials does not and will not violate, infringe and/or misappropriate the patent, trademark, intellectual property or other rights of any third party. The Author represents and warrants that it has and will continue to comply with all government, institutional and other regulations, including, without limitation all institutional, laboratory, hospital, ethical, human and animal treatment, privacy, and all other rules, regulations, laws, procedures or guidelines, applicable to the Materials, and that all research involving human and animal subjects has been approved by the Author's relevant institutional review board.

11. **JoVE Discretion.** If the Author requests the assistance of JoVE in producing the Video in the Author's facility, the Author shall ensure that the presence of JoVE employees, agents or independent contractors is in accordance with the relevant regulations of the Author's institution. If more than one author is listed at the beginning of this Agreement, JoVE may, in its sole

ARTICLE AND VIDEO LICENSE AGREEMENT

discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JoVE reserves the right, in its absolute and sole discretion and without giving any reason therefore, to accept or decline any work submitted to JoVE. JoVE and its employees, agents and independent contractors shall have full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JoVE has sole discretion as to the method of making and publishing the Materials, including, without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

12. **Indemnification.** The Author agrees to indemnify JoVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JoVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JoVE, making of videos by JoVE, or publication in JoVE or elsewhere by JoVE. The Author shall be responsible for, and shall hold JoVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to

the making of a video by JoVE its employees, agents or independent contractors. All sterilization, cleanliness or decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's expense. All indemnifications provided herein shall include JoVE's attorney's fees and costs related to said losses or damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JoVE, its employees, agents or independent contractors.

13. **Fees.** To cover the cost incurred for publication, JoVE must receive payment before production and publication the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JoVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.

14. **Transfer, Governing Law.** This Agreement may be assigned by JoVE and shall inure to the benefits of any of JoVE's successors and assignees. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to be one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement is required per submission.

CORRESPONDING AUTHOR

Name:

Timothy Jarome

Department:

Animal and Poultry Sciences and School of Neuroscience

Institution:

Virginia Polytechnic Institute and State University

Title:

Assistant Professor

Signature:



Date:

1/11/2019

Please submit a signed and dated copy of this license by one of the following three methods:

1. Upload an electronic version on the JoVE submission site
2. Fax the document to +1.866.381.2236
3. Mail the document to JoVE / Attn: JoVE Editorial / 1 Alewife Center #200 / Cambridge, MA 02140

612542.6 For questions, please contact us at submissions@jove.com or +1.617.945.9051.

Response to Editors

1. Please take this opportunity to thoroughly proofread the manuscript to ensure that there are no spelling or grammar issues.

Response: We have thoroughly reviewed the manuscript to ensure no spelling or grammar issues are present

2. JoVE cannot publish manuscripts containing commercial language. This includes company names before an instrument or reagent. Please remove all commercial language from your manuscript and use generic terms instead. All commercial products should be sufficiently referenced in the Table of Materials and Reagents.

Response: All references to commercial products have been removed and are instead listed in the Table of Materials and Reagents

3. Please split some long steps into two or more sub-steps so that each step contains 2-3 action and is less than 4 lines.

Response: We have separated steps into smaller parts. A few steps are still 4 lines, which is needed because no logical division can be made. However, most steps are 3 lines or less now. Notes that are present in some sections remain more than 4 lines and were included to address reviewer comments.

4. Please revise the text in Protocol to avoid the use of any personal pronouns (e.g., "we", "you", "our" etc.).

Response: All personal pronouns have been removed from the Protocol section

5. Step 4.1: Please ensure that all text is written in the imperative tense.

Response: This has been corrected

6. 5.1: Please ensure that all text is written in the imperative tense.

Response: This has been corrected

7. 5.1.1: Please ensure that all text is written in the imperative tense.

Response: This has been corrected

8. 5.1.2: Please ensure that all text is written in the imperative tense.

Response: This has been corrected

9. 5.1.3: Please ensure that all text is written in the imperative tense.

Response: This has been corrected

10. 5.2: Please ensure that all text is written in the imperative tense.

Response: This has been corrected

11. Please do not abbreviate journal titles for all references.

Response: This mistake was due to an error with Endnote and has been corrected

Response to Reviewers

We would like to thank the reviewers for their insightful comments and suggestions on our original manuscript. With these comments in mind, we have made a number of changes that we feel have strengthened our manuscript. Below is a point-by-point response to each concern raised by the reviewers for the previous version of our manuscript. Additionally, all changes are marked in red in the manuscript text.

Reviewer #1

1) In figure 3, it is difficult to observe Histone H3 in the nuclear fraction, and there isn't a blot to show lack of contamination from the cytoplasmic fraction (e.g. tubulin?) - CREB is both cytoplasmic and nuclear.

Response: Thank you for the suggestion of using tubulin as a cytoplasmic marker. We have included that image in revised Figure 2, which shows largely an absence of tubulin from the nuclear fraction. Additionally, we have replaced the Histone H3 image with a stronger development, making it easier to see the H3 protein.

2) In figures 4 and 5, the synaptic proteasome activity varies a lot. The authors do mention in the text that the ionic buffer used to collect the synaptic fraction could interfere with proteasome activity. The worry with this is how repeatable and reliable these assays will be if there are issues from the buffer used. Perhaps the buffer can be optimized or the discussion can be more thorough on this point?

Response: Thank you for pointing this out. The samples within the synaptic fraction (Fig. 3) contain a combination of samples from two groups, which have opposite patterns of proteasome activity; this can account for some of the variability observed. The synaptic fraction shown in Fig. 4 has the 2 groups separated out, which should help visualize how the variability can be attributed to this group effect. Furthermore, although the buffer used for the synaptic fraction contains an ionic detergent, we have used it during other experiments and can replicate results across experiments (e.g., Jarome et al., 2016, *Neurobiol Learn Mem*). However, we do acknowledge that the ionic detergents could be problematic in some preparations. To account for this, we have added a note within the revised manuscript in section 3.2.1 about the use of ionic detergents in our lysis buffer.

3) In figure 7: it is difficult to understand how those blots in figure 7A can actually be quantified. How are they normalized? For example, how did the authors ensure equal loads among different fractions? Also, are all of the images from the same blot, cropped and rearranged? I am assuming that the images from A do not really correspond to quantifications from B-E. Would it be possible to see blots that were used for actual quantifications? How were they quantified? How were these done, compared, loaded, probed? For figure 7A, kD ladders are necessary for the viewer to know where different bands are running, especially to also see if there is any mono-UB present.

Response: Thank you for pointing this out. All of our western blots were normalized to β -actin, which was used as a loading control, though we never explicitly stated that in the previous version of the manuscript. We have now included text in the manuscript to indicate this (Figure 6 caption). For the western blot images, originally we meant to show that the different ubiquitin tags/modifications could be visualized in all subcellular compartments using our protocol. However, we now realize that this may be confusing since these images do not reflect the graphs presented in the second part of the figure. In the revised manuscript we have corrected this by including representative western blot images for each of the bar graphs. Additionally, we have included ladder markers and actin loading control images.

4) Lastly, the proteasome activity assays described ONLY report 20S core activity, not ub-dependent activity. This specific and highly important differentiation needs to be made clear at multiple stages, otherwise those who are not well versed with the Ub field might misinterpret the use and outcome of the results. While it is beneficial to know the activity of the core, it is not, in any direct way, a read-out of the actual proteasome activity. It can only be used as a proxy of proteasome levels and if the core is active, but not whether the activity of the ub-proteasome system itself is compromised or enhanced. This is a significant limitation of the current methodology out there, and not the authors' fault. Nevertheless, I think that the authors should address this head-on from the very beginning, perhaps even in the abstract.

Response: This is an excellent point raised by Reviewer #1. This is an important distinction that needs to be made and was not clearly stated in the previous version of our manuscript. In the revised manuscript we have added information about the indirect nature of this assay where appropriate, including the abstract, introduction and methods (section 4.1.).

Reviewer #2

1) All buffers used were supplemented with a cocktail of protease inhibitors. Because this cocktail of inhibitors can also inhibit the proteasome, how can we assure that the proteasome activity measured is not underestimated at the end?

Response: This is an interesting question raised by Reviewer #2. While it is possible that our protease inhibitor cocktail could cause a slight reduction in proteasome activity, it is essential to preserve proteins and ensure that the data obtained from this assay can be attributed directly to proteasome activity as opposed to other proteases. As shown in Figure 3G-H, application of the highly specific proteasome inhibitor β lac completely blocks AMC release in our assay, suggesting that all AMC release is likely due specifically to the proteasome. However, you are right in stating that overall proteasome activity may be reduced some, providing an underestimation in the end. We have included a brief statement acknowledging this in the revised manuscript, Section 2.1.2.

2) The buffers used (e.g. Lysis buffer, Extraction buffer and TEVP buffer) lack glycerol and/or

ATP. This may induce proteasome disassembly, favoring the presence of only the catalytic core. Can you please comment on this?

Response: It is correct that our buffers generally lack glycerol and ATP, with the exception of the nuclear extraction buffer. While this may not result in the perfect conditions in which to extract proteasomes, we are able to consistently recover subcellular fractions that retain high levels of proteasome function, based on activity levels using our *in vitro* proteasome activity assay. However, the addition of glycerol and ATP could potentially improve proteasome recovery/activity on our assay. To ensure that the reader knows all options available, we have included a brief note to sections 2 and 3 indicating that the buffers could be supplemented with ATP and glycerol if desired.

3) Step 4.13.4: the way the results were calculated, as percentage of increase in fluorescence considering the last and the first scan, gives values that are above 100%, which may induce some misinterpretation. This may suggest that the proteasome is being stimulated under the conditions used. The authors may consider expressing the results using the slope of the line calculated based on the measurements performed.

Response: Thank you for this comment. Originally, we added the percentage of baseline calculation to account for any variance in baseline readings. However, we realize that this may be misleading in terms of describing proteasome activity. As the steps listed in 4.13.1 - 4.13.3 provide the necessary steps for calculation of proteasome activity against the AMC standard curve, we have removed section 4.13.4 in the revised manuscript and modified Figure 3 for the use of this analysis.

4) Figure 7A: In order to allow comparing the results obtained with the different subcellular fractions (for each antibody), the same membrane should be shown. Moreover, the results obtained with a loading control are missing.

Response: These images have been corrected and loading control information included in the revised manuscript. Please see response to Reviewer #1 comment 3 for more information.

5) Information should be provided i) regarding the methodology used in sample denaturation and ii) the reducing agent used for the immunoblotting experiments.

Response: For denaturing, normalized samples were mixed with an equal volume of Laemmli sample buffer (Bio-rad #1610737) supplemented with β -mercaptoethanol (Fisher Scientific #ICN19024280) to a percentage of 5% by volume, per the manufacturer's instruction. We have added this information to the manuscript in section 5.1 and the materials list.

6) Can you please provide information about the company and catalog number for the anti-H3, -PSD95 and -CREB antibody

Response: We have added the company and catalog information to the materials list. The image including the CREB antibody was not used in the revised manuscript, so we did not provide the information on that antibody.

7) Please provide information regarding the plate reader and the imaging system used to measure proteasome activity and in the immunoblot experiments, respectively.

Response: We have added this information to the materials list.

8) Page 11 (lines 375-379): These sentences suggest that proteasomes are more abundant in the cytoplasm than in the nuclear compartment, which contrasts with the results previously published (J Neurosci 16, 6331-6341 [1996]). This should be clarified.

Response: We apologize for incorrectly stating that proteasomes were most abundant in the cytoplasm. We have corrected this statement to indicate that proteasomes are abundant in the nucleus and cytoplasm.

9) The results of Figure 6 suggest that CaMKII activity regulates proteasome activity. This is surprising considering that the lysis and the assay buffers contain EDTA which chelates Ca²⁺. This should also be clarified.

Response: This is an interesting point. The AIP peptide acts as a selective inhibitor of CaMKII, targeting its autoinhibitory domain. In our brain lysates some CaMKII is in an autophosphorylated state, which likely occurred from calcium influx when the animal was alive (especially after learning). In the lysed samples, the CaMKII phosphorylation state is preserved but can be manipulated by selective inhibitors such as AIP. Without AIP, CaMKII that is already phosphorylated can target the proteasome, controlling its activity. When AIP is applied, the amount of phosphorylated CaMKII is reduced, leading to a reduction in proteasome activity. While EDTA chelates Ca²⁺, it would have no effect on CaMKII that is already in a phosphorylated state. As a result, we can still manipulate proteasome function in our lysates via targeting of CaMKII that is already in a phosphorylated state.

10) Page 5, line 109: for rapid flash freezing perhaps it would advisable to use liquid nitrogen.

Response: We have added this information to the revised manuscript (Section 1.1).

11) Page 5, line 118: should the microscope slide be chilled in order to preserve the tissue?

Response: Yes, the slide should be chilled. Thank you for pointing out that we did not include this important information. We have added this information to the revised manuscript (Section 1.2.2.).

12) Page 6 (line 149): the sentence should be corrected to '... 1.5 ml centrifuge microtube...'

Response: We made the correction to the revised manuscript.

13) Page 6 (sections 2.4 and 3.4): The supplier and reference of the glass Teflon homogenizer used should be provided.

Response: This information has been added to the material list, per journal instructions

14) Page 8 (line 268): the sentence should be corrected to '...substrate (Suc-LLVY-AMC) provided...'

Response: We made this correction in the revised manuscript.

15) Page 9, step 4.7: although it is understandable that the authors always use the same amount of protein for measurement of the proteasome activity, the changes in protein volume (directly influenced by the protein concentration in each sample) loaded in each well may increase pipetting errors. Therefore, it would be advisable to normalize protein concentration to the less concentrated sample, and from there using the same sample volume in all conditions.

Response: Thank you for suggesting this alternative. In the past we have used this method for protein normalization and have obtained similar results. In the revised manuscript we include this method of normalization as option for the reader/user

16) Page 9, step 4.9: does the plate need to be under constant agitation in order to keep the sample homogeneous? And during the measurement of proteasome activity?

Response: Thank you for asking this question. The plate does not need to be under constant agitation to get reliable signal on the proteasome activity assay, though can be should the user choose to. We have included this information in the revised manuscript (Section 4.12)

17) Page 10, step 5.1.2: incubation of the membrane with the secondary antibody after stripping should be recommended in order to confirm that membrane stripping was done properly.

Response: We have added the additional recommended information to the revised manuscript.

18) Page 10, step 5.2: please indicate which software was used to develop the immunoblots.

Response: The software used, Image Studio Ver 5.2, is provided by our imaging system, which is listed in the materials section. As the software is included with the imaging system, we did not include a separate listing for it in the revised manuscript.

19) The composition of the TBS solution should be indicated.

Response: We buy our TBS commercially. We have included the TBS product information in the materials list of the revised manuscript.

20) Figure 1: The units of temperature are missing.

Response: We apologize for this error, which occurred during conversion of the image to a tif format. The revised Figure 1 includes the temperature information

21) Figures 2 and 3 should be merged in a single figure.

Response: Thank you for this suggestion. In the revised manuscript we have merged Figures 2 and 3

22) Figure 7A: please indicate which protein standards were used in the Western blot experiments, the supplier and catalog number.

Response: We have added this information to the materials list of the revised manuscript.

23) Figure 2A-B: please indicate the molecular weight of the protein markers.

Response: All molecular weight markers are now included in the revised Figure 2

24) Can you please indicate the company and catalog number for Tween?

Response: We have added this information to the materials list of the revised manuscript

25) Clasto-lactacystin β -lactone is already included in the 20S Proteasome activity kit. It should be deleted from the list of reagents.

Response: The proteasome inhibitor lactacystin is included in the 20S Proteasome activity kit. While it is the parent compound to clasto lactacystin β -lactone, they have very different potencies. Since we choose to use the more potent clasto-lactacystin β -lactone, we have included it in the list of reagents because it has to be purchased separately from the proteasome kit.

Reviewer #3

1) The proteasome activity measured is that of the core of the proteasome (called the 20S proteasome) which does not degrade proteins tagged by ubiquitin (i.e. ubiquitinated proteins). While measuring this activity might be useful, it does not accurately reflect the activity of the 26S proteasome which is the type of proteasome that degrades ubiquitinated proteins. The

authors need to add a caveat about this. They could perhaps suggest that the subcellular fractions can be assayed for 26S proteasome activity by citing appropriate published papers.

Response: This is an excellent point raised by Reviewer #3. In the revised manuscript, we now explicitly state the limitation of this assay. Please see Reviewer #1 comment 4 response for more information.

2) In order for others to carry out the protocol successfully, the age, strain and the sex of the rat need to be specified.

Response: We have added this information to the revised manuscript.

3) For the synaptic fraction, data are shown for presence of postsynaptic components. Does the synaptic fraction have presynaptic components? Perhaps the authors can add a sentence about the need to characterize this fraction further. Or, if evidence for presences (or absence) of presynaptic components is available from published literature (that is applicable to the present study), such evidence should be provided.

Response: Thank you for asking this important question. The synaptic preparation we use includes presynaptic and postsynaptic proteins (Dunah and Standaert, 2001, *J Neurosci*). We now acknowledge this important information in the representative results, including this relevant citation.

4) For rigorous quantification of western blots (Fig. 7), the signal needs to be normalized to another protein which is known not to change under the experimental conditions. Just relying on protein quantification determined by Bradford assay to ensure equal loading is not adequate.

Response: Thank you for pointing this out. All of our western blots were normalized to β -actin, which was used as a loading control. This information has been added to the revised manuscript (Figure 6 caption). Please see response to Reviewer #1 comment 3 for more information regarding this.