

Multiwell-Screen

Application Commands Settings Help

Settings

Raw Data (Well A1)

21 31

Experimental Flow

Recording Phase

Recording Phase Timing: Timer

Recording Phase Duration: 2 min

Analysis Phase

Restrict Analysis to Time Segment

Analysis Phase Start: 0 s

Analysis Phase Duration: 2 min

Wash-In Phase

Use Wash-In Phase

Wash-In Timing: Timer

Wash-In Duration: 2 min

OK Cancel

Data Source Settings

Data Source

Sampling Rate: 20000 Hz

Filter Settings

High-Pass Filter: 2. order 0.1 Hz Make Permanent

Low-Pass Filter: 4. order 3500 Hz Make Permanent

OK Cancel

Multiwell Selection

1 2

A B C D

Temperature Control

Temperature [°C]

36 34 32 30 28 26 24 22

4.6 4.8 5 5.2 5.4

Disable Setpoint: 37.0 °C Range: 1 min

Raw Data Overview Averages (Well A1) Raw Data (Well A1)

0 0.5 1 1.5 2 2.5 3 3.5 4 4.5 5

Time [s]