

# Journal of Visualized Experiments

## Activated Cross-Linked Agarose for the Rapid Development of Affinity Chromatography Resins - Antibody Capture as a Case Study

--Manuscript Draft--

Article Type:	Invited Methods Article - JoVE Produced Video
Manuscript Number:	JoVE59933R2
Full Title:	Activated Cross-Linked Agarose for the Rapid Development of Affinity Chromatography Resins - Antibody Capture as a Case Study
Keywords:	Protein A replacements; chromatography resin; downstream processing (DSP); design of experiments (DoE); plant-made pharmaceuticals (PMPs); primary recovery
Corresponding Author:	Johannes Felix Buyel, Dr. Fraunhofer-Gesellschaft zur Förderung der angewandten Forschung e. V. Aachen, NRW GERMANY
Corresponding Author's Institution:	Fraunhofer-Gesellschaft zur Förderung der angewandten Forschung e. V.
Corresponding Author E-Mail:	johannes.buyel@rwth-aachen.de
Order of Authors:	Matthias Knödler Clemens Rühl Patrick Opdensteinen Johannes Felix Buyel
Additional Information:	
Question	Response
Please indicate whether this article will be Standard Access or Open Access.	Open Access (US\$4,200)
Please indicate the <b>city, state/province, and country</b> where this article will be <b>filmed</b> . Please do not use abbreviations.	Aachen, Germany

**TITLE:**

**Activated Cross-Linked Agarose for the Rapid Development of Affinity Chromatography Resins - Antibody Capture as a Case Study**

**AUTHORS:**

Knödler, Matthias<sup>1</sup>

[matthias.knoedler@ime.fraunhofer.de](mailto:matthias.knoedler@ime.fraunhofer.de)

Rühl, Clemens<sup>2</sup>

[clemens.ruehl@gmail.com](mailto:clemens.ruehl@gmail.com)

Opdensteinen, Patrick<sup>1</sup>

[patrick.opdensteinen@ime.fraunhofer.de](mailto:patrick.opdensteinen@ime.fraunhofer.de)

Buyel, Johannes F.<sup>1,3</sup>

[johannes.buyel@ime.fraunhofer.de](mailto:johannes.buyel@ime.fraunhofer.de)

1. Fraunhofer Institute for Molecular Biology and Applied Ecology IME, Fraunhofer-Gesellschaft zur Förderung der angewandten Forschung e. V., Aachen, Germany
2. Sanofi Deutschland GmbH, Industriepark Hoechst, Frankfurt am Main, Germany
3. Institute for Molecular Biotechnology, RWTH Aachen University, Aachen, Germany

**CORRESPONDING AUTHOR:**

Johannes F. Buyel

**KEYWORDS:**

Protein A replacements; chromatography resin; downstream processing (DSP); design of experiments (DoE); plant-made pharmaceuticals (PMPs); primary recovery

**SUMMARY:**

In this procedure, a DsRed-based epitope ligand is immobilized to produce a highly selective affinity resin for the capture of monoclonal antibodies from crude plant extracts or cell culture supernatants, as an alternative to Protein A.

**ABSTRACT:**

The purification of monoclonal antibodies (mAbs) is commonly achieved by Protein A affinity chromatography, which can account for up to 25% of the overall process costs. Alternative, cost-effective capture steps are therefore valuable for industrial-scale manufacturing, where large quantities of a single mAb are produced. Here we present a method for the immobilization of a DsRed-based epitope ligand to a cross-linked agarose resin allowing the selective capture of the HIV-neutralizing antibody 2F5 from crude plant extracts without using Protein A. The linear epitope ELDKWA was first genetically fused to the fluorescent protein DsRed and the fusion protein was expressed in transgenic tobacco (*Nicotiana tabacum*) plants before purification by immobilized metal-ion affinity chromatography. Furthermore, a method based on activated

cross-linked agarose was optimized for high ligand density, efficient coupling and low costs. The pH and buffer composition and the soluble ligand concentration were the most important parameters during the coupling procedure, which was improved using a design-of-experiments approach. The resulting affinity resin was tested for its ability to selectively bind the target mAb in a crude plant extract and the elution buffer was optimized for high mAb recovery, product activity and affinity resin stability. The method can easily be adapted to other antibodies with linear epitopes. The new resins allow gentler elution conditions than Protein A and could also reduce the costs of an initial capture step for mAb production.

## INTRODUCTION:

Biopharmaceutical products are important for the treatment of a broad spectrum of diseases in nearly every branch of medicine<sup>1</sup>. Monoclonal antibodies (mAbs) dominate the biopharmaceutical market, with worldwide sales expected to reach almost €110 billion in 2020<sup>2</sup>. The favored expression platform for mAbs are Chinese hamster ovary cell lines, which typically produce high mAb titers of up to 10 g·L<sup>-1</sup> in the culture supernatant<sup>3,4</sup>. However, the production of mAbs in mammalian cell cultures is expensive due to the high cost of the medium and the need for sterile fermentation<sup>5</sup>. Alternative expression platforms such as plants potentially offer a faster, simpler, less expensive and more scalable approach for industrial manufacturing<sup>6,7</sup>.

In addition to the costs associated with mammalian cell cultures, the widespread use of Protein A affinity chromatography for product capture is a major cost driver for the industrial production of mAbs. Protein A is naturally found on the surface of *Staphylococcus aureus* cells and it binds to the fragment crystallizable (Fc) region of certain murine and human antibodies, thereby acting as a defense mechanism to evade the humoral immune system<sup>8</sup>. Protein A has become the industry gold standard for the capture of mAbs from cell culture supernatants and is also widely used by the research community because it is highly selective, typically achieving mAb purities of ~95% in a single step<sup>8</sup>. Unsurprisingly, sales of Protein A over the last two decades have closely mirrored the sales of mAbs<sup>8</sup>. Depending on the production scale, the costs of Protein A can correspond to more than 25% of the total process costs and thereby affect the market price of therapeutic mAbs, which can be up to €2000 g<sup>-1</sup><sup>5,9</sup>. Therefore, alternative chromatography resins with a similar purification performance have the potential to substantially reduce production costs, making antibody-based therapies accessible for a larger number of patients<sup>10-12</sup>. Such alternatives may also circumvent the disadvantages of Protein A chromatography, including the harsh elution conditions at low pH (typically <3.5) that can potentially cause mAbs to undergo conformational changes that promote aggregation<sup>13</sup>. Importantly, Protein A is selective only for the Fc region of certain IgG subclasses, so non-functional molecules with truncated binding domains may co-purify with the intact product<sup>5</sup>, whereas mAb derivatives such as single-chain variable fragments do not bind to Protein A at all.

Here, we describe an alternative affinity chromatography resin for the capture of the HIV-neutralizing mAb 2F5 using its linear epitope ELDKWA (one letter amino acid code)<sup>5,14</sup>. We genetically fused the 2F5 epitope to the C-terminus of the fluorescent protein DsRed, which functioned as a carrier and reporter molecule, and produced the resulting protein DsRed-2F5-Epitope (DFE) in transgenic tobacco (*Nicotiana tabacum*) plants. DFE was purified by single-step

immobilized metal-ion affinity chromatography (IMAC). The immobilization of the purified DFE affinity ligand onto a cross-linked agarose resin was achieved by chemical coupling using N-hydroxysuccinimide (NHS)-activated cross-linked agarose columns. Statistical experimental designs were then used to optimize the immobilization procedure and coupling efficiency<sup>15</sup>. The purification strategy for mAb 2F5 was evaluated in terms of antibody purity, yield and ligand stability. In contrast to Protein A, which binds the Fc region, DFE bound to the complementarity-determining region of 2F5, ensuring the purification of molecules with an intact paratope. Our concept can easily be adapted to any mAb with a linear epitope or to other peptide-based affinity ligands which can be easily identified by microarray studies<sup>16</sup>.

[Place Figure 1 here]

## **PROTOCOL:**

### **1. Cultivate the transgenic tobacco plants**

NOTE: The design of the DFE fusion protein and the generation of transgenic plants are described elsewhere<sup>5,17</sup>.

1.1) Germinate the transgenic tobacco seedlings in soil. Irrigate the plants with a 1.0 g L<sup>-1</sup> fertilizer solution. Transfer the plants to single, 100 mm x 100 mm x 60 mm (length, width, height) pots when they grow to a diameter of ~0.04 m.

1.2) Cultivate the transgenic tobacco plants in a greenhouse with a 16-h photoperiod (25/22 °C light/dark temperature regime), with 70% relative humidity and automated fertilization with a 1.0 g L<sup>-1</sup> fertilizer solution for 15 min every hour.

1.3) After 7 weeks, harvest all leaves except the four cotyledon leaves at the base of the plant stem. Immediately process the harvested leaves as described below.

### **2. Heat precipitation of host cell proteins**

2.1) Set up a water bath with a 20 L working volume in a stainless-steel heated vessel (0.3 m x 0.3 m x 0.3 m). At the start, place a magnetic pump in the water bath to agitate the liquid permanently with a flow of 5 L min<sup>-1</sup>. Mount an adjustable thermostat on the water bath for temperature control (steps 2.4 and 2.8).

CAUTION: All the following steps up to 2.10 involve the handling of hot liquid. Wear appropriate personal protective equipment including thermally insulated gloves and goggles.

2.2) Add 15 L of deionized water to the water bath and heat the liquid to 70 °C.

2.3) Place a 10 L bucket filled with 5 L of deionized water on a magnetic stirrer. Add sodium phosphate to a final concentration of 120 mM. Adjust the pH to 8.14 using 10 M hydrochloric

acid. When all the components have dissolved, add the resulting concentrated blanching buffer (5 L) to the water bath from step 2.2.

2.4) Agitate the water bath until the temperature reaches 70 °C. Use 10 M hydrochloric acid to adjust the pH to 8.14 if necessary. Continue agitating for at least 15 min after the required temperature and pH are reached to ensure that the entire assembly is in thermal equilibrium.

2.5) Prepare a 20 L bucket (0.3 m x 0.3 m x 0.3 m) filled with deionized water at a temperature of ~17 °C (required for step 2.9).

2.6) Prepare 400 g aliquots of the harvested tobacco leaves from step 1.3. Place one aliquot in a perforated basket (0.2 m x 0.2 m x 0.2 m; pore width at least 0.02 m x 0.02 m). Avoid overfilling the basket with plant material or packing the leaves too tightly to avoid damaging the tissue.

2.7) Fully submerge the basket in the hot blanching buffer from step 2.4 and make sure all leaves remain under the liquid surface. Use a thermostable silicone spoon to hold the leaves under the surface if necessary.

2.8) Incubate the tobacco leaves for 1.5 min in the blanching fluid while the pump is still agitating the liquid. Monitor the liquid temperature during the entire incubation period. Avoid blocking the pump inlet with tobacco leaves.

2.9) Remove the basket of tobacco leaves from the blanching buffer and drain residual buffer from the leaves for 30 s. Transfer the basket to the bucket filled with cold water (from step 2.5) and immerse the leaves for 30 s. Remove the basket and drain residual water from the leaves for 30 s before homogenization (step 3).

2.10) Repeat steps 2.6 to 2.9 with fresh aliquots from 2.6 until the entire harvested biomass is processed. Constantly monitor and adjust the pH and the temperature of the blanching buffer in the water bath.

NOTE: Blanched plant material can be stored on ice for up to 30 min when several blanching cycles are required to process the entire biomass. Blanched leaves can also be stored in vacuum-sealed bags at -80 °C for at least 3 weeks. However, immediate processing is recommended because prolonged storage can reduce the DFE yield.

### **3. Protein extraction and clarification**

CAUTION: The following steps involve a blender with rotating blades. Do not work in the blender vessel when powered or mounted on the motor unit.

3.1) Place 150 g (wet mass) of blanched tobacco leaves (step 2.9) in the blender vessel and add 450 mL of extraction buffer (50 mM sodium phosphate, 500 mM sodium chloride, 10 mM

sodium disulfite, pH 7.5). Close the blender cap tightly to prevent the spilling of plant material or buffer.

3.2) Homogenize the leaves for  $3 \times 30$  s with 30 s breaks between each pulse. Ensure that all leaves are homogenized and that none stick to the top of the blender vessel. If necessary, open the blender during the breaks and push down leaves that stick at the upper part of the vessel, using a clean silicone spoon.

3.3) After homogenization, take a 1.0 mL sample of the homogenate for subsequent analysis (step 7).

3.4) Collect the homogenate in a vessel of adequate size (e.g. if working with a total biomass of 1.0 kg, use a vessel with a capacity of at least a 5 L). Repeat steps 3.1 to 3.3 until all of the blanched biomass is homogenized.

3.5) Mount a bag filter into a filter mount and place another adequately sized vessel (see 3.4) beneath the assembly. Apply the homogenate to the bag filter at a rate of  $\sim 0.15$  L  $\text{min}^{-1}$ . Take samples of the bag filtrate after each liter for subsequent analysis (step 7) and measure the turbidity of the bag filtrate pool as a 1:10 dilution in extraction buffer using a turbidimeter or similar device.

3.6) Further clarify the bag filtrate pool using 0.02  $\text{m}^2$  of a K700 (top) and KS50 (bottom) depth filter layer combination per liter of bag filtrate pool. Apply a flow rate of  $3.0$  L  $\text{min}^{-1}$   $\text{m}^{-2}$  up to a maximum pressure of 0.2 MPa. Then remove residual particles by passing the clarified filtrate through a 0.2- $\mu\text{m}$  sterile filter as previously described<sup>5</sup>.

#### **4. Purification of the DFE affinity ligand**

4.1) Prepare a fast protein liquid chromatography system by flushing with elution buffer (10 mM sodium phosphate, 300 mM imidazole, pH 7.4), wash buffer (10 mM sodium phosphate, pH 7.4) and equilibration buffer (20 mM sodium phosphate, 500 mM sodium chloride, pH 7.4). Mount a column containing  $\sim 50$  mL of chelating cross-linked agarose resin per kilogram of leaf biomass ( $\sim 2.5$  L depth filtrate).

4.2) Charge the column with nickel ions by applying 5 column volumes (CV) of 200 mM nickel sulfate solution and wash it with 5 CV of elution buffer. Use a flow rate of  $50$   $\text{cm h}^{-1}$  and monitor the UV signals at 260, 280 and 558 nm for all subsequent chromatography steps.

4.3) Equilibrate the column with 10 CV of equilibration buffer. Then load 50 CV of the clarified plant extract (from step 3.6) onto the conditioned column.

4.4) Wash the column with 10 CV of wash buffer. Elute the DFE fusion protein with 5 CV elution buffer and collect the product-containing fraction once the UV signals at 280 nm and 558 nm have increased to more than 5 mAU above the baseline.

4.5) Take a 0.2 mL sample from the elution fraction in order to measure the total soluble protein (TSP) concentration (step 7.1), DFE concentration (steps 7.2 and 7.3), and DFE purity (step 7.4).

4.6) Mount a column containing ~50 mL of cross-linked dextran resin on a chromatography system. Equilibrate the column with 5 CV of coupling buffer (200 mM HEPES, 500 mM sodium chloride, pH 8.5). Inject 10 mL of the DFE elution fraction (step 4.4) for buffer exchange and monitor the UV absorbance at 280 and 558 nm.

4.7) Collect the DFE-containing fraction once the UV signals at 280 nm and 558 nm have increased to more than 5 mAU above the baseline. Take a 0.2 mL sample and determine the TSP concentration, DFE concentration, and DFE purity (step 7).

4.8) Concentrate the purified DFE sample (from step 4.7) to 15 g L<sup>-1</sup> using a centrifugal concentrator tube at 3,000 × *g* at 4 °C for 30 min in a centrifuge. Continue with the coupling reaction (step 5).

NOTE: Store the DFE solution at 4 °C if the concentration or coupling steps cannot be carried out immediately.

## **5. Coupling DFE to the activated cross-linked agarose resin**

NOTE: Do not replace the isopropanol used for the storage of NHS-activated columns until all equipment and solutions for coupling are ready. Never let the columns run dry.

5.1) Setup a design of experiments (DoE) model to optimize the coupling of DFE to the activated resin, with pH, buffer composition and DFE concentration as factors. The details of the DoE method are described elsewhere<sup>15</sup>.

5.2) Prepare the affinity ligand solution (from step 4.8) in a concentration range from 0.5–15 g L<sup>-1</sup> or as defined in the DoE and store it on ice until the coupling reaction is ready (step 5.5). Fill at least ten 2 mL syringes with the DFE solution and prepare an adapter to mount the syringes on NHS-activated cross-linked agarose columns with a bed volume of 1.0 mL.

5.3) For every 10 columns used for coupling, prepare 30 mL of deactivation solution (0.5 M ethanolamine, 0.5 M sodium chloride, pH 8.3), 30 mL of low-pH solution (0.1 M sodium acetate, 0.5 M sodium chloride, pH 4.0) and 10 mL of storage solution (0.05 M disodium phosphate, 0.1% (m/v) sodium azide, pH 7.0).

5.4) Prepare 20 mL of 1 mM hydrochloric acid in a tube and incubate it on ice for at least 20 min. Prepare a precision scale to monitor the flow-through fractions for all steps during the coupling reaction (step 5.7).

5.5) Open a sealed NHS-activated cross-linked agarose column and mount the syringe adapter at the column inlet. Prevent any air entering the column by applying a drop of buffer onto the adapter inlet before connecting it to the syringe.

5.6) Wash the column with 6 mL of ice-cold 1 mM hydrochloric acid (from step 5.4) at a flow rate of  $<1 \text{ mL min}^{-1}$  and immediately proceed to step 5.7.

5.7) Inject 1.5 mL of DFE solution (step 5.2) using a 2 mL syringe at a flow rate of  $<1 \text{ mL min}^{-1}$  and collect the flow-through fraction on a precision scale (step 5.4) for subsequent analysis (step 7). Seal the column at both ends and incubate for 15–45 min at  $22^\circ\text{C}$ , depending on the DoE setup.

5.8) Inject 6 mL of deactivation solution followed by 6 mL of low-pH solution at a flow rate of  $<1 \text{ mL min}^{-1}$  to remove non-covalently bound ligands from the resin. Then inject 6 mL of deactivation solution and incubate the column for 15 min.

5.9) Inject 6 mL of low-pH solution into the column, followed by 6 mL of deactivation solution. Then inject another 6 mL of low-pH solution into the column.

5.10) Inject 2 mL of storage solution into the column and store at  $4^\circ\text{C}$ .

## **6. Testing the purification of mAbs from clarified plants extracts**

6.1) Prepare 100 mL of clarified plant extract containing 2F5<sup>5</sup> or the supernatant from the preferred cell-based expression system, also containing 2F5.

6.2) Prepare equilibration buffer (20 mM sodium phosphate, 500 mM sodium chloride, pH 7.4), low-pH elution buffer (0.05 M citrate, 0.05 M sodium chloride, pH 4.0–3.25), and high-ionic-strength elution buffer (1.0–4.0 M magnesium chloride, 0.1 M HEPES, pH 8.0)

6.3) Flush the chromatography system with the buffers. Mount a DFE affinity column (from step 5.10) on the chromatography system and equilibrate with 5 CV of equilibration buffer at a flow rate of  $1.0 \text{ mL min}^{-1}$ . Monitor the UV absorbance at 280 nm.

**NOTE:** Loading plant extract or cell culture supernatant onto the column can cause an increase in backpressure. Set a high-pressure alert at 0.2 MPa to avoid damage to the chromatography system or DFE column.

6.4) Load 80 mL of the clarified plant extract or supernatant (step 6.1) onto the column at a flow rate of  $0.5 \text{ mL min}^{-1}$  to guarantee a contact time of 2 min. Collect the flow-through samples in 2 mL fractions for breakthrough curve reconstruction (step 7.3). Store the flow-through samples at  $4^\circ\text{C}$  if immediate sample analysis is not possible.

6.5) Wash the column with 6 CV of equilibration buffer. Collect a sample of the wash at the



beginning, middle and end of this step.

6.6) Elute mAb 2F5 with 5 CV of low-pH elution buffer or high-ionic-strength elution buffer (0.1 M HEPES, 1.25 M magnesium chloride, pH 8.0). Collect the DFE fraction when the UV 280 nm signal has increased to 5 mAU above the baseline.

6.6.1) Optimize the elution buffer for each epitope–antibody pair. For 2F5, 1.25 M magnesium chloride achieved an optimal balance between product recovery and ligand stability.

NOTE: The magnesium chloride solution is prone to precipitation. Therefore, dissolve the magnesium chloride in ~700 mL of water. Separately dissolve the HEPES in 100 mL of water and adjust the pH to 8.0. Add the dissolved magnesium chloride solution to the HEPES solution and add water to a final volume of 1.0 L. Do not adjust the pH after dissolving the magnesium chloride because this will cause precipitation.

6.7) Analyze all samples taken during steps 6.4–6.6 using the Bradford method, lithium dodecylsulfate polyacrylamide gel electrophoresis (LDS-PAGE) and surface plasmon resonance (SPR) spectroscopy (step 7).

## 7. Sample analysis

7.1. Measure the TSP concentration using the Bradford method<sup>18,19</sup>.

7.1.1. In triplicate, pipette 2.5  $\mu$ L of each sample into a single well of a 96-well plate. Use eight bovine serum albumin (BSA) standards in triplicate covering the range 0–2,000 mg L<sup>-1</sup>.

7.1.2. Add 200  $\mu$ L of Bradford reagent to each well and mix by gently pipetting up and down. Keep the pipette level to avoid the formation of bubbles that distort the subsequent readout.

7.1.3. Incubate the plate for 10 min at 22 °C and measure the absorbance at 595 nm in a spectrophotometer. Calculate the TSP concentration in the samples based on a standard curve through the BSA reference points.

7.2. Quantify DFE by fluorimetry

7.2.1. In triplicate, pipette 50  $\mu$ L of each sample into single wells of a black 96-well half-area plate. Include six DsRed standards covering the range 0–225 mg L<sup>-1</sup>.

7.2.2. Measure the fluorescence twice using a 530  $\pm$  30 nm excitation filter and a 590  $\pm$  35 nm emission filter in a spectrophotometer. Calculate the DFE concentration in the samples based on a standard curve through the DsRed reference points.

7.3. Measure the 2F5 concentration by SPR spectroscopy<sup>20</sup>.

7.3.1. Prepare the HBS-EP running buffer (10 mM 4-(2-hydroxyethyl)-1-piperazineethanesulfonic acid (HEPES), 3 mM ethylenediaminetetraacetic acid (EDTA), 150 mM sodium chloride, 0.005% v/v Tween-20, pH 7.4) and filter sterilize it by passing it through a 0.2  $\mu\text{m}$  vacuum bottle-top filter. Degas the buffer for 15 min.

7.3.2. Connect the HBS-EP running buffer to the SPR instrument before docking a carboxymethylated dextran surface chip and equilibrate the system using the prime function. Start a manual run with a flow rate of 30  $\mu\text{L min}^{-1}$  and condition the chip surface by alternately injecting 30 mM hydrochloric acid and 25 mM sodium hydroxide (two injections each) over flow cells 1 and 2 at a flow rate of 30  $\mu\text{L min}^{-1}$  for 1 min.

7.3.3. Prepare 300  $\mu\text{L}$  of a 500  $\text{mg L}^{-1}$  Protein A solution in 10 mM sodium acetate (pH 4.0). Thaw vials containing 0.4 M 1-ethyl-3-(3-dimethylaminopropyl) carbodiimide hydrochloride (EDC) and 0.1 M NHS and centrifuge at 16,000  $\times g$  for 1 min before mixing 70  $\mu\text{L}$  of EDC and 70  $\mu\text{L}$  of NHS.

7.3.4. Transfer the EDC/NHS mixture and the Protein A solution to 7-mm plastic sample vials and place in the sample rack. Activate the carboxymethylated dextran chip surface by injecting the EDC/NHS mixture over flow cells 1 and 2 at a flow rate of 10  $\mu\text{L min}^{-1}$  for 10 min.

7.3.5. Couple Protein A to the activated carboxymethylated dextran surface by injecting the Protein A solution over flow cell 2 at a flow rate of 15  $\mu\text{L min}^{-1}$  for 15 min.

7.3.6. While coupling the Protein A, thaw a vial of 1.0 M ethanolamine hydrochloride (pH 8.5) and centrifuge at 16,000  $\times g$  for 1 min. Transfer 150  $\mu\text{L}$  to a 7-mm plastic vial and place into the instrument. When step 7.3.5 is complete, deactivate the chip surface by injecting the ethanolamine solution over flow cells 1 and 2 at a flow rate of 10  $\mu\text{L min}^{-1}$  for 7 min.

7.3.7. Condition the chip surface by alternately injecting 30 mM hydrochloric acid and 25 mM sodium hydroxide (two injections each) over flow cells 1 and 2 at a flow rate of 30  $\mu\text{L min}^{-1}$  for 0.5 min.

7.3.1) Prepare standards of antibody 2F5 in HBS-EP running buffer at a concentration of 500  $\mu\text{g L}^{-1}$  and dilute samples containing 2F5 in HBS-EP to a final concentration in the range 20–1,000  $\mu\text{g L}^{-1}$ . Inject samples and standards over flow cells 1 and 2 at a flow rate of 30  $\mu\text{L min}^{-1}$  for 3 min. Inject a 2F5 standard after every 15 samples.

7.3.2) Subtract the relative units (RU) signal of flow cell 1 (reference cell) from the signal of flow cell 2 (measurement cell) for each sample and calculate the antibody concentration based on the RU signal of 2F5 standard injections.

7.3.3) After every sample or standard injection, regenerate the chip surface by injecting 30 mM hydrochloric acid at a flow rate of 30  $\mu\text{L min}^{-1}$  for 0.5 min over flow cells 1 and 2.

7.3.4) Plot the 2F5 concentration for each flow-through sample from step 6.4 against the flow

through volume to obtain the 2F5 breakthrough curve. Calculate the amount of injected 2F5 when 10% of the loading 2F5 concentration was reached to obtain the 10% dynamic binding capacity (DBC).

#### 7.4. Analyze protein samples by LDS-PAGE.

7.4.1) Open a ready-to-use 4–12% BisTris LDS polyacrylamide gel and place it in an electrophoresis module. Transfer 800 mL of running buffer (50 mM MES, 50 mM Tris base, 0.1% (m/v) SDS, 1 mM EDTA, pH 7.3) into the module and add 0.5 mL of antioxidant solution.

7.4.1. In a 1.5-mL reaction tube, mix 10  $\mu$ L of loading buffer with 4  $\mu$ L of reducing agent and 26  $\mu$ L of the protein sample. Incubate the tube in a heat block for 10 min at 80 °C. Centrifuge the tube at  $500 \times g$  for 30 s.

7.4.2. Load the LDS polyacrylamide gel (10  $\mu$ L of sample per lane) and load 5  $\mu$ L of pre-stained protein standard (10–180 kDa) in a separate lane.

7.4.3. Close the electrophoresis chamber and connect the power supply. The electrophoresis should run for 40 min and at constant 200 V.

7.4.4. Remove the gel from the electrophoresis chamber. Open the gel sheathing and wash the gel for 15 min in water on a shaker at 19 rpm. Stain the gel for 1 h in staining solution on a shaker at 19 rpm.

7.4.5. Destain the gel for 1 h in water on a shaker at 19 rpm. Scan the gel using a film scanner.

### REPRESENTATIVE RESULTS:

#### *Expression and purification of the affinity ligand*

The fusion protein DFE was expressed in transgenic tobacco plants grown in a greenhouse. The yield was  $\sim 120 \text{ mg kg}^{-1}$  leaf mass with an average biomass of  $\sim 130 \text{ g}$  per plant. The DFE purity was  $<5\%$  of TSP in crude plant extracts before blanching but increased to  $\sim 40\%$  after heat treatment at 70 °C for 1.5 min, which precipitated  $>97\%$  of the host cell proteins. The blanching step was easily integrated into the harvesting and extraction routine (**Figure 1**) and took less than 2 h of extra time, including setting up the water bath. The overall recovery of DFE was  $23.5 \text{ mg kg}^{-1}$  with a purity of  $>90\%$ . The steps responsible for product loss were blanching, depth filtration and IMAC, with specific losses of 40%, 27% and 45%, respectively. The depth filter capacity was on average  $135 \pm 36 \text{ L m}^{-2}$  ( $\pm$ SD,  $n=3$ ) and thus in the upper range of values reported in the literature<sup>21</sup>. The DFE yield increased with plant age (**Figure 2**).

#### *Immobilization of the affinity ligand on NHS-activated chromatography columns*

During initial coupling tests, we found that HEPES buffer (pH 8.3) increased the coupling efficiency to  $89 \pm 6\%$  ( $\pm$ SD,  $n=3$ ) compared to  $78 \pm 9\%$  ( $\pm$ SD,  $n=3$ ) for the bicarbonate buffer recommended by the manufacturer. Therefore, HEPES was used for all subsequent coupling experiments. A DoE approach was selected to optimize the coupling efficiency of DFE to NHS-

activated cross-linked agarose resin. The absolute amount of DFE immobilized on the resin increased with the mass of DFE injected into the column and plateaued at  $\sim 15 \text{ g L}^{-1}$  whereas the coupling yield declined continuously as more DFE was injected (**Figure 2**). The coupling yield was also  $>50\%$  lower in an acidic buffer, indicating the need to screen for suitable coupling conditions for each ligand on a case-by-case basis. Ideal conditions in terms of coupling yield, absolute quantity of immobilized DFE and column costs were identified using the numerical optimization tool of the DoE software. The most desirable conditions (pH 9.0 and 7.0 mg of DFE per 1 mL of cross-linked agarose resin) were located on a large plateau and were therefore robust. The DFE molecules retained their red fluorescence even after coupling, and the color intensity corresponded to the total amount of immobilized DFE (**Figure 2**). Therefore, column color can be used as a simple quality control parameter to estimate the coupling efficiency and column quality. The fluorescence also confirmed that DFE fusion protein assembled in the tetrameric state of native DsRed.

[Place Figure 2 here]

#### *Testing 2F5 isolation using the DFE affinity resin*

The recombinant 2F5 antibody was transiently produced in *Nicotiana benthamiana* plants grown in a phytotron<sup>5</sup>. The capture of 2F5 from the crude plant extract was tested using affinity columns coupled with  $\sim 7.0 \text{ mg}$  purified DFE (step 6). Elution from Protein A resins usually involves an acidic buffer (pH  $\sim 3.3$ )<sup>13</sup>. Therefore, we initially evaluated different low-pH elution buffers (pH 6.0–3.25) for the DFE columns. The elution of 2F5 was successful at pH values below 4.5 with the highest recovery of  $\sim 35\%$  at pH 3.25. However, low-pH elution inactivated both the antibody (as confirmed by SPR spectrometry) and the DFE ligand (as indicated by the loss of color and the lower DBC, Figure 3). The latter was anticipated given that native DsRed denatures at pH  $< 4.0$ <sup>22,23</sup>. To avoid product and ligand denaturation, we tested magnesium chloride as an alternative elution agent because it has previously been used to elute mAbs from other affinity resins<sup>24</sup>. A magnesium chloride concentration of 1.25 M was sufficient to elute 2F5 from the DFE affinity resin with a recovery of  $105 \pm 11\%$  ( $\pm \text{SD}$ ,  $n=3$ ) and a purity of  $97 \pm 3\%$  ( $\pm \text{SD}$ ,  $n=3$ ). This performance was comparable to Protein A resins<sup>25,26</sup>. The equilibrium dissociation constant ( $K_D$ ) of DFE-purified 2F5 antibody and the synthetic ligand Fuzeon was 791 pM whereas that of a Protein A-purified counterpart was 763 pM<sup>5</sup>. Furthermore, no substantial color loss was observed in the resin over a total of 25 bind-and-elute cycles. The DBC of the DFE affinity resin at 10% 2F5 breakthrough declined linearly over the course of 25 cycles to  $\sim 15\%$  of the initial value (**Figure 3**).

[Place Figure 3 here]

#### **FIGURE LEGENDS:**

**Figure 1: Process flow scheme for the preparation of epitope affinity resins that can be used for the capture of mAbs from crude plant extracts or cell culture supernatants. (A)** The affinity ligand DFE was expressed in transgenic tobacco plants. A heat precipitation step **(B)** was included before harvested leaves were homogenized **(C)**. **(D)** The crude plant extract was clarified by bag filtration, depth filtration and  $0.2\text{-}\mu\text{m}$  sterile filtration. **(E)** DFE was then purified by IMAC. **(F, G)**

The purified DFE affinity ligand was immobilized on EDC/NHS-activated cross-linked agarose columns. (H) Bacterial cultures carrying T-DNA encoding antibody 2F5 were used for transient expression in *N. benthamiana* plants (I) grown in a phytotron. (J) *N. benthamiana* leaves were harvested and processed as described in D. (K) mAb 2F5 was purified from the clarified extract using the DFE affinity columns (L).

**Figure 2: Optimization of DFE immobilization to NHS-activated cross-linked agarose resin.**

(A) LDS-PAA gel with western blot overlay of homogenate and elution samples from unblanched and blanched transgenic DFE plant extracts. Harvest of plants was performed 38, 45 or 52 days after seeding. Western blots were performed using an anti-His<sub>6</sub>-antibodies<sup>5</sup>. (B) Total amount of coupled DFE affinity ligand in dependence of coupling pH and overall mass of purified DFE injected onto NHS-activated cross-linked agarose columns. Red dots indicate the actual experiments performed to build the response surface model. (C) DFE affinity columns after the coupling procedure. The numbers correspond to the coupling conditions highlighted in panel B. dps – days post seeding.

**Figure 3: Testing the isolation of 2F5 from clarified plant extracts using the DFE affinity resins.**

(A) DFE affinity resins after one elution cycle using buffers with different pH values in the range 4.0–3.0 and the same resins after a neutralization step at pH 6.0. (B) DFE affinity resins after one and six cycles of 2F5 purification using 1.25 M or 4.00 M magnesium chloride as eluent. (C) Chromatograms of frontal loading experiments (break-through curves) to determine the cycle-dependent dynamic binding capacity of DFE resin using magnesium chloride as eluent. The break-through curves were measured for 4.0 M and 1.25 M magnesium chloride elution buffers and various cycle numbers. (D) Cycle-dependent dynamic binding capacity of DFE using 1.25 M magnesium chloride as eluent.

**Table 1: Trouble-shooting**

**DISCUSSION:**

*Applications of the novel affinity resin*

We have shown that custom affinity chromatography resins for the capture of mAbs can be manufactured by immobilizing a ligand containing a mAb-specific epitope to NHS-activated cross-linked agarose. To design such a resin, it was necessary to know the epitope sequence and to use a linear epitope. The resulting resins are advantageous for the capture of mAbs because they could potentially replace expensive Protein A affinity chromatography steps. The interaction between 2F5 and DFE in our case study was mediated by epitope–paratope binding, so our ligand should be more selective than Protein A, which binds to the Fc region of most murine and human IgGs. Because individual ligands are needed for each mAb, our method may initially seem suitable mainly for antibodies that are produced on a very large scale. However, by combining our approach with rapid plant-based transient protein expression, new affinity ligands can be prepared in less than 2 weeks<sup>27</sup> with minimal effort<sup>28</sup>. Hence, the method is also suitable for small-scale mAb purification.

*Production and potential improvements of the affinity ligand*

Plants offer a fast and safe production platform for affinity ligands<sup>5,29,30</sup>, such as the DFE fusion protein featured in our case study. Blanching the plant material greatly reduced the quantity of host cell proteins in a single step and was easily integrated into a standard clarification routine. However, the recovery of the ligand was low in the current setup, probably due to its moderate thermal stability and some non-specific binding to the filter layers, as reported for other products<sup>31-33</sup>. Engineering the carrier to increase its thermal stability may therefore help to improve the ligand yield in the future, as described for the malaria vaccine candidate CCT, the antitumor enzyme PpADI or a mesophilic  $\beta$ -glucosidase<sup>34-36</sup>. The same holds true for the depth filtration step, where protein engineering may help to reduce non-specific binding to the filter material<sup>37</sup>. The production costs for DFE and similar ligands could also be reduced by improving the overall efficiency of clarification using flocculants or filter additives<sup>38,39</sup>.

When DsRed is used as a carrier, it forms a tetrameric complex. This is advantageous because it increases the number of epitopes per ligand, but it may also render the ligand more susceptible to disassembly or denaturation during affinity chromatography. A monomeric carrier protein such as mCherry may therefore be preferable, because it is stable at low pH<sup>40</sup>, and the inclusion of tandem repeats of the epitope would increase the avidity of the ligand and thus increase resin capacity<sup>5,26,41</sup>. For simple carrier-epitope proteins (i.e., those with no disulfide bonds or post-translational modifications) microbial production systems may reduce the manufacturing costs and make the ligands more competitive with Protein A. For example, green fluorescent protein has been expressed in bacterial cells with a yield of  $\sim 1 \text{ g kg}^{-1}$  biomass, which would significantly reduce ligand production costs<sup>42</sup>.

Regardless of the expression host, a purified affinity ligand was required during coupling to minimize the immobilization of host cell proteins or media components that can otherwise reduce resin selectivity and capacity. The inclusion of a poly-histidine tag for IMAC purification increased the purity to  $\sim 90\%$  in a single step, facilitating rapid and inexpensive ligand production<sup>5,43,44</sup>. However, the position of the fusion tag is important because it has the potential to sterically hinder epitope binding or to induce the cleavage of either the tag or the epitope from the carrier<sup>45,46</sup>.

#### *Immobilization of the affinity ligand on NHS-activated chromatography columns*

Immobilization was carried out manually or using a chromatography system. The small buffer volumes per column seemed to favor manual handling (e.g., due to the minimal waste volumes). However, if multiple/larger columns are needed, the chromatography system makes the coupling conditions easier to control (e.g., regulated flow rates) and is therefore more likely to achieve reproducible results in terms of DBC. Our data suggest that the coupling buffer and pH have an important effect on the coupling efficiency and overall column costs. Screening factors that influence the coupling reaction and adjusting them for each carrier protein (or even for each carrier–ligand fusion) could therefore improve coupling efficiency and resin performance, and we recommend this approach.

#### *Testing 2F5 isolation using the DFE affinity resin*

Product yield and purity are important aspects of resin performance, and in the case of DFE we

achieved a yield of  $105 \pm 11\%$  ( $\pm$ SD,  $n=3$ ) and a purity of  $97 \pm 3\%$  ( $\pm$ SD,  $n=3$ ), which is comparable with the performance of benchmark Protein A resins<sup>25,26</sup>. Another key performance indicator for resins in general (and particularly for those based on affinity ligands) is the DBC at 10% product breakthrough, because this parameter affects the amount of resin required for a specific process and thus the costs. For the DFE ligand, the initial DBC was  $\sim 4 \text{ g L}^{-1}$  resin, which is  $\sim 13\%$  of the corresponding value for Protein A under similar conditions (only 2 min contact time)<sup>25,47</sup> but about 15-fold higher compared to other custom affinity resins such as the anti-FSH-immunoaffinity ligand using the same residence time of 2 minutes<sup>48</sup>. The DBC of DFE declined to 15% of the starting value after 25 bind-and-elute cycles, whereas more than 50 cycles are required for the same loss of DBC in commercial Protein A resins<sup>49</sup>. However, it is important to note that our carrier has not yet been optimized to the same extent as Protein A, which has been comprehensively investigated and improved over the last four decades<sup>8</sup>.

Thus far we have improved the resin stability and product recovery by switching from a low-pH elution buffer to a high concentration of magnesium chloride (**Figure 3**), as recommended in earlier studies<sup>13</sup>. The characteristic red color of the affinity ligand did not fade substantially during the 25 bind-and-elute cycles, so we speculate that endogenous plant proteases in the clarified plant extracts<sup>31</sup> may have truncated and thus inactivated the epitope of the ligand. Therefore, designing protease-resistant linkers to connect the carrier and epitope may help to maintain the initial DBC over an extended number of cycles. Given the rapid and simple expression and purification of the DFE ligand, its straightforward coupling to commercial chromatography resins, and its excellent product yield and purity, we believe that our method offers a suitable alternative to Protein A for the purification of mAbs and antibody derivatives which do not bind to Protein A, especially if improvements to the carrier and linker can improve the DBC and ligand stability. This assumption was supported by the small difference in the dissociation constant of DFE-purified and Protein A-purified 2F5 antibody<sup>5</sup>, indicating that our new affinity ligand allows the recovery of high-quality mAbs.

#### *Benefits and current limitations of the method*

Producing the affinity ligand as a genetic fusion with a carrier protein increases solubility in aqueous buffers and thus compatibility with typical ligand coupling conditions. In contrast, blank peptides derived from solid phase peptide synthesis may have limited solubility under these conditions due to their sequence<sup>50</sup>, which cannot be changed because it is dictated by the amino acid sequence of epitope recognized by the mAb to be purified. Others have therefore used an on-resin synthesis of peptide ligands<sup>51</sup>. The static binding capacity of the resulting resin was high ( $\sim 80 \text{ g L}^{-1}$ ), but the process of resin preparation is lengthy, a dynamic binding capacity was not reported and the obtained purity and recovery were lower than in our approach. An additional advantage of a fusion protein ligand in laboratory scale is that the ligand and variants thereof can be rapidly produced, purified and tested with minimal effort in an easy-to-use high-through expression system<sup>52</sup>.

The two current limitations of the method presented here are the low dynamic binding capacity of  $3 \text{ g L}^{-1}$  and its 90% reduction over the course of 25 bind-and-elute cycles<sup>5</sup>. These limitations can be addressed in the future by applying less stringent loading conditions and replacing the

current DsRed carrier with an engineered, more stable variant respectively. For example, doubling the current contact time from 2 to 4 minutes has the potential to double the dynamic binding capacity as was shown for some Protein A resins<sup>26</sup>.

## Troubleshooting

The following table highlights potential problems that can be encountered during this protocol and provides hints on how to solve them (**Table 1**).

[place table 1 here]

## ACKNOWLEDGMENTS:

We would like to acknowledge Ibrahim Al Amedi for cultivating the transgenic tobacco plants and Dr. Thomas Rademacher for providing the tobacco expression vector. The authors wish to thank Dr. Richard M. Twyman for editorial assistance and Markus Sack for fruitful discussions about the DFE affinity ligand structure. This work was funded in part by the Fraunhofer-Gesellschaft Internal Programs under Grant No. Attract 125-600164 and the state of North-Rhine-Westphalia under the Leistungszentrum grant no. 423 “Networked, adaptive production”. This work was supported by the Deutsche Forschungsgemeinschaft (DFG) in the framework of the Research Training Group “Tumor-targeted Drug Delivery” grant 331065168. GE healthcare supported the open-access publication of this article.

## DISCLOSURES:

The authors have no conflicts of interest to disclose.

## REFERENCES:

- 1 Kesik-Brodacka, M. Progress in biopharmaceutical development. *Biotechnology and Applied Biochemistry*. **65** (3), 306-322, doi:10.1002/bab.1617 (2018).
- 2 Ecker, D.M., Jones, S.D., Levine, H.L. The therapeutic monoclonal antibody market. *MAbs*. **7** (1), 9-14, doi:10.4161/19420862.2015.989042 (2015).
- 3 Jayapal, K., Wlaschin, K. F., Hu, W. S., Yap, M. G. S. Recombinant Protein Therapeutics from CHO Cells - 20 Years and Counting. *Chemical Engineering Progress*. **103**, 40-47 (2007).
- 4 Kunert, R., Reinhart, D. Advances in recombinant antibody manufacturing. *Applied Microbiology and Biotechnology*. **100** (8), 3451-3461, doi:10.1007/s00253-016-7388-9 (2016).
- 5 Rühl, C., Knödler, M., Opdensteinen, P., Buyel, J. F. A linear epitope coupled to DsRed provides an affinity ligand for the capture of monoclonal antibodies. *Journal of Chromatography A*. **1571**, 55-64, doi:10.1016/j.chroma.2018.08.014 (2018).
- 6 Edgue, G., Twyman, R. M., Beiss, V., Fischer, R., Sack, M. Antibodies from plants for bionanomaterials. *Nanomedicine and Nanobiotechnology*. **9** (6), e1462, doi:10.1002/wnan.1462 (2017).
- 7 Buyel, J.F., Fischer, R. Very-large-scale production of antibodies in plants: The biologization of manufacturing. *Biotechnology Advances*. **35** (4), 458-465, doi:<https://doi.org/10.1016/j.biotechadv.2017.03.011> (2017).
- 8 Bolton, G.R., Mehta, K.K. The role of more than 40 years of improvement in protein A chromatography in the growth of the therapeutic antibody industry. *Biotechnology Progress*. **32**



661 (5), 1193-1202, doi:10.1002/btpr.2324 (2016).  
 662 9 Kelley, B. Industrialization of mAb production technology: the bioprocessing industry at a  
 663 crossroads. *MAbs*. **1** (5), 443-452 (2009).  
 664 10 Brochier, V., Ravault, V. High throughput development of a non protein A monoclonal  
 665 antibody purification process using mini-columns and bio-layer interferometry. *Engineering in*  
 666 *Life Sciences*. **16** (2), 152-159, doi:doi:10.1002/elsc.201400244 (2016).  
 667 11 Arakawa, T., Futatsumori-Sugai, M., Tsumoto, K., Kita, Y., Sato, H., Ejima, D. MEP HyperCel  
 668 chromatography II: binding, washing and elution. *Protein Expression and Purification*. **71** (2), 168-  
 669 173, doi:10.1016/j.pep.2009.11.004 (2010).  
 670 12 Barroso, T.B., R. J. Aguiar-Ricardo, A. Roque, A. C. Structural evaluation of an alternative  
 671 Protein A biomimetic ligand for antibody purification. *Journal of Computer-aided Molecular*  
 672 *Design*. **28** (1), 25-34, doi:10.1007/s10822-013-9703-1 (2014).  
 673 13 Mazzer, A.R., Perraud, X., Halley, J., O'Hara, J., Bracewell, D.G. Protein A chromatography  
 674 increases monoclonal antibody aggregation rate during subsequent low pH virus inactivation  
 675 hold. *Journal of Chromatography A*. **1415**, 83-90, doi:10.1016/j.chroma.2015.08.068 (2015).  
 676 14 Sack, M. et al. Functional analysis of the broadly neutralizing human anti-HIV-1 antibody  
 677 2F5 produced in transgenic BY-2 suspension cultures. *FASEB Journal*. **21** (8), 1655-1664,  
 678 doi:10.1096/fj.06-5863com (2007).  
 679 15 Buyel, J.F., Fischer, R. Characterization of Complex Systems Using the Design of  
 680 Experiments Approach: Transient Protein Expression in Tobacco as a Case Study. *Journal of*  
 681 *Visualized Experiments*. (83), 51216, doi:10.3791/51216 (2014).  
 682 16 Trasatti, J.P., Woo, J., Ladiwala, A., Cramer, S., Karande, P. Rational design of peptide  
 683 affinity ligands for the purification of therapeutic enzymes. *Biotechnology Progress*. **34** (4), 987-  
 684 998, doi:10.1002/btpr.2637 (2018).  
 685 17 Buyel, J.F., Gruchow, H. M. Boes, A., Fischer, R. Rational design of a host cell protein heat  
 686 precipitation step simplifies the subsequent purification of recombinant proteins from tobacco.  
 687 *Biochemical Engineering Journal*. **88**, 162-170, doi:<https://doi.org/10.1016/j.bej.2014.04.015>  
 688 (2014).  
 689 18 Simonian, M.H., Smith, J.A. Spectrophotometric and colorimetric determination of  
 690 protein concentration. *Current Protocols in Molecular Biology*. **76** (Chapter 10), 10.11.11-  
 691 10.11A.19., doi:10.1002/0471142727.mb1001as76 (2006).  
 692 19 Buyel, F.J., Kaefer, T., Buyel, J., Fischer, R. Predictive models for the accumulation of a  
 693 fluorescent marker protein in tobacco leaves according to the promoter/5'UTR combination.  
 694 *Biotchnology and Bioengineering*. **110** (2), 471-483, doi:10.1002/bit.24715 (2013).  
 695 20 Piliarik, M., Vaisocherova, H., Homola, J. Surface plasmon resonance biosensing. *Methods*  
 696 *in Molecular Biology*. **503**, 65-88, doi:10.1007/978-1-60327-567-5\_5 [doi] (2009).  
 697 21 Buyel, J.F., Fischer, R. Scale-down models to optimize a filter train for the downstream  
 698 purification of recombinant pharmaceutical proteins produced in tobacco leaves. *Biotechnology*  
 699 *Journal*. **9** (3), 415-425, doi:10.1002/biot.201300369 (2014).  
 700 22 Baird, G.S., Zacharias, D.A., Tsien, R.Y. Biochemistry, mutagenesis, and oligomerization of  
 701 DsRed, a red fluorescent protein from coral. *Proceedings of the National Academy of Sciences,*  
 702 *USA*. **97** (22), 11984-11989 (2000).  
 703 23 Vrzheschch, P.V., Akovbian, N.A., Varfolomeyev, S.D., Verkhusha, V.V. Denaturation and  
 704 partial renaturation of a tightly tetramerized DsRed protein under mildly acidic conditions. *FEBS*

705 *Letters*. **487** (2), 203-208, doi:10.1016/S0014-5793(00)02344-9 (2000).

706 24 Firer, M.A. Efficient elution of functional proteins in affinity chromatography. *Journal of*  
707 *Biochemical and Biophysical Methods*. **49** (1-3), 433-442, doi:10.1016/S0165-022X(01)00211-1  
708 (2001).

709 25 Pabst, T.M. et al. Engineering of novel Staphylococcal Protein A ligands to enable milder  
710 elution pH and high dynamic binding capacity. *Journal of Chromatography A*. **1362**, 180-185,  
711 doi:10.1016/j.chroma.2014.08.046 (2014).

712 26 Müller, E., Vajda, J. Routes to improve binding capacities of affinity resins demonstrated  
713 for Protein A chromatography. *Journal of Chromatography B*. **1021**, 159-168,  
714 doi:<https://doi.org/10.1016/j.ichromb.2016.01.036> (2016).

715 27 Shamloul, M., Trusa, Jason, Mett, Vadim, Yusibov, Vidadi. Optimization and Utilization of  
716 Agrobacterium-mediated Transient Protein Production in Nicotiana. *Journal of Visualized*  
717 *Experiments*. (86), 51204, doi:10.3791/51204 (2014).

718 28 Rademacher, T. et al. Plant cell packs: a scalable platform for recombinant protein  
719 production and metabolic engineering. *Plant Biotechnology Journal*. doi:10.1111/pbi.13081  
720 (2018).

721 29 Saxena, L., Iyer, B.K., Ananthanarayan, L. Purification of a bifunctional amylase/protease  
722 inhibitor from ragi (*Eleusine coracana*) by chromatography and its use as an affinity ligand.  
723 *Journal of Chromatography. B*. **878** (19), 1549-1554, doi:10.1016/j.jchromb.2010.04.009 (2010).

724 30 Kurppa, K., Reuter, L. J., Ritala, A., Linder, M. B., Joensuu, J. J. In-solution antibody  
725 harvesting with a plant-produced hydrophobin-Protein A fusion. *Plant Biotechnology Journal*. **16**  
726 (2), 404-414, doi:10.1111/pbi.12780 (2018).

727 31 Menzel, S. et al. Optimized Blanching Reduces the Host Cell Protein Content and  
728 Substantially Enhances the Recovery and Stability of Two Plant-Derived Malaria Vaccine  
729 Candidates. *Frontiers in Plant Science*. **7** (159), 1-7, doi:10.3389/fpls.2016.00159 (2016).

730 32 Yigzaw, Y., Piper, R., Tran, M., Shukla, A.A. Exploitation of the adsorptive properties of  
731 depth filters for host cell protein removal during monoclonal antibody purification. *Biotechnology*  
732 *Progress*. **22** (1), 288-296, doi:Doi 10.1021/Bp050274w (2006).

733 33 Menzel, S. et al. Downstream processing of a plant-derived malaria transmission-blocking  
734 vaccine candidate. *Protein Expression and Purification*. **152**, 122-130,  
735 doi:<https://doi.org/10.1016/j.pep.2018.07.012> (2018).

736 34 Vöpel, N., Boes, A., Edgü, G., Beiss, V. Malaria vaccine candidate antigen targeting the pre-  
737 erythrocytic stage of *Plasmodium falciparum* produced at high level in plants. *Biotechnology*  
738 *Journal*. **9** (11), 1435-1445, doi:10.1002/biot.201400350 (2014).

739 35 Lee, C.W., Wang, H. J., Hwang, J. K., Tseng, C. P. Protein thermal stability enhancement by  
740 designing salt bridges: a combined computational and experimental study. *PLoS ONE*. **9** (11),  
741 e112751, doi:10.1371/journal.pone.0112751 (2014).

742 36 Zhu, L., Cheng, F., Piatkowski, V., Schwaneberg, U. Protein engineering of the antitumor  
743 enzyme PpADI for improved thermal resistance. *Chembiochem*. **15** (2), 276-283,  
744 doi:10.1002/cbic.201300433 (2014).

745 37 Khanal, O. et al. Contributions of depth filter components to protein adsorption in  
746 bioprocessing. *Biotechnology and Bioengineering*. **115** (8), 1938-1948, doi:doi:10.1002/bit.26707  
747 (2018).

748 38 Buyel, J.F., Fischer, R. Downstream processing of biopharmaceutical proteins produced in

plants: the pros and cons of flocculants. *Bioengineered*. **5** (2), 138-142, doi:10.4161/bioe.28061 (2014).

39 Buyel, J.F. Procedure to Evaluate the Efficiency of Flocculants for the Removal of Dispersed Particles from Plant Extracts. *Journal of Visualized Experiments*. (110), e53940, doi:doi:10.3791/53940 (2016).

40 Fink, D. et al. Ubiquitous expression of the monomeric red fluorescent protein mCherry in transgenic mice. *Genesis*. **48** (12), 723-729, doi:10.1002/dvg.20677 (2010).

41 Gagnon, P. Technology trends in antibody purification. *Journal of Chromatography A*. **1221**, 57-70, doi:10.1016/j.chroma.2011.10.034 (2012).

42 Figueira, M., Laramée, L., Murrell, J. C., Groleau, D., Míguez, C. Production of green fluorescent protein by the methylotrophic bacterium *Methylobacterium extorquens*. *FEMS Microbiology Letters*. **193** (2), 195-200, doi:10.1016/S0378-1097(00)00463-8 (2001).

43 Bornhorst, J.A., Falke, J. J. Purification of proteins using polyhistidine affinity tags. *Methods in enzymology*. **326**, 245-254 (2000).

44 Sainsbury, F., Jutras, P. V., Vorster, J., Goulet, M., Michaud, D. A Chimeric Affinity Tag for Efficient Expression and Chromatographic Purification of Heterologous Proteins from Plants. *Frontiers in Plant Science*. **7**, 141-141, doi:10.3389/fpls.2016.00141 (2016).

45 Krupka, M. et al. The Position of His-Tag in Recombinant OspC and Application of Various Adjuvants Affects the Intensity and Quality of Specific Antibody Response after Immunization of Experimental Mice. *PLoS ONE*. **11** (2), e0148497, doi:10.1371/journal.pone.0148497 (2016).

46 Goel, A. et al. Relative position of the hexahistidine tag effects binding properties of a tumor-associated single-chain Fv construct. *Biochimica et Biophysica Acta*. **1523** (1), 13-20 (2000).

47 Tustian, A.D. et al. Development of a novel affinity chromatography resin for platform purification of bispecific antibodies with modified protein a binding avidity. *Biotechnology Progress*. doi:10.1002/btpr.2622 (2018).

48 Zandian, M., Jungbauer, A. An immunoaffinity column with a monoclonal antibody as ligand for human follicle stimulating hormone. *Journal of Separation Science*. **32** (10), 1585-1591, doi:10.1002/jssc.200900103 (2009).

49 Kelley, B. Very large scale monoclonal antibody purification: the case for conventional unit operations. *Biotechnology Progress*. **23** (5), 995-1008, doi:10.1021/bp070117s [doi] (2007).

50 Petrou, C., Sarigiannis, Y. In: *Peptide Applications in Biomedicine, Biotechnology and Bioengineering*. Koutsopoulos, S., ed., Ch. 1, Woodhead Publishing, Cambridge, 1-21 (2018).

51 Menegatti, S. et al. Design of protease-resistant peptide ligands for the purification of antibodies from human plasma. *Journal of Chromatography A*. **1445**, 93-104, doi:10.1016/j.chroma.2016.03.087 (2016).

52 Rademacher, T. et al. Plant cell packs: a scalable platform for recombinant protein production and metabolic engineering. *Plant Biotechnology Journal*. 1-7, doi:10.1111/pbi.13081 (2019).

Figure 1






Figure 2

[Click here to access/download;Figure;Figure2\\_DFE\\_coupling\\_optimization.png](#)




Figure 3

[Click here to access/download;Figure;Figure3\\_DFE\\_stability.png](#)




Table 1: Potential problems that can be encountered and possible fixes.

Protocol step	Problem
1	Plants do not grow
2 and 3	Large quantities of host cell proteins are present
2 and 3	No product found in the plant extract
3	Large stem or leaf parts remain after extraction
3	Rapid pressure increase during depth filtration
4	Little fusion protein during elution / a lot of fusion
4	Fusion protein lost during concentration
5	Low coupling yield
5 and 6	Low mAb yield
7	Very low/high protein concentrations in Bradford
7	Low mAb concentration during SPR measurement

Cause	Fix
Compromised growth conditions	Check the pH and conductivity of the fertilizer Check the temperature and light conditions
Incomplete precipitation	Check the temperature during blanching Check the agitation in the blanching bath
Blanching temperature too high	Check temperature and pH during blanching
pH in blanching buffer too low	
Incomplete mixing in blender	Make sure the plant material does not form a plug
Incorrect filter selection and/or orientation	Check the filter type and orientation
IMAC resin was not charged with metal ions	Check if the IMAC resin was correctly charged with metal ions
Fusion protein lost the affinity tag	Avoid intense sunlight and high temperatures during storage
Fusion protein bound to the membrane	Check the membrane type Make sure the concentration factor was not too high
Incorrect sequence of coupling reagent addition	Check the reagents labels and sequence of addition
Incorrect preparation of the columns before coupling	Check the conditions of column preparation
Low mAb expression in the plant biomass	Test mAb expression in biomass
Low ligand density	Check the purity of the fusion protein preparation
Bubble formation during pipetting	Check for bubbles in the 96-well plate
Compromised Protein A chip	Compare with results of standard mAb with known binding
Incorrect sample dilution	Check the dilution rate and buffer



g in the blender

th ions  
ing plant cultivation

igh  
ion

n concentration

Name of Material/ Equipment	Company	Catalog Number	Comments/Description
10L/20L Bucket	n/a	n/a	Blanching equipment
2100P Portable Turbidimeter	Hach	4650000	Turbidimeter
ÄKTApure	GE Helthcare	29018226	Chromatography system
Allegra 25R	Beckman Coulter	369434	Centrifuge
Amine Coupling Kit	GE Healthcare	BR100050	SPR chip coupling kit
Amine Coupling Kit	GE Healthcare	BR100050	SPR chip coupling kit
Antibody 2G12	Fraunhofer IME	n/a	Standard for SPR quantification
Blender	Waring	800EG	Blender
BP-410	Fuhr	2632410001	Bag filter
CanoScan 5600F	Canon	29258009	Scanner
Centrifuge tube 50 mL self-standing	Labomedic	1110504	Reaction tube
Chelating Sepharose FF	GE Helthcare	17-0575-01	Chromatography resin
Cond 3320	WTW	EKA 3338	Conductometer
Design-Expert(R) 8	Stat-Ease, Inc.	n/a	DoE software
Discovery Comfort	Gilson	F81029	Multichannel pipette
Sodium phosphate	Carl Roth GmbH	4984.3	Media component
Diverse bottles	Schott Duran	n/a	Glas bottles
Dri Block DB8	Techne	Z381373	Heat block
DsRed	Fraunhofer IME	n/a	Standart
EDTA	Carl Roth GmbH	8043.2	Buffer component
EnSpire	Perkin Elmer	2390-0000	Plate reader
ETHG-912	Oregon Scientific	086L001499-230	Thermometer
F9-C	GE Helthcare	29027743	Fraction collector
Ferty 2 Mega	Kammlott	5.220072	Fertilizer
Forma -86C ULT freezer	ThermoFisher	88400	Freezer
HEPES	Carl Roth GmbH	9105.3	Buffer component
Hettich Centrifuge Mikro 200	Hettich	2400	Centrifuge
HiPrep 26/10	GE Helthcare	GE17-5087-01	Chromatography column
HiTrap NHS-activated Sepharose HP, 1 mL	GE Helthcare	17-0716-01	Chromatography columns
Hydrochloric acid	Carl Roth GmbH	4625.1	Buffer component
Imidazole	Carl Roth GmbH	3899.2	Buffer component
K700	Pall	5302305	Depth filter layer
KM02 basic	IKA	n/a	Magnetic stirrer
KS50P 60D	Pall	B12486	Depth filter layer
L/S 24	Masterflex	SN-06508-24	Tubing
Lauda E300	Lauda Dr Wobser GmbH	Z90010	Immersion circulator
Magnesium chloride	Carl Roth GmbH	KK36.2	Buffer component
Masterflex L/S	Masterflex	HV-77921-75	Peristaltic pump
Minisart 0.2 µm	Sartorius	16534K	Filter unit
Nalgene Rapid-Flow PES bottle top filter	Thermo Fischer Scientific	595-4520	Vacuum filtration of SPR buffers
Nickel sulphate	Carl Roth GmbH	T111.1	Buffer component
Novex NuPAGE 4-12% BisTris LDS gels	Invitrogen	NP0336BOX	LDS-PAA gels
Novex X-cell Mini Cell	Invitrogen	EI0001	PAGE chamber
NuPAGE 20x running buffer	Invitrogen	NP0002	Buffer concentrate
NuPAGE antioxidant	Invitrogen	NP0005	Antioxidant
PageRuler protein ladder (10-180 kDa)	Invitrogen	26616	Protein standart
Perforated bucked	n/a	n/a	Blanching
PH 3110	WTW	2AA110	PH meter
PowerPac HC	Biorad	1645052	Electrophoresis module
Protein A from Staphylococcus aureus	Sigma-Aldrich	P7837-5MG	Coating of SPR chips
Sephadex G-25 fine, cross linked dextran	GE Helthcare	17003301	Chromatography resin
Silicone spoon	n/a	n/a	Spoon
Simply Blue SafeStain	Invitrogen	LC6060	Gel staining solution
Sodium acetate	Carl Roth GmbH	6773.1	Buffer component
Sodium acetate	Carl Roth GmbH	X891.1	Media component
Sodium azide	Sigma Aldrich	S2002-100G	Media component
Sodium chloride	Carl Roth GmbH	P029.2	Buffer component
Sodium citrate	Carl Roth GmbH	HN13.2	Buffer component
Sodium bisulfite	Carl Roth GmbH	243973-100G	Media component
Sodium phosophate	Carl Roth GmbH	T877.2	Media component
SPR Affinity Sensor - High Capacity Amine	Sierra Sensors GmbH/Bruker	SPR-AS-HCA	SPR chip
	Daltonics		
SPR-2/4 Surface Plasmon Resonance Analyzer	Sierra Sensors GmbH/Bruker	n/a	SPR device
	Daltonics		
SSM3	Stuart	10034264	Mini Gyro-rocker
Heated vessel, 20 L	Clatronic	n/a	Blanching chamber
Sterile syringes, 2 mL	B. Braun	4606027V	Syringes
Syringe adpter (Union Luer F)	GE Helthcare	181112-51	Syringe adapter
TE6101	Sartorius	TE6101	Precision scale
Tween-20 (Polysorbate)	Merck	8170721000	Buffer component
Unicorn 6.4	GE Helthcare	29056102	Chromatography software
Vacuum bags	Ikea	203.392.84	Plant storge
VelaPad 60	Pall	VP60G03KNH4	Filter housing
Vortex-Genie 2	Scientific industries	SI-0236	Vortex
XK-26/20 column housing	GE Helthcare	28-9889-48	Chromtography column

## ARTICLE AND VIDEO LICENSE AGREEMENT

Title of Article:

Activated Sepharose for the rapid development of affinity resins - antibo

Author(s):

Matthias Knödler, Clemens Rühl, Patrick Opdensteinen, Johannes Buye

Item 1: The Author elects to have the Materials be made available (as described at <http://www.jove.com/publish>) via:

☐

Standard Access

☒

Open Access

Item 2: Please select one of the following items:

☒

The Author is **NOT** a United States government employee.

☐

The Author is a United States government employee and the Materials were prepared in the course of his or her duties as a United States government employee.

☐

The Author is a United States government employee but the Materials were NOT prepared in the course of his or her duties as a United States government employee.

### ARTICLE AND VIDEO LICENSE AGREEMENT

1. **Defined Terms.** As used in this Article and Video License Agreement, the following terms shall have the following meanings: **"Agreement"** means this Article and Video License Agreement; **"Article"** means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein; **"Author"** means the author who is a signatory to this Agreement; **"Collective Work"** means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; **"CRC License"** means the Creative Commons Attribution-Non Commercial-No Derivs 3.0 Unported Agreement, the terms and conditions of which can be found at: <http://creativecommons.org/licenses/by-nc-nd/3.0/legalcode>; **"Derivative Work"** means a work based upon the Materials or upon the Materials and other pre-existing works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; **"Institution"** means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; **"JoVE"** means MyJove Corporation, a Massachusetts corporation and the publisher of The Journal of Visualized Experiments; **"Materials"** means the Article and / or the Video; **"Parties"** means the Author and JoVE; **"Video"** means any video(s) made by the Author, alone or in conjunction with any other parties, or by JoVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion

of the Article, and in which the Author may or may not appear.

2. **Background.** The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JoVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.

3. **Grant of Rights in Article.** In consideration of JoVE agreeing to publish the Article, the Author hereby grants to JoVE, subject to **Sections 4** and **7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the "Open Access" box has been checked in **Item 1** above, JoVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.

## ARTICLE AND VIDEO LICENSE AGREEMENT

4. **Retention of Rights in Article.** Notwithstanding the exclusive license granted to JoVE in **Section 3** above, the Author shall, with respect to the Article, retain the non-exclusive right to use all or part of the Article for the non-commercial purpose of giving lectures, presentations or teaching classes, and to post a copy of the Article on the Institution's website or the Author's personal website, in each case provided that a link to the Article on the JoVE website is provided and notice of JoVE's copyright in the Article is included. All non-copyright intellectual property rights in and to the Article, such as patent rights, shall remain with the Author.

5. **Grant of Rights in Video – Standard Access.** This **Section 5** applies if the "Standard Access" box has been checked in **Item 1** above or if no box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby acknowledges and agrees that, Subject to **Section 7** below, JoVE is and shall be the sole and exclusive owner of all rights of any nature, including, without limitation, all copyrights, in and to the Video. To the extent that, by law, the Author is deemed, now or at any time in the future, to have any rights of any nature in or to the Video, the Author hereby disclaims all such rights and transfers all such rights to JoVE.

6. **Grant of Rights in Video – Open Access.** This **Section 6** applies only if the "Open Access" box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby grants to JoVE, subject to **Section 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Video in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Video into other languages, create adaptations, summaries or extracts of the Video or other Derivative Works or Collective Works based on all or any portion of the Video and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. For any Video to which this **Section 6** is applicable, JoVE and the Author hereby grant to the public all such rights in the Video as provided in, but subject to all limitations and requirements set forth in, the CRC License.

7. **Government Employees.** If the Author is a United States government employee and the Article was prepared in the course of his or her duties as a United States government employee, as indicated in **Item 2** above, and any of the licenses or grants granted by the Author hereunder exceed the scope of the 17 U.S.C. 403, then the rights granted hereunder shall be limited to the maximum

rights permitted under such statute. In such case, all provisions contained herein that are not in conflict with such statute shall remain in full force and effect, and all provisions contained herein that do so conflict shall be deemed to be amended so as to provide to JoVE the maximum rights permissible within such statute.

8. **Protection of the Work.** The Author(s) authorize JoVE to take steps in the Author(s) name and on their behalf if JoVE believes some third party could be infringing or might infringe the copyright of either the Author's Article and/or Video.

9. **Likeness, Privacy, Personality.** The Author hereby grants JoVE the right to use the Author's name, voice, likeness, picture, photograph, image, biography and performance in any way, commercial or otherwise, in connection with the Materials and the sale, promotion and distribution thereof. The Author hereby waives any and all rights he or she may have, relating to his or her appearance in the Video or otherwise relating to the Materials, under all applicable privacy, likeness, personality or similar laws.

10. **Author Warranties.** The Author represents and warrants that the Article is original, that it has not been published, that the copyright interest is owned by the Author (or, if more than one author is listed at the beginning of this Agreement, by such authors collectively) and has not been assigned, licensed, or otherwise transferred to any other party. The Author represents and warrants that the author(s) listed at the top of this Agreement are the only authors of the Materials. If more than one author is listed at the top of this Agreement and if any such author has not entered into a separate Article and Video License Agreement with JoVE relating to the Materials, the Author represents and warrants that the Author has been authorized by each of the other such authors to execute this Agreement on his or her behalf and to bind him or her with respect to the terms of this Agreement as if each of them had been a party hereto as an Author. The Author warrants that the use, reproduction, distribution, public or private performance or display, and/or modification of all or any portion of the Materials does not and will not violate, infringe and/or misappropriate the patent, trademark, intellectual property or other rights of any third party. The Author represents and warrants that it has and will continue to comply with all government, institutional and other regulations, including, without limitation all institutional, laboratory, hospital, ethical, human and animal treatment, privacy, and all other rules, regulations, laws, procedures or guidelines, applicable to the Materials, and that all research involving human and animal subjects has been approved by the Author's relevant institutional review board.

11. **JoVE Discretion.** If the Author requests the assistance of JoVE in producing the Video in the Author's facility, the Author shall ensure that the presence of JoVE employees, agents or independent contractors is in accordance with the relevant regulations of the Author's institution. If more than one author is listed at the beginning of this Agreement, JoVE may, in its sole

## ARTICLE AND VIDEO LICENSE AGREEMENT

discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JoVE reserves the right, in its absolute and sole discretion and without giving any reason therefore, to accept or decline any work submitted to JoVE. JoVE and its employees, agents and independent contractors shall have full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JoVE has sole discretion as to the method of making and publishing the Materials, including, without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

**12. Indemnification.** The Author agrees to indemnify JoVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JoVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JoVE, making of videos by JoVE, or publication in JoVE or elsewhere by JoVE. The Author shall be responsible for, and shall hold JoVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to


the making of a video by JoVE its employees, agents or independent contractors. All sterilization, cleanliness or decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's expense. All indemnifications provided herein shall include JoVE's attorney's fees and costs related to said losses or damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JoVE, its employees, agents or independent contractors.

**13. Fees.** To cover the cost incurred for publication, JoVE must receive payment before production and publication of the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JoVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.

**14. Transfer, Governing Law.** This Agreement may be assigned by JoVE and shall inure to the benefits of any of JoVE's successors and assignees. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to be one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement is required per submission.

### CORRESPONDING AUTHOR

Name:	Johannes Buyel	
Department:	Integrated Production Platforms	
Institution:	Fraunhofer IME	
Title:	Head of Department	
Signature:	 <p>DN: cn=Johannes Buyel, o=Fraunhofer IME, ou=Integrated Production Platforms, email=johannes.buyel@rwth-aachen.de, c=DE Datum: 2019.02.25 09:00:11 +01'00'</p>	Date: 2019-02-25

Please submit a **signed** and **dated** copy of this license by one of the following three methods:

1. Upload an electronic version on the JoVE submission site
2. Fax the document to +1.866.381.2236
3. Mail the document to JoVE / Attn: JoVE Editorial / 1 Alewife Center #200 / Cambridge, MA 02140

## Response to review comments

Review Manuscript Number: 59933\_R0\_030119

Title: Activated Sepharose for the rapid development of affinity resins - antibody capture as a case study

### Editorial comments:

Editor comment 1: 4.2: How do you charge the column with nickel ions?

**Author response 1: We have added the information of how to charge the column with nickel ions.**

Editor comment 2: Figure 2A: Please define 'dps' in the legend.

**Author response 2: We added an definition of the abbreviation to the figure legend.**

Editor comment 3: Figure 2C: Most of the numbers here are not in Figure 2B.

**Author response 3: The figure legends states that the numbers in panel C correspond to those in panel B. In our understanding, this means: #3 in panel C = #3 in panel B. It does not imply completeness (i.e. for every number in C there is a number in B). We will gladly change the legend text if the Editor can suggest a better way to describe what we had in mind. Additionally, even though we did not claim this, most of the numbers (i.e. >50%) from panel C can be found in panel B (i.e. 5/8 =62.5%).**

### Reviewers' comments:

#### Reviewer #1:

Reviewer comment 1.1: The concept is good. However, the paper is written like a internal lab work protocol instead of publication

**Author response 1.1: The paper is written as a lab work protocol indeed, because this is the requirement of the journal.**

#### Major Concerns:

Reviewer comment 1.2: There are no results for product quality or comparison to product quality with historical methods of protein purification. It is difficult to assess efficiency or efficacy of product without any product quality impact. One of the reasons why other resins have failed to displace ProA is the efficiency of ProA resin to remove HCP. The author does not show any data.

**Author response 1.2: We agree that product quality is an important criterion when assessing a new purification method. Even though the manuscript at hand is meant to focus on the method of preparing the affinity ligand, rather than analyzing it, we mentioned in the results section "Testing 2F5 isolation using the DFE affinity resin" of the text that the recovery and purity (the latter being an indicator for the HCP removal mentioned by the reviewer) were similar to those achieved with Protein A. We have now expanded this information with data from our initial publication (Rühl et al. 2018 in J. Chrom. A) reporting that the equilibrium dissociation constant ( $K_D$ ) for DFE-purified and Protein A-purified 2F5 antibody were 791 pM and 763 pM respectively and thus very similar. We therefore believe that the 2F5 quality of both preparations was comparable.**

#### Reviewer #2:

Reviewer comment 2.1: My only point that presented a slight pause in my assessment of this paper was that the authors do not describe the reason that they choose to display their epitope on a protein, rather than using perhaps chemical synthesis to produce structurally similar compounds or even immobilize peptides on Sepharose. It is not an easy process to cultivate tobacco plants, produce and isolate the required protein and then immobilize it. IS it really significantly more advantageous than using Protein A,

because it does seem equally laborious and costly. Perhaps the authors should spend a paragraph commenting on this issue. However, I do believe that this paper is worthy of publication. A lot of work has been done and it will certainly prove useful.

Author response 2.1: The reviewer raises a valid point and we have added a paragraph in the discussion section to explain our perspective on the matter.

Reviewer #3:

Reviewer comment 3.1: Please add a Table to present the main faced problems during experiments and their possible solutions, for instance, a failure mode and effects analysis (FMEA) of the protocol proposed. This Table will contribute to troubleshooting when other people try to reproduce the experiment by following the protocol.

Author response 3.1: This is a good idea, we have added an according table at the end of the protocol section.

Reviewer comment 3.2: Please add some lines about the limitations and/or disadvantages and areas of improvement of the proposed method.

Author response 3.2: We have added a section highlighting the current limitations and challenges with the method at the end of the discussion.