

Journal of Visualized Experiments

Production of dynein and kinesin motor ensembles on DNA origami nanostructures for single molecule observation --Manuscript Draft--

Article Type:	Invited Methods Article - JoVE Produced Video
Manuscript Number:	JoVE60369R1
Full Title:	Production of dynein and kinesin motor ensembles on DNA origami nanostructures for single molecule observation
Keywords:	Cytoskeleton, dynein, kinesin, molecular motors, DNA origami, single molecule methods
Corresponding Author:	Nathan Derr, Ph.D. Smith College Northampton, MA UNITED STATES
Corresponding Author's Institution:	Smith College
Corresponding Author E-Mail:	nderr@smith.edu
Order of Authors:	Jingjie Hu Nathan Derr, Ph.D.
Additional Information:	
Question	Response
Please indicate whether this article will be Standard Access or Open Access.	Standard Access (US\$2,400)
Please indicate the city, state/province, and country where this article will be filmed . Please do not use abbreviations.	Northampton, MA, United States



Nathan Derr, Ph.D.

Department of Biological Sciences
44 College Lane
Smith College
Northampton, Massachusetts 01063
T (413) 585-3966
F (413) 585-3786
nderr@smith.edu

July 2, 2019

Dear Dr. Steindel,

We are pleased to submit our revised manuscript entitled “**Production of dynein and kinesin motor ensembles on DNA origami for single molecule observation.**”

We appreciate the investment of time and energy by the editors at *JoVE* and the reviewers who offered careful commentary and suggestions on our manuscript. We have addressed each point made by both editors and reviewers, and have included every suggestion in our revised manuscript. All major changes to the manuscript are listed in [blue text](#).

We have also included a Response to Reviewers, where we enumerate each change. Our responses are also listed in [blue text](#).

We appreciate the opportunity to submit our revised manuscript to *JoVE*!

Sincerely,

A handwritten signature in black ink, appearing to read 'Nathan Derr'.

Nathan D. Derr, Ph.D.
Assistant Professor of Biological Sciences
Director, Smith College Center for Microscopy and Imaging

TITLE:

Production of Dynein And Kinesin Motor Ensembles on DNA Origami Nanostructures for Single Molecule Observation

AUTHORS & AFFILIATIONS:

J. Hu¹ and N. D. Derr^{1,2}

¹Department of Biological Sciences, Smith College, Northampton MA

²Center for Microscopy and Imaging, Smith College, Northampton MA

Corresponding Author:

Nathan D. Derr (nderr@smith.edu)

Email Address of Co-Author:

Jingjie Hu (jhu60@smith.edu)

KEYWORDS:

cytoskeleton, dynein, kinesin, molecular motors, DNA origami, single molecule methods

SUMMARY:

The goal of this protocol is to form ensembles of molecular motors on DNA origami nanostructures and observe the ensemble motility using total internal reflection fluorescence microscopy.

ABSTRACT:

Cytoskeletal motors are responsible for a wide variety of functions in eukaryotic cells, including mitosis, cargo transport, cellular motility, and others. Many of these functions require motors to operate in ensembles. Despite a wealth of knowledge about the mechanisms of individual cytoskeletal motors, comparatively less is known about the mechanisms and emergent behaviors of motor ensembles, examples of which include changes to ensemble processivity and velocity with changing motor number, location, and configuration. Structural DNA nanotechnology, and the specific technique of DNA origami, enables the molecular construction of well-defined architectures of motor ensembles. The shape of cargo structures as well as the type, number and placement of motors on the structure can all be controlled. Here, we provide detailed protocols for producing these ensembles and observing them using total internal reflection fluorescence microscopy. Although these techniques have been specifically applied for cytoskeletal motors, the methods are generalizable to other proteins that assemble in complexes to accomplish their tasks. Overall, the DNA origami method for creating well-defined ensembles of motor proteins provides a powerful tool for dissecting the mechanisms that lead to emergent motile behavior.

INTRODUCTION:

Dynein and kinesin are cytoskeletal motor proteins responsible for myriad functions in eukaryotic cells¹. By converting the chemical energy of ATP hydrolysis into productive work, these motors translocate on microtubules to haul and distribute various intracellular cargos. They also

coordinate in the massive intracellular rearrangements associated with mitosis, where they exhibit orchestrated forces that contribute to the positioning and separation of chromosomes. Structural, biochemical, and biophysical assays, including single molecule observations, have revealed the mechanisms of these motors at the individual level (well-reviewed in previous works²⁻⁴). However, many of the motors' tasks require them to work in small ensembles of both similar and mixed motor types. Comparatively less is understood about the mechanisms that coordinate the activity and ultimate emergent motility of these ensembles^{5,6}. This knowledge gap is due, in part, to the difficulty in creating ensembles with controllable features, such as motor type and copy number. Over the past decade, the molecular construction techniques of DNA origami have been employed to solve this problem. For the microtubule based motors, some examples of these investigations include single molecule observations of ensembles of cytoplasmic dynein-1⁷⁻⁹, intraflagellar dynein¹¹, various kinesin motors^{12,13}, and mixtures of both dyneins and kinesins^{7,14,15}. Here, we provide details of the purification and oligonucleotide labeling of motors from yeast^{7,16-20}, the folding and purification of segmented DNA origami with tunable compliance⁸, and the imaging of the yeast motors propelling the chassis structures^{7,18}.

Constructing motor ensembles for in vitro single molecule observation requires three primary efforts. The first is the expression, purification and labeling of motor constructs suitable for attaching to DNA origami. The second is the production and purification of defined DNA origami structures (often termed "chassis"). And the third is the conjugation of the motors to the chassis structure followed by observation using total internal reflection fluorescence (TIRF) microscopy. Here, we provide established protocols for this process for recombinant microtubule-based motors purified from the yeast *Saccharomyces cerevisiae*^{7,16-19}. DNA origami-based motor ensembles have been investigated using both recombinant kinesin¹⁵ and dynein^{7,8,18} constructs produced in this yeast expression system¹⁶⁻¹⁹. This protocol is valid for these constructs, given that they are controlled by the galactose induced promoter, and fused to the same protein tags for purification (ZZ and TEV protease linker) and for DNA oligo conjugation (SNAPtag).

Specific yeast strains produce specific motor constructs. For example, the dynein used to study the role of cargo compliance was purified from strain RPY1084^{7,8}. In general, strains containing motor constructs with the appropriate genetic modifications for expression and purification can be requested from the laboratories having published the use of those motors. Constructs with novel attributes such as mutations or tags can be made using recombinant genetic techniques, such as lithium acetate transformation²¹ and commercial kits. Detailed protocols for creating modified motor proteins in yeast for single molecule studies have been published¹⁹. In addition to the motors being fused to the SNAPtag, the oligos used to label the motors must be conjugated to the SNAP substrate, benzylguanine (BG); previously published protocols describe the formation and purification of BG-oligo conjugates¹⁸. The overall strategy described here has also been employed for actin-based motors (see previous works for examples²²⁻²⁴), and motors purified from other organisms and expression systems (see previous works for examples^{7,9-14}).

Polymerized microtubules (MTs) are used in these experiments in two different procedures. MT affinity purification of functional motors requires MTs that are not labelled with other functional groups, while the motor-ensemble motility TIRF assay requires MTs labeled with biotin and

fluorophores. In all cases, MTs are stabilized with taxol to prevent denaturation. The MT affinity purification step is used to remove any non-motile motors with a high MT affinity, as these motors can alter ensemble motility if conjugated to a chassis. During this process, active motors unbind the MTs and remain in solution, while tight-binding motors spin down in the MT pellet. This helps ensure all motors on the chassis are from an active population.

A variety of DNA origami structures have been used to study cytoskeletal motor ensembles. As the mechanistic understanding of ensemble transport has increased, the DNA origami structures employed in experiments have grown in complexity. In principle, any structure could be adapted for this purpose provided it is modified to include single-stranded DNA attachment sites for motors and fluorophores. Specific chassis designs and attributes may be useful for probing particular questions about the emergent behavior of motors ensembles. For example, rigid rods have been used to develop foundational knowledge of how copy number affects transport by teams of dyneins and kinesins^{7,15,18}, and 2D platforms have been used to study myosin ensemble navigation of actin networks²². Structures with variable or tunable flexibility have been used to understand the roles of elastic coupling between motors and to probe how stepping synchronization affects motility^{8,24}. More recently, spherical structures are being used to gain insight into how geometrical constraints to motor-track binding affect the dynamics of motility²⁵.

In this protocol, we offer specific steps for ensemble experiments on segmented chassis with variable rigidity. Binding sites on the chassis are sometimes referred to as “handles”, while complementary DNA sequences that bind these handles are termed “antihandles”. The number of motors on these chassis is determined by which segments contain extended handle staples with complementarity to the antihandle oligo on the oligo-labeled motors. Using different handle sequences on different segments allows for binding of different types of motors to specific locations on the chassis. The chassis detailed here is composed of 7 sequential rigid segments, each comprised of 12 double-stranded DNA helices arranged in 2 concentric rings⁸. The rigid segments contain the motor handles and are connected through regions that can be either flexible single-stranded DNA or rigid double-stranded DNA, depending on the absence or presence, respectively, of “linker” staples. The compliance of the chassis structure is thus determined by the presence or absence of these “linker” staples. See previous reports for further details and specific DNA sequences⁸. In addition, multiple methods can be used to purify chassis²⁶. The rate-zonal glycerol gradient centrifugation method²⁷ is described here.

PROTOCOL:

1. Growth, expression and harvesting of motor proteins controlled by a galactose induced promoter

1.1. Using a yeast-peptone-dextrose (YPD) culture plate and a sterile inoculating stick, streak desired frozen yeast strain and incubate for 3–4 days at 30 °C.

1.2. Day 1 of culture growth: In the afternoon, add 10 mL of YP culture media with 2% dextrose to a 1" diameter glass culture tube and inoculate it with a single colony from the plate. Grow in a rapidly rotating roller drum at 30 °C overnight.

1.3. Day 2: In the afternoon, transfer the 10 mL of culture to a 250 mL flask containing 50 mL of YP culture media with 2% raffinose. Incubate overnight at 30 °C while shaking at 250 rpm on an orbital shaker.

1.4. Day 3: In the afternoon, transfer the 60 mL of culture to a 3 L flask containing 1 L of YP culture media with 2% galactose. Incubate overnight at 30 °C while shaking at 250 rpm on an orbital shaker.

1.5. Day 4: Starting in the mid-morning, monitor the optical density (OD) of the culture at 600 nm every 2 h. When the culture is between an OD 1.5 and 2, continue with the following steps to harvest the cells. Harvesting outside of this OD range may reduce yield due to low cell counts before OD 1.5, or cellular quiescence and protein degradation above OD 2.

1.6. Decant the cell culture into centrifuge bottles and spin at ~6200 x *g* at 4 °C for 8 min. Pour off the supernatant and discard.

1.7. Resuspend the cell pellet in the bottle with double-distilled H₂O (ddH₂O).

1.8. Spin the cells at ~6200 x *g* at 4 °C for 8 min. Pour off the supernatant and discard.

1.9. Depending on the viscosity of the cell pellet, add up to 2 mL of ddH₂O to create a solution fluid enough for pipetting.

1.10. Using a motorized pipette controller fitted with a 10 mL pipette, slowly dispense the cell slurry into liquid nitrogen one drop at a time. This will produce frozen pellets of yeast cells.

1.11. Store the frozen cells at -80 °C until ready to purify.

2. Purification of motor proteins from yeast cells

2.1. Cell lysis and soluble protein extraction

2.1.1. Prepare 1 mL of a fresh solution of 0.1 M PMSF in pure ethanol. Wait to add the PMSF to aqueous buffers as it is unstable in water; also, avoid exposure to water until within 20 min (the approximate half-life of PMSF in water) of use.

2.1.2. Prepare 50 mL of 4x lysis buffer on ice with supplements from the 5x stock of lysis buffer, but do not add PMSF until just before the buffer is applied to the yeast pellet.

2.1.3. Grind the liquid nitrogen-frozen yeast pellets into a fine powder with a blade-type coffee grinder that has been pre-chilled with liquid nitrogen.

Note: This protein extraction typically starts with the cell pellet harvested from 2 L of yeast culture. The starting amount of yeast culture can be adjusted up or down with commensurate adjustments to the volumes of IgG beads and DNA oligos in steps 2.2.2 and 2.3.1 below.

2.1.4. Aliquot 15 mL of the 4x lysis buffer with DTT and Mg-ATP on ice and add PMSF to the buffer to achieve a final concentration of 2 mM, completing the preparation of the 4x lysis buffer with supplements.

2.1.5. Collect the yeast powder into a pre-chilled ~100 mL glass beaker on ice. Add a small volume of the 4x lysis buffer with supplements into the powder so that the final concentration of the buffer does not exceed 1x. Typically, add ~1.5 mL of the buffer for every 10 mL of yeast powder.

2.1.6. Quickly thaw the powder to the liquid phase by placing the beaker in a 37 °C water bath with constant stirring using a spatula.

2.1.7. Place the beaker of lysate back on ice immediately after thawing, estimate the lysate volume using a 50 mL conical tube, and add more 4x lysis buffer with supplements so that the final concentration of the buffer is ~1x. Typically, 2 L of yeast culture yields a total of 25–35 mL of lysate containing 1x lysis buffer.

2.1.8. Evenly distribute the lysate to centrifuge bottles on ice. Carefully balance the bottles to a mass difference no more than 0.01 g between each pair, ensuring that each bottle is above the minimum volume required to prevent bottle collapse. Centrifuge at ~290,000 x *g* for 25 min at 4 °C.

2.1.9. Collect the supernatant containing soluble proteins in a 50 mL conical tube on ice, and discard the pellet containing cell debris and large organelles. Avoid collecting any cloudy portions of the supernatant, as they can clog the gravity-flow column used in subsequent steps. Save 10 µL of the supernatant for SDS-PAGE analysis.

2.2. IgG affinity purification

2.2.1. During the spin above, set up a gravity-flow chromatography column on ice or in a cold room.

2.2.2. Transfer 200 µL of 50% IgG affinity bead slurry to the column using a P-1000 pipette tip cut with a razor blade to enlarge the diameter of its entry port. If purifying more or less than 2 L of cell culture pellet, adjust the volume of bead slurry proportionately, with a minimum volume of 100 µL.

2.2.3. Aliquot another 15 mL of the 4x lysis buffer with DTT and Mg-ATP and add PMSF to the buffer to achieve a final concentration of 2mM. Dilute the 4x buffer on ice to 1x by adding ddH₂O.

2.2.4. Wash the beads 2x with 5 mL of 1x lysis buffer with supplements.

2.2.5. Resuspend the beads in 200 µL of 1x lysis buffer with supplements.

2.2.6. Add the bead suspension to the protein extract obtained from centrifugation and incubate the mixture at 4 °C for 1 h with gentle rotation.

2.2.7. During the incubation, prepare 25 mL of the wash buffer (recipe detailed in **Table 1**) on ice and 50 mL of 1x TEV buffer (recipe detailed in **Table 1**) on ice.

2.2.8. After the 1 h incubation, filter the lysate-bead mixture on ice or in the cold room through the same chromatography column used in step 2.2.2 above. Save 10 µL of the flow-through for SDS-PAGE analysis.

2.2.9. Wash the remaining motor-bound beads on ice 2x with 5 mL of wash buffer. Save 10 µL of the first wash for SDS-PAGE analysis.

2.2.10. Wash the beads on ice once with 5 mL of 1x TEV buffer. Allow the buffer to fully drain from the column.

2.3. Labeling with DNA oligonucleotides and TEV cleavage

2.3.1. Remove the chromatography column from the setup and cap the bottom of the column. Within the same column, incubate the motor-beads with 100 µL of 1x TEV buffer containing 10–20 µM of the purified BG-oligo at room temperature (RT) for 10–15 min. If purifying more or less than 2 L of cell culture pellet, adjust the volume of 1x TEV buffer and purified BG-oligos proportionately. Increase the incubation time according to manufacturer's instructions if needed to increase the yield and rate of motor labeling, but be aware that longer incubation times may also increase the proportion of dysfunctional motors due to protein denaturation at RT.

2.3.2. Gently resuspend the beads every minute of the incubation.

2.3.3. Wash the labeled motor-beads 4x with 4 mL of 1x TEV buffer using the same chromatography column and setup as before. Allow the final wash to fully drain from the column.

2.3.4. Cap the bottom of the column, resuspend the motor-beads in no more than 200 µL of 1x TEV Buffer, and transfer to a 2 mL round-bottom microcentrifuge tube.

2.3.5. Incubate the suspension with ~0.3 units of TEV protease per µL of motor-bead mixture at 16 °C for 1 h with gentle rotations. The tube should be mounted so as to minimize the total surface area of the tube with which the beads come into contact.

2.3.6. Centrifuge the tube in a microcentrifuge at 21,130 x *g* for 30 s in a cold room to concentrate the mixture at the bottom of the tube.

2.3.7. Still in the cold room, use a cut P-1000 pipette tip to transfer the mixture to a spin column and centrifuge at 21,130 x *g* for 30 s. Collect the filtrate containing the TEV-cleaved, purified motors. The TEV protease will be in the filtrate as well as the motors.

2.3.8. Save 10 µL of the filtrate for SDS-PAGE analysis. Aliquot the remaining filtrate in volumes of 2 µL for TIRF experiments or 50 µL for microtubule affinity purification. Flash freeze the aliquots in liquid nitrogen and store at -80 °C.

2.3.9. Resuspend the beads remaining in the filter in 1x TEV buffer for SDS-PAGE analysis. Use the same volume of 1x TEV buffer as in step 2.3.4 above.

3. Microtubule (MT) polymerization

3.1. Preparation of tubulin mixes for TIRF assays

3.1.1. In a cold room, separately dissolve the lyophilized tubulin of each type (unlabeled bovine tubulin, biotinylated tubulin, and fluorescent tubulin) in reconstitution buffer (recipe detailed in **Table 2**) to make a final concentration of 10 mg/mL. Let sit on ice for 10 min.

3.1.1.1. Mix the following components and let sit on ice for 10 min in the cold room (final concentrations are indicated parenthetically): 18 µL of bovine tubulin (~8.2 mg/mL), 2 µL of biotinylated tubulin (~0.9 mg/mL), and 2 µL of fluorescent tubulin (~0.9 mg/mL) .

3.1.2. Prepare 3 µL aliquots of the tubulin mixture and flash freeze in liquid nitrogen. Store the mixtures at -80 °C.

3.2. Preparation of tubulin for MT affinity purification of motors

3.2.1. In a cold room, dissolve lyophilized bovine tubulin in the reconstitution buffer to make a final concentration of 10 mg/mL. Let sit on ice for 10 min.

3.2.2. Prepare 3 µL aliquots of the tubulin mixture and flash freeze in liquid nitrogen. Store the mixtures at -80 °C.

3.3. Polymerization of tubulin into MTs

3.3.1. Remove a 3 µL aliquot of tubulin from the -80 °C freezer and very quickly and briefly thaw to the liquid phase by holding the bottom of the tube. Immediately place on ice and incubate for at least 3 min.

3.3.2. Gently layer 3 μ L of 2x polymerization mix (recipe detailed in **Table 2**) on top of the tubulin solution. Mix by gently flicking the tube, but do not mix by pipetting, as shear forces may disrupt MT nucleation and elongation.

3.3.3. Incubate the mixture in a water bath at 37 °C for 30 min.

3.3.4. Gently layer 6 μ L of RT 1x BRB80 with supplements (recipe detailed in **Table 2**) on top and mix by flicking. Do not mix by pipetting.

3.3.5. Incubate the mixture at 37 °C for at least 10 min.

3.3.6. Proceed to MT affinity purification, or if doing TIRF assays, incubate MTs in the dark at RT overnight to form longer MTs.

3.3.7. Store polymerized MTs for weeks in the dark at RT.

4. Microtubule (MT) affinity purification

4.1. Removal of unpolymerized tubulins by centrifugation

4.1.1. Mix the following components to make a 60% glycerol cushion: 1x BRB80 (without supplements; recipe detailed in **Table 2**), 20 μ M Taxol (dissolved in DMSO), 1 mM DTT, and 60% glycerol.

4.1.2. Transfer 60 μ L of the glycerol cushion to an ultracentrifuge tube. Gently layer 12 μ L of the unlabeled MTs polymerized previously on top of the cushion. To prevent shearing of MTs, transfer them using a cut pipette tip.

4.1.3. Centrifuge the mixture at $\sim 97,300 \times g$ at 22 °C for 15 min. After the spin, the excess tubulins remain in the top liquid layer, while the MTs form a pellet at the bottom. Mark the outer edge of the tube before the spin to help locate the pellet, as the pellet may not be visible with the naked eye.

4.1.4. Carefully remove the liquid layer ($\sim 12 \mu$ L) above the glycerol cushion and save 10 μ L for SDS-PAGE analysis.

4.1.5. Gently rinse the interface between the liquid layer and cushion with 20 μ L of 1x BRB80 with supplements. Remove and discard the 20 μ L rinse.

4.1.6. Carefully remove the glycerol cushion ($\sim 60 \mu$ L) and save 10 μ L for SDS-PAGE analysis. Be careful not to disturb the MT pellet.

4.1.7. Gently rinse the pellet with 60 μ L of 1x BRB80 with supplements (recipe detailed in **Table 2**). Be careful not to disturb the MT pellet. Discard the rinse solution.

4.1.8. Resuspend the MT pellet in 24 μL of Taxol-supplemented lysis buffer (recipe detailed in **Table 3**) using a cut pipette tip.

4.2. Purification of functional motors by MT-based affinity chromatography

4.3. Remove two 50 μL aliquots of motors purified from yeast from the -80°C freezer and quickly thaw to the liquid phase. Immediately place on ice.

4.4. Add the following components to a new ultracentrifuge tube in the indicated order and incubate the mixture at RT for 10 min: (1) 100 μL of the motors purified from yeast, (2) 29 μL of 5x ATP/Taxol mix (recipe detailed in **Table 3**), then (3) 12 μL of purified MTs transferred with a cut pipette tip.

4.5. Centrifuge the mixture at $\sim 97,300 \times g$ and 22°C for 15 min. After the spin, functional motors remain in the supernatant, and MTs form a pellet, along with the non-functional motors bound to them.

4.6. Collect the supernatant, flash freeze in 2 μL aliquots, and store at -80°C .

4.7. Save 10 μL of the supernatant for SDS-PAGE analysis. Resuspend the pellet in 141 μL of 1x BRB80 for SDS-PAGE analysis, too. Use protein standards, such as actin, to create a standard curve for quantifying the concentration of the purified motors as previously detailed¹⁷.

4.8. Use this concentration information when conjugating motors to the chassis in step 6.1.3.

5. Production of segmented DNA origami chassis

5.1. Formation of chassis

5.1.1. Order the staple oligonucleotides listed in previously published tables⁸ in 96 well plates wet at 250 μM in Tris buffer, or dry and then resuspend them to 250 μM with Tris buffer.

5.1.2. Create a pool of the core staples by mixing 5 μL of each core staple (see table S1 in Driller-Colangelo 2016⁸).

5.1.3. Create a pool of the linker staples by mixing 5 μL of each linker staple (see linker staple table in Driller-Colangelo 2016⁸).

5.1.4. Create a pool of the fluorophore binding staples by mixing 5 μL of each fluorophore handle staple (see fluorophore handle table in Driller-Colangelo 2016⁸).

5.1.5. Mix 50 μL folding reactions with the following components: 1x Folding buffer; 100 nM 8064 scaffold; 600 nM core staple pool; 600 nM fluorophore staple pool; 9 μM fluorophore strand (see

fluorophore antihandle table in Driller-Colangelo 2016⁸); for each motor binding site, either 4.2 μ M extended handle strands or 600 nM strands without handles as desired (see motor handle/antihandle staples table in Driller-Colangelo 2016⁸); 600 nM linkers as desired⁸; additional 6 mM $MgCl_2$; and water.

5.1.6. Fold in a thermal cycler using the following program: Rapid heating to 80 °C, cooling in single degree increments to 65 °C over 75 min, then additional cooling in single degree increments to 30 °C over 17.5 h.

5.1.7. Assay folding quality on a 2% agarose gel in 0.5x TBE buffer (see **Table 4** for recipe) supplemented with 11 mM $MgCl_2$ and DNA gel stain (see **Table of Materials**). Run gel in 0.5x TBE buffer supplemented with 11 mM $MgCl_2$ at 70 V for 60-90 min. Run gels in an ice water bath or in a cold room to prevent excessive heating and subsequent denaturation of chassis structures.

5.1.8. Image the gel using conditions suitable for the DNA gel stain used in step 5.1.7.

5.2. Purification of chassis

5.2.1. In the late afternoon the day before purification, create glycerol gradients by gently layering 80 μ L each of 45%, 40%, 35%, 30%, 25%, 20%, and 15% glycerol in 1x origami folding buffer in a centrifuge tube (see **Table 4** for recipe). Boundaries between layers should be slightly visible.

5.2.2. Incubate gradients at 4 °C overnight.

5.2.3. The next morning, add 45% glycerol in 1x origami folding buffer to the folded chassis solution to a final concentration of 10% glycerol. Mix gently and layer on the top of the gradient in the centrifuge tube.

5.2.4. Spin the gradient with the chassis at 243,000 x g for 130 min at 4 °C.

5.2.5. Collect 50 μ L fractions from the tube in a top to bottom direction.

5.2.6. Cast a 2% agarose gel in 0.5x TBE buffer supplemented with 11 mM $MgCl_2$ and DNA gel stain.

5.2.7. Load 5 μ L of each fraction on the gel and run gel in 0.5x TBE buffer supplemented with 11 mM $MgCl_2$ for 90–120 min at 70 V. Run gels in an ice water bath or in a cold room to prevent excessive heating and subsequent denaturation of chassis structures.

5.2.8. Image the gel using conditions suitable for the DNA gel stain used in step 5.2.6.

5.2.9. Select fractions for future experiments that exhibit well-folded monomeric structures free of unincorporated staples.

5.2.10. Quantified concentrations of selected fractions using appropriate spectroscopic methods, such as UV absorption at 260 nm. Use this concentration information when conjugating motors to the chassis in step 6.1.3.

6. Making slide assay chambers

6.1. Make a slide assay chamber by sticking two strips of double-sided tape to a glass slide and placing a coverslip on top. The chamber is the narrow space sandwiched between the coverslip and glass slide, flanked by the two strips of tape. **Figure 1** illustrates the assay chamber.

6.2. When preparing the slide with solutions, use a pipette to flow any fluid into the chamber from one side, and use a strip of filter paper to collect the flow of the fluid on the other side.

7. Motor-ensemble motility TIRF assay

7.1. Conjugation of motors to DNA chassis

7.1.1. On ice, prepare fresh 1 mL each of DTT-supplemented BRB80, Taxol-supplemented lysis buffer, and casein-Taxol-supplemented lysis buffer (recipes detailed in **Table 5**). Transfer 200 μ L of each buffer to a RT tube.

7.1.2. Still on ice, use the cold casein-Taxol-supplemented lysis buffer to prepare the 4x energy and 4x scavenger mixes (recipes detailed in **Table 5**).

7.1.3. Incubate 10 μ L of ~300 nM purified motors with 5 μ L of ~10 nM DNA chassis on ice for 15-30 min. These concentrations of motors have been shown to saturate the chassis' motor binding sites for this incubation time.

NOTE: motor occupancy is not 100%, likely due to stochastically missing handle staples in individual chassis structures^{8,28}.

7.1.4. During the incubation, dilute the biotin- and fluorophore-labeled MTs 100x with the RT Taxol-supplemented lysis buffer, and prepare a gel filtration resin appropriate for size exclusion column chromatography by following the procedure outlined below.

7.2. Removal of excess motors from motor-chassis conjugates by size exclusion column chromatography

7.2.1. In a 50 mL conical tube, wash 5 mL of the gel filtration resin 2x with 45 mL of ddH₂O. For each wash, mix the resin with water in the conical tube, and spin the mixture at 460 x g for 1 min. Discard the supernatant and save the resin.

7.2.2. Wash the resin 2x with 45 mL of 1x lysis buffer (recipe detailed in **Table 5**) using the same method described above.

7.2.3. Resuspend the washed resin in 1x lysis buffer in a 1:1 ratio, so that the resin becomes a ~50% slurry in buffer. The washed resin can be stored at 4 °C for at least 1 month.

7.2.4. Transfer 800 µL of the resin suspension to a spin column. Drain the excess buffer by gravity flow for 5 min. Remove any remaining buffer with a 2s spin at 1,000 x *g*. The final resin volume should be around 350–400 µL.

7.2.5. Dilute the motor-chassis mix with the casein-Taxol-supplemented lysis buffer to a final volume of 50 µL. Transfer the diluted mix to the spin column packed with resin.

7.2.6. Centrifuge the spin column at 1,000 x *g* for 6 s (this time is inclusive of the acceleration and deceleration time). Collect the pure motor-chassis conjugates by saving the filtrate and discard the column which retains the excess motors.

7.3. Preparation of slides for imaging

7.3.1. Flow 13 µL of 1 mg/mL biotinylated bovine serum albumin (biotin-BSA) into a slide assay chamber. Incubate for 2 min to allow binding of BSA to glass.

7.3.2. Wash the chamber 2x with 20 µL of the RT DTT-supplemented BRB80 by flowing in the buffer on one side of the chamber and wicking the excess fluid away from the other side using a strip of filter paper.

7.3.3. Flow in 20 µL of 0.5 mg/mL streptavidin. Incubate for 2 min to allow binding of streptavidin to the biotin on the BSA.

7.3.4. Wash the chamber 2x with 20 µL of the RT Taxol-supplemented lysis buffer.

7.3.5. Gently flow in 20 µL of the diluted MTs with a cut pipette tip. Incubate for 2 min to allow binding between the biotin on MTs and streptavidin.

7.3.6. Wash 2x with 20 µL of the RT casein-Taxol-supplemented lysis buffer. Incubate for 2 min to allow casein to permeate the entire chamber.

7.3.7. Dilute the purified motor-chassis conjugates 5x to 10x to single-molecule conditions (~10–100 pM) with the cold casein-Taxol-supplemented lysis buffer. Mix together the following components to produce a final motor-chassis mixture: 10 µL of the motor-chassis dilution, 5 µL of 4x energy mix, and 5 µL of 4x scavenger mix.

7.3.8. Flow in 20 µL of the final motor-chassis mixture to the assay slide chamber and proceed to TIRF microscopy.

7.4. Imaging and data acquisition

7.4.1. Image the slide immediately with a TIRF microscope. Typically, each slide remains useable for between 30 and 60 min.

7.4.2. Acquire a still image in the MT channel, and a movie in the chassis channel. For the motors in this protocol, 10 min movies with a frame rate of 0.5 fps and exposure time of 200 ms are appropriate.

7.4.3. From the chassis movie, generate one kymograph for each MT in ImageJ or a similar image processing software²⁹. Analyze the velocities and run lengths of motor-chassis ensembles by measuring the slopes and horizontal distances of the runs on the kymograph³⁰.

REPRESENTATIVE RESULTS:

Successful purifications of motors and chassis structures were assayed by gel electrophoresis. SDS-PAGE analysis confirmed the successful extraction of dynein from yeast (**Figure 2**), as the final filtrate collected in step 2.3.7 showed a clear, sharp band at the position of ~350 kDa. As expected, this dynein band was absent from the flowthrough and wash that removed unwanted proteins, and the beads from which dynein was cleaved. The observation suggests that the IgG affinity purification and TEV protease cleavage were both highly efficient. Additionally, TEV protease was also present in the final filtrate and formed a clear band at ~50 kDa.

The successful MT affinity purification of dynein and kinesin proteins was also confirmed with SDS-PAGE analysis (**Figure 3**). While dynein showed up as a clear single band at ~350 kDa, kinesin showed up as slightly smearing multiple bands at ~120 kDa, possibly due to the presence of both phosphorylated and dephosphorylated forms of the protein¹⁴ and variable yields in oligo-labeling. A comparison between the dynein and kinesin bands before and after this MT affinity purification revealed a reduction in the motor concentration, as indicated by the decrease in the band intensity, which was likely due to the removal of non-functional motors. Despite the reduction, the concentrations of motors retained were sufficient for effective TIRF assays. In addition to the TEV protease, a noticeable amount of tubulin (~51 kDa) was present in the final supernatant, most likely due to the gradual decomposition of MTs during the experiment, or incomplete removal of the excess tubulin through centrifugation. However, the consistent motility of motor-chassis ensembles shown in TIRF assays suggests that tubulin and TEV protease did not interfere with motor functions (see **Figure 6**).

Folding of DNA origami structures was assayed by agarose gel electrophoresis. **Figure 4** depicts the results of a gel analyzing the folding of a flexible segmented chassis with 2 motor binding sites. The shift in mobility between the pure unfolded scaffold strand (Lane 1) and the folding reaction (Lane 2) indicates origami folding. Additionally, the folding reaction in Lane 2 indicates the presence of some multimerization of chassis structures. Multimerization typically occurs, and requires subsequent purification of the well-folded monomeric structures. The unincorporated excess staples were also visible, displaying a high degree of mobility through the gel.

Folded origami reactions were purified to remove excess unincorporated staples and multimers of the chassis structure. **Figure 5** shows the results of a glycerol gradient purification of a flexible chassis with 7 motor binding sites. The early fractions correspond to the low glycerol density at the top of the tube. They contain the excess staples. The late fractions correspond to the high glycerol densities at the bottom of the tube and contain multimers and aggregates of folded structures. In this gel, fraction 7 indicates a suitable fraction containing well-folded monomeric chassis. Note that the well-folded structure is isolated from both excess staples and multimers. While this gel is representative, and fraction 7 often contains useable chassis, each purification experiment yields slightly different results and fractions should always be assayed to determine which fraction is best for motility assays.

The motility of motor-chassis ensembles is easily detectable and measurable on the kymographs generated from TIRF movies. For instance, the kymographs (**Figure 6**) of flexible chassis conjugated to seven dynein proteins (“7D” ensembles) show highly processive runs at relatively consistent velocities, demonstrating that the ensembles were active and motile in the reconstituted system, and that MT affinity purification successfully removed most of the non-functional, immobile dyneins that could slow or stall the ensembles. The same TIRF experiment has been successfully performed on other chassis types with different compliance and motor numbers to reveal the effects of these factors on dynein teamwork^{8,9}.

FIGURE AND TABLE LEGENDS:

Figure 1. Schematic of the slide assay chamber. The coverslip sits atop two strips of double-sided tape. Solutions are pipetted in one side, and extracted with filter paper on the other.

Figure 2. SDS-PAGE analysis of the purification of dynein from yeast. Yeast cells were lysed and centrifuged to collect the lysate (lane 1) containing total soluble proteins. The affinity between IgG on column beads and the ZZ tag on the recombinant dynein construct was exploited for this purification. The flow-through (lane 2) was collected from the mixture of the lysate and beads after it passed through a chromatography column. The dynein-bound beads were washed with buffers, and the first wash with the wash buffer (lane 3) was collected. Dynein was then conjugated to a DNA oligo. After TEV protease cleavage, centrifugation in a spin column separated the filtrate containing dynein (lane 4) and the residual beads (lane 5). All samples collected from the purification were denatured with 1x LDS Sample Buffer and 1x Reducing Agent, and loaded onto a 4–12% Bis-Tris gel. The gel was run in 1x MOPS buffer at 200 V for ~1 h, and stained with SYPRO Red Protein Gel Stain for imaging under UV light. Notably, the clear, sharp band at ~350 kDa in lane 4 indicates the presence of concentrated purified dynein in the final filtrate, while the band at ~50 kDa in the same lane indicates the co-presence of TEV protease.

Figure 3. SDS-PAGE analysis of the MT affinity purification of functional dynein and kinesin proteins. Kinesin and dynein purified from yeast (lanes 1 and 3) were mixed with polymerized MTs and ATP, and ultracentrifugation was performed to separate the functional motors in the

supernatants (lanes 2 and 4) from the non-functional motors that co-pellet with MTs. The motor samples before and after purification were denatured with 1x LDS Sample Buffer and 1x Reducing Agent, and loaded onto a 4–12% Bis-Tris gel. The gel was run in 1x MOPS buffer at 200 V for ~1 h, and stained with SYPRO Red Protein Gel Stain for imaging under UV light. The intensity of the motor bands (~120 kDa for kinesin, and ~350 kDa for dynein) appeared to decrease after the MT affinity purification, indicating a reduction in the motor concentrations. Noticeable concentrations of tubulin and TEV protease were present in the post-purification motor samples.

Figure 4. Agarose gel analysis of folded DNA origami structure. Folded DNA origami structures are assessed by gel electrophoresis. Lane 1 is the pure unfolded scaffold strand while lane 2 is the product of the folding reaction. Gel was run at 70 V in an ice water bath for 90 min in 0.5x TBE buffer supplemented with 11 mM MgCl₂. A shift in mobility indicates origami folding. Unincorporated staples and chassis multimers were also visible in Lane 2.

Figure 5. Agarose gel analysis of glycerol gradient purification of DNA origami chassis. The quality of the glycerol gradient purification can be determined by gel electrophoresis of the fractions from the centrifuge gradient. Low number fractions are from the top of the tube and correspond to low density of glycerol. High numbered fractions are from the bottom of the tube and correspond to high densities of glycerol. The excess staples, monomeric chassis, and multimeric chassis are all visible. Fraction 7 contains chassis suitable for TIRF microscopy as they are monomeric and excess staples are absent. Gel was run at 70 V in an ice water bath for 120 min in 0.5x TBE buffer supplemented with 11 mM MgCl₂.

Figure 6. Kymographs showing the motility of dynein-chassis ensembles on single MTs. Dynein proteins extracted from yeast and purified with MTs were conjugated to flexible chassis structures. Each chassis had seven binding sites for dynein and formed a “7D” ensemble. The movement of these ensembles on MTs was recorded in a 10 min movie (200 ms exposure, 0.5 fps) during a TIRF assay. Kymographs from two MTs were generated from this movie in ImageJ by tracing along a single MT. The vertical and horizontal red bars in the top left corner of each image indicate 2 min and 20 μm, respectively. Each bright line records the movement of one ensemble, with the inverse of the slope indicating the velocity and horizontal displacement indicating the run length. With TIRF assays and kymography, the motility of motor-chassis ensembles becomes easily detectable and directly measurable.

Table 1. Buffers for the IgG affinity purification of motor proteins from yeast cells (Protocol section 2).

Table 2. Buffers for the polymerization of microtubules (Protocol section 3).

Table 3. Buffers for the microtubule affinity purification of functional motor proteins (Protocol section 4).

Table 4. Buffers for the production of segmented DNA origami chassis (Protocol section 5).

Table 5. Buffers (a) and mixes (b) for the motor-ensemble motility TIRF assay (Protocol section 7). Details on the Oxygen Scavenger System used to make the Scavenger Mix can be found elsewhere¹⁷.

DISCUSSION:

The molecular construction techniques of DNA origami provide a unique way to construct motor ensembles with defined architectures, motor numbers, and types, enabling studies of how emergent behavior arises from specific motor configurations³¹. As structural and cellular studies continue to elucidate examples of cytoskeletal motors working in teams, techniques for isolating and investigating the biophysical and biochemical mechanisms of motors in ensembles are growing in utility. For example, cryo-EM has shown that dynactin can bind 2 individual dynein motors, and that such pairings have different motility than the individual motors¹⁰. In addition, DNA-based construction was used to determine if mammalian dynein, when activated by dynactin and bicD2, could match the force of kinesin-1 in a tug of war scenario¹⁴. In another mixed-motor study, DNA origami was used to decouple the effects of opposite polarity motors and their regulatory binding proteins by spatially separating them on an origami structure¹⁵. As more structural and regulatory determinants of motility are found, DNA-origami-based techniques should prove useful in determining the specific biochemical and biophysical contributors to the emergent motility of motor ensembles. The molecular construction techniques enabled by DNA origami are particularly useful because the emergent phenomenological outcomes of ensembles are difficult to predict. This is due in part to the myriad factors that contribute to the motility of the individual motors within the ensemble^{5,6,31}.

Examples of our previous structures include monolithic rods and segmented rods with variable rigidity. Current efforts explore spherical structures as well²⁵. Others have employed morphologies such as planar structures and rods^{22,24}. Likewise, by using motor handles with orthogonal DNA sequences, different types of motors can be bound to the same chassis structure^{7,14,15,18}. This approach enables studies of the opposing actions of dyneins and kinesins, minus- and plus-end directed kinesins, and minus- and plus-end directed myosins. It also enables the introduction of a mutant among otherwise wild-type ensembles, allowing the specific biochemical contributors to the emergent motility to be deciphered⁷. Because of the ability to bind multiple fluorophores to each individual structure, imaging in TIRF and subsequent analysis by kymography or particle tracking is possible. Previous reports show analysis of kymography data and statistical evaluation of chassis structures with variable compliance⁸. While cytoskeletal motors have proven to be an exciting early application of using DNA origami as a molecular breadboard³², other proteins and protein systems will also benefit from these methods.

ACKNOWLEDGEMENTS:

We thank K. Chau, J. Morgan, and A. Driller-Colangelo for contributing to the techniques of the segmented DNA origami chassis. We also thank former members of the Reck-Peterson and Shih laboratories for helpful discussions and contributions to the original development of these techniques. We thank J. Wopereis and the Smith College Center for Microscopy and Imaging and

L. Bierwert and the Smith College Center for Molecular Biology. We gratefully acknowledge the NSF MRI program for the acquisition of a TIRF microscope.

DISCLOSURES:


The authors have nothing to disclose.

REFERENCES:

1. Vale, R. D. The molecular motor toolbox for intracellular transport. *Cell*. **112** (4), 467–480 (2003).
2. Cianfrocco, M. A., DeSantis, M. E., Leschziner, A. E., Reck-Peterson, S. L. Mechanism and regulation of cytoplasmic dynein. *Annual Review of Cell and Developmental Biology*. **31** (1), 83–108, doi:10.1146/annurev-cellbio-100814-125438 (2015).
3. Roberts, A. J., Kon, T., Knight, P. J., Sutoh, K., Burgess, S. A. Functions and mechanics of dynein motor proteins. *Nature Reviews Molecular Cell Biology*. **14** (11), 713–726, doi:10.1038/nrm3667 (2013).
4. Hancock, W. O. The Kinesin-1 Chemomechanical Cycle: Stepping Toward a Consensus. *Biophysical Journal*. **110** (6), 1216–1225, doi:10.1016/j.bpj.2016.02.025 (2016).
5. McLaughlin, R. T., Diehl, M. R., Kolomeisky, A. B. Collective dynamics of processive cytoskeletal motors. *Soft Matter*. **12** (1), 14–21, doi:10.1039/c5sm01609f (2016).
6. Hancock W. O. Bidirectional cargo transport: moving beyond tug of war. *Nature Reviews Molecular Cell Biology*. **15** (9), 615–628 (2014).
7. Derr, N. D. et al. Tug-of-war in motor protein ensembles revealed with a programmable DNA origami scaffold. *Science*. **338** (6107), 662–665, doi:10.1126/science.1226734 (2012).
8. Driller-Colangelo, A. R., Chau, K. W. L., Morgan, J. M., Derr, N. D. Cargo rigidity affects the sensitivity of dynein ensembles to individual motor pausing. *Cytoskeleton*. **73** (12), 693–702, doi:10.1002/cm.21339 (2016).
9. Torisawa, T. et al. Autoinhibition and cooperative activation mechanisms of cytoplasmic dynein. *Nature Cell Biology*. **16** (11), 1118–1124, doi:10.1038/ncb3048 (2014).
10. Urnavicius, L. et al. Cryo-EM shows how dynactin recruits two dyneins for faster movement. *Nature*. **554** (7691), 202–206, doi:10.1038/nature25462 (2018).
11. Toropova, K., Mladenov, M., Roberts, A. J. Intraflagellar transport dynein is autoinhibited by trapping of its mechanical and track-binding elements. *Nature Structural & Molecular Biology*. **24** (5), 461–468, doi:10.1038/nsmb.3391 (2017).
12. Furuta, K., Furuta, A., Toyoshima, Y. Y., Amino, M., Oiwa, K., Kojima, H. Measuring collective transport by defined numbers of processive and nonprocessive kinesin motors. *Proceedings of the National Academy of Sciences of the United States of America*. **110** (2), 501–506, doi:10.1073/pnas.1201390110 (2013).
13. Rogers, A. R., Driver, J. W., Constantinou, P. E., Kenneth Jamison, D., Diehl, M. R. Negative interference dominates collective transport of kinesin motors in the absence of load. *Physical Chemistry Chemical Physics*. **11** (24), 4882–4889, doi:10.1039/b900964g (2009).

- 743 14. Belyy, V. et al. The mammalian dynein-dynactin complex is a strong opponent to
744 kinesin in a tug-of-war competition. *Nature Cell Biology*. **18** (9), 1018–1024,
745 doi:10.1038/ncb3393 (2016).
- 746 15. Roberts, A. J., Goodman, B. S., Reck-Peterson, S. L. Reconstitution of dynein transport
747 to the microtubule plus end by kinesin. *eLife*. **3**, e02641, doi:10.7554/eLife.02641
748 (2014).
- 749 16. Reck-Peterson, S. L. et al. Single-molecule analysis of dynein processivity and
750 stepping behavior. *Cell*. **126** (2), 335–348, doi:10.1016/j.cell.2006.05.046 (2006).
- 751 17. Gennerich, A., Reck-Peterson, S. L. Probing the force generation and stepping
752 behavior of cytoplasmic Dynein. *Methods in Molecular Biology*. **783** (Chapter 4), 63–
753 80, doi:10.1007/978-1-61779-282-3_4 (2011).
- 754 18. Goodman, B. S., Reck-Peterson, S. L. Engineering defined motor ensembles with DNA
755 origami. *Methods in Enzymology*. **540**, 169–188, doi:10.1016/B978-0-12-397924-
756 7.00010-8 (2014).
- 757 19. Rao L., Hülsemann M., Gennerich A. Combining Structure–Function and Single-
758 Molecule Studies on Cytoplasmic Dynein. In: Peterman E. (eds), *Single Molecule*
759 *Analysis. Methods in Molecular Biology*, vol. 1665. Humana Press, New York, NY.
760 (2018).
- 761 20. Qiu, W., Derr, N. D., et al. Dynein achieves processive motion using both stochastic
762 and coordinated stepping. *Nature Structural & Molecular Biology*. **19** (2), 193–200,
763 doi:10.1038/nsmb.2205 (2012).
- 764 21. Gietz R.D. Yeast Transformation by the LiAc/SS Carrier DNA/PEG Method. In: Smith
765 J., Burke D. (eds.), *Yeast Genetics. Methods in Molecular Biology (Methods and*
766 *Protocols)*, vol. 1205. Humana Press, New York, NY. (2014).
- 767 22. Hariadi, R. F., Cale, M., Sivaramakrishnan, S. Myosin lever arm directs collective
768 motion on cellular actin network. *Proceedings of the National Academy of Sciences*
769 *of the United States of America*. **111** (11), 4091–4096, doi:10.1073/pnas.1315923111
770 (2014).
- 771 23. Kremontsova, E. B., Furuta, K., Oiwa, K., Trybus, K. M., Ali, M. Y. Small teams of
772 myosin Vc motors coordinate their stepping for efficient cargo transport on actin
773 bundles. *The Journal of Biological Chemistry*. **292** (26), 10998–11008,
774 doi:10.1074/jbc.M117.780791 (2017).
- 775 24. Hariadi, R. F. et al. Mechanical coordination in motor ensembles revealed using
776 engineered artificial myosin filaments. *Nature Nanotechnology*. **10** (8), 696–700,
777 doi:10.1038/nnano.2015.132 (2015).
- 778 25. Hu, J. J., Morgan, J., Yang, Y., Lin, C., Derr, N. D. Spherical DNA Origami as a
779 Programmable Cargo Structure for Investigating the Emergent Motility of Dynein and
780 Kinesin Ensembles. *Biophysical Journal*. **116** (3), 408A,
781 doi:10.1016/j.bpj.2018.11.2200 (2019).
- 782 26. Wagenbauer, K. F. et al. How We Make DNA Origami. *ChemBiochem*. **18** (19), 1873–
783 1885, doi:10.1002/cbic.201700377 (2017).
- 784 27. Lin, C., Perrault, S. D., Kwak, M., Graf, F., Shih, W. M. Purification of DNA-origami
785 nanostructures by rate-zonal centrifugation. *Nucleic Acids Research*. **41** (2), e40,
786 doi:10.1093/nar/gks1070 (2013).


- 787 28. Ke, Y., Voigt, N. V., Gothelf, K. V., Shih, W. M. Multilayer DNA origami packed on
788 hexagonal and hybrid lattices. *Journal of the American Chemical Society*. **134**, 1770–
789 1774 (2012).
- 790 29. Schneider, C. A., Rasband, W. S., Eliceiri, K. W. NIH Image to ImageJ: 25 years of
791 image analysis. *Nature Methods*. **9** (7), 671–675 (2012).
- 792 30. Reck-Peterson, S. L., Derr, N. D., Stuurman, N. Imaging single molecules using total
793 internal reflection fluorescence microscopy (TIRFM). *Cold Spring Harbor Protocols*.
794 **2010** (3), pdb.top73–pdb.top73, doi:10.1101/pdb.top73 (2010).
- 795 31. Derr, N. D. Interactions of Multiple Dynein Motors Studied Using DNA Scaffolding. In
796 Hirose, K. (ed.), *Handbook of Dynein*, 2nd Edition. Stanford Publishing (2019).
- 797 32. Rothemund, P. W. K. Folding DNA to create nanoscale shapes and patterns. *Nature*.
798 **440** (7082), 297–302, doi:10.1038/nature04586 (2006).
799
800




Double-sided Tape

Coverslip

Glass Slide










Figure 6

Time



Position on the MT

Buffer Name	Composition	Step(s) Used
5x Lysis Buffer	150 mM HEPES (pH 7.4)	-
	250 mM KAcetate	
	10 mM MgAcetate	
	5 mM EGTA (pH 7.5)	
	50% glycerol	
4x Lysis Buffer With Supplements	4x Lysis Buffer	2.1-2.2
	4 mM DTT	
	0.4 mM Mg-ATP	
	2 mM PMSF	
Wash Buffer	1x Lysis Buffer With Supplements	2.2
	250 mM KCl	
	0.1% Triton X-100	
5x TEV Buffer	50 mM Tris-HCl (pH 8.0)	-
	150 mM KCl	
	10% Glycerol	
1x TEV Buffer	1x TEV Buffer	2.2-2.3
	1 mM DTT	
	0.1 mM Mg-ATP	
	0.1% Triton X-100	
	0.5 mM PMSF	

Comment

Filter sterilize the buffer. It can be stored at RT in a properly sealed container for a year.

Make this 4x buffer from the 5x lysis buffer above. Prepare a buffer without PMSF first, and add the compound (dissolved in pure ethanol) to a small aliquot of the buffer right before each step that requires it.

To make the buffer, add KCl, Triton X-100, and ddH₂O to the 4x lysis buffer with DTT and Mg-ATP, but do not add PMSF until right before use.

Filter sterilize the buffer. It can be stored at RT in a properly sealed container for a year.

To make the buffer, add DTT, Mg-ATP, Triton X-100, and ddH₂O to the 5x stock TEV buffer, but do not add PMSF until right before use.

Buffer Name	Composition	Step(s) Used
5x BRB80	400 mM PIPES	-
	10 mM MgCl ₂	
	5 mM EGTA	
	Adjust pH to 6.8 with KOH	
1x BRB80 With Supplements	1x BRB80	3.3 & 4.1-4.2
	20 M Taxol (dissolved in DMSO)	
	1 mM DTT	
2x Polymerization Mix	2x BRB80 (without supplements)	3.3
	2 mM DTT	
	2 mM Mg-GTP	
	20% DMSO	
Reconstitution Buffer	1x BRB80 (without supplements)	3.1
	1 mM DTT	
	1 mM Mg-GTP	

Comment

Filter sterilize the buffer. It can be stored at RT in a properly sealed container for a year.

Must be freshly made from the 5x BRB80 stock for every experiment.

Flash freeze the mix in small aliquots and store at -80 °C.

Must be freshly made for every experiment.

Buffer Name	Composition	Step(s) Used
Taxol-Supplemented Lysis Buffer	1x Lysis Buffer	4.1
	20 M Taxol (dissolved in DMSO)	
	1 mM DTT	
5x ATP/Taxol Mix	1x Lysis Buffer	4.2
	25 mM Mg-ATP	
	50 M Taxol (dissolved in DMSO)	
	5 mM DTT	

Comment

Must be freshly made for every purification.

Must be freshly made for every purification.

Buffer Name			Composition	Step(s) Used
20x Origami Buffer	Folding		100 mM Tris pH 8.0	-
			20 mM EDTA	
			200 mM MgCl ₂	
1x Origami Buffer	Folding		5 mM Tris pH 8.0	5.1-5.2
			1 mM EDTA	
			10 mM MgCl ₂	
0.5x TBE			45 mM Tris pH 8.0	5.1-5.2
			45 mM Boric Acid	
			1 mM EDTA	

Comment

Can be stored at RT in a properly sealed container for up to a year.

Make the buffer fresh by diluting the 20x stock with ddH₂O before every experiment.

Can be stored at RT in a properly sealed container for up to a year.

Buffer Name	Composition	Step(s) Used
DTT-Supplemented BRB80	1x BRB80 1 mM DTT	7.3
Taxol-Supplemented Lysis Buffer	1x Lysis Buffer 20 M Taxol (dissolved in DMSO) 1 mM DTT	7.1 & 7.3
Casein-Taxol-Supplemented Lysis Buffer	1x Lysis Buffer 20 M Taxol (dissolved in DMSO) 1 mM DTT ~2.5 mg/ml Casein (dissolved in Tris-HCl at pH 8.0)	7.1-7.3
1x Lysis Buffer	30 mM HEPES (pH 7.4) 50 mM KAcetate 2 mM MgAcetate 1 mM EGTA (pH 7.5) 10% glycerol	7.2
4x Energy Mix	22.5 µL 1x Casein-Taxol-Supplemented Lysis Buffer 1 µL 0.1 M Mg-ATP 1 µL 40% Glucose 0.5 µL β-Mercaptoethanol	7.3
4x Scavenger Mix	24 µL 1x Casein-Taxol-Supplemented Lysis Buffer 1 µL 1x Oxygen Scavenger System	7.3

Comment

Must be freshly made before every TIRF experiment.

Must be freshly made before every TIRF experiment.

Must be freshly made before every TIRF experiment.

Make this buffer fresh by diluting the 5x stock (recipe detailed in Table 1) with ddH₂O.

Must be freshly made before every TIRF experiment; volumes indicated are for a final volume of 25 µL.

Must be freshly made before every TIRF experiment; volumes indicated are for a final volume of 25 µL.

[illegible]

[illegible]

This piece of the submission is being sent via mail.

Response to review

We thank the editors and reviewers for their careful and considerate review of our manuscript. With the suggested revisions, we believe the manuscript will better serve the community and aid in the adoption of these techniques. We particularly appreciate the investment of time and energy to review, and we are grateful that the methods we have developed will be made more accessible through the production of these detailed protocols. We have addressed each point made by the editors and reviewers, and enumerate these changes on a point-by-point basis below. Both our responses to the reviewers and major changes in the manuscript are in blue text. By addressing all of these points, we believe we have made this protocol both stronger and more clear. (Please note that some changes we made based on the reviewers' comments have been moved out of the protocol itself and into the **Introduction** as per the style requirements of *JoVE*.)

Editorial comments:

General:

1. Please take this opportunity to thoroughly proofread the manuscript to ensure that there are no spelling or grammar issues.
2. Please ensure that the manuscript is formatted according to JoVE guidelines—letter (8.5" x 11") page size, 1-inch margins, 12 pt Calibri font throughout, all text aligned to the left margin, single spacing within paragraphs, and spaces between all paragraphs and protocol steps/substeps.
3. Please ensure that all references appear as numbered superscripts after the appropriate statement(s).
4. JoVE cannot publish manuscripts containing commercial language. This includes trademark symbols ([™]), registered symbols ([®]), and company names before an instrument or reagent. Please limit the use of commercial language from your manuscript and use generic terms instead. All commercial products should be sufficiently referenced in the Table of Materials and Reagents.

For example: Beckman Coulter, GE Healthcare, ProTEV, Millipore, 3M, Sephadryl,

We have made these needed edits.

Protocol:

1. Everything in the protocol should be in a numbered step in the imperative tense and of no more than 4 sentences, header, or 'Note'. Please move the introductory paragraphs of the protocol sections to the Introduction, Results, or Discussion (as appropriate) or break into steps. Please make recipes (2.4, 3.4, 4.3, 5.3, 7.5) into imperative steps or, alternatively, make into tables.

We have made these needed edits.

2. For each protocol step, please ensure you answer the "how" question, i.e., how is the step performed? Alternatively, add references to published material specifying how to perform the

protocol action. If revisions cause a step to have more than 2-3 actions and 4 sentences per step, please split into separate steps or substeps.

We believe we have made these needed edits. Please let us know if we have missed anything.

Figures, Tables, and Figure Legends:

1. Please remove the titles and Figure Legends from the uploaded figures. Please instead include all the Figure Legends together at the end of the Representative Results in the manuscript.

We have made these needed edits.

Reviewers' comments:

Reviewer #1:

Manuscript Summary:

In this methods article, Hu and Deer provide detailed protocols for the purification of the microtubule (MT)-associated motor proteins, kinesin and cytoplasmic dynein, and the generation of polymerized MTs and DNA origami. I believe the provided protocols are sufficient to reproduce the discussed experiments. The presented procedures are straightforward and relatively easy to implement in a lab. The protocol is written with care and I feel that it will be a timely and well-cited contribution. In conclusion, I believe the methods article is appropriate for JoVE and would recommend publication, subject to the following minor changes:

Minor Concerns:

1. The authors write in line 47 that "By converting the chemical energy of ATP into productive work and force generation, [...]". As work is the product of force times the distance over which the force is applied, I recommend to write either "By converting the chemical energy of ATP hydrolysis into productive work, [...]" or "By converting the chemical energy of ATP hydrolysis into force generation, [...]".

Thank you for catching this imprecise language. We have made the suggested edit.

2. In lines 623 to 626, the authors write that "The kymographs (Figure 6) of flexible chassis conjugated to seven dynein proteins ("7D" ensembles) show highly processive runs at relatively consistent and high velocities, which demonstrate the coordination and cooperation of multiple dynein motors in transporting one cargo." This statement is not correct. As the authors have shown in their Science paper (Fig. 2B in Derr et al. Science 2012), the velocity of the chassis with the 7 dynein motors is significantly lower than the velocity of a single dynein. This result suggests that yeast dynein motors, which are highly processive on their own, are not coordinating with each other and instead interfere. That the processivity of the multi-motor chassis is increased compared to the processivity of a single dynein, can be easily explained by the reduced likelihood of the chassis to detach from the MT as a result of the multiple MT-chassis dynein linkages. I recommend that the authors modify this sentence.

We agree and thank you for catching this. We have corrected it.

Reviewer #2:

Manuscript Summary:

This manuscript describes the procedure for assembling a complex hybrid nanostructure, consisting of a DNA-origami 'chassis' and purified molecular motors kinesin and dynein, and undertaking microscopy studies of the resulting nanostructures to elucidate motor mechanisms. It covers all the necessary steps, from yeast motor protein expression and purification, DNA origami assembly, loading of proteins onto the origami, and TIRF imaging.

This is an extremely comprehensive range of skills, in a very interdisciplinary project area, and so a JOVE video would be highly useful in broadening the applicability of this method. The efficacy of the protocol is well demonstrated, with reference to a number of recent high-impact publications. This protocol will be useful to scientists working in both motor protein and DNA origami areas. It is unlikely that many groups have sufficient background in both areas to develop such a protocol independently, so this video would act to bridge the expertise gap between the two areas. As many applications of DNA origami involve protein functionalization, this video would also serve as a useful example of such a protocol that could be adapted for other proteins.

Generally the discussion of background context and merits of the technique are unbiased and suitable. However, the protocol contains a large range of skills, and there is a high likelihood of scientists from diverse expertise backgrounds undertaking the protocol. Thus, some additional background discussion is required, as specified below. There are also a number of specific points that should be clarified.

Major Concerns:

None

Minor Concerns:

General discussion:

Abstract, Line 21: 'mechanisms the lead to emergent motile behavior'

Add a more specific example of the types of property that could be measured, eg 'how run length and velocity can vary with motor number'.

Thank you for this suggestion. We have added a phrase about motile properties to be measured.

Introduction, Line 78: 'segmented DNA origami with tunable compliance'

This sections needs more discussion on different DNA origami chassis/cargo designs and properties. For example, some examples of the different shapes, and a brief discussion of why

one would be chosen over another. As written it would not be clear to someone who is not already familiar with the referenced DNA origami papers.

We have added several new sentences to illustrate these ideas more fully.

5.1, Line 405: 'The compliance of the chassis structure is determined by the presence or absence of "linker" staples'

Similar to above comment, term 'Linkers' isn't very obvious in this context without looking up reference 7. Could describe it more specifically here, eg 'alternating regions of rigid DNA-origami helix bundles and regions of flexible ss or dsDNA, which act as linkers'

Thank you for noticing this. We have more fully explained these ideas as well.

1.1, Line 92: 'streak desired frozen yeast strain'

Additional detail needs to be added either here or before this on the yeast strain. This section would not be clear to those without detailed protein production experience, and could use some more references. For example, how is the strain selected? Does the yeast strain already contain the construct of interest? If so, where was this obtained from?

We have more thoroughly explained about the yeast strain and referenced a Methods paper that details how to create modified yeast proteins for single molecule assays.

4, Line 336: 'Microtubule (MT) affinity purification'

Would be good to give some background rationale for this step. Why do only non-functional motors bind to MTs under these conditions? 1-2 sentences would be suitable. Are there other potential affinity purifications that could be helpful?

Thank you for the advice. We have elaborated on the background rationale for this purification method. (We are unaware of other affinity purification methods that could be helpful at this step, so have not elaborated on that point.)

6.1.3, Line 504: 'A high absolute and relative concentration of motors enables saturation of the chassis' motor binding sites.'

Accurate control of motor stoichiometry is one of the key benefits of using DNA origami as a chassis for these experiments. Given this, it is necessary to give further discussion on how accurate binding site occupancy is at this point in the method, and in the background section. How are motor and chassis concentration estimated? What are appropriate ranges equivalent to 'high' and 'low' for good preparations? What is expected occupancy of motor binding sites, and how robust are these values?

We have added more details to address these concerns. We have added more details in sections 4.2.5 on PAGE analysis of motors and in section 5.2.10 on quantification of chassis concentration. We have also discussed the occupancy in section 6.1.3 as suggested.

6.4.3, line 557: 'generating kymographs from the chassis movies.'

Should provide suggestion of software or methods for this. Could provide references for methods papers covering this. This is mentioned in Figure caption for Figure 6, but should be mentioned earlier.

We have added more details on kymograph generation and analysis in Protocol section 6.4.3, suggested ImageJ as the software, and cited a method paper.

Specific points to clarify:

1.4, Line 103: 'When the culture is between an OD 1.5 and 2'

Can you suggest a time frequency for checking OD, eg. Every hour? For those not familiar with yeast culture, is shape of growth curve known? Is it better to under or overshoot the suggested OD range? Do you know how sensitive yield is to OD at this point?

We do not have data on all of these points ,but we have some more information here as suggested so as to facilitate adoption of this protocol by others.

2.1.4, Line 141: 'add PMSF to the buffer'

What would be an appropriate time between addition of PMSF to lysis buffer and addition of lysis buffer to yeast powder? 1 min? <5 min? Should clarify what time frame is meant by 'immediately' here.

We have clarified the time frame here.

2.1.5, Line 142: 'pre-chilled glass beaker on ice'

What is appropriate size for beaker?

We have included the beaker size that we typically use in our experiment.

2.1.7 Line 148: 'estimate the lysate volume, and add more 4x lysis buffer'

Add an example of typical/expected volumes for this point

We have added the expected volumes for this step.

2.1.8 Lines 153-156: 'The bottles must be carefully balanced in the rotor'

Balancing of centrifuge is key step, can you reference another JOVE video that shows this step if it is not shown directly in this video?

This will be in our video.

2.2.3, Line 171: '4x buffer prepared previously'

Does this require fresh PMSF?

We have clarified this point as suggested.

2.3.1, Line 193: '100 l of 1x TEV buffer containing ~15 uM of the purified BG-oligo'

You gives some indication of reaction yield dependence on reaction time, can you give some information on dependence on oligo concentration? How sensitive is it to oligo concentration, would be better to specify, eg. 10-20 uM instead of ~15 uM.

We agree 10-20 uM is a better instruction for this step and have changed it accordingly. We have added a comment on concentration as well.

2.3.2, Line 299: 'do not mix by pipetting'

Add a brief dicussion of why this is important, does pipetting/vortexing disrupt MT nucleation?

We have added this discussion.

5.1.1, Line 408: 'at 250 M concentration in Tris buffer'

Plate oligos can also be ordered dry and normalised, then resuspended, this is a more commonly available shipping option for many locations.

We have amended this step as suggested.

5.1.5, Line 418: '600 nM sore staple pool'

Typo, should be 'core' not 'sore'

This typo has been fixed.

5.1.7, Line 430 & 5.2.7, Line 456: 2% agarose gel in 0.5x TBE buffer supplemented with 11 mM MgCl₂'

Should specificy that gel running buffer also requires 11 mM MgCl₂.

Thank you for noticing these omissions. We have made the suggested edits.

5.2.10, Line 463: 'Concentrations of selected fractions can be quantified using appropriate spectroscopic methods.'

Give some examples here, eg. UV absorption at 260 nm.

We have added this suggestion.