

Journal of Visualized Experiments

Identification of host pathways targeted by bacterial effector proteins using yeast toxicity and suppressor screens --Manuscript Draft--

Article Type:	Invited Methods Article - JoVE Produced Video
Manuscript Number:	JoVE60488R1
Full Title:	Identification of host pathways targeted by bacterial effector proteins using yeast toxicity and suppressor screens
Section/Category:	JoVE Immunology and Infection
Keywords:	Bacterial effector protein, yeast suppressor, effector protein function, yeast screens, Chlamydia trachomatis, type III secreted effector protein, pathway identification
Corresponding Author:	Mary Weber UNITED STATES
Corresponding Author's Institution:	
Corresponding Author E-Mail:	mary-weber@uiowa.edu
Order of Authors:	Robert Faris Mary Weber
Additional Information:	
Question	Response
Please indicate whether this article will be Standard Access or Open Access.	Standard Access (US\$2,400)
Please indicate the city, state/province, and country where this article will be filmed . Please do not use abbreviations.	Iowa City



Mary M. Weber, Ph.D.
Assistant Professor
Department of Microbiology and Immunology
51 Newton Rd. BSB3-370
Iowa City, IA 52242
mary-weber@uiowa.edu

July 26th, 2019

Dear Editor,

We would like to submit a revised version of our manuscript (JoVE60488) "Identification of host pathways targeted by bacterial effector proteins using yeast toxicity and suppressor screens" for consideration by JoVE.

We have carefully reviewed the manuscript and have addressed all the concerns raised by the editor and peer-reviewers. We hope that our manuscript is found to be acceptable by JoVE, the editors, and peer-reviewers. If there is anything you need in relation to this resubmission, please don't hesitate to ask. We look forward to hearing back from you soon.

Sincerely,

A handwritten signature in black ink, appearing to read "Mary M. Weber". The signature is fluid and cursive, with a long horizontal stroke extending to the right.

Mary M. Weber, Ph.D.
Assistant Professor
Department of Microbiology and Immunology
University of Iowa

TITLE:

Identification of Host Pathways Targeted by Bacterial Effector Proteins using Yeast Toxicity and Suppressor Screens

AUTHORS AND AFFILIATIONS:

Robert Faris¹, Mary M. Weber¹

¹Department of Microbiology and Immunology, University of Iowa, Iowa City, IA, USA

Corresponding Author:

Mary M. Weber (mary-weber@uiowa.edu)

Email Addresses of Co-authors:

Robert Faris (robert-faris@uiowa.edu)

KEYWORDS:

bacterial effector protein, yeast suppressor, effector protein function, yeast screens, *Chlamydia trachomatis*, type III secreted effector protein, pathway identification

SUMMARY:

Bacterial pathogens secrete proteins into the host that target crucial biological processes. Identifying the host pathways targeted by bacterial effector proteins is key to addressing molecular pathogenesis. Here, a method using a modified yeast suppressor and toxicity screen to elucidate host pathways targeted by toxic bacterial effector proteins is described.

ABSTRACT:

Intracellular bacteria secrete virulence factors called effector proteins into the host cytosol that act to subvert host proteins and/or their associated biological pathways to the benefit of the bacterium. Identification of putative bacterial effector proteins has become more manageable due to advances in bacterial genome sequencing and the advent of algorithms that allow in silico identification of genes encoding secretion candidates and/or eukaryotic-like domains. However, identification of these important virulence factors is only an initial step. Naturally, the goal is to determine the molecular function of effector proteins and elucidate how they interact with the host. In recent years, techniques like the yeast two-hybrid screen and large-scale immunoprecipitations coupled with mass spectrometry have aided in the identification of protein-protein interactions. Although identification of a host binding partner is the crucial first step toward elucidating the molecular function of a bacterial effector protein, sometimes the host protein is found to have multiple biological functions (e.g., actin, clathrin, tubulin), or the bacterial protein may not physically bind host proteins, depriving the researcher of crucial information about the precise host pathway being manipulated. A modified yeast toxicity screen coupled with a suppressor screen has been adapted to identify host pathways impacted by bacterial effector proteins. The toxicity screen relies on a toxic effect in yeast caused by the effector protein interfering with the host biological pathways, which often manifests as a growth defect. Expression of a yeast genomic library is used to identify host factors that suppress the

toxicity of the bacterial effector protein and thus identify proteins in the pathway that the effector protein targets. This protocol contains detailed instructions for both the toxicity and suppressor screens. These techniques can be performed in any lab capable of molecular cloning and cultivation of yeast and *Escherichia coli*.

INTRODUCTION:

The first report of procedures similar to those presented here characterized the *Legionella pneumophila* type IV effector SidD, a deAMPylase that modifies Rab1¹. Comparable techniques were used for the characterization of several *L. pneumophila* effectors¹⁻³. The assay was adapted to characterize a *Coxiella burnetii* type IV effector protein⁴, and recently the utility of this technique was expanded for the characterization of *Chlamydia trachomatis* inclusion membrane proteins⁵.

This protocol can be broken into two major parts: 1) the yeast toxicity screen, in which the bacterial effector protein of interest is expressed in yeast and clones are screened for a toxic phenotype as evidenced by a growth defect, and 2) the yeast suppressor screen, in which the toxic phenotype is suppressed by expression of a yeast genomic library in the toxic strain. Thus, the toxicity screen is a screen for toxic phenotypes that manifest as growth defects when the bacterial effector of interest is overexpressed. Toxic clones, successfully transformed with and expressing the bacterial effector, are selected and saved for the next step. The second major step involves overexpressing a partially digested yeast genomic library in the toxic yeast clone. Plasmids making up the yeast genomic library suggested for the use in this protocol carry 5–20 kb inserts, usually corresponding to 3–13 yeast open reading frames (ORF) of an average gene size of ~1.5 kb across all plasmids, representing the whole yeast genome covered approximately 10x. This part of the assay is called the suppressor screen, as the goal is to suppress the toxicity of the bacterial effector protein. Potential suppressor plasmids are isolated from yeast, sequenced, and the suppressing ORFs identified. The rationale underlying the suppressor screen is that the effector protein binds, interacts with, and/or overwhelms components of the host pathway it targets, and that providing those host proteins back in excess can rescue the toxic effect on the pathway and thus, the growth defect. Thus, identified ORFs that suppress toxicity often represent multiple participants of a host pathway. Orthogonal experiments are then performed to verify that the bacterial effector indeed interacts with the implicated pathway. This is especially necessary if a binding partner such as clathrin or actin has been identified, because these proteins are involved in a multitude of host processes. Further experiments can then elucidate the physiological function of the effector protein during infection. The toxicity and suppressor screens are also powerful tools for deciphering the physiological function of bacterial effector proteins that do not physically bind host proteins with affinities sufficient to detect by immunoprecipitation or that interact with the host in enzymatic hit-and-run interactions that may not be detected by a yeast-two hybrid screen.

Although the suppressor screen can be a powerful method to reveal potential physiological interactions between bacterial effector proteins and host pathways, the bacterial effector protein must induce a growth defect in yeast, otherwise using it in the suppressor screen will be of little use. Furthermore, the toxic phenotype must result in at least a 2–3 log₁₀ deficit in growth

or it will be difficult to identify suppressors. If a laboratory is set up for cell culture, screening effector proteins for toxicity in common cell lines such as HeLa can often give insight as to whether it is worth the effort to proceed with the yeast toxicity screen. Ectopic expression of the effector protein in HeLa cells sometimes results in toxicity that correlates very strongly with toxicity in the yeast strain used for these screens⁴. Observable hallmarks of stress in HeLa cells include loss of stress fibers, cell detachment from the plate, and nuclear condensation indicating apoptosis. Any visual indication of stress in HeLa cells make the protein of interest a good candidate for inducing a growth defect in yeast, which replicate much more rapidly and are thus more responsive to perturbation of essential pathways.

It should be noted that the suppressor screen does not always identify host binding partners as suppressors, but it can still implicate critical components of the host pathway(s) targeted, yielding a holistic view of the biological processes being hijacked by the bacterial effector protein. On the surface, this seems counterintuitive, because providing the binding partner of the effector protein in excess would be expected to rescue the growth defect. In efforts to identify pathways targeted by the *C. trachomatis* effector protein CT229 (CpoS), which binds to at least 10 different Rab GTPases during infection⁵, none of the Rab binding partners suppressed the toxicity of CT229. However, numerous suppressors involved in clathrin-coated vesicle (CCV) trafficking were identified, which led to further work demonstrating that CT229 specifically subverts Rab-dependent CCV trafficking. Similarly, when investigating the *C. burnetii* effector protein Cbu0041 (CirA) several Rho GTPases that rescued the yeast growth defect were identified, and it was later found that CirA functions as a GTPase activating protein (GAP) for RhoA⁴.

The usefulness of the yeast suppressor screen for elucidating host pathways targeted by bacterial effector proteins cannot be overstated, and other researchers attempting to characterize intracellular bacterial effector proteins can greatly benefit from these techniques. These assays are of value if immunoprecipitations and/or yeast-two hybrid screens have failed to find a binding partner and can elucidate which pathways are targeted by the bacterial effector protein. Here, detailed protocols for the toxicity and suppressor screens to identify host biological pathways targeted by intracellular bacterial effector proteins are provided, as well as some of the common obstacles experienced when using these assays and their corresponding solutions.

PROTOCOL:

1. Preparation of media and reagents

NOTE: Plates should be prepared before the day of the assay and are good for 1 month. Media and reagents can be made at any point and are good for 1 month.

1.1. Prepare 1 L of the glucose solution (10% w/v) by dissolving 100 g of D-(+)-glucose in 800 mL of distilled water in a 1,000 mL beaker. Adjust the volume to 1 L with distilled water. Filter through a 0.2 µm sterile filter into a sterile 1 L media storage bottle.

1.2. Prepare 1 L of the galactose solution (10% w/v) by dissolving 100 g of D-(+)-galactose in 800

mL of distilled water in a 1 L beaker. Adjust the volume to 1 mL using distilled water and filter through a 0.2 µm sterile filter into a sterile 1,000 mL media storage bottle.

1.3. Prepare 1 L of yeast extract peptone dextrose (YPD) agar by dissolving 10 g of yeast extract, 20 g of peptone, 20 g of glucose, and 20 g of agar in 1,000 mL of distilled water. Autoclave for 20 min and cool in 56 °C water bath until cooled. Pour into 100 mm plates using ~20 mL of media per plate.

1.4. Prepare 1 L of YPD broth by dissolving 10 g of yeast extract, 20 g of peptone, and 20 g of glucose in 1,000 mL of distilled water. Autoclave for 20 min.

1.5. Prepare the synthetic dropout (SD) uracil (Ura⁻) glucose agar. For 1 L, dissolve 6.7 g of yeast nitrogen base without amino acids, 1.9 g of dropout supplement without uracil, and 15 g of agar in 800 mL of distilled water. Autoclave for 20 min and cool in a 56 °C water bath until the temperature is ~50–60 °C. Using a 50 mL serological pipette or a sterile graduated cylinder, add 200 mL of the glucose solution prepared in step 1.1. Mix well by gently swirling or place on stir plate. Pour into 100 mm plates using ~20 mL of media per plate.

1.6. Prepare the synthetic dropout (SD) uracil (Ura⁻) galactose agar. For 1 L, dissolve 6.7 g of yeast nitrogen base without amino acids, 1.9 g of dropout supplement without uracil, and 15 g of agar in 800 mL of distilled water. Autoclave for 20 min and cool in a 56 °C water bath until the temperature is ~50–60 °C. Using a 50 mL serological pipette or a sterile graduated cylinder, add 200 mL of the galactose solution prepared in step 1.2. Mix well by gently swirling or place on the stir plate. Pour into 100 mm plates using ~20 mL of media per plate.

1.7. Prepare the synthetic dropout (SD) uracil (Ura⁻) glucose broth. For 1 L, dissolve 6.7 g of yeast nitrogen base without amino acids and 1.9 g of dropout supplement without uracil in 800 mL of distilled water. Autoclave for 20 min and cool in a 56 °C water bath until the temperature is ~50–60 °C. Using a 50 mL serological pipette or a sterile graduated cylinder, add 200 ml of the glucose solution prepared in step 1.1.

1.8. Prepare the synthetic dropout (SD) uracil (Ura⁻) leucine (Leu⁻) glucose agar. For 1 L, dissolve 6.7 g of the yeast nitrogen base without amino acids; 1.9 g of dropout supplement without uracil, leucine, and tryptophan; 15 g of agar; and 0.076 g of tryptophan in 800 mL of distilled water. Autoclave for 20 min and cool in a 56 °C water bath until the temperature is ~50–60 °C. Using a 50 mL serological pipette or a sterile graduated cylinder, add 200 mL of the glucose solution prepared in step 1.1. Pour into 100 mm plates using ~20 mL of media per plate.

1.9. Prepare synthetic dropout (SD) uracil (Ura⁻) leucine (Leu⁻) galactose agar. For 1 L, dissolve 6.7 g of yeast nitrogen base without amino acids; 1.9 g of dropout supplement without uracil, leucine, and tryptophan; 15 g of agar; and 0.076 g of tryptophan in 800 mL of distilled water. Autoclave for 20 min and cool in a 56 °C water bath until the temperature is ~50–60 °C. Using a 50 mL serological pipette or a sterile graduated cylinder, add 200 mL of the galactose solution prepared in step 1.2. Mix well by gently swirling or placing on a stir plate. Pour into 100 mm plates

using ~20 mL of media per plate.

1.10. Prepare synthetic dropout (SD) uracil (Ura⁻) leucine (Leu⁻) glucose broth. For 1 L, dissolve 6.7 g of yeast nitrogen base without amino acids; 1.9 g of dropout supplement without uracil, leucine, tryptophan; and 0.076 g of tryptophan in 800 mL of distilled water. Autoclave for 20 min and cool in a 56–60 °C water bath until the temperature is ~50 °C. Using a 50 mL serological pipette or a sterile graduated cylinder, add 200 mL of the glucose solution prepared in step 1.1.

1.11. Prepare polyethylene glycol (PEG) solution. Add 50% w/v of polyethylene glycol 3350 in distilled water. Sterilize by autoclaving.

1.12. Prepare 1 M lithium acetate (LiAc) by dissolving 10.2 g of lithium acetate dehydrate in 200 mL of distilled water. Sterilize by autoclaving.

1.13. Prepare herring sperm DNA. Dilute 10 mg/mL herring sperm DNA to 2 mg/mL using distilled water. Heat at 100 °C for 5 min and immediately place on ice for 5 min.

2. Cloning the gene of interest into the yeast toxicity plasmid pYesNTA-Kan

NOTE: Currently, there are a variety yeast toxicity vectors available both commercially and academically. The yeast suppressor screen can be used in conjunction with many of these vectors provided the plasmid expressing the effector protein of interest uses a dropout selection other than uracil and an antibiotic marker other than Bla^R. This work used a modified pYesNTA vector^{4,5} that includes a kanamycin resistance cassette for easier screening for potential suppressors. The cloning scheme must allow for the effector protein to be in-frame with the Gal promoter and His-Tag. The *Chlamydia trachomatis* inclusion membrane protein CT229 was used as a proof of principle for these assays.

2.1. Use PCR to amplify CT229 from *C. trachomatis* genomic DNA following the manufacturer's instructions using CT229 +1 Kpn F (CCGGTACCAATGAGCTGTTCTAATGTTAATTCAGGT) and CT229 XhoI R (CCCTCGAGTTTTTTACGACGGGATGCC) primers. Use the following PCR conditions: (1) 98 °C for 30 s; (2) 98 °C for 10 s, 55 °C for 30 s, 72 °C for 2 min; (3) repeat step 2 for a total of 25x; (4) 72 °C for 10 min; and (5) 4 °C hold.

2.2. Analyze 5 µL of the PCR product on a 1% agarose gel (**Figure 1**).

2.3. Purify the remaining 45 µL of DNA using a PCR purification kit following the manufacturer's instructions.

2.4. Digest the 50 µL of purified PCR insert and pYesNTA-Kan (**Figure 2**) for 1 h at 37 °C in a water bath using KpnI-HF and XhoI-HF.

NOTE: For pYesNTA-Kan, digest 5 µg of plasmid using 5 µL of KpnI-HF, 5 µL of XhoI-HF, and 6 µL of buffer in a total volume of 60 µL. For the insert, digest the entire 50 µL of purified PCR product

using 2 μ L of KpnI-HF, 2 μ L of XhoI-HF, and 6 μ L of buffer.

2.5. Run the entire plasmid digest on a 1% agarose gel. Perform purification of the digested plasmid using a gel extraction kit following the manufacturer's instructions.

2.6. Purify the digested PCR insert using a PCR purification kit following the manufacturer's instructions.

2.7. Clone the insert into pYesNTA-Kan using 2 μ L of pYesNTA-Kan (step 2.5), 2 μ L of DNA ligase buffer, 1 μ L of T4 ligase, and 15 μ L of insert (step 2.6). Incubate at room temperature for 1 h.

2.8. Add the entire cloning reaction to 50 μ L of *E. coli* competent cells and perform the transformation reaction.

2.9. Plate the entire transformation reaction on an LB plate containing 100 μ g/mL of carbenicillin.

2.10. Inoculate 10 mL of the LB containing 100 μ g/mL of carbenicillin with three individual colonies. Incubate overnight at 37 °C with shaking at 150 rpm.

2.11. Isolate the plasmid using a plasmid miniprep kit following the manufacturer's instructions.

2.12. Sequence the isolated plasmids using gene-specific primers (step 2.1).

3. Test the protein of interest for toxicity in yeast

3.1 Streak *Saccharomyces cerevisiae* W303 on YPD agar (step 1.3) to obtain isolated colonies. Incubate at 30 °C for 24 h.

3.2 Inoculate 10 mL of YPD broth (step 1.4) with a single colony from the agar plate (step 3.1). Incubate at 30 °C with shaking at 150 rpm overnight.

3.3 Add 0.5 mL of the overnight culture to 10 mL of YPD broth (step 1.4) and incubate at 30 °C with shaking at 150 rpm for 4 h.

3.3.1 Pellet the 10 mL of culture at 3,000 x *g* for 10 min at 4 °C.

3.3.2 Resuspend the pellet in 1 mL of sterile water, transfer to the microcentrifuge tube, and pellet at 3,000 x *g* for 1 min at room temperature.

3.3.3 Resuspend the pellet in 1 mL of 1 mM lithium acetate (LiAc) and pellet at 3,000 x *g* for 1 min at room temperature.

3.3.4 Repeat the LiAc wash 2x.

3.3.5 Remove the wash. Resuspend the pellet in 2.4 mL of 50% PEG 3350. Add 360 μ L of 1M LiAc, 500 μ L of 2 mg/mL of herring sperm DNA, and 400 μ L of sterile water.

NOTE: This is sufficient for 20 transformations.

3.3.6 Add 180 μ L of the transformation mix from step 3.3.5 to a microcentrifuge tube containing 5 μ L of pYesNTA-Kan CT229 plasmid DNA (100–500 ng) from step 2.12.

3.3.7 Incubate in a water bath at 30 °C for 30 min.

3.3.8 Incubate in a water bath at 42 °C for 30 min.

3.3.9 Pellet at 3,000 x *g* for 1 min at room temperature.

3.3.10 Remove the transformation mix with the pipette. Resuspend the pellet in 100 μ L of sterile water. Plate the transformation on SD Ura⁻ agar with glucose (step 1.5).

3.3.11 Incubate plates at 30 °C for 48 h.

3.4 Inoculate 5 mL of SD Ura⁻ broth containing glucose (step 1.7) with a single colony from the plate. Incubate overnight at 30 °C with shaking at 150 rpm. Include yeast transformed with the vector alone as a negative control.

3.4.1 Add 180 μ L of sterile water to 5 wells of a 96 well plate (A2–A6).

3.4.2 Vortex the overnight culture to mix.

3.4.3 Add 180 μ L of yeast to the first well (A1). Serially dilute 1:10 (six samples in total including undiluted).

3.4.4 Using a multichannel pipette, spot 5 μ L of each dilution on SD Ura⁻ glucose (step 1.5) and Ura⁻ galactose (step 1.6) agar plates. Incubate at 30 °C for 48 h.

3.5 Assess the toxicity by comparing the growth of the yeast expressing the effector protein of interest grown on the galactose-containing media to the growth of yeast expressing the vector alone.

3.6 Confirm the expression of the His-tag fusion protein by Western blotting.

4. Transform toxic yeast with the yeast genomic library

4.1 Inoculate 100 mL of SD Ura⁻ glucose broth (step 1.7) using 1 mL of the stock from step 3.4. Incubate for 16–24 h at 30 °C with shaking at 150 rpm. Place 900 mL of SD Ura⁻ glucose broth at 30 °C overnight to prewarm the media.

4.2 Add the entire 100 mL of the overnight culture to the prewarmed 1 L flask. Incubate for 4–5 h at 30 °C with shaking at 150 rpm.

4.2.1 Pellet the culture at 6,000 x *g* for 10 min at 4 °C.

4.2.2 Discard the supernatant. Resuspend the pellet in 250 mL of sterile water. Pellet the culture at 6,000 x *g* for 5 min at 4 °C.

4.2.3 Discard the supernatant. Resuspend the pellet in 250 mL of 1 mM LiAc. Pellet the culture at 6,000 x *g* for 5 min at 4 °C.

4.2.4 Remove LiAc and resuspend the pellet in 9.6 mL of 50% PEG 3350. Add 1.44 mL of 1M LiAc, 2 mL of 2 mg/mL herring sperm DNA, 50 µg of the pYep13 genomic library (ATCC 37323). Adjust the volume to 15 mL with sterile water. Mix gently by inversion.

4.2.5 Incubate in a water bath at 30 °C for 30 min.

4.2.6 Add 750 µL of dimethyl sulfoxide (DMSO) to enhance the transformation efficiency. Incubate in a water bath at 42 °C for 30 min. Mix by gentle inversion every 10 min.

4.2.7 Pellet the yeast at 3,000 x *g* for 5 min at room temperature.

4.2.8 Discard the supernatant and resuspend the pellet in 10 mL of sterile water. Pellet at 3,000 x *g* for 5 min at room temperature.

4.2.9 Resuspend the pellet in 8 mL of sterile water.

4.2.10 To determine the transformation efficiency, dilute 1:10 and plate 100 µL of each dilution on SD Ura⁻ Leu⁻ glucose agar (step 1.8).

4.2.11 Plate 200 µL of the sample on SD Ura⁻ Leu⁻ galactose agar (step 1.9) plates. Use 50 plates in total.

4.2.12 Incubate at 30 °C for 48–96 h or until colonies appear.

NOTE: Colonies generally appear on glucose agar plates at 48–72 h and galactose agar plates at 72–96 h.

4.3 Patch colonies (potential rescues) on SD Ura⁻ Leu⁻ galactose agar (step 1.9) to expand and make stock. Incubate at 30 °C for 24–48h.

4.4 Inoculate 5 mL of SD Ura⁻ Leu⁻ glucose broth (step 1.10) using the part of the patch. Incubate overnight at 26 °C with shaking at 150 rpm.

4.4.1 Add 180 μ L of the sterile water to 5 wells of a 96 well plate (A2–A6).

4.4.2 Vortex the overnight culture to mix.

4.4.3 Add 180 μ L of yeast to the first well (A1). Serially dilute 1:10 (six samples in total including undiluted).

4.4.4 Using a multichannel pipette, spot 5 μ L of each dilution on SD Ura⁻ glucose and SD Ura⁻ galactose agar plates. Include toxic effector alone as a control.

4.4.5 Incubate at 30 °C for 48 h.

4.5 Compare the growth of the yeast expressing the toxic effector alone to yeast containing the potential suppressors. Only proceed with potential suppressors that have diminished toxicity compared to the yeast expressing the effector alone.

5. Identify and confirm suppressors

5.1 Inoculate 5 mL of SD Ura⁻ Leu⁻ glucose broth (step 1.10) with 100 μ L of yeast from step 4.4 and incubate overnight at 30 °C with shaking at 150 rpm.

5.1.1 Pellet the yeast at 3,000 $\times g$ for 2 min at room temperature. Gently discard the supernatant using a pipette.

5.1.2 Isolate the plasmid with a yeast plasmid miniprep kit following the manufacturer's instructions.

5.2 Transform the isolated plasmid into the *E. coli* competent cells and plate the entire transformation on LB agar with 100 μ g/mL of carbenicillin. Incubate at 37 °C for 24 h.

5.3 Inoculate 10 mL of LB broth containing 100 μ g/mL of carbenicillin, with three colonies from the plate and incubate at 37 °C overnight with shaking at 150 rpm.

5.3.1 Pellet cultures at 3,000 $\times g$ for 10 min at room temperature. Isolate plasmid using a miniprep kit following the manufacturer's instructions.

5.4 Retransform the toxic yeast with the isolated plasmids from step 5.3.1 following the steps in part 3 of the protocol.

5.4.1 Inoculate 5 mL of SD Ura⁻ Leu⁻ glucose broth with a colony from the transformation plate. Incubate overnight at 30 °C with shaking at 150 rpm.

5.4.2 Spot on Ura⁻ Leu⁻ glucose and galactose agar to confirm suppression of toxicity.

5.5 Sequence using pYep13 F (ACTACGCGATCATGGCGA) and pYep13 R (TGATGCCGCGCCACGATGC) primers to identify yeast ORFs.

REPRESENTATIVE RESULTS:

Before the actual yeast suppressor screen can be performed, the effector protein of interest must be tested for toxicity in yeast. This is accomplished by expressing the protein of interest in yeast under the control of a galactose-inducible promoter. Growth on glucose (noninducing conditions) should first be compared to ensure toxicity is specifically due to the expression of the protein of interest and is not a general defect. As shown in **Figure 3**, toxicity manifests as smaller colonies and/or reduced growth. A 2–3 log decrease in yeast growth is ideal for the yeast suppressor screens.

As shown in **Figure 4**, the toxicity and suppressor screen can be accomplished in eight steps using the methods described in the protocol section. The gene of interest is cloned into a suitable vector¹, and the resulting constructs are transformed into yeast to assess toxicity. This study used CT229 as the gene of interest and pYesNTA-Kan as the vector. Additionally, the expression of the His-tagged fusion protein was confirmed by Western blotting. Ideally, a 2–3 log reduction in yeast grown on galactose-containing media should be observed. If the protein of interest is toxic to yeast (**Figure 3** and **Figure 4**), the suppressor screen can be conducted by transforming the toxic strain with the yeast genomic library pYep13. Transformants are plated on glucose agar to determine transformation efficiency and galactose agar to identify potential suppressors. Using this protocol, a minimum transformation efficiency of 5×10^5 was achieved with a total of 10–250 colonies on galactose agar. Patching these colonies on SD Ura⁻ Leu⁻ galactose agar (**Figure 4**) aided in identification of true suppressors, because many false positives will not grow when patched. Here, 50 colonies were patched and eight clones that suppressed effector toxicity were obtained (**Figure 4**). Spotting suppressors on SD Ura⁻ Leu⁻ galactose agar was done to confirm suppression of toxicity. As shown in **Figure 4**, pSup1 and pSup 2 suppressed the toxicity of the effector protein, whereas pSup3 did not. Therefore, pSup3 was discarded. Plasmids were subsequently isolated from the suppressors and transformed into *E. coli* to increase the plasmid yield. The plasmid can then be retransformed into the toxic yeast to confirm that the isolated plasmid definitely suppresses effector toxicity (**Figure 4**). While most of the isolated plasmids should rescue toxicity, occasionally a plasmid that does not suppress toxicity when retransformed into the toxic yeast strain is obtained. Those plasmids are discarded, and only those that suppress toxicity following retransformation should be sequenced.

The isolated plasmid will contain multiple yeast ORFs⁴. To identify which is the true suppressor, each ORF can be individually cloned and expressed in the toxic yeast strain as previously described^{1–4}.

FIGURE AND TABLE LEGENDS:

Figure 1: PCR amplification of target gene. CT229 from *Chlamydia trachomatis* serovar L2 was PCR amplified from genomic DNA. PCR products were analyzed on a 1% agarose gel stained with

ethidium bromide.

Figure 2: pYesNTA-Kan plasmid map. The target gene of interest was cloned into the multiple cloning site (MCS) of pYesNTA-Kan. Kanamycin was used for selection in *E. coli* and uracil dropout was used for selection in *S. cerevisiae*. Expression of the His-tagged fusion protein was verified by Western blotting.

Figure 3: Representative yeast toxicity results. Effector proteins of interest were cloned into pYesNTA-Kan. Sequence-verified plasmids were transformed into *S. cerevisiae* and transformants were serially diluted and spotted on Ura⁻ glucose and galactose agar to assess toxicity.

Figure 4: Overview of the yeast suppressor screen. The target gene of interest was cloned into pYesNTA-Kan and plasmids were transformed into *S. cerevisiae*. Transformants were serially diluted and spotted on Ura⁻ glucose and galactose agar to assess toxicity. To identify the pathway targeted by the toxic effector protein, the yeast strain expressing the toxic protein was transformed with a yeast genomic library. Colonies were patched on Ura⁻ Leu⁻ galactose agar and spotted to assess toxicity. Plasmids were isolated from those that suppress the toxicity of the effector protein. Plasmids were retransformed into the toxic yeast strain to confirm that the isolated plasmid is a suppressor. Plasmids from suppressors were sequenced to identify the yeast ORFs present.

DISCUSSION:

This protocol outlines step-by-step procedures for identifying host biological pathways targeted by bacterial effector proteins using a modified yeast toxicity and suppressor screen. The yeast strain used, *S. cerevisiae* W303, is auxotrophic for both uracil and leucine. Uracil auxotrophy of the strain is used to select yeast carrying the protein of interest on the pYesNTA-Kan vector while leucine auxotrophy is used to select for the yeast genomic library vector pYep13. The yeast genomic library plasmids carry 3–13 ORFs, so each ORF should be individually cloned into the p415-ADH vector⁶ or a similar vector and retransformed into the toxic yeast to identify which is the true suppressor. It is typical to identify several suppressors representing a host biological pathway. Suppressors can sometimes be direct binding partners of the bacterial protein of interest, but not always.

Generating yeast clones carrying the bacterial protein of interest on pYesNTA-Kan is the first major step and great care should be taken to preserve the integrity of these clones as soon as possible because there is strong negative selection for yeast carrying toxic proteins, even without intentional galactose induction, and they may lose the plasmid. Occasionally there is plasmid loss even when the yeast is maintained under selection on dropout media. A method has been previously reported to decrease the occurrence of plasmid loss by linearizing the vector and integrating it within the yeast genome¹.

A major drawback of these techniques is that the suppressor screen only yields meaningful data if overexpression of the bacterial effector protein triggers an observable and consistent growth defect in yeast. The primary difficulty in this method is when no suppressors are identified. It is

hard to determine whether a failure to identify suppressors is due to biological reasons or technical issues. From a biological perspective, there could be several causes: the protein of interest could be so toxic that proteins from the genomic library are not able to significantly suppress toxicity, the targeted pathway may be incapable of rescue, or numerous co-expressed factors may be required to overcome the toxic effect. Approaching the biological explanation is challenging, as there are many undefined variables to dissect. Because there is an established protocol, exploring technical explanations is a much more tractable approach. From a technical standpoint, a failure to identify suppressors is likely the result of low transformation efficiency of the yeast genomic library. This protocol uses the ATCC *S. cerevisiae* AB320 partially digested genomic library in pYEp13. Other libraries may be sufficient as well as other vectors, but this protocol specifically uses the ATCC library and there are no studies currently published that report alternate library or vector sources. It is strongly recommended that coverage of the yeast genome should be at least 10x, or underrepresentation may hinder the identification of suppressors. Decreased transformation efficiency can bring these numbers under 1x. Thus, parts of the yeast genome may not be represented in the suppressor screen. Low transformation efficiency can be due to many issues, including bad reagents or poor library preparation. This method largely adheres to the DMSO/LiAc transformation method first put forth by Hill et al.⁷. It is recommended that new DMSO, 50% PEG, and fresh herring sperm DNA are purchased and prepared for the transformations and used exclusively for these assays to decrease contamination or degradation of the quality of the reagents. DMSO is highly hygroscopic and absorbs water from the atmosphere readily. Therefore, making many individual use aliquots is recommended to avoid repeated opening of a container's lid and exposure to the atmosphere. High quality herring sperm DNA can be stored at -20 °C and remains fit for use for several years, although a fresh preparation should be prepared from the stock for each transformation. Once boiled and cooled, the herring sperm DNA should not be reboiled or reused, so it is best to only remove what is needed and avoid making more than is necessary for the current transformation. Attention to detail can greatly improve the transformation efficiency of yeast and improve the chances of identifying a suppressor.

Considering that the yeast pYep13 vector carries 3–13 different yeast ORFs, reducing the number of false-positive plasmids is paramount to identifying suppressors. Yeast can take up and harbor multiple copies of similar vectors, whereas the general wisdom associated with plasmid incompatibility traditionally dictates that *E. coli* only maintains one type of plasmid with the same origin of replication (ORI) at a time, although this is not always the case. This can become a major issue when attempting to identify suppressor candidates. If suppressor yeast clones are identified, the plasmids will be extracted and retransformed into *E. coli* for propagation and subsequent sequencing for identification. Once the candidate suppressor plasmid has been isolated and retransformed into *E. coli*, it is strongly recommended that at least five *E. coli* colonies are picked for sequencing, because it is possible that individual clones may be carrying different plasmids. Since *E. coli* should only be able to propagate one plasmid at a time, if the yeast isolation yields multiple plasmids, different *E. coli* clones on the plate should carry only one of them at a time. Therefore picking five colonies based on normal transformation efficiency in *E. coli* should show whether the clones are carrying the same or different plasmids. If multiple plasmids arise subsequent to sequencing of the plasmids from *E. coli*, then each one should be

transformed into the toxic yeast clone and assessed for rescue of the growth defect. Furthermore, even if only a single suppressor plasmid is identified, for rigor and reproducibility retransformation of the plasmid into the toxic clone is encouraged.

The yeast suppressor assay is a very large experiment and appropriate time and materials should be planned accordingly. This protocol contains many details that must be adhered to strictly and it is recommended that the researcher conducting these assays carefully read and watch the protocol in full prior to starting an assay. This is a multipartite procedure including both bacteria and yeast that need to be cultivated at different temperatures. Isolation of plasmids from yeast can sometimes be difficult due to their tough cell wall and the use of DMSO helps achieve higher transformation efficiency.

One of the largest advantages of the yeast toxicity and suppressor screen is that the bacterial effector does not need to physically bind to the host substrates for any appreciable amount of time to detect the putative pathway targeted. Only a small handful of studies have reported using this or a similar technique¹⁻⁵. However, there are reports of techniques that identify host pathways and many techniques to screen for binding partners of bacterial effector proteins. Kramer et al. reported a novel systems biology approach to identify pathways targeted by bacterial effector proteins⁸. The technique used a haploid deletion strain collection of *S. cerevisiae* to screen for mutants sensitive to the *Shigella* effector Osp4. Using this technique, Osp4 was linked to regulation of MAPK signaling and attenuation of innate immunity in host cells. Array-based, high-throughput, automated yeast-two hybrid screens can now rapidly identify putative bacterial effector binding partners by screening many at the same time against millions of host proteins⁹. Similarly, high-throughput immunoprecipitations coupled with mass spectrometry are being used to mass identify putative binding partners of entire families of bacterial effector proteins¹⁰. While many of these studies can accelerate the timetable of effector protein characterization, they rarely provide insight as to what physiological processes are being manipulated and what the holistic outcome of these interactions are. Thus, techniques such as the yeast toxicity and suppressor screen are essential for illuminating the cellular processes being manipulated by these effector proteins. At some point, observations identifying binding partners or pathways targeted by bacterial effector proteins need to return to cellular infection models to determine whether these in vitro findings hold true in physiologically relevant scenarios. Genetic manipulation of obligate intracellular pathogens like *C. trachomatis* has been a major hurdle until recently¹¹⁻¹⁵ and generation of site-specific mutants to study effector proteins was not possible. In the last few years though, a revolution in generating chlamydial mutants, complementing them, and overexpressing target proteins has occurred. These recent advances have greatly enhanced the value of the information derived from studies using toxicity and suppressor screens as it is now possible to go from yeast directly into infection to probe the validity of the results.

ACKNOWLEDGMENTS:

We thank Shelby Andersen, Abby McCullough, and Laurel Woods for their assistance with these techniques. This study was funded by startup funds from the University of Iowa Department of Microbiology and Immunology to Mary M. Weber.

DISCLOSURES:

The authors declare that they have no competing financial interests.

REFERENCES:

1. Tan, Y., Luo, Z. Q. Legionella pneumophila SidD is a deAMPyase that modifies Rab1. *Nature*. **475** (7357), 506–509 (2011).
2. Guo, Z., Stephenson, R., Qiu, J., Zheng, S., Luo, Z. A Legionella effector modulates host cytoskeletal structure by inhibiting actin polymerization. *Microbes and Infection*. **16** (3), 225–236 (2014).
3. Tan, Y., Arnold, R. J., Luo, Z.-Q. Legionella pneumophila regulates the small GTPase Rab1 activity by reversible phosphorylation. *Proceedings of the National Academy of Sciences of the United States of America*. **108** (52), 21212–7 (2011).
4. Weber, M. M. et al. The type IV secreted effector protein CirA stimulates the GTPase activity of RhoA and is required for virulence in a mouse model of *Coxiella burnetii* infection. *Infection and Immunity*. **84** (June), IAI.01554-15 (2016).
5. Faris, R. et al. Chlamydia trachomatis CT229 subverts Rab GTPase-dependent CCV trafficking pathways to promote chlamydial infection. *Cell Reports*. **26**, 3380–3390 (2019).
6. Mumberg, D., Müller, R., Funk, M. Yeast vectors for the controlled expression of heterologous proteins in different genetic backgrounds. *Gene*. **156**, 119–122 (1995).
7. Hill, J., Donald, K. A. G., Griffiths, D. E. DMSO-enhanced whole cell yeast transformation. *Nucleic Acids Research*. **19** (20), 41070 (1991).
8. Kramer, R. W. et al. Yeast Functional Genomic Screens Lead to Identification of a Role for a Bacterial Effector in Innate Immunity Regulation. *PLoS Pathogens*. **3** (2), e21. (2007).
9. Häuser, R. T. S., Rajagopala, S. V., Uetz, P. Array-Based Yeast Two-Hybrid Screens: A Practical Guide. *Two Hybrid Technologies: Methods and Protocols, Methods in Molecular Biology*. 21–38 (2012).
10. Mirrashidi, K. M. et al. Global mapping of the inc-human interactome reveals that retromer restricts chlamydia infection. *Cell Host and Microbe*. **18** (1), 109–121 (2015).
11. Wang, Y. et al. Development of a transformation system for chlamydia trachomatis: Restoration of glycogen biosynthesis by acquisition of a plasmid shuttle vector. *PLoS Pathogens*. **7** (9), e1002258 (2011).
12. Mueller, K. E., Wolf, K., Fields, K. A. Gene deletion by fluorescence-reported allelic exchange mutagenesis in Chlamydia trachomatis. *mBio*. **7** (1), 1–9 (2016).
13. Johnson, C. M., Fisher, D. J. Site-Specific, Insertional Inactivation of incA in Chlamydia trachomatis Using a Group II Intron. *PLoS ONE*. **8** (12), e83989 (2013).
14. Weber, M. M. et al. Absence of specific Chlamydia trachomatis inclusion membrane proteins triggers premature inclusion membrane lysis and host cell death. *Cell Reports*. **19** (7), 1406–1417, 04.058 (2017).
15. Weber, M. M. et al. A functional core of IncA is required for Chlamydia trachomatis inclusion fusion. *Journal of Bacteriology*. **198** (8), 1347–1355 (2016).







Figure 2







Name of Material/ Equipment	Company	Catalog Number
Agar	Fisher Scientific	BP2641500
Galactose	MilliporeSigma	G0750-1KG
GeneJet Gel extraction kit	ThermoFisher Scientific	K0691
GeneJet PCR purification kit	ThermoFisher Scientific	K0701
GeneJet plasmid miniprep kit	Thermo	K0503
Glucose	MilliporeSigma	G8270-1KG
Herring Sperm DNA	Promega	D1811
KpnI-HF	New England Biolabs	R3142S
Lithium acetate dihydrate	MilliporeSigma	L6883-250G
Peptone	Fisher Scientific	
Phusion High-Fidelity DNA Polymerase	New England Biolabs	M0530
Poly(ethylene glycol) 3350	MilliporeSigma	1546547-1G
pYep13	ATCC	37323
T4 DNA ligase	New England Biolabs	M0202S
Tryptophan	MilliporeSigma	470031-1G
XhoI-HF	New England Biolabs	R0146S
Yeast extract	Fisher Scientific	BP1422-500
Yeast miniprep kit	Zymo	D2001
Yeast nitrogen base without amino acids	MilliporeSigma	Y0626-250G
Yeast Synthetic Drop-out Medium Supplements	MilliporeSigma	Y1501-20G
Yeast Synthetic Drop-out Medium Supplements	MilliporeSigma	Y1771-20G

Comments/Description

without uracil

without uracil, leucine, tryptophan



1 Alewife Center #200
Cambridge, MA 02140
tel. 617.945.9051
www.jove.com

ARTICLE AND VIDEO LICENSE AGREEMENT

Title of Article:	Identification of host pathways targeted by bacterial effector proteins using yeast toxicity and suppressor screens
Author(s):	Robert Faris and Mary M. Weber

Item 1: The Author elects to have the Materials be made available (as described at <http://www.jove.com/publish>) via:

☒ Standard Access ☐ Open Access

Item 2: Please select one of the following items:

- ☒ The Author is **NOT** a United States government employee.
- ☐ The Author is a United States government employee and the Materials were prepared in the course of his or her duties as a United States government employee.
- ☐ The Author is a United States government employee but the Materials were NOT prepared in the course of his or her duties as a United States government employee.

ARTICLE AND VIDEO LICENSE AGREEMENT

1. **Defined Terms.** As used in this Article and Video License Agreement, the following terms shall have the following meanings: **"Agreement"** means this Article and Video License Agreement; **"Article"** means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein; **"Author"** means the author who is a signatory to this Agreement; **"Collective Work"** means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; **"CRC License"** means the Creative Commons Attribution-Non Commercial-No Derivs 3.0 Unported Agreement, the terms and conditions of which can be found at: <http://creativecommons.org/licenses/by-nc-nd/3.0/legalcode>; **"Derivative Work"** means a work based upon the Materials or upon the Materials and other pre-existing works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; **"Institution"** means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; **"JoVE"** means MyJoVE Corporation, a Massachusetts corporation and the publisher of The Journal of Visualized Experiments; **"Materials"** means the Article and / or the Video; **"Parties"** means the Author and JoVE; **"Video"** means any video(s) made by the Author, alone or in conjunction with any other parties, or by JoVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion

of the Article, and in which the Author may or may not appear.

2. **Background.** The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JoVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.

3. **Grant of Rights in Article.** In consideration of JoVE agreeing to publish the Article, the Author hereby grants to JoVE, subject to **Sections 4** and **7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the "Open Access" box has been checked in **Item 1** above, JoVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.

ARTICLE AND VIDEO LICENSE AGREEMENT

4. **Retention of Rights in Article.** Notwithstanding the exclusive license granted to JoVE in **Section 3** above, the Author shall, with respect to the Article, retain the non-exclusive right to use all or part of the Article for the non-commercial purpose of giving lectures, presentations or teaching classes, and to post a copy of the Article on the Institution's website or the Author's personal website, in each case provided that a link to the Article on the JoVE website is provided and notice of JoVE's copyright in the Article is included. All non-copyright intellectual property rights in and to the Article, such as patent rights, shall remain with the Author.

5. **Grant of Rights in Video – Standard Access.** This **Section 5** applies if the "Standard Access" box has been checked in **Item 1** above or if no box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby acknowledges and agrees that, Subject to **Section 7** below, JoVE is and shall be the sole and exclusive owner of all rights of any nature, including, without limitation, all copyrights, in and to the Video. To the extent that, by law, the Author is deemed, now or at any time in the future, to have any rights of any nature in or to the Video, the Author hereby disclaims all such rights and transfers all such rights to JoVE.

6. **Grant of Rights in Video – Open Access.** This **Section 6** applies only if the "Open Access" box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby grants to JoVE, subject to **Section 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Video in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Video into other languages, create adaptations, summaries or extracts of the Video or other Derivative Works or Collective Works based on all or any portion of the Video and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. For any Video to which this **Section 6** is applicable, JoVE and the Author hereby grant to the public all such rights in the Video as provided in, but subject to all limitations and requirements set forth in, the CRC License.

7. **Government Employees.** If the Author is a United States government employee and the Article was prepared in the course of his or her duties as a United States government employee, as indicated in **Item 2** above, and any of the licenses or grants granted by the Author hereunder exceed the scope of the 17 U.S.C. 403, then the rights granted hereunder shall be limited to the maximum

rights permitted under such statute. In such case, all provisions contained herein that are not in conflict with such statute shall remain in full force and effect, and all provisions contained herein that do so conflict shall be deemed to be amended so as to provide to JoVE the maximum rights permissible within such statute.

8. **Protection of the Work.** The Author(s) authorize JoVE to take steps in the Author(s) name and on their behalf if JoVE believes some third party could be infringing or might infringe the copyright of either the Author's Article and/or Video.

9. **Likeness, Privacy, Personality.** The Author hereby grants JoVE the right to use the Author's name, voice, likeness, picture, photograph, image, biography and performance in any way, commercial or otherwise, in connection with the Materials and the sale, promotion and distribution thereof. The Author hereby waives any and all rights he or she may have, relating to his or her appearance in the Video or otherwise relating to the Materials, under all applicable privacy, likeness, personality or similar laws.

10. **Author Warranties.** The Author represents and warrants that the Article is original, that it has not been published, that the copyright interest is owned by the Author (or, if more than one author is listed at the beginning of this Agreement, by such authors collectively) and has not been assigned, licensed, or otherwise transferred to any other party. The Author represents and warrants that the author(s) listed at the top of this Agreement are the only authors of the Materials. If more than one author is listed at the top of this Agreement and if any such author has not entered into a separate Article and Video License Agreement with JoVE relating to the Materials, the Author represents and warrants that the Author has been authorized by each of the other such authors to execute this Agreement on his or her behalf and to bind him or her with respect to the terms of this Agreement as if each of them had been a party hereto as an Author. The Author warrants that the use, reproduction, distribution, public or private performance or display, and/or modification of all or any portion of the Materials does not and will not violate, infringe and/or misappropriate the patent, trademark, intellectual property or other rights of any third party. The Author represents and warrants that it has and will continue to comply with all government, institutional and other regulations, including, without limitation all institutional, laboratory, hospital, ethical, human and animal treatment, privacy, and all other rules, regulations, laws, procedures or guidelines, applicable to the Materials, and that all research involving human and animal subjects has been approved by the Author's relevant institutional review board.

11. **JoVE Discretion.** If the Author requests the assistance of JoVE in producing the Video in the Author's facility, the Author shall ensure that the presence of JoVE employees, agents or independent contractors is in accordance with the relevant regulations of the Author's institution. If more than one author is listed at the beginning of this Agreement, JoVE may, in its sole

ARTICLE AND VIDEO LICENSE AGREEMENT

discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JoVE reserves the right, in its absolute and sole discretion and without giving any reason therefore, to accept or decline any work submitted to JoVE. JoVE and its employees, agents and independent contractors shall have full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JoVE has sole discretion as to the method of making and publishing the Materials, including, without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

12. **Indemnification.** The Author agrees to indemnify JoVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JoVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JoVE, making of videos by JoVE, or publication in JoVE or elsewhere by JoVE. The Author shall be responsible for, and shall hold JoVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to

the making of a video by JoVE its employees, agents or independent contractors. All sterilization, cleanliness or decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's expense. All indemnifications provided herein shall include JoVE's attorney's fees and costs related to said losses or damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JoVE, its employees, agents or independent contractors.

13. **Fees.** To cover the cost incurred for publication, JoVE must receive payment before production and publication the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JoVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.

14. **Transfer, Governing Law.** This Agreement may be assigned by JoVE and shall inure to the benefits of any of JoVE's successors and assignees. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to be one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement is required per submission.

CORRESPONDING AUTHOR

Name:

Mary M. Weber

Department:

Microbiology and Immunology

Institution:

The University of Iowa, Carver College of Medicine

Title:

Assistant Professor

Signature:

Mary M. Weber

Date:

06-28-19

Please submit a **signed** and **dated** copy of this license by one of the following three methods:

1. Upload an electronic version on the JoVE submission site
2. Fax the document to +1.866.381.2236
3. Mail the document to JoVE / Attn: JoVE Editorial / 1 Alewife Center #200 / Cambridge, MA 02140

1. Please take this opportunity to thoroughly proofread the manuscript to ensure that there are no spelling or grammar issues. The JoVE editor will not copy-edit your manuscript and any errors in the submitted revision may be present in the published version.

We have proofread the manuscript.

2. Please ensure the Introduction includes all of the following with citations:

- a) A clear statement of the overall goal of this method
- b) The rationale behind the development and/or use of this technique
- c) The advantages over alternative techniques with applicable references to previous studies
- d) A description of the context of the technique in the wider body of literature
- e) Information to help readers to determine whether the method is appropriate for their application

The introduction contains all this information.

3. Please remove all commercial language from your manuscript and use generic terms instead. All commercial products should be sufficiently referenced in the Table of Materials and Reagents. For example: Zymo Research yeast plasmid miniprep, PEG 3350, etc.

We have removed the Zymo and Thermo references from the protocol.

4. Please ensure that all text in the protocol section is written in the imperative tense as if telling someone how to do the technique (e.g., “Do this,” “Ensure that,” etc.). The actions should be described in the imperative tense in complete sentences wherever possible. Avoid usage of phrases such as “could be,” “should be,” and “would be” throughout the Protocol. Any text that cannot be written in the imperative tense may be added as a “Note.” However, notes should be concise and used sparingly. Please include all safety procedures and use of hoods, etc.

We have confirmed the protocol is written in this format.

5. The Protocol should contain only action items that direct the reader to do something.

We have confirmed the protocol is written in this format.

6. Please revise the protocol text to avoid the use of any personal pronouns in the protocol (e.g., “we”, “you”, “our” etc.).

We have confirmed the protocol is written in this format.

7. Please ensure that individual steps of the protocol should only contain 2-3 actions per step.

We have confirmed the protocol is written in this format.

8. Please ensure that you include all specific details associated with your protocol in a step by step manner. Please use complete sentences throughout.

We have confirmed the protocol is written in this format.

9. Please bring out the link between different subheadings of the protocol.

We have brought out the link between different subheadings in the protocol.

10. Please ensure you answer the “how” question, i.e., how is the step performed?

We have confirmed the protocol is written in this format.

11. 1: Please include what gene was amplified in your case, where was the nucleic acid obtained for amplification, what are the primers used, reaction set up and program used for amplification. Please also include a step by step protocol to show how the cloning was performed in the MCS of the plasmid used. Please include a short note stating why this plasmid was used. After cloning, did you perform restriction digestion, colony pcr, etc.

We have updated the protocol to include this information.

12. 2.3: Do you check the O.D in this case?

No

13. 2.3.6: What is the transformation mix- the mix obtained in the step 2.3.4? The plasmid DNA is from step 1? Please include this for clarity. Do you add the culture medium at this stage? Do you perform the incubations in the water bath?

Additional details have been added.

14. 3.1.2: Do you measure the O.D in this case?

No

15. 3.1.6: importance of pYep13?

pYep13 is the yeast genomic library. This information has been added to the protocol.

16. 3.1.7: Why is DMSO added at this stage?

The addition of DMSO significantly increases the transformation efficiency. This information has been added to the protocol.

17. 4.1.3: If this step needs filming, please describe the actions associated with it in brief.

This step does not need filming.

18. There is a 10-page limit for the Protocol, but there is a 2.75-page limit for filmable content. Please highlight 2.75 pages or less of the Protocol (including headings and spacing) that identifies the essential steps of the protocol for the video, i.e., the steps that should be visualized to tell the most cohesive story of the Protocol. Please adjust the highlights accordingly.

Highlights have been adjusted accordingly.

19. Please include representative results for all the sections in the protocol. For example, PCR results, positive clone screen is not included.

Representative results have been added.

20. Please obtain explicit copyright permission to reuse any figures from a previous publication. Explicit permission can be expressed in the form of a letter from the editor or a link to the editorial policy that allows reprints. Please upload this information as a .doc or .docx file to your Editorial Manager account. The Figure must be cited appropriately in the Figure Legend, i.e. "This figure has been modified from [citation]."

All figures are new, no copyright permission is needed.

21. As we are a methods journal, please revise the Discussion to explicitly cover the following in detail in 3-6 paragraphs with citations:

- a) Critical steps within the protocol
- b) Any modifications and troubleshooting of the technique
- c) Any limitations of the technique
- d) The significance with respect to existing methods
- e) Any future applications of the technique

We have confirmed the discussion contains this information and have shortened it to 6 paragraphs.

22. Please sort the materials table in alphabetical order.

The materials are now listed in alphabetical order.

Reviewers' comments:

Reviewer #1:

Manuscript Summary:

This manuscript provides a detailed explanation of how to use a yeast toxicity-suppressor screen to elucidate eukaryotic pathways that may be targeted by bacterial effector proteins. The protocol is timely as it has recently been used in several high profile publications that have reported function of novel effectors.

Major Concerns:

No major concerns

Minor Concerns:

* the protocol makes no mention of checking for expression of the effector protein in yeast although this is mentioned in the discussion. This should be included in the protocol. Authors mention a His-tag - is it routine to include a tag to detect the effector? N- or C- terminal? Incorporated into the oligonucleotide for amplification of the effector?

The His-tag is on the plasmid and thus is not incorporated into the oligos for effector amplification. We have incorporated that effector expression should be confirmed in the protocol.

* The protocol refers to "agar" and "broth" in vague ways. I think the yeast media being used should be more clearly defined (can a recipe be included?). I assume it is a Yeast minimal media with other amino acids supplemented?

We have added a section that includes the recipes for all the media and reagents.

* what are the oligonucleotides used to sequence recovered pYep13 clones?

The oligonucleotides used for sequencing have been added.

Reviewer #2:

Manuscript Summary:

This manuscript by Faris and Weber describes a protocol to identify bacterial effectors toxic to eukaryotic cells via yeast viability and subsequently identify yeast proteins/pathways targeted by the toxic effector using an over-expression suppressor screen. Yeast are powerful model organisms to investigate bacterial effector protein function and this protocol will be of immense value to the scientific community. Minor comments are listed below.

Major Concerns:

None.

Minor Concerns:

1. Line 126: Is pYesNTA-Kan available commercially or from the authors upon request? Are there alternative plasmids that can be used?

We have added a note in the protocol. The pYesNTA-Kan is available from our laboratory and other commercially available plasmids can be used for the screen.

2. Are the restriction sites in frame with the 6xHis tag? Please revise Figure 1 so that the MCS shows the 6xHis codons.

The figure has been revised to include the 6X His-tag codons.

3. Step 2.3.1 describes a centrifugation step at 4°C. Are subsequent incubations/centrifugation steps also done at the this temperature? No temperature is specified.

Temperature for all centrifugation steps have been added.

4. Lines 158-169: fix numbering for steps and highlighting.

The numbering has been corrected.

5. Line 210: do you mean the growth of yeast harboring the toxic effector?

Yes, we have clarified this in the text.

6. Line 219: typo: manufacturers

The type has been corrected.

7. The use of idioms such as "pipe-dream" and "foolhardy" may decrease readability for non-native English readers. Consider revising verbiage.

We have revised the verbiage.