

Journal of Visualized Experiments

Combining Augmented Reality and 3D Printing to Display Patient Models on Your Smartphone

--Manuscript Draft--

Article Type:	Invited Methods Article - JoVE Produced Video
Manuscript Number:	JoVE60618R1
Full Title:	Combining Augmented Reality and 3D Printing to Display Patient Models on Your Smartphone
Section/Category:	JoVE Engineering
Keywords:	Augmented reality; 3D printing; Anatomical models; Clinical applications; Surgical navigation; Image guidance; Patient-to-model registration; education; Pre-operative planning
Corresponding Author:	Javier Pascau, Ph.D. Universidad Carlos III de Madrid Leganés, Madrid SPAIN
Corresponding Author's Institution:	Universidad Carlos III de Madrid
Corresponding Author E-Mail:	jpascau@ing.uc3m.es
Order of Authors:	Rafael Moreta-Martinez David García-Mato Mónica García-Sevilla Rubén Pérez-Mañanes José A. Calvo-Haro Javier Pascau, Ph.D.
Additional Information:	
Question	Response
Please indicate whether this article will be Standard Access or Open Access.	Standard Access (US\$2,400)
Please indicate the city, state/province, and country where this article will be filmed . Please do not use abbreviations.	Leganés, Madrid, Spain

TITLE:

Combining Augmented Reality and 3D Printing to Display Patient Models on a Smartphone

AUTHORS AND AFFILIATIONS:

Rafael Moreta-Martinez^{1,2}, David García-Mato^{1,2}, Mónica García-Sevilla^{1,2}, Rubén Pérez-Mañanes^{2,3,4}, José A. Calvo-Haro^{2,3,4}, Javier Pascau^{1,2}

¹Department of Bioengineering and Aerospace Engineering, Carlos III University of Madrid, Madrid, Spain

²Gregorio Marañón Health Research Institute, Madrid, Spain

³Department of Orthopedic Surgery and Traumatology, General University Hospital Gregorio Marañón, Madrid, Spain

⁴Department of Surgery, Faculty of Medicine, Complutense University of Madrid, Madrid, Spain

Corresponding Author:

Javier Pascau (jpascau@ing.uc3m.es)

Email Addresses of Co-authors:

Rafael Moreta-Martinez (rmoreta@pa.uc3m.es)

David García-Mato (dgmato@ing.uc3m.es)

Mónica García-Sevilla (mongarci@pa.uc3m.es)

Rubén Pérez-Mañanes (rubenperez.phd@gmail.com)

José A. Calvo-Haro (calvoharo@yahoo.es)

KEYWORDS:

augmented reality, 3D printing, anatomical models, clinical applications, surgical navigation, image guidance, patient-to-model registration, education, pre-operative planning

SUMMARY:

Presented here is a method to design an augmented reality smartphone application for the visualization of anatomical three-dimensional models of patients using a 3D-printed reference marker.

ABSTRACT:

Augmented reality (AR) has great potential in education, training, and surgical guidance in the medical field. Its combination with three-dimensional (3D) printing (3DP) opens new possibilities in clinical applications. Although these technologies have grown exponentially in recent years, their adoption by physicians is still limited, since they require extensive knowledge of engineering and software development. Therefore, the purpose of this protocol is to describe a step-by-step methodology enabling inexperienced users to create a smartphone app, which combines AR and 3DP for the visualization of anatomical 3D models of patients with a 3D-printed reference marker. The protocol describes how to create 3D virtual models of a patient's anatomy derived from 3D medical images. It then explains how to perform positioning of the 3D models with respect to marker references. Also provided are instructions for how to

3D print the required tools and models. Finally, steps to deploy the app are provided. The protocol is based on free and multi-platform software and can be applied to any medical imaging modality or patient. An alternative approach is described to provide automatic registration between a 3D-printed model created from a patient's anatomy and the projected holograms. As an example, a clinical case of a patient suffering from distal leg sarcoma is provided to illustrate the methodology. It is expected that this protocol will accelerate the adoption of AR and 3DP technologies by medical professionals.

INTRODUCTION:

AR and 3DP are technologies that provide increasing numbers of applications in the medical field. In the case of AR, its interaction with virtual 3D models and the real environment benefits physicians in regards to education and training¹⁻³, communication and interactions with other physicians⁴, and guidance during clinical interventions⁵⁻¹⁰. Likewise, 3DP has become a powerful solution for physicians when developing patient-specific customizable tools¹¹⁻¹³ or creating 3D models of a patient's anatomy, which can help improve preoperative planning and clinical interventions^{14,15}.

Both AR and 3DP technologies help to improve orientation, guidance, and spatial skills in medical procedures; thus, their combination is the next logical step. Previous work has demonstrated that their joint use can increase value in patient education¹⁶, facilitating explanations of medical conditions and proposed treatment, optimizing surgical workflow^{17,18} and improving patient-to-model registration¹⁹. Although these technologies have grown exponentially in recent years, their adoption by physicians is still limited, since they require extensive knowledge of engineering and software development. Therefore, the purpose of this work is to describe a step-by-step methodology that enables the use of AR and 3DP by inexperienced users without the need for broad technical knowledge.

This protocol describes how to develop an AR smartphone app that allows the superimposing of any patient-based 3D model onto a real-world environment using a 3D-printed marker tracked by the smartphone camera. In addition, an alternative approach is described to provide automatic registration between a 3D-printed biomodel (i.e., a 3D model created from a patient's anatomy) and the projected holograms. The protocol described is entirely based on free and multi-platform software.

In previous work, AR patient-to-image registration has been calculated manually⁵ with surface recognition algorithms¹⁰ or has been unavailable². These methods have been considered somewhat limited when an accurate registration is required¹⁹. To overcome these limitations, this work provides tools to perform accurate and simple patient-to-image registration in AR procedures by combining AR technology and 3DP.

The protocol is generic and can be applied to any medical imaging modality or patient. As an example, a real clinical case of a patient suffering from distal leg sarcoma is provided to illustrate the methodology. The first step describes how to easily segment the affected

anatomy from computed tomography (CT) medical images to generate 3D virtual models. Afterwards, positioning of the 3D models is performed, then the required tools and models are 3D-printed. Finally, the desired AR app is deployed. This app allows for the visualization of patient 3D models overlaid on a smartphone camera in real-time.

PROTOCOL:

This study was performed in accordance with the principles of the 1964 Declaration of Helsinki as revised in 2013. The anonymized patient data and pictures included in this paper are used after written informed consent was obtained from the participant and/or their legal representative, in which he/she approved the use of this data for dissemination activities including scientific publications.

1. Workstation set-up for segmentation, 3D models extraction, positioning, and AR app deployment

NOTE: This protocol has been tested with the specific software version indicated for each tool. It is likely to work with newer versions, although it is not guaranteed.

1.1. Use a computer with Microsoft Windows 10 or Mac OS as operating systems.

1.2. Install the following tools from the corresponding websites per the official instructions:

3D Slicer (v. 4.10.2): <<https://download.slicer.org/>>.

Meshmixer (v. 3.5): <<http://www.meshmixer.com/download.html>>.

Unity (v. 2019): <<https://unity3d.com/get-unity/download>>.

(For iOS deployment only) Xcode (last version): <<https://developer.apple.com/xcode/>>.

NOTE: All software tools required for completing the protocol can be freely downloaded for personal purposes. Software to be used in each step will be specifically indicated.

1.3. Download data from the following GitHub repository, found at <<https://github.com/BIIG-UC3M/OpenARHealth>>.

NOTE: The repository contains the following folders:

“/3DSlicerModule/”: 3D slicer module for positioning 3D models with respect to the 3D-printed marker. Used in section 3. Add the module into the 3D slicer following the instructions available at <<https://github.com/BIIG-UC3M/OpenARHealth>>.

“/Data/PatientData/Patient000_CT.nrrd”: CT of a patient suffering from distal leg sarcoma. The protocol is described using this image as an example.

“/Data/Biomodels/”: 3D models of the patient (bone and tumor).

“/Data/Markers/”: Markers that will be 3D-printed, which will be detected by the AR application to position the virtual 3D models. There are two markers available.

2. Biomodel creation

NOTE: The goal of this section is to create 3D models of the patient’s anatomy. They will be obtained by applying segmentation methods to a medical image (here, using a CT image). The process consists of three different steps: 1) loading the patient data into 3D slicer software, 2), segmentation of target anatomy volumes, and 3) exportation of segmentation as 3D models in OBJ format. The resulting 3D models will be visualized in the final AR application.

2.1. Load patient data (“/Data/PatientData/Patient000_CT.nrrd”) by dragging the medical image file into the 3D slicer software window. Click **OK**. The CT views (axial, sagittal, coronal) will appear on the corresponding windows.

NOTE: The data used here is found in “nearly raw raster data” (NRRD) format, but 3D Slicer allows for loading of medical image format (DICOM) files. Go to the following link for further instructions, found at <<https://www.slicer.org/wiki/Documentation/4.10/Training>>.

2.2. To segment the anatomy of the patient, go to the **Segment Editor** module in 3D slicer.

2.2.1. A “segmentation” item is created automatically when entering the module. Select the desired volume (a medical image of the patient) in the **Master Volume** section. Then, right-click below on the **Add** button to create a segment. A new segment will be created with the name “Segment_1”.

2.2.2. There is a panel called **Effects** that contains a variety of tools to properly segment the target area of the medical image. Select the most convenient tool for the target and segment onto the image windows area.

2.2.2.1. To segment the bone (tibia and fibula in this case), use the **Threshold** tool to set up minimum and maximum HU values from the CT image, which corresponds to bone tissue. By using this tool, other elements with HU outside these threshold values are removed, such as soft tissue.

2.2.2.2. Use the **Scissors** tool to remove undesired areas, such as the bed or other anatomical structures, from the segmented mask. Segment the sarcoma manually using the **Draw** and **Erase** tools, since the tumor is difficult to contour with automatic tools.

NOTE: To learn more details about the segmentation procedure, go to the link found at <<https://www.slicer.org/wiki/Documentation/4.10/Training#Segmentation>>.

2.2.3. Click on the **Show 3D** button to view a 3D representation of the segmentation.

2.3. Export the segmentation in a 3D model file format by going to the **Segmentations** module in 3D Slicer.

2.3.1. Go to **Export/import models and labelmaps**. Select **Export** in the operation section and **Models** in the output type section. Click **Export** to finish and create the 3D model from the segmented area.

2.3.2. Select **SAVE** (upper left) to save the model. Choose the elements to be saved. Then, change the file format of the 3D Model to “OBJ” within the **File format** column. Select the path where files will be stored and click on **Save**.

2.4. Repeat steps 2.2 and 2.3 to create additional 3D models of different anatomical regions.

NOTE: Pre-segmented models of the provided example can be found in the data previously downloaded in step 1.3 (“/Data/Biomodels/”).

3. **Biomodel positioning**

NOTE: In this section, the 3D models created in Section 2 will be positioned with respect to the marker for augmented reality visualization. **The ARHealth: Model Position module from 3D Slicer will be used for this task.** Follow the instructions provided in step 1.3 to add the module to 3D Slicer. There are two different alternatives to position the 3D models: “Visualization” mode and “Registration” mode.

3.1. **Visualization mode**

NOTE: **Visualization** mode allows positioning of the 3D patient models at any position with respect to the AR marker. With this option, the user are able to use the AR app to visualize biomodels using the 3D-printed AR marker as a reference. This mode may be used when precision is not required, and visualization of the virtual model can be displayed anywhere within the field-of-view of the smartphone camera and marker.

3.1.1. Go to the **ARHealth: Model Position** module, and (in the initialization section) select **Visualization** mode. Click on **Load Marker Model** to load the marker for this option.

3.1.2. Load the 3D models created in section 2 by clicking on the ... button to select the path of the saved models from section 2. Then, click on the **Load Model** button to load it in 3D Slicer. Models must be loaded one at a time. To delete any models previously loaded, click on that model followed by the **Remove Model** button, or click **Remove All** to delete all models loaded at once.

3.1.3. Click the **Finish and Center** button to center all models within the marker.

3.1.4. The position, orientation, and scaling of the 3D models can be modified with respect to the marker with different slider bars (i.e., translation, rotation, scale).

NOTE: There is an additional “Reset Position” button to reset the original position of the models before making any changes in the position.

3.1.5. Save the models at this position by choosing the path to store the files and clicking the **Save Models** button. The 3D models will be saved with the extension name “_registered.obj”.

3.2. Registration mode

NOTE: **Registration** mode allows combining of the AR marker with one 3D biomodel at any desired position. Then, any section of the combined 3D models (that includes the AR marker) can be extracted and 3D-printed. All biomodels will be displayed in the AR app using this combined 3D-printed biomodel as a reference. This mode allows the user to easily register the patient (here, a section of the patient’s bone) and virtual models using a reference marker.

3.2.1. Go to the **ARHealth: Model Position** module, and (in the initialization section) select **Registration** mode. Click on **Load Marker Model** to load the marker for this option.

3.2.2. Load the models as done in step 3.1.2.

3.2.3. Move the 3D models and ensure their intersection with the supporting structure of the cube marker, since these models will be combined and 3D-printed later. The height of the marker base can be modified. The position, orientation, and scaling of the 3D models can be modified with respect to the marker with different slider bars (i.e., translation, rotation, scale).

3.2.4. Save the models at this position by choosing the path to store the files and clicking the **Save Models** button. The 3D models will be saved with the extension name “_registered.obj”.

3.2.5. The anatomy model may be too large. If so, cut the 3D model around the marker adaptor and 3D-print only a section of the combination of both models using Meshmixer software.

3.2.6. Open Meshmixer and load the biomodel and supporting structure of the cube marker model saved in step 3.2.5. Combine these models by selecting both models in the **Object Browser** window. Click on the **Combine** option in the tool window that has just appeared in the upper left corner.

3.2.7. In Meshmixer, use the **Plane Cut** tool under the **Edit** menu to remove unwanted sections of the model that will not be 3D-printed.

3.2.8. To save the model to be 3D-printed, go to **File > Export** and select the desired format.

4. 3D printing

NOTE: The aim of this step is to 3D-print the physical models required for the final AR application. The marker to be detected by the application and the different objects needed depend on the mode selected in section 3. Any material can be used for 3D printing for the purpose of this work, when following the color material requirements requested at each step. Polylactic acid (PLA) or acrylonitrile butadiene styrene (ABS) are both sufficient choices.

4.1. Use a 3D printer to print the cubic marker. If a dual extruder 3D printer is not available, skip to step 4.2. Use a dual extruder 3D printer specifically to print the two-color marker provided in "Data/Markers/Marker1_TwoColorCubeMarker/". In the 3D printing software, select a white color material for the file "TwoColorCubeMarker_WHITE.obj" and black color material for "TwoColorCubeMarker_BLACK.obj".

NOTE: For better marker detection, print on high-quality mode with a small layer height.

4.2. If a dual extruder 3D printer is not available and step 4.1 was not performed, follow this step to print a 3D-printed marker with stickers as an alternative by doing the following:

4.2.1. Use a 3D printer to print the file "Data/Markers/ Marker2_StickerCubeMarker/ StickerCubeMarker_WHITE.obj" with white color material.

4.2.2. Use a conventional printer to print the file "Data/Markers/ Marker2_StickerCubeMarker/Stickers.pdf" on sticker paper. Then, use any cutting tool to precisely cut the images though the black frame by removing the black lines.

NOTE: It is recommended to use sticker paper to obtain a higher quality marker. However, the images can be printed on regular paper, and a common glue stick can be used to paste the images on the cube.

4.2.3. Place stickers in the 3D-printed cube obtained in step 4.2.1 in the corresponding order following instructions from the document "Data/Markers/ Marker2_StickerCubeMarker/Stickers.pdf".

NOTE: Stickers are smaller than the face of the cube. Leave a 1.5 mm frame between the sticker and edge of the face. "Data/Markers/Marker2_StickerCubeMarker/StickerPlacer.stl" can be 3D-printed to guide the sticker positioning and exactly match the center of the cube face.

4.3. 3D-print the adaptors, depending on the mode selected in section 3.

4.3.1. If **Visualization** mode (section 3.1), was selected, 3D-print "Data/3DPrinting/Option1/ MarkerBaseTable.obj", which is a base adaptor used to place the marker in vertical position on a horizontal surface.

4.3.2. If **Registration** mode (section 3.2) was selected, 3D-print the model created in step 3.2.8

with the marker adaptor attached.

NOTE: 3D printed objects from step 4.2 can be printed in any color material.

5. AR app deployment

NOTE: The goal of this section is to design a smartphone app in Unity engine that includes the 3D models created in the previous sections and deploy this app on a smartphone. A Vuforia Development License Key (free for personal use) is required for this step. The app can be deployed on Android or iOS devices.

5.1. Create a Vuforia Developer account to obtain a license key to use their libraries in Unity. Go to the link found at <<https://developer.vuforia.com/vui/auth/register>> and create an account.

5.1.1. Go to the link found at <<https://developer.vuforia.com/vui/develop/Licenses>> and select **Get Development Key**. Then, follow the instructions to add a free development license key into the user's account.

5.1.2. In the **License Manager** menu, select the key created in the previous step and copy the provided key, which will be used in step 5.3.3.

5.2. Set up the smartphone.

5.2.1. To get started with Unity and Android devices, go to the link found at <<https://docs.unity3d.com/Manual/android-GettingStarted.html>>.

5.2.2. To get started with Unity and iOS devices, go to the link found at <<https://docs.unity3d.com/Manual/iphone-GettingStarted.html>>.

5.3. Set up a Unity Project for the AR app by first opening Unity v.2019 and creating a new 3D project. Then, under **Build Settings** in the **File** menu, switch the platform to either an Android or iOS device.

5.3.1. Enable Vuforia into the project by selecting **Edit > Project Setting > Player Settings > XR Settings** and checking the box labeled **Vuforia Augmented Reality Support**.

5.3.2. Create an "ARCamera" under **Menubar > GameObject > Vuforia > ARCamera** and import Vuforia components when prompted.

5.3.3. Add the Vuforia License Key into Vuforia Configuration settings by selecting the **Resources** folder and clicking on **Vuforia Configuration**. Then, in the **App License Key** section, paste the key copied in section 5.1.2.

5.3.4. Import the Vuforia Target file provided in "/Data/Vuforia/ AR_Cube_3x3x3.unitypackage"

into Unity, which contains the files that Vuforia requires to detect the markers described in section 4.

5.3.5. Create a Vuforia MultiTarget under **Menubar > GameObject > Vuforia > MultiTarget**.

5.3.6. Select the marker type that will be used for detection by clicking on the MultiTarget created in the previous step. In the **Database** option under **Multi Target Behaviour**, select **ARHealth_3DPrintedCube_30x30x30**. In the **Multi Target** option under **Multi Target Behaviour**, select either **TwoColorCubeMarker** or **StickerCubeMarker**, depending on the marker created in section 4.

5.3.7. Load the 3D models created in section 3 into Unity Scene under **MultiTarget** by creating a new folder with the name “Models” under the “Resources” folder. Drag the 3D models into this folder. Once loaded in Unity, drag them under the “MultiTarget” item created in step 5.3.5. This will make them dependent on the marker.

NOTE: Models should be visible in the Unity 3D view scene.

5.3.8. Change the colors of the 3D models by creating a new material and assigning the new materials to the models.

5.3.8.1. Create a new folder named “Materials” under the “Resources” folder by going to **Menubar > Assets > Create > Material**. Select the material and change the color in the configuration section. Then, drag the file under the 3D model hierarchy.

5.3.9. Optional: if there is a webcam available, click on the play button located in the upper-middle portion to test its application on the computer. If the marker is visible to the webcam, it should be detected, and the 3D models should appear in the scene.

5.3.10. If an Android smartphone is used for app deployment, go to **File > Build Settings** in Unity, and select the plugged phone from the list. Select **Deploy and Run**. Save the file with extension .apk on the computer and allow the process to finish. Once deployment is done, the app should be on the phone and ready to run.

NOTE: This protocol has been tested on Android v.8.0 Oreo or above. Correct functionality is not guaranteed for older versions.

5.3.11. If the app will be deployed in an iOS device, go to **File > Build Settings** in Unity and select **Run**. Select the path to save the app files. Allow the process to finish. Go to the saved folder and open the file with the extension “.projectxcode”.

5.3.11.1. In Xcode, follow the instructions from step 5.2.2 to complete deployment.

NOTE: For more information about Vuforia in Unity, go to the link found at <https://Library.vuforia.com/articles/Training/getting-started-with-vuforia-in-unity.html>.

6. App visualization

6.1. Open the installed app, which will use the smartphone's camera. When running the app, look at the marker with the camera from a short distance away (40 cm minimum). Once the app detects the marker, the 3D models created in previous steps should appear exactly at the location defined during the procedure on the smartphone screen.

NOTE: Illumination can alter the precision of marker detection. It is recommended to use the app in environments with good lighting conditions.

REPRESENTATIVE RESULTS:

The protocol was applied to data from a patient suffering from distal leg sarcoma in order to visualize the affected anatomical region from a 3D perspective. Using the method described in section 2, the portion of the affected bone (here, the tibia and fibula) and tumor were segmented from the patient's CT scan. Then, using the segmentation tools from 3D Slicer, two biomodels were created: the bone (section of the tibia and fibula) (**Figure 1A**) and tumor (**Figure 1B**).

Next, the two 3D models were positioned virtually with respect to the marker for optimal visualization. Both modes described in section 3 were followed for this example. For visualization mode, the models were centered in the upper face of the marker (**Figure 2**). For registration mode, the marker adaptor was positioned in the bone (specifically, the tibia [**Figure 3**]). Then, a small section of the tibia was selected to be 3D-printed with a 3D marker adaptor (**Figure 4**). An Ultimaker 3 extended 3D printer with PLA material was used to create the 3D-printed markers (**Figure 5A,B**), marker holder base (**Figure 5C**) for the "visualization" mode, and section of the tibia for "registration" mode (**Figure 5D**). **Figure 5E** shows how the marker was attached to the "visualization" mode 3D-printed base. **Figure 5F** shows the attachment with the "registration" mode 3D-printed biomodel. Finally, Unity was used to create the app and deploy it in the smartphone.

Figure 6 shows how the app worked for "visualization" mode. The hologram was accurately located in the upper part of the cube as previously defined. **Figure 7** shows the application for "registration" mode, in which the app positioned the complete bone model on top of the 3D-printed section. The final visualization of the holograms was clear and realistic, maintained the real sizes of the biomodels, and positioned accurately. When using the smartphone application, the AR marker needs to be visible by the camera for the app to correctly display the holograms. In addition, the light conditions in the scene must be of good quality and constant for proper marker detection. Bad light conditions or reflections on the marker surface hinder the tracking of the AR marker and cause malfunctioning of the app.

The time required to create the app depends on several factors. The duration of section 1 is limited by the download speed. Regarding anatomy segmentation (section 2), factors affecting segmentation time include complexity of the region and medical imaging modality (i.e., CT is easily segmented, while MRI is more difficult). For the representative example of the tibia, approximately 10 min was required to generate both 3D models from the CT scan. Biomodel positioning (section 3) is simple and straightforward. Here, it took approximately 5 min to define the biomodel position with respect to the AR marker. For the 3D printing step, the duration is highly dependent on the selected mode. The “dual color marker” was manufactured at high quality in a period of 5 h and 20 min. The “sticker marker” was manufactured in a period of 1 h and 30 min, plus the time required to paste the stickers. The final step for app development can be time-consuming for those with no previous experience in Unity, but it can be easily completed following the protocol steps. Once the AR markers have been 3D-printed, the development of an entirely new AR app can be performed in less than 1 h. This duration can be further reduced with additional experience.

FIGURE AND TABLE LEGENDS:

Figure 1: Representation of 3D models created from CT image of a patient suffering from distal leg sarcoma. (A) Bone tissue represented in white (tibia and fibula). **(B)** Tumor represented in red.

Figure 2: Results showing how “visualization” mode in 3D Slicer positions the virtual 3D models of the bone and tumor with respect to 3D-printed marker reference. The patient 3D models **(A)** are positioned above the upper face of the marker cube **(B)**.

Figure 3: Results showing how “registration” mode in 3D Slicer positions the virtual 3D models of the bone and tumor (A) with respect to 3D-printed marker reference (B). The marker adaptor is attached to the bone tissue model.

Figure 4: Small section of the bone tissue and 3D marker adaptor. The two components are combined then 3D-printed.

Figure 5: 3D printed tools required for the final application. (A) “Two color cube marker” 3D-printed with two colors of materials. **(B)** “Sticker cube marker” 3D-printed, with stickers pasted. **(C)** Marker base cube adaptor. **(D)** Section of the patient’s bone tissue 3D model and marker cube adaptor. **(E)** “Sticker cube marker” placed in the marker base cube adaptor. **(F)** “Two color cube marker” placed in the marker adaptor attached to the patient’s anatomy.

Figure 6: App display when using “visualization” mode. The patient’s affected anatomy 3D models are positioned above the upper face of the 3D-printed cube.

Figure 7: AR visualization when using “registration” mode. The 3D-printed marker enables registration of the 3D-printed biomodel with the virtual 3D models.

DISCUSSION:

AR holds great potential in education, training, and surgical guidance in the medical field. Its combination with 3D printing opens many new possibilities in clinical applications. This protocol describes methodology that enables inexperienced users to create a smartphone app combining AR and 3DP for the visualization of anatomical 3D models of patients with 3D-printed reference markers.

In general, one of the most interesting clinical applications of AR and 3DP is to improve patient-to-physician communication by giving the patient a different perspective of the case, improving explanations of specific medical conditions or treatments. Another possible application includes surgical guidance for target localization, in which 3D-printed patient-specific tools (with a reference AR marker attached) can be placed on rigid structures (i.e., bone) and used as a reference for navigation. This application is especially useful for orthopedic and maxillofacial surgical procedures, in which bone tissue surface is easily accessed during surgery.

The protocol starts with section 1, describing the workstation set-up and software tools necessary. Section 2 describes how to use 3D Slicer software to easily segment target anatomies of the patient from any medical imaging modality to obtain 3D models. This step is crucial, as the virtual 3D models created are those displayed in the final AR application.

In section 3, 3D Slicer is used to register the 3D models created in the previous section with an AR marker. During this registration procedure, patient 3D models are efficiently and simply positioned with respect to the AR marker. The position defined in this section will determine the hologram relative position in the final app. It is believed that this solution reduces complexity and multiplies the possible applications. Section 3 describes two different options to define the spatial relationships between the models and AR markers: “visualization” and “registration” mode. The first option, “visualization” mode, allows the 3D models to be positioned anywhere with respect to the marker and displayed as the whole biomodel. This mode provides a realistic, 3D perspective of the patient’s anatomy and allows moving and rotating of the biomodels by moving the tracked AR marker. The second option, “registration” mode, allows attachment and combining of a marker adaptor to any part of the biomodel, offering an automatic registration process. With this option, a small section of the 3D model, including the marker adaptor, can be 3D-printed, and the app can display the rest of the model as a hologram.

Section 4 provides guidelines for the 3D printing process. First, the user can choose between two different markers: the “dual color marker” and “sticker marker”. The whole “dual color marker” can be 3D-printed but requires a dual extruder 3D printer. In case this printer is not available, the “sticker marker” is proposed. This is a simpler marker that can be obtained by 3D-printing the cubic structure, then pasting the images of the cube with sticker paper or sticker glue. Furthermore, both markers were designed with extensible sections to perfectly fit in a specific adaptor. Thus, the marker can be reused in several cases.

Section 5 describes the process to create a Unity project for AR using the Vuforia software development kit. This step may be the hardest portion for users with no programming experience, but with these guidelines, it should be easier to obtain the final application that is presented in section 6. The app displays the patient's virtual models over the smartphone screen when the camera recognizes the 3D-printed marker. In order for the app to detect the 3D marker, a minimum distance of approximately 40 cm or less from the phone to marker as well as good lighting conditions are required.

The final application of this protocol allows the user to choose the specific biomodels to visualize and in which positions. Additionally, the app can perform automatic patient-hologram registration using a 3D-printed marker and adaptor attached to the biomodel. This solves the challenge of registering virtual models with the environment in a direct and convenient manner. Moreover, this methodology does not require broad knowledge of medical imaging or software development, does not depend on complex hardware and expensive software, and can be implemented over a short time period. It is expected that this method will help accelerate the adoption of AR and 3DP technologies by medical professionals.

ACKNOWLEDGMENTS:

This report was supported by projects PI18/01625 and PI15/02121 (Ministerio de Ciencia, Innovación y Universidades, Instituto de Salud Carlos III and European Regional Development Fund "Una manera de hacer Europa") and IND2018/TIC-9753 (Comunidad de Madrid).


DISCLOSURES:


The authors have nothing to disclose.


REFERENCES:

1. Coles, T. R., John, N. W., Gould, D., Caldwell, D. G. Integrating Haptics with Augmented Reality in a Femoral Palpation and Needle Insertion Training Simulation. *IEEE Transactions on Haptics*. **4** (3), 199–209 (2011).
2. Pelargos, P. E. et al. Utilizing virtual and augmented reality for educational and clinical enhancements in neurosurgery. *Journal of Clinical Neuroscience*. **35**, 1–4 (2017).
3. Abhari, K. et al. Training for Planning Tumour Resection: Augmented Reality and Human Factors. *IEEE Transactions on Biomedical Engineering*. **62** (6), 1466–1477 (2015).
4. Uppot, R. et al. Implementing Virtual and Augmented Reality Tools for Radiology Education and Training, Communication, and Clinical Care. *Radiology*. **291**, 182210 (2019).
5. Pratt, P. et al. Through the HoloLens™ looking glass: augmented reality for extremity reconstruction surgery using 3D vascular models with perforating vessels. *European Radiology Experimental*. **2** (1), 2 (2018).
6. Rose, A. S., Kim, H., Fuchs, H., Frahm, J. -M. Development of augmented-reality applications in otolaryngology-head and neck surgery. *The Laryngoscope* (2019).
7. Zhou, C. et al. Robot-Assisted Surgery for Mandibular Angle Split Osteotomy Using Augmented Reality: Preliminary Results on Clinical Animal Experiment. *Aesthetic Plastic Surgery*. **41** (5), 1228–1236 (2017).



8. Heinrich, F., Joeres, F., Lawonn, K., Hansen, C. Comparison of Projective Augmented Reality Concepts to Support Medical Needle Insertion. *IEEE Transactions on Visualization and Computer Graphics*. **25** (6), 1 (2019).
9. Deng, W., Li, F., Wang, M., Song, Z. Easy-to-Use Augmented Reality Neuronavigation Using a Wireless Tablet PC. *Stereotactic and Functional Neurosurgery*. **92** (1), 17–24 (2014).
10. Fan, Z., Weng, Y., Chen, G., Liao, H. 3D interactive surgical visualization system using mobile spatial information acquisition and autostereoscopic display. *Journal of Biomedical Informatics*. **71**, 154–164 (2017).
11. Arnal-Burró, J., Pérez-Mañanes, R., Gallo-del-Valle, E., Igualada-Blazquez, C., Cuervas-Mons, M., Vaquero-Martín, J. Three dimensional-printed patient-specific cutting guides for femoral varization osteotomy: Do it yourself. *The Knee*. **24** (6), 1359–1368 (2017).
12. Vaquero, J., Arnal, J., Perez-Mañanes, R., Calvo-Haro, J., Chana, F. 3D patient-specific surgical printing cutting blocks guides and spacers for open- wedge high tibial osteotomy (HTO) - do it yourself. *Revue de Chirurgie Orthopédique et Traumatologique*. **102** (7, Supplement), S131 (2016).
13. De La Peña, A., De La Peña-Brambila, J., Pérez-De La Torre, J., Ochoa, M., Gallardo, G. Low-cost customized cranioplasty using a 3D digital printing model: a case report. *3D Printing in Medicine*. **4** (1), 1–9 (2018).
14. Kamio, T. et al. Utilizing a low-cost desktop 3D printer to develop a “one-stop 3D printing lab” for oral and maxillofacial surgery and dentistry fields. *3D Printing in Medicine*. **4** (1), 1–7 (2018).
15. Punyaratabandhu, T., Liacouras, P., Pairojboriboon, S. Using 3D models in orthopedic oncology: presenting personalized advantages in surgical planning and intraoperative outcomes. *3D Printing in Medicine*. **4** (1), 1–13 (2018).
16. Wake, N. et al. Patient-specific 3D printed and augmented reality kidney and prostate cancer models: impact on patient education. *3D Printing in Medicine*. **5** (1), 1–8 (2019).
17. Barber, S. R. et al. Augmented Reality, Surgical Navigation, and 3D Printing for Transcanal Endoscopic Approach to the Petrous Apex. *OTO Open: The Official Open Access Journal of the American Academy of Otolaryngology-Head and Neck Surgery Foundation*. **2** (4) (2018).
18. Witowski, J. et al. Augmented reality and three-dimensional printing in percutaneous interventions on pulmonary arteries. *Quantitative Imaging in Medicine and Surgery*. **9** (1) (2019).
19. Moreta-Martínez, R., García-Mato, D., García-Sevilla, M., Pérez-Mañanes, R., Calvo-Haro, J., Pascau, J. Augmented reality in computer-assisted interventions based on patient-specific 3D printed reference. *Healthcare Technology Letters*. (2018).







FRONT











Figure6

[Click here to access/download;Figure;Figure6.jpg](#)



Figure7

[Click here to access/download;Figure;Figure7.jpg](#)




Table of Materials		
Name of Material/Equipment	Company	Catalog Number
3D Printing material: Acrylonitrile Butadiene Styrene (ABS)		
3D Printing material: Polylactic Acid (PLA)		
3D Slicer		
Android	Alphabet, Inc.	
Autodesk Meshmixer	Autodesk, Inc.	
iPhone OS	Apple, Inc.	
Ultimaker 3 Extended	Ultimaker BV	
Unity	Unity Technologies	
Xcode	Apple, Inc.	

Comments/Description

Thermoplastic polymer material usually used in domestic 3D printers.

Bioplastic material usually used in domestic 3D printers.

Open-source software platform for medical image informatics, image processing, and three-dimensional visualization
Android is a mobile operating system developed by Google. It is based on a modified version of the Linux kernel and other open source software, and is designed primarily for touchscreen mobile devices such
Meshmixer is state-of-the-art software for working with triangle meshes. Free software.

iPhone OS is a mobile operating system created and developed by Apple Inc. exclusively for its hardware.

Fused deposition modeling 3D printer.
Unity is a real-time development platform to create 3D, 2D VR & AR visualizations for Games, Auto, Transportation, Film, Animation, Architecture, Engineering & more. Free software.
Xcode is a complete developer toolset for creating apps for Mac, iPhone, iPad, Apple Watch, and Apple TV. Free software.

ARTICLE AND VIDEO LICENSE AGREEMENT

Title of Article:

Combining Augmented Reality and 3D Printing to Display Patient Models on Your Smartphone

Author(s):

Rafael Moreta-Martinez, David García-Mato, Mónica García-Sevilla, Rubén Pérez-Mañanes,
José A. Calvo-Haro, Javier Pascau

Item 1: The Author elects to have the Materials be made available (as described at <http://www.jove.com/publish>) via:

☒

Standard Access

☐

Open Access

Item 2: Please select one of the following items:

☒

The Author is **NOT** a United States government employee.

☐

The Author is a United States government employee and the Materials were prepared in the course of his or her duties as a United States government employee.

☐

The Author is a United States government employee but the Materials were NOT prepared in the course of his or her duties as a United States government employee.

ARTICLE AND VIDEO LICENSE AGREEMENT

1. **Defined Terms.** As used in this Article and Video License Agreement, the following terms shall have the following meanings: **"Agreement"** means this Article and Video License Agreement; **"Article"** means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein; **"Author"** means the author who is a signatory to this Agreement; **"Collective Work"** means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; **"CRC License"** means the Creative Commons Attribution-Non Commercial-No Derivs 3.0 Unported Agreement, the terms and conditions of which can be found at: <http://creativecommons.org/licenses/by-nc-nd/3.0/legalcode>; **"Derivative Work"** means a work based upon the Materials or upon the Materials and other pre-existing works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; **"Institution"** means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; **"JoVE"** means MyJoVE Corporation, a Massachusetts corporation and the publisher of The Journal of Visualized Experiments; **"Materials"** means the Article and / or the Video; **"Parties"** means the Author and JoVE; **"Video"** means any video(s) made by the Author, alone or in conjunction with any other parties, or by JoVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion

of the Article, and in which the Author may or may not appear.

2. **Background.** The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JoVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.

3. **Grant of Rights in Article.** In consideration of JoVE agreeing to publish the Article, the Author hereby grants to JoVE, subject to **Sections 4** and **7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the "Open Access" box has been checked in **Item 1** above, JoVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.

ARTICLE AND VIDEO LICENSE AGREEMENT

4. **Retention of Rights in Article.** Notwithstanding the exclusive license granted to JoVE in **Section 3** above, the Author shall, with respect to the Article, retain the non-exclusive right to use all or part of the Article for the non-commercial purpose of giving lectures, presentations or teaching classes, and to post a copy of the Article on the Institution's website or the Author's personal website, in each case provided that a link to the Article on the JoVE website is provided and notice of JoVE's copyright in the Article is included. All non-copyright intellectual property rights in and to the Article, such as patent rights, shall remain with the Author.

5. **Grant of Rights in Video – Standard Access.** This **Section 5** applies if the "Standard Access" box has been checked in **Item 1** above or if no box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby acknowledges and agrees that, Subject to **Section 7** below, JoVE is and shall be the sole and exclusive owner of all rights of any nature, including, without limitation, all copyrights, in and to the Video. To the extent that, by law, the Author is deemed, now or at any time in the future, to have any rights of any nature in or to the Video, the Author hereby disclaims all such rights and transfers all such rights to JoVE.

6. **Grant of Rights in Video – Open Access.** This **Section 6** applies only if the "Open Access" box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby grants to JoVE, subject to **Section 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Video in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Video into other languages, create adaptations, summaries or extracts of the Video or other Derivative Works or Collective Works based on all or any portion of the Video and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. For any Video to which this **Section 6** is applicable, JoVE and the Author hereby grant to the public all such rights in the Video as provided in, but subject to all limitations and requirements set forth in, the CRC License.

7. **Government Employees.** If the Author is a United States government employee and the Article was prepared in the course of his or her duties as a United States government employee, as indicated in **Item 2** above, and any of the licenses or grants granted by the Author hereunder exceed the scope of the 17 U.S.C. 403, then the rights granted hereunder shall be limited to the maximum

rights permitted under such statute. In such case, all provisions contained herein that are not in conflict with such statute shall remain in full force and effect, and all provisions contained herein that do so conflict shall be deemed to be amended so as to provide to JoVE the maximum rights permissible within such statute.

8. **Protection of the Work.** The Author(s) authorize JoVE to take steps in the Author(s) name and on their behalf if JoVE believes some third party could be infringing or might infringe the copyright of either the Author's Article and/or Video.

9. **Likeness, Privacy, Personality.** The Author hereby grants JoVE the right to use the Author's name, voice, likeness, picture, photograph, image, biography and performance in any way, commercial or otherwise, in connection with the Materials and the sale, promotion and distribution thereof. The Author hereby waives any and all rights he or she may have, relating to his or her appearance in the Video or otherwise relating to the Materials, under all applicable privacy, likeness, personality or similar laws.

10. **Author Warranties.** The Author represents and warrants that the Article is original, that it has not been published, that the copyright interest is owned by the Author (or, if more than one author is listed at the beginning of this Agreement, by such authors collectively) and has not been assigned, licensed, or otherwise transferred to any other party. The Author represents and warrants that the author(s) listed at the top of this Agreement are the only authors of the Materials. If more than one author is listed at the top of this Agreement and if any such author has not entered into a separate Article and Video License Agreement with JoVE relating to the Materials, the Author represents and warrants that the Author has been authorized by each of the other such authors to execute this Agreement on his or her behalf and to bind him or her with respect to the terms of this Agreement as if each of them had been a party hereto as an Author. The Author warrants that the use, reproduction, distribution, public or private performance or display, and/or modification of all or any portion of the Materials does not and will not violate, infringe and/or misappropriate the patent, trademark, intellectual property or other rights of any third party. The Author represents and warrants that it has and will continue to comply with all government, institutional and other regulations, including, without limitation all institutional, laboratory, hospital, ethical, human and animal treatment, privacy, and all other rules, regulations, laws, procedures or guidelines, applicable to the Materials, and that all research involving human and animal subjects has been approved by the Author's relevant institutional review board.

11. **JoVE Discretion.** If the Author requests the assistance of JoVE in producing the Video in the Author's facility, the Author shall ensure that the presence of JoVE employees, agents or independent contractors is in accordance with the relevant regulations of the Author's institution. If more than one author is listed at the beginning of this Agreement, JoVE may, in its sole

ARTICLE AND VIDEO LICENSE AGREEMENT

discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JoVE reserves the right, in its absolute and sole discretion and without giving any reason therefore, to accept or decline any work submitted to JoVE. JoVE and its employees, agents and independent contractors shall have full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JoVE has sole discretion as to the method of making and publishing the Materials, including, without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

12. **Indemnification.** The Author agrees to indemnify JoVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JoVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JoVE, making of videos by JoVE, or publication in JoVE or elsewhere by JoVE. The Author shall be responsible for, and shall hold JoVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to

the making of a video by JoVE its employees, agents or independent contractors. All sterilization, cleanliness or decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's expense. All indemnifications provided herein shall include JoVE's attorney's fees and costs related to said losses or damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JoVE, its employees, agents or independent contractors.

13. **Fees.** To cover the cost incurred for publication, JoVE must receive payment before production and publication of the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JoVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.

14. **Transfer, Governing Law.** This Agreement may be assigned by JoVE and shall inure to the benefits of any of JoVE's successors and assignees. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to be one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement is required per submission.

CORRESPONDING AUTHOR

Name:

Javier Pascau

Department:

Bioingeniería e Ingeniería Aeroespacial

Institution:

Universidad Carlos III de Madrid

Title:

Associate Professor

Signature:



Date:

July 29 2019

Please submit a **signed** and **dated** copy of this license by one of the following three methods:

1. Upload an electronic version on the JoVE submission site
2. Fax the document to +1.866.381.2236
3. Mail the document to JoVE / Attn: JoVE Editorial / 1 Alewife Center #200 / Cambridge, MA 02140

Rebuttal Document

Manuscript JoVE60618

Title: Combining Augmented Reality and 3D Printing to Display Patient Models on Your Smartphone

We would like to thank the reviewers and the editorial for their insightful comments and their interest in our work.

First, we would like to summarize the major changes made for the paper revision:

- In the protocol section, the levels of subtitles have been modified to have a more appropriate format following the editorial instructions. Now sublevels 4 are sparingly used.
- Segmentation procedure in 2.2.2 is now more deeply explained.
- A new paragraph has been added in the introduction section to give a clearer vision of the proposed method in the context of the state-of-the-art.
- The introductory paragraph of step 3 of the protocol section has been modified to improve the explanation of “Visualization” and “Registration” modes.
- Possible applications of our method are discussed in a new paragraph in the discussion section.
- Software improvements in the repository provided for the 3D Slicer Module “AR Health: Model Position”. A new button has been added for restarting the process without having to restart 3D Slicer.
- Figures 2 and 3 have been updated, highlighting different elements in the image.

In the following pages of the document, you will find a separate answer to each reviewer’s comments.

Editorial comments:

We would like to thank the Editorial for their insightful review of the manuscript and for their comments.

General:

- 1. Please take this opportunity to thoroughly proofread the manuscript to ensure that there are no spelling or grammar issues.*
- 2. Please limit the use of personal pronouns (you, your*
- 3. Please include email addresses for all authors in the manuscript.*

Answer:

1. We apologize for any spelling or grammar issues. The paper has been reviewed again and corrected.
2. The use of personal pronouns has been reduced or avoided. Examples can be found in steps 2.2.1, 3.1.4, 3.2.3, 3.2.5, 4.2, 5.1.2, 5.2, 5.2.9.2, 5.3.10, 5.3.12.1.1., 6. All these steps are referenced to the first version of the manuscript.
3. Email addresses for all authors have been included in “Authors and Affiliations” section.

Protocol:

- 1. Please include an ethics statement before the numbered protocol steps, indicating that the protocol follows the guidelines of your institution’s human research ethics committee.*
- 2. If possible, please limit substeps to no more than 4 levels (e.g., 1.1.1.1), with a 4th level used sparingly.*
- 3. For each protocol step, please ensure you answer the “how” question, i.e., how is the step performed? Alternatively, add references to published material specifying how to perform the protocol action. If revisions cause a step to have more than 2-3 actions and 4 sentences per step, please split into separate steps or substeps.*

Answer:

1. The following ethics statement has been included before the numbered protocol steps: “This study was performed in accordance with the principles of the 1964 Declaration of Helsinki as revised in 2013. The anonymized patient data and pictures included in this paper are used after written informed consent was obtained from the participant and/or their legal representative, where he/she approved the use of this data for dissemination activities including scientific publications.”
2. Most of the level 4 steps have been removed to higher-level steps.
3. The manuscript steps and substeps have been modified to a more appropriated format following the instructions given by the Editorial.
 - Step 3.1.2 has been split, adding step 3.1.3, to avoid explaining too many actions in the same substep.
 - Step 3.2.5 and its substeps have been modified to level 3 substeps in order to avoid level 4 steps. “How” question has been answered specifically in some substeps.

- Steps and substeps from 4.1. have been extended to steps 4.1 and 4.2.
- Substeps in steps 5.3 have been adapted to avoid level 4 steps.

Table of Materials:

1. Please ensure the Table of Materials has information on all materials and equipment used, especially those mentioned in the Protocol.

Answer:

1. The Table of Materials has been modified by adding new materials that were not previously included.

Reviewers' comments

Reviewer #1

Manuscript Summary:

The authors describe a method for presenting segmented medical images in an augmented reality smartphone app with virtual content registered to a 3D printed model. Registration is performed by manually adjusting a virtual model to a Vuforia marker 3D model.

Major Concerns:

None.

Minor Concerns:

Item 2.2.2 summarizes the segmentation details with a link to the Slicer documentation on segmentation. Although this is an acceptable description, I would prefer to see 2-4 sentences on the specific effects the authors used to segment the distal leg sarcoma images. Or even one summary summarizing the effects used: "To segment the desired patient geometry, we used the threshold, etc... effects."

Answer:

We would like to thank Reviewer #1 for his time on reading and commenting on our paper. Now we will proceed to answer the reviewer's minor concerns:

Regarding the segmentation step in item 2.2.2, we agree with the reviewer that summarizing the segmentation procedure with a link to Slicer documentation tutorial may not be enough to show a user how to segment properly the proposed patient case. Therefore, a new substep 2.2.2.1 has been added explaining briefly which tools were used from the Slicer Segmentation module to finally obtain the given 3D models of the patient. The added step is the following: "2.2.2.1 To segment the bone (tibia and fibula in this case), the "Threshold" tool was used to set up a minimum and a maximum HU value from the CT image which correspond to bone tissue. By using this tool, other elements with HU outside these threshold values are removed, such as soft tissue. Then, "Scissors" tool was used to remove some undesired areas such as the bed or other anatomical structures from the segmented mask. The sarcoma region was segmented manually using "Draw" and "Erase" tools since the tumor is difficult to contour with automatic tools".

Reviewer #2

Manuscript Summary:

This article proposes a method to display AR scene based on markers on a smartphone. A step-by-step methodology was described to realize this work. 3D slicer was used to segment and reconstruct the anatomical 3D models of patients from CT or MR images. An AR Health module was also developed to design the location relationship between marker and visual 3D models through two modes of Visualization and Registration. Meshmixer was used to modify the combination of both models. Then, the physical models required were 3D printed for the final AR application. An AR smartphone app in Unity engine was designed and deployed on Android or iOS devices. Finally, a real clinical case of a patient suffering from a distal leg sarcoma was to describe the propose methodology.

This manuscript is well organized.

My major concerns are as follows:

- 1. The proposed method is not clear. Compared with other work about AR using smartphone, the differences and similarities among the proposed method and other AR methods should be described more clearly. Using markers to display AR scene has been proposed in previous articles.*
- 2. The accuracy of AR in the proposed method is not described. The author could show an in-situ superposition of 3D images on patients. Actually, it is very different to install the marker on bones and on deformed tissues. Tissue deformation affects registration accuracy. So, please give examples of some suitable clinical applications.*
- 3. In AR Health module, both Visualization and Registration modes obtain the location of the marker and the 3D models by translating and rotating them. It is hoped that the differences between the two operations can be further explained.*
- 4. The proposed method seems to be a combination of several software and the existing methods. For some diseases, medical image segmentation is difficult, especially for inexperienced doctors. Also, the process described in this manuscript is complicated. Authors may describe each process simply and express them in pictures.*
- 5. I have tried the AR Health module, which can only be executed once and needs further optimization. For examples, it cannot return to remove or add models after “Finish and Center”, and the registration process cannot be performed after the visualization process.*

Answer:

We would like to thank Reviewer #2 for his time on reading and commenting insightfully our paper. In the following paragraphs we answer all reviewer's major concerns:

1. The purpose of this paper is not to present a new method for showing AR scenes with data coming from medical images but to offer inexperienced users a complete protocol on how to use AR technology for that purpose. We include several references showing how this technology has been used previously in the biomedical field, providing context to the reader about the interest of these solutions. We have considered that a JOVE video and a full step-

by-step protocol will help novel users to try AR since our method is based on inexpensive tools. Markers-based AR has been proposed before and has shown its potential, particularly in one of our previous papers. With this protocol, we want to encourage more users to try this technology thanks to its simplicity. The novelty of our work is that we facilitate image-to-world registration, which has been previously considered a limiting factor. We have added a paragraph in the introduction section of the paper to better answer the reviewer's concerns.

2. We have previously evaluated the accuracy of the proposed protocol on a phantom in reference¹⁹, where we used HoloLens as the AR device. The protocol described in the paper is not intended to introduce any accuracy measure, but the reader can check our reference to better understand the limitations of the protocol. The limitations regarding the use of a marker attached to the bone were also identified in our previous paper. The referee is correct when pointing out that a marker attached to soft tissue would limit the utility of the protocol. We have added a paragraph in the discussion section with more clinical examples of the use of the proposed methodology: "In general, one of the most interesting clinical applications of AR and 3DP is to improve patient to physician communication by giving the patient a different perspective of her case, improving the explanation of a specific medical condition or treatment. Other possible application includes surgical guidance for target localization, where 3D printed patient-specific tools (with a reference AR marker attached) can be placed on rigid structures, such as bone, and used as a reference for the navigation. This application can be especially useful for orthopedic and maxillofacial surgical procedures in which bone tissue surface is easily accessed in the surgical field."
3. In order to describe the differences between Visualization and Registration methods, we have improved the description of the 3D Slicer module. In addition, the differences between both methods have been further explained in protocol sections.
4. Segmentation might be a difficult task for some diseases, even more, if the region of interest includes other tissue than bone. To answer the reviewer's concerns, we have added an explanation in section 2.2.2.1 of how to perform the segmentation step. We think that the whole protocol is now easier to follow and, although it may seem complex at some points, the document will be a good step-by-step guide for the reader when combined with the final video.
5. We apologize for the problems experienced when testing the version of the code available during the review. The 3D Slicer module has been updated and new functionalities have been added. One of the main improvements is the addition of a new button which allows the user to restart the process without the need of restarting 3D slicer when the button is pressed.

Minor Concerns:

1. Previous research proposed an AR method using a mobile device. You can refer to this article.
2. Fan, et al.: *3D Interactive Surgical Visualization System using Mobile Spatial Information Acquisition and Autostereoscopic Display*, *Journal of Biomedical Informatics*, Vol.71, pp.154-164, 2017.
2. You can mark some important components in figures for better understanding. For example, you can mark "3D CT model" and "3D printed marker reference".
3. A flowchart of the proposed method can be added for better understanding.

Answer:

The following paragraphs include our responses to the reviewer's minor concerns:

1. We would like to thank the reviewer for identifying that reference, that has been included in the introduction section of our paper (reference number 10).
2. Considering the reviewer comment, Figures 2 and 3 have been modified by remarking the different elements in the image.
3. A flowchart could help in understanding the method. However, we believe that the final video will serve as a flowchart and adding an extra figure illustrating the steps of the protocol will extend the paper length, so we did not include a flowchart.

References

19. Moreta-Martínez, R., García-Mato, D., García-Sevilla, M., Pérez-Mañanes, R., Calvo-Haro, J., Pascau, J. Augmented reality in computer-assisted interventions based on patient-specific 3D printed reference. *Healthcare Technology Letters*. (2018).