

HOW TO PREPARE ALIQUOTS FOR FIBROBLAST REPROGRAMMING

- STEP 1

PREPARE 9 ALIQUOTS OF **TOTAL NM-RNA-REPROGRAMMING COCKTAIL** MIXING IN EACH TUBE:

32 μL OSKMNL NM-RNA
+
24 μL EKB NM-RNA
+
5.6 μL NM-microRNAs
Tot. volume **61.6 μL**

x 9



- STEP 2

DIVIDE EACH OF THE NINE **61.6 μL** ALIQUOTS INTO FOUR TO OBTAIN A TOTAL OF **36** ALIQUOTS OF **15.4 μL**

