

Journal of Visualized Experiments

Controlled odor mimic permeation systems for olfactory training and field-testing --Manuscript Draft--

Article Type:	Invited Methods Article - JoVE Produced Video
Manuscript Number:	JoVE60846R1
Full Title:	Controlled odor mimic permeation systems for olfactory training and field-testing
Section/Category:	JoVE Chemistry
Keywords:	odor delivery, odor detection, olfaction, canine training, odor mimic, canine detection
Corresponding Author:	Lauryn tube DeGreeff, Ph.D. US Naval Research Laboratory Washington, DC UNITED STATES
Corresponding Author's Institution:	US Naval Research Laboratory
Corresponding Author E-Mail:	lauryn.degreeff@nrl.navy.mil
Order of Authors:	Lauryn E. DeGreeff, Ph.D. Alison G. Simon Michael S. Macias Howard K. Holness Kenneth G. Furton
Additional Information:	
Question	Response
Please indicate whether this article will be Standard Access or Open Access.	Standard Access (US\$2,400)
Please indicate the city, state/province, and country where this article will be filmed . Please do not use abbreviations.	Miami, Florida, United States



September 30, 2019

To whom it may concern,

Please find enclosed an original manuscript for consideration and peer review with *JoVE*. The manuscript, entitled ***Controlled odor mimic permeation systems for olfactory training and field-testing***, discusses a low-cost and field-portable option for the controlled delivery of odors for olfactory testing and training. The Controlled Odor Mimic Permeation System (COMPS) was developed to provide a field-appropriate method of odor delivery at controlled and reproducible rates. The manuscript will discuss the preparation of the COMPS and ways to control the odor quantity released, as well as methods of analysis including gravimetric analysis and gas chromatography / mass spectrometry. We thank you in advance for your consideration of this manuscript.

Sincerely,

Lauryn E. DeGreeff

Research Chemist
Chemistry Division, Code 6181
U.S. Naval Research Laboratory, 4555 Overlook Ave SW, Washington, DC 20375
Email: lauryn.degreeff@nrl.navy.mil
Phone: (202)767-4256

TITLE:**Controlled Odor Mimic Permeation Systems for Olfactory Training and Field Testing****AUTHORS AND AFFILIATIONS:**

Lauryn E. DeGreeff^{1,2}, Alison G. Simon^{3,4}, Michael S. Macias³, Howard K. Holness², Kenneth G. Furton²

¹Chemistry Division, US Naval Research Laboratory, Washington, DC, USA

²International Forensic Research Institute, Florida International University, Miami, FL, USA

³Formerly of International Forensic Research Institute, Florida International University, Miami, FL, USA

⁴Formerly of National Research Council Post-Doctoral Fellowship Program, US Naval Research Laboratory, Washington, DC, USA

Corresponding Author:

Lauryn E. DeGreeff (lauryn.degreeff@nrl.navy.mil or lauryn.degreeff@fiu.edu)

Email Addresses of Co-authors:

Alison G. Simon (alison.tcu@gmail.com)

Michael S. Macias (michael.s.macias@gmail.com)

Howard K. Holness (holnessh@fiu.edu)

Kenneth G. Furton (furtonk@fiu.edu)

KEYWORDS:

odor delivery, odor detection, olfaction, canine training, odor mimic, canine detection

SUMMARY:

The Controlled Odor Mimic Permeation System is a simple, field-portable, low-cost method of odor delivery for olfactory testing and training. It is constructed of an odorant retained on an adsorbent material and contained inside of a permeable polymer bag allowing controlled release of the odorant vapor over time.

ABSTRACT:

The Controlled Odor Mimic Permeation System (COMPS) was developed to provide a convenient field testing method of odor delivery at controlled and reproducible rates. COMPS are composed of an odorant of interest on an absorbent material sealed inside of a permeable polymer bag. The permeable layer allows for a constant release of the odorant over a given amount of time. The permeable bag is further stored in a secondary, impermeable bag. The double-containment procedure allows for equilibration of the odorant from the permeable bag but within the impermeable outer layer, resulting in an instant and reproducible source of odorant vapor upon removal from the outer packaging. COMPS are used in both olfactory testing for experimental scenarios and for olfactory detection training, such as with detection canines. COMPS can be used to contain a wide range of odorants (e.g., narcotics powders) and provide a controlled release of the associated odorants. Odor availability from COMPS is expressed in terms of

permeation rate (i.e., the rate of the odorant vapor released from a COMPS per unit time) and is typically measured by gravimetric means. The permeation rate for a given mass or volume of odorant can be adjusted as needed by varying the bag thickness, surface area, and/or polymer type. The available odor concentration from a COMPS can also be measured by headspace analysis techniques such as solid phase microextraction with gas chromatography/mass spectrometry (SPME-GC/MS).

INTRODUCTION:

Olfaction is a crucial, yet often overlooked, sensing mechanism used by most animals. For many it is the main mechanism for locating food, finding a mate, or sensing danger¹. Furthermore, the olfactory capabilities some animals, most notably canines, are regularly exploited by humans for the detection of contraband (e.g., narcotics or explosives), or other objects of interest, such as missing persons, invasive species, or diseases^{2,3}. For canine detection research or other olfaction research topics, investigators often study the process of olfaction and the strengths and limitations of the olfactory system. As such, it is generally desirable to control the release of an odorant vapor into the environment to reproducibly deliver known quantities of odorant during testing. Failure to account for variations in odor availability due to factors such as vapor pressure or environmental effects often complicates data interpretation and applicability⁴. It is similarly desirable to provide an established quantity of odor during training scenarios for detection canines. For example, studies by Hallowell et al.⁵ and by Papet⁶ have indicated the importance of odor intensity in odor perception, and that altering the intensity of an odorant can affect how it is perceived alone or in a mixture.

In laboratory settings, the use of analytical equipment such as permeation tubes with controllable ovens, vapor generators, or olfactometers may be used to control odor delivery. However, this type of equipment is impractical for use during field testing and training scenarios⁴. The Controlled Odor Mimic Permeation System (COMPS) was developed as a simple, low-cost, and disposable method for controlled odor delivery requiring no external power. Therefore, they can easily be incorporated into a variety of different testing and training scenarios⁷. COMPS units are simply composed of an odorant of interest on an absorbent material sealed inside of a permeable polymer bag, stored in a secondary containment system. The utilization of COMPS reduces variability between tests and improves consistency during training exercises⁸.

Odor delivery or availability from COMPS is measured in terms of permeation rate, as determined by gravimetric analysis in terms of mass of vapor released over time. Permeation rates can be controlled by a number of factors, including the thickness of the polymer bag, its available surface area, the type of absorbent material (substrate) used, and the amount of the odorant. Permeation rate is constant for a given period of time (hours or days) depending on the odorant being used. This allows for minimal variability in odor delivery during testing or training. During storage, COMPS come to equilibrium within the impermeable outer container, resulting in an instant source of odorant vapor at a known permeation rate.

COMPS were initially designed to contain odorants associated with explosive materials and to be used as odor mimics⁷. As defined by Macias et al., an odor mimic simulates a material of interest,

such as an explosive, by providing the dominant volatile compounds, or odorants, found in the headspace of that material without the presence of the parent material itself⁸. To create an odor mimic, the active odorants of the parent material must be determined. An active odorant, in this scenario, is described as a volatile compound that a trained explosive-detection canine detects, believing that there is an actual explosive material present. Having identified dominant volatile compounds in the headspace of several explosive materials, COMPS were prepared to release these individual odorants at a controlled rate for the duration of canine olfactory detection field trials and determine the active odorant associated with several explosive materials. COMPS were successfully used for this purpose^{7,9} and have since been used as odor mimics for further explosive detection training.

Macias et al. utilized COMPS containing piperonal, a pure chemical solid at room temperature that, in the vapor phase, has been shown to be the active odorant for MDMA (3,4-methylenedioxymethamphetamine), the psychoactive drug known as ecstasy. The researchers used varying thicknesses and surface areas of low-density polyethylene bags to adjust the permeation rate of piperonal vapor. This series of COMPS was then used to estimate piperonal detection threshold for trained narcotics-detection canines⁸. Conversely, in a separate study, COMPS bag thicknesses were adjusted to minimize the deviation of permeation rates between each compound in a homologous series though they possessed drastically varying vapor pressures. If a single bag thickness had been used in this study, those compounds with higher vapor pressures would have yielded much higher permeation rates. By increasing the bag thickness for the higher volatility compounds, the permeations rates were adjusted so that they were similar for all compounds⁴. Both studies demonstrate the utility and adaptability of the COMPS to control vapor release. Similar studies optimizing polymer bag thickness as well as absorbent material have been carried out in the creation of odor mimics for synthetic cathinones (i.e., bath salts)¹⁰, other narcotics (including heroin and marijuana¹¹), and human odor compounds^{12,13}. In a final example, Simon et al. investigated the active odorants associated with an invasive fungus species¹⁴. Whole pieces of infected tree bark, instead of the extracted odorants, were placed directly into the polymer bag to control release during canine olfaction testing¹⁴. COMPS can be utilized for a variety of scenarios, and the protocols discussed herein were chosen to demonstrate the diversity of this tool.

PROTOCOL:

1. Assembly of COMPS (Figure 1)

1.1. For neat (liquid) compound on a substrate (Figure 1A)

1.1.1. To impregnate the substrate with odorant, use a calibrated pipette to add 5 μ L of neat compound to a 2 x 2 inch cotton gauze pad or other substrate of choice (see **Table of Materials**).

1.1.2. Fold the gauze pad in half and place this (or an alternative substrate material) into a 2 x 3 inch low-density polyethylene permeable bag. The suggested bag thicknesses are between 1 MIL, for the fastest permeation rate, to 8 MIL, for a slower permeation rate.

NOTE: Variations in absorbent materials, permeable bag size, polymer chemistry, and thickness can be used, but these changes affect the permeation rate of the odorants (see further discussion in the Results Section).

1.1.3. Immediately seal the polymer bag closed with a heat sealer, eliminating as much air from within the bags as possible.

1.1.4. Store the bag in an outer nonpermeable bag, or if it will be used immediately, place it in a clean weigh boat in a fume hood (**Figure 1B**).

1.2. For solid material, no substrate necessary (**Figure 1C**)

1.2.1. Weigh the desired amount of solid material, which may be a pure compound or actual target material, and place in a 2 x 3 inch low-density polyethylene (LDPE) permeable bag. Again, the suggested bag thickness ranges from 1 MIL to 8 MIL.

1.2.2. Immediately heat-seal the polymer bag closed, eliminating as much air from within the bag as possible, and store or set aside in a weigh boat.

2. Gravimetric analysis to determine COMPS permeation rate

NOTE: A constant ambient temperature is important for accurate and reproducible measurements, both gravimetric and headspace. A constant temperature must be maintained during all testing. It is recommended to carry out all analytical measurements at desired temperatures during testing.

2.1. To determine the permeation rate of the odorants through the permeable bag, place a newly made COMPS in a weigh boat inside of a fume hood.

2.2. Place a clean, separate weigh boat on an analytical balance, and zero the balance.

2.3. Remove the COMPS from the fume hood and place on the balance. Record the mass and immediately return to the fume hood.

2.4. Continue recording the mass of the COMPS over regular time increments until the mass of the COMPS no longer changes ($\pm 5\%$). At this point the odor from the COMPS is depleted.

2.5. As a negative control, create an empty COMPS consisting of just the substrate material without odorants sealed into the permeable bag. Treat this negative control in the same manner as the COMPS with odorant to ensure minimal fluctuations in mass over time.

2.6. Calculate the permeation rate from the COMPS.

2.6.1. Plot the mass of the COMPS versus time on an X-Y plot in appropriate statistical analysis software.

2.6.2. Fit a linear trendline to just the linear portion of the graph and display an equation on a chart. The trendline should NOT be set to include the origin. The slope of the line (i.e., m in $y = m + b$) is the permeation rate in mass per unit time.

3. Headspace analysis by solid phase microextraction with gas chromatography/mass spectrometry (SPME-GC/MS) (optional)

3.1. Prepare a fresh COMPS following the instructions above and allow it to equilibrate in an open weigh boat inside of a fume hood for 30 min.

3.2. Remove the COMPS from the weigh boat, place it into a 1 pint epoxy-lined metal sample container without a lid, and put it into a 1 gallon epoxy-lined metal container. The containers should be kept in a fume hood for the duration of the experiment.

3.3. Allow at least 30 min for equilibration in the container prior to sampling.

3.4. For sampling after equilibration, place a lid with a previously drilled 1 cm hole atop the outer container. Insert an appropriate SPME fiber through the hole on the lid to extract the analyte of interest. When the SPME fiber is not used, cover the hole with paraffin film or the like. The extraction time and fiber coating will be dependent on the type and amount of analyte vapor present, as well as the size of the sampling vessel and environmental conditions¹⁵.

3.5. Remove the SPME fiber after the allotted extraction time and place in the heated inlet of a GC/MS for thermal desorption and analysis.

3.5.1. Run the GC/MS method appropriate for the compound used in the COMPS¹⁶.

3.6. For quantitation, compare the resulting peak area to an external calibration curve¹⁶ and/or internal standard¹⁷ as appropriate for the method and experimental design.

NOTE: 1) In this example, epoxy-lined metal sample containers were used, but other types of containers would also be suitable. To directly compare the odor availability to field olfactory evaluations, it would be best to use the same container, cleaned between each test, for both experiments; 2) For reproducible results, all aspects of the sampling procedure should be maintained in all replicate experiments, including but not limited to equilibration time, SPME extraction time, container type and size, and environmental conditions (i.e., temperature and humidity).

4. COMPS storage

4.1. Place a single COMPS in a metalized barrier bag (3.5 x 4.5 inch) and heat-seal to close,

removing as much air as possible from the bag prior to sealing (Figure 1B).

4.2. Store in cool ambient or refrigerated conditions but not below or close to freezing to avoid the formation of condensation as the COMPS thaws.

4.3. If testing multiple odorants or odorant delivery rates in a single experiment, secondary containment is recommended to eliminate any possible cross contamination during transportation and storage.

4.3.1. Place replicate multiple barrier bags each containing individual COMPS of the same analyte and permeation rate in an outer, larger metalized bag or glass jar for storage and transport.

5. Field olfactory testing

NOTE: Olfactory testing can be carried out in many different ways depending on the animal being tested, the goal of the experiment, and the environmental conditions. The protocol below describes one such manner of testing. All animal testing should first be reviewed and approved by an Institutional Animal Care and Use Committee (IACUC).

5.1. First, create blank or negative control COMPS as described above. Make enough so every container in the testing scenario will contain spare COMPS (2–3 depending on the number of animals involved in the experiment). Package all blank COMPS together in secondary containment (i.e., a larger metalized bag or glass jar with sealing lid).

5.2. Create fresh COMPS as needed for the intended field test protocol. Eliminate all possible sources of contamination between the COMPS and the metalized bag. This can be accomplished by regular changing of gloves and cleaning of the laboratory work surface.

5.3. Store the COMPS for at least 1 day prior to use to allow for equilibration. Store any replicates in the same secondary container. However, different COMPS should be in separate secondary containers.

5.4. To set up a basic canine olfactory test, lay out multiple lines of at least five identical containers (e.g., metal cans, boxes), with the number of lines dependent on the number of variables being tested.

5.4.1. Set up the trial so that each line contains one container with the target COMPS and four with blank COMPS. Positive control lines, prepared in the same manner but with a known target odor, may be used as appropriate for the experiment, training, or testing scenario. An additional negative control or blank line should contain five blank COMPS and no targets. Order this negative control line, positive control line (if using), and testing lines at random, and change order using a random number generator for each canine olfactory test as practical for the testing scenario.

5.4.1.1. Include one distractor odorant/material per line as well.

5.4.1.2. Randomize the order and location of the target and distractor odorants in each line for each canine being tested using a random number generator.

5.5. To prepare the containers, remove COMPS from the secondary and outer containers, placing only the permeable bag in the trial container.

5.5.1. Allow the COMPS to equilibrate in the container for a minimum of 30 min prior to testing.

5.5.2. Repeat steps for each COMPS being used in the test, starting with blanks, followed by positive controls (if using), and then testing odorants, changing gloves each time.

NOTE: Detailed examples of canine testing scenarios can be found in Simon et al.⁴ or Macias et al.⁸.

REPRESENTATIVE RESULTS:

The primary objective of using COMPS in olfactory testing/training is to control the release of the chosen odorants and deliver a controlled amount of the odorant over the duration of the test or training session. Odorant release is measured by gravimetric analysis in terms of mass loss per unit time. **Figure 2** gives an example of gravimetric results from the permeation of three identical COMPS prepared from 5 μ L of pentanoic acid on cotton gauze through a 3 MIL LDPE bag. A line of regression was added to the graph, and the slope of the line represents the permeation rate of 37 μ g/min for this set of COMPS.

It is often desirable to be able to adjust the amount of odorant released for a given test. This can be done in several ways, including adjusting the amount of material in the bag, the surface area of the permeable bag material, or the bag thickness. **Figure 3** shows how all three of these factors were used to control the release of piperonal. **Figure 3A** indicates a logarithmic relationship between the mass in piperonal in the permeable bag (3 x 3 inch, 2 MIL LDPE), where the permeation rate increased quickly at the lower masses, then slowed after 500 mg due to physical restriction in the amount of odorant that could be released from the given bag at a time. The data in **Figure 3B** depict a linear relationship between permeation rate and surface area of the permeable bag for 2 g of piperonal in a 2 MIL LDPE bag. Finally, permeation rate decreased linearly with increased bag thickness (2 g of piperonal in a 3 x 3 inch bag), as shown in **Figure 3C**, because the thicker bag restricts and slows emission.

In another example of the utility of controlling permeation rates, Simon et al.⁴ used bag thickness to standardize permeation rates for compounds of varying vapor pressures in order to present canines with similar odorant availability for each analyte during field testing. A volume of 5 μ Ls of each neat analyte was pipetted onto separate cotton gauze pads and placed into 2 x 3 inch LDPE permeable bags. The permeation rates were measured by gravimetric analysis. **Figure 4** shows the variation in vapor pressures (**Figure 4A**) across the groups of analytes (RSD = 138%) compared to the variation in permeation rate after adjusting the bag thickness (**Figure 4B**) to

control the rate and make them as similar as possible (RSD = 31.8%). Furthermore, adjusting the bag thickness allowed permeation rates to vary by three orders of magnitude (**Table 1**).

Headspace measurements can be used to better measure the amount of odorant available during a given testing or training scenario. Macias¹⁸ measured the amount of piperonal in the headspace of three COMPS with permeation rates of 1,000, 100, and 10 ng/min (**Figure 5**). The COMPS were placed in a 1 quart sampling can, and the headspace was extracted for 30 min using SPME. The resulting chromatograph in **Figure 5** shows the piperonal peak areas increasing with increasing permeation rate¹⁸.

Macias then used these three sets of piperonal COMPS in canine trials. Trained narcotics detection canines were tested on the 0 (blank), 10, 100, and 1,000 ng/s piperonal COMPS in a scent cage (**Table 2**). The results showed that as the permeation rate, and thus odor availability, increased the number of canines alerting to the appropriate COMPS increased¹⁸.

FIGURE AND TABLE LEGENDS:

Figure 1: Examples of COMPS. (A) A COMPS constructed from a cotton gauze pad in permeable polymer bag. Reproduced from Simon et al.⁴ (B) A COMPS inserted into an outer impermeable bag. (C) A COMPS containing infected wood as an odor source in a polymer bag. Figures B and C were reproduced with permission from Simon et al.¹⁹.

Figure 2: Example measurement of permeation rate by gravimetric analysis. The mass loss of the analyte (i.e., pentanoic acid) on gauze through a 3 MIL LDPE bag measured over time. A, B, and C indicate replicates of the same material, while "Average" is the averaged value of the three at each time point. The given equation describes the linear fit to the average data.

Figure 3: Examples of factors adjusting the permeation rate. Graphs of piperonal dissipation indicate experimentally measured permeation rates when changing (A) the mass of piperonal (3 x 3 inch, 2 MIL LDPE bag), (B) the surface area of the permeable bag (2 g of piperonal, 2 MIL LDPE), and (C) the bag thickness (2 g of piperonal, 3 x 3 inch bag). All error bars represent one standard deviation of the mean (some bars are within the size of the marker). These figures have been reproduced with permission from Macias et al.¹⁸.

Figure 4: Comparison of vapor pressure variation across a group of compounds versus permeation rate variation. (A) Vapor pressures for a selection of 12 compounds (RSD = 138%) compared to (B) permeation rates for the same compounds with selected COMPS thicknesses (RSD = 31.8%). Numbers in parentheses represent LDPE bag thickness in MIL. These figures have been reproduced with permission from Simon et al.⁴.

Figure 5: Headspace analysis of piperonal COMPS at three permeation rates. Overlaid chromatograms of the headspace components of piperonal COMPS adjusted to permeation rates of 1,000, 100, and 10 ng/s. Reproduced with permission from Macias et al.¹⁸.

Table 1: Permeation rate versus bag thickness. Permeation rate of methyl benzoate on cotton gauze in COMPS of varying bag thicknesses. Note that the lowest permeation rate was achieved

by placing a 4 MIL COMPS inside of a metalized bag with a 1/8 in hole. The R^2 value indicates the fit of the line to the gravimetric plot.

Table 2: Example of canine field trial results. Canine responses to piperonal COMPS with permeation rates ranging from 0–1,000 ng/s. Reproduced with permission from Macias et al.¹⁸.

DISCUSSION:

Controlled Odor Mimic Permeation Systems (COMPS) are easily created by sealing an odorant of interest into a permeable bag. This may be done by pipetting a neat liquid compound onto an absorbent material and then placing the absorbent material into the bag; by placing a pure, solid compound directly into the bag⁴, as was done in the case of piperonal⁸; or by placing the target material containing multiple or unknown odorants into a permeable bag, as was done with fungus-infected wood¹⁴. The permeable bag controls the release of the odorant so a known and reproducible amount can be delivered over a given training or testing period. Permeation rate is typically measured by gravimetric analysis, plotting mass loss over time, and can be adjusted by altering a number of parameters, including the absorbent material, the mass/volume of odorant, or parameters of the permeation bag (i.e. thickness, surface area, or polymer type). COMPS are stored in an outer nonpermeable envelope, which allows the COMPS to equilibrate prior to use, thus providing a known amount of odorant immediately upon use.

The greater the permeation rate of a COMPS, the greater concentration of odorant available during a training or testing scenario. To quantitate or compare the odor concentration emitted from a COMPS, headspace analysis of the COMPS in the testing/training container should be completed. This is most often done by extracting the odor using SPME with analysis by GC/MS. For quantitation or comparison purposes, it is recommended to use an internal standard and/or an external calibration curve.

COMPS serve as a low-cost, field-amenable devices for controlling the release of an odorant for olfactory training or testing, such as with canine detectors. COMPS can be used repeatedly until depletion, each time delivering the same odorant emission rate, although the length of time the emission rate is constant will change for each analyte and should be tested in the laboratory prior to use. This overcomes a widely recognized limitation of controlling odor delivery for field use and advances olfaction research and detection animal training.

ACKNOWLEDGMENTS:

This work was funded in part by the Office of Naval Research and the National Institute of Justice (2006-DN-BX-K027). The authors wish to thank the many “Furton Group” students that have participated in this project, as well as collaborators from the U.S. Naval Research Laboratory and the Naval Surface Warfare Center (Indian Head EOD Technology Division). Finally, the authors thank Peter Nunez of U.S. K-9 Academy, Tony Guzman of Metro-Dade K9 Services, and Miami-Dade area law enforcement canine teams.

DISCLOSURES:

No conflicts of interest to report.

REFERENCES:

1. Buck, L., Axel, R. A novel multigene family may encode odorant receptors: A molecular basis for odor recognition. *Cell*. **65**, 175–187 (1991).
2. Furton, K. G., Myers, L. J. Scientific foundation and efficacy of the use of canines as chemical detectors for explosives. *Talanta*. **54**, 487–500 (2001).
3. Leitch, O., Anderson, A., Kirkbride, K., Lennard, C. Biological organisms as volatile compound detectors: A review. *Forensic Science International*. **232**, 92–103 (2013).
4. Simon, A. G. et al. method for controlled odor delivery in canine olfactory testing. *Chemical Senses*. **44** (6), 399–408 (2019).
5. Hallowell, L. R. et al. Detection of hidden explosives: New challenges and progress (1998–2009). in *Forensic Investigation of Explosives*, 2nd Ed. Boca Raton, FL, CRC Press, 2012, pp.53–77.
6. Papet, L. Narcotic and explosive odors: Volatile organic compounds as training aids for olfactory detection. in *Canine Olfaction Science and Law*, Boca Raton, FL, CRC Press, 2016, pp. 265–278.
7. Furton, K., Harper, R. Controlled Odor Mimic Permeation System. U.S. Patent 20080295783A1, 7 July 2017.
8. Macias, M. S., Guerra-Diaz, P., Almirall, J. R., Furton, K. G. Detection of piperonal emitted from polymer controlled odor mimic permeation systems utilizing canis familiaris and solid phase microextract-ion mobility spectrometry. *Forensic Science International*. **195**, 132–138 (2010).
9. Harper, R., Almirall, J., Furton, K. Identification of dominant odor chemicals emanating from explosives for use in developing optimal training aid combinations and mimics for canine detection. *Talanta*. **67**, 313–327 (2005).
10. Francis, V. S. The identification of volatile organic compounds from synthetic cathinone derivatives for the development of odor mimic training aids. Florida International University, Miami, FL, 2017.
11. Huertas-Rivera, A. M. Identification of the active odors from illicit substances for the development of optimal canine training aids. Florida International University, Miami, FL, 2016.
12. DeGreeff, L. E., Furton, K. G. Collection and identification of human remains volatiles by non-contact, dynamic airflow sampling and SPME-GC/MS using various sorbent materials. *Analytical and Bioanalytical Chemistry*. **401**, 1295–1307 (2011).
13. DeGreeff, L. E., Curran, A. M., Furton, K. G. Evaluation of selected sorbent materials for the collection of volatile organic compounds related to human scent using non-contact sampling mode. *Forensic Science International*. **209** (1–3), 133–142 (2011).
14. Simon, A. G., Mills, D. K., Furton, K. G. Chemical and canine analysis as complimentary techniques for the identification of active odors of the invasive fungus, *Raffaelea lauicola*. *Talanta*. **168**, 320–328 (2017).
15. Penton, Z. Method development with solid phase microextraction. in *Solid Phase Microextraction: A Practical Guide*, New York, NY, Marcel Dekker, Inc., 1999, pp. 27–58.
16. Robards, K., Haddad, P. R., Jackson, P. E. *Principles and Practice of Modern Chromatographic Methods*, San Diego, CA: Elsevier Ltd., 2004.
17. MacCrehan, W., Moore, S., Schantz M. Evaluating headspace component vapor-time profiles by solid-phase microextraction with external sampling of an internal standard. *Analytical Chemistry*. **83**, 8560–8565 (2011).

- 436 18. Macias, M. S. The Development of an Optimized System of Narcotic and Explosive
437 Contraband Mimics for Calibration and Training of Biological Detectors. Dissertation, Florida
438 International University, Miami, FL, 2009.
- 439 19. Simon, A. G. The Detection of an Invasive Pathogen through Chemical and Biological Means
440 for the Protection of Commercial Crops. Dissertation, Florida International University, Miami, FL,
441 2017.
- 442

Bag thickness	Permeation	R-squared
	rate (mg/min)	
1 MIL	0.108	0.974
4 MIL	0.042	0.991
8 MIL	0.00499	0.99
4 MIL in metal b	0.000179	0.972

Piperonal COMPS permeation rate	Number of alerts	% alert
0 ng/s (blank)	0	0%
10 ng/s	4	25%
100 ng/s	7	44%
1000 ng/s	12	75%

Name of Material/ Equipment	Company	Catalog Number	Comments/Description
16 oz economy jars (70-450 finish)	Fillmore container	A16-08C-Case 12	
7890A gas chromatograph / 5975 mass selective detector	Agilent		
Analytical balance	Mettler Toledo	01-911-005	
Ball regular bands and dome lids	Fillmore container	J30000	
Cotton gauze (2" x 2")	Dukal		
Disposable weighing boats	VWR	10803-148	
Epoxy-lined sample containers, 1 gallon	TriTech Forensics	CANG-E	
Epoxy-lined sample containers, 1 pint	TriTech Forensics	CANPT-E	
Low density polyethylene bag	Uline	S-5373	
Rtx-Volatiles (30 m x 0.32 mmID) column	Restek	10901	
Silver metalized mylar barrier bag (3.5" x 4.5")	ESP Packaging	95509993779	
Silver metalized mylar barrier bag (5" x 8.5" x 3")	ESP Packaging	95509993793	
Solid phase microextraction fiber assembly (PDMS/DVB/CAR)	Sigma-Aldrich	57328-U	
Solid phase microextraction holder	Sigma-Aldrich	57330-U	
Tabletop Impulse Sealer	Uline	H-190	Heat sealer

Editorial comments:

The manuscript has been modified and the updated manuscript, 60846_R0.docx, is attached and located in your Editorial Manager account. Please use the updated version to make your revisions.

1. Please take this opportunity to thoroughly proofread the manuscript to ensure that there are no spelling or grammar issues.

2. Please obtain explicit copyright permission to reuse any figures from a previous publication. Explicit permission can be expressed in the form of a letter from the editor or a link to the editorial policy that allows re-prints. Please upload this information as a .doc or .docx file to your Editorial Manager account. The Figure must be cited appropriately in the Figure Legend, i.e. "This figure has been modified from [citation]." **Complete**

3. Please provide at least 6 keywords or phrases. There are only 5. **Complete**

4. Please use superscript arabic numerals to cite references in text. The superscript number is inserted immediately next to the word/group of words it applies to but before any punctuation. Please remove all brackets before and after the reference numbers. **Complete**

5. Step 5.4.1: Please ensure that all text is written in the imperative tense. **Complete**

6. Step 5.4.2.1: Please ensure that all text is written in the imperative tense. **Complete**

7. Step 5.4.2.2: Please ensure that all text is written in the imperative tense. **Complete**

8. Step 5.6: Please ensure that all text is written in the imperative tense. **Complete**

9. Please revise the Discussion to explicitly cover the following in detail in 3-6 paragraphs with citations:

a) Critical steps within the protocol **Already included in first and second paragraphs.**

b) Any modifications and troubleshooting of the technique **Options for COMPS creation included in first paragraph.**

c) Any limitations of the technique **Additional sentence added the last paragraph.**

d) The significance with respect to existing methods **Included in final paragraph.**

e) Any future applications of the technique **Not relevant.**

10. All tables should be uploaded separately to your Editorial Manager account in the form of an .xls or .xlsx file. **Complete**

11. Please sort the items in alphabetical order according to the name of material/equipment. **Complete**

Reviewers' comments:**Reviewer #1:**

Overall, I found this manuscript a good overview and description of the use of COMPS for canine testing. I think it will have several uses for a variety of scientists and canine handlers as well. I am happy to see this move forward.

I only have a few minor comments:

Line 107- not a clear description of the study. I think a little too brief for those not familiar with the recent paper. **An additional sentences describing the findings were added.**

Line 157- should there be any control for temperature? Temperature effects are generally not noted throughout. **This is a good point. A note was added in regards to handling temperature.**

Line 188- the description here is confusing. The video should correct that, but a little more clarity as to whether you are replacing the lid, or whether you are using a lid that can be pierced or plugged when the hole is not needed? **Additional information was added.**

Line 255- I wouldn't recommend making the second search always an all blank search. Over time, dogs will readily learn the second search is blank and will stop attending to that search. It should be random as to whether the first or second search doesn't contain the target. **The authors agree. This was misstated in the text and not the authors' intention. This was rephrased for clarity.**

Reviewer #2: No changes requested.

Manuscript Summary:

This manuscript describes the process used by "COMPS" to make controlled odor release systems to be used in experiments and training. It discusses this system for odor as a liquid placed on material, as well as solids. Authors discussed clearly how to measure the rate of the dissipation of the odor and how to vary the rate for different targets, if wanted or needed. The manuscript specifically discusses its use with canines, however, any odor detection research could use this system. The methods appear to be sufficient for someone to repeat and use this method in their odor detection research. Their COMPS system is backed up by other research, cited and briefly discussed in the manuscript.

Major Concerns:

None

Minor Concerns:

None

Figure 1.




Figure 2.




Figure 3.




Figure 4.




Figure 5.




Table 1.

Bag thickness	Permeation rate (mg/min)	R-squared
1 MIL	0.108	0.974
4 MIL	0.042	0.991
8 MIL	0.00499	0.990
4 MIL in metal bag w/ 1/8" hole	0.000179	0.972

Table 2.

Piperonal COMPS permeation rate	Number of alerts	% alert
0 ng/s (blank)	0	0%
10 ng/s	4	25%
100 ng/s	7	44%
1000 ng/s	12	75%

OXFORD UNIVERSITY PRESS LICENSE TERMS AND CONDITIONS

Nov 25, 2019

This Agreement between U.S. Naval Research Laboratory -- Lauryn DeGreeff ("You") and Oxford University Press ("Oxford University Press") consists of your license details and the terms and conditions provided by Oxford University Press and Copyright Clearance Center.

License Number	4716060797398
License date	Nov 25, 2019
Licensed Content Publisher	Oxford University Press
Licensed Content Publication	Chemical Senses
Licensed Content Title	A Method for Controlled Odor Delivery in Olfactory Field-Testing
Licensed Content Author	Simon, Alison G; DeGreeff, Lauryn E
Licensed Content Date	May 21, 2019
Licensed Content Volume	44
Licensed Content Issue	6
Type of Use	Journal
Requestor type	Author of this OUP content
Pharmaceutical support or sponsorship for this project	No

Format	Print and electronic
Portion	Figure/table
Number of figures/tables	1
Will you be translating?	No
Circulation/distribution	1
Title of new article	controlled odor mimic permeation systems for olfactory training and field-testing
Lead author	Lauryn DeGreeff
Title of targeted journal	Journal of Visualized Experiments
Publisher	JoVE
Expected publication date	Mar 2020
Portions	Figure 4, page 404
Requestor Location	U.S. Naval Research Laboratory U.S. Naval Research Laboratory 4555 Overlook Ave. SW Washington, D.C. 20375, USA WASHINGTON, DC 20002 United States Attn: U.S. Naval Research Laboratory
Publisher Tax ID	GB125506730
Total	0.00 USD
Terms and Conditions	

STANDARD TERMS AND CONDITIONS FOR REPRODUCTION OF MATERIAL FROM AN OXFORD UNIVERSITY PRESS JOURNAL

1. Use of the material is restricted to the type of use specified in your order details.
2. This permission covers the use of the material in the English language in the following territory: world. If you have requested additional permission to translate this material, the terms and conditions of this reuse will be set out in clause 12.
3. This permission is limited to the particular use authorized in (1) above and does not allow you to sanction its use elsewhere in any other format other than specified above, nor does it apply to quotations, images, artistic works etc that have been reproduced from other sources which may be part of the material to be used.
4. No alteration, omission or addition is made to the material without our written consent. Permission must be re-cleared with Oxford University Press if/when you decide to reprint.
5. The following credit line appears wherever the material is used: author, title, journal, year, volume, issue number, pagination, by permission of Oxford University Press or the sponsoring society if the journal is a society journal. Where a journal is being published on behalf of a learned society, the details of that society must be included in the credit line.
6. For the reproduction of a full article from an Oxford University Press journal for whatever purpose, the corresponding author of the material concerned should be informed of the proposed use. Contact details for the corresponding authors of all Oxford University Press journal contact can be found alongside either the abstract or full text of the article concerned, accessible from www.oxfordjournals.org Should there be a problem clearing these rights, please contact journals.permissions@oup.com
7. If the credit line or acknowledgement in our publication indicates that any of the figures, images or photos was reproduced, drawn or modified from an earlier source it will be necessary for you to clear this permission with the original publisher as well. If this permission has not been obtained, please note that this material cannot be included in your publication/photocopies.
8. While you may exercise the rights licensed immediately upon issuance of the license at the end of the licensing process for the transaction, provided that you have disclosed complete and accurate details of your proposed use, no license is finally effective unless and until full payment is received from you (either by Oxford University Press or by Copyright Clearance Center (CCC)) as provided in CCC's Billing and Payment terms and conditions. If full payment is not received on a timely basis, then any license preliminarily granted shall be deemed automatically revoked and shall be void as if never granted. Further, in the event that you breach any of these terms and conditions or any of CCC's Billing and Payment terms and conditions, the license is automatically revoked and shall be void as if never granted. Use of materials as described in a revoked license, as well as any use of the materials beyond the scope of an unrevoked license, may constitute copyright infringement and Oxford University Press reserves the right to take any and all action to protect its copyright in the materials.
9. This license is personal to you and may not be sublicensed, assigned or transferred by you to any other person without Oxford University Press's written permission.
10. Oxford University Press reserves all rights not specifically granted in the combination of (i) the license details provided by you and accepted in the course of this licensing

transaction, (ii) these terms and conditions and (iii) CCC's Billing and Payment terms and conditions.

11. You hereby indemnify and agree to hold harmless Oxford University Press and CCC, and their respective officers, directors, employs and agents, from and against any and all claims arising out of your use of the licensed material other than as specifically authorized pursuant to this license.

12. Other Terms and Conditions:

v1.4

Questions? customercare@copyright.com or +1-855-239-3415 (toll free in the US) or +1-978-646-2777.

From: [Jill Krefft](#)
To: [Lauryn E. Degreeff-silk, Code 6181](#)
Cc: "[mmacias22@hotmail.com](#)"; [Alison Simon](#)
Subject: RE: reproduction permission for publication from dissertation data
Date: Tuesday, November 19, 2019 11:26:42 AM

Hi Lauryn,

If Dr. Michael Macias and Dr. Alison Simon would like to use parts of their dissertation in a new publication they are co-authoring with others they are free to do so.

Authors who have a completed and published a dissertation at FIU retain all copyrights to their dissertation.

FIU is simply licensed to distribute the dissertation through our repository. If the publisher requires a letter, the letter should be from Dr. Macias and Dr. Smith since they are the copyright holders.

There is one caveat. If the data used in the dissertation is not their original data or used under some other license or agreement with a third party, they would need to honor that agreement. For example, we often have students who are working within a larger research project who may be restricted by the terms of the research project to share or use that data in other ways. Or another example would be if they retrieved the data from other published sources and used it with permission from another researcher they would need to request permission to use it in a different publication. Again, it is up to the authors of the dissertation to comply with those agreements.

Please let me know if you have additional questions.

Jill

-----Original Message-----

From: Lauryn E. Degreeff-silk, Code 6181 [<mailto:lauryn.degreeff@nrl.navy.mil>]
Sent: Monday, November 18, 2019 4:51 PM
To: Jill Krefft <jkrefft@fiu.edu>
Cc: 'mmacias22@hotmail.com' <mmacias22@hotmail.com>; Alison Simon <alison.tcu@gmail.com>
Subject: RE: reproduction permission for publication from dissertation data

Hello,

Thank you for getting back to me so quickly. We would like to derive the data from Dr. Michael Macias' dissertation as well as Dr. Alison Simon. They will also both be authored on the new manuscript in question. They are cc'd here for your convenience. After you receive their permission, will I be able to have something in writing that I can provide as proof to the publisher?

Thank you for your time.

Best,
Lauryn

Lauryn DeGreeff, Ph.D.
Research Chemist
Chemistry Division, Code 6181
U.S. Naval Research Laboratory
T 202.767.4256
www.nrl.navy.mil

-----Original Message-----

From: Jill Krefft <jkrefft@fiu.edu>

Sent: Monday, November 18, 2019 9:59 AM

To: Lauryn E. Degreeff-silk, Code 6181 <lauryn.degreeff@nrl.navy.mil>

Subject: RE: reproduction permission for publication from dissertation data

Hi Lauryn,

Which dissertation are you looking to derive the data from?

You will need permission from the authors to use the data, unless otherwise stated. I can forward your request on to the authors along with your contact information.

Best,

Jill

Institutional Repository Coordinator

FIU Libraries

<http://digitalcommons.fiu.edu>

-----Original Message-----

From: Lauryn E. Degreeff-silk, Code 6181 [<mailto:lauryn.degreeff@nrl.navy.mil>]

Sent: Friday, November 15, 2019 3:03 PM

To: Jill Krefft <jkrefft@fiu.edu>

Subject: reproduction permission for publication from dissertation data

Good afternoon,

The FIU librarian referred me to you. I am an FIU graduate school alumni (and current courtesy faculty member of the chemistry department). I am writing a new research paper co-authored with two other former FIU graduate school alumni - Dr. Alison Simon and Dr. Michael Macias. We would like to use data from Dr. Simon's and Dr. Macias' dissertation to make new figures for the new publication. Do I require permissions from FIU and/or the authors? If so, how can I go about obtaining permission?

Thank you very much.

Best regards,

Lauryn DeGreeff, Ph.D.

Research Chemist

Chemistry Division, Code 6181

U.S. Naval Research Laboratory

T 202.767.4256

www.nrl.navy.mil




Click here to access/download
Supplemental Coding Files
1.svg



Click here to access/download
Supplemental Coding Files
2.svg



Click here to access/download
Supplemental Coding Files
3.svg



Click here to access/download
Supplemental Coding Files
4.svg



Click here to access/download
Supplemental Coding Files
5.svg