Journal of Visualized Experiments

Engineering intracellular protein sensors in mammalian cells --Manuscript Draft--

Article Type:	Invited Methods Article - JoVE Produced Video		
Manuscript Number:	JoVE60878R3		
Full Title:	Engineering intracellular protein sensors in mammalian cells		
Section/Category:	JoVE Bioengineering		
Keywords:	Mammalian synthetic biology; protein sensors; sensing-actuating devices; HIV detection		
Corresponding Author:	Velia Siciliano Istituto Italiano di Tecnologia Center for Advanced Biomaterials for Healthcare Naples, Naples ITALY		
Corresponding Author's Institution:	Istituto Italiano di Tecnologia Center for Advanced Biomaterials for Healthcare		
Corresponding Author E-Mail:	veliasiciliano@gmail.com		
Order of Authors:	Marina Duhkinova		
	Claudia Crina		
	Ron Weiss		
	Velia Siciliano		
Additional Information:			
Question	Response		
Please indicate whether this article will be Standard Access or Open Access.	Standard Access (US\$2,400)		
Please indicate the city, state/province, and country where this article will be filmed . Please do not use abbreviations.	Naples, Italy		

1 TITLE:

Engineering Intracellular Protein Sensors in Mammalian Cells

2

- 4 AUTHORS:
- 5 Marina Dukhinova¹, Claudia Crina¹, Ron Weiss², Velia Siciliano¹

6

- 7 ¹Istituto Italiano di Tecnologia, Largo Barsanti e Matteucci, Naples, Italy
- 8 ²Synthetic Biology Center, Department of Biological Engineering, Massachusetts Institute of
- 9 Technology, Cambridge, MA, USA

10

- 11 Corresponding Author:
- 12 Velia Siciliano
- 13 velia.siciliano@iit.it

14

- 15 Email addresses of co-authors:
- 16 Marina Dukhinova (marina.dukhinova@gmail.com)
- 17 Claudia Crina (claudia.crina@iit.it)
- 18 Ron Weiss (rweiss@mit.edu)

19 20

KEYWORDS:

Mammalian synthetic biology, intracellular protein sensors, genetic sensor-actuators, HIV sensing, protein sensor platform, modular intracellular protein-sensing devices

2324

25

26

27

28

SUMMARY:

Here, we present a protocol for engineering genetically-encoded intracellular protein sensor-actuator(s). The device specifically detects target proteins through intracellular antibodies (intrabodies) and responds by switching on gene transcriptional output. A general framework is built to rapidly replace intrabodies, enabling rapid detection of any desired protein, without altering the general architecture.

293031

32

33

34 35

36

ABSTRACT:

Proteins can function as biomarkers of pathological conditions, such as neurodegenerative diseases, infections or metabolic syndromes. Engineering cells to sense and respond to these biomarkers may help the understanding of molecular mechanisms underlying pathologies, as well as to develop new cell-based therapies. While several systems that detect extracellular proteins have been developed, a modular framework that can be easily re-engineered to sense different intracellular proteins was missing.

373839

40

41

42

43

44 45

46

47

Here, we describe a protocol to implement a modular genetic platform that senses intracellular proteins and activates a specific cellular response. The device operates on intracellular antibodies or small peptides to sense with high specificity the protein of interest, triggering the transcriptional activation of output genes, through a TEV protease (TEVp)-based actuation module. TEVp is a viral protease that selectively cleaves short cognate peptides and is widely used in biotechnology and synthetic biology for its high orthogonality to the cleavage site. Specifically, we engineered devices that recognize and respond to protein-biomarkers of viral infections and genetic diseases, including mutated huntingtin, NS3 serine-protease, Tat and Nef proteins to detect Huntington's disease, hepatitis C virus (HCV)

and human immunodeficiency virus (HIV) infections, respectively. Importantly, the system can be hand tailored for the desired input-output functional outcome, such as fluorescent readouts for biosensors, stimulation of antigen presentation for immune response, or initiation of apoptosis to eliminate unhealthy cells.

INTRODUCTION:

The study and modulations of cellular responses via controllable engineered gene circuits are among the major challenges in current synthetic biology^{1–3}. These circuits represent prospective tools for relevant biological or medical applications in cancer⁴, infections⁵, or metabolic diseases⁶.

Reprogramming cell functions in response to specific signals requires the design of smart interfaces that link sensing of extracellular or intracellular dynamic changes (input) to downstream processing, triggering specific output either for diagnostic purposes (i.e., reporter genes) or to rewire cell response (therapeutics). The inputs detected by the sensing module can be small analytes⁷, proteins^{8–10} or microRNAs^{10–12}, specific for the onset or progression of a disease. Moreover, complex circuit regulation can be achieved by multiple input information processing, increasing the tight control over transgene expression in response to the defined conditions^{13–15}. For example, microRNA-based sensors can identify specific cell types, such as cancer cells, inducing their clearance with the expression of an apoptotic gene¹². Since microRNAs are easily implementable in synthetic circuits in a modular manner, they represent a widely used input for genetically encoded biosensors^{11,12,16}. Proteins are also a valid biomarker for genetic mutations, cancer and infections, and indeed a number of extracellular protein-sensing devices have been reported^{17,18}.

Many of the circuits that detect extracellular proteins rely on the use of engineered receptors, which tether a transcription factor (TF) to the membrane, fused to a TEVp-responsive cleavage site (TCS). A major advantage of the TEVp is the specificity of the cleavage and lack of interference with endogenous protein processing. In these systems, TEVp is fused to a second peptide that interacts with the engineered receptor upon binding of the extracellular molecules. Thus, the external inputs induce TEVp-mediated cleavage and TF release. Systems that function with this mechanism are Tango/TEVp¹⁷, light-induced¹⁹ and Modular Extracellular Sensor Architecture (MESA)¹⁸. Despite progress in detecting extracellular proteins, the technology for sensing intracellular proteins in a modular fashion was never realized before, with the limitation of going through many build-and-test-iterations for single devices responsive to a specific protein.

Our system is the first platform for intracellular protein sensing²⁰. The modularity is guaranteed by the use of intrabodies that define the specificity to the target, whereas the cell-reprogramming is TEVp-mediated. Specifically, one intrabody is membrane bound and fused to the C-terminal to a fluorescent protein mKate, a TCS and a TF (fusion protein 1); the second intrabody is fused to the TEVp and located in the cytosol (fusion protein 2) (**Figure 1**). Thus, the interaction between two intrabodies and the target protein occurs in the cytoplasm and leads to TCS cleavage by TEVp, resulting in TF translocation into the nucleus to activate functional output. The sensing—actuating device was successfully tested for four intracellular disease-specific proteins: NS3 serine protease expressed by the HCV virus²¹, Tat and Nef proteins from HIV infection^{22,23}, and mutated huntingtin (HTT) of the Huntington's disease²⁴.

Output expression includes fluorescent reporters, apoptotic gene (hBax)²⁵ and 95 immunomodulators (XCL-1)²⁶. We demonstrate that the system can also impair pathological 96 97 functionality of its targets. For instance, the Nef-responsive device interferes with the 98 spreading of viral infection by sequestering the target protein and reverting the downmodulation of HLA-I receptor on infected T cells²³. The described sensing-actuating 99 platform is the first of this kind for the detection of intracellular proteins and can be 100 potentially implemented to sense abnormal protein expression, post-translational or 101 epigenetic modifications, for diagnostic and therapeutic purposes²⁰. 102

_

PROTOCOL:

105106

103104

1. Design principles for construction and test the sensor-actuator device

107108

1.1. Select a protein of interest.

109

NOTE: We designed a system for proteins located in the cytoplasm or shuttling between the cytoplasm and other compartments.

112

1.2. Select two intrabodies binding different epitopes of the target protein. In our study we selected proteins for which the intrabodies were already developed and tested ^{27–30}. When intrabodies are not available, they should be newly developed.

116

NOTE: An alternative to intrabodies can be other molecules that interact with the protein of interest. For example, the Nef sensing device include one intrabody (sdAb19) and the SH3 domain of a p59fyn protein tyrosine kinase³¹.

120

121 1.3. Design in silico the plasmids with DNA sequences encoding the fusion proteins.

122

NOTE: Software for plasmid design are available for free or can be purchased.

123124

1.3.1. Generate sequences that include the following modules: membrane tag, fluorescent marker (e.g., mKate), intrabody, TCS and TF for fusion protein 1, and intrabody and TEVp for fusion protein 2.

128

NOTE: Building a new architecture based on the interaction of several components requires attention to parameters that may maximize the signal activation when the target protein is detected. Specifically, we focused on TEVp activity and flexibility of the chimeric proteins to facilitate protein-protein interaction (**Table 1**). All these features are specified in the following points.

134

1.3.2. Insert a flexible glycine-serine linker domain (LD) between the intrabody–TCS and/or TEVp–intrabody. The LD may be of variable length (0, 10 or 15 amino acids-AA).

- 138 1.3.3. Fusion protein 1: Design variants that include the high affinity TCS (S) or the low
- affinity TCS (L) to regulate TEVp cleavage activity. The expression of this construct is driven by a constitutive promoter. Include the GAL4VP16 transcription factor at the C-terminal of TCS
- to allow translocation to the nucleus upon TEVp cleavage.

142

NOTE: We used the hEF1 α promoter, but other constitutive promoters can be chosen.

144

1.3.4. Fusion protein 2: Use a constitutive (hEF1 α) or doxyclycine-inducible (pTET) promoter to tune TEVp transcription. Use a degradation domain (DD) regulated by the small molecule shield to modulate protein stability.

148

NOTE: Other inducible systems that tune protein expression can potentially be chosen.

150

1.3.5. Output: Place a reporter gene downstream UAS promoter that respond to GAL4VP16 transcription factor. Choose the output according to the application of the device (i.e., fluorescent protein, cytokine, etc.).

154

1.4. Perform cloning with any preferred strategy. Here, we built plasmids with golden gate or gateway technologies.

157158

159

160

1.5. Test the various configuration of the device in the desired cell lines by transient transfection. Co-transfect the cells with the obtained plasmids for fusion proteins 1 and 2 and fluorescent reporter gene, such as UAS-EYFP. Fusion protein 1 include a red fluorescent protein (mkate) used as transfection marker.

161162

NOTE: Start with simple-to-transfect cell lines such as HEK293FT cells to test the devices.

163164

1.5.1. Transfection in HEK293FT cells. Seed 2 x 10⁵ cells into a 24-well plate with 500 μL of
 DMEM. Transfect cells with 300 ng of total DNA using transfection reagent, vortex the mix
 and incubate at RT for 20 minutes.

168

169 1.5.1.1. Place cells in a 37 °C/5% CO₂ incubator. Add fresh culture medium 24 h after transfection.

171

1.5.2. Transfection in Jurkat cells. Use 3×10^5 cells in 500 μ L of RPMI cell culture medium for 24-well plate. Perform electroporation following manufacturer's instruction. Use 2-4 μ g of total DNA for each sample. Add 500 μ L of fresh cell culture media 24 h after transfection.

175

1.5.2.1. During electroporation, place cells and DNA mix on ice. Then, place the cells in the 37 $^{\circ}$ C/5% CO₂ incubator.

178

NOTE: Different cell lines and transfection methods may be chosen according to the specific application of the device.

181 182

183

184

1.6. Acquire microscope images to obtain a qualitative understanding of the functionality of the device. Confirm membrane localization of the fusion protein 1 by detecting the fluorescent signal (mkate) in proximity to the cell membrane. Visualize the fluorescent output (EYFP) in the presence or absence (negative control-NC) of the target protein.

185 186

187 1.6.1. Use the cell imaging system equipped with Tx Red and GFP light cubes to detect mKate and EYFP fluorescent proteins. Use the 10x or 20x objective.

NOTE:	Use non-transfected cells to ensure that there is no autofluorescence.
1.7.	Analyze the transfection with flow cytometry (FACS).

194 NOTE: See Protocol 2 for sample preparation. 195

196 197

198

200

203

208

211

215

216 217

220

223 224

225

230

234

1.7.1. Evaluate the devices in the presence of the target protein. Create an experimental matrix that combine variants of fusion protein 1 and 2, modulating TEVp with (i) constitutive hEf1 α promoter (ii) doxycycline-inducible TET promoter and (iii) shield for protein stability. Add doxycycline and shield at concentrations (0-1000 nM).

199

201 1.7.2. Perform FACS analysis to compare the fold induction of EYFP between samples with 202 or without protein of interest.

204 NOTE: The highest fold induction corresponds to the greatest input sensitivity in ON (presence 205 of the target protein) as compared to the OFF (absence of the target protein) mode.

- 206 207 1.8. Test the functionality of the device with respect to the desired applications.
- 209 1.8.1. Run FACS analysis to measure the levels of output proteins or apoptosis, or other relevant assays. 210
- 1.8.2. Perform RT PCR analysis if the output part of the device regulates expression levels of 212 target genes. Calculate 2-ddCT in respect to the house-keeping gene (GAPDH or equivalent) and 213 214 compare the fold change of output expression in On and Off modes.

2. Sample preparation for flow cytometry analysis

- 218 2.1. Perform quantitative analysis of protein devices with flow cytometer (FACS) 48 h after transfection. 219
- 2.1.1. Aspirate the DMEM media and wash HEK293FT cells with 300 μL of PBS. Add 100 μL 221 of trypsin and leave in the 37 °C/5% CO₂ incubator for 2-4 minutes. 222
 - NOTE: Jurkat cells do not require trypsin.
- 226 2.1.2. Add 400 µL of DMEM media without phenol red and resuspend the cells to avoid 227 clumps formation. Transfer in FACS tubes. 228
- 2.1.3. Use a flow cytometer equipped with 405, 488, and 561 nm lasers. 229
- 2.1.4. Collect 30,000 100,000 events selected from live cell population using the following 231 cytometer settings: 488 nm laser and 530/30 nm bandpass filter for EYFP/EGFP, 561 nm laser 232 233 and 610/20 nm filter for mKate, and 405 nm laser, 450/50 filter for EBFP.

NOTE: To set up the flow cytometer, include non-transfected cells and cells transfected with single fluorescent proteins.

237

238 2.1.5. Perform analysis of collected samples using preferred software (e.g., DIVA software, Flowjo, or TASBE³² method).

240241

3. Characterization of the HCV sensor-actuator device

242243

3.1. Select intrabodies to detect NS3 protein. scFv35 and scFv162 bind to two different epitopes and were thus chosen.

244245246

247

248

3.2. Design and test variants of the NS3-responsive device following the rational described in step 1.3 (e.g., Fusion protein 1: membrane tag, mKate, single chain fragment intrabody scFv35, LD0, TCS-L and TF GAL4VP16; Fusion protein 2: single chain fragment intrabody scFv162, LD0 and TEVp; Output: UAS-EYFP or UAS-hBax (apoptosis); NS3 (see **Table 2**)).

249250251

3.3. Test the constructs in HEK293FT cells.

252253

NOTE: For each experiment, set a negative control that is the delivery of the device and output quantification in the absence of the protein of interest.

254255

3.3.1. Test the devices for TEVp under a constitutive promoter in the presence or absence of NS3. Perform transfection of fusion protein 1 (step 1.3.3), fusion protein 2 (step 1.3.4) and EYFP output step (1.3.5) in HEK293FT cells as described in step 1.5.1. Analyze the samples 48 h post-transfection by measuring the fluorescent reporter (EYFP) by flow cytometry (Protocol 2).

261

3.3.2. Assess the interaction intrabody-NS3 by mKate and BFP colocalization using the BFP-NS3 fusion protein. Transfect cells with fusion protein 1 that includes mKate, and with a BFP-NS3. 48 h post-transfection, observe cells at confocal microscope (63x objective) to test for co-localization of fluorescent proteins.

266267

NOTE: Use a BFP (not fused to NS3) as negative control. In this condition BFP should show diffuse cellular localization and should not co-localize with the membrane-bound mKate.

268269

- 3.3.3. Test the devices for TEVp under the pTET promoter in the presence or absence of NS3.
 Compare EYFP in absence and presence of NS3. Co-transfect cells with fusion protein 1 (step
- 272 1.3.3), fusion protein 2 (step 1.3.4) and EYFP output (step 1.3.5) using the protocol described
- 273 in step 1.5.1. Transfection are performed in the presence or absence of NS3. Analyze the
- samples 48 h post-transfection by measuring the fluorescent reporter (EYFP) by flow
- 275 cytometry (Protocol 2).

276

277 3.4. Test cell killing driven by NS3 device.

278

3.4.1. Replace EYFP output with the apoptotic gene hBax.

- 3.4.2. Repeat step 3.3.3. Collect the cells 48 h after transfection and wash with PBS before
- 282 starting the assay.

283

3.4.3. Perform apoptosis assay with Annexin V conjugated to Pacific Blue. Stain cells with apoptotic marker Annexin V (conjugated with 2.5 μL of Pacific Blue) diluted in binding buffer (1:20 volume ratio) for 10 min at room temperature.

287

3.4.4. Analyze the cells by flow cytometry. Calculate the apoptosis level within population as number of Pacific-Blue (Annexin V) positive cells.

290

4. Characterization of Nef-responsive HIV sensor-actuator device

291292

293 4.1. Select intrabodies that bind to Nef. Single domain intrabody sdAb19 and SH3 were chosen as they bind to different epitopes.

295

NOTE: SH3 is not an intrabody but a domain of p59^{fyn} protein tyrosine kinase that physiologically interacts with Nef.

298

4.2. Design and test variants of Nef-responsive device following the rational described in step 1.3 (e.g., Fusion protein 1: membrane tag, mKate, sdAb19, LD0, TCS-L and GAL4-VP16; Fusion protein 2: SH3, LD0 and TEVp; Output: UAS-EYFP; Nef (see **Table 3**)).

302

303 4.3. Test the devices in HEK293FT and Jurkat T cells.

304

NOTE: For each experiment we set a negative control that is delivery of the device and output quantification in the absence of the protein of interest

307

4.3.1. Transfection of HEK293FT: Follow the procedure described in step 1.5.1 in presence or absence of Nef.

310

4.3.2. Transfection of Jurkat T cells: Follow the procedure described in step 1.5.2 in presence or absence of Nef.

313

4.3.3. 48h post-transfection, analyze the samples with flow cytometer as described in Protocol 2.

316

317 4.4. Infect Jurkat cells with HIV LAI strains for HIV detection.

318

- 4.4.1. Prepare viral stocks at least 48 h before the experiment. Use virus obtained from transfection of HEK-293T cells with the infectious molecular clones (NIH AIDS reagents program) and commercial reagents using manufacturer instructions. Concentrate virus after
- 40 h by ultracentrifugation for 1 h, 64,074 x g, 4°C on 20% sucrose to avoid viral particle-free
- 323 proteins. Titer viral stocks with HIV-1 p24 ELISA.

- 4.4.2. Infect cells with HIV-1 strains with 500 ng of p24 inoculum. 40 h post-infection quantify
- 326 the percentage of infected cells by detecting the viral protein p24 with KC57-FITC conjugated
- antibody and performing flow-cytometry analysis.

328329

4.4.3. Stain the cells with AlexaFluor 647 conjugated antibodies clone W6/32 against HLA-A, B, C. Fix the cells and analyze on a flow cytometer. Nef mediated downmodulation of HLA-I in Jurkat T cells in the presence or absence of the protein sensor device.

331332333

334

335

336337

338339

340

341

342

330

REPRESENTATIVE RESULTS:

An architecture for modular intracellular protein detection

As shown in **Figure 1**, the device is composed of: 1) intrabody 1 connected to the membrane-tethered fluorescent marker (mKate) and TEVp cleavage site (TCS), followed by a transcription activator GAL4VP16 (TF); 2) intrabody 2 fused to TEV protease (TEVp), free in the cytosol; 3) a synthetic promoter responsive to GAL4VP16, driving the expression of a reporter gene. The modularity is guaranteed by intrabodies that can be changed in a plug-and-play manner to target any desired protein. Also, the output can be chosen depending on the desired application. A critical aspect in order to optimize sensing-actuation performance is to avoid background activation of the system (in absence of the protein of interest). We thus tuned several features of the device (**Table 1**).

343344345

346

347

348

349

350

351 352

353

354

355356

357

358

359

360

HCV, Huntingtin and HIV sensor-actuator devices

The first device was built to recognize NS3 protein, expressed from the epatitis C virus HCV²⁸. The study of this sensor set the rules for the following ones. In specific, we co-transfected fusion protein 1 variants of the NS3-responsive devices along with constitutively expressed (hEF1 α) TEVp-scFv162 in HEK293FT cells. We observed NS3-independent output activation both with the TCS(S) and TCS(L) (Figure 2a), indicating that the system requires fine regulation of TEVp activity. To achieve input specific signal activation, we modulated TEVp expression transcriptionally and post-translationally. The first was mediated by a doxycycline (dox) responsive promoter (pTET), whereas the second was obtained by adding to TEVp a degradation domain tag (DD degron) regulated by the small molecule shield. The device was able to specifically induce gene expression following NS3 detection, in the absence of dox and shield when the system included the TCS(L) (Figure 2b). These results indicate that minimal TEVp concentration is required for efficient devices. On the contrary maximal activation of TEVp by dox and shield in the TCS(L) configuration, or the TCS(S) exhibits high background (Figure 2b, c). Altogether, the results demonstrate that the simple leakiness of the promoter (absence of doxycycline or the activator rtTA3) enables the best ON/OFF condition, highlighting the high sensitivity of the TEVp actuator.

361362363

364

365366

Based on the results we designed further NS3 device variants by changing parameters previously mentioned (**Table 1**) and demonstrated up to 24 EYFP fold induction (**Figure 2d,f**), and specific cell killing (data in²⁰). Further, we fused the NS3 to a blue fluorescent protein and demonstrated colocalization with fusion protein 1, indicating the interaction between the component of the system (**Figure 2e**).

367368369

370

371372

373

374

Nef sensor-actuator device

Differently from the previous example, the sensing module of Nef device is composed of a single domain antibody (sdAb19)³⁰, and SH3 domain from p59fyn protein tyrosine kinase that is highly expressed in T lymphocytes and implicated in antigen induced T-cell activation³³. Since Nef is involved in HLA-I downmodulation, we hypothesized that by sequestering Nef, the device could impair its pathological activity (depicted in **Figure 3**).

We first tested the Nef-responsive device in 293FT HEK and Jurkat T cells. The sensor exhibited up to 5-fold induction of EYFP in the presence of the target protein (**Figure 4a**). Different from the NS3 case, here the activation of TEVp by dox (D100) was required to obtain a functional device, whereas LD0 demonstrated lower background in the absence of the target protein. Moreover, we observed EYFP expression in Jurkat T cells infected with X4 LAI strain (**Figure 4b**). We also demonstrate that the device impairs Nef-mediated downregulation of HLA-I (**Figure 4c**). Finally, we rewired the output to induce XCL-1, a chemokine normally secreted by CD8+ T cells to induce a localized CTL response, demonstrating the potential towards selective activation of the immunomodulator (**Figure 4d**).

FIGURE AND TABLE LEGENDS:

Figure 1: Architecture of the sensor-actuator device for intracellular protein detection. The device is composed of intrabodies (or small peptides) to recognize with high specificity the protein of interest (POI). Interaction of intrabodies with POI brings TEVp in the proximity of the cleavage site (TCS). The membrane-tethered transcription factor (TF) fused to TCS translocates to the nucleus upon TEVp cleavage, activating programmed output. The figure has been modified from²⁰.

Figure 2: Test of the NS3 device. (a) TCS(S) and (L) variants of NS3 device were tested in HEK293FT cells in presence or absence on NS3. TEVp is driven by a constitutive promoter. Data show that the system is not sensitive to the presence of NS3 due to an overactivation of TEVp that leads to non-specific cleavage of the TCS. (**b-c**) NS3 induced EYFP reporter activation with TCS(L) but not TCS(S) in absence of Dox and Shield. (**d**) Test of the variants of NS3 sensor in HEK293FT cells show up to 24x induction of reporter gene when the HVC protein was detected. Data shows fold induction and standard deviation. (**e**) Co-localization in HEK293FT cells of NS3 fused to a BFP (BFP-nNS3) and fusion protein 1 which include the membrane tethered mKate. As control we used a BFP alone, that shows diffuse cellular localization. Confocal images (scale bar = 25 μm) (**f**) fluorescent reporter induction as result of NS3 detection in HEK293FT cells. BF: bright field; mKate: fluorescent marker included in fusion protein 1; EYFP: output of the device. Evos images (scale bar = 200 μm). The figure has been modified from 20 .

Figure 3: Model of Nef-responsive device. (a) In the absence of the engineered device, HIV escapes the immune response via Nef-mediated HLA downregulation. (b) Nef is detected and sequestered by the sdAb19 intrabody and SH3 domain to activate Xcl1 expression. Xcl1 stimulates anti-HIV immune response attracting dendritic cells and cytotoxic lymphocytes. In addition, Nef sequestration impairs HLA downmodulation such that CD4+T cells infected by HIV are more responsive to immune system attack. The figure has been modified from²⁰.

Figure 4: Nef sensor-actuator device detects HIV infection and impair HLA downmodulation.(a) HEK293FT cells induce output activation in the presence of Nef, at intermediate concentration of Dox. (b) Jurkat T cells infected with HIV-LAI strain show reporter activation.
(c) Nef device interferes with HLA downmodulation as compared to HLA levels in absence of sensor. (d) XCL-1 immunomodulator is triggered by the device upon Nef detection. The figure

has been modified from²⁰.

Table 1: Components for modular device design. Construction of efficient modular devices requires several elements, which can be adjusted to improve the system sensitivity, and specificity. TEVp activity can be adjusted in order to achieve the best signal-to-noise ratio in the system: (1) via transcriptional modulation of TEVp expression using a constitutive (hEF1α) or Doxycycline inducible promoter (pTET), or via posttranslational modulation, using a degradation domain tag (DD) controlled by the small molecule Shield; (2) via high or low affinity TEV cleavage sites (TCS) that include a serine (TCS-S) or a leucine (TCS-L) in P1 position; (3) via repositioning of TEVp to the N or C terminus of the intrabody. To take into account the intrabodies-protein interaction and TEVp accessibility to TCS we (4) included Glycine-Serine (G4S) Linker Domains (LD) of variable amino acid length between intrabody and TCS and between TEVp and intrabody. (5) Selection of intrabodies or their analogues to the target protein. Of note, *SH3 of the Nef-responsive device is not an intrabody but a domain of p59fyn protein tyrosine kinase that bind a proline-rich motif in Nef ³¹. These components were assembled in a plug-and-play fashion to build a modular system towards optimal input-output devices.

Table 2: Modules of the best NS3 device. Best performing sensory-actuator device is reported for NS3. Precisely, we specify the length of LD domain (LD0, LD15), affinity of the protease binding site TCS (S, L), type of TEVp regulation (hEF1 α or pTET, or DD degron).

Table 3: Modules of the best Nef device. Best performing sensory-actuator device is reported for Nef. Precisely, we specify the length of LD domain (LD0, LD15), affinity of the protease binding site TCS (S, L), type of TEVp regulation (hEF1 α or pTET).

DISCUSSION:

Until recently, interrogating cells based on intracellular environment was performed with systems developed de novo for specific targets. The present protocol describes an example of the most recent, cell engineering approach for protein sensing and actuating in one device, that can be rapidly adapted to new desired biomarkers.

This pioneering system sense intracellular proteins and provide a specific output to detect or neutralize the disease. The advantage of this class of genetic circuits is the modularity provided by intrabodies and the possibility to use other sensing "parts" to target the protein as shown for the Nef sensor. To our knowledge, a system that can be potentially applied to several intracellular proteins was never demonstrated before.

The modular architecture allows the researchers to adjust the components to maximize ON/OFF ratio of the device for any desired application. This suggests that the device must be carefully tested to achieve the best functional performance. In general, we observed that absence of linker domain (LD0) resulted in more effective device functionality, which perhaps may be the first design strategy for new protein-sensing devices. This is likely due to optimal physical interactions between target protein, intrabodies and TEVp in this device configuration³⁴.

With respect to the actuator module, we observed that TEVp levels, affinity to its cleavage sites, and location in the fusion protein are critical. Constitutive promoters driving protease expression resulted in high background in all devices. In addition, in this experimental set-up

TEVp cleavage sites with lower affinity (TCS-L) worked better for NS3, Tat and Nef sensing devices, whereas for HTT detection (reported in²⁰) the TCS-S showed better performances.

Last, we observed higher ON/OFF signal ratio with TEVp fused at N-terminal of the intrabody;

472 however, for Tat sensing device TEVp activity was equal at the N- or C-terminal of the chimeric 473 protein.

474 475

476 477

478

The described framework allows induction of desired functional response when the protein of interest is detected within the cell. One critical aspect is the sensitivity of the system to the target. We have shown for example that the output expression is a function of NS3 concentration (data shown in²⁰) and that Jurkat T cells were sensitive to HIV infection at titers commonly used in the immunology labs.

479 480 481

482

483

484

485

486

Applications of these circuits cover a broad range spanning from dynamic monitoring of cellular state to programmed therapeutic expression upon detection of a disease condition. In perspective, the rapid development of intrabodies with yeast and phage display libraries, that are smaller in size and stable in the cell (i.e. nanobodies) could lead to the design of devices specific for any wanted protein. In addition, protein design tools, and computational models that can predict circuitry performance help in designing effective devices minimizing the build-and-test reiterations.

487 488 489

Our framework provides a novel strategy to connect aberrant protein expression, post-translational modification or viral infection into an information that is translated in transcriptional activation of gene expression.

491 492 493

490

ACKNOWLEDGMENTS:

This work was supported by the Istituto Italiano di Tecnologia.

494 495 496

DISCLOSURES:

The authors declare no competing financial interests.

497 498 499

515

REFERENCES:

- 500 1. Bernardo, D., Marucci, L., Menolascina, F. & Siciliano, V. Predicting Synthetic Gene 501 Networks. in *Synthetic Gene Networks: Methods and Protocols* **813**, 57–81 (2012).
- 502 2. Krams, R. *et al.* Mammalian synthetic biology: emerging medical applications. *J. R. Soc. Interface* **12**, (2015).
- MacDonald, J. T. & Siciliano, V. Computational Sequence Design with R2oDNA
 Designer. in *Methods in molecular biology (Clifton, N.J.)* 1651, 249–262 (2017).
- Zah, E., Lin, M.-Y., Silva-Benedict, A., Jensen, M. C. & Chen, Y. Y. T cells expressing
 CD19/CD20 bi-specific chimeric antigen receptors prevent antigen escape by
 malignant B cells. *Cancer Immunology Research* 4, 498–508 (2016).
- 509 5. Wu, F., Bethke, J. H., Wang, M. & You, L. Quantitative and synthetic biology
 510 approaches to combat bacterial pathogens. *Current Opinion in Biomedical Engineering* 4, 116–126 (2017).
- 512 6. Ye, H. *et al.* Pharmaceutically controlled designer circuit for the treatment of the metabolic syndrome. *Proceedings of the National Academy of Sciences of the United States of America* **110**, 141–6 (2013).
 - 7. Thaker, M. N. & Wright, G. D. Opportunities for synthetic biology in antibiotics:


- Expanding glycopeptide chemical diversity. ACS Synthetic Biology 4, 195–206 (2015).
- 517 8. Culler, S. J., Hoff, K. G. & Smolke, C. D. Reprogramming cellular behavior with RNA
- 518 controllers responsive to endogenous proteins. *Science (New York, N.Y.)* **330**, 1251–5 (2010).
- 520 9. Ausländer, S. *et al.* A general design strategy for protein-responsive riboswitches in mammalian cells. *Nature Methods* **11**, 1154–1160 (2014).
- 522 10. Cella, F. & Siciliano, V. Protein-based parts and devices that respond to intracellular 523 and extracellular signals in mammalian cells. *Current Opinion in Chemical Biology* **52**, 524 47–53 (2019).
- 525 11. Miki, K. *et al.* Efficient Detection and Purification of Cell Populations Using Synthetic 526 MicroRNA Switches. *Cell Stem Cell* **16**, 699–711 (2015).
- Xie, Z., Wroblewska, L., Prochazka, L., Weiss, R. & Benenson, Y. Multi-input RNAi-based logic circuit for identification of specific cancer cells. *Science (New York, N.Y.)* 333, 1307–11 (2011).
- 530 13. Sun, J. *et al.* Engineered proteins with sensing and activating modules for automated reprogramming of cellular functions. *Nature Communications* **8**, 477 (2017).
- 532 14. Shigeto, H. *et al.* Insulin sensor cells for the analysis of insulin secretion responses in single living pancreatic β cells. *Analyst* **144**, 3765–3772 (2019).
- 534 15. Cella, F., Wroblewska, L., Weiss, R. & Siciliano, V. Engineering protein-protein devices 535 for multilayered regulation of mRNA translation using orthogonal proteases in 536 mammalian cells. *Nature Communications* **9**, (2018).
- 537 16. Siciliano, V. *et al.* MiRNAs confer phenotypic robustness to gene networks by suppressing biological noise. *Nature communications* **4**, 2364 (2013).
- Barnea, G. et al. The genetic design of signaling cascades to record receptor
 activation. Proceedings of the National Academy of Sciences of the United States of
 America 105, 64–9 (2008).
- 542 18. Schwarz, K. A., Daringer, N. M., Dolberg, T. B. & Leonard, J. N. Rewiring human 543 cellular input—output using modular extracellular sensors. *Nature Chemical Biology* 544 **13**, 202–209 (2016).
- 545 19. Wieland, M. & Fussenegger, M. Engineering molecular circuits using synthetic biology 546 in mammalian cells. *Annual review of chemical and biomolecular engineering* **3**, 209– 547 34 (2012).
- 548 20. Siciliano, V. *et al.* Engineering modular intracellular protein sensor-actuator devices.


 549 *Nature Communications* **9**, 1881 (2018).
- 550 21. Lin, C. HCV NS3-4A Serine Protease. (2006).
- Romani, B., Engelbrecht, S. & Glashoff, R. H. Functions of Tat: the versatile protein of human immunodeficiency virus type 1. *Journal of General Virology* **91**, 1–12 (2010).
- 553 23. Basmaciogullari, S. & Pizzato, M. The activity of Nef on HIV-1 infectivity. *Frontiers in microbiology* **5**, 232 (2014).
- Ross, C. A. & Tabrizi, S. J. Huntington's disease: from molecular pathogenesis to clinical treatment. *Lancet neurology* **10**, 83–98 (2011).
- 557 25. Pfeffer, C. M. & Singh, A. T. K. Apoptosis: A Target for Anticancer Therapy.


 International journal of molecular sciences **19**, 448 (2018).
- 559 26. Guzzo, C. *et al.* The CD8-derived chemokine XCL1/lymphotactin is a conformation-560 dependent, broad-spectrum inhibitor of HIV-1. *PLoS pathogens* **9**, e1003852 (2013).
- 561 27. Marasco, W. A., LaVecchio, J. & Winkler, A. Human anti-HIV-1 tat sFv intrabodies for gene therapy of advanced HIV-1-infection and AIDS. *Journal of Immunological*


563 *Methods* **231**, 223–238 (1999).

- 564 28. Gal-Tanamy, M. *et al.* HCV NS3 serine protease-neutralizing single-chain antibodies 565 isolated by a novel genetic screen. *Journal of molecular biology* **347**, 991–1003 566 (2005).
- 567 29. Southwell, A. L. *et al.* Intrabodies binding the proline-rich domains of mutant 568 huntingtin increase its turnover and reduce neurotoxicity. *The Journal of* 569 *neuroscience : the official journal of the Society for Neuroscience* **28**, 9013–9020 570 (2008).
- 571 30. Bouchet, J. *et al.* Inhibition of the Nef regulatory protein of HIV-1 by a single-domain antibody. *Blood* **117**, 3559–3568 (2011).
- 573 31. Arold, S. *et al.* The crystal structure of HIV-1 Nef protein bound to the Fyn kinase SH3
 574 domain suggests a role for this complex in altered T cell receptor signaling. *Structure*575 **5**, 1361–1372 (1997).
- Beal, J. Signal-to-Noise Ratio Measures Efficacy of Biological Computing Devices and Circuits. *Frontiers in Bioengineering and Biotechnology* **3**, 93 (2015).
- Arold, S. et al. The crystal structure of HIV-1 Nef protein bound to the Fyn kinase SH3
 domain suggests a role for this complex in altered T cell receptor signaling. Structure
 (London, England: 1993) 5, 1361–1372 (1997).
- Tyshchuk, O. *et al.* Detection of a phosphorylated glycine-serine linker in an IgG-based fusion protein. *mAbs* **9**, 94–103 (2017).









1. Regulation of TEVp expression						
Transcriptional		Constitutive (hEF1 α)				
		Inducible (pTET)				
Post-translational		Degradation Domain tag (DD degron) + Shield				
2. TEV Cleavage Site (TCS)						
Aa substitutions in P1			Affinity			
TCS-S	Serine		High			
TCS-L	Leucine			Low		
	3. Position of TEVp					
N-terminus	N-terminus of intrabodies					
C-terminus of intrabodies						
	4. Glycine-Serine (G₄S) Linker Domains (LD)					
LD0	0 Amino acids					
LD15	15 Amino acids					
5. Intrabodies and analogues						
		NS3	HTT	Tat	Nef	
Intrabody 1		scFv35	Happ1	scFv2	sdAb19	
Intrabody 2		scFv162	V _L 12.3	scFv3	SH3*	

Parameters of the protein sensor-actuator devices

TEVp affinity to the cleavage site

TCS (S) ENLYFQ S
TCS (L) ENLYFQ L

TEVp activity modulation by small molecules

Transcriptional Doxycycline
Translational Shield

TEVp access to the TCS

Fusion protein design N- or C- terminus fusion o TEVp to the intrabody

Protein flexibility

Flexible Linker (LD) no linker

15Aa (G4S)3

	Principal engin	eered modules
	Components	Short name
Nef	Protein of interest	Nef
	Fusion protein 1	phEF1a-sdAb19-LD0-TCS(L)
	Fusion protein 2	pTET-TEVp-LD0-SH3
	Output (cytokine)	UAS-XCL1
	Output (fluorescence)	UAS-EYFP

Name Company

AlexaFluor 647 mouse anti-human HLA-A, B, C antibody clone W6/32

Annexin V Attractene

BD Falcon Round-Bottom Tubes

Doxycycline

Dulbecco's modified Eagle medium Evos Cell Imaging System

FACSDiva8 software

Fast SYBR Green Master Mix

FBS (Fetal Bovin Serum)
Gateway System

Golden Gate System

HEK 293FT

Infusion Cloning System

JetPRIME reagent

Jurkat Cells

L-Glutamine

Lipofectamine LTX with Plus Reagent

LSR Fortessa flow cytometer (405, 488, and 561 nm lasers)

MicroAmp Fast Optical 96-Well Reaction Plate (0.1 mL)

Neon Transfection System

Non-essential amino acids

Opti-MEM I reduced serum medium

Penicillin/Streptomycin

QuantiTect Reverse Transcription Kit

RNeasy Mini Kit RPMI-1640

Shield

SpheroTech RCP-30-5-A beads

StepOnePlus 7500 Fast machine

Biolegend

LifeTechnologies

Qiagen

BD Biosciences

Clonetech

Cellgro

Life Technology

BD Biosciences

ThermoFisher Scientific

Atlanta BIO

Life Technologies

in-house

Invitrogen

Clonetech

Polyplus transfection

ATCC

Sigma-Aldrich

Thermo Fisher Scientific

BD Biosciences

ThermoFisher Scientific

Life Technologies

HyClone

Life Technologies

Sigma-Aldrich

Qiagen

Qiagen

ATCC

Clonetech

Spherotech

Applied Biosystems

Catalog Number	Comments	
311414	Antibodies	
A35122	Apoptosis marker	
301005	Transfection reagent	
352053	FACS tubes	
-	Cell Culture: Drugs	
10-013-CM	Cell Culture: Medium	
EVOS M5000	${\mathfrak z}$ systems; Infectious molecular clones	
659523	FACS software	
4385612	qPCR reaction	
S11050	Cell Culture: Medium	
-	Plasmid Construction	
-	Plasmid Construction	
R70007	Cell Culture: Cells	
638920	Plasmid Construction	
114-15	Transcfection reagent	
TIB-152	Cell Culture: Cells	
G7513-100ML	Cell Culture: Medium	
15338030	Transfection reagent	
649225	Flow cytometer	
4346907	qPCR reaction	
MPK10025	Transfection reagent	
SH3023801	Cell Culture: Medium	
31985070	Transfection medium	
P4458-100ML	Cell Culture: Medium	
205313	Rev Transcriptase kit	
74106	RNA extraction kit	
ATCC 30-2001	Cell Culture: Medium	
632189	Cell Culture: Drugs	
RCP-30- 5A-2	Compensation set up	
4351106	qPCR reaction	

Reviewers' comments:

Reviewer #1:

Manuscript Summary:

The manuscript describes a method for sensing intracellular proteins and linking this to a transcriptional output.

Major Concerns:

-A mathematical model is mentioned, but not really used in optimization or presented. It might be better to just omit this statement since it currently reads as if it is an afterthought. Otherwise, please clarify from where it is available and how to use it to guide the development of new devices.

We thank the reviewer for the feedback. Indeed, a mathematical model was used and described in the original manuscript but including it in this protocol would probably be out of the scope of this paper. We then decided to remove the reference to it to avoid confusion.

Minor Concerns:

-It may be necessary to explain what TEVp is in the abstract so that those who read it to decide if the protocol is of use can understand it.

We have now included the description of TEVp in the abstract (lines 40-42)

-Not clear what 'an interactor' is, perhaps a concrete example would help, though I am not sure if this is a real word or not.

We have replaced the word "interactor" with a deeper description of "molecules that interact with" throughout the manuscript

-Lines 75-77. Not clear if no solutions were proposed or if there were solutions proposed but these required bespoke engineering

The main novelty of this system is that it is a unique framework that can be easily adapted for theoretically any cytosolic protein. This is different from before where devices were built for a specific protein. We have provided a better description of this uniqueness in the revised version of the manuscript.

-The protocol describes transfection of HEK and Jurkat cells. This may lead readers to assume they must use one of these types of cells, but this is not the case. Perhaps clarify by including a Note. Also, do the authors feel that these exact protocols for transfection are necessary or could one transfect by another means or using a different optimized protocol?

We thank the reviewer for pointing this out. We have now included notes to explain that both cell lines and transfection methods can be chosen according to applications and preferences of the researchers.

-Lines 354-6, given that LD0 gives higher functionality, would you still recommend trying linkers? Or can this step be omitted?

Indeed, we observed better performances with LDO (along with other conditions mentioned in the discussion). We would suggest to try it (but not limiting it to it) as first option, and included it in

the revised discussion.

Reviewer #2:

Manuscript Summary:

This manuscript is interesting for the researches that are working on monitoring protein expression intracellularly. However, it is a very complex co-expression system. Every system needs at least two fusion proteins for co-transfection and each fusion protein 1 needs two to four components. The efficiency and sensitivity can be effect by many factors in this system, such as the efficiency of the TEVp, the introbody epitope overlap, the efficiency of the TF and the downstream output protein transcription and translation.

The overall of this manuscript is not well written, the quality of the figures and tables should be improved.

We thank the reviewer for the feedback. We have worked to substantially increase the quality of language as well as of the figures. The latter were generated with Adobe Illustrator and converted in PDF with the maximal quality in the set-up.

Major Concerns:

1. How sensitive is this system? The authors should suggest the expression level of the POI, which fits into this system.

We thank the reviewer for the comment. Indeed, sensitivity is a critical parameter to indicate the potential use of our devices. In our original manuscript we have tested as proof of concept the response of the device to increasing NS3 concentration. Although for format constraints reasons we have not added the data in this manuscript, we have mentioned their importance in the discussion (line 361-364) providing the reference to the original paper. In addition, we have tested the response to HIV infection which were concentrated to titers typically used in the immunology labs (performed by our collaborator at Ragon Institute).

2. The authors suggested that to optimize the sensing performance was to avoid background activation of this system. Therefore, the author should include 'in the absence of the protein of interest-POI' controls such as POI knock out cell line or blocking antibody treatment to set the background.

Our experiments and output quantification were performed using a Negative Control condition that is output quantification in the absence of the POI. We clarified this method including a note for each device we tested (Lines: 209-210,252-253)

- 3. How to control the expression efficiency of these two fusion proteins?

 One fusion protein (fusion protein 1) included a fluorescent protein (mKate), which allowed membrane visualization and to quantify transfection efficiency). The second fusion protein (fusion protein 2) did not include fluorescent marker, but for each test we ran a control in which we maximize TEVp expression to induce reporter expression. This control is a proxy of transfection efficiency of our systems.
- 4. Is the membrane tag is needed for the fusion protein1? The author should show member localized imaging data in the representative result.

We performed membrane localization imaging and showed it in the original manuscript. We now include these images in Figure 3d.

5. Why fusion protein 2 does not include a fluorescent marker to test the fusion protein2 expression and cell transfection.

We thank the reviewer for pointing this out. As we explain in the previous comment, we followed TEVp activity by output activation. Therefore, the reporter gene functions as a proxy of TEVp expression.

- 6. The author should show representative results for FACS and microscopy Flow cytometry data are presented in Figure 2 and 4 and plotted as ratio of reporter fold induction, calculated as the normalized fluorescence in the presence of the POI divided by the normalized florescence in the absence of the POI. For the normalization we used the mKate as transfection marker. We now included in Fig.2 microscope images of the reporter activation in the presence of NS3.
- 7. The author did not explain to the reader how to select between LD0 or LD15, TCS(S) or TCS(L), and who to choose the position of TEVp.

We thank the reviewer for the useful comment. The selection of all parameters was made to increase the likelihood to achieve a protein-specific induction of readout. Because of the complex architecture, we sought to build variants that allow the chimeric proteins to be more or less flexible, and to render the TCS accessible to TEVp upon interaction of the intrabodies with the target protein. All the features are included in Table 1 and mentioned in a note of the revised manuscript (lines 119-122). In the discussion we comment on the best results achieved as setting rules for future protein sensing-actuating devices.

- 8. The author did not maintain the affinity of the intrabodies. We are not sure to understand the comment. Since the intrabodies were previously published, the affinity to the proteins were already measured.
- 9. Only Figure 4 c has a statistic, but not other figures? We cannot get the conclusion for the effective range for the readout.

We could calculate the statistics for all the experiments performed, for which instead we plot the fold induction and relative standard deviation. However, we believe that since the scope of this manuscript is to describe clear protocols, we focused our effort to provide more useful information to the readers.

Minor Concerns:

1. Table 2 should be better organized and proved more information to make it correlated to the different domains of figure 1

Figure 1 provides the general structure of the framework. In Table 2 we opted to specify the single parts of the architecture that resulted in better sensing devices. Also, to avoid confusion we removed the HTT and Tat devices that are not described in this manuscript as well as the extended name of the plasmid that do not give additional information about the characteristics of the system. We hope that these modifications make table 2 more suitable to the readers.

2. Line 155. "Optional:.." the On and OFF models should be better defined. Same for line 176.

We specified what On and Off modes refer to.

3. The FACS method and the microscopy method to observe the output should write separate and clearer.

We have now created a protocol section of the FACS method (Protocol 2). We also explained better the microscopy method used for observation (Protocol 1.6) although since it was purely used for observations, we did not create a new section for it.

Reviewer #3:

Manuscript Summary:

The authors have described a protocol for engineering a modular genetic platform that sense intracellular proteins and activate a specific cellular response. The device operates on intracellular antibodies or small peptides to sense with high specificity the protein of interest. They have described the detection of proteins such as mutated huntingtin, NS3 serine-protease, Tat and Nef proteins to detect Huntington's disease, hepatitis C virus (HCV) and human immunodeficiency virus (HIV) infections respectively.

Major Concerns:

I did not notice any major issues in this submission.

Minor Concerns:

There are no images of the fluorescence detection of these proteins. If the authors could have taken microscopic images of the cells with fluorescence responses and the control cells without fluorescence responses, it would have helped the readers of the article immensely. Then, the readers would know what they should expect to see, when they try to repeat this reported protocol in their own labs. I am sure the apoptotic cells would have been very interesting to see in this study. If these images are provided in the video version of this article, then, this minor comment can be ignored.

In view of my comments above, I would recommend a minor revision of the manuscript.

We thank the reviewer for the positive comments. We have now included the microscopy images in the new figure 2e-f of the revised manuscript.