

Journal of Visualized Experiments

Multicellular human alveolar model composed of epithelial cells and primary immune cells for hazard assessment --Manuscript Draft--

Article Type:	Methods Article - JoVE Produced Video
Manuscript Number:	JoVE61090R2
Full Title:	Multicellular human alveolar model composed of epithelial cells and primary immune cells for hazard assessment
Section/Category:	JoVE Biology
Keywords:	in vitro; co-culture model; three dimensional human model; alveolar; primary macrophages; primary dendritic cells; human blood monocyte isolation; thawed monocytes
Corresponding Author:	Barbara Rothen-Rutishauser Adolphe Merkle Institute Fribourg, Fribourg SWITZERLAND
Corresponding Author's Institution:	Adolphe Merkle Institute
Corresponding Author E-Mail:	barbara.rothen@unifr.ch
Order of Authors:	Hana Barosova Barbara Drasler Alke Petri-Fink Barbara Rothen-Rutishauser
Additional Information:	
Question	Response
Please indicate whether this article will be Standard Access or Open Access.	Open Access (US\$4,200)
Please indicate the city, state/province, and country where this article will be filmed . Please do not use abbreviations.	Fribourg, Switzerland

TITLE:

Multicellular Human Alveolar Model Composed of Epithelial Cells and Primary Immune Cells for Hazard Assessment

AUTHORS AND AFFILIATIONS:

Hana Barosova^{1*}, Barbara Drasler^{1*}, Alke Petri-Fink¹, Barbara Rothen-Rutishauser¹

¹Adolphe Merkle Institute, University of Fribourg, Fribourg, Switzerland

*These authors contributed equally.

Corresponding Author:

Barbara Rothen-Rutishauser (barbara.rothen@unifr.ch)

Email Addresses of Co-Authors:

Hana Barosova (hana.barosova@unifr.ch)

Barbara Drasler (barbara.drasler@unifr.ch)

Alke Petri-Fink (alke.fink@unifr.ch)

KEYWORDS:

human lung cocultures, 3D human cocultures, air-liquid interface, isolation of human peripheral blood monocytes, primary macrophages and dendritic cells, cryopreserved primary immune cells

SUMMARY:

Presented here is a protocol for primary human blood monocyte isolation as well as their differentiation into macrophages and dendritic cells and assembly with epithelial cells into a multicellular human lung model. Biological responses of cocultures composed of immune cells differentiated from either freshly isolated or thawed monocytes, upon exposure to proinflammatory stimuli, are compared.

ABSTRACT:

A human alveolar cell coculture model is described here for simulation of the alveolar epithelial tissue barrier composed of alveolar epithelial type II cells and two types of immune cells (i.e., human monocyte-derived macrophages [MDMs] and dendritic cells [MDDCs]). A protocol for assembling the multicellular model is provided. Alveolar epithelial cells (A549 cell line) are grown and differentiated under submerged conditions on permeable inserts in two-chamber wells, then combined with differentiated MDMs and MDDCs. Finally, the cells are exposed to an air-liquid interface for several days. As human primary immune cells need to be isolated from human buffy coats, immune cells differentiated from either fresh or thawed monocytes are compared in order to tailor the method based on experimental needs. The three-dimensional models, composed of alveolar cells with either freshly isolated or thawed monocyte-derived immune cells, show a statistically significant increase in cytokine (interleukins 6 and 8) release upon exposure to proinflammatory stimuli (lipopolysaccharide and tumor necrosis factor α) compared to untreated cells. On the other hand, there is no statistically significant difference between the cytokine release observed in the cocultures. This shows that the presented model is responsive to inflammatory stimuli in

the presence of MDMs and MDDCs differentiated from fresh or thawed peripheral blood monocytes (PBMs). Thus, it is a powerful tool for investigations of acute biological response to different substances, including aerosolized drugs or nanomaterials.

INTRODUCTION:

In vitro lung cell cultures offer cost-effective, robust, and well-controlled platforms to assess the hazards of aerosols¹. As a model cell system for human alveolar pneumocytes, the epithelial A549 cell line isolated from a pulmonary adenocarcinoma is often used². These cells represent squamous type II epithelial cells from the alveolar region³ and are a widely used lung cell line for hazard and toxicity assessment^{1,4-10}. The A549 cell line possesses relevant features of alveolar epithelial type II cells, such as the presence of characteristic lamellar bodies containing densely packed phospholipids³.

It has been shown that when the cells are cultured at an air-liquid interface (ALI), a surfactant is released on the apical side of air-exposed epithelial cells, reducing surface tension¹¹⁻¹³. This feature is particularly important in nanomaterial respiratory hazard and toxicity investigations. Once inhaled nanomaterials/toxicants are deposited in the alveolar region, they first interact with pulmonary surfactant and are displaced by wetting forces into the aqueous hypophase, where the interaction with pulmonary cells takes place^{14,15}. Even though A549 cells form a monolayer (which can overgrow into multilayers at later timepoints when cultivated at ALI) and produce surfactant, a drawback is their insufficient tight junction formation, resulting in low transepithelial electrical resistance values¹⁶⁻¹⁸.

In the lungs, there is a variety of immune cell populations, including phagocytic and professional antigen-presenting cells (i.e., macrophages and dendritic cells) that directly communicate through cell-cell contact or intercellular signaling to control and maintain homeostasis. Macrophages and dendritic cells are critical innate immune effectors and initiators of the adaptive immune response¹⁹. Dendritic cells residing within or beneath the epithelium can form protrusions across the epithelium to the lumen to catch antigens. Alveolar macrophages are located on the apical surface of the epithelium and act as sentinel cells, representing the first cellular defense against foreign material as well as bacterial, viral and fungal infections. Their phenotypic plasticity allows rapid induction of proinflammatory reactions in response to such stimuli as well as shifting into triggering anti-inflammatory (i.e., inhibitory) reactions²⁰.

To simulate the human alveolar epithelial tissue barrier, we established a triple coculture model with A549 cells supplemented with human blood monocyte-derived macrophages (MDMs) and dendritic cells (MDDCs) on the apical and basal sides, respectively¹⁷. Cultivation of this model at ALI has been previously reported¹⁶, even up to 72 h post-exposure²¹. Acute immune responses to carbon nanotubes exposures was significantly enhanced in cell culture exposed to the ALI as compared to submerged conditions²². The coculture model, cultured and exposed to different materials at ALI, has been previously used to investigate cytotoxicity, oxidative stress and inflammatory responses upon exposures to zinc oxide,²³ graphene-related materials²⁴, gold nanoparticles^{25,26}, carbon nanotubes²¹, and volcanic ash and diesel exhaust particles²⁷.

Furthermore, the important role of macrophages and dendritic cells as immune effector cells in an in vitro human lung model has been confirmed. In particular, an increased proinflammatory response in the model was observed only in the presence of immune cells in comparison to monoculture systems⁷. Potential drawbacks of using primary monocyte-derived immune cells are the limited accessibility of PBMs as well as donor-to-donor variation. As a solution to these potential drawbacks, presented here is a protocol introducing cryopreservation of freshly isolated PBMs²⁸ for the cell culture model assembly. The aim of this study is to demonstrate the 3D human alveolar epithelial tissue model assembly, including the isolation of PBMs from human buffy coats. The responsiveness to proinflammatory stimuli is compared to the model composed of MDMs and MDDCs differentiated from fresh PBMs or differentiated from frozen/thawed PBMs.

Working with untested human blood samples involves specific care to prevent the potential transmission of infectious diseases, such as HIV (human immunodeficiency virus), hepatitis B, and hepatitis C. Therefore, the use of personal protective measures such as gloves, gowns, masks, and eye protection are crucial and must be in accordance with good laboratory practice principles. These protections reduce the risk of exposing the skin or mucous membranes to potentially infectious fluids. Additionally, for those involved in handling buffy coats and PBMs, vaccination against the hepatitis B virus is mandatory, and blood titer levels of antibodies anti-hepatitis B must be above 100 IU/L (country-specific legislative requirements need to be addressed). In addition, all work must be performed in biosafety level 2 laboratories (country-specific legislative requirements need to be addressed). Standard health and safety precautions associated with working in a laboratory environment and handling mammalian cell culture, including handling of waste, should be adopted when conducting the entire protocol.

PROTOCOL:

The work involving primary monocytes isolated from human blood was approved by the committee of the Federal Office for Public Health Switzerland (reference number: 611-1, Meldung A110635/2) for the Adolphe Merkle Institute.

1. Isolation of peripheral blood monocytes (PBMs) from human buffy coats

NOTE: The following section describes the isolation of immune cells from one 50 mL bag of a buffy coat, purchased from the Swiss Transfusion Centre in Bern, Switzerland.

1.1. Preparation of reagents

1.1.1. Prepare 100 mL of magnetic separation buffer per buffy coat: 0.5% (w/v) bovine serum albumin (BSA; in phosphate-buffered saline [PBS]) with 2 mM ethylenediaminetetraacetic acid (EDTA; adjust to pH = 7.2, sterile filter with 0.22 µm pore size). Keep at 4 °C throughout the procedure.

1.1.2. Prepare the cell culture medium (CCM): RPMI 1640 with 10% fetal bovine serum (FBS), 1% L-glutamine (here, 2 mM L-glutamine), and 1% penicillin-streptomycin (here, 100 units/mL penicillin and 100 µg/mL streptomycin).

NOTE: The required amount of each reagent depends on the number of cells to be seeded in the following steps.

1.2. Isolation of PBMs

NOTE: All glass and plasticware needs to be sterilized before use. For safety reasons, the use of plasticware is recommended when handling human blood samples to reduce the risk of injury with glassware.

1.2.1. Use scissors to cut open the hose end of the bag containing the buffy coat.

1.2.2. Distribute the buffy coat by pouring the bag's contents through the bag's duct directly into two conical centrifuge 50 mL tubes (~25 mL each).

1.2.3. Gently pour or pipette PBS into the tubes to reach 50 mL volumes. Mix the contents by turning the tube gently upside down 3x.

1.2.4. Divide the buffy coat-PBS mixture into four new 50 mL conical centrifuge tubes by pipetting 25 mL of the mixture in each fresh tube.

1.2.5. Slowly lay 13 mL of density gradient medium underneath the buffy coat-PBS mixture using a 10 mL serological pipette. Detach the filled pipette from the pipette holder and immediately plug the upper opening of the pipette with a thumb to prevent any additional leaking of the density gradient medium into the buffy coat-PBS mixture.

NOTE: Holding the upper opening with the thumb, place the filled pipette at the bottom of the conical centrifuge tube so that the density gradient medium slowly flows beneath the buffy coat-PBS mixture, leaving approximately 1 mL of the density gradient medium inside the pipette.

1.2.6. Repeat steps 1.2.4–1.2.6 with the other three tubes containing the buffy coat-PBS mixture.

1.2.7. Centrifuge all four tubes containing the mixtures for 20 min at $1,000 \times g$ and 25°C at a slow braking mode (here, maximum brake = 4). Use holders for centrifugation.

1.2.8. Open the lid of each tube, remove the upper layer containing plasma and platelets using a serological pipette, and dispose in a biohazard liquid waste container.

1.2.9. Use a serological pipette for collecting the peripheral blood mononuclear cell layer, which appears as a whitish turbid small fraction (~2–3 mm in thickness) between the plasma and the density gradient medium layers. The pellet contains red blood cells at the bottom. Avoid transferring erythrocytes that form the bottom-most layer. Repeat this for all four tubes.

NOTE: Peripheral blood mononuclear cells consist of PBMs and lymphocytes. PBMs will be separated from lymphocytes later during magnetic CD14+ separation.

1.2.10. Pool peripheral blood mononuclear cells from two tubes into one 50 mL tube.

1.2.11. Fill up the two tubes with PBS to 50 mL and cover with a lid.

1.2.12. Discard the leftover erythrocytes and plasma from the original four tubes in a biohazard liquid waste container.

1.2.13. Centrifuge the two tubes for 8 min at 500 x g and 18–20 °C at a regular centrifuge speed.

1.2.14. After centrifugation, remove the supernatant with a serological pipette and discard it in a biohazard liquid waste container.

1.2.15. Resuspend the cells with 5 mL of PBS using a serological pipette.

1.2.16. Pool the cell suspensions into one 50 mL conical centrifuge tube and fill to mL with PBS.

1.2.17. Use 5 µL of the cell suspension to count the cells with a cell counter using the trypan blue (45 µL) exclusion method.

1.2.18. Pipette 10 µL of the trypan blue-PBMs solution into a cell counter chamber and count the number of cells per the standard counting protocol. Use Equation 1 to calculate the total cell number, C_T .

$$C_T = \frac{\text{counted cells}}{\text{mL}} * 50 \quad (1)$$

1.2.19. After counting the cells, centrifuge the 50 mL tubes as done in step 1.2.13.

1.3. CD14 positive selection

1.3.1. Gently open the lid of each tube, then remove and discard the supernatant using a serological pipette without disturbing the pellet.

1.3.2. Add the calculated amount (Equation 2) of magnetic separation buffer (here, 80 µL of buffer per 1×10^7 total cells) and resuspend the cell pellet by pipetting the solution up and down.

$$\text{buffer volume} = \frac{C_T * 80}{10^7} = \mu\text{L} \quad (2)$$

1.3.3. Calculate using Equation 3 (here, 10 µL per 1×10^7 total cells) the corresponding volume of CD14+ magnetic beads and pipette the appropriate volume.

$$beads\ volume = \frac{c_T * 10}{10^7} = \mu L \quad (3)$$

1.3.4. Mix well by pipetting up and down, close the lid, and incubate the solution at 4 °C for 15 min.

1.3.5. Upon incubation, fill the tube up to 50 mL with magnetic separation buffer.

1.3.6. Centrifuge as done in step 1.2.13.

1.3.7. Aspirate and discard the supernatant using a serological pipette without disturbing the cell pellet.

1.3.8. Pipette the corresponding amount of magnetic separation buffer (Equation 4; here, 500 μ L of buffer per 1×10^8 cells) and gently mix by pipetting up and down 3x.

$$buffer\ volume = \frac{c_T * 50}{10^7} = \mu L \quad (4)$$

1.3.9. Disinfect the magnetic separation station by spraying and wiping it with a sterilizing agent. Place in the laminar flow hood along with the column for magnetic separation.

1.3.10. Place the magnetic separation column in the magnetic field and place an empty 50 mL conical centrifuge tube directly under the column for collecting the wash and unlabeled cells (i.e., waste).

1.3.11. Rinse the magnetic separation column by pipetting 3 mL of magnetic separation buffer into the column. Do not let the column dry out throughout the procedure.

1.3.12. Prepare a 15 mL conical centrifuge tube and pipette 1 mL of magnetic separation buffer.

1.3.13. Apply the cell suspension (prepared in step 1.3.8) to the magnetic separation column. In the 50 mL conical centrifuge tube under the filter, collect unlabeled cells that passed through.

NOTE: Do not exceed 2×10^9 cells per column to avoid blocking of the column.

1.3.14. As soon as the column reservoir is empty (i.e., when the cells have passed through the column), apply 3 mL of magnetic separation buffer using a serological pipette and let it pass through the column. Repeat this 3x.

1.3.15. Remove the magnetic separation column from the magnetic separator by gently pulling with hands, then place it into a 15 mL tube containing pre-pipetted 1 mL of magnetic separation buffer (prepared in step 1.3.12).

1.3.16. Add 5 mL of magnetic separation buffer to the column and flush out the magnetically labeled cells by firmly pushing the plunger into the column.

1.4. Reagent preparation for MDM and MDDC differentiation

1.4.1. Count the cells with a cell counter using the trypan blue exclusion method as done in steps 1.2.17 and 1.2.18.

1.4.2. Calculate the required volumes of CCM or FBS for further steps as follows: either the volume of CCM corresponding to a cell density of 1×10^6 cells/mL (step 1.4.4.1), or the volume of FBS corresponding to a cell density of 6×10^6 cells per 0.9 mL of FBS (step 1.4.4.2).

$$\text{Total cell number (} c_T \text{)} \quad c_T = \frac{\text{counted cells}}{\text{mL}} * 6$$

$$\text{CCM volume (step 1.5.1)} \quad CCM \text{ volume} = \frac{c_T}{10^6} = \mu L$$

$$\text{FBS volume (step 1.6.3)} \quad FBS \text{ volume} = \frac{c_T}{6 * 10^6} * 0.9 = \mu L$$

1.4.3. Close the lid, place the tube in the centrifuge, and centrifuge as done in step 1.2.13. Remove and discard the supernatant without disturbing the cell pellet. Proceed to step 1.5 for cell seeding or step 1.6 for cell freezing.

1.5. PBM seeding and differentiation into MDMs and MDDCs

1.5.1. Resuspend the cell pellet in the calculated volume of cell CCM as calculated in step 1.4.2 (here, a final concentration of 1×10^6 cells/mL) by pipetting up and down 3x.

1.5.2. Pipette the amount of cells intended to differentiate into MDMs and MDDCs in separate conical centrifuge tubes using a serological pipette.

1.5.3. Pipette differentiating factors to the CCM with PBMs and mix well by pipetting up and down. Differentiating factors are applied as follows:

1.5.3.1. For MDDCs: final concentration of 10 ng/mL interleukin-4 (IL-4) and 10 ng/mL granulocyte-macrophage colony-stimulating factor (GM-CSF).

1.5.3.2. For MDMs: final concentration of 10 ng/mL macrophage colony-stimulating factor (M-CSF).

1.5.4. Pipette the cell suspensions in CCM with the added differentiating factors into 6 well plates by distributing 3 mL of the suspension per well (corresponds to 3×10^6 cells/mL).

1.5.5. Place the 6 well plates in a cell culture incubator (37°C , 5% CO_2) and let them differentiate for 6 days without refreshing the CCM.

NOTE: Differentiation ranges from 5–8 days depending on the local availability of buffy coats and experimental setup, provided that the differentiation efficiency is determined using the appropriate techniques (refer to the discussion section).

1.6. PBM freezing

1.6.1. Resuspend the cell pellet in cryoprotective medium (here, FBS and dimethyl sulfoxide [DMSO; cytotoxic]) at a ratio of 9:1 (v/v) by pipetting a volume of prewarmed FBS. This corresponds to a final cell concentration of 6×10^6 cells/mL, considering a further addition of 10% DMSO (v/v).

1.6.2. Mark the desired number of cryovials in the laminar flow hood (i.e., record the date, isolation code, and number of cells).

1.6.3. Pipette 0.9 mL of cell suspension in pure FBS (here, 6×10^6 cells in 0.9 mL of FBS) to each cryovial. Subsequently, slowly pipette 0.1 mL of DMSO and mix the suspension well by turning the cryovials up and down 3x.

1.6.4. Transfer the cryovials to a cell-freezing container and immediately set it to -80°C for 24 h.

1.6.5. After 24 h, remove the cryovials from the -80°C freezer and container and place them into the liquid nitrogen tank suitable for cell storage.

1.7. PBM thawing and differentiation into MDMs and MDDCs

1.7.1. Warm all the required reagents to 37°C in a water bath (~20–30 min).

1.7.2. Prepare the appropriate number of 6 well plates corresponding to the number of thawed cells (here, one plate per 1.8×10^7 cells). Pipette 2 mL of CCM to each well under aseptic conditions. Place plates in the incubator (5 % CO_2 , 37°C) for 15 min to allow the pH to equilibrate.

1.7.3. Take the required amount of cryovials with frozen cells from a liquid nitrogen tank and gently swirl them in a 37°C water bath (1–2 min) to ensure uniform thawing of the cell suspension.

1.7.4. Remove the cryovial from the water bath and decontaminate with a sterilizing agent, ensuring that the agent does not interact with the lid and O-ring.

NOTE: Hereon, all steps must be completed under aseptic conditions.

1.7.5. Prepare the appropriate number of 15 mL conical centrifuge tubes corresponding to the number of cryovials to be thawed (here, 6×10^6 cells/tube). Pipette 9 mL of prewarmed CCM into each tube.

1.7.6. Pipette slowly (drop-by-drop) the contents of the cryovial into a tube containing CCM. Close the lid, repeat for each tube, and centrifuge at $200 \times g$ for 5 min at a regular centrifuge speed.

1.7.7. Discard the supernatant without disturbing the pellet.

1.7.8. Resuspend the pellet from each tube (which contains cells from one cryovial) in 2 mL of prewarmed CCM by pipetting up and down using a serological pipette (cell density corresponding to 3×10^6 cells/mL).

1.7.9. From each tube, pipette the resuspended cells into two wells of a 6 well plate containing 2 mL of previously prepared CCM (1 mL per well) to reach a cell density of 3×10^6 cells/well (corresponding to a final concentration of 1×10^6 cells/mL). Repeat this for all other tubes.

1.7.10. Proceed with differentiation as described in section 1.5.3.

1.7.11. Place the 6 well plates in a cell culture incubator (37°C , 5% CO_2) and let them differentiate for 6 days without refreshing the CCM.

2. Triple cell coculture model of human alveolar epithelial tissue

NOTE: This section provides instructions on volumes and cell numbers corresponding to 12 well plate inserts. **Figure 1** summarizes a proposed timeline for model assembly.

2.1. Epithelial cell (A549 cell line) seeding

2.1.1. Culture the epithelial cells according to the recommendations provided by the supplier (ATTC). Briefly, subculture the cells in CCM at 80% cell confluency (here, 2x–3x per week).

NOTE: Subculture A549 for at least four passages before the coculture model composition, using A549 cells in a passage range of 5–25.

2.1.2. Pipette 1.5 mL of prewarmed CCM into 12 well plates (the number of wells corresponds to the desired number of models).

2.1.3. Place individual 12 well cell culture inserts into wells of a 12 well plate using sterilized tweezers.

2.1.4. Detach the cells from a flask according to the subcultivation protocol (i.e., using a detaching agent, remove the agent by centrifugation as done in step 1.2.13). Resuspend in the appropriate volume of CCM corresponding to the final cell concentration of A549 (here, 50×10^4 cells/mL in 0.5 mL per insert and seeding density of 2.78×10^5 cells/cm²).

2.1.5. Pipette 0.5 mL of the cell suspension (i.e., 5.0×10^5 cells/mL) into the apical side of the insert using a 1 mL pipette.

2.1.6. Cover the plates with lids and place them in a cell culture incubator (37°C , 5 % CO_2) for 4 days.

NOTE: Regularly check the confluency of A549 cells under a phase-contrast microscope.

2.2. MDCC seeding

2.2.1. Aspirate CCM with unattached cells in the MDCCs containing 6 well plates.

2.2.2. Add 1 mL of fresh prewarmed CCM to each well.

2.2.3. Use a cell scraper, detach (scrape) adherent MDCCs from each well, gently wash the wells with the existing 1 mL of CCM 3x, and combine them into one conical centrifuge tube.

2.2.4. Count the cells with a cell counter using the trypan blue exclusion method using 10 μ L of cell suspension and 10 μ L of trypan blue solution.

2.2.5. Centrifuge the cell suspension as done in step 1.2.13.

2.2.6. Calculate the CCM volume required for resuspension (Equation 5):

Cell density (CD) = _____ cells/mL
Cell suspension volume (V_s) = _____ mL

$$C_T = \frac{CD}{mL} * V_s \quad (5)$$

Required MDCC density = 4.2×10^5 cells/mL; each insert requires 6.3×10^4 cells, which corresponds to a seeded cell density of 7×10^4 cell/cm².

CCM volume for cell resuspension (V_m ; Equation 6):

$$CCM \text{ volume} = \frac{CD}{42 \times 10^4} * V_s = mL \quad (6)$$

2.2.6. Gently aspirate and discard the CCM from the upper chamber of 12 well plates with growing A549 on the inserts.

2.2.7. Place the inserts with A549 cells in an upside-down position in a sterile Petri dish using sterilized tweezers. Prepare a conical centrifuge tube (50 mL) with PBS and pre-moisten a cell scraper.

2.2.8. Scrape off A549 cells from the basal surface of the insert (i.e., the upper part at the upside-down position), which should grow through the pores of the inserts.

NOTE: Rinse the scraper with PBS (prepared in a tube) between scraping individual samples and keep it wet throughout the procedure.

2.2.9. Upon centrifugation (step 2.2.5), aspirate and discard the supernatant, then redisperse the MDCC pellet in the calculated amount of CCM (step 2.2.6) and pipette up and down 3x.

2.2.10. Pipette 150 μ L of the cell suspension on top of each insert so that the entire basal surface of the insert is equally covered with the liquid and not containing bubbles.

2.2.11. Cover the dish with a lid and place in a cell culture incubator for 70 min. Aspirate the CCM from the cell culture plates (where the inserts have initially been placed), discard it in a biohazard liquid waste, and pipette 1.5 mL of fresh CCM into each well. Cover the plate with a lid and place in the cell incubator (37 °C, 5% CO₂).

NOTE: Do not exceed the time period mentioned above to avoid drying out the cells.

2.2.12. After the incubation, carefully hold each insert with sterilized tweezers and place them to the plates containing CCM in an ordinary position. Cover the plate with a lid and return it to the cell incubator (37°C, 5% CO₂).

2.3. Macrophage (MDM) seeding

2.3.1. Take 6 well plates containing predifferentiated MDMs. Place them from the cell incubator and place them to a laminar flow hood.

2.3.2. Aspirate and discard CCM with unattached MDMs grown in 6 well plates, and pipette 1 mL of fresh prewarmed CCM in each well.

2.3.3. Using a cell scraper, gently remove adherent MDMs from individual wells (as done in step 2.2.3 for the MDDCs).

2.3.4. Pipette 10 µL of trypan blue into a well or tube and add 10 µL of the MDM suspension to achieve the final dilution of 1:1 (v/v). Count the number of MDMs using the appropriate counting protocol.

2.3.5. Centrifuge the cell suspension as done in step 1.2.13.

2.3.6. Calculate the required volume (Equation 7):

Cell density (CD) = _____ cells/mL

Cell suspension volume (Vs) = _____ mL

Required MDM density = 2.5 x 10⁴ cells/mL in CCM (here, each insert requires 1.25 x 10⁴ cells, which corresponds to a seeded cell density of 2.5 x 10⁴ cell/cm²).

$$CCM\ volume = \frac{CD}{2.5 \times 10^4} * Vs = mL \quad (7)$$

2.3.7. Upon centrifugation, aspirate and discard the supernatant, redisperse the MDM pellet in the calculated amount of CCM (step 2.3.6), and pipette up and down 3x.

2.3.8. Carefully pipette 0.5 mL of the MDM suspension (prepared in step 2.3.7) on the wall of the cell culture inserts with A549 and MDDCs (not directly on epithelial cells) using a 1 mL pipette. Cover plates with lids and place in a cell culture incubator (37 °C, 5 % CO₂) for 24 h.

2.4. Transfer of coculture model to air-liquid interface (ALI)

2.4.1. Upon the end of the 24 h (± 2 h) incubation period of the assembled model in a cell culture incubator, aspirate and discard CCM from both apical and basal parts of cell culture inserts and from the wells.

2.4.2. Using sterilized tweezers, lift individual inserts from the wells and pipette 0.6 mL of fresh prewarmed CCM to each well using a 1 mL pipette. Do not add CCM to the apical side of the insert.

3. Exposure to selected positive controls (known stimuli for inducing proinflammatory response)

NOTE: Exposure of the coculture models to a known proinflammatory stimulus endotoxin lipopolysaccharide (LPS)⁷ and proinflammatory cytokine tumor necrosis factor α (TNF- α)⁷ is used to illustrate responsiveness of the model. Furthermore, exposure to a detergent (Triton X-100) is used to confirm the sensitivity of a lactate dehydrogenase (LDH) assay.

3.1. Prepare the positive control solutions: LPS stock (1 mg/mL in distilled water), TNF- α stock (100 μ g/mL in distilled water) and Triton-X 100 (2% [v/v] in PBS).

3.2. Upon 24 h incubation of the coculture model at ALI conditions, aspirate and discard the supernatant from the basal compartment. Using sterilized tweezers, lift individual inserts from the wells and pipette 0.6 mL of fresh prewarmed CCM into each well.

3.3. Prepare individual positive controls working solutions by diluting the stocks in CCM in conical centrifuge tubes as follows: 1 μ g/mL LPS, 1 μ g/mL TNF- α , and 0.2% Triton-X 100. The volumes correspond to the number of tested inserts (here, 100 μ L/insert). Mix the solutions well by pipetting up and down 3x.

3.4. Apply 100 μ L of each positive control solution by slowly pipetting onto the wall of the cell culture insert. Cover the well plate with a lid and place in the cell culture incubator (37 $^{\circ}$ C, 5% CO₂) for 24 h. Upon incubation, aspirate and discard the liquid at the apical side of the insert by holding individual inserts using tweezers.

3.5. Collect CCM in the basal compartments and store at 1) 4 $^{\circ}$ C for further LDH analysis, denoting cell membrane rupture-mediated cytotoxicity, and/or 2) store at -80 $^{\circ}$ C for further analysis of protein release via enzyme-linked immunosorbent assay (ELISA). Run the assays according to the kit supplier recommendations.

3.6. Upon the removal of CCM, wash the inserts with PBS 3x and fix the cells on cell culture inserts in 4% paraformaldehyde (in PBS, 15 min at room temperature [RT]) by ensuring both insert sides are well-covered with PFA solution. Subsequently wash 3x with PBS to remove PFA. Store the samples submerged in PBS at 4 $^{\circ}$ C for further immunostaining (an example of this method has been described previously¹⁷).

REPRESENTATIVE RESULTS:

Human lung coculture models, composed of alveolar epithelial cells and immune cells, were assembled either from fresh or frozen MDDCs and MDMs progenitors (here, human

peripheral blood-derived monocytes). As presented in **Figure 1**, A549 cells were seeded 3 days after the first section involving monocyte isolation/thawing. After 6 days of differentiation, the differentiated MDMs appeared round-shaped, whereas MDCCs formed a more elongated shape with observable protrusions. They also appeared as agglomerates, especially when differentiated from fresh monocytes (**Figure 2, Figure 3**). Epithelial cells formed a dense cell layer of cells after 3 days of growth on membrane inserts (**Figure 4**), when the cocultures were assembled. After 24 h of assembling and an additional 24 h of being subjected to ALI conditions, the cocultures were prepared for exposures.

Responsiveness of the 3D cell culture models was investigated upon exposure to known proinflammatory stimuli using a pseudo-ALI approach, as described previously²⁹. The proinflammatory stimuli, LPS and TNF- α , were added in low volumes (100 μ L) onto the apical surface of the air-exposed cell model. In parallel, the absence of membrane rupture as a measure of cytotoxicity was assessed via LDH assay. A significant increase in LDH release in CCM of the basal compartment was observed upon exposure to the positive control for membrane rupture, a detergent Triton-X 100 (**Figure 5**). These results proved responsiveness of the model to a cytotoxic substance, whereas no increase in LDH release was observed upon apical stimulation with TNF- α or LPS.

A possible reason for the different measured values of LDH in the samples assembled with either fresh or previously frozen PBMs can be attributed to the sample storage. Samples from fresh PBMs were stored for a longer time at -80 °C; therefore, the activity of LDH enzyme may drop. Notably, LDH is stable for only up to 4 days in CCM; thus, it is recommended to perform the assay the latest 2 days after collecting the supernatants. Alternatively, it is possible to freeze down the supernatants directly after collection. However, it is important to consider that freezing can decrease the enzymatic activity of LDH.

Secretion of proinflammatory mediators (here, TNF- α and interleukins 6 [IL-6] and 8 [IL-8]) into the basal CCM was quantified via ELISA. Statistically significant ($p < 0.05$, one-way ANOVA) increases in the release of IL-6 and IL-8 were observed in both LPS- and TNF- α -treated samples compared to respective untreated cells, as well as in the cell culture models assembled from either PBMs source (**Figure 6**). Although concentrations (pg/mL) of all the tested cytokines in the basal CCM were higher in the cocultures composed of fresh PBMs, differences among the two cocultures and monocultures were not statistically significant ($p > 0.05$) (**Figure 6**). To confirm the added value of coculture models with respect to a 2D epithelial cell culture, A549 monocultures were also exposed to LPS or TNF- α . As expected, the release of all investigated mediators from A549 monocultures was lower compared to both coculture models; although, the difference between them was not statistically significant ($p > 0.05$, one-way ANOVA).

Cellular morphology of the 3D human alveolar epithelial tissue barrier was assessed via confocal laser scanning microscopy (LSM). To visualize the composition of each model, macrophages within the coculture models (MDMs) were stained with mature macrophage marker 25F9 (also known as cluster of differentiation CD163). MDCCs were stained with CD83, which is an important marker for activated dendritic cells³⁰. Regarding cellular morphology, no difference was observed between the coculture models using MDMs and

MDDCs from fresh PBMs compared to those using thawed PBMs. In LPS- and TNF- α -exposed cocultures, both composed of fresh and frozen immune cells, a disrupted epithelial layer in LSM images was observed, which was not the case in untreated cells (**Figure 7, Figure 8**).

FIGURE LEGENDS:

Figure 1: Schematic timeline of the protocol. Presentation of the 3D coculture model preparation, assembly, and application (exposure to a tested substance). ALI = air-liquid interface, MDDCs = monocyte-derived dendritic cells, MDMs = monocyte-derived macrophages, PBMs = peripheral blood monocytes.

Figure 2: MDMs and MDDCs differentiated from fresh PBMs. Phase-contrast microscopy image of differentiated (A) MDMs and (B) MDDCs from fresh PBMs (6 days after cell isolation). MDMs are round-shaped, while MDDCs are often observed as agglomerates.

Figure 3: MDMs and MDDCs differentiated from frozen PBMs. Phase-contrast microscopy image of differentiated (A) MDMs and (B) MDDCs from thawed PBMs (6 days after thawing). MDMs are round-shaped, but some elongated cells can be observed. MDDCs also appear round-shaped with protrusions.

Figure 4: Epithelial cells growth on membrane inserts. Phase-contrast microscopy image of confluent A549 growing on a membrane insert 4 days after seeding, forming a dense layer of cells.

Figure 5: Cytotoxicity results investigated via membrane rupture-based (LDH) assay. The data is presented as a fold increase over untreated cells (mean \pm SD, n = 3, asterisk denotes statistically significant increase compared to untreated cells, **p < 0.01, ****p < 0.0001). In green models MDMs and MDDCs from fresh PBMs are represented, and in purple models those assembled from thawed PBMs are represented.

Figure 6: Proinflammatory reactions in the cocultures and monocultures. The proinflammatory mediators (TNF- α , IL-6, and IL-8) release in the cocultures upon 24 h challenge with LPS or TNF- α . The data is presented as relative to untreated cells (mean \pm SD, n = 3, **p < 0.01, ***p < 0.001, ****p < 0.0001). In green models MDMs and MDDCs from fresh PBMs are represented, and in purple models those assembled from thawed PBMs are represented. Grey represents A549 monocultures.

Figure 7: Morphology of cocultures composed of fresh immune cells. LSM images of apical and basal sides of the coculture model with xz projections of apical sides of the model using MDMs and MDDCs from fresh PBMs. Cyan represents nuclei (DAPI), magenta represents cytoskeleton (rhodamine-phalloidin), white represents MDMs (25F9, also known as CD163), and green represents MDDCs (CD 83). The white arrow denotes MDM, whereas the green arrow denotes MDDC.

Figure 8: Morphology of cocultures composed of frozen immune cells. LSM images of the apical side of the coculture model with corresponding xz projections, and basal side of the model using MDMs and MDDCs from thawed PBMs. Cyan represents nuclei (DAPI), magenta

represents cytoskeleton (rhodamine-phalloidin), white represents MDMs (25F9, also known as CD163), and green represents MDDCs (CD 83). The white arrow denotes MDM, whereas the green arrow denotes MDDC.

DISCUSSION:

Emerging production of novel materials, including chemicals and drugs, gradually increases the need for predictive in vitro models. To comply with the three principles of replacement, reduction, and refinement of animal testing³², in vitro cell models have become powerful tools for elucidating mechanisms of a drug's or material's action⁸⁻¹¹. Presented here is a detailed protocol of assembling the multicellular model using immune cells that are either freshly isolated or thawed from previously frozen monocytes. Also described is the cultivation of the model at ALI. Finally, the protocol illustrates an example of exposure to proinflammatory stimuli and compares the response of the two models containing either fresh or frozen monocytes.

Various studies have been performed to confirm and justify the added value of the enhanced complexity of the models grown and exposed under ALI conditions compared to conventional submerged exposure^{7,22,31}. Observation of the higher proinflammatory response in cocultures compared to monocultures of epithelial cells confirm a previous study. The study used the presented coculture model (stimulated with LPS) and showed a higher response at gene expression levels of *TNF* and *IL1B* compared to the A549 monoculture equivalent model⁷. On the other hand, both models showed higher variations within measured proinflammatory mediator release values compared to A549 monocultures. This can be explained by the use of immune cells from different donors (buffy coats) within biological replicates (i.e., one repetition, one donor), as previously shown⁷. If desired, the variations among the replicates can be overcome by 1) using thawed PBMs from the same donor or 2) pooling PBMs from different donors prior to freezing the cells, then subsequent usage of the same pool in each repetition. Including more biological replicates is also recommended.

The cell-freezing technique can be considered as a critical step; however, it is a common laboratory procedure for preserving cells for phenotypic and functional analysis. Various studies have demonstrated that the quality of frozen PBMs is vital to their survival, and an appropriate freezing technique is key to the success of subsequent assays with the same cells^{28,32}. Modification of the protocol can be performed by freezing PBMs, which provides flexibility in the experimental setup, as the availability of buffy coats is usually limited. Another advantage of using frozen PBMs (in several vials) over freshly isolated ones is that they can be used in subsequent experiments even after 1 year. This decreases the potential issue of donor-to-donor variability if this is a desired or required parameter in an experimental.

Outcomes of an interlaboratory comparison performed after up to 13 months show that PBMs, when properly stored in a liquid nitrogen tank, can be used over a long period without any effect on cell viability or cell recovery³³. Longer storage times (over 1 year) may be possible upon careful validation of cell viability and cell responsiveness before performing an experiment. Also, the temperature in the liquid nitrogen tank must remain stable at all

times. The main factor affecting the viability of cryopreserved PBMs was found to be DMSO concentration, with an optimal concentration of 10%–20 % (v/v)²⁸. To minimize potentially harmful effects of freezing, different sources of proteins, FBS or BSA (with a broad range of concentration from 40% up to 100 %³⁴) are often added to the freezing medium as natural protective components that can increase cell survival.

Due to the high cytotoxic potential of DMSO, it is recommended first to disperse PBMs in FBS, then add DMSO to PBMs already dispersed in FBS. Notably, although higher FBS concentrations (>40%) did not show any improvement in cell viability. However, at the same time, they did not cause harm to the cells²⁸. Nevertheless, freezing monocytes is a possible approach to overcoming issues of limited buffy coat availability. However, if the use of MDDCs and MDMs from fresh PBMs is desired, the immune cells can be differentiated and used 5–8 days after isolation^{7,16,17,35–37}. If experimental planning allows, at least 6 days of differentiation in both MDDCs and MDMs is recommended. However, consistency among different repetitions in the same experiment, along with routine inspections of their specific surface marker expressions, are crucial. The responsiveness to a proinflammatory stimulus, such as LPS, after the differentiation time should also be regularly checked.

Many investigations using the A549 cell line have been performed at ALI, either as a monoculture or combined with other cell types (macrophages, dendritic cells, or fibroblasts) into 3D coculture model^{22,24,29,38}. Using this 3D coculture model, the cytotoxicity, oxidative stress, or proinflammatory effect of (nano-)materials have been investigated for up to 72 h^{1,17,21,24,29}. The model's resemblance to in vivo tissue has previously been investigated based on confocal laser scanning imaging of the model¹⁶. When assembling the model, it is important to consider both cell proliferation (which can affect A549 in the model presented here) as well as performance of the primary (not proliferating) immune cells (here, MDDCs and MDMs). It is also important to consider that not all the CD14 positive monocytes differentiate into MDDCs and MDMs, and that the cells can be present both attached and suspended forms. Based on the nature of the coculture assembly (here, both cell types need to attach to the existing epithelial layer), it is recommended to use only the adherent subpopulations of both immune cell types. Additionally, routine analyses of monocytes, MDDC and MDM monoculture responsiveness to LPS, and expression of specific surface markers (CD14, CD163, CD86, CD93, or CD206, data not shown) have suggested that 6 and 7 days of differentiation are the optimum timepoints.

Although a realistic number of alveolar epithelial cells in human lungs corresponds to ~160,000 cells/cm², the number of A549 cells counted in the model is ~1,000,000 cells/cm² after 9 days cultured on the insert^{16,18}. Thus, this in vitro model's limitations need to be considered. First, the density of epithelial cells was established based on their ability to form a confluent layer on the growing membrane. It is also important to mention that the A549 represents an epithelial type II cell with a cuboidal form, contrary to epithelial type I cells, which are flat and outspread. On the other hand, the required number of immune cells was established based on the literature and presented in this protocol as cell number/surface area^{39–41}. The cell density of MDDCs in the range of 400 cells/mm² (4 cells/cm²)¹⁶ is comparable to the steady-state cell density of 500–750 cells/mm² (5– 7 cells/cm²) reported from in vivo studies³⁹. The density of MDMs in this model is within the same range of in vivo situation in the human alveolar region⁴⁰.

Mature macrophage staining (25F9) was observed both in the apical side (where MDMs are present) as well as basal side (i.e., at the site of dendritic cells). Translocation of immune cells through the membrane inserts pores is possible and has also been observed using this model¹⁶, which may explain the observed differences in staining intensities. However, another possible explanation is that the mature macrophage marker can also be expressed on dendritic cells, but the expression is highly donor-specific⁴². Also, the intensity of the 25F9/CD163 expression is much higher in MDMs (**Figure 7, Figure 8**). Both proinflammatory stimuli (LPS and TNF- α) affected the integrity of the pulmonary epithelial barrier in both cocultures (**Figure 7, Figure 8**). This was expected based on previous publications^{43,44} showing that proinflammatory cytokines and bacterial products disrupt the integrity of epithelial barriers.

The 3D multicellular model of the human alveolar epithelium, established and characterized previously¹⁷, has served as a powerful and useful tool for assessing biological responses (i.e., acute proinflammatory reactions, oxidative stress response, particle distribution, and cellular communication) in vitro^{21,24,25,45}. The results confirm the responsiveness of coculture models to proinflammatory stimuli (here, LPS and TNF- α). The response was slightly increased when using immune cells from fresh PBMs; however, there was no statistically significant difference between cocultures using fresh vs. thawed PBMs. Furthermore, the proinflammatory reactions of both coculture models were higher than those of epithelial cell monocultures cultivated under the same (ALI) conditions. In summary, the protocol describes the assembly of a 3D human alveolar epithelial tissue coculture model using either fresh or thawed PBMs for differentiation into MDMs and MDDCs. It is shown that both models are highly responsive to proinflammatory stimuli; therefore, they can serve as powerful tools for potential hazard and toxicity assessments.

DISCLOSURES:

The authors have nothing to disclose.

ACKNOWLEDGMENTS:

The authors would like to thank Dr. Miguel Spuch-Calvar for the coculture scheme in Figure 3 and to Dr. Bedia Begum Karakocak for critical reading. This study was supported by the PATROLS project, European Union's Horizon 2020 Research and Innovation Programme under grant agreement No 760813, and by the Adolphe Merkle Foundation. B.D. thanks the Peter und Traudl Engelhorn foundation for financial support.

REFERENCES:

- 1 Rothen-Rutishauser, B., Blank, F., Mühlfeld, C., Gehr, P. In vitro models of the human epithelial airway barrier to study the toxic potential of particulate matter. *Expert Opinion on Drug Metabolism and Toxicology*. **4** (8), 1075-1089 (2008).
- 2 Giard, D. et al. In vitro cultivation of human tumors: establishment of cell lines derived from a series of solid tumors. *Journal of National Cancer Institute*. **51** (5), 1417-1423 (1973).
- 3 Ochs, M., Weibel, E. R. in *Fishman's Pulmonary Diseases and Disorders, 5e.* eds MA Grippi et al. Ch. 2: Functional Design of the Human Lung for Gas Exchange (McGraw-Hill Education, 2008).

798 4 Foster, K. A., Oster, C. G., Mayer, M. M., Avery, M. L., Audus, K. L. Characterization of
799 the A549 Cell Line as a Type II Pulmonary Epithelial Cell Model for Drug Metabolism.
800 *Experimental Cell Research*. **243** (2), 359-366 (1998).

801 5 Guo, X.-Y., Lu, M., Chen, X.-Q., He, F.-D., Li, A. Correlation study of biological
802 characteristics of non-small cell lung cancer A549 cells after transfecting plasmid by
803 microbubble ultrasound contrast agent. *Asian Pacific Journal of Tropical Medicine*. **9** (6), 582-
804 586 (2016).

805 6 Cooper, J. R. et al. Long Term Culture of the A549 Cancer Cell Line Promotes
806 Multilamellar Body Formation and Differentiation towards an Alveolar Type II Pneumocyte
807 Phenotype. *PLoS ONE*. **11** (10), e0164438-e0164438 (2016).

808 7 Bisig, C., Voss, C., Petri-Fink, A., Rothen-Rutishauser, B. The crux of positive controls -
809 Proinflammatory responses in lung cell models. *Toxicology In Vitro*. **54**, 189-193 (2019).

810 8 Rothen-Rutishauser, B. et al. A newly developed in vitro model of the human
811 epithelial airway barrier to study the toxic potential of nanoparticles. *ALTEX*. **25**, (2008).

812 9 Braakhuis, H. M. et al. Progress and future of in vitro models to study translocation of
813 nanoparticles. *Archives of Toxicology*. **89** (9), 1469-1495 (2015).

814 10 Thai, P., Chen, Y., Dolganov, G., Wu, R. Differential regulation of MUC5AC/Muc5ac
815 and hCLCA-1/mGob-5 expression in airway epithelium. *American Journal of Respiratory Cell
816 and Molecular Biology*. **33** (6), 523-530 (2005).

817 11 Wu, J. et al. Characterization of air-liquid interface culture of A549 alveolar epithelial
818 cells. *Brazilian Journal of Medical and Biological Research*. **51** (2), e6950-e6950 (2017).

819 12 Shapiro, D. I., Nardone, L. L., Rooney, S. A., Motoyama, E. K., Munoz, J. L.
820 Phospholipid biosynthesis and secretion by a cell line (A549) which resembles type II
821 alveolar epithelial cells. *Biochimica and Biophysica Acta*. **530** (2), 197-207 (1978).

822 13 Balis, J., Bumgarner, S. D., Paciga, J. E., Paterson, J. F., Shelley, S. A. Synthesis of lung
823 surfactant-associated glycoproteins by A549 cells: description of an in vitro model for human
824 type II cell dysfunction. *Experimental Lung Research*. **6** (3-4), 197-213 (1984).

825 14 Schurch, S., Gehr, P., Im Hof, V., Geiser, M., Green, F. Surfactant displaces particles
826 toward the epithelium in airways and alveoli. *Respiration Physiology*. **80** (1), 17-32 (1990).

827 15 Gehr, P., Schurch, S., Berthiaume, Y., Hof, V. I., Geiser, M. Particle Retention in
828 Airways by Surfactant. *Journal of Aerosol Medicine*. **3** (1), 27-43 (2009).

829 16 Blank, F., Rothen-Rutishauser, B., Gehr, P. Dendritic Cells and Macrophages Form a
830 Transepithelial Network against Foreign Particulate Antigens. *American Journal of
831 Respiratory Cell and Molecular Biology*. **36** (6), (2007).

832 17 Rothen-Rutishauser, B. M., Kiama, S. G., Gehr, P. A three-dimensional cellular model
833 of the human respiratory tract to study the interaction with particles. *American Journal of
834 Respiratory Cell and Molecular Biology*. **32**, (2005).

835 18 Blank, F., Rothen-Rutishauser, B. M., Schurch, S., Gehr, P. An optimized in vitro model
836 of the respiratory tract wall to study particle cell interactions. *Journal of Aerosol Medicine*.
837 **19**, (2006).


838 19 Jardine, L. et al. Lipopolysaccharide inhalation recruits monocytes and dendritic cell
839 subsets to the alveolar airspace. *Nature Communications*. **10** (1), 1999, (2019).


840 20 Kopf, M., Schneider, C., Nobs, S. P. The development and function of lung-resident
841 macrophages and dendritic cells. *Nature Immunology*. **16** (1), 36-44 (2015).

842 21 Chortarea, S. et al. Repeated exposure to carbon nanotube-based aerosols does not
843 affect the functional properties of a 3D human epithelial airway model. *Nanotoxicology*. **9**
844 (8), 983-993 (2015).

845 22 Hilton, G., Barosova, H., Petri-Fink, A., Rothen-Rutishauser, B., Bereman, M.
 846 Leveraging proteomics to compare submerged versus air-liquid interface carbon nanotube
 847 exposure to a 3D lung cell model. *Toxicology In Vitro*. **54**, 58-66 (2019).
 848 23 Brandenberger, C. et al. Effects and uptake of gold nanoparticles deposited at the air-
 849 liquid interface of a human epithelial airway model. *Toxicology and Applied Pharmacology*.
 850 **242**, (2010).
 851 24 Drasler, B. et al. Single exposure to aerosolized graphene oxide and graphene
 852 nanoplatelets did not initiate an acute biological response in a 3D human lung model.
 853 *Carbon*. **137**, 125-135 (2018).
 854 25 Durantie, E. et al. Carbon nanodots: Opportunities and limitations to study their
 855 biodistribution at the human lung epithelial tissue barrier. *Biointerphases*. **13**, (2018).
 856 26 Brandenberger, C. et al. Quantitative evaluation of cellular uptake and trafficking of
 857 plain and polyethylene glycol-coated gold nanoparticles. *Small*. **6** (15), 1669-1678 (2010).
 858 27 Tomašek, I. et al. Combined exposure of diesel exhaust particles and respirable
 859 Soufrière Hills volcanic ash causes a (pro-)inflammatory response in an in vitro multicellular
 860 epithelial tissue barrier model. *Particle and Fibre Toxicology*. **13** (1), 67 (2016).
 861 28 Nazarpour, R. et al. Optimization of Human Peripheral Blood Mononuclear Cells
 862 (PBMCs) Cryopreservation. *International Journal of Molecular and Cellular Medicine*. **1** (2),
 863 88-93 (2012).
 864 29 Endes, C. et al. An in vitro testing strategy towards mimicking the inhalation of high
 865 aspect ratio nanoparticles. *Particle and Fibre Toxicology*. **11** (1), (2014).
 866 30 Ju, X. et al. The Analysis of CD83 Expression on Human Immune Cells Identifies a
 867 Unique CD83+-Activated T Cell Population. *Journal of Immunology*. **197** (12), 4613-4625
 868 (2016).
 869 31 Lenz, A. G. et al. Inflammatory and Oxidative Stress Responses of an Alveolar
 870 Epithelial Cell Line to Airborne Zinc Oxide Nanoparticles at the Air-Liquid Interface: A
 871 Comparison with Conventional, Submerged Cell-Culture Conditions. *BioMed Research*
 872 *International*. **12** (2013).
 873 32 Germann, A., Schulz, J. C., Kemp-Kamke, B., Zimmermann, H., von Briesen, H.
 874 Standardized serum-free cryomedia maintain peripheral blood mononuclear cell viability,
 875 recovery, and antigen-specific T-cell response compared to fetal calf serum-based medium.
 876 *Biopreservation and Biobanking*. **9** (3), 229-236 (2011).
 877 33 Weinberg, A. et al. Optimization and Limitations of Use of Cryopreserved Peripheral
 878 Blood Mononuclear Cells for Functional and Phenotypic T-Cell Characterization. *Clinical and*
 879 *Vaccine Immunology*. **16** (8), 1176 (2009).
 880 34 Freshney, R. I. in *Culture of animal cells: a manual of basic technique*. ed. R Ian
 881 Freshney. 321-334 (Wiley-Liss Inc., 2005).
 882 35 Lehmann, A. B. C., Blank, F., Gehr, P., Rothen-Rutishauser, B. in *Alternatives to animal*
 883 *testing*. eds ML Yarmush & RS Langer. 239-260 (Artech House, 2010).
 884 36 Steiner, S. et al. Reduction in (pro-)inflammatory responses of lung cells exposed
 885 in vitro to diesel exhaust treated with a non-catalyzed diesel particle filter. *Atmospheric*
 886 *Environment*. **81**, 117-124 (2013).
 887 37 Martinez, F. O., Gordon, S., Locati, M., Mantovani, A. Transcriptional Profiling of the
 888 Human Monocyte-to-Macrophage Differentiation and Polarization: New Molecules and
 889 Patterns of Gene Expression. *The Journal of Immunology*. **177** (10), 7303 (2006).

- 38 Chortarea, S. et al. Profibrotic activity of multi-walled carbon nanotubes upon prolonged exposures in different human lung cell types *Applied In Vitro Toxicology*. **5** (1), (2019).
- 39 Holt, P. G. Pulmonary Dendritic Cells in Local Immunity to Inert and Pathogenic Antigens in the Respiratory Tract. *Proceedings of the American Thoracic Society*. **2** (2), 116-120 (2005).
- 40 Pinkerton, K. E., Gehr, P., Castañeda, A., Crapo, J. D. in *Comparative Biology of the Normal Lung (Second Edition)*. ed Richard A. Parent. 105-117 (Academic Press, 2015).
- 41 Crapo, J., Barry, B., Gehr, P., Bachofen, M., Weibel, E. R. Cell number and cell characteristics of the normal human lung. *American Review of Respiratory Disease*. **126** (2), 332-337 (1982).
- 42 Maniecki, M. B., Møller, H. J., Moestrup, S. K., Møller, B. K. CD163 positive subsets of blood dendritic cells: The scavenging macrophage receptors CD163 and CD91 are coexpressed on human dendritic cells and monocytes. *Immunobiology*. **211** (6), 407-417 (2006).
- 43 Chignard, M., Balloy, V. Neutrophil recruitment and increased permeability during acute lung injury induced by lipopolysaccharide. *American Journal of Physiology-Lung Cellular and Molecular Physiology*. **279** (6), L1083-L1090 (2000).
- 44 Coyne, C. B. et al. Regulation of Airway Tight Junctions by Proinflammatory Cytokines. *Molecular Biology of the Cell*. **13** (9), 3218-3234 (2002).
- 45 Durantie, E. et al. Biodistribution of single and aggregated gold nanoparticles exposed to the human lung epithelial tissue barrier at the air-liquid interface. *Particle and Fibre Toxicology*. **14** (49) (2017).













Figure 4











Laboratory equipment

Benchmark microplate reader
Cell culture Incubator
Cell freezing container (for example Mr. Frosty)
Centrifuge
Confocal laser scanning microscope
Heamatocytometer, or automatic cell counter
Laminar bio-safety hood class II
MultiStand Macs (Macs Cell Separator)
pH meter
Phase contrast inverted light microscope
Pipette boy, pipettors (different volumes)
Scissors
Vacuum pump
Water bath

Disposable small equipment/glassware

15 mL and 50 mL conical centrifuge tubes
6- and 12-well cell culture plates, flat bottom, low evaporation lid, sterile
Cell culture inserts, transparent PET membrane, 12-well, 3 μm pore size
Cell scrapper
Cryovials
Glass autoclaved Petri Dishes
LS Columns
Sterile filtration cup for vacuum filtration, 0.2 μm pore size
Sterile Lab Bottle compatible with Filtration cup (min. 100 mL)
Sterile pipettes

Chemicals

Bovine serum albumine (BSA)
CD14+ MicroBeads human - magnetic beads
Deattachig agent Trypsin-EDTA, 0.05%, phenol red
Density gradient medium Lymphoprep
Dimethyl Sulfoxide (DMSO)
Ethylenediaminetetraacetic acid (EDTA)
Fetal bovine serum (heat inactivated)

Human granulocyte-macrophage colony-stimulating factor (GM-CSF), premium grade

Human Interleukin 4 (IL-4), premium grade
Human macrophage colony-stimulating factor (M-CSF), premium grade
L-glutamine
Lipopolysaccharid (LPS) from Escherichia coli
Paraformaldehyde (PFA)
Penicilin-Streptomycin
Phosphate Buffer Saline (PBS)
Roswell Park Memorial Institute-1640 Medium (RPMI)
Triton X-100
Trypan blue solution (0.4%)

Tumor necrosis factor alpha (TNF- α)

Assays used for cytotoxicity, (pro-)inflammatory response

Cytotoxicity Detection Kit (LDH)

Human IL-6 DuoSet ELISA

Human IL-8/CXCL8 DuoSet ELISA

Immunostaining

4',6-diamidino-2-phenylindole (DAPI), concentration 2 $\mu\text{g/mL}$

Goat anti-mouse IgG (H+L) Alexa Fluor 647 conjugated, concentration 20 $\mu\text{g/mL}$

Goat anti-rabbit IgG antibody (H+L) Dylight 488 conjugated, concentration 10 $\mu\text{g/mL}$

Mature Macrophage Marker Monoclonal Antibody, concentration 50 $\mu\text{g/mL}$

Phalloidin rhodamine, concentration 0.264 μM

Recombinant Anti-CD83 antibody, 1:50 dilution

Company

does not have to be specific, for example BioRad, Cressier, Switzerland

does not have to be specific

does not have to be specific

does not have to be specific

does not have to be specific, for example

does not have to be specific

does not have to be specific

Miltenyi, Germany

does not have to be specific

does not have to be specific

do not have to be specific

do not have to be specific

does not have to be specific

does not have to be specific

does not have to be specific

Falcon, Switzerland

Falcon, Switzerland

does not have to be specific, for example VWR, Switzerland

do not have to be specific

do not have to be specific

Miltenyi, Germany

does not have to be specific, for example VWR, Switzerland

does not have to be specific

do not have to be specific

Sigma-Aldrich, Switzerland

Miltenyi, Germany

Gibco, Switzerland

Alere Technologies AS, Norway

Sigma Aldrich, Switzerland

Sigma-Aldrich, Switzerland

Gibco, Switzerland

Miltenyi, Germany

Miltenyi, Germany

Miltenyi, Germany

Gibco, Switzerland

Sigma-Aldrich, Switzerland

Sigma-Aldrich, Switzerland

Gibco, Switzerland

Gibco, Switzerland

Gibco, Switzerland

Sigma-Aldrich, Switzerland

Sigma Aldrich, Switzerland

Immunotools

Roche, Switzerland

R&D, Biotechne, Switzerland

R&D, Biotechne, Switzerland

Sigma-Aldrich, Switzerland

Abcam, UK

Agrisera, Sweden

eBioScience, Thermo Fischer, Switzerland

Molecular Probes, Life Technologies, Switzerland

Abcam, UK

Type or catalogue number

Zeiss LSM 710 meta

130-042-303

Catalogue Number

353046 and 353043

353181

353085

130-042-401

10040-446

A7030-100g

130-097-052

25300054

1114547

D5879_1L

E6758-100g

10270-106

130-093-864

130-095-373

130-096-485

25030-024

4524-5mg

158127

31870-025

14190-094

11875093

T8787

11343015

11644793001

DY206

DY208

28718-90-3

ab150115

AS09 633

14-0115-82

R415

ab244204



Prof. Dr. Barbara Rothen-Rutishauser
Co-Chair BioNanomaterials
Adolphe Merkle Institute
University of Fribourg
Chemin des Verdiers 4
1700 Fribourg
Switzerland

Dr. Vineeta Bajaj
Review Editor

JoVe

T +41 26 300 95 02
www.am-institute.ch
barbara.rothen@unifr.ch

Fribourg, 27.2.2020

Manuscript Revision JoVE61090 - [EMID:9a197a938bf3fff0]

Dear Dr. Bajaj,

Please find attached author's letter to the reviewers' comments and the revised manuscript entitled "Multicellular human alveolar model composed of epithelial cells and primary immune cells for hazard assessment" by Hana Barosova, Barbara Drasler, Alke Petri-Fink, and myself.

We have carefully addressed all the Editorial comments in the revised version. We hope that you will find these changes adequate and further consider our article for publication in the highly regarded journal JoVE. Please do not hesitate to contact us if you have any queries or concerns regarding this rebuttal letter or changes to the manuscript.

We look forward to hearing from you at your earliest convenience.

On behalf of all authors, our kindest and sincerest regards,

Prof. Barbara Rothen-Rutishauser
Co-Chair BioNanomaterials

ARTICLE AND VIDEO LICENSE AGREEMENT

Title of Article:	Multicellular human alveolar model composed of epithelial cells and primary immune cells for hazard assessment
Author(s):	Hana Barosova, Barbara Drasler, Alke Petri-Fink, Barbara Rothen-Rutishauser

Item 1: The Author elects to have the Materials be made available (as described at <http://www.jove.com/publish>) via:



Standard Access



Open Access

Item 2: Please select one of the following items:



The Author is **NOT** a United States government employee.



The Author is a United States government employee and the Materials were prepared in the course of his or her duties as a United States government employee.



The Author is a United States government employee but the Materials were NOT prepared in the course of his or her duties as a United States government employee.

ARTICLE AND VIDEO LICENSE AGREEMENT

1. **Defined Terms.** As used in this Article and Video License Agreement, the following terms shall have the following meanings: “**Agreement**” means this Article and Video License Agreement; “**Article**” means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein; “**Author**” means the author who is a signatory to this Agreement; “**Collective Work**” means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; “**CRC License**” means the Creative Commons Attribution-Non Commercial-No Derivs 3.0 Unported Agreement, the terms and conditions of which can be found at: <http://creativecommons.org/licenses/by-nc-nd/3.0/legalcode>; “**Derivative Work**” means a work based upon the Materials or upon the Materials and other pre-existing works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; “**Institution**” means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; “**JoVE**” means MyJoVE Corporation, a Massachusetts corporation and the publisher of The Journal of Visualized Experiments; “**Materials**” means the Article and / or the Video; “**Parties**” means the Author and JoVE; “**Video**” means any video(s) made by the Author, alone or in conjunction with any other parties, or by JoVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion

of the Article, and in which the Author may or may not appear.

2. **Background.** The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JoVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.

3. **Grant of Rights in Article.** In consideration of JoVE agreeing to publish the Article, the Author hereby grants to JoVE, subject to **Sections 4** and **7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the “Open Access” box has been checked in **Item 1** above, JoVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.

ARTICLE AND VIDEO LICENSE AGREEMENT

4. **Retention of Rights in Article.** Notwithstanding the exclusive license granted to JoVE in **Section 3** above, the Author shall, with respect to the Article, retain the non-exclusive right to use all or part of the Article for the non-commercial purpose of giving lectures, presentations or teaching classes, and to post a copy of the Article on the Institution's website or the Author's personal website, in each case provided that a link to the Article on the JoVE website is provided and notice of JoVE's copyright in the Article is included. All non-copyright intellectual property rights in and to the Article, such as patent rights, shall remain with the Author.

5. **Grant of Rights in Video – Standard Access.** This **Section 5** applies if the "Standard Access" box has been checked in **Item 1** above or if no box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby acknowledges and agrees that, Subject to **Section 7** below, JoVE is and shall be the sole and exclusive owner of all rights of any nature, including, without limitation, all copyrights, in and to the Video. To the extent that, by law, the Author is deemed, now or at any time in the future, to have any rights of any nature in or to the Video, the Author hereby disclaims all such rights and transfers all such rights to JoVE.

6. **Grant of Rights in Video – Open Access.** This **Section 6** applies only if the "Open Access" box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby grants to JoVE, subject to **Section 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Video in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Video into other languages, create adaptations, summaries or extracts of the Video or other Derivative Works or Collective Works based on all or any portion of the Video and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. For any Video to which this **Section 6** is applicable, JoVE and the Author hereby grant to the public all such rights in the Video as provided in, but subject to all limitations and requirements set forth in, the CRC License.

7. **Government Employees.** If the Author is a United States government employee and the Article was prepared in the course of his or her duties as a United States government employee, as indicated in **Item 2** above, and any of the licenses or grants granted by the Author hereunder exceed the scope of the 17 U.S.C. 403, then the rights granted hereunder shall be limited to the maximum

rights permitted under such statute. In such case, all provisions contained herein that are not in conflict with such statute shall remain in full force and effect, and all provisions contained herein that do so conflict shall be deemed to be amended so as to provide to JoVE the maximum rights permissible within such statute.

8. **Protection of the Work.** The Author(s) authorize JoVE to take steps in the Author(s) name and on their behalf if JoVE believes some third party could be infringing or might infringe the copyright of either the Author's Article and/or Video.

9. **Likeness, Privacy, Personality.** The Author hereby grants JoVE the right to use the Author's name, voice, likeness, picture, photograph, image, biography and performance in any way, commercial or otherwise, in connection with the Materials and the sale, promotion and distribution thereof. The Author hereby waives any and all rights he or she may have, relating to his or her appearance in the Video or otherwise relating to the Materials, under all applicable privacy, likeness, personality or similar laws.

10. **Author Warranties.** The Author represents and warrants that the Article is original, that it has not been published, that the copyright interest is owned by the Author (or, if more than one author is listed at the beginning of this Agreement, by such authors collectively) and has not been assigned, licensed, or otherwise transferred to any other party. The Author represents and warrants that the author(s) listed at the top of this Agreement are the only authors of the Materials. If more than one author is listed at the top of this Agreement and if any such author has not entered into a separate Article and Video License Agreement with JoVE relating to the Materials, the Author represents and warrants that the Author has been authorized by each of the other such authors to execute this Agreement on his or her behalf and to bind him or her with respect to the terms of this Agreement as if each of them had been a party hereto as an Author. The Author warrants that the use, reproduction, distribution, public or private performance or display, and/or modification of all or any portion of the Materials does not and will not violate, infringe and/or misappropriate the patent, trademark, intellectual property or other rights of any third party. The Author represents and warrants that it has and will continue to comply with all government, institutional and other regulations, including, without limitation all institutional, laboratory, hospital, ethical, human and animal treatment, privacy, and all other rules, regulations, laws, procedures or guidelines, applicable to the Materials, and that all research involving human and animal subjects has been approved by the Author's relevant institutional review board.

11. **JoVE Discretion.** If the Author requests the assistance of JoVE in producing the Video in the Author's facility, the Author shall ensure that the presence of JoVE employees, agents or independent contractors is in accordance with the relevant regulations of the Author's institution. If more than one author is listed at the beginning of this Agreement, JoVE may, in its sole

ARTICLE AND VIDEO LICENSE AGREEMENT

discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JoVE reserves the right, in its absolute and sole discretion and without giving any reason therefore, to accept or decline any work submitted to JoVE. JoVE and its employees, agents and independent contractors shall have full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JoVE has sole discretion as to the method of making and publishing the Materials, including, without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

12. **Indemnification.** The Author agrees to indemnify JoVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JoVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JoVE, making of videos by JoVE, or publication in JoVE or elsewhere by JoVE. The Author shall be responsible for, and shall hold JoVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to


the making of a video by JoVE its employees, agents or independent contractors. All sterilization, cleanliness or decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's expense. All indemnifications provided herein shall include JoVE's attorney's fees and costs related to said losses or damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JoVE, its employees, agents or independent contractors.

13. **Fees.** To cover the cost incurred for publication, JoVE must receive payment before production and publication of the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JoVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.

14. **Transfer, Governing Law.** This Agreement may be assigned by JoVE and shall inure to the benefits of any of JoVE's successors and assignees. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to be one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement is required per submission.

CORRESPONDING AUTHOR

Name:	Barbara Rothen-Rutishauser		
Department:	Adolphe Merkle Institute		
Institution:	University Fribourg		
Title:	Prof.		
Signature:		Date:	12/04/2019

Please submit a **signed** and **dated** copy of this license by one of the following three methods:

1. Upload an electronic version on the JoVE submission site
2. Fax the document to +1.866.381.2236
3. Mail the document to JoVE / Attn: JoVE Editorial / 1 Alewife Center #200 / Cambridge, MA 02140

612542.6 For questions, please contact us at submissions@jove.com or +1.617.945.9051.

Signature Certificate

Document Ref.: NKNUZ-EZNFU-ELH6P-NS3BY

Document signed by:

	<p>Barbara Rothen-Rutishauser Verified E-mail: barbara.rothen@unifr.ch</p>	 
<p>IP: 77.58.160.171 Date: 04 Dec 2019 16:02:39 UTC</p>		

Document completed by all parties on:
04 Dec 2019 16:02:39 UTC

Page 1 of 1



Signed with PandaDoc.com

PandaDoc is the document platform that boosts your company's revenue by accelerating the way it transacts.

