

Table 1. Precision of the method in extraction and quantitation of HA and HFA from liquid commercial samples^a.

Humic substances, %							
Material	L16		L17		L2		
	HFA	HA	HFA	HA	HFA	HA	
Rep 1	1.44	17.00	6.59	7.76	0.36		4.46
Rep 2	1.39	16.03	6.25	7.79	0.42		4.93
Rep 3	1.34	16.44	6.02	7.55	0.40		4.46
Rep 4	1.54	16.75	6.20	7.69	0.33		4.53
Mean	1.43	16.56	6.27	7.70	0.38		4.60
SD	0.09	0.42	0. 24	0.11	0.04		0.23
RSD, %	6.29	2.53	3.80	1.39	10.4		4.91
Hor Rat(r)	1.58	0.72	1.25	0.47	2.31		1.55

^aExtraction conditions were 1 g in 1 L 0.1 M NaOH.