Fig3

-2 weeks	*	Culture Glia cells (change glia culture medium every 3-4days)
-4 days	*	Split and seed 10 ⁵ glia cells on 60mm culture dish
-3 days	*	Change glia cell dishes with neuronal culture medium Acid/basic treatment of coverslips
-2 days		Coat coverslips with poly-L-lysine
-1 days		Setup plating dishes with coated coverslips
Day 0	***	Dissecting hippocampi from PND0 AG22/23-f/f mice (seed 200K/60mm plating dish)
Day 3	茶木	Transfecting
Day 7-14	** 4	Fixing and staining of neurons
		Imaging and data analysis
`	▼	

