

Journal of Visualized Experiments

Open-Source Miniature Fluorimeter to Monitor Real-Time Isothermal Nucleic Acid Amplification Reactions in Resource-Limited Settings

--Manuscript Draft--

Article Type:	Methods Article - JoVE Produced Video
Manuscript Number:	JoVE62148R1
Full Title:	Open-Source Miniature Fluorimeter to Monitor Real-Time Isothermal Nucleic Acid Amplification Reactions in Resource-Limited Settings
Corresponding Author:	Rebecca Richards-Kortum UNITED STATES
Corresponding Author's Institution:	
Corresponding Author E-Mail:	rkortum@rice.edu
Order of Authors:	Jackson Coole Alex Kortum Yubo Tang Imran Vohra Yajur Maker Kathryn Kundrod Mary Natoli Rebecca Richards-Kortum
Additional Information:	
Question	Response
Please specify the section of the submitted manuscript.	Bioengineering
Please indicate whether this article will be Standard Access or Open Access.	Open Access (US\$4,200)
Please indicate the city, state/province, and country where this article will be filmed . Please do not use abbreviations.	Houston, TX, USA
Please confirm that you have read and agree to the terms and conditions of the author license agreement that applies below:	I agree to the Author License Agreement
Please provide any comments to the journal here.	It was unclear how to designate Alex Kortum as a co-first author in the "Manuscript Data" section of the submission, but it is designated in the actual manuscript itself.

TITLE:

Open-Source Miniature Fluorimeter to Monitor Real-Time Isothermal Nucleic Acid Amplification Reactions in Resource-Limited Settings

AUTHORS AND AFFILIATIONS:

Jackson Coole^{1*}, Alex Kortum^{1*}, Yubo Tang¹, Imran Vohra¹, Yajur Maker¹, Kathryn Kundrod¹, Mary Natoli¹, Rebecca Richards-Kortum^{1**}

¹Department of Bioengineering, Rice University, Houston, Texas, USA

*Authors contributed equally to the following work

** Corresponding author

jbc7@rice.edu

ask5@rice.edu

yt9@rice.edu

imran.vohra@rice.edu

ynm2@rice.edu

kk46@rice.edu

men2@rice.edu

**Corresponding Author

rkortum@rice.edu

KEYWORDS:

open-source, fluorimeter, isothermal amplification, nucleic acid amplification test, low-resource setting, SARS-CoV-2 detection

SUMMARY:

Detailed instructions are provided to build an open-source, modular fluorimeter that is compatible with many low-cost heaters to perform real-time, quantitative isothermal nucleic acid amplification.

ABSTRACT:

Traditional methods to detect and quantify nucleic acids rely on polymerase chain reaction (PCR) and require the use of expensive thermocyclers with integrated fluorescence detection of amplicons. Isothermal nucleic acid amplification technologies eliminate the need for thermal cycling; however, fluorescence-based detection of products is still required for real-time, quantitative results. Several portable isothermal heaters with integrated fluorescence detection are now commercially available; however, the cost of these devices remains a significant barrier to widespread adoption in resource-limited settings. Described here is a protocol for the design and assembly of a modular, low-cost fluorimeter constructed from off-the-shelf components. Enclosed in a compact 3D printed housing, the fluorimeter is designed to be placed atop a commercially available heat block holding a PCR tube. The fluorimeter described here was optimized to detect fluorescein isothiocyanate (FITC) dye, but the system can be modified for use

with dyes commonly used as reporters in real-time nucleic acid amplification reactions. Clinical applicability of the system is demonstrated by performing real-time nucleic acid detection with two isothermal amplification technologies: recombinase polymerase amplification (RPA) for detection of positive control DNA provided in a commercial kit and reverse transcription loop-mediated isothermal amplification (RT-LAMP) for detection of clinically meaningful levels of SARS-CoV-2 RNA.

INTRODUCTION:

Isothermal amplification technologies are widely used for detection of nucleic acids. Compared to traditional PCR approaches that require thermocycling, isothermal amplification allows nucleic acid amplification to occur at a single temperature, thus enabling faster time-to-results and better tolerance of inhibitors^{1,2}. Another key benefit of isothermal amplification is the reduced instrumentation complexity. Most isothermal amplification reactions only require a heat block and a detection modality—either real-time detection via fluorescence monitoring or endpoint detection, for example by lateral flow or gel electrophoresis^{3,4}. Real-time fluorescence detection is accomplished through the detection of fluorescence produced by intercalating dyes that activate in the presence of double-stranded DNA or quenched fluorescent probes that activate in the presence of specific double-stranded DNA sequences.

While commercially available benchtop isothermal fluorimeters exist, many lack customization for assay implementation. For example, many devices require specific or company-provided consumables, recommend preferred vendors, or utilize proprietary software to obtain advertised results. Most of these systems cost over \$5,000 USD, representing a significant barrier for widespread use in resource-limited settings. In addition, users in low-resource settings face challenges to maintain equipment designed for high-resource settings due to harsh environmental conditions, weak supply chains for spare parts, and specialized tools required for maintenance and repair⁵. To meet this need, described here is the design and assembly of a modular and low-cost fluorimeter constructed from off-the-shelf components enclosed in a compact 3D printed housing (**Figures 1A–C**) with two optional configurations. The first configuration of this device uses commercially available glass filters and a dichroic mirror to block excess background light and has a total cost of assembly of \$830 USD. While these filters are commonly used in fluorescence-based imaging systems, replacing expensive high-grade optical filter foils has previously been shown to allow for nucleic acid detection⁶. The second configuration of the fluorimeter incorporates these inexpensive filters and replaces the dichroic mirrors with $\phi\frac{1}{2}$ " beam splitters, reducing the overall cost of the system from \$830 to \$450 USD.

Representative images of the assembly are shown for the first configuration in **Figure 1** and **Figure 2**, but analogous images for the second configuration can be found in **Supplemental File 6**. To avoid the need for specialized optical alignment, the optical system has designated areas to place each optical component and can be made with a relatively low-end 3D printer, allowing for widespread use of the design. The only differences in construction and assembly for the two configurations are the files used for 3D printing and the optical components placed into the enclosure. The external dimensions of the 3D printed enclosure for both systems are the same. A cost comparison of the two systems is shown in **Table 1**.

As shown in **Figure 1A**, to maintain a small form factor, the fluorimeter consists of $\phi\frac{1}{2}$ " (~12.5 mm) optics, coupled with compact illumination and detection that are placed to measure signal through the top of the PCR tube. The system in **Figure 1** is designed to detect dyes with peak excitation and emission wavelengths near 490 nm and 525 nm, respectively, including FITC and closely related dyes such as SYBR and SYTO-9, which are commonly used as reporters in real-time nucleic acid amplification reactions^{7,8}. The excitation source, optical filters, and detector can readily be substituted for components compatible with different fluorescent dyes as desired. Nucleic acid amplification reactions are typically performed in PCR tubes, and the fluorimeter is designed to be placed atop any commercially available heat block that holds PCR tubes (**Figure 1D**) allowing for real time monitoring of isothermal reactions. Appropriate heat blocks are available in most biomedical laboratories and can be purchased for less than \$500 USD.

The use of single-board computers to provide a low-cost, point of care alternative for controlling imaging technologies has been previously demonstrated⁹. Building off that work, in this protocol a single-board computer-powered graphical user interface (**Figure 1D**) is used to facilitate real-time data logging and display of results at the point of care, eliminating the need for a laptop computer to process or visualize data. Fluorescence measurements were transferred through I²C protocol from the light sensors to a microcontroller, and then made available to the single-board computer through serial communication. Electrical connections for illumination and data transfer were provided through simplified wiring and soldering on miniaturized breadboards, negating the need for specialized printed circuit boards (PCBs). The software required to run the fluorimeter is available through open-source software frameworks and the code required to run the device is provided in the **supplemental coding files**. The complete fluorimeter can be assembled for between \$430 to \$830 USD, and results show that it provides accurate and reliable fluorescence measurements to monitor real-time isothermal amplification of nucleic acids.

PROTOCOL:

1. Preparation steps: 3D printing and soldering

NOTE: The optical system described in this protocol is designed for a standard dry block heater.

1.1. To create the first configuration, 3D print the CAD files provided as **Supplemental Files 1, 2, and 3** respectively:

1.2. To create the second configuration, 3D print the CAD files provided as **Supplemental Files 3, 4, and 5**, respectively:

NOTE: These parts are designed to be printed with supports. In this guide, a black polycarbonate filament is used that can maintain its form after being subjected to temperatures up to 110 °C. In general, any material that can be heated to the temperature of the desired isothermal reaction without significant deformation can be used. To minimize the effect of internal reflections and ambient light interference, a material that is black or another dark color is recommended.

1.3. Prepare the two light-to-digital sensor evaluation modules for parallel monitoring of two samples. On one of the sensor test boards, remove the R4 resistor, and solder a jumper wire from the right pad of the R4 area on the PCB to the top pad in the R1 area on the PCB. This will change the I²C address of the sensor, thus allowing simultaneous measurement of both sensors.

NOTE: The sensor used consists of two PCBs: a USB adapter board and a sensor test board containing the light sensor; only the sensor test board is needed for this device.

1.4. Solder wires to each of the two Light Emitting Diodes (LEDs). Connect a red wire (positive) to the pad labeled “1” on the LED and a black wire (negative) on the pad labeled “2” on the LED. Apply a thin layer of thermal adhesive on the back of the LED, place the LED on the top of an end cap, and wait until the thermal adhesive cures. On the other side of the end cap, add a heat sink.

NOTE: When testing the LEDs before they are sealed in the enclosure, ensure to wear a proper blue light blocking eye protection.

1.5. To create the second configuration, cut two 1/4-inch diameter circles from a blue excitation foil sheet and four 1/4-inch diameter circles from a yellow emission foil sheet with scissors or a razor blade.

1.6. Press an M2.5 hex shaped insert into each of the four holes on the slanted portion of the ‘LCD_Screen_Holder.stl’ part.

2. Optical assembly

2.1. Place a 3/16-inch-long 4-40 threaded insert into the hole on the top of the ‘Optics_Enclosure_Bottom.stl’ part. Place a 1/4-inch-long 4-40 threaded insert into all other holes of the 3D printed part as shown in **Figure 2A**.

2.2. Insert the sensor test board into the top cavity of the housing, with the five pins facing towards the top and closest to the center axis of the device. Secure with a 3/16-inch-long 4-40 screw through the hole on the sensor test board (**Figure 2B**).

2.3. Place one of the 20 mm focal length lenses in the section below the sensor test board with the convex side facing towards the bottom of the device and away from the test board (**Figure 2C**).

2.4. To create the first configuration, place the long pass filter into the next section below the 20 mm focal length lens (**Figure 2D**) placed in the previous step. To create the second configuration, place two yellow emission filter foils into the section below the lens.

2.5. To create the first configuration, place the dichroic mirror into the diagonal section near the center of the encasement while observing the filter orientation specified by the manufacturer

(**Figure 2E**). To create the second configuration, place the beam splitter into the diagonal section. No specific orientation is needed for the beam splitter.

2.6. Place a second 20 mm focal length lens into the section below the dichroic mirror (or beam splitter, depending on configuration) with the convex side pointing towards the top of the device (**Figure 2F**).

2.7. To create the first configuration, place the excitation filter in the section to the right of the dichroic mirror, making sure the arrow points towards the dichroic mirror (**Figure 2G**). To create the second configuration, place one blue excitation filter foil into the section to the right of the beam splitter.

2.8. Place the 15 mm focal length lens to the right of the excitation filter with the convex side facing the dichroic mirror (**Figure 2H**).

2.9. Place an LED into the remaining section of the print, with the LED facing towards the dichroic mirror (or beam splitter, depending on configuration). Ensure the two wires leading from the LED are inserted into the recessed channels so that the print will close tightly.

2.10. Repeat steps 2.3–2.9 for the other side of the 3D printed part (**Figure 2I**).

2.11. Close the empty side of the print on top of the print with the optical components by placing the extruded portions of the top half of the encasement into the recessed grooves of the bottom half of the encasement. Secure the two printed parts together with 3/8-inch-long 4-40 screws (**Figure 2J**).

NOTE: If the two printed parts are not completely closed, stray excitation light can escape from the optical housing. Ensure proper blue light blocking eye protection is worn until a proper seal is achieved. Reseal the enclosure until no excess light escapes.

3. Electronics and touchscreen assembly

3.1. Connect the two mini breadboards together, and then place the microcontroller onto one of the breadboards. Ensure that the microUSB port of the microcontroller faces outwards.

3.2. To connect LED modulation, connect the CTL pin of the LED (+) driver to a digital pin of the microcontroller and the LED (-) pin of the LED driver to a GND pin of the microcontroller.

3.3. Remove the plastic covers on the back of the breadboards. Press the adhesive backing of the breadboards to the 3D printed part to attach the combined breadboards to the inside of the back portion of the 'LCD_Screen_Holder.stl' printed part.

3.4. Secure the liquid crystal display (LCD) screen holder with the assembled breadboards inside to the optical enclosure assembled in section 2 with one-inch-long 4-40 screws.

3.5. To connect the LED power supply, connect the LED (+) pin of the LED driver to the positive wire of the first LED. Connect the negative wire of the first LED to the positive wire of the second LED on the breadboard. Connect the negative wire of the second LED to the LED (-) pin of the LED driver.

NOTE: The order of first or second LED is arbitrary.

3.6. To connect the LED driver power supply, connect the positive and negative wires of the 10 V power supply to the VIN+ and VIN- pins of the LED driver, respectively. (A barrel jack to two-pin adapter was used.)

3.7. Connect the sensor test board power supply and data transfer. Only four pins on the sensor test board are used: the SCK, SDA, VDUT, and GND. Take a 4-pin female to male jumper wire and connect those pins on the light-to-digital sensor test boards to the mini breadboard through the gap in the top right of the LCD Holder print.

3.8. On the breadboard, ensure connections between the following are in place: the 3.3 V pin of the microcontroller and the VDUT pin of both test boards; the GND pin of the microcontroller and the GND pin of both test boards; the analog 4 (A4) pin of the microcontroller and the SDA pin of both test boards; and the analog 5 (A5) pin of the microcontroller and the SCK pin of both test boards.

NOTE: Because I²C communication is used for the light sensors, the SCK and SDA pins of both sensors can both be routed to the same pins of the microcontroller.

3.9. Secure the single-board computer onto the LCD Screen holder with four M2.5 screws. Ensure that the HDMI and power adapter ports of the single-board computer face upwards and the single-board computer is centered on the 3D printed part.

3.10. Connect the touchscreen display to the single-board computer as per the touchscreen instructions, and then connect the HDMI port of the single-board computer to the touchscreen's HDMI port.

4. Software installation

4.1. Install and use the web editor to upload the custom sketch "MiniFluorimeter_2Diode.ino" provided in **Supplemental Coding File 1** onto the microcontroller. Ensure the "ClosedCube OPT3002" library is installed using the Library Manager.

4.2. Change the variable led_A_pin to the number of the digital pin used in step 3.3 (Electronics and touchscreen assembly section).

4.3. Adjust the number of milliseconds the LED is turned on, change the value of the variable

ExposureTime to acquire fluorescence measurements. To adjust the number of milliseconds between LED exposures by change the value of the variable led_A_Interval.

4.4. Change the variable led_Power to a number between zero and one to adjust the brightness of the LEDs during exposures. Zero gives the maximum amount of brightness and one gives the lowest amount of brightness.

4.5. Turn on the ability to control the display via touchscreen by following the manufacturer's instructions provided with the 3.5-inch display.

NOTE: If desired, the 3.5-inch screen can be used as a monitor without touchscreen capabilities, and a keyboard and mouse can be attached to the USB ports of the single-board computer for control of the single-board computer.

4.6. Download the "MiniFluorimeter_2Diode_GUI.py" file from **Supplemental Coding File 2** to a desired location on the single-board computer.

4.7. Ensure a working version of Python is installed on the single-board computer . Python 3.7 was used in the provided Python module, but any stable Python version could be used with appropriate changes to the provided script. Install the libraries needed for the Python program onto the single-board computer.

4.8. Change the variable measurement_time to the length of time desired to take measurements. The program ends acquisition and closes after the desired time is elapsed. The GUI also allows to end acquisition through a button on the user interface.

4.9. Change the variable serialPort to the serial address of the connected microcontroller.

5. Recording real-time fluorescence data

5.1. Turn on the commercial heat block and allow it to reach the desired temperature.

5.2. Power the single-board computer with a standard 5 V power supply provided with most single-board computer purchases. Connect the single-board computer to the microcontroller with a microUSB to USB cable.

5.3. Using the touchscreen, open the provided Python script. Change the measurement_time and serialPort variables to the desired values. Change the variable outputFilepath to the name of the data file the program generates. Ensure the filename ends in '.xlsx'.

5.4. Place two PCR tubes containing the reactions to be monitored into the heat block. Ensure that the placement of the PCR tubes aligns with the optical channels of the fluorimeter once it is placed onto the heat block.

5.5. Place the fluorimeter atop the heat block with the PCR tubes centered between the four pegs extruding from each optical channel of the fluorimeter. For optimal measurements, ensure that the 3D printed fluorimeter is attached securely into the wells of the heat block.

5.6. Securely attach the fluorimeter, plug in the power supply adapter for the LEDs.

5.7. Use the touchscreen to start the Python program. A graphical user interface (GUI) appears on the LCD screen and measures the real-time fluorescence.

5.8. Observe the real-time fluorescence measurements over time for both PCR tubes that are shown to the user on the GUI.

5.9. After the user-determined experiment time signal, acquisition ceases. View the measurements in the output data file saved in the user-defined location. To end the measurements early, click the button labeled “Stop Acquisition” on the user interface.

REPRESENTATIVE RESULTS:

Once assembled, fluorimeter performance can be validated by measuring fluorescence from a dilution series of FITC dye. In **Figure 3A**, measurements of FITC dye at concentrations of 0, 20, 40, 60, and 80 pg/ μ L prepared in 1x PBS on both channels of the first configuration of the fluorimeter are shown. Each sample was measured three times with an LED exposure of 1.5 s at 20 s intervals. Both channels of the fluorimeter show a linear response across the desired range.

Clinical applicability of the fluorimeter was further demonstrated by using the system together with a commercially available dry heat block to perform amplification with two isothermal amplification technologies: RPA and RT-LAMP.

Figure 3B demonstrates the baseline-subtracted time course of fluorescence measured during 39 °C amplification of 50 μ L of real-time RPA positive and negative control reactions for kit positive control DNA provided in a standard commercial kit and prepared per manufacturer’s instructions. RPA reactions, which produce a relatively low level of fluorescence, were measured using the first configuration of the fluorimeter that achieves better suppression of excitation light.

Figure 3C demonstrates the time course measurement of a custom RT-LAMP assay at 65 °C utilizing the N2, E1, and As1e primer sets described by Zhang et al¹⁰, and Rabe and Cepko¹¹. RT-LAMP reactions produce a greater amount of fluorescence and were measured using the second, lower-cost fluorimeter configuration. Oligonucleotides were purchased and resuspended in 2x TE buffer at a 1 mM concentration. Forward inner primer (FIP) and backward inner primer (BIP) oligos were ordered with high performance liquid chromatography purification. Each primer set (N2, E1, and As1e) was combined to make 1000 μ L of a 25x mix as follows: 40 μ L of FIP, 40 μ L of BIP, 5 μ L of F3, 5 μ L of B3, 10 μ L of LF, 10 μ L of LB, and 890 μ L of 1x TE buffer. To assemble each RT-LAMP reaction, 1 μ L of each primer set was added to 0.5 μ L of 50x fluorescent dye and 12.5 μ L of 2x master mix and the reaction volume was brought up to 20 μ L of with nuclease-free water

per manufacturer's instructions. The SARS-CoV-2 RNA Control was serially diluted in nuclease-free water to concentrations of 10, 100, or 1,000 copies per μL , and 5 μL were added for a total reaction volume of 25 μL . The no target control (NTC) used in all experiments was nuclease-free water. RT-LAMP reactions were overlaid with 25 μL of molecular biology-grade mineral oil.

RPA and RT-LAMP reactions were assembled in two wells of a 0.2 mL of low-profile 8-tube PCR strip and capped with ultraclear flat caps. Each RPA and RT-LAMP reaction was run in triplicate. In all tests, the mini-fluorimeter successfully quantified the temporal increase of fluorescence levels associated with DNA amplification.

FIGURE AND TABLE LEGENDS:

Figure 1: Optical housing and assembled miniature fluorimeter atop heat block. (A) Diagram of optical housing showing optical components placed in one detection channel. (B) Diagram of the first configuration of the miniature fluorimeter following assembly. (C) Photograph of the optical housing with optical components placed in one detection channel. (D) Photograph of assembled miniature fluorimeter placed on top of a commercially available heat block.

Figure 2: Assembly and electrical control diagram of the miniature fluorimeter. A-J) Step-by-step placement of the optical components in the 3D printed optical housing for the first system configuration. (K) Electrical diagram of the miniature fluorimeter for both configurations.

Figure 3: Representative measurements obtained with the miniature fluorimeter. (A) Measured fluorescence vs. FITC dye concentration in both channels shows linear response across the desired dynamic range. (B) Real-time fluorescence vs time for isothermal amplification of positive and negative controls of a commercially available kit. Amplification occurs as expected for the positive control. (C) Real-time fluorescence vs time for isothermal amplification of 50 copies of SARS-CoV-2 RNA and an NTC sample from a custom RT-LAMP assay. Amplification occurs as expected near the limit of detection of the assay.

Table 1: Cost comparison of the two configurations of the miniature fluorimeter.

Supplemental File 1: System1_Optics_Enclosure_Top.stl

Supplemental File 2: System1_Optics_Enclosure_Bottom.stl, and

Supplemental File 3: LCD_Screen_Holder.stl

Supplemental File 4: System2_Optics_Enclosure_Top.stl

Supplemental File 5: System2_Optics_Enclosure_Bottom.stl, and

Supplemental File 6: System2_BuildInstructions.pdf

Supplemental Coding File 1: MiniFluorimeter_2Diode.ino

Supplemental Coding File 2: MiniFluorimeter_2Diode_GUI.py

DISCUSSION:

Described here is an open-source, low-cost, modular, portable fluorimeter for quantitative fluorescence detection of isothermal amplification reactions. Open-source projects facilitate rapid and inexpensive maintenance with readily available replacement parts and allow users the flexibility to adapt the system to their needs based on modular design. This protocol describes the process of assembling mechanical, optical and electrical components, and validating optical performance. Furthermore, the flexibility of the fluorimeter to monitor two different types of isothermal amplification assays with significantly different temperature, volume, and fluorescence requirements, RPA exo and RT-LAMP, was demonstrated. RPA is performed at 39 °C in 50 µL reactions that utilize a sequence-specific FAM-tagged probe for fluorescence generation, while RT-LAMP is performed at 65 °C in a 25 µL reaction volume and uses an intercalating dye to report the presence of the amplified DNA. Because fluorescence measurements are made through the top of PCR tubes with flat caps, the fluorimeter is able to detect fluorescence from both assay volumes, and heat requirements are only limited by the commercial heat block chosen. Furthermore, the fluorescence intensity produced in RT-LAMP is almost in the order of magnitude greater than that produced in RPA, due to the dye- versus probe-based methods of fluorescence signal generation. However, the dynamic range of the optical sensor chosen can detect and quantify both the signals, and baseline subtraction algorithms account for these differences to produce reliable fluorescence readings.

To facilitate technology dissemination and minimize potential maintenance cost, a modular design that is compatible with heaters that are widely available in different settings was employed. In the current protocol, a common dry block heater was used; the same optical and electrical design can be readily adapted for other commercially available heaters. If another dry block heater is to be used, minimal changes in the 3D housing design will be required. Specifically, the bottom pegs of the optical enclosure STL files must be modified to ensure proper alignment with the wells of other commercial heat blocks. While the enclosures shown in the examples were printed on a relatively low-end 3D printer (see **Table of Materials**), care should be taken to ensure the printer resolution and/or print tolerances are adequate to accommodate the optical components and threaded inserts. In the STL files provided, a tolerance of 0.01–0.02 inches was added on either side of the optical components in the radial and axial directions based on the dimensions specified by the manufacturer. This ensures all optical components fit securely within the print and that the enclosure completely blocks excess light from entering or escaping. To ensure a proper press fit for the threaded inserts, a similar tolerance of 0.01–0.02 inches was subtracted from the manufacturer-provided diameter in the CAD file.

RPA reactions were successfully monitored using the first fluorimeter configuration, while RT-LAMP reactions could be monitored using either configuration. The improved stray light rejection of the first configuration was necessary to monitor the low levels of fluorescence produced by the fluorogenic probe in RPA reactions. In contrast, RT-LAMP utilizes an intercalating dye for signal generation, resulting in a higher fluorescence intensity that is compatible with the lower

dynamic range of the second configuration using photographic filter foils. Users should select the fluorimeter configuration that matches the fluorescence signal generating element—intercalating dye or fluorogenic probe—of their assay.

One limitation of this system is that heating is provided by a commercially available heat block powered through a standard wall outlet. This system could be further developed for use in areas lacking reliable access to electricity by incorporating portable and rechargeable battery packs as shown by other groups¹². Another limitation is the relatively low throughput of the system, which allows for simultaneous fluorescence measurement of only two samples at a time. Multiple prints of the enclosure can be placed atop the same heat block to increase throughput; however, the light sensor used only has four unique I2C addresses. This restricts the maximum number of samples that can be simultaneously measured to four. A different light sensor with a larger number of unique I²C addresses is needed to further increase throughput.

ACKNOWLEDGMENTS:

Special thanks to Chelsey Smith, Megan Chang, Emilie Newsham, Sai Paul, and Christopher Goh for their assistance with sample preparation. The authors thank Caroline Noxon for manuscript revision. Funding for this work was provided from the American people by USAID through an IAVI research grant CCID 9204 under award AID-OAA-A16-00032 between IAVI and USAID.

DISCLOSURES:

The authors declare no conflicts of interest.

REFERENCES:

1. Daher, R.K., Stewart, G., Boissinot, M., Bergeron, M.G. Recombinase polymerase amplification for diagnostic applications. *Clinical Chemistry*. **62** (7), 947–958 (2016).
2. Yan, L. et al. Isothermal Amplified Detection of DNA and RNA. *Molecular BioSystems*. **10** (5), 970–1003 (2014).
3. Giuffrida, M.C., Spoto, G. Integration of Isothermal Amplification Methods in Microfluidic Devices: Recent Advances. *Biosensors and Bioelectronics*. **90**, 174–186 (2017).
4. Gill, P., Ghaemi, A. Nucleic Acid Isothermal Amplification Technologies—A Review. *Nucleosides, Nucleotides and Nucleic Acids*. **27** (3), 224–243 (2008).
5. Richards-Kortum, R., Oden, M. Devices for Low-Resource Health Care. *Science*. **342** (6162), 1055–1057 (2013).
6. Katzmeier, F. et al. A Low-Cost Fluorescence Reader for in vitro Transcription and Nucleic Acid Detection with Cas13a. *PLOS One*. **14** (12): e0220091 (2019).
7. Safavieh, M. et al. Emerging Loop-Mediated Isothermal Amplification-Based Microchip and Microdevice Technologies for Nucleic Acid Detection. *ACS Biomaterials Science and Engineering*. **2** (3), 278–294 (2016).
8. Monis, P.T., Giglio, S., Saint, C.P. Comparison of SYTO9 and SYBR Green I for Real-Time Polymerase Chain Reaction and Investigation of the Effect of Dye concentration on Amplification and DNA Melting Curve Analysis. *Analytical Biochemistry*. **340** (1), 24–34 (2005).
9. Parra, S. et al. Development of Low-Cost Point-of-Care Technologies for Cervical Cancer Prevention Based on a Single-Board Computer. *IEEE Journal of Translational Engineering in*

Health and Medicine. **8** (4300210) (2020).

10. Zhang, Y. et al. Enhancing Colorimetric Loop-Mediated Isothermal Amplification Speed and Sensitivity with Guanidine Chloride. *Biotechniques*. **69** (3), 178–185 (2020).

11. Rabe, B.A., Cepko, C. SARS-CoV-2 Detection Using an Isothermal Amplification Reaction and a Rapid, Inexpensive Protocol for Sample Inactivation and Purification. *Proceedings of the National Academy of Sciences of the United States of America*. **117** (39), 24450–24458 (2020).

12. Snodgrass, R. et al. A Portable Device for Nucleic Acid Quantification Powered by Sunlight, a Flame or Electricity. *Nature Biomedical Engineering*. **2** (9), 657–665 (2018).




Figure 2






Figure 3



	System 1		System 2	
Item	Quantity	Total Price (USD)	Quantity	Total Price (USD)
Optical Components				
Lenses	6	158.14	6	158.14
Mirrors	2	244.56	2	60
Optical Filters	4	200	6	5
	Subtotal	602.7	Subtotal	223.14
Illumination and Detection				
LEDs	2	72.62	2	72.62
LED Driver	1	11.49	1	11.49
Photodiode	2	50	2	50
	Subtotal	134.11	Subtotal	134.11
Electronics and Display				
Arduino Nano	1	20.7	1	20.7
Raspberry Pi	1	35	1	35
LCD Screen	1	25	1	25
Mini Breadboard	1	4	1	4
10V Power Supply	1	8.6	1	8.6
	Subtotal	93.3	Subtotal	93.3
Total Cost (USD)		830.11		450.55

Name of Material/Equipment	Company	Catalog Number
1/4-inch-long 4-40 threaded insert	McMaster-Carr	90742A116
10v power supply	GlobTek, Inc.	WR9HU1800CCP-F(R6B)
15 mm focal length lens	Thorlabs	LA 1074
1-inch-long 4-40 screws	McMaster-Carr	
20 mm focal length lens	Thorlabs	LA 1540
2x WarmStart LAMP Master Mix	New England Biolabs, Inc	E1700
3.5" Touch Screen	Uctronics	BO10601
3/16-inch-long 4-40 screw	McMaster-Carr	90128A105
3/16-inch-long 4-40 threaded insert	McMaster-Carr	90742A115
3/8-inch-long 4-40 screws	McMaster-Carr	90128A108
3D printer filament	3D Universe	UMNFC-PC285-BLACK
3D printer used	Ultimaker	Ultimaker 2+
8-tube PCR strips	BioRad	#TLS0801
Advanced Mini Dry Block Heater	VWR International	10153-320
barrel jack to two-pin adapter	SparkFun Electronics	1568-1238-ND
Blue Excitation Filter Foil	LEE	LE071S
Blue LED - 460 nm	Mouser	LZ1-30DB00-0100
Dichroic Mirror	Thorlabs	DMLP505T
Emission Filter	Edmunds Optics	OG-515
Excitation Filter	Omega Filters	490AESP
LED Driver	LEDdynamics	3021-D-I-700
M2.5 Hex Shaped insert	McMaster-Carr	91292A009
Microcontroller	Arduino	Nano
Mini Breadboard	Adafruit	65
Molecular biology-grade mineral oil	Sigma Aldrich	69794
OPT3002EVM - Light-to-Digital Sensor	Texas Instruments	OPT3002EVM:
Purchased oligonucleotides	Integrated DNA Technologies	
RPA kit positive control DNA	TwistDx Limited	CONTROL01DNAE

SARS-CoV-2 RNA Control	Twist Biosciences	MN908947.3
Single board computer	Raspberry Pi	Raspberry Pi 3
TwistAmp RPA exo kit	TwistDx Limited	TAEXO02KIT
Ultraclear flat caps	BioRad	#TCS0803
Yellow Emmission Filter Foil	LEE	LE767S

Comments/Description

Used to secure the two sides of the 3D printed optical enclosure together.

AC/DC Wall Mount Adapter 10V 18W

Two total are used for the fluorimeter. This lens is used to focus the LED illumination.

Four total are used for the fluorimeter.

Master mix was used to create the LAMP reactions shown in Figure 3C

Used to secure the OPT3002 test board onto the 3D printed enclosure

Used to secure the two sides of the 3D printed optical enclosure together.

Black or another dark color preferred

The following heat blocks are acceptable substitutes without the need for redesigning the optical assembly: 949VWMNLUS, 949VWMHLUS, and 949VWMHLEU

Selected for use with FITC - other fluorescent dyes may require different filters.

Selected for use with FITC - other fluorescent dyes may require different parts

Selected for use with FITC - other fluorescent dyes may require different parts

Selected for use with FITC - other fluorescent dyes may require different parts. The arrow on the part points away from the illumination source.

Selected for use with FITC - other fluorescent dyes may require different parts

Used to secure the Raspberry Pi to the 3D printed LCD Screen Holder

Light-to-digital sensor used. Consists of two PCBs: a SM-USB_DIG board and the OPT3002 test board; only the OPT3002 test board is needed for this device.

Selected for use with FITC - other fluorescent dyes may require different parts

Author Response

Article Title: "Open-Source Miniature Fluorimeter to Monitor Real-Time Isothermal Nucleic Acid Amplification Reactions in Resource-Limited Settings"

Editorial comments:

Changes to be made by the Author(s):

1. Please take this opportunity to thoroughly proofread the manuscript to ensure that there are no spelling or grammar issues. Please define all abbreviations at first use.

We carefully reviewed the manuscript to correct spelling and grammar issues. We defined all abbreviations at first use.

2. Please provide an email address for each author.

Email addresses are provided under the list of authors and affiliations.

3. JoVE cannot publish manuscripts containing commercial language. This includes trademark symbols (™), registered symbols (®), and company names before an instrument or reagent. Please remove all commercial language from your manuscript and use generic terms instead. All commercial products should be sufficiently referenced in the Table of Materials and Reagents. For example: 949VWMNLEU Advanced 100 Mini Dry Block Heater (VWR #10153-320), 949VWMNLUS, 102 949VWMHLUS, and 949VWMHLEU; Texas Instruments; OPT3002; LZ1-30DB00-0100 LED Engines (Mouser); Adafruit; TwistDx Limited; TwistAmp RPA exo kit; Integrated DNA Technologies (IDT), Coralville IA; WarmStart LAMP Master Mix (New England BioLabs Inc, E1700); Twist Biosciences, MN908947.3; BioRad, #XYZ; VWR Advanced Dry Block heater etc

We replaced commercial language with generic terms throughout the manuscript. All commercial products are now referenced in the Table of Materials and Reagents.

4. Please revise the text, especially in the protocol, to avoid the use of any personal pronouns (e.g., "we", "you", "our" etc.).

We removed personal pronouns from the text of the revised manuscript.

5. Please note that your protocol will be used to generate the script for the video and must contain everything that you would like shown in the video. Please add more details to your protocol steps. Please ensure you answer the "how" question, i.e., how is the step performed? Alternatively, add references to published material specifying how to perform the protocol action. Please add more specific details (e.g., button clicks for software actions, numerical values for settings, etc) to your protocol steps. There should be enough detail in each step to supplement the actions seen in the video so that viewers can easily replicate the protocol.

We expanded the protocol sections to more fully address how each step is performed. We also expanded the "Software Installation" and "Recording real-time fluorescence data" sections to provide more specific detail on how to perform the protocol.

6. Please format the manuscript as: paragraph Indentation: 0 for both left and right and special: none, Line spacings: single. Please include a single line space between each step, substep and note in the protocol section. Please use Calibri 12 points and one-inch margins on all the side. Please include a one line space between each protocol step and then highlight up to 3 pages of protocol text for inclusion in the protocol section of the video.

We formatted the manuscript as directed.

7. Please remove the embedded Table of Materials from the manuscript. All tables should be uploaded separately to your Editorial Manager account in the form of an .xls or .xlsx file.

The embedded Table of Materials has been removed from the manuscript and attached to the submission as a .xlsx file.

8. As we are a methods journal, please revise the Discussion to explicitly cover the following in detail in 3-6 paragraphs with citations:

- a) Critical steps within the protocol
- b) Any modifications and troubleshooting of the technique
- c) Any limitations of the technique

We expanded the Discussion section to explain critical and troubleshooting steps, including:

- 1) how to adapt the 3D printed enclosure to work with a variety of commercial heat blocks;
- 2) how and why to ensure the system is tightly sealed; and
- 3) how to further reduce the cost of the fluorimeter based on the type of isothermal assay being performed.

Finally, we added a section discussing the current limitations of our technique.

9. Please ensure that the references appear as the following: [Lastname, F.I., LastName, F.I., LastName, F.I. Article Title. Source. Volume (Issue), FirstPage–LastPage (YEAR).] For more than 6 authors, list only the first author then et al. Please include volume and issue numbers for all references, and do not abbreviate journal names.

We reformatted the references using the requested format.

Reviewer #1:**Manuscript Summary:**

A nice piece of work. Clearly presented and reproducible build instruction are included. Needs a little more data to prove the system works.

Major Concerns:

Only one reaction for RPA and one reaction for LAMP is presented. This is a nice system and additional evidence to show efficacy and reproducibility would make it valuable to the reader. To assess the validity, utility and reproducibility of the system described, more representative results need to be presented. At a minimum 3 repeats of the positive and negative controls should be presented, ideally a serial dilution of the template would offer the ability to assess the performance of the system with different template starting quantities.

We thank the reviewer for their feedback, which helped us improve our manuscript. We added results to show data from a clinically relevant range of input copy numbers for RT-LAMP in triplicate. For RPA, we added results to show the positive and negative controls in triplicate; it was not possible to provide data with a range of input copy numbers due to the low concentration at which the RPA exo kit's positive control DNA is supplied.

In addition, we further reduced the cost of the system (from \$830 to \$450) by developing a configuration that uses lower-cost optical filters. We show that the lower-cost system can be used to monitor fluorescence generated in LAMP reactions, which use an intercalating dye to generate signal. RPA reactions, which rely on a fluorogenic probe, are best monitored using the original system.

Minor Concerns:

Figure 3a lists an absolute fluorescence intensity on the y axis as an irradiance quantity but I can not see what this irradiance quantity refers to. Is this the irradiance on the sensor surface? Can you please clarify where these absolute measures are made?

Our initial submission referred to the fluorescence intensity on the sensor surface (which is calculated using a conversion factor supplied by the manufacturer). For generality, we revised the manuscript to provide fluorescence in relative fluorescence units (RFU).

Reviewer #2:**Manuscript Summary:**

The paper describes design and assembly of a modular and low-cost fluorimeter for isothermal diagnostics. The paper demonstrates clinical applicability of the method by performing real-time nucleic acid detection with two different isothermal amplification technologies. These methods are needed for simple diagnostics.

Minor Concerns:

Why fluorescence of the positive control, shown on Fig. 3B, shows significantly lower fluorescence level than 50 copies (Fig. 3C)?

The difference in fluorescence intensity arises because RPA reactions (Fig 3B of initial submission) use a fluorogenic probe and produce less fluorescence than RT-LAMP reactions (Figure 3C) which utilize an intercalating fluorescent dye.

Fig 3C shows signal amplification from 50 copies. Would be useful to demonstrate signal amplification from other concentrations also, i.e., 100, 500 copies.

We thank the reviewer for this feedback, which helped us improve our manuscript. As described above, we added data (measured in triplicate) describing results across a clinically relevant range of copy numbers for RT-LAMP.

Reviewer #3:

Manuscript Summary:

This is a detailed description of home made real-time PCR system made by 3d printing and off the shelf components.

Major Concerns:

none

We thank the reviewer for their feedback.



October 2nd, 2020

Editor, *Journal of Visualized Experiments*
1 Alewife Center, Suite 200
Cambridge, MA 02140 USA
Phone: 617-945-9051
Email: submissions@jove.com

Dear Editor,

We are pleased to submit our manuscript titled “Open-Source Miniature Fluorimeter to Monitor Real-Time Isothermal Nucleic Acid Amplification Reactions in Resource-Limited Settings” for consideration to be published in *Journal of Visualized Experiments*. This manuscript provides detailed instructions describing how to build an open-source, modular fluorimeter that is compatible with many low-cost heaters to perform real-time, quantitative isothermal nucleic acid amplification. We believe many of your readers will find the design of our system insightful, especially researchers in the field of point of care diagnostics. We would like to request accelerated review due to the role the described technology could play in COVID-19 response across the world. We are aware that video production in affected areas is postponed and can film the video ourselves. Should the article be of interest, we would welcome a discussion of the best way to proceed with filming.

Please let me know if you require any additional information. We appreciate your consideration of our manuscript.

Sincerely,

A handwritten signature in black ink, reading "Rebecca RK", followed by a long horizontal line.

Rebecca Richards-Kortum, PhD
Malcolm Gillis University Professor
Department of Bioengineering
Rice University
Houston, Texas, USA
Phone: 713-348-3823
Email: rkortum@rice.edu



1 Alewife Center #200
Cambridge, MA 02140
tel. 617.945.9051
www.jove.com

ARTICLE AND VIDEO LICENSE AGREEMENT

Title of Article:	Open-Source Miniature Fluorimeter to Monitor Real-Time Isothermal Nucleic Acid Amplification Reactions in Resource-Limited Settings
Author(s):	Jackson Coole, Alex Kortum, Yubo Tang, Imran Vohra, Kathryn Kundrod, Mary Natoli, and Rebecca Richards-Kortum

Item 1: The Author elects to have the Materials be made available (as described at <http://www.jove.com/publish>) via:

☐

Standard Access

☒

Open Access

Item 2: Please select one of the following items:

☒

The Author is **NOT** a United States government employee.

☐

The Author is a United States government employee and the Materials were prepared in the course of his or her duties as a United States government employee.

☐

The Author is a United States government employee but the Materials were NOT prepared in the course of his or her duties as a United States government employee.

ARTICLE AND VIDEO LICENSE AGREEMENT

1. **Defined Terms.** As used in this Article and Video License Agreement, the following terms shall have the following meanings: **"Agreement"** means this Article and Video License Agreement; **"Article"** means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein; **"Author"** means the author who is a signatory to this Agreement; **"Collective Work"** means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; **"CRC License"** means the Creative Commons Attribution-Non Commercial-No Derivs 3.0 Unported Agreement, the terms and conditions of which can be found at: <http://creativecommons.org/licenses/by-nc-nd/3.0/legalcode>; **"Derivative Work"** means a work based upon the Materials or upon the Materials and other pre-existing works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; **"Institution"** means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; **"JoVE"** means MyJoVE Corporation, a Massachusetts corporation and the publisher of The Journal of Visualized Experiments; **"Materials"** means the Article and / or the Video; **"Parties"** means the Author and JoVE; **"Video"** means any video(s) made by the Author, alone or in conjunction with any other parties, or by JoVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion

of the Article, and in which the Author may or may not appear.

2. **Background.** The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JoVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.

3. **Grant of Rights in Article.** In consideration of JoVE agreeing to publish the Article, the Author hereby grants to JoVE, subject to **Sections 4** and **7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the "Open Access" box has been checked in **Item 1** above, JoVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.

ARTICLE AND VIDEO LICENSE AGREEMENT

4. **Retention of Rights in Article.** Notwithstanding the exclusive license granted to JoVE in **Section 3** above, the Author shall, with respect to the Article, retain the non-exclusive right to use all or part of the Article for the non-commercial purpose of giving lectures, presentations or teaching classes, and to post a copy of the Article on the Institution's website or the Author's personal website, in each case provided that a link to the Article on the JoVE website is provided and notice of JoVE's copyright in the Article is included. All non-copyright intellectual property rights in and to the Article, such as patent rights, shall remain with the Author.

5. **Grant of Rights in Video – Standard Access.** This **Section 5** applies if the "Standard Access" box has been checked in **Item 1** above or if no box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby acknowledges and agrees that, Subject to **Section 7** below, JoVE is and shall be the sole and exclusive owner of all rights of any nature, including, without limitation, all copyrights, in and to the Video. To the extent that, by law, the Author is deemed, now or at any time in the future, to have any rights of any nature in or to the Video, the Author hereby disclaims all such rights and transfers all such rights to JoVE.

6. **Grant of Rights in Video – Open Access.** This **Section 6** applies only if the "Open Access" box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby grants to JoVE, subject to **Section 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Video in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Video into other languages, create adaptations, summaries or extracts of the Video or other Derivative Works or Collective Works based on all or any portion of the Video and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. For any Video to which this **Section 6** is applicable, JoVE and the Author hereby grant to the public all such rights in the Video as provided in, but subject to all limitations and requirements set forth in, the CRC License.

7. **Government Employees.** If the Author is a United States government employee and the Article was prepared in the course of his or her duties as a United States government employee, as indicated in **Item 2** above, and any of the licenses or grants granted by the Author hereunder exceed the scope of the 17 U.S.C. 403, then the rights granted hereunder shall be limited to the maximum

rights permitted under such statute. In such case, all provisions contained herein that are not in conflict with such statute shall remain in full force and effect, and all provisions contained herein that do so conflict shall be deemed to be amended so as to provide to JoVE the maximum rights permissible within such statute.

8. **Protection of the Work.** The Author(s) authorize JoVE to take steps in the Author(s) name and on their behalf if JoVE believes some third party could be infringing or might infringe the copyright of either the Author's Article and/or Video.

9. **Likeness, Privacy, Personality.** The Author hereby grants JoVE the right to use the Author's name, voice, likeness, picture, photograph, image, biography and performance in any way, commercial or otherwise, in connection with the Materials and the sale, promotion and distribution thereof. The Author hereby waives any and all rights he or she may have, relating to his or her appearance in the Video or otherwise relating to the Materials, under all applicable privacy, likeness, personality or similar laws.

10. **Author Warranties.** The Author represents and warrants that the Article is original, that it has not been published, that the copyright interest is owned by the Author (or, if more than one author is listed at the beginning of this Agreement, by such authors collectively) and has not been assigned, licensed, or otherwise transferred to any other party. The Author represents and warrants that the author(s) listed at the top of this Agreement are the only authors of the Materials. If more than one author is listed at the top of this Agreement and if any such author has not entered into a separate Article and Video License Agreement with JoVE relating to the Materials, the Author represents and warrants that the Author has been authorized by each of the other such authors to execute this Agreement on his or her behalf and to bind him or her with respect to the terms of this Agreement as if each of them had been a party hereto as an Author. The Author warrants that the use, reproduction, distribution, public or private performance or display, and/or modification of all or any portion of the Materials does not and will not violate, infringe and/or misappropriate the patent, trademark, intellectual property or other rights of any third party. The Author represents and warrants that it has and will continue to comply with all government, institutional and other regulations, including, without limitation all institutional, laboratory, hospital, ethical, human and animal treatment, privacy, and all other rules, regulations, laws, procedures or guidelines, applicable to the Materials, and that all research involving human and animal subjects has been approved by the Author's relevant institutional review board.

11. **JoVE Discretion.** If the Author requests the assistance of JoVE in producing the Video in the Author's facility, the Author shall ensure that the presence of JoVE employees, agents or independent contractors is in accordance with the relevant regulations of the Author's institution. If more than one author is listed at the beginning of this Agreement, JoVE may, in its sole

ARTICLE AND VIDEO LICENSE AGREEMENT

discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JoVE reserves the right, in its absolute and sole discretion and without giving any reason therefore, to accept or decline any work submitted to JoVE. JoVE and its employees, agents and independent contractors shall have full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JoVE has sole discretion as to the method of making and publishing the Materials, including, without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

12. **Indemnification.** The Author agrees to indemnify JoVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JoVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JoVE, making of videos by JoVE, or publication in JoVE or elsewhere by JoVE. The Author shall be responsible for, and shall hold JoVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to

the making of a video by JoVE its employees, agents or independent contractors. All sterilization, cleanliness or decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's expense. All indemnifications provided herein shall include JoVE's attorney's fees and costs related to said losses or damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JoVE, its employees, agents or independent contractors.

13. **Fees.** To cover the cost incurred for publication, JoVE must receive payment before production and publication the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JoVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.

14. **Transfer, Governing Law.** This Agreement may be assigned by JoVE and shall inure to the benefits of any of JoVE's successors and assignees. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to be one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement is required per submission.

CORRESPONDING AUTHOR

Name:

Rebecca Richards-Kortum

Department:

Bioengineering

Institution:

Rice University

Title:

Malcolm Gillis University Professor

Signature:





Date:

10/2/2020

Please submit a **signed** and **dated** copy of this license by one of the following three methods:

1. Upload an electronic version on the JoVE submission site
2. Fax the document to +1.866.381.2236
3. Mail the document to JoVE / Attn: JoVE Editorial / 1 Alewife Center #200 / Cambridge, MA 02140





[Click here to access/download](#)

Supplemental Coding Files

System1_Optics_Enclosure_Top.STL





[Click here to access/download](#)


Supplemental Coding Files

System1_Optics_Enclosure_Bottom.STL







[Click here to access/download](#)
Supplemental Coding Files
LCD_Screen_Holder.STL



[Click here to access/download](#)

Supplemental Coding Files

System2_Optics_Enclosure_Top.STL






[Click here to access/download](#)

Supplemental Coding Files

System2_Optics_Enclosure_Bottom.STL





Click here to access/download
Supplemental Coding Files
MiniFluorimeter_2Diode.ino

