

A

Script file

Sequence

Step	Pump	Infuse	Volume	Rate	Pause	Tubing Diam.	Ch. Mode	No Movement	Stage X	Stage Y	Speed X	Speed Y
1	5	☒	0 uL	0 uL/min	0	3.17	Product	☐	0	0	0	0
2	1	☒	4000 uL	3500 uL/min	0	3.17	Product	☐	0	0	0	0
3	2	☒	4000 uL	3500 uL/min	0	3.17	Product	☐	0	0	0	0
4	3	☒	4000 uL	3500 uL/min	0	3.17	Product	☐	0	0	0	0
5	4	☒	4000 uL	3500 uL/min	0	3.17	Product	☐	0	0	0	0
6	5	☒	0 uL	0 uL/min	0	3.17	Waste	☐	0	0	0	0
7	5	☒	15 uL	300 uL/min	0	3.17	Waste	☒	2	2	2	2

+
Insert below
1
-

B

stop COMPORT for Pumps

SCRIPT

Step	Pump No	Direction	Rate	Volume	Dispensed	Pause	Elapsed	Diameter	CH Mode	Stage Moving	Stage X	Stage Y	Speed X	Speed Y
1	Pump 5	Infusing	0.00 ul/min	0.00 ul	0.00 ul	0 s	0 s	3.17	enabled	no stage	0.00 units	0.00 units	0.00 units/sec	0.00 units/sec
2	Pump 1	Infusing	3500 ul/min	3000 ul	0.00 ul	0 s	0 s	3.17	enabled	no stage	0.00 units	0.00 units	0.00 units/sec	0.00 units/sec
3	Pump 2	Infusing	3500 ul/min	3000 ul	0.00 ul	0 s	0 s	3.17	enabled	no stage	0.00 units	0.00 units	0.00 units/sec	0.00 units/sec
4	Pump 3	Infusing	3500 ul/min	3000 ul	0.00 ul	0 s	0 s	3.17	enabled	no stage	0.00 units	0.00 units	0.00 units/sec	0.00 units/sec
5	Pump 4	Infusing	3500 ul/min	3000 ul	0.00 ul	0 s	0 s	3.17	enabled	no stage	0.00 units	0.00 units	0.00 units/sec	0.00 units/sec
6	Pump 5	Infusing	0.00 ul/min	0.00 ul	0.00 ul	0 s	0 s	3.17	disabled	no stage	0.00 units	0.00 units	0.00 units/sec	0.00 units/sec
7	Pump 5	Infusing	300.00 ul/min	15.00 ul	0.00 ul	0 s	0 s	3.17	disabled	moving	2.00 units	2.00 units	2.00 units/sec	2.00 units/sec

X Axis Encoder (mm) ☒ X Moving ☒ X Direction ☒ X Negative Limit ☒ X Positive Limit

Y Axis Encoder (mm) ☒ Y Moving ☒ Y Direction ☒ Y Negative Limit ☒ Y Positive Limit

Device Handle