

# Journal of Visualized Experiments

## Visualizing and Quantifying Endonuclease-Based Site-Specific DNA Damage

--Manuscript Draft--

<b>Article Type:</b>	Methods Article - Author Produced Video
<b>Manuscript Number:</b>	JoVE62175R2
<b>Full Title:</b>	Visualizing and Quantifying Endonuclease-Based Site-Specific DNA Damage
<b>Corresponding Author:</b>	Pankotai Tibor, Habil, Ph.d. University of Szeged, Faculty of Medicine, Institute of Pathology Szeged, HUNGARY
<b>Corresponding Author's Institution:</b>	University of Szeged, Faculty of Medicine, Institute of Pathology
<b>Corresponding Author E-Mail:</b>	pankotai.tibor@szte.hu
<b>Order of Authors:</b>	Vasiliki Pantazi Ivett Berzsényi Barbara N Borsos Pankotai Tibor, Habil, Ph.d.
<b>Additional Information:</b>	
<b>Question</b>	<b>Response</b>
Please indicate whether this article will be Standard Access or Open Access.	
Please specify the section of the submitted manuscript.	Biology
Please confirm that you have read and agree to the terms and conditions of the author license agreement that applies below:	I agree to the <a href="#">Author License Agreement</a>
Please provide any comments to the journal here.	Nothing
Please indicate whether this article will be Standard Access or Open Access.	Standard Access (\$1400)
Please confirm that you have read and agree to the terms and conditions of the video release that applies below:	I agree to the <a href="#">Video Release</a>

**TITLE:**

Visualizing and Quantifying Endonuclease-Based Site-Specific DNA Damage

**AUTHORS AND AFFILIATIONS:**

Vasiliki Pantazi\*, Ivett Berzsenyi\*, Barbara N Borsos, Tibor Pankotai

Institute of Pathology, Faculty of Medicine, University of Szeged, Szeged, Hungary

\* The authors have contributed equally.

**Email Addresses of Co-Authors:**

Vasiliki Pantazi (vasiliki.pantazi@med.u-szeged.hu)

Ivett Berzsenyi (berzsenyi.ivett@szte.hu)

Barbara N Borsos (borsos.barbara.nikolett@szte.hu)

Tibor Pankotai (pankotai.tibor@szte.hu)

**Corresponding Author:**

Tibor Pankotai (pankotai.tibor@szte.hu)

**KEYWORDS**

DNA repair, chromatin immunoprecipitation, fluorescent microscopy, site-specific DSBs, immunostaining

**SUMMARY:**

This article introduces essential steps of immunostaining and chromatin immunoprecipitation. These protocols are commonly used to study DNA damage-related cellular processes and to visualize and quantify the recruitment of proteins implicated in DNA repair.

**ABSTRACT:**

Cells are continuously exposed to various DNA damaging agents, inducing different cellular responses. Applying biochemical and genetic approaches is essential in revealing cellular events associated with the recruitment and assembly of DNA repair complexes at the site of DNA damage. In the last few years, several powerful tools have been developed to induce site-specific DNA damage. Moreover, novel seminal techniques allow us to study these processes at the single-cell resolution level using both fixed and living cells. Although these techniques have been used to study various biological processes, herein we present the most widely used protocols in the field of DNA repair, Fluorescence Immunostaining (IF) and Chromatin Immunoprecipitation (ChIP), which in combination with endonuclease-based site-specific DNA damage make it possible to visualize and quantify the genomic occupancy of DNA repair factor in a directed and regulated fashion, respectively. These techniques provide powerful tools for the researchers to identify novel proteins bound to the damaged genomic locus as well as their post-translational modifications necessary for their fine-tune regulation during DNA repair.

**INTRODUCTION**

Our genome is constantly being challenged by various DNA damaging agents. These assaults can derive from environmental sources, such as UV light or irradiation, as well as from endogenous sources, such as metabolic by-products caused by oxidative stress or replication errors<sup>1,2</sup>. These lesions can affect the integrity of either one or both DNA strands, and if the generated errors become persistent, it frequently leads to translocations and genome instability, which may result in tumorigenesis<sup>3,4</sup>. To maintain genome integrity, multiple repair systems have been developed during evolution. According to the chemical and physical properties of specific types of DNA damage, multiple repair mechanisms can be activated. Mismatches, abasic sites, single-strand breaks, and 8-oxoguanine (8-oxoG) can be removed either by mismatch repair or base-excision repair pathway<sup>5,6</sup>. Lesions caused by UV-induced photoproducts and bulky adducts can be repaired either by nucleotide-excision repair (NER) or DNA double-strand break repair (DSBR) process<sup>7,8</sup>. NER consists of two main sub-pathways: transcription-coupled NER (TC-NER) and global genomic NER (GG-NER). Regarding the cell cycle phase, following DNA double-strand break induction, two sub-pathways can be activated non-homologous end joining (NHEJ) and homologous recombination (HR)<sup>1,9</sup>. NHEJ, which is the dominant pathway in resting cells, can be activated in all cell cycle phases, representing a faster but error-prone pathway<sup>10</sup>. On the other hand, HR is an error-free pathway, in which the DSBs are repaired based on sequence-homology search of the sister chromatids, therefore it is mainly present in S and G2 cell cycle phases<sup>11</sup>. Furthermore, microhomology-mediated end joining (MMEJ) is another DSB repair mechanism, distinct from the aforementioned ones, based on a KU70/80- and RAD51-independent way of re-ligation of previously resected microhomologous sequences flanking the broken DNA ends. Therefore, MMEJ is considered to be error-prone and highly mutagenic<sup>12</sup>. During DNA repair, DSBs can induce the DNA damage response (DDR), which results in the activation of checkpoint kinases that halt the cell cycle during repair<sup>13-15</sup>. The DDR is activated as a response to the recruitment and extensive spreading of initiator key players of the repair process around the lesions, contributing to the formation of a repair focus. In this early signaling cascade, the ATM (Ataxia Telangiectasia Mutated) kinase plays a pivotal role by catalyzing the phosphorylation of the histone variant H2AX at Ser-139 (referred to as  $\gamma$ H2AX) around the lesion<sup>16</sup>. This early event is responsible for the recruitment of additional repair factors and the initiation of downstream repair processes. Although the exact function of the recruited proteins at the repair focus has not yet been fully characterized, the formation and the dynamics of repair foci have been investigated by several laboratories. These markers are extensively used to follow the repair kinetics, but their precise role during the repair process remains elusive. Due to the great importance yet poor understanding of DNA repair-related cellular processes, several methods have been developed so far to induce and visualize the DDR.

Various methods and systems have been established to induce the desired type of DNA damage. For instance, some agents [such as neocarzinostatin (NCS), phleomycin, bleomycin,  $\gamma$ -irradiation, UV] can induce large numbers of random DNA breaks at non-predictive genomic positions, while others (endonucleases, such as AsiSI, I-PpoI or I-SceI, as well as laser striping) can induce DNA breaks at known genomic loci<sup>17-21</sup>. Here, we focus on the endonuclease-based techniques currently used to study the DDR in mammalian and yeast cells. Aside from highlighting the principles of these techniques, we emphasize both their advantages and disadvantages.

## PROTOCOL

### 1. Immunodetection of specific proteins

#### 1.1. Preparation of cell culture and experimental setup

1.1.1. Maintain U2OS cells in monolayers in DMEM culture medium supplemented with 10% fetal calf serum, 2 mM glutamine, and 1% antibiotic-antimycotic solution.

NOTE: For endonuclease-based DNA damage induction, use charcoal-treated or steroid-free medium to avoid system leakiness.

1.1.2. Grow cells in a humidified 5% CO<sub>2</sub> environment at 37 °C until 80% confluency, renewing medium every 2–3 days.

1.1.3. Aspirate the medium and wash the cells with 1x PBS. Detach cells with Trypsin-EDTA solution. When the cells detach, stop the trypsin activity by adding culture medium to the cells, yielding a cell suspension.

1.1.4. Count the cells using a cell counting chamber. Plate 2 x 10<sup>4</sup> cells/mL/well on a 24-well plate, with sterile 12 mm round coverslips in each well.

1.1.5. Incubate cells for 24 h at 37 °C in a humidified 5% CO<sub>2</sub> atmosphere to allow attachment onto the coverslips.

1.1.6. Treat the cells with 10 ng/mL of neocarzinostatin (NCS) by directly pipetting the damaging agent to the cultured medium. Incubate the cells with the NCS-containing medium for 15 min, then wash them with 1x PBS and add fresh, supplemented culture medium to the cells. Otherwise, use appropriate agent (i.e., 4-OHT) to induce DSBs via endonuclease-based systems without refreshing the medium<sup>22</sup>.

NOTE: Alternatively, use irradiation to induce DNA damage, ranging from 30 min up to 8 h of recovery time by using neutron flux between 2–20 Gy<sup>23</sup>.

1.1.7. Incubate cells for 1–8 h at 37 °C in a humidified 5% CO<sub>2</sub> atmosphere to follow the kinetics of DNA repair.

#### 1.2. Fixation of cells

NOTE: 300–500 µL of solutions/well should be used in the following steps (steps 1.1.2-1.1.5) to cover all cells adequately. Each incubation and washing step (except the antibody incubation) should be performed on an orbital shaker with gentle agitation.

1.2.1. Following DSB induction and incubation of cells, remove the medium from the attached

cells and wash the cells once with 1x PBS.

1.2.2. Fix cells with 4% formaldehyde-PBS solution for 20 min at 25 °C.

### **1.3. Permeabilization of cells**

1.3.1. Remove the fixing solution and wash cells three times with 1x PBS for 5 min each.

1.3.2. Remove the PBS and add 0.2% Triton X-100 dissolved in PBS. Incubate the samples for 20 min.

### **1.4. Blocking of non-specific binding sites**

1.4.1. Wash the cells three times with 1x PBS.

1.4.2. Block non-specific binding sites with 5% BSA (Bovine Serum Fraction V albumin) diluted in PBST (1x PBS supplemented with 0.1% Tween-20), and incubate the permeabilized samples for at least 20 min.

### **1.5. Immunofluorescence staining**

1.5.1. Add the proper amount of primary antibody (i.e., anti-γH2AX, anti-DNA-PKcs) diluted in 1% BSA-PBST solution. Place each coverslip upside down onto a paraffin film over a 10 μL droplet of the diluted anti-γH2AX antibody.

NOTE: In case of co-immunostaining dilute appropriately both antibodies in the same 1% BSA-PBST solution.

1.5.2. Incubate the samples in a humidity chamber for 1.5 h at 4 °C.

NOTE: Incubation can also be performed at 4 °C overnight.

1.5.3. Place the coverslips back being side up into the 24-well plate and wash three times for 5 min with 1x PBS.

1.5.4. Add the proper amount of secondary antibody diluted in 1% BSA-PBST. Place each coverslip upside down onto a paraffin film over a 10 μL droplet of the diluted antibody.

1.5.5. Incubate the samples in a humidity chamber at 4 °C for 1 h.

1.5.6. Place the coverslips back being side up into the 24-well plate and wash three times for 5 min with 1x PBS.

1.5.7. Before removing the last PBS washing solution, gently take out the coverslips using a

tweezer and needle and then place them upside down onto glass slides with droplets of mounting medium (supplemented with DAPI).

NOTE: Avoid the formation of air bubbles. When the mounting medium dries, it is recommended to seal the edges of the coverslips with nail polish to prevent shriveling of the samples.

## **2. Chromatin immunoprecipitation**

### **2.1. Cell collection, crosslinking, cell and nuclear lysis, and DNA fragmentation**

2.1.1. Culture approximately  $2 \times 10^7$  cells/mL in a 150 mm dish for each sample.

2.1.2. Remove the culture medium and wash the cells twice with ice-cold 1x PBS.

2.1.3. Fix the cells with 1% formaldehyde-PBS solution, place the plates on an orbital shaker, and agitate gently for 20 min.

NOTE: Formaldehyde is volatile; always prepare a fresh working solution. In some cases, the formaldehyde solution contains methanol to stabilize it, but it is better to use a methanol-free solution to avoid interference with downstream reactions.

2.1.4. Stop the fixation with 125 mM glycine and incubate on an orbital shaker with gentle agitation for 5 min at 25 °C.

2.1.5. Place the plates on ice and wash twice with ice-cold 1x PBS.

2.1.6. Scrape the cells in ice-cold 1x PBS and transfer them into 15 mL conical tubes.

2.1.7. Centrifuge the cells at  $2,500 \times g$  for 5 min at 4 °C.

2.1.8. Carefully aspirate the supernatant and resuspend the pellet in 2 mL of cell lysis buffer [5 mM PIPES pH 8.0, 85 mM KCl, 0.5% NP-40, 1x PIC (protease inhibitor cocktail)] and incubate on ice for 10 min.

2.1.9. Centrifuge the cell suspension at  $2,500 \times g$  for 5 min at 4 °C.

2.1.10. Carefully discard the supernatant and resuspend the pellet in 500–1,500 µL of nuclear lysis buffer (50 mM Tris-HCl pH 8.0, 10 mM EDTA pH 8.0, 0.8% SDS, 1x PIC) and incubate on ice for 30–60 min. Transfer the lysate into a polystyrene conical tube suitable for sonication.

NOTE: Since the nuclear lysis buffer contains SDS, it will precipitate on ice, and the solution will turn white. The solution should turn transparent following sonication.

2.1.11. Sonicate the lysate to shear DNA to an average fragment size of 300–1,000 bp.

NOTE: The appropriate sonication cycles and conditions should be set according to the cell type and the sonication equipment. Fragments smaller than 200 bp are not suitable for ChIP, because the nucleosome-DNA interactions can be disrupted.

## **2.2. Reversal of crosslinking, determination of sonicated fragment sizes**

2.2.1. Take out 100  $\mu$ L of the sonicated sample to verify the fragment size of the sonicated chromatin. The remaining chromatin should be stored at  $-80^{\circ}\text{C}$ .

2.2.2. Add 0.5 mg/mL RNase A to each 100  $\mu$ L of the sample and incubate them at  $37^{\circ}\text{C}$  for 20 min to activate the RNase.

2.2.3. Incubate the samples at  $65^{\circ}\text{C}$  overnight.

2.2.4. The next day, add 500  $\mu\text{g/mL}$  of Proteinase K and 0.5% SDS, and incubate the samples at  $50^{\circ}\text{C}$  for 3 h.

2.2.5. Add 0.5 volume of phenol and 0.5 volume chloroform-isoamyl alcohol mix (24:1) to each sample.

2.2.6. Vortex for 1 min.

2.2.7. Centrifuge at  $13,000 \times g$  for 10 min.

2.2.8. Transfer the upper aqueous phase to a new microcentrifuge tube.

2.2.9. Add 1 volume chloroform-isoamyl alcohol mix (24:1) to each sample.

2.2.10. Vortex for 1 min.

2.2.11. Centrifuge at  $13,000 \times g$  for 10 min.

2.2.12. Transfer the upper aqueous phase to a new microcentrifuge tube.

2.2.13. Add 2.5 volumes of 96% ethanol and 0.1 volume of 3 M Na-Acetate pH 5.2.

2.2.14. Incubate for at least 20 min at  $-80^{\circ}\text{C}$ .

2.2.15. Centrifuge the samples at  $13,000 \times g$  for 10 min at  $4^{\circ}\text{C}$ .

2.2.16. Remove the ethanol and wash the pellet with 400  $\mu$ L of 70% ethanol.

2.2.17. Centrifuge the samples at  $13,000 \times g$  for 10 min at  $4^{\circ}\text{C}$ .

2.2.18. Remove the ethanol and air dry the pellet.

2.2.19. Resuspend the pellet in 10  $\mu$ L of TE.

2.2.20. Run the samples on a 0.8% agarose gel. The sonicated chromatin size should be around 500 bp.

NOTE: Use bromophenol blue-free loading buffer because the size of this dye is approximately 500 bp, which can disturb the proper detection of chromatin fragments. Instead, it is recommended to use loading buffer complemented with xylene-cyanole, which is approximately 3,000 bp.

2.2.21. If the chromatin size is acceptable, dilute the frozen chromatin samples from step 2.1. in 3 volumes of dilution buffer (10 mM Tris-HCl pH 8.0, 0.5 mM EGTA pH 8.0, 1% Triton X-100, 140 mM NaCl, 1x PIC) and mix the samples via rotation for 10 min at 4 °C.

NOTE: This step is necessary to dilute the SDS present in the nuclear lysis buffer to avoid interference with downstream reactions, including the measurement of chromatin concentration.

2.2.22. Measure the DNA concentration of the chromatin samples at 260/280 nm using a spectrophotometer.

### **2.3. Preparation of beads, pre-clearing, and immunoprecipitation**

2.3.1. Prepare beads (sheep anti-rabbit or mouse IgG) for the pre-clearing and immunoprecipitation steps. Wash beads twice for 10 min at 4 °C with RIPA buffer (50 mM Tris-HCl pH 8.0, 1 mM EDTA pH 8.0, 1% Triton X-100, 0.1% Na-DOC, 0.1% SDS, 150 mM NaCl and 1X PIC).

2.3.2. Resuspend the beads in the same volume of RIPA buffer as in step 2.3.1.

2.3.3. Pre-clear 25–30  $\mu$ g chromatin of each sample with 4  $\mu$ L of the beads via rotation for 1–2 h at 4 °C.

NOTE: Add RIPA buffer to each chromatin sample up to 500  $\mu$ L of final volume to let the samples mix properly under rotation. Do not forget to take out chromatin for NAC (No Antibody Control) and TIC (Total Input Control) in the case of each sample set. TICs only require a final volume of up to 200  $\mu$ L.

2.3.4. Precipitate the beads with a magnet and transfer the supernatant to a new microcentrifuge tube.

2.3.5. Add the appropriate amount of antibody to each chromatin sample (except NAC and TIC) and rotate overnight at 4 °C.

2.3.6. Next day, add 40 µL of washed beads to each sample (except TIC) and incubate them overnight, rotating at 4 °C.

## **2.4. Washing**

2.4.1. Wash once with 300 µL of Low Salt buffer (20 mM Tris-HCl pH 8.0, 150 mM NaCl, 2 mM EDTA pH 8.0, 1% Triton X-100, 0.1% SDS, 1x PIC) for 10 min via rotation at 4 °C.

2.4.2. Wash once with 300 µL of High Salt buffer (20 mM Tris-HCl pH 8.0, 300 mM NaCl, 2 mM EDTA pH 8.0, 1% Triton X-100, 0.1% SDS, 1x PIC) for 10 min via rotation at 4 °C.

2.4.3. Wash once with 300 µL of LiCl buffer (250 mM LiCl, 1% NP-40, 1% Na-DOC, 1 mM EDTA pH 8.0, 10 mM Tris-HCl pH 8.0, 1x PIC) for 10 min via rotation at 4 °C.

2.4.4. Wash twice with 300 µL of TE (10 mM Tris-HCl pH 8.0, 1 mM EDTA pH 8.0) for 10 min via rotation, for the first wash at 4 °C, and the second wash at 25 °C.

## **2.5. Elution**

2.5.1. Add 200 µL of the elution buffer (1% SDS and 100 mM NaHCO<sub>3</sub>) to the beads and incubate at 65 °C in a thermo-shaker for 15 min with continuous shaking (approx. 400 RPM). Transfer the supernatant to a new tube and elute beads again in 200 µL of Elution buffer. Combine eluates (400 µL final volume).

2.5.2. Add NaCl to a final concentration of 200 mM in each sample. Supplement the TIC samples with 200 µL of Elution buffer and add NaCl as well.

NOTE: From this step, TIC should be handled under the same conditions as the other samples.

2.5.3. Incubate the samples at 65 °C (without shaking) for at least 6 h.

5.4. Add 1 mL of cold 100% ethanol to each sample, rotate the tubes twice to mix, and precipitate DNA overnight at -80 °C.

2.5.5. The next day, centrifuge for 30 min at 13,000 x g at 4 °C.

2.5.6. Discard the supernatant and wash the pellet with 70% EtOH.

2.5.7. Centrifuge the samples at 13,000 x g for 10 min at 4 °C.

2.5.8. Discard the supernatant and air dry the pellets.

2.5.9. Resuspend the pellets in 100  $\mu$ L of TE and add 0.5 mg/mL of RNase A to each sample. Incubate at 37 °C for 20 min to activate the RNase.

## **2.6. Reversal of crosslinking**

2.6.1. Add 500  $\mu$ g/mL of Proteinase K and 0.5% SDS and incubate the samples at 50 °C for 2 h.

NOTE: If proceeding to ChIP-seq, avoid phenol-chloroform extraction, as it inhibits the downstream NGS process. Instead, it is recommended to use a commercially available kit (see **Table of Materials**).

2.6.2. Add 0.5 volume of phenol and 0.5 volume of chloroform-isoamyl alcohol mix (24:1) to each sample.

2.6.3. Vortex for 1 min.

2.6.4. Centrifuge at 13,000 x *g* for 10 min.

2.6.5. Transfer the upper aqueous phase to a new microcentrifuge tube.

2.6.6. Add 1 volume of chloroform-isoamyl alcohol mix (24:1) to each sample.

2.6.7. Vortex for 1 min.

2.6.8. Centrifuge at 13,000 x *g* for 10 min.

2.6.9. Transfer the upper aqueous phase to a new microcentrifuge tube.

## **2.7. DNA extraction**

2.7.1. Add 2.5 volumes of 96% ethanol and 0.1 volume of 3 M Na-Acetate pH 5.2.

2.7.2. Incubate for at least 20 min at –80 °C.

2.7.3. Centrifuge the samples at 13,000 x *g* for 10 min at 4 °C.

2.7.4. Remove the ethanol and wash the pellet with 400  $\mu$ L of 70% ethanol.

2.7.5. Centrifuge the samples at 13,000 x *g* for 10 min at 4 °C.

2.7.6. Remove the ethanol and air dry the pellet.

2.7.7. Resuspend the pellet in 50  $\mu$ L of TE.

## REPRESENTATIVE RESULTS

Studying site-directed DSB-induced repair processes in cells is performed via either stable or transient transfection. However, it should be noted that stable transfection ensures a homogenous cell population, which gives a unified and thus more reliable cellular response. In the case of transient transfection, only a small proportion of the cell population takes up and maintains the plasmid, which introduces diversity into the experiment. Establishing ER-I-Ppol or ER-AsiSI endonuclease-based cell systems require a 50% confluent cell population, which is more effectively transfected with plasmids encoding the endonuclease. For transfection, commercially available transfection reagents or viral infection-based methods can also be used. If a microscopic visualization technique is to be applied and transient transfection is required, the directed DSBs can be induced by 4-OHT addition 24 h after transfection, which binds to the ER-fused endonucleases and allows the nuclear translocation and DSB induction. To determine the most appropriate time-points, immunofluorescence-based microscopy and western blot detection of  $\gamma$ H2AX at different time-points following 4-OHT treatment can be performed. Under physiological conditions, a maximum of 10–15  $\gamma$ H2AX foci per cell can be detected, and the formation of strong repair foci can be triggered by endonucleases (or various other techniques not discussed here e.g., laser microirradiation). A typical I-Ppol endonuclease leads to the formation of elevated  $\gamma$ H2AX signals around the nucleolus by inducing DSBs at ribosomal DNA (rDNA). If the breaks are repaired by NHEJ or HR, the number of repair foci decreases over time. For this reason, representative time-points at 0 h, 30 min, 1, 2, 4, and 8 h following 4-OHT treatments are recommended. To track DNA repair processes, the most commonly used cell line is U2OS, as all known repair pathways are fully functional in these cells. When investigating several proteins in the same cells, co-localization can be studied by combining antibodies conjugated with different fluorophores with different emission wavelengths raised in different animal species as shown in **Figure 1**. Therein, the induction of DSBs via an inducible stable cell line is represented which is based on the ER-AsiSI restriction endonuclease fused with hemagglutinin tag (HA). Doxycycline can induce the expression and sequestration of HA-ER-AsiSI in the cytoplasm which can be tracked using an antibody against HA (**Figure 1**. third column, second row). Incubation for 4 h with 4-OHT, 24 h after doxycycline addition, can induce high number of DSBs since the endonuclease has been translocated to the nucleus (**Figure 1**. third column, third row and second column, third row). DSBs can be visualized by using an antibody recognizing  $\gamma$ H2AX.

[Place **Figure 1** here].

Co-localization of repair proteins at the damage site indicates that they are recruited to the same DNA lesion site, but they do not necessarily interact with each other. The resolution of the confocal microscopy is approximately 300 nm; to determine the binding pattern of specific repair proteins at the break site, super-resolution microscopy (STORM) is instead recommended<sup>24</sup>. However, this method requires expensive microscopic equipment and an expert researcher. Alternatively, the binding pattern of the repair proteins can be examined by chromatin immunoprecipitation using DlvA or U2OS-pEP15 stable cell lines, which can express AsiSI and I-Ppol endonucleases in a regulated fashion, respectively<sup>17,21</sup>. Upon 4-OHT addition, both endonucleases can cut the DNA in a sequence specific manner which provides us the opportunity

to design locus specific primers to the expected break sites and their surrounding genomic regions. By applying  $\gamma$ H2AX antibody in the immunoprecipitation part of the ChIP, we can temporally follow the DNA repair kinetics upon different conditions (such as silencing or inhibition of certain repair factors of interest, i.e., DNA-PKcs). A typical experimental result obtained using ChIP-qPCR is represented in **Figure 2**. Therein, the temporal enrichment of  $\gamma$ H2AX is demonstrated as a response to I-Ppol-induced DNA damage. On the left part of the image, the timely detected  $\gamma$ H2AX signal is shown at the break site while on the right part, the  $\gamma$ H2AX distribution is represented at a control gene region at which DSBs have not been induced.

[Place **Figure 2** here]

## FIGURE LEGENDS

**Figure 1: Immunofluorescence microscopy is used to detect  $\gamma$ H2AX in cultured cells expressing HA-ER-AsiSI.** DOX (doxycycline) addition activates the cytoplasmic expression of HA-ER-AsiSI and 4-OHT (4-hydroxytamoxifen) (4 h) induces the nuclear translocation of the fusion protein, which leads to the induction of DSBs at known genomic positions. HA (hemagglutinin) staining (anti-HA antibody) represents the HA-ER-AsiSI fusion protein in green, and the induction of breaks is verified by  $\gamma$ H2AX staining (anti- $\gamma$ H2AX antibody) in red. Scale bars represent 20  $\mu$ m.

**Figure 2: Temporal enrichment of  $\gamma$ H2AX as determined by chromatin immunoprecipitation in response to I-Ppol-induced DNA damage.** Left,  $\gamma$ H2AX signal at the break site; right,  $\gamma$ H2AX distribution at a control region where DSBs were not induced (anti- $\gamma$ H2AX antibody). The represented results are derived from one biological experiment, and error bars indicate the variations of the corresponding sample replicates. N=3.

## DISCUSSION

Although DNA repair is a relatively recent research field, our knowledge is rapidly expanding with the help of various biochemical and microscopic methods. Preserving genetic information is crucial for cells since mutations occurring in genes involved in repair processes are among the leading causes of tumorigenesis and therefore elucidating the key steps of DNA repair pathways is essential.

Biochemical techniques (i.e., western blot, immunoprecipitation, mass-spectrometry, etc.) require large number of cells and the studied repair processes represent a snapshot of the desired cell population. Performing ChIP experiments is laboring, troublesome and many considerations have to be taken into account when designing specific experiment to study the process of DSB repair. The following steps are some examples: (I) cells should be properly lysed; it is highly recommended to apply a two-step lysis method by using separate cell- and nuclear lysis buffer to ensure higher accessibility to the chromatin fraction (II) chromatin should not be over- or undersonicated; the appropriate conditions of sonication should be optimized to each cell type beforehand (III) the appropriate amount of purified antibody should be determined since antibodies against the same protein from different companies exhibit non-identical characteristics (IV) for the efficient immunoprecipitation, 25-30  $\mu$ g initial chromatin should be used in each condition (V) the appropriate time of fixation should be optimized, as overfixation

can result in false positive results by cross-linking the distant protein complexes and underfixation can prevent the proper cross-linking between DNA and the desired protein (VI) based on the applied antibody, the type of the beads (protein A or G) should be carefully determined (VII) during washing steps, the order of the washing buffers should be thoroughly kept to avoid the release of the antibody from the beads (VIII) during DNA extraction, phenol traces should be properly eliminated to avoid decreasing the efficiency of the further downstream reactions. Since this method reduces the yield of recovered DNA, our personal advice is to rather use specific DNA-purification kit. When all critical points have been properly addressed, ChIP can provide valuable data for the occupancy of the desired protein at different genomic loci and can unravel critical steps in DNA repair processes.

However, ChIP combined with qPCR is an indirect approach to study the protein distribution at selected genomic regions and does not allow to specifically recognize the DNA binding site nor directly examine the function of the protein. Mono- or polyclonal antibodies used to capture the protein-DNA complexes can also cross-react with other proteins leading to false-positive data and therefore, the antibodies used in this technique should be ChIP-grade and highly specific against the protein of interest. However, ChIP is a widely used technique and further approaches based on this have been developed, such as ChIP-on-chip and ChIP-seq. The former one relies on hybridizing the immunoprecipitated and purified DNA fragments on a microarray with a large variety of small random DNA sequences which further amplifies the annealed sequences providing valuable information about protein binding sites. However, ChIP-seq has emerged as an attractive alternative approach since it provides genome-wide mapping of protein-DNA complexes with higher resolution than ChIP-on-chip and high-throughput genome sequencing. ChIP-seq has revolutionized the field of DNA repair by disclosing DNA binding sites of various transcription factors providing insights into the gene regulation and unravelling chromatin landscapes in a genome-wide scale<sup>25</sup>. According to this, the field of DNA repair has tremendously profited from ChIP-seq since these data play a crucial role in various diseases and biological pathways, such as cancer progression. Nonetheless, various modifications of the ChIP method have been developed such as HaloChIP, which does not require specific antibody against the protein of interest but rather uses sequences encoding DNA-binding proteins fused with HaloTag, which are transfected to cells and subsequently to crosslinking, the desired protein-DNA complexes can be captured by using HaloLink Resin. However, this technique relies on overexpression and not the endogenous level of the desired proteins which can result in misinterpreted data<sup>26</sup>.

Furthermore, microscopic techniques provide valuable information about the spatiotemporal tracking of DNA damage repair, even in a single-cell level. The rapid improvement of antibodies raised against specific repair proteins has led to a deeper understanding of the mechanism of NER and DSB repair sub-pathways, as well as their post-translational regulation. The microscopic field has been revolutionized by high-resolution techniques, such as super-resolution microscopy, which allows the visualization of DNA damage-induced cellular processes at the nucleosomal level, as well as ensuring the accurate mapping of protein co-localization<sup>24</sup>. However, it should be noted that the variance in cell lineages must be considered during the experiment as the rate of the repair can diverge, which can make the results difficult to interpret. Considering the rapid

evolution of fluorescence imaging methodology and the deliberate design of the experimental setup, a precious opportunity to investigate the DNA damage-induced cellular and molecular responses at a single protein resolution in a single-cell level is on its way to perfection.

In conclusion, the combination of super-resolution microscopy and single-cell sequencing methodology can significantly improve our understanding of the DNA repair field.

#### ACKNOWLEDGMENT:

This research was funded by the National Research, Development and Innovation Office grant GINOP-2.3.2-15-2016-00020, GINOP-2.3.2-15-2016-00036, GINOP-2.2.1-15-2017-00052, EFOP 3.6.3-VEKOP-16-2017-00009, NKFI-FK 132080, the János Bolyai Research Scholarship of the Hungarian Academy of Sciences BO/27/20, ÚNKP-20-5-SZTE-265, EMBO short-term fellowship 8513, and the Tempus Foundation.

#### DISCLOSURES:

None

#### REFERENCES

- 1 Borsos, B. N., Majoros, H. & Pankotai, T. Ubiquitylation-Mediated Fine-Tuning of DNA Double-Strand Break Repair. *Cancers (Basel)*. **12** (6), doi:10.3390/cancers12061617 (2020).
- 2 Borsos, B. N., Majoros, H. & Pankotai, T. Emerging Roles of Post-Translational Modifications in Nucleotide Excision Repair. *Cells*. **9** (6), doi:10.3390/cells9061466 (2020).
- 3 Stephens, P. J. *et al.* The landscape of cancer genes and mutational processes in breast cancer. *Nature*. **486** (7403), 400-404, doi:10.1038/nature11017 (2012).
- 4 Turnbull, C. *et al.* Gene-gene interactions in breast cancer susceptibility. *Human Molecular Genetics*. **21** (4), 958-962, doi:10.1093/hmg/ddr525 (2012).
- 5 Saxowsky, T. T., Meadows, K. L., Klungland, A. & Doetsch, P. W. 8-Oxoguanine-mediated transcriptional mutagenesis causes Ras activation in mammalian cells. *Proceedings of the National Academy of Sciences of the United States of America*. **105** (48), 18877-18882, doi:10.1073/pnas.0806464105 (2008).
- 6 Jiricny, J. The multifaceted mismatch-repair system. *Nature Reviews Molecular Cell Biology*. **7** (5), 335-346, doi:10.1038/nrm1907 (2006).
- 7 Hanawalt, P. C. & Spivak, G. Transcription-coupled DNA repair: two decades of progress and surprises. *Nature Reviews Molecular Cell Biology*. **9** (12), 958-970, doi:10.1038/nrm2549 (2008).
- 8 Kanaar, R., Wyman, C. & Rothstein, R. Quality control of DNA break metabolism: in the 'end', it's a good thing. *EMBO Journal*. **27** (4), 581-588, doi:10.1038/emboj.2008.11 (2008).
- 9 Lemaitre, C. *et al.* Nuclear position dictates DNA repair pathway choice. *Genes & Development*. **28** (22), 2450-2463, doi:10.1101/gad.248369.114 (2014).
- 10 Lieber, M. R. The mechanism of double-strand DNA break repair by the nonhomologous DNA end-joining pathway. *Annual Review of Biochemistry*. **79** 181-211, doi:10.1146/annurev.biochem.052308.093131 (2010).

573 11 Lambert, S. & Lopez, B. S. Characterization of mammalian RAD51 double strand break  
574 repair using non-lethal dominant-negative forms. *EMBO Journal*. **19** (12), 3090-3099,  
575 doi:10.1093/emboj/19.12.3090 (2000).

576 12 Xu, S. *et al.* p300-mediated acetylation of histone demethylase JMJD1A prevents its  
577 degradation by ubiquitin ligase STUB1 and enhances its activity in prostate cancer. *Cancer*  
578 *Research*. doi:10.1158/0008-5472.CAN-20-0233 (2020).

579 13 Kastan, M. B. & Bartek, J. Cell-cycle checkpoints and cancer. *Nature*. **432** (7015), 316-323,  
580 doi:10.1038/nature03097 (2004).

581 14 Roy, R., Chun, J. & Powell, S. N. BRCA1 and BRCA2: different roles in a common pathway  
582 of genome protection. *Nature Reviews Cancer*. **12** (1), 68-78, doi:10.1038/nrc3181 (2011).

583 15 Krenning, L., van den Berg, J. & Medema, R. H. Life or Death after a Break: What  
584 Determines the Choice? *Molecular Cell*. **76** (2), 346-358,  
585 doi:10.1016/j.molcel.2019.08.023 (2019).

586 16 Rogakou, E. P., Pilch, D. R., Orr, A. H., Ivanova, V. S. & Bonner, W. M. DNA double-stranded  
587 breaks induce histone H2AX phosphorylation on serine 139. *Journal of Biological*  
588 *Chemistry*. **273** (10), 5858-5868, doi:10.1074/jbc.273.10.5858 (1998).

589 17 Caron, P. *et al.* WWP2 ubiquitylates RNA polymerase II for DNA-PK-dependent  
590 transcription arrest and repair at DNA breaks. *Genes & Development*. **33** (11-12), 684-704,  
591 doi:10.1101/gad.321943.118 (2019).

592 18 Caron, P. *et al.* Cohesin protects genes against gammaH2AX Induced by DNA double-  
593 strand breaks. *PLoS Genetics*. **8** (1), e1002460, doi:10.1371/Journal.pgen.1002460 (2012).

594 19 Berkovich, E., Monnat, R. J., Jr. & Kastan, M. B. Assessment of protein dynamics and DNA  
595 repair following generation of DNA double-strand breaks at defined genomic sites. *Nature*  
596 *Protocols*. **3** (5), 915-922, doi:10.1038/nprot.2008.54 (2008).

597 20 Paques, F. & Duchateau, P. Meganucleases and DNA double-strand break-induced  
598 recombination: perspectives for gene therapy. *Current Gene Therapy*. **7** (1), 49-66,  
599 doi:10.2174/156652307779940216 (2007).

600 21 Pankotai, T., Bonhomme, C., Chen, D. & Soutoglou, E. DNAPKcs-dependent arrest of RNA  
601 polymerase II transcription in the presence of DNA breaks. *Nature Structural & Molecular*  
602 *Biology*. **19** (3), 276-282, doi:10.1038/nsmb.2224 (2012).

603 22 Iacovoni, J. S. *et al.* High-resolution profiling of gammaH2AX around DNA double strand  
604 breaks in the mammalian genome. *EMBO Journal*. **29** (8), 1446-1457,  
605 doi:10.1038/emboj.2010.38 (2010).


606 23 Poinsignon, C. *et al.* Phosphorylation of Artemis following irradiation-induced DNA  
607 damage. *European Journal of Immunology*. **34** (11), 3146-3155,  
608 doi:10.1002/eji.200425455 (2004).



609 24 Varga, D., Majoros, H., Ujfaludi, Z., Erdelyi, M. & Pankotai, T. Quantification of DNA  
610 damage induced repair focus formation via super-resolution dSTORM localization  
611 microscopy. *Nanoscale*. **11** (30), 14226-14236, doi:10.1039/c9nr03696b (2019).

612 25 Kim, J. A., Kruhlak, M., Dotiwala, F., Nussenzweig, A. & Haber, J. E. Heterochromatin is  
613 refractory to gamma-H2AX modification in yeast and mammals. *Journal of Cell Biology*.  
614 **178** (2), 209-218, doi:10.1083/jcb.200612031 (2007).

615 26 Daniels, D. L. & Urh, M. Isolation of intracellular protein--DNA complexes using HaloCHIP,  
616 an antibody-free alternative to chromatin immunoprecipitation. *Methods in Molecular*  
617 *Biology*. **977** 111-124, doi:10.1007/978-1-62703-284-1\_9 (2013).

618







[Click here to access/download](#)

**Table of Materials**

JoVE\_Materials\_Pantazi\_Berzsenyi\_Final.xlsx





Date: 17/05/2021

Dear Dr. Nam Nguyen,

Please find here a revised version of the manuscript entitled ‘**Visualising and quantifying endonuclease-based site-specific DNA damage**’, which was evaluated by three expert reviewers, and based on their opinion, you advised us to revise and resubmit it for publication in *JoVe*.

We thank the critiques and advices that we mostly accepted and modified the interpretation of our data accordingly. We sincerely believe that, as a result of these changes, the quality of our manuscript has been significantly improved, and the present form of the MS will be suitable for publication in *JoVe*. Our responses to the Editorial comments are provided here.

In short, the changes we introduced into the new version of the manuscript are as follows:

- We have made modifications in the text based on the recommendations of the Editor.
- We have modified the video and the narration accordingly to increase the homogeneity between the video file and the written manuscript. In the current form of the video and the written manuscript reflecting of each other.
- We have inserted chapter title cards and we applied American English spelling.
- We have removed commercials from the text and from the narration.
- We have added the requested chapters to the video file (result and conclusion. We also re-organized the chapter according to the request of the Editor.
- We have changed the narration accordingly.

We declare that our manuscript has not been submitted for publication elsewhere.

We thank all the Reviewer’s comments and your editorial contribution to the publication of this work.

We are looking forward to hearing your positive response.

Sincerely,

Tibor Pankotai

Tibor Pankotai, Ph.D.  
Research fellow and Principal Investigator  
Institute of Pathology, University of Szeged  
1, Allomas street Szeged H6725, HUNGARY  
e-mail: pankotai.tibor@szte.hu