

Journal of Visualized Experiments

Natural Killer (NK) and CAR-NK cell expansion method using membrane bound-IL-21-modified B cell line --Manuscript Draft--

| | |
|--|--|
| Article Type: | Methods Article - JoVE Produced Video |
| Manuscript Number: | JoVE62336R3 |
| Full Title: | Natural Killer (NK) and CAR-NK cell expansion method using membrane bound-IL-21-modified B cell line |
| Corresponding Author: | Minh Ma, BS Rutgers New Jersey Medical School Newark, NEW JERSEY UNITED STATES |
| Corresponding Author's Institution: | Rutgers New Jersey Medical School |
| Corresponding Author E-Mail: | mtm265@gsbs.rutgers.edu |
| Order of Authors: | Minh Ma, BS Saiaditya Badeti James Kim Dongfang Liu |
| Additional Information: | |
| Question | Response |
| Please indicate whether this article will be Standard Access or Open Access. | Open Access (US\$4,200) |
| Please specify the section of the submitted manuscript. | Cancer Research |
| Please indicate the city, state/province, and country where this article will be filmed . Please do not use abbreviations. | Newark, NJ, 07103 |
| Please confirm that you have read and agree to the terms and conditions of the author license agreement that applies below: | I agree to the Author License Agreement |
| Please provide any comments to the journal here. | |
| Please indicate whether this article will be Standard Access or Open Access. | Open Access (\$3900) |
| Please confirm that you have read and agree to the terms and conditions of the video release that applies below: | I agree to the Video Release |

TITLE:

Natural Killer (NK) and CAR-NK Cell Expansion Method using Membrane Bound-IL-21-Modified B Cell Line

AUTHORS AND AFFILIATIONS:

Minh Ma^{1,2}, Saiaditya Badeti¹, James K. Kim¹, Dongfang Liu^{1,3,4}

¹Department of Pathology, Immunology and Laboratory Medicine, 185 South Orange Avenue, Newark, NJ 07103, USA.

²Department of Microbiology, Biochemistry & Molecular Genetics, Public Health Research Institute Center, New Jersey Medical School, Rutgers University, Newark, New Jersey, USA

³Center for Immunity and Inflammation, New Jersey Medical School, Rutgers-The State University of New Jersey, Newark, NJ 07103, USA

Email addresses of co-authors:

Minh Ma (mtm265@gsbs.rutgers.edu)

Saiaditya Badeti (sb1511@njms.rutgers.edu)

James Kim (jkk111@gsbs.rutgers.edu)

Corresponding authors:

Dongfang Liu (dongfang.liu@rutgers.edu)

SUMMARY

Here, we present a method to expand peripheral blood natural killer (PBNK), NK cells from liver tissues, and chimeric antigen receptor (CAR)-NK cells derived from peripheral blood mononuclear cells (PBMCs) or cord blood (CB). This protocol demonstrates the expansion of NK and CAR-NK cells using 221-mIL-21 feeder cells in addition to the optimized purity of expanded NK cells.

ABSTRACT

Chimeric antigen receptor (CAR)-modified immune cell therapy has become an emerging treatment for cancers and infectious diseases. NK-based immunotherapy, particularly CAR-NK cell, is one of the most promising 'off-the-shelf' development without severe life-threatening toxicity. However, the bottleneck for developing a successful CAR-NK therapy is achieving sufficient numbers of non-exhaustive, long-lived, 'off-the-shelf' CAR-NK cells from a third party. Here, we developed a new CAR-NK expansion method using an Epstein-Barr virus- (EBV) transformed B cell line expressing a genetically modified membrane form of interleukin-21 (IL-21). In this protocol, step-by-step procedures are provided to expand NK and CAR-NK cells from cord blood and peripheral blood, as well as solid organ tissues. This work will significantly enhance the clinical development of CAR-NK immunotherapy.

INTRODUCTION

Natural killer (NK) cells are an important subset of lymphocytes that express CD56 and lack expression of the T cell marker, CD3^{1,2}. Conventional NK cells are classified as innate immune cells responsible for immunosurveillance of virally infected cells and cancerous cells. Unlike T cells, NK cells recognize infected or malignant cells using CD16 or other activating receptors and do not require the formation of a complex between antigen peptides and major histocompatibility complex (MHC) class I molecules^{3,4}. Recent clinical investigations using chimeric antigen receptor (CAR)-NK cells to treat relapsed or refractory CD19-positive cancers (non-Hodgkin's lymphoma or chronic lymphocytic leukemia [CLL]) showed the outstanding safety advantages of CAR-NK cells⁵. For instance, CAR-NK cell infusion is associated with minimal or negligible graft versus host disease (GVHD), neurotoxicity, cardiotoxicity, and cytokine release syndrome (CRS)^{6–10}. However, conventional methods to expand human NK cells showed exhaustive phenotypes with strong fratricidal killing and telomere shortage, which presents a major challenge in obtaining an adequate number of functional NK cells for adoptive immunotherapy¹¹.

To overcome these challenges, a method was developed to expand primary NK cells directly from unfractionated peripheral blood mononuclear cells (PBMCs) or cord blood (CB) using an irradiated and genetically engineered 721.221 (hereinafter, 221) cell line, a human B-lymphoblastoid cell line with low expression of MHC class I molecules³. Previous studies showed the importance of IL-21 in NK cell expansion; therefore, a genetically engineered membrane-bound IL-21 expressing a version of the 721.221 cell line (starting now, 221-mIL-21) was developed^{11–15}. The results showed that 221-mIL-21 feeder-cell-expanded primary NK cells were expanded to an average of >40,000-fold with persistent high NK cell purity for approximately 2–3 weeks. Additional information regarding the application of this protocol can be found in Yang et al.¹⁶.

This protocol aims to demonstrate the step-by-step procedure of the novel expansion of PBNK, CBNK, tissue-derived NK, and CAR-NK cells *ex vivo*.

PROTOCOL

Human tissues and blood-related work in this protocol follows the guidelines of the Rutgers University Institutional Review Board (IRB)

1. NK cell expansion from liver tissues (Day 0), as shown in Figure 1.

NOTE: Initial cell number and viability are strongly correlated with the time since organ removal and the initial tissue sample amount. However, if tissues are placed in 30 mL of Hank's Balanced Salt Solution (HBSS) and kept on ice or in the fridge at 4 °C overnight, NK cells can still be expanded at high purity and viability up to 24 h later.

1.1. Identify viable tissue areas to obtain lymphocytes from tissues and sections using sterile surgical equipment.

1.2. Place tissues in 30 mL of HBSS (w/o calcium or magnesium) and keep on ice until ready to prepare for isolation.

1.3. Mince the tissue into <0.5 cm cubes using sterile razor blades or scissors and forceps inside a biosafety cabinet.

1.4. Prepare a 1x collagenase IV solution (1 mg/mL) by diluting a 10x stock in HBSS (10x Collagenase IV: 10 mg/mL or ~200 U/mL).

1.5. Place the minced tissue pieces in the tissue dissociator tubes. Fill the tubes with no more than 4 g of tissue and immerse the tissue pieces in ~10 mL of 1x collagenase IV.

NOTE: Using DNase I is not recommended as it may slightly decrease NK viability and yield. Please refer to the **Table of Materials** for the specific tissue dissociator tubes used.

1.6. Place the tissue dissociator tubes into a tissue dissociator and blend at 37 °C to mince the tissue thoroughly.

NOTE: For liver tissue, this may take over 30 min. For more friable tissue, around 15 min may be sufficient. Please refer to the **Table of Materials** for specific tissue dissociator tubes and the tissue dissociator used.

1.7. Remove the tissue dissociator tubes and triturate through 40 µm nylon cell strainer using the backend of a 5 mL syringe. Collect the eluent and discard the large undigested fragments.

1.8. Spin down the eluent at 400 x g for 5 min at room temperature. Aspirate the supernatant.

1.9. Resuspend the cell pellets in 30% polyvinylpyrrolidone (PVP)-coated silica to remove fat cells that will otherwise contaminate the final lymphocyte fraction.

1.9.1. To prepare 1x PVP-coated silica, use 9:1 dilution of PVP-coated silica in 10x PBS.

NOTE: Refer to the **Table of Materials** for the specific PVP-coated silica used.

1.9.2. To prepare 30% PVP-coated silica, dilute 1x PVP-coated silica with PBS/HBSS.

1.10. Spin down the cell pellet at 400 x g for 5 min at room temperature. Aspirate the supernatant.

1.11. Resuspend the cell pellet in 9 mL of R-10 media.

NOTE: Refer to the **Table of Materials** for the composition of R-10 media

1.12. Carefully layer the cell suspension over 4 mL of Ficoll or lymphocyte separation media to separate lymphocytes from red blood cells and polymorphonuclear cells.

1.13. Separate the layers by centrifuging at $400 \times g$ for 23 min at room temperature with the acceleration and brakes off or at the lowest setting. Carefully decant the upper-medium layer and harvest the interphase containing tissue-infiltrating lymphocytes.

1.14. Rinse the cells with 10 mL of media and proceed to cell counting, flow cytometry, aliquoting and freezing of cells, or primary NK expansion protocol.

2. Primary NK cell expansion from PBMCs (or CB or organ tissues) (Day 0), as shown in Figure 2.

2.1. Thaw the frozen PBMC and the frozen, irradiated feeder cells in a 37°C water bath.

2.2. Wash the PBMC and 100 Gy-irradiated (Gy-irradiated) 221-mIL-21 cells by centrifugation at $400 \times g$ for 5 min with 10 mL of R-10 media separately.

2.3. Save 1×10^6 cells of PBMC for flow cytometry.

NOTE: The initial NK cell purity is an important factor in calculating NK cell expansion rate.

2.4. Resuspend the cells in 1 mL of R-10 media. Count the cells using Trypan Blue.

2.5. Mix 5×10^6 cells of PBMC with 10×10^6 cells of 100 Gy-irradiated 221-mIL-21 cells in a special 6-well plate (see **Table of Materials**).

2.6. Add 30 mL of R-10 media supplemented with human IL-2 (200 U/mL) and human IL-15 (5 ng/mL) (see **Table of Materials**).

2.7. Incubate the special 6-well plate at 37°C with 5% CO_2 .

2.8. Replace the media with R-10 media supplemented with human IL-2 (200 U/mL) and human IL-15 (5 ng/mL) to maintain NK cells every 3–4 days.

NOTE: Maintain less than 20×10^6 cells per well for further expansion at each media change. For the best viability, ensure that the total cell number in each well does not exceed 100×10^6 cells.

2.9. Record the total cell number, viability and perform flow cytometry every 3–4 days to calculate the NK cell expansion rate.

3. Attachment of 293T cells (Day 2), as shown in Figure 2.

3.1. Split 1.8×10^6 293T cells in 11 mL of D-10 media per treated 100 mm plate.

NOTE: Refer to the **Table of Materials** for the composition of D-10 media

3.2. Incubate 293T cells at 37 °C with 5% CO₂.

4. Retrovirus transfection (Day 3)

4.1. In a 1.7 mL tube, mix 470 µL of reduced serum media with 30 µL of transfection reagent.

NOTE: Refer to the **Table of Materials** for the specific reduced serum media and transfection reagent used.

4.2. In a separate 1.7 mL tube, add 2.5 µg of pRDF plasmid, 3.75 µg of Pegpam3 plasmid, and 2.5 µg of CAR construct in SFG vector into the reduced serum media so that the final volume is 500 µL.

4.3. Mix the solutions in steps 4.1 and 4.2 dropwise.

4.4. Incubate the tube at room temperature for 15 min.

4.5. Add 1 mL of the mixture from step 4.4 to 293T cell plate on Day 1 in a dropwise manner.

4.6. Incubate the plate(s) at 37 °C with 5% CO₂ for 48–72 h.

5. Retronectin plate-coating (Day 3)

5.1. Dilute retronectin protein with Phosphate Buffered Saline (PBS) to a final concentration of 50–100 µg/mL.

5.2. Add 500 µL of the diluted retronectin into each well of an untreated 24-well plate (5 wells per CAR construct). Seal the plate using parafilm and incubate the plate at 4 °C overnight.

6. Transduction (Day 4)

6.1. Centrifuge the retronectin plate at 2103 x g for 30 min at 4 °C. Discard the supernatant.

6.2. Block each well of the 24-well plate with 1 mL of R-10 medium.

6.3. Incubate the plate at 37 °C with 5% CO₂ for 1 h.

6.4. Pre-warm the centrifuge to 32 °C while the retronectin plate is being blocked.

6.5. Collect the retrovirus supernatant by filtering the transfected 293T cells using a 0.45 µm filter.

216
217 6.6. Aliquot 2 mL of the filtered retrovirus supernatant into each well.

218
219 6.7. Centrifuge the 24-well plate at $2103 \times g$ for 2 h at 32 °C.

220
221 6.8. During plate centrifugation, collect the expanded PBNK cells from Day 0 and count the
222 cells using Trypan Blue.

223
224 NOTE: Continue expanding PBNK cells by adding R-10 media supplemented with IL-2 (200 U/mL)
225 and IL-15 (5 ng/mL).

226
227 6.9. Dilute the expanded PBNK cells with R-10 media supplemented with IL-2 (200 U/mL) and
228 IL-15 (5 ng/mL) to 2.5×10^5 – 5×10^5 cells/mL (0.5×10^6 – 1×10^6 cells per well).

229
230 NOTE: Record the total cell number, viability, and save 5×10^5 expanded PBNK cells for flow
231 cytometry as these values are important in determining the NK cell expansion rate.

232
233 6.10. After centrifugation, partially aspirate the retrovirus supernatant from each well.

234
235 NOTE: Do not completely aspirate, i.e., leave approximately 100 µL of retrovirus supernatant per
236 well, as this will decrease the transduction efficiency.

237
238 6.11. Aliquot 2 mL of the diluted expanded PBNK cells from step 6.8 to each well.

239
240 6.12. Centrifuge the plate at $600 \times g$ for 10 min at 32 °C. Incubate the plate at 37 °C with 5%
241 CO₂ for 48–72 h.

242
243 NOTE: Do not parafilm the plate.

244 245 **7. CAR-NK cells collection (Day 6 or 7), as shown in Figure 2.**

246
247 7.1. Gently collect the cells from the 24-well plate and transfer the cells to a 50 mL centrifuge
248 tube

249
250 NOTE: Try not to generate bubbles, as this will result in a decrease in cell viability.

251
252 7.2. Centrifuge the tube at $400 \times g$ for 5 min.

253
254 7.3. Resuspend the pellet with 1 mL of R-10 media and count the cells using Trypan Blue.

255
256 NOTE: Save 5×10^5 cells for flow cytometry to determine the transduction efficiency.

257
258 7.4. Transfer the resuspended cells to a special 6-well plate containing 30 mL of R-10 media
259 supplemented with IL-2 (200 U/mL) and IL-15 (5 ng/mL).

7.5. Incubate the special 6-well plate at 37 °C with 5% CO₂.

7.6. Replace R-10 media supplemented with IL-2 (200 U/mL) and IL-15 (5 ng/mL) to maintain NK cells every 3 – 4 days.

NOTE: Maintain less than 20 x 10⁶ cells per well for further expansion at each change. For the best viability, ensure that the total cell number in each well does not exceed 100 x 10⁶ cells.

7.7. Record the total cell number, viability, and perform flow cytometry every 3–4 days to calculate the NK cell expansion rate.

7.8. Use the cells for appropriate *in vitro* or *in vivo* assays.

NOTE: The *ex vivo* expanded PBNK, and CAR-NK cells can be cultured in a 37 °C incubator for approximately 4 weeks.

7.9. Examine the NK cell number and purity at day 7, day 11, day 14, day 18, and day 21 by flow cytometry.

REPRESENTATIVE RESULTS

A schematic workflow of tissue-infiltrating NK cell isolation and PBNK cell expansion using the 221-mIL-21 feeder cell methodology is shown in **Figure 1** and **Figure 2**. Expanded PBNK cells were collected every 3 or 4 days for flow cytometry to determine the NK cell purity by staining cells with anti-human CD56 and anti-human CD3. The experiment was repeated using two different donors to show the reproducibility of the expansion system (**Figure 3**). PBNK cells expanded by 221-mIL-21 were shown to expand nearly 5 × 10⁴ folds (**Figure 3A**). Furthermore, the NK cell purity was highly maintained, around 85% throughout the 21-day expansion (**Figure 3B**). Using the 221-mIL-21 feeder cell expansion system, the NK cell purity consistently ranged between 85%–95%, independent of the donors (data not shown). To demonstrate the robustness of the 221-mIL-21 expansion system, PBMCs were stained for anti-CD56 and anti-CD3 prior to the expansion, which showed a cell purity of 7.09% for NK cells and a high percentage of T cells (**Figure 4A**). PBMCs were cocultured with 221-mIL-21 to expand NK cells; the NK purity was checked prior to CAR-NK transduction on Day 4 (**Figure 4A**). CAR-NK cells were collected and stained for anti-CD56, anti-CD3, and anti-hIgG(H+L) F(ab')₂, which showed a high NK cell population (86.9% on Day 7) and a high CAR transduction efficiency of approximately 70% (**Figure 4**). Higher transduction efficiencies (up to 95%) were also observed using the retrovirus packaging system. Altogether, these data show that the 221-mIL-21 feeder cells could successfully expand NK cells and preserve the NK cell purity *ex vivo*.

FIGURE AND TABLE LEGENDS

Figure 1: Diagram of NK cell expansion from solid human organ samples. Briefly, the obtained human liver samples are minced into small cubes for mechanical digestion. Dissociated cells are

then isolated using PVP-coated silica and Lymphocyte Separation Media. Further, the NK cells are expanded using the expansion protocol described in **Figure 2**.

Figure 2: Schematic workflow of CAR-NK cell generation from PBMCs. Briefly, 221-mIL21 feeder cells were irradiated at 100 Gy prior to coculturing with PBMCs supplemented with IL-2 and IL-15 on Day 0. In parallel, 293T cells were transfected with the retrovirus packaging system to produce CAR retrovirus that was then transduced into the expanded PBNK cells in the presence of IL-2 and IL-15. Primary CAR-NK cells were harvested on Day 7 and continued expansion for 21 days. This figure has been modified from Yang et al.¹⁶.

Figure 3: Dynamic time-lapsed expansion of PBNK cells. (A) Fold expansion of PBNK during a 22-day time course. Cells were stained with anti-CD56 and anti-CD3 at indicated days for flow cytometry. The total number of NK cells was determined by multiplying NK cell purity to the total number of cells. Expansion rate was generated as follows: $(\text{Number of NK cells})_{T_n} / (\text{Number of NK cells})_{T_0}$, where Number of NK cells = (percentage of NK cell purity) \times (total number of cells), T_0 is the number of NK cells at time day 0, and T_n is the number of NK cells at time day n . **(B)** NK cell purity during a 22-day time course. The NK cell expansion was repeated two times with two different donors. Error bars represent \pm SEM.

Figure 4: Representative flow cytometric analysis of CAR-NK cells. (A) Representative dot plots showing the dynamic time-lapse of NK cell purity of the CAR-NK cells during an 18-day course. Flow cytometry analysis was assessed by staining the cells with anti-CD56 and anti-CD3 at indicated timepoints. Day 0 indicates pre-expansion of PBNK. Day 4 indicates post-expansion of PBNK and pre-transduction of CAR-NK cells. Day 7 indicates the post-transduction of CAR-NK cells. **(B)** Representative dot plots showing the transduction efficiency of CAR-NK cells using the retrovirus packaging system. Cells were stained with anti-CD56, anti-CD3, and anti-hIgG(H+L) F(ab')₂ for flow cytometry. **(C)** Representative dot plots showing the CAR expression in various subsets, including CD56⁺CD3⁻, CD56⁻CD3⁺, CD56⁺CD3⁺, and CD56⁻CD3⁻ on Day 18. Cells were stained with anti-CD56, anti-CD3, and anti-hIgG(H+L) F(ab')₂ (indicating CAR expression) for flow cytometry.

DISCUSSION

Most of the current CAR-NK products in clinical trials utilize NK cell lines¹⁷, such as NK-92, a cell line isolated from a non-Hodgkin's lymphoma patient¹⁸, NK-92MI, IL-2 independent NK-92 cell line¹⁹, and NKL, isolated from a large granular lymphocyte patient²⁰, as these cell lines are easily proliferative for 'off-the-shelf' products. However, these cell lines, e.g., NK-92 cells, have marginal clinical efficacies and *in vivo* expansion, as they require irradiation prior to infusion, thus limiting their proliferation and cytotoxicity *in vivo*²¹. Given these reasons, various strategies are currently being explored to expand primary NK cells from several sources, including peripheral blood, CB, bone marrow (BM), human embryonic stem cells (HSCs), induced pluripotent stem cells (iPSCs), and tumor tissues²¹⁻²³. For instance, NK cells can be expanded *ex vivo* using interleukins including IL-15, IL-18, and IL-21. Lymphoblastoid cell lines such as K562 cells or Epstein-Barr Virus-transformed lymphoblastoid cell lines such as 721.221 cells, are also used for NK cell expansion¹⁶. However, the aforementioned strategies often generate insufficient number

of NK cells for an adoptive transfer of CAR-NK immunotherapy^{22,24}. To help solve the problem, the study here shows a protocol for an *ex vivo* NK cell expansion using a genetically modified EBV-transformed cell line, 221-mIL-21 feeder cells.

The expansion methodology using 221-mIL-21 feeder cells shown in this protocol is optimized to expand NK cells with an expansion rate of at least 10 to 100-fold higher than other leukemia cell lines, including HL-60 and OCI-AML3 expressing membrane IL-21, K562, and K652-mIL21 expressing OX40 ligand^{22,24,25}. The CAR expression is also evaluated for approximately 2 weeks *ex vivo*. More significantly, the 221-mIL-21 feeder cell expansion strategy can be applied to expand NK cells from various sources, including PBMCs, CB, and solid organs such as the liver, without an initial NK enrichment step. Although the 221-mIL-21 feeder system is not as donor-dependent as the aforementioned feeder cell lines, it is not entirely independent of donors. On average, the 221-mIL-21 expansion system can achieve 90% of NK cell purity with a high NK cell number, with approximately <5% of T cell contamination on day 14 post-expansion. Therefore, to eliminate the possibilities of T cell contamination, it is necessary to isolate NK cells from obtained samples prior to the *ex vivo* expansion or use a CD3⁺ selection system to eliminate T cells after an *ex vivo* expansion.

One of the criticisms in using an NK cell expansion system is that the feeder cells may not have been fully eradicated after the expansion or prior to a transfusion, which may possess significant regulatory concerns; therefore, complete eradication of feeder cells before a transfusion is crucial. However, recent CAR-NK clinical trials in which K562-mIL21-4-1BBL feeder cells were used for the *ex vivo* CBNK cell expansion^{24,25} showed no concerning complications. Furthermore, our preliminary data showed a gradual decrease of the irradiated 221-mIL-21 population as the expansion progressed (data not shown). However, more extensive studies are required for this expansion method to be implemented in a clinical setting. Collectively, the 221-mIL-21 expansion system helps solve the challenge of expanding primary CAR-NK cells, and therefore will significantly contribute to the broader use of CAR-NK cell-based immunotherapy in the near future.

ACKNOWLEDGMENTS

We would like to thank the members of the Liu laboratory (Dr. Hsiang-chi Tseng, Dr. Xuening Wang, and Dr. Chih-Hsiung Chen) for their comments on the manuscripts. We would like to thank Dr. Gianpietro Dotti for the SFG vectors and Dr. Eric Long for the 721.221 cells. This work was supported in part from HL125018 (D. Liu), AI124769 (D. Liu), AI129594 (D. Liu), AI130197 (D. Liu), and Rutgers-Health Advance Funding (NIH REACH program), U01HL150852 (R. Panettieri, S. Libutti, and R. Pasqualini), S10OD025182 (D. Liu), and Rutgers University-New Jersey Medical School Startup funding for D. Liu Laboratory.

DISCLOSURES


The authors declare no competing interests.


REFERENCES


1. Van Acker, H. H. et al. CD56 in the immune system: More than a marker for cytotoxicity? *Frontiers in Immunology*. **8**, 892 (2017)
2. Caligiuri, M. A. Human natural killer cells. *Blood*. **112** (3), 461–469 (2008).
3. Shimizu, Y. et al. Transfer and expression of three cloned human non-HLA-A,B,C class I major histocompatibility complex genes in mutant lymphoblastoid cells. *Proceedings of the National Academy of Sciences of the United States of America*. **85** (1), 227–231 (1988).
4. Wu, J., Lanier, L. L. Natural killer cells and cancer. *Advances in Cancer Research*. **90**, 127–156 (2003)
5. Liu, E. et al. Use of CAR-transduced natural killer cells in CD19-positive lymphoid tumors. *The New England Journal of Medicine*. **382** (6), 545–553 (2020).
6. Alonso-Camino, V. et al. Efficacy and toxicity management of CAR-T-cell immunotherapy: a matter of responsiveness control or tumour-specificity? *Biochemical Society Transactions*. **44** (2), 406–411 (2016).
7. Bonifant, C. L., Jackson, H. J., Brentjens, R. J., Curran, K. J. Toxicity and management in CAR T-cell therapy. *Molecular Therapy Oncolytics*. **3**, 16011 (2016).
8. Kalaitidou, M., Kueberuwa, G., Schitt, A., Gilham, D. E. CAR T-cell therapy: toxicity and the relevance of preclinical models. *Immunotherapy*. **7** (5), 487–97 (2015).
9. Gust, J. et al. Endothelial activation and blood-brain barrier disruption in neurotoxicity after adoptive immunotherapy with CD19 CAR-T cells. *Cancer Discovery*. **7** (12), 1404–1419 (2017).
10. Hay, K. A. et al. Kinetics and biomarkers of severe cytokine release syndrome after CD19 chimeric antigen receptor-modified T-cell therapy. *Blood*. **130** (21), 2295–2306 (2017).
11. Zhang, Y. et al. In vivo kinetics of human natural killer cells: the effects of ageing and acute and chronic viral infection. *Immunology*. **121** (2), 258–265 (2007).
12. Vidard, L. et al. CD137 (4-1BB) Engagement fine-tunes synergistic IL-15- and IL-21-driven nk cell proliferation. *Journal of Immunology (Baltimore, Md.:1950)*. **203** (3), 676–685 (2019).
13. Venkatasubramanian, S. et al. IL-21-dependent expansion of memory-like NK cells enhances protective immune responses against *Mycobacterium tuberculosis*. *Mucosal Immunology*. **10** (4), 1031–1042 (2017).
14. Ojo, E. O. et al. Membrane bound IL-21 based NK cell feeder cells drive robust expansion and metabolic activation of NK cells. *Scientific Reports*. **9** (1), 14916 (2019).
15. Denman, C. J. et al. Membrane-bound IL-21 promotes sustained ex vivo proliferation of human natural killer cells. *PLoS One*. **7** (1), e30264 (2012).
16. Yang, Y. et al. Superior expansion and cytotoxicity of human primary NK and CAR-NK cells from various sources via enriched metabolic pathways. *Molecular Therapy. Methods & Clinical Development*. **18**, 428–445 (2020).
17. Liu, S. et al. NK cell-based cancer immunotherapy: from basic biology to clinical development. *Journal of Hematology & Oncology*. **14** (1), 7 (2021).
18. Gong, J. H., G. Maki, Klingemann, H. G. Characterization of a human cell-line (Nk-92) with phenotypical and functional-characteristics of activated natural-killer-cells. *Leukemia*. **8** (4), 652–658 (1994).
19. Tam, Y. K. et al. Characterization of genetically altered, interleukin 2-independent natural killer cell lines suitable for adoptive cellular immunotherapy. *Human Gene Therapy*. **10** (8), 1359–1373 (1999).

- 434 20. Robertson, M. J. et al. Characterization of a cell line, NKL, derived from an aggressive
435 human natural killer cell leukemia. *Experimental Hematology*. **24** (3), 406–415 (1996).
- 436 21. Hu, Y., Tian, Z. G., Zhang, C. Chimeric antigen receptor (CAR)-transduced natural killer cells
437 in tumor immunotherapy. *Acta Pharmacologica Sinica*. **39** (2), 167–176 (2018).
- 438 22. Tseng, H. C. et al. Efficacy of anti-CD147 chimeric antigen receptors targeting
439 hepatocellular carcinoma. *Nature Communications*. **11** (1), 4810 (2020).
- 440 23. Easom, N. J. W. et al. IL-15 overcomes hepatocellular carcinoma-induced NK cell
441 dysfunction. *Frontiers in Immunology*. **9**, 1009 (2018).
- 442 24. Granzin, M. et al. Highly efficient IL-21 and feeder cell-driven ex vivo expansion of human
443 NK cells with therapeutic activity in a xenograft mouse model of melanoma. *Oncoimmunology*. **5**
444 (9), e1219007 (2016).
- 445 25. Liu, E. L. et al. Cord blood derived natural killer cells engineered with a chimeric antigen
446 receptor targeting CD19 and expressing IL-15 have long term persistence and exert potent anti-
447 leukemia activity. *Blood*. **126** (23) (2015).
- 448


Figure 1









A



B



C





[Click here to access/download](#)

Table of Materials

Materials list_03-09-2021_FINAL.xlsx



Editorial comments:

1. Line 39-42: Please move the lines to the introduction and cite the appropriate reference number.

Response: Lines 39-42 have been moved to the introduction and cited appropriately (highlighted in teal).

2. Please note that some of the references are duplicated (18 and 27, 6 and 30). Please replace/remove one of the references.

Response: We apologize for our carelessness. The citations have been removed.

3. Please ensure that the responses (Reviewer #3) in the rebuttal are included in the manuscript.

Response: Thank you for the editor's reminder. The questions asked by reviewer #3 had been addressed in the last 2 paragraphs in the discussion (now highlighted in teal).

Reviewers' comments:

Reviewer #4:

Major Concerns:

In my opinion this protocol does not give any additional value to the known results.

Response: Thank you for the reviewer's comment. We believe that this protocol will help advance the field of adoptive NK cell transfer therapy and CAR-NK therapy as follows:

1. The expansion system used by other groups require an additional NK cell isolation step[1]. As NK cells only represent 5-10% of the total lymphocytes[2], this step will result in a significant loss in starting NK cells, and increase the cost of manufacturing process, which results in lower expansion rate with high cost. Our current system **does not** require NK cell isolation step, which (1) simplifies the protocol and (2) reserves the limited amount of NK cells.
2. We and other groups (industry and academia) have compared our 221-mIL21 feeder cells with the currently well-known K562-mIL21-41BBL feeder cells and the results showed that the 221-mIL21 feeder cells out-perform K562-mIL21-41BBL feeder cells in both NK cell purity and expansion rate from in vitro and in vivo studies. This will significantly contribute to the field in making adoptive NK cell transfer therapy and CAR-NK therapy more affordable.
3. K562 is used to expand NK cells ex vivo because it does not express HLA molecules[1, 3]. However, the interaction between HLA molecules and NK cell receptors is necessary for NK cell "licensing", i.e., training NK cells to appropriately activate[4]. Therefore, NK cells expanded with K562 feeder cells may remain unlicensed and could result in hyporesponsiveness[5]. 221 cell line, however, expresses low amount of HLA ligands for

NK cell “licensing”[6], and therefore may prevent activation-induced NK cell death, as well as licensed NK and CAR-NK cells.

4. A current bottleneck in the field is that expanded NK cells by cytokine addition are phenotypically exhausted[7]. By using the 221-mIL21 feeder cell system, we not only demonstrated that the expanded NK cells are not exhausted but also have higher cytolytic functions compared to that of NK cells expanded by K562-mIL21[8].
5. To our knowledge, we are the first group to have demonstrated the antitumor responses of expanded NK cells *in vivo* using different feeder cell lines, and the results showed that NK cells expanded with 221-mIL21 have better anti-tumor responses compared to that of NK cells expanded with K562-mIL21[8].

Collectively, the expansion system described in this protocol adds substantial value to the adoptive NK cell transfer and CAR-NK therapies by providing (1) more simplified protocol, (2) robust and reproducible protocol, (3) higher NK cell purity and expansion rate, (4) less exhaustive NK cells, (5) more affordable adoptive NK cell transfer and CAR-NK therapies in the near future.

Reviewer #5:

Manuscript Summary:

Ma et al accurately explain the protocol of expanding CAR NK cells from liver tissue. The protocol is nice, clean and easy.

Major Concerns:

None

Minor Concerns:

None

Reviewer #6:

Manuscript Summary:

The protocol to expand NK and CAR-NK using 221-mil-21 feeder cells is described clearly and in details.

Major Concerns:

None

Minor Concerns:

please make sure of the formular listed in line 325. I think the expansion rate should be

calculated by dividing the number of NK cells after T days of culture by the respective number on day 0.

Response: Thank you for the reviewer's comment. We have changed the equation in the main text (highlighted in teal).

References

1. Liu EL, Tong YJ, Dotti G, Savoldo B, Muftuoglu M, Kondo K, et al. Cord Blood Derived Natural Killer Cells Engineered with a Chimeric Antigen Receptor Targeting CD19 and Expressing IL-15 Have Long Term Persistence and Exert Potent Anti-Leukemia Activity. *Blood*. 2015;126(23). doi: DOI 10.1182/blood.V126.23.3091.3091. PubMed PMID: WOS:000368020103286.
2. Koehl U, Brehm C, Huenecke S, Zimmermann SY, Kloess S, Bremm M, et al. Clinical grade purification and expansion of NK cell products for an optimized manufacturing protocol. *Front Oncol*. 2013;3:118. Epub 2013/06/05. doi: 10.3389/fonc.2013.00118. PubMed PMID: 23730623; PubMed Central PMCID: PMC3656406.
3. Denman CJ, Senyukov VV, Somanchi SS, Phatarpekar PV, Kopp LM, Johnson JL, et al. Membrane-Bound IL-21 Promotes Sustained Ex Vivo Proliferation of Human Natural Killer Cells. *Plos One*. 2012;7(1). doi: ARTN e30264 10.1371/journal.pone.0030264. PubMed PMID: WOS:000299771900053.
4. Kim S, Poursine-Laurent J, Truscott SM, Lybarger L, Song YJ, Yang L, et al. Licensing of natural killer cells by host major histocompatibility complex class I molecules. *Nature*. 2005;436(7051):709-13. Epub 2005/08/05. doi: 10.1038/nature03847. PubMed PMID: 16079848.
5. Anfossi N, Andre P, Guia S, Falk CS, Roetynck S, Stewart CA, et al. Human NK cell education by inhibitory receptors for MHC class I. *Immunity*. 2006;25(2):331-42. Epub 2006/08/12. doi: 10.1016/j.immuni.2006.06.013. PubMed PMID: 16901727.
6. Shimizu Y, Geraghty DE, Koller BH, Orr HT, DeMars R. Transfer and expression of three cloned human non-HLA-A,B,C class I major histocompatibility complex genes in mutant lymphoblastoid cells. *Proc Natl Acad Sci U S A*. 1988;85(1):227-31. Epub 1988/01/01. doi: 10.1073/pnas.85.1.227. PubMed PMID: 3257565; PubMed Central PMCID: PMC279517.
7. Keir ME, Butte MJ, Freeman GJ, Sharpe AH. PD-1 and its ligands in tolerance and immunity. *Annu Rev Immunol*. 2008;26:677-704. Epub 2008/01/05. doi: 10.1146/annurev.immunol.26.021607.090331. PubMed PMID: 18173375.
8. Yang Y, Badeti S, Tseng HC, Ma MT, Liu T, Jiang JG, et al. Superior Expansion and Cytotoxicity of Human Primary NK and CAR-NK Cells from Various Sources via Enriched Metabolic Pathways. *Mol Ther-Meth Clin D*. 2020;18:428-45. doi: 10.1016/j.omtm.2020.06.014. PubMed PMID: WOS:000568405100041.

The recent publication by Yang et. al. is an open access article under the CC BYNC-ND license (<http://creativecommons.org/licenses/by-nc-nd/4.0/>). We also cited the manuscript accordingly. In the Figure legend, we also added the sentence –“ This figure has been modified from (citation).”