

Journal of Visualized Experiments

A Rapid, Multiplex Dual Reporter IgG and IgM SARS-CoV-2 Neutralization Assay for a Multiplexed Bead-based Flow Analysis System --Manuscript Draft--

Article Type:	Invited Methods Article - JoVE Produced Video
Manuscript Number:	JoVE62487R2
Full Title:	A Rapid, Multiplex Dual Reporter IgG and IgM SARS-CoV-2 Neutralization Assay for a Multiplexed Bead-based Flow Analysis System
Corresponding Author:	Sherry Dunbar Luminex Corp Austin, Texas UNITED STATES
Corresponding Author's Institution:	Luminex Corp
Corresponding Author E-Mail:	sdunbar@luminexcorp.com
Order of Authors:	Stephen Angeloni Andrew Cameron Nicole Pecora Sherry Dunbar
Additional Information:	
Question	Response
Please indicate whether this article will be Standard Access or Open Access.	Open Access (US\$4,200)
Please specify the section of the submitted manuscript.	Immunology and Infection
Please indicate the city, state/province, and country where this article will be filmed . Please do not use abbreviations.	Austin, Texas, USA
Please confirm that you have read and agree to the terms and conditions of the author license agreement that applies below:	I agree to the Author License Agreement
Please provide any comments to the journal here.	Our editor is Lyndsay Troyer
Please indicate whether this article will be Standard Access or Open Access.	Open Access (\$3900)

TITLE:

A Rapid, Multiplex Dual Reporter IgG and IgM SARS-CoV-2 Neutralization Assay for a Multiplexed Bead-based Flow Analysis System

AUTHORS AND AFFILIATIONS:

Stephen Angeloni¹, Andrew Cameron², Nicole D. Pecora^{2,3}, Sherry Dunbar¹

¹Luminex Corporation, Austin, Texas, USA

²Departments of Pathology and Laboratory Medicine, Clinical Microbiology, University of Rochester Medical Center, Rochester, New York, USA

³Department of Microbiology and Immunology, University of Rochester Medical Center, Rochester, New York, USA

Email address of corresponding Author:

Sherry Dunbar (sdunbar@luminexcorp.com)

E-mail addresses of Co-Authors:

Stephen Angeloni (sangeloni@luminexcorp.com)

Andrew Cameron (Andrew_Cameron@URMC.Rochester.edu)

Nicole D. Pecora (Nicole_Pecora@URMC.Rochester.edu)

Sherry Dunbar (sdunbar@luminexcorp.com)

KEYWORDS:

multiplex, microsphere, bead, SARS-CoV-2, serological assay, neutralization assay

SUMMARY:

A flow analysis system for bead-based multiplexed assays which provides a two-reporter readout was used for the development of multiplex serological and antibody neutralization assays that can simultaneously measure neutralizing IgG and IgM antibodies for SARS-CoV-2.

ABSTRACT:

The COVID-19 pandemic has underscored the need for rapid high-throughput methods for sensitive and specific serological detection of infection with novel pathogens, such as SARS-CoV-2. Multiplex serological testing can be particularly useful because it can simultaneously analyze antibodies to multiple antigens that optimizes pathogen coverage, and controls for variability in the organism and the individual host response. Here we describe a SARS-CoV-2 IgG 3-plex fluorescent microsphere-based assay that can detect both IgM and IgG antibodies to three major SARS-CoV-2 antigens—the spike (S) protein, spike angiotensin-converting enzyme-2 (ACE2) receptor-binding domain (RBD), and nucleocapsid (Nc). The assay was shown to have comparable performance to a SARS-CoV-2 reference assay for IgG in serum obtained at ≥ 21 days from symptom onset but had higher sensitivity with samples collected at ≤ 5 days from symptom onset. Further, using soluble ACE2 in a neutralization assay format, inhibition of antibody binding was demonstrated for S and RBD.

INTRODUCTION:

The COVID-19 pandemic has emphasized the importance of being able to rapidly develop and implement fast, high-throughput, high-performing serological tests that can be used for diagnosis and surveillance of novel pathogens such as SARS-CoV-2. Multiplex serological testing can be extremely useful in a pandemic, as it can simultaneously analyze antibodies to multiple antigens, which provides for broad pathogen coverage and controls for variability in the organism and the host response to infection. The development of a multiplex fluorescent bead-based assay, the SARS-CoV-2 IgG 3Flex (3Flex), that can detect antibodies to the spike (S) protein, the spike angiotensin-converting enzyme-2 (ACE2) receptor-binding domain (RBD), and nucleocapsid (Nc)1 of SARS-CoV-2 was recently reported¹. 3Flex was shown to have comparable performance to a reference SARS-CoV-2 IgG assay for serum obtained at ≥ 21 days from symptom onset but had higher sensitivity with samples collected at ≤ 5 days from symptom onset. Additionally, inhibition of antibody binding was demonstrated for S and RBD antigens using soluble ACE2 in a neutralization assay format. Multiplex bead-based technology has been widely adopted as a proven technology for multiplex serological testing and has distinct advantages for large-scale screening, including short time-to-results, small sample requirements, and reduced labor²⁻⁴.

Recently, a new flow analysis instrument has become available that provides the same high-plex, high-throughput capability of the system demonstrated in previous work¹, but with the added benefit of two reporter channels. The ability to measure two fluorescent reporter signals in a single reaction allows the assessment of two results per analyte for each sample and is ideal for isotyping applications, which typically must be performed in separate reaction wells⁵. The dual reporter capability provides twice the data for each sample in half as many reaction wells with half the sample volume requirements, and thus has a clear advantage over single reporter systems. This study details the standard serological assay protocol and describes the adaptation to a neutralization assay. Additionally, this report describes the conversion of the single reporter assay to a dual reporter serological assay, and subsequently, a dual reporter neutralization assay, to simultaneously measure neutralizing antibodies of both IgG and IgM isotypes against the SARS-CoV-2 S, RBD, and Nc antigens. A visual overview of the assay protocol is provided in **Figure 1**.

PROTOCOL:

Residual human serum samples previously tested for SARS-CoV-2 diagnostic analyses were used in this study which was approved by the University of Rochester Institutional Review Board.

1. Reagents and equipment

1.1 Obtain coronavirus antigens and human IgG and IgM from commercial sources.

1.2 Obtain assay control (AC) bead sets (45 and 220) from the instrument vendor. The AC bead set 45 monitors Median Fluorescence Intensity (MFI) collection on single reporter channel (RP1) instruments. The AC bead set 220 monitors MFI collection on the second reporter channel (RP2).

1.3 Obtain detection antibodies and antigen-specific rabbit sera and antibodies used for coupling confirmation.

1.4 Perform the detection with biotin-labeled antibodies with either a streptavidin, R-phycoerythrin conjugate (SAPE) or a streptavidin conjugate of fluorophore Super Bright 436 (SASB; excitation at 405 nanometers (nm) and emission at 436 nm).

1.5 Obtain Human ACE2 protein (i.e., the SARS Receptor). Prior to use, rehydrate lyophilized ACE2 protein by dissolving to a 1 mg/mL concentration in nuclease-free water (Not DEPC-treated, autoclaved and 0.2- μ m sterile filtered) and store at 4 °C.

1.6 Perform all reactions in 96-well nonbinding surface (NBS) black flat-bottom plates. Seal the plates with 96 well microplate aluminum sealing tape for assay incubation steps. Use a magnetic plate separator for the assay wash steps.

2. Antigen-coupled and antibody-coupled control beads preparation

2.1 Covalently couple magnetic, polystyrene, and color-coded microsphere sets (6.5- μ m diameter) with S, RBD, and Nc antigens as described in Cameron et al.¹. Coronavirus antigens are coupled at 10 pmol/ 10^6 beads (S) and 100 pmol/ 10^6 beads (RBD and Nc). Control beads coupled with human IgG or IgM are generated as previously described¹. IgG and IgM are coupled at 5 μ g/ 10^6 beads.

2.2 After coupling, determine the bead concentrations and dilute each stock to 1×10^6 beads/mL using described procedures⁶.

2.3 Perform confirmation of antigen coupling with antigen-specific antibodies from rabbit or with positive human sera as described in Cameron et al.¹.

2.3.1 Confirm human IgG coupling with the assay's biotinylated anti-human IgG/SAPE detection reagents, and IgM coupling with a phycoerythrin (PE)-labeled goat anti-human IgM.

NOTE: Coupling confirmation is performed using a titration of protein-specific serum or antibody, because the optimal concentration to use for these reagents is not known. For serum, serial fold dilutions are used whereas for antibodies supplied as mg/mL stocks, serial μ g/mL dilution series are used.

2.4 Dilute assay control (AC) beads that monitor MFI collection in the two reporter channels to a concentration of 1×10^6 beads/mL prior to using in multiplex bead mixes.

3. Fresh multiplex bead mix preparation

3.1 Prepare multiplex beads mixes freshly each day from the individual coupled bead and control bead stocks at 1×10^6 beads/mL.

NOTE: Multiplex bead mixes are prepared to deliver 2,500 beads for each bead region in a 50- μ L volume. Calculations for how to prepare multiplex bead mixes for any number of reactions are described in the blog published by Angeloni (2020)⁶.

4. Sample dilution

4.1 Prepare a 1:100 serum sample dilution by mixing 10 μ L of serum into 990 μ L of phosphate-buffered saline (PBS) containing 0.05% Tween 20, 0.02% sodium azide, and 1% of bovine serum albumin (BSA) (PBS-TBN buffer).

4.2 Dilute samples again 10-fold by adding 20 μ L of the 1:100 dilution to 180 μ L of PBS-TBN.

NOTE: Serum samples consist of negative sera, from specimens collected in 2019 prior to the COVID-19 pandemic, and positive sera from SARS-CoV-2 PCR-positive patients, representing low and high antibody titers. Positive samples were collected between \approx 3 to >60 days from symptom onset (DFSO). All serum samples and IgG-stripped negative control sera are diluted to 1:1000 as described below. For neutralization assays, sera are diluted 1:500 as described in the neutralization protocols below (sections 6 and 8).

5. Single reporter serological assay protocol

5.1 Add 50 μ L of multiplex bead mix to each assigned well of a 96-well non-binding microtiter plate.

5.2 Pipette 50 μ L of 1:1000 serum or PBS-TBN into the appropriate wells; the final dilution of serum in the well is 1:2000.

5.3 Cover the plate with a microplate foil seal and incubate on a heated plate shaker at 37 $^{\circ}$ C for 15 min with shaking at 600 rpm.

5.4 Remove the plate from the plate shaker and place onto the magnetic plate separator for 2 min to separate the beads from the reaction mixture.

5.5 While retaining the plate on the magnet, remove the foil seal, carefully invert the plate over a waste container, and gently flick the supernatant out of the wells.

5.6 Blot the plate on absorbent paper while still holding the plate on the magnet.

5.7 Wash the reaction wells as described below (first post-sample incubation wash).

5.7.1 Remove the plate from the plate magnet and add 150 μ L of PBS-TBN to each well.

5.7.2 Cover the plate with a fresh foil seal and shake at 37 $^{\circ}$ C for 2 min at 600 rpm.

177
178 5.7.3 Remove the plate from the plate shaker and place onto the plate magnet for 2 min to
179 separate the beads from the reaction mixture.
180
181 5.7.4 While retaining the plate on the magnet, remove the foil seal, carefully invert the plate
182 over a waste container, and gently flick the supernatant out of the wells.
183
184 5.7.5 Blot the plate on absorbent paper while still holding the plate on the magnet.
185
186 5.8 Wash the reaction wells as described below (second post-sample incubation wash).
187
188 5.8.1 Remove the plate from the plate magnet and add 150 μ L of PBS-TBN to each well.
189
190 5.8.2 Cover the plate with a fresh foil seal and shake at 37 °C for 2 min at 600 rpm.
191
192 5.8.3 Remove the plate from the plate shaker and place onto the plate magnet for 2 min to
193 separate the beads from the reaction mixture.
194
195 5.8.4 While retaining the plate on the magnet, remove the foil seal, carefully invert the plate
196 over a waste container, and gently flick the supernatant out of the wells.
197
198 5.8.5 Remove the plate from the magnet.
199
200 5.9 Make a fresh mix of biotin-anti-human IgG with SAPE by diluting the antibody to 0.62
201 μ g/mL and the SAPE to 1 μ g/mL. Add 100 μ L to each well.
202
203 5.10 Cover the plate with a microplate foil seal and incubate on a heated plate shaker at 37 °C
204 for 15 min with shaking at 600 rpm.
205
206 5.11 Remove the plate from the plate shaker and place onto the magnetic plate separator for
207 2 min to separate the beads from the reaction mixture.
208
209 5.12 While retaining the plate on the magnet, remove the foil seal, carefully invert the plate
210 over a waste container, and gently flick the supernatant out of the wells.
211
212 5.13 Blot the plate on absorbent paper while still holding the plate on the magnet.
213
214 5.14 Wash the reaction wells as described below (first post detection antibody incubation
215 wash).
216
217 5.14.1 Remove the plate from the plate magnet and add 150 μ L of PBS-TBN to each well.
218
219 5.14.2 Cover the plate with a fresh foil seal and shake at 37 °C for 2 min at 600 rpm.
220

221 5.14.3 Remove the plate from the plate shaker and place onto the plate magnet for 2 min to
222 separate the beads from the reaction mixture.

223
224 5.14.4 While retaining the plate on the magnet, remove the foil seal, carefully invert the plate
225 over a waste container, and gently flick the supernatant out of the wells.

226
227 5.14.5 Blot the plate on absorbent paper while still holding the plate on the magnet.

228
229 5.15 Wash the reaction wells as described below (second post-detection antibody incubation
230 wash).

231
232 5.15.1 Remove the plate from the plate magnet and add 150 μ L of PBS-TBN to each well.

233
234 5.15.2 Cover the plate with a fresh foil seal and shake at 37 °C for 2 min at 600 rpm.

235
236 5.15.3 Remove the plate from the plate shaker and place onto the magnetic plate separator for
237 2 minutes to separate the beads from the reaction mixture.

238
239 5.15.4 While retaining the plate on the magnet, remove the foil seal, carefully invert the plate
240 over a waste container, and gently flick the supernatant out of the wells.

241
242 5.15.5 Remove the plate from the magnet.

243
244 5.16 Add 100 μ L of PBS-TBN to each well.

245
246 5.17 Cover with foil seal and shake for 2 min at 37 °C.

247
248 5.18 Remove the plate from the plate shaker. Then remove the foil seal and read 50 μ L of each
249 well on the flow analyzer according to the User Manual. A brief description is provided below.

250
251 5.18.1 Select the drop-down menu in the upper left corner of the screen and navigate to **Plate**
252 **Configuration**.

253
254 5.18.2 Select **Run Plate**.

255
256 5.18.3 Select the **Eject** icon to eject the plate carrier.

257
258 5.18.4 Load the plate onto the plate carrier and select the **Retract** icon to retract the plate
259 carrier.

260
261 5.18.5 Select the **Run** icon to run the plate.

262
263 **6. Single reporter neutralization assay protocol**

264

- 265 6.1 Prepare a fresh multiplex bead mix as described in section 3 above.
266
267 6.2 Dilute patient and control sera 1:500 as follows.
268
269 6.2.1 Dilute serum 1:100 by mixing 10 μ L of serum into 990 μ L of PBS-TBN.
270
271 6.2.2 Dilute serum a further 5-fold by adding 40 μ L of 1:100 serum into 160 μ L of PBS-TBN.
272
273 6.3 Add 50 μ L of multiplex bead mix to each assigned well of a 96-well non-binding microtiter
274 plate.
275
276 6.4 Add 25 μ L of 2 μ g/mL ACE2 dilution to each well.
277
278 6.5 Cover the plate with a microplate foil seal and incubate on a heated plate shaker at 37 °C
279 for 2 min with shaking at 600 rpm.
280
281 6.6 Remove the plate from the plate shaker and remove the foil seal.
282
283 6.7 Add 25 μ L of 1:500 serum or PBS-TBN to the assigned wells.
284
285 6.8 Cover the plate with a microplate foil seal and incubate on a heated plate shaker at 37 °C
286 for 15 min with shaking at 600 rpm.
287
288 6.9 For the remainder of the assay procedure, follow steps 5.4 to 5.18.5 as described above.
289

290 **7. Conversion to a dual reporter IgG and IgM serological assay**

- 291
292 7.1 Prepare antigen-coupled beads as described in Section 2 above.
293
294 NOTE: An additional bead set, coupled with human IgM, is included to monitor for the addition
295 of anti-IgM detection reagents, as well as an additional AC bead set (220) to monitor collection
296 of MFI for RP2. All bead stocks are at 1×10^6 beads/mL for inclusion in multiplex bead mixes as
297 described below.
298
299 7.2 Prepare fresh multiplex bead mixes as described in Section 3 above.
300
301 7.3 Dilute serum samples and IgG-stripped negative control sera are diluted 1:1000 in Section
302 4 above.
303
304 7.4 Add 50 μ L of multiplex bead mix to each assigned well of a 96-well non-binding microtiter
305 plate.
306
307 7.5 Pipette 50 μ L of 1:1000 serum or PBS-TBN to the appropriate wells. The final dilution of
308 serum in the well is 1:2000.

309
310 7.6 Cover the plate with a microplate foil seal and incubate on a heated plate shaker at 37 °C
311 for 15 min with shaking at 600 rpm.

312
313 7.7 Remove the plate from the plate shaker and place onto the magnetic plate separator for
314 2 min to separate the beads from the reaction mixture.

315
316 7.8 While retaining the plate on the magnet, remove the foil seal, carefully invert the plate
317 over a waste container, and gently flick the supernatant out of the wells.

318
319 7.9 Blot the plate on absorbent paper while still holding the plate on the magnet.

320
321 7.10 Wash the reaction wells as described below (first post-sample incubation wash).

322
323 7.10.1 Remove the plate from the plate magnet and add 150 µL of PBS-TBN to each well.

324
325 7.10.2 Cover the plate with a fresh foil seal and shake at 37 °C for 2 min at 600 rpm.

326
327 7.10.3 Remove the plate from the plate shaker and place onto the plate magnet for 2 min to
328 separate the beads from the reaction mixture.

329
330 7.10.4 While retaining the plate on the magnet, remove the foil seal, carefully invert the plate
331 over a waste container, and gently flick the supernatant out of the wells.

332
333 7.10.5 Blot the plate on absorbent paper while still holding the plate on the magnet.

334
335 7.11 Wash the reaction wells again with PBS-TBN (second post-sample incubation wash),
336 exactly as described above in steps 7.10.1–7.10.5.

337
338 7.12 Remove the plate from the magnet.

339
340 7.13 Make a fresh Detection Reagent mix with the required combination of reagents and add
341 50 µL or 100 µL to each well.

342
343 NOTE: Detection Reagent mixes containing the anti-Human IgG conjugated to the Brilliant Violet
344 421 reporter dye (BV; excitation at 405 nm and emission at 421 nm) were used at 50 µL/well due
345 to limiting amounts of the stock reagents available. All other Detection Reagent mixes are used
346 at 100 µL/well.

347
348 7.14 Cover the plate with a microplate foil seal and incubate on a heated plate shaker at 37 °C
349 for 15 min with shaking at 600 rpm.

350
351 7.15 Remove the plate from the plate shaker and place onto the magnetic plate separator for
352 2 min to separate the beads from the reaction mixture.

7.16 While retaining the plate on the magnet, remove the foil seal, carefully invert the plate over a waste container, and gently flick the supernatant out of the wells.

7.17 Blot the plate on absorbent paper while still holding the plate on the magnet.

7.18 Wash the reaction wells as described below with PBS-TBN (first wash after detection antibody incubation).

7.18.1 Remove the plate from the plate magnet and add 150 μ L of PBS-TBN to each well.

7.18.2 Cover the plate with a fresh foil seal and shake at 37 °C for 2 min at 600 rpm.

7.18.3 Remove the plate from the plate shaker and place onto the magnetic plate separator for 2 min to separate the beads from the reaction mixture.

7.18.4 While retaining the plate on the magnet, remove the foil seal, carefully invert the plate over a waste container, and gently flick the supernatant out of the wells.

7.18.5 Blot the plate on absorbent paper while still holding the plate on the magnet.

7.19 Wash the reaction wells again with PBS-TBN (second wash after detection antibody incubation), exactly as described above in steps 7.17.1–7.17.5.

7.20 Remove the plate from the magnet.

7.21 Add 100 μ L of PBS-TBN to each well.

7.22 Cover with foil seal and shake for 2 min at 37 °C.

7.23 Remove the plate from the plate shaker. Then remove the foil seal and read 50 μ L of each well on the flow analyzer as described above in steps 5.18.1–5.18-5.

8. Conversion to dual reporter neutralization assay

8.1 Use the coupled beads as described in step 7.1 above.

8.2 Prepare a fresh multiplex bead mix as described in step 7.2.

8.3 Serum samples and IgG-stripped negative control sera are diluted 1:500 as follows.

8.3.1 Prepare a 1:100 dilution by mixing 10 μ L of serum into 990 μ L of PBS-TBN.

8.3.2 Dilute the 1:100 samples a further 5-fold by adding 40 μ L of the 1:100 dilutions to 160 μ L

of PBS-TBN.

8.4 Add 50 μ L of multiplex bead mix to each assigned well of a 96-well non-binding microtiter plate.

8.5 Add 25 μ L of 2 μ g/mL ACE2 dilution to each well.

8.6 Cover the plate with a microplate foil seal and incubate on a heated plate shaker at 37 °C for 2 min with shaking at 600 rpm.

8.7 Remove the plate from the plate shaker and remove the foil seal.

8.8 Add 25 μ L of 1:500 serum or PBS-TBN to the assigned wells.

8.9 Cover the plate with a foil seal and incubate on a heated plate shaker at 37 °C for 15 min with shaking at 600 rpm.

8.10 Remove the plate from the plate shaker and place onto the magnetic plate separator for 2 min to separate the beads from the reaction mixture.

8.11 With the plate retained on the magnet, remove the foil seal, carefully invert the plate over a waste container, and gently flick the supernatant out of the wells.

8.12 Blot the plate on absorbent paper while still holding the plate on the magnet.

8.13 Wash the reaction wells as described below (first post sample incubation wash).

8.13.1 Remove the plate from the plate magnet and add 150 μ L of PBS-TBN to each well.

8.13.2 Cover the plate with a fresh foil seal and shake at 37 °C for 2 min at 600 rpm.

8.13.3 Remove the plate from the plate shaker and place onto the magnetic plate separator for 2 min to separate the beads from the reaction mixture.

8.13.4 With the plate retained on the magnet, remove the foil seal, carefully invert the plate over a waste container, and gently flick the supernatant out of the wells.

8.13.5 Blot the plate on absorbent paper while still holding the plate on the magnet.

8.14 Wash the reaction wells again with PBS-TBN (second post-sample incubation wash), exactly as described above in steps 8.13.1–8.13.5.

8.15 Remove the plate from the magnet.

8.16 Make a fresh Detection Reagent mix with the required combination of reagents and add 50 μ L or 100 μ L to each well.

NOTE: Detection Reagent mixes containing the anti-Human IgG conjugated to the BV reporter dye were used at 50 μ L/well due to limiting amounts of the stock reagents available (see **Table of Materials**). All other Detection Reagent mixes are used at 100 μ L/well.

8.17 Cover the plate with a microplate foil seal and incubate on a heated plate shaker at 37 °C for 15 min with shaking at 600 rpm.

8.18 Remove the plate from the plate shaker and place onto the magnetic plate separator for 2 min to separate the beads from the reaction mixture.

8.19 With the plate retained on the magnet, remove the foil seal, carefully invert the plate over a waste container, and gently flick the supernatant out of the wells.

8.20 Blot the plate on absorbent paper while still holding the plate on the magnet.

8.21 Wash the reaction wells as described below (first was after detection antibody incubation).

8.21.1 Remove the plate from the plate magnet and add 150 μ L of PBS-TBN to each well.

8.21.2 Cover the plate with a fresh foil seal and shake at 37 °C for 2 min at 600 rpm.

8.21.3 Remove the plate from the plate shaker and place onto the magnetic plate separator for 2 min to separate the beads from the reaction mixture.

8.21.4 With the plate retained on the magnet, remove the foil seal, carefully invert the plate over a waste container, and gently flick the supernatant out of the wells.

8.21.5 Blot the plate on absorbent paper while still holding the plate on the magnet.

8.22 Wash the reaction wells again with PBS-TBN (second wash after detection antibody incubation), exactly as detailed above in steps 8.21.1–8.21.5.

8.22.1 Remove the plate from the plate magnet and add 150 μ L of PBS-TBN to each well.

8.22.2 Cover the plate with a fresh foil seal and shake at 37 °C for 2 min at 600 rpm.

8.22.3 Remove the plate from the plate shaker and place onto the magnetic plate separator for 2 minutes to separate the beads from the reaction mixture.

8.22.4 With the plate retained on the magnet, remove the foil seal, carefully invert the plate

over a waste container, and gently flick the supernatant out of the wells.

8.22.5 Remove the plate from the magnet.

8.23 Add 100 μ L of PBS-TBN to each well.

8.24 Cover with foil seal and shake for 2 min at 37 °C.

8.25 Remove the plate from the plate shaker. Then remove the foil seal and read 50 μ L of each well on the flow analyzer as described above in steps 5.18.1–5.18-5.

REPRESENTATIVE RESULTS:

Optimization of the Single Reporter Serological Assay.

The coupling strategy for all SARS-CoV-2 antigens is described in Cameron, et al., 2020¹ and in step 2.1 above. The strategy for coupling confirmation is described in step 2.3 above. For an efficient coupling, the median fluorescence intensity (MFI) values should be significantly higher than the no serum/no antibody negative control reaction(s) and should demonstrate a dose response of decreasing MFI signal with decreasing amounts of detection reagent. MFIs of approximately 1,000 MFI units or greater over background MFI generally indicate efficient protein coupling and are dependent on the concentration and specificity of the detection reagent used. Using this method, confirmation of S, RBD and Nc antigen coupling were tested for two lots of beads (**Figure 2**). The MFI signals for all antigens were significantly higher than the background MFI and showed a dose response with titration of the detection reagents. The data also showed good reproducibility of MFI levels between the two different lots of coupled beads.

Optimization of the single reporter assay was done as described in Cameron, et al., 2020¹. Each antigen was tested at different coupling concentrations with different dilutions of several serum samples representing high, medium, low, and negative samples. In addition, because the assay measures IgG titers, an IgG-stripped serum was used as an additional negative sample. Results of a representative sample dilution series with different antigen coupling levels are shown in **Figure 3**. This initial dilution series was tested with a fixed concentration of detection reagents to get a potential range of representative antigen-specific and background MFI levels. From these data, additional sample dilution ranges and antigen coupling levels were determined and tested further with additional samples (data not shown).

Next, the concentration of the detection reagent was optimized for use with the different serum dilutions and antigen coupling levels. Representative data using 6 positive and 2 negative samples are shown in **Figure 4**. Following additional testing with several hundred negative pre-COVID-19 serum samples, the final optimal antigen coupling levels and detection reagent concentrations were selected for the final assay protocol, as described in section 5.

Conversion to a single reporter neutralization assay

Conversion of the single reporter serological assay to a neutralization assay was achieved by

adding an incubation step using ACE2 (as a competing reagent) with the beads prior to adding the serum sample. By titrating the amount of ACE2 added, inhibition of IgG binding to the S and RBD antigens was observed, as shown in **Figure 5**. While the 11 patient sera had different titers of IgG, each showed a significant decrease in MFI signal with increasing ACE2 concentration. As expected, ACE2 only impacts the MFI signals of S and RBD, but not the Nc antigen.

The percent residual MFI signal relative to 0 µg/mL ACE2 for all samples is shown in **Figure 6**. For most samples, a signal loss >30% was observed for S and RBD with increasing ACE2 concentration. As expected, there was no effect of ACE2 on the Nc signal.

Conversion to a dual reporter IgG and IgM serological assay

For the dual reporter IgG and IgM assay, the standard single reporter assay conditions were used to test different antibody and reporter combinations for the two reporter channels. The RP1 channel is similar to previous flow analysis instruments with detection of “orange” fluorescence for reporter dyes such as PE. The second channel, RP2, employs a violet excitation laser that can be used for detecting the “blue” fluorescence of dyes excited at 402 nm with emission spectra peaks in the 421–441 nm range.

Several different anti-human IgG and IgM antibody combinations were tested at various concentrations with the standard 2,000-fold sample dilution and incubation conditions detailed in section 6. An initial test of anti-IgM conjugated to the DyLight 405 reporter dye (DL 405; excitation at 400 nm and emission at 420 nm) for the RP2 channel was performed using a set of 11 PCR-positive samples with DFSO ranging from 5 to >60 days, and an IgG-stripped serum sample. This detection reagent did not produce a high signal in the RP2 channel as shown in **Figure 7A,B,C**. For samples with the high IgM RP2 MFI, the highest signals were seen for RBD and Nc (**Figure 7B,C**). While IgM titers should be elevated at this time in some samples, the observed MFI levels did not exceed 140 MFI units with the IgM detection reagent. Furthermore, the control bead for IgM lacked a significant dynamic range for MFI when using DL 405 conjugated to anti-IgM compared to a PE-labeled anti-IgM RP1 reporter used at the same concentrations (**Figure 7D**).

Because a suitable RP2 reporter-labeled detection reagent for IgM was not available, the original biotin anti-IgG with SASB was tested for use in the RP2 channel. The signal for IgG in the RP2 channel with SASB did not have the same dynamic range as with SAPE, but it still had a suitable range for measuring IgG titers for S, RBD, and Nc (**Figure 8**). These data led to testing different RP1 IgM and RP2 IgG detection reagent combinations, as described below.

Conversion to dual reporter neutralization assay

The dual reporter serological assay was converted to a neutralization assay by adding an incubation step with ACE2 and beads, prior to adding the serum sample. A 2-fold dilution series of ACE2 starting at 16 µg/ml was used, and then tested with sets of 11 samples and stripped sera as described above. In addition, 3 different combinations of RP1 IgM and RP2 IgG detection antibody mixes were tested on a set of 11 samples and stripped serum controls. The final combinations and reagent concentrations determined to be optimal are detailed in **Table 1**.

For IgG detection, the biotin anti-human IgG/SASB used in the original single reporter assay was tested along with another anti-IgG conjugated to the BV reporter dye. While the BV conjugate had high RP2 MFI signals, the MFI signal intensity for IgG titer varied across the ACE2 titrations (**Figure 9A,C,E**). The biotin anti-human IgG/SASB combination had more consistent MFI signal decrease across the ACE2 titration compared to the BV conjugate, demonstrating a dose response, although the MFI levels were lower. The signal fluctuation by the BV conjugate across ACE2 concentrations also interfered with determining the percent inhibition by ACE2 across the range of IgG titers compared to the biotin anti-human IgG/SASB combination (**Figure 9B,D,F**). Samples with low MFI signals, such as the DFSO 3 sample, do not have high enough titers to evaluate the performance of these reagents or measure IgG neutralizing capacity.

For determining IgM titers in the RP1 channel, a PE-conjugated anti-IgM was compared to an anti-human IgM conjugated to the DyLight 549 reporter dye (DL 549; excitation at 562 nm and emission at 576 nm) at the concentrations listed in **Table 1**. Of these two reagents, the PE anti-IgM displayed higher MFIs than those generated by DL 549 anti-human IgM for the samples tested (**Figure 9A,C,E**). The IgM control bead measuring the addition of these reagents had different average MFIs of $\approx 17,000$ for the PE anti-IgM and 2,723 for the DL 549 conjugate. Even with these differences, the impact of the various ACE2 concentrations on MFI was measurable with the DL 549 conjugate. For determining ACE2 percent inhibition of IgM binding, there was a slight but insignificant difference between the two IgM detection reagents (**Figure 9B,D,F**).

FIGURE AND TABLE LEGENDS:

Table 1: List of RP1 and RP2 reporter dye combinations tested. Several different RP1 and RP2 reporter dyes were evaluated for measuring IgG and IgM titers. The three combinations shown above were tested at different RP1 detection modes and for optimization of the neutralization assay. For RP2 combinations using SASB, concentrations for the biotinylated detection antibody and SASB are shown. *Final Concentration represents the final concentration in the reaction well.

Figure 1: Flowchart highlighting main assay protocol steps.

Figure 2: Confirmation of antigen coupling. Antigens were coupled at 100, 10 and 1 pmol/ 10^6 beads. For S and RBD antigens, rabbit anti-S sera was used and detected with PE anti-rabbit IgG. For Nc, a Protein G-purified rabbit polyclonal anti-Nc was used. The titration curves for beads coupled with 10 pmol S (**A**), 100 pmol RBD (**B**) and 100 pmol Nc (**C**) are shown.

Figure 3: Optimization of serum dilution with different antigen coupling levels. Serum samples representing high, medium, low, and negative samples were tested with the different pmol antigen/ 10^6 bead couplings. A 4-fold serum titration was tested with a multiplex mix of antigen-coupled beads using the assay protocol described in section 5. The MFI values of the titration curves for S (**A**), RBD (**B**), and Nc (**C**) are shown.

Figure 4: Determination of optimal serum dilution and detection reagent concentration. Eight

serum samples, including negative patients (2 samples) and low- to high-positive (6 samples), were tested at additional serum dilutions with three different detection reagent concentrations. The MFI results are shown for S (A), RBD (B), and Nc (C). Based on these and other data (not shown), a sample dilution of 1:2,000 and a biotin detection antibody concentration of 0.62 µg/mL was selected.

Figure 5: Optimization of single reporter neutralization assay. Modification of the single reporter serological assay to a neutralizing antibody detection assay was performed by adding an incubation step with ACE2 as a competitor, as described in Section 6. A set of 11 samples, ranging from low- to high-positive titer sera, and an IgG-stripped serum sample were tested with a two-fold dilution series of ACE2 starting at 4 µg/mL, using standard assay conditions. Across this range of ACE2 concentrations, a decrease in MFI can be seen with increased ACE2 for S (A) and RBD (B), but not for Nc (C).

Figure 6: Percent inhibition of signal with ACE2. The percent residual signals for the set of 11 samples (from Figure 4) are shown. For each sample, regardless of its IgG titer, a maximum decrease of ≥30% MFI signal was achieved for S (A) and RBD (B) at the highest ACE2 concentration. For the Nc antigen (C) there was no significant inhibition of signal with any amount of ACE2 tested.

Figure 7: Test of DyLight 405 anti-IgM as RP2 IgM detection reagent. A set of 11 samples and IgG-stripped sera were used to test DL 405-conjugated anti-IgM as an RP2 detection reagent. The RP2 IgM MFI signals for S (A), RBD (B), and Nc (C) are shown. The highest MFI signal for any antigen with any sample was approximately 140 MFI units for RBD with sample H (B). Even the signal on the IgM control bead set was less than 1,000 MFI in the RP2 channel across the whole antibody titer range and did not show the increased dynamic range observed with the PE anti-IgM detection antibody in the RP1 channel (D).

Figure 8: Optimization of SASB as RP2 reporter with Biotin anti-IgG. A set of 11 samples and IgG-stripped sera were used to test a dilution series of SASB with the standard 0.62 µg/mL biotin anti-IgG. The resulting IgG MFIs in the RP2 channel for S (A), RBD (B), and Nc (C) are shown.

Figure 9: Comparison of three different combinations of RP1 IgM and RP1 IgG detection mixes. Three different detection antibody mixes were tested on 11 samples and IgG stripped sera. The combinations and final concentrations used are described in Table 1. All samples were tested with a 2-fold ACE2 dilutions, starting at 16 µg/mL. Data for samples at DFSO of 3, 12, and 30 days are shown. The MFI signals for RP1 IgM (orange) and RP2 IgG (blue) are shown in A (DFSO 3), C (DFSO 12), and E (DFSO 30). The ACE2 percent residual MFIs are shown in B (DFSO 3), D (DFSO 12), and F (DFSO 30). Signals representing a ≥30% decrease as compared to No ACE2 controls are highlighted light green.

DISCUSSION:

This study expanded the functionality of a multiplex fluorescent microsphere assay, 3Flex, that qualitatively assesses the IgG response to infection by SARS-CoV-2¹. While many serological

assays for COVID-19 detect a single antigen, this assay simultaneously evaluates S, RBD, and Nc antigens, and includes an internal antibody isotype control⁷. In addition, ACE2 is used as an indicator to detect neutralizing antibody titers to SARS-CoV-2.

Multi-antigen assays may be particularly useful in the future for discriminating immune responses between infected and vaccinated persons. With SARS-CoV-2, if only S and RBD antibodies are evaluated, it would be impossible to discern whether this was due to a past infection or to a vaccination antibody response. None of the vaccines that are currently in use target the Nc antigen, so by evaluating S/RBD and Nc antigens, it is possible to distinguish past viral infection (Nc- and S/RBD-positive) from a vaccination response (S/RBD-positive only; Nc-negative).

In the current study, the 3Flex serological and neutralization assays (sections 5 and 6) were transferred onto a new dual reporter flow analysis instrument (sections 7 and 8). Typical immunoassay protocols for bead-based multiplex assays are similar to standard ELISA protocols, following comparable steps for sample incubation, detection antibody/reporter incubation, and analysis of results⁸. Previous studies have also demonstrated how standard ELISA assays can be transferred to a bead-based multiplexing platform⁹.

The assay protocols could be seamlessly transferred from the predecessor single reporter instrument to the new dual reporter system, as shown by the results obtained for the test samples and controls (**Figure 4** and **Figure 5**). In the serological assay, all positive serum samples demonstrated a characteristic dose-responsive decrease in signal corresponding to both a higher sample dilution and a lower detection antibody concentration, whereas the signals for the negative samples remained at background levels. Similarly, for the neutralization assay, decreasing signals for S and RBD, but not Nc, corresponded to increasing concentrations of the ACE2 competitor, whereas the IgG stripped serum was much less affected by the addition of ACE2. Pre-incubating the beads with ACE2 for 2 min prior to addition of the serum sample (as described in sections 6 and 8), was also compared to combining the beads with ACE2 and serum at the same time and produced little difference in the results (data not shown). However, because the results obtained with the beads plus ACE2 pre-incubation step were more consistent, it was the final procedure adopted for the neutralization assay protocol (sections 6 and 8).

Because the dual reporter system incorporates a second fluorescent reporter channel, both assays were converted into isotyping assays to simultaneously measure both IgG and IgM responses. The dual reporter functionality of the flow analyzer allows collection of fluorescent signals from two reporter detection reagents for each analyte in the multiplex for each sample, effectively doubling the data generated per sample in an assay. For development and optimization of the dual reporter assay protocols, several commercially available reporter dyes and detection antibodies were evaluated for the new RP2 channel (**Figure 7** and **Figure 8**) to determine the best option(s) available. The anti-IgM antibody conjugated with the DL 405 reporter dye did not perform well in the RP2 channel (**Figure 7**), so biotin anti-IgG with SASB was subsequently evaluated for detection in the RP2 channel (**Figure 8**). Three combinations of RP1

and RP2 detection reagents were compared in the dual reporter neutralization assay (**Table 1** and **Figure 9**) to determine the best performing combinations for simultaneous detection of IgG and IgM neutralizing antibodies. While there were significant differences in the signal intensity between the RP1 channel using PE-anti IgM and the RP2 channel using DL 549 anti-IgM, there was no significant difference in the measurement of the percent binding inhibition by ACE2.

Several limitations to the current method and this study are acknowledged. First, the samples tested and used in the method development and optimization were limited to sets of 8-11 samples. More samples will be tested in expanded studies to verify the method is fully optimized. Second, a limited number of reporter fluorophores and reporter pairs were tested. Further testing of all available reporters in all possible combinations will determine which reagents can be used successfully in this method. The anti-IgG antibody conjugated to biotin was different than the anti-IgG antibody conjugated to the BV 421 reporter. So, although variability in the signal intensity for the BV 421 anti-IgG existed, it could be due to the specificity of the antibody and not related to the reporter dye. Testing other available BV 421-conjugated antibodies may identify another antibody that would perform well in the RP2 channel. Future studies will also evaluate the combination of PE anti-IgM with BV 421 anti-human IgG for detection in RP1 and RP2, respectively. Last, even with the dual reporter capability of the instrument, assays are limited to two isotypes (two parameters) per reaction. If additional isotypes are to be measured or full isotyping is desired, then additional reactions using the same two reporter channels would need to be run in separate wells.

Bead-based multiplexing technology allows for rapid and flexible assay design with easy implementation into routine laboratory workflows^{3,4,10}. This study demonstrated the ease in transferring the previously described 3Flex serological and neutralization assays for SARS-CoV-2 onto a flow analysis system for measuring IgG with a single reporter, or with slight modification, simultaneously measuring both IgG and IgM responses using two reporters. The dual reporter option has clear advantages over single reporter platforms because it can provide twice as much data per sample, using half the sample volume, and in half the number of reaction wells. With the appropriate reporter dyes and detection antibodies, isotyping assays can be easily developed and performed on the dual reporter flow analysis instrument.

ACKNOWLEDGMENTS:

This report was funded by Luminex Corporation (Austin, TX). The authors thank Matthew Silverman PhD (Biomedical Publishing Solutions, Delray Beach, FL) for scientific editing assistance.

DISCLOSURES:

Stephen Angeloni and Sherry Dunbar are employees of Luminex Corporation that produces some of the reagents and instruments that were used in this study.

REFERENCES:

1. Cameron, A. et al. A multiplex microsphere IgG assay for SARS-CoV-2 using ACE2-mediated inhibition as a surrogate for neutralization. *Journal of Clinical Microbiology*. **59** (2),

e02489-20 (2020).

2. Waterboer, T. et al. Multiplex human papillomavirus serology based on in situ-purified glutathione s-transferase fusion proteins. *Clinical Chemistry*. **51** (10), 1845–1853 (2005).

3. Pickering, J. W. et al. A multiplexed fluorescent microsphere immunoassay for antibodies to pneumococcal capsular polysaccharides. *American Journal of Clinical Pathology*. **117** (4), 589–596 (2002).

4. Ayoub, A. et al. Development of a sensitive and specific serological assay based on Luminex technology for detection of antibodies to Zaire Ebola virus. *Journal of Clinical Microbiology*. **55** (1), 165–176 (2017).

5. Hornbeck, P., Fleisher, T.A., Papadopoulos, N.M. Isotype determination of antibodies. In: *Current Protocols in Immunology*. Vastl, J., ed. **116**, 2.2.1–2.2.7 (2017).

6. Angeloni, S. *Multiplex Bead Mixes Made Easy with an Excel-Based Bead Calculator*, <<https://www.luminexcorp.com/blog/multiplex-bead-mixes-made-easy-with-an-excel-based-bead-calculator/>> (2020).

7. FDA. *EUA Authorized Serology Test Performance*. <<https://www.fda.gov/medical-devices/coronavirus-disease-2019-covid-19-emergency-use-authorizations-medical-devices/eua-authorized-serology-test-performance>> (2020).

8. Abcam. *Indirect ELISA protocol*. <<https://www.abcam.com/ps/pdf/protocols/Indirect%20ELISA%20protocol.pdf>> (2021).

9. Baker, H.N., Murphy, R., Lopez, E., Garcia, C. Conversion of a capture ELISA to a Luminex xMAP assay using a multiplex antibody screening method. *Journal of Visualized Experiments*. (65), e4084 (2012).

10. Pickering, J.W., Martins, T.B., Schroder, M. C., Hill, H. R. Comparison of a multiplex flow cytometric assay with enzyme-linked immunosorbent assay for quantitation of antibodies to tetanus, diphtheria, and Haemophilus influenzae Type b. *Clinical and Diagnostic Laboratory Immunology*. **9** (4), 872–876 (2002).




Figure 2

[Click here to access/download;Figure;Angeloni Figure 2.jpg](#)




Figure 3

[Click here to access/download;Figure;Angeloni Figure 3.jpg](#)








Figure 4

[Click here to access/download;Figure;Angeloni Figure 4.jpg](#)













Figure 9

[Click here to access/download;Figure;Angeloni Figure 9.jpg](#)

Combination	Reporter
1	PE anti-IgM
	Biotin-Ig
	SASB
2	DyLight 549 Anti-Human IgM
	Brilliant Violet 421 Anti-Human IgG
3	DyLight 549 Anti-Human IgM
	Biotin-Ig
	SASB

Final Concentration In Well (µg/mL)	
RP1	RP2
2	-
-	0.62
-	4
1	-
-	1
1	-
-	0.62
-	4

Name of Material/ Equipment

ACE2 (human) (rec.)

Biotin-SP (long spacer) AffiniPure Goat Anti-Human IgG, Fcγ fragment specific

Bovine Serum Albumin

Branson 1800 Cleaner/water bath sonicator

Brilliant Violet 421 AffiniPure Donkey Anti-Human IgG (H+L)

Corning 96-well Black Flat Bottom Polystyrene NBS Microplate

Corning 96-well Microplate Aluminum Sealing Tape, Nonsterile

DyLight 405 AffiniPure Goat Anti-Human IgM, Fc5μ fragment specific

DyLight 549 AffiniPure F(ab')₂ Goat Anti-Human IgM Fc5μ specific

eBioscience Streptavidin Super Bright 436 Conjugate

Fisherbrand Incubating Microplate Shaker

Goat anti-Rabbit IgG (H+L) Cross-Adsorbed Secondary Antibody, PE

Luminex Magnetic Plate Separator

MagPlex beads

PE AffiniPure F(ab')₂ Goat Anti-Human IgM Fc5μ specific

Phosphate Buffered Saline

SARS Nucleocapsid Protein Antibody

SARS Nucleocapsid Protein Antibody

SARS Spike Protein Antibody

SARS Spike Protein Antibody

SARS-CoV-2 (2019-nCoV) Spike S1 Antibody, Rabbit MAb

Sodium Azide

Streptavidin, R-Phycoerythrin Conjugate (SAPE)

Tube Revolver/ Rotator

TWEEN 20

xMAP Antibody Coupling Kit

xMAP INTELLIFLEX DR-SE

Company	Catalog Number
AdipoGen Life Sciences	AG-40B-0192-3050
Jackson ImmunoResearch Laboratories	109-065-098
MilliporeSigma	A7888
Various	
Jackson ImmunoResearch Laboratories	709-675-149
Corning	3650
Corning	6570
Jackson ImmunoResearch Laboratories	109-475-043
Jackson ImmunoResearch Laboratories	109-506-129
ThermoFisher Scientific	62-4317-82
ThermoFisher Scientific	02-217-757
ThermoFisher Scientific	P-2771MP
Luminex Corporation	CN-0269-01
Luminex Corporation	MC100XX-ID
Jackson ImmunoResearch Laboratories	109-116-129
MilliporeSigma	P3813
Novos Biologicals	NB100-56683
Novos Biologicals	NB100-56049

Novos Biologicals	NB100-56047
Novos Biologicals	NB100-56048
SinoBiological	40150-R007
MilliporeSigma	S2002
ThermoFisher Scientific	S21388
ThermoFisher Scientific	88881001
MilliporeSigma	P9416
Luminex Corporation	40-50016
Luminex Corporation	INTELLIFLEX-DRSE-RUO

Comments/Description

Sigma-Aldrich

Discontinued but available from multiple vendors

Fisher Scientific

Sigma-Aldrich

Unconjugated, rabbit, polyclonal, Protein G purified

Unconjugated, rabbit, polyclonal, Protein A purified

Unconjugated, rabbit, polyclonal, Protein G purified
Unconjugated, rabbit, polyclonal, Protein G purified

Sigma-Aldrich
premium grade

Sigma-Aldrich

Sherry Dunbar, PhD
Senior Director, Global Scientific Affairs
Luminex Corporation
Austin, Texas, USA
sdunbar@luminexcorp.com

Vineeta Bajaj, PhD
Review Editor
JoVE
vineeta.bajaj@jove.com

March 16, 2021

Dear Dr. Bajaj,

We thank you, your editorial team, and the reviewers for your ongoing support in developing our manuscript, JoVE62487R1 "A Rapid, Multiplex Dual Reporter IgG and IgM SARS-CoV-2 Neutralization Assay for a Multiplexed Bead-based Flow Analysis System". We are confident that the changes resulting from your guidance have greatly enhanced the quality of our report.

We have strived to make all requested amendments to our paper, and outline all changes made below on a point-by-point basis, provided below. Our responses are shown in italic font immediately after each specific editorial comment. We have also provided two versions of our revised manuscript, one with changes showing/tracked and a second "Clean Copy" with changes accepted for easier reading.

Thank you again for your ongoing assistance.

Sincerely,

Sherry Dunbar, PhD

Editorial comments:

Changes to be made by the Author(s):

1. Please take this opportunity to thoroughly proofread the manuscript to ensure that there are no spelling or grammar issues. Please define all abbreviations at first use.

1. We have thoroughly reviewed our manuscript and have repaired any spelling or grammatical inconsistencies that may have existed. We also have made certain that all abbreviations used in the report are defined at first mention.

2. Please make the title concise by deleting "Development of".

2. *As requested, we have deleted “Development of” from the beginning of our paper’s title.*

3. Please provide an email address for each author.

3. *Per editorial reminder and JoVE Author Instructions, we now include on the Title Page the email addresses for all coauthors.*

4. Please revise the following lines to avoid overlap with previously published work: 32-34; 313-314; 343-344; 362-363; 378-379; 394-395; 712-714

4. *We have reviewed and revised all editor-identified lines mentioned above to avoid overlap with previously published work to the best of our ability. The lines indicated at 313-395 are protocol steps but we have also condensed several sections of these protocol steps themselves, cited in the above pagination, to avoid unnecessary redundancy. This is discussed again in more detail in Editorial Comment #11, below.*

5. JoVE cannot publish manuscripts containing commercial language. This includes trademark symbols (™), registered symbols (®), and company names before an instrument or reagent. Please remove all commercial language from your manuscript (text, figure legends, figures, tables) and use generic terms instead. All commercial products should be sufficiently referenced in the Table of Materials and Reagents.

For example: Streptavidin Super Bright 436 Conjugate (SASB); Tween®20; AC220; Brilliant Violet 421; DyLight 549 etc

5. *We have deleted all trademark, registry and company names, and any commercial language from all sections of the manuscript. We have by necessity retained the commercial names for several of the fluorophores used, simply because there is no generic term available to adequately describe them without excessive verbiage—in these few instances, we were certain not to include company sources or any commercial language or symbols and we created abbreviations for these names for any subsequent mentions in the report.*

6. Please revise the text, especially in the protocol, to avoid the use of any personal pronouns (e.g., "we", "you", "our" etc.).

6. *We have reworded all sentences in which personal pronouns had been used to eliminate their presence; these revisions are clearly shown in the tracked changes.*

7. Please ensure that all text in the protocol section is written in the imperative tense as if telling someone how to do the technique (e.g., “Do this,” “Ensure that,” etc.). The actions should be described in the imperative tense in complete sentences wherever possible. Avoid usage of phrases such as “could be,” “should be,” and “would be” throughout the Protocol. Any text that cannot be written in the imperative tense may be added as a “Note.” However, notes should be concise and used sparingly. Please include all safety procedures and use of hoods, etc.

7. *We have reviewed all manuscript text and have modified several sections and phrases to be more instructive and written in the imperative tense. We have occasionally added explanatory*

notes where descriptions are ancillary to protocol performance, but provide essential background information. We have not included any safety warnings because there are really no specific safety precautions needed beyond usual laboratory practices.

8. The Protocol should contain only action items that direct the reader to do something. Please move the discussion about the protocol to the Discussion.

8. We have rewritten the protocol steps so that they only contain “action items”.

9. Please note that your protocol will be used to generate the script for the video and must contain everything that you would like shown in the video. Please add more details to your protocol steps. Please ensure you answer the “how” question, i.e., how is the step performed? Alternatively, add references to published material specifying how to perform the protocol action. Please add more specific details (e.g., button clicks for software actions, numerical values for settings, etc) to your protocol steps. There should be enough detail in each step to supplement the actions seen in the video so that viewers can easily replicate the protocol.

9. We are confident that our protocol is now sufficiently detailed to allow easy and uncomplicated performance by any competent laboratory technician. There really are not many areas where settings and button clicks are used in the protocol, besides the final step of actually reading results on the flow cytometry unit. The User Manual is now referenced in these steps, and the key steps for loading and running the instrument are now briefly described. The start-up, calibration, programming, etc., for the instrument is conducted prior to running an assay and requires dozens of pages to describe in detail every step. Several chapters in the User Manual are devoted to these topics. The laboratory will also have their own preferred instrument settings, and procedures which may vary from lab-to-lab and run-to-run. Thus, we feel that detailed instructions for operating the instrument are beyond the scope of this report. Any technician performing this assay will already have received training and have experience with their instrument and will possess the equipment User Manual. We are certain that our protocol as provided is sufficient for assay reproduction that produces consistent and reliable results. We can also add the User Manual to the list of references if desired; however, it is not available for free download and is only provided with the instrument.

10. 2.3.1 (as an example): Please describe in detail how coupling was confirmed here or cite references that describe this so that readers know how to do this.

10. The procedures for coupling and coupling confirmation are described and cited in the first reference (Cameron, et al.). We have added a notation to coupling confirmation as a “Note” to this step, and have made certain that all other steps have sufficient explanatory detail. Coupling confirmation is a standard technique that is not unique to this protocol and is described in the cited prior work and its references.

11. Please make sure that you have highlighted only up to 3 pages of the protocol text including headings and spacings to indicate what is to be filmed.

11. We fully appreciate that our highlighted text was very close to the 3-page limit of what is to be filmed. We have condensed several steps throughout the protocol, including in the highlighted regions, to avoid redundancy. This has resulted in our highlighted steps for filming occupying less than the 3-page maximum allowance. Thank you for bringing this to our attention.

12. Please consider combining Tables 1 and 2 with the Table of Materials because they also contain information about materials used in the protocol.

12. As requested we have merged these sections, and appreciate that this recommended change greatly streamlines our materials presentation.

13. Please include an Acknowledgements section, containing any acknowledgments and all funding sources for this work.

13. We have now included an Acknowledgements section that provides funding information for this project, and also recognizes our enlistment of a professional medical writer to assist in making the editorial-requested changes to this report.

14. Please ensure that the references appear as the following: [Lastname, F.I., LastName, F.I., LastName, F.I. Article Title. Source (italics). Volume (bold) (Issue), FirstPage–LastPage (YEAR).] For more than 6 authors, list only the first author then et al. Please include volume and issue numbers for all references. Please do not abbreviate journal names; use title case for all journal names.

14. We have revised our Bibliography to precisely conform with the formatting specifications provided by our Editor, and in line with published JoVE Author Instructions.

15. Miscellaneous modifications performed during revision.

a. We have included a new Figure 1 Flowchart that highlights primary assay steps. This was a minor point suggested by Reviewer #2, which we were happy to accommodate.

b. We have performed minor edits specified by JoVE Author Instructions such as italicizing all Latin abbreviations, abbreviating “minutes” to “min”, etc.

Author Response to Reviewer Comments:

We thank the reviewers for their positive assessment of our report. Reviewer #1 had no modification request, neither minor nor major. Reviewer #2 made requested the inclusion of “summarized scheme of the whole protocol...”, which we now do with the creation and addition of new Figure 1, a Flowchart that clearly delineates the major steps associated with our assay. We appreciate this critique and agree that it enhances reader/operator comprehension. We have also reviewed all other included Tables and Figures for easy readability, including seeking impartial outside assessments, and have identified no inordinate issues that require modification. All Figures are at 300 dpi and we have enlarged them to ensure that they are at 11 inches at their tallest or widest dimension.

Reviewers' Comments:

Reviewer #1:

Manuscript Summary:

The manuscript "Development of a Rapid, Multiplex Dual Reporter IgG and IgM SARS-CoV-2 Neutralization Assay for a Multiplexed Bead-based Flow Analysis System" describe a flow analysis system for bead-based multiplexed assays which provides a two-reporter readout was used for the development of a multiplex antibody neutralization assay that can simultaneously measure neutralizing IgG and IgM antibodies for SARS-CoV-2. The author describe the limitations to the current method and the small number of samples.

The method described bead-based multiplexing technology can be useful for rapid and flexible assay design with easy implementation into routine laboratory workflows and can be useful tool to monitor the waves of infection.

Major Concerns:

none

Minor Concerns:

none

Reviewer #2:

Manuscript Summary:

The authors have developed a multiplex bead-based neutralizing assay for SARS-CoV-2 antibodies against 3 antigens, using flow cytometry. Steps of optimization are described to show the relevance of this method.

Major Concerns:

There are not major concerns about this manuscript, except the lack of comparison with an other neutralizing method, to confirm the results obtained. Thus, this manuscript should be considered as a protocol development study and may be confirm latter in clinical situations.

Minor Concerns:

Figures are not very easy to read.

Authors should provide a summarized scheme of the whole protocole to improve the understanding of the method developed.