

# Journal of Visualized Experiments

## Selection of transporter-targeted inhibitory nanobodies by SSM-based electrophysiology

--Manuscript Draft--

<b>Article Type:</b>	Methods Article - JoVE Produced Video
<b>Manuscript Number:</b>	JoVE62578R2
<b>Full Title:</b>	Selection of transporter-targeted inhibitory nanobodies by SSM-based electrophysiology
<b>Corresponding Author:</b>	Camilo Perez, Ph.D Universitat Basel Basel, Basel-Stadt SWITZERLAND
<b>Corresponding Author's Institution:</b>	Universitat Basel
<b>Corresponding Author E-Mail:</b>	camilo.perez@unibas.ch
<b>Order of Authors:</b>	Natalie Bärland Camilo Perez, Ph.D
<b>Additional Information:</b>	
<b>Question</b>	<b>Response</b>
Please indicate whether this article will be Standard Access or Open Access.	Standard Access (US\$2,400)
Please specify the section of the submitted manuscript.	Biochemistry
Please indicate the <b>city, state/province, and country</b> where this article will be <b>filmed</b> . Please do not use abbreviations.	Basel, Switzerland
Please confirm that you have read and agree to the terms and conditions of the author license agreement that applies below:	I agree to the <a href="#">Author License Agreement</a>
Please provide any comments to the journal here.	
Please indicate whether this article will be Standard Access or Open Access.	Standard Access (\$1400)

**TITLE:**

Selection of Transporter-Targeted Inhibitory Nanobodies by Solid-Supported-Membrane (SSM)-Based Electrophysiology

**AUTHORS AND AFFILIATIONS:**

Natalie Bärländ<sup>1</sup>, Camilo Perez<sup>1†</sup>

<sup>1</sup>Biozentrum, University of Basel, 4056 Basel, Switzerland

**Email Addresses of Co-Authors:**

Natalie Bärländ (natalie.baerland@unibas.ch)

Camilo Perez (camilo.perez@unibas.ch)

**Email Address of Corresponding Author:**

Camilo Perez (camilo.perez@unibas.ch)

**KEYWORDS:**

membrane transporter, membrane proteins, solid-supported-membranes electrophysiology, nanobodies, inhibitory nanobodies, nanobodies selection

**SUMMARY:**

Nanobodies are important tools in structural biology and pose a great potential for the development of therapies. However, the selection of nanobodies with inhibitory properties can be challenging. Here we demonstrate the use of solid-supported-membrane (SSM)-based electrophysiology for the classification of inhibitory and non-inhibitory nanobodies targeting electrogenic membrane transporters.

**ABSTRACT:**

Single domain antibodies (nanobodies) have been extensively used in mechanistic and structural studies of proteins and they pose an enormous potential as tools for developing clinical therapies, many of which depend on the inhibition of membrane proteins such as transporters. However, most of the methods used to determine the inhibition of transport activity are difficult to perform in high-throughput routines and depend on labeled substrates availability thereby complicating the screening of large nanobody libraries. Solid-supported membrane (SSM) electrophysiology is a high-throughput method, used for characterizing electrogenic transporters and measuring their transport kinetics and inhibition. Here we show the implementation of SSM-based electrophysiology to select inhibitory and non-inhibitory nanobodies targeting an electrogenic secondary transporter and to calculate nanobodies inhibitory constants. This technique may be especially useful for selecting inhibitory nanobodies targeting transporters for which labeled substrates are not available.

**INTRODUCTION**

Antibodies are composed of two identical heavy chains and two light chains that are responsible for the antigen binding. Camelids have heavy-chain only antibodies that exhibit similar affinity for their cognate antigen compared to conventional antibodies<sup>1,2</sup>. The single variable domain (VHH) of heavy-chain only antibodies retain the full antigen-binding potential and has been shown to be very stable<sup>1,2</sup>. These isolated VHH molecules or “nanobodies” have been

implemented in studies related to membrane proteins biochemistry as tools for stabilizing conformations<sup>3,4</sup>, as inhibitors<sup>5,6</sup>, as stabilization agents<sup>7</sup>, and as gadgets for the structure determination<sup>8-10</sup>. Nanobodies can be generated by the immunization of camelids for the pre-enrichment of B-cells that encode target-specific nanobodies and subsequent isolation of B cells, followed by cloning of the nanobody library and selection by phage display<sup>11-13</sup>. An alternate way to generate nanobodies is based on *in vitro* selection methods that rely on the construction of libraries and selection by phage display, ribosome display, or yeast display<sup>14,15-20</sup>. These *in vitro* methods require large library sizes but benefit from avoiding animal immunization and favor the selection of nanobodies targeting proteins with relatively low stability.

The small size of nanobodies, their high stability and solubility, strong antigen affinity, low immunogenicity, and relatively easy production, make them strong candidates for the development of therapeutics<sup>21,22,23</sup>. In particular, nanobodies inhibiting the activity of multiple membrane proteins are potential assets for clinical applications<sup>5,24-26</sup>. In the case of membrane transporters, to evaluate whether a nanobody has inhibitory activity, it is necessary to develop an assay that allows the detection of transported substrates and/or co-substrates. Such assays usually involve labeled molecules or the design of substrate-specific detection methods, which may lack a universal application. Furthermore, the identification of inhibitory nanobodies generally requires the screening of large numbers of binders. Thus, a method that can be used in a high-throughput mode and that does not rely on labeled substrates is essential for this selection.

SSM-based electrophysiology is an extremely sensitive, highly time-resolved technique that allows the detection of movement of charges across membranes (e.g., ion binding/transport)<sup>27,28</sup>. This technique has been applied to characterize electrogenic transporters, which are difficult to study using other electrophysiology techniques due to the relative low turnover of these proteins<sup>29-35</sup>. SSM electrophysiology does not require the use of labeled substrates, it is suitable for high-throughput screening, and either proteoliposomes or membrane vesicles containing the transporter of interest can be used. Here, we demonstrate that SSM-based electrophysiology can be used to classify transporter-targeted nanobodies with inhibitory and non-inhibitory properties. As a proof-of-principle, we describe the reconstitution of a bacterial choline transporter into liposomes, followed by detail steps for immobilization of proteoliposomes on the SSM sensors. We next describe how to perform SSM-based electrophysiology measurements of choline transport and how to determine the half-maximal effective concentration (EC<sub>50</sub>). We then show how to use SSM-based electrophysiology to screen multiple nanobodies and to identify inhibitors of choline transport. Finally, we describe how to determine the half maximal inhibitory concentrations (IC<sub>50</sub>) of selected inhibitory nanobodies.

## PROTOCOL:

### 1. Membrane protein reconstitution

1.1. Mix 3 mL of *E. coli* polar lipids with 1 mL of phosphatidylcholine in a round bottom flask under a ventilated hood.

1.2. Dry the lipid mixture for 20 min under vacuum using a rotary evaporator and a water

bath at 37 °C to remove chloroform. If needed, dry further under nitrogen or argon gas.

1.3. Using TS buffer (20 mM Tris-HCl pH 8.0, 150 mM NaCl) containing 2 mM β-mercaptoethanol, resuspend lipids to 25 mg/mL.

1.4. Resuspend lipids in 500 μL aliquots, flash freeze in liquid nitrogen, and store at -80 °C.

1.5. Thaw one 500 μL aliquot of lipids and dilute 1:1 using TS buffer containing 2 mM β-mercaptoethanol.

1.6. Extrude the lipid suspension 15 times using a 400 nm membrane.

1.7. Dilute the lipid suspension to have a final lipid concentration of 4.4 mg/mL.

1.8. Add n-dodecyl-β-D-maltoside (DDM) to have a final concentration of 0.2% and leave it rotating for 1 h at 200 rpm at room temperature (RT).

1.9. Add the purified protein to the lipids using a lipid-to-protein ratio between 1:10 and 1:100 (w:w).

NOTE: The ratio needs to be adjusted depending on the strength of the signal detected in SSM-electrophysiology measurements of substrate transport (see below). To obtain larger signals, use smaller lipid-to-protein ratios.

1.10. Incubate the mixture rotating at 200 rpm for 1 h at RT.

1.11. Add 30 mg/mL of polystyrene adsorbent beads, pre-washed in TS buffer.

NOTE: Add polystyrene beads stepwise.

1.12. Incubate the beads-lipids mixture for 30 min at RT under slow stirring.

1.13. To remove the beads, let the beads-lipids mixture stand so that the beads settle down. Transfer the solution to a new tube and leave the beads behind. Add 30 mg/mL of fresh polystyrene adsorbent beads to the separated lipid mixture.

1.14. Incubate the mixture for 1 h at 4 °C under slow stirring.

1.15. Separate the beads from the mixture as described in step 1.13 and add 30 mg/mL of fresh polystyrene adsorbent beads.

1.16. Incubate the mixture for 16 h at 4 °C.

1.17. Separate the beads from the mixture as described in step 13 and add 30 mg/mL of fresh polystyrene adsorbent beads.

1.18. Incubate the mixture for 2 h at 4 °C for a fourth and final wash.

142  
143 1.19. Centrifuge at 110,000 x *g* for 30 min at 4 °C.

144  
145 1.20. Wash the pellet with 500 µL of TS buffer containing 2 mM β-mercaptoethanol.

146  
147 1.21. Centrifuge again at 110,000 x *g* for 30 min at 4 °C.

148  
149 1.22. Resuspend the pellet to a final lipid concentration of 25 mg/mL in TS buffer with 2 mM  
150 β-mercaptoethanol.

151  
152 1.23. Estimate the protein concentration using an in-gel or amido black assay<sup>36</sup>.

153  
154 1.24. Aliquot the proteoliposomes, flash freeze in liquid nitrogen, and store at -80 °C.

## 155 156 **2. Chip preparation**

157  
158 2.1. Fill a single sensor chip with 50-100 µL of 0.5 mM 1-octadecanethiol solution (resus-  
159 pended in isopropanol).

160  
161 2.2. Incubate the chip with the solution for 30 min at RT.

162  
163 2.3. Remove the thiol solution by tapping the chip on a tissue.

164  
165 2.4. Rinse the sensor 3 times with 5 mL of pure isopropanol.

166  
167 2.5. Rinse the sensor 3 times with 5 mL of double distilled water.

168  
169 2.6. Dry the sensor by tapping on a tissue paper.

170  
171 2.7. Apply 1.5 µL of 7.5 µg/µL 1,2-diphytanoyl-sn-glycero-3-phosphocholine (lipids dried in a  
172 rotatory evaporator and resuspended in n-decane).

173  
174 2.8. Immediately after, fill the sensor with 50 µL of non-activating SSM buffer, which does not  
175 contain the substrate. This will lead to a spontaneous formation of the SSM layer.

176  
177 NOTE: The SSM buffer should be optimized beforehand to determine optimal conditions that  
178 have low background noise. A general buffer containing 30 mM HEPES pH 7.4, 5 mM MgCl<sub>2</sub>,  
179 140 mM NaCl, can be used as a starting point. The SSM buffer without the substrate is used  
180 for washing before and after the measurement (non-activating buffer). To avoid buffer mis-  
181 match, use the non-activating buffer to prepare the buffer containing the substrate (activat-  
182 ing buffer). The substrate can be added either directly as a powder or in a small volume from  
183 a high concentration stock to avoid dilutions.

184  
185 2.9. Thaw proteoliposomes from step 1.24 at RT.

186  
187 2.10. Dilute proteoliposomes between 1:5 and 1:100 (proteoliposomes:buffer, (v:v)) in the  
188 non-activating SSM buffer (here 1:20).

2.11. Sonicate proteoliposomes for 20-30 s or 3 times for 10 s, placing on ice in between sonication, if necessary. Here a water bath sonicator at 45 kHz was used.

2.12. Apply 5-10  $\mu\text{L}$  of the diluted sonicated proteoliposomes sample on the surface of the sensor without touching it.

2.13. Centrifuge the chips with the solution at RT for 30 min using a speed between 2,000 and 3,000  $\times g$ .

NOTE: Use 50 mL tubes with a flat bottom. Carefully place the sensor chips upright using tweezers. 6-well plates and a centrifuge with a plate holder can also be used.

2.14. Use the sensor chips on the same day.

### 3. Measuring the solute transportation: determination of saturation conditions

NOTE: As proof-of-principle, these experiments were performed using a bacterial choline transporter reconstituted in liposomes following the protocol described above. The step-by-step process of determining saturating conditions of the substrate choline prior to the measurement of inhibition by nanobodies is shown here.

3.1. Prepare 1-2 L of the non-activating SSM buffer.

NOTE: Prepare and use the same SSM buffer stock for all activating and non-activating buffers throughout all measurements.

3.2. Take 10 clean tubes and transfer 10 mL of the non-activating SSM buffer into each.

3.3. Add the substrate into the tubes from step 3.2 using a series of concentrations around the expected half maximum concentration (here 15, 10, 5, 1, 0.5, 0.1, 0.05, 0.01, 0.005, 0.001 mM of choline) to prepare the activating SSM buffers. Use a high concentration stock to avoid dilutions.

3.4. Switch on the SSM machine.

3.5. Start the SSM software and let the machine initialize automatically. Set the saving path for data and confirm by hitting the **OK** button. Select the standard **Initial Cleaning Protocol** in the **workflow** options and click **Run**.

3.6. Mount the proteoliposome coated chip on the socket, move the arm to lock the chip, and close the mounted chip with the cap.

3.7. Select the program **CapCom** in the workflow and let it **Run** to determine the conductivity and capacitance. Confirm that the conductivity is below 5 nS and the capacitance is between 15 and 35 nF before using it for the measurement.

NOTE: A capacitance value of 15-35 nF and conductance below 5 nS are recommended by the manufacturer when using a 3 mm chip.

3.8. Transfer the activating solutions into vials and position the buffers in the probe sampler.

3.9. Transfer the non-activating buffer into a reservoir and position it next to the chip holder at the reservoir position on the right.

3.10. Create a protocol for the workflow using a sequence of non-activating (B), activating (A), and non-activating (B) solutions (B-A-B sequence) and a loop that performs three measurements and moves to the next activating buffer for all 10 buffers prepared in step 3.3. Use the default flow rate at 200  $\mu$ L/s using 1 s – 1 s – 1 s flow times for the B-A-B sequence. Click **Play** to start the measurement.

NOTE: A typical experiment consists of the sequential flow of non-activating (B), activating (A), and non-activating (B) solutions (written as B-A-B; see **Figure 1A**). The immobilized proteoliposomes on the sensor will be washed with the solutions. Therefore, the B-A solution exchange generates a substrate concentration gradient, which drives the electrogenic transport reaction.

3.11. Save the protocol and let the workflow run by clicking on the **Play** button. Perform the same type of experiment but using protein-free liposomes. This is highly important as it would show the intensity of background currents. This should be considered when analyzing data of electrogenic transport measured with proteoliposomes (**Figure 1C**).

3.12. Use any preferred software for data analysis to plot the measured current versus time. Read out the peak current manually, or if using the software, use the function for peak height estimation in the range of the addition of the activating buffer.

3.13. Plot the peak current against the substrate concentration to determine the EC<sub>50</sub> of the substrate via the nonlinear regression (**Figure 1B,C**). Read out the lowest concentration at which the peak current reaches a maximum value, this concentration corresponds to saturating conditions.

NOTE: It is important to consider that the number of proteoliposomes that remain immobilized vary from chip-to-chip. This variation is evident as the peak currents at identical measurement conditions will show different heights. Therefore, it is necessary to normalize the current amplitudes of measurements performed on each chip separately before comparing measurements among different chip preparations.

#### 4. Serial classification of inhibitory and non-inhibitory nanobodies

NOTE: This section shows how to measure choline transport in the presence of nanobodies that bind specifically to the bacterial choline transporter. Smaller peak currents in the presence of nanobodies indicate transport inhibition. Non-inhibitory nanobodies will not impact substrate transport, i.e., no decrease of the peak current signal.

283 4.1. Prepare 1-2 L of non-activating SSM buffer.

284  
285 4.2. Transfer 50 mL of the non-activating SSM buffer into a clean tube. Add the substrate  
286 choline to a final concentration of 5 mM (saturating conditions). Use this for a positive control  
287 measurement.

288  
289 4.3. Transfer 10 mL of the non-activating SSM buffer into a clean tube. Add the substrate  
290 choline to a final concentration of 5 mM (saturating conditions) and add nanobody to a final  
291 concentration of 500 nM.

292  
293 4.4. Repeat step 4.3 for each nanobody to prepare the activating solutions.

294  
295 NOTE: If the purified nanobodies are resuspended in a different buffer than the SSM buffer,  
296 their addition to the activating and non-activating buffers will lead to a buffer mismatch. A  
297 buffer mismatch should be avoided since it can lead to high noise. Exchanging the buffer of  
298 the purified nanobodies with the SSM buffer can help to avoid this issue. Furthermore, using  
299 nanobody concentrations that allow to reach saturating conditions is recommended. Consid-  
300 ering that the binding constants of nanobodies are generally below 100 nM, a nanobody con-  
301 centration of 500 nM is recommended for this experiment. However, it is important to pre-  
302 screen for optimal concentrations.

303  
304 4.5. Start the SSM machine and measure the capacitance and conductivity of the proteolipo-  
305 some coated chip as described in steps 3.4-3.7.

306  
307 4.6. Transfer the activating solution without a nanobody into a vial and place the buffer in the  
308 probe sampler. Transfer the non-activating buffer without a nanobody into a reservoir and  
309 position it in the probe sampler.

310  
311 4.7. Transfer the activating solutions containing nanobodies into vials and position the buffers  
312 in the probe sampler. Transfer the non-activating buffers containing nanobodies into vials  
313 and position the buffers in the probe sampler.

314  
315 4.8. Create a protocol for the workflow using a sequence of non-activating (B), activating (A),  
316 and non-activating (B) solutions (B-A-B sequence).

317  
318 4.9. Create a loop that performs the following: three measurements of the B-A-B sequence  
319 using buffers without nanobody, two measurements of the B-A-B sequence with buffers con-  
320 taining a nanobody, 120 s delay time for incubation with the nanobody, then 3 measurements  
321 of the B-A-B sequence with buffers containing the nanobody.

322  
323 4.10. Save the workflow and let it run by clicking the **Play** button.

324  
325 NOTE: This workflow will measure the initial conditions of the transport without inhibition by  
326 running the B-A-B protocol 3 times using non-activating (B) and activating buffers (A) without  
327 the nanobody (**Figure 1B,C**), followed by the measurement of the nanobody effect on the  
328 transport by running the B-A-B protocol 5 times with the non-activating and activating buffers  
329 containing nanobodies (**Figure 2A**). The second order binding kinetics dictate the interaction



of nanobodies and their target proteins. Therefore, it is important to use a time delay in the B-A step in order to give enough time for nanobodies to bind to transporters in proteoliposomes on the chip. Optimal times depend on the nanobody concentration i.e., at lower concentrations longer times are required. The first two measurements are required to adapt the system to the new conditions and the second run should be performed after a delay time of 120 s. Only the following three measurements should be used for data analysis.

4.11. Create a new protocol for the workflow using a sequence of non-activating (B), activating (A), and non-activating (B) solutions (B-A-B sequence) and a loop of 5 measurements to wash out the reversibly bound nanobody.

NOTE: Optionally, include an incubation step to allow dissociation of nanobodies with slow kinetics.

4.12. Save the workflow and let it run by clicking the **Play** button.

4.13. Compare the last peak current of the measurements with the initial substrate-only measurement in step 4.10. The nanobody has been successfully washed out and the initial conditions have been reestablished if the peak current reaches the initial value, otherwise repeat steps 4.12-4.13 or change to a new chip.

4.14. Repeat steps 4.6-4.13 and use individual chips for each nanobody screen (**Figure 2C**) or repeat with multiple nanobodies using the same chip (**Figure 2D**).

4.15. Use any preferred software for data analysis to plot the measured current versus time. Read out the peak current manually, or if available, in the used software, automatically select the function for peak height estimation in the range of the addition of the activation buffer.

4.16. Normalize the peak current in presence of the nanobody, based on the preceding substrate-only measurement. Plot the peak currents in a histogram and compare the peak currents of the substrate only measurements to the peak currents measured in the presence of nanobodies (**Figure 2C,D**) to identify inhibitory nanobodies.

NOTE: Normalization of the determined peak currents of each individual run with a nanobody should be performed considering the peak current in the absence of the nanobody from the preceding measurement. Also, since the number of proteoliposomes that remain immobilized vary from chip-to-chip, it is important to normalize the current amplitudes of measurements performed on each chip separately before comparing measurements among different chip preparations.

## 5. IC<sub>50</sub> measurement with inhibitory nanobodies

NOTE: After identifying inhibitory nanobodies, it is possible to determine their half maximal inhibitory concentration (IC<sub>50</sub>). This is done by measuring the transport of choline at constant concentration, while varying concentrations of the inhibitory nanobody.

5.1. Prepare 1-2 L of the non-activating SSM buffer.

5.2. Transfer 50 mL of the non-activating SSM buffer into a clean tube. Add the substrate choline to a final concentration of 5 mM (saturating conditions). Used this as the activating solution for positive control.

5.3. Take 8 clean tubes and add 5 mL of the non-activating solution into each. Add the substrate choline to a final concentration of 5 mM (saturating conditions). Add the inhibitory nanobody to the tubes at concentrations in the expected  $IC_{50}$  range (here 500 nM – 1 nM).

5.4. Take 8 clean tubes and add 10 mL of the non-activating solution into each. Add inhibitory nanobody to each tube individually at the same concentration as in step 5.3. This corresponds to the non-activating buffer.

NOTE: This will generate a series of activating and non-activating buffer pairs at different concentrations of the same inhibitory nanobody.

5.5. Start the SSM setup and measure the capacitance and conductivity of the proteoliposome coated chip as described in steps 3.4-3.7.

5.6. Transfer the activating solution without nanobody into a vial and place it in the probe sampler. Transfer the non-activating buffer without nanobody into a reservoir and position it at the reservoir position next to the chip holder on the right.

5.7. Transfer the activating solutions containing nanobodies into vials and position the buffers in the probe sampler. Transfer the non-activating buffers containing nanobodies into vials and position the buffers in the probe sampler.

5.8. Create a protocol for the workflow using a sequence of non-activating (B), activating (A), and non-activating (B) solutions (B-A-B sequence). Include a loop to measure each concentration 2 times, incubate for 120 s and measure 3 more times. The workflow will start with the positive control of substrate-only followed by the lowest concentration of the nanobody. Each nanobody measurement will be followed by a subsequent measurement of the positive control with substrate-only to restore the initial peak amplitude, before moving to the next higher nanobody concentration.

NOTE: The second order kinetics dictate the binding of nanobodies. Therefore, it is important to use a time delay in the B-A step. Optimal times depend on the nanobody concentration i.e., at lower concentrations longer times are required; here 120 s was used with satisfactory results.

5.9. Use any preferred software for data analysis to plot the measured current versus time. Read out the peak current manually, or if using software, select the function for the peak height estimation in the range of the addition of the activation buffer (**Figure 3A**). Plot the peak currents against the nanobody concentration to determine the  $IC_{50}$  via non-linear regression (**Figure 3B**).

NOTE: Normalize the current amplitudes of measurements performed for each individual chip

before comparing measurements among different chip preparations.

## **6. Cleaning of sensors**

6.1. Rinse the single sensor chips after use with 10 mL of distilled water.

6.2. Dry the chip by tapping it on a tissue paper.

6.3. Fill the sensor cavity of the chip with 100  $\mu$ L of pure isopropanol and incubate for 10 min at RT.

6.4. Place a cotton swab in pure isopropanol and incubate for 1-3 min.

6.5. Use the presoaked cotton swabs and gently rotate on the sensor's surface without pressure to remove residues.

6.6. Rinse the sensor with 5 mL of pure isopropanol.

6.7. Rinse the sensor with 10 mL of distilled water.

6.8. Dry the sensor by tapping the chip on a tissue.

6.9. Let the sensor dry overnight at RT, and store at RT under dry conditions.

NOTE: Sensors can be re-used up to 4-5 times when cleaned and stored properly.

## **REPRESENTATIVE RESULTS**

SSM-based electrophysiology has been extensively used for the characterization of electrogenic transporters. In the protocol presented here, we show how to use SSM-based electrophysiology to classify nanobodies targeting a secondary transporter (here a bacterial choline symporter) based on their inhibitory and non-inhibitory properties. One of the most useful features of this technique is that it allows for the high-throughput screening of multiple buffer conditions. This particular characteristic is beneficial for the analysis of nanobody libraries, which after the selection of binders can be constituted from a few to dozens of nanobodies. In a standard experiment, a stable lipid monolayer is assembled on a sensor chip. After applying the proteoliposomes preparation containing the choline transporter, a check for the good conductivity and capacitance is performed as this is essential for the success of the experiment. In case that the integrity of the membrane is compromised during an experiment, which is easily observed due to the high noise background currents, changing to a new chip is recommended as recovering low noise conditions is rather difficult. In general, we have observed very good reproducibility among measurements of transport and inhibition by nanobodies when using different chips.

To decide about the substrate concentration to be used during a screening of nanobodies, electrogenic transport was first measured under different substrate concentrations to determine EC<sub>50</sub> (**Figure 1B,C**). A substrate concentration that corresponds to saturating conditions

was selected (**Figure 1C**). This substrate concentration was then kept constant in all activating buffers. For this particular example, we selected 5 mM choline.

For the screening of nanobodies, the nanobody must be added to both non-activating and activating buffers. When nanobodies were added to only the activating buffer, it was not possible to observe the inhibition of the electrogenic transport. We speculate that this is due to an incomplete occupation of all nanobody binding sites in the transporter population on the chip, thereby revealing the importance of pre-incubation with nanobodies in non-activating conditions. To ensure that all sites are likely to be occupied, a time delay step was included during the application of the first non-activating buffer step to allow the saturation of nanobody binding sites on the transporter population. Incubation times ranging from 2-60 min have been tested with reproducible results. Keep in mind that optimal times of incubation depend on the nature of the nanobody binder and its concentration during the experiment (as well as the concentration of transporter in proteoliposomes on the chip). Therefore, it is recommended to try different incubation times. In any case, as a rule of thumb, the lower the nanobody concentration, the longer the incubation time required. We tested incubation times of 2 min, 20 min, 30 min, and 60 min for different nanobodies but did not detect further transport inhibition.

The effect of inhibitory nanobodies on electrogenic transport is visualized from the decrease of peak currents amplitudes (**Figure 2A,C,D**). Non-inhibitory nanobodies, on the other hand, do not affect peak currents. After running the washing protocol to allow nanobodies unbinding, a recovery of 80 to 95% of the initial peak current amplitude was observed (**Figure 2A,C,D**). We have performed a similar experiment but in the presence of liposomes without the transporter protein. When changing from non-activating to activating conditions, no significant artifact currents was introduced by nanobodies present in these buffers (**Figure 2B**). Running this control experiment is recommended as it is important to know whether changes in peak currents arise from artifacts or not.

After the selection of nanobodies with inhibitory properties, we determined  $IC_{50}$  values for individual nanobodies (**Figure 3A,B**). For this particular experiment, it is recommended to start with a low concentration of nanobodies and then move towards the high concentration during the assay. The calculation of the inhibition for each concentration was then performed by comparing peak currents measured before and after the application of nanobody. To avoid unspecific binding of nanobodies to surfaces, which can be particularly problematic when using low nanobody concentrations, it is advised to follow a similar protocol to that described by Kermani et al.<sup>37</sup>, where 50  $\mu$ g/mL of bovine serum albumin was added to the buffers, preventing this deleterious effect. Adding detergents such as Tween or Triton for this purpose should be avoided as these would dissolve lipid membranes.

## FIGURES AND TABLES LEGENDS

**Figure 1: SSM-based electrophysiology.** (A) Protocol for transient currents measurement. A non-activating solution is replaced by an activating solution followed by the flow of non-activating solution to restore initial conditions. During the first step, nanobodies bind to the transporter. When switching to the activating solution, the substrate gradient drives the electrogenic transport (orange curve). In the presence of an inhibitory nanobody, the peak current

shows a smaller amplitude (blue curve). After finishing the protocol and running solutions without nanobody (wash), unbinding of nanobodies occurs. In the schematic, proteoliposomes with reconstituted protein (blue) are immobilized on the SSM sensor. Triangles and red circles represent nanobodies and substrate, respectively. (B) Electrogenic choline transport in the absence of nanobodies. Peak currents measured during activating conditions are shown for different substrate concentrations. (C) Representative measurement of currents during activating conditions in the absence of transporter protein at different substrate concentrations. (D) Plot of substrate concentration versus peak currents amplitude. The  $EC_{50}$  determined was  $95 \pm 11 \mu M$  choline. Error bars indicate standard deviation ( $n=3$  biological replicates,  $n=3$  technical replicates).

**Figure 2: Screening and classification of inhibitory and non-inhibitory nanobodies.** (A) Electrogenic choline transport in the presence of a nanobody. Peak currents measured during activating conditions are shown in the absence of nanobody (blue), in the presence of an inhibitory nanobody (red), and after nanobody unbinding (green). (B) Measurement of currents during activating conditions in the absence of transporter protein. Traces show recordings in the absence of nanobody (blue), in the presence of an inhibitory nanobody (green), and in the presence of a non-inhibitory nanobody (red). (C,D). Histograms showing peak currents measured during activating conditions in the presence of nanobodies and after nanobody unbinding (recovery). Panel C shows the results of measurements using individual chips per nanobody. Panel D shows the results from a serial measurement using one chip. Nanobodies are indicated as Nb. Error bars indicate the standard deviation ( $n=3$  biological replicates,  $n=2$  technical replicates).

**Figure 3: Determination of  $IC_{50}$  of an inhibitory nanobody.** (A) Electrogenic choline transport and inhibition by a nanobody. Peak currents measured during activating conditions are shown for different nanobody concentrations. (B) Plot of peak currents amplitude vs nanobody concentration from a serial measurement with an inhibitory nanobody. The  $IC_{50}$  determined was  $18 \pm 2 nM$ . Error bars indicate the standard deviation ( $n=3$  biological replicates,  $n=3$  technical replicates).

## DISCUSSION:

The technique presented here classifies nanobodies with inhibitory and non-inhibitory properties targeting electrogenic transporters. Assessing the substrate transport is possible due to the detection of the movement of charges through the transporter embedded in the membrane of proteoliposomes. Some of the critical steps during the setup of an experiment are reconstitution of active protein in liposomes, preparation of stable monolayers on SSM chips, and recovering of initial conditions after the application of the wash protocol to remove bound nanobody molecules. Once the membrane protein is reconstituted at an appropriate lipid-to-protein ratio, a general SSM protocol can be established using the native substrate. It is crucial to perform control experiments using protein-free liposomes to reveal noise currents that would need to be subtracted from the currents measured using proteoliposomes. However, if noise currents are too large, we recommend trying protein reconstitution using different lipids, or screen for buffer conditions that minimize these deleterious signals. After successfully establishing conditions for an SSM assay, screening of nanobodies can be performed. A very useful option that may help in faster screening of nanobodies is to perform high-throughput assays using a single SSM sensor chip. This reduces the time of manipulation

of chips and buffers and reduces the costs. However, because during this type of assay multiple nanobodies are applied sequentially, it is important to ensure that the applied nanobody can be washed away after the measurement. A stringent washing cycle may need to be implemented to unbind some nanobodies in case that reduced peak current amplitudes are detected in the absence of nanobody. We recommend using, as the starting point, the washing conditions described here. If increasing the washing volume or the number of cycles does not help, individual chips would need to be used to screen each nanobody separately. In all the cases examined here, the binding of nanobodies was reversible and a high-throughput protocol could be applied. In our experimental setup, we could not recover the full amplitude of the initial peak current after measurements with nanobodies (**Figure 1A; Figure 2C,D**). However, in most cases, the magnitude of the peak currents recovered ranged between 80 and 95% of the amplitude measured before applying nanobodies (**Fig. 2C,D**). We speculate that this could be a consequence of washing away a fraction of the proteoliposomes adhered to the chip, or due to slow kinetics of unbinding of some inhibitory nanobodies, or a combination of both. In either case, it was still possible to continue assaying further nanobodies as electrogenic transport was measurable. This is shown in **Figure 2D**, where we screened six nanobodies using a single chip preparation.

The high-throughput characteristic is one of the most significant advances of the method presented. In addition, in contrast to other approaches, this method allows for the selection of inhibitory nanobodies targeting electrogenic transporters for which labeled substrates are not available. A fast identification of inhibitory nanobodies can help to speed up the research aiming to identify novel applications of nanobodies as drugs. Their apparent advantages compared to similar therapies such as antibody treatments are numerous, starting with a smaller size which helps them propagate further into tissues or cells, to their low production costs and high stability.

SSM-based electrophysiology has been used in the past for the characterization of electrogenic transporters in membrane vesicles<sup>29,38</sup>. These types of experiments are advantageous as they do not depend on the protein purification and reconstitution protocols. We speculate that performing the selection of inhibitory nanobodies using membrane vesicles is feasible. This would help to reduce costs and avoid manipulations of purified proteins.

SSM-based electrophysiology is a strong technique to screen for nanobody inhibitors of multiple membrane proteins that exhibit electrogenic transport. We envision that SSM-based electrophysiology will become an important tool for the selection of inhibitory nanobodies and other antibodies with potential clinical applications.

## ACKNOWLEDGEMENT

We thank Cedric A. J. Hutter and Markus A. Seeger from the Institute of Medical Microbiology at the University of Zurich, and Gonzalo Cebrero from Biozentrum of the University of Basel for collaboration in the generation of synthetic nanobodies (sybodies). We thank Maria Barthmes and Andre Bazzone from NANION Technologies for technical assistance. This work was supported by the Swiss National Science Foundation (SNSF) (PP00P3\_170607 and NANION Research Grant Initiative to C.P.).

## DISCLOSURES

The authors declare no competing financial interests.





## REFERENCES





1. Braden, B. C., Goldman, E. R., Mariuzza, R. A., Poljak, R. J. Anatomy of an antibody molecule: structure, kinetics, thermodynamics, and mutational studies of the antilysozyme antibody D1.3. *Immunology Reviews*. **163**, 45-57 (1998).
2. Hamers-Casterman, C. et al. Naturally occurring antibodies devoid of light chains. *Nature*. **363**, 446-448 (1993).
3. Perez, C. et al. Structural basis of inhibition of lipid-linked oligosaccharide flippase PglK by a conformational nanobody. *Science Reports*. **7**, 46641 (2017).
4. Grahl, A., Abiko, L.A., Isogai, S., Sharpe, T., Grzesiek, S. A high-resolution description of beta1-adrenergic receptor functional dynamics and allosteric coupling from backbone NMR. *Nature Communication*. **11**, 2216 (2020).
5. Schenck, S. et al. Generation and characterization of anti-VGLUT nanobodies acting as inhibitors of transport. *Biochemistry*. **56**, 3962-3971 (2017).
6. Mireku, S. A., Sauer, M. M., Glockshuber, R., Locher, K.P. Structural basis of nanobody-mediated blocking of BtuF, the cognate substrate-binding protein of the Escherichia coli vitamin B12 transporter BtuCD. *Science Reports*. **7**, 14296 (2017).
7. Manglik, A., Kobilka, B. K., Steyaert, J. Nanobodies to study G protein-coupled receptor structure and function. *Annual Reviews of Pharmacology Toxicology*. **57**, 19-37 (2017).
8. Rasmussen, S. G. et al. Structure of a nanobody-stabilized active state of the beta(2) adrenoceptor. *Nature*. **469**, 175-80 (2011).
9. Jiang, X. et al. Crystal structure of a LacY-nanobody complex in a periplasmic-open conformation. *Proceeding of the National Academy of Science U. S. A.* **113**, 12420-12425 (2016).
10. Geertsma, E. R. et al. Structure of a prokaryotic fumarate transporter reveals the architecture of the SLC26 family. *Nature Structural Molecular Biology*. **22**, 803-808 (2015).
11. Harmsen, M. M., De Haard, H. J. Properties, production, and applications of camelid single-domain antibody fragments. *Applied Microbiology and Biotechnology*. **77**, 13-22 (2007).
12. Pardon, E. et al. A general protocol for the generation of nanobodies for structural biology. *Nature Protocols*. **9**, 674-93 (2014).
13. Nguyen, V. K., Desmyter, A., Muyldermans, S. Functional heavy-chain antibodies in Camelidae. *Advances in Immunology*. **79**, 261-296 (2001).
14. Zimmermann, I. et al. Synthetic single domain antibodies for the conformational trapping of membrane proteins. *Elife* **7**, e34317(2018).
15. McMahon, C. et al. Yeast surface display platform for rapid discovery of conformationally selective nanobodies. *Nature Structural Molecular Biology*. **25**, 289-296 (2018).
16. Olichon, A., de Marco, A. Preparation of a naive library of camelid single domain antibodies. *Methods in Molecular Biology*. **911**, 65-78 (2012).
17. Moutel, S. et al. NaLi-H1: A universal synthetic library of humanized nanobodies providing highly functional antibodies and intrabodies. *Elife* **5**, e16228 (2016).
18. Yan, J., Li, G., Hu, Y., Ou, W., Wan, Y. Construction of a synthetic phage-displayed nanobody library with CDR3 regions randomized by trinucleotide cassettes for diagnostic applications. *Journal of Translational Medicine*. **12**, 343 (2014).



19. Sabir, J. S. et al. Construction of naive camelids VHH repertoire in phage display-based library. *Comptes Rendus Biologies*. **337**, 244-249 (2014).
20. Yau, K. Y. et al. Selection of hapten-specific single-domain antibodies from a non-immunized llama ribosome display library. *Journal of Immunology Methods*. **281**, 161-175 (2003).
21. van der Linden, R.H. et al. Comparison of physical chemical properties of llama VHH antibody fragments and mouse monoclonal antibodies. *Biochimica et Biophysica Acta*. **1431**, 37-46 (1999).
22. Dumoulin, M. et al. Single-domain antibody fragments with high conformational stability. *Protein Science*. **11**, 500-515 (2002).
23. Iezzi, M. E., Policastro, L., Werbach, S., Podhajcer, O., Canziani, G. A. Single-domain antibodies and the promise of modular targeting in cancer imaging and treatment. *Frontiers in Immunology*. **9**, 273 (2018).
24. Yu, X. et al. Nanobodies derived from camelids represent versatile biomolecules for biomedical applications. *Biomaterials Science*. **8**, 3559-3573 (2020).
25. Jahnichen, S. et al. CXCR4 nanobodies (VHH-based single variable domains) potentially inhibit chemotaxis and HIV-1 replication and mobilize stem cells. *Proceedings of the National Academy of Science U. S. A.* **107**, 20565-70 (2010).
26. Nguyen, V. S. et al. Inhibition of type VI secretion by an anti-TssM llama nanobody. *PLoS One*. **10**, e0122187 (2015).
27. Bazzone, A., Barthmes, M., Fendler, K. SSM-Based Electrophysiology for Transporter Research. *Methods in Enzymology*. **594**, 31-83 (2017).
28. Schulz, P., Garcia-Celma, J. J., Fendler, K. SSM-based electrophysiology. *Methods*. **46**, 97-103 (2008).
29. Barthmes, M., Liao, J., Jiang, Y., Bruggemann, A., Wahl-Schott, C. Electrophysiological characterization of the archaeal transporter NCX\_Mj using solid supported membrane technology. *Journal of General Physiology*. **147**, 485-96 (2016).
30. Watzke, N., Diekert, K., Obrdlik, P. Electrophysiology of respiratory chain complexes and the ADP-ATP exchanger in native mitochondrial membranes. *Biochemistry*. **49**, 10308-10318 (2010).
31. Zuber, D. et al. Kinetics of charge translocation in the passive downhill uptake mode of the Na<sup>+</sup>/H<sup>+</sup> antiporter NhaA of Escherichia coli. *Biochim Biophys Acta*. **1709**, 240-250 (2005).
32. Garcia-Celma, J. J., Smirnova, I. N., Kaback, H. R., Fendler, K. Electrophysiological characterization of LacY. *Proceedings of the National Academy of Science U. S. A.* **106**, 7373-7378 (2009).
33. Bazzone, A., Madej, M. G., Kaback, H. R., Fendler, K. pH regulation of electrogenic sugar/H<sup>+</sup> symport in MFS sugar permeases. *PLoS One*. **11**, e0156392 (2016).
34. Williamson, G. et al. A two-lane mechanism for selective biological ammonium transport. *Elife*. **9**, e57183 (2020).
35. Mirandela, G. D., Tamburrino, G., Hoskisson, P. A., Zachariae, U., Javelle, A. The lipid environment determines the activity of the Escherichia coli ammonium transporter AmtB. *FASEB Journal*. **33**, 1989-1999 (2019).
36. Kaplan, R. S., Pedersen, P. L. Determination of microgram quantities of protein in the presence of milligram levels of lipid with amido black 10B. *Annals of Biochemistry*. **150**, 97-104 (1985).




- 703 37. Kermani, A. A. et al. The structural basis of promiscuity in small multidrug resistance  
704 transporters. *Nature Communication*. **11**, 6064 (2020).
- 705 38. Weitz, D. et al. Functional and structural characterization of a prokaryotic peptide  
706 transporter with features similar to mammalian PEPT1. *Journal of Biological Chemistry*. **282**,  
707 2832-9 (2007).







Name of Material/Equipment	Company	Catalog Number	Comments/Description
1-octadecanethiol solution	Sigma Aldrich	O1858-25ML	
1,2-diphytanoyl-sn-glycero-3-phosphocholine	Avanti Polar Lipids	850356C-25mg	
Bio-Beads SM-2 Adsorbent (Polystyrene adsorbent beads)	BioRad	#152-3920	
PD 10 Desalting Columns	GE Healthcare	GE17-0851-01	
Filter 200 nm membrane	Whatman	WHA800282	
2-Propanol	Nucleopore		
n-Decane	Merck	33539-1L-R	
	Sigma Aldrich	8034051000	
n-dodecyl- $\beta$ -D-maltoside (DDM)	Avanti Polar Lipids	850520P-25g	
Sodium Chloride	AppliChem	131659.1211	
(SSM setup) SURFE2R N1	Nanion	-----	
SURFE2R N1 Single Sensor Chips	Nanion	# 161001	
Trizma Base	Sigma Aldrich	T1503	



**University  
of Basel**

Department  
Biozentrum

**BIOZENTRUM**

The Center for  
Molecular Life Sciences

**Manuscript Title:**

Selection of transporter-targeted inhibitory nanobodies by SSM-based electrophysiology

JoVE62578


Dear Dr. Amit Krishnan,  
Review Editor  
JoVE

Thank you for your Email of March 9th, 2021, in which you invited us to revise our manuscript. We appreciate the insightful comments of the reviewers. We have responded to all the comments made by the editor and reviewers. The manuscript file contains track changes of the edits and a point by point response is found here below. We hope you find our manuscript acceptable after these modifications. Please do not hesitate to contact me in case of further comments or questions.

I look forward to hearing from you and I thank you for your efforts and consideration.

Yours sincerely,

Camilo Perez



### **Editorial comments:**


1. Please take this opportunity to thoroughly proofread the manuscript to ensure that there are no spelling or grammar issues.
2. Please revise the text to avoid the use of any personal pronouns (e.g., "we", "you", "our" etc.).
3. Please define all abbreviations before use (TS, DDM, etc.)
4. JoVE cannot publish manuscripts containing commercial language. This includes trademark symbols (™), registered symbols (®), and company names before an instrument or reagent. Please remove all commercial language from your manuscript and use generic terms instead. All commercial products should be sufficiently referenced in the Table of Materials. For example: Bio-Beads, SURFE2R N1, Nanion, etc.
5. Please adjust the numbering of the Protocol to follow the JoVE Instructions for Authors. For example, 1 should be followed by 1.1 and then 1.1.1 and 1.1.2 if necessary. Please refrain from using bullets or dashes.
6. Line 70: Please elaborate the process of drying lipids using Rotary evaporator. Please include specific requirements to perform this drying technique.
7. For time units, please use abbreviated forms for durations of less than one day when the unit is preceded by a numeral. Do not abbreviate day, week, month, and year. Examples: 5 h, 10 min, 100 s, 8 days, 10 weeks (Line 78, 80, 86, 118, 201, etc.)
8. For SI units, please use standard abbreviations when the unit is preceded by a numeral. Abbreviate liters to L to avoid confusion. Examples: 10 mL, 8 µL, 7 cm<sup>2</sup> (Line 81, 83, 85, 90, etc.)
9. Line 91: Please mention the conditions for centrifugation. Are they the same as in line 89?
10. Line 103/104: Please mention if there is any specific volume required for the solution used.
11. Line 114: Please mention if there is any specific frequency used for sonication.
12. Please include a one-line space between each protocol step and highlight up to 3 pages of the Protocol (including headings and spacing) that identifies the essential steps of the protocol for the video, i.e., the steps that should be visualized to tell the most cohesive story of the Protocol. Remember that non-highlighted Protocol steps will remain in the manuscript, and therefore will still be available to the reader.
13. As we are a methods journal, please revise the Discussion to include modifications/troubleshooting steps of the technique and limitations of the protocol.
14. Line 278-307: Please move the Figure Legends section to the end of the Representative Results section.
15. Please do not use the &-sign or the word "and" when listing authors. Please title case and italicize journal titles and book titles. Do not use any abbreviations. Article titles should start with a capital letter and end with a period and should appear exactly as they were published in the original work, without any abbreviations or truncations.

[Thanks for pointing this out. All these comments were addressed](#)

### **Reviewer #1:**

#### ***Major Concerns:***

Introduction is well written. I just note that there 2 impactful papers that have revised an old



transport system using SSME and lead to ground-breaking results. These papers are cited in the Nanion website, and in my opinion should be cited in the introduction.

1 - Williamson G, Tamburrino G, Bizior A, Boeckstaens M, Dias Mirandela G, Bage MG, Pislakov A, Ives CM, Terras E, Hoskisson PA, Marini AM, Zachariae U, Javelle A. A two-lane mechanism for selective biological ammonium transport. Elife. 2020 9:e57183. doi: 10.7554/eLife.57183

This paper should be cited after "This technique has been applied (...) turnover of these proteins".

2 - Mirandela GD, Tamburrino G, Hoskisson PA, Zachariae U, Javelle A. The lipid environment determines the activity of the Escherichia coli ammonium transporter AmtB. FASEB J. 2019 33:1989-1999. doi: 10.1096/fj.201800782R.

This paper should be cited as it presents a method to measure the effect of lipids on transport system, which have never been done before. It goes shows how to use SSME to measure external parameters on transport systems.

Both references have been added to the manuscript.

I am not sure what is the policy of the journal, but I find difficult to really assess and criticised the result as the transport system is not describe. It seems that the correct controls have been done, but again without knowing what transport system is use, it is difficult to evaluate the results.

Thanks for pointing this out. Initially, we didn't want to distract the readership from the main aim of the manuscript. We now mention that this is a "bacterial choline symporter" and specify that the substrate used in these experiments is "choline".

#### **Minor Concerns:**

Fig1 A. Should add "time (s)" on X axis.

This has been added

#### **Reviewer #2:**

We thank the reviewer for the very insightful comments.


#### **Major Comments:**

None. The overall structure and aim of this protocol is good. Please provide additional clarification/elaboration on points raised in the minor comments below.

#### **Minor Comments:**

General - is there a reason not to say what the transported substrate and transporter are? Even though the protocol is generally applicable, it seems unnecessarily evasive not to tell the reader what you're measuring in the example data shown.





Thanks for pointing this out. Initially, we didn't want to distract the readership from the main aim of the manuscript. We now mention that this is a "bacterial choline symporter" and specify that the substrate used in these experiments is "choline".

### Section 3.

1) The end of this section would be a good place to mention the negative controls that are essential for any SSM experiment (doing the experiment with empty liposomes). You might want to make a big deal out of it, for example "It is essential to measure a negative control when designing these assays..."

Thanks for pointing this out. A note on this important point was added.

### Section 4.

1) Step 1 - How much buffer is enough? For example, how much to go through the series of experiments with one antibody?

The required volume is now specified

2) Step 3 - why did you choose 500 nM for the antibody concentration? You may want to indicate a general range of concentrations typically seen for antibody binders.

This is now clarified

3) Step 9 -- This issue of the time delay before assessing binding is so important. I am wondering if you might want to elaborate on it a bit more. It would be more precise to describe the association kinetics (rather than "to interact better with proteoliposomes...") And mention that since nanobody association is a second order rate constant (rate is  $M^{-1}s^{-1}$ ), at very low (nM) concentrations of nanobody this delay might need to be extended.


This is now clarified in section 4.3 "...Note: It is recommended to use nanobody concentrations that allow to reach saturating conditions. Considering that the binding constants of nanobodies are generally below 100 nM, we suggest to use 500nM for this experiment. However, it is recommended to pre-screen for optimal concentrations."

4) Step 10 -- After the wash step, in some cases you might need to wait again before measuring the recovery currents, especially for tight binders with slow dissociation kinetics.

A note mentioning this was added

5) Step 13 -- Although you do elude to this, it would be helpful for people that are not familiar with the technique if you were explicit about the fact that you can't compare absolute current amplitudes from chip-to-chip since the vesicle fusions are not well-controlled.

Notes mentioning this were added



6) Frequently, when measuring protein-protein interactions, the proteins are "sticky" and adhere to surfaces, and this becomes especially problematic at lower antibody concentrations. How would an experimenter deal with non-specific binding? If we're assuming non-specialists, you may want to warn away from detergent additives people commonly use to prevent non-specific binding, like Tween/Triton. In my lab, we do similar SSM/antibody binding assays with addition of 50 ug/mL bovine serum albumin to combat this, and it works ok. (ie Kermani...Stockbridge Nat Comm, 2020. Please forgive the obnoxious insertion of my own work, and don't feel compelled to include the reference).

Thanks for this very helpful comment. Indeed, this is a very important point to take into account. We now mention this in lines 334-338 in the representative results section.

Section 5.

Step 5. Again, I think it's important to be precise that the binding kinetics are what dictate the length of the time delay. At low nM concentrations, longer time delays may be required since binding is a second order process (for example, at 1 nM, an antibody with a reasonable "on" rate of  $1 \times 10^6 \text{ M}^{-1}\text{s}^{-1}$  is going to have a time constant for association of >10 minutes).

A note mentioning this was added

Section 6. Cleaning of sensors

1) Can you re-use sensors after this type of experiment? How many times, typically?

A note clarifying this was added at the end of this section

Discussion.

Page 6/Line 217 - What do you mean by "nanobody libraries are generally large"? Do you mean after selection of binders to a target? Is large on the order of 10 or 100 candidates?

This has been now clarified in line 283

Page 6/Line 222 - Reproducibility of what, specifically?


This has been now clarified in line 289-290

Page 6/Line 229 - I would make a stronger statement than "highly recommend" - I would say the antibodies must be added to both buffers.

This chance was done

Page 6/Line 232 - Instead of "inefficient" occupation, maybe use "incomplete" or mention the kinetics again.

The changes were made and a note clarifying the importance of taking into account the kinetics was introduced in lines 315-321



Page 6/Line 237 - What do you mean by satisfactory results? Stable transporter fusions? (Also you might want to mention that you tried varying the length of the delay in the actual protocol, in section 4/step 3.)

This is now clarified in this part and in the notes introduced in the protocol

Page 7/Line 264 - Change "Inefficient unbinding" to "slow kinetics of unbinding"

This change was incorporated

Figures:

You show the negative control (empty liposomes) for the transported substrate in Figure 2, but you might want to put it in figure 1 since it's really important to do this at the early stages of the protocol.

Thanks for noticing this. Indeed, this is a very important point. A new figure 1C was added.