

# Journal of Visualized Experiments

## Evaluation of motor impairment in C. elegans models of amyotrophic lateral sclerosis (ALS)

--Manuscript Draft--

Article Type:	Invited Methods Collection - JoVE Produced Video
Manuscript Number:	JoVE62699R2
Full Title:	Evaluation of motor impairment in C. elegans models of amyotrophic lateral sclerosis (ALS)
Corresponding Author:	Nicole Liachko UNITED STATES
Corresponding Author's Institution:	
Corresponding Author E-Mail:	nliachko@uw.edu
Order of Authors:	Heather Currey Nicole Liachko
Additional Information:	
Question	Response
Please specify the section of the submitted manuscript.	Neuroscience
Please indicate whether this article will be Standard Access or Open Access.	Open Access (\$3900)
Please indicate the <b>city, state/province, and country</b> where this article will be <b>filmed</b> . Please do not use abbreviations.	Seattle, Washington, United States of America
Please confirm that you have read and agree to the terms and conditions of the author license agreement that applies below:	I agree to the <a href="#">Author License Agreement</a>
Please provide any comments to the journal here.	

**TITLE:**

Evaluation of Motor Impairment in *C. elegans* Models of Amyotrophic Lateral Sclerosis

**AUTHORS AND AFFILIATIONS:**

Heather N. Currey<sup>1</sup>, Nicole F. Liachko<sup>1,2</sup>

<sup>1</sup>Geriatrics Research Education and Clinical Center, Veterans Affairs Puget Sound Health Care System, Seattle, WA 98108, USA.

<sup>2</sup>Division of Gerontology and Geriatric Medicine, Department of Medicine, University of Washington, Seattle, WA 98104, USA.

**Author emails:**

Heather N. Currey (heather.currey@va.gov)

Nicole F. Liachko (nliachko@uw.edu)

**Corresponding Author:**

Nicole F. Liachko (nliachko@uw.edu)

**SUMMARY:**

This protocol describes two sensitive assays for discriminating among mild, moderate, and severe motor impairment in *C. elegans* models of amyotrophic lateral sclerosis, with general utility for *C. elegans* strains, with altered motility.

**ABSTRACT:**

The neurodegenerative disease amyotrophic lateral sclerosis (ALS) features progressive loss of motor neurons accompanied by muscle weakness and motor impairment that worsens with time. While considerable advances have been made in determining genetic drivers of ALS for a subset of patients, the majority of cases have an unknown etiology. Further, the mechanisms underlying motor neuron dysfunction and degeneration are not well understood; therefore, there is an ongoing need to develop and characterize representative models to study these processes. *Caenorhabditis elegans* can adapt their movement to the physical constraints of their surroundings, with two primary movement paradigms studied in a laboratory environment—crawling on a solid surface and swimming in liquid. These represent a complex interplay between sensation, motor neurons, and muscles. *C. elegans* models of ALS can exhibit impairment in one or both of these movement paradigms. This protocol describes two sensitive assays for evaluating motility in *C. elegans*: an optimized radial locomotion assay measuring crawling on a solid surface and an automated method for tracking and analyzing swimming in liquid (thrashing). In addition to the characterization of baseline motor impairment of ALS models, these assays can detect suppression or enhancement of the phenotypes from genetic or small molecule interventions. Thus, these methods have utility for studying ALS models and any *C. elegans* strain that exhibits altered motility.

**INTRODUCTION:**

Amyotrophic Lateral Sclerosis (ALS) is a debilitating, aging-related neurodegenerative disease

with a particular impact on motor neurons. The disease features loss of motor neurons in the brain and spinal cord and progressive motor impairment. This results in major functional disability and premature death, typically within 3-5 years after diagnosis<sup>1</sup>. Mutations in at least 38 genes can cause ALS; however, most patients with ALS accumulate ubiquitinated inclusions of the protein TDP-43 as their primary pathology in neurons and glial cells<sup>2-4</sup>. A number of animal models have been developed to study the underlying mechanisms causing or contributing to ALS *in vivo* (reviewed in <sup>5</sup>). In *C. elegans*, these models include genetic loss-of-function mutations in homologs of ALS-causing genes or transgenic expression of human ALS genes. There are numerous advantages to modeling ALS in *C. elegans*. *C. elegans* are a tractable simple animal with a differentiated nervous system, well-characterized behavioral paradigms, and considerable genetic homology to humans<sup>6,7</sup>. Many tools exist for working with *C. elegans*, including robust genome editing capabilities, *in vivo* fluorescent reporters of neurodegeneration, RNAi screening paradigms, tractable genetics, and established behavioral and phenotypic assays. *C. elegans* models of ALS recapitulate aspects of human disease, including accumulation of insoluble protein, neurodegeneration, and early death<sup>8,9</sup>. Further, motor dysfunction featuring disturbance of both crawling and swimming behaviors is present in many *C. elegans* ALS models.

This protocol describes two methods to characterize *C. elegans* motor phenotypes: the radial locomotion assay for evaluating crawling on a solid surface and the assessment of swimming in liquid (thrashing) using the WormLab automated tracking and analysis. These sensitive methods for characterizing motor deficits allow comparisons of severity and offer tools for measuring suppression and enhancements of motor phenotypes. The radial locomotion assay quantifies differences in crawling motility (sinusoidal movement on a solid surface) among populations of worms. This assay takes advantage of *C. elegans* natural unstimulated exploration behavior by placing worms in a single location on a plate and marking their final location after a given period of time<sup>10</sup>. Alternatively, swimming in liquid (thrashing) assays count body bends of individual worms over a set period of time. The manual counting of body bends by the human eye is time-intensive and typically exhibits considerable variability between experimenters. The use of computer-assisted automated tracking and analysis can eliminate much of that variability. In addition to the characterization of baseline motor impairment of ALS models, both radial locomotion and swimming assays can detect modulation of distinct locomotor phenotypes from genetic or small molecule interventions. These methods have utility for studying ALS models and any *C. elegans* strain that exhibits altered motility.

## PROTOCOL:

### 1. Radial locomotion assay

1.1. Prepare 100 mm or 150 mm diameter Petri dishes about half full of nematode growth media (NGM) and seeded with a uniform lawn of OP50 bacteria to cover the entire surface of the plate.

1.2. Flip the NGM assay plates upside-down and label the bottom with an identifier for the *C. elegans* strains to be assayed, 2 plates per strain (Figure 1A).

89  
90 1.3. Make a small dot with the marker in the center of the upside-down plate (Figure 1A).

91  
92 1.4. While working with a dissecting microscope, transfer 15-20 staged matched<sup>11</sup> worms to  
93 the center of the assay plate directly over the center dot (visible through the NGM) using a  
94 platinum-wire pick (Figure 1B).

95  
96 NOTE: Try to transfer all the worms at the same time or as quickly as possible.

97  
98 1.4.1. Set a timer for 30 min after worms are placed in the center of the plate. Put the lid back  
99 on the plate and set it aside (Figure 1C).

100  
101 1.5. Continue transferring worms until all strains are on the designated assay plates, making  
102 a note of approximately how long it takes to transfer between plates (aim for ~ 1 min between  
103 plates). Keep all plates in the same order that they were set up.

104  
105 1.6. After 30 min, begin scoring the first plate.

106  
107 1.6.1. Remove the lid and place the plate face down under the dissecting microscope, such  
108 that the labeled back of the plate is facing up, and the NGM agar is between the line of sight  
109 and the worms. The plate will be upside-down from the normal usage.

110  
111 1.6.2. Adjust the microscope focus until worms are visible through the agar.

112  
113 1.6.3. Using a different colored felt tip pen from the center point, put a small dot at the  
114 location of each worm - follow the worm tracks through the bacterial lawn to find the most  
115 distal animals.

116  
117 1.6.4. Check the edge of the plate as some worms may end up there. Also, count and record  
118 how many worms did not move from the center point (Figure 1D).

119  
120 1.7. Continue marking all assay plates in order, such that all plates have 30 min of activity  
121 time, being careful to account for the time it took to set up the plate. Then, perform manual or  
122 digital measurements as described in steps 1.8 or 1.9.

123  
124 1.8. Perform manual measurements using a ruler.

125  
126 1.8.1. Use a ruler to measure in mm, the distance from the center point to the final location  
127 markings for each worm and record the distance. To help keep track of progress during scoring,  
128 label the first scored point with a small marking line.

129  
130 1.8.2. Record length data for each dot consecutively by rotating the plate in a clockwise  
131 fashion (Figure 1E).

132



1.9. Perform digital measurement using a scanner and ImageJ.

1.9.1. Using a flatbed scanner, scan the back of the assay plate(s) next to a ruler (numbers are facing the scanning surface) (**Figure 2A**).

1.9.2. Open the scanned plate image in ImageJ or a similar program.

NOTE: ImageJ is a free Java-based image processing software hosted by the NIH.

1.9.3. Click on the **Straight-Line Tool** and then draw a line connecting the 1 cm and 2 cm points on the scanned image of the ruler (**Figure 2B**).

NOTE: Holding the **Shift** computer key will force the line to the nearest 45° angle.

1.9.4. Navigate to the **Analyze** menu and use it to set the scale (**Analyze | Set Scale**) (**Figure 2C**).

1.9.5. Enter the **Known Distance** (10) and **Units of Length** (mm) in the appropriate boxes, do not adjust the distance in pixels number. Check **Global** to apply the same scale to all images being analyzed in each ImageJ session; otherwise, set the scale for each image upon opening (**Figure 2D**). Click **Ok**.

NOTE: Steps 1.9.2—1.9.5 are performed to set the measurement scale. Measurements for this image will now automatically relate the length of lines to the defined pixel: length relationship. Images scanned at 300 dpi will have an approximate relationship of 11.8 pixels per mm.

1.9.6. Use the **Paintbrush Tool** to mark just above the first point to be scored.

NOTE: This is a visual indicator for the first worm scored.

1.9.7. Use the **Straight-Line Tool** to draw a line from the center point to the first worm data point.

1.9.8. Click the **M** key on the keyboard to measure the distance. Alternatively, navigate to the **Analyze** menu and choose **Measure (Analyze | Measure)**. In both cases, the measurement will appear in a new window; this window will collect all measurements per image.

1.9.9. Keeping the center point constant, move the end of the line touching the first worm point to the next worm point, moving clockwise, and click the **M** key to measure.

1.9.10. Continue measuring each point in a clockwise fashion until all points are scored.

1.9.11. **Copy** and **Paste** length measurement data into a spreadsheet to save it. Then, repeat the process for each scanned plate image.

NOTE: Steps 1.9.6 – 1.9.11 measure the worm displacement values (**Figure 2E-F**).

1.10. Perform statistical analysis.

1.10.1. Combine replicates within a single experiment so that a total number scored is between 30-40 worms.

1.10.2. Perform independent experimental replicates on different days and with independent populations of worms to end up with 3 independent experimental replicates.

1.10.3. Analyze data using appropriate statistical analyses such as Student's t-test for 2 strains or 1-way ANOVA for three or more strains.

## 2. Computer-analyzed swimming assay

NOTE: This protocol contains detailed instructions for the commercially available WormLab hardware and software system (see **Table of Materials**). However, the workflow can be applied to other computer-analyzed swimming assay systems.

### 2.1. Setting up and recording videos.

2.2. Raise the shielding on the imaging hardware (**Figure 3**) and verify that the camera height is set to the expected location.

NOTE: The optimal height for capturing worms on a 35 mm plate is such that the third from the bottom threaded screw bore is visible below the adjustable camera mounting bracket, but no additional mounting hole is seen above the screw bore (approximately 30.5 cm above the stage) (**Figure 3C**). This will result in a capture ratio of 11.43  $\mu\text{m}/\text{pixel}$ . Consider adding a "correct height" indicator so that this is easier to confirm prior to beginning the assay. See protocol step 2.18 to account for height adjustments.

2.3. Open the associated software (**Supplementary Figure 1A**).

2.4. Click the **Video Capture** icon (film with a right-pointing triangle in the middle) to open a new **Video Capture** window. The camera will be in the **Live View**, but the display will be black (**Supplementary Figure 1B**).

2.5. Press down on the dimmer knob and turn it clockwise to adjust the light.

NOTE: The dimmer knob is on the stage of the imaging system near the bright field light source. The **Video Capture** live view display will turn from black to illuminated.

2.6. Back in the **Video Capture** window, click on the **Settings** tab, adjust as follows: **Video mode**: 2456 x 2052\_Mono8, **Frame rate**: 14, **Output**: monochrome, **Exposure**: 0.00300 s, **Gain**: 1 dB, **Gamma**: 1, **Rotate 180** (check on).

2.7. Move back to the **Capture** tab. Select the recording folder (**Supplementary Figure 1B**).

NOTE: This is where all recordings and project files will be stored. Create unique folders for individual experiments to keep data organized.

2.8. Assign a file name by entering it into the **File prefix** text box.

NOTE: The software will automatically append a progressive number to the end of the prefix in the format "0001". It is advisable to have a consistent naming scheme for the files such as: "YYYYMMDD\_treatmentName\_rep#\_".

2.9. Set other **Capture** settings as follows: **Buffer**: 128 frames (default), **Duration** (check on): 00 h:01min:00 s. Ignore or turn off other settings in this tab.

2.10. Place the assay plate, lid up, on the device stage, and center it in the **Video Capture** screen; remove the lid.

NOTE: Prepare assay plates ahead of time. Each plate consists of a 35 mm Petri dish approximately half full of unseeded (no bacterial lawn) NGM.

2.11. Use a micropipette to wash the stage synchronized<sup>11</sup> worms in 1 mL of M9 onto the assay plate – aim for about 50 worms (**Figure 3D**). Gently swirl the plate to bring the animals to the center or use a micropipette to add a few drops of M9 to separate animals (**Figure 3E**). Set a timer for 60 s to allow worms to acclimate to swimming.

NOTE: The plate will not be 100% visible on the screen. Animals are not accurately tracked while overlapping.

2.12. While the timer counts down, manually adjust the camera focus by turning the focusing ring on the lens body of the camera while watching the display (**Figure 3F**). It is possible to zoom in on the worms using a scroll wheel mouse allowing for more precise focusing.

2.13. After setting the focus, adjust the light knob such that the display is as bright as possible without overexposure (shows up as red pixels on display).

2.14. After 60 s, press the **Record** button.

NOTE: The video will go from the **Live View** to **Recording** and a timer will be visible. When the capture is complete, the **Video Capture** display will return to the **Live View** (**Supplementary Figure 1B**).

2.15. Remove the plate after recording and set up the next plate. Rename the prefix as needed. Wash worms into the assay plate, incubate for 60 s, and record each treatment.

2.16. After completing all recordings for an experiment, close any open video tabs, and turn off the light source by depressing the dimmer knob.

2.17. Close the **Video Capture** window by clicking the red **x** button in the window's upper right corner; alternatively, close the entire software package and reopen.

2.18. Prepare videos for tracking.

2.19. From the **Workflow** sidebar menu, select the first button on this menu, **Import Image Sequence**, and find and double click on a video. The video will open in the middle window (Supplementary Figure 2A).

NOTE: The default display for the software is the **Workflow** sidebar menu; this menu can also be accessed through the view top dropdown menu under **View [View | Workflow]**. The **Workflow** menu has several view options; the **Track** view is appropriate at this point in the protocol. If the **Set Sequence Info** button is not visible in the menu, then this is not the correct view. First, look at the bottom of the **Menu** and check that the **Workflow** tab is active, next check that the **Track** view is active by clicking on the **Track** icon. Where menus appear in the window is dependent on monitor size and the location of other opened menus. Menus can be moved by drag and drop, resized, and closed. Most of the menus can also be opened through the **View** dropdown menu.

2.20. From the Workflow menu, select **Set Sequence Info** (Supplementary Figure 2B). In the new menu window, check the naming scheme, add notes, and validate metadata.

2.21. Check that the scaling is set correctly; **Scale**: 11.43  $\mu\text{m}/\text{pixel}$  and **Measure** will auto fill as 1000  $\mu\text{m}$  is 87 pixels. If the camera height has been adjusted differently than prescribed here, follow the instructions in this menu to update the scale. Click the **Save** button.

2.22. In the Workflow menu, select **Adjust Image**, and a new popup menu called **Image Adjustments** will open (Supplementary Figure 2C). Select **Dark worms on light background (bright field)**.

2.23. Adjust additional settings. Recommended **Image Processing Settings** are as follows: **Background smoothing**: 10, **Gaussian smoothing**: 5, **Fill holes**: 2, **Small object filter**: 0, and skip **Integral Derivative Segmentation**.

2.24. Adjust the **Threshold** level such that the animals are completely filled with green but still distinct from the background. Click **Apply**.

NOTE: It is permissible to pick up some objects that are not worms.

2.25. From the **Workflow** menu, select **Detect and Track** (Supplementary Figure 2D).

NOTE: A menu will appear to the side of the video with 3 tabs: **Detecting**, **Tracking**, and **Repair**. The cursor will look like a pointing finger when over the video display. Alternatively, click the **Select** icon (cursor arrow) in the icon menu in the upper left of the software.

2.25.1. Use the finger selector tool to select 3-7 worms and click the **Detect worms** button in the **Detection** tab.

NOTE: This selection helps the software accurately identify worms for given conditions and must be done for each video. After a few seconds, most of the worms will be highlighted in green and have a number in yellow. Worms that are touching or curled up will likely not be highlighted; however, the program will most likely find them during tracking. There is an option to adjust the detection parameters; however, the defaults work well under most conditions. If adjustments to the detection parameters are made, make sure to be consistent between samples and replicates.

2.26. Still in the **Detect and Track** window, move to the **Tracking** tab (Supplementary Figure 2E). In the **Tracking Parameters** check **Use backtracking**, uncheck **Track worms at the edge of the image**.

2.27. Set the **Max tracked hypothesis** to 5. Set the **Tracking mode** to **Swimming**.

NOTE: The Max tracked hypothesis setting adjusts how many potential worm-tracks a single tracked worm is allowed to have and gives the program some leeway in how accurately to track a single worm without losing track of it. 5 is a good setting for swimming.

2.28. Move to the **Advanced Settings** section and set the following: **Frames worms can touch boundary**: 50, **Frames worms can overlap**: 500, **Position tolerance**: 0.50, **Shape tolerance**: 0.50. Save these settings and deploy them for all videos in an experiment (and beyond) by saving them as a configuration.

NOTE: **Track filtering** is not recommended and can easily be done post-processing, leave off.

2.29. Navigate to the **Configuration Manager** (gear icon) in the upper-left icon menu. When clicked, the **Configuration Manager** sidebar menu will open. Next, click the **Save (gear +) icon**, give this configuration a name and description, then click the **OK** button.

2.29.1. To access this configuration in the future, open the **Configuration** menu, select the desired configuration, and click the **Load** (gear up arrow) icon. Then add sequence information, possibly fine-tune the threshold level, and select worms for detection. All other settings will be the same.

2.30. In the **Workflow** menu, click **Save Project**, verify the save location and add a project file name (matching the project file name to the video name helps with continuity). Close the video.

NOTE: These projects are saved with the .wpr filetype.

2.31. Repeat for all videos: Import Image Sequence (2.18), Load saved configuration (2.28.1), adjust the threshold (2.23), set sequence info (2.19), detect worms (2.24.1), save the project (2.29), close video.

## 2.32. Track videos

2.32.1. To track the project files, go to the **Workflow** menu and click the **Batch** (folder containing worm pages) icon (**Supplementary Figure 2F**).

NOTE: The **Batch Processing** panel will be displayed (**Supplementary Figure 2G**). Sometimes it is hidden but can be accessed at the very bottom of the **Batch Workflow** menu; the window may need to be expanded to see it.

2.32.2. Click the **Add** button under the **File Selection** section of this **Batch Processing** menu. Navigate to and select all the project files (.wpr) to be processed, click **Open**.

NOTE: The **Batch Processing** menu will fill with the selected file path, a status indicator and a gray progress bar for each file.

2.32.3. Click the **Start** (play button icon). Note the green progress indicator on the first file. Allow all the files to be processed. When all the files read as finished (green), the software can be closed (**Supplementary Figure 2G**).

2.32.4. If a project file gets stuck, consider running the hung-up file by itself, this usually fixes the problem. If this does not fix the problem, open the project file, delete the selected worms (**Detect** and **Track** menu), and reselect from a new training set. Then reprocess the project file.

## 2.33. Repairing tracks (recommended)

NOTE: Repairing tracks can increase the accuracy of the tracks. Common fixes are combining tracks that clearly belong to the same worm as viewed by the human eye as, deleting tracks that are not accurate or do not capture real worm data, and splitting tracks that are misrepresented because worms are touching (remove the bad section and combine tracks from before and after the poor section if practical).

2.33.1. Open the software and open a project file [**File | Open Project**] of interest. Click **Open** and click **yes** to use this file. Wait for the temporal filter to load. The processed video will be on

the screen; worms will be highlighted in green and outlined in yellow. Use the video player to scrub through or play the video.

2.33.2. From the **Workflow** menu, select **Detect and Track (Supplementary Figure 2D)**. Switch to the **Repair** tab. Activate the **Selector** cursor by clicking on the **Selector tool** icon (cursor arrow). Scrub through the video, finding and selecting a worm track, and then use the repair menu to create the desired change.

NOTE: To strategically fix most problems quickly, zoom in such that about 1/6<sup>th</sup> of the plate is visible, start with the area surrounding track #1. Quickly scrub through the video checking for worms that mysteriously appear out of nowhere, tracks where the numbers change, or worms bumping into each other, editing at problems arise. Follow worms in numerical order from the start position. Scrub back and forth to follow each worm as necessary until the whole plate is repaired. Pay careful attention to the edges of the plate; often, fake worms appear in the shadows of plate edges.

## 2.34. Data analysis

2.34.1. From the **Workflow** menu, select **Analyze Data (Supplementary Figure 3A)**.

NOTE: A new window named **Data Analysis and Plotting** will appear (**Supplementary Figure 3B**). This window has many ways to assess and view the movement data of a specific project file's worm population. Most of these analyses allow for the selection of specific, multiple, or all worms from a given video allowing for track information to be plotted in different ways.

2.34.2. To analyze turn count, navigate to the **Track Summary [Position & Speed | Track Summary]** (**Supplementary Figure 3C**). Use the **Export** button in the bottom right to export data in a spreadsheet readable format (.csv), choose a path to the desired location of the spreadsheet, including renaming it to something memorable. The default settings are a good choice (**Supplementary Figure 3D**).

2.34.3. Paste this data, as well as any other relevant project datasets into a consolidated spreadsheet workbook (**Supplementary Figure 3E**).

2.34.4. Use the **Turn Count** and **Track Duration** to calculate turns per min for each track in using spreadsheeting functions ( $\text{=turns/track duration} * 60$ ) or use a different metric as desired (**Supplementary Figure 3F**).

NOTE: Turns per minute captures very similar information as a manual thrashing assay which is a widely used metric of activity in liquid, but there are many options, and it is quite possible that a different metric may identify other interesting data for a given treatment or comparison.



2.34.5. If desired, screen for worm size, track duration, and/or some other chosen metric(s) using conditional formatting tools to highlight tracks for removal if they are potentially confounding. Maintain consistent data processing between treatments and replicates.

2.34.6. Analyze data using appropriate statistical analyses such as Student's t-test for 2 strains, or 1-way ANOVA for 3 or more strains. Perform independent experimental replicates on different days and with independent populations of worms for at least 3 independent experimental replicates.

#### **REPRESENTATIVE RESULTS:**

Both radial locomotion and swimming assays offer sensitive detection of motility impairment (**Figure 4** and **Figure 5**). To investigate the mechanisms underlying pathological TDP-43 in ALS, *C. elegans* models have been developed that express wild-type or ALS-mutant human TDP-43 pan-neuronally. These animals display molecular and cellular characteristics reminiscent of ALS, including motor dysfunction<sup>9</sup>. Importantly, they exhibit moderate motility impairment with the expression of wild-type human TDP-43 and more severe motility impairment in animals expressing ALS-mutant TDP-43 using both radial locomotion and swimming assays. Some mutant or transgenic animals will have greater impairment in crawling than swimming, or vice versa. By using two different motility assays, a clearer picture of phenotypic differences between strains is obtained.

#### **Radial locomotion**

When placed on a seeded agar plate, *C. elegans* explore their surroundings, including seeking out the boundaries of their food source. Radial locomotion assays are one way to take advantage of that behavior as a metric for physical fitness. By analyzing the locomotor behavior (crawling on a solid surface) in a controlled and quantifiable way, the radial locomotor assay offers a simple and effective tool for assessing the severity of motor deficits and other motor-related phenotypes. Radial locomotion assays capture differences in motility in moderately or severely impaired ALS model worms and offer a baseline to compare modulation of motility phenotypes or changes in motility with age (**Figure 6**). This strategy can be applied to quantify the crawling behavior of any strain that has altered movement from wild type (N2) or control worms. However, this method may not be a good choice to assess animals unable to crawl normally, such as roller mutants or animals that are paralyzed. Typically, wild-type worms will present an average displacement between 200-300  $\mu\text{m}/\text{min}$  when raised and tested at 20 °C. The example data presented in **Figure 4** show expected results comparing N2, two different transgenic strains expressing wild-type human TDP-43 with a mild phenotype [TDP-43(WT-mild), CK402 (bklIs402[*Psnb-1::TDP-43*; *Pmyo-3::GFP*])] or stronger phenotype [TDP-43(WT-moderate), CK410 (bklIs402[*Psnb-1::TDP-43*; *Pmyo-2::GFP*])], a transgenic strain expressing mutant human TDP-43 with a severe phenotype [TDP-43(M337V), CK423 (bklIs423[*Psnb-1::TDP-43*; *Pmyo-2::GFP*])], and a transgenic strain expressing another neurodegenerative-disease associated protein, wild-type human tau [tau(WT), CK144 (bklIs144[*Paex-3::tau(4R1N)*; *Pmyo-2::GFP*])]. The two strains expressing wild-type TDP-43 have different degrees of impairment as detected by radial locomotion. TDP-43(WT-low) is not significantly different from N2, while TDP-43(WT-high) does exhibit clear differences in motility. TDP-43(M337V) and tau (WT) strains have more severe impairments in motility.



## Swimming assay

*C. elegans* engage in stereotypical swimming (thrashing) movement when immersed in liquid. Immediately upon submersion, worms begin to characteristically bend head and tail towards each other with a bend angle of approximately 45°, with the angle vertex being the midpoint of the worm. Worms alternate bending in the ventral and dorsal directions. One thrash, as measured by the software, represents the movement of going from straight to a body bend angle of 20° or greater, irrespective of directionality (angle thresholding can be adjusted post-processing within the **Analysis** and **Plotting** window [**Workflow | Analyze Data | Body Shape | Bending Angle | Mid-Point | Amplitude Threshold: # degrees**]). The swimming assay described here uses automated computer-based tracking and analysis to provide unbiased scoring of swimming activity. It is expected that wild-type (N2) worms will average between 150-200 thrashes per min when raised and recorded at 20 °C. The example data presented in **Figure 5**, shows expected results comparing N2, two different transgenic strains expressing wild-type human TDP-43 with a mild phenotype [TDP-43(WT-mild), CK402 (bkl402[*Psnb-1::TDP-43*; *Pmyo-3::GFP*])] or stronger phenotype [TDP-43(WT-moderate), CK410 (bkl402[*Psnb-1::TDP-43*; *Pmyo-2::GFP*])], and a transgenic strain expressing another neurodegenerative-disease associated protein, wild-type human tau [tau(WT), CK144 (bkl144[*Paex-3::tau(4R1N)*; *Pmyo-2::GFP*])]. The transgenic strain expressing ALS-mutant TDP-43 [TDP-43(M337V)] does not thrash in liquid, and so it is denoted on the graph as ND (no data). This assay can discriminate different phenotypes than the radial locomotion assay shown in **Figure 4**. For example, in the radial locomotion assay (**Figure 4**), TDP-43(WT-low) was not significantly different from N2. However, in the swimming assay (**Figure 5**), both TDP-43(WT-low) and TDP-43(WT-high) are significantly different from N2, as well as significantly different from one another. Also, despite both TDP-43(M337V) and tau(WT) having severe crawling impairments by radial locomotion (**Figure 4**), only tau(WT) is able to thrash enough to be tracked by the software (**Figure 5**). TDP-43(M337V) animals are unable to thrash, and software-based analysis does not detect or track these worms accurately. Thus, data on these worms were not collected (ND, no data).

## FIGURE AND TABLE LEGENDS:

**Figure 1: Radial locomotion assay workflow.** Panels A-E show the generalized steps of the radial locomotion assay. The steps are as follows: (A) NGM Plates, with OP50 seeded to the edges, are prepared by labeling and marking central dot, (B) worms are placed in the center of the agar as marked by the central dot, (C) worms are allowed to move freely for a set amount of time, (D) plate is flipped bottom up and the final location of each worm is marked in a different color than the central dot, (E) The distance from the center dot to each final worm location is measured either by hand or digitally.

**Figure 2: Digital measurement of final location using ImageJ.** (A) Distance from center measurements can be measured by hand or digitally, to measure digitally using ImageJ scan backside of the plate with a ruler in the frame. (B) Draw a known length using the line tool. (C) Use the known length drawn in (B) to set scale [**Analyze | Set Scale...**]. (D) Use the line tool to draw a line from the center point to a final location mark, repeat for each mark. (E) A

paintbrush line marking the first mark measured aids in scoring (squiggly black line near the 1 mark). Final location marks are numbered in yellow to illustrate the directionality of scoring. (F) Measurements are recorded in the Results window. Save results elsewhere for statistical analysis.

**Figure 3: Materials and hardware adjustments for using the software.** (A) Materials used for the software-assisted swimming assay. (B) General set up of the software and needed materials, note that the housing is in the raised position. (C) The lens mounting bracket height can be adjusted. This image shows the adjustable track and a tape indicator marking the preferred height for performing swimming assays. It is necessary to adjust the brightfield light and camera focus before recording. (D) A 35 mm assay plate is placed on the stage, animals in M9 are added to the plate, and a 1 min timer is started. (E) It is sometimes advantageous to gently swirl the plate to move worms closer to the center – observe location on the live capture screen; alternatively, a few drops of M9 from a pipettor can be used to separate worms. (F) Adjust the focus ring on the camera lens, observe worms on the screen to determine optimal focus.

**Figure 4: Radial locomotion assays detect differences in crawling speed.** Unstimulated dispersion of developmentally staged L4 larvae was measured using the radial locomotion assay described above and graphed as  $\mu\text{m}/\text{min}$  traveled (A-B). The same data plotted in (A) and (B) demonstrate two different possible graphical presentations. In (A), data are displayed as a bar graph, making relative differences between strains clearer. In (B), the final displacement of each worm is plotted within the graph, allowing the variation within the population to be better visualized. To evaluate significance among strains tested, one-way analysis of variance (ANOVA) with Tukey's multiple comparison test was used.  $**p=0.0022$ ,  $****p<0.0001$ , ns=not significant. TDP-43(M337V) and tau(WT) are also significantly different from N2,  $p<0.0001$ . Error bars in (A) are standard error of the mean (SEM) and in (B) are standard deviation.

**Figure 5: Swimming assays detect differences of thrashing in liquid.** Rates of swimming (thrashing or undulation frequency in liquid) were measured using unbiased computer-assisted scoring and analysis described above and graphed as thrashes/min (A-B). Same data plotted in (A) and (B). In (A), data are graphed as a bar graph, which makes relative differences between strains easier to see. In (B), data from each individual worm scored are plotted within the graph, allowing the variation within the population to be better visualized. To evaluate significance among strains tested, a one-way analysis of variance (ANOVA) with Tukey's multiple comparison test was used.  $****p<0.0001$ , ns=not significant. TDP-43(WT-moderate) and tau(WT) are also significantly different from N2,  $p<0.0001$ . Error bars in (A) are standard error of the mean (SEM) and in (B) are standard deviation.

**Figure 6: Mutant TDP-43 worms travel less than wild-type worms.** Plates showing a representative difference in the final location of staged L4 N2 (wild-type) and TDP-43(M337V) (CK423 (bkl423[*Psnb-1::TDP-43*; *Pmyo-2::GFP*]), a strain expressing mutant human TDP-43 with severely impaired motor function after 1 hour of unstimulated crawling at room temperature.

Plates are marked with a red dot for the center point and blue dots for the final location of the animals.

**Supplementary Figure 1: Recording swimming behavior using the imaging and tracking software.** (A) The tracking workflow and the location of the video capture icon. (B) The **Video Capture** window is shown in **Live View**, and the most important aspects are highlighted. The green “Live view” words will change to a red “Recording” when the record button is activated.

**Supplementary Figure 2: Preparing and tracking swimming behavior.** This figure shows a series of screenshots to aid in setting up a video image sequence for tracking. The general protocol for preparing a video sequence is to open each workflow menu in this order: **Import Image Sequence | Set Sequence Info | Adjust Image | Detect and Track | Save Project**. The **Analyze Data** window can only be used after a project file is tracked. (A) shows an opened .avi (video) file after importing the sequence into the software from within the **Track** workflow. (B) The **Set Sequence Info** window provides a place for notes, set/change the scale, and to check the metadata for a video sequence. Instructions for changing the scale are found in this window and should be followed when the camera is lowered or raised (B). Image (C) shows the **Adjust Image** window settings. (D) Screenshot showing the detection tab of the **Detect and Track** window. Detecting worms is used to train the software and must be set for each video sequence. Additional detection parameters can be set here if desired. (E) Screenshot of **Tracking** tab of the **Detect and Track** window with the recommended settings for tracking swimming behavior. This is the last step in setting up a video sequence. Save the sequence as a project using the **Save Project** window. Repeat these steps for each video sequence. Settings can be saved as a configuration to reduce workload and ensure that all sequences are treated the same (not shown). To track project files for analysis, navigate to the **Batch** Workflow. (F) shows the icon used to navigate to the **Batch** workflow, and (G) shows the batch processing window, highlights the add and start buttons, and shows the expected appearance of a project file that has been tracked.

**Supplementary Figure 3: Swimming behavior analysis.** After a project file has been tracked, it can be opened using [**File | Open Project**]. Worms will appear in green as shown in (A). Scrubbing through the video offers a quick check that the video processed as expect, green highlighting will disappear while scrubbing. The Data Analysis and Plotting window is opened from the **Analyze Data** workflow item. (B) shows the Data Analysis and Plotting window with Position view opened (default), all tracks are highlighted. Below the datapoints is a plot showing the recorded track for each highlighted worm. The track summary analysis is used to calculated turns per minute. (C-D) show the track summary report and export details. (E-F) show the track summary in a spreadsheet and how to calculate turns per minute, the output measured in this assay.

## **DISCUSSION:**

### **Radial locomotion:**

The resolution of this assay is easily controlled by changing the time variable. Increasing the

length of time makes it easier to observe differences between animals with severe phenotypes, thereby identifying subtle differences. However, because this assay measures displacement, if the assay time is extended too long, animals with normal motility, such as N2, will travel to the edges of the plate, and foraging behavior will lead to backtracking. This will artificially decrease the measurement of distance traveled. Time periods that are too long may result in the disappearance of differences between strains, particularly between animals with less severe motor phenotypes, as animals become evenly dispersed across the plate. Shortening the time variable will prevent more active worms from finding the edges of the plate. This method does not track the total distance traveled for each worm but compresses the distance traveled for each worm into a linear distance from the center of the plate. As such, it is inherently less robust than a method that records the total track length of individual worms. However, the radial locomotion assay requires very little researcher training, uses relatively affordable reagents that are commonly available in most worm labs, and is sensitive enough to produce significant and reproducible results. For labs that prefer automated video tracking, several methods have been previously established to track and analyze crawling movements<sup>12</sup>, or the software parameters used for the swimming assay in this paper could be altered to permit crawling detection and analysis.

This experiment is usually done in triplicate independent replicates, with a set of 30-40 worms per replicate. Each replicate is split onto two different 100 mm or 150 mm plates, with 15-20 worms per plate. Using more worms than recommended per plate can make it hard to score efficiently. A total number scored of 90+ is sufficiently powered to establish significance for mild, moderate, or severe motility impairment. Being consistent with timing between strains scored is essential for accuracy and reproducibility. 30 minutes is generally long enough to establish differences between moderate to severe phenotypes such as transgenic strains expressing human mutant TDP-43 compared to wild-type worms (**Figure 4**). If the time variable will be extended, it is advisable to also increase the size of the plate from 100 mm to 150 mm. Environmental factors such as temperature and humidity can affect this assay, which is typically conducted at ambient room temperature, so it is important to always use a wild-type (N2) control when comparing across replicates. In addition, this assay can measure the motility of some strains that do not exhibit normal swimming behavior in liquid (thrashing), making it a useful complement to the swimming assay.

#### **Swimming assay:**

The use of the imaging and acquisition system to automate tracking and analysis of worm swimming allows for rigorous and unbiased data. However, there are several factors during initial setup of the experiment that needs to be controlled between samples. These include the time to acclimate to liquid before beginning recording, ambient conditions (i.e., temperature, humidity), and consistent light and recording settings. On the recording stage, there are several features that help reduce variability between plates. These include an integrated track-mounted camera and bright-field stage that makes recording video consistent between plates, shielding around the stage that prevents reflections, glare, and air movements while recording, and a robust software package that reliably detects worms and allows for manual correction of tracks in video post-processing. In this protocol, videos of a 35 mm plate with worms are recorded for 1 minute

and then processed using the software package. After processing, manual correction of tracks ensures that worm behaviors are recorded accurately without confounding tracking errors. Turn count and track duration data are used to determine turns per minute as a final readout. To ensure reproducibility, data is collected over a minimum of 3 independent replicate experiments, each with 40-50 animals scored, to achieve a combined final number of animals 120-150. This number is sufficient to discriminate small differences in swimming behavior from control worms. Some worms have motor deficits too severe to be captured by swimming assays. For example, if the animals placed in a liquid medium curl instead of performing the expected thrashing response, this assay will not record those movements accurately and another movement assay, such as radial locomotion, may better capture those motility defects. The provided protocol uses a commercially available imaging system (see **Table of Materials** for more details), but other worm tracking systems may provide a similar utility, with some being open source<sup>12</sup>. Previously published methods describe manual scoring of worm thrashing<sup>13</sup>. While the automated analysis produces a number of metrics for each individual worm, detection of body bends which is measured in thrashes per minute, provides consistent results between experiments and tracks well with conventional scoring of worm thrashes by eye.

#### **ACKNOWLEDGMENTS:**

We thank the reviewers for helpful comments and suggestions. We thank Aleen Saxton, Brandon Henderson, and Jade Stair for outstanding technical assistance. We thank Brian Kraemer and Rebecca Kow for assistance developing these assays. This material is the result of work supported with resources and the use of facilities at the VA Puget Sound Health Care System. This work was supported by a grant from the United States (U.S.) Department of Veterans Affairs (VA) Biomedical Laboratory Research and Development Service [Merit Review Grant #I01BX004044 to N.F.L.]

#### **DISCLOSURES:**

The authors declare that they have no disclosures.

#### **REFERENCES:**

- 1 Wijesekera, L. C., Leigh, P. N. Amyotrophic lateral sclerosis. *Orphanet Journal of Rare Diseases*. **4**, 3 (2009).
- 2 Arai, T. et al. TDP-43 is a component of ubiquitin-positive tau-negative inclusions in frontotemporal lobar degeneration and amyotrophic lateral sclerosis. *Biochemical and Biophysical Research Communication*. **351** (3), 602-611 (2006).
- 3 Neumann, M. et al. Ubiquitinated TDP-43 in frontotemporal lobar degeneration and amyotrophic lateral sclerosis. *Science*. **314** (5796), 130-133 (2006).
- 4 Mezzini, R. et al. ALS genetics, mechanisms, and therapeutics: Where are we now? *Frontiers in Neuroscience*. **13**, 1310 (2019).
- 5 Tan, R. H., Ke, Y. D., Ittner, L. M., Halliday, G. M. ALS/FTLD: Experimental models and reality. *Acta Neuropathologica*. **133** (2), 177-196 (2017).

699 6 Kim, W., Underwood, R. S., Greenwald, I., Shaye, D. D. OrthoList 2: A new comparative  
700 genomic analysis of human and *Caenorhabditis elegans* genes. *Genetics*. **210** (2), 445-461  
701 (2018).

702 7 White, J. G., Southgate, E., Thomson, J. N., Brenner, S. The structure of the nervous  
703 system of the nematode *Caenorhabditis elegans*. *Philosophical Transaction of the Royal Society  
704 of London B: Biological Sciences*. **314** (1165), 1-340 (1986).

705 8 Caldwell, K. A., Willicott, C. W., Caldwell, G. A. Modeling neurodegeneration in  
706 *Caenorhabditis elegans*. *Disease Models and Mechanisms*. **13** (10), dmm.046110, (2020).



707 9 Liachko, N. F., Guthrie, C. R., Kraemer, B. C. Phosphorylation promotes neurotoxicity in a  
708 *Caenorhabditis elegans* model of TDP-43 proteinopathy. *Journal of Neuroscience*. **30** (48),  
709 16208-16219 (2010).

710 10 Robatzek, M., Thomas, J. H. Calcium/calmodulin-dependent protein kinase II regulates  
711 *Caenorhabditis elegans* locomotion in concert with a G(o)/G(q) signaling network. *Genetics*. **156**  
712 (3), 1069-1082 (2000).

713 11 Porta-de-la-Riva, M., Fontrodona, L., Villanueva, A., Ceron, J. Basic *Caenorhabditis*  
714 *elegans* methods: synchronization and observation. *Journal of Visualized Experiments*. (64),  
715 e4019 (2012).


716 12 Husson, S. J., Costa, W. S., Schmitt, C., Gottschalk, A. Keeping track of worm trackers.  
717 *WormBook*. 1-17 (2013).

718 13 Ringstad, N., Horvitz, B., and Koelle, M. *WormBook*. Editor: A. C. Hart, 1551-8507 (2006).  
719


**A****Prepare Plate****B****Place Worms****C****Set Timer****D****Mark Final Location****E****Measure Distance**




A




B




C




D



E




F





A




Materials

B




WormLab Setup

C



Camera Height Indicator


D


Washed Worms Added  
to Assay Plate


E

Align Plate, Swirl Gently  
to Center Worms (optional)


F

Adjust Focus Ring on  
Camera Lens as Needed






WT



TDP-43(M337V)



Name of Material/ Equipment	Company	Catalog Number	Comments/Description
C. elegans	Caenorhabditis Genetics Center (CGC)	-	Acquire your strains as desired, N2 is a
Disposable pasteur pipets, borosilicate glass	VWR	14673-010	Glass pipet used to create worm pick
Disposable petri dishes, 35x10mm	VWR	10799-192	Assay plates for WormLab Imaging Sys
Disposable petri dishes, 60x15mm	VWR	25384-090	Stock plates for worms
Disposable petri dishes, 100x15mm	VWR	25384-302	Standard radial locomotion assay plat
Disposable petri dishes, 150x15mm	VWR	25384-326	Longer time frame radial locomotion
Dissecting microscope	Leica	M80	Scope for maintaining worms and sett
Fine-tipped markers	VWR	52877-810	Need at least 2 colors for radial locom
Flatbed Scanner	Amazon	V850	Optional for radial locomotion assay. I
ImageJ	NIH	-	Optional free software provided by th
M9 buffer	VWR	IC113037012	Medium used for swimming assay. Ca
NGM (Nematode Growth Medium)	VWR	76347-412	Medium used to cultivate C. elegans. I
OP50 bacteria	Caenorhabditis Genetics Center (CGC)	OP50	Primary food source for C. elegans
p1000 pipettor	VWR	76207-552	Pipettor, used in swimming assay
p1000 tips	VWR	83007-384	Tips for pipettor, used in swimming as
Platinum wire, 0.2032mm diameter	VWR	BT136585-5M	Fine gauge platinum wire used to crea
Ruler	VWR	56510-001	Need to score radial locomotion assay
WormLab Imaging System	MBF Bioscience	WormLab	The Imaging System includes WormLa

useful control strain

- hold glass pipette in one hand and ~1" of platinum wire (held by pliers) in the other over a flame to join.

stem

e

assay plate

setting up radial locomotion assays

locomotion assays. Fine tips required for accuracy.

Protocol assumes a resolution of 300dpi, most scanners would work fine

at NIH - <https://imagej.nih.gov/ij/>

can be made from scratch, see WormBook: Maintenance of *C. elegans*

Can be made from scratch, see WormBook: Maintenance of *C. elegans*

assay

make worm pick - hold glass pipette in one hand and ~1" of platinum wire (held by pliers) in the other over a flame to join.

15

lab hardware (bright field stage, camera, and housing) and WormLab software. <https://www.mbfbioscience.com/wormlab-imaging-system>



Dear Editorial Staff,

Thank you for your edits of our article now titled “Evaluate motor impairment in *C. elegans* models of amyotrophic lateral sclerosis (ALS)”. We appreciate the helpful comments and suggestions. As a result we have updated accompanying text according to the editor recommendations and moved several figures to the supplement.

best regards,  
Nicki

---


Nicki Liachko, Ph.D.

GRECC Core Investigator  
Geriatrics Research Education and Clinical Center  
VA Puget Sound Health Care System, S182  
1660 South Columbian Way  
Seattle, WA 98108  
[206-277-1587](tel:206-277-1587)

Research Assistant Professor  
Division of Gerontology and Geriatric Medicine  
University of Washington  
Seattle, WA 98195  
[206-931-3392](tel:206-931-3392)




A



B

Video Capture



**Camera**

**Live view**

Basler - acA2440-35um  
Time 00:09:05.629  
Frame 7639  
Rate 14.000  
Buffer

Display grid Display clipping

**Capture** **Settings**

Folder Tracker/TDP-43 Worm Test

File prefix 20210503\_CK1441;CK1943

Buffer 128 frames

☒ Duration 0 : 1 : 0

☐ Time-lapse 1.00 seconds

☒ Open video after recording

☐ Show Stimulus Controls

Record Stop

Snapshot

Live View – not recording

File location


File name

Duration (set to 1 min)


Recording Controls




A



B



C



Set threshold such that worms are green

Image Processing Settings

D

The screenshot displays the WormLab software interface. The main window shows a video frame with several worms. Some worms are highlighted in green, indicating they have been detected. The interface includes a menu bar (File, View, Tools, Help), a toolbar, and a workflow panel on the left with buttons for Track, Live, Batch, Import Image Sequence, Set Sequence Info, Adjust Image, Detect and Track, Analyze Data, and Save Project. The 'Detection and Tracking' panel on the right has tabs for Detection, Tracking, and Repair. The 'Detection' tab is active, showing a 'Detect worms' button (highlighted with a red box and a red arrow), a 'Delete worms' section with radio buttons for 'This frame' (selected) and 'All frames', and a 'Detection Parameters' section with various sliders and checkboxes. A text box in the bottom right corner provides instructions for the detection process.

1. Use cursor to select several worms  
2. Click “Detect Worms”  
3. Wait for most worms to be highlighted in green

E

WormLab : 20200316\_Rep\_N2\_0001.avi (unsaved changes)

File View Tools Help

Workflow

Track Live Batch

Import Image Sequence

Set Sequence Info

Adjust Image

Detect and Track

Analyze Data

Save Project

20200316\_Re...001.avi X

Detection and Tracking

Detection Tracking Repair

From frame 1 to 840

☒ Restrict tracking to labelled region

Start Pause Stop

Tracking Parameters

☒ Use back tracking

☐ Track worms at the edge of the image

Max tracked hypotheses 5

☐ Exhaustive hypotheses search

Tracking Mode

☐ Crawling

☒ Swimming

Advanced

Frames worms can touch boundary 50

Frames worms can overlap 500

Position tolerance 0.50

Shape tolerance 0.50

Track Filtering

☒ Minimum track duration (frames) 420

Copy parameters to Clipboard

254 256 258 260 262 264 266 268 270 272

263 / 840

00:00:18.714


14.00 fps

Ready

Save as project when done!



Tracking Settings

F






G




Add Project Files,  
Click “Start” to Run

A



B



C

Data Analysis and Plotting

Select an Analysis

- Position & Speed
  - Position (um)
  - Center Points (um)
  - Direction
  - Speed (um/s)
  - Moving Average Speed (um/s)
  - Smoothed Speed (um/s)
  - Reversal
  - Mobility
  - Validation
  - Position Speed Statistics
  - Pirouette
- Track Summary
- Body Shape
- Histograms
- Label Analysis
- Chemotaxis
- Swimming

Filtering option

☒ Kalman Smoothing

Export selected

Track Summary

Refresh Data

Data Report

Track Number

	1	2	3	4	5	6	7	8	9
<b>Mean Worm Length (um)</b>	674.26	686.52	732.80	720.40	718.74	788.12	862.28	722.57	743.37
<b>Mean Width (um)</b>	58.21	61.04	61.53	59.10	61.42	60.78	70.28	73.29	60.90
<b>Mean Area (um<sup>2</sup>)</b>	39924.05	42604.36	45874.72	43307.35	44904.44	48731.65	61642.91	54192.86	46048.25
<b>Track Length (um) [Center point trajectory]</b>	23794.35	11736.97	26822.62	28227.75	32378.65	23988.81	15957.33	12211.07	7806.49
<b>Center point Speed (um/s) [Center point trajectory/time]</b>	397.05	195.85	447.58	471.02	540.29	465.16	521.97	348.18	481.46
<b>Absolute Peristaltic Track Length (um) [forward + reversals]</b>	14875.82	5508.27	17183.21	15938.67	21243.16	17197.74	10625.01	7215.72	5305.72
<b>Peristaltic Speed (um/s) [Absolute peristaltic track length/time]</b>	248.23	91.91	286.73	265.96	354.47	333.47	347.55	205.74	327.23
<b>Peristaltic Track Length (um) [forward - reversals]</b>	6257.98	-3624.59	-13734.06	-15767.12	-21048.29	-13297.44	-7427.13	3110.29	-5267.94
<b>Peristaltic Speed (um/s) [peristaltic track length/time]</b>	104.42	-60.48	-229.17	-263.10	-351.22	-257.85	-242.94	88.68	-324.90
<b>Straight-line Distance (um)</b>	2687.76	80.24	1822.21	2160.97	1346.03	6120.54	1371.07	1455.47	1116.95
<b>Wavelength (um)</b>	388.19	404.05	429.81	422.49	449.54	494.14	532.99	403.69	465.05
<b>Mean Amplitude (um)</b>	65.38	48.61	66.61	71.57	70.03	61.67	83.92	63.14	71.47
<b>Max Amplitude (um)</b>	131.10	117.79	153.45	132.01	131.96	140.60	182.18	162.88	129.18
<b>Turn Count</b>	183.00	69.00	207.00	196.00	183.00	154.00	90.00	78.00	60.00
<b>Track duration (s)</b>	59.93	59.93	59.93	59.93	59.93	51.57	30.57	35.07	16.21
<b>Distance Traveled Forward(um)</b>	10566.90	941.84	1724.58	85.77	97.44	1950.15	1598.94	5163.01	18.89
<b>Distance Traveled in Reverse(um)</b>	4308.92	4566.43	15458.64	15852.89	21145.73	15247.59	9026.07	2052.72	5286.83
<b>Number of Reversals</b>	7.00	17.00	2.00	1.00	1.00	3.00	1.00	5.00	1.00
<b>Cumulative Reversal Time (s)</b>	24.07	45.00	46.50	60.00	60.00	48.00	20.50	11.57	16.21
<b>Cumulative Forward Time (s)</b>	35.93	15.00	13.50	0.00	0.00	3.64	10.14	23.57	0.07
<b>Number of Pirouette</b>	1.00	0.00	0.00	0.00	0.00	0.00	0.00	1.00	0.00
<b>Pirouette Time (s)</b>	0.86	0.00	0.00	0.00	0.00	0.00	0.00	1.57	0.00
<b>Mobility reversal time (s)</b>	25.50	47.29	46.21	60.00	60.00	47.93	20.93	13.07	16.29
<b>Mobility idle time (s)</b>	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00
<b>Mobility forward time (s)</b>	34.50	12.71	13.79	0.00	0.00	3.71	9.71	22.07	0.00

Copy Selected Data Export...

Export Data (all tracks)

D

WormLab: 20200316\_Rep\_N2\_0001.avi

File View Tools Help

Workflow

Track Live Batch

Import Image Sequence

Set Sequence Info

Adjust Image

Detect and Track

Analyze Data

Save Project

20200316\_Re...001.avi X

Data Analysis and Plotting

Select an Analysis

- Position & Speed
  - Position (um)
  - Center Points (um)
  - Direction
  - Speed (um/s)
  - Moving Average Speed (um/s)
  - Smoothed Speed (um/s)
  - Reversal
  - Mobility
  - Validation
  - Position Speed Statistics
  - Pirouette
  - Track Summary
- Body Shape
- Histograms
- Label Analysis
- Chemotaxis
- Swimming

Filtering option

☒ Kalman Smoothing

Export selected

Track Summary

Refresh Data

Data Report

Track Number

	84	85	86	87	88	89
Mean Worm Length (um)	77.16	751.47	741.08	557.65	691.05	657.40
Mean Width (um)	0.30	59.10	61.45	59.33	62.66	58.79
Mean Area (um^2)	7670.16	45166.27	46365.22	33656.74	44049.09	39310.60
Track Length (um) [Center point trajectory]	906.12	3125.10	5159.47	485.75	2739.22	1931.68
Center point Speed (um/s) [Center point trajectory/time]	67.78	514.72	472.11	2266.85	314.34	443.34
Absolute Period (s)	3.22	3475.31	484.35	1561.29	1240.59	
Peristaltic Speed (um/s)	7.71	318.00	2260.28	179.16	284.73	
Peristaltic Trajectory (um)	99.04	-3475.31	-484.35	-1365.62	-1240.59	
Peristaltic Speed (um/s)	6.31	-318.00	-2260.28	-156.71	-284.73	
Straight-line Distance (um)	37.06	1587.77	476.93	421.70	900.12	
Wavelength (um)	3.50	462.76	225.32	411.94	394.35	
Mean Amplitude (um)	63	74.46	31.83	65.39	49.98	
Max Amplitude (um)	9.26	140.97	132.36	57.99	117.87	95.70
Turn Count	4.00	19.00	39.00	1.00	16.00	15.00
Track duration (s)	660.93					
Distance Traveled Forward (um)	.00					
Distance Traveled in Reverse (um)	.00					
Number of Reversals	.00	1.00	1.00	0.00	1.00	1.00
Cumulative Reversal Time (s)	.00	4.57	11.00	0.00	7.71	4.43
Cumulative Forward Time (s)	.43	1.57	0.00	0.29	1.07	0.00
Number of Pirouette	.00	0.00	0.00	0.00	0.00	0.00
Pirouette Time (s)	.00	0.00	0.00	0.00	0.00	0.00

Export Data

From frame 1 through 840

File 200316\_TAUXDP43xSUT2\_rep2/TrackSummary.csv Choose...

☒ Include column names

☒ Include sequence information at top of file

☐ Use tab separator instead of comma

☒ Open file in external program when finished

Export Cancel

Copy Selected Data Export...

Ready Vertical

Project loaded 94.246 Mb

Check location and dataset name, export

A1		Worm Data																					
A		B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q	R	S	T	U	V	W
1	Worm Data	Track Summary																					
2	Sequence	20200316_Rep_N2_0001.avi																					
3																							
4		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
5	Mean Worm Length (um)	674.26	686.52	732.8	720.4	718.74	788.12	862.28	722.57	743.37	703.52	716.65	664.35	627.96	679.29	781.25	641.31	686.53	753.79	765.03	644.05	731.55	789.04
6	Mean Width (um)	58.21	61.04	61.53	59.1	61.42	60.78	70.28	73.29	60.9	64.31	59.97	58.42	55.76	65.7	64	61.01	59.89	58.56	62.05	62.74	60.46	65
7	Mean Area (um^2)	39924.05	42604.36	45874.72	43307.35	44904.44	48731.65	61642.91	54192.86	46048.25	45902.93	43715.97	39474.1	35630.87	45201.95	50858.77	39645.85	41818.33	44901.64	48286.08	41056.7	44965.99	52171.41
8	Track Length (um) [Center point trajectory]	23794.35	11736.97	26822.62	28227.75	32378.65	23988.81	15957.33	12211.07	7806.49	26370.85	29161.31	9494.03	11274.16	27159.58	29200.59	12764.09	5027.64	1158.44	22862.38	4290.85	20636.43	20624.62
9	Center point Speed (um/s) [Center point trajectory/time]	397.05	195.85	447.58	471.02	540.29	465.16	521.97	348.18	481.46	440.04	486.6	406.47	482.69	453.2	487.26	412.7	502.76	523.17	540.66	349.26	488.02	449.76
10	Absolute Peristaltic Track Length (um) [forward + reversals]	14875.82	5508.27	17183.21	15938.67	21243.16	17197.74	10625.01	7215.72	5305.72	17114.4	18793.41	6043.12	7578.37	17615.92	17503.27	8893.55	3118.43	770.73	15592.45	2338.46	13501.54	13008.22
11	Peristaltic Speed (um/s) [Absolute peristaltic track length/time]	248.23	91.91	286.73	265.96	354.47	333.47	347.55	205.74	327.23	285.58	313.6	258.73	324.46	293.95	292.07	287.55	311.84	348.07	368.74	190.34	319.29	283.67
12	Peristaltic Track Length (um) [forward - reversals]	6257.98	-3624.59	-13734.1	-15767.1	-21048.3	-13297.4	-7427.13	3110.29	-5267.94	7920.92	-11933.1	5882.6	-4039.66	8360.17	-8749.2	8089.91	1117.3	-766.9	-12441.3	1816.84	-7043.99	-6448.61
13	Peristaltic Speed (um/s) [peristaltic track length/time]	104.42	-60.48	-229.17	-263.1	-351.22	-257.85	-242.94	88.68	-324.9	132.17	-199.12	251.85	-172.95	139.5	-145.99	261.57	111.73	-346.34	-294.22	147.88	-166.58	-140.62
14	Straight-line Distance (um)	2687.76	80.24	1822.21	2160.97	1346.03	6120.54	1371.07	1455.47	1116.95	419.03	4530.41	1208.91	217.43	2064.61	3224.21	647.04	487.33	549.42	836.19	488.29	1712.89	1175.84
15	Wavelength (um)	388.19	404.05	429.81	422.49	449.54	494.14	532.99	403.69	465.05	399.77	434.94	411.48	361.89	386.15	479.97	376.61	378.01	514.33	476.36	316.62	460.29	492.08
16	Mean Amplitude (um)	65.38	48.61	66.61	71.57	70.03	61.67	83.92	63.14	71.47	73.87	68.13	61.3	62.68	60.38	68.34	43.43	65.78	67.85	74.86	60.9	68.82	72.55
17	Max Amplitude (um)	131.1	117.79	153.45	132.01	131.96	140.6	182.18	162.88	129.18	180.05	150.18	118.8	124.33	144.95	147.3							



Turn Count / Duration \* 60

[TrackSummary](#) [+](#)

EDIT