Journal of Visualized Experiments

Development of a lateral-flow immunochromatographic strip for Rapid and Quantitative Detection of small molecule compounds --Manuscript Draft--

| Article Type: | Invited Methods Article - Author Produced Video | | |
|---|--|--|--|
| Manuscript Number: | JoVE62754R3 | | |
| Full Title: | Development of a lateral-flow immunochromatographic strip for Rapid and Quantitative Detection of small molecule compounds | | |
| Corresponding Author: | Huihua Qu Beijing University of Chinese Medicine Beijing, Beijing CHINA | | |
| Corresponding Author's Institution: | Beijing University of Chinese Medicine | | |
| Corresponding Author E-Mail: | quhuihuadr@163.com | | |
| Order of Authors: | Yue Zhang | | |
| | Peng Cao | | |
| | Fang Lu | | |
| | Jinjun Cheng | | |
| | Huihua Qu | | |
| Additional Information: | | | |
| Question | Response | | |
| Please specify the section of the submitted manuscript. | Biochemistry | | |
| Please indicate whether this article will be Standard Access or Open Access. | Standard Access (\$1400) | | |
| Please confirm that you have read and agree to the terms and conditions of the author license agreement that applies below: | I agree to the Author License Agreement | | |
| Please provide any comments to the journal here. | none | | |
| Please confirm that you have read and agree to the terms and conditions of the video release that applies below: | I agree to the Video Release | | |

TITLE:

- 2 Development of a Lateral Flow Immunochromatographic Strip for Rapid and Quantitative
- 3 Detection of Small Molecule Compounds

4 5

1

AUTHORS AND AFFILIATIONS:

6 Yue Zhang^{1#}, Peng Cao^{2#}, Fang Lu¹, Jinjun Cheng³, Huihua Qu⁴

7

- 8 1. School of Life Science, Beijing University of Chinese Medicine, Beijing, China
- 9 2. Third Affiliated Hospital, Beijing University of Chinese Medicine, Beijing, China
- 10 3. National Institute of TCM Constitution and Preventive Medicine, Beijing University of
- 11 Chinese Medicine, Beijing, China
- 4. Center of Scientific Experiment, Beijing University of Chinese Medicine, Beijing, China

13 14

KEYWORDS:

lateral flow immunoassay; colloidal gold; monoclonal antibody; artificial antigen; small molecular compound; strip

17 18

19

20

21 22

23 24

25

26

ABSTRACT:

Membrane-based lateral flow immunochromatographic strips (ICSs) are useful tools for low-cost self-diagnosis and have been efficiently applied to toxin, physiological index and clinical biomarker detection. In this protocol, we provide a detailed description of the steps to develop a rapid, sensitive and quantitative lateral-flow immunoassay (using AuNPs as a marker and mAbs as a probe). The procedure describes the preparation and characterization of colloidal gold, synthesis of the AuNP-mAb conjugate, assembly of the immunochromatographic strip, and methodological investigation of the assay. The results showed that the final strips can be further utilized for the rapid and convenient self-diagnosis of a small molecule, which may provide an alternative tool in the rapid and precise analysis of physiological and biological indices.

272829

30

31

32

33

34

35 36

INTRODUCTION:

Membrane-based lateral flow immunochromatographic strips (ICSs) are useful tools for low-cost and rapid detection. The nitrocellulose membrane as the carrier and colloidal gold as markers of immune chromatography rapid diagnostic reagents are the most commonly used POCT (point of care testing) method, and the testing scope of the project is wider. From their original application in monitoring during pregnancy, their use has been extended to monitor blood coagulation state¹⁻², myocardial injury³, veterinary medicine⁴, pesticide residues⁵, infectious diseases⁶ and drug concentrations. More types of samples can be assessed, including urine, saliva, whole blood, serum and other body fluids⁷⁻⁸.

373839

40

41

In recent years, numerous novel assays have been developed for detecting biomarkers in the diagnosis of disorders, including HPLC, UPLC, LC–MS and ELISA, which are sensitive and accurate, credible and specific. However, these methods require sophisticated instrumentation, complex

preprocessing and time-consuming treatments ⁹. Hence, developing a more rapid and convenient point-of-care diagnostic strategy for the self- and real-time detection of medicinal active compounds is urgent ¹⁰⁻¹¹.

The popularity of ICSs, especially for common tests, is driven by their ease of use, as they do not require professionals or elaborate instrumental setups¹². In other words, people who do not have special training can operate strips or self-tests¹³. The results of the test can be obtained in 5 minutes, which means it can be used for site inspections¹⁴. Moreover, according to our calculations, the cost of strips could be lower than 1 RMB¹⁵, which means that the tests are inexpensive to promote¹⁶. Hence, the ICS is a relatively accurate, simple, and inexpensive disposable device. ICSs based on colloidal gold ¹⁷⁻¹⁸ are also applied in rapid COVID-19 detection.

The principle of ICS can be divided into sandwich ICS and competitive ICS. **Figure 1A** is a schematic diagram of the sandwich ICS, which is mainly used for detecting macromolecular substances such as proteins, including tumor markers, inflammatory factors, and human chorionic gonadotropin (HCG, early pregnancy antigen). In this method, paired antibodies targeted at different epitopes of the antigen are used, and the capture antibody is dried on the NC membrane as a test line. Labeled antibody is dried on the conjugate pad, and secondary antibody is used as the control line.

Figure 1B is a schematic diagram of competitive ICS, which is mainly used to detect small molecule substances (MWCO < 2000 Da). The coating antigen is fixed on the NC membrane as a test line, and the labeled antibody is dried on the conjugate pad. During detection, the sample and labeled antibody flow through the detection line under capillary action, and the coated antigen competitively binds free antigen in the sample and develops a red color on the detection line.

Recently, we described the procedure of monoclonal antibody generation against natural products¹⁹. In this work, we developed a novel lateral flow immunoassay based on the prepared anti-SSD mA ²⁰ for rapid, on-site detection. The results indicate that the immunochromatography assay is an indispensable and convenient tool for detecting natural product-derived compounds.

[Figure 1 here]

PROTOCOL:

All of the procedures performed in this study were approved by the Ethics Review Committee at the Beijing University of Chinese Medicine (approval number 2017BZYYL00120).

1. Preparation and Characterization of Colloidal Gold

83 84 85

86

87

88

NOTE: For colloidal gold synthesis, as colloidal gold is easily adsorbed on the inner wall of the vessel and is prone to precipitation by impurities, the vessel for synthesis and storage of colloidal gold should be thoroughly cleaned and soaked in acid (40 mL of distilled water, 360 mL of concentrated sulfuric acid, 20 g of potassium dichromate) or subjected to surface passivation treatment. A citric acid reduction method was used to synthesize colloidal gold.

89 90

91 1.1 Turn on the magnetic stirrer and place the flask (250 mL) on the mixer.

92

93 1.2 Prepare 4% gold chloride acid solution and 1% sodium citrate solution, respectively.

94

1.3 Add 120 mL of distilled water to a round bottom flask and heat it to boiling on a thermostatic magnetic stirrer.

97

98 1.4 Keep boiling, and quickly add 0.5 mL of 4% chloroauric acid and 5 mL of 1% sodium citrate.

99

1.5 Observe the color of the solution. The pale yellow chloroauric acid solution turns wine red within a few minutes.

102

1.6 Continue heating for 10 minutes, until the solution changes from colorless to transparent wine red.

105

1.7 Turn off the power of the thermostatic magnetic mixer, cool to room temperature, and move the mixture to a clean bottle. Store it at 4 °C.

108

1.8 Determine the size and morphology of the AuNPs by ultraviolet spectroscopy and TEM imaging.

111

NOTE: Different sizes of colloid gold particles for various applications can be obtained by changing the proportion of citrate sodium and chloroauric acid.

114115

2. Synthesis of AuNPs-mAb Conjugate

116

NOTE: Since antibodies bind to colloidal gold by electrostatic adsorption, charges on the surface of proteins and colloidal gold directly affect the binding intensity; therefore, the buffer pH value is an important factor affecting the stability of the antibody-colloidal gold conjugate. SSD and anti-SSD mAbs are used as examples in this protocol.

121

122 2.1 Evaluation of coupling pH

2.1.1 Add 100 μL of NaCl solution (10% m/v) into eight tubes.

2.1.2 Adjust the AuNP solutions at pH 5, 6, 7, 8, 9, 10, 11 and 12 with 0.1 M K₂CO₃.

2.1.3 Add 100 μL of colloidal gold solution (pH adjusted 5–12) to eight tubes containing NaCl.

2.1.4 Allow the solutions to stand for several minutes after blending. Observe the color change of each tube solution, and record the tube that remains red.

- 2.1.5 Choose the pH value with the least addition of K₂CO₃ and stable solution color as the
- optimum pH value for preparing AuNP-mAb conjugates.

NOTE: Do not use a pH meter because the probe may be destroyed by AuNPs.

2.2 Evaluation of antibody amount

2.2.1 Add 100 μL of NaCl solution (10% m/v) to eight tubes.

2.2.2 Add 100 μL of colloidal gold solution with optimum pH to each tube.

2.2.3 Add the monoclonal antibody solutions (protein concentration 0.1 mg/mL-3.2 mg/mL) to the abovementioned eight tubes.

2.2.4 Allow the solutions to stand for several minutes after blending. Observe the color change of each tube solution, and record the tube that remains red.

2.2.5 Choose the antibody amount with the lowest concentration of mAb and stable solution color as the optimum mAb amount for preparing AuNP-mAb conjugates.

2.3 Prepare the resuspension buffer: add 1 M Tris-HCl (pH 8.8), 1% (w/v) BSA, 0.5% (v/v) Tween 20 and 1% (v/v) PEG 20000 and blend well.

2.4 Synthesis of AuNP-mAb conjugate

2.4.1 Take 10 mL of colloidal gold solution and use 0.1 M K₂CO₃ to adjust the solution to the optimal pH value.

2.4.2 Slowly add SSD mAbs at appropriate concentrations and shake at room temperature for 30 min.

2.4.3 Centrifuge the mixture at 83 x g (1,000 rpm) for 10 min at 4 °C. Remove the precipitate, which contains impurities or precipitated colloidal gold.

166

2.4.4 Centrifuge the mixture at 8,330 x g (10,000 rpm) for 30 min at 4 °C. Discard the supernatant, and the precipitate is the colloidal gold-mAb conjugate.

169

2.4.5 Add the resuspension buffer to dissociate the precipitation.

171

3. Assembly of the Strip

172173

NOTE: For later flow immunoassays, the selection and pretreatment of membrane material will directly affect the test, which should be investigated. The immunochromatographic strip consists of a sample pad, a conjugate pad, a nitrocellulose (NC) membrane, an absorbent pad and a PVC board (Figure 1). The membrane material should be checked and evaluated by stereomicroscopy to eliminate inhomogeneity.

179

3.1. Paste the NC membrane on a PVC board 2 cm apart from the edge of the suction end of the board.

182

3.2 Add SSD-BSA (2 mg/mL) dropwise to the NC membrane (2 cm apart from the upper edge) as a test line (1 mm wide), and add goat anti-mouse IgG (1.5 mg/mL) dropwise on the NC membrane (2 cm apart from the lower edge) as a control line (1 mm wide). Control the amount of protein added.

187

3.3 Attach the absorbent pad to the PVC sheet above the NC membrane and overlap it with theNC membrane by 2 mm.

190

3.4 Submerge the glass fiber membrane into the AuNPs-mAb conjugate solution. Dry the wet membrane in an incubator at 37 °C.

193

3.5 Trim the glass fiber membrane to 5 cm long and 2 cm wide and use it as a conjugate pad.

195

3.6 Paste the pretreated conjugate pad under the NC membrane. The length of overlap with the NC membrane should be 0.1 cm.

198

NOTE: The glass fiber membrane has a strong ability to bind and release proteins.

200

3.7 Trim the fusion 3 membrane to 1.8 cm long and 3.5 cm wide and use it as the sample pad.

202

3.8 Paste the sample pad on the PVC board and overlap it with the conjugate pad by 2 mm.

3.9 Cut the assembled paper board into 3.5-mm-wide strips using a cutting machine and compact it using a batch lamination system.

207

3.10 Finally, place the test strips into the shell, seal them in an aluminum foil bag containing desiccant, and store them away from light. The ICSs are now assembled.

210

NOTE: The above is the laboratory procedure. In production, gold spraying equipment and cross-membrane instruments are used to spray gold and make the T and C lines.

213214

4. Quantitative Detection

215

216 4.1 Drop 50 μ L of sample solution onto the sample hole to observe the chromatography 217 process.

218

NOTE: As a result of the capillary action driven by the absorbent pad, the sample solution migrates to the other end of the strip. When the sample solution reaches the conjugate pad, the SSD (antigen) in the sample reacts with the AuNPs-mAb preloaded on the pad. When the solution migrates and reaches the T-line, the AuNPs-mAb without SSD can be selectively captured by SSD-BSA (antigen-carrier protein conjugate), showing as a red color on the T-line. Then, the solution migrates to the C line, where the AuNPs-mAb are captured by goat anti-mouse IgG in the region, thus showing as a red color on the C line.

226

4.2 Analyze the strips with a portable strip reader. The machine can provide the ratio of the test line to the control line (T/C).

229

4.3 Evaluate the specificity, sensitivity, repeatability and stability of the ICS test.

231

NOTE: Under qualitative detection, one red line indicates a positive result (control line). Two red lines indicate a negative result (test and control lines). If the control line is not present, the test is considered invalid.

235236

REPRESENTATIVE RESULTS:

237 Characterization of colloidal gold

- The prepared colloidal gold solutions were claret red. TEM analyses were used to determine the morphology and shape of AuNPs (Figure 2A-D). Figure 2A and Figure 2B reveal that the particles are polyhedral in shape and uniformly distributed. The average diameter of AuNPs was found to be approximately 14 nm (Figure 2C). A high-resolution TEM (HRTEM) image of AuNPs is shown in Figure 2D and Figure 2E. The HRTEM image taken from an individual AuNP shows a continuous fringe pattern with a spacing of 0.117 nm. The UV-vis absorption peaks of the five gold colloidal solutions were 518 nm, 521 nm, 524 nm, 534 nm and 540 nm, which showed that
- 245 the particle size increased slightly with decreasing sodium citrate (Figure 2F).

246247 [Figure 2 here]

Evaluation of the dipsticks

Sensitivity of the ICS

To determine the sensitivity of ICS, immunochromatographic strips were used to detect a variety of different concentrations of SSD (150,000, 60,000, 12,000, 2400, 480, and 96 ng/mL) standard samples (**Figure 3A**). Double-distilled water was used as the control. In the quantitative experiment, the strips were scanned by a JY1502GS portable ICS reader (**Figure 3B-C**). The linear regression equation was y = -0.113ln(x) + 1.5451, with a correlation coefficient (R2) of 0.983 (**Figure 3D**). It displayed good linearity at 96 ng/mL to 150 µg/mL. The IC50 value was 10.39 µg/mL. In the quantitative experiment, the optimal test time was suggested to be 10 min.

[Figure 3 here]

Specificity of the ICS

The specificity of the ICS was tested by evaluating the cross-reactivity with SSD-similar compounds. **Table 1** shows the cross-reactivity of ICS with SSD-related compounds. It can be clearly seen that the ICS only had low cross-reactivity with SSa and had no cross-reaction with other compounds, including SSc, SSb1, or SSb2 (**Table 1**), indicating that the prepared ICS has high specificity.

[Table 1 here]

Recovery rate

As shown in **Table 2**, the average recovery rate was 102.05% (mean \pm SD, n = 3). On the basis of its accuracy and consistency, this assay was sufficiently reliable for the determination of SSD in biological samples.

[Table 2 here]

Stability analysis of the ICS assay

To evaluate the stability of the ICS, the test strips stored for 8 and 16 weeks were tested and compared with the newly prepared strips. **Table 3** shows that the results of the stored and newly prepared strips for the negative and positive samples were essentially unchanged, indicating that the ICS can be stored at room temperature for at least several months and that it is suitable for large-scale promotion and application.

[Table 3 here]

Figure 1 Schematic diagram of the immunochromatography assay (**A**) Sandwich immunochromatographic test strips. (**B**) Indirect competitive immune chromatographic test strips. This figure has been modified from Zhang *et al.*, 2018²¹.

Figure 2 Characterization of colloidal gold. (A) TEM image of AuNPs (100,000× magnification). **(B)** TEM image of AuNPs (500,000× magnification). **(C)** The size distribution of AuNPs. As seen in the TEM image, the diameters of the AuNPs ranged from 10 nm to 18 nm, with an average diameter of approximately 14 nm. **(D, E)** High-resolution TEM (HRTEM) image of AuNPs. The HRTEM image taken from an individual AuNP shows a continuous fringe pattern with a spacing of 0.117 nm. **(F)** UV-vis absorption spectra of different AuNPs ranging from 10 nm to 18 nm. This figure has been modified from Zhang *et al.*, 2018²¹.

Figure 3 Characterization of the lateral-flow immunoassay for SSD. (A) Photographs of results for standard solutions containing different concentrations of SSD assayed using the ICS. (B) Photograph of the matched colloidal gold scan. (C) The intensity patterns of the test and control lines scanned by the colloidal gold quantitative instrument. (D) Standard curve of icELISA for SSD determination using the ICS. The regression equation is $y = -0.113\ln(x) + 1.5451$, with a correlation coefficient (R²) of 0.983. This figure has been modified from Zhang *et al.*, 2018²¹.

Table 1. Cross-reactivity of compounds measured by ICS and icELISA

The specificity of the ICS was evaluated by ICS (a) and icELISA (b). This table has been modified from Zhang et al., 2018^{21} .

Table 2. Recovery rate of SSD

Data are the mean \pm SD from triplicate samples at each spiked concentration of SSD. The percentage of recovery was calculated as follows: recovery (%) = measured amount/amount × 100%. This table has been modified from Zhang et al., 2018²¹.

Table 3. Variations among ICS used for the analysis of SSD

^aThe values indicate coefficients of variance for triplicate samples on 3 different strips used 1 day after manufacture. ^bThe values indicate coefficients of variance for triplicate samples on 3 different strips used after being stored for 4 weeks. ^cThe values indicate coefficients of variance for triplicate samples on 3 different strips used after being stored for 8 weeks. This table has been modified from Zhang *et al.*, 2018²¹.

DISCUSSION

In this work, we present a protocol for the preparation of mAbs against natural product-derived small molecules. The essential steps and the matters needing attention in the procedure have been outlined, and we have demonstrated the utility of this protocol using the small molecule SSD as an example. Example spectra, TEM images, quantitative results and methodological investigations are shown in representative data. Hence, we have demonstrated that the colloidal gold production, AuNP-mAb conjugation and strip assembly strategy presented here result in the creation of a novel immunoassay that can be used to rapidly detect the target molecule.

Immunochromatographic strips are easy to operate, and the reaction is very rapid. Only 50 µL of sample is needed, and detection is complete within 5 minutes. In contrast, the ELISA method, which is also an immunoassay, requires several hours, including multiple washes and long incubations. Multiple washes and long incubation periods can effectively remove the nonspecific reaction of the sample matrix to improve the signal-to-noise ratio and the sensitivity by increasing the signal value. This means that the anti-interference ability and optimization technology of tomographic paper are the key problems to be solved in this study. Technological optimizations can be achieved in several ways. Select suitable materials for the sample pad and conjugate pad. The length of the sample pad was increased appropriately to effectively improve the pretreatment ability of the sample. High pH, high ionic strength and high buffering capacity of material pretreatment can be used to avoid influencing factors such as the antigen—antibody response of the pH value of the biological sample. As the conditional setting of "three high" can interfere with the binding of antibodies, it should be considered comprehensively.

Moreover, it is also important to note that a monoclonal antibody (mAb) is the probe on the strip. Therefore, this procedure is associated with the sensitivity of the mAb. The mAb we applied was prepared in our previous work¹⁷, and mAb production can be referenced in our previous publication¹⁹.

Colloidal gold-based immunochromatographic strips are widely used strips because they are inexpensive and offer convenient detection. However, this method has its own limitations and stability issues. The signal strength of AuNPs is lower than those of fluorescent dyes, colloidal carbon and other markers, meaning that the sensitivity is limited. When exposed to air, the AuNPs change from wine red to purple, which will cause a change in the gray value and affect the quantitative detection. Hence, various modified immunochromatographic strips have emerged, and we also prepared quantum dot-labeled fluorescent strips²² and IgY-based immunochromatographic strips²³, which have similar procedures for production.

Membrane-based immunoassays have been widely used in clinical indicators, infectious diseases, pesticide residues and drug concentration monitoring. In this protocol, we focused on the preparation of an AuNP-based method for small molecule natural products. This procedure acts as a template for the development of various target immunoassays that can be effectively utilized as analytical tools for the evaluation and quality control of natural products.

This work provides a detailed protocol for the development of lateral-flow immunoassays based on immune chromatographic strips for use in the rapid, sensitive and quantitative detection of

small molecules. This procedure includes preparation of colloidal gold, synthesis of AuNP-mAb conjugates, assembly of the strip, and methodological investigation.

370 371

ACKNOWLEDGEMENTS:

- 372 This work was supported by the National Natural Science Foundation of China and the Classical
- 373 Prescription Basic Research Team at the Beijing University of Chinese Medicine.

374375

DISCLOSURES:


376 The authors have nothing to disclose.


377378


REFERENCES:


- 1. Huang, X. et al. Membrane-based lateral flow immunochromatographic strip with
- nanoparticles as reporters for detection: A review. *Biosensors and Bioelectronics*. **75**, 166–180
- 381 (2016).
- 382 2. Chang, H-F., Wang, J-Q., Wang, B., Deng, A-P. An immune chromatographic assay for rapid
- and simultaneous detection of levonorgestrel and methylprednisolone in water samples.
- 384 *Chinese Chemical Letters.* **24** (10), 937-940 (2013).
- 385 3. Lai, J. J., Stayton, P. S. Improving lateral-flow immunoassay (LFIA) diagnostics via biomarker
- enrichment for mHealth. *Methods in Molecular Biology*. **1256**, 71-84 (2015).
- 4. Zhang, M. Z. et al. Development of a colloidal gold-based lateral-flow immunoassay for the
- 388 rapid simultaneous detection of clenbuterol and ractopamine in swine urine. Analytical &
- 389 *Bioanalytical Chemistry*. **395** (8), 2591-9 (2009).
- 5. Kranthi, K. R. et al. Development of a colloidal-gold based lateral-flow immunoassay kit for
- 'quality-control' assessment of pyrethroid and endosulfan formulations in a novel single strip
- 392 format. *Crop Protection*. **28** (5), 428-434 (2009).
- 393 6. Qian, K. et al. Development and evaluation of an immunochromatographic strip for rapid
- detection of capsid protein antigen p27 of avian leukosis virus. Journal of Virological Methods.
- 395 **(221)** 115-118 **(2015)**.
- 7. Guo, H. et al. Lateral flow immunoassay devices for testing saliva and other liquid samples
- and methods of use of same. US Patent, US20040184954A1 (2003).
- 398 8. Miočević, O. et al. Quantitative Lateral Flow Assays for Salivary Biomarker Assessment: A
- 399 Review. *Frontiers in Public Health.* **5**, 1-13 (2017).
- 400 9. Lisa, M. et al. Gold nanoparticles based dipstick immunoassay for the rapid detection of
- dichlorodiphenyltrichloroethane: an organochlorine pesticide. *Biosensors and Bioelectronics*. **25**
- 402 (1), 224-227 (2009).
- 403 10. Zhang, Z. et al. Monoclonal Antibody-Europium Conjugate-Based Lateral Flow
- Time-Resolved Fluoroimmunoassay for Quantitative Determination of T-2 Toxin in Cereals and
- 405 Feed. *Analytical Methods*. **7** (6), 2822-2829 (2015).
- 406 11. Shen, H. et al. Facile synthesis of high-quality CuInZnxS2+x core/shell nanocrystals and their
- 407 application for detection of C-reactive protein. Journal of Materials Chemistry. 22 (35),
- 408 18623-18630 (2012).
- 409 12. Xiang, T. et al. A novel double antibody sandwich-lateral flow immunoassay for the rapid and
- 410 simple detection of hepatitis C virus. International Journal of Molecular Medicine. 30 (5),

- 411 1041-1047 (2012).
- 412 13. Yang, Q. et al. Quantum dot-based immunochromatography test strip for rapid, quantitative
- and sensitive detection of alpha fetoprotein. *Biosensors & Bioelectronics*. **30** (1), 145 (2011).
- 14. Song, L.W. et al. Rapid fluorescent lateral-flow immunoassay for hepatitis B virus genotyping.
- 415 *Analytical Chemistry*. **87**, 5173-5180 (2015).
- 416 15. Zhang, Y. et al. Quantum dot-based lateral-flow immunoassay for rapid detection of rhein
- using specific egg yolk antibodies. Artificial Cells, Nanomedicine, and Biotechnology. 1 (2017).
- 418 16. Qu, H. et al. Rapid Lateral-Flow Immunoassay for the Quantum Dot-based Detection of
- 419 GsRerarin. *Biosensors and Bioelectronics*. **81**, 358-362 (2016).
- 420 17. Li, Z. et al. Development and Clinical Application of a Rapid IgM-IgG Combined Antibody Test
- for SARS-CoV-2 Infection Diagnosis. Journal of Medical Virology. 92 (9) (2020).
- 422 18. Xiaomei, L., Jing., W., Ya., Z. The clinical application value analysis of the 2019-coronary virus
- disease was analyzed by the whole blood Sars-COV 2 specific antibody detection. *Natural*
- 424 Science Edition. **42** (2020).
- 425 19. Zhang, Y. et al. Generation of Monoclonal Antibodies Against Natural Products, Journal of
- 426 *Visualized Experiments.* e57116 (2019).
- 427 20. Sai, J. et al. Development of an Enzyme-Linked Immunosorbent Assay and Immunoaffinity
- 428 Column Chromatography for Saikosaponin d Using an Anti-Saikosaponin d Monoclonal Antibody.
- 429 *Planta Medica*. **82**, 432–439 (2016).
- 430 21. Yue, Z. et al. A Highly Sensitive Immunochromatographic Strip Test for Rapid and
- 431 Quantitative Detection of Saikosaponin d. *Molecules*. **23** (2), 338 (2018).
- 432 22. Qu, H. et al. Rapid Lateral-Flow Immunoassay for the Quantum Dot-based Detection of
- 433 Puerarin. *Biosensors and Bioelectronics*. **81**, 358-362 (2016).
- 434 23. Zhang, Y. et al. Quantum dot-based lateral-flow immunoassay for rapid detection of rhein
- using specific egg yolk antibodies. *Artificial Cells, Nanomedicine, and Biotechnology*. **1** (2017).









| Samples | ^a ICS | ^b icELISA | |
|---------|------------------|----------------------|--|
| | (%) | (%) | |
| SSd | 100 | 100 | |
| SSa | 4.30% | 4.97% | |
| SSb1 | < 0.09 | <0.09 | |
| SSb2 | < 0.09 | <0.09 | |
| SSc | < 0.09 | < 0.09 | |

| SSd | SSd concentration | | | |
|---------------|---------------------|----------------|--|--|
| concentration | established by the | Recovery (%) | | |
| (ng/mL) | test system (ng/mL) | | | |
| 100 | 135.72 ± 61.97 | 135.72 ± 61.97 | | |
| 1000 | 926.59 ± 114.24 | 92.66 ± 11.42 | | |
| 10000 | 11128.16 ± 745.75 | 111.28 ± 14.12 | | |

| _ | RSD % | | |
|-------------|--------------------|----------------------|----------------------|
| SSd (ng/mL) | 1 day ^a | 4 weeks ^b | 8 weeks ^c |
| 125 | 2.41 | 3.11 | 3.51 |
| 250 | 2.52 | 4 | 3.52 |
| 500 | 2.44 | 3 | 5.2 |
| 1000 | 3.12 | 2.71 | 4.5 |

Table of Materials

Click here to access/download **Table of Materials**Materials.xlsx

Manuscript No.:JoVE62754R1

Title: Development of lateral-flow immunochromatographic strip for Rapid

and Quantitative Detection of small molecule compounds

Dear editor,

Thank you very much for your reviewing and replying to our manuscript

(JoVE62754R1). We would like to express our sincere gratitude to the editor and

reviewers for their constructive and positive comments. We replied point by point to

the comments, revised the manuscript and video accordingly and the changes were

highlighted in red in the revised manuscript. We sincerely hope this revised

manuscript is more acceptable for publication.

Thank you again and all the reviewers for the kind advice and look forward to

hearing from you soon.

Yours sincerely,

Huihua Qu

Editorial comments:

Changes to be made by the Author(s):

1. Please employ professional copyediting services as the language in the manuscript is not publication grade.

Thank you for your comments, the manuscript has been edited by Elsevier language editing service. The language editing certificate have been attached.

2. Please provide reprint permissions for the figures and tables.

Thank you for your kind reminder, the figures is from *molecules* which is an open-access journal. According to MDPI open access information and policy, "Everyone is free to re-use the published material if proper citation of the original publication is given". As we have cited the original articles, the figures are not need reprint permissions.

Changes to be made by the Author(s) regarding the video:

1. Please change the title in the video to match the written manuscript: Development of a Lateral Flow Immunochromatographic Strip for Rapid and Quantitative Detection

of Small Molecule Compounds

Thank you for your kind reminder, we have modified them.

2. Please ensure that the protocol subheaders in the video match the protocol subheaders in the written manuscript. Use American English (characterization instead of characterisation)

Thank you for your kind reminder, we have modified them.

- 3. Video Editing Content:
- 03:50 Please remove a glitch happening between the shots.
- 05:08 05:13 Please remove dead air.
- Please use some image or video for conclusion section. Do not complete the section just by showing the chapter title card. Using of image or video on conclusion will be more effective for the audience.
- Please use cross dissolve on hard cut for making the video look smooth."

Thank you for your kind reminder, we have modified them.

- 4. Audio Editing and Pacing:
- Audio levels is low. Please increase the audio level and ensure audio level peaks average around -9 dB.

Thank you for your kind reminder, we have modified them.

Once done please ensure that the revised video is no more than 15 min in length. Please upload the revised file to https://www.jove.com/account/file-uploader?src=19140018

We have uploaded the revised video named JOVE62754R2.