# **Journal of Visualized Experiments**

# Adsorptive bioprocess improves yield of melanin from Pseudomonas stutzeri -- Manuscript Draft--

Article Type:	Invited Methods Collection - JoVE Produced Video		
Manuscript Number:	JoVE63339R3		
Full Title:	Adsorptive bioprocess improves yield of melanin from Pseudomonas stutzeri		
Corresponding Author:	Keyur Raval National Institute of Technology Karnataka Mangalore, Karnataka INDIA		
Corresponding Author's Institution:	National Institute of Technology Karnataka		
Corresponding Author E-Mail:	keyurnraval@nitk.edu.in		
Order of Authors:	Harsha Thaira		
	Keyur Raval		
	Ritu Raval		
Additional Information:			
Question	Response		
Please specify the section of the submitted manuscript.	Engineering		
Please indicate whether this article will be Standard Access or Open Access.	Standard Access (\$1400)		
Please indicate the city, state/province, and country where this article will be filmed. Please do not use abbreviations.	Mangalore, Karnataka, India		
Please confirm that you have read and agree to the terms and conditions of the author license agreement that applies below:	I agree to the Author License Agreement		
Please confirm that you have read and agree to the terms and conditions of the video release that applies below:	I agree to the <u>Video Release</u>		
Please provide any comments to the journal here.			

# TITLE:

Adsorptive Bioprocess Improves Yield of Melanin from *Pseudomonas stutzeri* 

# **AUTHORS AND AFFILIATIONS:**

5 Harsha Thaira<sup>1</sup>, Ritu Raval<sup>2</sup>, Keyur Raval<sup>1</sup>

<sup>1</sup>Department of Chemical Engineering, National Institute of Technology Karnataka 575025, India <sup>2</sup>Department of Biotechnology, Manipal Institute of Technology, MAHE, Manipal – 576104, India

- 10 Email addresses of co-authors:
- Harsha Thaira (harshathaira@gmail.com)
   Ritu Raval (ritu.raval@manipal.edu)
   Keyur Raval (keyurnraval@nitk.edu.in)

- 15 Corresponding author:
- 16 Keyur Raval (keyurnraval@nitk.edu.in)

### SUMMARY:

Here we present a protocol on adsorptive bioprocess to produce a biological product, especially a pigment that inhibits its own biosynthesis and inhibits the growth of the microorganism that produces it.

# **ABSTRACT:**

Melanins are natural pigments, and the presence of indole ring and numerous functional groups makes melanin an ideal choice for many applications such as UV protective agents, skincare, cosmetics etc. A marine *Pseudomonas stutzeri* produces melanin without the addition of tyrosine. The feedback inhibition was observed by melanin in the culture of a melanin-producing marine bacterium, *Pseudomonas stutzeri*. Melanin also demonstrated microbial growth inhibition. The Han-Levenspiel model-based analysis identified uncompetitive type product inhibition of melanin on the cell growth. Tyrosinase enzyme, which produces melanin, was inhibited by melanin. The double reciprocal plot of the enzymatic reaction in the presence of different melanin concentrations revealed uncompetitive product inhibition. An adsorbent-based adsorptive bioprocess was developed to reduce the feedback inhibition by melanin. Different adsorbents were screened to select the best adsorbent for melanin adsorption. Dosage amount and time were optimized to develop the adsorptive bioprocess, which resulted in an 8.8-fold enhancement in melanin production by the marine bacteria *Pseudomonas stutzeri* (153 mg/L to 1349 mg/L) without supplementation of tyrosine and yeast extract.

# **INTRODUCTION:**

Melanins are polyphenol compounds whose monomeric unit is the indole ring. Melanins are responsible for most of the black, brown, and grey colorations of plants, animals, and microorganisms<sup>1</sup>. The most common melanins are eumelanins, which are dark brown or black pigments widely distributed in vertebrates and invertebrates<sup>2</sup>. Cephalopods, plants, and microbes are the major sources of melanin<sup>3</sup>. Melanin obtained from microorganisms has several

advantages over melanin from cephalopods and plants. They are fast in growth and do not cause any problems of seasonal variations. Also, they modify themselves according to the growth medium and the operating conditions provided. Considering all these facts, the microbial sources are exploited as a major source of melanin with potential commercial applications<sup>4–6</sup>. Guo et al. used recombinant *Streptomyces kathirae* SC-1 to produce melanin up to 28.8 g/L in an optimized medium<sup>7</sup>. Lagunas-Muñoz et al. Obtained about 6 g/L of melanin in recombinant *E. coli* <sup>8</sup>. Wang et al. obtained about 4.2 g/L of soluble melanin after medium optimization<sup>9</sup>. Researchers also investigated marine cultures to produce melanin. Kiran et al. produced melanin from marine *Pseudomonas* species<sup>10</sup>. Kumar et al. isolated *Pseudomonas stutzeri* HMGM-7 from marine seaweed, which produced melanin in nutrient broth. All the reported melanin producers require medium supplements, yeast extract and/or peptone for melanin production, which may hinder the scale-up of the bioprocess. Raman et al. investigated melanin production from *Aspergillus fumigatus* AFGRD105, which did not require costly additives<sup>11</sup>.

A marine bacterium known as *Pseudomonas stutzeri*, isolated from seaweed, was selected for melanin production. This microorganism does not require additional supplementation of tyrosine and yeast extract for the production of soluble melanin and grows in seawater-based media, thereby reducing the potable water footprint and medium cost. Thaira et al. reported an optimized bioprocess using coconut cake meal for melanin nanopigment production using this marine bacterium in a previous publication<sup>12</sup>. The produced melanin was used for the highly efficient removal of heavy metals from the synthetic groundwater. However, during medium optimization, plausible feedback inhibition of melanin on the growth and product formation was identified. Therefore, in this protocol, feedback inhibition of melanin is investigated and quantified using the generalized Han-Levenspiel model. A simple and effective adsorptive bioprocess protocol is developed to reduce feedback inhibition to the maximum extent.

# **PROTOCOL**:

1. Growth and maintenance of *Pseudomonas stutzeri* culture.

1.1. Obtain the bacteria *Pseudomonas stutzeri* from MTCC Chandigarh (supplied as freezedried ampoule).

1.1.1. Disinfect the surface of the ampoule with alcohol (70%). Mark the ampoule near the middle of the gauge and break at the marked area.

NOTE: Care should be taken in opening the ampoule as the contents are in vacuum. Take precautions while breaking the ampoule, as a snap opening will draw the cotton plug to one end and a hasty opening will release the lyophilized fine particles of the microbe into the air.

1.1.2. Carefully remove the cotton plug with sterilized tweezers and add 0.3 mL of Luria Bertani (LB) broth to the ampoule. Let it stand for 30 min.

- 1.1.3. Inoculate 100 µL of the suspension to 10 mL of LB broth (in 100 mL conical flask) and 88 place it in the incubator shaker at 180 rpm maintained at 37 °C. Simultaneously transfer 100 μL 89 90 of the suspension to the solid LB agar and spread to obtain isolated colonies. Incubate it at 37 °C 91 for overnight.
- 1.2. Once the culture has grown, store them at 4 °C. 93
- 95 Do subculture after every 4 weeks by streaking a colony onto a fresh LB agar plate and 96 incubating the plate at 37 °C for overnight for the colonies to appear.
  - 2. Inhibition studies of melanin on the growth and enzyme activity
- 100 Prepare the culture medium by adding 51.4 g/L glucose, 49.8 g/L coconut cake meal, 0.09 g/L, 0.03 g/L and 0.05 g/L of CuSO<sub>4</sub>, FeSO<sub>4</sub> and MgSO<sub>4</sub>, respectively, in artificial seawater (**Table** 101 102 of Materials). Perform all the experiments in a 250 mL shake flask with 15 mL filling volume at 103 200 rpm, 37 °C, and initial pH of 7.
- 2.2. Add different concentrations of melanin (100 mg/L, 200 mg/L, 500 mg/L, and 1000 mg/L) 105 106 into each culture medium at the start of the bioprocess.
- Measure specific growth rate, the enzyme activity of the cells, and melanin at different 108 109 time intervals of the bioprocess for each added melanin concentration (section 3.2).
  - 3. Modeling melanin inhibition kinetics
- Use the following generalized equation of Han and Levenspiel<sup>13</sup> to model the inhibition 113 3.1. 114 kinetics of melanin production,

115 
$$\mu = \mu_{max} \cdot \left(1 - \frac{p}{p^*}\right)^n \cdot \frac{s}{\left(s + K_s \cdot \left(1 - \frac{p}{p^*}\right)^m\right)}$$
 (1)

116 Where,

92

94

97 98

99

104

107

110 111

112

- 117  $\mu$  = apparent specific growth rate (1/h)
- 118  $\mu_{max}$  = maximum specific growth rate (1/h)
- 119 p = product concentration (g/L)
- $p^*$  = critical product concentration at which growth ceases (g/L) 120
- 121 s = substrate concentration (g/L)
- 122 m,n = parameters
- Above equation can be rewritten as 123

124 
$$\mu = \mu_{obs} \cdot \frac{s}{(s + K_{sobs})}$$
 (2)

126 
$$\mu_{obs} = \mu_{max} \left( 1 - \frac{p}{p^*} \right)^n$$
 (3)  
127 and  $K_{obs} = K_s \cdot \left( 1 - \frac{p}{p^*} \right)^m$  (4)

127 and 
$$K_{obs} = K_s \cdot \left(1 - \frac{p}{p^*}\right)^m$$
 (4)

- 128 The above equation can be rearranged to a double reciprocal plot of  $1/\mu$  versus 1/s, as given
- 129 below
- 130  $\frac{1}{\mu} = \frac{1}{\mu_{obs}} + \frac{K_{obs}}{\mu_{obs}} \cdot \frac{1}{s}$  (5)
- 131
- 132 3.2. Draw a straight line using the data of different specific growth rates and substrate
- concentrations using the above equation. The intercept on the Y- and the X-axis, respectively,
- give  $\mu_{max}$  and  $K_s$ , which can be substituted in eqution 3 and 4 respectively to obtain  $\mu_{obs}$  and
- 135  $K_{obs}$ .
- 136
- 137 Equation (3) can be written as
- 138  $\ln(\mu_{obs}) = \ln(\mu_{max}) + n \cdot \ln\left(1 \frac{p}{p^*}\right)$  (6)
- 139
- 140 3.3. Calculate the specific growth rate of culture grown in step 2.2 and name it as  $\mu_{obs}$ .
- 141
- 142 3.4. Use the value of the specific growth rate obtained in step 3.3 in equation (6). Fit the data
- to a linear equation and calculate the slope and intercept.
- 144
- 145 3.5. Use equation (3) to calculate the maximum specific growth rate in the absence of melanin.
- 146147
- 4. Development of the adsorptive bioprocess
- 148
- 4.1. Perform all the experiments in a 250 mL shake flask with 15 mL filling volume at 200 rpm, 37 °C, and initial pH of 7.
- 151
- 4.2. Conduct the final optimized bioprocess in a 5 L stirred tank bioreactor with a 3 L filing
- volume. Maintain the pH of the culture between 7–7.5 and dissolved oxygen between 50%– 60%
- of the saturation in the bioreactor studies.
- 155
- 156 **5.** Selection of the best adsorbent from different available adsorbents
- 157
- 158 5.1. Add different adsorbents such as alumina, zeolite, celite, fuller's earth, and activated
- carbon in different flasks containing the culture broth.
- 160
- 161 5.2. Add each adsorbent separately at a concentration of 1 g/L to the culture broth and keep it in a rotary shaker at 150 rpm for 12.
- 163
- 5.3. Separate the adsorbents by centrifugation at 5000 x g for 10 min and measure the
- absorbance of the supernatant at 400 nm using a spectrophotometer.
- 166
- 5.4. Select the adsorbent which adsorbs the maximum amount of melanin from the culture broth.
- 169
- 170 6. Effect of dosage time of the adsorbent

171

172 6.1. Start the bioprocess and add the adsorbent yielding the highest adsorption of melanin 173 from the culture (selected after section 5). The concentration selected at the initial step is 1 g/L.

174

175 6.2. Measure the maximum melanin production according to the previously published protocol<sup>12</sup>.

177

178 6.3. Repeat step 6.1 by adding the adsorbent after 12 h of the start of the bioprocess and measure the melanin production at different times of the bioprocess.

180

181 6.4. Repeat step 6.3 by adding the adsorbent after 24 h, 36 h, and 48 h of the start of the bioprocess.

183

184 6.5. Compare the melanin produced in each bioprocess run and identify the best dosage time for the maximum melanin production and recovery.

186

7. Effect of amount of adsorbent dosage

187 188

7.1. Add different concentrations of the best adsorbent (0.5, 1, 2, 4, 6, 8, and 10 g/L) to the cultures at the time identified in section 6.

191

192 7.2. Harvest the bioprocess by centrifugation at 5000 x g for 10 min.

193

194 7.3. Collect the supernatant and calculate the amount of melanin by measuring absorbance at 400 nm.

196

7.4. Wash the pellet containing adsorbent two times each with distilled water and ethanol to remove cells and other compounds attached to the adsorbent. Perform the washing step by centrifugation at 5000 x g for 10 min.

200

7.5. Add a 0.1 M NaOH solution to the washed adsorbent to desorb the melanin.

201202203

7.6. Centrifuge the mixture at  $5000 \times g$  for 10 min and measure the melanin amount in the supernatant at 400 nm by a spectrophotometer.

205

204

206 7.7. Repeat the desorption process thrice or until the adsorbent is back to its normal color.

207

7.8. Reduce the pH of this concentrated desorbed melanin solution to 3 or less by adding 1 N
 HCl using a stirrer with magnetic beads to precipitate melanin. Store the acidic solution in the
 refrigerator for 1 h to increase precipitation, if required.

211

7.9. Wash the precipitated melanin thus obtained three times with distilled water by centrifugation at  $5000 \times g$  for 10 min, dry and store at 4 °C.

# 8. Mathematical modeling of the bioprocess

217 8.1. Biomass growth kinetics

8.1.1. Quantify the biomass growth rate using the following modified Verhulset model developed for product inhibition, as described by the equation given below<sup>14–17</sup>,

$$\frac{dx}{dt} = \mu_{obs} \cdot x \cdot \left(\frac{x_m - x}{x_m}\right) \tag{7}$$

- 222 Where,
- x = biomass concentration (g/L)
- $x_m$  = maximum biomass concentration (g/L)
- $\mu_{obs}$  = maximum specific growth rate (1/h) as given by equation (2)
- 226 Vazques et al. (2008) further derived the following equation for biomass concentration from
- 227 equation (7)

$$x = \frac{x_m}{1 + exp\left[2 + \frac{4 \cdot v_x}{x_m} \cdot (\tau_x - t)\right]} \tag{8}$$

- 229 Where,
- $v_x$  = biomass growth rate (g/L/h) =  $\mu_{obs}$ ·x
- $\mu_{obs}$  = observed specific growth rate (1/h)
- $\tau_x = \text{lag time (h)}$
- t = time (h)

8.1.2. Use Equation (8) to calculate the biomass concentration at different times of the bioprocess and simulate the growth profile of *Pseudomonas stutzeri*.

238 8.2. Substrate consumption kinetics

NOTE: The substrate is used for biomass growth, product formation, and maintenance in a bioprocess.

8.2.1. Use the following equation to balance substrate utilization in the bioprocess,

$$-\frac{ds}{dt} = \frac{1}{Y_{P/S}} \cdot \frac{dp}{dt} + \frac{1}{Y_{X/S}} \cdot \frac{dx}{dt} + m_S \cdot x \tag{9}$$

The majority of the glucose is used for biomass generation; therefore, for simplification, neglecting substrate consumption for product formation would reduce equation (9) to,

$$-\frac{ds}{dt} = \frac{1}{Y_{Y/S}} \cdot \frac{dx}{dt} + m_S \cdot x \tag{10}$$

- 249 Where,
- S = substrate concentration (g/L)
- $Y_{P/S}$  = product yield coefficient, i.e., melanin yield coefficient (g/g)
- P = product concentration (g/L), i.e., melanin concentration (g/L)
- $Y_{X/S}$  = biomass yield coefficient (g/g)
- $m_s$  = maintenance coefficient (g/g/h)

Integrating the above equation from time t = 0 to time t = t, and substituting the value of x from equation (8) gives the following equation for substrate concentration profile<sup>17</sup>,

258 
$$s = s_0 - \frac{1}{Y_{X/S}} [x - x_0] - \left(\frac{m_S \cdot x_m^2}{4 \cdot v_X}\right) ln \left[\frac{x_0 \cdot \left(exp^{\frac{4 \cdot v_X}{xm} \cdot t} - 1\right) + x_m}{x_m}\right]$$
 (11)

- Where,  $s_0$  = initial substrate concentration (g/L) and  $x_0$  = initial biomass concentration (g/L)
- 8.2.2. Use equation (11) to calculate the substrate concentration at different time intervals of the bioprocess or if substrate concentration is known, then use the same equation to estimate the parameters given in the equation.
- 265 8.3. Melanin production kinetics
- 8.3.1. Use Ludeking-Pirate equation to simulate melanin production. However, melanin is growth associated and therefore, the non-growth associated product formation term is neglected, and the simplified equation can be written as,
- $\frac{dp}{dt} = \alpha \cdot \frac{dx}{dt} \tag{12}$
- 271 Where,

260

264

266

276

278

280

283

290

- 272  $\alpha$  = Ludeking-Pirate constant linked to growth associated product formation (g/g)
- Integrating equation (12) from time t = 0 to t = t result in the following equation for product concentration profile

$$p = \alpha \cdot (x - x_0) \tag{13}$$

- 277 8.3.2. Use equation (13) to calculate  $\alpha$ .
- 279 8.4. Kinetic parameter estimation
- 8.4.1. Estimate all parameters mentioned above by the nonlinear least square method of parameter estimation using any suitable software tool (e.g., Microsoft Excel solver tool).
- 284 8.4.2. Use batch experiment data for parameter estimation and minimize the sum of squared error (SSE) using the Solver tool, as given by the following equation,
- 286 Minimize,

287 
$$SSE = \sum_{i=1}^{n} (x_{sim} - x_{exp})^{2}$$
 (14)

- 288 Where,
- 289  $x_{sim}$  = simulated value and  $x_{exp}$  = experimental value
- 291 9. Analytical techniques292
- 293 9.1. Melanin extraction
- 9.1.1. Centrifuge the growth culture at 5000 x g for 10 min to remove the biomass.

297 9.1.2. Add an equal volume of 1 M NaOH to the supernatant.

298

299 9.1.3. Acidify the alkaline supernatant to pH 2 by adding 1 N HCl to precipitate the melanin.

300

301 9.1.4. Centrifuge the acidified supernatant at 12000 x q for 20 min to collect melanin.

302

9.1.5. Wash the collected melanin thrice with 5 mL of distilled water, dry overnight at room temperature (RT) and store at 4 °C<sup>18</sup>.

305

306 9.2. Melanin quantification

307

9.2.1. Prepare different concentrations of melanin solutions (0–500 mg/L) by dissolving the appropriate amount of melanin in artificial seawater (pH 8).

310

311 9.2.2. Measure the absorbance at 400 nm using spectrophotometer<sup>18</sup>.

312

313 9.2.3. Prepare calibration curve using above data.

314

9.2.4. Centrifuge the growth culture at 5000 x *g* for 10 min to remove the biomass and measure the absorbance of the supernatant at 400 nm. Calculate the melanin concentration in the growth culture using the calibration curve prepared in step 9.2.3

318

319 9.3. Tyrosinase assay

320

- 9.3.1. Follow the protocol given below for the tyrosinase enzyme assay developed by Ren et al. (2013)<sup>19</sup> to determine the enzyme activity of tyrosinase extracted from *Pseudomonas stutzeri*
- 323 culture broth.

324

325 9.3.2. Centrifuge the bacterial culture at 5,000 x g to collect the cells.

326

9.3.3. Sonicate the cells in cold 100 mM phosphate buffer of pH 6.5 and centrifuge the homogenate at 9000 x *g* for 15 min and collect the supernatant, a source of tyrosinase enzyme.

329

9.3.4. Add 0.2 mL of the enzyme solution obtained in the above step to 2.8 mL of 0.1 M sodium phosphate buffer (pH 6.5) containing 1 mM L-tyrosine (substrate).

332

9.3.5. Observe the formation of dopachrome by measuring the absorbance at 475 nm.

334

- 9.3.6. Calculate the tyrosinase activity as follows. One international unit (IU) of tyrosinase
   activity is defined as the amount of enzyme required to oxidize 1 μmol of L-tyrosine to
   dopachrome per minute. Calculate this by using the molar extinction coefficient of dopachrome
- 338 (3600 L·M<sup>-1</sup>·cm<sup>-1</sup>) by the following equation,

$$340 \quad \frac{U}{mL} = \frac{\left(\frac{Absorbance}{minute}\right) \cdot \left(Assay\ volume(mL)\right) \cdot (dilutionfactor) \cdot (1000)}{(3600\ L\ M^{-1}.cm^{-1}) \cdot (1\ cm) \cdot (enzyme\ volume\ (mL))}$$
(15)

341342

9.4. Estimation of biomass

343

9.4.1. Collect 5 mL of the samples at an interval of 4 h in the Biosafety cabinet, maintaining sterility. Flame the cotton plug of the flask before and after taking the samples to reduce the chances of contamination. Measure the optical density (OD) of the samples at 660 nm for all these samples.

348

NOTE: The uninoculated media serves as the blank. These readings will be used to study the growth of the microorganisms.

351

9.4.2. Centrifuge the culture broth at  $5000 \times g$  for 10 min and 4 °C, collect the cell pellet and dry the pellet for 8 h at 60 °C in a hot air oven to obtain the dry cell mass.

354

9.4.3. Prepare calibration curve using the OD measured in step 9.4.1 and the corresponding biomass obtained in step 9.4.2

357

9.4.4. Use the calibration curve (OD = 0.82@Biomass (g/L)) to measure the unknown quantity of biomass by measuring the optical density of the culture broth at any given time.

360

361 9.5. Estimation of DNA

362

363 9.5.1. Estimate the amount of DNA in the sample by measuring the optical density at 260 nm.

364

365 9.5.2. Calculate the DNA using the following formula,

366 
$$Total DNA (\mu g) = \frac{(A_{260}) \cdot \left(50 \frac{\mu g}{mL}\right) \cdot (dilution factor) \cdot (0.1 mL)}{1}$$
 (16)

367 Where.

- 368  $A_{260}$  = Absorbance at 260 nm
- 369 0.1 = Volume of the sample
- 1 = Absorbance of DNA standard solution of concentration 50 μg/mL, at 260 nm

371

372 9.6. Estimation of protein

373

9.6.1. Estimate the amount of protein in the samples by using a Bradford protein estimation kit (Table of Materials).

376

9.6.2. Add 5 mL of the Bradford reagent to 0.1 mL of the supernatant and keep it for 5 min. The proteins present in the supernatant reacts with the Bradford reagent to form a colored complex.

379

9.6.3. Measure the absorbance at  $595~\text{nm}^{20}$  using the spectrophotometer.

382 9.7. Estimation of total reducing sugar (TRS)

383

384 9.7.1. Measure the TRS by dinitrosalicylic (DNS) method<sup>21</sup>.

385

9.7.2. Prepare the DNS solution by dissolving 1 g of DNS in 50 mL of distilled water. Add 30 g of sodium potassium tartarate to this solution in small amounts till the color of the solution becomes milky yellow, followed by transparent orange-yellow with the addition of 20 mL of 2 N NaOH. Bring the volume up to 100 mL using distilled water and store it in the dark.

390

9.7.3. Prepare glucose working solution (0.01 M) by dissolving 180 mg of glucose in 100 mL of distilled water in a standard flask.

393

9.7.4. Pipette out standard Glucose solution (0.2, 0.4, 0.6, 0.8 and 1 mL) into 5 separate test tubes. Use a test tube containing distilled water(1mL) as a blank solution.

396

397 9.7.5. Bring the volume up to 2 mL in each test tube by using distilled water.

398 399

399 9.7.6. Add 1 mL of DNS reagent to each test tube and cover with aluminum foil.

400

9.7.7. Heat the test tubes in a boiling water bath for 5 min and cool to RT.

402

9.7.8. Add 9 mL of distilled water to each test tube and mix well. Measure the OD of the content
 at 575 nm using a spectrophotometer.

405

9.7.9. Measure the concentration using a calibration chart prepared with the amount of glucose on the X-axis vs. absorbance at 575 nm (A<sub>575</sub>) on the Y-axis.

408

409 9.8. Particle size analysis

410

9.8.1. Disperse the biosynthesized melanin in distilled water (1 mL) using a sonicator.

412

9.8.2. Add the dispersed sample solution (1 mL) into the nanoparticle analyzer and measure the particle size per the manufacturer's protocol.

415

416 9.9. Transmission Electron Microscopy (TEM) analysis

417

9.9.1. Dry the melanin powder obtained after extraction by keeping it at 60 °C for 8 h.

419

9.9.2. Load the melanin powder in the TEM sample matrix and measure the size using the Transmission Electron Microscope following the operator protocol.

422 423

# REPRESENTATIVE RESULTS:

- Figure 1A,B demonstrates growth and melanin production under the influence of different initial
- 425 melanin concentrations. The cell growth and melanin production decreased substantially with

the increasing initial concentration of melanin. At 500 mg/L initial melanin concentration in the medium, the biomass and melanin reduced to almost 50% of the maximum value obtained in control. At 1000 mg/L of initial melanin concentration, there was marginal melanin production. The biomass growth rate reduced significantly at 1000 mg/L initial melanin concentration. **Figure 1C** shows intracellular tyrosinase activity in culture broth at different initial melanin concentrations. As mentioned in the figure, the tyrosinase activity decreased by a factor of 4 in the presence of 1000 mg/L initial melanin concentration in the broth compared to the control. These results clearly show that the higher melanin concentration in the culture adversely affects the tyrosinase enzyme activity and cell growth.

It is necessary to identify the type of inhibition so that a solution can be obtained accordingly. Therefore, the feedback inhibition of the melanin on the tyrosinase enzyme was identified using the Han-Levenspiel model (protocol section 3.1, equation (3)). **Figure 1D** quantifies the effect of melanin inhibition on the maximum specific growth rate of microorganisms. The maximum specific growth rate of *P. stutzeri* in the optimized medium decreased from about 0.25  $h^{-1}$  to about 0.01  $h^{-1}$ . The line drawn in the figure 1D is explained by the following equation,

$$\mu_{obs} = 0.25 \left( 1 - \frac{p}{p^*} \right)^{0.91}$$
for optimized medium. (17)

The values of coefficients m and n are greater than zero, confirming uncompetitive inhibition by melanin. The critical inhibitory concentration of melanin was about 1050 mg/L. The Lineweaver-Burke plot (**Figure 2A**) obtained from the tyrosinase activity also confirmed the uncompetitive inhibition of melanin onto crude tyrosinase. The values of Michelies constant and the maximum enzyme velocity changed at different inhibitor concentrations, as shown in **Figure 2A**.

 The above steps quantified the melanin inhibition onto the growth of the microorganisms and inhibition of its production. Therefore, it was necessary to remove melanin from the broth to reduce feedback inhibition and increase the productivity of the bioprocess. The present protocol describes the development of this adsorptive bioprocess. Different commercially viable and easily available adsorbents such as zeolite, celite, alumina and activated carbon were used to study their melanin adsorption potential from the broth. The method used to identify the potential adsorbent is explained in protocol sections 4 and 5. As Figure 2B indicates, alumina was the best adsorbent, which adsorbed almost 60% of the melanin produced from the broth at pH 8. Melanin is ionized in the anionic form in the basic pH of the broth and hence solubilizes in water. Further, it was necessary to determine the adsorption capacity of alumina for melanin. Adsorption isotherms describe the adsorption capacity of a given adsorbent. Melanin adsorption onto alumina was better explained by Langmuir isotherm, and therefore, this isotherm was used to find the binding capacity of alumina. Langmuir adsorption isotherm assumes monolayer adsorption, non-interaction between the adsorbed ions, the uniform energy of adsorption in the plane of the adsorbent surface and availability of the identical binding sites for monolayer adsorption. The linear form of Langmuir adsorption isotherm can be written as follows <sup>22</sup>:

$$\frac{C_e}{q_e} = \left(\frac{1}{q_{max}}\right) \cdot \left(\frac{1}{b}\right) + \frac{C_e}{q_{max}} \tag{18}$$

Where, Ce (mg/L) is the equilibrium concentration of melanin in the solution, qmax is the maximum adsorption capacity, which is the amount of adsorbate required to form a monolayer, qe (mg/g) is the amount of melanin adsorbed per specific amount of adsorbent at equilibrium and b is a constant, which represents the affinity of adsorbent to adsorbate.

**Figure 2C** shows experimental values fitted to the Langmuir adsorption isotherm equation (18). The maximum adsorption capacity of melanin onto alumina from pure melanin solution and the broth at pH 8 was 127.4 mg melanin/g of alumina and 113.8 mg melanin/g of alumina, respectively. The amount of proteins, sugars, and cells/DNA adsorbed onto alumina was negligible (>2 mg/g) and therefore is not reported.

**Figure 2D** shows the effect of dosage time of alumina on melanin production. Dosage time of 12 h demonstrated better melanin yield. At this stage, cells were in the exponential growth phase with a biomass concentration of 1.7 g/L. The pH of the medium was 7.3 and the melanin concentration was about  $130 \pm 15$  mg/L. There was about a 50% increase in melanin production when adsorbent was added to the culture after 12 h than the normal culture grown in the optimized medium. The effect of adsorbent concentration on melanin production was investigated in the next step of adsorptive bioprocess development.

Figure 3A shows that melanin production increased as the adsorbent dosage increased. Tyrosinase enzyme activity also improved by increasing the adsorbent dose, which is evident from the figure. At 10 g/L adsorbent concentration, the reducing sugars were completely consumed by P. stutzeri by 48 h. After identifying suitable alumina dosage time and amount, an adsorptive bioprocess was carried out in a 5-L stirred tank bioreactor with the addition of alumina (10 g/L) at the time when melanin concentration reached about 130 mg/L. The melanin production increased to about 1349 mg/L (Figure 3B), which is about an 8.8-fold increase as compared to the batch process in the unoptimized medium (Figure 3C). The solid lines in Figure **3B** represent simulation results of the adsorptive bioprocess. **Table 1** compares the growth and product formation parameters for a normal batch bioprocess and adsorptive bioprocess. The table is self-explanatory and proves that adsorptive bioprocess increased the specific growth rate, yield coefficients and hence productivity of melanin biosynthesis. Further, Figure 3C compares the melanin production at different time intervals in a normal batch bioprocess and an adsorptive bioprocess in a stirred tank bioreactor under optimized conditions. It also shows the amount of melanin remaining in the supernatant during the adsorptive bioprocess. The melanin concentration remained nearly constant in the supernatant, in the range of 95-110 mg/L in the presence of alumina, thus reducing the inhibitory effect. From the figure, it is evident that there is an enhancement in melanin production due to the addition of adsorbents in a batch process. Figure 3D demonstrates the adsorption of melanin onto alumina in the adsorptive bioprocess. The alumina adsorbed the melanin produced in the broth. The growth medium containing free melanin and melanin bound on alumina particles was collected at the end of the bioprocess. A 0.1 M NaOH solution was added to the alumina particles obtained to desorb melanin from alumina particles, which is exhibited in Figure 3D. Particle size analysis of the dispersed solution of biosynthesized melanin in water shows that particles are nanoparticles with a size of 32 nm ± 0.98 nm (**Figure 4A**). **Figure 4B** shows an aggregate of melanin molecules of size in the nano range. Thus, particle size analysis and transmission electron microscopy show that the melanin molecules are of nano size.

# FIGURE AND TABLE LEGENDS:

**Figure 1: Profiling of melanin with respect to growth.** (A) The growth of the cells. (B) The melanin production at different melanin concentrations. (C) Decrease in the enzyme activity with an increase in the inhibitor concentration. (D) Maximum specific growth rate at different melanin concentrations in the optimized medium.

**Figure 2: Adsorption profile of melanin**. **(A)** The Lineweaver-Burke plot for tyrosinase enzyme from Pseudomonas Stutzeri. **(B)** adsorption of melanin by different adsorbents and their zeta potential values. **(C)** Langmuir adsorption isotherm for adsorption from melanin solution. **(D)** Effect of time of adsorbent dosage on melanin production.

**Figure 3: Bioprocess sketch and melanin adsorption**. (A) Effect of adsorbent (alumina) on melanin production and enzyme activity of cells. (B) Melanin production in stirred tank bioreactor with the addition of adsorbent alumina. (— simulation results). (C) Production of melanin in normal batch bioprocess and adsorptive bioprocess. (D) Images of melanin adsorption onto alumina and its desorption

**Figure 4: Particle size and TEM analysis. (A)** Particle size analysis of the biosynthesized melanin. **(B)** An aggregate of melanin particles seen in TEM.

**Table 1: Growth parameters.** Values of growth parameters obtained through equations (8), (13) and (17) are the values of  $v_x$ ,  $\alpha$  and  $\mu_{max}$ , respectively.  $Y_{X/S}$  is calculated from grams of biomass produced per gram of carbon substrate consumed, and  $Y_{P/S}$  is calculated from grams of product formed per gram of carbon substrate consumed. The last column shows a percentage increase in values of growth parameters of adsorptive bioprocess with respect to the normal batch process.

# **DISCUSSION:**

Product inhibition is a major bottleneck in bioprocessing which leads to reduced productivity. Several methods exist to reduce product inhibition, such as continuous bioprocess and *in situ* product removal techniques. However, these options require a complete overhaul of the existing bioprocessing facility<sup>23–28</sup>. Adsorptive bioprocess does not require a significant overhaul in the existing facility and can be implemented very easily. The melanin exhibited feedback inhibition on tyrosinase as well as reduced the growth of the *Pseudomonas stutzeri* culture. The critical concentration of melanin to completely cease the growth was about 1050 mg/L. The extracellular melanin isolated from *Schizophyllum commune* showed significant antifungal activity against dermatophytic fungi, *Trichophyton simii*, and *Trichophyton rubrum* and antibacterial activity against *Escherichia coli*, *Proteus sp.*<sup>29</sup>. Melanin produced by *Auricularia auricula* displayed inhibitory activity on three bacterial strains, *Pseudomonas aeruginosa* PAO1, *Pseudomonas fluorescens* P-3 and *E. coli* K-12, and there was a significant reduction in biomass with the increase in pigment concentration<sup>30</sup>. Melanin produced by *Streptomyces lusitanus* DMZ-3 showed

cytotoxic activity on brine shrimps<sup>31</sup>. Although microbial melanin inhibition is reported in the literature, the finding that melanin also inhibits the source organism itself and exhibits feedback inhibition is not reported to date. Therefore, an adsorptive bioprocess was developed as given in the protocol to reduce/nullify melanin inhibition. Alumina was selected as the best adsorbent for melanin. Zeolite and fuller's earth are usually cation exchangers, due to which they did not bind large amounts of melanin. Celite and activated carbon are good physical adsorbents; however, they adsorb non-polar compounds.

563564565

566567

568

569

570

571

572

573

574

575

576

577

578

579

580

581

582

583

584

585

556

557

558

559

560

561

562

On the contrary, alumina is a polar adsorbent and hence was able to adsorb polar and anionic melanin from the fermentation broth. The maximum adsorption using alumina coincided with the maximum difference in zeta potential (Figure 2B), confirming the suitability of alumina as the adsorbent for melanin in basic pH. Therefore, alumina was chosen for the development of adsorptive bioprocess. The dosage time was optimized to 12 h after the start of the bioprocess (Figure 2D). The melanin production decreased when alumina was added just after the inoculation. Alumina is known to partially inhibit microbial growth. Since melanin was not biosynthesized in the medium just after the inoculation, free alumina might have inhibited the growth and hence production of melanin. Alumina addition at 36 hours after inoculation also did not show improvement in melanin production. By this time, melanin was already produced in quantities, which would have initiated inhibitory effects. These findings indicate that time of addition of adsorbent is crucial for successful adsorptive bioprocess, especially if the adsorbent is inhibiting the microbial growth. The adsorptive bioprocess using alumina as an adsorbent reduced melanin inhibition. The adsorptive bioprocess in a 5 L bioreactor increased melanin production from 520 mg/L to 1349 mg/L, using the same optimized medium. Thus, adsorptive bioprocess further increases the productivity by 2.6 times and 8.8 times as compared to an optimized batch process and an unoptimized process, respectively. Yield coefficient of melanin over substrate increased from 8.23 mg/g to 22.5 mg/g by increasing adsorbent dosage. Moreover, the total sugars were not detectable in the culture broth at the end of the bioprocess at an adsorbent dosage of 10 g/L, thus reducing the loading on the effluent treatment plant. The process described here can be easily incorporated to reduce/nullify product inhibition in any bioproduction facility without further modification in the existing plant, thus providing operational flexibility.

586 587 588

589

590

591

# **ACKNOWLEDGMENTS:**

We thank the Department of Science and Technology (DST/TSG/WP/2014/58), India, and the National Institute of Technology Karnataka for providing funding for the development of the above protocol.

592 593

# **DISCLOSURES:**

The authors do not have any conflict of interest

595596

597

598

599


594


# **REFERENCES:**


1. Langfelder, K., Streibel, M., Jahn, B., Haase, G., Brakhage, A. A. Biosynthesis of fungal melanins and their importance for human pathogenic fungi. *Fungal Genetics and Biology*. **38** (2), 143–158 (2003).


- 600 2. Riley, P. A. Melanin. *International Journal of Biochemistry and Cell Biology*. **29** (11), 1235–601 1239 (1997).
- Madaras, F., Gerber, J. P., Peddie, F., Kokkinn, M. J. The effect of sampling methods on
- the apparent constituents of ink from the squid sepioteuthis australis. *Journal of Chemical*
- 604 *Ecology*. **36** (11), 1171–1179 (2010).
- 4. Tran-Ly, A. N. *et al.* Microbial production of melanin and its various applications. *World Journal of Microbiology & Biotechnology.* **36** (11), 170 (2020).
- 607 5. Pavan, M. E., López, N. I., Pettinari, M. J. Melanin biosynthesis in bacteria, regulation and
- 608 production perspectives. *Applied Microbiology and Biotechnology*. 104 (4), 1357–1370 (2019).
- 609 6. Pombeiro-Sponchiado, S. R., Sousa, G. S., Andrade, J. C. R, Lisboa, H. F., Gonçalves, R. C.
- R. Production of melanin pigment by fungi and its biotechnological applications. *Melanin*. 47–77,
- 611 IntechOpen Limited, London, UK, (2017).
- 612 7. Guo, J. et al. Cloning and identification of a novel tyrosinase and its overexpression in
- Streptomyces kathirae SC-1 for enhancing melanin production. FEMS Microbiology Letters. 362,
- 614 41 (2015).
- 615 8. Lagunas-Muñoz, V. H., Cabrera-Valladares, N., Bolívar, F., Gosset, G., Martínez, A.
- 616 Optimum melanin production using recombinant Escherichia coli. Journal of Applied
- 617 *Microbiology*. **101** (5), 1002–1008 (2006).
- 618 9. Wang, L., Li, Y., Li, Y. Metal ions driven production, characterization and bioactivity of
- 619 extracellular melanin from Streptomyces sp. ZL-24. International Journal of Biological
- 620 *Macromolecules*. **123**, 521–530 (2019).
- 621 10. Kiran, G. S., Jackson, S. A., Priyadharsini, S., Dobson, A. D. W., Selvin, J. Synthesis of Nm-
- 622 PHB (nanomelanin-polyhydroxy butyrate) nanocomposite film and its protective effect against
- biofilm-forming multi drug resistant Staphylococcus aureus. *Scientific Reports.* **7**, 9167 (2017).
- 624 11. Raman, N. M., Shah, P. H., Mohan, M., Ramasamy, S. Improved production of melanin
- from Aspergillus fumigatus AFGRD105 by optimization of media factors. *AMB Express.* **5** (1), 72
- 626 (2015).
- 12. Thaira, H., Raval, K., Manirethan, V., Balakrishnan, R. M. Melanin nano-pigments for heavy
- metal remediation from water. Separation Science and Technology (Philadelphia). 54 (2), 265–
- 629 274 (2018).
- Han, K., Levenspiel, O. Extended monod kinetics for substrate, product, and cell inhibition.
- 631 *Biotechnology and Bioengineering.* **32** (4), 430–447 (1988).
- 632 14. Vázquez, J. A., Murado, M. A. Unstructured mathematical model for biomass, lactic acid
- and bacteriocin production by lactic acid bacteria in batchfermentation. Journal of Chemical
- 634 *Technology and Biotechnology.* **83** (1), 91–96 (2008).
- 635 15. Rohit, S. G., Jyoti, P. K., Subbi, R. R. T., Naresh, M., Senthilkumar, S. Kinetic modeling of
- 636 hyaluronic acid production in palmyra palm (Borassus flabellifer) based medium by Streptococcus
- cooepidemicus MTCC 3523. Biochemical Engineering Journal. 137, 284–293 (2018).
- 638 16. Nair, A. V, Gummadi, S. N., Doble, M. Process optimization and kinetic modelling of cyclic
- 639 (1 $\rightarrow$ 3, 1 $\rightarrow$ 6)-β-glucans production from Bradyrhizobium japonicum MTCC120. *Journal of*
- 640 *Biotechnology*. **226**, 35–43 (2016).
- 641 17. Reddy Tadi, S. R., E. V. R., Limaye, A. M., Sivaprakasam, S. Enhanced production of
- optically pure D (–) lactic acid from nutritionally rich Borassus flabellifer sugar and whey protein
- 643 hydrolysate based–fermentation medium. Biotechnology and Applied Biochemistry. 64 (2), 279–


- 644 289 (2017).
- 645 18. Ganesh Kumar, C., Sahu, N., Narender Reddy, G., Prasad, R. B. N., Nagesh, N., Kamal, A.
- Production of melanin pigment from Pseudomonas stutzeri isolated from red seaweed Hypnea
- musciformis. Letters in Applied Microbiology. **57** (4), 295–302 (2013).
- 648 19. Ren, Q., Henes, B., Fairhead, M., Thöny-Meyer, L. High level production of tyrosinase in
- recombinant Escherichia coli. *BMC Biotechnology*. **13**, 13–18 (2013).
- 650 20. Kruger, N. J. The Bradford method for protein quantitation. The Protein Protocols
- 651 Handbook. Springer Protocols Handbooks. Humana Press (2002).
- 652 21. Miller, G. L. Use of Dinitrosalicylic acid reagent for determination of reducing sugar.
- 653 *Analytical Chemistry.* **31** (3), 426–428 (1959).
- Dutta, A., Diao, Y., Jain, R., Rene, E. R., Dutta, S. Adsorption of cadmium from aqueous
- 655 solution onto coffee grounds and wheat straw: Equilibrium and kinetic study. Journal of
- 656 Environmental Engineering. **142** (9), 1–6 (2015).
- 657 23. Dhillon, G. S., Brar, S. K., Verma, M., Tyagi, R. D. Recent advances in citric acid bio-
- 658 production and recovery. *Food and Bioprocess Technology*. **4**, 505–529 (2011).
- 659 24. Roffler, S. R., Blanch, H. W., Wilke, C. R. In situ recovery of fermentation products. *Trends*
- 660 in Biotechnology. **2** (5), 129–136 (1984).
- 661 25. Maiti, S. et al. A re-look at the biochemical strategies to enhance butanol production.
- 662 Biomass and Bioenergy. **94**, 187–200 (2016).
- 663 26. Prakash, G., Srivastava, A. K. Integrated yield and productivity enhancement strategy for
- 664 biotechnological production of Azadirachtin by suspension culture of Azadirachta indica. Asia-
- 665 *Pacific Journal of Chemical Engineering*. **6** (1), 129–137 (2011).
- 666 27. Stark, D., von Stockar, U. In situ product removal (ISPR) in whole cell biotechnology during
- the last twenty years. Advances in Biochemical Engineering/Biotechnology. **80**, 149–175 (2003).
- 668 28. Kaur, G., Srivastava, A. K., Chand, S. Debottlenecking product inhibition in 1,3-propanediol
- 669 fermentation by In-Situ Product Recovery. *Bioresource Technology*. **197**, 451–457 (2015).
- 670 29. Arun, G., Eyini, M., Gunasekaran, P. Characterization and biological activities of
- 671 extracellular melanin produced by Schizophyllum commune (Fries). Indian Journal of
- 672 Experimental Biology. **53**, 380–387 (2015).
- 673 30. Bin, L., Wei, L., Xiaohong, C., Mei, J., Mingsheng, D. In vitro antibiofilm activity of the
- 674 melanin from Auricularia auricula, an edible jelly mushroom. Annals of Microbiology. 62 (4),
- 675 1523–1530 (2012).
- 676 31. Madhusudhan, D. N., Mazhari, B. B. Z., Dastager, S. G., Agsar, D. Production and
- 677 cytotoxicity of extracellular insoluble and droplets of soluble melanin by Streptomyces lusitanus
- 678 DMZ-3c. *BioMed Research International*. **2014**, 1–11 (2014).











Parameter	Normal batch process	Adsorptive bioprocess	% increase wrt batch process
$\mu_{max}$ (1/h)	0.145	0.25	72.4
$v_x$ (g/l/h)	0.15	0.28	86.6
$\alpha$ (g/g)	0.081	0.105	29.6
$Y_{X/S}$ (g/g)	0.14	0.21	50
$Y_{P/S}$ (mg/g)	8.3	22.5	177

Table of Materials

Click here to access/download **Table of Materials**Table of Materials-63339R3.xls

#### Dear Editors,

Thank you very much for reviewing the manuscript. We changed the manuscript as per the queries raised by you. Following are our replies to the queries raised.

Kindly note that all the changes have been made and are in the "TRACK CHANGES" mode in the word file

**Query 1:** Please note that the section numbers and step numbers have been formatted to fit the Journal standard. Please review.

**Reply:** As suggested by the editors, the manuscript was reviewed for the formatting changes.

**Query 2:** Howe is this done? Please add more details regarding retrieving and maintenance of cultures. **Reply:** As suggested by the editors, relevant details have been added for the retrieval and maintenance of the cultures.

**Query 3:** How is subculturing done. Please add the steps for subculturing. A citation would suffice. **Reply:** As suggested by the editors, the methodology for the subculturing has been added.

**Query 4:** Artificial seawater? Is this prepared? Any specific composition? **Reply:** The artificial seawater broth was purchased from Himedia, Bangalore, India. The specifications have been added in the text.

**Query 5:** If the details for each parameter mentioned in this step is detailed below, please reference the section numbers of the specific methods in this step?

**Reply:** The details of the section has been added in the text.

**Query 6:** The sentence s incomplete. Please revise.

Reply: The suggested modifications have been incorporated.

Query 7: Revised. Please check.

**Reply:** The corrections were made in the revised sentence.

**Query 8:** How?

Reply: The reference used for melanin estimation has been incorporated in the test.

**Query 9:** How is the bioprocess harvested?

**Reply:** The sentence has been modified to mention the harvesting step.

**Query 10:** How is the absorbent washed. Is any centrifugation or vortexing involved? Any specific volume of water and ethanol used? Any specific time duration for each washing step?

Reply: The centrifugation details required for washing have been incorporated in the text.

**Query 11:** What is added to reduce the pH? Please specify the added constituent and the volume added to achieve the mentioned pH.

**Reply:** The details of the HCL used for the acidification have been added. The volume of the HCl will depend on the melanin yield, hence the endpoint pH has been mentioned.

Query 12: Any centrifugation or voretxing involved? What is the time duration of each wash?

**Reply:** The details of the centrifugation have been incorporated in the manuscript.

**Query 13:** Is this stimulate or simulate?

**Reply:** The conditions were simulated for computer modeling.

Query 14: Please cite the reference number

**Reply:** The citation has been added in the text.

**Query 15:** How is the culture collected? How is the sterility ensured? What is the volume collected each time to measure the OD?

**Reply:** The details of the query raised has been addressed and incorporated in the manuscript.

**Query 16:** After how long of culturing is the culture centrifuged?

**Reply:** The samples were collected at an interval of 4h. This is mentioned in the preceding step.

Query 17: Sample? Please specify.

**Reply:** It is a commercial DNA sample and has been mentioned in the manuscript.

**Query 18:** What is the dilution?

**Reply:** The dilution factor has been included in the text.

**Query 19:** Is this commercial or prepared in the lab? If its prepared, please include a citation. In case it is purchased, ensure that it is included in the Table of Materials.

**Reply:** The estimation was done using the commercial kit from Himedia and has been mentioned in the manuscript.

Query 20: Volume?

**Reply:** The volume utilized has been added in the text.

Query 21: What is the volume of water used? How much of melanin is dispersed?

**Reply:** The volume has been included in the text.

**Query 22:** Volume?

**Reply:** The volume has been included in the text.

Query 23: Is this performed according to the manufacturer's protocol? Please cite.

**Reply:** The changes have been incorporated as suggested.

Query 24: How? Using an oven or dry at RT? For how long is the melanin dried?

**Reply:** The melanin is kept in the oven for 8h for drying. The conditions have been included.

**Query 25:** How is the sample loaded to the TEM sample matrix

**Reply:** The protocol as specified by the manufacturer was followed. This has been mentioned.

**Query 26:** Please check if the section number is correct.

**Reply:**The section no. as suggested has been corrected.

Query 27: Is this correct?

**Reply:** The notation of the units have been changed from 1/h to h<sup>-1</sup>

Query 27: Please check if the section numbers are correct.

**Reply:** The section number mentioned are correct.

**Query 28:** Please include a title for each figure.

**Reply:** The title has been included as suggested.

Query 29: Please include a brief description for the Table

**Reply:** The changes as suggested have been incorporated in the main text.

Query 30: As we are a methods journal, please ensure that the Discussion explicitly covers the following in detail in 3-6 paragraphs with citations: a) Critical steps within the protocol, b) Any

modifications and troubleshooting of the technique, c) Any limitations of the technique, d) The significance with respect to existing methods, e) Any future applications of the technique. **Reply:** The text explicitly speaks about the suggested points.