

Screenshot Summary:

- 67129_screenshot_1.mp4

- 2.12.1. SCREEN: Resonator tuning process with tuning wheel adjustments on-screen (00:04-00:57)
- 2.13.1. SCREEN: Microwave power adjustment and attenuation level to show peak at 60 ns (01:43-03:13).
- 2.13.2. SCREEN: FID centering and phasing display in progress (00:58-01:42).
- 2.15.1. SCREEN: Programmed sequence showing T1, T2 (03:13-06:14).
- 2.15.2. SCREEN: Cursor hovering 3D Electron Spin Echo data (19:19-24:03).
- 2.15.3. SCREEN: 4D Inversion Recovery Electron Spin Echo (IRESE) acquisition in progress (07:04-19:05).

- 67129_screenshot_2.mp4

- 2.16.1. SCREEN: Display of pO2 map with tumor mask visualization (00:04-00:37; 01:18-01:46).
- 2.16.2. SCREEN: Exporting voxel data (00:38;01:10).