

Journal of Visualized Experiments

Light-Induced Molecular Adsorption of Proteins using the PRIMO system for micro-patterning to study cell responses to extracellular matrix proteins --Manuscript Draft--

Article Type:	Invited Methods Article - JoVE Produced Video
Manuscript Number:	JoVE60092R4
Full Title:	Light-Induced Molecular Adsorption of Proteins using the PRIMO system for micro-patterning to study cell responses to extracellular matrix proteins
Keywords:	Micro-pattern, patterning, micro-patterning, LIMAP, extracellular matrix (ECM), neuron, cell behavior, migration, axonal growth, pathfinding, decision-making
Corresponding Author:	Christoph Ballestrem The University of Manchester Faculty of Biology Medicine and Health Manchester, UNITED KINGDOM
Corresponding Author's Institution:	The University of Manchester Faculty of Biology Medicine and Health
Corresponding Author E-Mail:	Christoph.Ballestrem@manchester.ac.uk
Order of Authors:	Cristina Melero Aljona Kolmogorova Paul Atherton Brian Derby Adam Reid Karin Jansen Christoph Ballestrem
Additional Information:	
Question	Response
Please indicate whether this article will be Standard Access or Open Access.	Open Access (US\$4,200)
Please indicate the city, state/province, and country where this article will be filmed . Please do not use abbreviations.	Manchester, Lancashire, UK (M13 9PT)

TITLE:

Light-Induced Molecular Adsorption of Proteins using the PRIMO System for Micro-Patterning to Study Cell Responses to Extracellular Matrix Proteins

AUTHORS AND AFFILIATIONS:

Cristina Melero^{*1}, Aljona Kolmogorova^{*1}, Paul Atherton¹, Brian Derby², Adam Reid^{3,4}, Karin Jansen⁵, Christoph Ballestrem¹

¹Faculty of Biology, Medicine and Health. Division of Cell Matrix, Biology and Regenerative Medicine, Wellcome Trust Centre for Cell-Matrix Research, The University of Manchester, Manchester, UK

²School of Materials, The University of Manchester, Manchester, UK

³Blond McIndoe Laboratories, School of Biological Sciences, Faculty of Biology, Medicine and Health, Division of Cell Matrix Biology and Regenerative Medicine, The University of Manchester, Manchester Academic Health Science Centre, Manchester, UK

⁴Department of Plastic Surgery & Burns, Wythenshawe Hospital, Manchester University NHS Foundation Trust, Manchester Academic Health Science Centre, Manchester, UK.

⁵UMC Utrecht, Department of Pathology, Utrecht, The Netherlands.

*These authors contributed equally

Email Address:

Cristina Melero: cristina.melero@manchester.ac.uk

Aljona Kolmogorova: aljona.kolmogorova@postgrad.manchester.ac.uk

Paul Atherton: paul.atherton@manchester.ac.uk

Brian Derby: brian.derby@manchester.ac.uk

Adam Reid: adam.reid@manchester.ac.uk

Corresponding Authors:

Christoph Ballestrem: christoph.ballestrem@manchester.ac.uk

Karin Jansen: kjansen6@umcutrecht.nl

KEYWORDS:

micro-pattern, patterning, micro-patterning, LIMAP, extracellular matrix (ECM), neuron, cell behavior, migration, axonal growth, pathfinding, decision-making

SUMMARY:

Our overall aim is to understand how cells sense extracellular cues that lead to directed axonal growth. Here, we describe the methodology of Light-Induced Molecular Adsorption of Proteins, used to produce defined micro-patterns of extracellular matrix components in order to study specific events that govern axon outgrowth and pathfinding.

ABSTRACT:

Cells sense a variety of extracellular cues, including the composition and geometry of the extracellular matrix, which is synthesized and remodeled by the cells themselves. Here, we present the method of Light-Induced Molecular Adsorption of Proteins (LIMAP) using the PRIMO

system as a patterning technique to produce micro-patterned extracellular matrix (ECM) substrates using a single or combination of proteins. The method enables printing of ECM patterns in micron resolution with excellent reproducibility. We provide a step-by-step protocol and demonstrate how this can be applied to study the processes of neuronal pathfinding. LIMAP has significant advantages over existing micro-printing methods in terms of the ease of patterning more than one component and the ability to generate a pattern with any geometry or gradient. The protocol can easily be adapted to study the contribution of almost any chemical component towards cell fate and cell behavior. Finally, we discuss common issues that can arise and how these can be avoided.

INTRODUCTION:

In recent years, the biological sciences have increasingly made use of the advances provided by material sciences. One prominent example is the micro-patterning of substrates, which can be used to study cellular responses such as cell proliferation^{1,2}, differentiation³⁻⁶, cell migration⁷⁻⁹ and pathfinding^{10,11}. There are a number of techniques available that enable the micro-patterning of substrates, such as multiphoton excited photochemistry¹², AFM dip-pen nanolithography¹³, pin and inkjet direct printing¹⁴, electron beam lithography¹⁵ or microfluidics¹⁶. However, two techniques that are widely used in the biological field are microcontact printing¹⁷⁻¹⁹ or laser-assisted patterning³ (**Figure 1**). Laser-assisted patterning is considered to deliver more reliable results in terms of protein and PEG stability and cell confinement on the patterns, compared to microcontact printing²⁰. A more novel approach for micro-patterning described here is the use of Light-Induced Molecular Adsorption of Proteins²¹ (LIMAP, **Figure 1D**) using a commercially available system (PRIMO, **Table of Materials**). Each of the methods has advantages and limitations that are briefly described below.

Microcontact printing uses PDMS molds (stamps) with desired micro-features that are generated from lithographed masters. The stamps are incubated with a chosen protein which is then transferred (stamped) onto the cell culture substrate¹⁸ (**Figure 1A**). Laser-assisted patterning uses UV light to cleave an anti-fouling film²²⁻²⁵, exposing regions that can subsequently be coated with the protein of interest (**Figure 1B**). While the resolution achieved with photo-patterning approaches is in the micron range^{25,26}, most of these techniques require a photo-mask, either in contact with the sample, or situated in the object plane of the microscope objective^{23,27,28}. The requirements for masks in both microcontact printing and photo-patterning can be a limitation; specific masks are required for every geometric pattern and size, which can be expensive and time consuming to generate. In contrast to these techniques, LIMAP does not require a mask (**Figure 1D**). Using the PRIMO system for LIMAP can be cost intensive at the beginning because it requires purchasing of equipment. However, open-source software is used to design patterns of any desired geometry, giving much more freedom and allowing more complex experiments including the use of protein concentration gradients.

The PRIMO laser is controlled and directed by a digitally-controlled micromirror device (DMD) to create patterns in any number of user-defined geometries. LIMAP requires the culture surface to be coated with molecules that prevent cell attachment. Polyethylene glycol (PEG) is most commonly used as such an “antifouling” reagent; it forms a dense anti-adhesive film on the glass

or plastic surface²⁹. Subsequently, a photo-initiator is added that allows the PEG film to be removed with high precision through a photoscission mechanism³⁰ by local exposure to UV light under the control of the DMD. These PEG-free regions can be coated with proteins that adsorb to the laser-etched surface, generating a micro-pattern. By varying the laser power, different quantities of PEG can be removed from the surface allowing the user to generate protein gradients. PEG removal and the coating procedure can be repeated to create patterns with two or more distinct proteins in the same micro-well²¹.

The generated micro-patterns provide adhesive surfaces for cells, allowing the study of cell behavior. In our studies, we use micro-patterning to study neurite or axon pathfinding of a neuronal cell line (CAD (Catecholaminergic-a differentiated) cells³¹) or primary rat dorsal-root ganglion (DRG) neurons, respectively.

Here, we outline a step-by-step protocol for LIMAP (**Figure 2**) using the commercially available PRIMO system and accompanying Leonardo software. We demonstrate how it can be used for the generation of patterns with defined geometries and multiple proteins, which we use to study axonal pathfinding. We discuss common issues that can arise and how these can be avoided.

PROTOCOL:

1. Design of pattern templates

NOTE: Templates for patterning are generated with digital drawing software (**Table of Materials**). Drawing in different grey levels will determine laser intensities. Using software to design pattern templates allows rapid generation of patterns with any desired geometry and gradients (**Figure 3**).

1.1 Digitally draw the desired pattern template using drawing software. Select an image size of 1824 pixels length and 1140 pixels width (which in this study corresponds to 415 μ m length and 260 μ m width). Save the pattern template as an 8-bit Tiff file.

NOTE: a step-by-step protocol to generate templates is available upon request to the brand who commercializes the micro-patterning equipment (**Table of Materials**).

2. Plasma cleaning

NOTE: Optimal results require plasma cleaning of the surfaces prior to patterning, which will remove all organic matter and activate the surface. In the present case, ambient air is sufficient for surface activation. A plasma cleaner (**Table of Materials**) was used with a process pressure of 1000-1300 mTorr and a power of 29.6 W for 1-5 min.

2.1 Use glass bottom dish/es with a 20 mm inner well size and glass thickness of 0.16-0.19 mm. For testing multiple conditions, use a 6-well glass bottom dish. Otherwise, use a single well glass bottom dish (**Table of Materials**).

NOTE: A step-by-step protocol for plasma cleaning is available upon request to the brand who commercializes the plasma cleaning equipment (**Table of Materials**).

3. Passivation

NOTE: This step generates an antifouling film that prevents protein adsorption to the glass surface. PEG offers high resistance to protein adsorption²⁹ as an antifouling agent. LIMAP uses a photo-initiator to locally remove PEG through UV photocission. The protein/s of interest will then adsorb to these PEG-free surfaces²¹, generating micro-patterns.

3.1 Passivation with PLL-PEG

3.1.1 Under sterile conditions, cut the PDMS stencils (see **Figure 4A,B** and **Table of Materials**) to ensure they fit in the inner glass bottom well, and remove the inner micro-well fillings with sterile forceps. Stick stencils on the glass well using forceps.

3.1.2 Ensure that the stencils stick tightly on the glass well, preventing the formation of air bubbles that can cause leakages during the passivation process.

NOTE: PDMS stencils can also be fabricated in-house using published protocols^{18,32}.

3.1.3 Prepare PLL-PEG solution (**Table of Materials**, 0.1 mg/mL) in phosphate-buffered saline (PBS). Add 20 μ L of PLL-PEG solution to each micro-well and incubate at room temperature, for 1 h.

3.1.4 Remove 15 μ L of PLL-PEG from the micro-wells and wash them five times with 20 μ L of PBS (**Table of Materials**) without letting them dry.

NOTE: Always leave approximately 5 μ L of PBS in between washes. Given their small volume, the micro-wells are particularly susceptible to drying. Drying will result in micro-patterns of poor quality.

3.1.5 Either keep the culture dish in PBS (3 mL per well) at 4 °C for up to 3 days or continue with the next step (step 3.1.5).

3.1.6 Remove 18 μ L of PBS from a single micro-well (for example, the top left micro-well, see **Figure 4D**), add 5 μ L of photo-initiator (PLPP, **Table of Materials**) and leave 20 μ L of PBS in the remaining micro-wells. This micro-well with PLPP will be used to create a reference pattern (see **Figure 4D,E**) during the system calibration step (step 4). Keep PLPP protected from light.

3.2 Passivation with long-lasting PEG-SVA

NOTE: To generate the anti-adhesive surface one can use: 1) a PEG linked to Poly-L-Lysine (PLL-PEG, step 3.1) or 2) a PEG-succinate N-hydroxysuccinimide (PEG-SVA). The decision to choose one or the other passivation option depends on the storage options (see step 10). Culture dishes incubated with PEG-SVA require double amount of passivation and photopatterning time.

3.2.1 Prepare stencils as described in step 3.1.1.

3.2.2 Add 20 μ L of 0.01% Poly-L-Lysine (PLL, **Table of Materials**) to each micro-well and incubate at room temperature for 30 min to pre-coat with PLL.

3.2.3 Remove 15 μ L of PLL from the micro-wells and wash them three times with 20 μ L of 1 M HEPES buffer (**Table of Materials**) without letting the wells dry out.

NOTE: Always leave approximately 5 μ L of HEPES or PBS in between washes. Given their small volume, the micro-wells are particularly susceptible to drying. Drying will result in micro-patterns of poor quality.

3.2.4 Prepare PEG-SVA solution. PEG-SVA solution (50 mg/mL in HEPES buffer 1M) should be prepared fresh each time immediately before use. Prepare 20 μ L of PEG-SVA per micro-well.

3.2.5 HEPES buffer must have a pH between 8-8.5. Test HEPES pH prior to preparing the PEG-SVA solution with pH paper. Weigh PEG-SVA in a centrifuge tube using a precision scale. Add HEPES buffer and vortex 30 s until dissolved. PEG-SVA is completely dissolved when the solution is transparent.

NOTE: SVA is the ester that allows binding of PEG to the previously coated PLL. Once the HEPES buffer is added to the PEG-SVA, the SVA has a half-life of 15 min and should be used immediately.

3.2.6 Add 20 μ L of PEG-SVA solution to each micro-well and incubate at room temperature, for 1 h. Remove 15 μ L of PEG-SVA from the micro-wells and wash five times with 20 μ L of PBS (**Table of Materials**) without letting the wells dry out.

3.2.7 Either prepare the culture dish for long-term storage (up to 1 month, see step 10.2) or proceed to the next step (step 3.2.8).

3.2.8 Remove 18 μ L of PBS from a single micro-well (for example, the top left micro-well, see **Figure 4D**), add 5 μ L of PLPP (**Table of Materials**) and leave 20 μ L of PBS in the remaining micro-wells. This micro-well will be used to create a reference pattern (see **Figure 4D, E**). Ensure that PLPP is homogenous on the whole surface of the micro-well. Keep PLPP protected from light.

4. System calibration

NOTE: In these steps, the focus of the laser will be adjusted to the particular type of culture dish (step 4.1). A reference pattern will be generated in only one micro-well (step 4.2) followed by

incubation with a protein solution (step 4.3) to ensure the optimal focus conditions of the laser (step 4.4), necessary to obtain sharp and defined patterns.

4.1 Laser calibration

NOTE: During the calibration process, a calibration laser image will be projected onto a fluorescently-highlighted glass surface (calibration well, marked with a fluorescent highlighter, **Figure 4C**), which needs to be placed in focus on the microscope.

4.1.1 Use a fluorescent highlighter (**Table of Materials**) to mark an empty inner glass well.

NOTE: The calibration well must have the same glass thickness (0.16-0.19 mm) as the culture dish on which micro-patterns will be generated. If a 6-well glass bottom dish is being used, an empty well can be used for calibration and should be marked with fluorescent highlighter under sterile conditions.

4.1.2 Switch on the microscope, the stage and the computer. Switch on the PRIMO micro-patterning equipment, open Micro-manager and Leonardo software. Leonardo software is operated through Micro-manager under **Plugins**. Check the brands/catalogue numbers of equipment and software in **Table of Materials**

4.1.3 In the initial menu of Leonardo, select **Calibrate**. Check that the dedicated PRIMO filter cube is in the correct position (optical path) in the filter turret; select the 20X objective (0.75 DIC S Plan Fluor, no phase ring) both on the microscope and on Leonardo software.

NOTE: This protocol is adjusted to Leonardo version 4.4. The protocol may need adjusting for other versions.

4.1.4 Position the previously highlighted calibration well (**Figure 4C**) above the objective. Select camera path. Adjust the objective focus until the laser projection of both the **PRIMO** logo and the tag line **Take care of your cells** are in focus.

4.1.5 Leave default camera exposure time at 25 ms. Adjust the laser intensity to see the letter **I** from the **PRIMO** logo projection in grey color and the rest of the letters in white.

4.1.6 Record the Z-position of the focal plane, (height of the objective to the sample) later referred to as Z-position of calibration. This will be an approximation to the optimal focus obtained after reference pattern generation (see step 4.2-4.4).

4.2 Reference pattern

4.2.1 Position the culture dish with the micro-well containing the photo-initiator (reference pattern micro-well, **Figure 4D**) above the objective and select **Pattern now** in the Leonardo software.

4.2.2 Visualize the edge of the micro-well with transmitted light through the camera and choose the ROI symbol from the right-hand menu. Set the diameter of the ROI circle to 4000 μm and align the edge of the digital ROI with the edge of the current micro-well.

4.2.3 Ensure that there is an accurate overlap between the digital ROI and the current micro-well by moving the stage around the edges of the micro-well. The ROI position will be coupled to the stage movements.

4.2.4 Select **Lock** in Leonardo software to lock the ROI in the desired position. Switch off the transmitted light.

4.2.5 Select **PRIMO** to upload the desired pattern template, which will be projected on the ROI as a design unit (see **Figure 5**). Patterns will appear in a dropdown list referred to as **Actions** and shown in the **Actions menu** on the software.

NOTE: Pattern templates need to be designed previously (see step 1) and saved as an 8-bit Tiff file before loading the template in the software.

4.2.6 Only a small pattern is needed for the reference pattern; for example, 3 lines, 1 column (see **Figure 4E** and **Figure 5B**). In the **Replication** menu, set the desired number of columns and rows (**Lines** in Leonardo software). Click **Refresh** to observe a digital preview of the pattern design.

4.2.7 Set the laser dose in the **Replication** menu. The optimal laser dose in this setting and using PLL-PEG is 1390 mJ/mm^2 .

NOTE: Laser power can vary between 5-7.5 mW/mm^2 . In this case, it is 7.5 mW/mm^2 , which takes approximately 30 s to pattern each design unit, using a laser dose of 1390 mJ/mm^2 . Higher laser doses may be required if the culture dish surface is passivated with PEG-SVA (approximately double laser dose compared to PLL-PEG). This needs to be tested beforehand.

4.2.8 Navigate to a peripheral region of the micro-well, (for example, top part) away from the main region of interest for pattern generation (central region) of the micro-well (see **Figure 4E**) and select **Lock**. Wait until the pattern is displayed.

4.2.9 Adjust the focus to the Z-position of calibration (see step 4.1.6).

NOTE: It is recommended to do an additional system calibration step if another user is using the same microscope between photopatterning rounds.

4.2.10 Select the **Play** symbol to start patterning. Ensure the laser is on in the software. Wait until the patterning process is finished. Patterning duration will be displayed in the **Estimated Time** panel. On Leonardo software version 4.4 patterning is completed when all **Actions** appear

blue in the **Visualization** menu.

4.3 Protein incubation on reference pattern

4.3.1 Under sterile conditions, wash the reference pattern micro-well three times with 20 μ L of PBS to remove the PLPP.

4.3.2 Add 20 μ L of fluorescently-labeled ECM protein (10 μ g/mL laminin, fibronectin or fibrinogen in PBS, see **Table of Materials** and step 6) to the reference pattern micro-well. Incubate at room temperature for 10-20 min (depending on the protein) protected from light (wrap the dish in aluminum foil).

4.3.3 After incubation, remove 18 μ L of protein solution and wash three times with 20 μ L of PBS; keep the micro-wells that are not used to create a reference pattern (see **Figure 4D,E**) in 20 μ L of PBS.

4.4 Visualization and setting optimal laser focus

4.4.1 Visualize the reference pattern using an epifluorescence microscope, 20X objective and suitable software (check software used in this study in **Table of Materials**). The reference pattern should be visible in the peripheral region (for example, top well region), in which the reference pattern was generated (see **Figure 4E**).

4.4.2 Adjust focus on the pattern edges through the camera path. Record the Z-position according to the best focus on the reference pattern. This adjusted Z-position will be the optimal laser focus used for subsequent patterning.

5. Software set-up and photopatterning

NOTE: Once the calibration of the system has been achieved (step 4), the user will upload the desired pattern templates (template configuration, **Figure 5**) for photopatterning, with the option to generate patterns for one or multiple proteins in each micro-well. The micro-patterning process involves photopatterning and protein incubation steps (see **Figure 2**).

5.1 Under sterile conditions, remove 18 μ L of PBS from all the micro-wells and add 5 μ L of PLPP to each micro-well. Ensure that PLPP is homogenous on the whole surface of the micro-wells.

5.2 Position the culture dish with the reference pattern micro-well (**Figure 4D,E**) above the objective and select **Pattern now** in the Leonardo software.

5.3 Visualize the micro-well with transmitted light through the camera path and choose the ROI symbol. Set the diameter of the ROI to 4,000 μ m and overlap the edge of the digital ROI to the edge of the current micro-well. Select **Lock**.

NOTE: the shape and diameter of the ROI depends on the design and size of the PDMS stencil being used. For example, if using 5,000 x 5,000 μm PDMS squared stencils, use a 5,000 x 5,000 μm squared ROI.

5.4 Repeat the overlapping step (step 5.3) for each micro-well of the dish. Upon completion, switch off the transmitted light.

5.5 Navigate to the center of the reference pattern micro-well, away from the reference pattern region, and select **PRIMO** to upload the desired pattern template, which will be projected on the ROI as a design unit (see **Figure 5**). Patterns will appear in a dropdown list referred to as **Actions** and shown in the **Actions menu** on the software.

NOTE: Pattern templates need to be designed ahead of the experiment and saved as an 8-bit Tiff file before loading the template in the software.

5.6 In order to create a pattern across the whole micro-well, the design unit needs to be replicated. A design unit covers around 0.1 mm^2 of the micro-well area. In the **Replication** menu, set the desired number of columns and rows (**Lines** in the software) (see **Figure 5**).

5.7 To generate a continuous pattern, adjust the spacing between columns and lines; in this study, patterns of continuous stripes are obtained using -20 to -35 μm spacing (negative spacing) between columns. This negative spacing creates an overlap among design units (**Figure 5B,C**).

5.8 Set the laser dose in the **Replication** menu. The optimal laser dose in this setting and using PLL-PEG is 1390 mJ/mm^2 . The patterning duration will be displayed in the **Estimated Time** panel.

NOTE: In this case laser power is 7.5 mW/mm^2 , taking approximately 30 s to pattern each design unit, using a dose of 1390 mJ/mm^2 . For example, a design unit (0.1 mm^2) replicated in 4 columns and 4 lines (around 1.6 mm^2), would take 8 min to be patterned. Higher laser doses may be required if the culture dish surface is passivated with PEG-SVA (approximately double laser dose compared to PLL-PEG).

5.9 Select **Lock** and wait until the pattern is virtually displayed.

5.10 To update the parameters of a pattern, click on the related **Action**, then unlock and update the parameters. Select **Lock** again when the pattern update is completed.

5.11 The patterning of multiple proteins in the same micro-well (sequential micro-patterning rounds) requires an accurate alignment of the patterns. To achieve such alignment, upload all sets of desired pattern templates simultaneously (patterns of first and second photopatterning rounds).

5.12 Set the replication and dose parameters of the pattern templates. After parameters are

set, patterns will appear as **Actions** in the **Actions** list. Save this template configuration as a file in software (see **Figure 5D**).

5.13 On the **Actions** list, select only the specific actions to be patterned during the first patterning round and deselect the actions that will be patterned during the second patterning round (see **Figure 5D**).

5.14 Navigate to the area where the reference pattern was produced in step 4.2 (for example, top region of the reference pattern micro-well, see **Figure 4E**). Adjust the focus to the optimal Z-position (obtained in step 4.4).

NOTE: It is strongly recommended to select perfect focus system on the microscope being used (if available), which ensures that the optimal Z-position for patterning will be maintained throughout the entire photopatterning process.

5.15 Select the **Play** symbol to start patterning. Patterning duration will be displayed in the **Estimated Time** panel. On Leonardo software version 4.4 patterning is completed when all **Actions** appear blue in the **Visualization** menu.

6. Protein incubation

NOTE: Micro-wells are incubated with ECM proteins (preferably fluorescently-labeled). These will only bind to the areas where PEG was cleaved through the photopatterning process described in step 5. Each well contains a PDMS stencil with 4 micro-wells, which will allow testing 4 different conditions simultaneously, for example, incubation of a different protein in each micro-well (see **Figure 4D**).

6.1 Use fluorescently-labeled proteins (for example, laminin, fibronectin or fibrinogen conjugated to red or green fluorophores) in order to visualize the micro-patterns (see step 6.7 and 9.4). Alternatively, adsorbed unlabeled proteins can be visualized at later stages using immunofluorescence.

NOTE: ECM proteins (for example, fibronectin) can be labeled using existing protocols³³ and commercially available fluorescence labeling kits (i.e., Alexa 488 labeling kit), or purchased readily labeled (for example, conjugated fibrinogen-488 or laminin-red fluorescent rhodamine, see **Table of Materials**).

6.2 Under sterile conditions, prepare desired concentration of ECM proteins (10 µg/mL laminin, fibronectin or fibrinogen in PBS, see **Table of Materials** and step 4.3).

NOTE: Fluorescently-labeled ECM proteins are light sensitive and should be protected from light (wrap dish in aluminum foil) and should be kept on ice at all times.

6.3 Under sterile conditions, wash the micro-wells three times with 20 µL of PBS to remove

the PLPP.

6.4 Remove 18 μ L of PBS from all the micro-wells and add 20 μ L of ECM protein solution to each micro-well. Incubate at room temperature, for 20-30 min and protect from light.

NOTE: Optimal coating incubation times may vary depending on the type and concentration of proteins.

6.5 After incubation, remove 15 μ L of ECM protein solution and wash the micro-wells three times with 20 μ L of PBS.

NOTE: Always leave approximately 5 μ L of PBS in between washes.

6.6 If patterning with only one protein (one round of micro-patterning), proceed either to storage of the culture dish (step 10) or to cell plating (step 11, and see **Figure 2**). If performing a second round of micro-patterning with a different protein in the same micro-well, proceed either to storage of the culture dish (step 10) or to blocking unspecific binding sites (step 7, and see **Figure 2**).

6.7 Optional quality control step: visualize and image printed patterns using an epifluorescence microscope prior to plating cells. Select appropriate fluorescent channels and adjust exposure times accordingly.

NOTE: To compare the fluorescence intensity of patterns between experiments, it is essential to use the same exposure times for the same proteins.

7. Blocking unspecific binding sites (only for coating with more than protein)

NOTE: Micro-patterning with multiple proteins in the same micro-well involves sequential patterning steps (see **Figure 2**). A blocking agent (PLL-PEG or BSA) is added to the micro-wells to prevent cross-binding, which occurs when the second incubated protein (step 9) binds to the first incubated protein (step 6), thus avoiding a mixture of proteins within the patterns.

7.1 Under sterile conditions, add 20 μ L of PLL-PEG (0.1 mg/mL in PBS) or BSA (1% BSA in PBS) as a blocking step to prevent cross-binding.

NOTE: Blocking efficiency may vary depending on the nature of proteins being used and the affinity among them. It is recommended to test both PLL-PEG and BSA as blocking agents beforehand (**Figure 10D-I**).

7.2 Incubate blocking agent at room temperature for 1 h and protect from light. Remove 15 μ L of blocking agent and wash all micro-wells three times with 20 μ L of PBS. Always leave approximately 5 μ L of PBS in between washes.

7.3 Remove 18 μ L of PBS from all the micro-wells and add 5 μ L of PLPP to each micro-well, ensuring that PLPP is homogenous on the whole surface of the micro-wells.

8. Second round of photopatterning (only for multiple protein patterns)

NOTE: After the first round of photopatterning and protein incubation, a micro-pattern is generated. During the second round of photopatterning, the pattern for the second protein will be generated in the same micro-well (see **Figure 2C** and **Figure 5D**). In the software, select the correct pattern templates (actions), which will be patterned during this round (see **Figure 5D**).

8.1 Navigate to the first micro-well (**Figure 4D**). Ensure appropriate overlap between the digital ROI and the micro-well.

8.2 Load the template configuration saved previously (step 5.12) and select the actions to be patterned during the second photopatterning round (deselect actions from first round, see **Figure 5D**).

8.3 Select the **Play** symbol to start patterning. Ensure the laser is on in the software.

9. Second round of protein incubation (only for multiple protein patterns)

NOTE: In this part of the protocol, fluorescently-labeled protein/s will be incubated on the culture dish after the second round of photopatterning.

9.1 Under sterile conditions, wash the micro-well three times with 20 μ L of PBS to remove the PLPP. Remove 18 μ L of PBS and add 20 μ L of ECM protein solution to each micro-well. Incubate at room temperature for 20-30 min and protect from light.

NOTE: Optimal coating incubation times may vary depending on the type and concentration of proteins.

9.2 After incubation, remove 15 μ L of ECM protein solution and wash the micro-wells three times with 20 μ L of PBS. Always leave approximately 5 μ L of PBS in between washes.

9.3 Either proceed to storage of culture dish (step 10) or to cell plating (step 11, and see **Figure 2**).

9.4 Optional quality control step: using an epifluorescence microscope, visualize and image printed patterns before plating cells. Select appropriate fluorescent channels and adjust exposure time.

NOTE: To compare the fluorescence intensity of patterns between experiments, it is essential to use the same exposure times for the same proteins.

10. Storage of micro-patterns

NOTE: Micro-patterns with adsorbed proteins can be stored during different steps of the protocol (see **Figure 2**). If micro-patterning with multiple proteins, micro-patterns can be stored after the first round of micro-patterning or after the two sequential rounds of micro-patterning have been completed (see **Figure 2B**).

10.1 When passivated with PLL-PEG, store the micro-patterns in PBS (3 mL per well) at 4 °C for up to 3 days.

10.2 If the culture dish was passivated with PEG-SVA, micro-patterns can be stored for up to 1 month. In order to do so, rinse micro-patterns intensively with double distilled deionized water and dry with a sterile air gun of argon or nitrogen, although normal air can also be used. After drying, micro-patterns can be stored at 4 °C for up to 1 month (PRIMO system support team, personal communication).

11. Plating cells

NOTE: During the next steps, cells will be plated on the prepared micro-patterned culture dish/es. In these studies, a neuronal cell line (CAD cells) is used³¹. However this protocol can be adjusted to study other cell types of interest (adjust cell plating protocol as required).

11.1 Differentiate CAD cells for 48 h using differentiation medium (DMEM supplemented with 1% Glutamine, without serum, and with 1% Pen/Strep, see **Table of Materials**).

11.2 Plate 1 mL of medium with cells in the inner glass well, covering all micro-wells and place the culture dish in a 37 °C, 5% CO₂ incubator.

REPRESENTATIVE RESULTS:

Following the above protocol results in micro-patterned surfaces, coated with ECM protein/s of interest. We are using these patterns to track neuronal pathfinding.

Generated patterns should be a precise representation of the template. An example is shown in **Figure 6** where a digital pattern template (**Figure 6A**) representing one design unit (**Figure 5B**), resulted in defined micro-patterns, ranging from 20 to 2 µm width, coated with labeled Fibrinogen (**Figure 6B**). Using ImageJ, fluorescence intensity measurements were obtained both vertically (**Figure 6C**) and horizontally (**Figure 6E**) along the stripe and from a corresponding background region 15 µm above each stripe. The background measurements were subtracted from the pattern measurements for each stripe width.

One limitation of the system is that an edge effect can be observed (**Figure 6B**, top stripe) when printing features ≥20 µm, with a higher intensity signal at the pattern edges compared to the center (**Figure 6D**, first peak of fluorescent intensity profile). In our experiments the resolution limit was approximately 2 µm; at this width we observed a significant decrease (by approximately

50%) in fibrinogen fluorescence intensity compared to the intensity of the wider stripes (**Figure 6F,G**). Patterning using the PRIMO system and the protocol outlined here produced reproducible patterns, with the highest standard deviation of the mean fluorescent intensity measured for the 2 μm width from four individual replicated design units (**Figure 6G**). Variation within the patterned stripes was also found to be low; the coefficient of variation ranged from 3 to 10%, with the 20 μm and 2 μm stripes having the largest internal variation. This is likely to be a result of the edge effect and the resolution limit of the system, respectively. Note that for these measurements we only measured the intensity at the center of the stripes, to avoid the uneven illumination resulting from the objective used to acquire these images (vignetting, **Figure 6E**).

Certain experiments may involve questions that require defined protein concentrations, which can be achieved in two ways: 1) varying the protein concentration (**Figure 7A,B**). Incubation with different concentrations of laminin, results in significantly different fluorescence intensities, increasing with higher protein concentrations (**Figure 7B**). 2) The laser dose that is used to cleave the anti-adhesive film (PEG) can be varied. Higher laser doses will remove the antifouling film to a greater extent, generating more binding sites for the proteins of interest (**Figure 7C,D**) resulting in significantly different fluorescence intensities, increasing with higher laser doses (**Figure 7D**).

Varying the laser dose allows the generation of protein gradients within the same pattern. This is displayed in **Figure 8A**, where a gradient template was designed using different greyscale levels, from black (no laser power) to white (maximum laser power).

Laser intensity is proportional to the greyscale level of the template (ranging from 0 to 255 in an 8-bit image), generating gradients of UV illumination. The measurement of the fluorescence intensity profile along the gradient stripe is linear in the pattern template (**Figure 8B**) and in the generated gradient pattern (**Figure 8C,D**). This is reproducible among all gradient stripes within the same template and gradient pattern (**Figure 8B,D**). Generation of such gradients is extremely useful and helps to mimic in vivo environments where cells often respond to gradients of bioactive proteins³⁴⁻³⁷.

Cells sense changing extracellular environments but assays that enable the study of cell behavior when cells encounter such changes are limited. LIMAP can be used to micro-pattern with multiple proteins in the same micro-well. Examples are shown in **Figure 9** where cross-patterns were generated with stripes of fibronectin (horizontal) and laminin (vertical). When creating patterns with multiple proteins, it is crucial to use a blocking step between the first and second protein incubation, to prevent cross-binding of proteins (see step 7). The blocking efficiency may vary depending on the biochemical characteristics of the proteins that are used for coating and we advise testing several blocking buffers including PLL-PEG (0.1 mg/mL) and BSA (1%). To evaluate this cross-binding effect, we performed fluorescence intensity measurements using ImageJ (**Figure 9**) and we showed that cross-binding can be reduced dramatically, using PLL-PEG buffer (0.1 mg/mL) for fibronectin and laminin cross-patterns (**Figure 9D,H**).

The generated cross-patterns were used for cellular assays with CAD cells (**Figure 10A,B**) or rat dorsal-root ganglion (DRG) neurons (**Figure 10C**). Their neurites (CAD) and axons (DRG) grow

along different lines. CAD cells are used as a neuronal model since they show a similar integrin expression profile compared to primary neurons and they still display actin-rich growth cones after 48 h in culture (**Figure 10B**), making them suitable for pathfinding studies.

In order to investigate the possible cytotoxic effects of the generated micro-patterns towards primary neurons, DRG neurons were isolated and cultured on micro-patterns following a previously published protocol³⁸. The results demonstrate that primary neurons tolerate the micro-patterns environment (**Figure 10C**). We are currently studying how a variety of ECM proteins influence axonal (neurite) pathfinding. Preliminary proof of concepts found in CAD cells will be further investigated using DRG neurons.

In order to validate the quality of generated micro-patterns, it is desirable to image patterns by fluorescence microscopy to ensure the pattern edges are well-defined before proceeding to cell plating. During the imaging process, it is important to ensure the optical adjustment between the microscope and the camera to avoid the peripheral darkening effect (vignetting) which affects the posterior analysis and interpretation of data. Additionally, acquire an image of a pattern-free region using the same exposure times that will be used to image the patterns and subtract this image from the pattern image.

In summary, for good quality micro-pattern generation, it is advisable to assess the protein concentration (**Figure 7A,B**), laser doses (**Figure 7C,D**), protein background levels (**Figure 6E,F,H**) and an efficient blocking step (**Figure 9**) when using multiple proteins. Conclusively, the quality of the micro-patterns generated with LIMAP is essential in order to obtain reliable and reproducible data from cellular assays.

FIGURE LEGENDS:

Figure 1. Scheme of micro-patterning techniques: microcontact printing and laser-assisted patterning. (A) Microcontact printing uses a lithographed master with defined micro-features to generate a PDMS stamp which is incubated with the protein of interest. This protein is then transferred (stamped) onto a glass surface, generating protein micro-patterns. (B) Laser-assisted patterning techniques include photopatterning and direct laser patterning. (C) Most photopatterning approaches use a UV light source and a photomask (either in contact with the substrate surface or in the focal plane of the objective) with desired geometries in order to cleave the PEG antifouling surface in specific positions, creating a defined pattern. A subsequent protein incubation step results in protein adsorption only to the laser-cleaved regions. (D) LIMAP is a photopatterning technique which does not require a photomask in contact with the substrate (i.e., a maskless and contactless approach). LIMAP uses a photo-initiator, which is activated by low doses of a laser, cleaving light-exposed regions of PEG. This creates attachment sites for sequential protein adsorption. (E) Direct laser patterning uses high energy light to directly etch the PEG film, allowing protein binding in those etched regions.

Figure 2. Scheme showing a summary of the steps in the micro-patterning protocol. (A) Micro-patterning with one protein involves only one round of micro-patterning (photopatterning and protein incubation) and can be performed in under 8 h. (B) Micro-patterning with multiple proteins requires two sequential rounds of micro-patterning and can be completed in 1-2 days,

depending on the number of micro-patterns being prepared. It is possible to go through the B version of the protocol in 1 day of work. Continuous arrows indicate direct flow of steps in the protocol. Discontinuous arrows indicate that there is a significant time gap between one step and the other (see step 6.6 and 9.3). (C) Schematic view of example patterns obtained after one round of micro-patterning (red stripes) or two sequential rounds of micro patterning (red and green stripes).

Figure 3. Pattern template generation is versatile with LIMAP. (A,B) Examples of pattern templates designed with ImageJ (A crossbows, B letters). Shapes drawn in white were projected at maximum laser power and shapes drawn in black were not projected. (C,D) Micro-patterns obtained with LIMAP from templates after incubation with 10 µg/mL fibrinogen (green). (C) Crossbows are 50 µm width and 50 µm height spaced by 75 µm horizontally and 50 µm vertically. (D) Letters are 80 µm width and 85 µm height. Scale bars in C and D represents 50 µm.

Figure 4. Essential materials for the LIMAP protocol. (A) Stencils used in this protocol are 20 mm diameter, thin circular-shape PDMS pieces (250 µm thickness) containing 4 micro-wells (4 mm diameter each). The volumes used in the micro-wells range from 5 to 20 µL, considerably reducing the amount of reagents and proteins needed for each experiment. (B) 6 well glass bottom dish where stencils have been already placed in each well. Micro-wells contain 20 µL of PBS to make them visible. (C) Calibration dish in which the inner glass well has been marked with a green highlighter, which will be used to calibrate the laser focus. (D) Schematic view of top left well from the 6-well glass bottom dish in B (outlined with dashed red circle). The inner glass bottom well is represented in white and the stencil is shown in grey. The stencil contains 4 micro-wells (numbered 1-4), for the testing of 4 different experimental conditions (e.g., different protein concentrations, pattern geometries, combinations of proteins, etc.). The asterisk represents the micro-well containing the reference pattern. (E) Schematic view of micro-well where a reference pattern has been generated in the top part (arrow with filled arrowhead). This reference pattern is required to obtain the optimal laser focus for patterning (see step 4). Arrow with empty arrowhead indicates the central area of the micro-well, which will be used for subsequent patterning after system calibration.

Figure 5. Software set-up for micro-patterning. (A) Pattern template with parallel stripes designed with ImageJ and saved as an 8-bit Tiff file. (B-D) Schematic view of digital ROIs (regions of interest) that will overlap with the current micro-wells where micro-patterns will be generated. (B) The pattern template to be used (length 1824 pixel=415 µm, width 1140 pixel=260 µm) is selected on Leonardo and is projected on the ROI as a design unit (red stripes in the black dashed rectangle), which will cover approximately 0.1 mm² of the micro-well area. The design unit is replicated in 4 columns and 4 lines in the **Replication** menu (template configuration), creating a pattern across the micro-well. NOTE the space among the columns. (C) To pattern continuous stripes, the spacing between columns has to be adjusted. In this case, to achieve an overlap among design units the spacing between columns is set in the replication menu as **Negative spacing**, -20 µm. (D) In order to pattern multiple proteins in the same micro-well, an accurate alignment of the patterns is required. During the software set-up step (step 5), upload all desired pattern templates simultaneously. On the **Actions** list, select only the specific actions

to be patterned during each patterning round and deselect the rest of the actions (step 5.12, 5.13 and 8.2).

Figure 6. Analysis of pattern variability using LIMAP. (A) Pattern template designed with ImageJ used to micro-pattern four stripes of varying width (20, 10, 5, 2 μm , from top to bottom). (B) Micro-pattern obtained after incubation with 10 $\mu\text{g}/\text{mL}$ of fluorescently-labeled fibrinogen (green). (C) Intensity measurements along a vertical line crossing the stripes of the micro-pattern. (D) Vertical fluorescence intensity profile obtained from measurement in (C). NOTE that at larger widths (20 μm) there is a variation in the vertical profile caused by the accumulation of protein at the edges of the stripe, resulting in two distinct fluorescence intensity peaks (edge effect). This effect is only seen in stripe widths ≥ 20 μm . (E) Intensity measurements along the depicted horizontal lines (fluorescence and background). (F) Horizontal fluorescence intensity profiles obtained from measurements in (E). (G) Graph showing the mean intensity for each stripe width, measured from four individual replicated design units (inter-pattern variation). NOTE the reduced protein adsorption to patterns of 2 μm stripe width. (H) The variation within the patterned stripes (coefficient of variation) was low for all stripe widths, ranging from 3 to 10%. Data in G and H shown as mean \pm SD. Statistical analysis in G and H was performed using one-way ANOVA (Kruskal-Wallis) non-parametric test with multiple comparisons. P value is <0.001 for ** significance.

Figure 7. The effect of variations in laser power and protein concentration for protein adsorption efficiency. (A) The PLL-PEG surface was laser-cleaved with a constant laser dose (1390 mJ/mm^2) and incubated with the indicated concentrations of fluorescently-labeled laminin (magenta). (B) Quantification of fluorescence intensity of the laminin stripes in (A). (C) The different indicated laser doses were applied followed by incubation with the same concentration (10 $\mu\text{g}/\text{mL}$) of fluorescently-labeled fibronectin (green). (D) Quantification of fluorescence intensity of the fibronectin stripes in (C) showing that higher laser doses correlate to higher levels of adsorbed protein. All measurements are background subtracted. Sample numbers are indicated at the bottom of columns; data is shown as mean \pm SEM. Statistical analysis was performed using non-parametric Mann-Whitney test with two-tailed calculation. P-value is <0.0001 for **** significance.

Figure 8. Generation of a protein concentration gradient within a micro-pattern. (A) Gradient pattern template in greyscale. (B) Fluorescence intensity profile measured from (A). (C) Pattern obtained with LIMAP from pattern template in (A) after incubation with 10 $\mu\text{g}/\text{mL}$ of fluorescently-labeled fibronectin (green). (D) Fluorescence intensity profile of $n = 3$ stripes and background represented as mean \pm SEM, showing the linear increase in intensity of the protein gradient.

Figure 9. The cross-binding effect when patterning several proteins sequentially. (A-C, E-G) Cross-patterns with 10 μm stripes of fluorescently-labeled fibronectin (cyan, horizontal) and fluorescently-labeled laminin (magenta, vertical). (A-C) Samples treated with BSA blocking buffer. (E-F) Samples treated with PLL-PEG for blocking unspecific binding sites (step 7). (A,E) Merged fluorescence channels showing both fibronectin and laminin. (B,F) Image showing fibronectin

only. (C,G) Image showing laminin only. For C, note the presence of laminin also on horizontal fibronectin positive stripes which is due to the ineffective blocking of unoccupied binding sites with BSA. For G, note that blocking with PLL-PEG prevents efficiently binding of laminin to the fibronectin stripes. (D,H) Fluorescence intensity profiles obtained from indicated measurements (diagonal yellow line) in A and E, respectively.

Figure 10. Cross-patterns to investigate neurite/axon pathfinding. (A-C) Cross-patterns with 10 μm stripes of fluorescently-labeled fibronectin (cyan, horizontal) and fluorescently-labeled laminin (magenta, vertical). (A,B) Fluorescent images of CAD cells with neurites growing along the micro-patterns. To visualize neurites, cells were cultured for 48 h, fixed with 4% PFA and stained for tubulin (A) or tubulin and actin (B). (C) Rat dorsal-root ganglion (DRG) neurons with axons growing along the micro-patterns. To visualize axons, DRG neurons were cultured for 72 h, fixed with 4% PFA and stained for tubulin.

Figure 11. Examples of common negative results obtained when generating micro-patterns with LIMAP. (A-C) Sub-optimal patterned stripes of 10 $\mu\text{g/mL}$ fluorescently-labeled fibrinogen (green) obtained under different circumstances. (A) The micro-well dried out during pattern generation. Note the high levels of fluorescence in the background (arrow) and the presence of PBS crystals (asterisks). (B) The stitching between the stripes was not properly adjusted during software set-up resulting in discontinuous stripes with gaps (arrow) among design units (see step 3.4.7). (C) The laser focus was sub-optimal causing diffused stripes (arrows) which do not represent the actual stripe widths of the pattern template, which should be 20, 10, 5, 2 μm from top to bottom, as in (B).

DISCUSSION:

LIMAP (PRIMO) micro-patterning advantages and comparison with microcontact printing

While microcontact printing is possibly the most commonly used micro-patterning technique in the biological field³⁹, there seems to be an increasing number of researchers using LIMAP technology⁴⁰⁻⁴⁴. Here, we presented a protocol using PRIMO, a commercially available system for LIMAP. Below we briefly discuss potential advantages and limitations of microcontact printing and LIMAP photo-patterning.

Microcontact printing requires lithographed masters produced by spin-coating a photo-mask (generally SU-8) onto a glass or silicon wafer, which is then laser-etched with the desired micro-features. These masters are used as templates to create a PDMS stamp⁴⁵. The stamp is incubated with a chosen protein that adsorbs to it, and is then transferred (stamped) onto the cell culture dish. The process of adsorption of the protein to the PDMS stamp is dependent on protein concentration, buffer and incubation time. These parameters need to be tested beforehand for optimal results⁴⁶.

Masters can be used in a substantial number of experiments, lasting for months or even years, if correctly preserved. However, a limiting factor of this technology is the necessity to re-design new lithographed masters for every desired modification. Changes in experimental designs may result in time-consuming production of new masters (up to several weeks) thus delaying

experiments. In comparison, LIMAP photopatterning does not require a physical master; it uses software-generated pattern templates that can be used to flexibly adapt desired geometries of micro-patterns to changing research questions. LIMAP can be also used to generate protein gradients within the same micro-pattern (**Figure 8**), which is harder to obtain in a reproducible manner using microcontact printing⁴⁷.

Furthermore, the micro-pattern resolution achieved with LIMAP, in our case, is 2 μm (**Figure 6B**). Approaching this resolution increased intra- and inter- pattern variability. Generating patterns around or above 10 μm width was highly reproducible (**Figure 6G,H**). On the contrary, with microcontact printing it is difficult to consistently obtain resolutions below 10 μm and it is common to find artefacts when stamping small features (data not shown).

We have shown that LIMAP can be used to micro-pattern multiple proteins (**Figure 9**) within the same micro-well, allowing further levels of complexity to be added to experiments. Although this could be achieved with microcontact printing, aligning different proteins with high level of precision can be technically rather demanding. Whilst patterning multiple proteins using LIMAP seems straight forward, it is important to mention that cross-binding of proteins through sequential coating procedures can be reduced through blocking reagents but not entirely eliminated (**Figure 9**).

Regarding the cost of one or the other technique, LIMAP as described here requires the purchase of micro-patterning equipment (PRIMO) that can be installed on different fluorescence microscopes and requires a motorized stage. Although this investment is initially cost intensive, there are no additional purchases other than consumable items (stencils, PEG and PLPP) on the long run associated with LIMAP. Alternatively, the PDMS stencils can also be produced in the lab by the own experimenter following published protocols^{18,32}. The largest costs for microcontact printing may be associated with the production of new masters, which can become substantial if experiments require new patterns.

One drawback of LIMAP is the relatively low throughput approach of this technique. Microcontact printing can produce a large number of micro-patterns quickly and efficiently in a simultaneous stamping step, compared to the required sequential laser micro-patterning with LIMAP. For example, it is possible to produce 6 stamped glass coverslips in about 2 h with microcontact printing using PDMS stamps (excluding stamp preparation); patterning a similar area (6-well dish) with LIMAP would take around 4 h, excluding the procedure of surface passivation (considering the pattern template configuration described in step 5.12 and see **Figure 5B**).

Another rate limiting factor of LIMAP technology is the long illumination time required for patterning large areas (30 s per design unit with a 7.5 mW/mm² laser). In these cases, microcontact printing might be a preferred option. A newly available photo-initiator (PLPP gel, **Table of Materials**) should considerably reduce the time taken for the patterning, allowing the generation of hundreds of micro-patterns in large areas (up to 8 mm²) in just a few minutes.

Another important factor to take into account when micro-patterning surfaces for cell culture is the reproducibility of the micro-patterns among different experimental repeats, in comparison to the variability obtained with microcontact printing. For example, the graphs shown in **Figure 7B,D** are representative data of three independent experimental repeats with very similar results (data not shown). Based on our experience and previous publications, this level of reproducibility is difficult to achieve with microcontact printing⁴⁸⁻⁵².

In contrast to other photo-patterning techniques that require either dedicated chemistry to engineer photosensitive materials or the use of photo-sensitizers, which are generally not very biocompatible³, the photosensitive component of LIMAP (PLPP) is biocompatible and well-tolerated by cells²¹; in our hands we have not experienced any cytotoxicity across a variety of cells, including CAD, DRG neurons (**Figure 10**), fibroblasts, epithelial cells, and melanoma cells (data not shown). Another advantage of LIMAP using PRIMO compared to other photo-patterning techniques is that no photomask is required. Similar to microcontact printing, new photomasks would need to be designed and generated for every desired pattern.

All the limitations mentioned above for microcontact printing, refer to the manual approach of the technique. However, it is possible to enhance the throughput and the reproducibility of microcontact printing using an automated device with stamp load and pressure control⁵³.

Key steps of the protocol and problems solving for LIMAP using PRIMO

One of the most common problems found during this protocol is having high levels of background fluorescence within the micro-patterns. This can be due to the drying out of micro-wells which often occurs due to their small volume. When this occurs, PBS crystals often appear surrounding the ECM patterns (**Figure 11A**).

Insufficient or inefficient washing steps after protein incubation can also result in high levels of background fluorescence. This can be observed particularly upon using protein concentrations of 10 µg/mL (**Figure 11B**) or higher. The excess of protein in the background can be reduced by including additional washing steps with PBS.

The presence of protein background needs to be measured and characterized in each experiment, calculating the background fluorescence intensity (**Figure 6E**) and subtracting it from the micro-patterns intensity (**Figure 6F-H** and **Figure 7B,D**). High protein background may have an impact in the attachment and sprouting of CAD cells, compromising the interpretation of results.

Having gaps between design units is a common problem when users have limited experience (**Figure 11B**), which occurs as a result of insufficient overlap between patterns. Two parameters in the Leonardo software can be adjusted to overcome this: 1) a negative spacing between columns may be required, depending on the design of the pattern (step 5.7 and see **Figure 5B,C**). Alternatively, 2) use the gradient option in the **Expert** menu to stitch the columns. A quick test to determine the optimal spacing parameters can be performed using UV adhesive (**Table of Materials**). A small drop of this adhesive is applied to a glass slide, which is then covered with a

glass coverslip, making a film. The embedded UV adhesive is photopatterned with the pattern template of interest using a low laser dose (30 mJ/mm²). The UV- exposed regions of the embedded adhesive will be cured, becoming visible under bright-field microscopy. The test results are visualized to evaluate the obtained spacing within the pattern. In our neuronal experiments, a gap between stripes may adversely affect cell behavior, producing variations in growth dynamics (either reduced speed or abandonment of the path).

In the latest update of Leonardo software (at the time of publication, Leonardo 4.11), it is possible to upload previously designed bigger pattern templates that cover a much larger area (up to 8 mm² using the 20X objective) of the micro-well surface compared to the current 0.1 mm² per design unit, eliminating the need to stitch together the smaller design units. Undefined edges can result from a lack of laser focus adjustment during pattern generation (**Figure 11C**). It is therefore critical to calibrate the laser and perform reference pattern steps (see step 4) prior to patterning. Poorly defined stripes result in variations in stripe width, making the correlation between axon growth dynamics and stripe width difficult. Axons also tend to abandon stripes that have diffused edges. Additionally, variability in edges can also be found when printing stripes of 10-20 µm width or higher, resulting in a higher protein content at the edges compared to the central regions of the pattern (**Figure 6B,D**). This edge effect is produced by a non-homogeneous diffusion of the photo-initiator during the photopatterning process. The photoscission reaction is oxygen dependent, which diffuses more at the edges. This edge effect can be minimized homogenizing the photo-initiator with a pipette in the micro-well during the photopatterning process. Furthermore, a new commercialized photo-initiator (PLPP gel), can also reduce the edge effect (PRIMO system support team, personal communication).

Micro-printing of more than one protein can result in cross-binding (**Figure 9A-D**). This can be minimized by increasing the blocking efficiency that is used to occupy unspecific binding sites between the incubation steps for the two different proteins. Cross-binding of proteins can perturb reproducibility of experimental outcomes and may lead to misinterpretation of data, since it is difficult to determine the contribution of each protein to axon growth dynamics and to other cell behaviors.

Conclusion

We hope that the provided protocol using LIMAP facilitates the generation of protein micro-patterns through the use of the PRIMO system. Whilst our protocol focuses on how to reliably produce micro-patterns in 2D glass surfaces, others have shown that it is possible to use LIMAP for micro-patterning of soft substrates⁵⁴, and microstructured surfaces for 3D cultures⁴². These micro-patterns can be a versatile tool to study cellular responses to changes in their micro-environment.

ACKNOWLEDGMENTS:

This work is supported by the BBSRC, EPSRC, MRC and Wellcome Trust. The C.B. laboratory is part of the Wellcome Trust Centre for Cell-Matrix research, University of Manchester, which is supported by core funding from the Wellcome Trust (grant number 088785/Z/09/Z). The authors wish to acknowledge the funding provided by the Biotechnology and Biological Sciences Research

Council (BBSRC) to C.M., K.J. (BB/M020630/1) and P.A. (BB/P000681/1) and by the Engineering and Physical Sciences Research Council (EPSRC) and Medical Research Council (MRC) Centre for Doctoral Training in Regenerative Medicine to A.K. (EP/L014904/1). The authors thank Alvéole for their correspondence and their after-sales support team. The authors thank Peter March and Roger Meadows from the Bioimaging Facility, University of Manchester for their help with the microscopy. The Bioimaging Facility microscopes used in this study were purchased with grants from BBSRC, Wellcome Trust and the University of Manchester Strategic Fund.

DISCLOSURES:

The authors have nothing to disclose.

REFERENCES:

- 1 Alamdari, O. G., Seyedjafari, E., Soleimani, M. & Ghaemi, N. Micropatterning of ECM Proteins on Glass Substrates to Regulate Cell Attachment and Proliferation. *Avicenna Journal of Medical Biotechnology*. **5** (4), 234-240 (2013).
- 2 Sunami, H., Yokota, I. & Igarashi, Y. Influence of the pattern size of micropatterned scaffolds on cell morphology, proliferation, migration and F-actin expression. *Biomaterials Science*. **2** (3), 399-409, doi:10.1039/c3bm60237k, (2014).
- 3 Thery, M. Micropatterning as a tool to decipher cell morphogenesis and functions. *Journal of Cell Science*. **123** (Pt 24), 4201-4213, doi:10.1242/jcs.075150, (2010).
- 4 Marino, A. *et al.* Two-photon polymerization of sub-micrometric patterned surfaces: investigation of cell-substrate interactions and improved differentiation of neuron-like cells. *ACS Applied Materials & Interfaces*. **5** (24), 13012-13021, doi:10.1021/am403895k, (2013).
- 5 Joo, S. *et al.* Effects of ECM protein micropatterns on the migration and differentiation of adult neural stem cells. *Scientific Reports*. **5**, doi:ARTN 1304310.1038/srep13043, (2015).
- 6 Morgani, S. M., Metzger, J. J., Nichols, J., Siggia, E. D. & Hadjantonakis, A. K. Micropattern differentiation of mouse pluripotent stem cells recapitulates embryo regionalized cell fate patterning. *Elife*. **7**, doi:10.7554/eLife.32839, (2018).
- 7 Javaherian, S., O'Donnell, K. A. & McGuigan, A. P. A Fast and Accessible Methodology for Micro-Patterning Cells on Standard Culture Substrates Using Parafilm (TM) Inserts. *Plos One*. **6** (6), doi:ARTN e20909 10.1371/journal.pone.0020909, (2011).
- 8 Smirnov, M. S., Cabral, K. A., Geller, H. M. & Urbach, J. S. The effects of confinement on neuronal growth cone morphology and velocity. *Biomaterials*. **35** (25), 6750-6757, doi:10.1016/j.biomaterials.2014.04.097, (2014).
- 9 Albert, P. J. & Schwarz, U. S. Dynamics of Cell Ensembles on Adhesive Micropatterns: Bridging the Gap between Single Cell Spreading and Collective Cell Migration. *PLOS Computational Biology*. **12** (4), e1004863, doi:10.1371/journal.pcbi.1004863, (2016).
- 10 Evans, A. R. *et al.* Laminin and fibronectin modulate inner ear spiral ganglion neurite outgrowth in an in vitro alternate choice assay. *Developmental Neurobiology*. **67** (13), 1721-1730, doi:10.1002/dneu.20540, (2007).
- 11 Nichol, R. H. t., Hagen, K. M., Lombard, D. C., Dent, E. W. & Gomez, T. M. Guidance of Axons by Local Coupling of Retrograde Flow to Point Contact Adhesions. *Journal of*

968 *Neuroscience*. **36** (7), 2267-2282, doi:10.1523/JNEUROSCI.2645-15.2016, (2016).

969 12 Burdick, J. A., Khademhosseini, A. & Langer, R. Fabrication of gradient hydrogels using a
970 microfluidics/photopolymerization process. *Langmuir*. **20** (13), 5153-5156,
971 doi:10.1021/la049298n, (2004).

972 13 Schwartz, P. V. Molecular transport from an atomic force microscope tip: A comparative
973 study of dip-pen nanolithography. *Langmuir*. **18** (10), 4041-4046, doi:UNSP
974 LA011652J10.1021/la011652j, (2002).

975 14 Barbulovic-Nad, I. *et al.* Bio-microarray fabrication techniques--a review. *Critical Reviews*
976 *in Biotechnology*. **26** (4), 237-259, doi:10.1080/07388550600978358, (2006).

977 15 Shafagh, R. Z., Vastesson, A., Guo, W. J., van der Wijngaart, W. & Haraldsson, T. E-Beam
978 Nanostructuring and Direct Click Biofunctionalization of Thiol-Ene Resist. *Acs Nano*. **12**
979 (10), 9940-9946, doi:10.1021/acsnano.8b03709, (2018).

980 16 Kobayashi, J., Yamato, M., Itoga, K., Kikuchi, A. & Okano, T. Preparation of microfluidic
981 devices using micropatterning of a photosensitive material by a maskless, liquid-crystal-
982 display projection method. *Advanced Materials*. **16** (22), 1997-+,
983 doi:10.1002/adma.200400312, (2004).

984 17 Bernard, A. *et al.* Printing patterns of proteins. *Langmuir*. **14** (9), 2225-2229, doi:DOI
985 10.1021/la980037l, (1998).

986 18 Ruiz, S. A. & Chen, C. S. Microcontact printing: A tool to pattern. *Soft Matter*. **3** (2), 168-
987 177, doi:10.1039/b613349e, (2007).

988 19 Qin, D., Xia, Y. & Whitesides, G. M. Soft lithography for micro- and nanoscale patterning.
989 *Nature Protocols*. **5** (3), 491-502, doi:10.1038/nprot.2009.234, (2010).

990 20 Fink, J. *et al.* Comparative study and improvement of current cell micro-patterning
991 techniques. *Lab Chip*. **7** (6), 672-680, doi:10.1039/b618545b, (2007).

992 21 Strale, P. O. *et al.* Multiprotein Printing by Light-Induced Molecular Adsorption. *Advanced*
993 *Materials*. **28** (10), 2024-2029, doi:10.1002/adma.201504154, (2016).

994 22 Belisle, J. M., Correia, J. P., Wiseman, P. W., Kennedy, T. E. & Costantino, S. Patterning
995 protein concentration using laser-assisted adsorption by photobleaching, LAPAP. *Lab*
996 *Chip*. **8** (12), 2164-2167, doi:10.1039/b813897d, (2008).

997 23 Belisle, J. M., Kunik, D. & Costantino, S. Rapid multicomponent optical protein patterning.
998 *Lab Chip*. **9** (24), 3580-3585, doi:10.1039/b911967a, (2009).

999 24 Heinz, W. F., Hoh, M. & Hoh, J. H. Laser inactivation protein patterning of cell culture
1000 microenvironments. *Lab Chip*. **11** (19), 3336-3346, doi:10.1039/c1lc20204a, (2011).

1001 25 Azioune, A., Carpi, N., Tseng, Q., Thery, M. & Piel, M. Protein Micropatterns: A Direct
1002 Printing Protocol Using Deep UVs. *Microtubules: In Vivo*. **97** 133-146, doi:10.1016/S0091-
1003 679x(10)97008-8, (2010).

1004 26 Vignaud, T., Ennomani, H. & Thery, M. Polyacrylamide hydrogel micropatterning.
1005 *Methods in Cell Biology*. **120** 93-116, doi:10.1016/B978-0-12-417136-7.00006-9, (2014).

1006 27 Waldbaur, A., Waterkotte, B., Schmitz, K. & Rapp, B. E. Maskless projection lithography
1007 for the fast and flexible generation of grayscale protein patterns. *Small*. **8** (10), 1570-1578,
1008 doi:10.1002/smll.201102163, (2012).

1009 28 Kang, J., Choi, J. C., Kim, M., Jung, H. R. & Doh, J. Photopatterning with a printed
1010 transparency mask and a protein-friendly photoresist. *Methods in Cell Biology*. **119** 55-72,
1011 doi:10.1016/B978-0-12-416742-1.00004-4, (2014).

1012 29 Falconnet, D., Csucs, G., Grandin, H. M. & Textor, M. Surface engineering approaches to
1013 micropattern surfaces for cell-based assays. *Biomaterials*. **27** (16), 3044-3063,
1014 doi:10.1016/j.biomaterials.2005.12.024, (2006).

1015 30 Morlat, S. & Gardette, J. L. Phototransformation of water-soluble polymers. Part II:
1016 photooxidation of poly(ethylene oxide) in aqueous solution. *Polymer*. **44** (26), 7891-7897,
1017 doi:10.1016/j.polymer.2003.10.022, (2003).

1018 31 Qi, Y., Wang, J. K., McMillian, M. & Chikaraishi, D. M. Characterization of a CNS cell line,
1019 CAD, in which morphological differentiation is initiated by serum deprivation. *Journal of*
1020 *Neuroscience*. **17** (4), 1217-1225 (1997).

1021 32 Shrirao, A. B. *et al.* A Versatile Method of Patterning Proteins and Cells. *Journal of*
1022 *Visualized Experiments*. (120), doi:10.3791/55513, (2017).

1023 33 Pankov, R. & Momchilova, A. Fluorescent labeling techniques for investigation of
1024 fibronectin fibrillogenesis (labeling fibronectin fibrillogenesis). *Methods in Molecular*
1025 *Biology*. **522** 261-274, doi:10.1007/978-1-59745-413-1_18, (2009).

1026 34 Dertinger, S. K., Jiang, X., Li, Z., Murthy, V. N. & Whitesides, G. M. Gradients of substrate-
1027 bound laminin orient axonal specification of neurons. *Proceedings of the National*
1028 *Academy of Sciences of the United States of America*. **99** (20), 12542-12547,
1029 doi:10.1073/pnas.192457199, (2002).

1030 35 Chelli, B. *et al.* Neural cell alignment by patterning gradients of the extracellular matrix
1031 protein laminin. *Interface Focus*. **4** (1), 20130041, doi:10.1098/rsfs.2013.0041, (2014).

1032 36 Tang, Y., Qiu, Q. F., Zhang, F. L., Xie, M. & Huang, W. H. Quantifying orientational
1033 regeneration of injured neurons by natural product concentration gradients in a 3D
1034 microfluidic device. *Lab Chip*. **18** (6), 971-978, doi:10.1039/c7lc01143a, (2018).

1035 37 Srinivasan, P., Zervantonakis, I. K. & Kothapalli, C. R. Synergistic effects of 3D ECM and
1036 chemogradients on neurite outgrowth and guidance: a simple modeling and microfluidic
1037 framework. *PLoS One*. **9** (6), e99640, doi:10.1371/journal.pone.0099640, (2014).

1038 38 de Luca, A. C., Faroni, A. & Reid, A. J. Dorsal root ganglia neurons and differentiated
1039 adipose-derived stem cells: an in vitro co-culture model to study peripheral nerve
1040 regeneration. *Journal of Visualized Experiments*. (96), doi:10.3791/52543, (2015).

1041 39 Khadpekar, A. J., Khan, M., Sose, A. & Majumder, A. Low Cost and Lithography-free Stamp
1042 fabrication for Microcontact Printing. *Scientific Reports*. **9** (1), 1024, doi:10.1038/s41598-
1043 018-36521-x, (2019).

1044 40 Delepine, C. *et al.* Altered microtubule dynamics and vesicular transport in mouse and
1045 human MeCP2-deficient astrocytes. *Human Molecular Genetics*. **25** (1), 146-157,
1046 doi:10.1093/hmg/ddv464, (2016).

1047 41 Decock, J., Schlenk, M. & Salmon, J. B. In situ photo-patterning of pressure-resistant
1048 hydrogel membranes with controlled permeabilities in PEGDA microfluidic channels. *Lab*
1049 *Chip*. **18** (7), 1075-1083, doi:10.1039/c7lc01342f, (2018).

1050 42 Stoecklin, C. *et al.* A New Approach to Design Artificial 3D Microniches with Combined
1051 Chemical, Topographical, and Rheological Cues. *Advanced Biosystems*. **2** (7), doi:ARTN
1052 170023710.1002/adbi.201700237, (2018).

1053 43 Toraille, L. *et al.* Optical Magnetometry of Single Biocompatible Micromagnets for
1054 Quantitative Magnetogenetic and Magnetomechanical Assays. *Nano Letters*.
1055 doi:10.1021/acs.nanolett.8b03222, (2018).

1056 44 Theodoly, O. *et al.* Live nanoscopic to mesoscopic topography reconstruction with an
1057 optical microscope for chemical and biological samples. *PLoS One*. **13** (12), e0207881,
1058 doi:10.1371/journal.pone.0207881, (2018).

1059 45 Ermis, M., Antmen, E. & Hasirci, V. Micro and Nanofabrication methods to control cell-
1060 substrate interactions and cell behavior: A review from the tissue engineering
1061 perspective. *Bioactive Materials*. **3** (3), 355-369, doi:10.1016/j.bioactmat.2018.05.005,
1062 (2018).

1063 46 von Philipsborn, A. C. *et al.* Microcontact printing of axon guidance molecules for
1064 generation of graded patterns. *Nature Protocols*. **1** (3), 1322-1328,
1065 doi:10.1038/nprot.2006.251, (2006).

1066 47 Ricoult, S. G., Kennedy, T. E. & Juncker, D. Substrate-bound protein gradients to study
1067 haptotaxis. *Frontiers in Bioengineering and Biotechnology*. **3** 40,
1068 doi:10.3389/fbioe.2015.00040, (2015).

1069 48 Bietsch, A. & Michel, B. Conformal contact and pattern stability of stamps used for soft
1070 lithography. *Journal of Applied Physics*. **88** (7), 4310-4318, doi:Doi 10.1063/1.1289816,
1071 (2000).

1072 49 Hui, C. Y., Jagota, A., Lin, Y. Y. & Kramer, E. J. Constraints on microcontact printing imposed
1073 by stamp deformation. *Langmuir*. **18** (4), 1394-1407, doi:10.1021/la0113567, (2002).

1074 50 Sharp, K. G., Blackman, G. S., Glassmaker, N. J., Jagota, A. & Hui, C. Y. Effect of stamp
1075 deformation on the quality of microcontact printing: theory and experiment. *Langmuir*.
1076 **20** (15), 6430-6438, doi:10.1021/la036332+, (2004).


1077 51 Delamarche, E., Schmid, H., Michel, B. & Biebuyck, H. Stability of molded
1078 polydimethylsiloxane microstructures. *Advanced Materials*. **9** (9), 741-746, doi:DOI
1079 10.1002/adma.19970090914, (1997).


1080 52 Perl, A., Reinhoudt, D. N. & Huskens, J. Microcontact Printing: Limitations and
1081 Achievements. *Advanced Materials*. **21** (22), 2257-2268, doi:10.1002/adma.200801864,
1082 (2009).


1083 53 Chakra, E. B., Hannes, B., Dilosquer, G., Mansfield, C. D. & Cabrera, M. A new instrument
1084 for automated microcontact printing with stamp load adjustment. *Review of Scientific
1085 Instruments*. **79** (6), doi:Artn 06410210.1063/1.2936259, (2008).

1086 54 Pasturel, S., Studer. Tailoring 3D cell culture templates with common hydrogels. *bioRxiv*.
1087 doi:10.1101/370882, (2019).


1088








A




B




C





D





E




















TABLE OF MATERIALS

Material/Equipment	Brand	Catalogue number (if applicable)
Alexa 488 protein labeling kit	Invitrogen	A10235
Alexa 647 protein labeling kit	Invitrogen	A20173
CAD cells	ECACC	8100805
Conjugated fibrinogen-488	Molecular Probes	F13191
DMEM culture medium	Gibco	11320033
Epifluorescence Microscope**	Nikon	Eclipse Ti inverted
Fibronectin	Sigma	F4759
Fiji-Image J	www.imagej.nih.gov	Version 2.0.0-rc-54/1.51f
Fluorescent highlighter	Stabilo	Stabilo Boss Original
HEPES	Gibco	15630080
Inkscape software	Inkscape	Check last update
Laminin-red fluorescent rhodamine	Cytoskeleton, Inc.	LMN01
Leonardo software	Alvéole	version 4.11
L-Glutamine	Sigma	G7513
Micro-manager software	Open imaging	Check last update
Motorized x/y stage	PRIOR Scientific	Proscan II
NIS Elements Software	Nikon	NIS Elements AR 4.60.00 64-bit (With Nikon jobs)

PBS (without Ca^{2+} , Mg^{2+})	Sigma	D8537
PDMS Stencils	Alvéole	visit www.alveolelab.com
PEG-SVA	Laysan bio, Inc.	MPEG-SVA-5000-1g
Phalloidin 405	Abcam	ab176752
Photo-initiator (PLPP)	Alvéole	Classic PLPP
Photo-initiator (PLPP gel)	Alvéole	PLPP gel
Plasma cleaner	Harrick Plasma	PDC-32G (115V) PDC-32G-2 (230V)
PLL-PEG	SuSoS (also distributed by Alvéole)	visit www.alveolelab.com
Poly-L-Lysine	Sigma	P4707
Primo equipment	Alvéole	visit www.alveolelab.com
Pen/Strep	Thermo Fisher	15140122
Tubulin anti-alpha antibody	Abcam	DM1A
Tubulin anti-beta 3 antibody	Sigma	T8660
UV adhesive	Norland Products	NOA81
1 well glass bottom dish	Cellvis	D35-20-1.5-N
6 well glass bottom dish	Cellvis	P06-20-1.5-N
20x objective**	Nikon	no phase ring (check updated catalogue)

Working concentration (if applicable)
N.A.
N.A.
N.A.
10 µg/ml
N.A.
N.A.
10 µg/ml (after labelling with Alexa 488 protein labeling kit, see above) (diluted in PBS)
N.A.
N.A.
1M
N.A.
10 µg/ml (diluted in PBS)
N.A.
1%
N.A.
N.A.
N.A.

1X
N.A.
50 mg/ml
1:1000
14.5 mg/ml
4.76% diluted in ethanol
N.A.
0.1 mg/ml (diluted in PBS)
0.01%
N.A.
1%
1:1000
1:500
N.A.
N.A.
N.A.
N.A.

Comments

ARTICLE AND VIDEO LICENSE AGREEMENT - UK

Title of Article:	Light-Induced Molecular Adsorption of Proteins using the PRIMO system for micro-patterning to study cell responses to extracellular matrix proteins
Author(s):	Cristina Melero*1, Aljona Kolmogorova*1, Paul Atherton1, Brian Derby2, Adam Reid3,4, Karin Jansen#5, Christoph Ballestrem#1

Item 1: The Author elects to have the Materials be made available (as described at <http://www.jove.com/publish>) via:

☐

Standard Access

☒

Open Access

Item 2: Please select one of the following items:

☒

The Author is **NOT** a United States government employee.

☐

The Author is a United States government employee and the Materials were prepared in the course of his or her duties as a United States government employee.

☐

The Author is a United States government employee but the Materials were NOT prepared in the course of his or her duties as a United States government employee.

ARTICLE AND VIDEO LICENSE AGREEMENT

1. **Defined Terms.** As used in this Article and Video License Agreement, the following terms shall have the following meanings: "**Agreement**" means this Article and Video License Agreement; "**Article**" means the article specified on the last page of this Agreement, including any associated materials such as texts, figures, tables, artwork, abstracts, or summaries contained therein; "**Author**" means the author who is a signatory to this Agreement; "**Collective Work**" means a work, such as a periodical issue, anthology or encyclopedia, in which the Materials in their entirety in unmodified form, along with a number of other contributions, constituting separate and independent works in themselves, are assembled into a collective whole; "**CRC License**" means the Creative Commons Attribution 3.0 Agreement (also known as CC-BY), the terms and conditions of which can be found at: <http://creativecommons.org/licenses/by/3.0/us/legalcode>; "**Derivative Work**" means a work based upon the Materials or upon the Materials and other pre-existing works, such as a translation, musical arrangement, dramatization, fictionalization, motion picture version, sound recording, art reproduction, abridgment, condensation, or any other form in which the Materials may be recast, transformed, or adapted; "**Institution**" means the institution, listed on the last page of this Agreement, by which the Author was employed at the time of the creation of the Materials; "**JoVE**" means MyJoVE Corporation, a Massachusetts corporation and the publisher of The Journal of Visualized Experiments; "**Materials**" means the Article and / or the Video; "**Parties**" means the Author and JoVE; "**Video**" means any video(s) made by the Author, alone or in conjunction with any other parties, or by JoVE or its affiliates or agents, individually or in collaboration with the Author or any other parties, incorporating all or any portion

of the Article, and in which the Author may or may not appear.

2. **Background.** The Author, who is the author of the Article, in order to ensure the dissemination and protection of the Article, desires to have the JoVE publish the Article and create and transmit videos based on the Article. In furtherance of such goals, the Parties desire to memorialize in this Agreement the respective rights of each Party in and to the Article and the Video.

3. **Grant of Rights in Article.** In consideration of JoVE agreeing to publish the Article, the Author hereby grants to JoVE, subject to **Sections 4 and 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Article in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Article into other languages, create adaptations, summaries or extracts of the Article or other Derivative Works (including, without limitation, the Video) or Collective Works based on all or any portion of the Article and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. If the "Open Access" box has been checked in **Item 1** above, JoVE and the Author hereby grant to the public all such rights in the Article as provided in, but subject to all limitations and requirements set forth in, the CRC License.

4. **Retention of Rights in Article.** Notwithstanding the exclusive license granted to JoVE in **Section 3** above, the Author shall, with respect to the Article, retain the non-exclusive right to use all or part of the Article for the non-commercial purpose of giving lectures, presentations or teaching classes, and to post a copy of the Article on the Institution's website or the Author's personal website, in each case provided that a link to the Article on the JoVE website is provided and notice of JoVE's copyright in the Article is included. All non-copyright intellectual property rights in and to the Article, such as patent rights, shall remain with the Author.

5. **Grant of Rights in Video - Standard Access.** This **Section 5** applies if the "Standard Access" box has been checked in **Item 1** above or if no box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby acknowledges and agrees that, subject to **Section 7** below, JoVE is and shall be the sole and exclusive owner of all rights of any nature, including, without limitation, all copyrights, in and to the Video. To the extent that, by law, the Author is deemed, now or at any time in the future, to have any rights of any nature in or to the Video, the Author hereby disclaims all such rights and transfers all such rights to JoVE.

6. **Grant of Rights in Video - Open Access.** This **Section 6** applies only if the "Open Access" box has been checked in **Item 1** above. In consideration of JoVE agreeing to produce, display or otherwise assist with the Video, the Author hereby grants to JoVE, subject to **Section 7** below, the exclusive, royalty-free, perpetual (for the full term of copyright in the Article, including any extensions thereto) license (a) to publish, reproduce, distribute, display and store the Video in all forms, formats and media whether now known or hereafter developed (including without limitation in print, digital and electronic form) throughout the world, (b) to translate the Video into other languages, create adaptations, summaries or extracts of the Video or other Derivative Works or Collective Works based on all or any portion of the Video and exercise all of the rights set forth in (a) above in such translations, adaptations, summaries, extracts, Derivative Works or Collective Works and (c) to license others to do any or all of the above. The foregoing rights may be exercised in all media and formats, whether now known or hereafter devised, and include the right to make such modifications as are technically necessary to exercise the rights in other media and formats.

7. **Government Employees.** If the Author is a United States government employee and the Article was prepared in the course of his or her duties as a United States government employee, as indicated in **Item 2** above, and any of the licenses or grants granted by the Author hereunder exceed the scope of the 17 U.S.C. 403, then the rights granted hereunder shall be limited to the maximum rights permitted under such statute. In such case, all provisions contained herein that are not in conflict with such statute shall remain in full force and effect, and all provisions contained herein that do so conflict shall be

deemed to be amended so as to provide to JoVE the maximum rights permissible within such statute.

8. **Protection of the work.** The Author(s) authorize JoVE to take steps in the Author(s) name and on their behalf if JoVE believes some third party could be infringing or might infringe the copyright of either the Author's Article and/or Video.

9. **Likeness, Privacy, Personality.** The Author hereby grants JoVE the right to use the Author's name, voice, likeness, picture, photograph, image, biography and performance in any way, commercial or otherwise, in connection with the Materials and the sale, promotion and distribution thereof. The Author hereby waives any and all rights he or she may have, relating to his or her appearance in the Video or otherwise relating to the Materials, under all applicable privacy, likeness, personality or similar laws.

10. **Author Warranties.** The Author represents and warrants that the Article is original, that it has not been published, that the copyright interest is owned by the Author (or, if more than one author is listed at the beginning of this Agreement, by such authors collectively) and has not been assigned, licensed, or otherwise transferred to any other party. The Author represents and warrants that the author(s) listed at the top of this Agreement are the only authors of the Materials. If more than one author is listed at the top of this Agreement and if any such author has not entered into a separate Article and Video License Agreement with JoVE relating to the Materials, the Author represents and warrants that the Author has been authorized by each of the other such authors to execute this Agreement on his or her behalf and to bind him or her with respect to the terms of this Agreement as if each of them had been a party hereto as an Author. The Author warrants that the use, reproduction, distribution, public or private performance or display, and/or modification of all or any portion of the Materials does not and will not violate, infringe and/or misappropriate the patent, trademark, intellectual property or other rights of any third party. The Author represents and warrants that it has and will continue to comply with all government, institutional and other regulations, including, without limitation all institutional, laboratory, hospital, ethical, human and animal treatment, privacy, and all other rules, regulations, laws, procedures or guidelines, applicable to the Materials, and that all research involving human and animal subjects has been approved by the Author's relevant institutional review board.

11. **JoVE Discretion.** If the Author requests the assistance of JoVE in producing the Video in the Author's facility, the Author shall ensure that the presence of JoVE employees, agents or independent contractors is in accordance with the relevant regulations of the Author's institution. If more than one author is listed at the beginning of this Agreement, JoVE may, in its sole discretion, elect not take any action with respect to the Article until such time as it has received complete, executed Article and Video License Agreements from each such author. JoVE reserves the right, in its absolute and sole

ARTICLE AND VIDEO LICENSE AGREEMENT - UK

discretion and without giving any reason therefore, to accept or decline any work submitted to JoVE. JoVE and its employees, agents and independent contractors shall have full, unfettered access to the facilities of the Author or of the Author's institution as necessary to make the Video, whether actually published or not. JoVE has sole discretion as to the method of making and publishing the Materials, including, without limitation, to all decisions regarding editing, lighting, filming, timing of publication, if any, length, quality, content and the like.

12. **Indemnification.** The Author agrees to indemnify JoVE and/or its successors and assigns from and against any and all claims, costs, and expenses, including attorney's fees, arising out of any breach of any warranty or other representations contained herein. The Author further agrees to indemnify and hold harmless JoVE from and against any and all claims, costs, and expenses, including attorney's fees, resulting from the breach by the Author of any representation or warranty contained herein or from allegations or instances of violation of intellectual property rights, damage to the Author's or the Author's institution's facilities, fraud, libel, defamation, research, equipment, experiments, property damage, personal injury, violations of institutional, laboratory, hospital, ethical, human and animal treatment, privacy or other rules, regulations, laws, procedures or guidelines, liabilities and other losses or damages related in any way to the submission of work to JoVE, making of videos by JoVE, or publication in JoVE or elsewhere by JoVE. The Author shall be responsible for, and shall hold JoVE harmless from, damages caused by lack of sterilization, lack of cleanliness or by contamination due to the making of a video by JoVE its employees, agents or independent contractors. All sterilization, cleanliness or


decontamination procedures shall be solely the responsibility of the Author and shall be undertaken at the Author's expense. All indemnifications provided herein shall include JoVE's attorney's fees and costs related to said losses or damages. Such indemnification and holding harmless shall include such losses or damages incurred by, or in connection with, acts or omissions of JoVE, its employees, agents or independent contractors.

13. **Fees.** To cover the cost incurred for publication, JoVE must receive payment before production and publication the Materials. Payment is due in 21 days of invoice. Should the Materials not be published due to an editorial or production decision, these funds will be returned to the Author. Withdrawal by the Author of any submitted Materials after final peer review approval will result in a US\$1,200 fee to cover pre-production expenses incurred by JoVE. If payment is not received by the completion of filming, production and publication of the Materials will be suspended until payment is received.

14. **Transfer, Governing Law.** This Agreement may be assigned by JoVE and shall inure to the benefits of any of JoVE's successors and assignees. This Agreement shall be governed and construed by the internal laws of the Commonwealth of Massachusetts without giving effect to any conflict of law provision thereunder. This Agreement may be executed in counterparts, each of which shall be deemed an original, but all of which together shall be deemed to be one and the same agreement. A signed copy of this Agreement delivered by facsimile, e-mail or other means of electronic transmission shall be deemed to have the same legal effect as delivery of an original signed copy of this Agreement.

A signed copy of this document must be sent with all new submissions. Only one Agreement is required per submission.

CORRESPONDING AUTHOR

Name:	Christoph Ballestrem	
Department:	Division of Cell Matrix Biology and Regenerative Medicine	
Institution:	University of Manchester	
Title:	Dr	
Signature:		Date: 12/06/2019

Please submit a **signed** and **dated** copy of this license by one of the following three methods:

1. Upload an electronic version on the JoVE submission site
2. Fax the document to +1.866.381.2236
3. Mail the document to JoVE / Attn: JoVE Editorial / 1 Alewife Center #200 / Cambridge, MA 02140

Rebuttal letter:***To the attention of JoVe Editors,***

We thank the reviewers for their constructive comments. Please find a detailed point-by-point response to the referee's comments.

Editor's comments:

All the changes suggested by the editor were done and we accompany our manuscript with a word document comparing the first submission version and the revised version.

Reviewers' comments:**Reviewer #1 and Reviewer #2:**

We thank the referees for the positive feedback.

Reviewer #3:

Manuscript Summary:

In this manuscript the authors describe a protocol for sequentially generating custom patterns of multiple proteins on glass substrates using the LIMAP technique, which is based on laser-assisted removal of the protective material (typically a PEG layer) such that proteins can bind on the glass substrate to form the desired pattern. The step-by-step protocol is split into several sections and described in detail.

Major Concerns:

1. The main issue with this manuscript is that it is relevant only to researchers who have access to a particular commercial equipment. As such, the title should be modified to include "using the PRIMO system". The authors fail to provide a comprehensive evaluation of other photopatterning techniques. In fact, parts of the manuscript sound like promotional material (e.g., lines 725-758), which is quite inappropriate for a scientific publication in my opinion. No matter how happy the authors are with the product they are using, they should avoid promoting it directly. To provide a fair and balanced review of the literature and other available techniques, the authors should critically evaluate these questions: (i) What other photopatterning techniques are available? Which of these are commercially available? (ii) What are their major technical and other limitations? (iii) How does the PRIMO system compare to existing systems/methods in terms of resolution, homogeneity, reproducibility, time, and cost? In the Introduction and Discussion sections, the authors briefly discuss some advantages and disadvantages of PRIMO over micro-contact printing, such as the facility of patterning multiple proteins, increased spatial resolution and increased process times. This analysis should be more elaborate and extend into other laser-assisted patterning techniques. One related issue is the use of PDMS stencils. These or similar stencils can be easily fabricated in any laboratory equipped with a scale, a vacuum desiccator and an oven. Using commercial stencils would be acceptable, but the authors should first explain the inexpensive/home-made option and then mention commercial options if necessary/preferred. It is particularly worrying that the stencils were

also bought from the same company. Same concern for the PLL-PEG solution, which should be available from multiple suppliers.

We are very grateful for the detailed comments of this referee to the manuscript. They were very constructive and helped us improving the protocol.

We have now taken many of the above suggestions on board.

- *As suggested we have now included “using the PRIMO system” in the title.*
- *Following suggestions, in the introduction, we have now mentioned other photopatterning techniques in the field (refs 12-15). However our main descriptive focus stays on LIMAP, and we have kept the previous comparison between microcontact printing and laser-assisted patterning in the discussion. We agree with the referee, that a more extensive evaluation between different micro-patterning techniques would be beneficial but we consider that type of broader analysis more appropriate for a review publication.*
- *As suggested by the referee, we have included the option to fabricate PDMS stencils in-house in order to show other options than the ones through Alveole. We had already included another PLL-PEG supplier in the Table of Materials in the manuscript.*

2. The technical limits of the system are poorly described. Specifically:

- In Fig. 3, the authors show an array of crossbows, patterned using LIMAP. First of all, these features are quite large (please provide minimum feature size) and it is not clear if the fidelity would be retained if the pattern was scaled down. This should be experimentally tested and reported. Second, there are obvious (visible to human eye) variations within and among different elements of the pattern (as evident from fluorescently-tagged fibrinogen). The authors should provide intra-pattern and inter-pattern variation (in terms of coefficient of variation) for this pattern, and do this as a function of minimum feature size.
- In Fig. 6, there is a clear decrease in signal intensity towards the ends of the stripes. Is this due to a failure in patterning or imaging? The authors should provide line profiles as in Fig. 8D and repeat this measurement (i.e., report SEM) to show variation within stripes and variation from process to process.
- In Fig. 7B and 7D, the data should be reported after background subtraction (by the way, the background signal is considerably high, can the authors estimate protein concentration outside the patterns and explain its cause). Also, line profiles should be given as in Fig. 8D.

We are grateful for the referee’s valid points. In the revised version we now provide:

- *The minimum feature sizes of the crossbows in Fig. 3. We have scaled down these patterns but only comment on the limitations in the discussion 915-916. Reproducibility is limited to the usage of a minimum stripe width of 2 μm .*
- *Towards Fig. 6, we now provide a detailed analysis of fidelity and stability of patterns. Fidelity of stripe patterns was retained upon scaling down; we present quantification. The intra- and inter-pattern variations were measured and the co-efficient of variation calculated. Quantifications are presented in Fig. 6E-H.*
- *In Fig. 6F line profiles were provided as requested. As recognised by the referee, a slight peripheral darkening effect was observed; this is due to the vignetting effect with the lenses we have used. We comment on this in lines 681-684 and 732-734 in the manuscript.*
- *7B and 7D data were presented after background subtraction. The patterning process was repeated three times and data of one representative experiment is shown in Fig 7B and D.*

We comment on potential background signals, which may vary depending on the type and concentration of proteins used (lines 973-984).

- Line profiles were not shown in Fig. 7B and D, to avoid repetition; line profiles are shown in Fig.6 and Fig.8.

3. The authors use CAD mouse neuroblastoma cells, which is a cell line that differentiates into neuron-like phenotypes and expresses certain neuron-specific proteins. However, to my knowledge, this cell line is not often (if ever) used to study axon elongation or pathfinding. The authors should avoid using the terms "neurons/neuronal" and "axons/axonal" which would be misleading (so is the "neuronal decision making assay"). The rationale for choosing this cell line (over more commonly used neuron-like cells, such as PC12 and SH-SY5Y) should be discussed in the manuscript. Figure 9 shows a few cells with elongated processes devoid of growth cones, thereby suggesting that the growth behavior of these cells may be drastically different from primary neurons. There is evidence in the literature that these cells can form growth cones (see 10.1083/jcb.201306067 for example); the lack of such structures in this study should also be discussed. It is also worrying that the authors mention biocompatibility of the PLPP (lines 681-688). In general, neuron-like cells much better tolerate subtle toxic insults compared to primary neurons. On the other hand, the authors state that they have not observed cytotoxicity in neurons. They need to provide evidence to support this claim by culturing neurons on the LMAP-generated pattern.

- Whilst CAD cells so far rarely appear in publications, we have now extensively characterised and compared them with other cell lines (including neuron-like PC12 and SH-SY5Y). We found that CAD cells (that derive from cortical neurons) are by far the most representative cell line in terms of similarity to primary cortical neurons (manuscript in preparation). Therefore they are an ideal cell for the purpose of pathfinding assays. As for primary cells, growth cones of CAD cells are enriched in F-actin and the shaft is enriched in microtubules (Fig. 10B).

- We have tested both primary and cell lines on our micropatterns; no toxic effect was seen in either of them. To demonstrate this we introduced a figure with primary cells in Fig. 10C.

- Following the referees suggestion, we now use the term "neurites" only for cell lines (as CAD cells) and the term "axons" for primary neurons.

Minor Concerns:

1. The first four sentences of the abstract are very general and hardly relevant to the manuscript.

We have revised the abstract to align it to the manuscript contents.

2. In Fig. 4, arrows are unnecessary and confusing. In panel D, the glass bottom can be referred to as "white" and the stencil as "grey". In panel E, the patterns can be shown in different colors. Also, please avoid using "highlighted" for two different things in the legend.

We agree and eliminated the arrows in Fig. 4 D and changed the wording in the legend. Regarding the colours used in Fig. 4E, we decided to keep the colours consistent throughout the Figures schemes with red for representation of the first patterned protein (horizontal) and green for the second patterned protein (vertical) (Fig. 1, 2, 4 and 5).

3. Protocol 1 is too detailed for a simple plasma treatment operation. For example, mentioning gas valve position is irrelevant for users of plasma systems other than the one used in this study. It would be sufficient (and necessary) to provide pressure, gas composition, power and plasma duration.

We shortened the plasma-cleaning protocol, as suggested and provided the requested information for pressure, gas composition, power and plasma duration.

4. It is not clear why two different passivation options are provided. Long-term storage is better than three-day storage; therefore PLL-PEG becomes redundant. Is there an advantage of using PLL-PEG? If so, please discuss.

We have outlined in the manuscript that using the long lasting PLL-PEG requires a longer passivation time and requires longer photopatterning steps. This can lead to a dramatic increase of time for larger experiments. We have added a more detailed explanation in lines 193-200 to make this clearer.

5. It is not intuitive to use "z-number"; "z-position" would be more precise/correct.

We are grateful for the comment and changed this according to the referee's suggestion.

6. A screenshot would be very helpful to accompany step 3.2.1.

We agree, but we have not included a screenshot because we believe this will all be shown in the JoVe video that will accompany the manuscript.

7. Parts of the protocol are repeated without an obvious reason. For example, step 4.8 is step 3.2.7 and step 8.4 is step 5.6. There may be repetitions in the process, but there is no point repeating the notes/warnings/caveats.

We have eliminated most repetitions following the referee's suggestions, when the information provided was redundant. However in some cases we have kept that information as a "Note" since we consider it key for the success of the protocol step (for example Protocol 4.2.7 and 5.8).

8. In step 6.1, the authors suggest testing different blocking reagents, but do not say how.

To clarify this point we have introduced a new Figure (Fig.9) and we have described more in detail how to test the blocking reagents in the text (709-714).

9. Last sentence of step 10.4 is not clear. Please rewrite.

We rephrased this sentence.

10. As a stylistic point, the authors can remove the last column of the Table of Materials and provide the information in the first column.

-Thank you very much. We have pointed out this issue to the Editor and left it to his choice whether to include that information in a column or as an independent paragraph below the Table of Materials.

11. The paragraph between lines 695 and 706 is quite vague. Which parameters in the Leonardo software? What is the gradient option? How would the UV adhesive be used? Please also provide dimensions instead of using superlatives such as "bigger" and "much larger".

We included now all requested information in lines 986-996.

12. What is the width of the stripes in Fig. 3A and their spacing?

We have stated now the width of the pattern and spacing in the Figure legend.

Kind regards,

The authors