

Materials List for:

Live Dissection of *Drosophila* Embryos: Streamlined Methods for Screening Mutant Collections by Antibody Staining

Hyung-Kook (Peter) Lee¹, Ashley P. Wright¹, Kai Zinn¹

¹Division of Biology, California Institute of Technology

Correspondence to: Hyung-Kook (Peter) Lee at hlee@caltech.edu, Ashley P. Wright at ashley@caltech.edu, Kai Zinn at zinnk@caltech.edu

URL: <https://www.jove.com/video/1647>

DOI: [doi:10.3791/1647](https://doi.org/10.3791/1647)

Materials

Name	Company	Catalog Number	Comments
Borosilicate Tubing With Filament	Sutter Instrument Co.	BF120-60-10	
Narishige model PC-10 needle puller	Narishige International		Tubes pulled using a heater level of 58.9