

Materials List for:

Determining the Reactivity and Titre of Serum using a Haemagglutination Assay

Maurizio Costabile¹

¹School of Pharmacy and Medical Sciences, University of South Australia

Correspondence to: Maurizio Costabile at maurizio.costabile@unisa.edu.au

URL: <https://www.jove.com/video/1752>

DOI: [doi:10.3791/1752](https://doi.org/10.3791/1752)

Materials

Name	Company	Catalog Number	Comments
Revercells	CSL	02490201	Use at 1% final
Anti-A serum	CSL	02611305	
96 well U bottom plate	TRP	92097	Must be U-bottom
Veronal buffered saline			Use at 1x final
37°C room/incubator			
Disposable gloves			